

线粒体应激研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	4E-BP1 Rabbit pAb	ER62532	Mouse, Rat, Monkey	WB,IHC-P
	4E-BP1 Mouse mAb [3G1]	RT1002	Human, Mouse, Rat	WB,IP,IF,IHC-P
Relative expression, Cell treatment	ATF3 Recombinant Rabbit mAb [PSH17-99]	HA723964	Human, Mouse	WB
Cell treatment	ATF-4 Recombinant Rabbit mAb [PSH07-08]	HA722796	Human, Mouse, Rat	WB,IP
Cell treatment	ATF-4 Recombinant Rabbit mAb [PSH07-08] - BSA and Azide free	HA751130	Human, Mouse, Rat	WB,IP
	ATF4 Recombinant Rabbit mAb [SD20-92]	ET1612-37	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	ATF4 Recombinant Rabbit mAb [SD20-92] - BSA and Azide free	HA750292	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	ATF5 Recombinant Rabbit mAb [SD2099]	ET1612-38	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	ATF5 Recombinant Rabbit mAb [SD2099] - BSA and Azide free	HA750293	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cell treatment	ATF6 Recombinant Rabbit mAb [PSH06-59] - BSA and Azide free	HA751087	Human, Mouse, Rat	WB,IHC-P,ChIP
Cell treatment	ATF6 Recombinant Rabbit mAb [PSH06-59]	HA722722	Human, Mouse, Rat	WB,IHC-P,ChIP
	ATF6 Mouse mAb [8D3]	EM1701-94	Human, Mouse, Rat	WB,IF-Cell
KD	ATF6 Recombinant Mouse mAb [8D3-R]	HA601321	Human, Mouse, Rat	WB,IF-Cell
KD	ATF6 Recombinant Mouse mAb [8D3-R] - BSA and Azide free	HA610193	Human, Mouse, Rat	WB,IF-Cell
	ATF6 Rabbit pAb	ER1706-34	Human, Rat	WB,IF-Cell,IHC-P,FC
	ATF6 beta Rabbit pAb	ER62756	Human	WB,IHC-P
	BACE1 Recombinant Rabbit mAb [PSH05-19] - BSA and Azide free	HA750990	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	BACE1 Recombinant Rabbit mAb [JE30-40]	HA721421	Human, Mouse, Rat	WB
	BACE1 Mouse mAb [7H4]	RT1066	Human, Mouse, Rat	WB,IP,IF,IHC-P
	BACE1 Recombinant Rabbit mAb [PSH05-19]	HA722244	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KO	Bax Recombinant Rabbit mAb [SZ3-07] - BSA and Azide free	HA750065	Human, Mouse, Rat, Goat	WB,IF-Cell,FC
KO	Bax Recombinant Rabbit mAb [SZ3-07]	ET1603-34	Human, Mouse, Rat, Goat	WB,IF-Cell,FC
KD	BAX Rabbit pAb	ER0907	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	BAX Recombinant Mouse mAb [PSH0-20] - BSA and Azide free	HA610038	Human	WB,IHC-P
	BAX Recombinant Mouse mAb [PSH0-20]	HA601123	Human	WB,IHC-P
	BAX Recombinant Rabbit mAb [PD01-44]	HA721324	Human	IHC-P,IF-Tissue
	BAX Mouse mAb [3-D2]	EM1203	Human, Mouse	WB,IF-Cell,IHC-P,FC
	BDNF Rabbit pAb	ER130915	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	BDNF Recombinant Rabbit mAb [JE59-35]	HA722912	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-P,IHC-Fr,FC
KD	BDNF Recombinant Rabbit mAb [SJ12-09]	ET1606-42	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC
KD	BDNF Recombinant Rabbit mAb [SJ12-09] - BSA and Azide free	HA750104	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC
KD	CLPP Recombinant Rabbit mAb [JE40-00]	HA721503	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Cpn10 Recombinant Rabbit mAb [JG82-34]	ET7108-67	Human, Mouse, Rat	WB,IHC-P,IP
	DNM1L Mouse mAb [C7-A4]	EM1708-45	Human, Mouse, Rat	WB,IF-Cell,FC
	DRP1 Rabbit pAb	HA500487	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
KD	eIF2A Recombinant Rabbit mAb [JE58-18]	ET7111-34	Human, Mouse, Rat	WB,IHC-P
	eIF2 beta Rabbit pAb	ER63330	Human, Mouse	WB,IHC-P
	EIF2S1 Rabbit pAb	HA500385	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell,IP
	EIF2S1 Recombinant Rabbit mAb [PSH04-29]	HA722112	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP
	EIF2S1 Recombinant Rabbit mAb [PSH04-29] - BSA and Azide free	HA750938	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP
	eIF4EBP1 Rabbit pAb	ER1706-64	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	eIF4EBP1 Recombinant Rabbit mAb [JJ09-25]	ET1701-83	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
KD	ERK5 Recombinant Rabbit mAb [SD2084] - BSA and Azide free	HA750299	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IP,FC
KD	ERK5 Recombinant Rabbit mAb [SD2084]	ET1612-7	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	ERO1L Recombinant Rabbit mAb [PSH07-30]	HA722818	Human, Mouse, Rat	WB,IHC-P
	ERO1L Rabbit pAb	ER1908-32	Human	WB
	ERO1L Recombinant Rabbit mAb [PSH07-30] - BSA and Azide free	HA751151	Human, Mouse, Rat	WB,IHC-P
	FGF21 Recombinant Rabbit mAb [JA10-97] - BSA and Azide free	HA750399	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	FGF21 Recombinant Rabbit mAb [JA10-97]	ET1704-04	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	FH Rabbit pAb	ER1901-10	Human, Mouse, Rat	WB,IHC-P
	FH Mouse mAb [7F1]	EM60064	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	FH Rabbit pAb	ER1902-01	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	FoxO1/3/4 Rabbit pAb	ER63465	Human, Mouse, Rat	WB,IHC-P
KD	FOXO1A Recombinant Rabbit mAb [SU33-01] - BSA and Azide free	HA750140	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
KD	FOXO1A Recombinant Rabbit mAb [SU33-01]	ET1608-25	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	iFluor™ 488 Conjugated FOXO1A Recombinant Rabbit mAb [SU33-01]	HA720161F	Human, Mouse, Rat	IF-Cell,FC
	FOXO3A Rabbit pAb	ER1706-79	Human, Mouse	WB,IF-Cell,IHC-P,FC
	FOXO3A Recombinant Rabbit mAb [SP0054]	ET1604-11	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue
	FOXO4 Recombinant Rabbit mAb [SD0817]	ET1612-66	Human, Rat	WB,IHC-P
	FOXO4 Recombinant Rabbit mAb [JJ09-11]	ET1701-67	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC,IP
	Fumarate hydratase Recombinant Rabbit mAb [PSH08-80]	HA723032	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	GADD34 Rabbit pAb	ER65755	Human	WB
	GADD34 Rabbit pAb	ER2001-56	Human	WB,IHC-P,FC
KD	GADD34 Mouse mAb [A4D7]	EM1902-26	Human	WB,IHC-P
Relative expression	GDF15 Recombinant Rabbit mAb [PSH10-09] - BSA and Azide free	HA751339	Human	WB,IHC-P
Relative expression	GDF15 Recombinant Rabbit mAb [PSH10-09]	HA723178	Human	WB,IHC-P
	GDF15 Recombinant Rabbit mAb [JE39-90]	HA722551	Human	WB
	HSP60 Recombinant Rabbit mAb [ST48-04] - BSA and Azide free	HA750175	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP
	HSP60 Rabbit pAb	1007-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	HSP60 Mouse mAb [9-49]	M1007-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HSP60 Mouse mAb [1-80]	EM00704	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	HSP60 Rabbit pAb	1007-1	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	HSP60 Recombinant Rabbit mAb [ST48-04]	ET1609-45	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP
KD	Hsp70 Recombinant Rabbit mAb [SA0379]	ET1601-11	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Hsp70 Recombinant Rabbit mAb [PD01-65]	HA721497	Human	WB,IHC-P,IF-Tissue
	Hsp70 Rabbit pAb	ER50802	Human, Mouse, Rat, Cow	WB,IF-Cell,IHC-P,FC
	HSP70 Mouse mAb [3G10]	EM60062	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	HSP70 (Acetyl K246) Rabbit pAb	ER60056	Human	WB
KD	HtrA2 / Omi Recombinant Rabbit mAb [ST54-05]	ET1609-73	Human, Mouse, Rat	WB,IHC-P,IP
	Human BACE1 Recombinant Rabbit mAb [PSH16-23] - BSA and Azide free (Detector)	HA725172	Human	ELISA(Det)
	Biotin Conjugated Human BACE1 Recombinant Rabbit mAb [PSH16-23] - Detector	HA725172B	Human	ELISA(Det),ELISA
	Human BACE1 Recombinant Rabbit mAb [PSH16-22] - BSA and Azide free (Capture)	HA725171	Human	ELISA(Cap)
	Human GDF-15 Recombinant Rabbit mAb [PSH16-00] - BSA and Azide free (Capture)	HA723805	Human	ELISA(Cap)
	Human GDF-15 Recombinant Rabbit mAb [PSH16-01] - BSA and Azide free (Detector)	HA723806	Human	ELISA(Det)
	Biotin Conjugated Human GDF-15 Recombinant Rabbit mAb [PSH16-02] - Detector	HA723807B	Human	ELISA(Det),ELISA
	Human GDF-15 Recombinant Rabbit mAb [PSH16-02] - BSA and Azide free (Detector)	HA723807	Human	ELISA(Det)
	Biotin Conjugated Human GDF-15 Recombinant Rabbit mAb [PSH16-01] - Detector	HA723806B	Human	ELISA(Det),ELISA
	Human IL-18 Recombinant Rabbit mAb [PSH01-70] - BSA and Azide free (Detector)	HA721723	Human	ELISA(Det)
Relative expression	Human IL-18 Recombinant Rabbit mAb [PSH02-09]	HA721768	Human	WB
	Human IL-18 Recombinant Rabbit mAb [PSH01-71] - BSA and Azide free (Capture)	HA721722	Human	ELISA(Cap)
Relative expression	Human IL-18 Recombinant Rabbit mAb [PSH02-09] - BSA and Azide free	HA750779	Human	WB
	Human IL-18 Recombinant Rabbit mAb [PSH02-08] - BSA and Azide free	HA750778	Human	WB,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-18 Recombinant Rabbit mAb [PSH02-08]	HA721767	Human	WB,IHC-P,IF-Tissue,IP
	Biotin Conjugated Human IL-18 Recombinant Rabbit mAb [PSH01-70]	HA722644B	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA721713	Human	WB
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - Detector	HA723882B	Human	ELISA(Det),ELISA
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-87] - BSA and Azide free (Capture)	HA723883	Human	ELISA(Cap)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - BSA and Azide free (Detector)	HA723882	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39]	HA721712	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free (Detector)	HA721279	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free	HA750745	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free (Capture)	HA721278	Human	ELISA(Cap)
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA722649B	Human	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-85] - BSA and Azide free (Capture)	HA723881	Human	ELISA(Cap)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free	HA750744	Human	WB
	Biotin Conjugated Human NLRP3 Recombinant Rabbit mAb [PSH11-09] - Detector	HA723279B	Human	ELISA(Det),ELISA
	Human NLRP3 Recombinant Rabbit mAb [PSH11-07] - BSA and Azide free (Capture)	HA723275	Human	ELISA(Cap)
	Biotin Conjugated Human NLRP3 Recombinant Rabbit mAb [PSH11-08] - Detector	HA723277B	Human	ELISA(Det),ELISA
	Human NLRP3 Recombinant Rabbit mAb [PSH11-08] - BSA and Azide free (Detector)	HA723276	Human	ELISA(Det)
	Human NLRP3 Recombinant Rabbit mAb [PSH11-09] - BSA and Azide free (Detector)	HA723278	Human	ELISA(Det)
	IKB beta Recombinant Rabbit mAb [SN0648]	ET1611-1	Human, Monkey	WB,IF-Cell,IHC-P,IP,FC
	IKK alpha Rabbit pAb	ER30911	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	IKK alpha Recombinant Rabbit mAb [SN06-99] - BSA and Azide free	HA750246	Human, Mouse, Rat	WB,IHC-P,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	IKK alpha Recombinant Rabbit mAb [SN06-99]	ET1611-15	Human, Mouse, Rat	WB,IHC-P,IP,FC
KD	IKK alpha + IKK beta Recombinant Rabbit mAb [SN63-02]	ET1611-23	Human, Mouse, Rat	WB,IF-Cell,IP
KD	IKK alpha + IKK beta Recombinant Rabbit mAb [SN63-02] - BSA and Azide free	HA750248	Human, Mouse, Rat	WB,IF-Cell,IP
	IKK beta Rabbit pAb	R1706-13	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	IKK beta Recombinant Rabbit mAb [JE32-81]	HA721715	Human	WB,FC
Relative expression	IL-18 Recombinant Rabbit mAb [PSH18-15]	HA723978	Human, Mouse, Rat	WB,IHC-P,IP
	IL-18 Recombinant Rabbit mAb [PS01-06]	HA721537	Human	WB,FC
	IL-18 Recombinant Rabbit mAb [PS00-41] - BSA and Azide free	HA750676	Human	WB,IHC-P
	IL-18 Recombinant Rabbit mAb [PS01-06] - BSA and Azide free	HA750677	Human	WB,FC
	IL-18 Recombinant Rabbit mAb [PS00-41]	HA721536	Human	WB,IHC-P
	IL-1 beta Recombinant Rabbit mAb [PSH18-00]	HA723965	Mouse, Rat	WB
Cell treatment	IL-1 beta Mouse mAb [A7F9]	HA601002	Human, Mouse, Rat	WB,ELISA,IF-Cell
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R]	HA601036	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Recombinant Rabbit mAb [JJ087-3]	ET1701-39	Human	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	IL-1 beta Mouse mAb [A7A10]	HA601035	Human	WB,ELISA,IF-Cell
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R] - BSA and Azide free	HA610004	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
	IRE1 Rabbit pAb	ER1902-90	Human, Mouse, Rat	WB,IHC-P
	IRE1 Recombinant Rabbit mAb [PSH10-60] - BSA and Azide free	HA751363	Human, Mouse, Rat, Monkey	WB
	IRE1 Recombinant Rabbit mAb [PSH10-60]	HA723225	Human, Mouse, Rat, Monkey	WB
	IRE1 Mouse mAb [C9-F2]	EM1708-88	Human	WB,IHC-P
	JNK1 Recombinant Rabbit mAb [JE63-04]	HA721899	Human, Mouse, Rat	WB,IHC-P
	JNK1 Rabbit pAb	R1309-1	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	JNK1+JNK2+JNK3 Recombinant Rabbit mAb [SA43-06] - BSA and Azide free	HA750018	Human, Mouse, Rat, Cow, Monkey, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	JNK1+JNK2+JNK3 Recombinant Rabbit mAb [SA43-06]	ET1601-28	Human, Mouse, Rat, Cow, Monkey, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP
	JNK1+JNK3 Recombinant Rabbit mAb [JJ08-86]	ET1701-47	Human, Mouse, Rat	WB,IP,IHC-P
	JNK2 Rabbit pAb	HA500110	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	JNK2 Recombinant Rabbit mAb [SC05-85]	ET1610-11	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	JNK3 Recombinant Rabbit mAb [SD082-09]	ET1612-68	Human, Mouse, Rat	WB,IF-Tissue,IHC-P
KD	Keap1 Recombinant Rabbit mAb [PSH0-90]	HA721525	Human, Mouse, Rat	WB,IHC-P,FC
KD	Keap1 Recombinant Rabbit mAb [PSH0-90] - BSA and Azide free	HA750670	Human, Mouse, Rat	WB,IHC-P,FC
	MK13 Rabbit pAb	ER60146	Human, Mouse, Rat	WB
	Mouse FGF-21 Recombinant Rabbit mAb [PSH15-26] - BSA and Azide free (Capture)	HA723731	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse FGF-21 Recombinant Rabbit mAb [PSH15-27] - Detector	HA723733B	Mouse	ELISA(Det),ELISA
	Mouse FGF-21 Recombinant Rabbit mAb [PSH15-27] - BSA and Azide free (Detector)	HA723732	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - BSA and Azide free (Detector)	HA721977	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH09-59] - BSA and Azide free	HA723121	Mouse	ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-31] - BSA and Azide free (Capture)	HA721976	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - Detector	HA722625B	Mouse	ELISA(Det),ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-66] - BSA and Azide free (Capture)	HA722166	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67]	HA722608B	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH09-58] - BSA and Azide free	HA723120	Mouse	ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67] - BSA and Azide free (Detector)	HA722167	Mouse	ELISA(Det)
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27]	HA721307	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04] - BSA and Azide free	HA750058	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	NF-kB p65 Rabbit pAb	HA500402	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27] - BSA and Azide free	HA750601	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	NF-kB p65 Rabbit pAb	ER0815	Human, Mouse, Rat, Zebrafish	WB,IHC-P,FC
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04]	ET1603-12	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cell treatment	NLRP3 Recombinant Rabbit mAb [SC06-23]	ET1610-93	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	NLRP3 Recombinant Rabbit mAb [PSH19-43]	HA724075	Human, Mouse	WB
Cell treatment	NLRP3 Recombinant Rabbit mAb [SC06-23] - BSA and Azide free	HA750236	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	NLRP3 Rabbit pAb	ER1706-72	Human, Mouse	WB,IF-Cell,IHC-P,FC
	NRF1 Recombinant Rabbit mAb [JM89-63]	ET1705-86	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	NRF1 Recombinant Rabbit mAb [JM89-63] - BSA and Azide free	HA750444	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Cell treatment	Nrf2 Recombinant Rabbit mAb [PSH11-20] - BSA and Azide free	HA751401	Human, Mouse, Rat	WB,IF-Cell,ChIP
Cell treatment	Nrf2 Recombinant Rabbit mAb [PSH11-20]	HA723302	Human, Mouse, Rat	WB,IF-Cell,ChIP
Cell treatment	Nrf2 Rabbit pAb	ER1706-41	Human, Mouse, Rat	IHC-P,IF-Cell,WB
	Nrf2 Rabbit pAb	R1312-8	Human, Mouse	WB,IHC-P,IF-Tissue
Cell treatment	Nrf2 Recombinant Rabbit mAb [JE30-14]	HA721432	Human, Mouse	WB,IF-Cell
	OMA1 Rabbit pAb	ER1914-36	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Cell treatment	OPA1 Recombinant Rabbit mAb [PSH06-39] - BSA and Azide free	HA751078	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
Cell treatment	OPA1 Recombinant Rabbit mAb [PSH06-39]	HA722673	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	OPA1 Recombinant Rabbit mAb [JM81-35]	ET1705-9	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC
	OPA1 Recombinant Rabbit mAb [JM81-35] - BSA and Azide free	HA750446	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC
KO	p38 alpha / MAPK14 Recombinant Rabbit mAb [JF55-07]	ET1702-65	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD, Cell treatment	p38 alpha / MAPK14 Recombinant Rabbit mAb [SR43-04]	ET1602-26	Human, Mouse, Rat	WB,IF-Cell,FC
	p38 alpha / MAPK14 Mouse mAb [9-B9]	EM40307	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	p38 alpha / MAPK14 Rabbit pAb	R1308-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD, Cell treatment	p38 alpha / MAPK14 Recombinant Rabbit mAb [SR43-04] - BSA and Azide free	HA750040	Human, Mouse, Rat	WB,IF-Cell,FC
	p38 alpha/MAPK14 + p38 beta/MAPK11 + p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-60]	HA722019	Human, Mouse, Rat, Monkey	WB
	p38 alpha/MAPK14 + p38 beta/MAPK11 + p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-60] - BSA and Azide free	HA750894	Human, Mouse, Rat, Monkey	WB
	p38 beta/MAPK11 + p38 gamma/MAPK12 + p38 delta/MAPK13 Recombinant Rabbit mAb [PSH03-59] - BSA and Azide free	HA750893	Human, Mouse, Rat, Monkey	WB
	p38 beta/MAPK11 + p38 gamma/MAPK12 + p38 delta/MAPK13 Recombinant Rabbit mAb [PSH03-59]	HA722018	Human, Mouse, Rat, Monkey	WB
	p38 delta / MAPK13 Rabbit pAb	HA500582	Human	WB
	p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-58] - BSA and Azide free	HA750892	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-58]	HA722017	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	p38 MAPK Recombinant Rabbit Multiclonal Antibody [PSH03-72]	HA722036	Human, Mouse, Rat, Monkey	WB
	p38 MAPK Recombinant Rabbit Multiclonal Antibody [PSH03-72] - BSA and Azide free	HA750906	Human, Mouse, Rat, Monkey	WB
	PERK Rabbit pAb	HA500326	Human	WB
	PERK Recombinant Rabbit mAb [PSH0-81]	HA721510	Human	WB
	PERK Rabbit pAb	ER64553	Human, Mouse, Rat	IF-Cell,WB,IHC-P,IF-Tissue
	PERK Recombinant Rabbit mAb [PSH0-81] - BSA and Azide free	HA750661	Human	WB
KD	Peroxiredoxin 3 Recombinant Rabbit mAb [JA53-21]	ET1704-92	Human	WB,IF-Cell,IHC-P,FC,IP
	PGC1 alpha+beta Recombinant Rabbit mAb [JF09-71]	ET1702-96	Human, Mouse, Rat	WB
	PGC1 beta Recombinant Rabbit mAb [JM33-24]	ET1705-65	Human, Rat	WB,IP,FC
	Phospho-4E-BP1/2/3 (T45) Mouse mAb [3G2]	RT1004	Human, Mouse, Rat	WB,IP,IF
Cell treatment	Phospho-D RP1 (S616) Recombinant Rabbit mAb [PSH06-77]	HA722760	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC
Cell treatment	Phospho-D RP1 (S616) Recombinant Rabbit mAb [PSH06-77] - BSA and Azide free	HA751109	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment, Relative expression	Phospho-EIF2S1 (S51) Recombinant Rabbit mAb [SZ01-06] - BSA and Azide free	HA750059	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cell treatment, Relative expression	Phospho-EIF2S1 (S51) Recombinant Rabbit mAb [SZ01-06]	ET1603-14	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Phospho-EIF2S1 (S52) Rabbit pAb	RT1198	Human, Mouse, Rat	WB,IP,IF,IHC-P
Cell treatment	Phospho-eIF4EBP1 (S65) Recombinant Rabbit mAb [PSH09-44] - BSA and Azide free	HA751291	Human	WB,IF-Cell,IP
Cell treatment	Phospho-eIF4EBP1 (S65) Recombinant Rabbit mAb [PSH09-44]	HA723106	Human	WB,IF-Cell,IP
Cell treatment	Phospho-eIF4EBP1 (T70) Recombinant Rabbit mAb [PSH09-91]	HA723222	Human, Mouse	WB,IHC-P,IP
Cell treatment	Phospho-eIF4EBP1 (T70) Recombinant Rabbit mAb [PSH09-91] - BSA and Azide free	HA751360	Human, Mouse	WB,IHC-P,IP
Cell treatment	Phospho-FOXO1A (T24)+ FOXO3A (T32)+ FOXO4 (T28)+ FOXO6 (T26) Recombinant Rabbit mAb [PSH06-00]	HA722499	Human, Mouse	WB,IF-Cell,FC
Cell treatment	Phospho-FOXO1A (T24)+ FOXO3A (T32)+ FOXO4 (T28)+ FOXO6 (T26) Recombinant Rabbit mAb [PSH06-00] - BSA and Azide free	HA751040	Human, Mouse	WB,IF-Cell,FC
Cell treatment	Phospho-FOXO3a (S253) Recombinant Rabbit mAb [ST49-01]	ET1609-49	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-FOXO3a (S253) Recombinant Rabbit mAb [ST49-01] - BSA and Azide free	HA750178	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-IRE1 (S724) Recombinant Rabbit mAb [PSH03-35]	HA721980	Human, Mouse	WB,IF-Cell,IHC-P,FC
Cell treatment	Phospho-IRE1 (S724) Recombinant Rabbit mAb [PSH03-35] - BSA and Azide free	HA750869	Human, Mouse	WB,IF-Cell,IHC-P,FC
Cell treatment	Phospho-JNK1/2/3 (T183 + T183 + T221) Recombinant Rabbit mAb [ST500]	ET1609-42	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
Cell treatment	Phospho-JNK1/2/3 (T183 + T183 + T221) Recombinant Rabbit mAb [ST500] - BSA and Azide free	HA750173	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
Cell treatment	Phospho-NF-kB p65 (S529) Recombinant Rabbit mAb [SP07-00]	ET1604-27	Human	WB,IHC-P,FC,IP
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58] - BSA and Azide free	HA751361	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58]	HA723223	Human, Mouse, Rat	WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment, Relative expression	Phospho-Nrf2 (S40) Recombinant Rabbit mAb [SU0334] - BSA and Azide free	HA750142	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Cell treatment, Relative expression	Phospho-Nrf2 (S40) Recombinant Rabbit mAb [SU0334]	ET1608-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Cell treatment	Phospho-p38 alpha (T180) Recombinant Rabbit mAb [PSH06-25]	HA722657	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180) Recombinant Rabbit mAb [PSH06-25] - BSA and Azide free	HA751062	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180 + Y182) Recombinant Rabbit mAb [PSH06-35] - BSA and Azide free	HA751074	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180 + Y182) Recombinant Rabbit mAb [PSH06-35]	HA722669	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-p38 (T180) Rabbit pAb	ER2001-51	Human, Mouse, Rat	WB
Cell treatment	Phospho-p38 (T180) Recombinant Rabbit mAb [JE77-37]	HA722150	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-p38 (T180 + Y182) Rabbit pAb	ER2001-52	Human, Mouse, Rat	WB
	Phospho-p38 (T180 + Y182) Rabbit pAb	ER1903-01	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
Cell treatment	Phospho-SIRT1 (S47) Recombinant Rabbit mAb [SN06-40]	ET1611-31	Human	WB,IF-Cell,IHC-P
	Phospho-SIRT1 (T530) Recombinant Rabbit mAb [JJ206-6]	ET1701-27	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	PUMA Recombinant Rabbit mAb [SR42-09] - BSA and Azide free	HA750038	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	PUMA Recombinant Rabbit mAb [SR42-09]	ET1602-24	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	PUMA Rabbit pAb	ER31215	Human, Mouse	WB,IF-Cell,IHC-P,FC
	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-59] - BSA and Azide free (Capture)	HA721836	Rat	ELISA(Cap)
	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60] - BSA and Azide free (Detector)	HA721837	Rat	ELISA(Det)
	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60]	HA722634B	Rat	ELISA(Det)
	Rat IL-1 beta Recombinant Rabbit mAb [PS01-86] - BSA and Azide free (Capture)	HA721609	Rat	ELISA(Cap)
	Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - BSA and Azide free (Detector)	HA721610	Rat	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - Detector	HA721611	Rat	ELISA(Det),ELISA
	S2P Rabbit pAb	ER64902	Human, Mouse, Rat	IF-Cell
KD	SDHA Recombinant Rabbit mAb [JM10-83]	ET1703-40	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC
	SEC61B Rabbit pAb	ER61055	Human, Mouse	WB
	SEC61B Recombinant Rabbit mAb [PSH09-48] - BSA and Azide free	HA751295	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,IF-Tissue
	SEC61B Recombinant Rabbit mAb [PSH09-48]	HA723110	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,IF-Tissue
KD	SIRT1 Recombinant Rabbit mAb [JE11-04]	ET7110-49	Human	WB,IF-Cell
KD	SIRT1 Recombinant Rabbit mAb [SZ04-01] - BSA and Azide free	HA750064	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP
KD	SIRT1 Rabbit pAb	ER130811	Human, Mouse	WB,IHC-P,FC,IF-Cell
KD	SIRT1 Recombinant Rabbit mAb [SZ04-01]	ET1603-3	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP
KD	SIRT1 Mouse mAb [7-C5-B2]	M1506-3	Human, Mouse	WB,IF-Cell,IHC-P
KD	SIRT2 Recombinant Rabbit mAb [JJ08-83]	ET1701-43	Human	WB
	SIRT2 Recombinant Rabbit mAb [SN70-04]	ET1611-72	Human, Cynomolgus monkey	WB,IP,IHC-P
	SIRT2 Recombinant Rabbit mAb [PSH17-48]	HA723923	Human, Mouse, Rat, Cynomolgus monkey	IHC-Fr,IHC-P,WB
KD	SIRT3 Recombinant Rabbit mAb [PSH05-22] - BSA and Azide free	HA750993	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue,IP,ChIP
KD	SIRT3 Recombinant Rabbit mAb [PSH05-22]	HA722251	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue,IP,ChIP
	SIRT3 Rabbit pAb	R1511-3	Human, Mouse, Rat	WB
	TRIB3 Recombinant Rabbit mAb [PSH02-13] - BSA and Azide free	HA750783	Human	WB
	TRIB3 Rabbit pAb	ER61692	Human	WB
	TRIB3 Recombinant Rabbit mAb [PSH02-13]	HA721772	Human	WB
	TXNRD2 Recombinant Rabbit mAb [PSH02-27] - BSA and Azide free	HA750795	Human, Mouse, Rat	WB,IHC-P
KD	TXNRD2 Recombinant Rabbit mAb [PSH02-26]	HA721806	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
KD	TXNRD2 Recombinant Rabbit mAb [PSH02-26] - BSA and Azide free	HA750794	Human, Mouse, Rat	WB,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	TXNRD2 Recombinant Rabbit mAb [PSH02-27]	HA721807	Human, Mouse, Rat	WB,IHC-P
	UCP-1 Rabbit pAb	ER1918-37	Mouse	WB,FC
Relative expression	UCP1 Recombinant Rabbit mAb [PSH15-73] - BSA and Azide free	HA751580	Mouse, Rat	WB,IHC-Fr,IHC-P,IP
	UCP1 Rabbit pAb	ER62425	Human, Mouse, Rat	WB
Relative expression	UCP1 Recombinant Rabbit mAb [PSH15-73]	HA723776	Mouse, Rat	WB,IHC-Fr,IHC-P,IP
	UCP1 Rabbit pAb	HA500118	Human, Mouse, Rat	WB,IHC-P