

细胞焦亡研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
Cell treatment, KD	Active+Pro Caspase-9 Recombinant Rabbit mAb [JJ08-05]	ET1701-22	Human	WB,IHC-P,FC
Cell treatment, KD	Active+Pro Caspase-9 Recombinant Rabbit mAb [SC65-05]	ET1610-95	Human, Mouse	WB,IP
Cell treatment, KD	Active+Pro Caspase-9 Recombinant Rabbit mAb [SC65-05] - BSA and Azide free	HA750238	Human, Mouse	WB,IP
Cell treatment, KD	Active+Pro Caspase-9 Rabbit pAb	R1308-12	Human	WB,IHC-P
	AIM2 Rabbit pAb	HA500270	Human	WB,IF-Cell,FC
	AKT1 Recombinant Rabbit mAb [PSH0-50] - BSA and Azide free	HA750627	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
KD	AKT1 Recombinant Rabbit mAb [ST05-09]	ET1609-47	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
	AKT1 Recombinant Rabbit mAb [PSH0-50]	HA721364	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
KD	AKT1 Recombinant Rabbit mAb [ST05-09] - BSA and Azide free	HA750177	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
	AKT1 Mouse mAb [D9-9-C9]	EM40507	Human, Mouse, Rat	WB,IHC-P
	AKT1/2/3 Rabbit pAb	ER1901-33	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	AKT1/2/3 Recombinant Rabbit mAb [ST48-09]	ET1609-51	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75-09] - BSA and Azide free	HA750840	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	AKT1/2/3 Recombinant Rabbit mAb [ST48-09] - BSA and Azide free	HA750179	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75-09]	HA721870	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	Akt1/3 Rabbit pAb	ER62638	Human, Mouse, Rat	WB,IHC-P
	Akt2 Mouse mAb [2A1]	RT1027	Human, Mouse, Rat	WB,IP,IF,IHC-P
	AKT2 Rabbit pAb	HA500091	Human, Mouse, Rat	WB,IHC-P
	Annexin A1 Rabbit pAb	ER1803-70	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Annexin A1 Recombinant Rabbit mAb [PSH04-13]	HA722094	Human	IHC-P,IF-Tissue
KO	ATG5 Recombinant Rabbit mAb [SN73-07] - BSA and Azide free	HA750251	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO	ATG5 Recombinant Rabbit mAb [SN73-07]	ET1611-38	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	Bak Recombinant Rabbit mAb [SU32-07]	ET1608-21	Human	WB,IF-Tissue,IHC-P,IP
	Beta III Tubulin Rabbit pAb	0807-2	Human, Mouse, Rat	WB,IF-Cell,FC
Relative expression	Beta III Tubulin Recombinant Rabbit mAb [PSH15-32] - BSA and Azide free	HA751554	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP
	iFluor™ 488 Conjugated Beta III Tubulin Mouse mAb [A8-D10]	HA600106F	Human, Mouse, Rat	IF-Cell,FC,IF-Tissue
	Beta III Tubulin Rabbit pAb	R1401-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	Beta III Tubulin Recombinant Rabbit mAb [PSH15-32]	HA723738	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP
	iFluor™ 594 Conjugated Beta III Tubulin Mouse mAb [A8-D10]	HA600133F	Human, Rat	IF,ICC,FC
	Beta III Tubulin Recombinant Rabbit mAb [SP06-00]	ET1604-17	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
	Beta III Tubulin Mouse mAb [A8-D10]	M0805-8	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
	Beta III Tubulin Rabbit pAb	0805-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
	Beta III Tubulin Recombinant Rabbit mAb [SP06-00] - BSA and Azide free	HA750071	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
	Carbonic anhydrase 9 Recombinant Rabbit mAb [JJ088-9]	ET1701-51	Human	WB,IHC-P,IF-Tissue,mIHC
	Caspase-1 Recombinant Rabbit mAb [SU40-07]	ET1608-69	Human, Mouse, Rat	WB,FC,IF-Cell
	Caspase-1 Recombinant Rabbit mAb [SU40-07] - BSA and Azide free	HA750161	Human, Mouse, Rat	WB,FC,IF-Cell
	Caspase-1 Rabbit pAb	R1510-23	Human	IF-Cell,IHC-P,FC
	Caspase-1 Rabbit pAb	R1510-22	Human	IHC-P,FC
Relative expression	Caspase-1 + p10 + p12 Recombinant Rabbit mAb [JE56-35]	HA722222	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Caspase-1 p20 Rabbit pAb	ER1905-47	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
Cell treatment, KD	Caspase-3 p12 Recombinant Rabbit mAb [JE07-51]	HA722235	Human, Mouse, Rat	WB,IHC-P
	Caspase-4 Mouse mAb [A9G2]	HA601431	Human, Cynomolgus monkey	WB
	Caspase-4 Rabbit pAb	ER62908	Human	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Caspase-5 Recombinant Rabbit mAb [SD203-2]	ET1612-29	Human	WB,IF-Tissue,IHC-P,IP
Cell treatment	Caspase-7 Recombinant Rabbit mAb [PSH13-69] - BSA and Azide free	HA751482	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Caspase-7 Recombinant Rabbit mAb [PSH13-69]	HA723558	Human, Mouse, Rat	WB,IHC-P
Cell treatment, KD	Caspase-8 Recombinant Rabbit mAb [PSH04-77] - BSA and Azide free	HA750971	Human	WB,IF-Cell,IHC-P,FC,IP
Relative expression, KD	Caspase-8 Recombinant Rabbit mAb [SD08-06]	ET1612-70	Human	WB,IF-Tissue,IHC-P
	Caspase-8 Mouse mAb [6G1]	RT1099	Human	WB,IP,IF,IHC-P
Cell treatment, KD	Caspase-8 Recombinant Rabbit mAb [PSH04-77]	HA722181	Human	WB,IF-Cell,IHC-P,FC,IP
Relative expression, KD	Caspase-8 Recombinant Rabbit mAb [SZ01-08]	ET1603-16	Human	WB,IF-Tissue,IHC-P
	Caspase-9 Recombinant Rabbit mAb [SZ29-01]	ET1603-27	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Caspase-9 Recombinant Rabbit mAb [SZ29-01] - BSA and Azide free	HA750062	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	CCL2 Mouse mAb [E11-D8]	EM1710-22	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
	CCL2 / MCP-1 Rabbit pAb	HA500042	Human	WB,IHC-P,IF-Cell
	CCL2 / MCP-1 Rabbit pAb	HA500267	Mouse	WB,IF-Cell,FC
Cell treatment	CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-16]	HA723979	Mouse, Rat	WB,IF-Cell,FC
	CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-44]	HA723691	Mouse	WB,IF-Cell,FC
	CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-44] - BSA and Azide free	HA751533	Mouse	WB,IF-Cell,FC
	Cdk9 Recombinant Rabbit mAb [SD204-07]	ET1612-78	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Cleaved Caspase-3 Rabbit pAb	IRS227RB	Mouse	mIHC
Cell treatment	Cleaved Caspase-3 Recombinant Rabbit mAb [SR01-02] - BSA and Azide free	HA750053	Human, Pig	WB,IF-Cell,IHC-P
Cell treatment	Cleaved Caspase-3 Recombinant Rabbit mAb [SR01-02]	ET1602-47	Human, Pig	WB,IF-Cell,IHC-P
	Cleaved Caspase-3 Recombinant Rabbit mAb	IRS136RB	Human	mIHC
	Cleaved-Caspase-4/5 p20 (D270/D311) Rabbit pAb	ER60012	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Cleaved-Caspase-4 p20 (Q81) Rabbit pAb	ER60011	Human, Mouse	WB,IHC-P
	Cleaved Caspase-7 (Ser199) Rabbit pAb	ER60002	Human, Mouse	WB,IHC-P,IF-Cell
	Cleaved-Caspase-9 (D353) Rabbit pAb	ER60008	Mouse, Rat, Human	IF-Cell,WB,IHC-P
Cell treatment, KD	Cleaved PARP Recombinant Rabbit mAb [SU0314]	ET1608-10	Human	WB,IF-Cell,IP
Cell treatment, KD	Cleaved PARP Recombinant Rabbit mAb [SU0314] - BSA and Azide free	HA750136	Human	WB,IF-Cell,IP
Cell treatment	Cleaved PARP1 Recombinant Rabbit mAb [JE00-15]	HA722218	Human, Mouse, Rat, Monkey	WB
Cell treatment, KO	Cleaved+pro Caspase-3 Recombinant Rabbit mAb [SU38-04] - BSA and Azide free	HA750159	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment, KD	Cleaved+pro Caspase-3 Recombinant Rabbit mAb [JE75-05]	HA722367	Human, Mouse	WB
Cell treatment, KO	Cleaved+pro Caspase-3 Recombinant Rabbit mAb [SU38-04]	ET1608-64	Human, Mouse, Rat	WB,IHC-P,IP
	COX2 Mouse mAb [A3F8]	EM1902-11	Human, Rat	WB,IHC-P,FC
KD	COX2 Mouse mAb [A3F7]	EM1902-12	Human, Rat	WB,IHC-P,FC
	COX2 Recombinant Rabbit mAb [6J8]	RE6064	Human, Mouse	WB,ICC
	COX2 Mouse mAb	RT1159	Human, Mouse, Rat	WB,IP,IF,IHC-P,ELISA
KD	COX2/Cyclooxygenase 2 Recombinant Rabbit mAb [SC56-06] - BSA and Azide free	HA750207	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
KD	COX2/Cyclooxygenase 2 Recombinant Rabbit mAb [SC56-06]	ET1610-23	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	DDX3 Rabbit pAb	HA500093	Human, Mouse, Rat	WB,IF-Cell
	DDX3X Mouse mAb	EM1707-67	Human, Mouse, Monkey, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	DFNA5 / GSDME Mouse mAb [A1H2]	EM1901-49	Human	WB
Cell treatment, Relative expression	DFNA5 / GSDME Recombinant Rabbit mAb [PSH10-83]	HA723250	Human	WB
	DFNA5 / GSDME Recombinant Rabbit mAb [PSH07-48] - BSA and Azide free	HA751169	Human, Mouse, Rat	WB
	DFNA5 / GSDME Recombinant Rabbit mAb [PSH07-48]	HA722836	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	DFNA5 / GSDME Mouse mAb [A1H3]	EM1901-48	Human	WB,IF-Cell,IHC-P,FC
Cell treatment, Relative expression	DFNA5 / GSDME Recombinant Rabbit mAb [PSH09-51]	HA723113	Human, Mouse	WB,IP
Cell treatment, Relative expression	DFNA5 / GSDME Recombinant Rabbit mAb [PSH10-84]	HA723251	Human, Mouse	WB,IP
	DFNA5 / GSDME Rabbit pAb	ER1901-12	Human	WB,IHC-P,FC
Cell treatment, Relative expression	DFNA5 / GSDME Recombinant Rabbit mAb [PSH10-83] - BSA and Azide free	HA751375	Human	WB
Cell treatment, Relative expression	DFNA5 / GSDME Recombinant Rabbit mAb [PSH09-51] - BSA and Azide free	HA751298	Human, Mouse	WB,IP
Cell treatment, Relative expression	DFNA5 / GSDME Recombinant Rabbit mAb [PSH10-84] - BSA and Azide free	HA751376	Human, Mouse	WB,IP
	Dnmt1 Recombinant Rabbit mAb [JF09-89] - BSA and Azide free	HA750359	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Dnmt1 Recombinant Rabbit mAb [JF09-89]	ET1702-77	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	DNMT1 Mouse mAb [A2-C11]	EM1707-16	Human, Monkey	WB,IHC-P
	DNMT1 Recombinant Mouse mAb [A3A8-R]	HA601281	Human	WB,IF-Cell,IHC-P
	DNMT1 Mouse mAb [A3A9]	EM1901-84	Human	WB,IHC-P,FC
	DNMT1 Mouse mAb [A3A8]	EM1901-83	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	DNMT1 Recombinant Mouse mAb [A3A8-R] - BSA and Azide free	HA610156	Human	WB,IF-Cell,IHC-P
	FoxO1/3/4 Rabbit pAb	ER63465	Human, Mouse, Rat	WB,IHC-P
KD	FOXO1A Recombinant Rabbit mAb [SU33-01]	ET1608-25	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
KD	FOXO1A Recombinant Rabbit mAb [SU33-01] - BSA and Azide free	HA750140	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	iFluor™ 488 Conjugated FOXO1A Recombinant Rabbit mAb [SU33-01]	HA720161F	Human, Mouse, Rat	IF-Cell,FC
	FOXO3A Rabbit pAb	ER1706-79	Human, Mouse	WB,IF-Cell,IHC-P,FC
	FOXO3A Recombinant Rabbit mAb [SP0054]	ET1604-11	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue
	FOXO4 Recombinant Rabbit mAb [JJ09-11]	ET1701-67	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	FOXO4 Recombinant Rabbit mAb [SD0817]	ET1612-66	Human, Rat	WB,IHC-P
Relative expression, KD	GAP43 Recombinant Rabbit mAb [SC60-06]	ET1610-94	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Relative expression	GAP43 Rabbit pAb	ER40201	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	Gasdermin D (N terminal) Recombinant Rabbit mAb [PD00-18] - BSA and Azide free	HA750486	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Gasdermin D (N terminal) Rabbit pAb	ER1901-37	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
Cell treatment	Gasdermin D (N terminal) Recombinant Rabbit mAb [PD00-18]	HA721144	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Gasdermin D (N terminal) Recombinant Rabbit mAb [PSH10-87] - BSA and Azide free	HA751379	Human	WB,IP
Cell treatment	Gasdermin D (N terminal) Recombinant Rabbit mAb [PSH10-87]	HA723254	Human	WB,IP
Relative expression	Granzyme A Recombinant Rabbit mAb [PSH14-39]	HA723633	Human	WB
Relative expression	Granzyme A Recombinant Rabbit mAb [PSH14-39] - BSA and Azide free	HA751521	Human	WB
	Granzyme A Rabbit pAb	ER63626	Human	IHC-P,IF-Cell
	Granzyme B Recombinant Rabbit mAb [PSH14-20] - BSA and Azide free	HA751513	Mouse	WB
	Granzyme B Rabbit pAb	IRS032	Human	mIHC
	Granzyme B Rabbit pAb	HA500252	Human	IHC-P,mIHC
	Granzyme B Rabbit pAb	IRS054	Mouse	mIHC
	Granzyme B Recombinant Rabbit mAb [PSH14-20]	HA723609	Mouse	WB
Relative expression	Granzyme B Recombinant Rabbit mAb [PSH14-66] - BSA and Azide free	HA751527	Human, Mouse	WB,IHC-P
	Granzyme B Mouse mAb [4G9]	RT1259	Human, Mouse, Rat	WB,IP,IF,IHC-P
Relative expression	Granzyme B Recombinant Rabbit mAb [PSH14-66]	HA723664	Human, Mouse	WB,IHC-P
KD	GRIM19 Recombinant Rabbit mAb [JG82-33]	ET7108-66	Human	WB,IHC-P,IF-Cell
	GSDMD Mouse mAb [A8D7]	HA601046	Human, Mouse	WB
	GSDMD Mouse mAb [A8D8]	HA601047	Human	WB
KD	GSK3 alpha Rabbit pAb	R1312-1	Human	WB,IF-Cell,IHC-P,FC
	GSK3 alpha Recombinant Rabbit mAb [JE33-27]	HA722217	Human, Mouse, Rat, Monkey	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	GSK3 alpha/beta Rabbit pAb	ER63648	Human, Mouse, Rat	WB,IHC-P
KD	GSK3 beta Recombinant Rabbit mAb [SY28-03]	ET1607-71	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
KD	GSK3 beta Recombinant Rabbit mAb [SY28-03] - BSA and Azide free	HA750131	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
	GSK3 beta Rabbit pAb	R1510-4	Human, Mouse, Rat	IF-Cell,IHC-P,FC
KO	HDAC2 Recombinant Rabbit mAb [SY29-02] - BSA and Azide free	HA750134	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO	HDAC2 Recombinant Rabbit mAb [SY29-02]	ET1607-78	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO	HDAC2 Recombinant Rabbit mAb [SD0816]	ET1612-50	Human, Mouse, Rat	WB,IP,IHC-P,FC
KO	HDAC2 Rabbit pAb	R1601-7	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KO	HMGB1 Recombinant Rabbit mAb [SA39-03] - BSA and Azide free	HA750011	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
KO	HMGB1 Recombinant Rabbit mAb [SA39-03]	ET1601-2	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
KO	HMGB1 Rabbit pAb	R1309-12	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	HMGB1 Rabbit pAb	ER0913	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Human Caspase-9 Recombinant Rabbit mAb [PSH17-20] - BSA and Azide free (Detector)	HA723907	Human	ELISA(Det)
	Biotin Conjugated Human Caspase-9 Recombinant Rabbit mAb [PSH17-20] - Detector	HA723907B	Human	ELISA(Det),ELISA
	Human Caspase-9 Recombinant Rabbit mAb [PSH17-19] - BSA and Azide free (Capture)	HA723906	Human	ELISA(Cap)
	Biotin Conjugated Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-22] - Detector	HA725196B	Human	ELISA(Det)
	Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-22] - BSA and Azide free (Detector)	HA725196	Human	ELISA(Det)
	Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-21] - BSA and Azide free (Capture)	HA725195	Human	ELISA(Cap)
	Biotin Conjugated Human CCL5 Recombinant Rabbit mAb [PSH13-17] - Detector	HA725057B	Human	ELISA(Det),ELISA
	Human CCL5 Recombinant Rabbit mAb [PSH13-17] - BSA and Azide free (Detector)	HA725056	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HRP Conjugated Human CCL5 Recombinant Rabbit mAb [PSH13-17] - Detector	HA725057H	Human	ELISA(Det),ELISA
	Human CCL5 Recombinant Rabbit mAb [PSH13-16] - BSA and Azide free (Capture)	HA725055	Human	ELISA(Cap)
	Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-37] - BSA and Azide free (Detector)	HA725008	Human	ELISA(Det)
	Biotin Conjugated Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-37] - Detector	HA725009B	Human	ELISA(Det),ELISA
	Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-36] - BSA and Azide free (Capture)	HA725007	Human	ELISA(Cap)
	Human Granzyme A Recombinant Rabbit mAb [PSH09-34] - BSA and Azide free (Capture)	HA723092	Human	ELISA(Cap)
	Human Granzyme A Recombinant Rabbit mAb [PSH09-37] - BSA and Azide free (Detector)	HA723096	Human	ELISA(Det)
	Human Granzyme A Recombinant Rabbit mAb [PSH09-35] - BSA and Azide free (Detector)	HA723093	Human	ELISA(Det)
	Biotin Conjugated Human Granzyme A Recombinant Rabbit mAb [PSH09-37] - Detector	HA723097B	Human	ELISA(Det),ELISA
	Human Granzyme A Recombinant Rabbit mAb [PSH09-36] - BSA and Azide free (Capture)	HA723095	Human	ELISA(Cap)
	Biotin Conjugated Human Granzyme A Recombinant Rabbit mAb [PSH09-35] - Detector	HA723094B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human Granzyme B Recombinant Rabbit mAb [PSH15-99] - Detector	HA723804B	Human	ELISA(Det),ELISA
	Human Granzyme B Recombinant Rabbit mAb [PSH15-98] - BSA and Azide free (Capture)	HA723803	Human	ELISA(Cap)
	Human Granzyme B Recombinant Rabbit mAb [PSH15-99] - BSA and Azide free (Detector)	HA723804	Human	ELISA(Det)
	Human IL-18 Recombinant Rabbit mAb [PSH01-71] - BSA and Azide free (Capture)	HA721722	Human	ELISA(Cap)
Relative expression	Human IL-18 Recombinant Rabbit mAb [PSH02-09]	HA721768	Human	WB
Relative expression	Human IL-18 Recombinant Rabbit mAb [PSH02-09] - BSA and Azide free	HA750779	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human IL-18 Recombinant Rabbit mAb [PSH01-70]	HA722644B	Human	ELISA(Det)
	Human IL-18 Recombinant Rabbit mAb [PSH02-08] - BSA and Azide free	HA750778	Human	WB,IHC-P,IF-Tissue
	Human IL-18 Recombinant Rabbit mAb [PSH01-70] - BSA and Azide free (Detector)	HA721723	Human	ELISA(Det)
	Human IL-18 Recombinant Rabbit mAb [PSH02-08]	HA721767	Human	WB,IHC-P,IF-Tissue,IP
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA722649B	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA721713	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-87] - BSA and Azide free (Capture)	HA723883	Human	ELISA(Cap)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39]	HA721712	Human	WB
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free	HA750745	Human	WB
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - Detector	HA723882B	Human	ELISA(Det),ELISA
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-85] - BSA and Azide free (Capture)	HA723881	Human	ELISA(Cap)
	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free (Detector)	HA721279	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free	HA750744	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free (Capture)	HA721278	Human	ELISA(Cap)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - BSA and Azide free (Detector)	HA723882	Human	ELISA(Det)
	Human IL-6 Recombinant Rabbit mAb [PS00-15] - BSA and Azide free (Detector)	HA721168	Human	ELISA(Det)
	Human IL-6 Mouse mAb [A6D1]	HA601051	Human	WB,ELISA
	Human IL-6 Recombinant Rabbit mAb [PS00-21] - BSA and Azide free (Capture)	HA721167	Human	ELISA(Cap)
	Biotin Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA722233B	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HRP Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA721293	Human	ELISA(Det)
	Human IL-6 Mouse mAb [A6D2]	HA601050	Human	WB,ELISA
	Biotin Conjugated Human IL-8 Recombinant Rabbit mAb [PS01-22]	HA722611B	Human	ELISA(Det)
	Human IL-8 Recombinant Rabbit mAb [PS00-83] - BSA and Azide free (Capture)	HA722098	Human	ELISA(Cap)
	Human IL-8 Recombinant Rabbit mAb [PS01-22] - BSA and Azide free (Detector)	HA722099	Human	ELISA(Det)
	Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-68] - BSA and Azide free (Detector)	HA725110	Human	ELISA(Det)
	Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-67] - BSA and Azide free (Capture)	HA725109	Human	ELISA(Cap)
	HRP Conjugated Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-68] - Detector	HA725111H	Human	ELISA(Det),ELISA
	Biotin Conjugated Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-68] - Detector	HA725111B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human MMP-3 Recombinant Rabbit mAb [PSH07-62]	HA722879B	Human	ELISA(Det)
	Human MMP-3 Recombinant Rabbit mAb [PSH07-61] - BSA and Azide free (Capture)	HA722877	Human	ELISA(Cap)
	Human MMP-3 Recombinant Rabbit mAb [PSH07-62] - BSA and Azide free (Detector)	HA722878	Human	ELISA(Det)
	Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-42] - BSA and Azide free (Detector)	HA723331	Human, Mouse, Rat	ELISA(Det)
	Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-41] - BSA and Azide free (Capture)	HA723330	Human, Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-42] - Detector	HA723332B	Human, Mouse, Rat	ELISA(Det),ELISA
	Biotin Conjugated Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-44] - Detector	HA723335B	Human, Mouse, Rat	ELISA(Det),ELISA
	Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-43] - BSA and Azide free (Capture)	HA723333	Human, Mouse, Rat	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-44] - BSA and Azide free (Detector)	HA723334	Human, Mouse, Rat	ELISA(Det)
	Biotin Conjugated Human NLRP3 Recombinant Rabbit mAb [PSH11-08] - Detector	HA723277B	Human	ELISA(Det),ELISA
	Human NLRP3 Recombinant Rabbit mAb [PSH11-08] - BSA and Azide free (Detector)	HA723276	Human	ELISA(Det)
	Biotin Conjugated Human NLRP3 Recombinant Rabbit mAb [PSH11-09] - Detector	HA723279B	Human	ELISA(Det),ELISA
	Human NLRP3 Recombinant Rabbit mAb [PSH11-09] - BSA and Azide free (Detector)	HA723278	Human	ELISA(Det)
	Human NLRP3 Recombinant Rabbit mAb [PSH11-07] - BSA and Azide free (Capture)	HA723275	Human	ELISA(Cap)
	Human PAI1 Recombinant Rabbit mAb [PSH12-20] - BSA and Azide free (Capture)	HA723438	Human	ELISA(Cap)
	Biotin Conjugated Human PAI1 Recombinant Rabbit mAb [PSH12-19] - Detector	HA723437B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human PAI1 Recombinant Rabbit mAb [PSH12-18] - Detector	HA723439B	Human	ELISA(Det),ELISA
	Human PAI1 Recombinant Rabbit mAb [PSH12-19] - BSA and Azide free (Detector)	HA723436	Human	ELISA(Det)
	Human PAI1 Recombinant Rabbit mAb [PSH12-18] - BSA and Azide free (Capture/Detector)	HA723435	Human	ELISA(Cap),ELISA(Det)
	HRP Conjugated Human PAI1 Recombinant Rabbit mAb [PSH12-18] - Detector	HA723439H	Human	ELISA(Det),ELISA
	Biotin Conjugated Human PD-L1 Recombinant Rabbit mAb [PSH15-37] - Detector	HA723742B	Human	ELISA(Det),ELISA
	Human PD-L1 Recombinant Rabbit mAb [PSH15-37] - BSA and Azide free (Detector)	HA723741	Human	ELISA(Det)
	Human PD-L1 Recombinant Rabbit mAb [PSH15-36] - BSA and Azide free (Capture)	HA723740	Human	ELISA(Cap)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40]	HA721286	Human	IHC-P,mlHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HRP Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721294	Human	ELISA(Det)
	Biotin Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721559	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free (Detector)	HA721266	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free	HA750593	Human	IHC-P
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free (Capture)	HA721265	Human	ELISA(Cap)
	Human TNF alpha Recombinant Rabbit mAb	IRS040	Human	mIHC
	Human TRAIL Recombinant Rabbit mAb [PSH08-05] - BSA and Azide free (Capture)	HA722943	Human	ELISA(Cap)
	Human TRAIL Recombinant Rabbit mAb [PSH08-06] - BSA and Azide free (Detector)	HA722944	Human	ELISA(Det)
	Biotin Conjugated Human TRAIL Recombinant Rabbit mAb [PSH08-06]	HA722945B	Human	ELISA(Det)
	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-93] - Detector	HA723497B	Human	ELISA(Det),ELISA
	Human VEGFA Recombinant Rabbit mAb [PSH12-91] - BSA and Azide free (Capture)	HA723493	Human	ELISA(Cap)
	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-92] - Detector	HA723495B	Human	ELISA(Det),ELISA
	Human VEGFA Recombinant Rabbit mAb [PSH12-92] - BSA and Azide free (Detector)	HA723494	Human	ELISA(Det)
	Human VEGFA Recombinant Rabbit mAb [PSH12-93] - BSA and Azide free (Detector)	HA723496	Human	ELISA(Det)
	IKB alpha Recombinant Rabbit mAb [SZ00-07]	ET1603-6	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	IKB alpha Recombinant Rabbit mAb [SZ00-07] - BSA and Azide free	HA750068	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	IKBKE Recombinant Rabbit mAb [JU06-72]	ET1706-20	Human	WB,IF-Cell,FC,IHC-P
	IL-18 Recombinant Rabbit mAb [PS01-06] - BSA and Azide free	HA750677	Human	WB,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	IL-18 Recombinant Rabbit mAb [PS00-41] - BSA and Azide free	HA750676	Human	WB,IHC-P
	IL-18 Recombinant Rabbit mAb [PS00-41]	HA721536	Human	WB,IHC-P
Relative expression	IL-18 Recombinant Rabbit mAb [PSH18-15]	HA723978	Human, Mouse, Rat	WB,IHC-P,IP
	IL-18 Recombinant Rabbit mAb [PS01-06]	HA721537	Human	WB,FC
Cell treatment	IL-1 beta Mouse mAb [A7F9]	HA601002	Human, Mouse, Rat	WB,ELISA,IF-Cell
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R] - BSA and Azide free	HA610004	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Mouse mAb [A7A10]	HA601035	Human	WB,ELISA,IF-Cell
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R]	HA601036	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Recombinant Rabbit mAb [JJ087-3]	ET1701-39	Human	WB,IF-Cell,IF-Tissue,IHC-P
	IL-6 Mouse mAb [3-A4]	M1303-3	Mouse	WB,IF-Cell,FC,ELISA
	IL-6 Rabbit pAb	R1203-2	Mouse	WB,ELISA
	IL-6 Recombinant Rabbit mAb	IRS160RB	Human	mIHC
	IL-6 Mouse mAb	EM30301	Human	WB
	IL-6 Mouse mAb [1-6]	EM1701-45	Mouse, Rat	WB,IF-Cell
	IL-6 Rabbit pAb	R1412-2	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-96]	HA723263	Human	WB,FC
Cell treatment	IL-8 Rabbit pAb	R1511-15	Human	WB,IHC-P,IF-Cell,FC
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-97] - BSA and Azide free	HA751387	Human	WB
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-96] - BSA and Azide free	HA751386	Human	WB,FC
Cell treatment	IL-8 Rabbit pAb	ER1901-61	Human	WB,IHC-P,IF-Cell,FC
	IL-8 Rabbit pAb	ER1706-67	Human	ELISA,IF-Cell,IHC-P
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-97]	HA723264	Human	WB
Cell treatment	IP10 Recombinant Rabbit mAb [JA10-82] - BSA and Azide free	HA750411	Human	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	IP10 Recombinant Rabbit mAb [JA10-82]	ET1704-27	Human	WB,IF-Cell,IF-Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	JAK2 Mouse mAb [6-D3]	M1501-8	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	JAK2 Recombinant Rabbit mAb [SY0245]	ET1607-35	Human, Mouse, Rat	WB,IF-Cell
KD	JAK2 Recombinant Rabbit mAb [SY0245] - BSA and Azide free	HA750118	Human, Mouse, Rat	WB,IF-Cell
	JAK2 Rabbit pAb	ER1706-58	Human, Mouse, Rat	IF-Cell,IHC-P
	JAK2 Rabbit pAb	RT1343	Human, Mouse, Rat	WB,IP,IF,IHC-P
	KLF4 Rabbit pAb	R1308-1	Human, Mouse, Rat	WB,IHC-P,FC
	KLF4 Recombinant Rabbit mAb [JF98-08] - BSA and Azide free	HA750356	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	KLF4 Recombinant Rabbit mAb [JF98-08]	ET1702-71	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Lipocalin-2 Recombinant Rabbit mAb [JM10-67]	ET1703-39	Human	WB,IF-Tissue,IHC-P
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Chicken IgY (Chimeric) - BSA and Azide free	HA610301	Mouse, Rat	IHC-Fr,WB
Relative expression	MAP2 Recombinant Rabbit mAb [PSH08-73]	HA723025	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IHC-Fr,IF-Cell
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610292	Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell
	MAP2 Recombinant Rabbit mAb [PSH08-73]	IRS071	Mouse	mIHC
Relative expression	MAP2 Recombinant Rabbit mAb [PSH08-73] - BSA and Azide free	HA751238	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IHC-Fr,IF-Cell
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Guinea pig IgG2 (Chimeric)	HA601457	Mouse, Rat	IHC-Fr,WB
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Chicken IgY (Chimeric)	HA601467	Mouse, Rat	IHC-Fr,WB
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Rat IgG1 (Chimeric)	HA601459	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Mouse IgG1 (Chimeric)	HA601458	Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610291	Mouse, Rat	IHC-Fr,WB
	MAP2 Rabbit pAb	HA500177	Human, Mouse, Rat	WB,IF-Cell,IHC-P,mIHC
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610293	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell
KD	MCL1 Recombinant Rabbit mAb [SI16-04] - BSA and Azide free	HA750093	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
	MCL1 Rabbit pAb	R1510-26	Human, Mouse	IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	MCL1 Recombinant Rabbit mAb [SI16-04]	ET1606-14	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
KD	METTL14 Recombinant Rabbit mAb [PSH07-85]	HA722910	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP
KD	METTL14 Recombinant Rabbit mAb [PSH07-85] - BSA and Azide free	HA751192	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP
	METTL3 Recombinant Rabbit mAb [JE58-20]	HA720002	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	METTL3 Recombinant Rabbit mAb [JE58-20] - BSA and Azide free	HA750470	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	MK13 Rabbit pAb	ER60146	Human, Mouse, Rat	WB
	MMP-3 Rabbit pAb	ER1706-77	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	MMP-3 Recombinant Rabbit mAb [JM46-22]	ET1705-98	Human, Mouse, Rat	WB,IHC-P
	Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-45] - BSA and Azide free (Detector)	HA725041	Mouse	ELISA(Det)
	Biotin Conjugated Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-45] - Detector	HA725042B	Mouse	ELISA(Det),ELISA
	Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-44] - BSA and Azide free (Capture)	HA725040	Mouse	ELISA(Cap)
	Mouse CCL5 Recombinant Rat mAb	HA601395	Human, Mouse	ELISA
	Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-15] - BSA and Azide free (Capture)	HA725046	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-16] - Detector	HA725048B	Mouse	ELISA(Det),ELISA
	Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-16] - BSA and Azide free (Detector)	HA725047	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PS01-74] - Detector	HA723643B	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PS01-83] - BSA and Azide free (Detector)	HA721602	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PS01-83] - BSA and Azide free (Detector)	HA721601	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PSH0-59] - Detector	HA722646B	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IL-6 Recombinant Rabbit mAb [PSH0-60] - BSA and Azide free (Detector)	HA721413	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-58] - BSA and Azide free (Capture)	HA721411	Mouse	ELISA(Cap)
	Mouse IL-6 Recombinant Rabbit mAb [PSH14-47] - BSA and Azide free (Capture)	HA723642	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PSH0-60] - Detector	HA722647B	Mouse	ELISA(Det)
	HRP Conjugated Mouse IL-6 Mouse mAb [3-A4]	M1303-2	Mouse	WB,ELISA
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-59] - BSA and Azide free (Detector)	HA721412	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PS01-74] - BSA and Azide free (Capture/Detector)	HA721600	Mouse	ELISA(Cap),ELISA(Det)
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-29] - BSA and Azide free (Capture)	HA721974	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30]	HA722626B	Mouse, Rat	ELISA(Det)
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30] - BSA and Azide free (Detector)	HA721975	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77]	HA722623B	Mouse	ELISA(Det)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77] - BSA and Azide free (Detector)	HA722041	Mouse	ELISA(Det)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-76] - BSA and Azide free (Capture)	HA722040	Mouse	ELISA(Cap)
	Mouse VEGF Recombinant Rabbit mAb [PSH02-72] - BSA and Azide free (Capture)	HA721850	Mouse	ELISA(Cap)
	Mouse VEGF Recombinant Rabbit mAb [PSH02-73] - BSA and Azide free (Detector)	HA721851	Mouse	ELISA(Det)
	Biotin Conjugated Mouse VEGF Recombinant Rabbit mAb [PSH02-73]	HA722636B	Mouse	ELISA(Det)
	Mre11 Recombinant Rabbit mAb [JM11-18]	ET1703-99	Human, Mouse, Rat	WB,IHC-P
Relative expression	MTAP Recombinant Mouse mAb [PSH13-98]	HA601445	Human, Mouse, Rat	WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	MTAP Recombinant Rabbit mAb [JE63-74]	HA721955	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,FC
	MTAP Rabbit pAb	ER1901-47	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	mtTFA Recombinant Rabbit mAb [JM32-44]	ET1705-64	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Muc1 Rabbit pAb	ER1902-10	Human	WB,IF-Cell,IHC-P
	iFluor™ 647 Conjugated MUC1 Mouse mAb [PSH0-49]	HA600108F	Human	FC
Relative expression, KD	MUC1 Recombinant Rabbit mAb [SN06-80]	ET1611-14	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Relative expression	MUC1 Mouse mAb [PSH0-47]	HA601140	Human	IHC-P,WB
Relative expression, KD	MUC1 Recombinant Rabbit mAb [SN06-80] - BSA and Azide free	HA750245	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Relative expression	MUC1 Mouse mAb [A4D9]	EM1902-32	Human	IHC-P,FC,IF-Cell
Relative expression	MUC1 Mouse mAb [A9G4]	HA601139	Human	IHC-P,IF-Cell
Relative expression	MUC1 Mouse mAb [PSH0-48]	HA601141	Human	IHC-P,FC,IF-Cell
	iFluor™ 647 Conjugated MUC1 Mouse mAb [PSH0-48]	HA600107F	Human	FC
Relative expression	MUC1 Mouse mAb [PSH0-49]	HA601142	Human	IHC-P,FC,IF-Cell
	MUC1 Mouse mAb	IRS174MS	Human	mIHC
	NAL12 Rabbit pAb	ER65953	Human, Mouse	WB,IHC
	NeuN Recombinant Antibody [SR45-07] - Rat IgG1 (Chimeric)	HA601397	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
	NeuN Recombinant Rabbit mAb	IRS178RB	Mouse	mIHC
	NeuN Recombinant Mouse mAb [PD01-45] - BSA and Azide free	HA610031	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IHC-Fr
Relative expression	NeuN Recombinant Antibody [SR45-07] - Guinea pig IgG2 (Chimeric)	HA601481	Mouse, Rat	IHC-Fr,WB
Relative expression	NeuN Recombinant Antibody [SR45-07] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610315	Mouse, Rat	IHC-Fr,WB
	NeuN Recombinant Antibody [SR45-07] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610237	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
	NeuN Recombinant Antibody [SR45-07] - Chicken IgY (Chimeric) - BSA and Azide free	HA610314	Mouse	IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	NeuN Recombinant Rabbit mAb [SR45-07]	ET1602-12	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr,mIHC
	iFluor™ 594 Conjugated NeuN Recombinant Rabbit mAb [SR45-07]	HA720167F	Human, Mouse, Rat	IF-Tissue,IHC-Fr
Relative expression	NeuN Recombinant Antibody [SR45-07] - Mouse IgG1 (Chimeric)	HA601482	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	NeuN Recombinant Antibody [SR45-07] - Chicken IgY (Chimeric)	HA601480	Mouse	IHC-Fr
	NeuN Recombinant Antibody - Rat IgG1 (Chimeric)	IRS179RT	Mouse	mIHC
Relative expression	NeuN Recombinant Antibody [SR45-07] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610316	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	iFluor™ 488 Conjugated NeuN Recombinant Rabbit mAb [SR45-07]	HA720168F	Human, Mouse, Rat	IF-Tissue
	NeuN Recombinant Mouse mAb [PD01-45]	HA601111	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IHC-Fr
KD	NF-kB p105 / p50 Recombinant Rabbit mAb [SZ20-01]	ET1603-18	Human, Mouse, Rat	WB,IHC-P,ICC,FC,IF-Tissue
KD	NF-kB p105 / p50 Recombinant Rabbit mAb [SZ20-01] - BSA and Azide free	HA750060	Human, Mouse, Rat	WB,IHC-P,ICC,FC,IF-Tissue
	NF-kB p105 / p50 Rabbit pAb	R1309-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	NF-kB p105 / p50 Mouse mAb [A2-7]	M1505-8	Human, Mouse	WB,IF-Cell,IHC-P
Relative expression	NF-L Recombinant Rabbit mAb [PSH0-77]	HA721498	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
Relative expression	NF-L Recombinant Rabbit mAb [PS02-10] - BSA and Azide free	HA750678	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,IF-Tissue,IHC-Fr
Relative expression	NF-L Recombinant Rabbit mAb [PS02-10]	HA721538	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,IF-Tissue,IHC-Fr
Relative expression	NF-L Recombinant Antibody [PS02-10] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610282	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	NF-L Recombinant Antibody [PS02-10] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610297	Mouse, Rat	IHC-Fr
Relative expression	NF-L Recombinant Antibody [PS02-10] - Chicken IgY (Chimeric) - BSA and Azide free	HA610306	Mouse, Rat	IHC-Fr,WB
Relative expression	NF-L Recombinant Antibody [PS02-10] - Chicken IgY (Chimeric)	HA601472	Mouse, Rat	IHC-Fr,WB
Relative expression	NF-L Recombinant Antibody [PS02-10] - Mouse IgG1 (Chimeric)	HA601448	Human, Mouse, Rat	IHC-Fr,IHC-P,WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	NF-L Recombinant Antibody [PS02-10] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610249	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
Relative expression	NF-L Recombinant Rabbit mAb [PSH0-77] - BSA and Azide free	HA750657	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
Relative expression	NF-L Recombinant Antibody [PS02-10] - Rat IgG1 (Chimeric)	HA601409	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
	NF-L Recombinant Antibody [PS02-10] - Guinea pig IgG2 (Chimeric)	HA601463	Mouse, Rat	IHC-Fr
	NF-L Rabbit pAb	ER65439	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	NLRP14 Rabbit pAb	ER1913-88	Human, Mouse, Rat	WB
Relative expression	NLRP3 Recombinant Rabbit mAb [PSH19-43]	HA724075	Human, Mouse	WB
Cell treatment	NLRP3 Recombinant Rabbit mAb [SC06-23]	ET1610-93	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	NLRP3 Recombinant Rabbit mAb [SC06-23] - BSA and Azide free	HA750236	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	NLRP3 Rabbit pAb	ER1706-72	Human, Mouse	WB,IF-Cell,IHC-P,FC
	NLRP9 Rabbit pAb	ER1913-89	Mouse	WB
	NLRX1 Rabbit pAb	ER65729	Human, Mouse, Rat	WB,IHC-P
	Nrf2 Rabbit pAb	R1312-8	Human, Mouse	WB,IHC-P,IF-Tissue
Cell treatment	Nrf2 Recombinant Rabbit mAb [PSH11-20]	HA723302	Human, Mouse, Rat	WB,IF-Cell,ChIP
Cell treatment	Nrf2 Rabbit pAb	ER1706-41	Human, Mouse, Rat	IHC-P,IF-Cell,WB
Cell treatment	Nrf2 Recombinant Rabbit mAb [JE30-14]	HA721432	Human, Mouse	WB,IF-Cell
Cell treatment	Nrf2 Recombinant Rabbit mAb [PSH11-20] - BSA and Azide free	HA751401	Human, Mouse, Rat	WB,IF-Cell,ChIP
KD, Cell treatment	p38 alpha / MAPK14 Recombinant Rabbit mAb [SR43-04] - BSA and Azide free	HA750040	Human, Mouse, Rat	WB,IF-Cell,FC
KO	p38 alpha / MAPK14 Recombinant Rabbit mAb [JF55-07]	ET1702-65	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	p38 alpha / MAPK14 Rabbit pAb	R1308-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	p38 alpha / MAPK14 Mouse mAb [9-B9]	EM40307	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD, Cell treatment	p38 alpha / MAPK14 Recombinant Rabbit mAb [SR43-04]	ET1602-26	Human, Mouse, Rat	WB,IF-Cell,FC
	p38 alpha/MAPK14 + p38 beta/MAPK11 + p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-60] - BSA and Azide free	HA750894	Human, Mouse, Rat, Monkey	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	p38 alpha/MAPK14 + p38 beta/MAPK11 + p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-60]	HA722019	Human, Mouse, Rat, Monkey	WB
	p38 beta/MAPK11 + p38 gamma/MAPK12 + p38 delta/MAPK13 Recombinant Rabbit mAb [PSH03-59]	HA722018	Human, Mouse, Rat, Monkey	WB
	p38 beta/MAPK11 + p38 gamma/MAPK12 + p38 delta/MAPK13 Recombinant Rabbit mAb [PSH03-59] - BSA and Azide free	HA750893	Human, Mouse, Rat, Monkey	WB
	p38 delta / MAPK13 Rabbit pAb	HA500582	Human	WB
	p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-58] - BSA and Azide free	HA750892	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-58]	HA722017	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	p38 MAPK Recombinant Rabbit Multiclonal Antibody [PSH03-72]	HA722036	Human, Mouse, Rat, Monkey	WB
	p38 MAPK Recombinant Rabbit Multiclonal Antibody [PSH03-72] - BSA and Azide free	HA750906	Human, Mouse, Rat, Monkey	WB
Relative expression	p40 Recombinant Rabbit mAb [PDH0-02]	HA721222	Human, Mouse, Rat	IHC-P,WB,IF-Cell
	p40 Recombinant Rabbit mAb [PDH0-06]	HA721255	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment, KD, Relative expression	p53 Recombinant Mouse mAb [6B5-R]	HA601315	Human, Mouse, Monkey	WB,IF-Cell
KD	p53 Rabbit pAb	R1308-7	Human	WB,IF-Cell
KD	p53 Recombinant Rabbit mAb [JB42-26]	ET7107-33	Human	WB,IHC-P,FC,IP,IF-Cell,IF-Tissue
Cell treatment	p53 Recombinant Rabbit mAb [PSH03-95]	HA722074	Mouse	WB,IF-Cell,IP
KD	p53 Recombinant Mouse mAb [6B2-R]	HA601322	Human, Monkey	WB,IF-Cell
KD	p53 Recombinant Mouse mAb [6B2-R] - BSA and Azide free	HA610194	Human, Monkey	WB,IF-Cell
Relative expression, KD, Cell treatment	p53 Recombinant Rabbit mAb [SA39-07] - BSA and Azide free	HA750005	Human	WB,IF-Cell,ChIP
	p53 Mouse mAb [C8-A11]	M1312-2	Human	WB,IF-Cell,IHC-P
Cell treatment	p53 Recombinant Rabbit mAb [PSH03-95] - BSA and Azide free	HA750920	Mouse	WB,IF-Cell,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment, Relative expression, KD	p53 Recombinant Mouse mAb [6B5-R] - BSA and Azide free	HA610187	Human, Mouse, Monkey	WB,IF-Cell
	p53 Recombinant Rabbit mAb [PD01-32] - BSA and Azide free	HA750545	Human	IHC-P
Relative expression	p53 Mouse mAb [6B2]	EM1701-91	Human, Mouse	WB,IF-Cell
	p53 Mouse mAb [6B5]	EM1701-92	Human	WB
KO	p53 Mouse mAb [7-5]	EM20603	Human	WB,IF-Cell,IHC-P,FC
	p53 Recombinant Rabbit mAb [PD01-32]	HA721213	Human	IHC-P
Relative expression, KD	p53 Recombinant Rabbit mAb [SY010-6]	ET1605-16	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Relative expression, KD, Cell treatment	p53 Recombinant Rabbit mAb [SA39-07]	ET1601-13	Human	WB,IF-Cell,ChIP
Cell treatment	p53 (acetyl K370) Recombinant Rabbit mAb [SR40-09]	ET1602-38	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC
Cell treatment	p53 (acetyl K370) Recombinant Rabbit mAb [SR40-09] - BSA and Azide free	HA750047	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [PSH06-87] - BSA and Azide free	HA751111	Human, Mouse	WB,IF-Cell
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [PSH06-88]	HA722763	Human	WB,IF-Cell
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [PSH06-87]	HA722762	Human, Mouse	WB,IF-Cell
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [PSH06-88] - BSA and Azide free	HA751112	Human	WB,IF-Cell
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [SD0801]	ET1612-71	Human, Mouse	WB,IF-Cell,IHC-P
	p63 Rabbit pAb	R1311-3	Human, Mouse, Rat	WB,IHC-P
Relative expression	p63 Recombinant Rabbit mAb [JE53-53]	ET7110-47	Human, Mouse	WB,IF-Cell,IHC-P
Relative expression	p63 Recombinant Rabbit mAb [SC06-31] - BSA and Azide free	HA750218	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	p63 Recombinant Rabbit mAb [SC06-31]	ET1610-44	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	PAI1 Rabbit pAb	HA500124	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	PAI1 Mouse mAb [A9G8]	HA601152	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment, KD	PARP1 Recombinant Rabbit mAb [SU03-68]	ET1608-56	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	PARP1 Rabbit pAb	ER1802-67	Human, Mouse	WB,IF-Cell,FC
Cell treatment, KD	PARP1 Recombinant Rabbit mAb [SU03-68] - BSA and Azide free	HA750154	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	PDIA3 / ERp57 Recombinant Rabbit mAb [PSH08-79]	HA723031	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
KD	PDIA3 / ERp57 Recombinant Rabbit mAb [JM11-44]	ET1703-68	Human	WB,IP,IHC-P
KD	PDIA3 / ERp57 Recombinant Rabbit mAb [PSH08-79] - BSA and Azide free	HA751244	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	PD-L1 Recombinant Rabbit Multiclonal Antibody [PSH13-23] - BSA and Azide free	HA751459	Human, Mouse	WB
Cell treatment, Relative expression, KD	PD-L1 Recombinant Rabbit mAb [PD01-02]	HA721176	Human	IHC-P,mIHC,WB
	PD-L1 Recombinant Rabbit mAb [JJ08-95]	ET1701-41	Human	WB,IHC-P
	PD-L1 Recombinant Rabbit mAb [PD00-97]	HA722300	Human	IHC-P
	PD-L1 Recombinant Rabbit mAb	IRS009	Human	mIHC
	PD-L1 Recombinant Rabbit mAb [JJ08-95] - BSA and Azide free	HA750309	Human	WB,IHC-P
	PD-L1 Recombinant Rabbit mAb [PD00-97] - BSA and Azide free	HA751015	Human	IHC-P
	PD-L1 Recombinant Rabbit Multiclonal Antibody [PSH13-23]	HA723510	Human, Mouse	WB
KO	PFKFB3 Recombinant Rabbit mAb [JM43-43]	ET1705-66	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP,IF-Tissue
Relative expression	PGP9.5 Recombinant Antibody [JM10-59] - Rat IgG1 (Chimeric)	HA601515	Human, Mouse, Rat	IHC-Fr,IF-Tissue,IHC-P,WB,IF-Cell
KD	PGP9.5 Recombinant Rabbit mAb [JF93-08]	ET1702-83	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Relative expression, KD	PGP9.5 Recombinant Rabbit mAb [JM10-59]	ET1703-22	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
Relative expression	PGP9.5 Recombinant Antibody [JM10-59] - Mouse IgG1 (Chimeric)	HA601514	Human, Mouse, Rat	IHC-Fr,IF-Tissue,IHC-P,WB,IF-Cell
Relative expression, KD	PGP9.5 Mouse mAb [B1-5-6]	EM1701-86	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression, KD	PGP9.5 Recombinant Rabbit mAb [JM10-59] - BSA and Azide free	HA750375	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
Relative expression	PGP9.5 Recombinant Antibody [JM10-59] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610344	Human, Mouse, Rat	IHC-Fr,IF-Tissue,IHC-P,WB,IF-Cell
Relative expression	PGP9.5 Recombinant Rabbit mAb [PSH13-92]	HA723585	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue
Relative expression	PGP9.5 Recombinant Antibody [JM10-59] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610345	Human, Mouse, Rat	IHC-Fr,IF-Tissue,IHC-P,WB,IF-Cell
KD	PGP9.5 Rabbit pAb	R1511-8	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	Phospho-AKT1 (S124) Recombinant Rabbit mAb [JJ08-46]	ET1701-36	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-AKT1 (S129) Recombinant Rabbit mAb [JE59-50]	HA721677	Human, Mouse, Rat	WB
Cell treatment	Phospho-AKT1 (T450) Recombinant Rabbit mAb [SD08-12]	ET1612-73	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell,FC
Cell treatment	Phospho-AKT2 (S474) Recombinant Rabbit mAb [PSH09-46]	HA723108	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-AKT2 (S474) Rabbit pAb	HA500116	Human, Mouse	WB,IHC-P
Cell treatment	Phospho-AKT2 (S474) Recombinant Rabbit mAb [PSH09-46] - BSA and Azide free	HA751293	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05]	ET1607-73	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44]	HA722129	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05] - BSA and Azide free	HA750132	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44] - BSA and Azide free	HA750953	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10]	HA722951	Human	WB,IF-Cell,FC
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10] - BSA and Azide free	HA751205	Human	WB,IF-Cell,FC
Cell treatment	Phospho-FOXO1A (T24)+ FOXO3A (T32)+ FOXO4 (T28)+ FOXO6 (T26) Recombinant Rabbit mAb [PSH06-00] - BSA and Azide free	HA751040	Human, Mouse	WB,IF-Cell,FC
Cell treatment	Phospho-FOXO1A (T24)+ FOXO3A (T32)+ FOXO4 (T28)+ FOXO6 (T26) Recombinant Rabbit mAb [PSH06-00]	HA722499	Human, Mouse	WB,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-FOXO3a (S253) Recombinant Rabbit mAb [ST49-01] - BSA and Azide free	HA750178	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-FOXO3a (S253) Recombinant Rabbit mAb [ST49-01]	ET1609-49	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-GSK3 (alpha + beta) (Y216 + Y279) Recombinant Rabbit mAb [SY26-05]	ET1607-54	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,IF-Tissue,FC
Cell treatment	Phospho-GSK3 (alpha + beta) (Y216 + Y279) Recombinant Rabbit mAb [SY26-05] - BSA and Azide free	HA750127	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,IF-Tissue,FC
Cell treatment	Phospho-GSK3 alpha (S21) Recombinant Rabbit mAb [PSH09-11] - BSA and Azide free	HA751265	Human, Mouse	WB,IHC-P,Dot Blot
Cell treatment	Phospho-GSK3 alpha (S21) Recombinant Rabbit mAb [PSH09-11]	HA723069	Human, Mouse	WB,IHC-P,Dot Blot
Cell treatment	Phospho-GSK3 beta (S9) Recombinant Rabbit mAb [SY02-71]	ET1607-60	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-GSK3 beta (S9) Recombinant Rabbit mAb [SY02-71] - BSA and Azide free	HA750128	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-IKB alpha (S32) Recombinant Rabbit mAb [ST53-05] - BSA and Azide free	HA750196	Human, Mouse	WB,IF-Cell,IP
Cell treatment	Phospho-IKB alpha (S32) Recombinant Rabbit mAb [ST53-05]	ET1609-78	Human, Mouse	WB,IF-Cell,IP
Cell treatment	Phospho-IKB alpha (S32/S36) Recombinant Rabbit mAb [PSH06-76]	HA722770	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Phospho-IKB alpha (S32/S36) Recombinant Rabbit mAb [PSH06-76] - BSA and Azide free	HA751119	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Phospho-IKB alpha (S36) Recombinant Rabbit mAb [JE56-65]	HA721802	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-JAK2 (Y1007 + Y1008) Recombinant Rabbit mAb [SY24-03] - BSA and Azide free	HA750117	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-JAK2 (Y1007 + Y1008) Recombinant Rabbit mAb [SY24-03]	ET1607-34	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-MAP2 (S136) Recombinant Rabbit mAb [JE44-25]	HA722351	Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Tissue
	Phospho-NLRC4 (S533) Rat mAb [B7-B7]	EM1707-24	Human, Mouse	WB,IF-Cell
Cell treatment, Relative expression	Phospho-Nrf2 (S40) Recombinant Rabbit mAb [SU0334] - BSA and Azide free	HA750142	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment, Relative expression	Phospho-Nrf2 (S40) Recombinant Rabbit mAb [SU0334]	ET1608-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Cell treatment	Phospho-p38 alpha (T180) Recombinant Rabbit mAb [PSH06-25] - BSA and Azide free	HA751062	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180) Recombinant Rabbit mAb [PSH06-25]	HA722657	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180 + Y182) Recombinant Rabbit mAb [PSH06-35]	HA722669	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180 + Y182) Recombinant Rabbit mAb [PSH06-35] - BSA and Azide free	HA751074	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-p38 (T180) Rabbit pAb	ER2001-51	Human, Mouse, Rat	WB
Cell treatment	Phospho-p38 (T180) Recombinant Rabbit mAb [JE77-37]	HA722150	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-p38 (T180 + Y182) Rabbit pAb	ER2001-52	Human, Mouse, Rat	WB
	Phospho-p38 (T180 + Y182) Rabbit pAb	ER1903-01	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
Cell treatment	Phospho-p53 (S15) Recombinant Rabbit mAb [PSH01-98] - BSA and Azide free	HA750772	Human	WB,IF-Cell
Cell treatment	Phospho-p53 (S15) Recombinant Rabbit mAb [PSH01-98]	HA721756	Human	WB,IF-Cell
Cell treatment	Phospho-p53 (S20) Recombinant Rabbit mAb [PSH06-86] - BSA and Azide free	HA751110	Human	WB
Cell treatment	Phospho-p53 (S20) Recombinant Rabbit mAb [PSH06-86]	HA722761	Human	WB
Cell treatment	Phospho-p53 (S33) Recombinant Rabbit mAb [JE41-97]	HA722703	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-p53 (S376) Recombinant Rabbit mAb [ST0440]	ET1609-14	Human, Mouse	WB
Cell treatment	Phospho-p53 (S392) Recombinant Rabbit mAb [SI17-04] - BSA and Azide free	HA750097	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-p53 (S392) Recombinant Rabbit mAb [SI17-04]	ET1606-24	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-p53 (T55) Recombinant Rabbit mAb [ST0439]	ET1609-13	Human	WB,IHC-P
	Phospho-PTEN (S380) Recombinant Rabbit mAb [JJ08-81]	ET1701-46	Human, Mouse, Rat	WB,IF-Cell
	Phospho-SHIP (Y1020) Recombinant Rabbit mAb [JJ082-03]	ET1701-5	Human	WB,IF-Cell,IP
Cell treatment	Phospho-SIRT1 (S47) Recombinant Rabbit mAb [SN06-40]	ET1611-31	Human	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Phospho-SIRT1 (T530) Recombinant Rabbit mAb [JJ206-6]	ET1701-27	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-Smad2 (S250) Recombinant Rabbit mAb [SD207-1]	ET1612-32	Human, Mouse, Rat	WB,FC,IP,IHC-P
Cell treatment	Phospho-Smad2 (S250) Recombinant Rabbit mAb [SD207-1] - BSA and Azide free	HA750291	Human, Mouse, Rat	WB,FC,IP,IHC-P
Cell treatment	Phospho-Smad2 (S255) Recombinant Rabbit mAb [JF0882]	ET1702-34	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-SMAD2 (S465+467) Recombinant Rabbit mAb [PSH05-76] - BSA and Azide free	HA751026	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-SMAD2 (S465+467) Recombinant Rabbit mAb [PSH05-76]	HA722443	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-Smad3 (S423 + S425) Recombinant Rabbit mAb [ST0493] - BSA and Azide free	HA750172	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-Smad3 (S423 + S425) Recombinant Rabbit mAb [ST0493]	ET1609-41	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment, Relative expression	Phospho-STAT1 (S727) Recombinant Rabbit mAb [SN67-04] - BSA and Azide free	HA750247	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment, Relative expression	Phospho-STAT1 (S727) Recombinant Rabbit mAb [SN67-04]	ET1611-20	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-STAT1 (Y701) Recombinant Rabbit mAb [PSH04-02]	HA722083	Human	WB,IF-Cell
Cell treatment	Phospho-STAT1 (Y701) Recombinant Rabbit mAb [PSH04-02] - BSA and Azide free	HA750925	Human	WB,IF-Cell
Cell treatment, KD	Phospho-STAT3 (S727) Recombinant Rabbit mAb [SY24-09] - BSA and Azide free	HA750121	Human, Rat, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment, KD	Phospho-STAT3 (S727) Recombinant Rabbit mAb [SY24-09]	ET1607-39	Human, Rat, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-STAT3 (Y705) Recombinant Rabbit mAb [SZ43-01]	ET1603-40	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-STAT3 (Y705) Recombinant Rabbit mAb [SZ43-01] - BSA and Azide free	HA750067	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-STAT5 (Y694) Recombinant Rabbit mAb [SC05-31]	ET1610-48	Human, Rat, Mouse	WB,IHC-P,IF-Cell
Cell treatment	Phospho-STAT5 (Y694) Recombinant Rabbit mAb [SC05-31] - BSA and Azide free	HA750220	Human, Rat, Mouse	WB,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-STING (S366) Recombinant Rabbit mAb [PSH09-72]	HA723137	Human	WB, Dot Blot, IP
Cell treatment	Phospho-STING (S366) Recombinant Rabbit mAb [PSH09-72] - BSA and Azide free	HA751314	Human	WB, Dot Blot, IP
Cell treatment	Phospho-Vimentin (S39) Recombinant Rabbit mAb [JE43-26]	HA721442	Human, Mouse, Rat	WB, IF-Cell
Cell treatment	PINK1 Recombinant Rabbit mAb [PSH08-69] - BSA and Azide free	HA751234	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB, IHC-P, IHC-Fr
Cell treatment	PINK1 Recombinant Rabbit mAb [PSH08-69]	HA723021	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB, IHC-P, IHC-Fr
Cell treatment	PINK1 Rabbit pAb	ER1706-27	Human	WB
Relative expression	PKM1 Recombinant Rabbit mAb [PSH04-37] - BSA and Azide free	HA750946	Human, Mouse, Rat	WB, IF-Cell, IHC-P, FC
Relative expression	PKM1 Recombinant Rabbit mAb [PSH04-37]	HA722120	Human, Mouse, Rat	WB, IF-Cell, IHC-P, FC
	PKM1/2 Recombinant Rabbit mAb [JE07-34]	HA722524	Human, Mouse, Rat	WB, IP
Relative expression	PKM2 Recombinant Rabbit mAb [PSH11-74] - BSA and Azide free	HA751429	Human, Mouse, Rat, Monkey	WB, IF-Cell, IHC-P, FC, IP
Relative expression	PKM2 Rabbit pAb	R1603-5	Human, Mouse, Rat	WB, IF-Cell, IHC-P, FC
KO	PKM2 Rabbit pAb	ER1802-70	Human	WB, IF-Cell, IHC-P
Relative expression	PKM2 Recombinant Rabbit mAb [PSH11-74]	HA723370	Human, Mouse, Rat, Monkey	WB, IF-Cell, IHC-P, FC, IP
KO	PKM2 Rabbit pAb	ER1901-90	Human, Mouse	WB, FC
Cell treatment, Relative expression, KO	Pro Caspase-3 Recombinant Rabbit mAb [SR03-01] - BSA and Azide free	HA750048	Human, Monkey	WB, IF-Cell, IF-Tissue, IHC-P, IP, FC
Cell treatment, Relative expression, KO	Pro Caspase-3 Recombinant Rabbit mAb [SR03-01]	ET1602-39	Human, Monkey	WB, IF-Cell, IF-Tissue, IHC-P, IP, FC
KD	Pro Caspase-3 Recombinant Rabbit mAb [SZ02-08]	ET1603-26	Human, Monkey, Mouse	WB, IF-Tissue, IHC-P
	Pro Caspase-3 Rabbit pAb	ER30804	Human, Mouse	WB, IF-Cell, IHC-P, FC, IHC-Fr
	pro Caspase-7 Recombinant Rabbit mAb [SD20-57]	ET1612-28	Human	WB, IF-Cell, IF-Tissue, IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PSD95 Recombinant Rabbit mAb [SR38-09] - BSA and Azide free	HA750035	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,IP,IHC-Fr
Relative expression	PSD95 Recombinant Antibody [SR38-09] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610256	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,WB
	PSD95 Recombinant Rabbit mAb [SR38-09]	ET1602-20	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,IHC-Fr
	PSD95 Rabbit pAb	R1508-4	Human, Mouse, Rat	WB,IHC-P,FC
Relative expression	PSD95 Recombinant Antibody [SR38-09] - Rat IgG1 (Chimeric)	HA601416	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,WB
	PSD95 Mouse mAb [B5-E6]	M1511-4	Human, Mouse, Rat	WB,IF-Cell,FC
Relative expression, KD	PTEN Recombinant Rabbit mAb [JE51-96]	HA721962	Human, Mouse, Rat	WB,IHC-P,FC
Relative expression, KD	PTEN Recombinant Rabbit mAb [SJ19-03] - BSA and Azide free	HA750105	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	PTEN Mouse mAb [4D8]	RT1519	Human, Mouse, Rat	WB,IP,IF
	PTEN Rabbit pAb	ER1902-27	Human, Mouse, Rat	WB,IHC-P,FC
KD	PTEN Rabbit pAb	ER1902-33	Human, Mouse, Rat	WB,IHC-P,FC
Relative expression, KD	PTEN Recombinant Rabbit mAb [SJ19-03]	ET1606-43	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	PTEN Rabbit pAb	ER1706-51	Human, Mouse	WB,IF-Cell,IHC-P,FC
Cell treatment	Pyrin Recombinant Rabbit mAb [PSH06-12] - BSA and Azide free	HA751050	Human	WB,FC
Cell treatment	Pyrin Recombinant Rabbit mAb [PSH06-12]	HA722509	Human	WB,FC
	RANTES Recombinant Rabbit mAb [JM03-45]	ET1705-70	Human, Mouse	WB,IHC-P,IF-Cell,IF-Tissue
Cell treatment	S100A4 Recombinant Rabbit mAb [SD200-08]	ET1612-13	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cell treatment	S100A4 Recombinant Rabbit mAb [SD200-08] - BSA and Azide free	HA750282	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO	SDHB Rabbit pAb	ER1803-63	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	SDHB Recombinant Rabbit mAb [JE77-70]	HA722901	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell
KO	SDHB Recombinant Rabbit mAb [JU02-33]	ET1706-30	Human, Mouse, Rat	WB,IHC-P
	SERPINE1 Mouse mAb [H9-D5]	EM1709-36	Human, Mouse, Rat, Rabbit	IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	SET7 Recombinant Rabbit mAb [PSH02-90]	HA721909	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
KD	SET7 Recombinant Rabbit mAb [PSH02-90] - BSA and Azide free	HA750845	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
	SETD7 Mouse mAb [1B3B8]	EM1712-11	Human, Rat, Monkey	WB,FC
	SGK1 Recombinant Rabbit mAb [SC05-71]	ET1610-19	Human, Mouse, Rat	WB,IHC-P,IP,FC,IF-Tissue
	SHIP Recombinant Rabbit mAb [SY11-08]	ET1605-30	Human	WB,IHC-P,IP,FC
KD	SIRT1 Recombinant Rabbit mAb [SZ04-01]	ET1603-3	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP
KD	SIRT1 Rabbit pAb	ER130811	Human, Mouse	WB,IHC-P,FC,IF-Cell
KD	SIRT1 Recombinant Rabbit mAb [JE11-04]	ET7110-49	Human	WB,IF-Cell
KD	SIRT1 Recombinant Rabbit mAb [SZ04-01] - BSA and Azide free	HA750064	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP
KD	SIRT1 Mouse mAb [7-C5-B2]	M1506-3	Human, Mouse	WB,IF-Cell,IHC-P
KD	SIRT3 Recombinant Rabbit mAb [PSH05-22] - BSA and Azide free	HA750993	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue,IP,ChIP
KD	SIRT3 Recombinant Rabbit mAb [PSH05-22]	HA722251	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue,IP,ChIP
	SIRT3 Rabbit pAb	R1511-3	Human, Mouse, Rat	WB
KO	Smad2 Recombinant Rabbit mAb [SP06-05]	ET1604-22	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	Smad2 Mouse mAb [C9-B0]	EM1701-44	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Smad2 Rabbit pAb	R1511-29	Human, Mouse, Rat	IF-Cell,IHC-P
KO	Smad2 Recombinant Rabbit mAb [SP06-05] - BSA and Azide free	HA750073	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	Smad2/3 Rabbit pAb	ER64981	Human, Mouse, Rat, Monkey	WB,IHC-P
	Smad2/3 Mouse mAb [4D3]	RT1566	Human, Mouse, Rat	WB,IP,IF,IHC-P
	Smad2/3 (Acetyl K19) Rabbit pAb	ER60061	Human, Mouse, Rat	IHC-P,WB
	Smad3 Rabbit pAb	R1510-24	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
KO	Smad3 Recombinant Rabbit mAb [SY25-01] - BSA and Azide free	HA750123	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Smad3 Rabbit pAb	R1510-25	Human	WB,IF-Cell,IHC-P,FC
KO	Smad3 Recombinant Rabbit mAb [SY25-01]	ET1607-41	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Smad3 Mouse mAb [8E1]	RT1567	Human, Mouse, Rat	WB,IP,IF,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	STAT1 Recombinant Rabbit mAb [SD20-75]	ET1612-22	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	STAT1 Recombinant Mouse mAb [G3-B11-R] - BSA and Azide free	HA610181	Human	WB,IHC-P
	STAT1 Recombinant Mouse mAb [G3-B11-R]	HA601309	Human	WB,IHC-P
	STAT1 Mouse mAb [G3-B11]	M1407-1	Human	WB,IF-Cell,IHC-P,IF-Tissue
	STAT1 Rabbit pAb	R1408-2	Human, Mouse, Rat, Zebrafish	WB,IF-Cell
KD	STAT1 alpha Recombinant Rabbit mAb [SJ01-89] - BSA and Azide free	HA750102	Human, Mouse	WB,IP,FC,IHC-P,IF-Cell,IF-Tissue
KD	STAT1 alpha Recombinant Rabbit mAb [SJ01-89]	ET1606-39	Human, Mouse	WB,IP,FC,IHC-P,IF-Cell,IF-Tissue
	STAT1 alpha/beta Rabbit pAb	ER31221	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	STAT3 Recombinant Rabbit mAb [SY34-01] - BSA and Azide free	HA750087	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	STAT3 Recombinant Rabbit mAb [SY24-08]	ET1607-38	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
KD	STAT3 Recombinant Rabbit mAb [SY24-08] - BSA and Azide free	HA750120	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	STAT3 Rabbit pAb	RT1591	Human, Mouse, Rat	WB,IP,IF,IHC-P
KD	STAT3 Recombinant Rabbit mAb [SY34-01]	ET1605-45	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	STAT5a Recombinant Rabbit mAb [SC66-09]	ET1610-58	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	STAT5a Recombinant Rabbit mAb [SC66-09] - BSA and Azide free	HA750225	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	STAT 5A+B Recombinant Rabbit mAb [JJ08-78]	ET1701-45	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	STAT 5A+B Recombinant Rabbit mAb [JJ08-78] - BSA and Azide free	HA750310	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	STAT5b Recombinant Rabbit mAb [SD08-08]	ET1612-63	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	STING Recombinant Rabbit mAb [JM03-47]	ET1705-68	Human, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue,IP
	STING Recombinant Rabbit mAb [JM03-47] - BSA and Azide free	HA750437	Human, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue,IP
	STING Recombinant Rabbit mAb [PSH07-44] - BSA and Azide free	HA751165	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,IP
	STING Recombinant Rabbit mAb [PSH07-44]	HA722832	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,IP
	Synaptophysin Mouse mAb [A1-7]	M0407-30	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Synaptophysin Rabbit pAb	0407-2	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	Synaptophysin Recombinant Rabbit mAb [SJ26-85]	ET1606-56	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	Synaptophysin Recombinant Rabbit mAb [JE42-51]	HA721599	Human, Mouse, Rat	WB,IHC
KO	TLR2 Recombinant Rabbit mAb [JM22-41]	ET1705-92	Human, Mouse, Rat	WB,FC,IF-Cell
KO	TLR2 Recombinant Rabbit mAb [JM22-41] - BSA and Azide free	HA750447	Human, Mouse, Rat	WB,FC,IF-Cell
	TLR4 Rabbit pAb	ER1706-43	Human, Mouse, Rat	IF-Cell,IHC-P,FC
	TLR4 Mouse mAb	RT1666	Human, Mouse, Rat	WB
	TMS1 Recombinant Rabbit mAb [SN07-10]	ET1611-62	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	TMS1 / ASC Recombinant Rabbit mAb [PSH17-34] - BSA and Azide free	HA751622	Human, Mouse	WB,IHC-Fr,IHC-P
Relative expression	TMS1 / ASC Recombinant Rabbit mAb [PSH17-34]	HA723917	Human, Mouse	WB,IHC-Fr,IHC-P
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free	HA750897	Human, Mouse	WB
	TNF alpha Rabbit pAb	ER65189	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	TNF alpha Rabbit pAb	ER1919-22	Human, Mouse	WB
Cell treatment	TNF alpha Rabbit pAb	R1203-1	Mouse	WB,IF-Cell,ELISA
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07]	HA722022	Human, Mouse	WB
Relative expression, KD	TOMM20 Recombinant Rabbit mAb [ST04-72] - BSA and Azide free	HA750166	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP,IHC-Fr
	TOMM20 Recombinant Antibody [ST04-72] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610288	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression, KD	TOMM20 Recombinant Rabbit mAb [ST04-72]	ET1609-25	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP,IHC-Fr
	TOMM20 Recombinant Antibody [ST04-72] - Mouse IgG1 (Chimeric)	HA601454	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Total OXPHOS Antibody Cocktail Recombinant Rabbit Multiclonal Antibody [PSH15-50]	HA723756	Human	WB,IHC-P
	TRAIL Rabbit pAb	ER1917-99	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	TRAIL Recombinant Rabbit mAb [PSH14-40] - BSA and Azide free	HA751522	Human	WB,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	TRAIL Recombinant Rabbit mAb [PSH14-40]	HA723634	Human	WB,FC
	UCP-1 Rabbit pAb	ER1918-37	Mouse	WB,FC
	UCP1 Rabbit pAb	HA500118	Human, Mouse, Rat	WB,IHC-P
Relative expression	UCP1 Recombinant Rabbit mAb [PSH15-73]	HA723776	Mouse, Rat	WB,IHC-Fr,IHC-P,IP
	UCP1 Rabbit pAb	ER62425	Human, Mouse, Rat	WB
Relative expression	UCP1 Recombinant Rabbit mAb [PSH15-73] - BSA and Azide free	HA751580	Mouse, Rat	WB,IHC-Fr,IHC-P,IP
	VEGF Recombinant Rabbit mAb [SP07-01]	ET1604-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	VEGF Rabbit pAb	ER30607	Human	WB,IF-Cell,IHC-P,FC
	VEGF Recombinant Rabbit mAb [SP07-01] - BSA and Azide free	HA750074	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	VEGFA Recombinant Rabbit mAb [PSH19-17]	HA724051	Human, Mouse, Rat, Monkey	WB
	Vimentin Recombinant Rabbit mAb	IRS060	Mouse	mIHC
	Vimentin Recombinant Rabbit mAb [PDH0-01]	HA721174	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue,mIHC
	iFluor™ 647 Conjugated Vimentin Recombinant Rabbit mAb [SC60-05]	HA720165F	Human, Mouse, Rat	IF-Cell,IF-Tissue,FC
KO, Relative expression	Vimentin Recombinant Rabbit mAb [SC60-05]	ET1610-39	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC,IP
	Vimentin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS214RT	Human	mIHC
KO, Relative expression	Vimentin Recombinant Rabbit mAb [SC60-05] - BSA and Azide free	HA750215	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC,IP
	Vimentin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS263RT	Mouse	mIHC
KO	Vimentin Mouse mAb [D4-B11]	EM0401	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,IHC-Fr,IP
	Vimentin Recombinant Antibody - Rat IgG1 (Chimeric)	HA601564	Human, Mouse, Rat	IHC,ICC
	iFluor™ 488 Conjugated Vimentin Recombinant Rabbit mAb [SC60-05]	HA720166F	Human, Mouse, Rat	IF-Cell,IF-Tissue,FC
Relative expression, KO	Vimentin Mouse mAb [A6-C1]	M1412-1	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	Vimentin Recombinant Mouse mAb [A6-C1-R] - BSA and Azide free	HA610126	Human, Mouse, Rat	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	Vimentin Recombinant Mouse mAb [A6-C1-R]	HA601251	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KO	Vimentin Rabbit pAb	R1308-6	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Vimentin Recombinant Rabbit mAb	IRS044	Human	mIHC
KD	XBP1 Recombinant Rabbit mAb [JM10-62]	ET1703-23	Human, Mouse	WB,IHC-P,IF-Cell,IF-Tissue,FC