

胆固醇代谢研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	ABCA1 Rabbit pAb	R1510-42	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	ABCA1 Mouse mAb [4G2]	RT1009	Human, Mouse	WB,IP,IF,IHC-P
	ABCA6 Rabbit pAb	ER62544	Human	WB
	ABCG1 Recombinant Rabbit mAb [SU03-26]	ET1608-20	Human, Rat, Mouse	WB,IHC-P
	ABCG5 Mouse mAb [E10-B5]	EM1708-58	Human	WB,IHC-P,FC
	ABCG8 Rabbit pAb	ER1903-17	Human	WB,IHC-P,IF-Tissue
KD	ACAA2 Recombinant Rabbit mAb [JE56-67]	HA720006	Human, Mouse, Rat	WB,IHC-P
KD	ACAT2 Recombinant Rabbit mAb [JE58-83]	HA720032	Human, Mouse, Rat	WB,IHC-P,FC
	AMPK beta 2 Rabbit pAb	ER60091	Human, Mouse, Rat	WB
	AMPK gamma 1 Recombinant Rabbit mAb [JE07-47]	HA722691	Human, Mouse, Rat	WB,IP
	AMPK gamma 1/2/3 Rabbit pAb	ER62661	Human, Mouse, Rat	WB,IF-Cell
	AMPK gamma 2 Rabbit pAb	ER62662	Human, Mouse	WB,IHC-P,IF-Cell
	Angiopoietin-like 4 Rabbit pAb	R1310-5	Mouse	WB
	Angiopoietin-like 4 Rabbit pAb	ER1903-68	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Angiotensinogen Recombinant Rabbit mAb [SD201-02]	ET1612-48	Human	WB,IF-Cell,IF-Tissue,IP
	ANGPTL3 Recombinant Rabbit mAb [JE40-17]	HA722388	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	ANGPTL3 Rabbit pAb	ER1903-70	Human	WB
KD	ANGPTL4 Recombinant Rabbit mAb [PSH07-35]	HA722823	Human, Monkey	WB
KD	ANGPTL4 Recombinant Rabbit mAb [PSH07-35] - BSA and Azide free	HA751156	Human, Monkey	WB
	ANGPTL8 Mouse mAb [22G2]	EM2001-08	Human	WB
	ApoB Mouse mAb [B1-G6]	EM1710-01	Human	IF-Cell,FC
	APOC2 Rabbit pAb	ER65921	Human	IHC-P
	APOC2 Rabbit pAb	ER61698	Human	WB
	APOC3 Rabbit pAb	ER1802-47	Human	WB,IHC-P
Cell treatment, KD	APOE4 Recombinant Rabbit mAb [PSH02-85] - BSA and Azide free	HA750838	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment, KD	APOE4 Recombinant Rabbit mAb [PSH02-85]	HA721867	Human	WB
	ApoF Rabbit pAb	ER62688	Human	WB
KD	Apo-H Recombinant Rabbit mAb [JE37-02]	HA721649	Human	WB
	Apolipoprotein A1 Rabbit pAb	ER1802-29	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Apolipoprotein A1 Rabbit pAb	RT1046	Human, Mouse, Rat	WB,IP,IF
	Apolipoprotein A1 Mouse mAb	M0809-13	Human	WB,IF-Cell,IHC-P
	Apolipoprotein A1 Recombinant Rabbit mAb [JF0548]	ET1702-23	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Apolipoprotein A II Recombinant Rabbit mAb [JB70-35]	ET7107-47	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	Apolipoprotein B Recombinant Rabbit mAb [PSH19-76]	HA724100	Human, Mouse, Rat	WB
Cell treatment	Apolipoprotein CI / Apo-CI Recombinant Rabbit mAb [JE36-30]	HA722783	Human	WB
	Apolipoprotein E Rabbit pAb	0809-6	Human	WB
	Apolipoprotein E Recombinant Rabbit mAb [SC0536] - BSA and Azide free	HA750206	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Apolipoprotein E Recombinant Rabbit mAb [SC0536]	ET1610-22	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Relative expression	Apolipoprotein E Recombinant Antibody [PSH10-52] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610233	Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Apolipoprotein E Recombinant Antibody [PSH10-52] - Rat IgG1 (Chimeric)	HA601391	Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Apolipoprotein E Recombinant Rabbit mAb [PSH10-52]	HA723216	Mouse, Rat	WB,IHC-Fr,IHC-P
Relative expression	Apolipoprotein E Recombinant Rabbit mAb [PSH10-52] - BSA and Azide free	HA751355	Mouse, Rat	WB,IHC-Fr,IHC-P
	AR beta 3 Rabbit pAb	ER62746	Human, Mouse, Rat	WB
	BMAL1 Recombinant Rabbit mAb [JM17-34]	ET1705-5	Human, Mouse, Rat	WB,IF-Cell
	BMAL1 Rabbit pAb	ER1706-71	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	BRCA1 Mouse mAb [A7E3]	HA601009	Human, Mouse, Rat	IHC-P
	BRCA1 Recombinant Mouse mAb [A7E3-R] - BSA and Azide free	HA610177	Human, Mouse, Rat	IHC-P,IF-Cell
	BRCA1 Rabbit pAb	HA500015	Human, Mouse, Rat	WB,IHC-P
	BRCA1 Mouse mAb [D2-C6]	EM1710-25	Human	IHC-P
	BRCA1 Recombinant Mouse mAb [A7E3-R]	HA601304	Human, Mouse, Rat	IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	Cannabinoid receptor 1 Recombinant Rabbit mAb [PSH02-43] - BSA and Azide free	HA750808	Human, Mouse, Rat	WB,IHC-P
Relative expression	Cannabinoid receptor 1 Recombinant Rabbit mAb [PSH02-43]	HA721820	Human, Mouse, Rat	WB,IHC-P
	Caveolin-1 Mouse mAb [A3E4]	EM1902-01	Human, Mouse	WB,IHC-P,FC
	Caveolin-1 Recombinant Rabbit mAb [SZ02-01]	ET1603-1	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Caveolin-1 Rabbit pAb	0407-6	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Caveolin-1 Recombinant Rabbit mAb	IRS100RB	Mouse	mIHC
	Caveolin-1 Mouse mAb [A3E8]	EM1901-88	Human	WB,IF-Cell,IHC-P,FC
	Caveolin-1 Mouse mAb [A3E7]	EM1901-87	Human	WB,IF-Cell,IHC-P,FC
	Caveolin-1 Mouse mAb [19-D6]	EM40728	Human, Mouse, Rat	WB,IHC-P
	Caveolin-1 Mouse mAb [A3E6]	HA601005	Human	WB,IHC-P
Relative expression	CD36 Recombinant Rabbit mAb [PSH18-47]	HA724000	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	CD36 Recombinant Rabbit mAb	IRS118RB	Human	mIHC
	CD36 Rabbit pAb	ER1802-7	Human, Mouse, Rat	IF-Cell,IHC-P,FC
Cell treatment	CD36 Recombinant Rabbit mAb [JJ2005]	ET1701-24	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	CD36 Recombinant Rabbit mAb [JJ2005] - BSA and Azide free	HA750307	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	CETP Recombinant Rabbit mAb [JE02-61]	HA722752	Human, Mouse, Rat	WB,IHC-P
	CIDE-B Rabbit pAb	ER63015	Human, Mouse	WB,IHC-P,IF-Cell
	Clusterin Rabbit pAb	ER63048	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Clusterin Recombinant Rabbit mAb [PSH04-33]	HA722116	Human	WB,IHC-P,IF-Tissue
	Clusterin Rabbit pAb	ER2001-58	Human	WB,IHC-P,FC
	CNR1 Rabbit pAb	ER1906-46	Human, Rat	WB,IHC-P,IF-Tissue
KD	COX2 Mouse mAb [A3F7]	EM1902-12	Human, Rat	WB,IHC-P,FC
	COX2 Mouse mAb	RT1159	Human, Mouse, Rat	WB,IP,IF,IHC-P,ELISA
	COX2 Mouse mAb [A3F8]	EM1902-11	Human, Rat	WB,IHC-P,FC
	COX2 Recombinant Rabbit mAb [6J8]	RE6064	Human, Mouse	WB,ICC
KD	COX2/Cyclooxygenase 2 Recombinant Rabbit mAb [SC56-06] - BSA and Azide free	HA750207	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
KD	COX2/Cyclooxygenase 2 Recombinant Rabbit mAb [SC56-06]	ET1610-23	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CYB5R1 Rabbit pAb	ER63152	Human, Mouse, Rat	WB
KD	CYB5R3 Recombinant Rabbit mAb [JE63-51]	HA721008	Human	WB,IHC-P
	CYP11A1 Recombinant Rabbit mAb [PSH07-26] - BSA and Azide free	HA751147	Human, Mouse, Rat	WB
	CYP11A1 Recombinant Rabbit mAb [PSH07-26]	HA722814	Human, Mouse, Rat	WB
	CYP11A1 Rabbit pAb	ER1906-98	Rat	WB
	CYP11B1 Rabbit pAb	ER1906-99	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	CYP11B1/2 Rabbit pAb	ER63158	Human	WB
	CYP1B1 Rabbit pAb	HA500494	Human, Rat	WB,IHC-P
	CYP27A1 Recombinant Rabbit mAb [JG40-36]	ET7109-05	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	CYP7A1 Rabbit pAb	HA500207	Human, Rat	WB,IF-Cell,FC
	CYP7B1 Rabbit pAb	ER63187	Human	WB,IHC-P
	DHCR24 Rabbit pAb	ER1919-32	Human, Mouse, Rat	WB,IHC
	DNAJA4 Recombinant Rabbit mAb [JE65-20]	HA722730	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	EHD1 Recombinant Rabbit mAb [JE38-59]	HA722340	Human, Mouse, Rat	WB,IHC-P,IP
	Endothelial lipase Rabbit pAb	ER1908-11	Rat	WB,IHC-P,IF-Tissue
	EPHX2 Recombinant Rabbit mAb [JE64-07]	HA721056	Human	WB
KD	FDFT1 Recombinant Rabbit mAb [JE60-69]	HA720048	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KD	FDPS Recombinant Rabbit mAb [JE63-33]	HA721031	Human	WB,FC
	HDLBP Rabbit pAb	ER2001-48	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Hepatic Lipase Rabbit pAb	ER63675	Human	WB,IHC-P
	HMGCR Recombinant Rabbit mAb [JF0981]	ET1702-41	Human, Mouse, Rat	WB,IP,IHC-P
	HMGCR Recombinant Rabbit mAb [JF0981] - BSA and Azide free	HA750345	Human, Mouse, Rat	WB,IP,IHC-P
Relative expression	HMGCS2 Recombinant Rabbit mAb [JE54-76]	HA720083	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	HNF-4-alpha Recombinant Rabbit mAb	IRS055	Mouse	mIHC
	HNF-4-alpha Recombinant Rabbit mAb	IRS155RB	Human	mIHC
	HNF-4-alpha Recombinant Rabbit mAb [JE63-17]	HA721006	Human, Mouse, Rat	WB,IHC-P,mIHC,IF-Tissue
	HNF-4-alpha Recombinant Rabbit mAb [SN72-03]	ET1611-43	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	HNF-4-alpha Recombinant Rabbit mAb [SN72-03] - BSA and Azide free	HA750254	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	HNF-4 alpha/gamma Rabbit pAb	ER63693	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	Hormone sensitive lipase / HSL Recombinant Rabbit mAb [PSH07-22] - BSA and Azide free	HA751143	Human, Mouse, Rat	WB,IHC-P,IP
Relative expression	Hormone sensitive lipase / HSL Recombinant Rabbit mAb [PSH07-22]	HA722810	Human, Mouse, Rat	WB,IHC-P,IP
	Hormone sensitive lipase / HSL Recombinant Rabbit mAb	IRS157RB	Human	mIHC
	Human ANGPTL3 Recombinant Rabbit mAb [PSH17-40] - BSA and Azide free (Detector)	HA725080	Human	ELISA(Det)
	Biotin Conjugated Human ANGPTL3 Recombinant Rabbit mAb [PSH17-40] - Detector	HA725081B	Human	ELISA(Det),ELISA
	Human ANGPTL3 Recombinant Rabbit mAb [PSH17-39] - BSA and Azide free (Capture)	HA725079	Human	ELISA(Cap)
	Human ANGPTL4 Recombinant Rabbit mAb [PSH20-15] - BSA and Azide free (Detector)	HA725390	Human	ELISA(Det)
	Biotin Conjugated Human ANGPTL4 Recombinant Rabbit mAb [PSH20-15] - Detector	HA725390B	Human	ELISA(Det),ELISA
	Human ANGPTL4 Recombinant Rabbit mAb [PSH20-14] - BSA and Azide free (Capture)	HA725389	Human	ELISA(Cap)
	Human Apolipoprotein A1 Recombinant Rabbit mAb [PSH18-91] - BSA and Azide free (Detector)	HA724028	Human	ELISA(Det)
	Human Apolipoprotein A1 Recombinant Rabbit mAb [PSH18-90] - BSA and Azide free (Detector)	HA724027	Human	ELISA(Det)
	Biotin Conjugated Human Apolipoprotein A1 Recombinant Rabbit mAb [PSH18-90] - Detector	HA724027B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human Apolipoprotein A1 Recombinant Rabbit mAb [PSH18-91] - Detector	HA724028B	Human	ELISA(Det),ELISA
	Human Apolipoprotein A1 Recombinant Rabbit mAb [PSH18-89] - BSA and Azide free (Capture)	HA724026	Human	ELISA(Cap)
	Human Clusterin Recombinant Rabbit mAb [PSH15-55] - BSA and Azide free (Capture)	HA725091	Human	ELISA(Cap)
	Human Clusterin Recombinant Rabbit mAb [PSH15-56] - BSA and Azide free (Detector)	HA725092	Human	ELISA(Det)
	Biotin Conjugated Human Clusterin Recombinant Rabbit mAb [PSH15-56] - Detector	HA725093B	Human	ELISA(Det),ELISA
	HRP Conjugated Human Clusterin Recombinant Rabbit mAb [PSH15-56] - Detector	HA725093H	Human	ELISA(Det),ELISA
	Biotin Conjugated Human Lipoprotein A Recombinant Rabbit mAb [PSH13-82] - Detector	HA723572B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Lipoprotein A Recombinant Rabbit mAb [PSH13-82] - BSA and Azide free (Detector)	HA723571	Human	ELISA(Det)
	Human Lipoprotein A Recombinant Rabbit mAb [PSH13-83] - BSA and Azide free (Detector)	HA723573	Human	ELISA(Det)
	Biotin Conjugated Human Lipoprotein A Recombinant Rabbit mAb [PSH13-83] - Detector	HA723574B	Human	ELISA(Det),ELISA
	Human Lipoprotein A Recombinant Rabbit mAb [PSH13-81] - BSA and Azide free (Capture)	HA723570	Human	ELISA(Cap)
	HRP Conjugated Human Lipoprotein A Recombinant Rabbit mAb [PSH13-83] - Detector	HA723574H	Human	ELISA(Det),ELISA
	Biotin Conjugated Human PCSK9 Recombinant Rabbit mAb [PSH12-86] - Detector	HA725126B	Human	ELISA(Det),ELISA
	Human PCSK9 Recombinant Rabbit mAb [PSH12-85] - BSA and Azide free (Capture)	HA725124	Human	ELISA(Cap)
	HRP Conjugated Human PCSK9 Recombinant Rabbit mAb [PSH12-86] - Detector	HA725126H	Human	ELISA(Det),ELISA
	Human PCSK9 Recombinant Rabbit mAb [PSH12-86] - BSA and Azide free (Detector)	HA725125	Human	ELISA(Det)
	KLOTB Rabbit pAb	ER61816	Human, Mouse	WB
	LCAT Recombinant Rabbit mAb [JE35-16]	HA722868	Mouse, Rat	WB,IHC-P
	LDL Receptor Recombinant Rabbit mAb [SJ0197]	ET1606-47	Human, Rat	WB,IF-Cell,FC
	LDL Receptor Recombinant Rabbit mAb [SJ0197] - BSA and Azide free	HA750107	Human, Rat	WB,IF-Cell,FC
	LDL Receptor Rabbit pAb	R1512-12	Human, Mouse	WB,IF-Cell,FC
KD	LIMPII Recombinant Rabbit mAb [PSH03-25]	HA721967	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell
KD	LIMPII Recombinant Rabbit mAb [PSH03-25] - BSA and Azide free	HA750865	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell
	Lipin 1 Recombinant Rabbit mAb [JG37-69]	ET7108-51	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Relative expression	Lipoprotein a Recombinant Rabbit mAb [PSH11-29] - BSA and Azide free	HA751408	Human	WB,IHC-P,IF-Tissue
Relative expression	Lipoprotein a Recombinant Rabbit mAb [PSH11-29]	HA723311	Human	WB,IHC-P,IF-Tissue
	Lipoprotein a Recombinant Rabbit mAb [JE39-47]	HA722312	Human	WB,IHC-P
	Lipoprotein lipase Recombinant Rabbit mAb [JA22-02]	ET1704-48	Human	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Liver Carboxylesterase 1 Recombinant Rabbit mAb [JE33-04]	HA721416	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	LOX 1 Rabbit pAb	ER1912-42	Human, Mouse, Rat, Rabbit	WB,IHC-P,IF-Tissue
	LOX-1 Rabbit pAb	ER1706-53	Human, Mouse	WB,IHC-P,IF-Cell,FC
	LRP-1 Mouse mAb [A4-F5]	M1211-4	Mouse	WB,IHC-P
KD	LRP1 Recombinant Rabbit mAb [SA0290]	ET1601-1	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP,IHC-Fr
Relative expression	Lrp2 / Megalin Recombinant Rabbit mAb [PSH04-52] - BSA and Azide free	HA750959	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	Lrp2 / Megalin Recombinant Rabbit mAb [PSH02-91] - BSA and Azide free	HA750846	Human	IHC-P,IF-Tissue,WB
Relative expression	Lrp2 / Megalin Recombinant Rabbit mAb [PSH04-52]	HA722151	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	Lrp2 / Megalin Recombinant Rabbit mAb [PSH02-91]	HA721910	Human	IHC-P,IF-Tissue,WB
	LXR alpha Recombinant Rabbit mAb [JA20-38]	ET1704-51	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	LXR alpha Recombinant Rabbit mAb [JA20-38] - BSA and Azide free	HA750414	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	LXR alpha + LXR beta Rabbit pAb	ER1912-55	Human, Mouse	WB,FC
KD	LXR beta Recombinant Rabbit mAb [PSH0-56] - BSA and Azide free	HA750633	Human, Mouse, Rat	WB,IF-Cell
KD	LXR beta Recombinant Rabbit mAb [PSH0-56]	HA721403	Human, Mouse, Rat	WB,IF-Cell
	LXR beta Rabbit pAb	ER1912-56	Mouse	WB
	MC2-R Rabbit pAb	ER63925	Human, Mouse	WB,IF-Cell
	MC4 Receptor Rabbit pAb	ER1912-84	Human, Mouse, Rat	WB,IF-Cell
	Mouse Apolipoprotein E Recombinant Rabbit mAb [PSH16-83] - BSA and Azide free (Capture)	HA725221	Mouse	ELISA(Cap)
	Mouse Apolipoprotein E Recombinant Rabbit mAb [PSH16-84] - BSA and Azide free (Detector)	HA725222	Mouse	ELISA(Det)
	Biotin Conjugated Mouse Apolipoprotein E Recombinant Rabbit mAb [PSH16-84] - Detector	HA725222B	Mouse	ELISA(Det),ELISA
	MRAP Rabbit pAb	ER1913-16	Human	WB
	MTCO2 Recombinant Rabbit mAb [SC06-88]	ET1610-72	Human	WB,IF-Cell,IF-Tissue,IHC-P
	MVD Rabbit pAb	ER1913-40	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	MVD Rabbit pAb	HA500275	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	MVK Mouse mAb [A9G12]	HA601151	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	MVK Rabbit pAb	ER2001-64	Human	WB,FC
	MVK Rabbit pAb	ER64069	Human, Monkey	WB,IHC-P,IF-Cell
	MYLIP Rabbit pAb	ER64081	Human, Mouse, Rat	WB,IHC-P
KO	Niemann Pick C1 Recombinant Rabbit mAb [JB87-33]	ET7107-57	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Niemann Pick C2 Recombinant Rabbit mAb [PSH05-52] - BSA and Azide free	HA751011	Human, Mouse, Rat	WB,IHC-P
	Niemann Pick C2 Recombinant Rabbit mAb [PSH05-52]	HA722294	Human, Mouse, Rat	WB,IHC-P
	NPC2 Rabbit pAb	ER1914-03	Human	WB
KD	NSDHL Recombinant Rabbit mAb [PSH02-37]	HA721814	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
KD	NSDHL Recombinant Rabbit mAb [PSH02-36]	HA721813	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
KD	NSDHL Recombinant Rabbit mAb [PSH02-37] - BSA and Azide free	HA750802	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
KD	NSDHL Recombinant Rabbit mAb [JE64-86]	HA721088	Human	WB,IF-Cell,IHC-P
KD	NSDHL Recombinant Rabbit mAb [PSH02-36] - BSA and Azide free	HA750801	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
KD	PBR Recombinant Rabbit mAb [SA90-03]	ET1601-19	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	PBR Recombinant Rabbit mAb [SA90-03] - BSA and Azide free	HA750010	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	PCK1 Recombinant Rabbit mAb [JU84-39]	ET7106-81	Human	IF-Cell,IHC-P,WB
KD	PCK2 Recombinant Rabbit mAb [JB52-39]	ET7107-29	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	PCSK9 Recombinant Rabbit mAb [JE30-37]	HA722528	Human	WB,IHC-P,IP
Relative expression	PCSK9 Mouse mAb [2F1]	EM1701-26	Human, Mouse	WB,IHC-P,IF-Cell,FC
Relative expression	PCSK9 Mouse mAb [2F1] - BSA and Azide free	HA610371	Human, Mouse	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-Hormone sensitive lipase (S660) Recombinant Rabbit mAb [PSH10-28] - BSA and Azide free	HA751345	Mouse, Rat, Human	WB
Cell treatment	Phospho-Hormone sensitive lipase (S660) Recombinant Rabbit mAb [PSH10-28]	HA723195	Mouse, Rat, Human	WB
Cell treatment	Phospho-Hormone sensitive lipase (S853) Recombinant Rabbit mAb [SN06-39]	ET1611-19	Human, Mouse, Rat	WB,IHC-P
	PRKAG3 Mouse mAb [A7-C7]	EM1708-52	Human	IHC,IF
	S2P Rabbit pAb	ER64902	Human, Mouse, Rat	IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Scavenging Receptor SR-BI Recombinant Rabbit mAb [SR37-06]	ET1602-32	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	SOAT1 Rabbit pAb	ER1916-93	Mouse, Rat	WB
	SOAT2 Rabbit pAb	ER1916-94	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Sortilin Recombinant Rabbit mAb [JE60-00]	HA722070	Human	WB,IHC-P
	Sortilin / NT3 Recombinant Rabbit mAb [PSH06-16]	HA722513	Human, Mouse, Rat	WB,IHC-P
	Sortilin / NT3 Recombinant Rabbit mAb [PSH06-16] - BSA and Azide free	HA751054	Human, Mouse, Rat	WB,IHC-P
	SPTLC1 Rabbit pAb	ER65501	Human, Mouse, Rat	WB,IHC-P
	SRBP2 Rabbit pAb	ER60160	Human, Mouse, Rat	WB
KD	SREBP1 Recombinant Rabbit mAb [PSH04-61]	HA722160	Human	WB,IHC-P,IF-Cell
	SREBP1 Rabbit pAb	HA500210	Human, Mouse	WB,IHC-P
KD	SREBP1 Recombinant Rabbit mAb [PSH04-61] - BSA and Azide free	HA750968	Human	WB,IHC-P,IF-Cell
	SREBP1 Rabbit pAb	ER1917-19	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	StAR Recombinant Rabbit mAb [JE32-23]	HA721469	Human	WB,IHC-P
	StAR Rabbit pAb	ER1917-23	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	STAR3 Rabbit pAb	ER62478	Human, Mouse	WB
	Total OXPHOS Antibody Cocktail Recombinant Rabbit Multiclonal Antibody [PSH15-50]	HA723756	Human	WB,IHC-P
	VAPA Recombinant Rabbit mAb [JE77-48]	HA722372	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	VDAC1 Recombinant Rabbit mAb [SA93-03] - BSA and Azide free	HA750012	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
KD	VDAC1 Rabbit pAb	R1307-1	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P
KD	VDAC1 Recombinant Rabbit mAb [SA93-03]	ET1601-20	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
	VDAC1 / 2 / 3 Recombinant Rabbit mAb [PSH14-90] - BSA and Azide free	HA751534	Human, Mouse, Rat	WB
	VDAC1 / 2 / 3 Recombinant Rabbit mAb [PSH14-90]	HA723692	Human, Mouse, Rat	WB
	VDAC3 Rabbit pAb	ER1918-49	Mouse, Rat	WB,IHC-P,IF-Tissue
	VLDL Receptor Rabbit pAb	HA500194	Human, Mouse	WB,IF-Cell,IHC-P,FC
	VLDL Receptor Rabbit pAb	ER1918-54	Human	WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

