

成骨 & 软骨研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	ALCAM Mouse mAb [C7-E12]	EM1708-92	Human	WB,IF-Cell,IHC-P,FC
	ALCAM/CD166 Mouse mAb [7C2]	RT1029	Human, Mouse, Rat	WB,IP,IF,IHC-P
	Alkaline Phosphatase Recombinant Rabbit mAb [SA40-00] - BSA and Azide free	HA750013	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Alkaline Phosphatase Mouse mAb [6-10]	M0806-10	Human, Mouse	WB,IHC-P,FC
	Alkaline Phosphatase Recombinant Rabbit mAb [SA40-00]	ET1601-21	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Alkaline Phosphatase(For 4 isoform) Rabbit pAb	0806-3	Human	WB
	Alkaline Phosphatase(PLAP-1/PLAP-like) Rabbit pAb	0806-4	Human	WB,IHC-P,FC
KD	Cathepsin H Recombinant Rabbit mAb [JE56-25]	HA721352	Human	WB,IF-Cell
	Cathepsin H Rabbit pAb	ER60014	Human	WB
	Cathepsin K Rabbit pAb	HA500469	Human	WB,IF-Cell
KD	Cathepsin L Recombinant Rabbit mAb [PSH03-91]	HA722063	Human, Mouse, Rat, Monkey	WB,IHC-P,FC
KD	Cathepsin L Recombinant Rabbit mAb [PSH03-91] - BSA and Azide free	HA750916	Human, Mouse, Rat, Monkey	WB,IHC-P,FC
KD	Cathepsin L Recombinant Rabbit mAb [PSH03-92] - BSA and Azide free	HA750917	Human, Mouse, Rat, Monkey	WB,IHC-P
	Cathepsin L Rabbit pAb	ER1905-49	Human, Rat	WB,IHC-P,IF-Tissue,FC
KD	Cathepsin L Recombinant Rabbit mAb [PSH03-92]	HA722064	Human, Mouse, Rat, Monkey	WB,IHC-P
KD	Cathepsin L/V/K/H Recombinant Rabbit mAb [JM10-78]	ET1703-44	Human, Mouse, Rat	WB,IF-Tissue,IHC-P
	CCL2 Mouse mAb [E11-D8]	EM1710-22	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
Cell treatment	CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-16]	HA723979	Mouse, Rat	WB,IF-Cell,FC
	CCL2 / MCP-1 Rabbit pAb	HA500042	Human	WB,IHC-P,IF-Cell
	CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-44] - BSA and Azide free	HA751533	Mouse	WB,IF-Cell,FC
	CCL2 / MCP-1 Rabbit pAb	HA500267	Mouse	WB,IF-Cell,FC

	CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-44]	HA723691	Mouse	WB,IF-Cell,FC
Cell treatment, Relative expression	CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-27]	HA723309	Human	WB,IF-Cell,FC
Cell treatment, Relative expression	CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-28]	HA723310	Human, Mouse	WB,IF-Cell,IHC-P
Cell treatment, Relative expression	CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-27] - BSA and Azide free	HA751406	Human	WB,IF-Cell,FC
	CCL3 / MIP-1 alpha Rabbit pAb	R1510-9	Human, Mouse, Rat	IF-Cell,IHC-P,FC,ELISA
Cell treatment	CCL3 / MIP-1 alpha Rabbit pAb	R1510-10	Human	WB
Cell treatment, Relative expression	CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-28] - BSA and Azide free	HA751407	Human, Mouse	WB,IF-Cell,IHC-P
Cell treatment	CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH19-79]	HA724103	Mouse, Rat	WB,IF-Cell,FC,IP
	CD105 Mouse mAb [H4-A9]	EM1709-99	Human	WB,IF-Cell,IHC-P,FC
Relative expression	CD105 Recombinant Rabbit mAb [PSH15-96]	HA723801	Human	FC
	CD105 Rabbit pAb	HA500266	Human	WB,IHC-P
	CD105 Mouse mAb [A5E8]	HA600077	Human	WB
Cell treatment	CD105 Recombinant Rabbit mAb [JE60-59]	HA720072	Human	WB,IHC-P,IF-Tissue
Relative expression	CD105 Recombinant Rabbit mAb [PSH15-97]	HA723802	Human	FC
Relative expression	CD105 Recombinant Rabbit mAb [PSH15-96] - BSA and Azide free	HA751587	Human	FC
Relative expression	CD105 Recombinant Rabbit mAb [PSH15-97] - BSA and Azide free	HA751588	Human	FC
Relative expression	CD11b Mouse mAb [A2-9]	M1211-6	Human	WB,IHC-P,FC,IF-Tissue
	CD11b Recombinant Rabbit mAb	IRS014	Human	mlIHC
Relative expression	CD11b Recombinant Rabbit mAb [PSH03-96]	HA722075	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	CD11b Recombinant Rabbit mAb [JU93-81] - BSA and Azide free	HA750451	Human	WB,IP,IHC-P,IF-Cell,IF-Tissue
	CD11b Mouse mAb [A2-B5-D2]	EM1902-35	Human, Mouse	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD11b Recombinant Rabbit mAb [JU93-81]	ET1706-04	Human	WB,IP,IHC-P,IF-Cell,IF-Tissue,mIHC
	CD11b Mouse mAb [4-6]	EM1701-42	Human	IHC-P,IF-Cell,FC,WB,IF-Tissue
Relative expression	CD11b Recombinant Rabbit mAb [PSH03-96] - BSA and Azide free	HA750921	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	CD11b Recombinant Rabbit mAb [JE01-17]	HA722856	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue
	CD11b/ITAM Mouse mAb [4G2]	RT1107	Human, Mouse	WB,IP,IF,FC
KD	CD14 Recombinant Antibody [SC69-02] - Rat IgG1 (Chimeric)	HA601420		
	CD14 Rabbit pAb	R1107-2	Human, Mouse	WB,IF-Cell,FC
	CD14 Recombinant Rabbit mAb	IRS015	Human	mIHC
KD	CD14 Recombinant Antibody [SC69-02] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610260	0	
	CD14 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS106RT	Human	mIHC
KD	CD14 Recombinant Rabbit mAb [PD01-63]	HA721321	Human, Mouse	WB,IHC-P,IF-Tissue
KD	CD14 Recombinant Rabbit mAb [SC69-02]	ET1610-85	Human, Mouse	WB,IF-Cell,IHC-P,mIHC,FC
KD	CD14 Recombinant Rabbit mAb [SC69-02] - BSA and Azide free	HA750232	Human, Mouse	WB,IF-Cell,IHC-P,FC
	CD166 Recombinant Rabbit mAb [JE57-47]	HA722705	Human, Mouse, Rat	WB,IHC-P,FC,IHC-Fr
KD	CD166 Mouse mAb [A6A8]	HA600083	Human	IF-Cell,IHC-P,WB
KD	CD166 Mouse mAb [11-6]	M1012-7	Human, Mouse, Rat	WB,IHC-P,FC
KD	CD166 Mouse mAb [A6A9]	HA600082	Human	IF-Cell,IHC-P,WB
KD	CD166 Mouse mAb [5-5]	EM1902-36	Human	WB,IHC-P,IF-Cell,FC
	CD19 Recombinant mAb	IRS110RB	Mouse	mIHC
Relative expression	CD19 Recombinant Rabbit mAb [PSH13-44] - BSA and Azide free	HA751463	Human	WB,IF-Cell,FC
Relative expression	CD19 Recombinant Rabbit mAb [PSH13-44]	HA723533	Human	WB,IF-Cell,FC
	PE Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723763F	Human	FC
Relative expression	CD19 Recombinant Rabbit mAb [PSH03-94]	HA722073	Mouse	WB,IHC-P,IF-Tissue,IP,IHC-Fr,FC,IF-Cell,mIHC
	CD19 Recombinant Rabbit mAb	IRS018	Human	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD19 Recombinant Rabbit mAb [PSH03-94] - BSA and Azide free	HA750919	Mouse	WB,IHC-P,IF-Tissue,IP,IHC-Fr,FC,IF-Cell
	CD19 Mouse mAb [A3-A8]	EM40308	Human, Mouse	WB,IHC-P,FC
	CD19 Rabbit pAb	R1107-7	Human, Mouse	WB
Relative expression	CD19 Recombinant Rabbit mAb [JF100-06]	ET1702-93	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,mIHC
	FITC Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723762F	Human	FC
	CD19 Recombinant Rabbit mAb [JF099-9]	ET1702-74	Human, Mouse, Rat	WB,FC
	FITC Conjugated CD19 Recombinant Mouse mAb [PSH14-97]	HA601518F	Human	FC
Relative expression	CD19 Recombinant Rabbit mAb [JF100-06] - BSA and Azide free	HA750365	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	APC Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723761F	Human	FC
	CD191 Mouse mAb [8E10F9]	EM1710-71	Human	WB,FC
	CD34 Recombinant Mouse mAb [PDM0-12]	HA601109	Human	IHC-P,IF-Tissue
	CD34 Mouse mAb [15H1]	EM1901-01	Human, Mouse, Rat	IHC-P,FC
	CD34 Rabbit pAb	0108-2	Human, Mouse, Zebrafish	WB,IF-Cell,IHC-P
	CD34 Rabbit pAb	ER0802	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	CD34 Recombinant Rabbit mAb	IRS050	Mouse	mIHC
	CD34 Rabbit pAb	ER1706-66	Human, Mouse	WB,IHC-P
	CD34 Recombinant Rabbit mAb [SI16-01] - BSA and Azide free	HA750092	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	CD34 Recombinant Rabbit mAb [SI16-01]	ET1606-11	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,mIHC,FC
	CD34 Recombinant Rabbit mAb	IRS024	Human	mIHC
	CD40L Recombinant Rabbit mAb [JM11-34] - BSA and Azide free	HA750387	Human, Mouse, Rat	WB
	CD40L Recombinant Rabbit mAb [JM11-34]	ET1703-64	Human, Mouse, Rat	WB
	CD44 Mouse mAb [A7-1]	EM010805	Human	WB,IHC-P,FC
	CD44 Recombinant mAb	IRS122RB	Human	mIHC
Relative expression, KD	CD44 Recombinant Rabbit mAb [JE64-01]	HA722845	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD44 Recombinant Rabbit mAb [ST57-03] - BSA and Azide free	HA750192	Human	WB,IHC-P
	CD44 Recombinant Rabbit mAb	IRS123RB	Mouse	mlHC
Relative expression	CD44 Recombinant Rabbit mAb [ST57-03]	ET1609-74	Human	WB,IHC-P
Relative expression	CD45 Recombinant Rabbit Multiclonal Antibody [PSH14-19] - BSA and Azide free	HA751510	Human, Mouse	WB,IHC-P,IF-Cell
	CD45 Rabbit pAb	ER1901-29	Human	WB,IHC-P,FC
	CD45 Rabbit pAb	ER1905-93	Human	FC
	CD45 Recombinant Rabbit mAb	IRS025	Human	mlHC
Relative expression	CD45 Recombinant Rabbit mAb [PSH13-02]	HA723507	Mouse, Rat	WB,IHC-Fr,IF-Tissue
	CD45 Recombinant Rabbit mAb [JE03-05]	ET7111-03	Human	WB,IHC-P,FC,IF-Tissue,mlHC
Relative expression	CD45 Recombinant Rabbit mAb [PSH13-02] - BSA and Azide free	HA751457	Mouse, Rat	WB,IHC-Fr,IF-Tissue
Relative expression	CD45 Recombinant Rabbit Multiclonal Antibody [PSH14-19]	HA723606	Human, Mouse	WB,IHC-P,IF-Cell
	CD45 Rabbit pAb	ER1902-02	Human	WB,FC
	PE Conjugated CD45RO Recombinant Mouse mAb [PSH13-94]	HA601531F	Human	FC
	CD45RO Recombinant Mouse mAb [PSH13-94]	HA601441	Human	IHC-P,FC
	CD45RO Recombinant Mouse mAb	IRS124MS	Human	mlHC
	CD73 Mouse mAb [6E-G3]	EM1706-12	Human	WB,IHC-P,FC
Relative expression	CD73 Mouse mAb [A6F7]	HA601010	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
Relative expression	CD73 Recombinant Rabbit mAb [PSH09-57]	HA723119	Human	WB,IF-Cell,FC
Relative expression	CD73 Mouse mAb [A6F9]	HA601011	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC,IF-Cell
Relative expression	CD73 Recombinant Rabbit mAb [PSH09-57] - BSA and Azide free	HA751304	Human	WB,IF-Cell,FC
Relative expression	CD73 Mouse mAb [A6G2]	HA601004	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Relative expression	CD73 Recombinant Mouse mAb [A6G2-R] - BSA and Azide free	HA610172	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Relative expression	CD73 Mouse mAb [A6G1]	HA601012	Human	WB,IHC-P,IF-Tissue,FC,IF-Cell
Relative expression	CD73 Recombinant Mouse mAb [A6G2-R]	HA601299	Human, Mouse, Rat	WB,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD73 Recombinant Rabbit mAb [JM11-40]	ET1703-41	Human	WB,IHC-P,IF-Tissue
	CD73 Rabbit pAb	R1107-6	Human, Mouse	WB,IF-Cell,FC
KD	CD90 / THY1 Recombinant Rabbit mAb [JF10-09]	ET1702-92	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
	CD90 / THY1 Rabbit pAb	R1107-9	Human, Mouse, Rat	WB
	CD90 / THY1 Recombinant Rabbit mAb [SU35-07]	ET1608-46	Human, Rat	WB,IHC-P
KD	CD90 / THY1 Recombinant Rabbit mAb [JF10-09] - BSA and Azide free	HA750364	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
	CD90 / THY1 Recombinant Rabbit mAb [SU35-07] - BSA and Azide free	HA750149	Human, Rat	WB,IHC-P
	C/EBP Beta Recombinant Rabbit mAb [SC55-05]	ET1610-9	Human, Rat, Mouse	WB,FC,IP
	C/EBP Beta Recombinant Rabbit mAb [SC55-05] - BSA and Azide free	HA750235	Human, Rat, Mouse	WB,FC,IP
	CKR-1 Rabbit pAb	ER65775	Human	IHC-P
	Cleaved-Cathepsin L2 (L114) Rabbit pAb	ER60017	Human, Monkey	WB
	COL1A1 Recombinant Rabbit mAb	IRS228RB	Mouse	mIHC
	COL1A1 Recombinant Rabbit mAb [ST58-04] - BSA and Azide free	HA750188	Human, Cow	WB,IHC-P
	COL1A1 Recombinant Rabbit mAb [PSH06-20] - BSA and Azide free	HA751058	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP,IHC-Fr
	COL1A1 Recombinant Rabbit mAb [PSH06-20]	HA722517	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP,mIHC,IHC-Fr
	COL1A1 Recombinant mAb	IRS137RB	Human	mIHC
	COL1A1 Recombinant Rabbit mAb [ST58-04]	ET1609-68	Human, Cow	WB,IHC-P
	CXCL12/SDF1 Recombinant Rabbit mAb [JJ09-20]	ET1701-77	Human	WB,IF-Cell,IF-Tissue,IHC-P
	CXCL12/SDF1 Rabbit pAb	ER1802-6	Human	WB,IHC-P,FC
	CXCL12/SDF1 Rabbit pAb	ER1902-35	Human, Mouse, Rat	WB,IHC-P
	DAP12 Recombinant Rabbit mAb [JG38-70]	ET7108-76	Human	WB,IHC-P,FC
	DKK1 Rabbit pAb	ER1512-20	Human, Mouse	WB,IF-Cell,IHC-P,FC
Relative expression	DKK1 Recombinant Rabbit mAb [SC06-86] - BSA and Azide free	HA750228	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	DKK1 Recombinant Rabbit mAb [SC06-86]	ET1610-63	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Dlx-5 Rabbit pAb	ER63256	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KD	FAK Recombinant Rabbit mAb [SR46-04]	ET1602-25	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	FAK Rabbit pAb	RT1215	Human, Mouse, Rat	WB,IP,ICC,IHC-P
KD	FAK Recombinant Rabbit mAb [SR46-04] - BSA and Azide free	HA750039	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	Gli-1 Rabbit pAb	HA500465	Human	IHC-P,IF-Cell
	Gli1 Recombinant Rabbit mAb [JF09-08]	ET1702-85	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Gli1 Recombinant Rabbit mAb [JF09-08] - BSA and Azide free	HA750361	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	GLI1 Mouse mAb [F8-B9]	EM1708-22	Human	IF-Cell,IHC-P,FC
	Gli3 Recombinant Rabbit mAb [JG83-39]	ET7108-68	Human	WB,IF-Cell
	GLI4 Rabbit pAb	ER1909-64	Human	WB,IHC-P,IF-Tissue
	Glucose Transporter GLUT1 Recombinant Antibody [SA0377] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610348	Human, Mouse, Rat	IHC-Fr,IHC-P,IF-Tissue,WB
	iFluor™ 647 Conjugated Glucose Transporter GLUT1 Recombinant Rabbit mAb [SA0377]	HA720170F	Human, Mouse	IF-Tissue
	Glucose Transporter GLUT1 Recombinant Antibody [SA0377] - Rat IgG1 (Chimeric)	HA601521	Human, Mouse, Rat	IHC-Fr,IHC-P,IF-Tissue,WB
KD	Glucose Transporter GLUT1 Recombinant Rabbit mAb [SA0377] - BSA and Azide free	HA750002	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Glucose Transporter GLUT1 Recombinant Rabbit mAb	IRS151RB	Human	mIHC
	Glucose Transporter GLUT1 Rabbit pAb	ER1510-11	Human, Mouse	WB,IF-Cell,IHC-P,FC
KD	Glucose Transporter GLUT1 Recombinant Rabbit mAb [SA0377]	ET1601-10	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	iFluor™ 488 Conjugated Glucose Transporter GLUT1 Recombinant Rabbit mAb [SA0377]	HA720154F	Human, Mouse, Rat	IF-Cell,IF-Tissue
	Gremlin Rabbit pAb	ER1910-23	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	HAND2 Rabbit pAb	ER61828	Human, Rat, Mouse	WB
	Hex Recombinant Rabbit mAb [JE37-45]	HA722753	Human	WB,IP
	Histamine H2 receptor Rabbit pAb	HA500154	Human	WB,IHC-P
	HLA-DR Recombinant Rabbit mAb [SC06-78] - BSA and Azide free	HA750230	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	HLA-DR Rabbit pAb	R1510-34	Human, Mouse	WB,IF-Cell,IHC-P,FC
	HLA-DR Recombinant Mouse mAb [10-D8-R]	HA601252	Human	WB,IF-Cell,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HLA-DR Recombinant Rabbit mAb [SC06-78]	ET1610-66	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	HLA-DR Mouse mAb [10-D8]	M1701-5	Human	WB,IHC-P,IF-Cell
	HLA-DR Rabbit pAb	R1510-35	Human	WB,IF-Cell,IHC-P,FC
	HLA-DR Recombinant Mouse mAb	IRS154MS	Human	mIHC
	HLA-DR Recombinant Mouse mAb [10-D8-R] - BSA and Azide free	HA610127	Human	WB,IF-Cell,IHC-P,IF-Tissue
	HoxA11/D11 Rabbit pAb	ER63705	Human, Mouse	WB,IHC-P,IF-Cell
	Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-22] - BSA and Azide free (Detector)	HA725196	Human	ELISA(Det)
	Biotin Conjugated Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-22] - Detector	HA725196B	Human	ELISA(Det)
	Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-21] - BSA and Azide free (Capture)	HA725195	Human	ELISA(Cap)
	Human CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-77] - BSA and Azide free (Detector)	HA723378	Human	ELISA(Det)
	Biotin Conjugated Human CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-77] - Detector	HA723379B	Human	ELISA(Det)
	Human CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-76] - BSA and Azide free (Capture)	HA723377	Human	ELISA(Cap)
	Human CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-78] - BSA and Azide free (Capture)	HA723380	Human	ELISA(Cap)
	HRP Conjugated Human CCL5 Recombinant Rabbit mAb [PSH13-17] - Detector	HA725057H	Human	ELISA(Det),ELISA
	Human CCL5 Recombinant Rabbit mAb [PSH13-16] - BSA and Azide free (Capture)	HA725055	Human	ELISA(Cap)
	Human CCL5 Recombinant Rabbit mAb [PSH13-17] - BSA and Azide free (Detector)	HA725056	Human	ELISA(Det)
	Biotin Conjugated Human CCL5 Recombinant Rabbit mAb [PSH13-17] - Detector	HA725057B	Human	ELISA(Det),ELISA
	FITC Conjugated Human CD11b Recombinant Mouse mAb [A2-9-R]	HA600130F	Human	FC
	FITC Conjugated Human CD14 Recombinant Rabbit mAb [SC69-02]	HA720194F	Human	FC
	iFluor™ 647 Conjugated Human CD14 Recombinant Rabbit mAb [SC69-02]	HA720184F	Human	FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	FITC Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600125F	Human	FC
	iFluor™ 555 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA601556F	Human	FC
	iFluor™ 488 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600109F	Human	FC
	Human CD45 Recombinant Mouse mAb [PSH03-64]	HA601285	Human	IF-Cell,FC
	Human CD45 Recombinant Mouse mAb [PSH03-64] - BSA and Azide free	HA610158	Human	IF-Cell,FC
	FITC Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600112F	Human	FC
	Percp-Cy5.5 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600128F	Human	FC
	APC Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600127F	Human	FC
	PE Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600126F	Human	FC
	FITC Conjugated Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA600131F	Human	FC
	PE Conjugated Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA601525F	Human	FC
	Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA601341	Human	FC,IF-Cell,IHC-P
	Human CD73 Recombinant Rabbit mAb [PSH08-15] - BSA and Azide free (Capture)	HA722957	Human	ELISA(Cap)
	Biotin Conjugated Human CD73 Recombinant Rabbit mAb [PSH08-16] - Detector	HA722959B	Human	ELISA(Det)
	Human CD73 Recombinant Rabbit mAb [PSH08-16] - BSA and Azide free (Detector)	HA722958	Human	ELISA(Det)
	Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-36] - BSA and Azide free (Capture)	HA725007	Human	ELISA(Cap)
	Biotin Conjugated Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-37] - Detector	HA725009B	Human	ELISA(Det),ELISA
	Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-37] - BSA and Azide free (Detector)	HA725008	Human	ELISA(Det)
	Human Endoglin / CD105 Recombinant Rabbit mAb [PSH14-60] - BSA and Azide free (Capture)	HA723657	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Endoglin / CD105 Recombinant Rabbit mAb [PSH14-58] - BSA and Azide free (Capture)	HA723654	Human	ELISA(Cap)
	Biotin Conjugated Human Endoglin / CD105 Recombinant Rabbit mAb [PSH14-59] - Detector	HA723656B	Human	ELISA(Det),ELISA
	Human Endoglin / CD105 Recombinant Rabbit mAb [PSH14-59] - BSA and Azide free (Detector)	HA723655	Human	ELISA(Det)
	Human IL-4R Recombinant Rabbit mAb [PSH11-46] - BSA and Azide free (Detector)	HA723337	Human	ELISA(Det)
	Biotin Conjugated Human IL-4R Recombinant Rabbit mAb [PSH11-46] - Detector	HA723338B	Human	ELISA(Det),ELISA
	Human IL-4R Recombinant Rabbit mAb [PSH11-45] - BSA and Azide free (Capture)	HA723336	Human	ELISA(Cap)
	Human M-CSF Recombinant Rabbit mAb [PSH11-36] - BSA and Azide free (Capture)	HA723323	Human	ELISA(Cap)
	Biotin Conjugated Human M-CSF Recombinant Rabbit mAb [PSH11-38] - Detector	HA723326B	Human	ELISA(Det),ELISA
	Human M-CSF Recombinant Rabbit mAb [PSH11-37] - BSA and Azide free (Capture)	HA723324	Human	ELISA(Cap)
	Human M-CSF Recombinant Rabbit mAb [PSH11-38] - BSA and Azide free (Detector)	HA723325	Human	ELISA(Det)
	Human Osteopontin Recombinant Rabbit mAb [PSH08-58] - BSA and Azide free (Capture)	HA723004	Human	ELISA(Cap)
	Biotin Conjugated Human Osteopontin Recombinant Rabbit mAb [PSH08-59] - Detector	HA723007B	Human	ELISA(Det),ELISA
	Human Osteopontin Recombinant Rabbit mAb [PSH08-59] - BSA and Azide free (Detector)	HA723005	Human	ELISA(Det)
	Biotin Conjugated Human Osteopontin Recombinant Rabbit mAb [PSH08-60] - Detector	HA723008B	Human	ELISA(Det),ELISA
	Human Osteopontin Recombinant Rabbit mAb [PSH08-60] - BSA and Azide free (Detector)	HA723006	Human	ELISA(Det)
	Human Osteoprotegerin Recombinant Rabbit mAb [PSH13-31] - BSA and Azide free (Capture)	HA725067	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Osteoprotegerin Recombinant Rabbit mAb [PSH13-32] - BSA and Azide free (Detector)	HA725068	Human	ELISA(Det)
	Biotin Conjugated Human Osteoprotegerin Recombinant Rabbit mAb [PSH13-32] - Detector	HA725069B	Human	ELISA(Det),ELISA
	Human RANKL Recombinant Rabbit mAb [PSH16-59] - BSA and Azide free (Detector)	HA723861	Human	ELISA(Det)
	Biotin Conjugated Human RANKL Recombinant Rabbit mAb [PSH16-59] - Detector	HA723861B	Human	ELISA(Det),ELISA
	Human RANKL Recombinant Rabbit mAb [PSH16-60] - BSA and Azide free (Capture)	HA723862	Human	ELISA(Cap)
	Human RANKL Recombinant Rabbit mAb [PSH16-58] - BSA and Azide free (Capture)	HA723860	Human	ELISA(Cap)
	Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - BSA and Azide free (Detector)	HA723159	Human	ELISA(Det)
	Human sTNF RI Recombinant Rabbit mAb [PSH09-89] - BSA and Azide free (Capture)	HA723158	Human	ELISA(Cap)
	HRP Conjugated Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - Detector	HA723374H	Human	ELISA(Det)
	Biotin Conjugated Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - BSA and Azide free (Detector)	HA723160B	Human	ELISA(Det),ELISA
	Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - BSA and Azide free (Detector)	HA725065	Human	ELISA(Det)
	Biotin Conjugated Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - Detector	HA725066B	Human	ELISA(Det),ELISA
	Human VCAM-1 Recombinant Rabbit mAb [PSH12-46] - BSA and Azide free (Capture)	HA725064	Human	ELISA(Cap)
	HRP Conjugated Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - Detector	HA725066H	Human	ELISA(Det),ELISA
	HXA11 Rabbit pAb	ER60762	Human, Mouse	WB
	HXD11 Rabbit pAb	ER60770	Human, Mouse	WB
Cell treatment	IGF1 Recombinant Rabbit mAb [JE30-01]	HA722527	Human	WB
	IGF1 (E peptide) Rabbit pAb	HA500107	Human	WB,IHC-P,FC
	IGF I Rabbit pAb	ER1911-10	Human, Mouse, Rat, Dog	WB,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	IL-4R Rabbit pAb	ER1911-37	Human	WB
	IL-4R Recombinant Rabbit mAb [PSH10-86] - BSA and Azide free	HA751378	Human, Mouse, Rat	WB,IP,IF-Cell
	IL-4R Rabbit pAb	HA500453	Human	WB,IF-Cell
	IL-4R Recombinant Rabbit mAb [PSH10-86]	HA723253	Human, Mouse, Rat	WB,IP,IF-Cell
	IL-4R Rabbit pAb	HA500454	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
KD	Integrin beta 1 Recombinant Rabbit mAb [SR30-03]	ET1602-10	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Integrin beta 1 Mouse mAb [B10-A5]	EM1901-77	Human, Mouse	WB,IHC-P,FC
KD	Integrin beta 1 Rabbit pAb	ER31001	Human, Mouse	WB,IHC-P,FC
KD	Integrin beta 1 Recombinant Rabbit mAb [SA40-08] - BSA and Azide free	HA750008	Human, Mouse, Rat, Zebrafish	WB,IHC-P
KD	Integrin beta 1 Recombinant Rabbit mAb [SA40-08]	ET1601-17	Human, Mouse, Rat, Zebrafish	WB,IHC-P
Cell treatment	IP10 Recombinant Rabbit mAb [JA10-82]	ET1704-27	Human	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	IP10 Recombinant Rabbit mAb [JA10-82] - BSA and Azide free	HA750411	Human	WB,IF-Cell,IF-Tissue,IHC-P
	M-CSF Recombinant Rabbit mAb [SU0413]	ET1609-1	Human, Mouse	WB,IHC-P,IF-Cell,IP
	Biotin Conjugated Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-45] - Detector	HA725042B	Mouse	ELISA(Det),ELISA
	Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-44] - BSA and Azide free (Capture)	HA725040	Mouse	ELISA(Cap)
	Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-45] - BSA and Azide free (Detector)	HA725041	Mouse	ELISA(Det)
	Biotin Conjugated Mouse CCL3 Recombinant Rabbit mAb [PSH08-40] - Detector	HA722984B	Mouse	ELISA(Det),ELISA
	Mouse CCL3 Recombinant Rabbit mAb [PSH08-40] - BSA and Azide free (Detector)	HA722983	Mouse	ELISA(Det)
	Mouse CCL3 Recombinant Rabbit mAb [PSH08-39] - BSA and Azide free (Capture)	HA722982	Mouse	ELISA(Cap)
	Mouse CCL5 Recombinant Rat mAb	HA601395	Human, Mouse	ELISA
	Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-16] - BSA and Azide free (Detector)	HA725047	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-15] - BSA and Azide free (Capture)	HA725046	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-16] - Detector	HA725048B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse DKK1 Recombinant Rabbit mAb [PSH17-70] - Detector	HA723939B	Mouse	ELISA(Det),ELISA
	Mouse DKK1 Recombinant Rabbit mAb [PSH17-70] - BSA and Azide free (Detector)	HA723939	Mouse	ELISA(Det)
	Mouse DKK1 Recombinant Rabbit mAb [PSH17-69] - BSA and Azide free (Capture)	HA723938	Mouse	ELISA(Cap)
	Mouse M-CSF/CSF1 Recombinant Rabbit mAb [PSH15-84] - BSA and Azide free (Capture)	HA723792	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse M-CSF/CSF1 Recombinant Rabbit mAb [PSH15-85] - Detector	HA723794B	Mouse	ELISA(Det),ELISA
	Mouse M-CSF/CSF1 Recombinant Rabbit mAb [PSH15-85] - BSA and Azide free (Detector)	HA723793	Mouse	ELISA(Det)
	Biotin Conjugated Mouse Osteopontin Recombinant Rabbit mAb [PSH17-50] - Detector	HA723925B	Mouse	ELISA(Det),ELISA
	HRP Conjugated Mouse Osteopontin Recombinant Rabbit mAb [PSH17-50] - Detector	HA723925H	Mouse	ELISA(Det),ELISA
	Mouse Osteopontin Recombinant Rabbit mAb [PSH17-50] - BSA and Azide free (Detector)	HA723925	Mouse	ELISA(Det)
	Mouse Osteopontin Recombinant Rabbit mAb [PSH17-49] - BSA and Azide free (Capture)	HA723924	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse Pro-Collagen I alpha 1 Recombinant Rabbit mAb [PSH07-99]	HA722936B	Mouse	ELISA(Det)
	Mouse Pro-Collagen I alpha 1 Recombinant Rabbit mAb [PSH07-99] - BSA and Azide free (Detector)	HA722935	Mouse	ELISA(Det)
	Mouse Pro-Collagen I alpha 1 Recombinant Rabbit mAb [PSH07-98] - BSA and Azide free (Capture)	HA722934	Mouse	ELISA(Cap)
	HRP Conjugated Mouse Pro-Collagen I alpha 1 Recombinant Rabbit mAb [PSH07-99] - Detector	HA722936H	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Mouse RANKL Recombinant Rabbit mAb [PSH07-91] - BSA and Azide free (Capture)	HA722923	Mouse	ELISA(Cap)
Cell treatment	Mouse RANKL Recombinant Rabbit mAb [PSH07-92] - BSA and Azide free (Detector)	HA722924	Mouse	ELISA(Det)
Cell treatment	Mouse RANKL Recombinant Rabbit mAb [PSH07-90] - BSA and Azide free (Capture)	HA722922	Mouse	ELISA(Cap)
Cell treatment	Biotin Conjugated Mouse RANKL Recombinant Rabbit mAb [PSH07-92]	HA722925B	Mouse	ELISA(Det)
	HRP Conjugated Mouse VCAM1 Recombinant Rabbit mAb [PSH16-32] - Detector	HA723830H	Mouse	ELISA(Det),ELISA
	Mouse VCAM1 Recombinant Rabbit mAb [PSH16-33] - BSA and Azide free (Capture)	HA723831	Mouse	ELISA(Cap)
	Mouse VCAM1 Recombinant Rabbit mAb [PSH16-32] - BSA and Azide free (Detector)	HA723830	Mouse	ELISA(Det)
	Mouse VCAM1 Recombinant Rabbit mAb [PSH16-31] - BSA and Azide free (Capture)	HA723829	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse VCAM1 Recombinant Rabbit mAb [PSH16-32] - Detector	HA723830B	Mouse	ELISA(Det),ELISA
	Msx-2 Rabbit pAb	ER60005	Human, Mouse	WB,IHC-P
	ODF2 Recombinant Rabbit mAb [JE35-01]	HA721425	Human	WB
	OSTCN Rabbit pAb	ER61846	Human, Rat, Mouse	WB
	Osteocalcin Rabbit pAb	ER1919-20	Human, Mouse	WB,IHC-P,IF-Tissue
	Osteopontin Rabbit pAb	HA500464	Human, Mouse	WB,FC
Cell treatment	Osteopontin Recombinant Rabbit mAb [PSH09-24] - BSA and Azide free	HA751278	Human, Mouse	WB,IF-Cell,FC
	Osteopontin Recombinant Rabbit mAb [PSH14-11] - BSA and Azide free	HA751506	Mouse, Rat	IHC-P,IHC-Fr,WB,IF-Cell
Cell treatment	Osteopontin Recombinant Rabbit mAb [PSH09-24]	HA723082	Human, Mouse	WB,IF-Cell,FC
Cell treatment	Osteopontin Rabbit pAb	ER1802-16	Human, Mouse	IF-Cell,IHC-P,ELISA
	Osteopontin Rabbit pAb	ER1802-15	Human, Mouse	IHC-P,ELISA
	Osteopontin Recombinant Rabbit mAb	IRS200RB	Human	mIHC
	Osteopontin Recombinant Rabbit mAb [PSH14-11]	HA723601	Mouse, Rat	IHC-P,IHC-Fr,WB,IF-Cell
	Osteopontin Rabbit pAb	0806-8	Human, Zebrafish	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Osteopontin Rabbit pAb	0806-6	Mouse, Human, Rat	WB
	Osteopontin Recombinant Rabbit mAb [PSH09-23] - BSA and Azide free	HA751277	Human	IHC-Fr,IHC-P,IF-Tissue
	Osteopontin Recombinant Rabbit mAb [PSH09-23]	HA723081	Human	IHC-Fr,IHC-P,IF-Tissue
	Osteoprotegerin Mouse mAb [9C2]	EM1701-99	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Osteoprotegerin Rabbit pAb	R1608-4	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Osteoprotegerin Mouse mAb [9C1]	EM1701-98	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Osterix Rabbit pAb	ER1914-47	Human, Mouse	WB
	Phospho-FAK (S722) Mouse mAb [4G1]	RT1216	Human, Mouse, Rat	WB,IP,IF,IHC-P
Cell treatment	Phospho-FAK (Y397) Recombinant Rabbit mAb [SC54-07] - BSA and Azide free	HA750212	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-FAK (Y397) Recombinant Rabbit mAb [SC54-07]	ET1610-34	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-Syk (Y352) Recombinant Rabbit mAb [PSH11-66] - BSA and Azide free	HA751422	Human	WB,IF-Cell,FC
Cell treatment	Phospho-Syk (Y352) Recombinant Rabbit mAb [PSH11-66]	HA723362	Human	WB,IF-Cell,FC
Cell treatment	Phospho-Syk (Y525 + Y526) Recombinant Rabbit mAb [PSH03-70]	HA722034	Human	WB,IF-Cell
Cell treatment	Phospho-Syk (Y525 + Y526) Recombinant Rabbit mAb [PSH03-70] - BSA and Azide free	HA750904	Human	WB,IF-Cell
	Placental alkaline phosphatase Recombinant Rabbit mAb [JM22-53]	ET1705-21	Human	WB,IF-Tissue,IHC-P
	Placental alkaline phosphatase (PLAP) Recombinant Rabbit mAb [PD00-19]	HA721151	Human	WB
	PRX1 Rabbit pAb	ER64686	Human, Mouse, Rat	WB
	PTH-rP Rabbit pAb	ER65726	Human	WB
	PTPRC Mouse mAb [C11-C10]	EM1708-69	Human	WB,IHC-P
	PU.1/Spi1 Recombinant Rabbit mAb [JB72-35]	ET7107-49	Human	WB,IHC-P,FC
	PU.1 / Spi1 Recombinant Rabbit mAb [PD00-14]	HA721162	Human	WB,IHC-P
	RANKL Rabbit pAb	HA500369	Human, Mouse, Rat	WB,IF-Cell
	RANTES Recombinant Rabbit mAb [JM03-45]	ET1705-70	Human, Mouse	WB,IHC-P,IF-Cell,IF-Tissue
	RUNX1 Recombinant Rabbit mAb [PSH20-30]	HA724141	Human	IHC-P,ChIP
	RUNX1+RUNX2+RUNX3 Recombinant Rabbit mAb [SD0803]	ET1612-49	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	RUNX1+RUNX2+RUNX3 Recombinant Rabbit mAb [SD0803] - BSA and Azide free	HA750296	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	RUNX2 Rabbit pAb	ER1802-78	Human	WB,IHC-P,IF-Cell
	RUNX2 Recombinant Rabbit mAb	IRS238RB	Mouse	mIHC
	RUNX2 Recombinant Rabbit mAb	IRS047	Human	mIHC
Relative expression, KD	RUNX2 Recombinant Rabbit mAb [SD208-0] - BSA and Azide free	HA750295	Human, Mouse, Rat	IF-Cell,IF-Tissue,IHC-P,WB,FC,IP
Relative expression, KD	RUNX2 Recombinant Rabbit mAb [SD208-0]	ET1612-47	Human, Mouse, Rat	IF-Cell,IF-Tissue,IHC-P,WB,FC,IP,mIHC
Relative expression	RUNX3 Mouse mAb [A7B7]	HA600058	Human, Rat	WB,IHC-P
	RUNX3 Mouse mAb [B8-E2]	EM1707-92	Human	IF-Cell,IHC-P,FC
	Sclerostin Rabbit pAb	ER1916-29	Human, Rat	WB,IHC-P,IF-Tissue
	SIAL Rabbit pAb	ER62043	Human, Rat, Mouse	WB
Relative expression	SOX9 Recombinant Rabbit mAb [PSH13-59] - BSA and Azide free	HA751478	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell,IP,FC,ChIP
	SOX9 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS280RT	Human	mIHC
Relative expression, KD	SOX9 Recombinant Rabbit mAb [SN74-09]	ET1611-56	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-P,IP,FC,CUT&Tag-seq,mIHC
	SOX9 Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS279MS	Human	mIHC
Relative expression	SOX9 Recombinant Rabbit mAb [PSH13-59]	HA723548	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell,IP,FC,ChIP
	SOX9 Recombinant Rabbit mAb	IRS241RB	Mouse	mIHC
Relative expression, KD	SOX9 Recombinant Rabbit mAb [SN74-09] - BSA and Azide free	HA750262	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-P,IP,FC,CUT&Tag-seq
	SOX9 Recombinant Rabbit mAb	IRS199RB	Human	mIHC
	Sp7 / Osterix Recombinant Rabbit mAb [PSH07-29] - BSA and Azide free	HA751150	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
	Sp7 / Osterix Recombinant Rabbit mAb [PSH07-29]	HA722817	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
KD	Syk Recombinant Rabbit mAb [PSH0-76] - BSA and Azide free	HA750655	Human	WB,IF-Cell
	Syk Rabbit pAb	HA500080	Human, Rat	WB,IF-Cell,IHC-P
KD	Syk Recombinant Rabbit mAb [PSH0-76]	HA721494	Human	WB,IF-Cell
	Tafazzin Rabbit pAb	ER65093	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	TGF beta 3 Recombinant Rabbit mAb [PSH09-26]	HA723084	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	TGF beta 3 Recombinant Rabbit mAb [PSH09-26] - BSA and Azide free	HA751280	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	TGF beta 3 Rabbit pAb	ER1917-64	Human, Mouse, Rat, Pig	WB,IHC-P,IF-Tissue,FC
	Thy-1 (OX7) Mouse mAb [1G2]	RT1615	Mouse, Rat	WB,IP,IF,IHC-P,FC
	TNFR1 Mouse mAb [2G2]	RT1624	Human, Mouse, Rat	WB,IP,IF,IHC-P
	TNF Receptor I Rabbit pAb	HA500140	Human	WB,IHC-P,FC
	TNFSF11 Mouse mAb [8A7B9]	EM1710-42	Human, Monkey	WB,IHC-P,IF-Cell,FC
	Transferrin Receptor (CD71) Rabbit pAb	ER1706-35	Human	IHC-P
KD	Transferrin Receptor (CD71) Recombinant Rabbit mAb [JF0956] - BSA and Azide free	HA750333	Human, Mouse, Rat	WB,IP,FC
KD	Transferrin Receptor (CD71) Recombinant Rabbit mAb [JF0956]	ET1702-06	Human, Mouse, Rat	WB,IP,FC
Relative expression	Transferrin Receptor (CD71) Recombinant Rabbit mAb [PSH13-56]	HA723545	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP,IF-Tissue
	Twist Mouse mAb [1G3]	RT1635	Human, Mouse, Rat	WB,IP,IF,FC
	Twist1 Recombinant Rabbit mAb [PSH11-98]	HA723411	Human, Mouse, Rat	WB
	Twist1 Recombinant Rabbit mAb [PSH11-98] - BSA and Azide free	HA751431	Human, Mouse, Rat	WB
	TWIST1 Mouse mAb [E6-F4]	EM1710-24	Human, Mouse	WB,IF-Cell,IHC-P,FC
Relative expression, KD	VCAM1 Recombinant Rabbit mAb [SA05-04] - BSA and Azide free	HA750009	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	VCAM1 Rabbit pAb	R1512-13	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression, KD	VCAM1 Recombinant Rabbit mAb [SA05-04]	ET1601-18	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层