

神经标志物研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	5-HTR3 Rabbit pAb	ER1903-14	Human, Mouse, Rat	WB
KD	ABCG2 Recombinant Rabbit mAb [JF0994]	ET1702-40	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,FC
	ACL6B Rabbit pAb	ER60371	Human, Mouse	WB
KD	ALDH1A1 Recombinant Rabbit mAb [SY11-02]	ET1605-24	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IP
	ALDH1A1 Rabbit pAb	R1706-2	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	ALDH1A1 Mouse mAb [8F1]	EM1801-07	Human, Rat, Mouse	WB,IHC-P,IF-Cell,FC,IHC-Fr
	ALDH1A1 Rabbit pAb	R1706-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	ALDH1A1 Recombinant Mouse mAb [8F1-R]	HA601215	Human	WB,IF-Cell,FC
KD	ALDH1A1 Recombinant Mouse mAb [8F1-R] - BSA and Azide free	HA610093	Human	WB,IF-Cell,FC
	ALDH1L1 Recombinant Rabbit mAb [JU53-54] - BSA and Azide free	HA750457	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
	ALDH1L1 Recombinant Rabbit mAb [JU53-54]	ET7106-64	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
	Beta III Tubulin Rabbit pAb	0807-2	Human, Mouse, Rat	WB,IF-Cell,FC
	iFluor™ 594 Conjugated Beta III Tubulin Mouse mAb [A8-D10]	HA600133F	Human, Rat	IF,ICC,FC
	Beta III Tubulin Recombinant Rabbit mAb [SP06-00]	ET1604-17	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
	Beta III Tubulin Rabbit pAb	R1401-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	Beta III Tubulin Recombinant Rabbit mAb [PSH15-32] - BSA and Azide free	HA751554	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP
	Beta III Tubulin Rabbit pAb	0805-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
	Beta III Tubulin Mouse mAb [A8-D10]	M0805-8	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
	iFluor™ 488 Conjugated Beta III Tubulin Mouse mAb [A8-D10]	HA600106F	Human, Mouse, Rat	IF-Cell,FC,IF-Tissue
	Beta III Tubulin Recombinant Rabbit mAb [SP06-00] - BSA and Azide free	HA750071	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
Relative expression	Beta III Tubulin Recombinant Rabbit mAb [PSH15-32]	HA723738	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	bFGF Rabbit pAb	ER1904-22	Human, Mouse	WB,FC
Relative expression	BLBP Recombinant Rabbit mAb [PSH0-89] - BSA and Azide free	HA750669	Human	WB,IHC-P
Relative expression	BLBP Recombinant Rabbit mAb [PSH0-88] - BSA and Azide free	HA750668	Human	WB
Relative expression	BLBP Recombinant Rabbit mAb [PSH0-88]	HA721523	Human	WB
Relative expression	BLBP Recombinant Rabbit mAb [PSH0-89]	HA721524	Human	WB,IHC-P
	BMP2 Recombinant Rabbit mAb [JE10-29]	HA722693	Human, Mouse, Rat	WB
	BMP2 Rabbit pAb	ER80602	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Brain-derived tau/BD-tau Recombinant Rabbit mAb [PSH16-47] - BSA and Azide free (Detector)	HA723848	Human	ELISA(Det)
	Biotin Conjugated Brain-derived tau/BD-tau Recombinant Rabbit mAb [PSH16-47] - Detector	HA723849B	Human	ELISA(Det),ELISA
	Brain-derived tau/BD-tau Recombinant Rabbit mAb [PSH16-46] - BSA and Azide free (Capture)	HA723847	Human, Mouse, Rat	ELISA(Cap)
Relative expression	Calbindin Recombinant Rabbit mAb [JF05-01] - BSA and Azide free	HA750351	Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,IHC-Fr
	Calbindin Recombinant Rabbit mAb	IRS098RB	Human	mIHC
Relative expression	Calbindin Recombinant Rabbit mAb [JF05-01]	ET1702-54	Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,mIHC,IHC-Fr
Relative expression	Calbindin Recombinant Antibody [JF05-01] - Rat IgG1 (Chimeric)	HA601563	Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Calbindin Recombinant Antibody [JF05-01] - Guinea pig IgG2 (Chimeric)	HA601562	Mouse, Rat	IHC-Fr,WB
	Calbindin Recombinant Rabbit mAb	IRS221RB	Mouse	mIHC
Relative expression	CD11b Recombinant Rabbit mAb [JU93-81] - BSA and Azide free	HA750451	Human	WB,IP,IHC-P,IF-Cell,IF-Tissue
	CD11b Mouse mAb [A2-B5-D2]	EM1902-35	Human, Mouse	WB
Relative expression	CD11b Recombinant Rabbit mAb [PSH03-96]	HA722075	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	CD11b Recombinant Rabbit mAb [JU93-81]	ET1706-04	Human	WB,IP,IHC-P,IF-Cell,IF-Tissue,mIHC
Relative expression	CD11b Mouse mAb [A2-9]	M1211-6	Human	WB,IHC-P,FC,IF-Tissue
	CD11b Mouse mAb [4-6]	EM1701-42	Human	IHC-P,IF-Cell,FC,WB,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD11b Recombinant Rabbit mAb [JE01-17]	HA722856	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue
	CD11b Recombinant Rabbit mAb	IRS014	Human	mIHC
Relative expression	CD11b Recombinant Rabbit mAb [PSH03-96] - BSA and Azide free	HA750921	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	CD11b/ITAM Mouse mAb [4G2]	RT1107	Human, Mouse	WB,IP,IF,FC
	CD133 Rabbit pAb	ER1901-63	Human	WB,FC
Cell treatment	CD133 Recombinant Rabbit mAb [PSH06-48] - BSA and Azide free	HA751072	Mouse	WB,IHC-P,IHC-Fr
	CD133 Rabbit pAb	0804-5	Human, Mouse, Rat	WB,FC
Relative expression	CD133 Recombinant Mouse mAb [A8C3-R] - BSA and Azide free	HA610009	Human	WB,IHC-P
	CD133 Mouse mAb [A8C2]	HA601025	Human	WB,IHC-P,FC
	CD133 Rabbit pAb	ER31008	Human	WB,IHC-P,FC
Relative expression	CD133 Recombinant Mouse mAb [A8C3-R]	HA601062	Human	WB,IHC-P
Relative expression, Cell treatment	CD133 Mouse mAb [A8C3]	HA601024	Human	WB,IHC-P
Cell treatment	CD133 Recombinant Rabbit mAb [PSH06-48]	HA722667	Mouse	WB,IHC-P,IHC-Fr
	CD271 Rabbit pAb	1003-9	Mouse, Rat	WB,IF-Cell,IHC-P
	CD40 Recombinant Rabbit mAb [PSH10-01] - BSA and Azide free	HA751331	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue
	CD40 Mouse mAb [A1E5]	EM1901-41	Human	WB,IHC-P
Relative expression	CD40 Recombinant Rabbit mAb [PSH05-90] - BSA and Azide free	HA751037	Human, Cynomolgus monkey, Pig	WB,IF-Cell,FC
	CD40 Recombinant Rabbit mAb [PSH10-01]	HA723170	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue
	CD40 Rabbit pAb	ER1803-54	Human, Rat, Mouse	WB,IF-Cell,IHC-P,FC
Relative expression	CD40 Mouse mAb [A1E6]	EM1901-40	Human	WB,IF-Cell,IHC-P
Relative expression	CD40 Recombinant Rabbit mAb [PSH05-90]	HA722488	Human, Cynomolgus monkey, Pig	WB,IF-Cell,FC
Relative expression	CD68 Recombinant Mouse mAb [PDM0-13] - BSA and Azide free	HA610034	Human	IHC-P,IF-Cell,FC,IF-Tissue
Relative expression	CD68 Recombinant Mouse mAb [A3C3-R]	HA601290	Human	WB,IF-Cell,IHC-P,FC
	CD68 Recombinant Antibody [PDM0-13] - Rabbit IgG (Chimeric)	HA723708	Human	IF-Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD68 Mouse mAb [A3C4]	EM1901-95	Human	WB,IHC-P,FC,mIHC
Cell treatment	CD68 Recombinant Rabbit mAb [PSH05-47] - BSA and Azide free	HA751006	Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IF-Tissue
Relative expression	CD68 Recombinant Mouse mAb [A3C2-R] - BSA and Azide free	HA610164	Human	WB,IF-Cell,IHC-P
Relative expression	CD68 Recombinant Mouse mAb [PDM0-13]	HA601115	Human	IHC-P,IF-Cell,FC,mIHC,IF-Tissue
Relative expression, Cell treatment	CD68 Recombinant Rabbit mAb [PSH06-65] - BSA and Azide free	HA751094	Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,IHC-Fr
Relative expression	CD68 Recombinant Mouse mAb [A3C4-R] - BSA and Azide free	HA610165	Human	WB,IF-Cell,IHC-P
	CD68 Recombinant Antibody - Rabbit IgG (Chimeric)	IRS266RB	Human	mIHC
	CD68 Recombinant Antibody [PDM0-13] - Rabbit IgG (Chimeric) - BSA and Azide free	HA751544	Human	IF-Tissue,IHC-P
Relative expression	CD68 Recombinant Mouse mAb [A3C2-R]	HA601291	Human	WB,IF-Cell,IHC-P
Relative expression	CD68 Recombinant Mouse mAb [A3C4-R]	HA601292	Human	WB,IF-Cell,IHC-P
	CD68 Recombinant Mouse mAb	IRS004	Human	mIHC
	CD68 Recombinant Rabbit mAb [SN07-27]	ET1611-53	Human	WB
Cell treatment	CD68 Recombinant Rabbit mAb [PSH05-47]	HA722285	Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IF-Tissue
Relative expression, Cell treatment	CD68 Recombinant Rabbit mAb [PSH06-65]	HA722737	Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,IHC-Fr
	CD68 Rabbit pAb	ER1901-32	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	CD68 Mouse mAb [A3C3]	EM1901-94	Human	WB,IHC-P,FC
	CD68 Recombinant Antibody [PDM0-13] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610347	Human	IF-Tissue,IHC-P
Relative expression	CD68 Mouse mAb [A3C2]	EM1901-93	Human	WB,IHC-P,FC
	CD68 Recombinant Antibody [PDM0-13] - Rat IgG1 (Chimeric)	HA601520	Human	IF-Tissue,IHC-P
	CD68 Mouse mAb [F3-D7]	EM1706-11	Human	WB,IHC-P,IF-Cell,FC
	CD68 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS264RT	Human	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	CD90 / THY1 Recombinant Rabbit mAb [JF10-09]	ET1702-92	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
	CD90 / THY1 Recombinant Rabbit mAb [SU35-07]	ET1608-46	Human, Rat	WB,IHC-P
	CD90 / THY1 Recombinant Rabbit mAb [SU35-07] - BSA and Azide free	HA750149	Human, Rat	WB,IHC-P
	CD90 / THY1 Rabbit pAb	R1107-9	Human, Mouse, Rat	WB
KD	CD90 / THY1 Recombinant Rabbit mAb [JF10-09] - BSA and Azide free	HA750364	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
	c-Fos Recombinant Rabbit mAb [PSH07-51]	HA722839	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IF-Cell,IHC-P,ChIP,IF-Tissue,IP
Cell treatment	c-Fos Recombinant Antibody [PSH05-77] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610220	Mouse, Rat, Human, Cynomolgus monkey, Pig	IHC-Fr,WB
	c-Fos Mouse mAb [A10-G2]	EM1710-29	Human	WB,IHC-P,FC
Cell treatment	c-Fos Recombinant Rabbit mAb [JJ0938] - BSA and Azide free	HA750329	Human, Rat	WB,IF-Cell,IHC-P
	c-Fos Recombinant Rabbit mAb [PSH07-51] - BSA and Azide free	HA751172	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IF-Cell,IHC-P,ChIP,IF-Tissue,IP
Cell treatment	c-Fos Recombinant Antibody [PSH05-77] - Guinea pig IgG2 (Chimeric)	HA601374	Mouse, Rat, Human, Cynomolgus monkey, Pig	IHC-Fr,WB
Cell treatment	c-Fos Recombinant Antibody [PSH07-51] - Mouse IgG1 (Chimeric)	HA601440	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Cell,WB
Cell treatment	c-Fos Recombinant Antibody [PSH05-77] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610219	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Tissue,IF-Cell,WB
Cell treatment	c-Fos Recombinant Rabbit mAb [JJ0938]	ET1701-95	Human, Rat	WB,IF-Cell,IHC-P
Cell treatment	c-Fos Recombinant Rabbit mAb [PSH05-77] - BSA and Azide free	HA751071	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,IHC-Fr,IF-Tissue
Cell treatment	c-Fos Recombinant Antibody [PSH05-77] - Rat IgG1 (Chimeric)	HA601373	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Tissue,IF-Cell,WB
	c-Fos Recombinant Rabbit mAb	IRS131RB	Mouse	mIHC
Cell treatment	c-Fos Recombinant Rabbit mAb [PSH05-77]	HA722666	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,IHC-Fr,mIHC,IF-Tissue
Cell treatment	c-Fos Recombinant Antibody [PSH07-51] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610275	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Cell,WB
	Choline Acetyltransferase Recombinant Rabbit mAb [JA67-11] - BSA and Azide free	HA750406	Human, Mouse, Rat	WB,IP,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Choline Acetyltransferase Recombinant Rabbit mAb [PSH15-05] - BSA and Azide free	HA751543	Human, Mouse, Rat	WB,IHC-Fr
	Choline Acetyltransferase Recombinant Rabbit mAb [JA67-11]	ET1704-16	Human, Mouse, Rat	WB,IP,IF-Cell,FC
	Choline Acetyltransferase Recombinant Rabbit mAb [PSH15-05]	HA723705	Human, Mouse, Rat	WB,IHC-Fr
	Choline Acetyltransferase Recombinant Rabbit mAb [PSH15-49] - BSA and Azide free	HA751564	Human, Mouse, Rat	IHC-P
	Choline Acetyltransferase Recombinant Rabbit mAb [PSH15-49]	HA723755	Human, Mouse, Rat	IHC-P
	Cleaved-Notch 2 (A1734) Rabbit pAb	ER60031	Human, Mouse, Rat	WB
	Cleaved-Notch 2 (D1733) Rabbit pAb	ER60032	Human, Mouse, Rat	IF-Cell,WB,IHC-P
	Cleaved-Notch 2 (V1697) Rabbit pAb	ER60033	Human, Mouse, Rat	WB
Relative expression	CNPase Recombinant Rabbit mAb [JF10-25]	ET1702-46	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IF-Tissue
	CNTF Rabbit pAb	ER60106	Human, Rat, Mouse	WB
	CREB Recombinant Rabbit mAb [PSH15-48] - BSA and Azide free	HA751563	Human, Mouse, Rat, Monkey	WB,IHC-P,IHC-Fr,IF-Cell,FC
	CREB Recombinant Rabbit mAb [PSH15-48]	HA723754	Human, Mouse, Rat, Monkey	WB,IHC-P,IHC-Fr,IF-Cell,FC
	CREB Recombinant Rabbit mAb [SA04-04]	ET1601-15	Human, Mouse, Zebrafish, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	CREB Recombinant Rabbit mAb [SA04-04] - BSA and Azide free	HA750007	Human, Mouse, Zebrafish, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	CREB1 Rabbit pAb	R1309-14	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	CXCR4 Recombinant Rabbit mAb [JE03-35]	HA722304	Human	WB,IF-Cell,IHC-P
	CXCR4 Rabbit pAb	ER1802-28	Human, Rat	WB,IF-Cell,FC
KD	Cystatin C Recombinant Rabbit mAb [JJ09-16]	ET1701-72	Human, Mouse, Rat	WB,IHC-P
	Cystatin C Rabbit pAb	0808-1	Human, Zebrafish	WB,IF-Cell
Relative expression	Cystatin C Recombinant Rabbit mAb [PSH18-60]	HA724013	Human, Mouse, Rat	WB,IHC-P
Relative expression	Delta-catenin Recombinant Rabbit mAb [JE75-34]	HA722438	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IP,IHC-Fr
	Delta-catenin Rabbit pAb	HA500428	Human, Mouse, Rat	WB
	Dlx-2 Rabbit pAb	ER65853	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	Doublecortin Recombinant Antibody [JJ0959] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610326	Human, Mouse, Rat	IHC-Fr,WB
	Doublecortin Rabbit pAb	IRS252RB	Mouse	mIHC
Relative expression	Doublecortin Recombinant Antibody [JJ0959] - Rat IgG1 (Chimeric)	HA601398	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
Relative expression	Doublecortin Recombinant Antibody [JJ0959] - Guinea pig IgG2 (Chimeric)	HA601496	Human, Mouse, Rat	IHC-Fr,WB
	Doublecortin Recombinant Rabbit mAb [JJ0959]	ET1701-98	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
Relative expression	Doublecortin Recombinant Antibody [JJ0959] - Mouse IgG1 (Chimeric)	HA601497	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Doublecortin Recombinant Antibody [JJ0959] - Chicken IgY (Chimeric)	HA601495	Human, Mouse, Rat	IHC-Fr,WB
Relative expression	Doublecortin Recombinant Antibody [JJ0959] - Chicken IgY (Chimeric) - BSA and Azide free	HA610325	Human, Mouse, Rat	IHC-Fr,WB
Relative expression	Doublecortin Recombinant Antibody [JJ0959] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610238	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
Relative expression	Doublecortin Recombinant Antibody [JJ0959] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610327	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Cell treatment	EAAT1 Recombinant Rabbit mAb [JA30-35] - BSA and Azide free	HA750415	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
Cell treatment	EAAT1 Recombinant Rabbit mAb [JA30-35]	ET1704-54	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression	EAAT2 Recombinant Antibody [PS01-62] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610334	Mouse, Rat	IHC-Fr,WB
Relative expression	EAAT2 Recombinant Rabbit mAb [PS01-62] - BSA and Azide free	HA750609	Mouse, Rat, Human	WB,IHC-P,IHC-Fr
Relative expression	EAAT2 Recombinant Rabbit mAb [PS01-62]	HA721323	Mouse, Rat, Human	WB,IHC-P,IHC-Fr
	EAAT2 Recombinant Antibody [PS01-62] - Chicken IgY (Chimeric) - BSA and Azide free	HA610333	Mouse	IHC-Fr
	EAAT2 Rabbit pAb	ER60109	Human, Rat, Mouse	WB
Relative expression	EAAT2 Recombinant Rabbit mAb [PSH07-43]	HA722831	Human, Mouse, Rat	WB,IHC-P,IHC-Fr
Relative expression	EAAT2 Recombinant Antibody [PS01-62] - Guinea pig IgG2 (Chimeric)	HA601504	Mouse, Rat	IHC-Fr,WB
	EAAT2 Recombinant Antibody [PS01-62] - Chicken IgY (Chimeric)	HA601503	Mouse	IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	EAAT2 Recombinant Antibody [PS01-62] - Rat IgG1 (Chimeric)	HA601403	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
Relative expression	EAAT2 Recombinant Rabbit mAb [PSH07-43] - BSA and Azide free	HA751164	Human, Mouse, Rat	WB,IHC-P,IHC-Fr
Relative expression	EAAT2 Recombinant Antibody [PS01-62] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610243	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
	E-Cadherin Rabbit pAb	0407-25	Human, Mouse	WB,IF-Cell,FC
KD	E-Cadherin Mouse mAb [A0-G11-2]	EM0502	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	E-Cadherin Recombinant Rabbit mAb	IRS030	Human	mIHC
	E-Cadherin Recombinant Mouse mAb [A0-G11-2-R] - BSA and Azide free	HA610047	Human, Mouse, Rat	WB,IHC
Relative expression, KD	E-Cadherin Recombinant Rabbit mAb [SY0287]	ET1607-75	Human, Cynomolgus monkey	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP,mIHC
Relative expression	E-Cadherin Rabbit pAb	ER63312	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IHC-Fr,FC
	E-Cadherin Recombinant Mouse mAb	IRS229MS	Mouse	mIHC
	E-Cadherin Mouse mAb [3-F9]	M1405-3	Human, Mouse	IHC-P
Relative expression, KD	E-Cadherin Recombinant Rabbit mAb [SY0287] - BSA and Azide free	HA750133	Human, Cynomolgus monkey	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	iFluor™ 488 Conjugated E-Cadherin Recombinant Rabbit mAb [SY0287]	HA720159F	Human	IF-Cell,FC
	iFluor™ 594 Conjugated E-Cadherin Recombinant Rabbit mAb [SY0287]	HA720164F	Human	IF-Tissue
Relative expression	E-Cadherin Recombinant Rabbit mAb [PSH13-75] - BSA and Azide free	HA751488	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Cell,FC,IP
	E-Cadherin Recombinant Mouse mAb [A0-G11-2-R]	HA601143	Human, Mouse, Rat	WB,IHC,mIHC
	E-Cadherin Recombinant Mouse mAb	IRS031	Human	mIHC
Relative expression	E-Cadherin Recombinant Rabbit mAb [PSH13-75]	HA723564	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Cell,FC,IP
	E-Cadherin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS284RT	Mouse	mIHC
	Egr1 Rabbit pAb	ER1907-98	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	EGR1 Recombinant Rabbit mAb [JE57-62]	ET7111-41	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P
	Egr-2 Rabbit pAb	ER63323	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	EGR2 Recombinant Rabbit mAb [JG78-39]	ET7108-57	Human, Mouse, Rat	WB
	EOMES Rabbit pAb	ER63369	Human, Mouse	WB
	F4/80 Recombinant Rabbit mAb	IRS052	Mouse	mIHC
	F4/80 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS143RT	Mouse	mIHC
Relative expression	F4/80 Recombinant Antibody [PSH01-87] - Mouse IgG1 (Chimeric)	HA601537	Mouse	WB,IHC-P
	F4/80 Recombinant Antibody [PSH0-85] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610319	Mouse, Rat	IHC-P
Relative expression	F4/80 Recombinant Antibody [PSH01-87] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610356	Mouse	WB,IHC-P
	F4/80 Recombinant Antibody [PSH01-87] - Rat IgG1 (Chimeric)	HA601538	Mouse	WB,IHC-Fr,IHC-P
	F4/80 Recombinant Antibody [PSH0-85] - Rat IgG1 (Chimeric)	HA601485	Mouse, Rat	IHC-P
Relative expression	F4/80 Recombinant Rabbit mAb [PSH01-87] - BSA and Azide free	HA750761	Mouse	WB,IHC-P,IF-Tissue,IHC-Fr
	F4/80 Recombinant Antibody [PSH01-87] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610357	Mouse	WB,IHC-Fr,IHC-P
Relative expression	F4/80 Recombinant Rabbit mAb [PSH0-85] - BSA and Azide free	HA750665	Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression	F4/80 Recombinant Antibody [PSH01-87] - Guinea pig IgG2 (Chimeric)	HA601567	Mouse	WB,IHC-Fr
Relative expression	F4/80 Recombinant Rabbit mAb [PSH0-85]	HA721520	Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr,mIHC
Relative expression	F4/80 Recombinant Rabbit mAb [PSH01-87]	HA721745	Mouse	WB,IHC-P,IF-Tissue,IHC-Fr
	FGF2 Recombinant Rabbit mAb [JM28-10]	ET1703-18	Human, Mouse	WB,IP,FC,IHC-P
KD	FGFR3 Recombinant Rabbit mAb [JM81-10]	ET1703-36	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	FGFR3 Recombinant Rabbit mAb [JM81-10] - BSA and Azide free	HA750379	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	FGFR3 Recombinant Rabbit mAb [JM110-33] - BSA and Azide free	HA750388	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	FGFR3 Recombinant Rabbit mAb [JM110-33]	ET1703-65	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	FGFR-4 Rabbit pAb	ER65688	Human	WB
	FGFR4 Recombinant Rabbit mAb [JE59-96]	HA721578	Human	WB,FC
KD	Fibronectin Recombinant Rabbit mAb [JF0582] - BSA and Azide free	HA750340	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
	Fibronectin Recombinant Rabbit mAb [PSH14-21]	HA723610	Human, Mouse	WB,IHC-P,IF-Cell,FC
	Fibronectin Recombinant Rabbit mAb [PSH14-22] - BSA and Azide free	HA751515	Human, Mouse	WB,IF-Cell,FC
	Fibronectin Mouse mAb [3G4]	RT1224	Human, Mouse, Rat	WB,IP,IF,IHC-P
	Fibronectin Recombinant Rabbit mAb [PSH14-22]	HA723611	Human, Mouse	WB,IF-Cell,FC
KD	Fibronectin Recombinant Rabbit mAb [JF0582]	ET1702-25	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
	Fibronectin Recombinant Rabbit mAb [PSH14-21] - BSA and Azide free	HA751514	Human, Mouse	WB,IHC-P,IF-Cell,FC
	Fibronectin Recombinant mAb	IRS144RB	Mouse	mIHC
	FOXA2 Recombinant mAb	IRS146RB	Mouse	mIHC
	FOXA2 Recombinant Rabbit mAb [JM10-64]	ET1703-76	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-P,IHC-Fr
	FoxD3 Rabbit pAb	ER63453	Human, Mouse	WB,IHC-P
	FOXD3 Rabbit pAb	R1107-4	Human	WB
	Frizzled-9 Rabbit pAb	ER63472	Human, Mouse, Monkey	WB,IF-Cell
Relative expression	GABA B Receptor 1 Recombinant Rabbit mAb [PSH06-21]	HA722653	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr
	GABA B Receptor 1 Recombinant Rabbit mAb [JE60-09]	HA721338	Human, Mouse	WB,IHC-P,IF-Tissue
Relative expression	GABA B Receptor 1 Recombinant Rabbit mAb [PSH03-87] - BSA and Azide free	HA750912	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P
Relative expression	GABA B Receptor 1 Recombinant Rabbit mAb [PSH06-21] - BSA and Azide free	HA751059	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr
Relative expression	GABA B Receptor 1 Recombinant Rabbit mAb [PSH03-87]	HA722053	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P
	GAD65 Recombinant Rabbit mAb [PSH06-80] - BSA and Azide free	HA751102	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
	GAD65 Recombinant Rabbit mAb [PSH06-80]	HA722750	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	GAD-65/67 Rabbit pAb	ER63491	Human, Mouse	WB,IHC-P,IF-Cell
	GAD67 Recombinant Rabbit mAb [JM11-11]	ET1703-71	Human, Mouse, Rat	WB,IP,IHC-P
	GAD67 Recombinant Mouse mAb [PSH08-30] - BSA and Azide free	HA610216	Mouse, Rat	WB,IHC-P,IHC-Fr
	GAD67 Recombinant Mouse mAb [PSH08-30]	HA601370	Mouse, Rat	WB,IHC-P,IHC-Fr
Relative expression	GAP43 Rabbit pAb	ER40201	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression, KD	GAP43 Recombinant Rabbit mAb [SC60-06]	ET1610-94	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	GDF-7 Rabbit pAb	ER65661	Human, Mouse, Rat	WB,IHC-P
	iFluor™ 594 Conjugated GFAP Mouse mAb [1-D4]	HA600102F	Human, Rat	IF-Tissue
	Biotin Conjugated GFAP Mouse mAb [PSH05-17]	HA601355B	Human, Mouse, Rat	ELISA(Det)
	GFAP Recombinant Mouse mAb [PSH05-17] - BSA and Azide free (Detector)	HA601342	Human, Mouse, Rat	ELISA(Det)
	iFluor™ 488 Conjugated GFAP Mouse mAb [1-D4]	HA600101F	Human, Rat	IF-Tissue,IHC-Fr
	GFAP Recombinant Rabbit mAb	IRS069	Mouse	mIHC
	GFAP Rabbit pAb	R1308-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	GFAP Recombinant Mouse mAb [PSH10-24] - BSA and Azide free	HA610225	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue
	GFAP Recombinant Rabbit mAb [PSH05-16] - BSA and Azide free (Capture)	HA722234	Human, Mouse, Rat	ELISA(Cap)
	GFAP Recombinant Antibody - Rat IgG1 (Chimeric)	IRS148RT	Mouse	mIHC
	GFAP Mouse mAb [1-D4]	EM140707	Human, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression	GFAP Recombinant Antibody [SA03-04] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610317	Mouse, Rat	IHC-Fr,WB
Relative expression	GFAP Recombinant Mouse mAb [PSH10-24]	HA601380	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue
Relative expression	GFAP Recombinant Antibody [SA03-04] - Rat IgG1 (Chimeric)	HA601483	Mouse, Rat	IHC-Fr,WB
	GFAP Rabbit pAb	R1308-10	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	iFluor™ 647 Conjugated GFAP Mouse mAb [1-D4]	HA600103F	Human, Rat	IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	GFAP Recombinant Rabbit mAb [SA03-04]	ET1601-23	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr,mIHC,IP
Relative expression	GIRK2 Recombinant Rabbit mAb [PSH05-91]	HA722478	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IP,IHC-Fr
Relative expression	GIRK2 Recombinant Rabbit mAb [PSH05-91] - BSA and Azide free	HA751027	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IP,IHC-Fr
	Glutamine Synthetase Recombinant Antibody [PSH10-48] - Rat IgG1 (Chimeric)	HA601428	Human, Mouse, Rat	IHC-Fr,WB
	Glutamine Synthetase Recombinant Rabbit mAb [PSH10-48] - BSA and Azide free	HA751351	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr
	Glutamine Synthetase Recombinant Mouse mAb	IRS231MS	Mouse	mIHC
	Glutamine Synthetase Recombinant Antibody [PSH10-48] - Chicken IgY (Chimeric)	HA601516	Mouse, Rat	IHC-Fr,WB
	Glutamine Synthetase Mouse mAb [A3G2]	EM1902-39	Human, Mouse, Rat	WB,IHC-P,FC,mIHC,IF-Tissue
	Glutamine Synthetase Recombinant Mouse mAb [A3G2-R] - BSA and Azide free	HA610100	Human, Mouse, Rat	WB,IHC-P
	Glutamine Synthetase Recombinant Rabbit mAb [PSH10-48]	HA723212	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr
	Glutamine Synthetase Recombinant Antibody [PSH10-48] - Chicken IgY (Chimeric) - BSA and Azide free	HA610346	Mouse, Rat	IHC-Fr,WB
	Glutamine Synthetase Recombinant Mouse mAb [A3G2-R]	HA601225	Human, Mouse, Rat	WB,IHC-P
	Glutamine Synthetase Recombinant Antibody [PSH10-48] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610268	Human, Mouse, Rat	IHC-Fr,WB
	Glutamine Synthetase Mouse mAb	IRS033	Human	mIHC
	HB9/HLXB9/MNX1 Recombinant Rabbit mAb [JE32-08]	HA721410	Human, Mouse, Rat	WB
	Hes1 Recombinant Rabbit mAb [SC06-21] - BSA and Azide free	HA750240	Human, Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC
	Hes1 Recombinant Rabbit mAb [SC06-21]	ET1610-97	Human, Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC
	HES5 Rabbit pAb	ER1910-51	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	HES5 Recombinant Rabbit mAb [JE35-96]	HA722699	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P
	HMGA2 Rabbit pAb	ER62068	Human, Mouse	WB
	HMGA2 Recombinant Rabbit mAb [PSH05-53]	HA722295	Human, Mouse	WB,IF-Cell,IHC-P,FC,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HMGA2 Recombinant Rabbit mAb [PSH05-53] - BSA and Azide free	HA751012	Human, Mouse	WB,IF-Cell,IHC-P,FC,IF-Tissue
	HTR3A Mouse mAb [3B11C10]	EM1710-35	Human	IHC-P,IF-Cell,FC
	FITC Conjugated Human CD11b Recombinant Mouse mAb [A2-9-R]	HA600130F	Human	FC
	FITC Conjugated Human CD68 Recombinant Mouse mAb [A3C3-R]	HA600132F	Human	FC
	Human CNTF Recombinant Rabbit mAb [PSH11-82] - BSA and Azide free (Detector)	HA723385	Human, Rat	ELISA(Det)
	Biotin Conjugated Human CNTF Recombinant Rabbit mAb [PSH11-80] - Detector	HA723383B	Human, Rat	ELISA(Det),ELISA
	Biotin Conjugated Human CNTF Recombinant Rabbit mAb [PSH11-82] - Detector	HA723386B	Human, Rat	ELISA(Det),ELISA
	Human CNTF Recombinant Rabbit mAb [PSH11-81] - BSA and Azide free (Capture)	HA723384	Human, Rat	ELISA(Cap)
	Human CNTF Recombinant Rabbit mAb [PSH11-80] - BSA and Azide free (Detector)	HA723382	Human, Rat	ELISA(Det)
	Human CNTF Recombinant Rabbit mAb [PSH11-79] - BSA and Azide free (Capture)	HA723381	Human, Rat	ELISA(Cap)
	Biotin Conjugated Human CSPG4 Recombinant Rabbit mAb [PSH10-79] - Detector	HA723247B	Human	ELISA(Det),ELISA
	HRP Conjugated Human CSPG4 Recombinant Rabbit mAb [PSH10-79] - Detector	HA723376H	Human	ELISA(Det)
	Human CSPG4 Recombinant Rabbit mAb [PSH10-78] - BSA and Azide free (Capture)	HA723245	Human	ELISA(Cap)
	Human CSPG4 Recombinant Rabbit mAb [PSH10-79] - BSA and Azide free (Detector)	HA723246	Human	ELISA(Det)
	Human CSPG4 Recombinant Rabbit mAb [PSH10-77] - BSA and Azide free (Capture)	HA723244	Human	ELISA(Cap)
	Biotin Conjugated Human E-Cadherin Recombinant Rabbit mAb [SY0287] - BSA and Azide free (Detector)	HA723057B	Human	ELISA(Det),ELISA
	Human E-Cadherin Recombinant Rabbit mAb [SY0287] - BSA and Azide free (Detector)	HA723056	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human E-Cadherin Recombinant Rabbit mAb [PSH09-02] - BSA and Azide free (Capture)	HA723055	Human	ELISA(Cap)
	Human FGF-2 Recombinant Rabbit mAb [PSH02-09] - BSA and Azide free (Detector)	HA721781	Human	ELISA(Det)
	Biotin Conjugated Human FGF-2 Recombinant Rabbit mAb [PSH02-09]	HA722639B	Human	ELISA(Det)
	Human FGF-2 Recombinant Rabbit mAb [PSH02-08] - BSA and Azide free (Capture)	HA721780	Human	ELISA(Cap)
	Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-42] - BSA and Azide free (Detector)	HA723331	Human, Mouse, Rat	ELISA(Det)
	Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-44] - BSA and Azide free (Detector)	HA723334	Human, Mouse, Rat	ELISA(Det)
	Biotin Conjugated Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-44] - Detector	HA723335B	Human, Mouse, Rat	ELISA(Det),ELISA
	Biotin Conjugated Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-42] - Detector	HA723332B	Human, Mouse, Rat	ELISA(Det),ELISA
	Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-41] - BSA and Azide free (Capture)	HA723330	Human, Mouse, Rat	ELISA(Cap)
	Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-43] - BSA and Azide free (Capture)	HA723333	Human, Mouse, Rat	ELISA(Cap)
	Human NCAM1 Recombinant Rabbit mAb [PSH08-08] - BSA and Azide free (Capture/Detector)	HA722947	Human	ELISA(Det),ELISA(Cap)
	Human NCAM1 Recombinant Rabbit mAb [PSH08-09] - BSA and Azide free (Detector)	HA722949	Human	ELISA(Det)
	Biotin Conjugated Human NCAM1 Recombinant Rabbit mAb [PSH08-09] - Detector	HA722950B	Human	ELISA(Det)
	Biotin Conjugated Human NCAM1 Recombinant Rabbit mAb [PSH08-08] - Detector	HA722948B	Human	ELISA(Det)
	Human NCAM1 Recombinant Rabbit mAb [PSH08-07] - BSA and Azide free (Capture)	HA722946	Human	ELISA(Cap)
	PE Conjugated Human NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53]	HA720215F	Human	FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	APC Conjugated Human NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53]	HA723757F	Human	FC
	Human PCSK9 Recombinant Rabbit mAb [PSH12-85] - BSA and Azide free (Capture)	HA725124	Human	ELISA(Cap)
	HRP Conjugated Human PCSK9 Recombinant Rabbit mAb [PSH12-86] - Detector	HA725126H	Human	ELISA(Det),ELISA
	Biotin Conjugated Human PCSK9 Recombinant Rabbit mAb [PSH12-86] - Detector	HA725126B	Human	ELISA(Det),ELISA
	Human PCSK9 Recombinant Rabbit mAb [PSH12-86] - BSA and Azide free (Detector)	HA725125	Human	ELISA(Det)
	Human VE Cadherin Recombinant Rabbit mAb [PSH02-61] - BSA and Azide free (Capture)	HA721838	Human	ELISA(Cap)
	Human VE Cadherin Recombinant Rabbit mAb [PSH02-62] - BSA and Azide free (Detector)	HA721839	Human	ELISA(Det)
	Biotin Conjugated Human VE Cadherin Recombinant Rabbit mAb [PSH02-62] - Detector	HA722635B	Human	ELISA(Det)
Relative expression	Iba1 Recombinant Antibody [JM36-62] - Chicken IgY (Chimeric) - BSA and Azide free	HA610221	Human, Mouse, Rat, Pig, Cynomolgus monkey	IHC-Fr,WB,FC
Relative expression	Iba1 Recombinant Rabbit mAb [JM36-62]	ET1705-78	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC,IP,mIHC
	Iba1 Recombinant Antibody [JM36-62] - Rat IgG1 (Chimeric)	HA601367	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
Relative expression	Iba1 Recombinant Rabbit mAb [JM36-62] - BSA and Azide free	HA750440	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC,IP
	Iba1 Rabbit pAb	ER1802-20	Human, Mouse	IHC-P
	iFluor™ 488 Conjugated Iba1 Recombinant Rabbit mAb [JM36-62]	HA720158F	Human, Mouse, Rat	IF-Cell,FC
	Iba1 Recombinant Mouse mAb [8G8]	RT1316	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,FC,IP,IHC-Fr,IF-Cell,IF-Tissue
	Iba1 Recombinant Rabbit mAb	IRS068	Mouse	mIHC
	Iba1 Recombinant Antibody [JM36-62] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610214	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,WB
	Iba1 Recombinant Antibody [JM36-62] - Guinea pig IgG2 (Chimeric)	HA601368	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,WB
Relative expression	Iba1 Recombinant Antibody [JM36-62] - Chicken IgY (Chimeric)	HA601376	Human, Mouse, Rat, Pig, Cynomolgus monkey	IHC-Fr,WB,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	iFluor™ 594 Conjugated Iba1 Recombinant Rabbit mAb [JM36-62]	HA723821F	Mouse	IHC-Fr
	Iba1 Recombinant Mouse mAb [8G8] - BSA and Azide free	HA610321	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,FC,IP,IHC-Fr,IF-Cell,IF-Tissue
	Iba1 Recombinant Antibody [JM36-62] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610213	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
	Iba1 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS158RT	Mouse	mIHC
Cell treatment	Integrin alpha 4 Recombinant Rabbit mAb [JA09-36]	ET1704-31	Human, Mouse, Rat	WB,IHC-P,FC
	Integrin alpha 4 Rabbit pAb	ER65508	Human, Mouse, Rat	WB,IHC-P
	IPMK Rabbit pAb	ER63775	Human, Mouse, Rat, Monkey	WB,IHC-P
	Islet 1 Recombinant Rabbit mAb [SC05-64] - BSA and Azide free	HA750211	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-P,FC,IP,IHC-Fr
	Islet 1 Recombinant Rabbit mAb [SC05-64]	ET1610-33	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-P,FC,IP,IHC-Fr
	Islet-2 Rabbit pAb	ER63782	Human, Mouse, Rat	WB,IHC
	Ki67 Recombinant Rabbit mAb	IRS036	Human	mIHC
	iFluor™ 647 Conjugated Ki67 Recombinant Rabbit mAb [SR00-02]	HA720163F	Human, Mouse, Rat	IF-Cell,IF-Tissue,FC,IHC-Fr
	Ki67 Recombinant Rabbit mAb [ST50-01]	ET1609-34	Human	WB,IF-Cell,IHC-P,FC
	Ki67 Rabbit pAb	ER1902-75	Human	WB,IF-Cell,FC
	Ki67 Rabbit pAb	RT1349	Human	WB,IP,IF,IHC-P
	iFluor™ 488 Conjugated Ki67 Recombinant Rabbit mAb [SR00-02]	HA720162F	Human	IF-Cell,IF-Tissue,FC
	Ki67 Recombinant Rabbit mAb [SR00-02]	HA721115	Human, Mouse, Rat, Cynomolgus monkey	IHC-P,IF-Tissue,mIHC,IF-Cell,FC
	Ki67 Recombinant Mouse mAb [PD00-10]	HA601053	Human, Mouse	IHC-P,IF-Tissue
	Ki67 Rabbit pAb	ER1706-46	Human	WB,IF-Cell,IHC-P,FC
	Ki67 Recombinant Rabbit mAb [PSH0-02]	HA721254	Human	IHC-P,FC
	Ki67 Recombinant Rabbit mAb [PSH0-02] - BSA and Azide free	HA750570	Human	IHC-P,FC
	Ki67 Recombinant Rabbit mAb	IRS057	Mouse	mIHC
	Lhx2 Recombinant Rabbit mAb [JE54-99]	HA721900	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	LMX1B Rabbit pAb	ER63864	Human, Mouse	WB,IHC-P
Relative expression	MAG Recombinant Rabbit mAb [PSH02-41]	HA721818	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IP,IHC-Fr
Relative expression	MAG Recombinant Rabbit mAb [PSH02-41] - BSA and Azide free	HA750806	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IP,IHC-Fr
Relative expression	MAG Recombinant Rabbit mAb [PSH02-42]	HA721819	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	MAG Recombinant Antibody [PSH02-41] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610254	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB,IF-Tissue
Relative expression	MAG Recombinant Rabbit mAb [PSH02-42] - BSA and Azide free	HA750807	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	MAG Rabbit pAb	ER1912-62	Human, Rat	WB
Relative expression	MAG Recombinant Antibody [PSH02-41] - Rat IgG1 (Chimeric)	HA601414	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB,IF-Tissue
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610293	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell
Relative expression	MAP2 Recombinant Rabbit mAb [PSH08-73] - BSA and Azide free	HA751238	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IHC-Fr,IF-Cell
Relative expression	MAP2 Recombinant Rabbit mAb [PSH08-73]	HA723025	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IHC-Fr,IF-Cell
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Chicken IgY (Chimeric) - BSA and Azide free	HA610301	Mouse, Rat	IHC-Fr,WB
	MAP2 Rabbit pAb	HA500177	Human, Mouse, Rat	WB,IF-Cell,IHC-P,mIHC
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610292	Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Chicken IgY (Chimeric)	HA601467	Mouse, Rat	IHC-Fr,WB
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Rat IgG1 (Chimeric)	HA601459	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Mouse IgG1 (Chimeric)	HA601458	Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610291	Mouse, Rat	IHC-Fr,WB
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Guinea pig IgG2 (Chimeric)	HA601457	Mouse, Rat	IHC-Fr,WB
	MAP2 Recombinant Rabbit mAb [PSH08-73]	IRS071	Mouse	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	MOG Recombinant Rabbit mAb [JM23-06]	ET1705-16	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
	Biotin Conjugated Mouse CD68 Recombinant Rabbit mAb [PSH12-25] - Detector	HA723448B	Mouse	ELISA(Det),ELISA
	Mouse CD68 Recombinant Rabbit mAb [PSH12-26] - BSA and Azide free (Capture)	HA723449	Mouse	ELISA(Cap)
	Mouse CD68 Recombinant Rabbit mAb [PSH12-25] - BSA and Azide free (Detector)	HA723447	Mouse	ELISA(Det)
	Biotin Conjugated Mouse CD68 Recombinant Rabbit mAb [PSH12-27] - Detector	HA723451B	Mouse	ELISA(Det),ELISA
	Mouse CD68 Recombinant Rabbit mAb [PSH12-24] - BSA and Azide free (Capture)	HA723446	Mouse	ELISA(Cap)
	Mouse CD68 Recombinant Rabbit mAb [PSH12-27] - BSA and Azide free (Detector)	HA723450	Mouse	ELISA(Det)
	Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-25] - BSA and Azide free (Detector)	HA723912	Mouse, Rat	ELISA(Det)
	Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-24] - BSA and Azide free (Capture)	HA723911	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-25] - Detector	HA723912B	Mouse, Rat	ELISA(Det),ELISA
	HRP Conjugated Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-25] - Detector	HA723912H	Mouse, Rat	ELISA(Det),ELISA
	MSX1 Mouse mAb [D11-A8]	EM1710-19	Human	WB
Relative expression	MTAP Recombinant Mouse mAb [PSH13-98]	HA601445	Human, Mouse, Rat	WB,IF-Cell
	MTAP Rabbit pAb	ER1901-47	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	MTAP Recombinant Rabbit mAb [JE63-74]	HA721955	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,FC
	Musashi 1 Recombinant Rabbit mAb [SJ201]	ET1606-51	Human, Mouse, Pig, Cynomolgus monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Relative expression	Myelin Basic Protein Recombinant Rabbit mAb [JF0943]	ET1702-15	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr,mIHC
Relative expression	Myelin Basic Protein Recombinant Rabbit mAb [JF0943] - BSA and Azide free	HA750337	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
	Myelin Basic Protein Recombinant Rabbit mAb	IRS064	Mouse	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression, KD	N Cadherin Recombinant Rabbit mAb [SY02-46]	ET1607-37	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
Relative expression, KD	N Cadherin Recombinant Rabbit mAb [SY02-46] - BSA and Azide free	HA750119	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	N Cadherin Recombinant Antibody [SY02-46] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610298	Human, Mouse, Rat	IHC-Fr,IHC-P
	N Cadherin Recombinant Antibody [SY02-46] - Mouse IgG1 (Chimeric)	HA601464	Human, Mouse, Rat	IHC-Fr,IHC-P
	N-Cadherin Mouse mAb [A10-4]	M1304-1	Human, Mouse	WB,IHC-P
	N-Cadherin Rabbit pAb	ER0503	Human, Mouse, Rat	WB,IHC,IF
	NCAM1 / CD56 Rabbit pAb	R1204-1	Human, Mouse, Rat	WB,IF-Cell
	NCAM1 / CD56 Recombinant Rabbit mAb [JF1021] - BSA and Azide free	HA750346	Human, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53] - BSA and Azide free	HA751105	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
Relative expression	NCAM1 / CD56 Recombinant Mouse mAb [PDH0-07]	HA601126	Human	IHC-P,WB
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH08-89]	HA723041	Human	WB,IHC-P,IF-Tissue
	NCAM1 / CD56 Recombinant Rabbit mAb	IRS026	Human	mIHC
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53]	HA722756	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-82] - BSA and Azide free	HA751104	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-82]	HA722755	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue
	NCAM1 / CD56 Recombinant Rabbit mAb [JF1021]	ET1702-43	Human, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP,mIHC
	NEDD1 Rabbit pAb	HA500294	Human, Mouse	WB,IHC-P,IF-Cell,FC
KD	Nestin Recombinant Rabbit mAb [SN06-27]	ET1611-21	Human	IF-Cell,IF-Tissue,IHC-P,WB
	Nestin Recombinant Antibody [PSH07-89] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610246	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell
Relative expression, KD	Nestin Recombinant Rabbit mAb [PSH07-89] - BSA and Azide free	HA751196	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IF-Tissue,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression, KD	Nestin Recombinant Rabbit mAb [PSH07-89]	HA722919	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IF-Tissue,mIHC,IP
Relative expression	Nestin Recombinant Antibody [PSH07-89] - Guinea pig IgG2 (Chimeric)	HA601468	Human, Mouse, Rat	IHC-Fr,WB
Relative expression	Nestin Recombinant Antibody [PSH07-89] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610302	Human, Mouse, Rat	IHC-Fr,WB
	Nestin Recombinant Rabbit mAb	IRS065	Mouse	mIHC
	Nestin Recombinant Antibody [PSH07-89] - Rat IgG1 (Chimeric)	HA601406	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell
	Nestin Rabbit pAb	R1510-19	Human, Mouse	IHC-P
	NeuN Recombinant Mouse mAb [PD01-45] - BSA and Azide free	HA610031	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IHC-Fr
	NeuN Recombinant Antibody [SR45-07] - Chicken IgY (Chimeric) - BSA and Azide free	HA610314	Mouse	IHC-Fr
Relative expression	NeuN Recombinant Rabbit mAb [SR45-07]	ET1602-12	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr,mIHC
	NeuN Recombinant Antibody [SR45-07] - Chicken IgY (Chimeric)	HA601480	Mouse	IHC-Fr
Relative expression	NeuN Recombinant Antibody [SR45-07] - Guinea pig IgG2 (Chimeric)	HA601481	Mouse, Rat	IHC-Fr,WB
Relative expression	NeuN Recombinant Antibody [SR45-07] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610316	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	NeuN Recombinant Antibody [SR45-07] - Rat IgG1 (Chimeric)	HA601397	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
Relative expression	NeuN Recombinant Antibody [SR45-07] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610315	Mouse, Rat	IHC-Fr,WB
	NeuN Recombinant Antibody - Rat IgG1 (Chimeric)	IRS179RT	Mouse	mIHC
Relative expression	NeuN Recombinant Antibody [SR45-07] - Mouse IgG1 (Chimeric)	HA601482	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	iFluor™ 488 Conjugated NeuN Recombinant Rabbit mAb [SR45-07]	HA720168F	Human, Mouse, Rat	IF-Tissue
	NeuN Recombinant Mouse mAb [PD01-45]	HA601111	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IHC-Fr
	NeuN Recombinant Rabbit mAb	IRS178RB	Mouse	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	NeuN Recombinant Antibody [SR45-07] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610237	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
	iFluor™ 594 Conjugated NeuN Recombinant Rabbit mAb [SR45-07]	HA720167F	Human, Mouse, Rat	IF-Tissue,IHC-Fr
	NeuroD1 Recombinant Rabbit mAb [JM11-10]	ET1703-73	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IP,IHC-P,IHC-Fr
	NeuroD2 Recombinant Rabbit mAb [JE50-57]	ET7110-17	Mouse, Rat, Human	IHC-P,WB
	Neurofilament Rabbit pAb	R1706-17	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Neuroglycan C Rabbit pAb	ER1913-72	Human	WB
	NF-H Mouse mAb [2G1]	RT1424	Human, Mouse, Rat	WB,IP,IF,IHC-P
Relative expression	NF-H Recombinant Rabbit mAb [JF99-07] - BSA and Azide free	HA750357	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
Relative expression	NF-H Recombinant Rabbit mAb [PSH18-14]	HA723977	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	NF-H Recombinant Rabbit mAb [JF99-07]	ET1702-72	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
Relative expression	NF-L Recombinant Antibody [PS02-10] - Chicken IgY (Chimeric)	HA601472	Mouse, Rat	IHC-Fr,WB
Relative expression	NF-L Recombinant Rabbit mAb [PS02-10]	HA721538	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,IF-Tissue,IHC-Fr
Relative expression	NF-L Recombinant Antibody [PS02-10] - Chicken IgY (Chimeric) - BSA and Azide free	HA610306	Mouse, Rat	IHC-Fr,WB
	NF-L Recombinant Antibody [PS02-10] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610297	Mouse, Rat	IHC-Fr
Relative expression	NF-L Recombinant Rabbit mAb [PSH0-77] - BSA and Azide free	HA750657	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
	NF-L Rabbit pAb	ER65439	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	NF-L Recombinant Antibody [PS02-10] - Guinea pig IgG2 (Chimeric)	HA601463	Mouse, Rat	IHC-Fr
Relative expression	NF-L Recombinant Rabbit mAb [PSH0-77]	HA721498	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
Relative expression	NF-L Recombinant Rabbit mAb [PS02-10] - BSA and Azide free	HA750678	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,IF-Tissue,IHC-Fr
Relative expression	NF-L Recombinant Antibody [PS02-10] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610249	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
Relative expression	NF-L Recombinant Antibody [PS02-10] - Rat IgG1 (Chimeric)	HA601409	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	NF-L Recombinant Antibody [PS02-10] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610282	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	NF-L Recombinant Antibody [PS02-10] - Mouse IgG1 (Chimeric)	HA601448	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	NF-M Recombinant Antibody [JM11-20] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610279	Mouse, Rat	IHC-Fr,WB
	NF-M Recombinant Antibody [JM11-20] - Chicken IgY (Chimeric)	HA601507	Mouse	IHC-Fr
	NF-M Recombinant Antibody [JM11-20] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610248	Human, Mouse, Rat, Cynomolgus monkey	IHC-Fr,IHC-P,WB
	NF-M Recombinant Antibody [JM11-20] - Mouse IgG1 (Chimeric)	HA601443	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	NF-M Recombinant Antibody [JM11-20] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610278	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	NF-M Recombinant Antibody [JM11-20] - Guinea pig IgG2 (Chimeric)	HA601444	Mouse, Rat	IHC-Fr,WB
Relative expression	NF-M Recombinant Rabbit mAb [JM11-20]	ET1703-57	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IHC-Fr,IF-Tissue
	NF-M Recombinant Antibody [JM11-20] - Chicken IgY (Chimeric) - BSA and Azide free	HA610337	Mouse	IHC-Fr
	NF-M Recombinant Antibody [JM11-20] - Rat IgG1 (Chimeric)	HA601408	Human, Mouse, Rat, Cynomolgus monkey	IHC-Fr,IHC-P,WB
Relative expression	NG2 Recombinant Rabbit mAb [PSH09-98]	HA723167	Human	WB,IF-Cell,FC
Relative expression	NG2 Recombinant Rabbit mAb [JM10-13] - BSA and Azide free	HA750372	Human, Mouse, Rat	WB,IHC-P
Relative expression	NG2 Recombinant Rabbit mAb [PSH09-98] - BSA and Azide free	HA751329	Human	WB,IF-Cell,FC
Relative expression	NG2 Recombinant Rabbit mAb [JM10-13]	ET1703-16	Human, Mouse, Rat	WB,IHC-P
	NGF Receptr Mouse mAb [2G2]	RT1429	Human, Mouse, Rat	WB,IP,IF
	NGN1 Rabbit pAb	ER60918	Human	WB
	NGN2 Rabbit pAb	ER61851	Human, Mouse	WB
Relative expression	NMDAR1 Recombinant Rabbit mAb [JM11-26]	ET1703-75	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,IHC-Fr
Relative expression	NMDAR1 Recombinant Rabbit mAb [JM11-26] - BSA and Azide free	HA750391	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	NMDAR2B Recombinant Rabbit mAb [PSH05-46] - BSA and Azide free	HA751005	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
	NMDAR2B Rabbit pAb	HA500208	Mouse, Rat, Human	WB,IF-Cell,IHC-P,FC
	NMDAR2B Rabbit pAb	ER1913-90	Human, Rat	WB,IHC-P,IF-Tissue
Relative expression	NMDAR2B Recombinant Rabbit mAb [PSH05-46]	HA722284	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
	Biotin Conjugated Non-Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-45] - Detector	HA723846B	Human	ELISA(Det),ELISA
	Non-Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-43] - BSA and Azide free (Capture)	HA723843	Human, Mouse, Rat	ELISA(Cap)
	Non-Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-45] - BSA and Azide free (Detector)	HA723845	Human	ELISA(Det)
	Non-Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-44] - BSA and Azide free (Capture)	HA723844	Human	ELISA(Cap)
KD	Notch 1 Recombinant Rabbit mAb [SJ205]	ET1606-55	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	Notch 1 Recombinant Rabbit mAb [SJ205] - BSA and Azide free	HA750111	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Notch2 Recombinant Rabbit mAb [PSH0-54] - BSA and Azide free	HA750631	Human, Mouse, Rat	WB
	Notch2 Recombinant Rabbit mAb [PSH0-55]	HA721402	Human, Mouse, Rat	WB
	Notch2 Recombinant Rabbit mAb [PSH0-54]	HA721401	Human, Mouse, Rat	WB
	Notch2 Recombinant Rabbit mAb [PSH0-55] - BSA and Azide free	HA750632	Human, Mouse, Rat	WB
	NOTCH3 Rabbit pAb	HA500519	Human	WB,FC
	NR2E1 Rabbit pAb	ER61909	Human, Mouse	WB
	NR2E1 Recombinant Rabbit mAb [JE38-26]	HA722548	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Nucleostemin Recombinant Rabbit mAb [JG37-82]	ET7108-55	Human	WB,IF-Cell,IHC-P,FC
	Olig2 Recombinant Rabbit mAb	IRS067	Mouse	mIHC
	Olig2 Recombinant Antibody [SP07-02] - Rat IgG1 (Chimeric)	HA601558	Mouse, Rat	IHC-Fr,IHC-P
Relative expression	Olig2 Recombinant Rabbit mAb [SP07-02]	ET1604-29	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,mIHC,IF-Cell
	p75 NGF Receptor Recombinant Antibody [PSH13-26] - Rabbit IgG (Chimeric)	HA723532	Mouse	IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	iFluor™ 488 Conjugated p75 NGF Receptor Recombinant Rabbit mAb [SA39-02]	HA720190F	Human, Mouse, Rat	IF-Cell
KD	p75 NGF Receptor Recombinant Rabbit mAb [SA39-02] - BSA and Azide free	HA750014	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	p75 NGF Receptor Recombinant Rabbit mAb [SA39-02]	ET1601-22	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
	p75 NGF Receptor Recombinant Antibody [PSH13-26] - Rabbit IgG (Chimeric) - BSA and Azide free	HA751462	Mouse	IHC-Fr
	p75 NGF Receptor Recombinant Rat mAb [PSH13-26] - BSA and Azide free	HA610274	Mouse, Rat	IHC-Fr,WB
	p75 NGF Receptor Recombinant Rat mAb [PSH13-26]	HA601436	Mouse, Rat	IHC-Fr,WB
	Pan Cadherin Recombinant Rabbit mAb	RE6045	Human, Mouse, Rat	WB,ICC,IHC-P,FC,IP
	Pan-Cadherin Recombinant Rabbit mAb	IRS185RB	Mouse	mIHC
	Pan-Cadherin Recombinant Rabbit mAb [ST54-01] - BSA and Azide free	HA750190	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,IHC-Fr
	Pan-Cadherin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS288RT	Mouse	mIHC
	Pan-Cadherin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS286MS	Mouse	mIHC
	Pan-Cadherin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS287RT	Human	mIHC
	Pan-Cadherin Recombinant Rabbit mAb [ST54-01]	ET1609-70	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,IHC-Fr
	Pan-Cadherin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS285MS	Human	mIHC
Relative expression	Parvalbumin Recombinant Rabbit mAb [JM100-08]	ET1703-15	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IP
Relative expression	Parvalbumin Recombinant Rabbit mAb [JM100-08] - BSA and Azide free	HA750371	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IP
	PAX3 Mouse mAb [G7-D8]	EM1709-18	Human, Mouse	WB
Relative expression	PAX6 Recombinant Rabbit mAb [SD08-31]	ET1612-89	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	PAX6 Recombinant Rabbit mAb	IRS066	Mouse	mIHC
Relative expression	PAX6 Recombinant Rabbit mAb [SD081-03]	ET1612-58	Human, Mouse, Rat	WB,IHC-P,FC,mIHC,IF-Cell,IF-Tissue,IHC-Fr
	PCNA Mouse mAb [4-E11]	M1405-8	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	PCNA Recombinant Mouse mAb [A6-G11-R]	HA601172	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,IP
KD	PCNA Recombinant Rabbit mAb [SY12-07]	ET1605-38	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	PCNA Mouse mAb [A6-G11]	EM111201	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
KD	PCNA Recombinant Mouse mAb [A6-G11-R] - BSA and Azide free	HA610056	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,IP
KD	PCNA Recombinant Rabbit mAb [SY12-07] - BSA and Azide free	HA750084	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	PCNA Rabbit pAb	R1306-5	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	PCSK9 Mouse mAb [2F1] - BSA and Azide free	HA610371	Human, Mouse	WB,IHC-P,IF-Cell,FC
Relative expression	PCSK9 Mouse mAb [2F1]	EM1701-26	Human, Mouse	WB,IHC-P,IF-Cell,FC
Relative expression	PCSK9 Recombinant Rabbit mAb [JE30-37]	HA722528	Human	WB,IHC-P,IP
	PENK Rabbit pAb	ER65744	Human	WB
Cell treatment	Pet1 Recombinant Rabbit mAb [PSH07-77] - BSA and Azide free	HA751189	Human	WB,IF-Cell,IHC-P
Cell treatment	Pet1 Recombinant Rabbit mAb [PSH07-77]	HA722893	Human	WB,IF-Cell,IHC-P
Cell treatment	Phospho-Creb (S133) Recombinant Rabbit mAb [JB25-40] - BSA and Azide free	HA750465	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,IF-Tissue
Cell treatment	Phospho-Creb (S133) Recombinant Rabbit mAb [JB25-40]	ET7107-93	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,IF-Tissue
Cell treatment	Phospho-MAP2 (S136) Recombinant Rabbit mAb [JE44-25]	HA722351	Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Tissue
Cell treatment	Phospho-Tau (S198) Recombinant Rabbit mAb [JE66-79]	HA721668	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S199) Recombinant Rabbit mAb [JE66-80]	HA721667	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S199) Recombinant Rabbit mAb [JE46-95]	HA721800	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S202) Recombinant Rabbit mAb [JE66-81]	HA721666	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S214) Recombinant Rabbit mAb [JE66-83]	HA721797	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S324) Recombinant Rabbit mAb [JE30-73]	HA721798	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S356) Recombinant Rabbit mAb [JE46-86]	HA721799	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Tau (S396) Recombinant Rabbit mAb [SN62-09] - BSA and Azide free	HA750268	Mouse, Rat, Human, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-Tau (S396) Recombinant Rabbit mAb [SN62-09]	ET1611-68	Mouse, Rat, Human, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-Tau (S404) Recombinant Rabbit mAb [JE47-88]	HA721801	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-Tau (S404) Recombinant Rabbit mAb [JE66-82]	HA721665	Human, Mouse, Rat	WB,IHC
Cell treatment	Phospho-Tau (T181) Recombinant Rabbit mAb [PSH05-42]	HA722271	Human, Mouse, Cynomolgus monkey	WB,IHC-P
Cell treatment	Phospho-Tau (T181) Recombinant Rabbit mAb [PSH05-42] - BSA and Azide free	HA751001	Human, Mouse, Cynomolgus monkey	WB,IHC-P
Cell treatment	Phospho-Tau (T217) Recombinant Rabbit mAb [PSH09-33]	HA723091	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr
	Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-41] - BSA and Azide free (Capture)	HA723841	Human, Mouse	ELISA(Cap)
Cell treatment	Phospho-Tau (T217) Recombinant Rabbit mAb [PSH09-33] - BSA and Azide free	HA751287	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr
	Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-42] - BSA and Azide free (Capture)	HA723842	Human	ELISA(Cap)
	Phospho-Tau (T217) Recombinant Rabbit mAb [PSH17-41] - BSA and Azide free (Capture)	HA723922	Human, Mouse	ELISA(Cap)
	Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-75] - BSA and Azide free (Capture)	HA723877	Human	ELISA(Cap)
	Phospho-Tau (T231) Recombinant Rabbit mAb [PSH01-04] - BSA and Azide free	HA750680	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Phospho-Tau (T231) Recombinant Rabbit mAb [PSH01-04]	HA721570	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Phospho-Tau (T231) Recombinant Rabbit mAb [PSH01-05] - BSA and Azide free	HA750679	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-Tau (T231) Recombinant Rabbit mAb [SC58-08]	ET1610-31	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IP,IF-Tissue,IHC-Fr
	Phospho-Tau (T231) Recombinant Rabbit mAb [PSH01-05]	HA721569	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-Tau (T231) Recombinant Rabbit mAb [SC58-08] - BSA and Azide free	HA750210	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IP,IF-Tissue,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Vimentin (S39) Recombinant Rabbit mAb [JE43-26]	HA721442	Human, Mouse, Rat	WB,IF-Cell
	PLAGL1 Mouse mAb [C5-D8]	EM1708-70	Human	IHC-P
Relative expression	PSD95 Recombinant Antibody [SR38-09] - Rat IgG1 (Chimeric)	HA601416	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,WB
	PSD95 Mouse mAb [B5-E6]	M1511-4	Human, Mouse, Rat	WB,IF-Cell,FC
	PSD95 Recombinant Rabbit mAb [SR38-09]	ET1602-20	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,IHC-Fr
	PSD95 Recombinant Rabbit mAb [SR38-09] - BSA and Azide free	HA750035	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,IP,IHC-Fr
	PSD95 Rabbit pAb	R1508-4	Human, Mouse, Rat	WB,IHC-P,FC
Relative expression	PSD95 Recombinant Antibody [SR38-09] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610256	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,WB
	PU.1/Spi1 Recombinant Rabbit mAb [JB72-35]	ET7107-49	Human	WB,IHC-P,FC
	PU.1 / Spi1 Recombinant Rabbit mAb [PD00-14]	HA721162	Human	WB,IHC-P
	R-cadherin Rabbit pAb	ER64773	Human, Mouse	IHC-P,IF-Cell,WB
	Biotin Conjugated S100B Recombinant Rabbit mAb [PSH15-25] - Detector	HA723730B	Human, Mouse, Rat	ELISA(Det),ELISA
	S100B Recombinant Rabbit mAb [PSH15-25] - BSA and Azide free (Detector)	HA723729	Human, Mouse, Rat	ELISA(Det)
	S100B Recombinant Rabbit mAb [SC57-02] - BSA and Azide free (Capture)	HA723728	Human, Mouse, Rat	ELISA(Cap)
	S100 beta Rabbit pAb	RT1547	Human, Mouse, Rat	WB,IP,IF
	S100 beta Recombinant Antibody [SC57-02] - Mouse IgG1 (Chimeric)	HA601551	Mouse, Rat	IHC-Fr,IHC-P
	S100 beta Recombinant Antibody [SC57-02] - Guinea pig IgG2 (Chimeric)	HA601550	Mouse	IHC-Fr
Relative expression	S100 beta Recombinant Rabbit mAb [PSH17-13]	HA723900	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Tissue
	iFluor™ 594 Conjugated S100 beta Recombinant Rabbit mAb [SC57-02]	HA720125F	Human	IF-Tissue
	S100 beta Recombinant Rabbit mAb [SC57-02] - BSA and Azide free	HA750209	Human, Mouse, Rat, Zebrafish, Goat, Cynomolgus monkey, Green monkey, Pig, Rabbit	WB,IF-Tissue,IHC-P,IHC-Fr,IF-Cell,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	S100 beta Recombinant Rabbit mAb [PSH17-13] - BSA and Azide free	HA751617	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Tissue
	S100 beta Recombinant Antibody [SC57-02] - Rat IgG1 (Chimeric)	HA601552	Mouse, Rat	IHC-Fr,IHC-P
	S100 beta Recombinant Antibody [SC57-02] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610366	Mouse, Rat	IHC-Fr,IHC-P
	S100 beta Recombinant Antibody [SC57-02] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610367	Mouse, Rat	IHC-Fr,IHC-P
Relative expression	S100 beta Recombinant Mouse mAb [PD00-11]	HA601098	Human, Mouse, Rat	IHC-P,WB,IF-Tissue
	S100 beta Recombinant Antibody [SC57-02] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610365	Mouse	IHC-Fr
	S100 beta Recombinant Rabbit mAb [SC57-02]	ET1610-3	Human, Mouse, Rat, Zebrafish, Goat, Cynomolgus monkey, Green monkey, Pig, Rabbit	WB,IF-Tissue,IHC-P,IHC-Fr,IF-Cell,FC,IP
Relative expression	Serotonin transporter Recombinant Rabbit mAb [PSH19-18]	HA724052	Mouse, Rat	IHC-Fr,IHC-P,WB
	Serotonin transporter Rabbit pAb	ER1916-44	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Serotonin transporter Recombinant Rabbit mAb [JE34-53]	HA721544	Human, Mouse, Rat	WB,FC
Relative expression	SLUG Recombinant Rabbit mAb [PSH07-40]	HA722828	Human	WB,IF-Cell
Relative expression	SLUG Recombinant Rabbit mAb [PSH07-40] - BSA and Azide free	HA751161	Human	WB,IF-Cell
	SNAI2 Mouse mAb [D5-B6]	EM1706-65	Human	WB,IF-Cell,IHC-P,FC
	Sonic Hedgehog Protein/SHH(C-Product) Recombinant Rabbit mAb [SC05-28]	ET1610-6	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	SOX1 Recombinant Rabbit mAb [JJ20-40]	ET1701-38	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	SOX1 Recombinant Rabbit mAb [PSH10-51]	HA723215	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IHC-P
	SOX1 Recombinant Rabbit mAb [PSH10-51] - BSA and Azide free	HA751354	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IHC-P
	SOX1 Recombinant Rabbit mAb [JJ20-40] - BSA and Azide free	HA750308	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	SOX10 Mouse mAb [1-11]	M1311-1	Mouse	WB,IHC-P
Relative expression	SOX10 Recombinant Rabbit mAb [PDH0-05]	HA721239	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	SOX10 Recombinant Rabbit mAb [PDH0-03]	HA721240	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,IF-Tissue,FC,IHC-Fr,IP
Relative expression	SOX10 Recombinant Rabbit mAb [PDH0-04]	HA721241	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
KD	SOX10 Recombinant Rabbit mAb [SD204-04]	ET1612-79	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	SOX10 Recombinant Antibody [PDH0-03] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610253	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Tissue
	SOX10 Recombinant Antibody [PDH0-03] - Rat IgG1 (Chimeric)	HA601413	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Tissue
	SOX10 Rabbit pAb	HA500503	Human, Mouse	WB,IHC-P
	SOX10 Recombinant Rabbit mAb [SN20-41]	ET1611-89	Human, Mouse, Rat	WB
	SOX11 Recombinant Rabbit mAb [SD203-08]	ET1612-76	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IP,IHC-P
	Sox-12 Rabbit pAb	ER64999	Human, Mouse	IHC-P
Relative expression	SOX2 Recombinant Rabbit mAb [PO00-28]	HA721155	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,WB,IF-Cell,ChIP,IHC-Fr,IF-Tissue
	iFluor™ 488 Conjugated SOX2 Recombinant Rabbit mAb [PO00-28]	HA720179F	Human, Mouse	IF-Cell,FC
	SOX2 Recombinant Rabbit mAb	IRS240RB	Mouse	mIHC
	SOX2 Recombinant Antibody [PO00-28] - Chicken IgY (Chimeric) - BSA and Azide free	HA610329	Human, Mouse	IHC-Fr,WB
	SOX2 Mouse mAb [C9-F10]	EM1708-84	Human	WB,IF-Cell,IHC-P
	SOX2 Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS281MS	Human	mIHC
	SOX2 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS273RT	Human	mIHC
	iFluor™ 594 Conjugated SOX2 Recombinant Rabbit mAb [PO00-28]	HA720180F	Human, Mouse	IF-Cell,FC
	SOX2 Recombinant Antibody [PO00-28] - Chicken IgY (Chimeric)	HA601499	Human, Mouse	IHC-Fr,WB
	SOX2 Recombinant Antibody [PO00-28] - Guinea pig IgG2 (Chimeric)	HA601500	Human, Mouse	IHC-Fr,WB
	SOX2 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS274RT	Mouse	mIHC
	SOX2 Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS282MS	Mouse	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	SOX2 Rabbit pAb	IRS198RB	Human	mIHC
	SOX2 Rabbit pAb	R1106-1	Human, Mouse, Rat	WB
	SOX2 Recombinant Antibody [PO00-28] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610330	Human, Mouse	IHC-Fr, WB
	SOX2 Recombinant Antibody [PO00-28] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610236	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr, IHC-P, WB
	SOX2 Recombinant Antibody [PO00-28] - Rat IgG1 (Chimeric)	HA601396	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr, IHC-P, WB
	SOX4 Rabbit pAb	ER1916-97	Human, Mouse, Rat	WB, IHC-P, IF-Tissue, FC
Relative expression	SOX9 Recombinant Rabbit mAb [PSH13-59]	HA723548	Human, Mouse, Rat	IHC-Fr, IHC-P, WB, IF-Cell, IP, FC, ChIP
	SOX9 Recombinant Rabbit mAb	IRS241RB	Mouse	mIHC
	SOX9 Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS279MS	Human	mIHC
	SOX9 Recombinant Rabbit mAb	IRS199RB	Human	mIHC
Relative expression	SOX9 Recombinant Rabbit mAb [PSH13-59] - BSA and Azide free	HA751478	Human, Mouse, Rat	IHC-Fr, IHC-P, WB, IF-Cell, IP, FC, ChIP
Relative expression, KD	SOX9 Recombinant Rabbit mAb [SN74-09] - BSA and Azide free	HA750262	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB, IF-Cell, IHC-P, IP, FC, CUT&Tag-seq
Relative expression, KD	SOX9 Recombinant Rabbit mAb [SN74-09]	ET1611-56	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB, IF-Cell, IHC-P, IP, FC, CUT&Tag-seq, mIHC
	SOX9 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS280RT	Human	mIHC
	SST Mouse mAb [E7-G5]	EM1709-27	Human	WB, IHC-P, FC
	SST/Somatostatin 28 Recombinant Rabbit mAb [JE43-52] - BSA and Azide free	HA750648	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P, IHC-Fr, IF-Tissue
	SST/Somatostatin 28 Recombinant Rabbit mAb	IRS239RB	Mouse	mIHC
	SST/Somatostatin 28 Recombinant Rabbit mAb [JE43-52]	HA721460	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P, IHC-Fr, IF-Tissue
	SST/Somatostatin 28 Recombinant Rabbit mAb	IRS197RB	Human	mIHC
	Synaptophysin Recombinant Rabbit mAb [JE42-51]	HA721599	Human, Mouse, Rat	WB, IHC
	Synaptophysin Rabbit pAb	0407-2	Human, Mouse, Rat	WB, IF-Cell, IHC-P, FC
Relative expression	Synaptophysin Recombinant Rabbit mAb [SJ26-85]	ET1606-56	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB, IHC-P, IHC-Fr, IF-Tissue
	Synaptophysin Mouse mAb [A1-7]	M0407-30	Human, Mouse, Rat	WB, IF-Cell, IHC-P, FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Tau Recombinant Rabbit mAb [PSH16-49] - Detector	HA723852B	Human	ELISA(Det),ELISA
	Tau Recombinant Rabbit mAb [PSH17-27] - BSA and Azide free	HA723914	Human	ELISA
	Tau Recombinant Rabbit mAb [PSH10-42] - BSA and Azide free (Detector)	HA723209	Human, Mouse, Rat	ELISA(Det)
Relative expression	Tau Recombinant Rabbit mAb [PSH05-50] - BSA and Azide free	HA751009	Mouse, Rat	WB,IHC-P,IF-Tissue,IP
	Tau Recombinant Rabbit mAb [SD205-09]	ET1612-44	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
Relative expression	Tau Recombinant Rabbit mAb [PSH05-50]	HA722288	Mouse, Rat	WB,IHC-P,IF-Tissue,IP
	Biotin Conjugated Tau Recombinant Rabbit mAb [PSH17-11] - Detector	HA723898B	Human	ELISA(Det),ELISA
	Tau Recombinant Rabbit mAb [PSH17-26] - BSA and Azide free	HA723913	Human	ELISA
	Tau Recombinant Rabbit mAb [PSH16-49] - BSA and Azide free (Detector)	HA723851	Human	ELISA(Det)
	Tau Recombinant Rabbit mAb [PSH17-11] - BSA and Azide free (Detector)	HA723898	Human	ELISA(Det)
	Tau Rabbit pAb	R1401-7	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Biotin Conjugated Tau Recombinant Rabbit mAb [PSH10-42] - Detector	HA723209B	Human, Mouse, Rat	ELISA(Det),ELISA
KD	Tau Recombinant Rabbit mAb [SZ03-03]	ET1603-2	Human, Mouse, Rat	WB,IHC-P,IP,IHC-Fr
	Tau Recombinant Rabbit mAb [PSH16-48] - BSA and Azide free (Capture)	HA723850	Human, Mouse, Rat	ELISA(Cap)
Relative expression	TBR1 Recombinant Antibody [JF10-00] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610262	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
	TBR1 Recombinant Rabbit mAb [JF10-00] - BSA and Azide free	HA750366	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression	TBR1 Recombinant Antibody [JF10-00] - Rat IgG1 (Chimeric)	HA601422	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
	TBR1 Recombinant Rabbit mAb	IRS070	Mouse	mIHC
	TBR1 Recombinant Rabbit mAb [JF10-00]	ET1702-97	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,mIHC,IHC-Fr
Relative expression	TBR2 Recombinant Rabbit mAb [PSH07-88]	HA722918	Human, Mouse, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	TBR2 Recombinant Rabbit mAb [PSH07-88] - BSA and Azide free	HA751195	Human, Mouse, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
	Tenascin C Recombinant Rabbit mAb [SU36-01]	ET1608-50	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	TGIF Rabbit pAb	ER65155	Human, Mouse, Rat	WB
	Thy-1 (OX7) Mouse mAb [1G2]	RT1615	Mouse, Rat	WB,IP,IF,IHC-P,FC
	TMEM119 Recombinant Rabbit mAb [PSH12-98]	HA723503	Human, Mouse, Rat	WB,IF-Cell,FC,IP
	TMEM119 Recombinant Rabbit mAb [PSH12-98] - BSA and Azide free	HA751453	Human, Mouse, Rat	WB,IF-Cell,FC,IP
	TMEM119 Rabbit pAb	ER2001-37	Human	WB
	TMEM119 Rabbit pAb	HA500260	Human, Mouse	WB,IF-Cell,FC
	TPH2 Rabbit pAb	ER65200	Human, Mouse, Rat, Monkey	WB
KD	Transferrin Receptor (CD71) Recombinant Rabbit mAb [JF0956] - BSA and Azide free	HA750333	Human, Mouse, Rat	WB,IP,FC
Relative expression	Transferrin Receptor (CD71) Recombinant Rabbit mAb [PSH13-56]	HA723545	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP,IF-Tissue
KD	Transferrin Receptor (CD71) Recombinant Rabbit mAb [JF0956]	ET1702-06	Human, Mouse, Rat	WB,IP,FC
	Transferrin Receptor (CD71) Rabbit pAb	ER1706-35	Human	IHC-P
KD	Tryptophan Hydroxylase 1 (TPH1) Recombinant Rabbit mAb [SC53-07] - BSA and Azide free	HA750214	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	Tryptophan Hydroxylase 1 (TPH1) Recombinant Rabbit mAb [SC53-07]	ET1610-37	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Tyrosine Hydroxylase Rabbit pAb	R1706-18	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Tyrosine Hydroxylase Recombinant Rabbit mAb	IRS206RB	Mouse	mIHC
	Tyrosine Hydroxylase Recombinant Rabbit mAb [SN59-03] - BSA and Azide free	HA750243	Human, Mouse, Rat, Cynomolgus monkey	WB,IF-Cell,IHC-P,FC,IF-Tissue,IHC-Fr
	Tyrosine Hydroxylase Recombinant Rabbit mAb [SD080-02]	ET1612-65	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IF-Tissue,IHC-Fr
	Tyrosine Hydroxylase Recombinant Rabbit mAb [SN59-03]	ET1611-12	Human, Mouse, Rat, Cynomolgus monkey	WB,IF-Cell,IHC-P,FC,IF-Tissue,mIHC,IHC-Fr
	Tyrosine Hydroxylase Recombinant Rabbit mAb [SD080-02] - BSA and Azide free	HA750298	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IF-Tissue,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	VE Cadherin Recombinant Rabbit mAb [PSH03-15] - BSA and Azide free	HA750856	Human	WB,IF-Cell,FC
KD	VE Cadherin Recombinant Rabbit mAb [PSH03-15]	HA721943	Human	WB,IF-Cell,FC
	VGLUT1 Recombinant Antibody [PSH10-53] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610284	Mouse, Rat	IHC-Fr,IHC-P,WB
	VGLUT1 Recombinant Antibody [PSH10-53] - Chicken IgY (Chimeric)	HA601478	Mouse	IHC-Fr
	VGLUT1 Recombinant Antibody [PSH10-53] - Mouse IgG1 (Chimeric)	HA601450	Mouse, Rat	IHC-Fr,IHC-P,WB
	VGLUT1 Recombinant Antibody [PSH10-53] - Guinea pig IgG2 (Chimeric)	HA601449	Mouse, Rat	IHC-Fr,WB
	VGLUT1 Mouse mAb [A7A3]	HA601014	Human, Mouse, Rat	WB
	VGLUT1 Rabbit pAb	HA500488	Human, Mouse, Rat	WB
	VGLUT1 Rabbit pAb	HA500083	Human, Mouse, Rat	WB,IHC-P
	VGLUT1 Recombinant Antibody [PSH10-53] - Chicken IgY (Chimeric) - BSA and Azide free	HA610312	Mouse	IHC-Fr
	VGLUT1 Recombinant Mouse mAb [PSH08-34]	HA601371	Mouse, Rat	WB,IHC-P
Relative expression	VGLUT1 Recombinant Rabbit mAb [PSH10-53] - BSA and Azide free	HA751356	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IHC-P,IF-Tissue
	VGLUT1 Recombinant Antibody [PSH10-53] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610283	Mouse, Rat	IHC-Fr,WB
	VGLUT1 Recombinant Mouse mAb [PSH08-34] - BSA and Azide free	HA610217	Mouse, Rat	WB,IHC-P
Relative expression	VGLUT1 Recombinant Antibody [PSH10-53] - Rat IgG1 (Chimeric)	HA601392	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB,IF-Tissue
Relative expression	VGLUT1 Recombinant Antibody [PSH10-53] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610234	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB,IF-Tissue
Relative expression	VGLUT1 Recombinant Rabbit mAb [PSH10-53]	HA723217	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IHC-P,IF-Tissue
KO	Vimentin Rabbit pAb	R1308-6	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KO, Relative expression	Vimentin Recombinant Rabbit mAb [SC60-05] - BSA and Azide free	HA750215	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC,IP
	Vimentin Recombinant Antibody - Rat IgG1 (Chimeric)	HA601564	Human, Mouse, Rat	IHC,ICC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KO	Vimentin Mouse mAb [D4-B11]	EM0401	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,IHC-Fr,IP
	iFluor™ 647 Conjugated Vimentin Recombinant Rabbit mAb [SC60-05]	HA720165F	Human, Mouse, Rat	IF-Cell,IF-Tissue,FC
	Vimentin Recombinant Rabbit mAb [PDH0-01]	HA721174	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue,mIHC
Relative expression, KO	Vimentin Mouse mAb [A6-C1]	M1412-1	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Vimentin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS214RT	Human	mIHC
Relative expression	Vimentin Recombinant Mouse mAb [A6-C1-R]	HA601251	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Vimentin Recombinant Rabbit mAb	IRS044	Human	mIHC
	iFluor™ 488 Conjugated Vimentin Recombinant Rabbit mAb [SC60-05]	HA720166F	Human, Mouse, Rat	IF-Cell,IF-Tissue,FC
	Vimentin Recombinant Rabbit mAb	IRS060	Mouse	mIHC
	Vimentin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS263RT	Mouse	mIHC
KO, Relative expression	Vimentin Recombinant Rabbit mAb [SC60-05]	ET1610-39	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC,IP
Relative expression	Vimentin Recombinant Mouse mAb [A6-C1-R] - BSA and Azide free	HA610126	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	VIP Recombinant Antibody [PSH10-56] - Rat IgG1 (Chimeric)	HA601401	Human, Mouse, Rat, Cynomolgus monkey	IHC-Fr,IHC-P
	VIP Recombinant Antibody [PSH10-56] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610241	Human, Mouse, Rat, Cynomolgus monkey	IHC-Fr,IHC-P
	VIP Recombinant Rabbit mAb [PSH10-56] - BSA and Azide free	HA751359	Human, Mouse, Rat, Cynomolgus monkey	IHC-Fr,IHC-P,IF-Tissue
	VIP Recombinant Rabbit mAb [PSH10-56]	HA723220	Human, Mouse, Rat, Cynomolgus monkey	IHC-Fr,IHC-P,IF-Tissue
	VIP Rabbit pAb	ER65300	Human	IHC-P,WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层