

破骨细胞分化研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	AKT1 Recombinant Rabbit mAb [PSH0-50]	HA721364	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
KD	AKT1 Recombinant Rabbit mAb [ST05-09]	ET1609-47	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
KD	AKT1 Recombinant Rabbit mAb [ST05-09] - BSA and Azide free	HA750177	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
	AKT1 Mouse mAb [D9-9-C9]	EM40507	Human, Mouse, Rat	WB,IHC-P
	AKT1 Recombinant Rabbit mAb [PSH0-50] - BSA and Azide free	HA750627	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75-09] - BSA and Azide free	HA750840	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	AKT1/2/3 Recombinant Rabbit mAb [ST48-09]	ET1609-51	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
	AKT1/2/3 Recombinant Rabbit mAb [ST48-09] - BSA and Azide free	HA750179	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
	AKT1/2/3 Rabbit pAb	ER1901-33	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75-09]	HA721870	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	Akt1/3 Rabbit pAb	ER62638	Human, Mouse, Rat	WB,IHC-P
	Akt2 Mouse mAb [2A1]	RT1027	Human, Mouse, Rat	WB,IP,IF,IHC-P
	AKT2 Rabbit pAb	HA500091	Human, Mouse, Rat	WB,IHC-P
	AP-1/Jun D Rabbit pAb	ER62677	Human, Mouse, Rat	WB,IHC-P
	ATP2A1 / SERCA1 Recombinant Rabbit mAb [JE45-75]	ET7109-43	Human	WB,IP,IHC-P
Relative expression	ATP2A1 / SERCA1 Rabbit pAb	ER1803-34	Human, Mouse, Rat	WB,IHC-P,FC
	BLNK Recombinant Rabbit mAb [JE63-50]	HA721028	Human, Mouse, Rat	WB
	BLNK Mouse mAb [G1-F9]	EM1709-94	Human, Mouse	WB,IF-Cell,IHC-P,FC
KO, Relative expression	BTK Recombinant Rabbit mAb [SN69-07] - BSA and Azide free	HA750252	Human, Mouse	WB,IHC-P,IP
KO, Relative expression	BTK Recombinant Rabbit mAb [SN69-07]	ET1611-4	Human, Mouse	WB,IHC-P,IP
	Calcineurin A Recombinant Rabbit mAb [JA64-11]	ET1704-02	Human, Mouse, Rat	WB,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Calcineurin B Rabbit pAb	HA500346	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	CAMKIV Recombinant Rabbit mAb [JB33-13]	ET7107-96	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	CANB2 Rabbit pAb	ER60537	Human, Rat	WB
KD	Cathepsin H Recombinant Rabbit mAb [JE56-25]	HA721352	Human	WB,IF-Cell
	Cathepsin H Rabbit pAb	ER60014	Human	WB
	Cathepsin K Rabbit pAb	HA500469	Human	WB,IF-Cell
KD	Cathepsin L Recombinant Rabbit mAb [PSH03-91] - BSA and Azide free	HA750916	Human, Mouse, Rat, Monkey	WB,IHC-P,FC
	Cathepsin L Rabbit pAb	ER1905-49	Human, Rat	WB,IHC-P,IF-Tissue,FC
KD	Cathepsin L Recombinant Rabbit mAb [PSH03-92]	HA722064	Human, Mouse, Rat, Monkey	WB,IHC-P
KD	Cathepsin L Recombinant Rabbit mAb [PSH03-91]	HA722063	Human, Mouse, Rat, Monkey	WB,IHC-P,FC
KD	Cathepsin L Recombinant Rabbit mAb [PSH03-92] - BSA and Azide free	HA750917	Human, Mouse, Rat, Monkey	WB,IHC-P
KD	Cathepsin L/V/K/H Recombinant Rabbit mAb [JM10-78]	ET1703-44	Human, Mouse, Rat	WB,IF-Tissue,IHC-P
Relative expression	CD16 Recombinant Rabbit mAb [PD00-12]	HA721149	Human	IHC-P,IF-Tissue,IF-Cell,mIHC,WB
	CD16 Recombinant Rabbit mAb [JE49-79]	ET7109-97	Human, Mouse, Rat	WB
	CD16 Recombinant Rabbit mAb	IRS017	Human	mIHC
	CD16 Mouse mAb [Q32]	EM60067	Human	WB,IHC-P
	CD172b Rabbit pAb	ER62937	Human	WB
	CD172g Rabbit pAb	ER62938	Human, Monkey	WB
	CD32 Rabbit pAb	ER62945	Human, Mouse, Rat	WB,IHC-P
Cell treatment	CD32 Recombinant Rabbit mAb [JM10-70] - BSA and Azide free	HA750382	Human	WB,IHC-P,IF-Cell,FC,IP
Cell treatment	CD32 Recombinant Rabbit mAb [JM10-70]	ET1703-43	Human	WB,IHC-P,IF-Cell,FC,IP
	CD32A + CD32B + CD32C Recombinant Rabbit mAb [JE64-17]	HA721930	Human	WB,IHC-P
	CD32B Mouse mAb [A9C11]	HA601121	Human	WB,IHC-P
	CD32-C Rabbit pAb	ER62946	Human	WB
	CD64 Mouse mAb [A7E4]	HA601003	Human	WB
	CD64 Mouse mAb [5F5F8]	EM1710-99	Human	WB,IHC-P,FC
	CD85c Rabbit pAb	ER65604	Human	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD85e/h Rabbit pAb	ER65514	Human	WB,IHC-P
	CD85f Rabbit pAb	ER65588	Human	WB
	CD85g Rabbit pAb	ER65515	Human	WB,IHC-P
	CDC42 Recombinant Rabbit mAb [JJ086-04] - BSA and Azide free	HA750320	Human, Mouse, Rat	WB,FC,IF-Cell
	CDC42 Recombinant Rabbit mAb [JJ086-04]	ET1701-7	Human, Mouse, Rat	WB,FC,IF-Cell
Cell treatment	c-Fos Recombinant Antibody [PSH05-77] - Guinea pig IgG2 (Chimeric)	HA601374	Mouse, Rat, Human, Cynomolgus monkey, Pig	IHC-Fr,WB
Cell treatment	c-Fos Recombinant Antibody [PSH07-51] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610275	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Cell,WB
Cell treatment	c-Fos Recombinant Antibody [PSH05-77] - Rat IgG1 (Chimeric)	HA601373	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Tissue,IF-Cell,WB
Cell treatment	c-Fos Recombinant Antibody [PSH05-77] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610219	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Tissue,IF-Cell,WB
	c-Fos Recombinant Rabbit mAb	IRS131RB	Mouse	mIHC
Cell treatment	c-Fos Recombinant Rabbit mAb [JJ0938] - BSA and Azide free	HA750329	Human, Rat	WB,IF-Cell,IHC-P
Cell treatment	c-Fos Recombinant Antibody [PSH07-51] - Mouse IgG1 (Chimeric)	HA601440	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Cell,WB
	c-Fos Recombinant Rabbit mAb [PSH07-51] - BSA and Azide free	HA751172	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IF-Cell,IHC-P,ChIP,IF-Tissue,IP
Cell treatment	c-Fos Recombinant Rabbit mAb [PSH05-77] - BSA and Azide free	HA751071	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,IHC-Fr,IF-Tissue
Cell treatment	c-Fos Recombinant Rabbit mAb [JJ0938]	ET1701-95	Human, Rat	WB,IF-Cell,IHC-P
	c-Fos Recombinant Rabbit mAb [PSH07-51]	HA722839	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IF-Cell,IHC-P,ChIP,IF-Tissue,IP
Cell treatment	c-Fos Recombinant Antibody [PSH05-77] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610220	Mouse, Rat, Human, Cynomolgus monkey, Pig	IHC-Fr,WB
	c-Fos Mouse mAb [A10-G2]	EM1710-29	Human	WB,IHC-P,FC
Cell treatment	c-Fos Recombinant Rabbit mAb [PSH05-77]	HA722666	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,IHC-Fr,mIHC,IF-Tissue
	c-Jun Recombinant Rabbit mAb [SY0298]	ET1608-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	c-Jun Recombinant Rabbit mAb [SY0298] - BSA and Azide free	HA750143	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	Cleaved-Cathepsin L2 (L114) Rabbit pAb	ER60017	Human, Monkey	WB
	CREB Recombinant Rabbit mAb [PSH15-48] - BSA and Azide free	HA751563	Human, Mouse, Rat, Monkey	WB,IHC-P,IHC-Fr,IF-Cell,FC
	CREB Recombinant Rabbit mAb [PSH15-48]	HA723754	Human, Mouse, Rat, Monkey	WB,IHC-P,IHC-Fr,IF-Cell,FC
	CREB Recombinant Rabbit mAb [SA04-04]	ET1601-15	Human, Mouse, Zebrafish, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	CREB Recombinant Rabbit mAb [SA04-04] - BSA and Azide free	HA750007	Human, Mouse, Zebrafish, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	CREB1 Rabbit pAb	R1309-14	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	CSF-1-R Recombinant Rabbit mAb [PSH06-70]	HA722742	Mouse	WB,IF-Cell
Relative expression	CSF-1-R Recombinant Rabbit mAb [PSH09-77] - BSA and Azide free	HA751319	Human, Mouse, Rat	WB,IHC-P
Relative expression	CSF-1-R Recombinant Rabbit mAb [PSH06-70] - BSA and Azide free	HA751099	Mouse	WB,IF-Cell
	CSF-1-R Recombinant Rabbit mAb [PD00-45]	HA721180	Human	IHC-P,FC
Relative expression	CSF-1-R Recombinant Rabbit mAb [PSH09-77]	HA723142	Human, Mouse, Rat	WB,IHC-P
	CSF1R Mouse mAb [B9-B9]	EM1708-56	Human	WB,IF-Cell,IHC-P
	CY24A Rabbit pAb	ER61573	Human	WB
	DAP12 Recombinant Rabbit mAb [JG38-70]	ET7108-76	Human	WB,IHC-P,FC
KD	ERK1 Recombinant Rabbit mAb [SP05-09]	ET1604-16	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	ERK1 Mouse mAb [7-D6-E5]	M1505-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	ERK1 Recombinant Rabbit mAb [SP00-74]	ET1604-32	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Erk1/2 Mouse mAb [2G1]	RT1484	Human, Mouse, Rat	WB,IP,IF
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03]	ET1601-29	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03] - BSA and Azide free	HA750019	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
	ERK1/2 Rabbit pAb	ER131011	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	ERK2 Mouse mAb [B10-8]	EM1901-53	Human, Mouse, Rat	WB,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	ERK2 Recombinant Rabbit mAb [SZ25-01]	ET1603-23	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
	ERK2 Mouse mAb [C3B12]	EM1901-55	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Mouse mAb [A10D2]	EM1901-54	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Recombinant Mouse mAb [C3B12-R]	HA601201	Human, Mouse, Rat	WB
	ERK2 Mouse mAb [5-D2]	M1407-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	ERK2 Recombinant Mouse mAb [5-D2-R]	HA601202	Human, Mouse, Rat	WB
KD	ERK2 Recombinant Mouse mAb [C3B12-R] - BSA and Azide free	HA610079	Human, Mouse, Rat	WB
KD	ERK2 Recombinant Mouse mAb [5-D2-R] - BSA and Azide free	HA610080	Human, Mouse, Rat	WB
KD	ERK2 Rabbit pAb	ER131218	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	FCG3B Rabbit pAb	ER62049	Human	WB
	FCGR1A/CD64 Rabbit pAb	ER1908-78	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
	FosB Recombinant Rabbit mAb [SN703]	ET1611-75	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	FOSL2 Rabbit pAb	ER63449	Human, Mouse, Rat	WB,IHC-P
	Fra-1 Rabbit pAb	ER63468	Human, Mouse, Rat	WB,IHC-P
	FRA1 Recombinant Rabbit mAb [PSH06-95] - BSA and Azide free	HA751120	Human	WB
Relative expression	FRA1 Recombinant Rabbit mAb [JE33-20]	HA721447	Human	WB,IHC-P,IF-Cell,FC,ChIP
	FRA1 Recombinant Rabbit mAb [PSH06-95]	HA722771	Human	WB
KD	Fyn Recombinant Rabbit mAb [JE42-22]	HA721577	Human, Mouse, Rat	WB
	Fyn Rabbit pAb	ER1706-74	Human, Mouse, Rat	WB,IF-Cell,FC
	FYN Mouse mAb [D2-A10]	EM1709-88	Human	WB,IF-Cell,IHC-P,FC
	GRB2 Recombinant Rabbit mAb [JE49-73]	HA721923	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
	GRB2 Rabbit pAb	ER31203	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	GRB2 Mouse mAb [11-B11]	EM1701-41	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	PE Conjugated Human CD16 Recombinant Mouse mAb [PSH04-95]	HA600129F	Human	FC
	Human CD16 Recombinant Mouse mAb [PSH04-95]	HA601340	Human	FC
	Biotin Conjugated Human CD85d Recombinant Rabbit mAb [PSH05-62]	HA722592B	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CD85d Recombinant Rabbit mAb [PSH05-62] - BSA and Azide free (Detector)	HA722382	Human	ELISA(Det)
	Human CD85d Recombinant Rabbit mAb [PSH05-61] - BSA and Azide free (Capture)	HA722381	Human	ELISA(Cap)
	Human CSF-1-R Recombinant Rabbit mAb [PSH12-66] - BSA and Azide free (Capture/Detector)	HA723474	Human	ELISA(Cap),ELISA(Det)
	Human CSF-1-R Recombinant Rabbit mAb [PSH12-68] - BSA and Azide free (Capture)	HA723477	Human	ELISA(Cap)
	Biotin Conjugated Human CSF-1-R Recombinant Rabbit mAb [PSH12-67] - Detector	HA723476B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human CSF-1-R Recombinant Rabbit mAb [PSH12-66] - Detector	HA723478B	Human	ELISA(Det),ELISA
	Human CSF-1-R Recombinant Rabbit mAb [PSH12-67] - BSA and Azide free (Detector)	HA723475	Human	ELISA(Det)
	HRP Conjugated Human CSF-1-R Recombinant Rabbit mAb [PSH12-66] - Detector	HA723478H	Human	ELISA(Det),ELISA
	Human IFN beta Recombinant Rabbit mAb [PSH17-30] - BSA and Azide free (Detector)	HA725053	Human	ELISA(Det)
	Human IFN beta Recombinant Rabbit mAb [PSH17-29] - BSA and Azide free (Capture)	HA725052	Human	ELISA(Cap)
	Biotin Conjugated Human IFN beta Recombinant Rabbit mAb [PSH17-30] - Detector	HA725054B	Human	ELISA(Det),ELISA
	Human IFNGR1 Recombinant Rabbit mAb [PSH06-43] - BSA and Azide free (Capture)	HA722679	Human	ELISA(Cap)
	Biotin Conjugated Human IFNGR1 Recombinant Rabbit mAb [PSH06-44]	HA722681B	Human	ELISA(Det)
	Human IFNGR1 Recombinant Rabbit mAb [PSH06-44] - BSA and Azide free (Detector)	HA722680	Human	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-87] - BSA and Azide free (Capture)	HA723883	Human	ELISA(Cap)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free	HA750745	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - BSA and Azide free (Detector)	HA723882	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-1 beta Recombinant Rabbit mAb [PSH16-85] - BSA and Azide free (Capture)	HA723881	Human	ELISA(Cap)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39]	HA721712	Human	WB
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA722649B	Human	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free (Detector)	HA721279	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free	HA750744	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free (Capture)	HA721278	Human	ELISA(Cap)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA721713	Human	WB
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - Detector	HA723882B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human Interferon gamma Recombinant Rabbit mAb [PS01-76] - Detector	HA721576	Human	ELISA(Det),ELISA
	Human Interferon gamma Recombinant Rabbit mAb [PS01-75] - BSA and Azide free (Capture)	HA721574	Human	ELISA(Cap)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-45] - BSA and Azide free (Capture)	HA721683	Human	ELISA(Cap)
	Human Interferon gamma Recombinant Rabbit mAb [PS01-76] - BSA and Azide free (Detector)	HA721575	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-47] - BSA and Azide free (Detector)	HA721685	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-46] - BSA and Azide free (Capture)	HA721684	Human	ELISA(Cap)
	Biotin Conjugated Human Interferon gamma Recombinant Rabbit mAb [PSH01-47]	HA722645B	Human	ELISA(Det)
	Biotin Conjugated Human M-CSF Recombinant Rabbit mAb [PSH11-38] - Detector	HA723326B	Human	ELISA(Det),ELISA
	Human M-CSF Recombinant Rabbit mAb [PSH11-38] - BSA and Azide free (Detector)	HA723325	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human M-CSF Recombinant Rabbit mAb [PSH11-37] - BSA and Azide free (Capture)	HA723324	Human	ELISA(Cap)
	Human M-CSF Recombinant Rabbit mAb [PSH11-36] - BSA and Azide free (Capture)	HA723323	Human	ELISA(Cap)
	Human Osteoprotegerin Recombinant Rabbit mAb [PSH13-31] - BSA and Azide free (Capture)	HA725067	Human	ELISA(Cap)
	Human Osteoprotegerin Recombinant Rabbit mAb [PSH13-32] - BSA and Azide free (Detector)	HA725068	Human	ELISA(Det)
	Biotin Conjugated Human Osteoprotegerin Recombinant Rabbit mAb [PSH13-32] - Detector	HA725069B	Human	ELISA(Det),ELISA
	Human RANKL Recombinant Rabbit mAb [PSH16-59] - BSA and Azide free (Detector)	HA723861	Human	ELISA(Det)
	Human RANKL Recombinant Rabbit mAb [PSH16-58] - BSA and Azide free (Capture)	HA723860	Human	ELISA(Cap)
	Human RANKL Recombinant Rabbit mAb [PSH16-60] - BSA and Azide free (Capture)	HA723862	Human	ELISA(Cap)
	Biotin Conjugated Human RANKL Recombinant Rabbit mAb [PSH16-59] - Detector	HA723861B	Human	ELISA(Det),ELISA
	Human SIRP alpha Recombinant Rabbit mAb [PSH17-63] - BSA and Azide free (Capture)	HA725347	Human	ELISA(Cap)
	Human SIRP alpha Recombinant Rabbit mAb [PSH17-64] - BSA and Azide free (Detector)	HA725348	Human	ELISA(Det)
	Biotin Conjugated Human SIRP alpha Recombinant Rabbit mAb [PSH17-64] - Detector	HA725348B	Human	ELISA(Det),ELISA
	Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - BSA and Azide free (Detector)	HA723159	Human	ELISA(Det)
	Biotin Conjugated Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - BSA and Azide free (Detector)	HA723160B	Human	ELISA(Det),ELISA
	Human sTNF RI Recombinant Rabbit mAb [PSH09-89] - BSA and Azide free (Capture)	HA723158	Human	ELISA(Cap)
	HRP Conjugated Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - Detector	HA723374H	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HRP Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721294	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb	IRS040	Human	mIHC
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free (Capture)	HA721265	Human	ELISA(Cap)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free	HA750593	Human	IHC-P
	Human TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free (Detector)	HA721266	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40]	HA721286	Human	IHC-P,mIHC
	Biotin Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721559	Human	ELISA(Det)
	Human TREM2 Recombinant Rabbit mAb [PSH09-85] - BSA and Azide free (Detector)	HA723151	Human	ELISA(Det)
	Biotin Conjugated Human TREM2 Recombinant Rabbit mAb [PSH09-83] - Detector	HA723153B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human TREM2 Recombinant Rabbit mAb [PSH09-85] - Detector	HA723152B	Human	ELISA(Det),ELISA
	Human TREM2 Recombinant Rabbit mAb [PSH09-84] - BSA and Azide free (Capture)	HA723150	Human	ELISA(Cap)
	HRP Conjugated Human TREM2 Recombinant Rabbit mAb [PSH09-83] - Detector	HA723372	Human	ELISA(Det)
	Human TREM2 Recombinant Rabbit mAb [PSH09-83] - BSA and Azide free (Capture/Detector)	HA723149	Human	ELISA(Cap),ELISA(Det)
	IFN alpha R1 Mouse mAb [2G1]	RT1319	Mouse	WB,IP,FC
	IFNAR1 Recombinant Rabbit mAb [SR45-08]	ET1602-37	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	IFNAR1 Rabbit pAb	R1511-30	Human	IF-Cell,IHC-P,FC,IF-Tissue
	IFNAR1 Recombinant Rabbit mAb [SR45-08] - BSA and Azide free	HA750046	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	IFNAR1 Mouse mAb [A5-A3]	EM1701-46	Human, Mouse	WB,IHC-P,IF-Cell
	IFNAR2 Rabbit pAb	HA500425	Human, Rat	WB
	IFNAR2 Rabbit pAb	HA500284	Human	IHC-P,IF-Cell,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	IFN beta Rabbit pAb	ER65794	Human	IHC-P
Relative expression, KD	IFNGR1 Recombinant Rabbit mAb [JE40-50]	HA721952	Human, Mouse, Rat	WB,IHC-P
Relative expression, KD	IFNGR1 Recombinant Rabbit mAb [PSH06-11]	HA722508	Human	WB,IF-Cell
Relative expression, KD	IFNGR1 Recombinant Rabbit mAb [PSH06-11] - BSA and Azide free	HA751049	Human	WB,IF-Cell
KD, Relative expression	IFNGR1 Rabbit pAb	ER63732	Human	WB,IF-Cell,FC
	IFNGR2 Rabbit pAb	HA500363	Human, Mouse, Rat	WB,IF-Cell
	IFNGR2 Rabbit pAb	HA500246	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	IKB alpha Recombinant Rabbit mAb [SZ00-07] - BSA and Azide free	HA750068	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	IKB alpha Recombinant Rabbit mAb [SZ00-07]	ET1603-6	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	IKK alpha Recombinant Rabbit mAb [SN06-99]	ET1611-15	Human, Mouse, Rat	WB,IHC-P,IP,FC
	IKK alpha Rabbit pAb	ER30911	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	IKK alpha Recombinant Rabbit mAb [SN06-99] - BSA and Azide free	HA750246	Human, Mouse, Rat	WB,IHC-P,IP,FC
KD	IKK alpha + IKK beta Recombinant Rabbit mAb [SN63-02] - BSA and Azide free	HA750248	Human, Mouse, Rat	WB,IF-Cell,IP
KD	IKK alpha + IKK beta Recombinant Rabbit mAb [SN63-02]	ET1611-23	Human, Mouse, Rat	WB,IF-Cell,IP
	IKK beta Rabbit pAb	R1706-13	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	IKK beta Recombinant Rabbit mAb [JE32-81]	HA721715	Human	WB,FC
KD	IKK gamma Recombinant Rabbit mAb [JA11-83]	ET1704-40	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,IP,FC
Cell treatment	IL-1 alpha Recombinant Rabbit mAb [PSH15-65] - BSA and Azide free	HA751576	Mouse	WB
	IL-1 alpha Recombinant Rabbit mAb [PSH15-66]	HA723773	Mouse	IF-Cell,FC
	IL-1 alpha Recombinant Rabbit mAb [PSH15-66] - BSA and Azide free	HA751577	Mouse	IF-Cell,FC
	IL-1 alpha Rabbit pAb	ER1802-36	Human, Mouse	IHC-P,FC,WB
Cell treatment	IL-1 alpha Recombinant Rabbit mAb [PSH15-65]	HA723772	Mouse	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R]	HA601036	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Mouse mAb [A7A10]	HA601035	Human	WB,ELISA,IF-Cell
Cell treatment	IL-1 beta Recombinant Rabbit mAb [JJ087-3]	ET1701-39	Human	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R] - BSA and Azide free	HA610004	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
	IL-1 beta Recombinant Rabbit mAb [PSH18-00]	HA723965	Mouse, Rat	WB
Cell treatment	IL-1 beta Mouse mAb [A7F9]	HA601002	Human, Mouse, Rat	WB,ELISA,IF-Cell
	IL1R1 Mouse mAb [D12-D2]	EM1708-33	Human	WB,FC
	ILT-3 Rabbit pAb	ER63760	Human	IF-Cell
KD	Integrin alpha V Recombinant Rabbit mAb [SC56-07] - BSA and Azide free	HA750201	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
KD	Integrin alpha V Recombinant Rabbit mAb [SC56-07]	ET1610-15	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
	Integrin Alpha V Rabbit pAb	ER1911-52	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,IP
Relative expression	Integrin beta 3 Recombinant Rabbit mAb [PSH09-27]	HA723085	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Integrin beta 3 Rabbit pAb	ER1601-4	Human, Mouse, Rat	WB,IHC-P
	Integrin beta 3 Recombinant Rabbit mAb [JE22-64]	ET7110-61	Human, Mouse	WB,IHC-P
KD	Integrin beta 3 Recombinant Rabbit mAb [SJ19-09]	ET1606-49	Human, Mouse	WB,IHC-P,IF-Tissue
	Interferon gamma Recombinant Rabbit mAb [JM10-10]	ET1703-17	Human	WB,FC
	Interferon gamma Recombinant Rabbit mAb [JM10-10] - BSA and Azide free	HA750373	Human	WB,FC
Relative expression	IP3 receptor Recombinant Mouse mAb [PSH14-18] - BSA and Azide free	HA610290	Mouse, Rat	WB,IHC-P
Relative expression	IP3 receptor Recombinant Mouse mAb [PSH14-18]	HA601456	Mouse, Rat	WB,IHC-P
	IP3 Receptor Recombinant Rabbit mAb [JA11-35]	ET1704-77	Human, Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Tissue
	IP3 Receptor Rabbit pAb	ER1803-17	Human, Mouse, Rat	WB,IHC-P
Cell treatment	IRF9 Recombinant Rabbit mAb [PSH08-20]	HA722963	Human	WB,IF-Cell,IP,ChIP
Cell treatment	IRF9 Recombinant Rabbit mAb [PSH08-20] - BSA and Azide free	HA751211	Human	WB,IF-Cell,IP,ChIP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	IRF9 Recombinant Rabbit mAb [PSH07-87] - BSA and Azide free	HA751193	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
Cell treatment	IRF9 Recombinant Rabbit mAb [PSH07-87]	HA722916	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	ITPR2 Rabbit pAb	ER1803-65	Human, Mouse, Rat	Dot Blot,IF-Cell,IHC-P,FC
	JAK1 Recombinant Rabbit mAb [JM75-03] - BSA and Azide free	HA750442	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell,IF-Tissue
	JAK1 Mouse mAb [1G1]	RT1341	Human, Mouse, Rat	WB,IP,IF,IHC-P
	JAK1 Recombinant Rabbit mAb [JM75-03]	ET1705-84	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell,IF-Tissue
	JNK1 Rabbit pAb	R1309-1	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	JNK1 Recombinant Rabbit mAb [JE63-04]	HA721899	Human, Mouse, Rat	WB,IHC-P
	JNK1+JNK2+JNK3 Recombinant Rabbit mAb [SA43-06]	ET1601-28	Human, Mouse, Rat, Cow, Monkey, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP
	JNK1+JNK2+JNK3 Recombinant Rabbit mAb [SA43-06] - BSA and Azide free	HA750018	Human, Mouse, Rat, Cow, Monkey, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP
	JNK1+JNK3 Recombinant Rabbit mAb [JJ08-86]	ET1701-47	Human, Mouse, Rat	WB,IP,IHC-P
	JNK2 Rabbit pAb	HA500110	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	JNK2 Recombinant Rabbit mAb [SC05-85]	ET1610-11	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	JNK3 Recombinant Rabbit mAb [SD082-09]	ET1612-68	Human, Mouse, Rat	WB,IF-Tissue,IHC-P
	JunB Recombinant Rabbit mAb [SD081-08]	ET1612-57	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	JunD Recombinant Rabbit mAb [SD0830]	ET1612-92	Human	WB,IF-Cell,IF-Tissue,IP,IHC-P,FC
	JunD Recombinant Rabbit mAb [JE59-51]	HA722913	Human, Mouse, Rat	WB,IF-Cell,IHC-P,ChIP,IP
	JunD Recombinant Rabbit mAb [PSH08-85]	HA723037	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP,ChIP
	JunD Recombinant Rabbit mAb [PSH08-85] - BSA and Azide free	HA751250	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP,ChIP
	Lck Rabbit pAb	ER63849	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Relative expression	Lck Recombinant Rabbit mAb [PSH09-25] - BSA and Azide free	HA751279	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	Lck Rabbit pAb	HA500280	Human, Rat	WB,FC
Relative expression	Lck Recombinant Rabbit mAb [PSH09-25]	HA723083	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	LILRB1 Recombinant Rabbit mAb [PSH04-59] - BSA and Azide free	HA750966	Human	WB
Relative expression	LILRB1 Recombinant Rabbit mAb [PSH04-59]	HA722158	Human	WB
	LILRB2 Mouse mAb [2G7E11]	EM1711-09	Human	WB,FC
	LILRB2 Rabbit pAb	ER1912-34	Mouse	WB
	LIR-1 Rabbit pAb	ER65632	Human	WB
	LIR-6 Rabbit pAb	ER65589	Human	WB,IHC-P
	MAP3K14 Mouse mAb [2G2B4]	EM1712-23	Human	WB,IHC-P,FC
	MAP3K7 Mouse mAb [1G5]	RT1377	Human, Mouse, Rat	WB,IP,IF
	M-CSF Recombinant Rabbit mAb [SU0413]	ET1609-1	Human, Mouse	WB,IHC-P,IF-Cell,IP
KD	MEK1 Recombinant Rabbit mAb [SZ22-01]	ET1603-20	Human, Mouse, Rat, Cow, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK1/2 Recombinant Rabbit mAb [SR13-07]	ET1602-3	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK1/2 Recombinant Rabbit mAb [SR13-07] - BSA and Azide free	HA750042	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK2 Recombinant Rabbit mAb [SD207-08]	ET1612-6	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	MEK3 Recombinant Rabbit mAb [SD20-93]	ET1612-2	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK3+MEK6 Recombinant Rabbit mAb [JJ087-09]	ET1701-48	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	MEK-6 Rabbit pAb	ER63945	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	MEK6 Recombinant Rabbit mAb [PSH07-86]	HA722917	Human, Mouse, Rat	WB,IP
	MEK6 Recombinant Rabbit mAb [PSH07-86] - BSA and Azide free	HA751194	Human, Mouse, Rat	WB,IP
KD	MEK7 Recombinant Rabbit mAb [SD20-87]	ET1612-35	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	MiTF Recombinant Rabbit mAb [JF100-01]	ET1702-86	Human, Mouse, Rat	WB,IF-Cell,FC
KD	MiTF Recombinant Rabbit mAb [JF100-01] - BSA and Azide free	HA750362	Human, Mouse, Rat	WB,IF-Cell,FC
	MK13 Rabbit pAb	ER60146	Human, Mouse, Rat	WB
	PE Conjugated Mouse CD16/32 Recombinant Rabbit mAb [PSH06-60]	HA720216F	Mouse	FC
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-41] - BSA and Azide free (Detector)	HA723527	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-40] - BSA and Azide free (Capture)	HA723526	Mouse	ELISA(Cap)
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-42] - BSA and Azide free (Capture)	HA723529	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-43] - Detector	HA723531B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-41] - Detector	HA723528B	Mouse	ELISA(Det),ELISA
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-43] - BSA and Azide free (Detector)	HA723530	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67]	HA722608B	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - Detector	HA722625B	Mouse	ELISA(Det),ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH09-58] - BSA and Azide free	HA723120	Mouse	ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-66] - BSA and Azide free (Capture)	HA722166	Mouse	ELISA(Cap)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67] - BSA and Azide free (Detector)	HA722167	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH09-59] - BSA and Azide free	HA723121	Mouse	ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - BSA and Azide free (Detector)	HA721977	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-31] - BSA and Azide free (Capture)	HA721976	Mouse	ELISA(Cap)
	Mouse M-CSF/CSF1 Recombinant Rabbit mAb [PSH15-85] - BSA and Azide free (Detector)	HA723793	Mouse	ELISA(Det)
	Biotin Conjugated Mouse M-CSF/CSF1 Recombinant Rabbit mAb [PSH15-85] - Detector	HA723794B	Mouse	ELISA(Det),ELISA
	Mouse M-CSF/CSF1 Recombinant Rabbit mAb [PSH15-84] - BSA and Azide free (Capture)	HA723792	Mouse	ELISA(Cap)
Cell treatment	Mouse RANKL Recombinant Rabbit mAb [PSH07-91] - BSA and Azide free (Capture)	HA722923	Mouse	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Biotin Conjugated Mouse RANKL Recombinant Rabbit mAb [PSH07-92]	HA722925B	Mouse	ELISA(Det)
Cell treatment	Mouse RANKL Recombinant Rabbit mAb [PSH07-90] - BSA and Azide free (Capture)	HA722922	Mouse	ELISA(Cap)
Cell treatment	Mouse RANKL Recombinant Rabbit mAb [PSH07-92] - BSA and Azide free (Detector)	HA722924	Mouse	ELISA(Det)
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30] - BSA and Azide free (Detector)	HA721975	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30]	HA722626B	Mouse, Rat	ELISA(Det)
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-29] - BSA and Azide free (Capture)	HA721974	Mouse, Rat	ELISA(Cap)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77] - BSA and Azide free (Detector)	HA722041	Mouse	ELISA(Det)
	Biotin Conjugated Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77]	HA722623B	Mouse	ELISA(Det)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-76] - BSA and Azide free (Capture)	HA722040	Mouse	ELISA(Cap)
Relative expression	NFAT1 Recombinant Rabbit mAb [JA11-08]	ET1704-14	Human	WB,IF-Tissue,IHC-P
Relative expression	NFAT1 Recombinant Rabbit mAb [JA11-08] - BSA and Azide free	HA750405	Human	WB,IF-Tissue,IHC-P
	NFAT2 Recombinant Rabbit mAb [JA81-03]	ET1704-45	Human	WB,IHC-P,FC
	NFkB Inducing Kinase NIK Rabbit pAb	ER1913-80	Human, Mouse	WB,IHC-P,IF-Tissue,FC
KD	NF-kB p100 / NFkB2 Mouse mAb [16B2]	EM1901-78	Human	WB,IHC-P,FC
KD	NF-kB p100 / NFkB2 Recombinant Rabbit mAb [JM82-03]	ET1705-45	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
KD	NF-kB p100 / NFkB2 Mouse mAb [16B3]	EM1901-79	Human, Mouse, Rat	WB,IHC-P,FC
	NF-kB p105 / p50 Mouse mAb [A2-7]	M1505-8	Human, Mouse	WB,IF-Cell,IHC-P
KD	NF-kB p105 / p50 Recombinant Rabbit mAb [SZ20-01] - BSA and Azide free	HA750060	Human, Mouse, Rat	WB,IHC-P,ICC,FC,IF-Tissue
KD	NF-kB p105 / p50 Recombinant Rabbit mAb [SZ20-01]	ET1603-18	Human, Mouse, Rat	WB,IHC-P,ICC,FC,IF-Tissue
	NF-kB p105 / p50 Rabbit pAb	R1309-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	NF-kB p65 Rabbit pAb	ER0815	Human, Mouse, Rat, Zebrafish	WB,IHC-P,FC
	NF-kB p65 Rabbit pAb	HA500402	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27] - BSA and Azide free	HA750601	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04] - BSA and Azide free	HA750058	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04]	ET1603-12	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27]	HA721307	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	Nox1 Rabbit pAb	ER1913-99	Human, Mouse	WB
	NOXA2 Recombinant Rabbit mAb [JG86-31]	ET7108-72	Human, Mouse, Rat	WB,IHC-P,IP
	OSCAR Rabbit pAb	ER61888	Human, Mouse	WB
	Osteoprotegerin Mouse mAb [9C1]	EM1701-98	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Osteoprotegerin Mouse mAb [9C2]	EM1701-99	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Osteoprotegerin Rabbit pAb	R1608-4	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	p38 alpha / MAPK14 Rabbit pAb	R1308-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD, Cell treatment	p38 alpha / MAPK14 Recombinant Rabbit mAb [SR43-04] - BSA and Azide free	HA750040	Human, Mouse, Rat	WB,IF-Cell,FC
KD, Cell treatment	p38 alpha / MAPK14 Recombinant Rabbit mAb [SR43-04]	ET1602-26	Human, Mouse, Rat	WB,IF-Cell,FC
	p38 alpha / MAPK14 Mouse mAb [9-B9]	EM40307	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KO	p38 alpha / MAPK14 Recombinant Rabbit mAb [JF55-07]	ET1702-65	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	p38 alpha/MAPK14 + p38 beta/MAPK11 + p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-60] - BSA and Azide free	HA750894	Human, Mouse, Rat, Monkey	WB
	p38 alpha/MAPK14 + p38 beta/MAPK11 + p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-60]	HA722019	Human, Mouse, Rat, Monkey	WB
	p38 beta/MAPK11 + p38 gamma/MAPK12 + p38 delta/MAPK13 Recombinant Rabbit mAb [PSH03-59] - BSA and Azide free	HA750893	Human, Mouse, Rat, Monkey	WB
	p38 beta/MAPK11 + p38 gamma/MAPK12 + p38 delta/MAPK13 Recombinant Rabbit mAb [PSH03-59]	HA722018	Human, Mouse, Rat, Monkey	WB
	p38 delta / MAPK13 Rabbit pAb	HA500582	Human	WB
	p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-58]	HA722017	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-58] - BSA and Azide free	HA750892	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	p38 MAPK Recombinant Rabbit Multiclonal Antibody [PSH03-72] - BSA and Azide free	HA750906	Human, Mouse, Rat, Monkey	WB
	p38 MAPK Recombinant Rabbit Multiclonal Antibody [PSH03-72]	HA722036	Human, Mouse, Rat, Monkey	WB
	p40-phox Recombinant Rabbit mAb [JE58-57]	HA720040	Human, Mouse, Rat	WB,IHC-P
	p47-phox Recombinant Mouse mAb [A6B10-R] - BSA and Azide free	HA610139	Human	WB
	p47-phox Recombinant Mouse mAb [A6B10-R]	HA601264	Human	WB
	p47-phox Mouse mAb [A6B10]	HA600050	Human, Mouse	WB,IF-Cell,IHC-P,FC
	p47-phox Rabbit pAb	ER64501	Human, Mouse, Rat	WB,IHC-P
	p47-phox Mouse mAb [A6B11]	HA600049	Human, Mouse	WB,IF-Cell,IHC-P,FC
Cell treatment	Phospho-AKT1 (S124) Recombinant Rabbit mAb [JJ08-46]	ET1701-36	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-AKT1 (S129) Recombinant Rabbit mAb [JE59-50]	HA721677	Human, Mouse, Rat	WB
Cell treatment	Phospho-AKT1 (T450) Recombinant Rabbit mAb [SD08-12]	ET1612-73	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell,FC
Cell treatment	Phospho-AKT2 (S474) Recombinant Rabbit mAb [PSH09-46] - BSA and Azide free	HA751293	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-AKT2 (S474) Rabbit pAb	HA500116	Human, Mouse	WB,IHC-P
Cell treatment	Phospho-AKT2 (S474) Recombinant Rabbit mAb [PSH09-46]	HA723108	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44]	HA722129	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44] - BSA and Azide free	HA750953	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05] - BSA and Azide free	HA750132	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05]	ET1607-73	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10]	HA722951	Human	WB,IF-Cell,FC
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10] - BSA and Azide free	HA751205	Human	WB,IF-Cell,FC
Cell treatment	Phospho-BTK (Y223) Recombinant Rabbit mAb [ST51-05]	ET1609-55	Human, Mouse	WB,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-c-Jun (S63) Recombinant Rabbit mAb [SY0297]	ET1608-4	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-c-Jun (S63) Recombinant Rabbit mAb [SY0297] - BSA and Azide free	HA750146	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit mAb [PSH07-12]	HA722800	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit mAb [PSH07-12] - BSA and Azide free	HA751133	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit mAb [JE65-46]	HA722475	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-c-Jun (T91) Recombinant Rabbit mAb [JJ080-9]	ET1701-32	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-Creb (S133) Recombinant Rabbit mAb [JB25-40] - BSA and Azide free	HA750465	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,IF-Tissue
Cell treatment	Phospho-Creb (S133) Recombinant Rabbit mAb [JB25-40]	ET7107-93	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4]	ET1603-22	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4] - BSA and Azide free	HA750061	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01]	ET1610-13	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01] - BSA and Azide free	HA750199	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-IKB alpha (S32) Recombinant Rabbit mAb [ST53-05]	ET1609-78	Human, Mouse	WB,IF-Cell,IP
Cell treatment	Phospho-IKB alpha (S32) Recombinant Rabbit mAb [ST53-05] - BSA and Azide free	HA750196	Human, Mouse	WB,IF-Cell,IP
Cell treatment	Phospho-IKB alpha (S32/S36) Recombinant Rabbit mAb [PSH06-76] - BSA and Azide free	HA751119	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Phospho-IKB alpha (S32/S36) Recombinant Rabbit mAb [PSH06-76]	HA722770	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Phospho-IKB alpha (S36) Recombinant Rabbit mAb [JE56-65]	HA721802	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-JNK1/2/3 (T183 + T183 + T221) Recombinant Rabbit mAb [ST500]	ET1609-42	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
Cell treatment	Phospho-JNK1/2/3 (T183 + T183 + T221) Recombinant Rabbit mAb [ST500] - BSA and Azide free	HA750173	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-JunD (S255) Recombinant Rabbit mAb [JJ08-21]	ET1701-35	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment, Relative expression	Phospho-LYN (Y397)/LCK (Y394)/HCK (Y411)/BLK (Y389) Recombinant Rabbit mAb [PSH09-43] - BSA and Azide free	HA751290	Human, Mouse, Rat	WB,IF-Cell
Cell treatment, Relative expression	Phospho-LYN (Y397)/LCK (Y394)/HCK (Y411)/BLK (Y389) Recombinant Rabbit mAb [PSH09-43]	HA723105	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-MEK1/2 (S218 + S222) Recombinant Rabbit mAb [ST0490]	ET1609-50	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Phospho-MEK1 (S298) Recombinant Rabbit mAb [SD206-7]	ET1612-40	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-MEK1 (T292) Recombinant Rabbit mAb [SD2088]	ET1612-42	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-MEK3 (S189)+MEK6 (S207) Recombinant Rabbit mAb [PSH11-17]	HA723299	Human, Mouse	WB
Cell treatment	Phospho-MEK3 (S189)+MEK6 (S207) Recombinant Rabbit mAb [PSH11-17] - BSA and Azide free	HA751398	Human, Mouse	WB
Cell treatment	Phospho-NF-kB p65 (S529) Recombinant Rabbit mAb [SP07-00]	ET1604-27	Human	WB,IHC-P,FC,IP
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58]	HA723223	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58] - BSA and Azide free	HA751361	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-p38 alpha (T180) Recombinant Rabbit mAb [PSH06-25] - BSA and Azide free	HA751062	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180) Recombinant Rabbit mAb [PSH06-25]	HA722657	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180 + Y182) Recombinant Rabbit mAb [PSH06-35]	HA722669	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180 + Y182) Recombinant Rabbit mAb [PSH06-35] - BSA and Azide free	HA751074	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-p38 (T180) Recombinant Rabbit mAb [JE77-37]	HA722150	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-p38 (T180) Rabbit pAb	ER2001-51	Human, Mouse, Rat	WB
Cell treatment	Phospho-p38 (T180 + Y182) Rabbit pAb	ER2001-52	Human, Mouse, Rat	WB
	Phospho-p38 (T180 + Y182) Rabbit pAb	ER1903-01	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-39] - BSA and Azide free	HA750717	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-38]	HA721672	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-39]	HA721673	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-38] - BSA and Azide free	HA750716	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-Rac1+Cdc42 (S71) Recombinant Rabbit mAb [JJ086-9]	ET1701-40	Human	WB,IF-Cell,IHC-P
Cell treatment, Relative expression	Phospho-STAT1 (S727) Recombinant Rabbit mAb [SN67-04] - BSA and Azide free	HA750247	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment, Relative expression	Phospho-STAT1 (S727) Recombinant Rabbit mAb [SN67-04]	ET1611-20	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-STAT1 (Y701) Recombinant Rabbit mAb [PSH04-02] - BSA and Azide free	HA750925	Human	WB,IF-Cell
Cell treatment	Phospho-STAT1 (Y701) Recombinant Rabbit mAb [PSH04-02]	HA722083	Human	WB,IF-Cell
Cell treatment	Phospho-Syk (Y352) Recombinant Rabbit mAb [PSH11-66]	HA723362	Human	WB,IF-Cell,FC
Cell treatment	Phospho-Syk (Y352) Recombinant Rabbit mAb [PSH11-66] - BSA and Azide free	HA751422	Human	WB,IF-Cell,FC
Cell treatment	Phospho-Syk (Y525 + Y526) Recombinant Rabbit mAb [PSH03-70] - BSA and Azide free	HA750904	Human	WB,IF-Cell
Cell treatment	Phospho-Syk (Y525 + Y526) Recombinant Rabbit mAb [PSH03-70]	HA722034	Human	WB,IF-Cell
Cell treatment	Phospho-TAK1 (S412) Recombinant Rabbit mAb [PSH14-64]	HA723662	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-TAK1 (S412) Recombinant Rabbit mAb [PSH14-64] - BSA and Azide free	HA751525	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-TAK1 (S439) Recombinant Rabbit mAb [JE44-93]	HA722242	Human, Mouse, Rat	WB
	PI3 Kinase p110 beta Recombinant Rabbit mAb [JE65-38]	HA722474	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	PI 3 Kinase p110 delta Recombinant Rabbit mAb [JM83-24]	ET1705-46	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	PI3-kinase p110 subunit beta Recombinant Rabbit mAb [SC0580]	ET1610-36	Human	WB,IP,IHC-P,FC
	PI3 Kinase p110 α Recombinant Rabbit mAb [SJ0186]	ET1606-36	Human, Mouse	WB,IF-Cell,IF-Tissue,IP,IHC-P
	PI 3 Kinase p55 gamma Recombinant Rabbit mAb [SR44-06]	ET1602-27	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
KD	PI 3 Kinase p85 alpha Mouse mAb [L3-D6]	EM1701-62	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
	PI 3 Kinase p85 alpha Recombinant Rabbit mAb [SU04-07]	ET1608-70	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	PI 3 Kinase p85 alpha Recombinant Rabbit mAb [SU04-07] - BSA and Azide free	HA750162	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	PI 3 Kinase p85 alpha Recombinant Mouse mAb [A3-D0-R] - BSA and Azide free	HA610084	Human	WB,IF-Cell
KD	PI 3 Kinase p85 alpha Recombinant Mouse mAb [A3-D0-R]	HA601206	Human	WB,IF-Cell
	PI 3 Kinase p85 alpha Rabbit pAb	R1312-6	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	PI 3 Kinase p85 alpha Mouse mAb [A3-D0]	M1510-2	Human	WB,IF-Cell
	PI3 Kinase p85 alpha Recombinant Rabbit mAb [JE04-11]	HA722522	Human, Mouse, Rat	WB
	PI 3 Kinase p85 alpha/gamma Rabbit pAb	ER64588	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	PI 3 Kinase p85 beta Recombinant Rabbit mAb [ST04-77]	ET1609-30	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	PLCG 2 Recombinant Rabbit mAb [JE39-12]	HA721477	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	PLCG2 Rabbit pAb	R1512-4	Human	WB,IF-Cell,IHC-P,FC
	PPP3CB Recombinant Rabbit mAb [PSH0-21] - BSA and Azide free	HA750595	Human, Mouse, Rat	WB,IHC-P,ICC,FC
	PPP3CB Rabbit pAb	HA500370	Human, Mouse, Rat	WB,IHC-P
	PPP3CB Recombinant Rabbit mAb [PSH0-21]	HA721292	Human, Mouse, Rat	WB,IHC-P,ICC,FC
	PPP3CC Rabbit pAb	HA500101	Human, Mouse	WB
	PPP3CC Rabbit pAb	ER1915-38	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	PU.1/Spi1 Recombinant Rabbit mAb [JB72-35]	ET7107-49	Human	WB,IHC-P,FC
	PU.1 / Spi1 Recombinant Rabbit mAb [PD00-14]	HA721162	Human	WB,IHC-P
	RAC1 Recombinant Rabbit mAb [PSH05-09] - BSA and Azide free	HA750982	Human, Rat	WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	RAC1 Recombinant Rabbit mAb [PSH05-09]	HA722225	Human, Rat	WB,IF-Cell
	Rac1-2-3 Recombinant Rabbit mAb [JM11-29]	ET1703-80	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	RAC1/2/3 Rabbit pAb	HA500424	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	RAC2 Rabbit pAb	ER61092	Human, Mouse	WB
KD	RAC3 Recombinant Rabbit mAb [JE56-48]	HA721890	Human, Mouse, Rat	WB
	RAC3 Rabbit pAb	ER61093	Human, Mouse	WB
	RANK Rabbit pAb	ER1915-70	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	RANKL Rabbit pAb	HA500369	Human, Mouse, Rat	WB,IF-Cell
	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-59] - BSA and Azide free (Capture)	HA721836	Rat	ELISA(Cap)
	Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - BSA and Azide free (Detector)	HA721610	Rat	ELISA(Det)
	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60] - BSA and Azide free (Detector)	HA721837	Rat	ELISA(Det)
	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60]	HA722634B	Rat	ELISA(Det)
	Rat IL-1 beta Recombinant Rabbit mAb [PS01-86] - BSA and Azide free (Capture)	HA721609	Rat	ELISA(Cap)
	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - Detector	HA721611	Rat	ELISA(Det),ELISA
	RelB Recombinant Rabbit mAb [SD07-39]	ET1612-18	Human	WB,IF-Tissue,IHC-P,IP
	SERCA2 Recombinant Rabbit mAb [JM10-20]	ET1703-01	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	SERCA3 Rabbit pAb	ER1916-41	Human	WB
KD	SIRP alpha Recombinant Rabbit mAb [JE36-20]	HA721316	Human, Mouse, Rat	WB,IHC-P,FC
	SIRP alpha Mouse mAb [A3E3]	EM1902-37	Human	WB,IF-Cell,IHC-P,FC
	SIRP Alpha Rabbit pAb	ER1916-61	Human, Mouse, Rat	WB,FC
	SLP76 Rabbit pAb	ER65421	Human, Mouse, Rat	WB,IHC-P
Relative expression	SLP76 Recombinant Rabbit mAb [PSH09-22]	HA723080	Human	WB,IHC-P,IF-Cell,IP
Relative expression	SLP76 Recombinant Rabbit mAb [PSH09-22] - BSA and Azide free	HA751276	Human	WB,IHC-P,IF-Cell,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	SOCS1 Recombinant Rabbit mAb [PSH09-68]	HA723133	Human, Mouse, Rat, Monkey	WB
Relative expression	SOCS1 Recombinant Rabbit mAb [PSH09-68] - BSA and Azide free	HA751310	Human, Mouse, Rat, Monkey	WB
	SOCS1 Rabbit pAb	ER64997	Human, Mouse, Rat	IF-Cell, WB, IHC-P
	SOCS3 Rabbit pAb	ER65765	Human, Mouse, Rat	WB, IHC-P, FC
	STAT1 Recombinant Mouse mAb [G3-B11-R] - BSA and Azide free	HA610181	Human	WB, IHC-P
	STAT1 Rabbit pAb	R1408-2	Human, Mouse, Rat, Zebrafish	WB, IF-Cell
	STAT1 Mouse mAb [G3-B11]	M1407-1	Human	WB, IF-Cell, IHC-P, IF-Tissue
KD	STAT1 Recombinant Rabbit mAb [SD20-75]	ET1612-22	Human, Mouse, Rat	WB, IF-Cell, IF-Tissue, IHC-P, FC, IP
	STAT1 Recombinant Mouse mAb [G3-B11-R]	HA601309	Human	WB, IHC-P
KD	STAT1 alpha Recombinant Rabbit mAb [SJ01-89]	ET1606-39	Human, Mouse	WB, IP, FC, IHC-P, IF-Cell, IF-Tissue
KD	STAT1 alpha Recombinant Rabbit mAb [SJ01-89] - BSA and Azide free	HA750102	Human, Mouse	WB, IP, FC, IHC-P, IF-Cell, IF-Tissue
	STAT1 alpha/beta Rabbit pAb	ER31221	Human, Mouse, Rat	WB, IF-Cell, IHC-P
	STAT2 Recombinant Rabbit mAb [JF0884]	ET1702-18	Human, Mouse, Rat	WB, IF-Cell, IHC-P, FC
KD	Syk Recombinant Rabbit mAb [PSH0-76]	HA721494	Human	WB, IF-Cell
KD	Syk Recombinant Rabbit mAb [PSH0-76] - BSA and Azide free	HA750655	Human	WB, IF-Cell
	Syk Rabbit pAb	HA500080	Human, Rat	WB, IF-Cell, IHC-P
KD	TAB1 Rabbit pAb	HA500400	Human, Mouse, Rat	WB, IHC-P, IF-Cell, FC
	TAB1 Mouse mAb [1G2]	RT1603	Human	WB, IP, IF
KD	TAK1 Recombinant Rabbit mAb [JM73-19]	ET1705-14	Human, Mouse, Rat, Monkey	WB
KD	TAK1 Recombinant Rabbit mAb [JM73-19] - BSA and Azide free	HA750424	Human, Mouse, Rat, Monkey	WB
	TAK1 Rabbit pAb	ER1802-90	Human, Mouse, Rat	IF-Cell, IHC-P, FC
	Tartrate-resistant acid phosphatase type 5 Recombinant Rabbit mAb [JE35-95]	HA721897	Human	WB, IHC-P
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17]	HA721143	Human, Mouse, Rat, Monkey	WB, IHC-P, IF-Tissue
	TGF beta 1 Mouse mAb [B2-F6]	EM1709-93	Human	WB, IHC-P, FC
	Biotin Conjugated TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - Detector	HA723727B	Human, Mouse, Rat	ELISA(Det), ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - BSA and Azide free (Detector)	HA723726	Human, Mouse, Rat	ELISA(Det)
	TGF beta 1 Recombinant Rabbit mAb [PSH15-21] - BSA and Azide free (Capture)	HA723725	Human, Mouse, Rat	ELISA(Cap)
	TGF beta 1 Recombinant Rabbit mAb [PSH18-59]	HA724012	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
	TGF beta 1 Rabbit pAb	HA500496	Human	WB
	TGF beta 1 Recombinant Rabbit mAb	IRS205RB	Human	mIHC
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17] - BSA and Azide free	HA750485	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	TGF beta 2 Rabbit pAb	ER1917-63	Human, Mouse	WB,IHC-P,IF-Tissue
	TGF beta 2 Rabbit pAb	ER1917-62	Human, Mouse, Rat	WB
	TGF beta Receptor I Recombinant Rabbit mAb [PSH14-14] - BSA and Azide free	HA751509	Human, Mouse, Rat	WB,IP
	TGF beta Receptor I Recombinant Rabbit mAb [PSH14-14]	HA723604	Human, Mouse, Rat	WB,IP
	TGF beta Receptor I Rabbit pAb	ER1917-65	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	TGF beta Receptor II Recombinant Rabbit mAb [JE59-90]	HA721480	Human, Rat	WB
	TGF beta Receptor II Recombinant Rabbit mAb [PSH01-50] - BSA and Azide free	HA750726	Human	WB
	TGF beta Receptor II Recombinant Rabbit mAb [PSH01-50]	HA721693	Human	WB
	TGF beta Receptor II Rabbit pAb	ER1917-66	Human, Mouse, Rat, Cow	IHC-P,IF-Tissue
	TGFβ RII Rabbit pAb	ER65957	Human, Mouse, Rat	IF-Cell,WB,IHC
	TNF alpha Rabbit pAb	ER65189	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	TNF alpha Rabbit pAb	ER1919-22	Human, Mouse	WB
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07]	HA722022	Human, Mouse	WB
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free	HA750897	Human, Mouse	WB
Cell treatment	TNF alpha Rabbit pAb	R1203-1	Mouse	WB,IF-Cell,ELISA
	TNFR1 Mouse mAb [2G2]	RT1624	Human, Mouse, Rat	WB,IP,IF,IHC-P
	TNF Receptor I Rabbit pAb	HA500140	Human	WB,IHC-P,FC
	TNFSF11 Mouse mAb [8A7B9]	EM1710-42	Human, Monkey	WB,IHC-P,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	TRAF2 Recombinant Rabbit mAb [SD205-06]	ET1612-5	Human, Mouse, Rat	WB,IHC-P,IP
	TRAF2 Recombinant Rabbit mAb [SD205-06] - BSA and Azide free	HA750297	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
	TRAF6 Rabbit pAb	R1311-2	Human, Mouse, Rat, Chicken, Zebrafish	IF-Cell,WB
	TRAF6 Recombinant Rabbit mAb [JE01-42]	HA722857	Human, Mouse, Rat	WB,IF-Cell,FC
	TREM2 Rabbit pAb	ER2001-05	Human	WB,IHC-P,FC
Relative expression	TREM2 Rabbit pAb	ER1902-96	Human	WB,IHC-P,FC
	TREM2 Rabbit pAb	ER1918-04	Human, Mouse	WB
	TREM2 Rabbit pAb	ER2001-59	Human, Rat	WB,FC
	TYK2 Rabbit pAb	ER1803-05	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	TYK2 Recombinant Rabbit mAb [PSH18-53]	HA724006	Human, Mouse, Rat, Monkey	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层