

心肌发育研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	Actin Recombinant Antibody [JM36-62] - Human IgG1 (Chimeric)	HA721878	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	HRP Conjugated Actin Recombinant Rabbit mAb [JJ09-29]	HA721186	Human, Mouse, Rat	WB
	Actin Recombinant Rabbit mAb [JF47-01]	ET1702-52	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IF-Cell,FC
	Actin Recombinant Antibody [JM36-62] - Human IgG1 (Chimeric) - BSA and Azide free	HA750844	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Actin Recombinant Rabbit mAb [JJ09-29] - BSA and Azide free	HA750323	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P
	iFluor™ 488 Conjugated Actin Recombinant Rabbit mAb [JJ09-29]	HA720160F	Human, Mouse, Rat	IF-Cell,IF-Tissue,FC,WB
	Actin Recombinant Rabbit mAb [JJ09-29]	ET1701-80	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P
	Actinin alpha 2/3 Rabbit pAb	ER62588	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Actinin alpha 3 Rabbit pAb	ER62589	Human, Mouse, Rat	WB,IF-Cell
	Activin Receptor Type IIB / ACVR2B Recombinant Rabbit mAb [JE33-52]	HA722725	Human, Mouse, Rat	WB
	ACVR2A Rabbit pAb	HA500084	Human, Mouse, Rat	WB
	ACVR2A Mouse mAb [A7C8]	HA601091	Human	WB
	ACVR2B Mouse mAb [A5H6]	HA600066	Human, Mouse, Rat, Zebrafish	WB
	ACVR2B Mouse mAb [A5H7]	HA600067	Human, Mouse, Rat, Zebrafish	WB
	ADORA2B Rabbit pAb	ER1903-44	Human, Mouse, Rat	WB
	ALK-1 Recombinant Rabbit mAb [JE30-65]	HA722529	Human, Mouse, Rat	WB,IHC-P
	ALK-1 Rabbit pAb	ER62643	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	alpha Actinin Rabbit pAb	ER1803-60	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	alpha Actinin Mouse mAb [A1G2]	EM1901-52	Human, Rat	WB,IF-Cell,IHC-P,FC
	alpha Actinin Recombinant Mouse mAb [A1G2-R]	HA601194	Human, Mouse, Rat	WB,IHC-P
	alpha Actinin Recombinant Mouse mAb [A1G2-R] - BSA and Azide free	HA610074	Human, Mouse, Rat	WB,IHC-P

	Alpha Actinin 2 Recombinant Rabbit mAb [JF0980]	ET1702-30	Human, Mouse, Rat	WB,IHC-P,IP
	Alpha-cardiac actin Recombinant Mouse mAb [33-32-R]	HA601208	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Alpha-cardiac actin Mouse mAb [33-32]	M1206-1	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,IF-Tissue,IHC-Fr
	Alpha-cardiac actin Rabbit pAb	0407-3	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P
	Alpha-cardiac actin Recombinant Mouse mAb [33-32-R] - BSA and Azide free	HA610086	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	alpha sarcoglycan Recombinant Rabbit mAb [JA51-81]	ET1704-25	Human, Mouse, Rat	WB,IP,IHC-P,IF-Tissue
KD	alpha smooth muscle Actin Recombinant Rabbit mAb [SY25-03]	ET1607-43	Human, Mouse, Rat	WB,IHC-P,FC,IP
	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610362	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	alpha smooth muscle Actin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS094RT	Human	mIHC
	alpha smooth muscle Actin Recombinant Rabbit mAb	IRS048	Mouse	mIHC
	alpha smooth muscle Actin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS093MS	Human	mIHC
KD, Relative expression	alpha smooth muscle Actin Recombinant Rabbit mAb [SY02-64] - BSA and Azide free	HA750126	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Mouse IgG1 (Chimeric)	HA601546	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
KD, Relative expression	alpha smooth muscle Actin Recombinant Rabbit mAb [SY02-64]	ET1607-53	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,mIHC
KD, Relative expression	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Rat IgG1 (Chimeric)	HA601427	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD, Relative expression	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610267	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	alpha smooth muscle Actin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS220RT	Mouse	mIHC
	alpha smooth muscle Actin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS219MS	Mouse	mIHC
	alpha smooth muscle Actin Recombinant Rabbit mAb	IRS012	Human	mIHC
	AMH Rabbit pAb	HA500137	Human, Mouse, Rat	WB,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	AMH Recombinant Rabbit mAb [PSH13-48]	HA723537	Mouse, Rat, Human	WB,IHC-Fr,IHC-P
Relative expression	AMH Recombinant Rabbit mAb [PSH13-48] - BSA and Azide free	HA751467	Mouse, Rat, Human	WB,IHC-Fr,IHC-P
	AMHR2 Rabbit pAb	HA500309	Human, Mouse, Rat	WB,IHC-P
	AMHR2 Rabbit pAb	ER1903-64	Human, Mouse, Rat	WB
	AP-1/Jun D Rabbit pAb	ER62677	Human, Mouse, Rat	WB,IHC-P
	ATF2 Recombinant Rabbit mAb [SZ17-01]	ET1603-15	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	ATP1A2 Recombinant Rabbit mAb [JE34-00]	HA721422	Human, Mouse, Rat	WB,IHC-P,FC
	ATP1A3 Recombinant Rabbit mAb [JE35-24]	HA721596	Human, Mouse, Rat	WB
Relative expression	ATP1B4 Rabbit pAb	HA500311	Human, Mouse, Rat	WB
	ATP2A1 / SERCA1 Recombinant Rabbit mAb [JE45-75]	ET7109-43	Human	WB,IP,IHC-P
Relative expression	ATP2A1 / SERCA1 Rabbit pAb	ER1803-34	Human, Mouse, Rat	WB,IHC-P,FC
	ATPase Na ⁺ / K ⁺ beta 2 Rabbit pAb	ER1904-02	Mouse, Rat	WB
	beta 1 Sodium Potassium ATPase Rabbit pAb	ER60078	Mouse, Rat	WB
	BMP10 Rabbit pAb	HA500312	Human	WB
	BMP2 Recombinant Rabbit mAb [JE10-29]	HA722693	Human, Mouse, Rat	WB
	BMP2 Rabbit pAb	ER80602	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	BMP-3 Rabbit pAb	ER62817	Human, Mouse, Rat	WB,IHC-P
	BMP4 Recombinant Rabbit mAb [JE42-44]	HA721953	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	BMP4 Mouse mAb [1C-17]	EM1706-23	Human, Rat	WB,IF-Cell,FC
	BMP7 Rabbit pAb	ER1904-26	Human, Mouse, Rat, Rabbit	WB,IHC-P,IF-Tissue
	BMP-9 Rabbit pAb	ER65787	Human, Mouse	IHC-P
	BMPR1A Rabbit pAb	R1510-1	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	BMPR2 Mouse mAb [F8-G6]	EM1709-92	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P
	CACNA1C Rabbit pAb	ER1803-49	Human, Mouse, Rat	Dot Blot,IF-Cell,IHC-P,FC
	Cacna2d1 Rabbit pAb	ER62862	Human, Mouse, Rat	WB
	Cacna2d1 Rabbit pAb	ER60097	Rat	WB,IHC-P
	Cacna2d4 Rabbit pAb	ER62863	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CACNB1 Rabbit pAb	ER1803-29	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	CACNB3 Rabbit pAb	ER1803-24	Human, Mouse, Rat	WB,IHC-P,FC
	CACNG1 Rabbit pAb	ER1803-27	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	CACNG2 Rabbit pAb	ER1901-89	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	CACNG3 Rabbit pAb	ER1803-38	Human, Rat	WB,IF-Cell,IHC-P,FC
	Cacng4 Rabbit pAb	ER1803-13	Human, Rat	WB,IF-Cell,IHC-P,FC
	CACNG5 Rabbit pAb	ER1803-28	Human, Rat	WB,IF-Cell,IHC-P,FC
	Calcineurin A Recombinant Rabbit mAb [JA64-11]	ET1704-02	Human, Mouse, Rat	WB,FC
Relative expression	Cardiac Troponin I Recombinant Rabbit mAb [JF0550]	ET1702-37	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Cardiac Troponin T Mouse mAb [5-C7]	M1302-1	Human, Mouse	WB,IF-Cell,IHC-P,FC
	Cardiac Troponin T Recombinant Rabbit mAb [SC64-03] - BSA and Azide free	HA750222	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Cardiac Troponin T Mouse mAb [3-G6]	EM1701-39	Human, Mouse	WB,IHC-P,IF-Cell,IF-Tissue
	Cardiac Troponin T Rabbit pAb	ER21102	Human, Mouse, Rat	WB,IHC-P
	Cardiac Troponin T Recombinant Rabbit mAb [SC64-03]	ET1610-51	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Cav1.3 Rabbit pAb	ER60100	Human, Rat, Mouse	IHC-P
	CaV alpha 2 delta 2 Rabbit pAb	ER60098	Human, Rat, Mouse	WB,IHC-P
	CaV alpha 2 delta 3 Rabbit pAb	ER60104	Human, Rat, Mouse	WB
	CBP Rabbit pAb	ER62918	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	CCG8 Rabbit pAb	ER61403	Human, Mouse, Rat	WB
	CD298 Rabbit pAb	ER62941	Human	WB
	Cer1 Rabbit pAb	HA500223	Human, Mouse, Rat	WB,IHC-P
Cell treatment	c-Fos Recombinant Antibody [PSH07-51] - Mouse IgG1 (Chimeric)	HA601440	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Cell,WB
	c-Fos Recombinant Rabbit mAb	IRS131RB	Mouse	mIHC
	c-Fos Recombinant Rabbit mAb [PSH07-51] - BSA and Azide free	HA751172	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IF-Cell,IHC-P,ChIP,IF-Tissue,IP
Cell treatment	c-Fos Recombinant Antibody [PSH07-51] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610275	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Cell,WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	c-Fos Recombinant Rabbit mAb [PSH07-51]	HA722839	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IF-Cell,IHC-P,ChIP,IF-Tissue,IP
Cell treatment	c-Fos Recombinant Antibody [PSH05-77] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610219	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Tissue,IF-Cell,WB
Cell treatment	c-Fos Recombinant Rabbit mAb [JJ0938] - BSA and Azide free	HA750329	Human, Rat	WB,IF-Cell,IHC-P
Cell treatment	c-Fos Recombinant Rabbit mAb [PSH05-77]	HA722666	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,IHC-Fr,mIHC,IF-Tissue
Cell treatment	c-Fos Recombinant Antibody [PSH05-77] - Rat IgG1 (Chimeric)	HA601373	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Tissue,IF-Cell,WB
Cell treatment	c-Fos Recombinant Rabbit mAb [PSH05-77] - BSA and Azide free	HA751071	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,IHC-Fr,IF-Tissue
Cell treatment	c-Fos Recombinant Rabbit mAb [JJ0938]	ET1701-95	Human, Rat	WB,IF-Cell,IHC-P
Cell treatment	c-Fos Recombinant Antibody [PSH05-77] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610220	Mouse, Rat, Human, Cynomolgus monkey, Pig	IHC-Fr,WB
	c-Fos Mouse mAb [A10-G2]	EM1710-29	Human	WB,IHC-P,FC
Cell treatment	c-Fos Recombinant Antibody [PSH05-77] - Guinea pig IgG2 (Chimeric)	HA601374	Mouse, Rat, Human, Cynomolgus monkey, Pig	IHC-Fr,WB
	c-Jun Recombinant Rabbit mAb [SY0298]	ET1608-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	c-Jun Recombinant Rabbit mAb [SY0298] - BSA and Azide free	HA750143	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	Connexin-37 Rabbit pAb	ER1901-85	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Connexin-40 Rabbit pAb	ER1906-57	Human, Mouse	WB,IHC-P,IF-Tissue
	Connexin 43 Recombinant Rabbit mAb [PSH0-32] - BSA and Azide free	HA750606	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
	Connexin 43 Recombinant Rabbit mAb [PSH0-32]	HA721312	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
	COX1 Rabbit pAb	ER60256	Human, Mouse	WB
	COX3 Rabbit pAb	ER1906-59	Mouse, Rat	WB
	COX4I2 Rabbit pAb	ER63102	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	COX5a Rabbit pAb	ER63103	Human, Mouse, Rat	WB,IHC-P
	COX5A Recombinant Rabbit mAb [PSH13-57]	HA723546	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	COX5A Recombinant Rabbit mAb [PSH13-57] - BSA and Azide free	HA751476	Human, Mouse, Rat	WB,IHC-P
	COX5B Rabbit pAb	ER1902-87	Human	WB,IHC-P,FC
	COX5B Rabbit pAb	HA500163	Human, Rat	WB,IHC-P
	COX6B2 Rabbit pAb	ER65927	Human	IHC-P
	COX6c Rabbit pAb	ER63105	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	COX7a2/3 Rabbit pAb	ER63106	Human, Mouse, Rat	WB,IHC-P
	COX7B Rabbit pAb	ER65924	Human	IHC-P
	COX7R Rabbit pAb	ER61571	Human, Mouse	WB
	COX82 Rabbit pAb	ER65925	Human	IHC-P
	COX IV Recombinant Rabbit mAb [JJ09-05] - BSA and Azide free	HA750317	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	COX IV Rabbit pAb	R1706-1	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	COX IV Recombinant Rabbit mAb [JJ09-05]	ET1701-63	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	COX IV Rabbit pAb	RT1158	Human, Mouse, Rat	WB,IP,IF,IHC-P
	CREB Recombinant Rabbit mAb [PSH15-48] - BSA and Azide free	HA751563	Human, Mouse, Rat, Monkey	WB,IHC-P,IHC-Fr,IF-Cell,FC
	CREB Recombinant Rabbit mAb [PSH15-48]	HA723754	Human, Mouse, Rat, Monkey	WB,IHC-P,IHC-Fr,IF-Cell,FC
	CREB Recombinant Rabbit mAb [SA04-04]	ET1601-15	Human, Mouse, Zebrafish, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	CREB Recombinant Rabbit mAb [SA04-04] - BSA and Azide free	HA750007	Human, Mouse, Zebrafish, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	CREB1 Rabbit pAb	R1309-14	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	CREBBP Recombinant Rabbit mAb [PS01-16] - BSA and Azide free	HA751014	Human, Mouse	IHC-P
	CREBBP Rabbit pAb	HA500056	Human, Mouse, Rat	WB,IHC-P
	CREBBP Rabbit pAb	HA500191	Human	WB,IHC-P
	CREBBP Recombinant Rabbit mAb [PS01-16]	HA722297	Human, Mouse	IHC-P
	Cripto1 Recombinant Rabbit mAb [JE33-31]	HA721587	Human	WB,IF-Cell,FC
	Cripto1 Recombinant Rabbit mAb [JE48-67]	HA721586	Human	WB
	Cx40 / GJA5 Recombinant Rabbit mAb [PSH10-82]	HA723249	Human, Mouse, Rat, Monkey	WB,IP
	Cx40 / GJA5 Recombinant Rabbit mAb [PSH10-82] - BSA and Azide free	HA751374	Human, Mouse, Rat, Monkey	WB,IP
	CX6A1 Rabbit pAb	ER60268	Human, Rat	WB
	CX6A2 Rabbit pAb	ER61572	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Cytochrome B Rabbit pAb	ER1907-11	Human, Mouse, Rat	WB,FC
	Cytochrome b-c1 complex subunit 9/UQCR10 Rabbit pAb	HA500222	Human	WB,IHC-P,IF-Cell,FC
	Cytochrome c1 Rabbit pAb	ER63200	Human, Mouse	WB
	Desmin Rabbit pAb	ER40113	Human, Mouse, Zebrafish	WB,IF-Cell,IHC-P
KO	Desmin Mouse mAb [3-F7]	M1501-10	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Desmin Recombinant Rabbit mAb [SI18-00]	ET1606-30	Human, Mouse, Rat, Zebrafish, Chicken	WB,IF-Cell,IF-Tissue,IHC-P,FC,mIHC
	Desmin Recombinant Rabbit mAb [SI18-00] - BSA and Azide free	HA750100	Human, Mouse, Rat, Zebrafish, Chicken	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Desmin Recombinant mAb	IRS140RB	Mouse	mIHC
	Endofin Rabbit pAb	HA500231	Human, Mouse	WB,IHC-P
Relative expression	Fast Myosin Skeletal Heavy chain Recombinant Rabbit mAb [PSH05-12]	HA722228	Human, Mouse, Rat	IHC-P,IF-Tissue
Relative expression	Fast Myosin Skeletal Heavy chain Recombinant Rabbit mAb [PSH05-12] - BSA and Azide free	HA750985	Human, Mouse, Rat	IHC-P,IF-Tissue
	Fast Myosin Skeletal Heavy chain Recombinant mAb	IRS145RB	Mouse	mIHC
Relative expression	GATA1 Recombinant Rabbit mAb [JA10-32]	ET1704-41	Human, Mouse	WB,IHC-P
	GATA4 Recombinant Rabbit mAb [JA11-41]	ET1704-32	Human, Mouse, Rat	WB,IF-Cell,FC
	GATA4 Mouse mAb [7-4]	M1306-2	Human	WB,IHC-P
	GATA4 Rabbit pAb	R1107-5	Human, Mouse	WB,FC
	GATA-6 Rabbit pAb	ER63509	Human, Mouse, Rat	WB
	HCN4 Rabbit pAb	ER1910-38	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
Relative expression	heavy chain Myosin / MYH3 Recombinant Rabbit mAb [PSH18-12]	HA723975	Human, Mouse, Rat	WB,IHC-P,IP
	HIF-1 alpha Mouse mAb [5H7]	RT1278	Human, Mouse, Rat	WB,IP,ICC,IHC-P
Cell treatment	HIF-1 alpha Recombinant Rabbit mAb [JE75-33]	HA721997	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	HIF-1 alpha Recombinant Rabbit mAb [PSH07-00]	HA722778	Human, Mouse	WB,IF-Cell,IHC-P
Cell treatment	HIF-1 alpha Recombinant Rabbit mAb [PSH07-00] - BSA and Azide free	HA751127	Human, Mouse	WB,IF-Cell,IHC-P
Cell treatment	HIF-1 alpha Rabbit pAb	ER1802-41	Human	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human SIRP alpha Recombinant Rabbit mAb [PSH17-64] - BSA and Azide free (Detector)	HA725348	Human	ELISA(Det)
	Biotin Conjugated Human SIRP alpha Recombinant Rabbit mAb [PSH17-64] - Detector	HA725348B	Human	ELISA(Det),ELISA
	Human SIRP alpha Recombinant Rabbit mAb [PSH17-63] - BSA and Azide free (Capture)	HA725347	Human	ELISA(Cap)
	JNK1 Recombinant Rabbit mAb [JE63-04]	HA721899	Human, Mouse, Rat	WB,IHC-P
	JNK1 Rabbit pAb	R1309-1	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	JNK1+JNK2+JNK3 Recombinant Rabbit mAb [SA43-06] - BSA and Azide free	HA750018	Human, Mouse, Rat, Cow, Monkey, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP
	JNK1+JNK2+JNK3 Recombinant Rabbit mAb [SA43-06]	ET1601-28	Human, Mouse, Rat, Cow, Monkey, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP
	JNK1+JNK3 Recombinant Rabbit mAb [JJ08-86]	ET1701-47	Human, Mouse, Rat	WB,IP,IHC-P
	JNK2 Rabbit pAb	HA500110	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	JNK2 Recombinant Rabbit mAb [SC05-85]	ET1610-11	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	JNK3 Recombinant Rabbit mAb [SD082-09]	ET1612-68	Human, Mouse, Rat	WB,IF-Tissue,IHC-P
	JunD Recombinant Rabbit mAb [PSH08-85]	HA723037	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP,ChIP
	JunD Recombinant Rabbit mAb [JE59-51]	HA722913	Human, Mouse, Rat	WB,IF-Cell,IHC-P,ChIP,IP
	JunD Recombinant Rabbit mAb [PSH08-85] - BSA and Azide free	HA751250	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP,ChIP
	JunD Recombinant Rabbit mAb [SD0830]	ET1612-92	Human	WB,IF-Cell,IF-Tissue,IP,IHC-P,FC
	KAT3B / p300 Recombinant Rabbit mAb [PSH15-31]	HA723737	Human, Mouse, Rat, Monkey	WB,IF-Cell
	KAT3B / p300 Recombinant Rabbit mAb [PSH15-30]	HA723736	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	KAT3B / p300 Recombinant Rabbit mAb [PSH15-31] - BSA and Azide free	HA751553	Human, Mouse, Rat, Monkey	WB,IF-Cell
	KAT3B / p300 Recombinant Rabbit mAb [PSH15-30] - BSA and Azide free	HA751552	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	KAT3B / p300 Rabbit pAb	ER1908-18	Human, Mouse, Rat, Monkey	WB
	Kir2.1 Rabbit pAb	ER1911-91	Mouse, Rat	WB,FC
	Kir2.1 Recombinant Rabbit mAb [JE54-56]	HA721227	Human, Mouse, Rat	WB,IHC-P
	LRG1 Rabbit pAb	HA500576	Human, Mouse	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	LRG1 Recombinant Rabbit mAb [JB34-31]	ET7107-11	Human	WB,IHC-P,IP
	L-type Ca++ CP gamma7 Rabbit pAb	ER63880	Human, Mouse, Rat	WB
	MADH7/Smad7 Rabbit pAb	ER1912-61	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
	MAP3K7 Mouse mAb [1G5]	RT1377	Human, Mouse, Rat	WB,IP,IF
	MEF2A Recombinant Rabbit mAb [JA33-04]	ET1704-58	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	MEF2A + MEF2C Recombinant Rabbit mAb [JE63-26]	HA721030	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
Relative expression	MEF2C Recombinant Rabbit mAb [JE35-45]	HA721298	Human, Rat	WB,IF-Cell,IHC-P,IP,FC
	MK13 Rabbit pAb	ER60146	Human, Mouse, Rat	WB
	MTCO1 Recombinant Rabbit mAb [PSH07-50] - BSA and Azide free	HA751171	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
	MTCO1 Rabbit pAb	HA500517	Human, Mouse	WB,IHC-P
	MTCO1 Recombinant Rabbit mAb [PSH07-50]	HA722838	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
	MTCO2 Recombinant Rabbit mAb [SC06-88]	ET1610-72	Human	WB,IF-Cell,IF-Tissue,IHC-P
Relative expression	MYH2 Recombinant Rabbit mAb [JE53-30]	HA721862	Human, Mouse, Rat	WB,IHC-P
	MYH3 Rabbit pAb	ER60864	Human, Rat, Mouse	WB
	MYH4 Rabbit pAb	ER64077	Human	IHC-P
	MYH7B Recombinant Rabbit mAb [JE54-44]	ET7110-69	Rat, Mouse, Human	WB
	MYH-pan Rabbit pAb	ER65575	Human, Mouse, Rat	WB,IHC-P
	MYL3 Rabbit pAb	ER64078	Human	WB
	MYL4 Rabbit pAb	HA500471	Human, Mouse, Rat	WB
	Myosin heavy chain Recombinant Rabbit mAb [JF097-7]	ET1702-88	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Myosin Light Chain 2 Rabbit pAb	R1512-19	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Myosin Light Chain 2 Recombinant Rabbit mAb [SN67-09] - BSA and Azide free	HA750244	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IP
	Myosin Light Chain 2 Recombinant Rabbit mAb [SN67-09]	ET1611-13	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IP
	Na+/K+ ATPase Mouse mAb [1G1]	RT1412	Human, Mouse, Rat	WB,IP,IF,IHC-P
	Na+/K+-ATPase alpha 1 Rabbit pAb	ER64103	Rat	WB,IF-Cell
	Na+/K+-ATPase alpha 2 Rabbit pAb	ER64104	Human, Mouse, Rat, Monkey	WB
	NCX1 Rabbit pAb	ER65446	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	NCX1 Recombinant Rabbit mAb [PSH20-37]	HA724148	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue
	NF-kB p105 / p50 Mouse mAb [A2-7]	M1505-8	Human, Mouse	WB,IF-Cell,IHC-P
KD	NF-kB p105 / p50 Recombinant Rabbit mAb [SZ20-01] - BSA and Azide free	HA750060	Human, Mouse, Rat	WB,IHC-P,ICC,FC,IF-Tissue
KD	NF-kB p105 / p50 Recombinant Rabbit mAb [SZ20-01]	ET1603-18	Human, Mouse, Rat	WB,IHC-P,ICC,FC,IF-Tissue
	NF-kB p105 / p50 Rabbit pAb	R1309-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	NF-kB p65 Rabbit pAb	HA500402	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	NF-kB p65 Rabbit pAb	ER0815	Human, Mouse, Rat, Zebrafish	WB,IHC-P,FC
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27] - BSA and Azide free	HA750601	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27]	HA721307	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04]	ET1603-12	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04] - BSA and Azide free	HA750058	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	NHE-1 Rabbit pAb	ER1902-51	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	NHE-6 Rabbit pAb	ER64186	Human, Mouse	WB
	NHE-7 Rabbit pAb	ER64187	Human, Mouse, Rat	WB
	Noggin Recombinant Rabbit mAb [PSH10-62]	HA723227	Human, Mouse, Rat	WB,IHC-P
	Noggin Rabbit pAb	0805-2	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Noggin Recombinant Rabbit mAb [PSH10-62] - BSA and Azide free	HA751365	Human, Mouse, Rat	WB,IHC-P
	Noggin Rabbit pAb	ER1802-30	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	p38 alpha / MAPK14 Rabbit pAb	R1308-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD, Cell treatment	p38 alpha / MAPK14 Recombinant Rabbit mAb [SR43-04] - BSA and Azide free	HA750040	Human, Mouse, Rat	WB,IF-Cell,FC
	p38 alpha / MAPK14 Mouse mAb [9-B9]	EM40307	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD, Cell treatment	p38 alpha / MAPK14 Recombinant Rabbit mAb [SR43-04]	ET1602-26	Human, Mouse, Rat	WB,IF-Cell,FC
KO	p38 alpha / MAPK14 Recombinant Rabbit mAb [JF55-07]	ET1702-65	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	p38 alpha/MAPK14 + p38 beta/MAPK11 + p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-60]	HA722019	Human, Mouse, Rat, Monkey	WB
	p38 alpha/MAPK14 + p38 beta/MAPK11 + p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-60] - BSA and Azide free	HA750894	Human, Mouse, Rat, Monkey	WB
	p38 beta/MAPK11 + p38 gamma/MAPK12 + p38 delta/MAPK13 Recombinant Rabbit mAb [PSH03-59] - BSA and Azide free	HA750893	Human, Mouse, Rat, Monkey	WB
	p38 beta/MAPK11 + p38 gamma/MAPK12 + p38 delta/MAPK13 Recombinant Rabbit mAb [PSH03-59]	HA722018	Human, Mouse, Rat, Monkey	WB
	p38 delta / MAPK13 Rabbit pAb	HA500582	Human	WB
	p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-58]	HA722017	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-58] - BSA and Azide free	HA750892	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	p38 MAPK Recombinant Rabbit Multiclonal Antibody [PSH03-72]	HA722036	Human, Mouse, Rat, Monkey	WB
	p38 MAPK Recombinant Rabbit Multiclonal Antibody [PSH03-72] - BSA and Azide free	HA750906	Human, Mouse, Rat, Monkey	WB
	Pan-Actin Recombinant Mouse mAb [A2E1-R]	HA601287	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	Pan-Actin Recombinant Mouse mAb [A2E1-R] - BSA and Azide free	HA610160	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	Pan-Actin Mouse mAb [A2E1]	HA600032	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	PGC1 alpha+beta Recombinant Rabbit mAb [JF09-71]	ET1702-96	Human, Mouse, Rat	WB
	PGC1 beta Recombinant Rabbit mAb [JM33-24]	ET1705-65	Human, Rat	WB,IP,FC
Cell treatment	Phospho-ATF2 (T71) Recombinant Rabbit mAb [SC05-90]	ET1610-30	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IP,IHC-P
Cell treatment	Phospho-c-Jun (S63) Recombinant Rabbit mAb [SY0297] - BSA and Azide free	HA750146	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-c-Jun (S63) Recombinant Rabbit mAb [SY0297]	ET1608-4	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit mAb [JE65-46]	HA722475	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit mAb [PSH07-12] - BSA and Azide free	HA751133	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit mAb [PSH07-12]	HA722800	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-c-Jun (T91) Recombinant Rabbit mAb [JJ080-9]	ET1701-32	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Creb (S133) Recombinant Rabbit mAb [JB25-40]	ET7107-93	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,IF-Tissue
Cell treatment	Phospho-Creb (S133) Recombinant Rabbit mAb [JB25-40] - BSA and Azide free	HA750465	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,IF-Tissue
Cell treatment	Phospho-JNK1/2/3 (T183 + T183 + T221) Recombinant Rabbit mAb [ST500]	ET1609-42	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
Cell treatment	Phospho-JNK1/2/3 (T183 + T183 + T221) Recombinant Rabbit mAb [ST500] - BSA and Azide free	HA750173	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
Cell treatment	Phospho-JunD (S255) Recombinant Rabbit mAb [JJ08-21]	ET1701-35	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-NF-kB p65 (S529) Recombinant Rabbit mAb [SP07-00]	ET1604-27	Human	WB,IHC-P,FC,IP
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58]	HA723223	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58] - BSA and Azide free	HA751361	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-p38 alpha (T180) Recombinant Rabbit mAb [PSH06-25]	HA722657	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180) Recombinant Rabbit mAb [PSH06-25] - BSA and Azide free	HA751062	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180 + Y182) Recombinant Rabbit mAb [PSH06-35]	HA722669	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180 + Y182) Recombinant Rabbit mAb [PSH06-35] - BSA and Azide free	HA751074	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-p38 (T180) Rabbit pAb	ER2001-51	Human, Mouse, Rat	WB
Cell treatment	Phospho-p38 (T180) Recombinant Rabbit mAb [JE77-37]	HA722150	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-p38 (T180 + Y182) Rabbit pAb	ER2001-52	Human, Mouse, Rat	WB
	Phospho-p38 (T180 + Y182) Rabbit pAb	ER1903-01	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
Cell treatment	Phospho-Smad1/5/9 (S463/S465/S467) Recombinant Rabbit mAb [JE59-46]	HA722566	Human, Mouse, Rat	WB,IHC,IF-Cell
Cell treatment	Phospho-Smad5 (S463 + S465) Recombinant Rabbit mAb [SY09-03]	ET1605-5	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-Smad5 (S463 + S465) Recombinant Rabbit mAb [SY09-03] - BSA and Azide free	HA750089	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment, KD	Phospho-STAT3 (S727) Recombinant Rabbit mAb [SY24-09] - BSA and Azide free	HA750121	Human, Rat, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment, KD	Phospho-STAT3 (S727) Recombinant Rabbit mAb [SY24-09]	ET1607-39	Human, Rat, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-STAT3 (Y705) Recombinant Rabbit mAb [SZ43-01] - BSA and Azide free	HA750067	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-STAT3 (Y705) Recombinant Rabbit mAb [SZ43-01]	ET1603-40	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-TAK1 (S412) Recombinant Rabbit mAb [PSH14-64]	HA723662	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-TAK1 (S412) Recombinant Rabbit mAb [PSH14-64] - BSA and Azide free	HA751525	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-TAK1 (S439) Recombinant Rabbit mAb [JE44-93]	HA722242	Human, Mouse, Rat	WB
	PPAR alpha Mouse mAb	EM1707-71	Human, Monkey	FC
	PRDC Rabbit pAb	HA500225	Human	WB,FC
	Ryanodine receptor 2 / RYR-2 Rabbit pAb	ER1803-22	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	Ryanodine receptor 2 / RYR-2 Recombinant Rabbit mAb [PSH15-54] - BSA and Azide free	HA751567	Mouse, Rat	WB,IHC-P
Relative expression	Ryanodine receptor 2 / RYR-2 Recombinant Rabbit mAb [PSH15-54]	HA723760	Mouse, Rat	WB,IHC-P
	Sarcoglycan beta Rabbit pAb	ER64909	Human, Mouse	WB,IHC-P
	Sarcomeric Alpha Actinin Recombinant Rabbit mAb [JE45-22]	ET7109-42	Human, Mouse, Rat	WB,IP,IHC-P
	SERCA2 Recombinant Rabbit mAb [JM10-20]	ET1703-01	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	SERCA3 Rabbit pAb	ER1916-41	Human	WB
	SFT Rabbit pAb	HA500229	Human, Mouse, Rat	WB,IHC-P
	SFT Recombinant Rabbit mAb [JE64-46]	HA722578	Human, Mouse, Rat	WB,IHC-P,IP
	SIRP alpha Mouse mAb [A3E3]	EM1902-37	Human	WB,IF-Cell,IHC-P,FC
KD	SIRP alpha Recombinant Rabbit mAb [JE36-20]	HA721316	Human, Mouse, Rat	WB,IHC-P,FC
	SIRP Alpha Rabbit pAb	ER1916-61	Human, Mouse, Rat	WB,FC
	SKI Rabbit pAb	HA500220	Human, Mouse, Rat	WB,IHC-P
Relative expression	Slow Skeletal Myosin Heavy chain Recombinant Rabbit mAb [PSH06-58]	HA722721	Human, Mouse, Rat	WB,IHC-P
Relative expression	Slow Skeletal Myosin Heavy chain Recombinant Rabbit mAb [PSH07-54]	HA722735	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression	Slow Skeletal Myosin Heavy chain Recombinant Rabbit mAb [PSH07-54] - BSA and Azide free	HA751092	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	Slow Skeletal Myosin Heavy chain Recombinant Rabbit mAb [PSH06-58] - BSA and Azide free	HA751086	Human, Mouse, Rat	WB,IHC-P
KO	Smad1 Recombinant Rabbit mAb [SY0254]	ET1607-42	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
	Smad1 Recombinant Rabbit mAb [SD2090]	ET1612-39	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
KO	Smad1 Recombinant Rabbit mAb [SY0254] - BSA and Azide free	HA750124	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
	Smad1/5/9 Rabbit pAb	ER64980	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Smad4 Recombinant Rabbit mAb [PD01-37]	HA721234	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KO	Smad4 Recombinant Rabbit mAb [SP05-05] - BSA and Azide free	HA750070	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Smad4 Rabbit pAb	ER1706-90	Human, Mouse, Rat	IF-Cell,IHC-P,FC
KO	Smad4 Recombinant Rabbit mAb [SP05-05]	ET1604-12	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Smad5 Recombinant Rabbit mAb [SI17-06]	ET1606-26	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	SMAD6 Mouse mAb [H6-H3]	EM1709-56	Human	WB
	SMURF1 Recombinant Rabbit mAb [PSH10-12]	HA723181	Human	WB
	SMURF1 Recombinant Rabbit mAb [PSH10-12] - BSA and Azide free	HA751341	Human	WB
	SMURF1 Rabbit pAb	HA500057	Human, Mouse, Rat	WB
	SMURF 2 Recombinant Rabbit mAb [JU32-34]	ET7106-76	Human	WB
	Sodium hydrogen exchanger 1 Rabbit pAb	ER1916-95	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
KD	Sodium Potassium ATPase Recombinant Rabbit mAb [ST0533]	ET1609-76	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
	iFluor™ 488 Conjugated Sodium Potassium ATPase Recombinant Rabbit mAb [ST0533]	HA720174F	Human, Mouse, Rat	IF-Tissue
KD	Sodium Potassium ATPase Recombinant Rabbit mAb [ST0533] - BSA and Azide free	HA750194	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
	iFluor™ 594 Conjugated Sodium Potassium ATPase Recombinant Rabbit mAb [ST0533]	HA720175F	Human, Mouse, Rat	IF-Tissue
	iFluor™ 647 Conjugated Sodium Potassium ATPase Recombinant Rabbit mAb [ST0533]	HA720176F	Human, Mouse, Rat	IF-Tissue,FC
	SRF Recombinant Rabbit mAb [PSH05-06] - BSA and Azide free	HA750979	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	SRF Recombinant Rabbit mAb [PSH05-06]	HA722220	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	STAT3 Recombinant Rabbit mAb [SY34-01]	ET1605-45	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	STAT3 Recombinant Rabbit mAb [SY24-08]	ET1607-38	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
KD	STAT3 Recombinant Rabbit mAb [SY24-08] - BSA and Azide free	HA750120	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
KD	STAT3 Recombinant Rabbit mAb [SY34-01] - BSA and Azide free	HA750087	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	STAT3 Rabbit pAb	RT1591	Human, Mouse, Rat	WB,IP,IF,IHC-P
KD	TAB1 Rabbit pAb	HA500400	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	TAB1 Mouse mAb [1G2]	RT1603	Human	WB,IP,IF
KD	TAK1 Recombinant Rabbit mAb [JM73-19]	ET1705-14	Human, Mouse, Rat, Monkey	WB
	TAK1 Rabbit pAb	ER1802-90	Human, Mouse, Rat	IF-Cell,IHC-P,FC
KD	TAK1 Recombinant Rabbit mAb [JM73-19] - BSA and Azide free	HA750424	Human, Mouse, Rat, Monkey	WB
	TNNC1 Rabbit pAb	ER1902-25	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	TNNC1 Rabbit pAb	HA500407	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	TNNT1 Rabbit pAb	ER61345	Human, Mouse, Rat	WB
	Total OXPHOS Antibody Cocktail Recombinant Rabbit Multiclonal Antibody [PSH15-50]	HA723756	Human	WB,IHC-P
	Tropomyosin 1 Recombinant Rabbit mAb [JG38-66]	ET7108-74	Human, Mouse, Rat	WB,IHC-P,FC
KD	Tropomyosin 2 Recombinant Rabbit mAb [PSH02-63] - BSA and Azide free	HA750820	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell
KD	Tropomyosin 2 Recombinant Rabbit mAb [PSH02-63]	HA721841	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell
	UBE2D3 Rabbit pAb	HA500230	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Ubiquinol-Cytochrome C Reductase Core Protein I Recombinant Rabbit mAb [PSH01-23] - BSA and Azide free	HA750699	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell
	Ubiquinol-Cytochrome C Reductase Core Protein I Recombinant Rabbit mAb [PSH01-23]	HA721640	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell
	UQCRC2 Rabbit pAb	ER1803-08	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	UQCRC2 Recombinant Rabbit mAb [JE32-56]	HA721872	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	UQCRFS1 Recombinant Rabbit mAb [JG62-31]	ET7108-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	UQCRH Recombinant Rabbit mAb [JE62-25]	HA720113	Human, Mouse, Rat	WB,IHC-P,FC
	ZNF144 Rabbit pAb	ER1918-94	Human, Rat	WB,IHC-P,IF-Tissue