

心肌细胞标志-产品清单

高级验证	产品名	货号	反应种属	验证应用
	HRP Conjugated Actin Recombinant Rabbit mAb [JJ09-29]	HA721186	Human, Mouse, Rat	WB
	Actin Recombinant Rabbit mAb [JJ09-29]	ET1701-80	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P
	iFluor™ 488 Conjugated Actin Recombinant Rabbit mAb [JJ09-29]	HA720160F	Human, Mouse, Rat	IF-Cell,IF-Tissue,FC,WB
	Actin Recombinant Rabbit mAb [JF47-01]	ET1702-52	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IF-Cell,FC
	Actin Recombinant Antibody [JM36-62] - Human IgG1 (Chimeric) - BSA and Azide free	HA750844	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Actin Recombinant Rabbit mAb [JJ09-29] - BSA and Azide free	HA750323	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P
	Actin Recombinant Antibody [JM36-62] - Human IgG1 (Chimeric)	HA721878	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Actinin alpha 2/3 Rabbit pAb	ER62588	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Actinin alpha 3 Rabbit pAb	ER62589	Human, Mouse, Rat	WB,IF-Cell
	ADORA2B Rabbit pAb	ER1903-44	Human, Mouse, Rat	WB
	alpha Actinin Recombinant Mouse mAb [A1G2-R] - BSA and Azide free	HA610074	Human, Mouse, Rat	WB,IHC-P
	alpha Actinin Mouse mAb [A1G2]	EM1901-52	Human, Rat	WB,IF-Cell,IHC-P,FC
	alpha Actinin Rabbit pAb	ER1803-60	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	alpha Actinin Recombinant Mouse mAb [A1G2-R]	HA601194	Human, Mouse, Rat	WB,IHC-P
	Alpha Actinin 2 Recombinant Rabbit mAb [JF0980]	ET1702-30	Human, Mouse, Rat	WB,IHC-P,IP
	Alpha-cardiac actin Rabbit pAb	0407-3	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P
	Alpha-cardiac actin Recombinant Mouse mAb [33-32-R] - BSA and Azide free	HA610086	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Alpha-cardiac actin Recombinant Mouse mAb [33-32-R]	HA601208	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Alpha-cardiac actin Mouse mAb [33-32]	M1206-1	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,IF-Tissue,IHC-Fr
	alpha sarcoglycan Recombinant Rabbit mAb [JA51-81]	ET1704-25	Human, Mouse, Rat	WB,IP,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	alpha smooth muscle Actin Recombinant Rabbit mAb [SY25-03]	ET1607-43	Human, Mouse, Rat	WB,IHC-P,FC,IP
KD, Relative expression	alpha smooth muscle Actin Recombinant Rabbit mAb [SY02-64] - BSA and Azide free	HA750126	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	alpha smooth muscle Actin Recombinant Rabbit mAb	IRS048	Mouse	mIHC
KD, Relative expression	alpha smooth muscle Actin Recombinant Rabbit mAb [SY02-64]	ET1607-53	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,mIHC
	alpha smooth muscle Actin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS094RT	Human	mIHC
	alpha smooth muscle Actin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS219MS	Mouse	mIHC
KD, Relative expression	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Rat IgG1 (Chimeric)	HA601427	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Mouse IgG1 (Chimeric)	HA601546	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	alpha smooth muscle Actin Recombinant Rabbit mAb	IRS012	Human	mIHC
KD, Relative expression	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610267	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610362	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	alpha smooth muscle Actin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS220RT	Mouse	mIHC
	alpha smooth muscle Actin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS093MS	Human	mIHC
	ANP Recombinant Rabbit mAb [JM08-36]	ET1705-75	Human	WB,IP
	BMP4 Recombinant Rabbit mAb [JE42-44]	HA721953	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	BMP4 Mouse mAb [1C-17]	EM1706-23	Human, Rat	WB,IF-Cell,FC
	BNP Recombinant Rabbit mAb [JE30-94]	HA721946	Human, Mouse, Rat	WB,IHC-P
Relative expression	Cardiac Troponin I Recombinant Rabbit mAb [JF0550]	ET1702-37	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Cardiac Troponin T Mouse mAb [5-C7]	M1302-1	Human, Mouse	WB,IF-Cell,IHC-P,FC
	Cardiac Troponin T Recombinant Rabbit mAb [SC64-03] - BSA and Azide free	HA750222	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Cardiac Troponin T Recombinant Rabbit mAb [SC64-03]	ET1610-51	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Cardiac Troponin T Rabbit pAb	ER21102	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Cardiac Troponin T Mouse mAb [3-G6]	EM1701-39	Human, Mouse	WB,IHC-P,IF-Cell,IF-Tissue
	CD3 Recombinant Rabbit mAb [JE80-02]	HA720082	Human	WB,IHC-P,IF-Tissue,IP,mIHC,IF-Cell
	CD3 Recombinant Rabbit mAb	IRS022	Human	mIHC
	CD3 Recombinant Mouse mAb [PD01-43]	HA601110	Human	WB,IHC-P,IF-Tissue
Relative expression	CD3 Recombinant Rabbit mAb [PSH19-70]	HA724094	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	CD31 Recombinant Rabbit mAb	IRS117RB	Mouse	mIHC
Relative expression	CD31 Mouse mAb [7-A1]	M1511-8	Human, Mouse	WB,IF-Cell,IHC-P,mIHC
	CD31 Mouse mAb	IRS023	Human	mIHC
	CD31 Recombinant Rabbit mAb [SU03-59] - BSA and Azide free	HA750151	Human	WB,IHC-P,FC,IP,IF-Cell,IF-Tissue
Relative expression	CD31 Recombinant Mouse mAb [7-A1-R]	HA601240	Human	WB,IF-Cell,IHC-P
	CD31 Recombinant Rabbit mAb [PSH20-35]	HA724146	Human, Mouse, Rat	IHC-P
Relative expression	CD31 Recombinant Antibody [7-A1-R] - Rabbit IgG (Chimeric)	HA723897	Human	WB,IF-Cell,IHC-P,FC
	CD31 Recombinant Rabbit mAb [SU03-59]	ET1608-48	Human	WB,IHC-P,FC,IP,IF-Cell,IF-Tissue
Relative expression	CD31 Recombinant Mouse mAb [7-A1-R] - BSA and Azide free	HA610115	Human	WB,IF-Cell,IHC-P
	CD31 Recombinant Antibody - Rabbit IgG (Chimeric)	IRS116RB	Human	mIHC
Relative expression	CD31 Recombinant Antibody [7-A1-R] - Rabbit IgG (Chimeric) - BSA and Azide free	HA751615	Human	WB,IF-Cell,IHC-P,FC
	CD31 Recombinant Antibody [7-A1-R] - Rat IgG1 (Chimeric)	HA601522	Human	IHC-P
Relative expression	CD31 Rabbit pAb	ER31219	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
	CD31 Recombinant Rabbit mAb	IRS075	Mouse	mIHC
	iFluor™ 488 Conjugated CD31 Recombinant Rabbit mAb [SU03-59]	HA720220F	Human	IF-Tissue
	CD31 Recombinant Antibody [7-A1-R] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610349	Human	IHC-P
	CD3D Recombinant Rabbit mAb [JJ08-97]	ET1701-57	Human	WB,IF-Cell,IHC-P,FC,IP,IF-Tissue
	CD3E Rabbit pAb	HA500344	Human, Mouse, Rat	WB
	CD3 epsilon Recombinant Rabbit mAb [SY0239]	ET1607-29	Human	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD3 epsilon Recombinant Rabbit mAb [JE17-74]	ET7110-63	Human	WB,IHC-P,FC
	CD3 epsilon Recombinant Rabbit mAb [SY0239] - BSA and Azide free	HA750115	Human	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
	CD3 epsilon Rabbit pAb	ER80501	Human, Mouse, Rat	WB
	CD3G Recombinant Rabbit mAb [JB38-29]	ET7107-55	Human, Mouse	WB,IHC-P,FC,IP,IF-Cell,IF-Tissue
Relative expression	CD3 zeta Recombinant Rabbit mAb [JE09-37]	HA722692	Human, Rat	WB,IHC-P,IF-Tissue
	CD3 zeta Rabbit pAb	0805-7	Human	WB,IF-Cell,FC
	COL1A1 Recombinant Rabbit mAb [PSH06-20] - BSA and Azide free	HA751058	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP,IHC-Fr
	COL1A1 Recombinant Rabbit mAb [ST58-04] - BSA and Azide free	HA750188	Human, Cow	WB,IHC-P
	COL1A1 Recombinant mAb	IRS137RB	Human	mIHC
	COL1A1 Recombinant Rabbit mAb [ST58-04]	ET1609-68	Human, Cow	WB,IHC-P
	COL1A1 Recombinant Rabbit mAb [PSH06-20]	HA722517	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP,mIHC,IHC-Fr
	Connexin-37 Rabbit pAb	ER1901-85	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Connexin-40 Rabbit pAb	ER1906-57	Human, Mouse	WB,IHC-P,IF-Tissue
	Connexin 43 Recombinant Rabbit mAb [PSH0-32] - BSA and Azide free	HA750606	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
	Connexin 43 Recombinant Rabbit mAb [PSH0-32]	HA721312	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
	Cripto1 Recombinant Rabbit mAb [JE48-67]	HA721586	Human	WB
	Cripto1 Recombinant Rabbit mAb [JE33-31]	HA721587	Human	WB,IF-Cell,FC
Relative expression	CSF-1-R Recombinant Rabbit mAb [PSH09-77]	HA723142	Human, Mouse, Rat	WB,IHC-P
Relative expression	CSF-1-R Recombinant Rabbit mAb [PSH09-77] - BSA and Azide free	HA751319	Human, Mouse, Rat	WB,IHC-P
	CSF-1-R Recombinant Rabbit mAb [PD00-45]	HA721180	Human	IHC-P,FC
	CSF1R Mouse mAb [B9-B9]	EM1708-56	Human	WB,IF-Cell,IHC-P
	Cx40 / GJA5 Recombinant Rabbit mAb [PSH10-82]	HA723249	Human, Mouse, Rat, Monkey	WB,IP
	Cx40 / GJA5 Recombinant Rabbit mAb [PSH10-82] - BSA and Azide free	HA751374	Human, Mouse, Rat, Monkey	WB,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Desmin Recombinant Rabbit mAb [SI18-00] - BSA and Azide free	HA750100	Human, Mouse, Rat, Zebrafish, Chicken	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Desmin Recombinant Rabbit mAb [SI18-00]	ET1606-30	Human, Mouse, Rat, Zebrafish, Chicken	WB,IF-Cell,IF-Tissue,IHC-P,FC,mIHC
	Desmin Recombinant mAb	IRS140RB	Mouse	mIHC
	Desmin Rabbit pAb	ER40113	Human, Mouse, Zebrafish	WB,IF-Cell,IHC-P
KO	Desmin Mouse mAb [3-F7]	M1501-10	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	FABP3 Recombinant Rabbit mAb [JU73-03]	ET1706-50	Mouse, Rat	WB,IHC-P
	FABP3 Recombinant Rabbit mAb [JE53-37]	HA721859	Human, Mouse, Rat	WB,IHC-P
Relative expression	FABP4 Recombinant Rabbit mAb [JM10-99] - BSA and Azide free	HA750398	Human, Mouse, Rat	WB,IF-Cell,IHC-Fr
Relative expression	FABP4 Recombinant Rabbit mAb [JM10-99]	ET1703-98	Human, Mouse, Rat	WB,IF-Cell,IHC-Fr
Relative expression	Fast Myosin Skeletal Heavy chain Recombinant Rabbit mAb [PSH05-12]	HA722228	Human, Mouse, Rat	IHC-P,IF-Tissue
Relative expression	Fast Myosin Skeletal Heavy chain Recombinant Rabbit mAb [PSH05-12] - BSA and Azide free	HA750985	Human, Mouse, Rat	IHC-P,IF-Tissue
	Fast Myosin Skeletal Heavy chain Recombinant mAb	IRS145RB	Mouse	mIHC
	GATA4 Rabbit pAb	R1107-5	Human, Mouse	WB,FC
	GATA4 Mouse mAb [7-4]	M1306-2	Human	WB,IHC-P
	GATA4 Recombinant Rabbit mAb [JA11-41]	ET1704-32	Human, Mouse, Rat	WB,IF-Cell,FC
	GATA-6 Rabbit pAb	ER63509	Human, Mouse, Rat	WB
	HADHA Recombinant Rabbit mAb [PSH01-29]	HA721653	Human, Mouse, Rat, Monkey	WB
	HADHA Recombinant Rabbit mAb [PSH01-28]	HA721652	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell
	HADHA Recombinant Rabbit mAb [PSH01-28] - BSA and Azide free	HA750706	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell
	HADHA Recombinant Rabbit mAb [PSH01-29] - BSA and Azide free	HA750707	Human, Mouse, Rat, Monkey	WB
	HCN4 Rabbit pAb	ER1910-38	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
Relative expression	heavy chain Myosin / MYH3 Recombinant Rabbit mAb [PSH18-12]	HA723975	Human, Mouse, Rat	WB,IHC-P,IP
	H-FABP Recombinant Rabbit mAb [PSH16-98] - BSA and Azide free (Capture)	HA725177	Human, Mouse	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	H-FABP Recombinant Rabbit mAb [PSH16-99] - BSA and Azide free (Detector)	HA725178	Human, Mouse	ELISA(Det)
	Biotin Conjugated H-FABP Recombinant Rabbit mAb [PSH16-99] - Detector	HA725178B	Human, Mouse	ELISA(Det),ELISA
	PE Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600117F	Human	FC
	Percp-Cy5.5 Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA601548F	Human	FC
	Human CD3 Recombinant Mouse mAb [PSH04-92]	HA601337	Human	FC
	FITC Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600116F	Human	FC
	APC Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600118F	Human	FC
	Human CD31 Recombinant Rabbit mAb [PSH07-83] - BSA and Azide free (Detector)	HA722907	Human	ELISA(Det)
	Human CD31 Recombinant Rabbit mAb [PSH07-82] - BSA and Azide free (Capture)	HA722906	Human	ELISA(Cap)
	Biotin Conjugated Human CD31 Recombinant Rabbit mAb [PSH07-83] - Detector	HA722908B	Human	ELISA(Det)
	Human CSF-1-R Recombinant Rabbit mAb [PSH12-66] - BSA and Azide free (Capture/Detector)	HA723474	Human	ELISA(Cap),ELISA(Det)
	HRP Conjugated Human CSF-1-R Recombinant Rabbit mAb [PSH12-66] - Detector	HA723478H	Human	ELISA(Det),ELISA
	Human CSF-1-R Recombinant Rabbit mAb [PSH12-67] - BSA and Azide free (Detector)	HA723475	Human	ELISA(Det)
	Human CSF-1-R Recombinant Rabbit mAb [PSH12-68] - BSA and Azide free (Capture)	HA723477	Human	ELISA(Cap)
	Biotin Conjugated Human CSF-1-R Recombinant Rabbit mAb [PSH12-67] - Detector	HA723476B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human CSF-1-R Recombinant Rabbit mAb [PSH12-66] - Detector	HA723478B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human SIRP alpha Recombinant Rabbit mAb [PSH17-64] - Detector	HA725348B	Human	ELISA(Det),ELISA
	Human SIRP alpha Recombinant Rabbit mAb [PSH17-64] - BSA and Azide free (Detector)	HA725348	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human SIRP alpha Recombinant Rabbit mAb [PSH17-63] - BSA and Azide free (Capture)	HA725347	Human	ELISA(Cap)
	Kir2.1 Recombinant Rabbit mAb [JE54-56]	HA721227	Human, Mouse, Rat	WB,IHC-P
	Kir2.1 Rabbit pAb	ER1911-91	Mouse, Rat	WB,FC
	LRG1 Recombinant Rabbit mAb [JB34-31]	ET7107-11	Human	WB,IHC-P,IP
	LRG1 Rabbit pAb	HA500576	Human, Mouse	WB
	Mast Cell Tryptase Recombinant Rabbit mAb [SC68-07] - BSA and Azide free	HA750229	Human, Mouse	WB,IHC-P,IP
	Mast Cell Tryptase Recombinant Rabbit mAb [SC68-07]	ET1610-64	Human, Mouse	WB,IHC-P,IP,mIHC
	Mast Cell Tryptase Recombinant Rabbit mAb	IRS041	Human	mIHC
	Mast Cell Tryptase Recombinant Antibody [SC68-07] - Rat IgG1 (Chimeric)	HA601419		
	Mast Cell Tryptase Recombinant Antibody - Rat IgG1 (Chimeric)	IRS209RT	Human	mIHC
	Mast Cell Tryptase Recombinant Antibody [SC68-07] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610259	0	
	Mast Cell Tryptase Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS210MS	Human	mIHC
	MEF2A Recombinant Rabbit mAb [JA33-04]	ET1704-58	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	MEF2A + MEF2C Recombinant Rabbit mAb [JE63-26]	HA721030	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
Relative expression	MEF2C Recombinant Rabbit mAb [JE35-45]	HA721298	Human, Rat	WB,IF-Cell,IHC-P,IP,FC
	MLC-2 Rabbit pAb	ER63968	Human, Mouse, Rat	IHC-P,IF-Cell
	MYBPC3 Rabbit pAb	ER65870	Human, Mouse, Rat	WB
Relative expression	MYH2 Recombinant Rabbit mAb [JE53-30]	HA721862	Human, Mouse, Rat	WB,IHC-P
	MYH3 Rabbit pAb	ER60864	Human, Rat, Mouse	WB
	MYH4 Rabbit pAb	ER64077	Human	IHC-P
	MYH7B Recombinant Rabbit mAb [JE54-44]	ET7110-69	Rat, Mouse, Human	WB
	MYH-pan Rabbit pAb	ER65575	Human, Mouse, Rat	WB,IHC-P
	MYL12B Recombinant Rabbit mAb [JE62-13]	HA722355	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	MYL4 Rabbit pAb	HA500471	Human, Mouse, Rat	WB
	MYL7 Rabbit pAb	ER64080	Human, Mouse, Rat, Monkey	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	MYL9 Recombinant Rabbit mAb [PSH02-07] - BSA and Azide free	HA750777	Human, Mouse, Rat	WB,FC
	MYL9 Recombinant Rabbit mAb [PSH02-07]	HA721766	Human, Mouse, Rat	WB,FC
	MYL9 Recombinant Rabbit mAb [JE34-63]	HA721361	Human, Mouse, Rat	WB,IHC-P
	Myosin heavy chain Recombinant Rabbit mAb [JF097-7]	ET1702-88	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Myosin Light Chain 2 Recombinant Rabbit mAb [SN67-09] - BSA and Azide free	HA750244	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IP
	Myosin Light Chain 2 Rabbit pAb	R1512-19	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Myosin Light Chain 2 Recombinant Rabbit mAb [SN67-09]	ET1611-13	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IP
	Na+ CP-pan Rabbit pAb	ER64102	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Na+ CP type II alpha Rabbit pAb	ER64098	Human, Mouse, Rat	WB,IHC-P
	Na+ CP type IV alpha Rabbit pAb	ER64100	Human, Mouse, Rat	WB
	NaV1.7 Rabbit pAb	ER1706-57	Human, Rat	IF-Cell,IHC-P
	NaV1.7 Mouse mAb [1D6]	EM1701-35	Human	IHC-P,IF-Cell,FC,WB
	NaV1.7 Rabbit pAb	ER1706-44	Human	IF-Cell,IHC-P,FC
	Pan-Actin Mouse mAb [A2E1]	HA600032	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Pan-Actin Recombinant Mouse mAb [A2E1-R] - BSA and Azide free	HA610160	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	Pan-Actin Recombinant Mouse mAb [A2E1-R]	HA601287	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	PCDH7 Rabbit pAb	ER1914-84	Human, Mouse	WB,IHC-P,IF-Tissue
	Periostin Recombinant Rabbit mAb [JA63-40]	ET1704-65	Human	WB
	Phospho-Myosin light chain (S20) Rabbit pAb	ER1913-46	Human, Mouse	WB,IHC-P,IF-Tissue
	RYR1 Rabbit pAb	ER1803-19	Human, Mouse	Dot Blot,IHC-P,FC
	Sarcoglycan beta Rabbit pAb	ER64909	Human, Mouse	WB,IHC-P
	Sarcomeric Alpha Actinin Recombinant Rabbit mAb [JE45-22]	ET7109-42	Human, Mouse, Rat	WB,IP,IHC-P
	SCN1A Rabbit pAb	ER65886	Human, Mouse, Rat	IHC-P,WB
KD	SIRP alpha Recombinant Rabbit mAb [JE36-20]	HA721316	Human, Mouse, Rat	WB,IHC-P,FC
	SIRP alpha Mouse mAb [A3E3]	EM1902-37	Human	WB,IF-Cell,IHC-P,FC
	SIRP Alpha Rabbit pAb	ER1916-61	Human, Mouse, Rat	WB,FC
Relative expression	Slow Skeletal Myosin Heavy chain Recombinant Rabbit mAb [PSH06-58] - BSA and Azide free	HA751086	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	Slow Skeletal Myosin Heavy chain Recombinant Rabbit mAb [PSH07-54]	HA722735	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression	Slow Skeletal Myosin Heavy chain Recombinant Rabbit mAb [PSH07-54] - BSA and Azide free	HA751092	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression	Slow Skeletal Myosin Heavy chain Recombinant Rabbit mAb [PSH06-58]	HA722721	Human, Mouse, Rat	WB,IHC-P
	smooth muscle Myosin heavy chain 11 Recombinant Rabbit mAb	IRS176RB	Human	mIHC
	smooth muscle Myosin heavy chain 11 Recombinant Rabbit mAb [JE60-21]	HA722863	Human, Mouse, Rat	WB,IHC-P
	smooth muscle Myosin heavy chain 11 Recombinant Rabbit mAb [JA03-35]	ET1704-61	Human, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
	smooth muscle Myosin heavy chain 11 Recombinant Rabbit mAb [JA03-35] - BSA and Azide free	HA750417	Human, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
	SOX9 Recombinant Rabbit mAb	IRS199RB	Human	mIHC
Relative expression, KD	SOX9 Recombinant Rabbit mAb [SN74-09] - BSA and Azide free	HA750262	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-P,IP,FC,CUT&Tag-seq
Relative expression	SOX9 Recombinant Rabbit mAb [PSH13-59]	HA723548	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell,IP,FC,ChIP
Relative expression, KD	SOX9 Recombinant Rabbit mAb [SN74-09]	ET1611-56	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-P,IP,FC,CUT&Tag-seq,mIHC
	SOX9 Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS279MS	Human	mIHC
	SOX9 Recombinant Rabbit mAb	IRS241RB	Mouse	mIHC
	SOX9 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS280RT	Human	mIHC
Relative expression	SOX9 Recombinant Rabbit mAb [PSH13-59] - BSA and Azide free	HA751478	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell,IP,FC,ChIP
	TGF beta 1 Recombinant Rabbit mAb [PSH18-59]	HA724012	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
	TGF beta 1 Recombinant Rabbit mAb [PSH15-21] - BSA and Azide free (Capture)	HA723725	Human, Mouse, Rat	ELISA(Cap)
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17]	HA721143	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17] - BSA and Azide free	HA750485	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	TGF beta 1 Recombinant Rabbit mAb	IRS205RB	Human	mIHC
	TGF beta 1 Mouse mAb [B2-F6]	EM1709-93	Human	WB,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - BSA and Azide free (Detector)	HA723726	Human, Mouse, Rat	ELISA(Det)
	Biotin Conjugated TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - Detector	HA723727B	Human, Mouse, Rat	ELISA(Det),ELISA
	TGF beta 1 Rabbit pAb	HA500496	Human	WB
	TNNC1 Rabbit pAb	HA500407	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	TNNC1 Rabbit pAb	ER1902-25	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	TNNT1 Rabbit pAb	ER61345	Human, Mouse, Rat	WB
	TTN Rabbit pAb	ER65816	Human, Mouse, Rat	IHC-P
	Wilms Tumor Protein Recombinant Rabbit mAb	IRS072	Mouse	mIHC
KD	Wilms Tumor Protein Recombinant Rabbit mAb [SC06-41] - BSA and Azide free	HA750219	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
KD	Wilms Tumor Protein Recombinant Rabbit mAb [SC06-41]	ET1610-45	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层