

心肌肥大研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	Actin Recombinant Rabbit mAb [JJ09-29]	ET1701-80	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P
	Actin Recombinant Antibody [JM36-62] - Human IgG1 (Chimeric) - BSA and Azide free	HA750844	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Actin Recombinant Rabbit mAb [JF47-01]	ET1702-52	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IF-Cell,FC
	iFluor™ 488 Conjugated Actin Recombinant Rabbit mAb [JJ09-29]	HA720160F	Human, Mouse, Rat	IF-Cell,IF-Tissue,FC,WB
	Actin Recombinant Rabbit mAb [JJ09-29] - BSA and Azide free	HA750323	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P
	HRP Conjugated Actin Recombinant Rabbit mAb [JJ09-29]	HA721186	Human, Mouse, Rat	WB
	Actin Recombinant Antibody [JM36-62] - Human IgG1 (Chimeric)	HA721878	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Actinin alpha 2/3 Rabbit pAb	ER62588	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Actinin alpha 3 Rabbit pAb	ER62589	Human, Mouse, Rat	WB,IF-Cell
	ADORA2B Rabbit pAb	ER1903-44	Human, Mouse, Rat	WB
	AKAP 13 Rabbit pAb	ER62629	Human	WB,IHC-P
	Alpha Actinin 2 Recombinant Rabbit mAb [JF0980]	ET1702-30	Human, Mouse, Rat	WB,IHC-P,IP
	Alpha-cardiac actin Rabbit pAb	0407-3	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P
	Alpha-cardiac actin Recombinant Mouse mAb [33-32-R] - BSA and Azide free	HA610086	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Alpha-cardiac actin Mouse mAb [33-32]	M1206-1	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,IF-Tissue,IHC-Fr
	Alpha-cardiac actin Recombinant Mouse mAb [33-32-R]	HA601208	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Mouse IgG1 (Chimeric)	HA601546	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
KD, Relative expression	alpha smooth muscle Actin Recombinant Rabbit mAb [SY02-64] - BSA and Azide free	HA750126	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	alpha smooth muscle Actin Recombinant Rabbit mAb [SY25-03]	ET1607-43	Human, Mouse, Rat	WB,IHC-P,FC,IP
	alpha smooth muscle Actin Recombinant Rabbit mAb	IRS048	Mouse	mIHC

KD, Relative expression	alpha smooth muscle Actin Recombinant Rabbit mAb [SY02-64]	ET1607-53	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,mIHC
	alpha smooth muscle Actin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS094RT	Human	mIHC
	alpha smooth muscle Actin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS219MS	Mouse	mIHC
	alpha smooth muscle Actin Recombinant Rabbit mAb	IRS012	Human	mIHC
KD, Relative expression	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610267	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	alpha smooth muscle Actin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS220RT	Mouse	mIHC
	alpha smooth muscle Actin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS093MS	Human	mIHC
	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610362	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
KD, Relative expression	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Rat IgG1 (Chimeric)	HA601427	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	ANP Recombinant Rabbit mAb [JM08-36]	ET1705-75	Human	WB,IP
	AP-1/Jun D Rabbit pAb	ER62677	Human, Mouse, Rat	WB,IHC-P
	AP-A Recombinant Rabbit mAb [JE46-03]	ET7109-45	Human	WB,IHC-P
	ASK1 Recombinant Rabbit mAb [SU36-06]	ET1608-54	Human, Mouse	WB,IF-Tissue,IHC-P
	ASK1 Rabbit pAb	ER1802-84	Human, Mouse	WB,IF-Cell,IHC-P
	BMP10 Rabbit pAb	HA500312	Human	WB
	BMPR2 Mouse mAb [F8-G6]	EM1709-92	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P
	BNP Recombinant Rabbit mAb [JE30-94]	HA721946	Human, Mouse, Rat	WB,IHC-P
	CACNA1C Rabbit pAb	ER1803-49	Human, Mouse, Rat	Dot Blot,IF-Cell,IHC-P,FC
	CAMK2D Rabbit pAb	R1305-2	Human	WB,IF-Cell,IHC-P
	CAMK2G Mouse mAb [10C-G1]	EM1706-16	Human, Rat	WB,IHC-P,FC
	CaMKI Recombinant Rabbit mAb [JE36-77]	HA722544	Human, Mouse, Rat	WB,IF-Cell
	CaMKI alpha Rabbit pAb	ER62889	Human, Mouse, Rat	WB,IHC-P
	CaMKII Recombinant Rabbit mAb [SU03-57] - BSA and Azide free	HA750150	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr,FC
	CaMKII Recombinant Rabbit mAb [SU03-57]	ET1608-47	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CaMKII alpha Recombinant Rabbit mAb [JM11-07]	ET1703-54	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr,IF-Cell,FC,IP
	CaMKII alpha/beta/delta Rabbit pAb	ER62885	Human, Mouse, Rat	IHC-P,WB
	CaMKII alpha/delta Rabbit pAb	ER62886	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	CaMKII beta/gamma Rabbit pAb	ER62887	Human, Mouse, Rat	WB,IF-Cell
	CaMKII beta/gamma/delta Rabbit pAb	ER62888	Human, Mouse, Rat	WB,IHC-P
KD	CaMKII delta Recombinant Rabbit mAb [JA30-07]	ET1704-37	Human, Mouse, Rat, Zebrafish	WB,IHC-P
Relative expression	Cardiac Troponin I Recombinant Rabbit mAb [JF0550]	ET1702-37	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Cardiac Troponin T Mouse mAb [5-C7]	M1302-1	Human, Mouse	WB,IF-Cell,IHC-P,FC
	Cardiac Troponin T Recombinant Rabbit mAb [SC64-03]	ET1610-51	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Cardiac Troponin T Mouse mAb [3-G6]	EM1701-39	Human, Mouse	WB,IHC-P,IF-Cell,IF-Tissue
	Cardiac Troponin T Recombinant Rabbit mAb [SC64-03] - BSA and Azide free	HA750222	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Cardiac Troponin T Rabbit pAb	ER21102	Human, Mouse, Rat	WB,IHC-P
	Cardiotrophin-1 Rabbit pAb	ER62897	Human, Mouse, Rat	WB,IHC-P
	Caveolin-3 Recombinant Rabbit mAb [SY22-06]	ET1607-16	Human, Mouse, Rat	WB,IP,IHC-P
	CD130 / GP130 Mouse mAb [B11-B9-A8]	EM1701-43	Human, Mouse	WB,IHC-P,IF-Cell
	C/EBP Beta Recombinant Rabbit mAb [SC55-05] - BSA and Azide free	HA750235	Human, Rat, Mouse	WB,FC,IP
	C/EBP Beta Recombinant Rabbit mAb [SC55-05]	ET1610-9	Human, Rat, Mouse	WB,FC,IP
	CIITA Rabbit pAb	ER63016	Human, Mouse, Rat	WB
	c-Jun Recombinant Rabbit mAb [SY0298] - BSA and Azide free	HA750143	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	c-Jun Recombinant Rabbit mAb [SY0298]	ET1608-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	CLC-6 Rabbit pAb	ER63029	Human, Mouse, Rat, Monkey	WB
Relative expression	CLP Recombinant Rabbit mAb [PSH01-73] - BSA and Azide free	HA750748	Human, Mouse, Rat	IHC,IF-Tissue
Relative expression	CLP Recombinant Rabbit mAb [PSH01-72]	HA721724	Human, Mouse, Rat	IHC,IF-Tissue
Relative expression	CLP Recombinant Rabbit mAb [PSH01-73]	HA721725	Human, Mouse, Rat	IHC,IF-Tissue
Relative expression	CLP Recombinant Rabbit mAb [PSH01-72] - BSA and Azide free	HA750747	Human, Mouse, Rat	IHC,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Desmin Recombinant Rabbit mAb [SI18-00] - BSA and Azide free	HA750100	Human, Mouse, Rat, Zebrafish, Chicken	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Desmin Recombinant mAb	IRS140RB	Mouse	mIHC
	Desmin Recombinant Rabbit mAb [SI18-00]	ET1606-30	Human, Mouse, Rat, Zebrafish, Chicken	WB,IF-Cell,IF-Tissue,IHC-P,FC,mIHC
	Desmin Rabbit pAb	ER40113	Human, Mouse, Zebrafish	WB,IF-Cell,IHC-P
KO	Desmin Mouse mAb [3-F7]	M1501-10	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	DGK-zeta Rabbit pAb	ER63239	Human	IHC-P
	EKLF/CKLF/UKLF Rabbit pAb	ER63340	Human, Mouse, Rat	WB,IF-Cell
	EPO-R Rabbit pAb	HA500529	Human	WB
	EPOR Rabbit pAb	ER1908-28	Human, Rat	WB,IHC-P,IF-Tissue
	EYA1/EYA4 Rabbit pAb	ER65840	Human, Mouse	IHC-P
Relative expression	Fast Myosin Skeletal Heavy chain Recombinant Rabbit mAb [PSH05-12]	HA722228	Human, Mouse, Rat	IHC-P,IF-Tissue
Relative expression	Fast Myosin Skeletal Heavy chain Recombinant Rabbit mAb [PSH05-12] - BSA and Azide free	HA750985	Human, Mouse, Rat	IHC-P,IF-Tissue
	Fbx32 Recombinant Rabbit mAb [JE41-27]	ET7109-25	Human, Mouse, Rat	WB,IF-Cell
	FHL1 Recombinant Rabbit mAb [JE40-35]	HA720075	Human, Mouse, Rat	WB,IHC-P
	GATA4 Recombinant Rabbit mAb [JA11-41]	ET1704-32	Human, Mouse, Rat	WB,IF-Cell,FC
	GATA4 Rabbit pAb	R1107-5	Human, Mouse	WB,FC
	GATA4 Mouse mAb [7-4]	M1306-2	Human	WB,IHC-P
	GATA-5 Rabbit pAb	ER63508	Human, Mouse	WB,IHC-P
Relative expression	GDF15 Recombinant Rabbit mAb [PSH10-09] - BSA and Azide free	HA751339	Human	WB,IHC-P
Relative expression	GDF15 Recombinant Rabbit mAb [PSH10-09]	HA723178	Human	WB,IHC-P
	GDF15 Recombinant Rabbit mAb [JE39-90]	HA722551	Human	WB
	GPR14 Rabbit pAb	ER1909-98	Human	WB
	GRK 6 Rabbit pAb	ER63642	Human, Mouse, Rat	WB,IF-Cell
KD	GRK6 Recombinant Rabbit mAb [JE63-06]	HA721884	Human	WB,FC
KD	GSK3 alpha Rabbit pAb	R1312-1	Human	WB,IF-Cell,IHC-P,FC
	GSK3 alpha Recombinant Rabbit mAb [JE33-27]	HA722217	Human, Mouse, Rat, Monkey	WB,IHC-P
	GSK3 alpha/beta Rabbit pAb	ER63648	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	GSK3 beta Recombinant Rabbit mAb [SY28-03]	ET1607-71	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
KD	GSK3 beta Recombinant Rabbit mAb [SY28-03] - BSA and Azide free	HA750131	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
	GSK3 beta Rabbit pAb	R1510-4	Human, Mouse, Rat	IF-Cell,IHC-P,FC
	HAND2 Rabbit pAb	ER61828	Human, Rat, Mouse	WB
KO	HDAC1 Recombinant Rabbit mAb [SY12-04] - BSA and Azide free	HA750083	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,IF-Tissue
KO	HDAC1 Recombinant Rabbit mAb [SY12-04]	ET1605-35	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,IF-Tissue
KO	HDAC2 Rabbit pAb	R1601-7	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KO	HDAC2 Recombinant Rabbit mAb [SY29-02]	ET1607-78	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO	HDAC2 Recombinant Rabbit mAb [SD0816]	ET1612-50	Human, Mouse, Rat	WB,IP,IHC-P,FC
KO	HDAC2 Recombinant Rabbit mAb [SY29-02] - BSA and Azide free	HA750134	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	HDAC3 Recombinant Rabbit mAb [SC55-02] - BSA and Azide free	HA750221	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
KD	HDAC3 Mouse mAb [B2-G9]	M1511-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	HDAC3 Recombinant Rabbit mAb [SC55-02]	ET1610-5	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
KD	HDAC4 Recombinant Rabbit mAb [SD08-35]	ET1612-51	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	HDAC5 Recombinant Rabbit mAb [PSH08-70]	HA723022	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,ChIP
	HDAC5 Rabbit pAb	ER64245	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	HDAC5 Recombinant Rabbit mAb [PSH08-70] - BSA and Azide free	HA751235	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,ChIP
	HDAC5/9 Rabbit pAb	ER63671	Human, Mouse	WB,IHC-P
	HDAC7 Recombinant Rabbit mAb [SD082-4]	ET1612-67	Human, Mouse, Rat	WB,IF-Cell,FC
	HDAC8 Recombinant Rabbit mAb [JJ0845]	ET1701-12	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	HDAC8 Recombinant Rabbit mAb [SD082-07]	ET1612-90	Human, Mouse, Rat	WB,FC,IP
KD	HDAC9 Recombinant Rabbit mAb [JU30-44]	ET1706-36	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP
Relative expression	heavy chain Myosin / MYH3 Recombinant Rabbit mAb [PSH18-12]	HA723975	Human, Mouse, Rat	WB,IHC-P,IP
	Hexim1 Rabbit pAb	ER63677	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Histone deacetylase 9 Rabbit pAb	ER63682	Human	WB,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human GDF-15 Recombinant Rabbit mAb [PSH16-02] - BSA and Azide free (Detector)	HA723807	Human	ELISA(Det)
	Human GDF-15 Recombinant Rabbit mAb [PSH16-00] - BSA and Azide free (Capture)	HA723805	Human	ELISA(Cap)
	Biotin Conjugated Human GDF-15 Recombinant Rabbit mAb [PSH16-01] - Detector	HA723806B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human GDF-15 Recombinant Rabbit mAb [PSH16-02] - Detector	HA723807B	Human	ELISA(Det),ELISA
	Human GDF-15 Recombinant Rabbit mAb [PSH16-01] - BSA and Azide free (Detector)	HA723806	Human	ELISA(Det)
	Biotin Conjugated Human IL1RL1 / ST2 Recombinant Rabbit mAb [PSH19-93] - Detector	HA725376B	Human	ELISA(Det),ELISA
	Human IL1RL1 / ST2 Recombinant Rabbit mAb [PSH19-92] - BSA and Azide free (Capture)	HA725375	Human	ELISA(Cap)
	Human IL1RL1 / ST2 Recombinant Rabbit mAb [PSH19-93] - BSA and Azide free (Detector)	HA725376	Human	ELISA(Det)
	Human IL-33 Recombinant Rabbit mAb [PSH03-05] - BSA and Azide free (Capture)	HA721935	Human	ELISA(Cap)
	Biotin Conjugated Human IL-33 Recombinant Rabbit mAb [PSH03-06]	HA722628B	Human	ELISA(Det)
	Human IL-33 Recombinant Rabbit mAb [PSH03-06] - BSA and Azide free (Detector)	HA721936	Human	ELISA(Det)
	IL33 Rabbit pAb	ER1911-36	Human, Rat	WB,IHC-P,IF-Tissue
	Junctophilin-2 Rabbit pAb	0407-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Junctophilin-2 Recombinant Mouse mAb [3-2-R] - BSA and Azide free	HA610170	Human, Mouse, Rat	WB,IHC-P
	Junctophilin-2 Recombinant Mouse mAb [3-2-R]	HA601297	Human, Mouse, Rat	WB,IHC-P
	Junctophilin-2 Mouse mAb [3-2]	M0407-20	Human, Mouse, Rat	WB,IHC-P
	JunD Recombinant Rabbit mAb [PSH08-85]	HA723037	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP,ChIP
	JunD Recombinant Rabbit mAb [JE59-51]	HA722913	Human, Mouse, Rat	WB,IF-Cell,IHC-P,ChIP,IP
	JunD Recombinant Rabbit mAb [SD0830]	ET1612-92	Human	WB,IF-Cell,IF-Tissue,IP,IHC-P,FC
	JunD Recombinant Rabbit mAb [PSH08-85] - BSA and Azide free	HA751250	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP,ChIP
	Kinesin 5B Rabbit pAb	ER1911-89	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	KLF10 Recombinant Rabbit mAb [JE34-09]	HA722324	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	KLF4 Rabbit pAb	R1308-1	Human, Mouse, Rat	WB,IHC-P,FC
	KLF4 Recombinant Rabbit mAb [JF98-08] - BSA and Azide free	HA750356	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	KLF4 Recombinant Rabbit mAb [JF98-08]	ET1702-71	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	KLF5 Recombinant Rabbit mAb [PSH15-42] - BSA and Azide free	HA751557	Human, Mouse, Rat	IF-Cell,FC
	KLF5 Recombinant Rabbit mAb [PSH15-42]	HA723748	Human, Mouse, Rat	IF-Cell,FC
	KLF7 Rabbit pAb	ER1911-99	Mouse	WB
	Krs-1/2 Rabbit pAb	ER63820	Human, Mouse	WB,IHC-P
	LDB3 Rabbit pAb	ER62089	Human, Mouse	WB
	MEF2A Recombinant Rabbit mAb [JA33-04]	ET1704-58	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	MEF2A + MEF2C Recombinant Rabbit mAb [JE63-26]	HA721030	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	MEF-2B Rabbit pAb	ER63939	Human, Mouse	WB,IHC-P
Relative expression	MEF2C Recombinant Rabbit mAb [JE35-45]	HA721298	Human, Rat	WB,IF-Cell,IHC-P,IP,FC
	MEF-2D Rabbit pAb	ER63940	Human, Mouse, Rat	IHC-P
	Mouse GP130 / IL6ST Recombinant Rabbit mAb [PSH19-82] - BSA and Azide free (Capture)	HA724106	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse GP130 / IL6ST Recombinant Rabbit mAb [PSH19-83] - Detector	HA724107B	Mouse	ELISA(Det),ELISA
	Mouse GP130 / IL6ST Recombinant Rabbit mAb [PSH19-83] - BSA and Azide free (Detector)	HA724107	Mouse	ELISA(Det)
	MST1 Rabbit pAb	ER1913-22	Human, Rat	IHC-P,IF-Cell,IF-Tissue,FC
KD	MST2 Recombinant Rabbit mAb [SC05-83]	ET1610-8	Human, Mouse, Rat	WB,IHC-P,IP
	MTPN Rabbit pAb	ER62458	Human, Mouse, Rat	WB
	MuRF1 Mouse mAb [B8-E1]	EM1706-67	Human	IF-Cell,IHC-P,FC
Relative expression	MURF1 Recombinant Rabbit mAb [PSH18-10]	HA723973	Human, Mouse, Rat	WB,FC,IP
	MYBPC3 Rabbit pAb	ER65870	Human, Mouse, Rat	WB
Relative expression	MYH2 Recombinant Rabbit mAb [JE53-30]	HA721862	Human, Mouse, Rat	WB,IHC-P
	MYH3 Rabbit pAb	ER60864	Human, Rat, Mouse	WB
	MYH4 Rabbit pAb	ER64077	Human	IHC-P
	MYH-pan Rabbit pAb	ER65575	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	MYL3 Rabbit pAb	ER64078	Human	WB
	Myosin heavy chain Recombinant Rabbit mAb [JF097-7]	ET1702-88	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Myosin Light Chain 2 Recombinant Rabbit mAb [SN67-09]	ET1611-13	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IP
	Myosin Light Chain 2 Rabbit pAb	R1512-19	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Myosin Light Chain 2 Recombinant Rabbit mAb [SN67-09] - BSA and Azide free	HA750244	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IP
Cell treatment	Neuropeptide Y Recombinant Rabbit mAb [PSH16-52]	HA723856	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	Biotin Conjugated Neuropeptide Y Recombinant Rabbit mAb [PSH16-89] - Detector	HA723885B	Human, Mouse, Rat	ELISA(Det),ELISA
	Neuropeptide Y Recombinant Rabbit mAb [PSH16-90] - BSA and Azide free (Detector)	HA723886	Human, Mouse, Rat	ELISA(Det)
	Neuropeptide Y Recombinant Rabbit mAb [PSH16-88] - BSA and Azide free (Capture)	HA723884	Human, Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Neuropeptide Y Recombinant Rabbit mAb [PSH16-90] - Detector	HA723886B	Human, Mouse, Rat	ELISA(Det),ELISA
Cell treatment	Neuropeptide Y Recombinant Rabbit mAb [PSH16-52] - BSA and Azide free	HA751601	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	Neuropeptide Y Recombinant Rabbit mAb [PSH16-89] - BSA and Azide free (Detector)	HA723885	Human, Mouse, Rat	ELISA(Det)
Relative expression	NFAT1 Recombinant Rabbit mAb [JA11-08]	ET1704-14	Human	WB,IF-Tissue,IHC-P
Relative expression	NFAT1 Recombinant Rabbit mAb [JA11-08] - BSA and Azide free	HA750405	Human	WB,IF-Tissue,IHC-P
	NPY Rabbit pAb	ER64220	Human, Mouse, Rat	IHC-P,WB,IF
	Pan-Actin Mouse mAb [A2E1]	HA600032	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Pan-Actin Recombinant Mouse mAb [A2E1-R] - BSA and Azide free	HA610160	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	Pan-Actin Recombinant Mouse mAb [A2E1-R]	HA601287	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
Cell treatment	Phospho-CaMKII (T286) Recombinant Rabbit mAb [JE51-81]	HA721794	Human, Mouse, Rat	WB
	Phospho-CaMKII (T286) Mouse mAb [2G1]	RT1092	Human, Mouse, Rat	WB,IP,IF
Cell treatment	Phospho-c-Jun (S63) Recombinant Rabbit mAb [SY0297]	ET1608-4	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-c-Jun (S63) Recombinant Rabbit mAb [SY0297] - BSA and Azide free	HA750146	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit mAb [JE65-46]	HA722475	Human, Mouse, Rat	WB,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit mAb [PSH07-12] - BSA and Azide free	HA751133	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit mAb [PSH07-12]	HA722800	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-c-Jun (T91) Recombinant Rabbit mAb [JJ080-9]	ET1701-32	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-GSK3 (alpha + beta) (Y216 + Y279) Recombinant Rabbit mAb [SY26-05] - BSA and Azide free	HA750127	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,IF-Tissue,FC
Cell treatment	Phospho-GSK3 (alpha + beta) (Y216 + Y279) Recombinant Rabbit mAb [SY26-05]	ET1607-54	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,IF-Tissue,FC
Cell treatment	Phospho-GSK3 alpha (S21) Recombinant Rabbit mAb [PSH09-11] - BSA and Azide free	HA751265	Human, Mouse	WB,IHC-P,Dot Blot
Cell treatment	Phospho-GSK3 alpha (S21) Recombinant Rabbit mAb [PSH09-11]	HA723069	Human, Mouse	WB,IHC-P,Dot Blot
Cell treatment	Phospho-GSK3 beta (S9) Recombinant Rabbit mAb [SY02-71] - BSA and Azide free	HA750128	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-GSK3 beta (S9) Recombinant Rabbit mAb [SY02-71]	ET1607-60	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-JunD (S255) Recombinant Rabbit mAb [JJ08-21]	ET1701-35	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Phospholamban Mouse mAb [1G3]	RT1491	Human, Mouse, Rat	WB,IP,IF
	Phospholamban Rabbit pAb	HA500103	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-MST1/MST2 (T180 + T183) Recombinant Rabbit mAb [JE46-64]	HA721737	Human, Mouse, Rat	WB
Cell treatment	Phospho-Phospholamban (S16) Recombinant Rabbit mAb [JE47-45]	HA722291	Mouse, Rat	WB
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH15-91] - BSA and Azide free	HA751582	Human, Mouse, Rat	WB,IF-Cell,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH15-91]	HA723796	Human, Mouse, Rat	WB,IF-Cell,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH16-27]	HA723825	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH16-27] - BSA and Azide free	HA751595	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,Dot Blot
Cell treatment	Phospho-PKC alpha (T497) Recombinant Rabbit mAb [JE44-81]	HA722278	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-PKC alpha (T638) Recombinant Rabbit mAb [JF0964]	ET1702-17	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,Dot Blot
Cell treatment	Phospho-PKC delta (Y311) Recombinant Rabbit mAb [JE44-20]	HA722458	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-PKC (S729) Recombinant Rabbit mAb [JE77-54]	HA722846	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-PKC zeta (T560) Recombinant Rabbit mAb [SN0753]	ET1611-78	Human, Mouse, Rat	WB,IHC-P
Cell treatment, KD	Phospho-STAT3 (S727) Recombinant Rabbit mAb [SY24-09] - BSA and Azide free	HA750121	Human, Rat, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment, KD	Phospho-STAT3 (S727) Recombinant Rabbit mAb [SY24-09]	ET1607-39	Human, Rat, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-STAT3 (Y705) Recombinant Rabbit mAb [SZ43-01]	ET1603-40	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-STAT3 (Y705) Recombinant Rabbit mAb [SZ43-01] - BSA and Azide free	HA750067	Human, Mouse, Rat	WB,IHC-P,IP
	PKC Rabbit pAb	ER64620	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
KD	PKC alpha Recombinant Rabbit mAb [SU31-08]	ET1608-15	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	PKC beta Rabbit pAb	ER1802-14	Human, Mouse	WB,IF-Cell,IHC-P
	PKC beta 1 Recombinant Rabbit mAb [JM43-26]	ET1705-30	Human, Mouse, Rat	WB,IHC-P,FC
	PKC beta 2 Recombinant Rabbit mAb [ST48-06]	ET1609-44	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	PKC delta Recombinant Rabbit mAb [JJ093-02] - BSA and Azide free	HA750325	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	PKC delta Recombinant Rabbit mAb [JJ093-02]	ET1701-85	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	PKC epsilon Recombinant Rabbit mAb [JA38-21]	ET1704-85	Human, Mouse, Rat	WB,IHC-P
	PKC eta Recombinant Rabbit mAb [JE32-41]	HA722389	Human, Mouse, Rat, Monkey	WB
	PKC gamma Recombinant Rabbit mAb [JG69-32]	ET7108-42	Human, Mouse, Rat	WB,IHC-P
	PKC theta Rabbit pAb	ER64621	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Relative expression	PKC zeta Recombinant Rabbit mAb [PSH06-74] - BSA and Azide free	HA751122	Human	WB
Relative expression	PKC zeta Recombinant Rabbit mAb [PSH06-74]	HA722773	Human	WB
	PNKD Rabbit pAb	ER61038	Human, Mouse	WB
	Rad GTPase Rabbit pAb	ER64742	Human	WB,IHC-P,IF-Cell
	RGS4 Rabbit pAb	ER61127	Human, Mouse, Rat	WB
	RGS5 Rabbit pAb	ER64808	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Sarcomeric Alpha Actinin Recombinant Rabbit mAb [JE45-22]	ET7109-42	Human, Mouse, Rat	WB,IP,IHC-P
	SEMA3A Rabbit pAb	ER64921	Human, Mouse, Rat	IHC-P
	Sigma1-receptor Rabbit pAb	HA500511	Human, Mouse	WB,IHC-P
Relative expression	Slow Skeletal Myosin Heavy chain Recombinant Rabbit mAb [PSH07-54] - BSA and Azide free	HA751092	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression	Slow Skeletal Myosin Heavy chain Recombinant Rabbit mAb [PSH06-58]	HA722721	Human, Mouse, Rat	WB,IHC-P
Relative expression	Slow Skeletal Myosin Heavy chain Recombinant Rabbit mAb [PSH07-54]	HA722735	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression	Slow Skeletal Myosin Heavy chain Recombinant Rabbit mAb [PSH06-58] - BSA and Azide free	HA751086	Human, Mouse, Rat	WB,IHC-P
KO	Smad4 Recombinant Rabbit mAb [SP05-05] - BSA and Azide free	HA750070	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KO	Smad4 Recombinant Rabbit mAb [SP05-05]	ET1604-12	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Smad4 Recombinant Rabbit mAb [PD01-37]	HA721234	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Smad4 Rabbit pAb	ER1706-90	Human, Mouse, Rat	IF-Cell,IHC-P,FC
	SMAD6 Mouse mAb [H6-H3]	EM1709-56	Human	WB
Cell treatment	ST2 Mouse mAb [PSH01-69]	HA601171	Human	WB,IF-Cell,FC
Cell treatment	ST2 Mouse mAb [PSH01-68]	HA601170	Human	WB,IF-Cell,FC
	ST2 Rabbit pAb	0902-2	Human, Mouse, Zebrafish	WB
	STAT3 Rabbit pAb	RT1591	Human, Mouse, Rat	WB,IP,IF,IHC-P
KD	STAT3 Recombinant Rabbit mAb [SY34-01]	ET1605-45	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	STAT3 Recombinant Rabbit mAb [SY24-08] - BSA and Azide free	HA750120	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
KD	STAT3 Recombinant Rabbit mAb [SY34-01] - BSA and Azide free	HA750087	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	STAT3 Recombinant Rabbit mAb [SY24-08]	ET1607-38	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	STK4 / MST-1 Rabbit pAb	HA500031	Human, Mouse	WB,IHC-P
	STK4 / MST-1 Rabbit pAb	HA500043	Human, Mouse, Rat	WB,IHC-P
	TELT Rabbit pAb	ER61299	Human	WB
	TNNC1 Rabbit pAb	ER1902-25	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	TNNC1 Rabbit pAb	HA500407	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	TP4A3 Rabbit pAb	ER62006	Human, Mouse	WB
	Tropomyosin 1 Recombinant Rabbit mAb [JG38-66]	ET7108-74	Human, Mouse, Rat	WB,IHC-P,FC
	Troponin I-C Rabbit pAb	ER65220	Mouse, Rat	WB,IHC-P,IF-Cell
	TRPC1 Recombinant Rabbit mAb [JE35-25]	HA721427	Human, Mouse, Rat	WB
	TTN Rabbit pAb	ER65816	Human, Mouse, Rat	IHC-P
	WISP2 Rabbit pAb	ER1918-66	Human, Mouse	WB,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层