

心血管研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	ACE2 Rabbit pAb	ER1902-40	Human	WB,IHC-P,FC
KD	ACE2 Recombinant Rabbit mAb [SN0754]	ET1611-58	Human, Mouse, Rat, Hamster	WB,IF-Cell,IHC-P,IF-Tissue,IP
KD	ACE2 Recombinant Rabbit mAb [SN0754] - BSA and Azide free	HA750264	Human, Mouse, Rat, Hamster	WB,IF-Cell,IHC-P,IF-Tissue,IP
KD	ADFP / PLIN2 Recombinant Rabbit mAb [JA31-81]	ET1704-17	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC,IHC-P
KD	ADFP / PLIN2 Recombinant Rabbit mAb [JA31-81] - BSA and Azide free	HA750407	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC,IHC-P
	Adiponectin receptor protein 1 Recombinant Rabbit mAb [SC69-04]	ET1610-86	Human, Mouse, Rat, Zebrafish	IF-Cell,IF-Tissue,IHC-P,FC
	Adiponectin receptor protein 1 Recombinant Rabbit mAb [SC69-04] - BSA and Azide free	HA750233	Human, Mouse, Rat, Zebrafish	IF-Cell,IF-Tissue,IHC-P,FC
	AKT1 Mouse mAb [D9-9-C9]	EM40507	Human, Mouse, Rat	WB,IHC-P
	AKT1 Recombinant Rabbit mAb [PSH0-50]	HA721364	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
KD	AKT1 Recombinant Rabbit mAb [ST05-09]	ET1609-47	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
KD	AKT1 Recombinant Rabbit mAb [ST05-09] - BSA and Azide free	HA750177	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
	AKT1 Recombinant Rabbit mAb [PSH0-50] - BSA and Azide free	HA750627	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75-09] - BSA and Azide free	HA750840	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	AKT1/2/3 Rabbit pAb	ER1901-33	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	AKT1/2/3 Recombinant Rabbit mAb [ST48-09]	ET1609-51	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75-09]	HA721870	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	AKT1/2/3 Recombinant Rabbit mAb [ST48-09] - BSA and Azide free	HA750179	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
	Akt1/3 Rabbit pAb	ER62638	Human, Mouse, Rat	WB,IHC-P
	Akt2 Mouse mAb [2A1]	RT1027	Human, Mouse, Rat	WB,IP,IF,IHC-P
	AKT2 Rabbit pAb	HA500091	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Albumin Recombinant Rabbit mAb [JF32-10] - BSA and Azide free	HA750352	Human, Mouse, Rat, Cow	WB,IHC-P,IF-Tissue
	Albumin Recombinant Rabbit mAb	IRS087RB	Mouse	mIHC
	Albumin Rabbit pAb	0806-9	Human, Mouse	WB,IF-Cell,IHC-P,FC
	Albumin Mouse mAb [C5-A8]	EM1901-80	Human, Rat	WB,IHC-P,FC,ELISA
	Albumin Recombinant Rabbit mAb [JF32-10]	ET1702-55	Human, Mouse, Rat, Cow	WB,IHC-P,mIHC,IF-Tissue
	Albumin Mouse mAb [1-G10]	EM1901-81	Human	WB,FC,ELISA
KO	AMPK alpha 1 Recombinant Rabbit mAb [SU03-48]	ET1608-40	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC,IP,IHC-Fr,IHC-P
KO	AMPK alpha 1 Recombinant Rabbit mAb [SU03-48] - BSA and Azide free	HA750147	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC,IP,IHC-Fr,IHC-P
	AMPK alpha 1 Mouse mAb [B1-E10]	M1409-1	Human	WB,IF-Cell
	AMPK alpha 1/2 Rabbit pAb	ER62660	Human, Mouse, Rat, Monkey	WB,IHC-P
	AMPK alpha 1/2 Recombinant Rabbit mAb [PSH08-86]	HA723038	Human, Mouse, Rat	WB,IHC-P
	AMPK alpha 1/2 Recombinant Rabbit mAb [PSH08-86] - BSA and Azide free	HA751251	Human, Mouse, Rat	WB,IHC-P
	AMPK alpha 2 Mouse mAb [A6A11]	HA600079	Human, Mouse, Rat	IF-Cell,IHC-P
	AMPK alpha 2 Mouse mAb [A6A12]	HA600080	Human, Mouse, Rat	WB
	AMPK alpha 2 Mouse mAb [A6A10]	HA600078	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	ANP Recombinant Rabbit mAb [JM08-36]	ET1705-75	Human	WB,IP
	Apolipoprotein A1 Recombinant Rabbit mAb [JF0548]	ET1702-23	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Apolipoprotein A1 Rabbit pAb	RT1046	Human, Mouse, Rat	WB,IP,IF
	Apolipoprotein A1 Mouse mAb	M0809-13	Human	WB,IF-Cell,IHC-P
	Apolipoprotein A1 Rabbit pAb	ER1802-29	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	BDNF Rabbit pAb	ER130915	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	BDNF Recombinant Rabbit mAb [SJ12-09]	ET1606-42	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC
	BDNF Recombinant Rabbit mAb [JE59-35]	HA722912	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-P,IHC-Fr,FC
KD	BDNF Recombinant Rabbit mAb [SJ12-09] - BSA and Azide free	HA750104	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC
	BNP Recombinant Rabbit mAb [JE30-94]	HA721946	Human, Mouse, Rat	WB,IHC-P
	BSA Mouse mAb [2-B8]	M0809-3	Bovine	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	BSA Mouse mAb [A7D10]	HA600085	Cow	ELISA,WB
	BSA Mouse mAb [A7D11]	HA600086	Cow	ELISA,WB
	BSA Mouse mAb [A7D12]	HA600087	Cow	ELISA
	Cardiac Troponin T Recombinant Rabbit mAb [SC64-03] - BSA and Azide free	HA750222	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Cardiac Troponin T Mouse mAb [5-C7]	M1302-1	Human, Mouse	WB,IF-Cell,IHC-P,FC
	Cardiac Troponin T Rabbit pAb	ER21102	Human, Mouse, Rat	WB,IHC-P
	Cardiac Troponin T Mouse mAb [3-G6]	EM1701-39	Human, Mouse	WB,IHC-P,IF-Cell,IF-Tissue
	Cardiac Troponin T Recombinant Rabbit mAb [SC64-03]	ET1610-51	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	CCR2 Recombinant Rabbit mAb [SN707]	ET1611-65	Human, Mouse	WB,IHC-P,IP,FC
	CCR2 Recombinant Rabbit mAb [SN707] - BSA and Azide free	HA750266	Human, Mouse	WB,IHC-P,IP,FC
	CD34 Mouse mAb [15H1]	EM1901-01	Human, Mouse, Rat	IHC-P,FC
	CD34 Rabbit pAb	ER1706-66	Human, Mouse	WB,IHC-P
	CD34 Recombinant Rabbit mAb [SI16-01]	ET1606-11	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,mIHC,FC
	CD34 Rabbit pAb	0108-2	Human, Mouse, Zebrafish	WB,IF-Cell,IHC-P
	CD34 Recombinant Rabbit mAb [SI16-01] - BSA and Azide free	HA750092	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	CD34 Recombinant Mouse mAb [PDM0-12]	HA601109	Human	IHC-P,IF-Tissue
	CD34 Recombinant Rabbit mAb	IRS024	Human	mIHC
	CD34 Recombinant Rabbit mAb	IRS050	Mouse	mIHC
	CD34 Rabbit pAb	ER0802	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
Cell treatment	CD36 Recombinant Rabbit mAb [JJ2005] - BSA and Azide free	HA750307	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Relative expression	CD36 Recombinant Rabbit mAb [PSH18-47]	HA724000	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	CD36 Rabbit pAb	ER1802-7	Human, Mouse, Rat	IF-Cell,IHC-P,FC
	CD36 Recombinant Rabbit mAb	IRS118RB	Human	mIHC
Cell treatment	CD36 Recombinant Rabbit mAb [JJ2005]	ET1701-24	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Relative expression	CD73 Mouse mAb [A6F7]	HA601010	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
Relative expression	CD73 Recombinant Rabbit mAb [PSH09-57] - BSA and Azide free	HA751304	Human	WB,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD73 Recombinant Mouse mAb [A6G2-R] - BSA and Azide free	HA610172	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Relative expression	CD73 Mouse mAb [A6G2]	HA601004	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	CD73 Rabbit pAb	R1107-6	Human, Mouse	WB,IF-Cell,FC
Relative expression	CD73 Mouse mAb [A6F9]	HA601011	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC,IF-Cell
Relative expression	CD73 Recombinant Rabbit mAb [PSH09-57]	HA723119	Human	WB,IF-Cell,FC
	CD73 Recombinant Rabbit mAb [JM11-40]	ET1703-41	Human	WB,IHC-P,IF-Tissue
Relative expression	CD73 Mouse mAb [A6G1]	HA601012	Human	WB,IHC-P,IF-Tissue,FC,IF-Cell
Relative expression	CD73 Recombinant Mouse mAb [A6G2-R]	HA601299	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	CD73 Mouse mAb [6E-G3]	EM1706-12	Human	WB,IHC-P,FC
Relative expression, KD	CD9 Recombinant Rabbit mAb [PSH0-95] - BSA and Azide free	HA750675	Human	WB,IHC-P
	CD9 Recombinant Rabbit mAb [SA35-08]	ET1601-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	CD9 Recombinant Rabbit mAb [SA35-08] - BSA and Azide free	HA750027	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
Relative expression, KD	CD9 Recombinant Rabbit mAb [PSH0-95]	HA721533	Human	WB,IHC-P
	CD9 Rabbit pAb	ER80402	Human, Mouse	WB
	COX2 Mouse mAb	RT1159	Human, Mouse, Rat	WB,IP,IF,IHC-P,ELISA
	COX2 Mouse mAb [A3F8]	EM1902-11	Human, Rat	WB,IHC-P,FC
KD	COX2 Mouse mAb [A3F7]	EM1902-12	Human, Rat	WB,IHC-P,FC
	COX2 Recombinant Rabbit mAb [6J8]	RE6064	Human, Mouse	WB,ICC
KD	COX2/Cyclooxygenase 2 Recombinant Rabbit mAb [SC56-06] - BSA and Azide free	HA750207	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
KD	COX2/Cyclooxygenase 2 Recombinant Rabbit mAb [SC56-06]	ET1610-23	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	FABP3 Recombinant Rabbit mAb [JU73-03]	ET1706-50	Mouse, Rat	WB,IHC-P
	FABP3 Recombinant Rabbit mAb [JE53-37]	HA721859	Human, Mouse, Rat	WB,IHC-P
Relative expression	FABP4 Recombinant Rabbit mAb [JM10-99] - BSA and Azide free	HA750398	Human, Mouse, Rat	WB,IF-Cell,IHC-Fr
Relative expression	FABP4 Recombinant Rabbit mAb [JM10-99]	ET1703-98	Human, Mouse, Rat	WB,IF-Cell,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KO	Fatty Acid Synthase Recombinant Rabbit mAb [JJ0939]	ET1701-91	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	Fatty Acid Synthase Rabbit pAb	R1706-8	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KO	Fatty Acid Synthase Recombinant Rabbit mAb [JJ0939] - BSA and Azide free	HA750327	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	FGFR3 Recombinant Rabbit mAb [JM110-33] - BSA and Azide free	HA750388	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
KD	FGFR3 Recombinant Rabbit mAb [JM81-10] - BSA and Azide free	HA750379	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	FGFR3 Recombinant Rabbit mAb [JM110-33]	ET1703-65	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
KD	FGFR3 Recombinant Rabbit mAb [JM81-10]	ET1703-36	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	Fibronectin Recombinant Rabbit mAb [PSH14-21] - BSA and Azide free	HA751514	Human, Mouse	WB,IHC-P,IF-Cell,FC
	Fibronectin Recombinant Rabbit mAb [PSH14-21]	HA723610	Human, Mouse	WB,IHC-P,IF-Cell,FC
	Fibronectin Recombinant Rabbit mAb [PSH14-22] - BSA and Azide free	HA751515	Human, Mouse	WB,IF-Cell,FC
	Fibronectin Recombinant mAb	IRS144RB	Mouse	mIHC
	Fibronectin Recombinant Rabbit mAb [PSH14-22]	HA723611	Human, Mouse	WB,IF-Cell,FC
	Fibronectin Mouse mAb [3G4]	RT1224	Human, Mouse, Rat	WB,IP,IF,IHC-P
KD	Fibronectin Recombinant Rabbit mAb [JF0582] - BSA and Azide free	HA750340	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
KD	Fibronectin Recombinant Rabbit mAb [JF0582]	ET1702-25	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
	Galectin 3 Rabbit pAb	ER1803-82	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression, KD	Galectin 3 Recombinant Rabbit mAb [JF39-10]	ET1702-48	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC,IHC-P
Relative expression	GDF15 Recombinant Rabbit mAb [PSH10-09]	HA723178	Human	WB,IHC-P
	GDF15 Recombinant Rabbit mAb [JE39-90]	HA722551	Human	WB
Relative expression	GDF15 Recombinant Rabbit mAb [PSH10-09] - BSA and Azide free	HA751339	Human	WB,IHC-P
KD	Glucose Transporter GLUT1 Recombinant Rabbit mAb [SA0377] - BSA and Azide free	HA750002	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	iFluor™ 647 Conjugated Glucose Transporter GLUT1 Recombinant Rabbit mAb [SA0377]	HA720170F	Human, Mouse	IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	Glucose Transporter GLUT1 Recombinant Rabbit mAb [SA0377]	ET1601-10	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Glucose Transporter GLUT1 Recombinant Rabbit mAb	IRS151RB	Human	mIHC
	iFluor™ 488 Conjugated Glucose Transporter GLUT1 Recombinant Rabbit mAb [SA0377]	HA720154F	Human, Mouse, Rat	IF-Cell,IF-Tissue
	Glucose Transporter GLUT1 Recombinant Antibody [SA0377] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610348	Human, Mouse, Rat	IHC-Fr,IHC-P,IF-Tissue,WB
	Glucose Transporter GLUT1 Recombinant Antibody [SA0377] - Rat IgG1 (Chimeric)	HA601521	Human, Mouse, Rat	IHC-Fr,IHC-P,IF-Tissue,WB
	Glucose Transporter GLUT1 Rabbit pAb	ER1510-11	Human, Mouse	WB,IF-Cell,IHC-P,FC
	Heme Oxygenase 1 (HO-1) Recombinant Rabbit mAb	IRS232RB	Mouse	mIHC
	Heme Oxygenase 1 (HO-1) Recombinant Rabbit mAb [SP08-07] - BSA and Azide free	HA750079	Human	WB,IHC-P,IP,FC,IF-Tissue
KD	Heme Oxygenase 1 (HO-1) Rabbit pAb	ER1802-73	Human, Mouse	WB,IF-Cell,IHC-P,FC
	iFluor™ 594 Conjugated Heme Oxygenase 1 (HO-1) Recombinant Rabbit mAb [SP08-07]	HA720169F	Human	IF-Tissue
KD	Heme Oxygenase 1 (HO-1) Recombinant Rabbit mAb [PSH02-76] - BSA and Azide free	HA750831	Human, Mouse	WB,IF-Cell,IHC-P
	Heme Oxygenase 1 (HO-1) Recombinant Rabbit mAb	IRS156RB	Human	mIHC
	Heme Oxygenase 1 (HO-1) Recombinant Rabbit mAb [SP08-07]	ET1604-45	Human	WB,IHC-P,IP,FC,IF-Tissue
KD	Heme Oxygenase 1 (HO-1) Recombinant Rabbit mAb [PSH02-76]	HA721854	Human, Mouse	WB,IF-Cell,IHC-P
	H-FABP Recombinant Rabbit mAb [PSH16-99] - BSA and Azide free (Detector)	HA725178	Human, Mouse	ELISA(Det)
	H-FABP Recombinant Rabbit mAb [PSH16-98] - BSA and Azide free (Capture)	HA725177	Human, Mouse	ELISA(Cap)
	Biotin Conjugated H-FABP Recombinant Rabbit mAb [PSH16-99] - Detector	HA725178B	Human, Mouse	ELISA(Det),ELISA
	HIF-1 alpha Mouse mAb [5H7]	RT1278	Human, Mouse, Rat	WB,IP,ICC,IHC-P
Cell treatment	HIF-1 alpha Rabbit pAb	ER1802-41	Human	WB,IHC-P
Cell treatment	HIF-1 alpha Recombinant Rabbit mAb [JE75-33]	HA721997	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	HIF-1 alpha Recombinant Rabbit mAb [PSH07-00] - BSA and Azide free	HA751127	Human, Mouse	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	HIF-1 alpha Recombinant Rabbit mAb [PSH07-00]	HA722778	Human, Mouse	WB,IF-Cell,IHC-P
KD	HIF-1 beta Recombinant Rabbit mAb [ST05-25] - BSA and Azide free	HA750187	Human, Mouse, Rat	WB,IHC-P
KD	HIF-1 beta Recombinant Rabbit mAb [ST05-25]	ET1609-67	Human, Mouse, Rat	WB,IHC-P
	HIF-2 alpha Recombinant Rabbit mAb [JB24-42]	ET7107-32	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	HMGCR Recombinant Rabbit mAb [JF0981]	ET1702-41	Human, Mouse, Rat	WB,IP,IHC-P
	HMGCR Recombinant Rabbit mAb [JF0981] - BSA and Azide free	HA750345	Human, Mouse, Rat	WB,IP,IHC-P
	Human ACE2 Recombinant Rabbit mAb [PSH17-17] - BSA and Azide free (Capture)	HA723904	Human	ELISA(Cap)
	Human ACE2 Recombinant Rabbit mAb [PSH17-18] - BSA and Azide free (Detector)	HA723905	Human	ELISA(Det)
	Biotin Conjugated Human ACE2 Recombinant Rabbit mAb [PSH17-18] - Detector	HA723905B	Human	ELISA(Det),ELISA
	Human Apolipoprotein A1 Recombinant Rabbit mAb [PSH18-89] - BSA and Azide free (Capture)	HA724026	Human	ELISA(Cap)
	Human Apolipoprotein A1 Recombinant Rabbit mAb [PSH18-91] - BSA and Azide free (Detector)	HA724028	Human	ELISA(Det)
	Biotin Conjugated Human Apolipoprotein A1 Recombinant Rabbit mAb [PSH18-91] - Detector	HA724028B	Human	ELISA(Det),ELISA
	Human Apolipoprotein A1 Recombinant Rabbit mAb [PSH18-90] - BSA and Azide free (Detector)	HA724027	Human	ELISA(Det)
	Biotin Conjugated Human Apolipoprotein A1 Recombinant Rabbit mAb [PSH18-90] - Detector	HA724027B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human CD73 Recombinant Rabbit mAb [PSH08-16] - Detector	HA722959B	Human	ELISA(Det)
	Human CD73 Recombinant Rabbit mAb [PSH08-16] - BSA and Azide free (Detector)	HA722958	Human	ELISA(Det)
	Human CD73 Recombinant Rabbit mAb [PSH08-15] - BSA and Azide free (Capture)	HA722957	Human	ELISA(Cap)
	Human GDF-15 Recombinant Rabbit mAb [PSH16-00] - BSA and Azide free (Capture)	HA723805	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human GDF-15 Recombinant Rabbit mAb [PSH16-02] - BSA and Azide free (Detector)	HA723807	Human	ELISA(Det)
	Biotin Conjugated Human GDF-15 Recombinant Rabbit mAb [PSH16-01] - Detector	HA723806B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human GDF-15 Recombinant Rabbit mAb [PSH16-02] - Detector	HA723807B	Human	ELISA(Det),ELISA
	Human GDF-15 Recombinant Rabbit mAb [PSH16-01] - BSA and Azide free (Detector)	HA723806	Human	ELISA(Det)
	Human HGF enzyme-linked immunoassay kit	EH0019	Human	ELISA
	Human IFN-gamma enzyme-linked immunoassay kit	EH0039	Human	ELISA
	Human IL-18 enzyme-linked immunoassay kit	EH0053	Human	ELISA
	Human IL-1 β enzyme-linked immunoassay kit	EH0006	Human	ELISA
	Human IL-2 enzyme-linked immunoassay kit	EH0013	Human	ELISA
	Human IL-6 enzyme-linked immunoassay kit	EH0001	Human	ELISA
	Human IL-8/CXCL8 enzyme-linked immunoassay kit	EH0014	Human	ELISA
	Human Leptin enzyme-linked immunoassay kit	EH0056	Human	ELISA
	Biotin Conjugated Human Leptin Receptor Recombinant Rabbit mAb [PSH14-72] - Detector	HA725060B	Human	ELISA(Det),ELISA
	Human Leptin Receptor Recombinant Rabbit mAb [PSH14-71] - BSA and Azide free (Capture)	HA725058	Human	ELISA(Cap)
	Human Leptin Receptor Recombinant Rabbit mAb [PSH14-72] - BSA and Azide free (Detector)	HA725059	Human	ELISA(Det)
	Human MMP-9 Recombinant Rabbit mAb [PSH04-24] - BSA and Azide free (Detector)	HA722107	Human	ELISA(Det)
	Human MMP-9 Recombinant Rabbit mAb [PSH04-63] - BSA and Azide free (Detector)	HA722163	Human	ELISA(Det)
	Biotin Conjugated Human MMP-9 Recombinant Rabbit mAb [PSH04-24]	HA722615B	Human	ELISA(Det)
	Human MMP-9 Recombinant Rabbit mAb [PSH04-23] - BSA and Azide free (Capture)	HA722106	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human MMP-9 Recombinant Rabbit mAb [PSH04-62] - BSA and Azide free (Capture)	HA722162	Human	ELISA(Cap)
	Human MMP-9 enzyme-linked immunoassay kit	EH0020	Human	ELISA
	Human/Monkey Albumin Recombinant Rabbit mAb [PSH20-16] - BSA and Azide free (Capture)	HA724127	Human, Monkey	ELISA(Cap)
	Human/Monkey Albumin Recombinant Rabbit mAb [PSH20-17] - BSA and Azide free (Detector)	HA724128	Human, Monkey	ELISA(Det)
	Biotin Conjugated Human/Monkey Albumin Recombinant Rabbit mAb [PSH20-17] - Detector	HA724128B	Human, Monkey	ELISA(Det),ELISA
	Human/Mouse/Rat/Porcine/Canine TGF- β 1 enzyme-linked immunoassay kit	EH0012	Human, Mouse, Rat, Porcine, Canine	ELISA
	Human Neuropilin-1 Recombinant Rabbit mAb [PSH14-73] - BSA and Azide free (Capture)	HA725112	Human	ELISA(Cap)
	Human Neuropilin-1 Recombinant Rabbit mAb [PSH14-74] - BSA and Azide free (Detector)	HA725113	Human	ELISA(Det)
	HRP Conjugated Human Neuropilin-1 Recombinant Rabbit mAb [PSH14-74] - Detector	HA725114H	Human	ELISA(Det),ELISA
	Human Osteopontin (OPN) enzyme-linked immunoassay kit	EH0057	Human	ELISA
	Human Osteoprotegerin/TNFRSF11B enzyme-linked immunoassay kit	EH0070	Human	ELISA
	Human PDGF-AA enzyme-linked immunoassay kit	EH0058	Human	ELISA
	HRP Conjugated Human Prealbumin Recombinant Rabbit mAb [PSH17-43] - Detector	HA725260H	Human	ELISA(Det),ELISA
	Human Prealbumin Recombinant Rabbit mAb [PSH17-42] - BSA and Azide free (Capture)	HA725259	Human	ELISA(Cap)
	Human Prealbumin Recombinant Rabbit mAb [PSH17-43] - BSA and Azide free (Detector)	HA725260	Human	ELISA(Det)
	Biotin Conjugated Human Prealbumin Recombinant Rabbit mAb [PSH17-43] - Detector	HA725260B	Human	ELISA(Det),ELISA
	Human Serum Albumin Mouse mAb [8-A3]	M0809-4	Human	WB,IHC-P
	Human TGF- α enzyme-linked immunoassay kit	EH0059	Human	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human TNF Receptor II Recombinant Rabbit mAb [PSH11-93] - Detector	HA723404B	Human	ELISA(Det),ELISA
	Human TNF Receptor II Recombinant Rabbit mAb [PSH11-93] - BSA and Azide free (Detector)	HA723403	Human	ELISA(Det)
	Human TNF Receptor II Recombinant Rabbit mAb [PSH11-94] - BSA and Azide free (Capture)	HA723405	Human	ELISA(Cap)
	Human TNF Receptor II Recombinant Rabbit mAb [PSH11-92] - BSA and Azide free (Capture)	HA723402	Human	ELISA(Cap)
	Human TNF- α enzyme-linked immunoassay kit	EH0002	Human	ELISA
	Human uPAR enzyme-linked immunoassay kit	EH0077	Human	ELISA
	Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - BSA and Azide free (Detector)	HA725065	Human	ELISA(Det)
	HRP Conjugated Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - Detector	HA725066H	Human	ELISA(Det),ELISA
	Human VCAM-1 Recombinant Rabbit mAb [PSH12-46] - BSA and Azide free (Capture)	HA725064	Human	ELISA(Cap)
	Biotin Conjugated Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - Detector	HA725066B	Human	ELISA(Det),ELISA
	Human VCAM-1 / CD106 enzyme-linked immunoassay kit	EH0067	Human	ELISA
	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-92] - Detector	HA723495B	Human	ELISA(Det),ELISA
	Human VEGFA Recombinant Rabbit mAb [PSH12-91] - BSA and Azide free (Capture)	HA723493	Human	ELISA(Cap)
	Human VEGFA Recombinant Rabbit mAb [PSH12-92] - BSA and Azide free (Detector)	HA723494	Human	ELISA(Det)
	Human VEGFA Recombinant Rabbit mAb [PSH12-93] - BSA and Azide free (Detector)	HA723496	Human	ELISA(Det)
	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-93] - Detector	HA723497B	Human	ELISA(Det),ELISA
	Human VEGF enzyme-linked immunoassay kit	EH0075	Human	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Von Willebrand Factor Recombinant Rabbit mAb [PSH09-04] - BSA and Azide free (Detector)	HA723059	Human	ELISA(Det)
	Biotin Conjugated Human Von Willebrand Factor Recombinant Rabbit mAb [PSH09-04] - BSA and Azide free (Detector)	HA723060B	Human	ELISA(Det),ELISA
	HRP Conjugated Human Von Willebrand Factor Recombinant Rabbit mAb [PSH09-04] - Detector	HA723431H	Human	ELISA(Det)
	Human Von Willebrand Factor Recombinant Rabbit mAb [PSH09-03] - BSA and Azide free (Capture)	HA723058	Human	ELISA(Cap)
KD	Integrin alpha V Recombinant Rabbit mAb [SC56-07]	ET1610-15	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
KD	Integrin alpha V Recombinant Rabbit mAb [SC56-07] - BSA and Azide free	HA750201	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
	Integrin Alpha V Rabbit pAb	ER1911-52	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,IP
Relative expression	Integrin beta 3 Recombinant Rabbit mAb [PSH09-27]	HA723085	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
KD	Integrin beta 3 Recombinant Rabbit mAb [SJ19-09]	ET1606-49	Human, Mouse	WB,IHC-P,IF-Tissue
	Integrin beta 3 Recombinant Rabbit mAb [JE22-64]	ET7110-61	Human, Mouse	WB,IHC-P
	Integrin beta 3 Rabbit pAb	ER1601-4	Human, Mouse, Rat	WB,IHC-P
	IREB2 Recombinant Mouse mAb [9D1-R]	HA601271	Human, Mouse, Rat	WB
	IREB2 Mouse mAb [1G3]	RT1335	Human	WB,IP,IF,IHC-P
	IREB2 Mouse mAb [9D1]	HA600015	Human, Mouse, Rat	WB,IHC-P
	IREB2 Rabbit pAb	R1706-12	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	IREB2 Recombinant Mouse mAb [9D1-R] - BSA and Azide free	HA610146	Human, Mouse, Rat	WB
	KCNA5 / KV1-5 Rabbit pAb	ER1803-23	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	KCNA5 / KV1-5 Recombinant Rabbit mAb [PSH12-30] - BSA and Azide free	HA751444	Human, Mouse, Rat	WB
Relative expression	KCNA5 / KV1-5 Recombinant Rabbit mAb [PSH12-30]	HA723455	Human, Mouse, Rat	WB
KD	Lactate Dehydrogenase Recombinant Rabbit mAb [SU39-06]	ET1608-57	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	Lactate Dehydrogenase Recombinant Rabbit mAb [SU39-06] - BSA and Azide free	HA750155	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	LDHA Recombinant Rabbit mAb [PSH13-58] - BSA and Azide free	HA751477	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	LDHA Mouse mAb [20-179]	M1007-7	Human, Mouse	IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	LDHA Rabbit pAb	ER00702	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	LDHA Recombinant Rabbit mAb [PSH13-58]	HA723547	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	LDHA + LDHB Recombinant Rabbit mAb [PSH12-75] - BSA and Azide free	HA751451	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	LDHA + LDHB Recombinant Rabbit mAb [PSH12-75]	HA723486	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
KD	LDHB Rabbit pAb	0807-1	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
	LDL Receptor Rabbit pAb	R1512-12	Human, Mouse	WB,IF-Cell,FC
	LDL Receptor Recombinant Rabbit mAb [SJ0197] - BSA and Azide free	HA750107	Human, Rat	WB,IF-Cell,FC
	LDL Receptor Recombinant Rabbit mAb [SJ0197]	ET1606-47	Human, Rat	WB,IF-Cell,FC
	Leptin Receptor Recombinant Rabbit mAb [JA73-01] - BSA and Azide free	HA750413	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Leptin Receptor Recombinant Rabbit mAb [JA73-01]	ET1704-44	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	LYVE1 Recombinant Rabbit mAb [PSH18-29]	HA723982	Human	IHC-P
KD	LYVE1 Recombinant Rabbit mAb [JF0979] - BSA and Azide free	HA750343	Human	WB
KD	LYVE1 Recombinant Rabbit mAb [JF0979]	ET1702-29	Human	WB
	LYVE1 Recombinant Rabbit mAb	IRS256RB	Human	mIHC
Cell treatment, KD	MMP-9 Recombinant Rabbit mAb [JA80-73]	ET1704-69	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC
Cell treatment	MMP-9 Recombinant Rabbit mAb [PSH06-33]	HA722665	Human	WB,IHC-P
	MMP-9 Recombinant Mouse mAb [10A1]	EM1801-22	Human	ELISA
	MMP-9 Recombinant Rabbit mAb	IRS172RB	Human	mIHC
	MMP-9 Mouse mAb [2G3]	RT1401	Human	WB,IP,IF,IHC-P
Cell treatment, KD	MMP-9 Recombinant Rabbit mAb [JA80-73] - BSA and Azide free	HA750420	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC
Cell treatment	MMP-9 Recombinant Rabbit mAb [PSH06-33] - BSA and Azide free	HA751070	Human	WB,IHC-P
	Mouse CCL2/MCP-1 enzyme-linked immunoassay kit	EM0025	Mouse	ELISA
	Mouse CRP/C-Reactive Protein enzyme-linked immunoassay kit	EM0035	Mouse	ELISA
	Mouse IFN-γ enzyme-linked immunoassay kit	EM0009	Mouse	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IL-18 enzyme-linked immunoassay kit	EM0034	Mouse	ELISA
	Mouse IL-2 enzyme-linked immunoassay kit	EM0002	Mouse	ELISA
	Mouse IL-6 enzyme-linked immunoassay kit	EM0004	Mouse	ELISA
	Mouse Leptin enzyme-linked immunoassay kit	EM0027	Mouse	ELISA
	HRP Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - Detector	HA723788H	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-80] - Detector	HA723786B	Mouse	ELISA(Det),ELISA
	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - BSA and Azide free (Detector)	HA723787	Mouse	ELISA(Det)
	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-79] - BSA and Azide free (Capture)	HA723784	Mouse	ELISA(Cap)
	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-80] - BSA and Azide free (Capture/Detector)	HA723785	Mouse	ELISA(Det),ELISA(Cap)
	Biotin Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - Detector	HA723788B	Mouse	ELISA(Det),ELISA
	Mouse TNF- α enzyme-linked immunoassay kit	EM0010	Mouse	ELISA
	Mouse VCAM1 Recombinant Rabbit mAb [PSH16-33] - BSA and Azide free (Capture)	HA723831	Mouse	ELISA(Cap)
	Mouse VCAM1 Recombinant Rabbit mAb [PSH16-31] - BSA and Azide free (Capture)	HA723829	Mouse	ELISA(Cap)
	Mouse VCAM1 Recombinant Rabbit mAb [PSH16-32] - BSA and Azide free (Detector)	HA723830	Mouse	ELISA(Det)
	HRP Conjugated Mouse VCAM1 Recombinant Rabbit mAb [PSH16-32] - Detector	HA723830H	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse VCAM1 Recombinant Rabbit mAb [PSH16-32] - Detector	HA723830B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse VEGF Recombinant Rabbit mAb [PSH02-73]	HA722636B	Mouse	ELISA(Det)
	Mouse VEGF Recombinant Rabbit mAb [PSH02-73] - BSA and Azide free (Detector)	HA721851	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse VEGF Recombinant Rabbit mAb [PSH02-72] - BSA and Azide free (Capture)	HA721850	Mouse	ELISA(Cap)
	Mouse VEGF enzyme-linked immunoassay kit	EM0013	Mouse	ELISA
KD	MYO1B Recombinant Rabbit mAb [JE65-01]	HA721119	Human, Mouse, Rat	WB,IHC-P
Relative expression, KD	Neuropilin 1 Recombinant Rabbit mAb [ST05-30]	ET1609-69	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Neuropilin 1 Mouse mAb [A3-D1]	M0910-5	Human, Mouse	WB,IF-Cell,IHC-P
Relative expression, KD	Neuropilin 1 Recombinant Rabbit mAb [ST05-30] - BSA and Azide free	HA750189	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04]	ET1603-12	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04] - BSA and Azide free	HA750058	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	NF-kB p65 Rabbit pAb	ER0815	Human, Mouse, Rat, Zebrafish	WB,IHC-P,FC
	NF-kB p65 Rabbit pAb	HA500402	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27]	HA721307	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27] - BSA and Azide free	HA750601	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
KD	OS9 Recombinant Rabbit mAb [JG38-07]	ET7108-61	Human, Mouse, Rat	WB,IHC-P,IP
	PDGFR beta Rabbit pAb	R1510-44	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	PDGFR beta Recombinant Rabbit mAb	IRS038	Human	mIHC
	PDGFR beta Recombinant Antibody - Rat IgG1 (Chimeric)	IRS283RT	Mouse	mIHC
	PDGFR beta Recombinant Rabbit mAb [SY10-08]	ET1605-20	Human, Mouse, Rat	WB,IHC-P,IP,FC,IF-Tissue,IF-Cell,mIHC,IHC-Fr
	PDGFR beta Recombinant Rabbit mAb [SY10-08] - BSA and Azide free	HA750082	Human, Mouse, Rat	WB,IHC-P,IP,FC,IF-Tissue,IF-Cell,IHC-Fr
	PDGFR beta Recombinant Rabbit mAb	IRS058	Mouse	mIHC
KO	PFKFB3 Recombinant Rabbit mAb [JM43-43]	ET1705-66	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP,IF-Tissue
	PHD1 / prolyl hydroxylase Recombinant Rabbit mAb [SP00-48]	ET1604-5	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-AKT1 (S124) Recombinant Rabbit mAb [JJ08-46]	ET1701-36	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-AKT1 (S129) Recombinant Rabbit mAb [JE59-50]	HA721677	Human, Mouse, Rat	WB
Cell treatment	Phospho-AKT1 (T450) Recombinant Rabbit mAb [SD08-12]	ET1612-73	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell,FC
Cell treatment	Phospho-AKT2 (S474) Recombinant Rabbit mAb [PSH09-46]	HA723108	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-AKT2 (S474) Recombinant Rabbit mAb [PSH09-46] - BSA and Azide free	HA751293	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-AKT2 (S474) Rabbit pAb	HA500116	Human, Mouse	WB,IHC-P
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44] - BSA and Azide free	HA750953	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05]	ET1607-73	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44]	HA722129	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05] - BSA and Azide free	HA750132	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10] - BSA and Azide free	HA751205	Human	WB,IF-Cell,FC
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10]	HA722951	Human	WB,IF-Cell,FC
	Phospho-AMPK alpha 1 (S485) Recombinant Rabbit mAb [PSH18-01]	HA723966	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-AMPK alpha 1 (S496) Recombinant Rabbit mAb [SD0810]	ET1612-72	Human, Mouse	WB,IP,IHC-Fr
Cell treatment	Phospho-AMPK alpha 1 (S496) Recombinant Rabbit mAb [SD0810] - BSA and Azide free	HA750300	Human, Mouse	WB,IP,IHC-Fr
Cell treatment	Phospho-AMPK alpha 1 (T183) +AMPK alpha 2 (T172) Recombinant Rabbit mAb [PSH14-13]	HA723603	Human, Mouse, Rat	WB,IHC-P
	Phospho-AMPK alpha 1 (T183) +AMPK alpha 2 (T172) Recombinant Rabbit mAb [PSH20-31]	HA724142	Human	WB
Cell treatment	Phospho-AMPK alpha 1 (T183) +AMPK alpha 2 (T172) Recombinant Rabbit mAb [PSH14-13] - BSA and Azide free	HA751508	Human, Mouse, Rat	WB,IHC-P
	Phospho-AMPK alpha 2 (S345) Recombinant Rabbit mAb [JJ08-19]	ET1701-37	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-AMPK alpha 2 (S491) Recombinant Rabbit mAb [JE59-49]	HA722567	Human	WB
Cell treatment	Phospho-NF-kB p65 (S529) Recombinant Rabbit mAb [SP07-00]	ET1604-27	Human	WB,IHC-P,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58]	HA723223	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58] - BSA and Azide free	HA751361	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-PDGFR beta (Y1009) Recombinant Rabbit mAb [JE46-72]	HA721736	Human, Mouse, Rat	WB
Cell treatment	Phospho-PDGFR beta (Y1021) Recombinant Rabbit mAb [JE48-54]	HA721833	Mouse, Rat, Human	WB
Cell treatment	Phospho-PDGFR beta (Y740) Recombinant Rabbit mAb [SN0646]	ET1611-29	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue
Cell treatment	Phospho-PDGFR beta (Y751) Recombinant Rabbit mAb [PSH13-61] - BSA and Azide free	HA751480	Mouse, Rat, Human	WB
Cell treatment	Phospho-PDGFR beta (Y751) Recombinant Rabbit mAb [PSH13-61]	HA723550	Mouse, Rat, Human	WB
	Prealbumin Recombinant Rabbit mAb [JM11-43] - BSA and Azide free	HA750392	Human	WB,IP,IF-Cell,IF-Tissue,IHC-P,FC
	Prealbumin Recombinant Rabbit mAb [JM11-43]	ET1703-83	Human	WB,IP,IF-Cell,IF-Tissue,IHC-P,FC
	Prealbumin Rabbit pAb	0809-10	Human	WB,IF-Cell
	Prealbumin Mouse mAb [C5-E10]	EM1707-12	Human	WB,IF-Cell,IHC-P,FC
	Prealbumin Mouse mAb [3-E11]	M0809-12	Human	WB
	Prealbumin Recombinant Rabbit mAb [JF1016]	ET1702-42	Human	WB,IHC-P,IP
	PRKAA2 Mouse mAb [H5-E11]	EM1708-49	Human, Monkey, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	Rat IFN- γ enzyme-linked immunoassay kit	ER0006	Rat	ELISA
	Rat IL-1 β enzyme-linked immunoassay kit	ER0001	Rat	ELISA
	Rat IL-2 enzyme-linked immunoassay kit	ER0002	Rat	ELISA
	Rat IL-6 enzyme-linked immunoassay kit	ER0004	Rat	ELISA
	Rat TNF- α enzyme-linked immunoassay kit	ER0007	Rat	ELISA
	Serum Albumin Rabbit pAb	R1212-2	Human	WB
	Somatostatin Receptor 2 Recombinant Rabbit mAb [JE03-37]	HA722521	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
	SSTR2 Rabbit pAb	ER65707	Human, Mouse, Rat	WB
	TGF beta 1 Rabbit pAb	HA500496	Human	WB
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17] - BSA and Azide free	HA750485	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	TGF beta 1 Mouse mAb [B2-F6]	EM1709-93	Human	WB,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	TGF beta 1 Recombinant Rabbit mAb	IRS205RB	Human	mIHC
	Biotin Conjugated TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - Detector	HA723727B	Human, Mouse, Rat	ELISA(Det),ELISA
	TGF beta 1 Recombinant Rabbit mAb [PSH18-59]	HA724012	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
	TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - BSA and Azide free (Detector)	HA723726	Human, Mouse, Rat	ELISA(Det)
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17]	HA721143	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	TGF beta 1 Recombinant Rabbit mAb [PSH15-21] - BSA and Azide free (Capture)	HA723725	Human, Mouse, Rat	ELISA(Cap)
	Tissue Factor Recombinant Rabbit mAb [SN20-16] - BSA and Azide free	HA750278	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Tissue Factor Recombinant Rabbit mAb [SN20-16]	ET1611-95	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	TNF Receptor II Rabbit pAb	HA500332	Human, Mouse, Rat	IHC-P,FC
	TNF Receptor II Recombinant Rabbit mAb [JM113-01]	ET1703-63	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP,IHC-Fr
	TNF Receptor II Recombinant Rabbit mAb [JM113-01] - BSA and Azide free	HA750386	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP,IHC-Fr
Relative expression	TNF Receptor II Recombinant Rabbit mAb [PSH12-12] - BSA and Azide free	HA751441	Human	FC
Relative expression	TNF Receptor II Recombinant Rabbit mAb [PSH12-12]	HA723428	Human	FC
KD	Transferrin Receptor (CD71) Recombinant Rabbit mAb [JF0956]	ET1702-06	Human, Mouse, Rat	WB,IP,FC
Relative expression	Transferrin Receptor (CD71) Recombinant Rabbit mAb [PSH13-56]	HA723545	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP,IF-Tissue
KD	Transferrin Receptor (CD71) Recombinant Rabbit mAb [JF0956] - BSA and Azide free	HA750333	Human, Mouse, Rat	WB,IP,FC
	Transferrin Receptor (CD71) Rabbit pAb	ER1706-35	Human	IHC-P
	VCAM1 Rabbit pAb	R1512-13	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression, KD	VCAM1 Recombinant Rabbit mAb [SA05-04] - BSA and Azide free	HA750009	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Relative expression, KD	VCAM1 Recombinant Rabbit mAb [SA05-04]	ET1601-18	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	VEGF Recombinant Rabbit mAb [SP07-01]	ET1604-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	VEGF Rabbit pAb	ER30607	Human	WB,IF-Cell,IHC-P,FC
	VEGF Recombinant Rabbit mAb [SP07-01] - BSA and Azide free	HA750074	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	VEGFA Recombinant Rabbit mAb [PSH19-17]	HA724051	Human, Mouse, Rat, Monkey	WB
Relative expression	Von Willebrand Factor Recombinant Rabbit mAb [PSH07-45]	HA722833	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP
	Von Willebrand Factor Recombinant Rabbit mAb	IRS244RB	Mouse	mIHC
	Von Willebrand Factor Recombinant Rabbit mAb [JJ09-34]	ET1701-87	Human, Mouse	WB
	Von Willebrand Factor Recombinant Rabbit mAb	IRS215RB	Human	mIHC
Relative expression	Von Willebrand Factor Recombinant Rabbit mAb [PSH07-45] - BSA and Azide free	HA751166	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层