

表观遗传与核信号研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
Cell treatment	3mC Recombinant Rabbit mAb [PSH03-21]	HA721957	Species independent	IF-Cell, Dot Blot
Cell treatment	3mC Recombinant Rabbit mAb [PSH03-22]	HA721958	Species independent	ELISA, IF-Cell, Dot Blot
	ALKBH5 Recombinant Rabbit mAb [PSH01-16]	HA721628	Human, Mouse, Rat	WB, IHC-P
KD	ATF6 Recombinant Mouse mAb [8D3-R]	HA601321	Human, Mouse, Rat	WB, IF-Cell
	BRCC36 Recombinant Rabbit mAb [PSH05-10]	HA722226	Human, Mouse, Rat	WB, IF-Cell, IHC-P, FC
	CARM1 Recombinant Rabbit mAb [PSH05-08]	HA722224	Human, Mouse, Rat	WB, IHC-P, IP
	CGBP Recombinant Rabbit mAb [PSH01-54]	HA721697	Human, Mouse, Rat	WB, IHC-P
	CNOT7 Recombinant Rabbit mAb [PSH02-81]	HA721861	Human, Mouse, Rat	WB, IF-Cell, FC
	DLX3 Recombinant Rabbit mAb [JE64-45]	HA722281	Human	WB, IF-Cell, IHC-P, FC, IF-Tissue
KD	Dnmt2 Recombinant Rabbit mAb [PSH03-43]	HA721988	Human, Monkey	WB, IF-Cell, IHC-P, FC
	DOT1L Recombinant Rabbit mAb [PSH01-57]	HA721700	Human, Mouse, Rat	WB, IHC-P
	E2F4 Recombinant Rabbit mAb [JE40-78]	HA721374	Human	WB, IHC-P, FC
	EIF2S1 Recombinant Rabbit mAb [PSH04-29]	HA722112	Human, Mouse, Rat, Monkey	WB, IF-Cell, IHC-P, FC, IP
	eIF3A Recombinant Rabbit mAb [JE58-60]	HA721375	Human, Mouse, Rat	WB, IHC-P, IF-Cell
	eIF4E Recombinant Rabbit mAb [JE60-70]	HA721192	Human, Mouse, Rat	WB, IHC-P, IF-Cell, FC
	Elp4 Recombinant Rabbit mAb [JE61-03]	HA720054	Human, Mouse	WB
Relative expression, KD	Emerin Recombinant Rabbit mAb [PSH01-86]	HA721744	Human	WB, IHC-P, IF-Cell, FC
	EZH1 Recombinant Rabbit mAb [PS01-09]	HA721187	Human, Mouse, Rat	WB, IHC-P, IF-Cell
	FBXO3 Recombinant Rabbit mAb [PSH04-09]	HA722090	Human, Mouse	WB
	FOXC1 Recombinant Rabbit mAb [PSH03-89]	HA722060	Human, Mouse, Rat	WB, IF-Cell, IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	GATA2 Recombinant Rabbit mAb [PSH04-53]	HA722152	Human, Rat	WB,IF-Cell
	Geminin Recombinant Rabbit mAb [JE35-51]	HA722239	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Histone H2A (acetyl K9) Recombinant Rabbit mAb [PS01-45]	HA721284	Human, Mouse, Rat	WB,IHC-P,IF-Cell,CUT&Tag-seq,ChIP
Cell treatment	Histone H3 (acetyl K27) Mouse mAb [A6D7]	HA600047	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,ChIP
Cell treatment	Histone H3 (acetyl K27) Mouse mAb [A6D6]	HA600048	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,ChIP
Cell treatment	Histone H3 (acetyl K4) Recombinant Rabbit mAb [PSH03-71]	HA722035	Human, Mouse, Rat	WB,IF-Cell,ChIP,Dot Blot
Cell treatment	Histone H3 (acetyl K9) Recombinant Rabbit mAb [PSH04-47] - ChIP Grade	HA722132	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,FC,ChIP,Dot Blot,IP
Cell treatment	Histone H3 (di methyl K27) Recombinant Rabbit mAb [PSH05-88]	HA722486	Human, Mouse, Rat	WB,IF-Cell,ChIP
	Histone H3 (di methyl K4) Recombinant Rabbit mAb [JE00-98]	HA722135	Human, Mouse, Rat, Monkey	WB,IF-Cell,ChIP,Dot Blot,IHC-P
	Histone H3 (mono methyl K14) Recombinant Rabbit mAb [JE43-29]	HA722232	Human, Mouse, Rat	WB,IF-Cell,ChIP,Dot Blot
	Histone H3 (mono methyl K27) Recombinant Rabbit mAb [PSH05-67]	HA722376	Human, Mouse, Rat	WB,IF-Cell,FC,ChIP
	Histone H3 (mono methyl K9) Recombinant Rabbit mAb [JE43-27]	HA722350	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,ChIP,IF-Tissue,Dot Blot
Cell treatment, Relative expression	Histone H3 (mutated K27 M) Recombinant Rabbit mAb [JE46-51]	HA722200	Human, Mouse	WB,IF-Cell,IHC-P,ChIP
Cell treatment	Histone H3 (tri methyl K27) Recombinant Rabbit mAb [PSH05-15]	HA722231	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,ChIP
	HMG2 Recombinant Rabbit mAb [PSH05-53]	HA722295	Human, Mouse	WB,IF-Cell,IHC-P,FC,IF-Tissue
KD	HNRNPA0 Recombinant Rabbit mAb [PSH02-77]	HA721855	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	hnRNP G-T Recombinant Rabbit mAb [JE64-33]	HA721104	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	HNRNPM Recombinant Rabbit mAb [JE64-62]	HA721103	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	hSET1/SET1 Recombinant Rabbit mAb [PH00-02]	HA721226	Human, Mouse, Rat	WB,FC
	hSET1/SET1 Recombinant Rabbit mAb [PH00-01]	HA721224	Human, Mouse, Rat	WB,IHC-P
	IKZF3 Recombinant Rabbit mAb [JE62-81]	HA720110	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Kaiso Recombinant Rabbit mAb [JE38-42]	HA722275	Human, Rat	WB,IF-Cell,FC
	KDM2B Recombinant Rabbit mAb [PSH05-38]	HA722267	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,IP,ChIP
	KDM3A / JHDM2A Recombinant Rabbit mAb [PSH01-77]	HA721729	Human, Mouse	WB,IHC
	KDM4A Recombinant Rabbit mAb [PSH03-67]	HA722026	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	KDM4B Recombinant Rabbit mAb [JE36-58]	HA721350	Human, Mouse	WB,IHC-P,IF-Cell
	KDM5B Recombinant Rabbit mAb [PSH03-97]	HA722076	Human	WB,IP
	KMT2C / MLL3 Recombinant Rabbit mAb [PSH02-11]	HA721770	Human, Rat, Mouse	WB
	KMT3B Recombinant Rabbit mAb [1B1]	HA721166	Human	IF-Cell
	KMT6 / EZH2 Recombinant Rabbit mAb [PSH04-14]	HA722095	Human, Mouse, Rat, Monkey	WB,IF-Cell
KD	LSD2 Recombinant Rabbit mAb [JE32-40]	HA721462	Human	WB
	MBD1 Recombinant Rabbit mAb [PSH01-78]	HA721730	Human, Mouse, Rat	WB
	MCM2 Recombinant Rabbit mAb [PD00-90]	HA721177	Human	WB,IHC-P,IF-Cell,FC
KD	Menin Recombinant Rabbit mAb [JE30-97]	HA721313	Human	WB,IHC-P
Relative expression	MyoD1 Recombinant Rabbit mAb [JE60-47]	HA722179	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	NSUN2 Recombinant Rabbit mAb [PSH05-05]	HA722219	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,IF-Tissue,IP
KD	NSUN3 Recombinant Rabbit mAb [PSH01-97]	HA721755	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	NSUN5 Rabbit pAb	HA500544	Human, Mouse, Rat	WB,IHC-P
	NSUN6 Recombinant Rabbit mAb [PSH02-80]	HA721858	Human, Mouse, Rat	WB
	PABP Recombinant Rabbit mAb [PSH05-51]	HA722289	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC,IF-Tissue
	PACT (PKR activating protein) / PRKRA Recombinant Rabbit mAb [PSH02-15]	HA721774	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
Cell treatment	PADI4 Recombinant Rabbit mAb [PSH01-33]	HA721657	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Pan-Lactyl-lysine Recombinant Rabbit mAb [PSH03-74]	HA722038	Species independent	WB,IHC-P,IF-Cell,Dot Blot

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Pan-Lactyl-lysine Recombinant Rabbit mAb [PSH03-73]	HA722037	Species independent	WB,IHC-P,IF-Cell,Dot Blot
Cell treatment	Phospho-eIF4G (S1108) Recombinant Rabbit mAb [PS01-08]	HA721276	Human	WB,IHC-P
Cell treatment	Phospho-Histone H3.3 (S31) Recombinant Rabbit mAb [JE62-92]	HA721034	Human, Rat	WB,IF-Cell
KD	PML Protein Recombinant Rabbit mAb [PSH02-89]	HA721877	Human	WB,IHC-P,IF-Cell,FC
	PRMT6 Recombinant Rabbit mAb [JE36-57]	HA721333	Human, Mouse, Rat	WB
Relative expression	PSIP1 Recombinant Rabbit mAb [PSH02-55]	HA721831	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	RING1 Recombinant Rabbit mAb [JE34-66]	HA721461	Human	WB,IHC-P
	ROR gamma T Recombinant Rabbit mAb [PSH04-38]	HA722121	Human, Rat	WB,IHC-P
KD	SET7 Recombinant Rabbit mAb [PSH02-90]	HA721909	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
KD	SIRT3 Recombinant Rabbit mAb [PSH05-22]	HA722251	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue,IP,ChIP
	SMCHD1 Recombinant Rabbit mAb [JE34-29]	HA722141	Human	WB,IHC-P
	SP110 Recombinant Rabbit mAb [PSH03-66]	HA722025	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	TEAD1 Recombinant Rabbit mAb [JE45-37]	HA722144	Human, Mouse, Rat	WB,IHC-P
	TFIIB Recombinant Rabbit mAb [JE38-37]	HA722213	Human, Mouse, Rat	WB
	Transcription factor AP-2-alpha Recombinant Rabbit mAb [PSH02-06]	HA721765	Human	WB,IHC-P,IF-Cell
	Transcription factor AP-2-alpha Recombinant Rabbit mAb [PSH02-05]	HA721764	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	USP39 Recombinant Rabbit mAb [JE59-07]	HA720034	Human, Mouse, Rat	WB,IHC-P
	WDR77 Recombinant Rabbit mAb [PSH0-07]	HA721260	Human	WB,IHC-P,FC
	YTHDF1 Recombinant Rabbit mAb [PSH0-23]	HA721302	Human, Mouse, Rat	WB,IHC-P,FC
	ZBTB7A Recombinant Rabbit mAb [JE34-74]	HA722205	Human	WB,IF-Cell,IHC-P,IP
	ZEB1 Recombinant Rabbit mAb [JE57-64]	HA721318	Human	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

