

组蛋白研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	Acetyl-Histone H3 Antibody Sampler Kit	HAK21081	Human, Mouse, Rat	WB,IHC-P,ChIP
	H2B1C Rabbit pAb	ER62473	Human, Mouse	WB
	HDAC5/9 Rabbit pAb	ER63671	Human, Mouse	WB,IHC-P
KD	HDAC9 Recombinant Rabbit mAb [JU30-44]	ET1706-36	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Histone deacetylase Rabbit pAb	HA500121	Rice	WB
	Histone deacetylase 9 Rabbit pAb	ER63682	Human	WB,IHC-P,IF-Cell
	Histone H1.0 Recombinant Rabbit mAb [SD206-04]	ET1612-24	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
KD	Histone H1.2 Recombinant Rabbit mAb [JU43-48]	ET1706-26	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP
	Histone H1.3 Rabbit pAb	ER1919-58	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Histone H1 (Acetyl K25) Rabbit pAb	ER60044	Human, Monkey	WB,IHC-P
	Histone H1 (di methyl K25) Rabbit pAb	ER60071	Human, Mouse, Rat	WB
	Histone H1 (tri methyl K25) Rabbit pAb	ER60076	Human	WB
	Histone H1.X Rabbit pAb	ER1910-65	Human, Mouse	WB,IHC-P,IF-Tissue
	Histone H2A Acetyl Rabbit pAb	RT1283	Broad mammalian, Zebrafish, Drosophila, Xenopus, C. elegans	WB,IP,IF
	Histone H2A (Acetyl K15) Rabbit pAb	ER60073	Human, Mouse, Rat	WB
	Histone H2A (Acetyl K5) Rabbit pAb	ER60045	Human, Mouse, Rat	IHC-P,WB
	Histone H2A (Acetyl K5) Rabbit pAb	ER60072	Human, Mouse, Rat	WB
Cell treatment	Histone H2A (acetyl K9) Recombinant Rabbit mAb [SR4-15]	ET1602-34	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP,CUT&Tag-seq
Cell treatment	Histone H2A (acetyl K9) Recombinant Rabbit mAb [PS01-45] - BSA and Azide free	HA750591	Human, Mouse, Rat	WB,IHC-P,IF-Cell,CUT&Tag-seq,ChIP
Cell treatment	Histone H2A (acetyl K9) Recombinant Rabbit mAb [PS01-45]	HA721284	Human, Mouse, Rat	WB,IHC-P,IF-Cell,CUT&Tag-seq,ChIP
	Histone H2A (hydroxyl Y39) Recombinant Rabbit mAb [SR4-17]	ET1602-35	Human, Mouse, Rat	WB,IHC-P
	Histone H2A.x Rabbit pAb	ER1901-70	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC

	Histone H2A.X Recombinant Rabbit mAb [JM06-42]	ET1705-97	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
	Histone H2A.Z Recombinant Rabbit mAb [ST46-09]	ET1609-23	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,ChIP
	Histone H2A.Z (Acetyl K4) Rabbit pAb	ER60075	Human, Mouse, Rat	WB
	Histone H2A.Z (Acetyl K7) Rabbit pAb	ER60074	Human, Mouse, Rat	WB
	Histone H2B Rabbit pAb	HA500381	Human, Mouse, Rat, Zebrafish	WB,IHC-P
	Histone H2B Recombinant Rabbit mAb [SD20-63] - BSA and Azide free	HA750287	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP,ChIP
	Histone H2B Recombinant Rabbit mAb [SD20-63]	ET1612-25	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP,ChIP
	Histone H2B (Acetyl K12) Rabbit pAb	ER60047	Human, Mouse, Monkey	WB,IHC-P,IF-Cell
	Histone H2B (Acetyl K15) Rabbit pAb	ER60048	Human, Mouse	WB,IHC-P,IF-Cell
Cell treatment	Histone H2B (acetyl K20) Recombinant Rabbit mAb [PS01-46]	HA721285	Human, Mouse, Rat	WB,IHC-P,IF-Cell,ChIP
Cell treatment	Histone H2B (acetyl K20) Recombinant Rabbit mAb [SN72-09]	ET1611-48	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP,ChIP,CUT&Tag-seq
Cell treatment	Histone H2B (acetyl K20) Recombinant Rabbit mAb [PS01-46] - BSA and Azide free	HA750592	Human, Mouse, Rat	WB,IHC-P,IF-Cell,ChIP
	Histone H2B (Acetyl K5) Rabbit pAb	ER60046	Human, Mouse	WB,IHC-P,IF-Cell
	Histone H2B (mono methyl R79) Recombinant Rabbit mAb [ST0430]	ET1609-9	Human, Mouse, Rat	IF-Cell,IHC-P
	Histone H3 Mouse mAb [6-A7]	M1306-4	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,ChIP
	Histone H3 Recombinant Rabbit mAb [JJ090-07]	ET1701-64	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP,IP
	Histone H3 Mouse mAb [3-C4]	EM30605	Human, Mouse, Rat	WB,IHC-P,mIHC,IF-Tissue,ChIP
	Histone H3 Recombinant Mouse mAb [A11-D7-R]	HA601278	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Histone H3 Recombinant Rabbit mAb [JJ090-07] - BSA and Azide free	HA750318	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP,IP
	Histone H3 Recombinant Mouse mAb [6-A7-R] - BSA and Azide free	HA610203	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,ChIP
	Histone H3 Recombinant Mouse mAb [A11-D7-R] - BSA and Azide free	HA610153	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Histone H3 Recombinant Mouse mAb [6-A7-R] - Loading control	HA601335	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,ChIP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Histone H3 Recombinant Rabbit mAb [JJ092-08]	ET1701-81	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,ChIP
	Histone H3 Mouse mAb [A11-D7]	M1309-1	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,IF-Tissue
	Histone H3 Rabbit pAb	R1105-1	Human, Mouse, Rat	WB
	Histone H3.1 Rabbit pAb	ER1802-79	Human, Mouse	WB,IF-Cell,IHC-P
	Histone H3.3 Recombinant Rabbit mAb [ST50-08] - BSA and Azide free	HA750184	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,ChIP
	Histone H3.3 Recombinant Rabbit mAb [ST50-08]	ET1609-59	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,ChIP
Cell treatment	Histone H3 (acetyl K14) Recombinant Rabbit mAb [JU43-26] - BSA and Azide free	HA750454	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,SNAP-ChIP,CUT&Tag-seq
Cell treatment	Histone H3 (acetyl K14) Recombinant Rabbit mAb [JU43-26]	ET1706-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,SNAP-ChIP,CUT&Tag-seq
Cell treatment	Histone H3 (acetyl K18) Mouse mAb [A7B5]	HA600090	Human, Mouse, Rat	WB,IHC-P,ChIP
Cell treatment	Histone H3 (acetyl K18) Mouse mAb [A7B6]	HA600091	Human	WB,IF-Cell
Cell treatment	Histone H3 (acetyl K18) Rabbit pAb	HA500047	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Histone H3 (acetyl K23) Recombinant Rabbit mAb [JE46-39]	HA721840	Human, Mouse	WB,IF-Cell
Cell treatment	Histone H3 (acetyl K27) Mouse mAb [A6D7]	HA600047	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,ChIP
Cell treatment	Histone H3 (acetyl K27) Rabbit pAb	HA500046	Human, Mouse, Rat	WB,IHC-P,ChIP
Cell treatment	Histone H3 (acetyl K27) Recombinant Rabbit mAb [PSH13-93]	HA723586	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Histone H3 (acetyl K27) Mouse mAb [A6D6]	HA600048	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,ChIP
Cell treatment	Histone H3 (acetyl K27) Recombinant Rabbit mAb [PSH13-93] - BSA and Azide free	HA751495	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Histone H3 (acetyl K36) Recombinant Rabbit mAb [JE77-32]	HA722415	Human, Mouse, Rat	WB,IHC-P,IF-Cell,ChIP
Cell treatment	Histone H3 (acetyl K4) Recombinant Rabbit mAb [PSH03-71] - BSA and Azide free	HA750905	Human, Mouse, Rat	WB,IF-Cell,ChIP,Dot Blot
Cell treatment	Histone H3 (acetyl K4) Recombinant Rabbit mAb [PSH03-71]	HA722035	Human, Mouse, Rat	WB,IF-Cell,ChIP,Dot Blot
Cell treatment	Histone H3 (acetyl K56) Recombinant Rabbit mAb [SU30-10]	ET1608-9	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP,CUT&Tag-seq

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Histone H3 (acetyl K56) Recombinant Rabbit mAb [SU30-10] - BSA and Azide free	HA750163	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP,CUT&Tag-seq
Cell treatment	Histone H3 (acetyl K9) Recombinant Rabbit mAb [PSH04-47] - ChIP Grade	HA722132	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,FC,ChIP,Dot Blot,IP
Cell treatment	Histone H3 (acetyl K9) Recombinant Rabbit mAb [PSH04-47] - BSA and Azide free	HA750956	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,FC,ChIP,Dot Blot,IP
Cell treatment	Histone H3 (di methyl K27) Recombinant Rabbit mAb [PSH05-88]	HA722486	Human, Mouse, Rat	WB,IF-Cell,ChIP
Cell treatment	Histone H3 (di methyl K27) Recombinant Rabbit mAb [PSH05-88] - BSA and Azide free	HA751035	Human, Mouse, Rat	WB,IF-Cell,ChIP
Cell treatment	Histone H3 (di methyl K4) Rabbit pAb	R1110-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,Dot Blot
	Histone H3 (di methyl K4) Recombinant Rabbit mAb [JE00-98]	HA722135	Human, Mouse, Rat, Monkey	WB,IF-Cell,ChIP,Dot Blot,IHC-P
	Histone H3 (di methyl K9) Recombinant Rabbit mAb [SN07-30] - BSA and Azide free	HA750259	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Histone H3 (di methyl K9) Recombinant Rabbit mAb [SN07-30]	ET1611-51	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Histone H3 (mono+di+tri methyl K36) Recombinant Rabbit mAb [PSH07-14] - BSA and Azide free	HA751135	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,ChIP
	Histone H3 (mono+di+tri methyl K36) Recombinant Rabbit mAb [PSH07-14]	HA722802	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,ChIP
	Histone H3 (mono+di+tri methyl K4) Recombinant Rabbit mAb [PSH06-30] - BSA and Azide free	HA751067	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,ChIP
	Histone H3 (mono+di+tri methyl K4) Recombinant Rabbit mAb [PSH06-30]	HA722662	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,ChIP
	Histone H3 (mono+di+tri methyl K79) Recombinant Rabbit mAb [SR42-06] - BSA and Azide free	HA750051	Human, Mouse, Rat	WB,IHC-P,ChIP
	Histone H3 (mono+di+tri methyl K79) Recombinant Rabbit mAb [SR42-06]	ET1602-41	Human, Mouse, Rat	WB,IHC-P,ChIP
	Histone H3 (mono+di+tri methyl K79) Recombinant Rabbit mAb [PSH09-10]	HA723068	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,ChIP
	Histone H3 (mono+di+tri methyl K79) Recombinant Rabbit mAb [PSH09-10] - BSA and Azide free	HA751264	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,ChIP
	Histone H3 (mono methyl K14) Recombinant Rabbit mAb [JE43-29]	HA722232	Human, Mouse, Rat	WB,IF-Cell,ChIP,Dot Blot
	Histone H3 (mono methyl K18) Recombinant Rabbit mAb [SA42-07]	ET1601-14	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Histone H3 (mono methyl K18) Recombinant Rabbit mAb [SA42-07] - BSA and Azide free	HA750006	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P
	Histone H3 (mono methyl K27) Recombinant Rabbit mAb [PSH05-67] - BSA and Azide free	HA751019	Human, Mouse, Rat	WB,IF-Cell,FC,ChIP
	Histone H3 (mono methyl K27) Recombinant Rabbit mAb [PSH05-67]	HA722376	Human, Mouse, Rat	WB,IF-Cell,FC,ChIP
	Histone H3 (mono methyl K36) Recombinant Rabbit mAb [SR04-20]	ET1602-48	Human, Mouse	WB,IF-Cell,IF- Tissue,ChIP
	Histone H3 (mono methyl K4) Rabbit pAb	R1110-1	Human, Mouse, Rat	WB,IF-Cell,IHC- P,FC,Dot Blot
	Histone H3 (mono methyl K4) Recombinant Rabbit mAb [PSH07-15]	HA722803	Human, Mouse, Rat	WB,IF-Cell,ChIP,Dot Blot
	Histone H3 (mono methyl K4) Recombinant Rabbit mAb [PSH07-15] - BSA and Azide free	HA751136	Human, Mouse, Rat	WB,IF-Cell,ChIP,Dot Blot
KD	Histone H3 (mono methyl K79) Recombinant Rabbit mAb [PSH07-66]	HA722848	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC- P,FC,ChIP
KD	Histone H3 (mono methyl K79) Recombinant Rabbit mAb [PSH07-66] - BSA and Azide free	HA751175	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC- P,FC,ChIP
	Histone H3 (mono methyl K9) Recombinant Rabbit mAb [JE43-27]	HA722350	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC- P,FC,ChIP,IF-Tissue,Dot Blot
	Histone H3 (mono methyl R2) Recombinant Rabbit mAb [ST0427] - BSA and Azide free	HA750197	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P,ChIP
	Histone H3 (mono methyl R2) Recombinant Rabbit mAb [ST0427]	ET1609-8	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P,ChIP
Cell treatment, Relative expression	Histone H3 (mutated K27 M) Recombinant Rabbit mAb [JE46-51]	HA722200	Human, Mouse	WB,IF-Cell,IHC-P,ChIP
Cell treatment	Histone H3 (tri methyl K27) Recombinant Rabbit mAb [PSH05-15] - BSA and Azide free	HA750988	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF- Tissue,ChIP
Cell treatment	Histone H3 (tri methyl K27) Recombinant Rabbit mAb [PSH05-15]	HA722231	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF- Tissue,ChIP
	Histone H3 (tri methyl K36) Recombinant Rabbit mAb [PSH15-01] - BSA and Azide free	HA751539	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,IP,Dot Blot,ChIP
	Histone H3 (tri methyl K36) Recombinant Rabbit mAb [PSH15-01]	HA723701	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,IP,Dot Blot,ChIP
	Histone H3 (tri methyl K9) Mouse mAb [2G1]	RT1287	Human	WB,IF
	Histone H3 (tri methyl K9) Mouse mAb [1-6]	M1112-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF- Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Histone H3 (tri methyl K9) Recombinant Rabbit mAb [PSH19-58]	HA724084	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,ChIP
	Histone H4 Recombinant Rabbit mAb [SD209-04] - BSA and Azide free	HA750294	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P,ChIP
	Histone H4 Mouse mAb [H2-E7]	EM1706-19	Human	WB,IF-Cell,IHC-P,FC
	Histone H4 Recombinant Rabbit mAb [SD209-04]	ET1612-43	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P,ChIP
Cell treatment	Histone H4 (acetyl K16) Recombinant Rabbit mAb [JB21-44] - BSA and Azide free	HA750464	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell,IF- Tissue,ChIP
Cell treatment	Histone H4 (acetyl K16) Recombinant Rabbit mAb [JB21-44]	ET7107-89	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell,IF- Tissue,ChIP
Cell treatment	Histone H4 (acetyl K5) Recombinant Rabbit mAb [SR31-07]	ET1602-40	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC- P,IP,ChIP,CUT&Tag-seq
Cell treatment	Histone H4 (acetyl K5) Recombinant Rabbit mAb [SR31-07] - BSA and Azide free	HA750050	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC- P,IP,ChIP,CUT&Tag-seq
	Histone H4 (acetyl K8) Recombinant Rabbit mAb [PS01-25]	HA721249	Human, Mouse, Rat	WB,IF-Cell,IHC-P,ChIP
	Histone H4 (acetyl K8) Recombinant Rabbit mAb [PS01-25] - BSA and Azide free	HA750568	Human, Mouse, Rat	WB,IF-Cell,IHC-P,ChIP
	Histone H4 (tri methyl K20) Recombinant Rabbit mAb [ST044-07]	ET1609-16	Human, Mouse, Rat	IHC-P
	macroH2A.1 Rabbit pAb	ER1803-04	Human, Mouse	WB,IHC-P,FC
	macroH2A.1 Recombinant Rabbit mAb [JU22-83]	ET1706-15	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P
	Mouse/Rat CXCL1/GRO alpha Recombinant Rabbit mAb [PSH13-15] - BSA and Azide free (Detector)	HA725029	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse/Rat CXCL1/GRO alpha Recombinant Rabbit mAb [PSH13-15] - Detector	HA725030B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat CXCL1/GRO alpha Recombinant Rabbit mAb [PSH13-14] - BSA and Azide free (Capture)	HA725028	Mouse, Rat	ELISA(Cap)
Cell treatment	Phospho-Histone H1.3 (T17)+Histone H1.4 (T17) Recombinant Rabbit mAb [SR38-03]	ET1602-11	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P
Cell treatment	Phospho-Histone H2A.X (S139) Recombinant Rabbit mAb [SR33-09] - BSA and Azide free	HA750034	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IF-Cell,IF- Tissue
Cell treatment	Phospho-Histone H2A.X (S139) Recombinant Rabbit mAb [SR33-09]	ET1602-2	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IF-Cell,IF- Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Histone H3.3 (S31) Recombinant Rabbit mAb [JE62-92]	HA721034	Human, Rat	WB,IF-Cell
Cell treatment	Phospho-Histone H3 (S10) Recombinant Rabbit mAb [SA31-01]	ET1601-30	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,ChIP
Cell treatment	Phospho-Histone H3 (S10) Recombinant Rabbit mAb [SA31-01] - BSA and Azide free	HA750020	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,ChIP
Cell treatment	Phospho-Histone H3 (S28) Recombinant Rabbit mAb [JE56-06]	HA721795	Human, Rat, Mouse	WB,IF-Cell,IHC-P,FC,ChIP
Cell treatment	Phospho-Histone H3 (T3) Recombinant Rabbit mAb [JE42-48]	HA721796	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC