

荧光二抗-产品清单

高级验证	产品名	货号	反应种属	验证应用
	iFluor™ 488 Conjugated Alpaca anti-Rabbit IgG FC, Recombinant VHH monoclonal Antibody	HA1129	Rabbit	IF-Cell, IHC-Fr
	AP Conjugated Goat anti-mouse IgG polyclonal Antibody	HA1020	Mouse	WB, ELISA
	cy3 Conjugated Goat anti-mouse IgG polyclonal Antibody	HA1109	Mouse	ICC, IF
	cy5 Conjugated Goat anti-mouse IgG polyclonal Antibody	HA1110	Mouse	IF-Cell, FC
	Biotin Conjugated Goat anti-mouse IgG polyclonal Antibody	HA1108	Mouse	WB, ELISA, IHC-P, IF-Tissue, FC
	Alexa Fluor® 555 Conjugated Goat anti-mouse IgG polyclonal Antibody	HA1118	Mouse	IF-Tissue
	FITC Conjugated Goat anti-Mouse IgG polyclonal Antibody	HA1003	Mouse	IF-Cell
	iFluor™ 647 Conjugated Goat anti-mouse IgG polyclonal Antibody	HA1127	Mouse	IF-Cell, IF-Tissue, IHC-Fr
	iFluor™ 594 Conjugated Goat anti-mouse IgG polyclonal Antibody	HA1126	Mouse	IF-Cell, IF-Tissue, FC, IHC-Fr
	PE Conjugated Goat anti-mouse IgG polyclonal Antibody	HA1015	Mouse	FC, IF
	iFluor™ 488 Conjugated Goat anti-mouse IgG polyclonal Antibody	HA1125	Mouse	IF-Cell, IF-Tissue, FC, IHC-Fr
	TRITC Conjugated Goat Anti-Mouse IgG H&L polyclonal Antibody	HA1017	Mouse	IF-Cell, IF-Tissue, FC
	FITC Conjugated Goat anti-mouse IgG polyclonal Antibody	HA1128	Mouse	IF-Cell, IF-Tissue, FC, IHC-Fr
	Alexa Fluor® 555 Conjugated Goat anti-rabbit IgG polyclonal Antibody	HA1117	Rabbit	IF-Tissue
	AP Conjugated Goat anti-rabbit IgG polyclonal Antibody	HA1019	Rabbit	WB, ELISA
	cy3 Conjugated Goat anti-rabbit IgG polyclonal Antibody	HA1102	Rabbit	IF-Cell, IF-Tissue
	cy5 Conjugated Goat anti-rabbit IgG polyclonal Antibody	HA1103	Rabbit	IF-Cell, IF-Tissue
	FITC Conjugated Goat anti-Rabbit IgG polyclonal Antibody	HA1004	Rabbit	IF-Cell
	Biotin Conjugated Goat anti-rabbit IgG polyclonal Antibody	HA1101	Rabbit	WB, ELISA, IHC-P, IF-Tissue, FC, IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	iFluor™ 488 Conjugated Goat anti-rabbit IgG polyclonal Antibody	HA1121	Rabbit	IF-Cell,IF-Tissue,FC,IHC-Fr,WB
	iFluor™ 594 Conjugated Goat anti-rabbit IgG polyclonal Antibody	HA1122	Rabbit	IF-Cell,IF-Tissue,FC,IHC-Fr
	iFluor™ 647 Conjugated Goat anti-rabbit IgG polyclonal Antibody	HA1123	Rabbit	IF-Cell,IF-Tissue,IHC-Fr,FC
	FITC Conjugated Goat anti-rabbit IgG polyclonal Antibody	HA1124	Rabbit	IF-Cell,IF-Tissue,FC,IHC-Fr
	TRITC Conjugated Goat Anti-Rabbit IgG H&L polyclonal Antibody	HA1016	Rabbit	IF-Cell,IF-Tissue,FC
	PE Conjugated Goat anti-rabbit IgG polyclonal Antibody	HA1014	Rabbit	FC,IF
	iFluor™ 555 Conjugated Goat Anti-Rabbit IgG (H+L) Goat pAb	HA1135	Rabbit	IF-Cell,FC
	iFluor™ 750 Conjugated Goat Anti-Rabbit IgG (H+L) Goat pAb	HA1134	Rabbit	IF-Cell,IF-Tissue,FC
	iFluor™ 594 Conjugated Goat anti-Rat IgG(H+L) Goat pAb	HA1136	Rat	IF-Tissue,IHC-Fr,IF-Cell,FC
	iFluor™ 488 Conjugated Goat anti-Rat IgG(H+L) Goat pAb	HA1133	Rat	IF-Cell,IF-Tissue,FC,IHC-Fr
	iFluor™ 488 Conjugated Rabbit anti-Goat IgG Rabbit pAb	HA1131	Goat	IF-Cell,IF-Tissue,FC
	iFluor™ 594 Conjugated Rabbit anti-Goat IgG Rabbit pAb	HA1132	Goat	IF-Cell,IF-Tissue,FC
	FITC Conjugated Rabbit Anti-Goat IgG polyclonal Antibody	HA1007	Goat	IF-Cell,IF-Tissue,FC
	TRITC Conjugated Rabbit Anti-Goat IgG(H+L) polyclonal Antibody	HA1116	Goat	IF-Cell,IF-Tissue,FC
	FITC Conjugated Rabbit Anti-Mouse IgG polyclonal Antibody	HA1009	Mouse	IF-Cell,IF-Tissue,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层