

HRP二抗-产品清单

高级验证	产品名	货号	反应种属	验证应用
	HRP Conjugated Alpaca anti-Rabbit IgG Fc, Recombinant VHH monoclonal Antibody	HA1031	Rabbit	IP,ELISA,IHC-P,WB
	HRP Conjugated Goat Anti-human IgG polyclonal Antibody	HA1018	Human	WB,ELISA
	HRP Conjugated Goat anti-Mouse IgG polyclonal Antibody	HA1006	Mouse	WB,ELISA,IHC-P
	HRP Conjugated Goat anti-Mouse IgG UltraPolymer Antibody	HA1120	Mouse	IHC-P
	HRP Conjugated Goat anti-Rabbit IgG polyclonal Antibody	HA1001	Rabbit	WB,ELISA,IHC-P
	HRP Conjugated Goat anti-Rabbit IgG UltraPolymer Antibody	HA1119	Rabbit	IHC-P
	HRP Conjugated Goat anti-Rat IgG(H+L) polyclonal Antibody	HA1023	Rat	WB,ELISA
	HRP Conjugated Mouse Anti-Human IgG heavy chain [54-76-R] monoclonal Antibody	HA601284	Human	WB
	HRP Conjugated Rabbit anti-Chicken IgY polyclonal Antibody	HA1026	Chicken	ELISA,WB
	HRP Conjugated Rabbit anti-Goat IgG polyclonal Antibody	HA1005	Goat	WB,ELISA
	HRP Conjugated Rabbit anti-Guinea pig IgG polyclonal Antibody	HA1021	Guinea pig	WB,ELISA
	HRP Conjugated Rabbit anti-Human IgG polyclonal Antibody	HA1024	Human	WB,ELISA
	HRP Conjugated Rabbit anti-Mouse IgG polyclonal Antibody	HA1008	Mouse	WB,ELISA
	HRP Conjugated Rabbit anti-Rat IgG(H+L) polyclonal Antibody	HA1022	Rat	WB,ELISA