

2025年1月-新上市产品清单

高级验证	产品名	货号	反应种属	验证应用
	ACOX1 Rabbit pAb	ER1919-52	Human, Mouse, Rat	WB
Relative expression	Apolipoprotein E Recombinant Antibody [PSH10-52] - Rat IgG1 (Chimeric)	HA601391	Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Bik Recombinant Rabbit mAb [PSH11-51]	HA723343	Human	WB,IF-Cell,IP
	C9orf72 Recombinant Rabbit mAb [PSH11-11]	HA723290	Human, Mouse, Rat	WB
Cell treatment, Relative expression	CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-27]	HA723309	Human	WB,IF-Cell,FC
Cell treatment, Relative expression	CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-28]	HA723310	Human, Mouse	WB,IF-Cell,IHC-P
	CEPT1 (N-term) Rabbit pAb	ER1919-53	Human, Mouse, Rat	WB,IHC
	CRLS1 Rabbit pAb	ER1919-51	Human	WB,IHC
Relative expression	Doublecortin Recombinant Antibody [JJ0959] - Rat IgG1 (Chimeric)	HA601398	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
Cell treatment	FAM38A / PIEZO1 Recombinant Rabbit mAb [PSH11-49]	HA723341	Mouse	WB,IF-Cell
Cell treatment	FAM38B / PIEZO2 Recombinant Rabbit mAb [PSH11-50]	HA723342	Mouse	WB,IF-Cell
	Gli2 Recombinant Rabbit mAb [PSH11-23]	HA723305	Human, Mouse, Rat	WB
Relative expression	GM2A Rabbit pAb	HA500560	Human, Mouse, Rat	WB,IHC-P,IP
	HIP2 / LIG Recombinant Rabbit mAb [PSH11-24]	HA723306	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	Human CKMB Recombinant Rabbit mAb [PSH11-56] - BSA and Azide free (Detector)	HA723349	Human, Rat, Mouse, Monkey	ELISA(Det)
	Biotin Conjugated Human CKMB Recombinant Rabbit mAb [PSH11-56] - Detector	HA723350B	Human, Rat, Mouse, Monkey	ELISA(Det),ELISA
	Human CKMB Recombinant Rabbit mAb [PSH11-55] - BSA and Azide free (Capture)	HA723348	Human, Rat, Mouse, Monkey	ELISA(Cap)
	Biotin Conjugated Human CTGF Recombinant Rabbit mAb [PSH11-16] - Detector	HA723298B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CTGF Recombinant Rabbit mAb [PSH11-16] - BSA and Azide free (Detector)	HA723297	Human	ELISA(Det)
	Human CTGF Recombinant Rabbit mAb [PSH11-15] - BSA and Azide free (Capture)	HA723296	Human	ELISA(Cap)
	Human CXCL6/GCP-2 Recombinant Rabbit mAb [PSH11-30] - BSA and Azide free (Capture)	HA723314	Human	ELISA(Cap)
	Human CXCL6/GCP-2 Recombinant Rabbit mAb [PSH11-31] - BSA and Azide free (Detector)	HA723315	Human	ELISA(Det)
	Biotin Conjugated Human CXCL6/GCP-2 Recombinant Rabbit mAb [PSH11-31] - Detector	HA723316B	Human	ELISA(Det),ELISA
	Human FOLR1 Recombinant Rabbit mAb [PSH11-60] - BSA and Azide free (Detector)	HA723355	Human	ELISA(Det)
	Human FOLR1 Recombinant Rabbit mAb [PSH11-59] - BSA and Azide free (Capture)	HA723354	Human	ELISA(Cap)
	Biotin Conjugated Human FOLR1 Recombinant Rabbit mAb [PSH11-60] - Detector	HA723356B	Human	ELISA(Det),ELISA
	Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-33] - BSA and Azide free (Detector)	HA723318	Human	ELISA(Det)
	Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-32] - BSA and Azide free (Capture)	HA723317	Human	ELISA(Cap)
	Biotin Conjugated Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-33] - Detector	HA723319B	Human	ELISA(Det),ELISA
	Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-34] - BSA and Azide free (Capture)	HA723320	Human	ELISA(Cap)
	Biotin Conjugated Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-35] - Detector	HA723322B	Human	ELISA(Det),ELISA
	Human ICOSLG/B7-H2/CD275 Recombinant Rabbit mAb [PSH11-35] - BSA and Azide free (Detector)	HA723321	Human	ELISA(Det)
	Biotin Conjugated Human IL3RA/CD123 Recombinant Rabbit mAb [PSH11-58] - Detector	HA723353B	Human	ELISA(Det),ELISA
	Human IL3RA/CD123 Recombinant Rabbit mAb [PSH11-57] - BSA and Azide free (Capture)	HA723351	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL3RA/CD123 Recombinant Rabbit mAb [PSH11-58] - BSA and Azide free (Detector)	HA723352	Human	ELISA(Det)
	Human IL-4R Recombinant Rabbit mAb [PSH11-45] - BSA and Azide free (Capture)	HA723336	Human	ELISA(Cap)
	Human IL-4R Recombinant Rabbit mAb [PSH11-46] - BSA and Azide free (Detector)	HA723337	Human	ELISA(Det)
	Biotin Conjugated Human IL-4R Recombinant Rabbit mAb [PSH11-46] - Detector	HA723338B	Human	ELISA(Det),ELISA
	Human M-CSF Recombinant Rabbit mAb [PSH11-36] - BSA and Azide free (Capture)	HA723323	Human	ELISA(Cap)
	Human M-CSF Recombinant Rabbit mAb [PSH11-37] - BSA and Azide free (Capture)	HA723324	Human	ELISA(Cap)
	Biotin Conjugated Human M-CSF Recombinant Rabbit mAb [PSH11-38] - Detector	HA723326B	Human	ELISA(Det),ELISA
	Human M-CSF Recombinant Rabbit mAb [PSH11-38] - BSA and Azide free (Detector)	HA723325	Human	ELISA(Det)
	Human MICA Recombinant Rabbit mAb [PSH11-40] - BSA and Azide free (Detector)	HA723328	Human	ELISA(Det)
	Human MICA Recombinant Rabbit mAb [PSH11-39] - BSA and Azide free (Capture)	HA723327	Human	ELISA(Cap)
	Biotin Conjugated Human MICA Recombinant Rabbit mAb [PSH11-40] - Detector	HA723329B	Human	ELISA(Det),ELISA
	Human MMP13 Recombinant Rabbit mAb [PSH11-13] - BSA and Azide free (Capture)	HA723293	Human	ELISA(Cap)
	Biotin Conjugated Human MMP13 Recombinant Rabbit mAb [PSH11-14] - Detector	HA723295B	Human	ELISA(Det),ELISA
	Human MMP13 Recombinant Rabbit mAb [PSH11-14] - BSA and Azide free (Detector)	HA723294	Human	ELISA(Det)
	Biotin Conjugated Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-44] - Detector	HA723335B	Human, Mouse, Rat	ELISA(Det),ELISA
	Biotin Conjugated Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-42] - Detector	HA723332B	Human, Mouse, Rat	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-44] - BSA and Azide free (Detector)	HA723334	Human, Mouse, Rat	ELISA(Det)
	Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-42] - BSA and Azide free (Detector)	HA723331	Human, Mouse, Rat	ELISA(Det)
	Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-41] - BSA and Azide free (Capture)	HA723330	Human, Mouse, Rat	ELISA(Cap)
	Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-43] - BSA and Azide free (Capture)	HA723333	Human, Mouse, Rat	ELISA(Cap)
	IFT88 Recombinant Rabbit mAb [PSH11-61]	HA723357	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell
	Integrin beta 5 Recombinant Rabbit mAb [PSH11-53]	HA723345	Human	IF-Cell,FC
Relative expression	Lipoprotein a Recombinant Rabbit mAb [PSH11-29]	HA723311	Human	WB,IHC-P,IF-Tissue
Relative expression	MMP-2 Recombinant Rabbit mAb [PSH11-26]	HA723308	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA723347	Mouse	FC
	MYOF Rabbit pAb	ER65966	Human, Mouse	WB
Relative expression	Nectin 2 Recombinant Rabbit mAb [PSH11-54]	HA723346	Human	WB,IF-Cell,FC
	NeuN Recombinant Antibody [SR45-07] - Rat IgG1 (Chimeric)	HA601397	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
Cell treatment	Nrf2 Recombinant Rabbit mAb [PSH11-20]	HA723302	Human, Mouse, Rat	WB,IF-Cell,ChIP
Cell treatment	Phospho-Bcl-2 (S70) Recombinant Rabbit mAb [PSH11-02]	HA723269	Human	WB,IF-Cell,FC
	Phospho-ELK1 (S383) Rabbit pAb	ER65967	Human, Mouse, Rat	WB,IHC
Cell treatment	Phospho-Lamin A + Lamin C (S22) Recombinant Rabbit mAb [PSH11-47]	HA723339	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-Lamin A + Lamin C (S22) Recombinant Rabbit mAb [PSH11-48]	HA723340	Human, Mouse, Rat	WB
Cell treatment	Phospho-MEK3 (S189)+MEK6 (S207) Recombinant Rabbit mAb [PSH11-17]	HA723299	Human, Mouse	WB
	PI 3-Kinase p110 γ Rabbit pAb	ER65965	Human, Mouse	WB,IHC,IF,ELISA
	PPAR gamma Mouse mAb [E-8]	RT1668	Human, Mouse, Rat	WB,IHC-P,IP
	Presenilin 1 / PS-1 Recombinant Rabbit mAb [PSH11-19]	HA723301	Human, Mouse, Rat	WB,IP
	Presenilin 2 / AD5 Recombinant Rabbit mAb [PSH11-18]	HA723300	Human, Mouse, Rat	WB,IF-Cell,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	SLC32A1 / VGAT Recombinant Antibody [PSH10-47] - Rabbit IgG (Chimeric)	HA723292	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
	SMOC2 Recombinant Rabbit mAb [PSH11-25]	HA723307	Human, Mouse, Rat	WB,IHC-P
	SOX2 Recombinant Antibody [PO00-28] - Rat IgG1 (Chimeric)	HA601396	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
	iFluor™ 594 Conjugated Strep-Tag II Recombinant Rabbit mAb [PSH05-13]	HA720228F	Species independent	IF-Cell
	FITC Conjugated Strep-Tag II Recombinant Rabbit mAb [PSH05-13]	HA720229F	Species independent	IF-Cell
Relative expression	TIMP2 Recombinant Rabbit mAb [PSH11-62]	HA723358	Human, Monkey	WB
Cell treatment	TMC1 Recombinant Rabbit mAb [PSH11-12]	HA723291	Human, Mouse	WB,IF-Cell
	iFluor™ 594 Conjugated V5 tag Recombinant Rabbit mAb [SY30-01]	HA720226F	Species independent	IF-Cell
	FITC Conjugated V5 tag Recombinant Rabbit mAb [SY30-01]	HA720227F	Species independent	IF-Cell
Relative expression	VGluT1 Recombinant Antibody [PSH10-53] - Rat IgG1 (Chimeric)	HA601392	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB,IF-Tissue
Relative expression	VGLUT2 Recombinant Antibody [PSH10-54] - Rat IgG1 (Chimeric)	HA601393	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,WB
Relative expression	ZAP70 Recombinant Rabbit mAb [PSH11-52]	HA723344	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层