

2025年2月—新上市产品清单

高级验证	产品名	货号	反应种属	验证应用
Relative expression	Acetylcholinesterase Recombinant Antibody [PSH10-55] - Rat IgG1 (Chimeric)	HA601400	Mouse, Rat	WB,IHC-Fr
Relative expression	AFP Recombinant Rabbit mAb [PSH11-75]	HA723371	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	Alx1 Rabbit pAb	HA500565	Human, Mouse, Rat	WB,IHC-P
Relative expression	Antithrombin III / ATIII Recombinant Rabbit mAb [PSH12-13]	HA723429	Human, Mouse, Rat	WB,IHC-P
	BAF170 Recombinant Rabbit mAb [PSH12-00]	HA723413	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,ChIP
Relative expression	beta 2 Microglobulin Recombinant Rabbit mAb [PSH12-14]	HA723430	Human	WB,IF-Cell,IHC-P,FC,IP
Relative expression	BLBP Recombinant Antibody [PSH10-68] - Rat IgG1 (Chimeric)	HA601411	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
	BRN3A Recombinant Rabbit mAb [PSH11-65]	HA723361	Mouse, Rat	IHC-Fr,IHC-P,WB
	BRN3A Recombinant Antibody [PSH11-65] - Rat IgG1 (Chimeric)	HA601412	Mouse	IHC-Fr,IHC-P
	CD11c Recombinant Antibody [SI19-06] - Rat IgG1 (Chimeric)	HA601418		
KD	CD14 Recombinant Antibody [SC69-02] - Rat IgG1 (Chimeric)	HA601420		
Relative expression	CD38 Recombinant Antibody [PD01-49] - Rat IgG1 (Chimeric)	HA601417		
Relative expression	CGRP Recombinant Antibody [PSH08-45] - Rat IgG1 (Chimeric)	HA601423	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	CHD8 Recombinant Rabbit mAb [PSH11-71]	HA723367	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,ChIP
Relative expression	Claudin 18.2 Recombinant Rabbit mAb [PSH11-63]	HA723359	Human, Mouse, Rat	IHC-P
	DNA helicase HEL308 Rabbit pAb	HA500562	Human, Mouse, Rat	WB
	Dopamine Transporter Recombinant Antibody [PSH10-49] - Rat IgG1 (Chimeric)	HA601404	Mouse, Rat, Human, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
	FITC Conjugated DYKDDDDK Tag (FLAG) Recombinant Mouse mAb [A2-A4-R]	HA600139F	Species independent	IF-Cell
	DYNC1H1 Rabbit pAb	ER1919-57	Human, Mouse, Rat	WB,IHC-P

Relative expression	EAAT2 Recombinant Antibody [PS01-62] - Rat IgG1 (Chimeric)	HA601403	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
KD	EEF2K Recombinant Rabbit mAb [PSH11-67]	HA723363	Human, Mouse, Rat	WB,IF-Cell,FC
	ELL3 Rabbit pAb	HA500564	Mouse, Rat	WB
	ELMO2 Rabbit pAb	ER1919-56	Human, Mouse	WB,IHC-P
	Exosome Essentials Antibody Sampler Kit	HAK21125	Human	WB
	Ferredoxin-2 Rabbit pAb	HA500563	Human, Mouse	WB
KD	GALC Recombinant Rabbit mAb [PSH04-32]	HA722115	Human	WB,IHC-P
	GLS2 Mouse mAb	EM2001-11	Human, Mouse, Rat	WB,IHC-P
	iFluor™ 594 Conjugated Goat anti-Rat IgG(H+L) Goat pAb	HA1136	Rat	IF-Tissue,IHC-Fr,IF-Cell,FC
	Histone H1.3 Rabbit pAb	ER1919-58	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	HLA-ABC (MHC I) Recombinant Rabbit mAb [PSH12-02]	HA723415	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	HLA Class I Recombinant Mouse mAb [PSH12-15]	HA601424	Human	IF-Cell,FC
Relative expression	Homer1 Recombinant Antibody [JE65-19] - Rat IgG1 (Chimeric)	HA601415	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
	Hsp40 Recombinant Rabbit mAb [PSH12-10]	HA723426	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	Human CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-78] - BSA and Azide free (Capture)	HA723380	Human	ELISA(Cap)
	Human CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-76] - BSA and Azide free (Capture)	HA723377	Human	ELISA(Cap)
	Biotin Conjugated Human CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-77] - Detector	HA723379B	Human	ELISA(Det)
	Human CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-77] - BSA and Azide free (Detector)	HA723378	Human	ELISA(Det)
	Human CD47, C-His Tag (ECD) Protein	HA211050	Human	
	Biotin Conjugated Human CNTF Recombinant Rabbit mAb [PSH11-80] - Detector	HA723383B	Human, Rat	ELISA(Det),ELISA
	Human CNTF Recombinant Rabbit mAb [PSH11-80] - BSA and Azide free (Detector)	HA723382	Human, Rat	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CNTF Recombinant Rabbit mAb [PSH11-82] - BSA and Azide free (Detector)	HA723385	Human, Rat	ELISA(Det)
	Human CNTF Recombinant Rabbit mAb [PSH11-81] - BSA and Azide free (Capture)	HA723384	Human, Rat	ELISA(Cap)
	Human CNTF Recombinant Rabbit mAb [PSH11-79] - BSA and Azide free (Capture)	HA723381	Human, Rat	ELISA(Cap)
	Biotin Conjugated Human CNTF Recombinant Rabbit mAb [PSH11-82] - Detector	HA723386B	Human, Rat	ELISA(Det),ELISA
	HRP Conjugated Human CSPG4 Recombinant Rabbit mAb [PSH10-79] - Detector	HA723376H	Human	ELISA(Det)
	HRP Conjugated Human CTGF Recombinant Rabbit mAb [PSH11-16] - Detector	HA723375H	Human	ELISA(Det)
	Human CTGF/CCN2, C-Flag Tag Protein	HA211051	Human	
	Human CXCL13 / BCA-1, Tag Free Protein	HA211065	Human	
	Human Fas/TNFRSF6/CD95, C-His Tag (ECD) Protein	HA211064	Human	
	Human FOLR1, Tag Free Protein	HA211055	Human	
	Human Fractalkine/CX3CL1, C-His Tag (ECD) Protein	HA211057	Human	
	Biotin Conjugated Human IGFBP3 Recombinant Rabbit mAb [PSH12-17] - Detector	HA723434B	Human	ELISA(Det),ELISA
	Human IGFBP3 Recombinant Rabbit mAb [PSH12-17] - BSA and Azide free (Detector)	HA723433	Human	ELISA(Det)
	Human IGFBP3 Recombinant Rabbit mAb [PSH12-16] - BSA and Azide free (Capture)	HA723432	Human	ELISA(Cap)
	Human Leptin Recombinant Rabbit mAb [PSH12-07] - BSA and Azide free (Detector)	HA723422	Human	ELISA(Det)
	Biotin Conjugated Human Leptin Recombinant Rabbit mAb [PSH12-07] - Detector	HA723423B	Human	ELISA(Det),ELISA
	Human Leptin Recombinant Rabbit mAb [PSH12-06] - BSA and Azide free (Capture)	HA723421	Human	ELISA(Cap)
	Human NG2/MCSP/CSPG4, Tag Free (ECD) Protein	HA211047	Human	
	HRP Conjugated Human PD-L2 Recombinant Rabbit mAb [PSH09-82] - Detector	HA723373H	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human Pentraxin 3 / PTX3 Recombinant Rabbit mAb [PSH11-85] - Detector	HA723391B	Human	ELISA(Det),ELISA
	Human Pentraxin 3 / PTX3 Recombinant Rabbit mAb [PSH11-84] - BSA and Azide free (Detector)	HA723388	Human	ELISA(Det)
	Human Pentraxin 3 / PTX3 Recombinant Rabbit mAb [PSH11-85] - BSA and Azide free (Detector)	HA723390	Human	ELISA(Det)
	Biotin Conjugated Human Pentraxin 3 / PTX3 Recombinant Rabbit mAb [PSH11-84] - Detector	HA723389B	Human	ELISA(Det),ELISA
	Human Pentraxin 3 / PTX3 Recombinant Rabbit mAb [PSH11-83] - BSA and Azide free (Capture)	HA723387	Human	ELISA(Cap)
	Human Siglec 7 Recombinant Rabbit mAb [PSH11-88] - BSA and Azide free (Detector)	HA723395	Human	ELISA(Det)
	Biotin Conjugated Human Siglec 7 Recombinant Rabbit mAb [PSH11-87] - Detector	HA723394B	Human	ELISA(Det),ELISA
	Human Siglec 7 Recombinant Rabbit mAb [PSH11-87] - BSA and Azide free (Capture/Detector)	HA723393	Human	ELISA(Cap),ELISA(Det)
	Biotin Conjugated Human Siglec 7 Recombinant Rabbit mAb [PSH11-88] - Detector	HA723396B	Human	ELISA(Det),ELISA
	Human Siglec 7 Recombinant Rabbit mAb [PSH11-86] - BSA and Azide free (Capture)	HA723392	Human	ELISA(Cap)
	Human Siglec 9 Recombinant Rabbit mAb [PSH11-89] - BSA and Azide free (Capture)	HA723397	Human	ELISA(Cap)
	Biotin Conjugated Human Siglec 9 Recombinant Rabbit mAb [PSH11-90] - Detector	HA723399B	Human	ELISA(Det),ELISA
	Human Siglec 9 Recombinant Rabbit mAb [PSH11-90] - BSA and Azide free (Detector)	HA723398	Human	ELISA(Det)
	Biotin Conjugated Human Siglec 9 Recombinant Rabbit mAb [PSH11-91] - Detector	HA723401B	Human	ELISA(Det),ELISA
	Human Siglec 9 Recombinant Rabbit mAb [PSH11-91] - BSA and Azide free (Detector)	HA723400	Human	ELISA(Det)
	HRP Conjugated Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - Detector	HA723374H	Human	ELISA(Det)
	Human TNF Receptor II Recombinant Rabbit mAb [PSH11-92] - BSA and Azide free (Capture)	HA723402	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human TNF Receptor II Recombinant Rabbit mAb [PSH11-94] - BSA and Azide free (Capture)	HA723405	Human	ELISA(Cap)
	Human TNF Receptor II Recombinant Rabbit mAb [PSH11-93] - BSA and Azide free (Detector)	HA723403	Human	ELISA(Det)
	Biotin Conjugated Human TNF Receptor II Recombinant Rabbit mAb [PSH11-93] - Detector	HA723404B	Human	ELISA(Det),ELISA
	HRP Conjugated Human TREM2 Recombinant Rabbit mAb [PSH09-83] - Detector	HA723372	Human	ELISA(Det)
	Human TREM2, Tag Free Protein	HA211052	Human	
	Human ULBP1 Recombinant Rabbit mAb [PSH11-96] - BSA and Azide free (Detector)	HA723407	Human	ELISA(Det)
	Biotin Conjugated Human ULBP1 Recombinant Rabbit mAb [PSH11-96] - Detector	HA723408B	Human	ELISA(Det),ELISA
	Human ULBP1 Recombinant Rabbit mAb [PSH11-95] - BSA and Azide free (Capture)	HA723406	Human	ELISA(Cap)
	Human ULBP1 Recombinant Rabbit mAb [PSH11-97] - BSA and Azide free (Detector)	HA723409	Human	ELISA(Det)
	Biotin Conjugated Human ULBP1 Recombinant Rabbit mAb [PSH11-97] - Detector	HA723410B	Human	ELISA(Det),ELISA
Relative expression	IGF1 Receptor beta Recombinant Rabbit mAb [PSH11-72]	HA723368	Human, Mouse, Rat	WB,IF-Cell
	LARP1 Recombinant Rabbit mAb [PSH12-09]	HA723425	Human, Monkey	WB,IF-Cell
Relative expression	MAG Recombinant Antibody [PSH02-41] - Rat IgG1 (Chimeric)	HA601414	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB,IF-Tissue
	Mast Cell Tryptase Recombinant Antibody [SC68-07] - Rat IgG1 (Chimeric)	HA601419		
	MBPHL Rabbit pAb	ER65968	Human, Mouse, Rat	WB
	MDH1B Rabbit pAb	HA500561	Mouse	WB
	MEK4 / MKK4 Recombinant Rabbit mAb [PSH11-99]	HA723412	Human, Mouse, Rat, Monkey	WB,IHC-P,FC
	Mouse CCL20/MIP-3, Tag Free Protein	HA211071	Mouse	
	Mouse CCL5 Recombinant Rat mAb	HA601395	Human, Mouse	ELISA
	Mouse CD30/TNFRSF8, C-His Tag Protein	HA211062	Mouse	
	Mouse CXCL16, C-His Tag Protein	HA211059	Mouse	
	Mouse CXCL1, Tag Free (HEK293) Protein	HA211049	Mouse	

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse Erythropoietin/EPO, C-His Tag Protein	HA211068	Mouse	
	Mouse ICAM-1/CD54, C-His Tag (ECD) Protein	HA211066	Mouse	
	Biotin Conjugated Mouse IL-18 Recombinant Rabbit mAb [PSH12-04] - Detector	HA723418B	Mouse	ELISA(Det),ELISA
	Mouse IL-18 Recombinant Rabbit mAb [PSH12-04] - BSA and Azide free (Detector)	HA723417	Mouse	ELISA(Det)
	Mouse IL-18 Recombinant Rabbit mAb [PSH12-03] - BSA and Azide free (Capture)	HA723416	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-21 Recombinant Rabbit mAb [PSH04-75] - Detector	HA723420B	Mouse	ELISA(Det)
	Mouse IL-23, C-His Tag Protein	HA211060	Mouse	
	Mouse MDC/CCL22, Tag Free Protein	HA211072	Mouse	
	Mouse MerTK Rabbit pAb	ER1919-54	Mouse, Rat	WB
	Mouse SCF, C-His Tag Protein	HA211067	Mouse	
	Mouse TNF-alpha, Tag Free Protein	HA211054	Mouse	
	Mouse TREM2, Tag Free (ECD) Protein	HA211061	Mouse	
	Mouse VEGFR1/FLT1, C-His Tag Protein	HA211053	Mouse	
Relative expression	Myelin PLP Recombinant Antibody [PSH08-21] - Rat IgG1 (Chimeric)	HA601407	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	NAK / TBK1 Recombinant Rabbit mAb [PSH11-69]	HA723365	Human, Mouse, Rat	WB,IF-Cell,FC,IP
	Nestin Recombinant Antibody [PSH07-89] - Rat IgG1 (Chimeric)	HA601406	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell
Relative expression	NF-L Recombinant Antibody [PS02-10] - Rat IgG1 (Chimeric)	HA601409	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
	NF-M Recombinant Antibody [JM11-20] - Rat IgG1 (Chimeric)	HA601408	Human, Mouse, Rat, Cynomolgus monkey	IHC-Fr,IHC-P,WB
	NSFL1C Rabbit pAb	ER1919-59	Human, Mouse, Rat	WB
	Otx2 Recombinant Antibody [PSH01-61] - Rat IgG1 (Chimeric)	HA601399	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
Relative expression	Pericyte Antibody Sampler Kit	HAK21124	Human	WB
	Peripherin Recombinant Antibody [PSH10-70] - Rat IgG1 (Chimeric)	HA601405	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB,IF-Cell
	PGAM5 Recombinant Rabbit mAb [PSH11-68]	HA723364	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Src (S17) Recombinant Rabbit mAb [PSH11-03]	HA723270	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP
	Phospho-STAT4 (Y693) Rabbit pAb	ER1919-55	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-Syk (Y352) Recombinant Rabbit mAb [PSH11-66]	HA723362	Human	WB,IF-Cell,FC
Relative expression	PKM2 Recombinant Rabbit mAb [PSH11-74]	HA723370	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP
	POLR2A Rabbit pAb	HA500566	Human, Mouse, Rat	WB
Relative expression	PSD95 Recombinant Antibody [SR38-09] - Rat IgG1 (Chimeric)	HA601416	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,WB
	Rat CCL3/MIP-1 alpha, Tag Free Protein	HA211056	Rat	
	Rat ICAM-1/CD54, C-His Tag (ECD) Protein	HA211063	Rat	
	Rat IL-1 alpha/IL1A, Tag Free Protein	HA211058	Rat	
Relative expression	SLAM / CD150 Recombinant Rabbit mAb [PSH12-08]	HA723424	Human	WB,IP
	SOX10 Recombinant Antibody [PDH0-03] - Rat IgG1 (Chimeric)	HA601413	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Tissue
Cell treatment	Stanniocalcin 2 / STC-2 Recombinant Rabbit mAb [PSH12-01]	HA723414	Human	WB,IP
	Survivin Recombinant Rabbit mAb [PSH11-73]	HA723369	Human, Mouse	WB,IHC-P
Relative expression	SV2A Recombinant Antibody [PSH10-69] - Rat IgG1 (Chimeric)	HA601421	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
Relative expression	Synaptopodin Recombinant Antibody [PSH10-45] - Rat IgG1 (Chimeric)	HA601402	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	TBR1 Recombinant Antibody [JF10-00] - Rat IgG1 (Chimeric)	HA601422	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
	TFEB Recombinant Rabbit mAb [PSH12-05]	HA723419	Human, Mouse, Rat	WB
Relative expression	TNF Receptor II Recombinant Rabbit mAb [PSH12-12]	HA723428	Human	FC
	Twist1 Recombinant Rabbit mAb [PSH11-98]	HA723411	Human, Mouse, Rat	WB
	VACht / SLC18A3 Recombinant Antibody [PSH10-50] - Rat IgG1 (Chimeric)	HA601410	Mouse, Species independent, Pig	IHC-Fr
	Vasopressin Recombinant Rabbit mAb [PSH11-64]	HA723360	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Tissue
	VIP Recombinant Antibody [PSH10-56] - Rat IgG1 (Chimeric)	HA601401	Human, Mouse, Rat, Cynomolgus monkey	IHC-Fr,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

