

炎症小体信号通路-产品清单

高级验证	产品名	货号	反应种属	验证应用
KD	14-3-3 epsilon Recombinant Rabbit mAb [JJ08-40]	ET1701-1	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P
	53BP1 Recombinant Rabbit mAb [JA66-11] - BSA and Azide free	HA750400	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	53BP1 Recombinant Rabbit mAb [JA66-11]	ET1704-05	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	AIM2 Rabbit pAb	HA500270	Human	WB,IF-Cell,FC
	Angiotensinogen Recombinant Rabbit mAb [SD201-02]	ET1612-48	Human	WB,IF-Cell,IF-Tissue,IP
	AP-1/Jun D Rabbit pAb	ER62677	Human, Mouse, Rat	WB,IHC-P
	ATG12 Rabbit pAb	ER1802-82	Human	WB,IHC-P
	ATG12 Recombinant Rabbit mAb [JE30-61]	HA721504	Human	WB,IF-Cell,FC
	ATG12 Rabbit pAb	R1404-1	Human, Mouse, Rat	WB,IHC-P
	ATG16L1 Recombinant Rabbit mAb [JU30-62]	ET7106-65	Human, Mouse, Rat	WB,IHC-P
	ATG16L1 Mouse mAb [A11-H9]	EM1708-14	Human, Mouse, Rat	WB,FC
	ATG4B Recombinant Rabbit mAb [JE59-97]	HA721334	Human, Mouse, Rat, Monkey	WB
KD	ATG4B Recombinant Rabbit mAb [PSH02-53] - BSA and Azide free	HA750816	Human, Mouse, Rat, Monkey	WB
KD	ATG4B Recombinant Rabbit mAb [PSH02-53]	HA721829	Human, Mouse, Rat, Monkey	WB
	ATG4B Rabbit pAb	ER62759	Human, Mouse, Rat, Monkey	WB,IHC-P
KO	ATG5 Recombinant Rabbit mAb [SN73-07] - BSA and Azide free	HA750251	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO	ATG5 Recombinant Rabbit mAb [SN73-07]	ET1611-38	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	AVPR1A Rabbit pAb	HA500002	Human, Mouse	WB
KD	Bcl-2 Recombinant Rabbit mAb [SZ10-03]	ET1603-11	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	FITC Conjugated Bcl-2 Recombinant Rabbit mAb [JF104-8]	HA720210F	Human, Mouse	FC
KD	Bcl-2 Recombinant Rabbit mAb [JF104-8] - BSA and Azide free	HA750350	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Bcl-2 Mouse mAb [F7-1-A1]	M1206-4	Human, Mouse, Rat	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Bcl-2 Mouse mAb [9F3]	EM1701-83	Human	WB,IHC-P,FC
	Bcl-2 Recombinant Mouse mAb [9F3-R]	HA601256	Human	WB,IHC-P,FC
	Bcl-2 Rabbit pAb	ER1802-97	Human, Mouse	WB,IF-Cell,FC
	Bcl-2 Rabbit pAb	ER0602	Human	WB,IF-Cell,IHC-P,FC
KD	Bcl-2 Recombinant Rabbit mAb [PD01-41]	HA721235	Human	WB,IHC-P,IF-Cell,FC
	Bcl-2 Recombinant Mouse mAb [9F3-R] - BSA and Azide free	HA610131	Human	WB,IHC-P,FC
	Bcl-2 Rabbit pAb	ER1706-47	Human	WB,IF-Cell,IHC-P,FC
	Bcl-2 Mouse mAb [9F1]	EM1701-82	Human	WB,IHC-P
KD	Bcl-2 Recombinant Rabbit mAb [SZ10-03] - BSA and Azide free	HA750057	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Bcl-2 Recombinant Rabbit mAb [PSH09-49]	HA723111	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	Bcl-2 Recombinant Rabbit mAb [PSH09-49] - BSA and Azide free	HA751296	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
KD	Bcl-2 Recombinant Rabbit mAb [JF104-8]	ET1702-53	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	BCL2L1 Rabbit pAb	R1308-16	Human, Mouse, Rat	WB,IHC-P
	Bcl-x Recombinant Rabbit mAb	RE6009	Human, Mouse, Rat	WB,ICC,IHC-P,FC,IP
	Bcl-XL Recombinant Rabbit mAb [SZ3-03] - BSA and Azide free	HA750063	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	Bcl-XL Recombinant Rabbit mAb [SZ3-03]	ET1603-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	BIRC2 Rabbit pAb	ER62503	Human, Mouse	WB
	BRCC36 Recombinant Rabbit mAb [PSH05-10] - BSA and Azide free	HA750983	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	BRCC36 Recombinant Rabbit mAb [PSH05-10]	HA722226	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	CARD9 Recombinant Rabbit mAb [PSH0-41]	HA721343	Human, Mouse, Rat	WB,IF-Cell,FC
	CARD9 Recombinant Rabbit mAb [PSH0-41] - BSA and Azide free	HA750619	Human, Mouse, Rat	WB,IF-Cell,FC
	Caspase-1 Rabbit pAb	R1510-22	Human	IHC-P,FC
	Caspase-1 Recombinant Rabbit mAb [SU40-07] - BSA and Azide free	HA750161	Human, Mouse, Rat	WB,FC,IF-Cell
	Caspase-1 Rabbit pAb	R1510-23	Human	IF-Cell,IHC-P,FC
	Caspase-1 Recombinant Rabbit mAb [SU40-07]	ET1608-69	Human, Mouse, Rat	WB,FC,IF-Cell
Relative expression	Caspase-12 Recombinant Rabbit mAb [PSH07-49] - BSA and Azide free	HA751170	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Caspase-12 Rabbit pAb	ER62907	Human	WB,IHC-P
Relative expression	Caspase-12 Recombinant Rabbit mAb [PSH07-49]	HA722837	Human, Mouse, Rat	WB,IHC-P
	Caspase-12 Rabbit pAb	HA500144	Human, Mouse, Rat	WB,IHC-P
Relative expression	Caspase-1 + p10 + p12 Recombinant Rabbit mAb [JE56-35]	HA722222	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Caspase-1 p20 Rabbit pAb	ER1905-47	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	Caspase-5 Recombinant Rabbit mAb [SD203-2]	ET1612-29	Human	WB,IF-Tissue,IHC-P,IP
	Caspase-8 Mouse mAb [6G1]	RT1099	Human	WB,IP,IF,IHC-P
Cell treatment, KD	Caspase-8 Recombinant Rabbit mAb [PSH04-77] - BSA and Azide free	HA750971	Human	WB,IF-Cell,IHC-P,FC,IP
Cell treatment, KD	Caspase-8 Recombinant Rabbit mAb [PSH04-77]	HA722181	Human	WB,IF-Cell,IHC-P,FC,IP
Relative expression, KD	Caspase-8 Recombinant Rabbit mAb [SD08-06]	ET1612-70	Human	WB,IF-Tissue,IHC-P
Relative expression, KD	Caspase-8 Recombinant Rabbit mAb [SZ01-08]	ET1603-16	Human	WB,IF-Tissue,IHC-P
Cell treatment	Caspase-8 (p18) Recombinant Rabbit mAb [PSH05-83]	HA722482	Mouse	WB,IF-Cell
Cell treatment	Caspase-8 (p18) Recombinant Rabbit mAb [PSH05-83] - BSA and Azide free	HA751031	Mouse	WB,IF-Cell
	CaSR Rabbit pAb	ER62910	Human, Mouse, Rat	WB,IHC-P
Relative expression	CaSR Recombinant Rabbit mAb [PSH03-47] - BSA and Azide free	HA750881	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	CaSR Recombinant Rabbit mAb [PSH03-47]	HA722006	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	Cathepsin B Recombinant Rabbit mAb [JA11-02]	ET1704-03	Human	WB,IF-Cell,IF-Tissue,IHC-P
	Cathepsin B Mouse mAb [J11-A11]	M1506-1	Human, Zebrafish, Mouse, Rat	WB,IF-Cell
	Cathepsin B Rabbit pAb	ER1905-48	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
Relative expression	Cathepsin B Rabbit pAb	R1407-1	Human	WB
	Cathepsin S Rabbit pAb	HA500422	Human, Mouse, Rat	WB,IF-Cell
	CCL2 Mouse mAb [E11-D8]	EM1710-22	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
	CCL2 / MCP-1 Rabbit pAb	HA500267	Mouse	WB,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-44] - BSA and Azide free	HA751533	Mouse	WB,IF-Cell,FC
	CCL2 / MCP-1 Rabbit pAb	HA500042	Human	WB,IHC-P,IF-Cell
	CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-44]	HA723691	Mouse	WB,IF-Cell,FC
Cell treatment	CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-16]	HA723979	Mouse, Rat	WB,IF-Cell,FC
Relative expression	CFTR Recombinant Rabbit mAb [PSH09-99] - BSA and Azide free	HA751330	Human	WB,IP
Relative expression	CFTR Recombinant Rabbit mAb [PSH09-99]	HA723168	Human	WB,IP
KD	cIAP1 Recombinant Rabbit mAb [PSH02-10]	HA721769	Human	WB
KD	cIAP1 Recombinant Rabbit mAb [PSH06-71]	HA722743	Human, Mouse, Rat	WB
	cIAP1 Recombinant Rabbit mAb [PSH02-10] - BSA and Azide free (Capture)	HA723281	Human	ELISA(Cap)
	cIAP1 Recombinant Rabbit mAb [PSH06-71] - BSA and Azide free (Detector)	HA723282	Human, Mouse, Rat	ELISA(Det)
KD	cIAP1 Recombinant Rabbit mAb [PSH02-10] - BSA and Azide free	HA750780	Human	WB
	Biotin Conjugated cIAP1 Recombinant Rabbit mAb [PSH06-71] - Detector	HA723283B	Human, Mouse, Rat	ELISA(Det),ELISA
KD	cIAP1 Recombinant Rabbit mAb [PSH06-71] - BSA and Azide free	HA751100	Human, Mouse, Rat	WB
	cIAP1 Recombinant Rabbit mAb [PSH06-72]	HA722744	Human, Mouse, Rat	WB,IF-Cell,FC
	cIAP1 Recombinant Rabbit mAb [PSH06-72] - BSA and Azide free (Capture)	HA723280	Human, Mouse, Rat	ELISA(Cap)
	cIAP1 Recombinant Rabbit mAb [PSH06-72] - BSA and Azide free	HA751101	Human, Mouse, Rat	WB,IF-Cell,FC
	cIAP1/2 Rabbit pAb	HA500535	Human, Mouse	WB,IF-Cell,IHC-P,FC
	cIAP2 Recombinant Rabbit mAb [PSH06-92]	HA722767	Human	WB
	cIAP2 Recombinant Rabbit mAb [PSH06-92] - BSA and Azide free	HA751116	Human	WB
	c-Jun Recombinant Rabbit mAb [SY0298]	ET1608-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	c-Jun Recombinant Rabbit mAb [SY0298] - BSA and Azide free	HA750143	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	Cleaved-Caspase-4/5 p20 (D270/D311) Rabbit pAb	ER60012	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CLP Recombinant Rabbit mAb [PSH01-72]	HA721724	Human, Mouse, Rat	IHC,IF-Tissue
Relative expression	CLP Recombinant Rabbit mAb [PSH01-73] - BSA and Azide free	HA750748	Human, Mouse, Rat	IHC,IF-Tissue
Relative expression	CLP Recombinant Rabbit mAb [PSH01-72] - BSA and Azide free	HA750747	Human, Mouse, Rat	IHC,IF-Tissue
Relative expression	CLP Recombinant Rabbit mAb [PSH01-73]	HA721725	Human, Mouse, Rat	IHC,IF-Tissue
Cell treatment	CXCL1 / GRO alpha Recombinant Rabbit mAb [PSH07-53] - BSA and Azide free	HA751174	Human	WB
Cell treatment	CXCL1 / GRO alpha Recombinant Rabbit mAb [PSH07-53]	HA722841	Human	WB
Cell treatment	CXCL2 Recombinant Rabbit mAb [PSH14-44] - BSA and Azide free	HA751523	Mouse	WB,IF-Cell,IHC-P,FC
Cell treatment	CXCL2 Recombinant Rabbit mAb [PSH14-44]	HA723640	Mouse	WB,IF-Cell,IHC-P,FC,IP
	DFNA5 / GSDME Recombinant Rabbit mAb [PSH07-48] - BSA and Azide free	HA751169	Human, Mouse, Rat	WB
Cell treatment, Relative expression	DFNA5 / GSDME Recombinant Rabbit mAb [PSH10-83]	HA723250	Human	WB
Cell treatment, Relative expression	DFNA5 / GSDME Recombinant Rabbit mAb [PSH10-84] - BSA and Azide free	HA751376	Human, Mouse	WB,IP
Cell treatment, Relative expression	DFNA5 / GSDME Recombinant Rabbit mAb [PSH09-51] - BSA and Azide free	HA751298	Human, Mouse	WB,IP
Cell treatment, Relative expression	DFNA5 / GSDME Recombinant Rabbit mAb [PSH10-83] - BSA and Azide free	HA751375	Human	WB
	DFNA5 / GSDME Recombinant Rabbit mAb [PSH07-48]	HA722836	Human, Mouse, Rat	WB
Cell treatment, Relative expression	DFNA5 / GSDME Recombinant Rabbit mAb [PSH09-51]	HA723113	Human, Mouse	WB,IP
Relative expression	DFNA5 / GSDME Mouse mAb [A1H2]	EM1901-49	Human	WB
	DFNA5 / GSDME Mouse mAb [A1H3]	EM1901-48	Human	WB,IF-Cell,IHC-P,FC
	DFNA5 / GSDME Rabbit pAb	ER1901-12	Human	WB,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment, Relative expression	DFNA5 / GSDME Recombinant Rabbit mAb [PSH10-84]	HA723251	Human, Mouse	WB,IP
	DNM1L Mouse mAb [C7-A4]	EM1708-45	Human, Mouse, Rat	WB,IF-Cell,FC
	DRP1 Rabbit pAb	HA500487	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
	ERK1 Mouse mAb [7-D6-E5]	M1505-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	ERK1 Recombinant Rabbit mAb [SP05-09]	ET1604-16	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	ERK1 Recombinant Rabbit mAb [SP00-74]	ET1604-32	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Erk1/2 Mouse mAb [2G1]	RT1484	Human, Mouse, Rat	WB,IP,IF
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03] - BSA and Azide free	HA750019	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
	ERK1/2 Rabbit pAb	ER131011	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03]	ET1601-29	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
	ERK2 Mouse mAb [5-D2]	M1407-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	ERK2 Recombinant Mouse mAb [5-D2-R]	HA601202	Human, Mouse, Rat	WB
KD	ERK2 Recombinant Mouse mAb [C3B12-R] - BSA and Azide free	HA610079	Human, Mouse, Rat	WB
KD	ERK2 Recombinant Mouse mAb [C3B12-R]	HA601201	Human, Mouse, Rat	WB
	ERK2 Mouse mAb [C3B12]	EM1901-55	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Mouse mAb [A10D2]	EM1901-54	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Recombinant Mouse mAb [5-D2-R] - BSA and Azide free	HA610080	Human, Mouse, Rat	WB
KD	ERK2 Rabbit pAb	ER131218	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	ERK2 Recombinant Rabbit mAb [SZ25-01]	ET1603-23	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
KD	ERK2 Mouse mAb [B10-8]	EM1901-53	Human, Mouse, Rat	WB,IHC-P,FC
	FADD Rabbit pAb	0810-1	Mouse	WB
	FADD Recombinant Rabbit mAb [PSH07-78] - BSA and Azide free	HA751190	Human, Mouse, Rat	WB,IHC-P
	FADD Recombinant Rabbit mAb [JE32-76]	HA722140	Mouse, Rat	WB
	FADD Recombinant Rabbit mAb [PSH07-78]	HA722894	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	FADD Recombinant Rabbit mAb [PSH01-13] - BSA and Azide free	HA750689	Human	WB,IF-Cell,FC
	FADD Rabbit pAb	ER1908-51	Human, Rat	WB,IHC-P,IF-Tissue
	FADD Recombinant Rabbit mAb [JE45-25]	HA721353	Human	WB,IHC-P,IF-Cell
	FADD Recombinant Rabbit mAb [PSH01-13]	HA721625	Human	WB,IF-Cell,FC
	GABARAP Recombinant Rabbit mAb [JM30-30]	ET1705-53	Human, Mouse, Rat	WB,IHC-P
	GABARAPL2 Rabbit pAb	ER1901-73	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Gasdermin D (N terminal) Recombinant Rabbit mAb [PSH10-87] - BSA and Azide free	HA751379	Human	WB,IP
Cell treatment	Gasdermin D (N terminal) Recombinant Rabbit mAb [PD00-18] - BSA and Azide free	HA750486	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Gasdermin D (N terminal) Recombinant Rabbit mAb [PD00-18]	HA721144	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Gasdermin D (N terminal) Rabbit pAb	ER1901-37	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
Cell treatment	Gasdermin D (N terminal) Recombinant Rabbit mAb [PSH10-87]	HA723254	Human	WB,IP
	GBP1 Recombinant Rabbit mAb [JE60-78]	HA720049	Human, Mouse, Rat	WB,IHC-P
	GBP2 Rabbit pAb	ER2001-12	Human	WB,IHC-P
	GRO alpha Rabbit pAb	ER63645	Human, Mouse, Rat	IHC-P
	GSDMD Mouse mAb [A8D8]	HA601047	Human	WB
	GSDMD Mouse mAb [A8D7]	HA601046	Human, Mouse	WB
Relative expression	HNP-1 Mouse mAb [2D2]	EM1701-52	Human, Mouse, Rat	WB,IHC-P
KD	HNRNPA0 Recombinant Rabbit mAb [PSH02-77]	HA721855	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
KD	HNRNPA0 Recombinant Rabbit mAb [PSH02-77] - BSA and Azide free	HA750832	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	Hsp90 Recombinant Rabbit mAb [JE05-06]	HA722689	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP
	Hsp90 alpha Rabbit pAb	R1510-29	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Hsp90 alpha Recombinant Rabbit mAb [SY14-06]	ET1605-57	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P
	Hsp90 alpha Recombinant Rabbit mAb [SY14-06] - BSA and Azide free	HA750091	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P
	Hsp90 alpha Mouse mAb [A10-A5-D2]	M1603-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	Hsp90 beta Recombinant Rabbit mAb [SY46-01]	ET1605-56	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	Hsp90 beta Recombinant Rabbit mAb [SY46-01] - BSA and Azide free	HA750090	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	HSP90 Beta Recombinant Mouse mAb [2-1-G3-R] - BSA and Azide free	HA610060	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	HSP90 Beta Mouse mAb [2-1-G3]	EM21103	Human, Mouse, Rat	WB,IHC-P
	HSP90 Beta Recombinant Mouse mAb [2-1-G3-R]	HA601176	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Biotin Conjugated Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-22] - Detector	HA725196B	Human	ELISA(Det)
	Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-22] - BSA and Azide free (Detector)	HA725196	Human	ELISA(Det)
	Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-21] - BSA and Azide free (Capture)	HA725195	Human	ELISA(Cap)
	Human CCL5 Recombinant Rabbit mAb [PSH13-16] - BSA and Azide free (Capture)	HA725055	Human	ELISA(Cap)
	HRP Conjugated Human CCL5 Recombinant Rabbit mAb [PSH13-17] - Detector	HA725057H	Human	ELISA(Det),ELISA
	Human CCL5 Recombinant Rabbit mAb [PSH13-17] - BSA and Azide free (Detector)	HA725056	Human	ELISA(Det)
	Biotin Conjugated Human CCL5 Recombinant Rabbit mAb [PSH13-17] - Detector	HA725057B	Human	ELISA(Det),ELISA
	Human cIAP2 Recombinant Rabbit mAb [PSH09-39] - BSA and Azide free (Capture/Detector)	HA723099	Human	ELISA(Det),ELISA(Cap)
	Human cIAP2 Recombinant Rabbit mAb [PSH09-38] - BSA and Azide free (Capture)	HA723098	Human	ELISA(Cap)
	Biotin Conjugated Human cIAP2 Recombinant Rabbit mAb [PSH09-39] - Detector	HA723101B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human cIAP2 Recombinant Rabbit mAb [PSH09-40] - Detector	HA723102B	Human	ELISA(Det),ELISA
	Human cIAP2 Recombinant Rabbit mAb [PSH09-40] - BSA and Azide free (Detector)	HA723100	Human	ELISA(Det)
	Human CXCL1/GRO alpha Recombinant Rabbit mAb [PSH05-94] - BSA and Azide free (Capture)	HA722493	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CXCL1/GRO alpha Recombinant Rabbit mAb [PSH05-95] - BSA and Azide free (Detector)	HA722494	Human	ELISA(Det)
	Biotin Conjugated Human CXCL1/GRO alpha Recombinant Rabbit mAb [PSH05-95]	HA722587B	Human	ELISA(Det)
	Human CXCL2/MIP-2 Recombinant Rabbit mAb [PSH17-89] - BSA and Azide free (Detector)	HA725248	Human	ELISA(Det)
	Biotin Conjugated Human CXCL2/MIP-2 Recombinant Rabbit mAb [PSH17-89] - Detector	HA725248B	Human	ELISA(Det),ELISA
	Human CXCL2/MIP-2 Recombinant Rabbit mAb [PSH17-88] - BSA and Azide free (Capture)	HA725247	Human	ELISA(Cap)
	Human IFN-alpha Recombinant Rabbit mAb [PS00-80] - BSA and Azide free (Capture)	HA721290	Human	ELISA(Cap)
	Human IFN-alpha Recombinant Rabbit mAb [PS00-79] - BSA and Azide free (Detector)	HA721291	Human	ELISA(Det)
	Biotin Conjugated Human IFN-alpha Recombinant Rabbit mAb [PS00-79] - BSA and Azide free (Detector)	HA721369	Human	ELISA(Det)
	Human IL-18 Recombinant Rabbit mAb [PSH01-70] - BSA and Azide free (Detector)	HA721723	Human	ELISA(Det)
Relative expression	Human IL-18 Recombinant Rabbit mAb [PSH02-09]	HA721768	Human	WB
	Human IL-18 Recombinant Rabbit mAb [PSH01-71] - BSA and Azide free (Capture)	HA721722	Human	ELISA(Cap)
	Human IL-18 Recombinant Rabbit mAb [PSH02-08]	HA721767	Human	WB,IHC-P,IF-Tissue,IP
	Biotin Conjugated Human IL-18 Recombinant Rabbit mAb [PSH01-70]	HA722644B	Human	ELISA(Det)
	Human IL-18 Recombinant Rabbit mAb [PSH02-08] - BSA and Azide free	HA750778	Human	WB,IHC-P,IF-Tissue
Relative expression	Human IL-18 Recombinant Rabbit mAb [PSH02-09] - BSA and Azide free	HA750779	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-87] - BSA and Azide free (Capture)	HA723883	Human	ELISA(Cap)
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA722649B	Human	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - BSA and Azide free (Detector)	HA723882	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - Detector	HA723882B	Human	ELISA(Det),ELISA
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free	HA750745	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free (Detector)	HA721279	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free	HA750744	Human	WB
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39]	HA721712	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free (Capture)	HA721278	Human	ELISA(Cap)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA721713	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-85] - BSA and Azide free (Capture)	HA723881	Human	ELISA(Cap)
	Biotin Conjugated Human IL-33 Recombinant Rabbit mAb [PSH03-06]	HA722628B	Human	ELISA(Det)
	Human IL-33 Recombinant Rabbit mAb [PSH03-05] - BSA and Azide free (Capture)	HA721935	Human	ELISA(Cap)
	Human IL-33 Recombinant Rabbit mAb [PSH03-06] - BSA and Azide free (Detector)	HA721936	Human	ELISA(Det)
	HRP Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA721293	Human	ELISA(Det)
	Human IL-6 Mouse mAb [A6D2]	HA601050	Human	WB,ELISA
	Biotin Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA722233B	Human	ELISA(Det)
	Human IL-6 Recombinant Rabbit mAb [PS00-21] - BSA and Azide free (Capture)	HA721167	Human	ELISA(Cap)
	Human IL-6 Recombinant Rabbit mAb [PS00-15] - BSA and Azide free (Detector)	HA721168	Human	ELISA(Det)
	Human IL-6 Mouse mAb [A6D1]	HA601051	Human	WB,ELISA
	Biotin Conjugated Human IL-8 Recombinant Rabbit mAb [PS01-22]	HA722611B	Human	ELISA(Det)
	Human IL-8 Recombinant Rabbit mAb [PS00-83] - BSA and Azide free (Capture)	HA722098	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-8 Recombinant Rabbit mAb [PS01-22] - BSA and Azide free (Detector)	HA722099	Human	ELISA(Det)
	Human NLRP3 Recombinant Rabbit mAb [PSH11-08] - BSA and Azide free (Detector)	HA723276	Human	ELISA(Det)
	Biotin Conjugated Human NLRP3 Recombinant Rabbit mAb [PSH11-08] - Detector	HA723277B	Human	ELISA(Det),ELISA
	Human NLRP3 Recombinant Rabbit mAb [PSH11-07] - BSA and Azide free (Capture)	HA723275	Human	ELISA(Cap)
	Human NLRP3 Recombinant Rabbit mAb [PSH11-09] - BSA and Azide free (Detector)	HA723278	Human	ELISA(Det)
	Biotin Conjugated Human NLRP3 Recombinant Rabbit mAb [PSH11-09] - Detector	HA723279B	Human	ELISA(Det),ELISA
	HRP Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721294	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40]	HA721286	Human	IHC-P,mIHC
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free	HA750593	Human	IHC-P
	Biotin Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721559	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free (Detector)	HA721266	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free (Capture)	HA721265	Human	ELISA(Cap)
	Human TNF alpha Recombinant Rabbit mAb	IRS040	Human	mIHC
	IFN alpha R1 Mouse mAb [2G1]	RT1319	Mouse	WB,IP,FC
	IFNAR1 Recombinant Rabbit mAb [SR45-08]	ET1602-37	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	IFNAR1 Recombinant Rabbit mAb [SR45-08] - BSA and Azide free	HA750046	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	IFNAR1 Rabbit pAb	R1511-30	Human	IF-Cell,IHC-P,FC,IF-Tissue
	IFNAR1 Mouse mAb [A5-A3]	EM1701-46	Human, Mouse	WB,IHC-P,IF-Cell
	IFNAR2 Rabbit pAb	HA500425	Human, Rat	WB
	IFNAR2 Rabbit pAb	HA500284	Human	IHC-P,IF-Cell,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	IKB alpha Recombinant Rabbit mAb [SZ00-07] - BSA and Azide free	HA750068	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	IKB alpha Recombinant Rabbit mAb [SZ00-07]	ET1603-6	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	IKB beta Recombinant Rabbit mAb [SN0648]	ET1611-1	Human, Monkey	WB,IF-Cell,IHC-P,IP,FC
	IKBKE Recombinant Rabbit mAb [JU06-72]	ET1706-20	Human	WB,IF-Cell,FC,IHC-P
	IKK alpha Recombinant Rabbit mAb [SN06-99] - BSA and Azide free	HA750246	Human, Mouse, Rat	WB,IHC-P,IP,FC
	IKK alpha Rabbit pAb	ER30911	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	IKK alpha Recombinant Rabbit mAb [SN06-99]	ET1611-15	Human, Mouse, Rat	WB,IHC-P,IP,FC
KD	IKK alpha + IKK beta Recombinant Rabbit mAb [SN63-02] - BSA and Azide free	HA750248	Human, Mouse, Rat	WB,IF-Cell,IP
KD	IKK alpha + IKK beta Recombinant Rabbit mAb [SN63-02]	ET1611-23	Human, Mouse, Rat	WB,IF-Cell,IP
	IKK beta Rabbit pAb	R1706-13	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	IKK beta Recombinant Rabbit mAb [JE32-81]	HA721715	Human	WB,FC
KD	IKK gamma Recombinant Rabbit mAb [JA11-83]	ET1704-40	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,IP,FC
	IL-18 Recombinant Rabbit mAb [PS01-06] - BSA and Azide free	HA750677	Human	WB,FC
Relative expression	IL-18 Recombinant Rabbit mAb [PSH12-04] - BSA and Azide free	HA751511	Mouse	WB,IF-Cell,FC
	IL-18 Recombinant Rabbit mAb [PS00-41] - BSA and Azide free	HA750676	Human	WB,IHC-P
Relative expression	IL-18 Recombinant Rabbit mAb [PSH18-15]	HA723978	Human, Mouse, Rat	WB,IHC-P,IP
	IL-18 Recombinant Rabbit mAb [PS01-06]	HA721537	Human	WB,FC
Relative expression	IL-18 Recombinant Rabbit mAb [PSH12-04]	HA723607	Mouse	WB,IF-Cell,FC
	IL-18 Recombinant Rabbit mAb [PS00-41]	HA721536	Human	WB,IHC-P
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R]	HA601036	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Recombinant Rabbit mAb [JJ087-3]	ET1701-39	Human	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R] - BSA and Azide free	HA610004	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	IL-1 beta Recombinant Rabbit mAb [PSH18-00]	HA723965	Mouse, Rat	WB
Cell treatment	IL-1 beta Mouse mAb [A7A10]	HA601035	Human	WB,ELISA,IF-Cell
Cell treatment	IL-1 beta Mouse mAb [A7F9]	HA601002	Human, Mouse, Rat	WB,ELISA,IF-Cell
	IL33 Rabbit pAb	ER1911-36	Human, Rat	WB,IHC-P,IF-Tissue
	IL-6 Mouse mAb	EM30301	Human	WB
	IL-6 Rabbit pAb	R1203-2	Mouse	WB,ELISA
	IL-6 Rabbit pAb	R1412-2	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	IL-6 Mouse mAb [3-A4]	M1303-3	Mouse	WB,IF-Cell,FC,ELISA
	IL-6 Mouse mAb [1-6]	EM1701-45	Mouse, Rat	WB,IF-Cell
	IL-6 Recombinant Rabbit mAb	IRS160RB	Human	mIHC
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-97]	HA723264	Human	WB
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-96] - BSA and Azide free	HA751386	Human	WB,FC
	IL-8 Rabbit pAb	ER1706-67	Human	ELISA,IF-Cell,IHC-P
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-96]	HA723263	Human	WB,FC
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-97] - BSA and Azide free	HA751387	Human	WB
Cell treatment	IL-8 Rabbit pAb	R1511-15	Human	WB,IHC-P,IF-Cell,FC
Cell treatment	IL-8 Rabbit pAb	ER1901-61	Human	WB,IHC-P,IF-Cell,FC
	Interferon alpha Rabbit pAb	ER1802-4	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	Interferon alpha 2 Recombinant Rabbit mAb [PSH08-23] - BSA and Azide free	HA751214	Human	WB,IF-Cell,IHC-P
Cell treatment	Interferon alpha 2 Recombinant Rabbit mAb [PSH08-23]	HA722967	Human	WB,IF-Cell,IHC-P
Relative expression	IP3 receptor Recombinant Mouse mAb [PSH14-18] - BSA and Azide free	HA610290	Mouse, Rat	WB,IHC-P
Relative expression	IP3 receptor Recombinant Mouse mAb [PSH14-18]	HA601456	Mouse, Rat	WB,IHC-P
	IP3 Receptor Rabbit pAb	ER1803-17	Human, Mouse, Rat	WB,IHC-P
	IP3 Receptor Recombinant Rabbit mAb [JA11-35]	ET1704-77	Human, Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Tissue
	IRAK4 Recombinant Rabbit mAb [JE41-77]	HA721347	Human, Mouse, Monkey	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	IRAK4 Mouse mAb [B9-H2]	EM1709-86	Human, Mouse, Monkey	WB,IHC-P,FC
KD	IRF3 Recombinant Rabbit mAb [SD2062]	ET1612-14	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	IRF3 Recombinant Rabbit mAb [SD2062] - BSA and Azide free	HA750283	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	IRF7 Recombinant Rabbit mAb [SC0617] - BSA and Azide free	HA750234	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	IRF7 Recombinant Rabbit mAb [SC0617]	ET1610-89	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Cell treatment	IRF9 Recombinant Rabbit mAb [PSH08-20]	HA722963	Human	WB,IF-Cell,IP,ChIP
Cell treatment	IRF9 Recombinant Rabbit mAb [PSH07-87]	HA722916	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
Cell treatment	IRF9 Recombinant Rabbit mAb [PSH08-20] - BSA and Azide free	HA751211	Human	WB,IF-Cell,IP,ChIP
Cell treatment	IRF9 Recombinant Rabbit mAb [PSH07-87] - BSA and Azide free	HA751193	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	ITPR2 Rabbit pAb	ER1803-65	Human, Mouse, Rat	Dot Blot,IF-Cell,IHC-P,FC
	JAK1 Mouse mAb [1G1]	RT1341	Human, Mouse, Rat	WB,IP,IF,IHC-P
	JAK1 Recombinant Rabbit mAb [JM75-03]	ET1705-84	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell,IF-Tissue
	JAK1 Recombinant Rabbit mAb [JM75-03] - BSA and Azide free	HA750442	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell,IF-Tissue
	JNK1 Rabbit pAb	R1309-1	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	JNK1 Recombinant Rabbit mAb [JE63-04]	HA721899	Human, Mouse, Rat	WB,IHC-P
	JNK1+JNK2+JNK3 Recombinant Rabbit mAb [SA43-06] - BSA and Azide free	HA750018	Human, Mouse, Rat, Cow, Monkey, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP
	JNK1+JNK2+JNK3 Recombinant Rabbit mAb [SA43-06]	ET1601-28	Human, Mouse, Rat, Cow, Monkey, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP
	JNK1+JNK3 Recombinant Rabbit mAb [JJ08-86]	ET1701-47	Human, Mouse, Rat	WB,IP,IHC-P
	JNK2 Rabbit pAb	HA500110	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	JNK2 Recombinant Rabbit mAb [SC05-85]	ET1610-11	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	JNK3 Recombinant Rabbit mAb [SD082-09]	ET1612-68	Human, Mouse, Rat	WB,IF-Tissue,IHC-P
	JunD Recombinant Rabbit mAb [SD0830]	ET1612-92	Human	WB,IF-Cell,IF-Tissue,IP,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	JunD Recombinant Rabbit mAb [JE59-51]	HA722913	Human, Mouse, Rat	WB,IF-Cell,IHC-P,ChIP,IP
	JunD Recombinant Rabbit mAb [PSH08-85]	HA723037	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP,ChIP
	JunD Recombinant Rabbit mAb [PSH08-85] - BSA and Azide free	HA751250	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP,ChIP
Cell treatment, KD	LC3B Recombinant Rabbit mAb [JJ090-6]	ET1701-65	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,IP,mIHC
Cell treatment	LC3B Recombinant Antibody [JJ090-6] - Mouse IgG1 (Chimeric)	HA601532	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment, KD	LC3B Recombinant Rabbit mAb [JJ090-6] - BSA and Azide free	HA750319	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,IP
Cell treatment	LC3B Recombinant Antibody [JJ090-6] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610352	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	MAP1LC3A Recombinant Rabbit mAb [SN07-40]	ET1611-37	Human, Mouse, Rat	WB,IHC-P
Cell treatment	MAP1LC3A Recombinant Rabbit mAb [ST47-03]	ET1609-26	Human, Mouse, Rat	WB,IF-Cell,FC,IP
	MAP3K7 Mouse mAb [1G5]	RT1377	Human, Mouse, Rat	WB,IP,IF
	MATER Rabbit pAb	ER50101	Human, Mouse	IF-Cell
	MAVS Recombinant Rabbit mAb [PSH01-65] - BSA and Azide free	HA750741	Human	WB,IHC-P,IF-Cell,IF-Tissue
	MAVS Rabbit pAb	ER1912-80	Mouse	WB
	MAVS Recombinant Rabbit mAb [PSH0-30] - BSA and Azide free	HA750604	Human, Mouse	WB,IHC-P,IF-Cell,IP
	MAVS Recombinant Rabbit mAb [PSH0-30]	HA721310	Human, Mouse	WB,IHC-P,IF-Cell,IP
	MAVS Recombinant Rabbit mAb [PSH01-65]	HA721708	Human	WB,IHC-P,IF-Cell,IF-Tissue
	MAVS Rabbit pAb	HA500530	Human	WB,IF-Cell,IHC-P,FC
	MCU Rabbit pAb	ER1803-57	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	MCU Recombinant Rabbit mAb [PSH13-74] - BSA and Azide free	HA751487	Human, Mouse, Rat, Monkey	WB,IHC-Fr,IHC-P,IF-Cell,IP
	MCU Mouse mAb [A2E11]	EM1901-85	Human, Mouse, Rat	WB
	MCU Recombinant Rabbit mAb [PSH13-74]	HA723563	Human, Mouse, Rat, Monkey	WB,IHC-Fr,IHC-P,IF-Cell,IP
	Mitofusin 1 Recombinant Rabbit mAb [JF0954]	ET1702-01	Human	WB,IP
	Mitofusin 2 Rabbit pAb	ER1802-23	Human, Rat	WB,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mitofusin 2 Recombinant Rabbit mAb [JE59-66]	HA720073	Human, Mouse, Rat	WB
	MK13 Rabbit pAb	ER60146	Human, Mouse, Rat	WB
	Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-44] - BSA and Azide free (Capture)	HA725040	Mouse	ELISA(Cap)
	Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-45] - BSA and Azide free (Detector)	HA725041	Mouse	ELISA(Det)
	Biotin Conjugated Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-45] - Detector	HA725042B	Mouse	ELISA(Det),ELISA
	Mouse CCL5 Recombinant Rat mAb	HA601395	Human, Mouse	ELISA
	Biotin Conjugated Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-16] - Detector	HA725048B	Mouse	ELISA(Det),ELISA
	Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-15] - BSA and Azide free (Capture)	HA725046	Mouse	ELISA(Cap)
	Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-16] - BSA and Azide free (Detector)	HA725047	Mouse	ELISA(Det)
	Mouse CXCL2 Recombinant Rabbit mAb [PSH12-60] - BSA and Azide free (Capture)	HA725070	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CXCL2 Recombinant Rabbit mAb [PSH12-61] - Detector	HA725072B	Mouse	ELISA(Det),ELISA
	Mouse CXCL2 Recombinant Rabbit mAb [PSH12-61] - BSA and Azide free (Detector)	HA725071	Mouse	ELISA(Det)
	Mouse IFN beta Recombinant Rabbit mAb [PSH04-72] - BSA and Azide free (Capture)	HA722172	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IFN beta Recombinant Rabbit mAb [PSH04-73]	HA722605B	Mouse	ELISA(Det),ELISA
	Mouse IFN beta Recombinant Rabbit mAb [PSH04-73] - BSA and Azide free (Detector)	HA722173	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-18 Recombinant Rabbit mAb [PSH12-04] - Detector	HA723418B	Mouse	ELISA(Det),ELISA
	Mouse IL-18 Recombinant Rabbit mAb [PSH12-04] - BSA and Azide free (Detector)	HA723417	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IL-18 Recombinant Rabbit mAb [PSH12-03] - BSA and Azide free (Capture)	HA723416	Mouse	ELISA(Cap)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH09-59] - BSA and Azide free	HA723121	Mouse	ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - BSA and Azide free (Detector)	HA721977	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH09-58] - BSA and Azide free	HA723120	Mouse	ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67] - BSA and Azide free (Detector)	HA722167	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-31] - BSA and Azide free (Capture)	HA721976	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67]	HA722608B	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-66] - BSA and Azide free (Capture)	HA722166	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - Detector	HA722625B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PS01-83] - BSA and Azide free (Detector)	HA721602	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PSH0-60] - Detector	HA722647B	Mouse	ELISA(Det)
	HRP Conjugated Mouse IL-6 Mouse mAb [3-A4]	M1303-2	Mouse	WB,ELISA
	Mouse IL-6 Recombinant Rabbit mAb [PSH14-47] - BSA and Azide free (Capture)	HA723642	Mouse	ELISA(Cap)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-60] - BSA and Azide free (Detector)	HA721413	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PS01-74] - Detector	HA723643B	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PS01-83] - BSA and Azide free (Detector)	HA721601	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PS01-74] - BSA and Azide free (Capture/Detector)	HA721600	Mouse	ELISA(Cap),ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IL-6 Recombinant Rabbit mAb [PSH0-59] - BSA and Azide free (Detector)	HA721412	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-58] - BSA and Azide free (Capture)	HA721411	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PSH0-59] - Detector	HA722646B	Mouse	ELISA(Det)
	Mouse Interferon alpha 1 Recombinant Rabbit mAb [PSH14-77] - BSA and Azide free (Capture)	HA723672	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse Interferon alpha 1 Recombinant Rabbit mAb [PSH14-76] - Detector	HA723671B	Mouse	ELISA(Det),ELISA
	Mouse Interferon alpha 1 Recombinant Rabbit mAb [PSH14-76] - BSA and Azide free (Detector)	HA723670	Mouse	ELISA(Det)
	Mouse Interferon alpha 1 Recombinant Rabbit mAb [PSH14-75] - BSA and Azide free (Capture)	HA723669	Mouse	ELISA(Cap)
	Mouse/Rat CXCL1/GRO alpha Recombinant Rabbit mAb [PSH13-14] - BSA and Azide free (Capture)	HA725028	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat CXCL1/GRO alpha Recombinant Rabbit mAb [PSH13-15] - Detector	HA725030B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat CXCL1/GRO alpha Recombinant Rabbit mAb [PSH13-15] - BSA and Azide free (Detector)	HA725029	Mouse, Rat	ELISA(Det)
	Mouse/Rat IL-18 Recombinant Rabbit mAb [PSH18-85] - BSA and Azide free (Detector)	HA724022	Mouse, Rat	ELISA(Det)
	Mouse/Rat IL-18 Recombinant Rabbit mAb [PSH18-84] - BSA and Azide free (Capture)	HA724021	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat IL-18 Recombinant Rabbit mAb [PSH18-85] - Detector	HA724022B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30] - BSA and Azide free (Detector)	HA721975	Mouse, Rat	ELISA(Det)
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-29] - BSA and Azide free (Capture)	HA721974	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30]	HA722626B	Mouse, Rat	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77] - BSA and Azide free (Detector)	HA722041	Mouse	ELISA(Det)
	Biotin Conjugated Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77]	HA722623B	Mouse	ELISA(Det)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-76] - BSA and Azide free (Capture)	HA722040	Mouse	ELISA(Cap)
	MyD88 Recombinant Rabbit mAb [PSH09-74] - BSA and Azide free	HA751316	Human, Mouse, Rat, Monkey	WB
	MyD88 Recombinant Rabbit mAb [PSH05-58]	HA722378	Human, Monkey	WB,IF-Cell,IHC-P,FC,IP
	MyD88 Recombinant Rabbit mAb [PSH09-74]	HA723139	Human, Mouse, Rat, Monkey	WB
	MyD88 Recombinant Rabbit mAb [SC65-04] - BSA and Azide free	HA750231	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	MyD88 Recombinant Rabbit mAb [PSH05-58] - BSA and Azide free	HA751021	Human, Monkey	WB,IF-Cell,IHC-P,FC,IP
	MyD88 Recombinant Rabbit mAb [SC65-04]	ET1610-81	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	NAIP Rabbit pAb	ER65842	Human, Mouse	IHC-P
	NAIP Recombinant Rabbit mAb [JE39-24]	HA722550	Human	WB
	NAK / TBK1 Recombinant Rabbit mAb [PSH11-69] - BSA and Azide free	HA751425	Human, Mouse, Rat	WB,IF-Cell,FC,IP
	NAK / TBK1 Mouse mAb [2G1]	RT1413	Human, Mouse, Rat	WB,IP
	NAK / TBK1 Recombinant Rabbit mAb [PSH11-69]	HA723365	Human, Mouse, Rat	WB,IF-Cell,FC,IP
KD	NAK / TBK1 Mouse mAb [A8D9]	HA601045	Human, Mouse, Rat	WB
KD	NAK / TBK1 Mouse mAb [A8D10]	HA601044	Human, Mouse, Rat	WB
	NEK7 Recombinant Rabbit mAb [JE51-70]	ET7110-16	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
KD	NF-kB p100 / NFKB2 Mouse mAb [16B2]	EM1901-78	Human	WB,IHC-P,FC
KD	NF-kB p100 / NFKB2 Recombinant Rabbit mAb [JM82-03]	ET1705-45	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
KD	NF-kB p100 / NFKB2 Mouse mAb [16B3]	EM1901-79	Human, Mouse, Rat	WB,IHC-P,FC
KD	NF-kB p105 / p50 Recombinant Rabbit mAb [SZ20-01] - BSA and Azide free	HA750060	Human, Mouse, Rat	WB,IHC-P,ICC,FC,IF-Tissue
	NF-kB p105 / p50 Rabbit pAb	R1309-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	NF-kB p105 / p50 Recombinant Rabbit mAb [SZ20-01]	ET1603-18	Human, Mouse, Rat	WB,IHC-P,ICC,FC,IF-Tissue
	NF-kB p105 / p50 Mouse mAb [A2-7]	M1505-8	Human, Mouse	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	NF-kB p65 Rabbit pAb	HA500402	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04] - BSA and Azide free	HA750058	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27]	HA721307	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04]	ET1603-12	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	NF-kB p65 Rabbit pAb	ER0815	Human, Mouse, Rat, Zebrafish	WB,IHC-P,FC
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27] - BSA and Azide free	HA750601	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	NLRP3 Rabbit pAb	ER1706-72	Human, Mouse	WB,IF-Cell,IHC-P,FC
Cell treatment	NLRP3 Recombinant Rabbit mAb [SC06-23] - BSA and Azide free	HA750236	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	NLRP3 Recombinant Rabbit mAb [PSH19-43]	HA724075	Human, Mouse	WB
Cell treatment	NLRP3 Recombinant Rabbit mAb [SC06-23]	ET1610-93	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	NOX2 / gp91phox Recombinant Rabbit mAb [SN07-16]	ET1611-44	Human, Mouse, Rat	WB
	NOX2 / gp91phox Recombinant Rabbit mAb [SN07-16] - BSA and Azide free	HA750255	Human, Mouse, Rat	WB
Relative expression	NOX2 / gp91phox Recombinant Rabbit mAb [PSH10-59] - BSA and Azide free	HA751362	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IP
Relative expression	NOX2 / gp91phox Rabbit pAb	ER1511-34	Human, Mouse	IF-Cell,IHC-P,FC,WB
Relative expression	NOX2 / gp91phox Recombinant Rabbit mAb [PSH10-59]	HA723224	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IP
	OAS3 Rabbit pAb	HA500039	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	P2RX7 Rabbit pAb	ER1919-21	Human, Mouse	WB,IHC-P,IF-Cell,IF-Tissue
	P2RX7 Mouse mAb [G1-D6]	EM1708-40	Human	WB,IHC-P,FC
	P2X7 Rabbit pAb	ER1901-50	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	P2X7 Rabbit pAb	ER1901-99	Human, Mouse, Rat	WB,IHC-P,FC
KD, Cell treatment	p38 alpha / MAPK14 Recombinant Rabbit mAb [SR43-04] - BSA and Azide free	HA750040	Human, Mouse, Rat	WB,IF-Cell,FC
KD, Cell treatment	p38 alpha / MAPK14 Recombinant Rabbit mAb [SR43-04]	ET1602-26	Human, Mouse, Rat	WB,IF-Cell,FC
KO	p38 alpha / MAPK14 Recombinant Rabbit mAb [JF55-07]	ET1702-65	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	p38 alpha / MAPK14 Rabbit pAb	R1308-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	p38 alpha / MAPK14 Mouse mAb [9-B9]	EM40307	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	p38 alpha/MAPK14 + p38 beta/MAPK11 + p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-60]	HA722019	Human, Mouse, Rat, Monkey	WB
	p38 alpha/MAPK14 + p38 beta/MAPK11 + p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-60] - BSA and Azide free	HA750894	Human, Mouse, Rat, Monkey	WB
	p38 beta/MAPK11 + p38 gamma/MAPK12 + p38 delta/MAPK13 Recombinant Rabbit mAb [PSH03-59] - BSA and Azide free	HA750893	Human, Mouse, Rat, Monkey	WB
	p38 beta/MAPK11 + p38 gamma/MAPK12 + p38 delta/MAPK13 Recombinant Rabbit mAb [PSH03-59]	HA722018	Human, Mouse, Rat, Monkey	WB
	p38 delta / MAPK13 Rabbit pAb	HA500582	Human	WB
	p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-58]	HA722017	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-58] - BSA and Azide free	HA750892	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	p38 MAPK Recombinant Rabbit Multiclonal Antibody [PSH03-72]	HA722036	Human, Mouse, Rat, Monkey	WB
	p38 MAPK Recombinant Rabbit Multiclonal Antibody [PSH03-72] - BSA and Azide free	HA750906	Human, Mouse, Rat, Monkey	WB
	Pannexin 1 Recombinant Rabbit mAb [JE38-97]	HA721629	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-Bcl-2 (S70) Recombinant Rabbit mAb [PSH11-02]	HA723269	Human	WB,IF-Cell,FC
Cell treatment	Phospho-Bcl-2 (S70) Recombinant Rabbit mAb [PSH11-02] - BSA and Azide free	HA751392	Human	WB,IF-Cell,FC
Cell treatment	Phospho-c-Jun (S63) Recombinant Rabbit mAb [SY0297] - BSA and Azide free	HA750146	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-c-Jun (S63) Recombinant Rabbit mAb [SY0297]	ET1608-4	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit mAb [JE65-46]	HA722475	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit mAb [PSH07-12] - BSA and Azide free	HA751133	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit mAb [PSH07-12]	HA722800	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-c-Jun (T91) Recombinant Rabbit mAb [JJ080-9]	ET1701-32	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-DRP1 (S616) Recombinant Rabbit mAb [PSH06-77] - BSA and Azide free	HA751109	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC
Cell treatment	Phospho-DRP1 (S616) Recombinant Rabbit mAb [PSH06-77]	HA722760	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4]	ET1603-22	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4] - BSA and Azide free	HA750061	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01]	ET1610-13	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01] - BSA and Azide free	HA750199	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-IKB alpha (S32) Recombinant Rabbit mAb [ST53-05]	ET1609-78	Human, Mouse	WB,IF-Cell,IP
Cell treatment	Phospho-IKB alpha (S32) Recombinant Rabbit mAb [ST53-05] - BSA and Azide free	HA750196	Human, Mouse	WB,IF-Cell,IP
Cell treatment	Phospho-IKB alpha (S32/S36) Recombinant Rabbit mAb [PSH06-76]	HA722770	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Phospho-IKB alpha (S32/S36) Recombinant Rabbit mAb [PSH06-76] - BSA and Azide free	HA751119	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Phospho-IKB alpha (S36) Recombinant Rabbit mAb [JE56-65]	HA721802	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-IRF3 (S386) Recombinant Rabbit mAb [SU03-28] - BSA and Azide free	HA750138	Human	WB,IF-Cell,IF-Tissue
Cell treatment	Phospho-IRF3 (S386) Recombinant Rabbit mAb [SU03-28]	ET1608-22	Human	WB,IF-Cell,IF-Tissue
Cell treatment	Phospho-IRF3 (S396) Recombinant Rabbit mAb [PSH06-75] - BSA and Azide free	HA751121	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-IRF3 (S396) Recombinant Rabbit mAb [PSH06-75]	HA722772	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-JNK1/2/3 (T183 + T183 + T221) Recombinant Rabbit mAb [ST500]	ET1609-42	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
Cell treatment	Phospho-JNK1/2/3 (T183 + T183 + T221) Recombinant Rabbit mAb [ST500] - BSA and Azide free	HA750173	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
Cell treatment	Phospho-JunD (S255) Recombinant Rabbit mAb [JJ08-21]	ET1701-35	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Phospholipase C beta 1 Rabbit pAb	HA500302	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Phospholipase C beta 2 Rabbit pAb	ER1915-09	Human, Mouse	WB,IHC-P,IF-Tissue,FC
	Phospholipase C beta 3/PLCB3 Recombinant Rabbit mAb [JE39-14]	HA721561	Human	WB,IHC-P
Cell treatment	Phospho-NF-kB p65 (S529) Recombinant Rabbit mAb [SP07-00]	ET1604-27	Human	WB,IHC-P,FC,IP
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58]	HA723223	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58] - BSA and Azide free	HA751361	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-p38 alpha (T180) Recombinant Rabbit mAb [PSH06-25]	HA722657	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180) Recombinant Rabbit mAb [PSH06-25] - BSA and Azide free	HA751062	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180 + Y182) Recombinant Rabbit mAb [PSH06-35] - BSA and Azide free	HA751074	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180 + Y182) Recombinant Rabbit mAb [PSH06-35]	HA722669	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-p38 (T180) Recombinant Rabbit mAb [JE77-37]	HA722150	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-p38 (T180) Rabbit pAb	ER2001-51	Human, Mouse, Rat	WB
Cell treatment	Phospho-p38 (T180 + Y182) Rabbit pAb	ER2001-52	Human, Mouse, Rat	WB
	Phospho-p38 (T180 + Y182) Rabbit pAb	ER1903-01	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH15-91]	HA723796	Human, Mouse, Rat	WB,IF-Cell,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH16-27]	HA723825	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH15-91] - BSA and Azide free	HA751582	Human, Mouse, Rat	WB,IF-Cell,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH16-27] - BSA and Azide free	HA751595	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,Dot Blot
Cell treatment	Phospho-PKC alpha (T497) Recombinant Rabbit mAb [JE44-81]	HA722278	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-PKC alpha (T638) Recombinant Rabbit mAb [JF0964]	ET1702-17	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,Dot Blot
Cell treatment	Phospho-PKC delta (Y311) Recombinant Rabbit mAb [JE44-20]	HA722458	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-PKC (S729) Recombinant Rabbit mAb [JE77-54]	HA722846	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-PKC zeta (T560) Recombinant Rabbit mAb [SN0753]	ET1611-78	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-RIP3 (S232) Recombinant Rabbit mAb [JE44-62]	HA721428	Mouse, Rat	WB,IF-Cell
Cell treatment, Relative expression	Phospho-STAT1 (S727) Recombinant Rabbit mAb [SN67-04]	ET1611-20	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment, Relative expression	Phospho-STAT1 (S727) Recombinant Rabbit mAb [SN67-04] - BSA and Azide free	HA750247	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-STAT1 (Y701) Recombinant Rabbit mAb [PSH04-02]	HA722083	Human	WB,IF-Cell
Cell treatment	Phospho-STAT1 (Y701) Recombinant Rabbit mAb [PSH04-02] - BSA and Azide free	HA750925	Human	WB,IF-Cell
Cell treatment	Phospho-STING (S366) Recombinant Rabbit mAb [PSH09-72]	HA723137	Human	WB,Dot Blot,IP
Cell treatment	Phospho-STING (S366) Recombinant Rabbit mAb [PSH09-72] - BSA and Azide free	HA751314	Human	WB,Dot Blot,IP
Cell treatment	Phospho-TAK1 (S412) Recombinant Rabbit mAb [PSH14-64] - BSA and Azide free	HA751525	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-TAK1 (S412) Recombinant Rabbit mAb [PSH14-64]	HA723662	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-TAK1 (S439) Recombinant Rabbit mAb [JE44-93]	HA722242	Human, Mouse, Rat	WB
	Pin1 Recombinant Rabbit mAb [PSH10-81]	HA723248	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	Pin1 Recombinant Rabbit mAb [PSH10-81] - BSA and Azide free	HA751373	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	PKC Rabbit pAb	ER64620	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
KD	PKC alpha Recombinant Rabbit mAb [SU31-08]	ET1608-15	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	PKC beta Rabbit pAb	ER1802-14	Human, Mouse	WB,IF-Cell,IHC-P
	PKC beta 1 Recombinant Rabbit mAb [JM43-26]	ET1705-30	Human, Mouse, Rat	WB,IHC-P,FC
	PKC beta 2 Recombinant Rabbit mAb [ST48-06]	ET1609-44	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	PKC delta Recombinant Rabbit mAb [JJ093-02] - BSA and Azide free	HA750325	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PKC delta Recombinant Rabbit mAb [JJ093-02]	ET1701-85	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	PKC epsilon Recombinant Rabbit mAb [JA38-21]	ET1704-85	Human, Mouse, Rat	WB,IHC-P
	PKC eta Recombinant Rabbit mAb [JE32-41]	HA722389	Human, Mouse, Rat, Monkey	WB
	PKC gamma Recombinant Rabbit mAb [JG69-32]	ET7108-42	Human, Mouse, Rat	WB,IHC-P
Relative expression	PKC zeta Recombinant Rabbit mAb [PSH06-74] - BSA and Azide free	HA751122	Human	WB
Relative expression	PKC zeta Recombinant Rabbit mAb [PSH06-74]	HA722773	Human	WB
KD	PKN2 Recombinant Rabbit mAb [JE62-44]	HA720102	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	PLCB2 Rabbit pAb	HA500399	Human	WB,IF-Cell
	PLCB3 Rabbit pAb	HA500391	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	PLC beta 3 Rabbit pAb	ER64629	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	PRMT1 Recombinant Rabbit mAb [JA11-17]	ET1704-19	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Pyrin Recombinant Rabbit mAb [PSH06-12]	HA722509	Human	WB,FC
Cell treatment	Pyrin Recombinant Rabbit mAb [PSH06-12] - BSA and Azide free	HA751050	Human	WB,FC
	Rab 2A Rabbit pAb	ER64714	Human, Mouse, Rat	WB,IHC-P
	RANTES Recombinant Rabbit mAb [JM03-45]	ET1705-70	Human, Mouse	WB,IHC-P,IF-Cell,IF-Tissue
	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60] - BSA and Azide free (Detector)	HA721837	Rat	ELISA(Det)
	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - Detector	HA721611	Rat	ELISA(Det),ELISA
	Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - BSA and Azide free (Detector)	HA721610	Rat	ELISA(Det)
	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60]	HA722634B	Rat	ELISA(Det)
	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-59] - BSA and Azide free (Capture)	HA721836	Rat	ELISA(Cap)
	Rat IL-1 beta Recombinant Rabbit mAb [PS01-86] - BSA and Azide free (Capture)	HA721609	Rat	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Rat IL-6 Recombinant Rabbit mAb [PS01-94] - BSA and Azide free (Detector)	HA721617	Rat	ELISA(Det)
	Rat IL-6 Recombinant Rabbit mAb [PS01-90] - BSA and Azide free (Capture)	HA721615	Rat	ELISA(Cap)
	Rat IL-6 Recombinant Rabbit mAb [PS01-94] - BSA and Azide free (Detector)	HA721616	Rat	ELISA(Det)
	Biotin Conjugated Rat TNF alpha Recombinant Rabbit mAb [PSH04-71]	HA722606B	Rat	ELISA(Det)
	Rat TNF alpha Recombinant Rabbit mAb [PSH04-70] - BSA and Azide free (Capture)	HA722170	Rat	ELISA(Cap)
	Rat TNF alpha Recombinant Rabbit mAb [PSH04-71] - BSA and Azide free (Detector)	HA722171	Rat	ELISA(Det)
	RhoA Recombinant Rabbit mAb [SN0612]	ET1611-10	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	RhoA Recombinant Rabbit mAb [SN0612] - BSA and Azide free	HA750242	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Rho A+B+C Recombinant Rabbit mAb [SR38-00]	ET1602-23	Human, Mouse, Rat	WB,IF-Cell,FC
	Rho B Rabbit pAb	ER64810	Human, Mouse, Rat	WB
	RIP Recombinant Rabbit mAb [JJ092-04]	ET1701-79	Human	WB,FC,IP
	RIP Recombinant Rabbit mAb [PSH19-57]	HA724083	Human, Mouse, Rat, Monkey	WB
	RIP2 Rabbit pAb	ER1915-87	Human, Mouse	WB,IHC-P,IF-Tissue,FC
KO	RIP3 Rabbit pAb	ER1901-27	Mouse, Rat	WB,IHC-P
	RIP3 Rabbit pAb	ER1902-67	Human	WB,IHC-P,FC
Relative expression	RIP3 Recombinant Rabbit mAb [PSH04-79]	HA722183	Human	WB,IHC-P
Relative expression	RIP3 Recombinant Rabbit mAb [PSH04-78] - BSA and Azide free	HA750972	Human	WB,IHC-P
Relative expression	RIP3 Recombinant Rabbit mAb [PSH04-78]	HA722182	Human	WB,IHC-P
Relative expression	RIP3 Recombinant Rabbit mAb [PSH04-79] - BSA and Azide free	HA750973	Human	WB,IHC-P
KD	RIPK2 Recombinant Rabbit mAb [JE51-90]	HA721300	Human	WB
Cell treatment	RNase L Recombinant Rabbit mAb [JE36-11]	HA722221	Human	WB
KD	STAT1 Recombinant Rabbit mAb [SD20-75]	ET1612-22	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	STAT1 Rabbit pAb	R1408-2	Human, Mouse, Rat, Zebrafish	WB,IF-Cell
	STAT1 Recombinant Mouse mAb [G3-B11-R]	HA601309	Human	WB,IHC-P
	STAT1 Mouse mAb [G3-B11]	M1407-1	Human	WB,IF-Cell,IHC-P,IF-Tissue
	STAT1 Recombinant Mouse mAb [G3-B11-R] - BSA and Azide free	HA610181	Human	WB,IHC-P
KD	STAT1 alpha Recombinant Rabbit mAb [SJ01-89]	ET1606-39	Human, Mouse	WB,IP,FC,IHC-P,IF-Cell,IF-Tissue
KD	STAT1 alpha Recombinant Rabbit mAb [SJ01-89] - BSA and Azide free	HA750102	Human, Mouse	WB,IP,FC,IHC-P,IF-Cell,IF-Tissue
	STAT1 alpha/beta Rabbit pAb	ER31221	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	STAT2 Recombinant Rabbit mAb [JF0884]	ET1702-18	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	STING Recombinant Rabbit mAb [JM03-47] - BSA and Azide free	HA750437	Human, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue,IP
	STING Recombinant Rabbit mAb [PSH07-44]	HA722832	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,IP
	STING Recombinant Rabbit mAb [JM03-47]	ET1705-68	Human, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue,IP
	STING Recombinant Rabbit mAb [PSH07-44] - BSA and Azide free	HA751165	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,IP
	TAB1 Mouse mAb [1G2]	RT1603	Human	WB,IP,IF
KD	TAB1 Rabbit pAb	HA500400	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	TAB3 Rabbit pAb	HA500421	Human, Mouse	WB,IHC-P
KD	TAK1 Recombinant Rabbit mAb [JM73-19]	ET1705-14	Human, Mouse, Rat, Monkey	WB
	TAK1 Rabbit pAb	ER1802-90	Human, Mouse, Rat	IF-Cell,IHC-P,FC
KD	TAK1 Recombinant Rabbit mAb [JM73-19] - BSA and Azide free	HA750424	Human, Mouse, Rat, Monkey	WB
	TBP-1 Rabbit pAb	HA500413	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
KO	TDG Recombinant Rabbit mAb [PSH10-11]	HA723180	Human, Mouse, Rat	WB,IHC-P
	TDG Rabbit pAb	ER65129	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KO	TDG Recombinant Rabbit mAb [PSH10-11] - BSA and Azide free	HA751531	Human, Mouse, Rat	WB,IHC-P
	Thioredoxin Rabbit pAb	R1511-16	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Thioredoxin Recombinant Rabbit mAb [JF101-8]	ET1702-90	Human, Mouse, Rat	WB,IHC-P,FC
	Thioredoxin Recombinant Rabbit mAb [JM10-019]	ET1703-03	Human	WB,IF-Cell,IF-Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Thioredoxin 2 Recombinant Rabbit mAb [JE35-76]	HA722142	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	TLR4 Rabbit pAb	ER1706-43	Human, Mouse, Rat	IF-Cell,IHC-P,FC
	TLR4 Mouse mAb	RT1666	Human, Mouse, Rat	WB
	TMS1 Recombinant Rabbit mAb [SN07-10]	ET1611-62	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	TMS1 / ASC Recombinant Rabbit mAb [PSH17-34]	HA723917	Human, Mouse	WB,IHC-Fr,IHC-P
Relative expression	TMS1 / ASC Recombinant Rabbit mAb [PSH17-34] - BSA and Azide free	HA751622	Human, Mouse	WB,IHC-Fr,IHC-P
	TNFAIP3 Recombinant Rabbit mAb [SN07-31] - BSA and Azide free	HA750253	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	TNFAIP3 Rabbit pAb	ER1706-32	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	TNFAIP3 Recombinant Rabbit mAb [SN07-31]	ET1611-40	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	TNF alpha Rabbit pAb	ER65189	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	TNF alpha Rabbit pAb	ER1919-22	Human, Mouse	WB
Cell treatment	TNF alpha Rabbit pAb	R1203-1	Mouse	WB,IF-Cell,ELISA
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free	HA750897	Human, Mouse	WB
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07]	HA722022	Human, Mouse	WB
Relative expression, KD	TOMM20 Recombinant Rabbit mAb [ST04-72]	ET1609-25	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP,IHC-Fr
Relative expression, KD	TOMM20 Recombinant Rabbit mAb [ST04-72] - BSA and Azide free	HA750166	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP,IHC-Fr
	TOMM20 Recombinant Antibody [ST04-72] - Mouse IgG1 (Chimeric)	HA601454	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	TOMM20 Recombinant Antibody [ST04-72] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610288	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	TRAF2 Recombinant Rabbit mAb [SD205-06] - BSA and Azide free	HA750297	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
	TRAF2 Recombinant Rabbit mAb [SD205-06]	ET1612-5	Human, Mouse, Rat	WB,IHC-P,IP
	TRAF3 Recombinant Rabbit mAb [PSH01-30] - BSA and Azide free	HA750708	Human, Mouse, Rat	WB
	TRAF3 Recombinant Rabbit mAb [PSH01-30]	HA721654	Human, Mouse, Rat	WB
	TRAF3 Rabbit pAb	ER65205	Human, Mouse	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	TRAF6 Rabbit pAb	R1311-2	Human, Mouse, Rat, Chicken, Zebrafish	IF-Cell, WB
	TRAF6 Recombinant Rabbit mAb [JE01-42]	HA722857	Human, Mouse, Rat	WB, IF-Cell, FC
	TRIF Recombinant Rabbit mAb [PSH05-39] - BSA and Azide free	HA750998	Human	WB
	TRIF Rabbit pAb	ER65721	Human	WB
	TRIF Recombinant Rabbit mAb [PSH05-39]	HA722268	Human	WB
	TRPM2 Rabbit pAb	HA500437	Human, Mouse, Rat	WB, FC
	TRPM7 Rabbit pAb	ER61334	Human, Rat, Mouse	WB
	TRPM7 Recombinant Rabbit mAb [JE38-36]	HA721501	Human, Mouse, Rat	WB, FC
KO	TXNIP Recombinant Rabbit mAb [JM60-35] - BSA and Azide free	HA750438	Human, Mouse, Rat	WB, IF-Cell, IF-Tissue, IHC-P, FC
KO	TXNIP Recombinant Rabbit mAb [JM60-35]	ET1705-72	Human, Mouse, Rat	WB, IF-Cell, IF-Tissue, IHC-P, FC
	TYK2 Rabbit pAb	ER1803-05	Human, Mouse, Rat	WB, IF-Cell, IHC-P, FC
	TYK2 Recombinant Rabbit mAb [PSH18-53]	HA724006	Human, Mouse, Rat, Monkey	WB, IHC-P
KD	VDAC1 Rabbit pAb	R1307-1	Human, Mouse, Rat, Zebrafish	WB, IF-Cell, IHC-P
KD	VDAC1 Recombinant Rabbit mAb [SA93-03]	ET1601-20	Human, Mouse, Rat	WB, IF-Cell, IHC-P, FC, IF-Tissue
KD	VDAC1 Recombinant Rabbit mAb [SA93-03] - BSA and Azide free	HA750012	Human, Mouse, Rat	WB, IF-Cell, IHC-P, FC, IF-Tissue
	VDAC1 / 2 / 3 Recombinant Rabbit mAb [PSH14-90]	HA723692	Human, Mouse, Rat	WB
	VDAC1 / 2 / 3 Recombinant Rabbit mAb [PSH14-90] - BSA and Azide free	HA751534	Human, Mouse, Rat	WB
KD	XIAP Recombinant Rabbit mAb [PSH04-30]	HA722113	Human, Mouse, Rat, Monkey	WB, IHC-P
	XIAP Rabbit pAb	ER65320	Human, Mouse, Rat	WB, IHC-P
KD	XIAP Recombinant Rabbit mAb [PSH04-30] - BSA and Azide free	HA750939	Human, Mouse, Rat, Monkey	WB, IHC-P
	XIAP/BIRC4 Rabbit pAb	ER1918-75	Human, Mouse, Rat	WB, IHC-P, IF-Cell, IF-Tissue, FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层