

细胞因子研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	Angiotensinogen Recombinant Rabbit mAb [SD201-02]	ET1612-48	Human	WB,IF-Cell,IF-Tissue,IP
	BCL9L Rabbit pAb	HA500459	Human, Mouse, Rat	WB,IHC-P
	Biotin Conjugated Beta Nerve Growth Factor Recombinant Rabbit mAb [PSH07-81]	HA722905B	Human	ELISA(Det)
	bFGF Rabbit pAb	ER1904-22	Human, Mouse	WB,FC
	C1s Recombinant Rabbit mAb [JE61-95]	HA722292	Human	WB,IHC-P
	CCL23 Rabbit pAb	ER61224	Human	WB
	CCL25 Rabbit pAb	ER61225	Human, Mouse	WB
Cell treatment, Relative expression	CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-28]	HA723310	Human, Mouse	WB,IF-Cell,IHC-P
Cell treatment	CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH19-79]	HA724103	Mouse, Rat	WB,IF-Cell,FC,IP
	CCL3 / MIP-1 alpha Rabbit pAb	R1510-9	Human, Mouse, Rat	IF-Cell,IHC-P,FC,ELISA
Cell treatment, Relative expression	CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-27] - BSA and Azide free	HA751406	Human	WB,IF-Cell,FC
Cell treatment	CCL3 / MIP-1 alpha Rabbit pAb	R1510-10	Human	WB
Cell treatment, Relative expression	CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-28] - BSA and Azide free	HA751407	Human, Mouse	WB,IF-Cell,IHC-P
Cell treatment, Relative expression	CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-27]	HA723309	Human	WB,IF-Cell,FC
	PE Conjugated CD123 / IL3RA Recombinant Mouse mAb [PSH14-51]	HA601528F	Human	FC
	CD123 / IL3RA Recombinant Mouse mAb [PSH14-51]	HA601490	Human	FC
	CD25 / IL-2 Receptor alpha Recombinant Mouse mAb [PSH14-46]	HA601486	Human, Mouse	FC,IHC-P,IF-Tissue
	APC Conjugated CD25 / IL-2 Receptor alpha Recombinant Mouse mAb [PSH14-46]	HA601530F	Human	FC
Relative expression	CFTR Recombinant Rabbit mAb [PSH09-99] - BSA and Azide free	HA751330	Human	WB,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CFTR Recombinant Rabbit mAb [PSH09-99]	HA723168	Human	WB,IP
Cell treatment, Relative expression	CXCL13 Recombinant Rabbit mAb [PSH04-34]	HA722117	Human	WB,IF-Cell,IHC-P
Cell treatment, Relative expression	CXCL13 Recombinant Rabbit mAb [PSH04-34] - BSA and Azide free	HA750943	Human	WB,IF-Cell,IHC-P
Cell treatment	CXCL13 Rabbit pAb	ER65489	Human	WB
Cell treatment	CXCL1 / GRO alpha Recombinant Rabbit mAb [PSH07-53] - BSA and Azide free	HA751174	Human	WB
Cell treatment	CXCL1 / GRO alpha Recombinant Rabbit mAb [PSH07-53]	HA722841	Human	WB
Relative expression	EAAT2 Recombinant Rabbit mAb [PS01-62]	HA721323	Mouse, Rat, Human	WB,IHC-P,IHC-Fr
Relative expression	EAAT2 Recombinant Rabbit mAb [PS01-62] - BSA and Azide free	HA750609	Mouse, Rat, Human	WB,IHC-P,IHC-Fr
Relative expression	EAAT2 Recombinant Antibody [PS01-62] - Guinea pig IgG2 (Chimeric)	HA601504	Mouse, Rat	IHC-Fr,WB
Relative expression	EAAT2 Recombinant Rabbit mAb [PSH07-43]	HA722831	Human, Mouse, Rat	WB,IHC-P,IHC-Fr
Relative expression	EAAT2 Recombinant Rabbit mAb [PSH07-43] - BSA and Azide free	HA751164	Human, Mouse, Rat	WB,IHC-P,IHC-Fr
Relative expression	EAAT2 Recombinant Antibody [PS01-62] - Rat IgG1 (Chimeric)	HA601403	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
Relative expression	EAAT2 Recombinant Antibody [PS01-62] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610334	Mouse, Rat	IHC-Fr,WB
Relative expression	EAAT2 Recombinant Antibody [PS01-62] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610243	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
	EAAT2 Recombinant Antibody [PS01-62] - Chicken IgY (Chimeric)	HA601503	Mouse	IHC-Fr
	EAAT2 Recombinant Antibody [PS01-62] - Chicken IgY (Chimeric) - BSA and Azide free	HA610333	Mouse	IHC-Fr
	Ebi3 Rabbit pAb	ER65656	Human, Mouse	WB
	EGF Recombinant Rabbit mAb [JA14-15]	ET1704-33	Human, Mouse, Rat	WB,IHC-P,IP
	FGF2 Recombinant Rabbit mAb [JM28-10]	ET1703-18	Human, Mouse	WB,IP,FC,IHC-P
	GM-CSF Mouse mAb [A6C6]	HA601032	Human	WB,ELISA
	GM-CSF Mouse mAb [A6C5]	HA601033	Human	WB,ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	GRO alpha Rabbit pAb	ER63645	Human, Mouse, Rat	IHC-P
	HGF Recombinant Rabbit mAb [PSH20-29]	HA724140	Human	IHC-P
	HGF Recombinant Rabbit mAb [PSH03-79]	HA723896	Human	WB
	HGF Rabbit pAb	ER65684	Human	WB,IHC-P
	HGF Recombinant Rabbit mAb [PSH03-79] - BSA and Azide free	HA751614	Human	WB
KD	HNRNPA0 Recombinant Rabbit mAb [PSH02-77]	HA721855	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
KD	HNRNPA0 Recombinant Rabbit mAb [PSH02-77] - BSA and Azide free	HA750832	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	Human Beta Nerve Growth Factor Recombinant Rabbit mAb [PSH07-80] - BSA and Azide free (Capture)	HA722903	Human	ELISA(Cap)
	Human Beta Nerve Growth Factor Recombinant Rabbit mAb [PSH07-79] - BSA and Azide free (Capture)	HA722902	Human	ELISA(Cap)
	Human Beta Nerve Growth Factor Recombinant Rabbit mAb [PSH07-81] - BSA and Azide free (Detector)	HA722904	Human	ELISA(Det)
	Human CCL23 Recombinant Rabbit mAb [PSH09-05] - BSA and Azide free (Capture)	HA723061	Human	ELISA(Cap)
	Human CCL23 Recombinant Rabbit mAb [PSH09-07] - BSA and Azide free (Detector)	HA723063	Human	ELISA(Det)
	Human CCL23 Recombinant Rabbit mAb [PSH09-06] - BSA and Azide free (Capture)	HA723062	Human	ELISA(Cap)
	Biotin Conjugated Human CCL23 Recombinant Rabbit mAb [PSH09-07] - BSA and Azide free (Detector)	HA723064B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human CCL25 Recombinant Rabbit mAb [PSH08-66] - Detector	HA723018B	Human	ELISA(Det),ELISA
	Human CCL25 Recombinant Rabbit mAb [PSH08-66] - BSA and Azide free (Detector)	HA723016	Human	ELISA(Det)
	Biotin Conjugated Human CCL25 Recombinant Rabbit mAb [PSH08-67] - Detector	HA723019B	Human	ELISA(Det),ELISA
	Human CCL25 Recombinant Rabbit mAb [PSH08-65] - BSA and Azide free (Capture)	HA723015	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CCL25 Recombinant Rabbit mAb [PSH08-67] - BSA and Azide free (Detector)	HA723017	Human	ELISA(Det)
	Human CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-78] - BSA and Azide free (Capture)	HA723380	Human	ELISA(Cap)
	Human CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-76] - BSA and Azide free (Capture)	HA723377	Human	ELISA(Cap)
	Biotin Conjugated Human CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-77] - Detector	HA723379B	Human	ELISA(Det)
	Human CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-77] - BSA and Azide free (Detector)	HA723378	Human	ELISA(Det)
Cell treatment	Human CXCL13 Recombinant Rabbit mAb [PSH07-05] - BSA and Azide free (Capture)	HA722791	Human	ELISA(Cap)
Cell treatment	Biotin Conjugated Human CXCL13 Recombinant Rabbit mAb [PSH07-06] - Detector	HA722793B	Human	ELISA(Det),ELISA
Cell treatment	Biotin Conjugated Human CXCL13 Recombinant Rabbit mAb [PSH07-07] - Detector	HA722795B	Human	ELISA(Det),ELISA
Cell treatment	Human CXCL13 Recombinant Rabbit mAb [PSH07-07] - BSA and Azide free (Detector)	HA722794	Human	ELISA(Det)
Cell treatment	Human CXCL13 Recombinant Rabbit mAb [PSH07-06] - BSA and Azide free (Detector)	HA722792	Human	ELISA(Det)
	Biotin Conjugated Human CXCL1/GRO alpha Recombinant Rabbit mAb [PSH05-95]	HA722587B	Human	ELISA(Det)
	Human CXCL1/GRO alpha Recombinant Rabbit mAb [PSH05-94] - BSA and Azide free (Capture)	HA722493	Human	ELISA(Cap)
	Human CXCL1/GRO alpha Recombinant Rabbit mAb [PSH05-95] - BSA and Azide free (Detector)	HA722494	Human	ELISA(Det)
	Human CXCL6/GCP-2 Recombinant Rabbit mAb [PSH11-30] - BSA and Azide free (Capture)	HA723314	Human	ELISA(Cap)
	Human CXCL6/GCP-2 Recombinant Rabbit mAb [PSH11-31] - BSA and Azide free (Detector)	HA723315	Human	ELISA(Det)
	Biotin Conjugated Human CXCL6/GCP-2 Recombinant Rabbit mAb [PSH11-31] - Detector	HA723316B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human EGF Recombinant Rabbit mAb [PSH02-20] - BSA and Azide free (Detector)	HA721783	Human	ELISA(Det)
	Human EGF Recombinant Rabbit mAb [PSH02-19] - BSA and Azide free (Capture)	HA721782	Human	ELISA(Cap)
	Biotin Conjugated Human EGF Recombinant Rabbit mAb [PSH02-20] - Detector	HA722638B	Human	ELISA(Det)
	Biotin Conjugated Human FGF-2 Recombinant Rabbit mAb [PSH02-09]	HA722639B	Human	ELISA(Det)
	Human FGF-2 Recombinant Rabbit mAb [PSH02-08] - BSA and Azide free (Capture)	HA721780	Human	ELISA(Cap)
	Human FGF-2 Recombinant Rabbit mAb [PSH02-09] - BSA and Azide free (Detector)	HA721781	Human	ELISA(Det)
	Human GM-CSF Recombinant Rabbit mAb [PS01-52] - BSA and Azide free (Capture)	HA721288	Human	ELISA(Cap)
	Biotin Conjugated Human GM-CSF Recombinant Rabbit mAb [PS01-38]	HA722648B	Human	ELISA(Det)
	Human GM-CSF Recombinant Rabbit mAb [PS01-38] - BSA and Azide free (Detector)	HA721289	Human	ELISA(Det)
	Human HGF Recombinant Rabbit mAb [PSH03-79] - BSA and Azide free (Detector)	HA722043	Human	ELISA(Det)
	Human HGF Recombinant Rabbit mAb [PSH02-07] - BSA and Azide free (Detector)	HA721779	Human	ELISA(Det)
	Biotin Conjugated Human HGF Recombinant Rabbit mAb [PSH03-79]	HA722619B	Human	ELISA(Det)
	Human HGF Recombinant Rabbit mAb [PSH02-06] - BSA and Azide free (Capture)	HA721778	Human	ELISA(Cap)
	Human HGF Recombinant Rabbit mAb [PSH03-78] - BSA and Azide free (Capture)	HA722042	Human	ELISA(Cap)
	Biotin Conjugated Human HGF Recombinant Rabbit mAb [PSH02-07]	HA722640B	Human	ELISA(Det)
	Human IFN-alpha Recombinant Rabbit mAb [PS00-80] - BSA and Azide free (Capture)	HA721290	Human	ELISA(Cap)
	Biotin Conjugated Human IFN-alpha Recombinant Rabbit mAb [PS00-79] - BSA and Azide free (Detector)	HA721369	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IFN-alpha Recombinant Rabbit mAb [PS00-79] - BSA and Azide free (Detector)	HA721291	Human	ELISA(Det)
	Human IFNGR1 Recombinant Rabbit mAb [PSH06-44] - BSA and Azide free (Detector)	HA722680	Human	ELISA(Det)
	Human IFNGR1 Recombinant Rabbit mAb [PSH06-43] - BSA and Azide free (Capture)	HA722679	Human	ELISA(Cap)
	Biotin Conjugated Human IFNGR1 Recombinant Rabbit mAb [PSH06-44]	HA722681B	Human	ELISA(Det)
	Human IL-10 Recombinant Rabbit mAb [PS00-81] - BSA and Azide free (Detector)	HA721237	Human	ELISA(Det)
	Biotin Conjugated Human IL-10 Recombinant Rabbit mAb [PS00-81]	HA722652B	Human	ELISA(Det)
	Human IL-10 Recombinant Rabbit mAb [PS00-98] - BSA and Azide free (Capture)	HA721236	Human	ELISA(Cap)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH02-57] - BSA and Azide free (Capture)	HA721834	Human	ELISA(Cap)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH02-58] - BSA and Azide free (Detector)	HA721835	Human	ELISA(Det)
	Biotin Conjugated Human IL-12B p40 Recombinant Rabbit mAb [PSH02-58]	HA722633B	Human	ELISA(Det)
	Biotin Conjugated Human IL-12B p40 Recombinant Rabbit mAb [PSH03-04]	HA722629B	Human	ELISA(Det)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH03-04] - BSA and Azide free (Detector)	HA721934	Human	ELISA(Det)
	Human IL-13R alpha 2 Recombinant Rabbit mAb [PSH07-63] - BSA and Azide free (Capture)	HA722880	Human	ELISA(Cap)
	Human IL-13R alpha 2 Recombinant Rabbit mAb [PSH07-64] - BSA and Azide free (Detector)	HA722881	Human	ELISA(Det)
	Biotin Conjugated Human IL-13R alpha 2 Recombinant Rabbit mAb [PSH07-64]	HA722882B	Human	ELISA(Det)
	Human IL-15 Recombinant Rabbit mAb [PSH04-20] - BSA and Azide free (Detector)	HA722103	Human	ELISA(Det)
	Human IL-15 Recombinant Rabbit mAb [PSH04-19] - BSA and Azide free (Capture)	HA722102	Human	ELISA(Cap)
	Biotin Conjugated Human IL-15 Recombinant Rabbit mAb [PSH04-20]	HA722612B	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-17A Recombinant Rabbit mAb [PS01-24] - BSA and Azide free (Detector)	HA721233	Human	ELISA(Det)
	Biotin Conjugated Human IL-17A Recombinant Rabbit mAb [PS01-24] - Detector	HA721631	Human	ELISA(Det)
	HRP Conjugated Human IL-17A Recombinant Rabbit mAb [PS01-24]	HA721296	Human	ELISA(Det)
	Human IL-17A Recombinant Rabbit mAb [PS01-20] - BSA and Azide free (Capture)	HA721232	Human	ELISA(Cap)
Relative expression	Human IL-18 Recombinant Rabbit mAb [PSH02-09] - BSA and Azide free	HA750779	Human	WB
	Human IL-18 Recombinant Rabbit mAb [PSH02-08] - BSA and Azide free	HA750778	Human	WB,IHC-P,IF-Tissue
	Biotin Conjugated Human IL-18 Recombinant Rabbit mAb [PSH01-70]	HA722644B	Human	ELISA(Det)
	Human IL-18 Recombinant Rabbit mAb [PSH01-71] - BSA and Azide free (Capture)	HA721722	Human	ELISA(Cap)
Relative expression	Human IL-18 Recombinant Rabbit mAb [PSH02-09]	HA721768	Human	WB
	Human IL-18 Recombinant Rabbit mAb [PSH01-70] - BSA and Azide free (Detector)	HA721723	Human	ELISA(Det)
	Human IL-18 Recombinant Rabbit mAb [PSH02-08]	HA721767	Human	WB,IHC-P,IF-Tissue,IP
	Human IL-18BP Recombinant Rabbit mAb [PSH10-72] - BSA and Azide free (Detector)	HA723236	Human	ELISA(Det)
	Biotin Conjugated Human IL-18BP Recombinant Rabbit mAb [PSH10-72] - Detector	HA723237B	Human	ELISA(Det),ELISA
	Human IL-18BP Recombinant Rabbit mAb [PSH10-71] - BSA and Azide free (Capture)	HA723235	Human	ELISA(Cap)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-87] - BSA and Azide free (Capture)	HA723883	Human	ELISA(Cap)
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - Detector	HA723882B	Human	ELISA(Det),ELISA
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39]	HA721712	Human	WB
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free	HA750744	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - BSA and Azide free (Detector)	HA723882	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free	HA750745	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-85] - BSA and Azide free (Capture)	HA723881	Human	ELISA(Cap)
	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free (Detector)	HA721279	Human	ELISA(Det)
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA722649B	Human	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free (Capture)	HA721278	Human	ELISA(Cap)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA721713	Human	WB
	Human IL-1RA Recombinant Rabbit mAb [PSH07-03] - BSA and Azide free (Capture)	HA722788	Human	ELISA(Cap)
	Human IL-1RA Recombinant Rabbit mAb [PSH07-04] - BSA and Azide free (Detector)	HA722789	Human	ELISA(Det)
	Biotin Conjugated Human IL-1RA Recombinant Rabbit mAb [PSH07-04]	HA722790B	Human	ELISA(Det)
	Human IL-2 Recombinant Rabbit mAb [PSH01-42] - BSA and Azide free (Capture)	HA721680	Human	ELISA(Cap)
	Human IL-2 Mouse mAb [A8B7] - BSA and Azide free (Detector)	HA601006	Human	WB,ELISA(Det)
	Biotin Conjugated Human IL-2 Mouse mAb [A8B7]	HA601361B	Human	ELISA(Det)
	Human IL-21 Recombinant Mouse mAb [PSH02-48] - BSA and Azide free (Capture)	HA601195	Human	ELISA(Cap)
	Human IL-21 Recombinant Mouse mAb [PSH02-49] - BSA and Azide free (Detector)	HA601196	Human	ELISA(Det)
	Biotin Conjugated Human IL-21 Recombinant Mouse mAb [PSH02-49]	HA601358B	Human	ELISA(Det)
	Biotin Conjugated Human IL-22 Recombinant Rabbit mAb [PSH03-12] - Detector	HA722630B	Human	ELISA(Det),ELISA
	Human IL-22 Recombinant Rabbit mAb [PSH03-11] - BSA and Azide free (Capture)	HA721939	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-22 Recombinant Rabbit mAb [PSH03-12] - BSA and Azide free (Detector)	HA721940	Human	ELISA(Det)
	Human IL-23A Recombinant Rabbit mAb [PSH03-03] - BSA and Azide free (Capture)	HA721933	Human	ELISA(Cap)
	Human IL-27 Recombinant Rabbit mAb [PSH03-81] - BSA and Azide free (Capture)	HA722048	Human	ELISA(Cap)
	Human IL-27 Recombinant Rabbit mAb [PSH03-82] - BSA and Azide free (Detector)	HA722049	Human	ELISA(Det)
	Biotin Conjugated Human IL-27 Recombinant Rabbit mAb [PSH03-82] - Detector	HA722622B	Human	ELISA(Det)
	Human IL-27 Recombinant Rabbit mAb [PSH03-83] - BSA and Azide free (Detector)	HA722047	Human	ELISA(Det)
	Biotin Conjugated Human IL-27 Recombinant Rabbit mAb [PSH03-83] - Detector	HA722621B	Human	ELISA(Det)
	Human IL-27 Recombinant Rabbit mAb [PSH03-82] - BSA and Azide free (Capture)	HA722046	Human	ELISA(Cap)
	Biotin Conjugated Human IL-27 p28 Recombinant Rabbit mAb [PSH03-81] - Detector	HA722620B	Human	ELISA(Det)
	Human IL-27 p28 Recombinant Rabbit mAb [PSH03-80] - BSA and Azide free (Capture)	HA722044	Human	ELISA(Cap)
	Human IL-27 p28 Recombinant Rabbit mAb [PSH03-81] - BSA and Azide free (Detector)	HA722045	Human	ELISA(Det)
	Biotin Conjugated Human IL-2 Receptor alpha Recombinant Rabbit mAb [PS01-03]	HA722650B	Human	ELISA(Det)
	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PSH05-66] - BSA and Azide free (Detector)	HA722386	Human	ELISA(Det)
	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PS01-03] - BSA and Azide free (Detector)	HA721281	Human	ELISA(Det)
	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PSH05-65] - BSA and Azide free (Capture)	HA722385	Human	ELISA(Cap)
	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PS00-52] - BSA and Azide free (Capture)	HA721280	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human IL-2 Receptor alpha Recombinant Rabbit mAb [PSH05-66]	HA722594B	Human	ELISA(Det)
	Biotin Conjugated Human IL-3 Mouse mAb [PSH0-98]	HA601359B	Human	ELISA(Det)
	Human IL-3 Mouse mAb [PSH0-98] - BSA and Azide free (Detector)	HA601155	Human	ELISA(Det)
	Human IL-3 Mouse mAb [PSH0-97] - BSA and Azide free (Capture)	HA601154	Human	ELISA(Cap)
	Human IL-3 Mouse mAb [PSH0-96] - BSA and Azide free (Capture)	HA601153	Human	ELISA(Cap)
	Biotin Conjugated Human IL-31 Recombinant Rabbit mAb [PSH03-02] - Detector	HA722631B	Human	ELISA(Det),ELISA
	Human IL-31 Recombinant Rabbit mAb [PSH03-02] - BSA and Azide free (Detector)	HA721932	Human	ELISA(Det)
	Human IL-31 Recombinant Rabbit mAb [PSH03-01] - BSA and Azide free (Capture)	HA721931	Human	ELISA(Cap)
	Human IL-33 Recombinant Rabbit mAb [PSH03-06] - BSA and Azide free (Detector)	HA721936	Human	ELISA(Det)
	Biotin Conjugated Human IL-33 Recombinant Rabbit mAb [PSH03-06]	HA722628B	Human	ELISA(Det)
	Human IL-33 Recombinant Rabbit mAb [PSH03-05] - BSA and Azide free (Capture)	HA721935	Human	ELISA(Cap)
	Human IL3RA/CD123 Recombinant Rabbit mAb [PSH11-57] - BSA and Azide free (Capture)	HA723351	Human	ELISA(Cap)
	Biotin Conjugated Human IL3RA/CD123 Recombinant Rabbit mAb [PSH11-58] - Detector	HA723353B	Human	ELISA(Det),ELISA
	Human IL3RA/CD123 Recombinant Rabbit mAb [PSH11-58] - BSA and Azide free (Detector)	HA723352	Human	ELISA(Det)
	Human IL-4 Recombinant Rabbit mAb [PS00-67] - BSA and Azide free (Capture)	HA721230	Human	ELISA(Cap)
	Human IL-4 Recombinant Rabbit mAb [PS00-66] - BSA and Azide free (Detector)	HA721231	Human	ELISA(Det)
	HRP Conjugated Human IL-4 Recombinant Rabbit mAb [PS00-66] - Detector	HA721295	Human	ELISA(Det)
	Biotin Conjugated Human IL-4 Recombinant Rabbit mAb [PS00-66] - Detector	HA721368	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human IL-4R Recombinant Rabbit mAb [PSH11-46] - Detector	HA723338B	Human	ELISA(Det),ELISA
	Human IL-4R Recombinant Rabbit mAb [PSH11-46] - BSA and Azide free (Detector)	HA723337	Human	ELISA(Det)
	Human IL-4R Recombinant Rabbit mAb [PSH11-45] - BSA and Azide free (Capture)	HA723336	Human	ELISA(Cap)
	Biotin Conjugated Human IL-5 Recombinant Rabbit mAb [PSH04-65] - Detector	HA722609B	Human	ELISA(Det)
	Human IL-5 Recombinant Rabbit mAb [PSH04-65] - BSA and Azide free (Detector)	HA722165	Human	ELISA(Det)
	Human IL-5 Recombinant Rabbit mAb [PSH04-64] - BSA and Azide free (Capture)	HA722164	Human	ELISA(Cap)
	HRP Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA721293	Human	ELISA(Det)
	Biotin Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA722233B	Human	ELISA(Det)
	Human IL-6 Recombinant Rabbit mAb [PS00-21] - BSA and Azide free (Capture)	HA721167	Human	ELISA(Cap)
	Human IL-6 Mouse mAb [A6D1]	HA601051	Human	WB,ELISA
	Human IL-6 Recombinant Rabbit mAb [PS00-15] - BSA and Azide free (Detector)	HA721168	Human	ELISA(Det)
	Human IL-6 Mouse mAb [A6D2]	HA601050	Human	WB,ELISA
	Human IL-7 Recombinant Rabbit mAb [PSH03-10] - BSA and Azide free (Detector)	HA721938	Human	ELISA(Det)
	Biotin Conjugated Human IL-7 Recombinant Rabbit mAb [PSH03-10]	HA722632B	Human	ELISA(Det)
	Human IL-7 Recombinant Rabbit mAb [PSH03-09] - BSA and Azide free (Capture)	HA721937	Human	ELISA(Cap)
	Human IL-8 Recombinant Rabbit mAb [PS00-83] - BSA and Azide free (Capture)	HA722098	Human	ELISA(Cap)
	Biotin Conjugated Human IL-8 Recombinant Rabbit mAb [PS01-22]	HA722611B	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-8 Recombinant Rabbit mAb [PS01-22] - BSA and Azide free (Detector)	HA722099	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PS01-75] - BSA and Azide free (Capture)	HA721574	Human	ELISA(Cap)
	Human Interferon gamma Recombinant Rabbit mAb [PS01-76] - BSA and Azide free (Detector)	HA721575	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-45] - BSA and Azide free (Capture)	HA721683	Human	ELISA(Cap)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-47] - BSA and Azide free (Detector)	HA721685	Human	ELISA(Det)
	Biotin Conjugated Human Interferon gamma Recombinant Rabbit mAb [PSH01-47]	HA722645B	Human	ELISA(Det)
	Biotin Conjugated Human Interferon gamma Recombinant Rabbit mAb [PS01-76] - Detector	HA721576	Human	ELISA(Det),ELISA
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-46] - BSA and Azide free (Capture)	HA721684	Human	ELISA(Cap)
	Human LAP3 Recombinant Rabbit mAb [PSH01-34] - BSA and Azide free (Detector)	HA723285	Human, Mouse, Rat	ELISA(Det)
	Human LAP3 Recombinant Rabbit mAb [PSH01-26] - BSA and Azide free (Capture)	HA723284	Human, Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Human LAP3 Recombinant Rabbit mAb [PSH01-34] - Detector	HA723286B	Human, Mouse, Rat	ELISA(Det),ELISA
	Biotin Conjugated Human Leptin Recombinant Rabbit mAb [PSH12-07] - Detector	HA723423B	Human	ELISA(Det),ELISA
	Human Leptin Recombinant Rabbit mAb [PSH12-07] - BSA and Azide free (Detector)	HA723422	Human	ELISA(Det)
	Human Leptin Recombinant Rabbit mAb [PSH12-06] - BSA and Azide free (Capture)	HA723421	Human	ELISA(Cap)
	Biotin Conjugated Human M-CSF Recombinant Rabbit mAb [PSH11-38] - Detector	HA723326B	Human	ELISA(Det),ELISA
	Human M-CSF Recombinant Rabbit mAb [PSH11-36] - BSA and Azide free (Capture)	HA723323	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human M-CSF Recombinant Rabbit mAb [PSH11-38] - BSA and Azide free (Detector)	HA723325	Human	ELISA(Det)
	Human M-CSF Recombinant Rabbit mAb [PSH11-37] - BSA and Azide free (Capture)	HA723324	Human	ELISA(Cap)
	Human/Mouse VEGFC Recombinant Rabbit mAb [PSH10-76] - BSA and Azide free (Detector)	HA723242	Human, Mouse	ELISA(Det)
	Human/Mouse VEGFC Recombinant Rabbit mAb [PSH10-75] - BSA and Azide free (Capture)	HA723241	Human, Mouse	ELISA(Cap)
	Biotin Conjugated Human/Mouse VEGFC Recombinant Rabbit mAb [PSH10-76] - Detector	HA723243B	Human, Mouse	ELISA(Det),ELISA
	Human Prolactin Recombinant Rabbit mAb [PSH03-14] - BSA and Azide free (Detector)	HA721942	Human	ELISA(Det)
	Human Prolactin Recombinant Rabbit mAb [PSH03-13] - BSA and Azide free (Capture)	HA721941	Human	ELISA(Cap)
	Biotin Conjugated Human Prolactin Recombinant Rabbit mAb [PSH03-14] - Detector	HA722627B	Human	ELISA(Det)
	Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - BSA and Azide free (Detector)	HA723159	Human	ELISA(Det)
	Biotin Conjugated Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - BSA and Azide free (Detector)	HA723160B	Human	ELISA(Det),ELISA
	Human sTNF RI Recombinant Rabbit mAb [PSH09-89] - BSA and Azide free (Capture)	HA723158	Human	ELISA(Cap)
	HRP Conjugated Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - Detector	HA723374H	Human	ELISA(Det)
	Biotin Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721559	Human	ELISA(Det)
	HRP Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721294	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40]	HA721286	Human	IHC-P, mIHC
	Human TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free (Detector)	HA721266	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb	IRS040	Human	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free (Capture)	HA721265	Human	ELISA(Cap)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free	HA750593	Human	IHC-P
	Biotin Conjugated Human TNF-beta Recombinant Rabbit mAb [PS01-19]	HA722651B	Human	ELISA(Det)
	Human TNF-beta Recombinant Rabbit mAb [PS01-19] - BSA and Azide free (Detector)	HA721243	Human	ELISA(Det)
	Human TNF-beta Recombinant Rabbit mAb [PS01-05] - BSA and Azide free (Capture)	HA721242	Human	ELISA(Cap)
	Biotin Conjugated Human TNF Receptor II Recombinant Rabbit mAb [PSH11-93] - Detector	HA723404B	Human	ELISA(Det),ELISA
	Human TNF Receptor II Recombinant Rabbit mAb [PSH11-94] - BSA and Azide free (Capture)	HA723405	Human	ELISA(Cap)
	Human TNF Receptor II Recombinant Rabbit mAb [PSH11-93] - BSA and Azide free (Detector)	HA723403	Human	ELISA(Det)
	Human TNF Receptor II Recombinant Rabbit mAb [PSH11-92] - BSA and Azide free (Capture)	HA723402	Human	ELISA(Cap)
	Human VEGFA Recombinant Rabbit mAb [PSH12-91] - BSA and Azide free (Capture)	HA723493	Human	ELISA(Cap)
	Human VEGFA Recombinant Rabbit mAb [PSH12-93] - BSA and Azide free (Detector)	HA723496	Human	ELISA(Det)
	Human VEGFA Recombinant Rabbit mAb [PSH12-92] - BSA and Azide free (Detector)	HA723494	Human	ELISA(Det)
	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-93] - Detector	HA723497B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-92] - Detector	HA723495B	Human	ELISA(Det),ELISA
Relative expression, KD	IFNGR1 Recombinant Rabbit mAb [PSH06-11]	HA722508	Human	WB,IF-Cell
Relative expression, KD	IFNGR1 Recombinant Rabbit mAb [PSH06-11] - BSA and Azide free	HA751049	Human	WB,IF-Cell
KD, Relative expression	IFNGR1 Rabbit pAb	ER63732	Human	WB,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression, KD	IFNGR1 Recombinant Rabbit mAb [JE40-50]	HA721952	Human, Mouse, Rat	WB,IHC-P
Cell treatment	IL-10 Recombinant Rabbit mAb [JE59-71]	HA722032	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	IL-10 Rabbit pAb	ER1911-19	Human, Mouse	WB,IHC-P,IF-Tissue
	IL-10 Rabbit pAb	ER65464	Human, Mouse, Rat	WB,IHC-P
	IL-12 Rabbit pAb	ER1911-22	Rat	WB,IHC-P,IF-Tissue
	IL-12A Recombinant Rabbit mAb [PSH16-10] - BSA and Azide free	HA723815	Human	ELISA
	IL12B Rabbit pAb	ER1803-31	Human, Mouse	WB
	IL-13 Rabbit pAb	ER1911-23	Human	IHC-P,IF-Tissue
Relative expression	IL-13R alpha 2 Recombinant Rabbit mAb [PSH03-85]	HA722051	Human	WB,IF-Cell,IHC-P,FC
Relative expression	IL-13R alpha 2 Recombinant Rabbit mAb [PSH03-85] - BSA and Azide free	HA750910	Human	WB,IF-Cell,IHC-P,FC
	IL-17 Rabbit pAb	ER1706-91	Human	WB,IHC-P,FC
	IL17 Mouse mAb [4G2]	RT1326	Human, Mouse, Rat	WB,IP,IF
	IL-17A Rabbit pAb	ER1902-37	Human, Mouse, Rat	WB,IHC-P,FC,mIHC
	IL-17A Rabbit pAb	IRS056	Mouse	mIHC
	IL-18 Recombinant Rabbit mAb [PS01-06]	HA721537	Human	WB,FC
Relative expression	IL-18 Recombinant Rabbit mAb [PSH18-15]	HA723978	Human, Mouse, Rat	WB,IHC-P,IP
	IL-18 Recombinant Rabbit mAb [PS01-06] - BSA and Azide free	HA750677	Human	WB,FC
	IL-18 Recombinant Rabbit mAb [PS00-41]	HA721536	Human	WB,IHC-P
	IL-18 Recombinant Rabbit mAb [PS00-41] - BSA and Azide free	HA750676	Human	WB,IHC-P
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R]	HA601036	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R] - BSA and Azide free	HA610004	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Mouse mAb [A7A10]	HA601035	Human	WB,ELISA,IF-Cell
Cell treatment	IL-1 beta Mouse mAb [A7F9]	HA601002	Human, Mouse, Rat	WB,ELISA,IF-Cell
Cell treatment	IL-1 beta Recombinant Rabbit mAb [JJ087-3]	ET1701-39	Human	WB,IF-Cell,IF-Tissue,IHC-P
	IL-1 beta Recombinant Rabbit mAb [PSH18-00]	HA723965	Mouse, Rat	WB
	IL-1RA Rabbit pAb	ER1911-28	Human	IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	IL-1RA Recombinant Rabbit mAb [JE07-58]	HA722446	Human, Mouse, Rat	WB
	IL-2 Rabbit pAb	ER63748	Human, Rat	WB,IHC-P
	IL-2 Mouse mAb [A8B8]	HA601015	Human	WB,ELISA
	IL-2 Mouse mAb [A8B9]	HA601016	Human	WB,ELISA
	IL-21 Rabbit pAb	ER63752	Human, Mouse	WB
	IL-22 Rabbit pAb	ER1803-85	Human, Mouse, Rat	ELISA,WB,IHC-P
	IL22 Rabbit pAb	ER1911-31	Human, Mouse	WB,IHC-P,IF-Tissue,FC
	IL-23 Recombinant Rabbit mAb [JE38-99]	HA721037	Human	WB,IHC-P,IF-Tissue
	IL-27A Rabbit pAb	ER65802	Human, Mouse, Rat	IHC-P
	IL-2 receptor alpha Rabbit pAb	HA500076	Human, Rat	WB
	IL-31 Mouse mAb [B7-B0]	M1305-3	Human	WB,ELISA
	IL-31 Mouse mAb [8-C3]	EM1701-79	Human	WB,IHC-P,IF-Cell
Relative expression, Cell treatment	IL-33 Recombinant Rabbit mAb [PSH17-37] - BSA and Azide free	HA751625	Mouse, Rat	WB,IHC-P
Relative expression, Cell treatment	IL-33 Recombinant Rabbit mAb [PSH17-37]	HA723920	Mouse, Rat	WB,IHC-P
	IL33 Rabbit pAb	ER1911-36	Human, Rat	WB,IHC-P,IF-Tissue
	IL3RA Mouse mAb [C5-E11]	EM1709-19	Human	WB,IF-Cell,IHC-P,FC
	IL3RA / CD123 Rabbit pAb	HA500426	Human	WB
	IL3RA / CD123 Rabbit pAb	HA500291	Human	WB,IF-Cell,IHC-P,FC
Relative expression	IL3RA / CD123 Recombinant Rabbit mAb [PSH09-18]	HA723076	Human	WB,IF-Cell,IHC-P,FC,IF-Tissue
	IL3RA / CD123 Recombinant mAb	IRS104RB	Human	mIHC
	IL-4 Mouse mAb [A7D7]	HA601034	Human	ELISA
	IL-4 Recombinant Rabbit mAb [JA86-10]	ET1704-28	Human	WB,IF-Cell,IF-Tissue,IHC-P
	IL-4 Rabbit pAb	ER1706-55	Human	WB,IF-Cell,IHC-P
	IL-4R Rabbit pAb	HA500453	Human	WB,IF-Cell
	IL-4R Rabbit pAb	HA500454	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	IL-4R Recombinant Rabbit mAb [PSH10-86]	HA723253	Human, Mouse, Rat	WB,IP,IF-Cell
	IL-4R Recombinant Rabbit mAb [PSH10-86] - BSA and Azide free	HA751378	Human, Mouse, Rat	WB,IP,IF-Cell
	IL-4R Rabbit pAb	ER1911-37	Human	WB
	IL-5 Rabbit pAb	HA500081	Human	WB,IF-Cell,ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	IL-6 Recombinant Rabbit mAb	IRS160RB	Human	mIHC
	IL-6 Rabbit pAb	R1203-2	Mouse	WB,ELISA
	IL-6 Mouse mAb [1-6]	EM1701-45	Mouse, Rat	WB,IF-Cell
	IL-6 Rabbit pAb	R1412-2	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	IL-6 Mouse mAb	EM30301	Human	WB
	IL-6 Mouse mAb [3-A4]	M1303-3	Mouse	WB,IF-Cell,FC,ELISA
	IL-7 Rabbit pAb	HA500171	Human	WB,IF-Cell,IHC-P
	IL-7 Mouse mAb [A8B11]	HA601048	Human	WB
	IL7 Mouse mAb [AH31-34]	EM170412	Human, Rat	ELISA,WB,ICC,IHC-P
Cell treatment	IL-8 Rabbit pAb	ER1901-61	Human	WB,IHC-P,IF-Cell,FC
Cell treatment	IL-8 Rabbit pAb	R1511-15	Human	WB,IHC-P,IF-Cell,FC
	IL-8 Rabbit pAb	ER1706-67	Human	ELISA,IF-Cell,IHC-P
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-96]	HA723263	Human	WB,FC
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-97]	HA723264	Human	WB
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-96] - BSA and Azide free	HA751386	Human	WB,FC
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-97] - BSA and Azide free	HA751387	Human	WB
	IL-9 Rabbit pAb	ER1911-41	Mouse, Rat	WB,IHC-P,IF-Tissue
	Insulin Mouse mAb [A6-6]	EM80714	Human, Mouse, Rat	IHC-P,ELISA,IF-Tissue,IHC-Fr
	Insulin Recombinant Rabbit mAb	IRS161RB	Human	mIHC
	Insulin Recombinant Rabbit mAb [PSH04-15] - BSA and Azide free (Capture)	HA722096	Human, Mouse, Rat	ELISA(Cap)
	Insulin Recombinant Rabbit mAb	IRS233RB	Mouse	mIHC
	Biotin Conjugated Insulin Recombinant Rabbit mAb [PSH04-16] - Detector	HA722610B	Human, Mouse, Rat	ELISA(Det),ELISA
	Insulin Recombinant Rabbit mAb [PSH04-16] - BSA and Azide free (Detector)	HA722097	Human, Mouse, Rat	ELISA(Det)
Relative expression	Insulin Recombinant Rabbit mAb [SA0410] - BSA and Azide free	HA750004	Human, Mouse, Rat	IF-Cell,IF-Tissue,IHC-P,IHC-Fr,WB
	Insulin Mouse mAb [5E1]	EM1701-33	Human	IHC-P,ELISA
Relative expression	Insulin Recombinant Rabbit mAb [SA0410]	ET1601-12	Human, Mouse, Rat	IF-Cell,IF-Tissue,IHC-P,mIHC,IHC-Fr,WB
	Insulin B Chain Rabbit pAb	0807-11	Human, Mouse, Rat	IHC-P,ELISA,IF-Tissue
	Interferon alpha Rabbit pAb	ER1802-4	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Interferon gamma Recombinant Rabbit mAb [JM10-10]	ET1703-17	Human	WB,FC
	Interferon gamma Recombinant Rabbit mAb [JM10-10] - BSA and Azide free	HA750373	Human	WB,FC
KD	LAP3 Recombinant Rabbit mAb [PSH01-26]	HA721650	Human, Mouse, Rat	WB,IHC-P
KD	LAP3 Recombinant Rabbit mAb [PSH01-34] - BSA and Azide free	HA750712	Human, Mouse, Rat	WB,IHC-P
KD	LAP3 Recombinant Rabbit mAb [PSH01-26] - BSA and Azide free	HA750704	Human, Mouse, Rat	WB,IHC-P
KD	LAP3 Recombinant Rabbit mAb [PSH01-34]	HA721658	Human, Mouse, Rat	WB,IHC-P
	Leptin Rabbit pAb	ER1912-27	Human, Mouse, Rat	IHC-P,IF-Tissue
	LIM1 Rabbit pAb	ER1902-83	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	LT-alpha Mouse mAb [6-C10-10]	M0807-15	Human, Rat	WB
	LT-alpha Rabbit pAb	0807-5	Human, Rat	WB
	LT-alpha Mouse mAb [C2-F2]	M0807-12	Human	WB
Relative expression	LT-alpha (C-term) Rabbit pAb	0807-6	Human, Mouse, Rat	WB,IF-Cell
	M-CSF Recombinant Rabbit mAb [SU0413]	ET1609-1	Human, Mouse	WB,IHC-P,IF-Cell,IP
	Mouse CCL3 Recombinant Rabbit mAb [PSH08-39] - BSA and Azide free (Capture)	HA722982	Mouse	ELISA(Cap)
	Mouse CCL3 Recombinant Rabbit mAb [PSH08-40] - BSA and Azide free (Detector)	HA722983	Mouse	ELISA(Det)
	Biotin Conjugated Mouse CCL3 Recombinant Rabbit mAb [PSH08-40] - Detector	HA722984B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse IFN beta Recombinant Rabbit mAb [PSH04-73]	HA722605B	Mouse	ELISA(Det),ELISA
	Mouse IFN beta Recombinant Rabbit mAb [PSH04-72] - BSA and Azide free (Capture)	HA722172	Mouse	ELISA(Cap)
	Mouse IFN beta Recombinant Rabbit mAb [PSH04-73] - BSA and Azide free (Detector)	HA722173	Mouse	ELISA(Det)
	Mouse IL-10 Recombinant Rabbit mAb [PS01-93] - BSA and Azide free (Detector)	HA721518	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-10 Recombinant Rabbit mAb [PS01-93] - Detector	HA721519	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IL-10 Recombinant Rabbit mAb [PS01-85] - BSA and Azide free (Capture)	HA721517	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL12 p40 Recombinant Rabbit mAb [PSH04-88]	HA722603B	Mouse	ELISA(Det)
	Mouse IL12 p40 Recombinant Rabbit mAb [PSH04-87] - BSA and Azide free (Capture)	HA722193	Mouse	ELISA(Cap)
	Mouse IL12 p40 Recombinant Rabbit mAb [PSH04-88] - BSA and Azide free (Detector)	HA722194	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-13 Recombinant Rabbit mAb [PSH02-05]	HA722641B	Mouse	ELISA(Det)
	Mouse IL-13 Recombinant Rabbit mAb [PSH02-04] - BSA and Azide free (Capture)	HA721784	Mouse	ELISA(Cap)
	Mouse IL-13 Recombinant Rabbit mAb [PSH02-05] - BSA and Azide free (Detector)	HA721785	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-17A Recombinant Rabbit mAb [PS01-92] - Detector	HA721608	Mouse	ELISA(Det),ELISA
	Mouse IL-17A Recombinant Rabbit mAb [PS01-91] - BSA and Azide free (Capture)	HA721606	Mouse	ELISA(Cap)
	Mouse IL-17A Recombinant Rabbit mAb [PS01-92] - BSA and Azide free (Detector)	HA721607	Mouse	ELISA(Det)
	Mouse IL-17A Recombinant Rabbit mAb [PSH05-35] - BSA and Azide free (Capture)	HA722264	Mouse, Rat	ELISA(Cap)
	Mouse IL-17A Recombinant Rabbit mAb [PSH05-36] - BSA and Azide free (Detector)	HA722265	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse IL-17A Recombinant Rabbit mAb [PSH05-36]	HA722600B	Mouse, Rat	ELISA(Det),ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH09-58] - BSA and Azide free	HA723120	Mouse	ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67] - BSA and Azide free (Detector)	HA722167	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-66] - BSA and Azide free (Capture)	HA722166	Mouse	ELISA(Cap)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH09-59] - BSA and Azide free	HA723121	Mouse	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - Detector	HA722625B	Mouse	ELISA(Det),ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-31] - BSA and Azide free (Capture)	HA721976	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67]	HA722608B	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - BSA and Azide free (Detector)	HA721977	Mouse	ELISA(Det)
	Mouse IL-1RA Recombinant Rabbit mAb [PSH15-00] - BSA and Azide free (Capture)	HA723700	Mouse	ELISA(Cap)
	Mouse IL-1RA Recombinant Rabbit mAb [PSH14-99] - BSA and Azide free (Detector)	HA723698	Mouse	ELISA(Det)
	Mouse IL-1RA Recombinant Rabbit mAb [PSH14-98] - BSA and Azide free (Capture)	HA723697	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-1RA Recombinant Rabbit mAb [PSH14-99] - Detector	HA723699B	Mouse	ELISA(Det)
	Mouse IL-2 Recombinant Rabbit mAb [PSH04-27] - BSA and Azide free (Capture)	HA722110	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-2 Recombinant Rabbit mAb [PSH04-28]	HA722617B	Mouse	ELISA(Det)
	Mouse IL-2 Recombinant Rabbit mAb [PSH04-28] - BSA and Azide free (Detector)	HA722111	Mouse	ELISA(Det)
	Mouse IL-21 Recombinant Rabbit mAb [PSH14-56] - BSA and Azide free (Detector)	HA723650	Mouse	ELISA(Det)
	Mouse IL-21 Recombinant Rabbit mAb [PSH14-57] - BSA and Azide free (Detector)	HA723652	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-21 Recombinant Rabbit mAb [PSH14-56] - Detector	HA723651B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse IL-21 Recombinant Rabbit mAb [PSH14-57] - Detector	HA723653B	Mouse	ELISA(Det),ELISA
	Mouse IL-21 Recombinant Rabbit mAb [PSH04-75] - BSA and Azide free (Detector)	HA722175	Mouse	ELISA(Det)
	Mouse IL-21 Recombinant Rabbit mAb [PSH04-74] - BSA and Azide free (Capture)	HA722174	Mouse	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IL-21 Recombinant Rabbit mAb [PSH04-75] - BSA and Azide free (Capture)	HA723649	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-21 Recombinant Rabbit mAb [PSH04-75] - Detector	HA723420B	Mouse	ELISA(Det)
	Mouse IL-22 Recombinant Rabbit mAb [PSH13-05] - BSA and Azide free (Detector)	HA725032	Mouse	ELISA(Det)
	Mouse IL-22 Recombinant Rabbit mAb [PSH13-04] - BSA and Azide free (Capture)	HA725031	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-22 Recombinant Rabbit mAb [PSH13-05] - Detector	HA725033B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse IL-33 Recombinant Rabbit mAb [PSH16-62] - Detector	HA723864B	Mouse	ELISA(Det),ELISA
	Mouse IL-33 Recombinant Rabbit mAb [PSH08-42] - BSA and Azide free (Detector)	HA722986	Mouse	ELISA(Det)
	Mouse IL-33 Recombinant Rabbit mAb [PSH16-62] - BSA and Azide free (Detector)	HA723864	Mouse	ELISA(Det)
	Mouse IL-33 Recombinant Rabbit mAb [PSH16-63] - BSA and Azide free (Capture)	HA723865	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-33 Recombinant Rabbit mAb [PSH08-42] - Detector	HA722987B	Mouse	ELISA(Det),ELISA
	Mouse IL-33 Recombinant Rabbit mAb [PSH08-41] - BSA and Azide free (Capture)	HA722985	Mouse	ELISA(Cap)
	Mouse IL-33 Recombinant Rabbit mAb [PSH16-61] - BSA and Azide free (Capture)	HA723863	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-4 Recombinant Rabbit mAb [PSH04-69]	HA722607B	Mouse	ELISA(Det)
	Mouse IL-4 Recombinant Rabbit mAb [PSH05-34] - BSA and Azide free (Detector)	HA722263	Mouse	ELISA(Det)
	Mouse IL-4 Recombinant Rabbit mAb [PSH04-68] - BSA and Azide free (Capture)	HA722168	Mouse	ELISA(Cap)
	Mouse IL-4 Recombinant Rabbit mAb [PSH04-69] - BSA and Azide free (Detector)	HA722169	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IL-4 Recombinant Rabbit mAb [PSH05-33] - BSA and Azide free (Capture)	HA722262	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-4 Recombinant Rabbit mAb [PSH05-34]	HA722597B	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH14-47] - BSA and Azide free (Capture)	HA723642	Mouse	ELISA(Cap)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-58] - BSA and Azide free (Capture)	HA721411	Mouse	ELISA(Cap)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-60] - BSA and Azide free (Detector)	HA721413	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PSH0-60] - Detector	HA722647B	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PSH0-59] - Detector	HA722646B	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PS01-74] - BSA and Azide free (Capture/Detector)	HA721600	Mouse	ELISA(Cap),ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PS01-83] - BSA and Azide free (Detector)	HA721602	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-59] - BSA and Azide free (Detector)	HA721412	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PS01-74] - Detector	HA723643B	Mouse	ELISA(Det)
	HRP Conjugated Mouse IL-6 Mouse mAb [3-A4]	M1303-2	Mouse	WB,ELISA
	Mouse IL-6 Recombinant Rabbit mAb [PS01-83] - BSA and Azide free (Detector)	HA721601	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-9 Recombinant Rabbit mAb [PS01-82] - Detector	HA721605	Mouse	ELISA(Det),ELISA
	Mouse IL-9 Recombinant Rabbit mAb [PS01-81] - BSA and Azide free (Capture)	HA721603	Mouse	ELISA(Cap)
	Mouse IL-9 Recombinant Rabbit mAb [PS01-82] - BSA and Azide free (Detector)	HA721604	Mouse	ELISA(Det)
	Biotin Conjugated Mouse Interferon alpha 1 Recombinant Rabbit mAb [PSH14-76] - Detector	HA723671B	Mouse	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse Interferon alpha 1 Recombinant Rabbit mAb [PSH14-76] - BSA and Azide free (Detector)	HA723670	Mouse	ELISA(Det)
	Mouse Interferon alpha 1 Recombinant Rabbit mAb [PSH14-77] - BSA and Azide free (Capture)	HA723672	Mouse	ELISA(Cap)
	Mouse Interferon alpha 1 Recombinant Rabbit mAb [PSH14-75] - BSA and Azide free (Capture)	HA723669	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse M-CSF/CSF1 Recombinant Rabbit mAb [PSH15-85] - Detector	HA723794B	Mouse	ELISA(Det),ELISA
	Mouse M-CSF/CSF1 Recombinant Rabbit mAb [PSH15-84] - BSA and Azide free (Capture)	HA723792	Mouse	ELISA(Cap)
	Mouse M-CSF/CSF1 Recombinant Rabbit mAb [PSH15-85] - BSA and Azide free (Detector)	HA723793	Mouse	ELISA(Det)
	Biotin Conjugated Mouse Prolactin Recombinant Rabbit mAb [PSH17-91] - Detector	HA723956B	Mouse	ELISA(Det),ELISA
	Mouse Prolactin Recombinant Rabbit mAb [PSH17-90] - BSA and Azide free (Capture)	HA723955	Mouse	ELISA(Cap)
	Mouse Prolactin Recombinant Rabbit mAb [PSH17-91] - BSA and Azide free (Detector)	HA723956	Mouse	ELISA(Det)
	Mouse/Rat Leptin Recombinant Rabbit mAb [PSH08-95] - BSA and Azide free (Detector)	HA723047	Mouse, Rat	ELISA(Det)
	Mouse/Rat Leptin Recombinant Rabbit mAb [PSH08-94] - BSA and Azide free (Capture)	HA723046	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat Leptin Recombinant Rabbit mAb [PSH08-95] - Detector	HA723048B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-29] - BSA and Azide free (Capture)	HA721974	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30]	HA722626B	Mouse, Rat	ELISA(Det)
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30] - BSA and Azide free (Detector)	HA721975	Mouse, Rat	ELISA(Det)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-76] - BSA and Azide free (Capture)	HA722040	Mouse	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77] - BSA and Azide free (Detector)	HA722041	Mouse	ELISA(Det)
	Biotin Conjugated Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77]	HA722623B	Mouse	ELISA(Det)
	Biotin Conjugated Mouse VEGF Recombinant Rabbit mAb [PSH02-73]	HA722636B	Mouse	ELISA(Det)
	Mouse VEGF Recombinant Rabbit mAb [PSH02-72] - BSA and Azide free (Capture)	HA721850	Mouse	ELISA(Cap)
	Mouse VEGF Recombinant Rabbit mAb [PSH02-73] - BSA and Azide free (Detector)	HA721851	Mouse	ELISA(Det)
	NGF Rabbit pAb	ER1802-59	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	NGF Recombinant Rabbit mAb [SI79-01]	ET1606-29	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
	NGF Recombinant Rabbit mAb [SI79-01] - BSA and Azide free	HA750099	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
	NGF beta Mouse mAb [1G4]	RT1428	Human, Mouse, Rat	WB,IP,IF
	Ob Rabbit pAb	ER65824	Human, Mouse, Rat	IHC-P
	Pin1 Recombinant Rabbit mAb [PSH10-81]	HA723248	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	Pin1 Recombinant Rabbit mAb [PSH10-81] - BSA and Azide free	HA751373	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	PRMT1 Recombinant Rabbit mAb [JA11-17]	ET1704-19	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Relative expression	Prolactin / PRL Recombinant Rabbit mAb [PSH19-33]	HA724065	Human, Mouse, Rat	IHC-P
Relative expression	Prolactin / PRL Recombinant Rabbit mAb [PSH18-58]	HA724011	Mouse	WB,IF-Tissue,IHC-P
	Prolactin / PRL Rabbit pAb	ER1803-72	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	Prolactin / PRL Recombinant Rabbit mAb [PSH18-57]	HA724010	Human, Mouse, Rat	WB,IF-Tissue,IHC-P
	Rab 2A Rabbit pAb	ER64714	Human, Mouse, Rat	WB,IHC-P
	Rat IFN gamma Recombinant Rabbit mAb [PSH05-69] - BSA and Azide free (Capture)	HA722421	Rat	ELISA(Cap)
	Biotin Conjugated Rat IFN gamma Recombinant Rabbit mAb [PS01-97] - Detector	HA721620	Rat	ELISA(Det)
	Biotin Conjugated Rat IFN gamma Recombinant Rabbit mAb [PSH05-70]	HA722589B	Rat	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Rat IFN gamma Recombinant Rabbit mAb [PS01-97] - BSA and Azide free (Detector)	HA721619	Rat	ELISA(Det)
	Rat IFN gamma Recombinant Rabbit mAb [PS01-96] - BSA and Azide free (Capture)	HA721618	Rat	ELISA(Cap)
	Rat IFN gamma Recombinant Rabbit mAb [PSH05-70] - BSA and Azide free (Detector)	HA722422	Rat	ELISA(Det)
	Rat IL-10 Recombinant Rabbit mAb [PSH02-22] - BSA and Azide free (Detector)	HA721787	Rat	ELISA(Det)
	Biotin Conjugated Rat IL-10 Recombinant Rabbit mAb [PSH02-22]	HA722637B	Rat	ELISA(Det)
	Rat IL-10 Recombinant Rabbit mAb [PSH02-21] - BSA and Azide free (Capture)	HA721786	Rat	ELISA(Cap)
	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60]	HA722634B	Rat	ELISA(Det)
	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-59] - BSA and Azide free (Capture)	HA721836	Rat	ELISA(Cap)
	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60] - BSA and Azide free (Detector)	HA721837	Rat	ELISA(Det)
	Rat IL-1 beta Recombinant Rabbit mAb [PS01-86] - BSA and Azide free (Capture)	HA721609	Rat	ELISA(Cap)
	Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - BSA and Azide free (Detector)	HA721610	Rat	ELISA(Det)
	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - Detector	HA721611	Rat	ELISA(Det),ELISA
	Rat IL-2 Recombinant Rabbit mAb [PS01-88] - BSA and Azide free (Capture)	HA721612	Rat	ELISA(Cap)
	Biotin Conjugated Rat IL-2 Recombinant Rabbit mAb [PSH05-72]	HA722590B	Rat	ELISA(Det)
	Rat IL-2 Recombinant Rabbit mAb [PS01-89] - BSA and Azide free (Detector)	HA721613	Rat	ELISA(Det)
	Rat IL-2 Recombinant Rabbit mAb [PSH05-71] - BSA and Azide free (Capture)	HA722423	Rat	ELISA(Cap)
	Rat IL-2 Recombinant Rabbit mAb [PSH05-72] - BSA and Azide free (Detector)	HA722424	Rat	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Rat IL-2 Recombinant Rabbit mAb [PS01-89] - Detector	HA721614	Rat	ELISA(Det)
	Biotin Conjugated Rat IL-4 Recombinant Rabbit mAb [PSH04-22]	HA722614B	Rat	ELISA(Det)
	Rat IL-4 Recombinant Rabbit mAb [PSH04-22] - BSA and Azide free (Detector)	HA722105	Rat	ELISA(Det)
	Rat IL-4 Recombinant Rabbit mAb [PSH04-21] - BSA and Azide free (Capture)	HA722104	Rat	ELISA(Cap)
	Rat IL-6 Recombinant Rabbit mAb [PS01-94] - BSA and Azide free (Detector)	HA721616	Rat	ELISA(Det)
	Biotin Conjugated Rat IL-6 Recombinant Rabbit mAb [PS01-94] - BSA and Azide free (Detector)	HA721617	Rat	ELISA(Det)
	Rat IL-6 Recombinant Rabbit mAb [PS01-90] - BSA and Azide free (Capture)	HA721615	Rat	ELISA(Cap)
	Rat TNF alpha Recombinant Rabbit mAb [PSH04-70] - BSA and Azide free (Capture)	HA722170	Rat	ELISA(Cap)
	Biotin Conjugated Rat TNF alpha Recombinant Rabbit mAb [PSH04-71]	HA722606B	Rat	ELISA(Det)
	Rat TNF alpha Recombinant Rabbit mAb [PSH04-71] - BSA and Azide free (Detector)	HA722171	Rat	ELISA(Det)
	TBP-1 Rabbit pAb	HA500413	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	TDG Rabbit pAb	ER65129	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KO	TDG Recombinant Rabbit mAb [PSH10-11] - BSA and Azide free	HA751531	Human, Mouse, Rat	WB,IHC-P
KO	TDG Recombinant Rabbit mAb [PSH10-11]	HA723180	Human, Mouse, Rat	WB,IHC-P
	TNF alpha Rabbit pAb	ER1919-22	Human, Mouse	WB
Cell treatment	TNF alpha Rabbit pAb	R1203-1	Mouse	WB,IF-Cell,ELISA
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free	HA750897	Human, Mouse	WB
	TNF alpha Rabbit pAb	ER65189	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07]	HA722022	Human, Mouse	WB
	TNFR1 Mouse mAb [2G2]	RT1624	Human, Mouse, Rat	WB,IP,IF,IHC-P
	TNF Receptor I Rabbit pAb	HA500140	Human	WB,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	TNF Receptor II Recombinant Rabbit mAb [PSH12-12]	HA723428	Human	FC
	TNF Receptor II Recombinant Rabbit mAb [JM113-01]	ET1703-63	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP,IHC-Fr
Relative expression	TNF Receptor II Recombinant Rabbit mAb [PSH12-12] - BSA and Azide free	HA751441	Human	FC
	TNF Receptor II Rabbit pAb	HA500332	Human, Mouse, Rat	IHC-P,FC
	TNF Receptor II Recombinant Rabbit mAb [JM113-01] - BSA and Azide free	HA750386	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP,IHC-Fr
	TRIM29 Recombinant Rabbit mAb [PSH02-99] - BSA and Azide free	HA750854	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	TRIM29 Recombinant Rabbit mAb [PSH02-99]	HA721918	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	UCHL5IP Recombinant Rabbit mAb [JE65-32]	HA721077	Human	WB,IHC-P
KD	VEGF Rabbit pAb	ER30607	Human	WB,IF-Cell,IHC-P,FC
	VEGF Recombinant Rabbit mAb [SP07-01] - BSA and Azide free	HA750074	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	VEGF Recombinant Rabbit mAb [SP07-01]	ET1604-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	VEGFA Recombinant Rabbit mAb [PSH19-17]	HA724051	Human, Mouse, Rat, Monkey	WB
	VEGFC Rabbit pAb	R1307-4	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层