

胶质细胞功能与调控研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	ABCA1 Mouse mAb [4G2]	RT1009	Human, Mouse	WB,IP,IF,IHC-P
	ABCA1 Rabbit pAb	R1510-42	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	ABCA7 Recombinant Rabbit mAb [PSH18-11]	HA723974	Mouse, Rat	WB,IHC-P
	ABCG1 Recombinant Rabbit mAb [SU03-26]	ET1608-20	Human, Rat, Mouse	WB,IHC-P
	ALDH1L1 Recombinant Rabbit mAb [JU53-54]	ET7106-64	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
	ALDH1L1 Recombinant Rabbit mAb [JU53-54] - BSA and Azide free	HA750457	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
KD, Relative expression	Aldolase A/B/C Recombinant Rabbit mAb [PSH01-81] - BSA and Azide free	HA750755	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell
Relative expression	Aldolase A/B/C Recombinant Antibody [PSH01-81] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610342	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Aldolase A/B/C Recombinant Antibody [PSH01-81] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610341	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Aldolase A/B/C Recombinant Antibody [PSH01-81] - Mouse IgG1 (Chimeric)	HA601511	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Aldolase A/B/C Recombinant Antibody [PSH01-81] - Rat IgG1 (Chimeric)	HA601512	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
KD, Relative expression	Aldolase A/B/C Recombinant Rabbit mAb [PSH01-81]	HA721733	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell
	Aldolase C Rabbit pAb	ER65479	Human, Mouse, Rat	WB,IHC-P
Relative expression	Aldolase C Recombinant Rabbit mAb [PSH01-82]	HA721734	Mouse, Rat, Human	WB,IHC-P
Relative expression	Aldolase C Recombinant Rabbit mAb [PSH01-82] - BSA and Azide free	HA750756	Mouse, Rat, Human	WB,IHC-P
	alpha+beta Synuclein Recombinant Rabbit mAb [ST05-21]	ET1609-66	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment, Relative expression	Alpha-Synuclein Recombinant Rabbit mAb [PSH08-83] - BSA and Azide free	HA751248	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue

Cell treatment, Relative expression	Alpha-Synuclein Recombinant Rabbit mAb [PSH08-83]	HA723035	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	Alpha-Synuclein Recombinant Rabbit mAb [JB42-33]	ET7107-31	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Alpha-Synuclein Mouse mAb [2G5]	RT1036	Human	WB,IP,IF
	Amyloid Beta A4 Precursor (APP) Rabbit pAb	1007-5	Human, Mouse, Rat	WB,IF-Cell
	APC Rabbit pAb	ER1802-39	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,IHC-Fr
	APC Recombinant Rabbit mAb [SC61-03]	ET1610-80	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IP,IHC-P,IHC-Fr
Cell treatment, Relative expression	APP/beta Amyloid Recombinant Rabbit mAb [PSH04-90]	HA722207	Human, Mouse	WB,IF-Cell,IHC-P
Relative expression	APP/beta Amyloid Recombinant Rabbit mAb [JE30-54]	HA722139	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
Cell treatment, Relative expression	APP/ β -Amyloid Recombinant Rabbit mAb [PSH04-90] - BSA and Azide free	HA750977	Human, Mouse	WB,IF-Cell,IHC-P
	AP repeat Rabbit pAb	HA500112	Human	WB
Relative expression	Aquaporin 4 Recombinant Rabbit mAb [PSH06-38] - BSA and Azide free	HA751077	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr,IP
Relative expression	Aquaporin 4 Recombinant Antibody [PSH06-38] - Guinea pig IgG2 (Chimeric)	HA601461	Mouse	IHC-Fr,WB
	Aquaporin 4 Recombinant Antibody [PSH06-38] - Chicken IgY (Chimeric)	HA601460	Mouse	IHC-Fr
	Aquaporin 4 Rabbit pAb	HA500125	Human, Mouse, Rat	WB,IHC-P
	Aquaporin 4 Recombinant Antibody [PSH06-38] - Chicken IgY (Chimeric) - BSA and Azide free	HA610294	Mouse	IHC-Fr
Relative expression	Aquaporin 4 Recombinant Antibody [PSH06-38] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610295	Mouse	IHC-Fr,WB
	Aquaporin 4 Rabbit pAb	ER1903-87	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
Relative expression	Aquaporin 4 Recombinant Antibody [PSH06-38] - Rat IgG1 (Chimeric)	HA601462	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Aquaporin 4 Recombinant Rabbit mAb [PSH06-38]	HA722672	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	Aquaporin 4 Recombinant Rabbit mAb [JE59-76]	HA721954	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	Aquaporin 4 Recombinant Antibody [PSH06-38] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610296	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
KD	BDNF Recombinant Rabbit mAb [SJ12-09]	ET1606-42	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC
	BDNF Recombinant Rabbit mAb [JE59-35]	HA722912	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-P,IHC-Fr,FC
KD	BDNF Recombinant Rabbit mAb [SJ12-09] - BSA and Azide free	HA750104	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC
	BDNF Rabbit pAb	ER130915	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	beta Amyloid Recombinant Rabbit mAb [PSH08-18]	HA722961	Human, Mouse	WB,Dot Blot
Relative expression	beta Amyloid Recombinant Rabbit mAb [PSH08-19]	HA722962	Human, Mouse	WB,IHC-P,IF-Tissue,Dot Blot
	beta Amyloid Recombinant Rabbit mAb	IRS251RB	Mouse	mIHC
Relative expression	beta Amyloid Recombinant Rabbit mAb [PSH08-19] - BSA and Azide free	HA751210	Human, Mouse	WB,IHC-P,IF-Tissue,Dot Blot
	beta Amyloid Recombinant Rabbit mAb [PSH08-18] - BSA and Azide free	HA751209	Human, Mouse	WB,Dot Blot
Relative expression	beta Amyloid 1-40 Recombinant Rabbit mAb [PSH05-86] - BSA and Azide free	HA751034	Human, Mouse	WB,IHC-P,IF-Tissue
Relative expression	beta Amyloid 1-40 Recombinant Rabbit mAb [PSH02-82] - BSA and Azide free	HA750789	Human, Mouse, Cynomolgus monkey, Pig	IHC-P,WB
	Biotin Conjugated beta Amyloid 1-40 Recombinant Rabbit mAb [PSH03-98]	HA722618B	Human, Mouse	ELISA(Det)
	beta Amyloid 1-40 Recombinant Rabbit mAb [PSH02-82] - BSA and Azide free (Capture)	HA722077	Human, Mouse	ELISA(Cap)
	beta Amyloid 1-40 Recombinant Rabbit mAb [PSH03-98] - BSA and Azide free (Detector)	HA722078	Human, Mouse	ELISA(Det)
Relative expression	beta Amyloid 1-40 Recombinant Rabbit mAb [PSH05-86]	HA722485	Human, Mouse	WB,IHC-P,IF-Tissue
Relative expression	beta Amyloid 1-40 Recombinant Rabbit mAb [PSH02-82]	HA721788	Human, Mouse, Cynomolgus monkey, Pig	IHC-P,WB
	Biotin Conjugated beta Amyloid 1-42 Recombinant Rabbit mAb [PSH03-99]	HA722293	Human	ELISA(Det)
Relative expression	beta Amyloid 1-42 Recombinant Rabbit mAb [PSH06-81] - BSA and Azide free	HA751103	Human, Mouse	WB,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	beta Amyloid 1-42 Recombinant Rabbit mAb [PSH02-83] - BSA and Azide free	HA750790	Human, Mouse, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue
	beta Amyloid 1-42 Recombinant Rabbit mAb [PSH02-83] - BSA and Azide free (Capture)	HA722079	Human	ELISA(Cap)
	beta Amyloid 1-42 Recombinant Rabbit mAb	IRS247RB	Mouse	mIHC
Relative expression	beta Amyloid 1-42 Recombinant Rabbit mAb [PSH02-83]	HA721789	Human, Mouse, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue
	beta Amyloid 1-42 Recombinant Rabbit mAb [PSH03-99] - BSA and Azide free (Detector)	HA722080	Human	ELISA(Det)
Relative expression	beta Amyloid 1-42 Recombinant Rabbit mAb [PSH06-81]	HA722754	Human, Mouse	WB,IHC-P,IF-Tissue
Relative expression, KD	Caspr Recombinant Rabbit mAb [JE31-54]	HA722531	Human, Mouse, Rat	WB,IHC-P,IHC-Fr
	Caspr Mouse mAb [A2D3]	EM1902-07	Human, Mouse	WB,IF-Cell,IHC-P,FC
	Caspr Mouse mAb [A2D2]	EM1902-10	Human, Mouse	IHC-P,FC
	Caspr2 Recombinant Rabbit mAb [JE41-03]	HA721598	Human, Mouse, Rat	WB,IHC-P
	CCL2 Mouse mAb [E11-D8]	EM1710-22	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
	CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-44] - BSA and Azide free	HA751533	Mouse	WB,IF-Cell,FC
	CCL2 / MCP-1 Rabbit pAb	HA500267	Mouse	WB,IF-Cell,FC
	CCL2 / MCP-1 Rabbit pAb	HA500042	Human	WB,IHC-P,IF-Cell
Cell treatment	CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-16]	HA723979	Mouse, Rat	WB,IF-Cell,FC
	CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-44]	HA723691	Mouse	WB,IF-Cell,FC
	CD11b Mouse mAb [A2-B5-D2]	EM1902-35	Human, Mouse	WB
Relative expression	CD11b Recombinant Rabbit mAb [PSH03-96]	HA722075	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	CD11b Recombinant Rabbit mAb [JE01-17]	HA722856	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue
Relative expression	CD11b Recombinant Rabbit mAb [JU93-81]	ET1706-04	Human	WB,IP,IHC-P,IF-Cell,IF-Tissue,mIHC
Relative expression	CD11b Recombinant Rabbit mAb [JU93-81] - BSA and Azide free	HA750451	Human	WB,IP,IHC-P,IF-Cell,IF-Tissue
Relative expression	CD11b Mouse mAb [A2-9]	M1211-6	Human	WB,IHC-P,FC,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD11b Mouse mAb [4-6]	EM1701-42	Human	IHC-P,IF-Cell,FC,WB,IF-Tissue
	CD11b Recombinant Rabbit mAb	IRS014	Human	mIHC
Relative expression	CD11b Recombinant Rabbit mAb [PSH03-96] - BSA and Azide free	HA750921	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	CD11b/ITAM Mouse mAb [4G2]	RT1107	Human, Mouse	WB,IP,IF,FC
Relative expression	CD45 Recombinant Rabbit mAb [PSH13-01] - BSA and Azide free	HA751456	Mouse	WB,IHC-P,IHC-Fr,FC,IF-Cell,IF-Tissue
	CD45 Rabbit pAb	HA500575	Mouse	WB,IF-Cell,FC
Relative expression	CD45 Recombinant Rabbit mAb [PSH13-02] - BSA and Azide free	HA751457	Mouse, Rat	WB,IHC-Fr,IF-Tissue
Relative expression	CD45 Recombinant Rabbit Multiclonal Antibody [PSH14-19]	HA723606	Human, Mouse	WB,IHC-P,IF-Cell
	CD45 Recombinant Rabbit mAb	IRS025	Human	mIHC
	CD45 Recombinant Rabbit mAb [JE03-05]	ET7111-03	Human	WB,IHC-P,FC,IF-Tissue,mIHC
	CD45 Rabbit pAb	ER1901-29	Human	WB,IHC-P,FC
	CD45 Rabbit pAb	ER1902-02	Human	WB,FC
Relative expression	CD45 Recombinant Rabbit Multiclonal Antibody [PSH14-19] - BSA and Azide free	HA751510	Human, Mouse	WB,IHC-P,IF-Cell
Relative expression	CD45 Recombinant Rabbit mAb [PSH13-01]	HA723506	Mouse	WB,IHC-P,IHC-Fr,FC,IF-Cell,IF-Tissue
Relative expression	CD45 Recombinant Rabbit mAb [PSH13-02]	HA723507	Mouse, Rat	WB,IHC-Fr,IF-Tissue
	CD45 Rabbit pAb	ER1905-93	Human	FC
	CD45R Rat mAb	IRS076	Mouse	mIHC
	PE Conjugated CD45RO Recombinant Mouse mAb [PSH13-94]	HA601531F	Human	FC
	CD45RO Recombinant Mouse mAb [PSH13-94]	HA601441	Human	IHC-P,FC
	CD45RO Recombinant Mouse mAb	IRS124MS	Human	mIHC
	Claudin 11 Rabbit pAb	ER65869	Human, Mouse, Rat	WB
	Claudin 5 Recombinant Rabbit mAb [JM11-22] - BSA and Azide free	HA750384	Human	WB,IHC-P
	Claudin 5 Recombinant Rabbit mAb [JM11-22]	ET1703-58	Human	WB,IHC-P
Relative expression	CNPase Recombinant Rabbit mAb [JF10-25]	ET1702-46	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Complement C3 Mouse mAb [6-A10]	M1510-39	Human	WB,IHC-P,ELISA
	Complement C3 Recombinant Rabbit mAb [JF10-30] - BSA and Azide free	HA750367	Human, Mouse, Rat	WB,IF-Cell,FC
	Complement C3 Recombinant Rabbit mAb [PSH18-31]	HA723984	Mouse, Rat	WB,IHC-Fr,IHC-P
	Complement C3 Mouse mAb [10-B3]	M1510-40	Human	WB,IHC-P,ELISA
	Complement C3 Recombinant Rabbit mAb [JF10-30]	ET1702-99	Human, Mouse, Rat	WB,IF-Cell,FC
	Complement C3 Recombinant Rabbit mAb [JU42-68]	ET7106-75	Human	WB,IHC-P,IP
	Complement C3 (beta chain) Mouse mAb [11-10]	M0809-14	Human	WB,IF-Cell,FC
	Complement C3 (beta chain) Rabbit pAb	0809-8	Human	WB
	Complement C5 Recombinant Rabbit mAb [PSH04-01]	HA722082	Human	WB,IF-Cell,IHC-P,FC
	Complement C5 Recombinant Rabbit mAb [PSH04-01] - BSA and Azide free	HA750924	Human	WB,IF-Cell,IHC-P,FC
	Connexin 43 Recombinant Rabbit mAb [PSH0-32]	HA721312	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
	Connexin 43 Recombinant Rabbit mAb [PSH0-32] - BSA and Azide free	HA750606	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
	CX3CR1 Mouse mAb [A5C7]	HA600074	Human	IF-Cell
	CX3CR1 Recombinant Rabbit mAb	IRS250RB	Mouse	mIHC
	Dnmt1 Recombinant Rabbit mAb [JF09-89]	ET1702-77	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Dnmt1 Recombinant Rabbit mAb [JF09-89] - BSA and Azide free	HA750359	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	DNMT1 Recombinant Mouse mAb [A3A8-R]	HA601281	Human	WB,IF-Cell,IHC-P
	DNMT1 Mouse mAb [A2-C11]	EM1707-16	Human, Monkey	WB,IHC-P
	DNMT1 Mouse mAb [A3A9]	EM1901-84	Human	WB,IHC-P,FC
	DNMT1 Mouse mAb [A3A8]	EM1901-83	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	DNMT1 Recombinant Mouse mAb [A3A8-R] - BSA and Azide free	HA610156	Human	WB,IF-Cell,IHC-P
	Dnmt3a Rabbit pAb	ER1902-62	Human, Mouse, Rat	WB,IHC-P,FC
KD	Dnmt3a Recombinant Rabbit mAb [ST04-78]	ET1609-31	Human, Mouse, Rat	WB,IF-Tissue,IHC-P
KD	Dnmt3a Recombinant Rabbit mAb [ST04-78] - BSA and Azide free	HA750169	Human, Mouse, Rat	WB,IF-Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	EAAT1 Recombinant Rabbit mAb [JA30-35]	ET1704-54	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
Cell treatment	EAAT1 Recombinant Rabbit mAb [JA30-35] - BSA and Azide free	HA750415	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
	Ermin Rabbit pAb	ER1908-31	Mouse, Rat	WB,IHC-P,IF-Tissue
	GDNF Recombinant Rabbit mAb [JA93-10]	ET1704-46	Human, Mouse, Rat	WB,IF-Cell
	iFluor™ 488 Conjugated GFAP Mouse mAb [1-D4]	HA600101F	Human, Rat	IF-Tissue,IHC-Fr
	GFAP Recombinant Rabbit mAb	IRS069	Mouse	mIHC
Relative expression	GFAP Recombinant Mouse mAb [PSH10-24]	HA601380	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue
	Biotin Conjugated GFAP Mouse mAb [PSH05-17]	HA601355B	Human, Mouse, Rat	ELISA(Det)
	GFAP Rabbit pAb	R1308-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	GFAP Recombinant Mouse mAb [PSH10-24] - BSA and Azide free	HA610225	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue
	iFluor™ 594 Conjugated GFAP Mouse mAb [1-D4]	HA600102F	Human, Rat	IF-Tissue
	GFAP Recombinant Rabbit mAb [PSH05-16] - BSA and Azide free (Capture)	HA722234	Human, Mouse, Rat	ELISA(Cap)
Relative expression	GFAP Recombinant Antibody [SA03-04] - Rat IgG1 (Chimeric)	HA601483	Mouse, Rat	IHC-Fr,WB
Relative expression	GFAP Recombinant Antibody [SA03-04] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610317	Mouse, Rat	IHC-Fr,WB
	GFAP Rabbit pAb	R1308-10	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	GFAP Recombinant Antibody - Rat IgG1 (Chimeric)	IRS148RT	Mouse	mIHC
	GFAP Recombinant Mouse mAb [PSH05-17] - BSA and Azide free (Detector)	HA601342	Human, Mouse, Rat	ELISA(Det)
	GFAP Recombinant Rabbit mAb [SA03-04]	ET1601-23	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr,mIHC,IP
	GFAP Mouse mAb [1-D4]	EM140707	Human, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
	iFluor™ 647 Conjugated GFAP Mouse mAb [1-D4]	HA600103F	Human, Rat	IF-Tissue
KD	Glucose Transporter GLUT1 Recombinant Rabbit mAb [SA0377] - BSA and Azide free	HA750002	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Glucose Transporter GLUT1 Recombinant Rabbit mAb	IRS151RB	Human	mIHC
	Glucose Transporter GLUT1 Recombinant Antibody [SA0377] - Rat IgG1 (Chimeric)	HA601521	Human, Mouse, Rat	IHC-Fr,IHC-P,IF-Tissue,WB
KD	Glucose Transporter GLUT1 Recombinant Rabbit mAb [SA0377]	ET1601-10	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	iFluor™ 647 Conjugated Glucose Transporter GLUT1 Recombinant Rabbit mAb [SA0377]	HA720170F	Human, Mouse	IF-Tissue
	Glucose Transporter GLUT1 Rabbit pAb	ER1510-11	Human, Mouse	WB,IF-Cell,IHC-P,FC
	iFluor™ 488 Conjugated Glucose Transporter GLUT1 Recombinant Rabbit mAb [SA0377]	HA720154F	Human, Mouse, Rat	IF-Cell,IF-Tissue
	Glucose Transporter GLUT1 Recombinant Antibody [SA0377] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610348	Human, Mouse, Rat	IHC-Fr,IHC-P,IF-Tissue,WB
KD	Glucose Transporter GLUT3 Recombinant Rabbit mAb [JA50-31]	ET1704-59	Human, Mouse, Rat	WB,FC
KD	Glucose Transporter GLUT3 Recombinant Rabbit mAb [JA50-31] - BSA and Azide free	HA750416	Human, Mouse, Rat	WB,FC
	Glucose Transporter GLUT4 Mouse mAb [3-A10]	M1505-6	Human, Mouse	WB,IF-Cell,IHC-P
Relative expression	Glucose Transporter GLUT4 Rabbit pAb	R1402-3	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	GLUT12 Rabbit pAb	ER1909-73	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
	GLUT2 Recombinant Rabbit mAb [JJ20-21]	ET1701-34	Human, Mouse, Rat	WB,IHC-P,FC
	Glut5 Rabbit pAb	ER65566	Human	WB
	GLUT8 Rabbit pAb	ER1909-74	Mouse, Rat	WB
	GPC1 Rabbit pAb	R1512-15	Human	WB,IF-Cell,IHC-P,FC
	GPR105 Rabbit pAb	ER63561	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	GRM3 Mouse mAb [G2-G9]	EM1706-64	Human	WB,IF-Cell,IHC-P,FC
KO	HDAC1 Recombinant Rabbit mAb [SY12-04] - BSA and Azide free	HA750083	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,IF-Tissue
KO	HDAC1 Recombinant Rabbit mAb [SY12-04]	ET1605-35	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,IF-Tissue
KO	HDAC2 Rabbit pAb	R1601-7	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KO	HDAC2 Recombinant Rabbit mAb [SY29-02]	ET1607-78	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO	HDAC2 Recombinant Rabbit mAb [SY29-02] - BSA and Azide free	HA750134	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO	HDAC2 Recombinant Rabbit mAb [SD0816]	ET1612-50	Human, Mouse, Rat	WB,IP,IHC-P,FC
KD	HDAC3 Recombinant Rabbit mAb [SC55-02] - BSA and Azide free	HA750221	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
KD	HDAC3 Mouse mAb [B2-G9]	M1511-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	HDAC3 Recombinant Rabbit mAb [SC55-02]	ET1610-5	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
KD	HDAC4 Recombinant Rabbit mAb [SD08-35]	ET1612-51	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	HDAC5 Rabbit pAb	ER64245	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	HDAC5 Recombinant Rabbit mAb [PSH08-70] - BSA and Azide free	HA751235	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,ChIP
	HDAC5 Recombinant Rabbit mAb [PSH08-70]	HA723022	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,ChIP
	HDAC5/9 Rabbit pAb	ER63671	Human, Mouse	WB,IHC-P
KD	HDAC6 Recombinant Rabbit mAb [JJ09-09]	ET1701-66	Human, Monkey	WB,IF-Tissue,IHC-P,IP
	HDAC7 Recombinant Rabbit mAb [SD082-4]	ET1612-67	Human, Mouse, Rat	WB,IF-Cell,FC
	HDAC8 Recombinant Rabbit mAb [SD082-07]	ET1612-90	Human, Mouse, Rat	WB,FC,IP
	HDAC8 Recombinant Rabbit mAb [JJ0845]	ET1701-12	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
KD	HDAC9 Recombinant Rabbit mAb [JU30-44]	ET1706-36	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Histone deacetylase 9 Rabbit pAb	ER63682	Human	WB,IHC-P,IF-Cell
	HRG1 Rabbit pAb	ER1910-81	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IF-Tissue
	HRG beta 1 Rabbit pAb	ER1910-80	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
	Human Alpha-Synuclein Recombinant Rabbit mAb [JB42-33] - BSA and Azide free (Capture)	HA723467	Human	ELISA(Cap)
	Human Alpha-Synuclein Recombinant Rabbit mAb [PSH12-56] - BSA and Azide free (Detector)	HA723468	Human	ELISA(Det)
	Biotin Conjugated Human Alpha-Synuclein Recombinant Rabbit mAb [PSH12-56] - Detector	HA723469B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-22] - Detector	HA725196B	Human	ELISA(Det)
	Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-22] - BSA and Azide free (Detector)	HA725196	Human	ELISA(Det)
	Human CCL2 / MCP-1 Recombinant Rabbit mAb [PSH18-21] - BSA and Azide free (Capture)	HA725195	Human	ELISA(Cap)
	FITC Conjugated Human CD11b Recombinant Mouse mAb [A2-9-R]	HA600130F	Human	FC
	FITC Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600125F	Human	FC
	iFluor™ 555 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA601556F	Human	FC
	Human CD45 Recombinant Mouse mAb [PSH03-64] - BSA and Azide free	HA610158	Human	IF-Cell,FC
	Human CD45 Recombinant Mouse mAb [PSH03-64]	HA601285	Human	IF-Cell,FC
	PerCP-Cy5.5 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600128F	Human	FC
	APC Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600127F	Human	FC
	PE Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600126F	Human	FC
	iFluor™ 488 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600109F	Human	FC
	FITC Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA600112F	Human	FC
	Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA601341	Human	FC,IF-Cell,IHC-P
	FITC Conjugated Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA600131F	Human	FC
	PE Conjugated Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA601525F	Human	FC
	Biotin Conjugated Human Complement C3 Recombinant Rabbit mAb [PSH15-08] - Detector	HA725142B	Human	ELISA(Det),ELISA
	Human Complement C3 Recombinant Rabbit mAb [PSH15-08] - BSA and Azide free (Detector)	HA725141	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Complement C3 Recombinant Rabbit mAb [PSH15-07] - BSA and Azide free (Capture)	HA725140	Human	ELISA(Cap)
	HRP Conjugated Human Complement C3 Recombinant Rabbit mAb [PSH15-08] - Detector	HA725142H	Human	ELISA(Det),ELISA
	Human Complement C5 Recombinant Rabbit mAb [PSH06-45] - BSA and Azide free (Capture)	HA722682	Human	ELISA(Cap)
	Human Complement C5 Recombinant Rabbit mAb [PSH06-46] - BSA and Azide free (Detector)	HA722683	Human	ELISA(Det)
	HRP Conjugated Human Complement C5 Recombinant Rabbit mAb [PSH06-46] - Detector	HA725127H	Human	ELISA(Det)
	Biotin Conjugated Human Complement C5 Recombinant Rabbit mAb [PSH06-46]	HA722684B	Human	ELISA(Det)
	Biotin Conjugated Human CSPG4 Recombinant Rabbit mAb [PSH10-79] - Detector	HA723247B	Human	ELISA(Det),ELISA
	Human CSPG4 Recombinant Rabbit mAb [PSH10-77] - BSA and Azide free (Capture)	HA723244	Human	ELISA(Cap)
	Human CSPG4 Recombinant Rabbit mAb [PSH10-78] - BSA and Azide free (Capture)	HA723245	Human	ELISA(Cap)
	HRP Conjugated Human CSPG4 Recombinant Rabbit mAb [PSH10-79] - Detector	HA723376H	Human	ELISA(Det)
	Human CSPG4 Recombinant Rabbit mAb [PSH10-79] - BSA and Azide free (Detector)	HA723246	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free	HA750745	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-85] - BSA and Azide free (Capture)	HA723881	Human	ELISA(Cap)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free	HA750744	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-87] - BSA and Azide free (Capture)	HA723883	Human	ELISA(Cap)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - BSA and Azide free (Detector)	HA723882	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39]	HA721712	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free (Detector)	HA721279	Human	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free (Capture)	HA721278	Human	ELISA(Cap)
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - Detector	HA723882B	Human	ELISA(Det),ELISA
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA721713	Human	WB
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA722649B	Human	ELISA(Det)
	Human NLRP3 Recombinant Rabbit mAb [PSH11-08] - BSA and Azide free (Detector)	HA723276	Human	ELISA(Det)
	Biotin Conjugated Human NLRP3 Recombinant Rabbit mAb [PSH11-09] - Detector	HA723279B	Human	ELISA(Det),ELISA
	Human NLRP3 Recombinant Rabbit mAb [PSH11-07] - BSA and Azide free (Capture)	HA723275	Human	ELISA(Cap)
	Human NLRP3 Recombinant Rabbit mAb [PSH11-09] - BSA and Azide free (Detector)	HA723278	Human	ELISA(Det)
	Biotin Conjugated Human NLRP3 Recombinant Rabbit mAb [PSH11-08] - Detector	HA723277B	Human	ELISA(Det),ELISA
	Human TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free (Detector)	HA721266	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb	IRS040	Human	mIHC
	Biotin Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721559	Human	ELISA(Det)
	HRP Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721294	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free (Capture)	HA721265	Human	ELISA(Cap)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40]	HA721286	Human	IHC-P,mIHC
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free	HA750593	Human	IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-92] - Detector	HA723495B	Human	ELISA(Det),ELISA
	Human VEGFA Recombinant Rabbit mAb [PSH12-91] - BSA and Azide free (Capture)	HA723493	Human	ELISA(Cap)
	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-93] - Detector	HA723497B	Human	ELISA(Det),ELISA
	Human VEGFA Recombinant Rabbit mAb [PSH12-92] - BSA and Azide free (Detector)	HA723494	Human	ELISA(Det)
	Human VEGFA Recombinant Rabbit mAb [PSH12-93] - BSA and Azide free (Detector)	HA723496	Human	ELISA(Det)
	Iba1 Recombinant Antibody [JM36-62] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610213	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
Relative expression	Iba1 Recombinant Rabbit mAb [JM36-62] - BSA and Azide free	HA750440	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC,IP
Relative expression	Iba1 Recombinant Rabbit mAb [JM36-62]	ET1705-78	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC,IP,mIHC
	iFluor™ 594 Conjugated Iba1 Recombinant Rabbit mAb [JM36-62]	HA723821F	Mouse	IHC-Fr
	Iba1 Recombinant Antibody [JM36-62] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610214	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,WB
	Iba1 Recombinant Antibody [JM36-62] - Rat IgG1 (Chimeric)	HA601367	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
	Iba1 Rabbit pAb	ER1802-20	Human, Mouse	IHC-P
	Iba1 Recombinant Mouse mAb [8G8] - BSA and Azide free	HA610321	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,FC,IP,IHC-Fr,IF-Cell,IF-Tissue
	iFluor™ 488 Conjugated Iba1 Recombinant Rabbit mAb [JM36-62]	HA720158F	Human, Mouse, Rat	IF-Cell,FC
	Iba1 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS158RT	Mouse	mIHC
	Iba1 Recombinant Rabbit mAb	IRS068	Mouse	mIHC
Relative expression	Iba1 Recombinant Antibody [JM36-62] - Chicken IgY (Chimeric)	HA601376	Human, Mouse, Rat, Pig, Cynomolgus monkey	IHC-Fr,WB,FC
	Iba1 Recombinant Antibody [JM36-62] - Guinea pig IgG2 (Chimeric)	HA601368	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,WB
	Iba1 Recombinant Mouse mAb [8G8]	RT1316	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,FC,IP,IHC-Fr,IF-Cell,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	Iba1 Recombinant Antibody [JM36-62] - Chicken IgY (Chimeric) - BSA and Azide free	HA610221	Human, Mouse, Rat, Pig, Cynomolgus monkey	IHC-Fr,WB,FC
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R]	HA601036	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Mouse mAb [A7A10]	HA601035	Human	WB,ELISA,IF-Cell
	IL-1 beta Recombinant Rabbit mAb [PSH18-00]	HA723965	Mouse, Rat	WB
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R] - BSA and Azide free	HA610004	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Mouse mAb [A7F9]	HA601002	Human, Mouse, Rat	WB,ELISA,IF-Cell
Cell treatment	IL-1 beta Recombinant Rabbit mAb [JJ087-3]	ET1701-39	Human	WB,IF-Cell,IF-Tissue,IHC-P
	Kir4.1 Rabbit pAb	ER60120	Human, Rat, Mouse	IHC-P
	mGluR2 Recombinant Rabbit mAb [JJ091-09]	ET1701-74	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	mGluR3 Recombinant Rabbit mAb [JE41-18]	HA722434	Human, Mouse, Rat	WB,IHC-P
Relative expression	MOG Recombinant Rabbit mAb [JM23-06]	ET1705-16	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
	Biotin Conjugated Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-45] - Detector	HA725042B	Mouse	ELISA(Det),ELISA
	Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-45] - BSA and Azide free (Detector)	HA725041	Mouse	ELISA(Det)
	Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-44] - BSA and Azide free (Capture)	HA725040	Mouse	ELISA(Cap)
	Mouse CD45R Recombinant Rat mAb [PSH10-80]	HA601385	Mouse	IHC-P,FC,IF-Tissue
	Mouse CD45R Recombinant Rat mAb [PSH10-80] - BSA and Azide free	HA610229	Mouse	IHC-P,FC,IF-Tissue
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67] - BSA and Azide free (Detector)	HA722167	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - BSA and Azide free (Detector)	HA721977	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH09-59] - BSA and Azide free	HA723121	Mouse	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IL-1 beta Recombinant Rabbit mAb [PSH09-58] - BSA and Azide free	HA723120	Mouse	ELISA
	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - Detector	HA722625B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67]	HA722608B	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-31] - BSA and Azide free (Capture)	HA721976	Mouse	ELISA(Cap)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-66] - BSA and Azide free (Capture)	HA722166	Mouse	ELISA(Cap)
	Mouse/Rat Complement C3 Recombinant Rabbit mAb [PSH16-81] - BSA and Azide free (Capture)	HA725143	Mouse, Rat	ELISA(Cap)
	Mouse/Rat Complement C3 Recombinant Rabbit mAb [PSH16-82] - BSA and Azide free (Detector)	HA725144	Mouse, Rat	ELISA(Det)
	HRP Conjugated Mouse/Rat Complement C3 Recombinant Rabbit mAb [PSH16-82] - Detector	HA725145H	Mouse, Rat	ELISA(Det),ELISA
	Biotin Conjugated Mouse/Rat Complement C3 Recombinant Rabbit mAb [PSH16-82] - Detector	HA725145B	Mouse, Rat	ELISA(Det),ELISA
	Biotin Conjugated Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30]	HA722626B	Mouse, Rat	ELISA(Det)
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30] - BSA and Azide free (Detector)	HA721975	Mouse, Rat	ELISA(Det)
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-29] - BSA and Azide free (Capture)	HA721974	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77]	HA722623B	Mouse	ELISA(Det)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77] - BSA and Azide free (Detector)	HA722041	Mouse	ELISA(Det)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-76] - BSA and Azide free (Capture)	HA722040	Mouse	ELISA(Cap)
	Mouse TREM2 Recombinant Rabbit mAb [PSH02-03] - BSA and Azide free (Detector)	HA721761	Mouse	ELISA(Det)
	Mouse TREM2 Recombinant Rabbit mAb [PSH02-02] - BSA and Azide free (Capture)	HA721760	Mouse	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse TREM2 Recombinant Rabbit mAb [PSH02-03] - Detector	HA722642B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse VEGF Recombinant Rabbit mAb [PSH02-73]	HA722636B	Mouse	ELISA(Det)
	Mouse VEGF Recombinant Rabbit mAb [PSH02-72] - BSA and Azide free (Capture)	HA721850	Mouse	ELISA(Cap)
	Mouse VEGF Recombinant Rabbit mAb [PSH02-73] - BSA and Azide free (Detector)	HA721851	Mouse	ELISA(Det)
	MRF/C11orf9 Rabbit pAb	ER1913-17	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	Myelin Basic Protein Recombinant Rabbit mAb [JF0943]	ET1702-15	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr,mIHC
	Myelin Basic Protein Recombinant Rabbit mAb	IRS064	Mouse	mIHC
Relative expression	Myelin Basic Protein Recombinant Rabbit mAb [JF0943] - BSA and Azide free	HA750337	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression	Myelin PLP Recombinant Antibody [PSH08-21] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610247	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Myelin PLP Recombinant Antibody [PSH08-21] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610304	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Myelin PLP Recombinant Antibody [PSH08-21] - Rat IgG1 (Chimeric)	HA601407	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Myelin PLP Recombinant Rabbit mAb [PSH08-21] - BSA and Azide free	HA751212	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression	Myelin PLP Chicken pAb	HA500573	Mouse, Rat	IHC-Fr,WB
Relative expression	Myelin PLP Recombinant Antibody [PSH08-21] - Guinea pig IgG2 (Chimeric)	HA601469	Mouse, Rat	IHC-Fr,WB
Relative expression	Myelin PLP Recombinant Antibody [PSH08-21] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610303	Mouse, Rat	IHC-Fr,WB
Relative expression	Myelin PLP Recombinant Antibody [PSH08-21] - Mouse IgG1 (Chimeric)	HA601470	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Myelin PLP Recombinant Rabbit mAb [PSH08-21]	HA722964	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression	Myelin PLP Rabbit pAb	HA500202	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Neuregulin-1 SMDf Rabbit pAb	ER64158	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KD	NF-kB p100 / NFkB2 Recombinant Rabbit mAb [JM82-03]	ET1705-45	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	NF-kB p100 / NFKB2 Mouse mAb [16B2]	EM1901-78	Human	WB,IHC-P,FC
KD	NF-kB p100 / NFKB2 Mouse mAb [16B3]	EM1901-79	Human, Mouse, Rat	WB,IHC-P,FC
KD	NF-kB p105 / p50 Recombinant Rabbit mAb [SZ20-01]	ET1603-18	Human, Mouse, Rat	WB,IHC-P,ICC,FC,IF-Tissue
KD	NF-kB p105 / p50 Recombinant Rabbit mAb [SZ20-01] - BSA and Azide free	HA750060	Human, Mouse, Rat	WB,IHC-P,ICC,FC,IF-Tissue
	NF-kB p105 / p50 Mouse mAb [A2-7]	M1505-8	Human, Mouse	WB,IF-Cell,IHC-P
	NF-kB p105 / p50 Rabbit pAb	R1309-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27] - BSA and Azide free	HA750601	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	NF-kB p65 Rabbit pAb	HA500402	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04]	ET1603-12	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	NF-kB p65 Rabbit pAb	ER0815	Human, Mouse, Rat, Zebrafish	WB,IHC-P,FC
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27]	HA721307	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04] - BSA and Azide free	HA750058	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Relative expression	NG2 Recombinant Rabbit mAb [JM10-13]	ET1703-16	Human, Mouse, Rat	WB,IHC-P
Relative expression	NG2 Recombinant Rabbit mAb [PSH09-98]	HA723167	Human	WB,IF-Cell,FC
Relative expression	NG2 Recombinant Rabbit mAb [PSH09-98] - BSA and Azide free	HA751329	Human	WB,IF-Cell,FC
Relative expression	NG2 Recombinant Rabbit mAb [JM10-13] - BSA and Azide free	HA750372	Human, Mouse, Rat	WB,IHC-P
	NLGN3 Rabbit pAb	ER60922	Human, Mouse, Rat	WB
	NLRP3 Rabbit pAb	ER1706-72	Human, Mouse	WB,IF-Cell,IHC-P,FC
Relative expression	NLRP3 Recombinant Rabbit mAb [PSH19-43]	HA724075	Human, Mouse	WB
Cell treatment	NLRP3 Recombinant Rabbit mAb [SC06-23]	ET1610-93	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	NLRP3 Recombinant Rabbit mAb [SC06-23] - BSA and Azide free	HA750236	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Nrf2 Rabbit pAb	R1312-8	Human, Mouse	WB,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Nrf2 Recombinant Rabbit mAb [JE30-14]	HA721432	Human, Mouse	WB,IF-Cell
Cell treatment	Nrf2 Rabbit pAb	ER1706-41	Human, Mouse, Rat	IHC-P,IF-Cell,WB
Cell treatment	Nrf2 Recombinant Rabbit mAb [PSH11-20]	HA723302	Human, Mouse, Rat	WB,IF-Cell,ChIP
Cell treatment	Nrf2 Recombinant Rabbit mAb [PSH11-20] - BSA and Azide free	HA751401	Human, Mouse, Rat	WB,IF-Cell,ChIP
	Occludin Recombinant Rabbit mAb	IRS183RB	Mouse	mIHC
	Occludin Rabbit pAb	R1510-33	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	Occludin Recombinant Rabbit mAb [JJ091-08]	ET1701-76	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IP,FC,IF-Cell,IHC-P,IHC-Fr,IF-Tissue
	OLIG1 Rabbit pAb	ER61887	Human, Rat, Mouse	WB
Relative expression	Olig2 Recombinant Rabbit mAb [SP07-02]	ET1604-29	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,mIHC,IF-Cell
	Olig2 Recombinant Antibody [SP07-02] - Rat IgG1 (Chimeric)	HA601558	Mouse, Rat	IHC-Fr,IHC-P
	Olig2 Recombinant Rabbit mAb	IRS067	Mouse	mIHC
	Omgp Rabbit pAb	ER1914-37	Human, Mouse, Rat	WB
	P2RX1 Rabbit pAb	ER1901-53	Human, Mouse, Rat	WB,IHC-P
	P2RX4 Rabbit pAb	ER60964	Human, Rat, Mouse	WB
	P2RX7 Mouse mAb [G1-D6]	EM1708-40	Human	WB,IHC-P,FC
	P2RX7 Rabbit pAb	ER1919-21	Human, Mouse	WB,IHC-P,IF-Cell,IF-Tissue
	P2RY1 Mouse mAb [4G5D6]	EM1710-48	Human, Rat	WB,IHC-P,FC
	P2RY14 Mouse mAb [5A3F1]	EM1710-54	Human	WB,IHC-P,IF-Cell,FC
	P2RY2 Rabbit pAb	ER62266	Human, Rat, Mouse	WB
	P2RY4 Mouse mAb [3B4C11]	EM1710-51	Human	WB,IHC-P
	P2RY8 Mouse mAb [1G5A11]	EM1710-43	Human	WB,IHC-P,IF-Cell,FC
	P2X2 Rabbit pAb	ER1902-59	Human, Mouse, Rat	WB,IHC-P,FC
	P2X3 Rabbit pAb	ER1902-88	Human, Mouse, Rat	WB,IHC-P
	P2X6 Rabbit pAb	ER1901-49	Human, Mouse, Rat	WB,IHC-P,FC
	P2X7 Rabbit pAb	ER1901-99	Human, Mouse, Rat	WB,IHC-P,FC
	P2X7 Rabbit pAb	ER1901-50	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	P2Y10 Rabbit pAb	ER1914-60	Human	WB,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	P2Y11 Rabbit pAb	ER64494	Human	WB
	P2Y12 Recombinant Antibody [PSH04-48] - Rat IgG1 (Chimeric)	HA601388	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
	P2Y12 Recombinant Rabbit mAb [PSH04-48]	HA722133	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr,IF-Tissue
	P2Y12 Recombinant Rabbit mAb [PSH04-48] - BSA and Azide free	HA750957	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr,IF-Tissue
	P2Y12 Recombinant Antibody [PSH04-48] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610231	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
	P2Y12 Recombinant Rabbit mAb [JE63-23]	HA721011	Human, Mouse	WB
	P2Y4 Rabbit pAb	ER64495	Human, Mouse, Rat	WB,IF-Cell
	P2Y6 Recombinant Rabbit mAb [JG37-77]	ET7108-52	Human	WB,IP,IHC-P,IF-Cell
	P2Y8 Rabbit pAb	ER1914-62	Human	WB,FC
	Phospho-alpha Synuclein (S129) Recombinant Rabbit mAb [JB22-44]	ET7107-30	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Phospho-Glucose Transporter GLUT4 (S488) Recombinant Rabbit mAb [JE43-23]	HA721866	Mouse, Human, Rat	WB
Cell treatment	Phospho-NF-kB p65 (S529) Recombinant Rabbit mAb [SP07-00]	ET1604-27	Human	WB,IHC-P,FC,IP
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58] - BSA and Azide free	HA751361	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58]	HA723223	Human, Mouse, Rat	WB,IF-Cell
Cell treatment, Relative expression	Phospho-Nrf2 (S40) Recombinant Rabbit mAb [SU0334] - BSA and Azide free	HA750142	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Cell treatment, Relative expression	Phospho-Nrf2 (S40) Recombinant Rabbit mAb [SU0334]	ET1608-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	PMP22 Recombinant Rabbit mAb [JM52-30]	ET1705-15	Human, Mouse, Rat	WB
	PMP22 Recombinant Rabbit mAb [JM52-30] - BSA and Azide free	HA750425	Human, Mouse, Rat	WB
	Pro-neuregulin-1, membrane-bound isoform Rabbit pAb	ER1915-45	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	PTPRC Mouse mAb [C11-C10]	EM1708-69	Human	WB,IHC-P
	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - Detector	HA721611	Rat	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-59] - BSA and Azide free (Capture)	HA721836	Rat	ELISA(Cap)
	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60] - BSA and Azide free (Detector)	HA721837	Rat	ELISA(Det)
	Rat IL-1 beta Recombinant Rabbit mAb [PS01-86] - BSA and Azide free (Capture)	HA721609	Rat	ELISA(Cap)
	Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - BSA and Azide free (Detector)	HA721610	Rat	ELISA(Det)
	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60]	HA722634B	Rat	ELISA(Det)
Relative expression	S100A8 Recombinant Rabbit mAb [PSH09-70] - BSA and Azide free	HA751312	Human	WB,IHC-P,IF-Cell,FC,IF-Tissue,IP
Relative expression	S100A8 Recombinant Rabbit mAb [PSH09-70]	HA723135	Human	WB,IHC-P,IF-Cell,FC,IF-Tissue,IP
	S100A8+S100A9 Recombinant Rabbit Multiclonal Antibody	IRS196RB	Mouse	mIHC
	S100 alpha Recombinant Rabbit mAb [JM51-31] - BSA and Azide free	HA750389	Human, Mouse, Rat	WB,IP,IHC-P,FC,IF-Cell
	S100 alpha Recombinant Rabbit mAb [JM51-31]	ET1703-72	Human, Mouse, Rat	WB,IP,IHC-P,FC,IF-Cell
	Biotin Conjugated S100B Recombinant Rabbit mAb [PSH15-25] - Detector	HA723730B	Human, Mouse, Rat	ELISA(Det),ELISA
	S100B Recombinant Rabbit mAb [SC57-02] - BSA and Azide free (Capture)	HA723728	Human, Mouse, Rat	ELISA(Cap)
	S100B Recombinant Rabbit mAb [PSH15-25] - BSA and Azide free (Detector)	HA723729	Human, Mouse, Rat	ELISA(Det)
	S100 beta Recombinant Rabbit mAb [SC57-02] - BSA and Azide free	HA750209	Human, Mouse, Rat, Zebrafish, Goat, Cynomolgus monkey, Green monkey, Pig, Rabbit	WB,IF-Tissue,IHC-P,IHC-Fr,IF-Cell,FC,IP
	S100 beta Recombinant Antibody [SC57-02] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610365	Mouse	IHC-Fr
	S100 beta Recombinant Rabbit mAb [PSH17-13] - BSA and Azide free	HA751617	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Tissue
Relative expression	S100 beta Recombinant Mouse mAb [PD00-11]	HA601098	Human, Mouse, Rat	IHC-P,WB,IF-Tissue
	iFluor™ 594 Conjugated S100 beta Recombinant Rabbit mAb [SC57-02]	HA720125F	Human	IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	S100 beta Recombinant Rabbit mAb [PSH17-13]	HA723900	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Tissue
	S100 beta Recombinant Antibody [SC57-02] - Guinea pig IgG2 (Chimeric)	HA601550	Mouse	IHC-Fr
	S100 beta Recombinant Antibody [SC57-02] - Mouse IgG1 (Chimeric)	HA601551	Mouse, Rat	IHC-Fr,IHC-P
	S100 beta Rabbit pAb	RT1547	Human, Mouse, Rat	WB,IP,IF
	S100 beta Recombinant Rabbit mAb [SC57-02]	ET1610-3	Human, Mouse, Rat, Zebrafish, Goat, Cynomolgus monkey, Green monkey, Pig, Rabbit	WB,IF-Tissue,IHC-P,IHC-Fr,IF-Cell,FC,IP
	S100 beta Recombinant Antibody [SC57-02] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610367	Mouse, Rat	IHC-Fr,IHC-P
	S100 beta Recombinant Antibody [SC57-02] - Rat IgG1 (Chimeric)	HA601552	Mouse, Rat	IHC-Fr,IHC-P
	S100 beta Recombinant Antibody [SC57-02] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610366	Mouse, Rat	IHC-Fr,IHC-P
	SLC12A1 / NKCC2 Recombinant Rabbit mAb	IRS236RB	Mouse	mIHC
	SLC12A1 / NKCC2 Recombinant Rabbit mAb	IRS181RB	Human	mIHC
	SLC12A1 / NKCC2 Recombinant Rabbit mAb [JE64-26]	HA721906	Human, Mouse, Rat	WB,IHC-P,IF-Cell,mIHC
	SLC12A4 Rabbit pAb	ER60119	Human, Rat, Mouse	WB,IHC-P
	Slc12a5 / KCC2 Recombinant Mouse mAb [PSH13-20]	HA601433	Mouse, Rat	IHC-Fr,WB
	Slc12a5 / KCC2 Recombinant Mouse mAb [PSH13-20] - BSA and Azide free	HA610271	Mouse, Rat	IHC-Fr,WB
	Survivin Recombinant Rabbit mAb [SR44-02]	ET1602-43	Mouse, Rat	WB,IHC-P,IP
	Survivin Recombinant Rabbit mAb [PSH11-73]	HA723369	Human, Mouse	WB,IHC-P
	Survivin Recombinant Rabbit mAb [SP07-06] - BSA and Azide free	HA750076	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Survivin Recombinant Rabbit mAb [SP07-06]	ET1604-34	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Survivin Recombinant Rabbit mAb [PSH11-73] - BSA and Azide free	HA751428	Human, Mouse	WB,IHC-P
	Thrombospondin 1 Recombinant Rabbit mAb [PSH02-97] - BSA and Azide free	HA750852	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Thrombospondin 1 Recombinant Rabbit mAb [PSH02-97]	HA721916	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
	TLR4 Rabbit pAb	ER1706-43	Human, Mouse, Rat	IF-Cell,IHC-P,FC
	TLR4 Mouse mAb	RT1666	Human, Mouse, Rat	WB
	TMEM163 Rabbit pAb	ER1901-87	Human, Mouse	WB,IF-Cell,FC
	TMEM163 Rabbit pAb	ER1901-91	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	TNF alpha Rabbit pAb	ER1919-22	Human, Mouse	WB
	TNF alpha Rabbit pAb	ER65189	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07]	HA722022	Human, Mouse	WB
Cell treatment	TNF alpha Rabbit pAb	R1203-1	Mouse	WB,IF-Cell,ELISA
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free	HA750897	Human, Mouse	WB
Relative expression	TREM2 Recombinant Rabbit mAb [PSH18-55]	HA724008	Mouse	WB,IP
	VEGF Recombinant Rabbit mAb [SP07-01] - BSA and Azide free	HA750074	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	VEGF Rabbit pAb	ER30607	Human	WB,IF-Cell,IHC-P,FC
	VEGF Recombinant Rabbit mAb [SP07-01]	ET1604-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	VEGFA Recombinant Rabbit mAb [PSH19-17]	HA724051	Human, Mouse, Rat, Monkey	WB
	VEGF-D Rabbit pAb	ER65298	Human, Mouse, Rat, Monkey	WB,IHC-P
	VEGFD Recombinant Rabbit mAb [JE54-75]	ET7110-81	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层