

肿瘤血管生成研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	AKT1 Recombinant Rabbit mAb [PSH0-50] - BSA and Azide free	HA750627	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
	AKT1 Mouse mAb [D9-9-C9]	EM40507	Human, Mouse, Rat	WB,IHC-P
	AKT1 Recombinant Rabbit mAb [PSH0-50]	HA721364	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
KD	AKT1 Recombinant Rabbit mAb [ST05-09] - BSA and Azide free	HA750177	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
KD	AKT1 Recombinant Rabbit mAb [ST05-09]	ET1609-47	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75-09]	HA721870	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	AKT1/2/3 Rabbit pAb	ER1901-33	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	AKT1/2/3 Recombinant Rabbit mAb [ST48-09]	ET1609-51	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
	AKT1/2/3 Recombinant Rabbit mAb [ST48-09] - BSA and Azide free	HA750179	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75-09] - BSA and Azide free	HA750840	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	Akt1/3 Rabbit pAb	ER62638	Human, Mouse, Rat	WB,IHC-P
	Akt2 Mouse mAb [2A1]	RT1027	Human, Mouse, Rat	WB,IP,IF,IHC-P
	AKT2 Rabbit pAb	HA500091	Human, Mouse, Rat	WB,IHC-P
	ANG I Rabbit pAb	ER65619	Human	WB,IHC-P
Relative expression	Angiomotin Recombinant Rabbit mAb [PSH07-46]	HA722834	Human	WB,IHC-P,IF-Cell,IP
Relative expression	Angiomotin Recombinant Rabbit mAb [PSH07-46] - BSA and Azide free	HA751167	Human	WB,IHC-P,IF-Cell,IP
	Angiopoietin 2 Recombinant Rabbit mAb [JM71-34]	ET1705-6	Human, Mouse, Rat	WB,IHC-P,FC
	Angiopoietin 2 Recombinant Rabbit mAb [JM71-34] - BSA and Azide free	HA750436	Human, Mouse, Rat	WB,IHC-P,FC
	Angiopoietin-2 Mouse mAb [4A2]	RT1042	Human, Mouse, Rat	WB,IP,IF,IHC-P
	Angiotensin II type 1A receptor Rabbit pAb	ER1903-69	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Angiotensin II Type 1 Receptor Mouse mAb [11B3]	HA600029	Human, Mouse, Rat	WB,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	ARP-1 Rabbit pAb	ER62729	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	bFGF Rabbit pAb	ER1904-22	Human, Mouse	WB,FC
Relative expression	Cathepsin B Recombinant Rabbit mAb [JA11-02]	ET1704-03	Human	WB,IF-Cell,IF-Tissue,IHC-P
Relative expression	Cathepsin B Rabbit pAb	R1407-1	Human	WB
	Cathepsin B Mouse mAb [J11-A11]	M1506-1	Human, Zebrafish, Mouse, Rat	WB,IF-Cell
	Cathepsin B Rabbit pAb	ER1905-48	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	Cathepsin D Recombinant Rabbit mAb [SU0360]	ET1608-49	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	Cathepsin D Recombinant Mouse mAb [13F4-R]	EM1901-03	Human	WB,IHC-P
	Cathepsin D Rabbit pAb	R1601-2	Human, Mouse, Rat	WB,IHC-P,FC
	Cathepsin D Recombinant Mouse mAb [13F3-R] - BSA and Azide free	HA610001	Human	WB,IHC-P
	Cathepsin D Recombinant Mouse mAb [13F3-R]	EM1901-02	Human	WB,IHC-P
	Cathepsin D Recombinant Mouse mAb [13F4-R] - BSA and Azide free	HA610002	Human	WB,IHC-P
	Cathepsin D Mouse mAb [13F3]	EM1901-16	Human	WB,IHC-P
	Cathepsin D Recombinant Rabbit mAb [JE24-31]	ET7110-58	Human	IF-Cell,IHC-P
	Cathepsin D Rabbit pAb	R1512-9	Human, Mouse	WB,IF-Cell,IHC-P,FC
	Cathepsin D Mouse mAb [13F2]	EM1901-15	Human	WB,IHC-P
	Cathepsin D Mouse mAb [13F4]	EM1901-17	Human	WB,IHC-P
	Cathepsin D Recombinant Rabbit mAb [SU0360] - BSA and Azide free	HA750152	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	Cathepsin D Recombinant mAb	IRS099RB	Human	mIHC
	CBP Rabbit pAb	ER62918	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Cleaved-COL4A3 (L1425) Rabbit pAb	ER60021	Human	WB
	Cleaved-MMP-10 (F99) Rabbit pAb	ER60028	Human	WB
	Cleaved-Notch 2 (A1734) Rabbit pAb	ER60031	Human, Mouse, Rat	WB
	Cleaved-Notch 2 (D1733) Rabbit pAb	ER60032	Human, Mouse, Rat	IF-Cell,WB,IHC-P
	Cleaved-Notch 2 (V1697) Rabbit pAb	ER60033	Human, Mouse, Rat	WB
	Cleaved-Tumstatin (P1426) Rabbit pAb	ER60036	Human, Monkey	WB
	COL4A3 Rabbit pAb	ER63072	Human	IHC-P,IF-Cell
	CRDL2 Rabbit pAb	ER61936	Human, Mouse	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CREBBP Recombinant Rabbit mAb [PS01-16] - BSA and Azide free	HA751014	Human, Mouse	IHC-P
	CREBBP Rabbit pAb	HA500191	Human	WB,IHC-P
	CREBBP Recombinant Rabbit mAb [PS01-16]	HA722297	Human, Mouse	IHC-P
	CREBBP Rabbit pAb	HA500056	Human, Mouse, Rat	WB,IHC-P
	CSF3 Mouse mAb [7F-4E]	EM1706-38	Human	IF-Cell,IHC-P
	EGF Recombinant Rabbit mAb [JA14-15]	ET1704-33	Human, Mouse, Rat	WB,IHC-P,IP
Relative expression, KD	EGFR Recombinant Rabbit mAb [SZ40-19]	ET1603-37	Human, Mouse, Monkey, Rat	WB,IF-Cell,IP,FC,mIHC
	EGFR Rabbit pAb	R1608-3	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	EGFR Recombinant Rabbit mAb	IRS046	Human	mIHC
Relative expression, KO	EGFR Recombinant Rabbit mAb [SP00-86]	ET1604-44	Human, Mouse, Rat, Monkey	WB,IF-Cell,IP,FC
Relative expression, KD	EGFR Recombinant Rabbit mAb [PD00-75]	HA721220	Human, Monkey	WB,IHC-P,IF-Cell,FC
	EGFR Rabbit pAb	ER1512-6	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
	EGFR Mouse mAb [B6-E5-D9]	EM1701-22	Human, Mouse, Monkey, Rat	WB,IHC-P,IF-Cell,FC
Relative expression, KO	EGFR Mouse mAb [7-F2-F8]	EM1901-67	Human, Monkey	WB,IHC-P,FC
	EGFR Rabbit pAb	R1511-18	Human, Monkey	WB,IF-Cell,IHC-P
Relative expression, KD	EGFR Mouse mAb [A8E12]	HA601073	Human, Monkey	WB,IHC-P
Relative expression, KO	EGFR Recombinant Rabbit mAb [SP00-86] - BSA and Azide free	HA750078	Human, Mouse, Rat, Monkey	WB,IF-Cell,IP,FC
	EGFR Rabbit pAb	0407-21	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P
Relative expression, KD	EGFR Recombinant Rabbit mAb [PD00-75] - BSA and Azide free	HA750550	Human, Monkey	WB,IHC-P,IF-Cell,FC
	eIF4E Recombinant Rabbit mAb [JE60-70]	HA721192	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Endostatin Rabbit pAb	ER1908-10	Human, Mouse	WB,IHC-P,IF-Tissue
	EphA1 Rabbit pAb	ER60108	Human, Rat, Mouse	WB,IHC-P
	EphA2/3/4 Rabbit pAb	ER63372	Human, Rat	WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	EphA2 (N-term) Rabbit pAb	ER2001-31	Human, Mouse	WB,IHC-P,FC
	EphA3 Rabbit pAb	ER1908-23	Mouse, Rat	WB,IHC-P,IF-Tissue
	EphA7 Rabbit pAb	ER1908-24	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	EphB2 Rabbit pAb	ER1908-25	Rat, Mouse, Human	WB,IHC-P,IF-Tissue
	EphB4 Rabbit pAb	ER63378	Human, Mouse	WB,IHC-P,IF-Cell
	Eph receptor A2 Recombinant Rabbit mAb [PSH09-21]	HA723079	Human	WB,IF-Cell,IP,FC
	Eph receptor A2 Recombinant Rabbit mAb [PSH09-21] - BSA and Azide free	HA751275	Human	WB,IF-Cell,IP,FC
Relative expression	Eph receptor B2 Recombinant Rabbit mAb [PSH14-91]	HA723693	Human	WB,IF-Cell,IP
Relative expression	Eph receptor B2 Recombinant Rabbit mAb [PSH14-91] - BSA and Azide free	HA751535	Human	WB,IF-Cell,IP
	Ephrin-B1 Rabbit pAb	ER63382	Human, Mouse, Rat, Monkey	WB
	Ephrin-B1/2 Rabbit pAb	ER63383	Human, Mouse, Rat	WB,IHC-P
	Ephrin B2 Recombinant Rabbit mAb [JM53-21] - BSA and Azide free	HA750428	Human, Mouse, Rat	WB,IF-Cell
	Ephrin B2 Recombinant Rabbit mAb [JM53-21]	ET1705-33	Human, Mouse, Rat	WB,IF-Cell
	Ephrin-B3 Rabbit pAb	ER63384	Human, Mouse, Rat	WB,IHC-P
KD	ERK1 Recombinant Rabbit mAb [SP00-74]	ET1604-32	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	ERK1 Recombinant Rabbit mAb [SP05-09]	ET1604-16	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	ERK1 Mouse mAb [7-D6-E5]	M1505-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Erk1/2 Mouse mAb [2G1]	RT1484	Human, Mouse, Rat	WB,IP,IF
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03]	ET1601-29	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03] - BSA and Azide free	HA750019	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
	ERK1/2 Rabbit pAb	ER131011	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	ERK2 Recombinant Mouse mAb [5-D2-R] - BSA and Azide free	HA610080	Human, Mouse, Rat	WB
	ERK2 Mouse mAb [C3B12]	EM1901-55	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Recombinant Mouse mAb [5-D2-R]	HA601202	Human, Mouse, Rat	WB
KD	ERK2 Recombinant Mouse mAb [C3B12-R]	HA601201	Human, Mouse, Rat	WB
	ERK2 Mouse mAb [5-D2]	M1407-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	ERK2 Recombinant Mouse mAb [C3B12-R] - BSA and Azide free	HA610079	Human, Mouse, Rat	WB
KD	ERK2 Recombinant Rabbit mAb [SZ25-01]	ET1603-23	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
KD	ERK2 Mouse mAb [B10-8]	EM1901-53	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Rabbit pAb	ER131218	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	ERK2 Mouse mAb [A10D2]	EM1901-54	Human, Mouse, Rat	WB,IHC-P,FC
	FGF-1 Rabbit pAb	ER63425	Human, Mouse, Rat	IHC-P
	FGF1 Recombinant Rabbit mAb [PSH18-61]	HA724014	Human, Mouse, Rat	WB,FC
	FGF1 Recombinant Rabbit mAb [JE33-93]	HA722695	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	FGF2 Recombinant Rabbit mAb [JM28-10]	ET1703-18	Human, Mouse	WB,IP,FC,IHC-P
	Flk-1/Flt-4 Rabbit pAb	ER63443	Human, Mouse, Rat	IHC-P,IF-Cell
	Flt-4 Rabbit pAb	ER65750	Human, Mouse, Rat	WB,IHC-P
	FRP-2 Rabbit pAb	ER63473	Human, Mouse, Rat	WB
	GPR124 Rabbit pAb	ER63568	Human, Mouse	WB,IF-Cell
	HBEGF Recombinant Rabbit mAb [JE34-79]	HA721470	Human, Rat	WB,IHC-P
	HBEGF Rabbit pAb	ER61976	Human, Rat, Mouse	WB
	HGF Recombinant Rabbit mAb [PSH20-29]	HA724140	Human	IHC-P
	HGF Recombinant Rabbit mAb [PSH03-79] - BSA and Azide free	HA751614	Human	WB
	HGF Recombinant Rabbit mAb [PSH03-79]	HA723896	Human	WB
	HGF Rabbit pAb	ER65684	Human	WB,IHC-P
Cell treatment	HIF-1 alpha Recombinant Rabbit mAb [JE75-33]	HA721997	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	HIF-1 alpha Recombinant Rabbit mAb [PSH07-00]	HA722778	Human, Mouse	WB,IF-Cell,IHC-P
Cell treatment	HIF-1 alpha Recombinant Rabbit mAb [PSH07-00] - BSA and Azide free	HA751127	Human, Mouse	WB,IF-Cell,IHC-P
	HIF-1 alpha Mouse mAb [5H7]	RT1278	Human, Mouse, Rat	WB,IP,ICC,IHC-P
Cell treatment	HIF-1 alpha Rabbit pAb	ER1802-41	Human	WB,IHC-P
KD	HIF-1 beta Recombinant Rabbit mAb [ST05-25]	ET1609-67	Human, Mouse, Rat	WB,IHC-P
KD	HIF-1 beta Recombinant Rabbit mAb [ST05-25] - BSA and Azide free	HA750187	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Angiogenin Recombinant Rabbit mAb [PSH12-64] - BSA and Azide free (Capture)	HA725076	Human	ELISA(Cap)
	Human Angiogenin Recombinant Rabbit mAb [PSH12-65] - BSA and Azide free (Detector)	HA725077	Human	ELISA(Det)
	HRP Conjugated Human Angiogenin Recombinant Rabbit mAb [PSH12-65] - Detector	HA725078H	Human	ELISA(Det),ELISA
	Biotin Conjugated Human Angiogenin Recombinant Rabbit mAb [PSH12-65] - Detector	HA725078B	Human	ELISA(Det),ELISA
	Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-37] - BSA and Azide free (Detector)	HA725008	Human	ELISA(Det)
	Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-36] - BSA and Azide free (Capture)	HA725007	Human	ELISA(Cap)
	Biotin Conjugated Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-37] - Detector	HA725009B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human EGF Recombinant Rabbit mAb [PSH02-20] - Detector	HA722638B	Human	ELISA(Det)
	Human EGF Recombinant Rabbit mAb [PSH02-20] - BSA and Azide free (Detector)	HA721783	Human	ELISA(Det)
	Human EGF Recombinant Rabbit mAb [PSH02-19] - BSA and Azide free (Capture)	HA721782	Human	ELISA(Cap)
	Human Endostatin Recombinant Rabbit mAb [PSH15-82] - BSA and Azide free (Capture)	HA723789	Human	ELISA(Cap)
	HRP Conjugated Human Endostatin Recombinant Rabbit mAb [PSH15-83] - Detector	HA723791H	Human	ELISA(Det),ELISA
	Human Endostatin Recombinant Rabbit mAb [PSH15-83] - BSA and Azide free (Detector)	HA723790	Human	ELISA(Det)
	Biotin Conjugated Human Endostatin Recombinant Rabbit mAb [PSH15-83] - Detector	HA723791B	Human	ELISA(Det),ELISA
	Human EphA2 Recombinant Rabbit mAb [PSH15-13] - BSA and Azide free (Detector)	HA723715	Human	ELISA(Det)
	Human EphA2 Recombinant Rabbit mAb [PSH15-12] - BSA and Azide free (Capture)	HA723714	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human EphA2 Recombinant Rabbit mAb [PSH15-13] - Detector	HA723716B	Human	ELISA(Det),ELISA
	Human FGF-2 Recombinant Rabbit mAb [PSH02-09] - BSA and Azide free (Detector)	HA721781	Human	ELISA(Det)
	Biotin Conjugated Human FGF-2 Recombinant Rabbit mAb [PSH02-09]	HA722639B	Human	ELISA(Det)
	Human FGF-2 Recombinant Rabbit mAb [PSH02-08] - BSA and Azide free (Capture)	HA721780	Human	ELISA(Cap)
	Human G-CSF Recombinant Rabbit mAb [PSH13-12] - BSA and Azide free (Capture)	HA725001	Human	ELISA(Cap)
	Biotin Conjugated Human G-CSF Recombinant Rabbit mAb [PSH13-13] - Detector	HA725003B	Human	ELISA(Det),ELISA
	Human G-CSF Recombinant Rabbit mAb [PSH13-13] - BSA and Azide free (Detector)	HA725002	Human	ELISA(Det)
	Human HGF Recombinant Rabbit mAb [PSH02-06] - BSA and Azide free (Capture)	HA721778	Human	ELISA(Cap)
	Human HGF Recombinant Rabbit mAb [PSH03-79] - BSA and Azide free (Detector)	HA722043	Human	ELISA(Det)
	Human HGF Recombinant Rabbit mAb [PSH02-07] - BSA and Azide free (Detector)	HA721779	Human	ELISA(Det)
	Biotin Conjugated Human HGF Recombinant Rabbit mAb [PSH02-07]	HA722640B	Human	ELISA(Det)
	Biotin Conjugated Human HGF Recombinant Rabbit mAb [PSH03-79]	HA722619B	Human	ELISA(Det)
	Human HGF Recombinant Rabbit mAb [PSH03-78] - BSA and Azide free (Capture)	HA722042	Human	ELISA(Cap)
	Biotin Conjugated Human IL-17A Recombinant Rabbit mAb [PS01-24] - Detector	HA721631	Human	ELISA(Det)
	Human IL-17A Recombinant Rabbit mAb [PS01-24] - BSA and Azide free (Detector)	HA721233	Human	ELISA(Det)
	Human IL-17A Recombinant Rabbit mAb [PS01-20] - BSA and Azide free (Capture)	HA721232	Human	ELISA(Cap)
	HRP Conjugated Human IL-17A Recombinant Rabbit mAb [PS01-24]	HA721296	Human	ELISA(Det)
	Human IL-17F Recombinant Rabbit mAb [PSH18-69] - BSA and Azide free (Detector)	HA725316	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human IL-17F Recombinant Rabbit mAb [PSH18-69] - Detector	HA725316B	Human	ELISA(Det),ELISA
	Human IL-17F Recombinant Rabbit mAb [PSH18-68] - BSA and Azide free (Capture)	HA725315	Human	ELISA(Cap)
	Human IL-18 Recombinant Rabbit mAb [PSH02-08]	HA721767	Human	WB,IHC-P,IF-Tissue,IP
	Human IL-18 Recombinant Rabbit mAb [PSH02-08] - BSA and Azide free	HA750778	Human	WB,IHC-P,IF-Tissue
Relative expression	Human IL-18 Recombinant Rabbit mAb [PSH02-09] - BSA and Azide free	HA750779	Human	WB
	Human IL-18 Recombinant Rabbit mAb [PSH01-71] - BSA and Azide free (Capture)	HA721722	Human	ELISA(Cap)
Relative expression	Human IL-18 Recombinant Rabbit mAb [PSH02-09]	HA721768	Human	WB
	Human IL-18 Recombinant Rabbit mAb [PSH01-70] - BSA and Azide free (Detector)	HA721723	Human	ELISA(Det)
	Biotin Conjugated Human IL-18 Recombinant Rabbit mAb [PSH01-70]	HA722644B	Human	ELISA(Det)
	Human IL-6 Mouse mAb [A6D1]	HA601051	Human	WB,ELISA
	Biotin Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA722233B	Human	ELISA(Det)
	Human IL-6 Mouse mAb [A6D2]	HA601050	Human	WB,ELISA
	Human IL-6 Recombinant Rabbit mAb [PS00-21] - BSA and Azide free (Capture)	HA721167	Human	ELISA(Cap)
	Human IL-6 Recombinant Rabbit mAb [PS00-15] - BSA and Azide free (Detector)	HA721168	Human	ELISA(Det)
	HRP Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA721293	Human	ELISA(Det)
	Human IL-8 Recombinant Rabbit mAb [PS00-83] - BSA and Azide free (Capture)	HA722098	Human	ELISA(Cap)
	Human IL-8 Recombinant Rabbit mAb [PS01-22] - BSA and Azide free (Detector)	HA722099	Human	ELISA(Det)
	Biotin Conjugated Human IL-8 Recombinant Rabbit mAb [PS01-22]	HA722611B	Human	ELISA(Det)
	Biotin Conjugated Human Interferon gamma Recombinant Rabbit mAb [PS01-76] - Detector	HA721576	Human	ELISA(Det),ELISA
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-46] - BSA and Azide free (Capture)	HA721684	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Interferon gamma Recombinant Rabbit mAb [PSH01-45] - BSA and Azide free (Capture)	HA721683	Human	ELISA(Cap)
	Human Interferon gamma Recombinant Rabbit mAb [PS01-76] - BSA and Azide free (Detector)	HA721575	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-47] - BSA and Azide free (Detector)	HA721685	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PS01-75] - BSA and Azide free (Capture)	HA721574	Human	ELISA(Cap)
	Biotin Conjugated Human Interferon gamma Recombinant Rabbit mAb [PSH01-47]	HA722645B	Human	ELISA(Det)
	Human MMP-1 Recombinant Rabbit mAb [PSH07-59] - BSA and Azide free (Capture)	HA722874	Human	ELISA(Cap)
	Biotin Conjugated Human MMP-1 Recombinant Rabbit mAb [PSH07-60]	HA722876B	Human	ELISA(Det)
	Human MMP-1 Recombinant Rabbit mAb [PSH07-60] - BSA and Azide free (Detector)	HA722875	Human	ELISA(Det)
	Human MMP13 Recombinant Rabbit mAb [PSH11-14] - BSA and Azide free (Detector)	HA723294	Human	ELISA(Det)
	Human MMP13 Recombinant Rabbit mAb [PSH11-13] - BSA and Azide free (Capture)	HA723293	Human	ELISA(Cap)
	Biotin Conjugated Human MMP13 Recombinant Rabbit mAb [PSH11-14] - Detector	HA723295B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human MMP-3 Recombinant Rabbit mAb [PSH07-62]	HA722879B	Human	ELISA(Det)
	Human MMP-3 Recombinant Rabbit mAb [PSH07-62] - BSA and Azide free (Detector)	HA722878	Human	ELISA(Det)
	Human MMP-3 Recombinant Rabbit mAb [PSH07-61] - BSA and Azide free (Capture)	HA722877	Human	ELISA(Cap)
	Human MMP-7 Recombinant Rabbit mAb [PSH07-94] - BSA and Azide free (Detector)	HA722927	Human	ELISA(Det)
	Biotin Conjugated Human MMP-7 Recombinant Rabbit mAb [PSH07-94] - Detector	HA722928B	Human	ELISA(Det)
	Human MMP-7 Recombinant Rabbit mAb [PSH07-93] - BSA and Azide free (Capture)	HA722926	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human MMP-9 Recombinant Rabbit mAb [PSH04-62] - BSA and Azide free (Capture)	HA722162	Human	ELISA(Cap)
	Human MMP-9 Recombinant Rabbit mAb [PSH04-63] - BSA and Azide free (Detector)	HA722163	Human	ELISA(Det)
	Human MMP-9 Recombinant Rabbit mAb [PSH04-24] - BSA and Azide free (Detector)	HA722107	Human	ELISA(Det)
	Human MMP-9 Recombinant Rabbit mAb [PSH04-23] - BSA and Azide free (Capture)	HA722106	Human	ELISA(Cap)
	Biotin Conjugated Human MMP-9 Recombinant Rabbit mAb [PSH04-24]	HA722615B	Human	ELISA(Det)
	Human/Mouse/Rat FGF1 Recombinant Rabbit mAb [PSH18-34] - BSA and Azide free (Capture)	HA723987	Human, Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Human/Mouse/Rat FGF1 Recombinant Rabbit mAb [PSH18-35] - Detector	HA723988B	Human, Mouse, Rat	ELISA(Det),ELISA
	Human/Mouse/Rat FGF1 Recombinant Rabbit mAb [PSH18-35] - BSA and Azide free (Detector)	HA723988	Human, Mouse, Rat	ELISA(Det)
	Human/Mouse VEGFC Recombinant Rabbit mAb [PSH10-75] - BSA and Azide free (Capture)	HA723241	Human, Mouse	ELISA(Cap)
	Human/Mouse VEGFC Recombinant Rabbit mAb [PSH10-76] - BSA and Azide free (Detector)	HA723242	Human, Mouse	ELISA(Det)
	Biotin Conjugated Human/Mouse VEGFC Recombinant Rabbit mAb [PSH10-76] - Detector	HA723243B	Human, Mouse	ELISA(Det),ELISA
	Human PDGF-AA Recombinant Rabbit mAb [PSH13-27] - BSA and Azide free (Capture)	HA725061	Human	ELISA(Cap)
	Human PDGF-AA Recombinant Rabbit mAb [PSH13-28] - BSA and Azide free (Detector)	HA725062	Human	ELISA(Det)
	Biotin Conjugated Human PDGF-AA Recombinant Rabbit mAb [PSH13-28] - Detector	HA725063B	Human	ELISA(Det),ELISA
	Human PDGF-BB Recombinant Rabbit mAb [PSH05-99] - BSA and Azide free (Detector)	HA722498	Human	ELISA(Det)
	Biotin Conjugated Human PDGF-BB Recombinant Rabbit mAb [PSH05-99] - Detector	HA722585B	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human PDGF-BB Recombinant Rabbit mAb [PSH05-98] - BSA and Azide free (Capture)	HA722497	Human	ELISA(Cap)
	Human Pentraxin 3 / PTX3 Recombinant Rabbit mAb [PSH11-83] - BSA and Azide free (Capture)	HA723387	Human	ELISA(Cap)
	Biotin Conjugated Human Pentraxin 3 / PTX3 Recombinant Rabbit mAb [PSH11-85] - Detector	HA723391B	Human	ELISA(Det),ELISA
	Human Pentraxin 3 / PTX3 Recombinant Rabbit mAb [PSH11-84] - BSA and Azide free (Detector)	HA723388	Human	ELISA(Det)
	Biotin Conjugated Human Pentraxin 3 / PTX3 Recombinant Rabbit mAb [PSH11-84] - Detector	HA723389B	Human	ELISA(Det),ELISA
	Human Pentraxin 3 / PTX3 Recombinant Rabbit mAb [PSH11-85] - BSA and Azide free (Detector)	HA723390	Human	ELISA(Det)
	Human Plasminogen Recombinant Rabbit mAb [PSH20-22] - BSA and Azide free (Detector)	HA724133	Human	ELISA(Det)
	Human Plasminogen Recombinant Rabbit mAb [PSH20-21] - BSA and Azide free (Capture)	HA724132	Human	ELISA(Cap)
	Biotin Conjugated Human Plasminogen Recombinant Rabbit mAb [PSH20-22] - Detector	HA724133B	Human	ELISA(Det),ELISA
	HRP Conjugated Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - Detector	HA723374H	Human	ELISA(Det)
	Biotin Conjugated Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - BSA and Azide free (Detector)	HA723160B	Human	ELISA(Det),ELISA
	Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - BSA and Azide free (Detector)	HA723159	Human	ELISA(Det)
	Human sTNF RI Recombinant Rabbit mAb [PSH09-89] - BSA and Azide free (Capture)	HA723158	Human	ELISA(Cap)
	Human TIMP1 Recombinant Rabbit mAb [PSH18-78] - BSA and Azide free (Capture)	HA725367	Human	ELISA(Cap)
	Biotin Conjugated Human TIMP1 Recombinant Rabbit mAb [PSH18-79] - Detector	HA725368B	Human	ELISA(Det),ELISA
	Human TIMP1 Recombinant Rabbit mAb [PSH18-79] - BSA and Azide free (Detector)	HA725368	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human TIMP2 Recombinant Rabbit mAb [PSH20-12] - BSA and Azide free (Capture)	HA725387	Human	ELISA(Cap)
	Biotin Conjugated Human TIMP2 Recombinant Rabbit mAb [PSH20-13] - Detector	HA725388B	Human	ELISA(Det),ELISA
	Human TIMP2 Recombinant Rabbit mAb [PSH20-13] - BSA and Azide free (Detector)	HA725388	Human	ELISA(Det)
	Human TNF-beta Recombinant Rabbit mAb [PS01-05] - BSA and Azide free (Capture)	HA721242	Human	ELISA(Cap)
	Human TNF-beta Recombinant Rabbit mAb [PS01-19] - BSA and Azide free (Detector)	HA721243	Human	ELISA(Det)
	Biotin Conjugated Human TNF-beta Recombinant Rabbit mAb [PS01-19]	HA722651B	Human	ELISA(Det)
	Human TNF Receptor II Recombinant Rabbit mAb [PSH11-93] - BSA and Azide free (Detector)	HA723403	Human	ELISA(Det)
	Human TNF Receptor II Recombinant Rabbit mAb [PSH11-94] - BSA and Azide free (Capture)	HA723405	Human	ELISA(Cap)
	Biotin Conjugated Human TNF Receptor II Recombinant Rabbit mAb [PSH11-93] - Detector	HA723404B	Human	ELISA(Det),ELISA
	Human TNF Receptor II Recombinant Rabbit mAb [PSH11-92] - BSA and Azide free (Capture)	HA723402	Human	ELISA(Cap)
	Human VEGFA Recombinant Rabbit mAb [PSH12-91] - BSA and Azide free (Capture)	HA723493	Human	ELISA(Cap)
	Human VEGFA Recombinant Rabbit mAb [PSH12-93] - BSA and Azide free (Detector)	HA723496	Human	ELISA(Det)
	Human VEGFA Recombinant Rabbit mAb [PSH12-92] - BSA and Azide free (Detector)	HA723494	Human	ELISA(Det)
	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-92] - Detector	HA723495B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-93] - Detector	HA723497B	Human	ELISA(Det),ELISA
	Human VEGF Receptor 1 Recombinant Rabbit mAb [PSH19-22] - BSA and Azide free (Capture)	HA725365	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human VEGF Receptor 1 Recombinant Rabbit mAb [PSH19-23] - BSA and Azide free (Detector)	HA725366	Human	ELISA(Det)
	Biotin Conjugated Human VEGF Receptor 1 Recombinant Rabbit mAb [PSH19-23] - BSA and Azide free (Detector)	HA725366B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human VEGF Receptor 2 Recombinant Rabbit mAb [PSH12-84] - Detector	HA725123B	Human	ELISA(Det),ELISA
	Human VEGF Receptor 2 Recombinant Rabbit mAb [PSH12-83] - BSA and Azide free (Capture)	HA725121	Human	ELISA(Cap)
	Human VEGF Receptor 2 Recombinant Rabbit mAb [PSH12-84] - BSA and Azide free (Detector)	HA725122	Human	ELISA(Det)
	Human Von Willebrand Factor Recombinant Rabbit mAb [PSH09-04] - BSA and Azide free (Detector)	HA723059	Human	ELISA(Det)
	HRP Conjugated Human Von Willebrand Factor Recombinant Rabbit mAb [PSH09-04] - Detector	HA723431H	Human	ELISA(Det)
	Human Von Willebrand Factor Recombinant Rabbit mAb [PSH09-03] - BSA and Azide free (Capture)	HA723058	Human	ELISA(Cap)
	Biotin Conjugated Human Von Willebrand Factor Recombinant Rabbit mAb [PSH09-04] - BSA and Azide free (Detector)	HA723060B	Human	ELISA(Det),ELISA
	Id1 Recombinant Rabbit mAb	IRS159RB	Human	mIHC
	Id1 Recombinant Rabbit mAb [JM92-13] - BSA and Azide free	HA750431	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Id1 Recombinant Rabbit mAb [JM92-13]	ET1705-49	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	IGF1R Mouse mAb [D5-G2]	HA600033	Human	WB,FC
	IGF1R Rabbit pAb	R1310-6	Human	WB,IHC-P
Relative expression	IGF1 Receptor beta Recombinant Rabbit mAb [PSH11-72]	HA723368	Human, Mouse, Rat	WB,IF-Cell
Relative expression	IGF1 Receptor beta Recombinant Rabbit mAb [PSH11-72] - BSA and Azide free	HA751427	Human, Mouse, Rat	WB,IF-Cell
	IGF-IR Rabbit pAb	ER63734	Human, Mouse, Rat	IF-Cell,WB,IHC-P
	IL-12 Rabbit pAb	ER1911-22	Rat	WB,IHC-P,IF-Tissue
	IL12 Rabbit pAb	ER1911-21	Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	IL-12A Recombinant Rabbit mAb [PSH16-10] - BSA and Azide free	HA723815	Human	ELISA
	IL-17 Rabbit pAb	ER1706-91	Human	WB,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	IL17 Mouse mAb [4G2]	RT1326	Human, Mouse, Rat	WB,IP,IF
	IL-17A Rabbit pAb	IRS056	Mouse	mIHC
	IL-17A Rabbit pAb	ER1902-37	Human, Mouse, Rat	WB,IHC-P,FC,mIHC
	IL-17F Rabbit pAb	ER65837	Human	IHC-P
Relative expression	IL-18 Recombinant Rabbit mAb [PSH12-04] - BSA and Azide free	HA751511	Mouse	WB,IF-Cell,FC
Relative expression	IL-18 Recombinant Rabbit mAb [PSH12-04]	HA723607	Mouse	WB,IF-Cell,FC
	IL-18 Recombinant Rabbit mAb [PS00-41] - BSA and Azide free	HA750676	Human	WB,IHC-P
	IL-18 Recombinant Rabbit mAb [PS01-06]	HA721537	Human	WB,FC
	IL-18 Recombinant Rabbit mAb [PS00-41]	HA721536	Human	WB,IHC-P
	IL-18 Recombinant Rabbit mAb [PS01-06] - BSA and Azide free	HA750677	Human	WB,FC
Relative expression	IL-18 Recombinant Rabbit mAb [PSH18-15]	HA723978	Human, Mouse, Rat	WB,IHC-P,IP
	IL-6 Rabbit pAb	R1412-2	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	IL-6 Recombinant Rabbit mAb	IRS160RB	Human	mIHC
	IL-6 Mouse mAb [3-A4]	M1303-3	Mouse	WB,IF-Cell,FC,ELISA
	IL-6 Rabbit pAb	R1203-2	Mouse	WB,ELISA
	IL-6 Mouse mAb [1-6]	EM1701-45	Mouse, Rat	WB,IF-Cell
	IL-6 Mouse mAb	EM30301	Human	WB
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-97]	HA723264	Human	WB
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-96] - BSA and Azide free	HA751386	Human	WB,FC
Cell treatment	IL-8 Rabbit pAb	R1511-15	Human	WB,IHC-P,IF-Cell,FC
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-97] - BSA and Azide free	HA751387	Human	WB
	IL-8 Rabbit pAb	ER1706-67	Human	ELISA,IF-Cell,IHC-P
Cell treatment	IL-8 Rabbit pAb	ER1901-61	Human	WB,IHC-P,IF-Cell,FC
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-96]	HA723263	Human	WB,FC
KD, Relative expression	Insulin Receptor Beta Recombinant Rabbit mAb [SN20-13] - BSA and Azide free	HA750276	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC
KD, Relative expression	Insulin Receptor Beta Recombinant Rabbit mAb [SN20-13]	ET1611-90	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	Integrin alpha 2 Recombinant Rabbit mAb [SN0752]	ET1611-57	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	Integrin alpha 2 Recombinant Rabbit mAb [SN0752] - BSA and Azide free	HA750263	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	Integrin alpha 5 Recombinant Rabbit mAb [JJ08-94]	ET1701-58	Human, Mouse, Monkey, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Integrin alpha 6 Recombinant Rabbit mAb [SR45-00]	ET1602-1	Human, Mouse, Rat	WB,IHC-P
	Integrin alpha 6 Recombinant Rabbit mAb [SR45-00] - BSA and Azide free	HA750028	Human, Mouse, Rat	WB,IHC-P
KD	Integrin beta 1 Recombinant Rabbit mAb [SA40-08] - BSA and Azide free	HA750008	Human, Mouse, Rat, Zebrafish	WB,IHC-P
KD	Integrin beta 1 Recombinant Rabbit mAb [SR30-03]	ET1602-10	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Integrin beta 1 Mouse mAb [B10-A5]	EM1901-77	Human, Mouse	WB,IHC-P,FC
KD	Integrin beta 1 Rabbit pAb	ER31001	Human, Mouse	WB,IHC-P,FC
KD	Integrin beta 1 Recombinant Rabbit mAb [SA40-08]	ET1601-17	Human, Mouse, Rat, Zebrafish	WB,IHC-P
	Interferon gamma Recombinant Rabbit mAb [JM10-10]	ET1703-17	Human	WB,FC
	Interferon gamma Recombinant Rabbit mAb [JM10-10] - BSA and Azide free	HA750373	Human	WB,FC
Cell treatment	IP10 Recombinant Rabbit mAb [JA10-82] - BSA and Azide free	HA750411	Human	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	IP10 Recombinant Rabbit mAb [JA10-82]	ET1704-27	Human	WB,IF-Cell,IF-Tissue,IHC-P
	Kallistatin Rabbit pAb	ER65492	Human	WB,IHC-P
	KDR / VEGFR-2 Mouse mAb [1G1]	RT1347	Human, Mouse, Rat	WB,IP,IF,IHC-P
	KLF2 Rabbit pAb	ER1911-98	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	KLK1 Rabbit pAb	ER65582	Human	WB
	LECT1 Rabbit pAb	ER62052	Human, Rat, Mouse	WB
	LOXL4 Rabbit pAb	ER1912-43	Human	WB
	LT-alpha Mouse mAb [C2-F2]	M0807-12	Human	WB
	LT-alpha Mouse mAb [6-C10-10]	M0807-15	Human, Rat	WB
	LT-alpha Rabbit pAb	0807-5	Human, Rat	WB
Relative expression	LT-alpha (C-term) Rabbit pAb	0807-6	Human, Mouse, Rat	WB,IF-Cell
	MFAP5 Rabbit pAb	ER62023	Human, Mouse	WB
	MMP-1 Rabbit pAb	ER31211	Human, Mouse, Monkey	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	MMP-1 Recombinant Rabbit mAb [PSH09-56]	HA723118	Human	WB,IHC-P
Relative expression	MMP-1 Recombinant Rabbit mAb [PSH09-56] - BSA and Azide free	HA751303	Human	WB,IHC-P
	MMP-10 Rabbit pAb	ER63978	Human, Mouse, Rat	IHC-P,IF-Cell,WB
	MMP-11 Recombinant Rabbit mAb [SN74-08] - BSA and Azide free	HA750250	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	MMP-11 Recombinant Rabbit mAb [SN74-08]	ET1611-33	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
KD	MMP-12 Recombinant Rabbit mAb [SR03-23] - BSA and Azide free	HA750052	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	MMP-12 Recombinant Rabbit mAb [SR03-23]	ET1602-42	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MMP-13 Recombinant Rabbit mAb [JF0893]	ET1702-14	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	MMP-13 Recombinant Rabbit mAb [PSH19-21]	HA724055	Mouse	WB,IHC-P
	MMP-14 Recombinant Rabbit mAb [SJ18-09] - BSA and Azide free	HA750108	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	MMP-14 Mouse mAb [11A-1D]	EM1706-30	Human	IHC-P,FC
	MMP-14 Recombinant Rabbit mAb [SJ18-09]	ET1606-48	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	MMP-15 Rabbit pAb	ER1913-06	Mouse	WB
	MMP-16 Rabbit pAb	ER64067	Human, Mouse, Rat	WB,IHC-P
	MMP-17 Recombinant Rabbit mAb [JU33-03]	ET7106-84	Human, Mouse, Rat	IF-Cell,IF-Tissue,IHC-P,WB
	MMP-19 Rabbit pAb	ER1913-07	Human, Mouse	WB
	MMP-3 Recombinant Rabbit mAb [JM46-22]	ET1705-98	Human, Mouse, Rat	WB,IHC-P
	MMP-3 Rabbit pAb	ER1706-77	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	MMP-7 Recombinant Rabbit mAb [PSH08-91]	HA723043	Human	WB,IF-Cell,FC,IHC-P
	MMP-7 Recombinant Rabbit mAb [PSH08-91] - BSA and Azide free	HA751256	Human	WB,IF-Cell,FC,IHC-P
	MMP-7 Rabbit pAb	ER1913-08	Human, Mouse	WB,IHC-P,IF-Tissue
	MMP-7 Rabbit pAb	HA500305	Human	WB,IF-Cell,IHC-P,FC
Cell treatment	MMP-9 Recombinant Rabbit mAb [PSH06-33] - BSA and Azide free	HA751070	Human	WB,IHC-P
	MMP-9 Mouse mAb [2G3]	RT1401	Human	WB,IP,IF,IHC-P
	MMP-9 Recombinant Mouse mAb [10A1]	EM1801-22	Human	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	MMP-9 Recombinant Rabbit mAb [PSH06-33]	HA722665	Human	WB,IHC-P
Cell treatment, KD	MMP-9 Recombinant Rabbit mAb [JA80-73]	ET1704-69	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC
	MMP-9 Recombinant Rabbit mAb	IRS172RB	Human	mIHC
Cell treatment, KD	MMP-9 Recombinant Rabbit mAb [JA80-73] - BSA and Azide free	HA750420	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC
	Mouse G-CSF Recombinant Rabbit mAb [PSH15-74] - BSA and Azide free (Capture)	HA723777	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse G-CSF Recombinant Rabbit mAb [PSH15-75] - Detector	HA723779B	Mouse	ELISA(Det),ELISA
	Mouse G-CSF Recombinant Rabbit mAb [PSH15-75] - BSA and Azide free (Detector)	HA723778	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IFN gamma Recombinant Rabbit mAb [PSH14-81] - Detector	HA723678B	Mouse	ELISA(Det),ELISA
	Mouse IFN gamma Recombinant Rabbit mAb [PSH03-20] - BSA and Azide free (Capture)	HA721956	Mouse	ELISA(Cap)
	Mouse IFN gamma Recombinant Rabbit mAb [PSH14-80] - BSA and Azide free (Capture)	HA723676	Mouse	ELISA(Cap)
	Mouse IFN gamma Recombinant Rabbit mAb [PSH14-81] - BSA and Azide free (Detector)	HA723677	Mouse	ELISA(Det)
	Mouse IL-17A Recombinant Rabbit mAb [PS01-92] - BSA and Azide free (Detector)	HA721607	Mouse	ELISA(Det)
	Mouse IL-17A Recombinant Rabbit mAb [PSH05-35] - BSA and Azide free (Capture)	HA722264	Mouse, Rat	ELISA(Cap)
	Mouse IL-17A Recombinant Rabbit mAb [PSH05-36] - BSA and Azide free (Detector)	HA722265	Mouse, Rat	ELISA(Det)
	Mouse IL-17A Recombinant Rabbit mAb [PS01-91] - BSA and Azide free (Capture)	HA721606	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-17A Recombinant Rabbit mAb [PS01-92] - Detector	HA721608	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse IL-17A Recombinant Rabbit mAb [PSH05-36]	HA722600B	Mouse, Rat	ELISA(Det),ELISA
	Mouse IL-17F Recombinant Rabbit mAb [PSH17-14] - BSA and Azide free (Capture)	HA723901	Mouse	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse IL-17F Recombinant Rabbit mAb [PSH17-15] - Detector	HA723902B	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-17F Recombinant Rabbit mAb [PSH17-16] - Detector	HA723903B	Mouse	ELISA(Det)
	Mouse IL-17F Recombinant Rabbit mAb [PSH17-16] - BSA and Azide free (Detector)	HA723903	Mouse	ELISA(Det)
	Mouse IL-17F Recombinant Rabbit mAb [PSH17-15] - BSA and Azide free (Capture/Detector)	HA723902	Mouse	ELISA(Det),ELISA(Cap)
	Mouse IL-18 Recombinant Rabbit mAb [PSH12-03] - BSA and Azide free (Capture)	HA723416	Mouse	ELISA(Cap)
	Mouse IL-18 Recombinant Rabbit mAb [PSH12-04] - BSA and Azide free (Detector)	HA723417	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-18 Recombinant Rabbit mAb [PSH12-04] - Detector	HA723418B	Mouse	ELISA(Det),ELISA
	HRP Conjugated Mouse IL-6 Mouse mAb [3-A4]	M1303-2	Mouse	WB,ELISA
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PSH0-60] - Detector	HA722647B	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PS01-83] - BSA and Azide free (Detector)	HA721601	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-59] - BSA and Azide free (Detector)	HA721412	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PS01-74] - BSA and Azide free (Capture/Detector)	HA721600	Mouse	ELISA(Cap),ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-58] - BSA and Azide free (Capture)	HA721411	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PS01-74] - Detector	HA723643B	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-60] - BSA and Azide free (Detector)	HA721413	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PS01-83] - BSA and Azide free (Detector)	HA721602	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PSH0-59] - Detector	HA722646B	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IL-6 Recombinant Rabbit mAb [PSH14-47] - BSA and Azide free (Capture)	HA723642	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-80] - Detector	HA723786B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - Detector	HA723788B	Mouse	ELISA(Det),ELISA
	HRP Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - Detector	HA723788H	Mouse	ELISA(Det),ELISA
	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-79] - BSA and Azide free (Capture)	HA723784	Mouse	ELISA(Cap)
	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - BSA and Azide free (Detector)	HA723787	Mouse	ELISA(Det)
	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-80] - BSA and Azide free (Capture/Detector)	HA723785	Mouse	ELISA(Det),ELISA(Cap)
	Mouse/Rat HGF Recombinant Rabbit mAb [PSH18-33] - BSA and Azide free (Detector)	HA723986	Mouse, Rat	ELISA(Det)
	Mouse/Rat HGF Recombinant Rabbit mAb [PSH18-32] - BSA and Azide free (Capture)	HA723985	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat HGF Recombinant Rabbit mAb [PSH18-33] - Detector	HA723986B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat IL-18 Recombinant Rabbit mAb [PSH18-85] - BSA and Azide free (Detector)	HA724022	Mouse, Rat	ELISA(Det)
	Mouse/Rat IL-18 Recombinant Rabbit mAb [PSH18-84] - BSA and Azide free (Capture)	HA724021	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat IL-18 Recombinant Rabbit mAb [PSH18-85] - Detector	HA724022B	Mouse, Rat	ELISA(Det),ELISA
	Biotin Conjugated Mouse/Rat PDGF-BB Recombinant Rabbit mAb [PSH16-06] - Detector	HA723811B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat PDGF-BB Recombinant Rabbit mAb [PSH16-06] - BSA and Azide free (Detector)	HA723811	Mouse, Rat	ELISA(Det)
	Mouse/Rat PDGF-BB Recombinant Rabbit mAb [PSH16-05] - BSA and Azide free (Capture)	HA723810	Mouse, Rat	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse/Rat PDGF-BB Recombinant Rabbit mAb [PSH16-07] - BSA and Azide free (Capture)	HA723812	Mouse, Rat	ELISA(Cap)
	Mouse/Rat VEGFR3 / FLT4 Recombinant Rabbit mAb [PSH19-86] - BSA and Azide free (Capture)	HA724110	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat VEGFR3 / FLT4 Recombinant Rabbit mAb [PSH19-87] - Detector	HA724111B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat VEGFR3 / FLT4 Recombinant Rabbit mAb [PSH19-87] - BSA and Azide free (Detector)	HA724111	Mouse, Rat	ELISA(Det)
	Mouse VEGF Recombinant Rabbit mAb [PSH02-73] - BSA and Azide free (Detector)	HA721851	Mouse	ELISA(Det)
	Biotin Conjugated Mouse VEGF Recombinant Rabbit mAb [PSH02-73]	HA722636B	Mouse	ELISA(Det)
	Mouse VEGF Recombinant Rabbit mAb [PSH02-72] - BSA and Azide free (Capture)	HA721850	Mouse	ELISA(Cap)
	NDRG1 Recombinant Rabbit mAb [JM32-02]	ET1705-20	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Relative expression	Nestin Recombinant Antibody [PSH07-89] - Guinea pig IgG2 (Chimeric)	HA601468	Human, Mouse, Rat	IHC-Fr,WB
KD	Nestin Recombinant Rabbit mAb [SN06-27]	ET1611-21	Human	IF-Cell,IF-Tissue,IHC-P,WB
	Nestin Recombinant Rabbit mAb	IRS065	Mouse	mIHC
Relative expression	Nestin Recombinant Antibody [PSH07-89] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610302	Human, Mouse, Rat	IHC-Fr,WB
Relative expression, KD	Nestin Recombinant Rabbit mAb [PSH07-89]	HA722919	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IF-Tissue,mIHC,IP
	Nestin Recombinant Antibody [PSH07-89] - Rat IgG1 (Chimeric)	HA601406	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell
	Nestin Rabbit pAb	R1510-19	Human, Mouse	IHC-P
Relative expression, KD	Nestin Recombinant Rabbit mAb [PSH07-89] - BSA and Azide free	HA751196	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IF-Tissue,IP
	Nestin Recombinant Antibody [PSH07-89] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610246	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell
KD	NF-kB p100 / NFkB2 Recombinant Rabbit mAb [JM82-03]	ET1705-45	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
KD	NF-kB p100 / NFkB2 Mouse mAb [16B3]	EM1901-79	Human, Mouse, Rat	WB,IHC-P,FC
KD	NF-kB p100 / NFkB2 Mouse mAb [16B2]	EM1901-78	Human	WB,IHC-P,FC
	NF-kB p105 / p50 Mouse mAb [A2-7]	M1505-8	Human, Mouse	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	NF-kB p105 / p50 Recombinant Rabbit mAb [SZ20-01]	ET1603-18	Human, Mouse, Rat	WB,IHC-P,ICC,FC,IF-Tissue
KD	NF-kB p105 / p50 Recombinant Rabbit mAb [SZ20-01] - BSA and Azide free	HA750060	Human, Mouse, Rat	WB,IHC-P,ICC,FC,IF-Tissue
	NF-kB p105 / p50 Rabbit pAb	R1309-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	Notch 1 Recombinant Rabbit mAb [SJ205]	ET1606-55	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	Notch 1 Recombinant Rabbit mAb [SJ205] - BSA and Azide free	HA750111	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Notch2 Recombinant Rabbit mAb [PSH0-54] - BSA and Azide free	HA750631	Human, Mouse, Rat	WB
	Notch2 Recombinant Rabbit mAb [PSH0-55] - BSA and Azide free	HA750632	Human, Mouse, Rat	WB
	Notch2 Recombinant Rabbit mAb [PSH0-54]	HA721401	Human, Mouse, Rat	WB
	Notch2 Recombinant Rabbit mAb [PSH0-55]	HA721402	Human, Mouse, Rat	WB
	NOTCH3 Rabbit pAb	HA500519	Human	WB,FC
	NOTCH4 Recombinant Rabbit mAb [JM101-05]	ET1703-20	Mouse	WB,IP
	PAI-3 / Serpin A5 Rabbit pAb	0812-5	Human	WB
	PDGF-A Rabbit pAb	ER1914-91	Human, Mouse	WB,IHC-P,IF-Tissue,FC
	PDGF B Recombinant Rabbit mAb [JE75-16]	HA722003	Human, Mouse, Rat	WB,IHC-P
	PDGF-B Rabbit pAb	ER1914-92	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	PDGF-D Rabbit pAb	ER65537	Human	WB,IHC-P
KD	PDGFR alpha Recombinant Rabbit mAb [JF104-6]	ET1702-49	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr
	PDGFR beta Rabbit pAb	R1510-44	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	PDGFR beta Recombinant Rabbit mAb [SY10-08]	ET1605-20	Human, Mouse, Rat	WB,IHC-P,IP,FC,IF-Tissue,IF-Cell,mIHC,IHC-Fr
	PDGFR beta Recombinant Antibody - Rat IgG1 (Chimeric)	IRS283RT	Mouse	mIHC
	PDGFR beta Recombinant Rabbit mAb	IRS038	Human	mIHC
	PDGFR beta Recombinant Rabbit mAb	IRS058	Mouse	mIHC
	PDGFR beta Recombinant Rabbit mAb [SY10-08] - BSA and Azide free	HA750082	Human, Mouse, Rat	WB,IHC-P,IP,FC,IF-Tissue,IF-Cell,IHC-Fr
	PEDF Rabbit pAb	ER1914-97	Human, Mouse	WB,IHC-P,IF-Tissue
	PEDF Recombinant Rabbit mAb [JE41-62]	HA721661	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-AKT1 (S124) Recombinant Rabbit mAb [JJ08-46]	ET1701-36	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-AKT1 (S129) Recombinant Rabbit mAb [JE59-50]	HA721677	Human, Mouse, Rat	WB
Cell treatment	Phospho-AKT1 (T450) Recombinant Rabbit mAb [SD08-12]	ET1612-73	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell,FC
Cell treatment	Phospho-AKT2 (S474) Recombinant Rabbit mAb [PSH09-46]	HA723108	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-AKT2 (S474) Rabbit pAb	HA500116	Human, Mouse	WB,IHC-P
Cell treatment	Phospho-AKT2 (S474) Recombinant Rabbit mAb [PSH09-46] - BSA and Azide free	HA751293	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44]	HA722129	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05]	ET1607-73	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44] - BSA and Azide free	HA750953	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05] - BSA and Azide free	HA750132	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10] - BSA and Azide free	HA751205	Human	WB,IF-Cell,FC
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10]	HA722951	Human	WB,IF-Cell,FC
	Phospho-EGFR (S695) Recombinant Rabbit mAb [SN07-36]	ET1611-50	Human	WB,IF-Cell,IF-Tissue
Cell treatment	Phospho-EGFR (T669) Recombinant Rabbit mAb [PSH10-03] - BSA and Azide free	HA751333	Human, Mouse, Rat	WB,IP
Cell treatment	Phospho-EGFR (T669) Recombinant Rabbit mAb [PSH10-03]	HA723172	Human, Mouse, Rat	WB,IP
Cell treatment	Phospho-EGFR (Y1045) Recombinant Rabbit mAb [PSH03-28] - BSA and Azide free	HA750868	Human	WB,IF-Cell,FC
Cell treatment	Phospho-EGFR (Y1045) Recombinant Rabbit mAb [PSH03-28]	HA721970	Human	WB,IF-Cell,FC
Cell treatment	Phospho-EGFR (Y1068) Recombinant Rabbit mAb [SD2055]	ET1612-30	Human	WB,IF-Cell,IF-Tissue,IP,FC
Cell treatment	Phospho-EGFR (Y1068) Recombinant Rabbit mAb [SD2055] - BSA and Azide free	HA750290	Human	WB,IF-Cell,IF-Tissue,IP,FC
Cell treatment	Phospho-EGFR (Y1092) Recombinant Rabbit mAb [SJ0194]	ET1606-44	Human, Mouse	WB,IF-Tissue,IHC-P
	Phospho-EGFR (Y1148) Rabbit pAb	HA500580	Human	WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-EGFR (Y1173) Recombinant Rabbit mAb [SC57-04]	ET1610-4	Human, Mouse	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-EGFR (Y845) Recombinant Rabbit mAb [JE46-77]	HA721488	Human	WB
Cell treatment	Phospho-eIF4E (S209) Recombinant Rabbit mAb [SU0396]	ET1608-66	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
Cell treatment	Phospho-eIF4E (S209) Recombinant Rabbit mAb [SU0396] - BSA and Azide free	HA750160	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4]	ET1603-22	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4] - BSA and Azide free	HA750061	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01]	ET1610-13	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01] - BSA and Azide free	HA750199	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-NDRG1 (T346) Recombinant Rabbit mAb [PSH07-17]	HA722805	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue
Cell treatment	Phospho-NDRG1 (T346) Recombinant Rabbit mAb [PSH07-17] - BSA and Azide free	HA751138	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue
Cell treatment	Phospho-PDGFR beta (Y1009) Recombinant Rabbit mAb [JE46-72]	HA721736	Human, Mouse, Rat	WB
Cell treatment	Phospho-PDGFR beta (Y1021) Recombinant Rabbit mAb [JE48-54]	HA721833	Mouse, Rat, Human	WB
Cell treatment	Phospho-PDGFR beta (Y740) Recombinant Rabbit mAb [SN0646]	ET1611-29	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue
Cell treatment	Phospho-PDGFR beta (Y751) Recombinant Rabbit mAb [PSH13-61]	HA723550	Mouse, Rat, Human	WB
Cell treatment	Phospho-PDGFR beta (Y751) Recombinant Rabbit mAb [PSH13-61] - BSA and Azide free	HA751480	Mouse, Rat, Human	WB
	Plasminogen Rabbit pAb	ER1915-13	Human	WB,IHC-P,IF-Tissue
	PLGF Mouse mAb [3G2]	EM1701-90	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	PLGF Mouse mAb [3G3]	EM1701-88	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	PLGF Recombinant Rabbit mAb [JA63-15]	ET1704-99	Human, Mouse, Rat	WB,IHC-P
	PLGF Rabbit pAb	ER1706-88	Human, Mouse	IHC-P,FC
	Ptx3 Rabbit pAb	ER64704	Human, Mouse	WB,IHC-P,IF-Cell
	Biotin Conjugated Rat IFN gamma Recombinant Rabbit mAb [PSH05-70]	HA722589B	Rat	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Rat IFN gamma Recombinant Rabbit mAb [PSH05-69] - BSA and Azide free (Capture)	HA722421	Rat	ELISA(Cap)
	Rat IFN gamma Recombinant Rabbit mAb [PS01-97] - BSA and Azide free (Detector)	HA721619	Rat	ELISA(Det)
	Rat IFN gamma Recombinant Rabbit mAb [PS01-96] - BSA and Azide free (Capture)	HA721618	Rat	ELISA(Cap)
	Biotin Conjugated Rat IFN gamma Recombinant Rabbit mAb [PS01-97] - Detector	HA721620	Rat	ELISA(Det)
	Rat IFN gamma Recombinant Rabbit mAb [PSH05-70] - BSA and Azide free (Detector)	HA722422	Rat	ELISA(Det)
	Rat IL-6 Recombinant Rabbit mAb [PS01-90] - BSA and Azide free (Capture)	HA721615	Rat	ELISA(Cap)
	Rat IL-6 Recombinant Rabbit mAb [PS01-94] - BSA and Azide free (Detector)	HA721616	Rat	ELISA(Det)
	Biotin Conjugated Rat IL-6 Recombinant Rabbit mAb [PS01-94] - BSA and Azide free (Detector)	HA721617	Rat	ELISA(Det)
	ROBO1 Recombinant Mouse mAb [10E2-R] - BSA and Azide free	HA610007	Human, Mouse, Rat	WB,IHC-P
	ROBO1 Rabbit pAb	ER1706-24	Human, Mouse	IF-Cell,IHC-P
	ROBO1 Mouse mAb [10E2]	EM1901-72	Human, Mouse, Rat	WB,IHC-P
	ROBO1 Recombinant Mouse mAb [10E2-R]	HA601056	Human, Mouse, Rat	WB,IHC-P
	ROBO4 Rabbit pAb	ER61895	Human, Rat, Mouse	WB
	Slit2 Recombinant Rabbit mAb [JE31-87]	HA721407	Human, Mouse, Rat	WB,IF-Cell,FC
	SLIT2 Rabbit pAb	ER62446	Human, Mouse	WB
	SPATA13 Rabbit pAb	ER1902-32	Human, Rat, Mouse	WB,IHC-P,FC
	TGF alpha Recombinant Rabbit mAb [JB58-36]	ET7107-40	Human	WB,IF-Tissue,IHC-P,IP
	TGF beta 1 Recombinant Rabbit mAb [PSH18-59]	HA724012	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
	Biotin Conjugated TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - Detector	HA723727B	Human, Mouse, Rat	ELISA(Det),ELISA
	TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - BSA and Azide free (Detector)	HA723726	Human, Mouse, Rat	ELISA(Det)
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17] - BSA and Azide free	HA750485	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	TGF beta 1 Recombinant Rabbit mAb [PSH15-21] - BSA and Azide free (Capture)	HA723725	Human, Mouse, Rat	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17]	HA721143	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	TGF beta 1 Recombinant Rabbit mAb	IRS205RB	Human	mIHC
	TGF beta 1 Rabbit pAb	HA500496	Human	WB
	TGF beta 1 Mouse mAb [B2-F6]	EM1709-93	Human	WB,IHC-P,FC
	TGF beta 2 Rabbit pAb	ER1917-63	Human, Mouse	WB,IHC-P,IF-Tissue
	TGF beta 2 Rabbit pAb	ER1917-62	Human, Mouse, Rat	WB
	TGF beta 3 Recombinant Rabbit mAb [PSH09-26] - BSA and Azide free	HA751280	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	TGF beta 3 Rabbit pAb	ER1917-64	Human, Mouse, Rat, Pig	WB,IHC-P,IF-Tissue,FC
	TGF beta 3 Recombinant Rabbit mAb [PSH09-26]	HA723084	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Thrombospondin 1 Recombinant Rabbit mAb [PSH02-97]	HA721916	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
	Thrombospondin 1 Recombinant Rabbit mAb [PSH02-97] - BSA and Azide free	HA750852	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
	Thrombospondin 2 Rabbit pAb	ER1917-73	Human	WB
	TIMP-1 Mouse mAb [4C2]	EM2001-03	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	TIMP1 Recombinant Rabbit mAb [JE59-72]	HA721741	Human	WB,IHC-P,FC
Relative expression	TIMP2 Recombinant Rabbit mAb [PSH11-62]	HA723358	Human, Monkey	WB
Relative expression	TIMP2 Recombinant Rabbit mAb [PSH11-62] - BSA and Azide free	HA751418	Human, Monkey	WB
Relative expression	TIMP3 Recombinant Rabbit mAb [PSH04-41] - BSA and Azide free	HA750950	Human, Rat	WB,IF-Cell
Relative expression	TIMP3 Recombinant Rabbit mAb [PSH04-41]	HA722126	Human, Rat	WB,IF-Cell
	TIMP4 Rabbit pAb	ER65170	Human, Mouse	WB
	TNFR1 Mouse mAb [2G2]	RT1624	Human, Mouse, Rat	WB,IP,IF,IHC-P
	TNF Receptor I Rabbit pAb	HA500140	Human	WB,IHC-P,FC
	TNF Receptor II Recombinant Rabbit mAb [JM113-01]	ET1703-63	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP,IHC-Fr
	TNF Receptor II Rabbit pAb	HA500332	Human, Mouse, Rat	IHC-P,FC
Relative expression	TNF Receptor II Recombinant Rabbit mAb [PSH12-12]	HA723428	Human	FC
	TNF Receptor II Recombinant Rabbit mAb [JM113-01] - BSA and Azide free	HA750386	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	TNF Receptor II Recombinant Rabbit mAb [PSH12-12] - BSA and Azide free	HA751441	Human	FC
	VEGF Recombinant Rabbit mAb [SP07-01] - BSA and Azide free	HA750074	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	VEGF Rabbit pAb	ER30607	Human	WB,IF-Cell,IHC-P,FC
	VEGF Recombinant Rabbit mAb [SP07-01]	ET1604-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	VEGFA Recombinant Rabbit mAb [PSH19-17]	HA724051	Human, Mouse, Rat, Monkey	WB
	VEGFB167 Rabbit pAb	ER1918-50	Rat	WB
	VEGFC Rabbit pAb	R1307-4	Human	WB
	VEGF-D Rabbit pAb	ER65298	Human, Mouse, Rat, Monkey	WB,IHC-P
	VEGFD Recombinant Rabbit mAb [JE54-75]	ET7110-81	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	VEGF Receptor 1 Recombinant Rabbit mAb [SY09-09]	ET1605-11	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	VEGF Receptor 1 Rabbit pAb	ER50803	Human, Mouse, Rat	WB,IHC-P
KD	VEGF Receptor 1 Recombinant Rabbit mAb [SY09-09] - BSA and Azide free	HA750081	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	VEGF Receptor 1 Rabbit pAb	RT1651	Human, Mouse, Rat	WB,IP,IF,IHC-P
Relative expression	VEGF Receptor 2 Recombinant Rabbit mAb [SU03-42]	ET1608-33	Human, Rat	WB,IF-Cell,IHC-P,IP
Relative expression	Von Willebrand Factor Recombinant Rabbit mAb [PSH07-45] - BSA and Azide free	HA751166	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP
Relative expression	Von Willebrand Factor Recombinant Rabbit mAb [PSH07-45]	HA722833	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP
	Von Willebrand Factor Recombinant Rabbit mAb	IRS244RB	Mouse	mIHC
	Von Willebrand Factor Recombinant Rabbit mAb [JJ09-34]	ET1701-87	Human, Mouse	WB
	Von Willebrand Factor Recombinant Rabbit mAb	IRS215RB	Human	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层