

2025年3月—新上市产品清单

高级验证	产品名	货号	反应种属	验证应用
KD, Relative expression	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Rat IgG1 (Chimeric)	HA601427	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Caspase-4 Mouse mAb [A9G2]	HA601431	Human, Cynomolgus monkey	WB
	CD3 Recombinant Rabbit mAb	IRS074	Mouse	mIHC
	CD31 Recombinant Rabbit mAb	IRS075	Mouse	mIHC
	CD45R Rat mAb	IRS076	Mouse	mIHC
	CD8a Mouse mAb	IRS007	Human	mIHC
Cell treatment	DPP4 / CD26 Recombinant Rabbit mAb [PSH12-59]	HA723473	Human	WB,IHC-P,IF-Cell,FC,IP
Relative expression	ENPP1 / PC1 Recombinant Rabbit mAb [PSH12-11]	HA723427	Human	WB,IHC
	FCN1 / M-Ficolin Rabbit pAb	HA500567	Human, Mouse	WB,IF-Cell
	GFP Nanoselector Agarose	019-101-002	Species independent	IP,ChIP,MS,Purification
	Glutamine Synthetase Recombinant Antibody [PSH10-48] - Rat IgG1 (Chimeric)	HA601428	Human, Mouse, Rat	IHC-Fr,WB
	HA tag Nanoselector Agarose	003-101-002	Species independent	IP,ChIP,MS,Purification
Relative expression	HDAC11 Rabbit pAb	HA500568	Human, Mouse, Rat	WB
	HSPB2 Rabbit pAb	ER1919-63	Human, Mouse, Rat	WB,IHC-P
	Human Adiponectin enzyme-linked immunoassay kit	EH0099	Human	ELISA
	Biotin Conjugated Human Alpha-Synuclein Recombinant Rabbit mAb [PSH12-56] - Detector	HA723469B	Human	ELISA(Det),ELISA
	Human Alpha-Synuclein Recombinant Rabbit mAb [JB42-33] - BSA and Azide free (Capture)	HA723467	Human	ELISA(Cap)
	Human Alpha-Synuclein Recombinant Rabbit mAb [PSH12-56] - BSA and Azide free (Detector)	HA723468	Human	ELISA(Det)
	Human Angiogenin Recombinant Rabbit mAb [PSH12-65] - BSA and Azide free (Detector)	HA725077	Human	ELISA(Det)

	Biotin Conjugated Human Angiogenin Recombinant Rabbit mAb [PSH12-65] - Detector	HA725078B	Human	ELISA(Det),ELISA
	Human Angiogenin Recombinant Rabbit mAb [PSH12-64] - BSA and Azide free (Capture)	HA725076	Human	ELISA(Cap)
	Human CCL20/MIP-3 alpha Recombinant Rabbit mAb [PSH12-33] - BSA and Azide free (Detector)	HA725011	Human	ELISA(Det)
	Biotin Conjugated Human CCL20/MIP-3 alpha Recombinant Rabbit mAb [PSH12-33] - Detector	HA725012B	Human	ELISA(Det),ELISA
	Human CCL20/MIP-3 alpha Recombinant Rabbit mAb [PSH12-32] - BSA and Azide free (Capture)	HA725010	Human	ELISA(Cap)
	Human CD28 Recombinant Rabbit mAb [PSH12-50] - BSA and Azide free (Capture)	HA723457	Human	ELISA(Cap)
	Biotin Conjugated Human CD28 Recombinant Rabbit mAb [PSH12-51] - Detector	HA723459B	Human	ELISA(Det),ELISA
	Human CD28 Recombinant Rabbit mAb [PSH12-51] - BSA and Azide free (Detector)	HA723458	Human	ELISA(Det)
	HRP Conjugated Human Complement C5 Recombinant Rabbit mAb [PSH06-46] - Detector	HA725127H	Human	ELISA(Det)
	Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-36] - BSA and Azide free (Capture)	HA725007	Human	ELISA(Cap)
	Biotin Conjugated Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-37] - Detector	HA725009B	Human	ELISA(Det),ELISA
	Human CXCL10 / IP-10 Recombinant Rabbit mAb [PSH12-37] - BSA and Azide free (Detector)	HA725008	Human	ELISA(Det)
	Human CXCL9 Recombinant Rabbit mAb [PSH12-35] - BSA and Azide free (Detector)	HA725005	Human	ELISA(Det)
	Biotin Conjugated Human CXCL9 Recombinant Rabbit mAb [PSH12-35] - Detector	HA725006B	Human	ELISA(Det),ELISA
	Human CXCL9 Recombinant Rabbit mAb [PSH12-34] - BSA and Azide free (Capture)	HA725004	Human	ELISA(Cap)
	Human GTR/TNFRSF18 enzyme-linked immunoassay kit	EH0095	Human	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HRP Conjugated Human ICAM1 Recombinant Rabbit mAb [PSH10-16] - Detector	HA723474H	Human	ELISA(Det),ELISA
	Human IFN-gamma R1 enzyme-linked immunoassay kit	EH0104	Human	ELISA
	Human IL-13 Recombinant Rabbit mAb [PSH12-39] - BSA and Azide free (Detector)	HA725014	Human	ELISA(Det)
	Biotin Conjugated Human IL-13 Recombinant Rabbit mAb [PSH12-39] - Detector	HA725015B	Human	ELISA(Det),ELISA
	Human IL-13 Recombinant Rabbit mAb [PSH12-38] - BSA and Azide free (Capture)	HA725013	Human	ELISA(Cap)
	Human IL-27 enzyme-linked immunoassay kit	EH0094	Human	ELISA
	Human LIF enzyme-linked immunoassay kit	EH0096	Human	ELISA
	Human Livin Recombinant Rabbit mAb [PSH12-96] - BSA and Azide free (Capture)	HA723501	Human	ELISA(Cap)
	Biotin Conjugated Human Livin Recombinant Rabbit mAb [PSH12-95] - Detector	HA723500B	Human	ELISA(Det),ELISA
	Human Livin Recombinant Rabbit mAb [PSH12-95] - BSA and Azide free (Detector)	HA723499	Human	ELISA(Det)
	Human Livin Recombinant Rabbit mAb [PSH12-94] - BSA and Azide free (Capture)	HA723498	Human	ELISA(Cap)
	Human MMP-8 enzyme-linked immunoassay kit	EH0103	Human	ELISA
	Human Nectin-2 Recombinant Rabbit mAb [PSH12-52] - BSA and Azide free (Detector)	HA723461	Human	ELISA(Det)
	Human Nectin-2 Recombinant Rabbit mAb [PSH12-53] - BSA and Azide free (Capture)	HA723463	Human	ELISA(Cap)
	Biotin Conjugated Human Nectin-2 Recombinant Rabbit mAb [PSH12-52] - Detector	HA723462B	Human	ELISA(Det),ELISA
	Human Nectin-2 Recombinant Rabbit mAb [PSH11-54] - BSA and Azide free (Capture)	HA723460	Human	ELISA(Cap)
	Biotin Conjugated Human Nectin-4 Recombinant Rabbit mAb [PSH12-58] - Detector	HA723472B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Nectin-4 Recombinant Rabbit mAb [PSH12-57] - BSA and Azide free (Capture)	HA723470	Human	ELISA(Cap)
	Human Nectin-4 Recombinant Rabbit mAb [PSH12-58] - BSA and Azide free (Detector)	HA723471	Human	ELISA(Det)
	Human Oxytocin-neurophysin 1 Recombinant Rabbit mAb [PSH12-29] - BSA and Azide free (Detector)	HA723453	Human	ELISA(Det)
	Human Oxytocin-neurophysin 1 Recombinant Rabbit mAb [PSH12-28] - BSA and Azide free (Capture)	HA723452	Human	ELISA(Cap)
	Biotin Conjugated Human Oxytocin-neurophysin 1 Recombinant Rabbit mAb [PSH12-29] - Detector	HA723454B	Human	ELISA(Det),ELISA
	Human PAI1 Recombinant Rabbit mAb [PSH12-20] - BSA and Azide free (Capture)	HA723438	Human	ELISA(Cap)
	Human PAI1 Recombinant Rabbit mAb [PSH12-19] - BSA and Azide free (Detector)	HA723436	Human	ELISA(Det)
	Biotin Conjugated Human PAI1 Recombinant Rabbit mAb [PSH12-19] - Detector	HA723437B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human PAI1 Recombinant Rabbit mAb [PSH12-18] - Detector	HA723439B	Human	ELISA(Det),ELISA
	Human PAI1 Recombinant Rabbit mAb [PSH12-18] - BSA and Azide free (Capture/Detector)	HA723435	Human	ELISA(Cap),ELISA(Det)
	Biotin Conjugated Human Poliovirus Receptor/PVR Recombinant Rabbit mAb [PSH12-55] - Detector	HA723466B	Human	ELISA(Det),ELISA
	Human Poliovirus Receptor/PVR Recombinant Rabbit mAb [PSH12-54] - BSA and Azide free (Capture)	HA723464	Human	ELISA(Cap)
	Human Poliovirus Receptor/PVR Recombinant Rabbit mAb [PSH12-55] - BSA and Azide free (Detector)	HA723465	Human	ELISA(Det)
	Human VCAM-1 Recombinant Rabbit mAb [PSH12-46] - BSA and Azide free (Capture)	HA725064	Human	ELISA(Cap)
	Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - BSA and Azide free (Detector)	HA725065	Human	ELISA(Det)
	Biotin Conjugated Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - Detector	HA725066B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human VE Cadherin enzyme-linked immunoassay kit	EH0091	Human	ELISA
	HRP Conjugated Human Von Willebrand Factor Recombinant Rabbit mAb [PSH09-04] - Detector	HA723431H	Human	ELISA(Det)
Cell treatment, Relative expression	IDO Recombinant Rabbit mAb [PSH12-31]	HA723456	Human	WB,IF-Cell,IHC-P,FC,IP
Relative expression	KCNA5 / KV1-5 Recombinant Rabbit mAb [PSH12-30]	HA723455	Human, Mouse, Rat	WB
	MBP Nanoselector Agarose	015-101-002	Species independent	IP,ChIP,MS,Purification
	Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-45] - BSA and Azide free (Detector)	HA725041	Mouse	ELISA(Det)
	Biotin Conjugated Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-45] - Detector	HA725042B	Mouse	ELISA(Det),ELISA
	Mouse CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-44] - BSA and Azide free (Capture)	HA725040	Mouse	ELISA(Cap)
	Mouse CD68 Recombinant Rabbit mAb [PSH12-26] - BSA and Azide free (Capture)	HA723449	Mouse	ELISA(Cap)
	Mouse CD68 Recombinant Rabbit mAb [PSH12-27] - BSA and Azide free (Detector)	HA723450	Mouse	ELISA(Det)
	Mouse CD68 Recombinant Rabbit mAb [PSH12-24] - BSA and Azide free (Capture)	HA723446	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CD68 Recombinant Rabbit mAb [PSH12-25] - Detector	HA723448B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse CD68 Recombinant Rabbit mAb [PSH12-27] - Detector	HA723451B	Mouse	ELISA(Det),ELISA
	Mouse CD68 Recombinant Rabbit mAb [PSH12-25] - BSA and Azide free (Detector)	HA723447	Mouse	ELISA(Det)
	Mouse CRP Recombinant Rabbit mAb [PSH12-48] - BSA and Azide free (Capture)	HA725073	Mouse	ELISA(Cap)
	Mouse CRP Recombinant Rabbit mAb [PSH12-49] - BSA and Azide free (Detector)	HA725074	Mouse	ELISA(Det)
	Biotin Conjugated Mouse CRP Recombinant Rabbit mAb [PSH12-49] - Detector	HA725075B	Mouse	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse CXCL2 Recombinant Rabbit mAb [PSH12-61] - Detector	HA725072B	Mouse	ELISA(Det),ELISA
	Mouse CXCL2 Recombinant Rabbit mAb [PSH12-60] - BSA and Azide free (Capture)	HA725070	Mouse	ELISA(Cap)
	Mouse CXCL2 Recombinant Rabbit mAb [PSH12-61] - BSA and Azide free (Detector)	HA725071	Mouse	ELISA(Det)
	Mouse CXCL9 Recombinant Rabbit mAb [PSH12-62] - BSA and Azide free (Capture)	HA725022	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CXCL9 Recombinant Rabbit mAb [PSH12-63]	HA725024B	Mouse	ELISA(Det),ELISA
	Mouse CXCL9 Recombinant Rabbit mAb [PSH12-63] - BSA and Azide free (Detector)	HA725023	Mouse	ELISA(Det)
	Mouse Factor XII Recombinant Rabbit mAb [PSH12-22] - BSA and Azide free (Detector)	HA723441	Mouse	ELISA(Det)
	Mouse Factor XII Recombinant Rabbit mAb [PSH12-21] - BSA and Azide free (Capture)	HA723440	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse Factor XII Recombinant Rabbit mAb [PSH12-22] - Detector	HA723442B	Mouse	ELISA(Det),ELISA
	Mouse GM-CSF Recombinant Rabbit mAb [PSH12-41] - BSA and Azide free (Detector)	HA725020	Mouse	ELISA(Det)
	Biotin Conjugated Mouse GM-CSF Recombinant Rabbit mAb [PSH12-41] - Detector	HA725021B	Mouse	ELISA(Det),ELISA
	Mouse GM-CSF Recombinant Rabbit mAb [PSH12-40] - BSA and Azide free (Capture)	HA725019	Mouse	ELISA(Cap)
	Mouse IL-1 β enzyme-linked immunoassay kit	EM0029	Mouse	ELISA
	Mouse/Rat Factor XII Recombinant Rabbit mAb [PSH12-21] - BSA and Azide free (Capture)	HA723443	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat Factor XII Recombinant Rabbit mAb [PSH12-23] - Detector	HA723445B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat Factor XII Recombinant Rabbit mAb [PSH12-23] - BSA and Azide free (Detector)	HA723444	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse/Rat IGF1 Recombinant Rabbit mAb [PSH12-43] - Detector	HA725039B	Mouse, Rat	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse/Rat IGF1 Recombinant Rabbit mAb [PSH12-43] - BSA and Azide free (Detector)	HA725038	Mouse, Rat	ELISA(Det)
	Mouse/Rat IGF1 Recombinant Rabbit mAb [PSH12-42] - BSA and Azide free (Capture)	HA725037	Mouse, Rat	ELISA(Cap)
	Nrf3 Rabbit pAb	HA500570	Human, Mouse	WB
	Oxytocin-neurophysin 1 / OXT Recombinant Antibody [PSH08-24] - Rat IgG1 (Chimeric)	HA601426	Mouse, Rat, Human, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Tissue
	Oxytocin-neurophysin 1 / OXT Recombinant Antibody [PSH08-24] - Guinea pig IgG2 (Chimeric)	HA601425	Mouse, Rat, Human, Cynomolgus monkey, Pig	IHC-Fr,IF-Tissue
	PD-L1 Recombinant Rabbit mAb	IRS009	Human	mIHC
	PPAR Rabbit pAb	HA500569	Mouse, Rat	WB
	PPP2R2A (N-term) Rabbit pAb	ER1919-61	Human, Mouse, Rat	WB
	TIM 3 Recombinant Rabbit mAb	IRS077	Human	mIHC
	TMEM97 (N-term) Rabbit pAb	ER1919-62	Human, Mouse, Rat	WB
	TNRC6A Rabbit pAb	ER1919-60	Human, Mouse	WB
	猪伪狂犬病毒gB(PRV gB)抗体ELISA检测试剂盒	AE0021		ELISA
	猪繁殖与呼吸综合征病毒抗体检测试剂盒(OD 450nm)	AES001		ELISA
	非洲猪瘟病毒(ASFV) P54抗体检测试剂盒	AE0012		ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层