

CAR 信号转导研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	AKT1 Recombinant Rabbit mAb [PSH0-50] - BSA and Azide free	HA750627	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
	AKT1 Recombinant Rabbit mAb [PSH0-50]	HA721364	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
KD	AKT1 Recombinant Rabbit mAb [ST05-09]	ET1609-47	Human, Mouse, Rat	WB,IHC- P,IP,FC,IHC-Fr,IF- Cell
	AKT1 Mouse mAb [D9-9-C9]	EM40507	Human, Mouse, Rat	WB,IHC-P
KD	AKT1 Recombinant Rabbit mAb [ST05-09] - BSA and Azide free	HA750177	Human, Mouse, Rat	WB,IHC- P,IP,FC,IHC-Fr,IF- Cell
	AKT1/2/3 Rabbit pAb	ER1901-33	Human, Mouse, Rat	WB,IF-Cell,IHC- P,FC
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75- 09] - BSA and Azide free	HA750840	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	AKT1/2/3 Recombinant Rabbit mAb [ST48- 09] - BSA and Azide free	HA750179	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF- Tissue,IHC- P,IP,FC,IHC-Fr
	AKT1/2/3 Recombinant Rabbit mAb [ST48- 09]	ET1609-51	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF- Tissue,IHC- P,IP,FC,IHC-Fr
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75- 09]	HA721870	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	Akt1/3 Rabbit pAb	ER62638	Human, Mouse, Rat	WB,IHC-P
	Akt2 Mouse mAb [2A1]	RT1027	Human, Mouse, Rat	WB,IP,IF,IHC-P
	AKT2 Rabbit pAb	HA500091	Human, Mouse, Rat	WB,IHC-P
	AP-1/Jun D Rabbit pAb	ER62677	Human, Mouse, Rat	WB,IHC-P
	ATF2 Recombinant Rabbit mAb [SZ17-01]	ET1603-15	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P,IP
	Bcl10 Recombinant Rabbit mAb [SN74-04]	ET1611-79	Human	WB,IF-Cell,IF- Tissue,IHC-P,IP
	Bcl10 Recombinant Rabbit mAb [PSH13- 49] - BSA and Azide free	HA751468	Human, Mouse, Rat	WB,IF-Cell,FC,IP
	Bcl10 Recombinant Rabbit mAb [PSH13- 49]	HA723538	Human, Mouse, Rat	WB,IF-Cell,FC,IP
	BCMA Recombinant Rabbit mAb [PSH16- 80] - BSA and Azide free	HA751607	Human	FC

	BCMA Recombinant Rabbit mAb [PSH16-80]	HA723880	Human	FC
Relative expression	BCMA Recombinant Rabbit mAb [PSH17-94]	HA723959	Human	FC
Relative expression	BCMA Recombinant Rabbit Multiclonal Antibody [PSH18-02]	HA723967	Human	IF-Cell,FC
	CAMK2D Rabbit pAb	R1305-2	Human	WB,IF-Cell,IHC-P
	CAMK2G Mouse mAb [10C-G1]	EM1706-16	Human, Rat	WB,IHC-P,FC
	CaMKII Recombinant Rabbit mAb [SU03-57] - BSA and Azide free	HA750150	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr,FC
	CaMKII Recombinant Rabbit mAb [SU03-57]	ET1608-47	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr,FC
	CaMKII alpha Recombinant Rabbit mAb [JM11-07]	ET1703-54	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr,IF-Cell,FC,IP
	CaMKII alpha/beta/delta Rabbit pAb	ER62885	Human, Mouse, Rat	IHC-P,WB
	CaMKII alpha/delta Rabbit pAb	ER62886	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	CaMKII beta/gamma Rabbit pAb	ER62887	Human, Mouse, Rat	WB,IF-Cell
	CaMKII beta/gamma/delta Rabbit pAb	ER62888	Human, Mouse, Rat	WB,IHC-P
KD	CaMKII delta Recombinant Rabbit mAb [JA30-07]	ET1704-37	Human, Mouse, Rat, Zebrafish	WB,IHC-P
KD	CARD11 Recombinant Rabbit mAb [JE37-28]	HA722309	Human, Mouse, Rat	WB,IHC-P
Relative expression	CD11b Recombinant Rabbit mAb [PSH03-96]	HA722075	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	CD11b Recombinant Rabbit mAb [JE01-17]	HA722856	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue
Relative expression	CD11b Recombinant Rabbit mAb [PSH03-96] - BSA and Azide free	HA750921	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	CD11b Mouse mAb [4-6]	EM1701-42	Human	IHC-P,IF-Cell,FC,WB,IF-Tissue
	CD11b Recombinant Rabbit mAb	IRS014	Human	mIHC
Relative expression	CD11b Mouse mAb [A2-9]	M1211-6	Human	WB,IHC-P,FC,IF-Tissue
	CD11b Mouse mAb [A2-B5-D2]	EM1902-35	Human, Mouse	WB
Relative expression	CD11b Recombinant Rabbit mAb [JU93-81] - BSA and Azide free	HA750451	Human	WB,IP,IHC-P,IF-Cell,IF-Tissue
Relative expression	CD11b Recombinant Rabbit mAb [JU93-81]	ET1706-04	Human	WB,IP,IHC-P,IF-Cell,IF-Tissue,mIHC
	CD11b/ITAM Mouse mAb [4G2]	RT1107	Human, Mouse	WB,IP,IF,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD134 Mouse mAb [2H1]	EM1701-24	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Relative expression	CD134 Recombinant Rabbit mAb [PSH10-14] - BSA and Azide free	HA751343	Human	WB,IF-Cell,FC,IP
Relative expression	CD134 Recombinant Rabbit mAb [JE36-48]	HA722700	Human	WB
Relative expression	CD134 Recombinant Rabbit mAb [PSH10-14]	HA723183	Human	WB,IF-Cell,FC,IP
	CD137 Mouse mAb [2G1]	M1701-11	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	CD19 Mouse mAb [A3-A8]	EM40308	Human, Mouse	WB,IHC-P,FC
Relative expression	CD19 Recombinant Rabbit mAb [PSH13-44]	HA723533	Human	WB,IF-Cell,FC
	CD19 Rabbit pAb	R1107-7	Human, Mouse	WB
	FITC Conjugated CD19 Recombinant Mouse mAb [PSH14-97]	HA601518F	Human	FC
Relative expression	APC Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723761F	Human	FC
Relative expression	CD19 Recombinant Rabbit mAb [PSH13-44] - BSA and Azide free	HA751463	Human	WB,IF-Cell,FC
	CD19 Recombinant Rabbit mAb [JF099-9]	ET1702-74	Human, Mouse, Rat	WB,FC
	CD19 Recombinant Rabbit mAb	IRS018	Human	mIHC
Relative expression	CD19 Recombinant Rabbit mAb [JF100-06]	ET1702-93	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,mIHC
	PE Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723763F	Human	FC
Relative expression	CD19 Recombinant Rabbit mAb [JF100-06] - BSA and Azide free	HA750365	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	FITC Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723762F	Human	FC
	CD223 Mouse mAb [4H6A3]	EM1710-97	Human	WB,FC
	APC Conjugated CD25 / IL-2 Receptor alpha Recombinant Mouse mAb [PSH14-46]	HA601530F	Human	FC
	CD25 / IL-2 Receptor alpha Recombinant Mouse mAb [PSH14-46]	HA601486	Human, Mouse	FC,IHC-P,IF-Tissue
	CD269 Rabbit pAb	ER65683	Human	WB
	CD28 Rabbit pAb	ER2001-42	Human, Mouse, Rat	WB,IHC-P,FC
	CD3 Recombinant Mouse mAb [PD01-43]	HA601110	Human	WB,IHC-P,IF-Tissue
Relative expression	CD3 Recombinant Rabbit mAb [PSH19-70]	HA724094	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD3 Recombinant Rabbit mAb [JE80-02]	HA720082	Human	WB,IHC-P,IF-Tissue,IP,mIHC,IF-Cell
	CD3 Recombinant Rabbit mAb	IRS022	Human	mIHC
Relative expression	CD30 Recombinant Rabbit mAb [PSH04-12]	HA722093	Human	WB,IF-Cell,IHC-P,IF-Tissue
Relative expression	CD30 Recombinant Rabbit mAb [PSH04-10]	HA722091	Human	WB,IHC-P,FC,IF-Tissue
Relative expression	CD30 Recombinant Rabbit mAb [PSH04-11]	HA722092	Human	WB,IF-Cell,IHC-P
Relative expression	CD30 Recombinant Rabbit mAb [JE57-71]	HA721579	Human	WB,IHC-P
	CD30-L Rabbit pAb	ER62943	Human, Mouse	WB,IF-Cell
	CD319 Mouse mAb [5G5A4]	EM1710-79	Human	WB,IHC-P,FC
	CD38 Recombinant Rabbit mAb [JE29-35]	ET7110-53	Human	WB,IHC-P,FC,IF-Tissue
	CD38 Recombinant Antibody [PD01-49] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610369	Human	WB,IHC-P,FC
Relative expression	CD38 Recombinant Antibody [PD01-49] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610257		
	CD38 Recombinant Rabbit mAb [SN07-20]	ET1611-66	Human	WB,IHC-P,FC,IF-Tissue
	CD38 Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS120MS	Human	mIHC
	CD38 Recombinant Antibody [PD01-49] - Mouse IgG1 (Chimeric)	HA601554	Human	WB,IHC-P,FC
Relative expression	CD38 Recombinant Rabbit mAb [PD01-49]	HA721268	Human	WB,IHC-P,IF-Cell,FC,mIHC,IF-Tissue
	CD38 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS119RT	Human	mIHC
	CD38 Recombinant Rabbit mAb	IRS006	Human	mIHC
Relative expression	CD38 Recombinant Antibody [PD01-49] - Rat IgG1 (Chimeric)	HA601417		
	CD3E Rabbit pAb	HA500344	Human, Mouse, Rat	WB
	CD3 epsilon Recombinant Rabbit mAb [JE17-74]	ET7110-63	Human	WB,IHC-P,FC
	CD3 epsilon Recombinant Rabbit mAb [SY0239] - BSA and Azide free	HA750115	Human	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
	CD3 epsilon Recombinant Rabbit mAb [SY0239]	ET1607-29	Human	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
	CD3 epsilon Rabbit pAb	ER80501	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD3 zeta Recombinant Rabbit mAb [JE09-37]	HA722692	Human, Rat	WB,IHC-P,IF-Tissue
	CD3 zeta Rabbit pAb	0805-7	Human	WB,IF-Cell,FC
Cell treatment	c-Fos Recombinant Antibody [PSH07-51] - Mouse IgG1 (Chimeric)	HA601440	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Cell,WB
	c-Fos Recombinant Rabbit mAb	IRS131RB	Mouse	mIHC
Cell treatment	c-Fos Recombinant Antibody [PSH05-77] - Guinea pig IgG2 (Chimeric)	HA601374	Mouse, Rat, Human, Cynomolgus monkey, Pig	IHC-Fr,WB
	c-Fos Mouse mAb [A10-G2]	EM1710-29	Human	WB,IHC-P,FC
Cell treatment	c-Fos Recombinant Antibody [PSH05-77] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610220	Mouse, Rat, Human, Cynomolgus monkey, Pig	IHC-Fr,WB
	c-Fos Recombinant Rabbit mAb [PSH07-51] - BSA and Azide free	HA751172	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IF-Cell,IHC-P,ChIP,IF-Tissue,IP
Cell treatment	c-Fos Recombinant Antibody [PSH07-51] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610275	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Cell,WB
Cell treatment	c-Fos Recombinant Rabbit mAb [JJ0938]	ET1701-95	Human, Rat	WB,IF-Cell,IHC-P
Cell treatment	c-Fos Recombinant Rabbit mAb [PSH05-77] - BSA and Azide free	HA751071	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,IHC-Fr,IF-Tissue
Cell treatment	c-Fos Recombinant Antibody [PSH05-77] - Rat IgG1 (Chimeric)	HA601373	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Tissue,IF-Cell,WB
Cell treatment	c-Fos Recombinant Rabbit mAb [PSH05-77]	HA722666	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,IHC-Fr,mIHC,IF-Tissue
	c-Fos Recombinant Rabbit mAb [PSH07-51]	HA722839	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IF-Cell,IHC-P,ChIP,IF-Tissue,IP
Cell treatment	c-Fos Recombinant Rabbit mAb [JJ0938] - BSA and Azide free	HA750329	Human, Rat	WB,IF-Cell,IHC-P
Cell treatment	c-Fos Recombinant Antibody [PSH05-77] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610219	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Tissue,IF-Cell,WB
	c-Jun Recombinant Rabbit mAb [SY0298]	ET1608-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	c-Jun Recombinant Rabbit mAb [SY0298] - BSA and Azide free	HA750143	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	CL12A Rabbit pAb	ER61745	Human, Mouse	WB
	CLL1 Mouse mAb [4D5A8]	EM1711-46	Human	WB,IF-Cell,FC
	CREB Recombinant Rabbit mAb [SA04-04]	ET1601-15	Human, Mouse, Zebrafish, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CREB Recombinant Rabbit mAb [PSH15-48]	HA723754	Human, Mouse, Rat, Monkey	WB,IHC-P,IHC-Fr,IF-Cell,FC
	CREB Recombinant Rabbit mAb [PSH15-48] - BSA and Azide free	HA751563	Human, Mouse, Rat, Monkey	WB,IHC-P,IHC-Fr,IF-Cell,FC
	CREB Recombinant Rabbit mAb [SA04-04] - BSA and Azide free	HA750007	Human, Mouse, Zebrafish, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	CREB1 Rabbit pAb	R1309-14	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	CTLA-4 Recombinant Rabbit mAb [PS01-34] - BSA and Azide free	HA750581	Human, Mouse	IHC-P
	CTLA-4 Recombinant Rabbit mAb	IRS016	Human	mIHC
	CTLA-4 Recombinant Rabbit mAb [PS01-34]	HA721269	Human, Mouse	IHC-P,mIHC
	CTLA4 Rabbit pAb	ER1906-81	Human, Mouse, Rat	WB,FC
	DLL3 Rabbit pAb	ER1907-51	Human, Mouse	WB
Relative expression, KO	EGFR Recombinant Rabbit mAb [SP00-86] - BSA and Azide free	HA750078	Human, Mouse, Rat, Monkey	WB,IF-Cell,IP,FC
	EGFR Rabbit pAb	R1608-3	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
Relative expression, KO	EGFR Recombinant Rabbit mAb [SP00-86]	ET1604-44	Human, Mouse, Rat, Monkey	WB,IF-Cell,IP,FC
Relative expression, KD	EGFR Mouse mAb [A8E12]	HA601073	Human, Monkey	WB,IHC-P
	EGFR Rabbit pAb	ER1512-6	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
Relative expression, KD	EGFR Recombinant Rabbit mAb [PD00-75]	HA721220	Human, Monkey	WB,IHC-P,IF-Cell,FC
	EGFR Recombinant Rabbit mAb	IRS046	Human	mIHC
	EGFR Rabbit pAb	R1511-18	Human, Monkey	WB,IF-Cell,IHC-P
	EGFR Mouse mAb [B6-E5-D9]	EM1701-22	Human, Mouse, Monkey, Rat	WB,IHC-P,IF-Cell,FC
Relative expression, KO	EGFR Mouse mAb [7-F2-F8]	EM1901-67	Human, Monkey	WB,IHC-P,FC
Relative expression, KD	EGFR Recombinant Rabbit mAb [SZ40-19]	ET1603-37	Human, Mouse, Monkey, Rat	WB,IF-Cell,IP,FC,mIHC
Relative expression, KD	EGFR Recombinant Rabbit mAb [PD00-75] - BSA and Azide free	HA750550	Human, Monkey	WB,IHC-P,IF-Cell,FC
	EGFR Rabbit pAb	0407-21	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	ERK1 Mouse mAb [7-D6-E5]	M1505-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	ERK1 Recombinant Rabbit mAb [SP00-74]	ET1604-32	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	ERK1 Recombinant Rabbit mAb [SP05-09]	ET1604-16	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Erk1/2 Mouse mAb [2G1]	RT1484	Human, Mouse, Rat	WB,IP,IF
	ERK1/2 Rabbit pAb	ER131011	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03]	ET1601-29	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03] - BSA and Azide free	HA750019	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
KD	ERK2 Rabbit pAb	ER131218	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	ERK2 Recombinant Rabbit mAb [SZ25-01]	ET1603-23	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
KD	ERK2 Mouse mAb [B10-8]	EM1901-53	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Recombinant Mouse mAb [5-D2-R] - BSA and Azide free	HA610080	Human, Mouse, Rat	WB
KD	ERK2 Recombinant Mouse mAb [C3B12-R]	HA601201	Human, Mouse, Rat	WB
	ERK2 Mouse mAb [5-D2]	M1407-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	ERK2 Mouse mAb [A10D2]	EM1901-54	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Recombinant Mouse mAb [5-D2-R]	HA601202	Human, Mouse, Rat	WB
KD	ERK2 Recombinant Mouse mAb [C3B12-R] - BSA and Azide free	HA610079	Human, Mouse, Rat	WB
	ERK2 Mouse mAb [C3B12]	EM1901-55	Human, Mouse, Rat	WB,IHC-P,FC
Relative expression	Granzyme B Recombinant Rabbit mAb [PSH14-66] - BSA and Azide free	HA751527	Human, Mouse	WB,IHC-P
	Granzyme B Mouse mAb [4G9]	RT1259	Human, Mouse, Rat	WB,IP,IF,IHC-P
	Granzyme B Rabbit pAb	IRS032	Human	mIHC
	Granzyme B Rabbit pAb	HA500252	Human	IHC-P,mIHC
Relative expression	Granzyme B Recombinant Rabbit mAb [PSH14-66]	HA723664	Human, Mouse	WB,IHC-P
	GRB2 Recombinant Rabbit mAb [JE49-73]	HA721923	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
	GRB2 Mouse mAb [11-B11]	EM1701-41	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	GRB2 Rabbit pAb	ER31203	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	GTPase HRAS Recombinant Rabbit mAb [JE42-15]	HA721355	Human, Mouse, Rat	WB,IHC-P
	HER2 / ErbB2 Recombinant Rabbit mAb [PD00-53]	HA721178	Human	IHC-P,WB,IF-Cell,FC,mIHC
	HER2 / ErbB2 Rabbit pAb	R1511-33	Human	WB,IF-Cell,IHC-P,FC
	HER2 / ErbB2 Rabbit pAb	ER0106	Human	WB,IHC-P
	HER2 / ErbB2 Recombinant Rabbit mAb [PD00-99]	HA721210	Human	WB,IHC-P,IF-Tissue
	HER2 / ErbB2 Recombinant Rabbit mAb	IRS034	Human	mIHC
	Biotin Conjugated Human BCMA Recombinant Rabbit mAb [PSH13-19] - Detector	HA725084B	Human	ELISA(Det),ELISA
	Human BCMA Recombinant Rabbit mAb [PSH13-18] - BSA and Azide free (Capture)	HA725082	Human	ELISA(Cap)
	Human BCMA Recombinant Rabbit mAb [PSH13-19] - BSA and Azide free (Detector)	HA725083	Human	ELISA(Det)
	FITC Conjugated Human CD11b Recombinant Mouse mAb [A2-9-R]	HA600130F	Human	FC
	Human CD137 Recombinant Rabbit mAb [PSH08-61] - BSA and Azide free (Capture)	HA723009	Human	ELISA(Cap)
	Human CD137 Recombinant Rabbit mAb [PSH08-63] - BSA and Azide free (Capture)	HA723012	Human	ELISA(Cap)
	Biotin Conjugated Human CD137 Recombinant Rabbit mAb [PSH08-64] - Detector	HA723014B	Human	ELISA(Det),ELISA
	Human CD137 Recombinant Rabbit mAb [PSH08-64] - BSA and Azide free (Detector)	HA723013	Human	ELISA(Det)
	Human CD137 Recombinant Rabbit mAb [PSH08-62] - BSA and Azide free (Detector)	HA723010	Human	ELISA(Det)
	Biotin Conjugated Human CD137 Recombinant Rabbit mAb [PSH08-62] - Detector	HA723011B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human CD28 Recombinant Rabbit mAb [PSH12-51] - Detector	HA723459B	Human	ELISA(Det),ELISA
	Human CD28 Recombinant Rabbit mAb [PSH12-50] - BSA and Azide free (Capture)	HA723457	Human	ELISA(Cap)
	Human CD28 Recombinant Rabbit mAb [PSH12-51] - BSA and Azide free (Detector)	HA723458	Human	ELISA(Det)
	APC Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600118F	Human	FC
	PE Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600117F	Human	FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CD3 Recombinant Mouse mAb [PSH04-92]	HA601337	Human	FC
	FITC Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600116F	Human	FC
	Percp-Cy5.5 Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA601548F	Human	FC
	Human CD30 Recombinant Rabbit mAb [PSH04-82] - BSA and Azide free (Detector)	HA722188	Human	ELISA(Det)
	Biotin Conjugated Human CD30 Recombinant Rabbit mAb [PSH04-82]	HA722601B	Human	ELISA(Det)
	Human CD30 Recombinant Rabbit mAb [PSH04-81] - BSA and Azide free (Capture)	HA722187	Human	ELISA(Cap)
	Human CD30-L Recombinant Rabbit mAb [PSH18-25] - BSA and Azide free (Capture)	HA725339	Human	ELISA(Cap)
	Biotin Conjugated Human CD30-L Recombinant Rabbit mAb [PSH18-26] - Detector	HA725340B	Human	ELISA(Det),ELISA
	Human CD30-L Recombinant Rabbit mAb [PSH18-26] - BSA and Azide free (Detector)	HA725340	Human	ELISA(Det)
	FITC Conjugated Human CD38 Recombinant Rabbit mAb [SN07-20]	HA720196F	Human	FC
	Biotin Conjugated Human CTLA-4 Recombinant Rabbit mAb [PSH08-14]	HA722956B	Human	ELISA(Det)
	Human CTLA-4 Recombinant Rabbit mAb [PSH08-13] - BSA and Azide free (Capture)	HA722954	Human	ELISA(Cap)
	Human CTLA-4 Recombinant Rabbit mAb [PSH08-14] - BSA and Azide free (Detector)	HA722955	Human	ELISA(Det)
	Human Granzyme B Recombinant Rabbit mAb [PSH15-99] - BSA and Azide free (Detector)	HA723804	Human	ELISA(Det)
	Biotin Conjugated Human Granzyme B Recombinant Rabbit mAb [PSH15-99] - Detector	HA723804B	Human	ELISA(Det),ELISA
	Human Granzyme B Recombinant Rabbit mAb [PSH15-98] - BSA and Azide free (Capture)	HA723803	Human	ELISA(Cap)
	Human IL-13R alpha 2 Recombinant Rabbit mAb [PSH07-64] - BSA and Azide free (Detector)	HA722881	Human	ELISA(Det)
	Human IL-13R alpha 2 Recombinant Rabbit mAb [PSH07-63] - BSA and Azide free (Capture)	HA722880	Human	ELISA(Cap)
	Biotin Conjugated Human IL-13R alpha 2 Recombinant Rabbit mAb [PSH07-64]	HA722882B	Human	ELISA(Det)
	Biotin Conjugated Human IL-2 Mouse mAb [A8B7]	HA601361B	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-2 Recombinant Rabbit mAb [PSH01-42] - BSA and Azide free (Capture)	HA721680	Human	ELISA(Cap)
	Human IL-2 Mouse mAb [A8B7] - BSA and Azide free (Detector)	HA601006	Human	WB,ELISA(Det)
	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PS01-03] - BSA and Azide free (Detector)	HA721281	Human	ELISA(Det)
	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PSH05-65] - BSA and Azide free (Capture)	HA722385	Human	ELISA(Cap)
	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PSH05-66] - BSA and Azide free (Detector)	HA722386	Human	ELISA(Det)
	Biotin Conjugated Human IL-2 Receptor alpha Recombinant Rabbit mAb [PSH05-66]	HA722594B	Human	ELISA(Det)
	Biotin Conjugated Human IL-2 Receptor alpha Recombinant Rabbit mAb [PS01-03]	HA722650B	Human	ELISA(Det)
	Human IL-2 Receptor alpha Recombinant Rabbit mAb [PS00-52] - BSA and Azide free (Capture)	HA721280	Human	ELISA(Cap)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-46] - BSA and Azide free (Capture)	HA721684	Human	ELISA(Cap)
	Biotin Conjugated Human Interferon gamma Recombinant Rabbit mAb [PS01-76] - Detector	HA721576	Human	ELISA(Det),ELISA
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-47] - BSA and Azide free (Detector)	HA721685	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PS01-76] - BSA and Azide free (Detector)	HA721575	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-45] - BSA and Azide free (Capture)	HA721683	Human	ELISA(Cap)
	Biotin Conjugated Human Interferon gamma Recombinant Rabbit mAb [PSH01-47]	HA722645B	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PS01-75] - BSA and Azide free (Capture)	HA721574	Human	ELISA(Cap)
	Human LAG-3 Recombinant Rabbit mAb [PSH05-30] - BSA and Azide free (Detector)	HA722259	Human	ELISA(Det)
	Biotin Conjugated Human LAG-3 Recombinant Rabbit mAb [PSH05-30]	HA722596B	Human	ELISA(Det)
	Human LAG-3 Recombinant Rabbit mAb [PSH05-29] - BSA and Azide free (Capture)	HA722258	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Mesothelin Recombinant Rabbit mAb [PSH02-30] - BSA and Azide free (Capture)	HA723287	Human	ELISA(Cap)
	Biotin Conjugated Human Mesothelin Recombinant Rabbit mAb [PSH0-57] - Detector	HA723289B	Human	ELISA(Det),ELISA
	HRP Conjugated Human Mesothelin Recombinant Rabbit mAb [PSH0-57] - Detector	HA723289H	Human	ELISA(Det),ELISA
	Human Mesothelin Recombinant Rabbit mAb [PSH0-57] - BSA and Azide free (Detector)	HA723288	Human	ELISA(Det)
	Human OX40/TNFRSF4 Recombinant Rabbit mAb [PSH10-17] - BSA and Azide free (Capture)	HA723187	Human	ELISA(Cap)
	Biotin Conjugated Human OX40/TNFRSF4 Recombinant Rabbit mAb [PSH10-18] - Detector	HA723189B	Human	ELISA(Det),ELISA
	Human OX40/TNFRSF4 Recombinant Rabbit mAb [PSH10-18] - BSA and Azide free (Detector)	HA723188	Human	ELISA(Det)
	Human PD1 Recombinant Rabbit mAb [PSH13-38] - BSA and Azide free (Detector)	HA723522	Human	ELISA(Det)
	Human PD1 Recombinant Rabbit mAb [SJ01-91] - BSA and Azide free (Capture)	HA723521	Human	ELISA(Cap)
	Human PD1 Recombinant Rabbit mAb [PSH13-39] - BSA and Azide free (Detector)	HA723524	Human	ELISA(Det)
	Biotin Conjugated Human PD1 Recombinant Rabbit mAb [PSH13-39] - Detector	HA723525B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human PD1 Recombinant Rabbit mAb [PSH13-38] - Detector	HA723523B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721559	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free (Capture)	HA721265	Human	ELISA(Cap)
	Human TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free (Detector)	HA721266	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free	HA750593	Human	IHC-P
	HRP Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721294	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40]	HA721286	Human	IHC-P,mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human TNF alpha Recombinant Rabbit mAb	IRS040	Human	mIHC
	ICOS Recombinant Rabbit mAb [JE39-99]	HA722247	Mouse, Rat	WB
	ICOS Rabbit pAb	ER65686	Human, Mouse, Rat	WB,IHC-P
	IKK alpha Recombinant Rabbit mAb [SN06-99]	ET1611-15	Human, Mouse, Rat	WB,IHC-P,IP,FC
	IKK alpha Recombinant Rabbit mAb [SN06-99] - BSA and Azide free	HA750246	Human, Mouse, Rat	WB,IHC-P,IP,FC
	IKK alpha Rabbit pAb	ER30911	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	IKK alpha + IKK beta Recombinant Rabbit mAb [SN63-02]	ET1611-23	Human, Mouse, Rat	WB,IF-Cell,IP
KD	IKK alpha + IKK beta Recombinant Rabbit mAb [SN63-02] - BSA and Azide free	HA750248	Human, Mouse, Rat	WB,IF-Cell,IP
KD	IKK beta Recombinant Rabbit mAb [JE32-81]	HA721715	Human	WB,FC
	IKK beta Rabbit pAb	R1706-13	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	IL-13R alpha 2 Recombinant Rabbit mAb [PSH03-85] - BSA and Azide free	HA750910	Human	WB,IF-Cell,IHC-P,FC
Relative expression	IL-13R alpha 2 Recombinant Rabbit mAb [PSH03-85]	HA722051	Human	WB,IF-Cell,IHC-P,FC
	IL-2 Rabbit pAb	ER63748	Human, Rat	WB,IHC-P
	IL-2 Mouse mAb [A8B8]	HA601015	Human	WB,ELISA
	IL-2 Mouse mAb [A8B9]	HA601016	Human	WB,ELISA
	IL-2 receptor alpha Rabbit pAb	HA500076	Human, Rat	WB
	Interferon gamma Recombinant Rabbit mAb [JM10-10]	ET1703-17	Human	WB,FC
	Interferon gamma Recombinant Rabbit mAb [JM10-10] - BSA and Azide free	HA750373	Human	WB,FC
	JNK1 Rabbit pAb	R1309-1	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	JNK1 Recombinant Rabbit mAb [JE63-04]	HA721899	Human, Mouse, Rat	WB,IHC-P
	JNK1+JNK2+JNK3 Recombinant Rabbit mAb [SA43-06]	ET1601-28	Human, Mouse, Rat, Cow, Monkey, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP
	JNK1+JNK2+JNK3 Recombinant Rabbit mAb [SA43-06] - BSA and Azide free	HA750018	Human, Mouse, Rat, Cow, Monkey, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP
	JNK1+JNK3 Recombinant Rabbit mAb [JJ08-86]	ET1701-47	Human, Mouse, Rat	WB,IP,IHC-P
	JNK2 Rabbit pAb	HA500110	Human, Mouse, Rat	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	JNK2 Recombinant Rabbit mAb [SC05-85]	ET1610-11	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	JNK3 Recombinant Rabbit mAb [SD082-09]	ET1612-68	Human, Mouse, Rat	WB,IF-Tissue,IHC-P
	K-Ras Rabbit pAb	ER65706	Human, Mouse, Rat	WB,IHC-P
	KRAS Mouse mAb [A8E4]	HA601058	Human, Mouse, Rat	WB,IF-Cell,FC
	KRAS Recombinant Mouse mAb [A8E6] - BSA and Azide free	HA610008	Human, Mouse, Rat	WB,IHC-P
	KRAS Recombinant Mouse mAb [A8E6]	HA601060	Human, Mouse, Rat	WB,IHC-P
	KRAS Mouse mAb [A8E5]	HA601059	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-43] - BSA and Azide free	HA750621	Human, Mouse	WB,IHC-P
	LAG-3 Recombinant Rabbit mAb	IRS020	Human	mIHC
Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-43]	HA721345	Human, Mouse	WB,IHC-P
Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-44]	HA721346	Human, Mouse	WB,IHC-P,IF-Cell
Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-44] - BSA and Azide free	HA750622	Human, Mouse	WB,IHC-P,IF-Cell
Cell treatment	LAG-3 Recombinant Rabbit mAb [PS01-33]	HA721358	Human	IHC-P,IF-Cell,WB,mIHC
Relative expression	LAT Recombinant Rabbit mAb [PSH09-76] - BSA and Azide free	HA751318	Human	WB,IF-Cell,FC
Relative expression	LAT Recombinant Rabbit mAb [PSH09-76]	HA723141	Human	WB,IF-Cell,FC
	LAT Rabbit pAb	ER63843	Human, Mouse, Rat	WB
	Lck Rabbit pAb	ER63849	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Relative expression	Lck Recombinant Rabbit mAb [PSH09-25]	HA723083	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
Relative expression	Lck Recombinant Rabbit mAb [PSH09-25] - BSA and Azide free	HA751279	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	Lck Rabbit pAb	HA500280	Human, Rat	WB,FC
	MALT1 Recombinant Rabbit mAb [PSH02-32]	HA721689	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	MALT1 Recombinant Rabbit mAb [PSH02-32] - BSA and Azide free	HA750722	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KD	MALT1 Recombinant Rabbit mAb [PSH02-31] - BSA and Azide free	HA750563	Human, Mouse, Rat	WB,IHC-P
	MALT1 Rabbit pAb	ER63902	Human, Mouse	WB,IHC-P,IF-Cell
KD	MALT1 Recombinant Rabbit mAb [PSH02-31]	HA721244	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	MAP3K7 Mouse mAb [1G5]	RT1377	Human, Mouse, Rat	WB,IP,IF
KD	MEK1 Recombinant Rabbit mAb [SZ22-01]	ET1603-20	Human, Mouse, Rat, Cow, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK1/2 Recombinant Rabbit mAb [SR13-07] - BSA and Azide free	HA750042	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK1/2 Recombinant Rabbit mAb [SR13-07]	ET1602-3	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK2 Recombinant Rabbit mAb [SD207-08]	ET1612-6	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	MEK7 Recombinant Rabbit mAb [SD20-87]	ET1612-35	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO	MEKK2 Recombinant Rabbit mAb [SC68-02]	ET1610-70	Human, Mouse, Rat	WB,IHC-P,IP,FC
KD	MEKK3 Recombinant Rabbit mAb [SD0839]	ET1612-98	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Relative expression	Mesothelin Recombinant Rabbit mAb [PD00-68]	HA721175	Human	WB,IHC-P,IF-Tissue
Relative expression	Mesothelin Recombinant Rabbit mAb [PSH02-30]	HA721810	Human	WB,IF-Cell,IF-Tissue,IHC-P
Relative expression	Mesothelin Recombinant Rabbit mAb [PSH0-57]	HA721404	Human	WB,IHC-P,IF-Tissue
Relative expression	Mesothelin Recombinant Rabbit mAb [PSH0-57] - BSA and Azide free	HA750634	Human	WB,IHC-P,IF-Tissue
	Mesothelin Rabbit pAb	ER1803-76	Human	WB,IF-Cell,IHC-P
	Biotin Conjugated Mouse IL-2 Recombinant Rabbit mAb [PSH04-28]	HA722617B	Mouse	ELISA(Det)
	Mouse IL-2 Recombinant Rabbit mAb [PSH04-27] - BSA and Azide free (Capture)	HA722110	Mouse	ELISA(Cap)
	Mouse IL-2 Recombinant Rabbit mAb [PSH04-28] - BSA and Azide free (Detector)	HA722111	Mouse	ELISA(Det)
Relative expression	NFAT1 Recombinant Rabbit mAb [JA11-08]	ET1704-14	Human	WB,IF-Tissue,IHC-P
Relative expression	NFAT1 Recombinant Rabbit mAb [JA11-08] - BSA and Azide free	HA750405	Human	WB,IF-Tissue,IHC-P
	NFAT2 Recombinant Rabbit mAb [JA81-03]	ET1704-45	Human	WB,IHC-P,FC
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04]	ET1603-12	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	NF-kB p65 Rabbit pAb	HA500402	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27]	HA721307	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04] - BSA and Azide free	HA750058	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27] - BSA and Azide free	HA750601	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	NF-kB p65 Rabbit pAb	ER0815	Human, Mouse, Rat, Zebrafish	WB,IHC-P,FC
	N/H/K-Ras Rabbit pAb	ER64095	Human, Mouse, Rat	WB,IHC-P
KO	p38 alpha / MAPK14 Recombinant Rabbit mAb [JF55-07]	ET1702-65	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD, Cell treatment	p38 alpha / MAPK14 Recombinant Rabbit mAb [SR43-04] - BSA and Azide free	HA750040	Human, Mouse, Rat	WB,IF-Cell,FC
KD, Cell treatment	p38 alpha / MAPK14 Recombinant Rabbit mAb [SR43-04]	ET1602-26	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	p38 alpha / MAPK14 Rabbit pAb	R1308-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	p38 alpha / MAPK14 Mouse mAb [9-B9]	EM40307	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	p38 alpha/MAPK14 + p38 beta/MAPK11 + p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-60]	HA722019	Human, Mouse, Rat, Monkey	WB
	p38 alpha/MAPK14 + p38 beta/MAPK11 + p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-60] - BSA and Azide free	HA750894	Human, Mouse, Rat, Monkey	WB
	p38 MAPK Recombinant Rabbit Multiclonal Antibody [PSH03-72]	HA722036	Human, Mouse, Rat, Monkey	WB
	p38 MAPK Recombinant Rabbit Multiclonal Antibody [PSH03-72] - BSA and Azide free	HA750906	Human, Mouse, Rat, Monkey	WB
	PD1 Recombinant Rabbit mAb [SJ01-91] - BSA and Azide free	HA750103	Human	WB,IHC-P
Cell treatment	PD1 Mouse mAb [1F2]	EM1707-60	Human	WB,IF-Cell,IHC-P,FC
	PD1 Mouse mAb	IRS190MS	Human	mIHC
	PD1 Recombinant Rabbit mAb [SJ01-91]	ET1606-41	Human	WB,IHC-P
	PD1 Blocking Mouse mAb [3A11]	EM1701-17	Human	ELISA,FC
	PDPK1 Recombinant Rabbit mAb [SD204-9]	ET1612-27	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-AKT1 (S124) Recombinant Rabbit mAb [JJ08-46]	ET1701-36	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-AKT1 (S129) Recombinant Rabbit mAb [JE59-50]	HA721677	Human, Mouse, Rat	WB
Cell treatment	Phospho-AKT1 (T450) Recombinant Rabbit mAb [SD08-12]	ET1612-73	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell,FC
Cell treatment	Phospho-AKT2 (S474) Rabbit pAb	HA500116	Human, Mouse	WB,IHC-P
Cell treatment	Phospho-AKT2 (S474) Recombinant Rabbit mAb [PSH09-46] - BSA and Azide free	HA751293	Human, Mouse, Rat	WB,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-AKT2 (S474) Recombinant Rabbit mAb [PSH09-46]	HA723108	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44]	HA722129	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05] - BSA and Azide free	HA750132	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05]	ET1607-73	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44] - BSA and Azide free	HA750953	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10]	HA722951	Human	WB,IF-Cell,FC
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10] - BSA and Azide free	HA751205	Human	WB,IF-Cell,FC
Cell treatment	Phospho-ATF2 (T71) Recombinant Rabbit mAb [SC05-90]	ET1610-30	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IP,IHC-P
Cell treatment	Phospho-CaMKII (T286) Recombinant Rabbit mAb [JE51-81]	HA721794	Human, Mouse, Rat	WB
	Phospho-CaMKII (T286) Mouse mAb [2G1]	RT1092	Human, Mouse, Rat	WB,IP,IF
Cell treatment	Phospho-c-Jun (S63) Recombinant Rabbit mAb [SY0297]	ET1608-4	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-c-Jun (S63) Recombinant Rabbit mAb [SY0297] - BSA and Azide free	HA750146	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit mAb [PSH07-12]	HA722800	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit mAb [PSH07-12] - BSA and Azide free	HA751133	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit mAb [JE65-46]	HA722475	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-c-Jun (T91) Recombinant Rabbit mAb [JJ080-9]	ET1701-32	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-Creb (S133) Recombinant Rabbit mAb [JB25-40]	ET7107-93	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,IF-Tissue
Cell treatment	Phospho-Creb (S133) Recombinant Rabbit mAb [JB25-40] - BSA and Azide free	HA750465	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,IF-Tissue
	Phospho-EGFR (S695) Recombinant Rabbit mAb [SN07-36]	ET1611-50	Human	WB,IF-Cell,IF-Tissue
Cell treatment	Phospho-EGFR (T669) Recombinant Rabbit mAb [PSH10-03]	HA723172	Human, Mouse, Rat	WB,IP
Cell treatment	Phospho-EGFR (T669) Recombinant Rabbit mAb [PSH10-03] - BSA and Azide free	HA751333	Human, Mouse, Rat	WB,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-EGFR (Y1045) Recombinant Rabbit mAb [PSH03-28] - BSA and Azide free	HA750868	Human	WB,IF-Cell,FC
Cell treatment	Phospho-EGFR (Y1045) Recombinant Rabbit mAb [PSH03-28]	HA721970	Human	WB,IF-Cell,FC
Cell treatment	Phospho-EGFR (Y1068) Recombinant Rabbit mAb [SD2055] - BSA and Azide free	HA750290	Human	WB,IF-Cell,IF-Tissue,IP,FC
Cell treatment	Phospho-EGFR (Y1068) Recombinant Rabbit mAb [SD2055]	ET1612-30	Human	WB,IF-Cell,IF-Tissue,IP,FC
Cell treatment	Phospho-EGFR (Y1092) Recombinant Rabbit mAb [SJ0194]	ET1606-44	Human, Mouse	WB,IF-Tissue,IHC-P
	Phospho-EGFR (Y1148) Rabbit pAb	HA500580	Human	WB,IF-Cell
Cell treatment	Phospho-EGFR (Y1173) Recombinant Rabbit mAb [SC57-04]	ET1610-4	Human, Mouse	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-EGFR (Y845) Recombinant Rabbit mAb [JE46-77]	HA721488	Human	WB
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4] - BSA and Azide free	HA750061	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4]	ET1603-22	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01]	ET1610-13	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01] - BSA and Azide free	HA750199	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
	Phospho-HER2 / ErbB2 (T686) Mouse mAb [4G3]	RT1275	Human, Mouse, Rat	WB,IP
Cell treatment	Phospho-HER2 / ErbB2 (Y1139) Recombinant Rabbit mAb [JE60-90]	HA722201	Human	WB,IF-Cell,FC
Cell treatment	Phospho-HER2 / ErbB2 (Y1221 + Y1222) Recombinant Rabbit mAb [JE44-12]	HA721433	Human	WB
Cell treatment	Phospho-HER2 / ErbB2 (Y1248) Recombinant Rabbit mAb [PSH04-03]	HA722084	Human	WB,IF-Cell
Cell treatment	Phospho-HER2 / ErbB2 (Y1248) Recombinant Rabbit mAb [PSH04-03] - BSA and Azide free	HA750926	Human	WB,IF-Cell
Cell treatment	Phospho-JNK1/2/3 (T183 + T183 + T221) Recombinant Rabbit mAb [ST500] - BSA and Azide free	HA750173	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
Cell treatment	Phospho-JNK1/2/3 (T183 + T183 + T221) Recombinant Rabbit mAb [ST500]	ET1609-42	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment, Relative expression	Phospho-LYN (Y397)/LCK (Y394)/HCK (Y411)/BLK (Y389) Recombinant Rabbit mAb [PSH09-43] - BSA and Azide free	HA751290	Human, Mouse, Rat	WB,IF-Cell
Cell treatment, Relative expression	Phospho-LYN (Y397)/LCK (Y394)/HCK (Y411)/BLK (Y389) Recombinant Rabbit mAb [PSH09-43]	HA723105	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-MEK1/2 (S218 + S222) Recombinant Rabbit mAb [ST0490]	ET1609-50	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Phospho-MEK1 (S298) Recombinant Rabbit mAb [SD206-7]	ET1612-40	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-MEK1 (T292) Recombinant Rabbit mAb [SD2088]	ET1612-42	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-NF-kB p65 (S529) Recombinant Rabbit mAb [SP07-00]	ET1604-27	Human	WB,IHC-P,FC,IP
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58]	HA723223	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58] - BSA and Azide free	HA751361	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-p38 alpha (T180) Recombinant Rabbit mAb [PSH06-25]	HA722657	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180) Recombinant Rabbit mAb [PSH06-25] - BSA and Azide free	HA751062	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180 + Y182) Recombinant Rabbit mAb [PSH06-35]	HA722669	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180 + Y182) Recombinant Rabbit mAb [PSH06-35] - BSA and Azide free	HA751074	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-p38 (T180) Recombinant Rabbit mAb [JE77-37]	HA722150	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-p38 (T180) Rabbit pAb	ER2001-51	Human, Mouse, Rat	WB
Cell treatment	Phospho-p38 (T180 + Y182) Rabbit pAb	ER2001-52	Human, Mouse, Rat	WB
	Phospho-p38 (T180 + Y182) Rabbit pAb	ER1903-01	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
Cell treatment	Phospho-PLC-gamma-1 (S1248) Recombinant Rabbit mAb [PSH13-62] - BSA and Azide free	HA751481	Human	WB,IP
Cell treatment	Phospho-PLC-gamma-1 (S1248) Recombinant Rabbit mAb [PSH13-62]	HA723551	Human	WB,IP
Cell treatment	Phospho-PLC-gamma-1 (Y783) Recombinant Rabbit mAb [PSH06-29]	HA722661	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-PLC-gamma-1 (Y783) Recombinant Rabbit mAb [PSH06-29] - BSA and Azide free	HA751066	Human, Mouse, Rat	WB
Cell treatment, KD	Phospho-STAT3 (S727) Recombinant Rabbit mAb [SY24-09]	ET1607-39	Human, Rat, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment, KD	Phospho-STAT3 (S727) Recombinant Rabbit mAb [SY24-09] - BSA and Azide free	HA750121	Human, Rat, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-STAT3 (Y705) Recombinant Rabbit mAb [SZ43-01] - BSA and Azide free	HA750067	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-STAT3 (Y705) Recombinant Rabbit mAb [SZ43-01]	ET1603-40	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-Syk (Y352) Recombinant Rabbit mAb [PSH11-66] - BSA and Azide free	HA751422	Human	WB,IF-Cell,FC
Cell treatment	Phospho-Syk (Y352) Recombinant Rabbit mAb [PSH11-66]	HA723362	Human	WB,IF-Cell,FC
Cell treatment	Phospho-Syk (Y525 + Y526) Recombinant Rabbit mAb [PSH03-70]	HA722034	Human	WB,IF-Cell
Cell treatment	Phospho-Syk (Y525 + Y526) Recombinant Rabbit mAb [PSH03-70] - BSA and Azide free	HA750904	Human	WB,IF-Cell
Cell treatment	Phospho-TAK1 (S412) Recombinant Rabbit mAb [PSH14-64] - BSA and Azide free	HA751525	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-TAK1 (S412) Recombinant Rabbit mAb [PSH14-64]	HA723662	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-TAK1 (S439) Recombinant Rabbit mAb [JE44-93]	HA722242	Human, Mouse, Rat	WB
	PI3 Kinase p110 α Recombinant Rabbit mAb [SJ0186]	ET1606-36	Human, Mouse	WB,IF-Cell,IF-Tissue,IP,IHC-P
	PI 3-Kinase p110 γ Rabbit pAb	ER65965	Human, Mouse	WB,IHC,IF,ELISA
	PLCG1 Mouse mAb [D11-C2]	EM1707-02	Human, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	PLCG 2 Recombinant Rabbit mAb [JE39-12]	HA721477	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	PLCG2 Rabbit pAb	R1512-4	Human	WB,IF-Cell,IHC-P,FC
	PLC-gamma-1 Recombinant Rabbit mAb [SN06-10]	ET1611-25	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	Ras Rabbit pAb	ER40115	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Ras Recombinant Rabbit mAb [JF10-11]	ET1702-94	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IP,FC
	Ras Recombinant Rabbit mAb [JE44-13]	HA721883	Human, Mouse, Rat	WB,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Ras Recombinant Rabbit mAb [SA39-05]	ET1601-16	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC,IP
	Biotin Conjugated Rat IL-2 Recombinant Rabbit mAb [PSH05-72]	HA722590B	Rat	ELISA(Det)
	Rat IL-2 Recombinant Rabbit mAb [PS01-88] - BSA and Azide free (Capture)	HA721612	Rat	ELISA(Cap)
	Rat IL-2 Recombinant Rabbit mAb [PSH05-71] - BSA and Azide free (Capture)	HA722423	Rat	ELISA(Cap)
	Rat IL-2 Recombinant Rabbit mAb [PS01-89] - BSA and Azide free (Detector)	HA721613	Rat	ELISA(Det)
	Rat IL-2 Recombinant Rabbit mAb [PSH05-72] - BSA and Azide free (Detector)	HA722424	Rat	ELISA(Det)
	Biotin Conjugated Rat IL-2 Recombinant Rabbit mAb [PS01-89] - Detector	HA721614	Rat	ELISA(Det)
Relative expression	SLP76 Recombinant Rabbit mAb [PSH09-22] - BSA and Azide free	HA751276	Human	WB,IHC-P,IF-Cell,IP
Relative expression	SLP76 Recombinant Rabbit mAb [PSH09-22]	HA723080	Human	WB,IHC-P,IF-Cell,IP
	SLP76 Rabbit pAb	ER65421	Human, Mouse, Rat	WB,IHC-P
KD	STAT3 Recombinant Rabbit mAb [SY24-08] - BSA and Azide free	HA750120	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
KD	STAT3 Recombinant Rabbit mAb [SY24-08]	ET1607-38	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
KD	STAT3 Recombinant Rabbit mAb [SY34-01] - BSA and Azide free	HA750087	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	STAT3 Recombinant Rabbit mAb [SY34-01]	ET1605-45	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	STAT3 Rabbit pAb	RT1591	Human, Mouse, Rat	WB,IP,IF,IHC-P
KD	Syk Recombinant Rabbit mAb [PSH0-76] - BSA and Azide free	HA750655	Human	WB,IF-Cell
	Syk Rabbit pAb	HA500080	Human, Rat	WB,IF-Cell,IHC-P
KD	Syk Recombinant Rabbit mAb [PSH0-76]	HA721494	Human	WB,IF-Cell
	TAK1 Rabbit pAb	ER1802-90	Human, Mouse, Rat	IF-Cell,IHC-P,FC
KD	TAK1 Recombinant Rabbit mAb [JM73-19]	ET1705-14	Human, Mouse, Rat, Monkey	WB
KD	TAK1 Recombinant Rabbit mAb [JM73-19] - BSA and Azide free	HA750424	Human, Mouse, Rat, Monkey	WB
Relative expression	TIM 3 Recombinant Rabbit mAb [PSH10-00]	HA723169	Human	WB,IF-Cell,FC,IP
Relative expression	TIM 3 Recombinant Rabbit mAb [PSH07-75]	HA722891	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	TIM 3 Recombinant Antibody [PSH13-95] - Mouse IgG1 (Chimeric)	HA601442	Human	WB,IF-Cell,FC
	TIM 3 Recombinant Rabbit mAb	IRS077	Human	mIHC
Relative expression	TIM 3 Recombinant Rabbit mAb [PSH16-50] - BSA and Azide free	HA751599	Human, Mouse, Rat	WB,IHC-P,IP
	TIM 3 Mouse mAb [1E3]	EM1701-19	Human	WB,IHC-P,IF-Cell,FC
Relative expression	TIM 3 Recombinant Antibody [PSH13-95] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610277	Human	WB,IF-Cell,FC
	TIM 3 Mouse mAb [1E6]	EM1701-18	Human	WB,IF-Cell,IHC-P,FC
Relative expression	TIM 3 Recombinant Rabbit mAb [PSH16-50]	HA723854	Human, Mouse, Rat	WB,IHC-P,IP
	TIM 3 Mouse mAb [1E5]	EM1701-20	Human	WB,IF-Cell,IHC-P,FC
Relative expression	TIM 3 Recombinant Rabbit mAb [PSH07-75] - BSA and Azide free	HA751187	Human	WB
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free	HA750897	Human, Mouse	WB
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07]	HA722022	Human, Mouse	WB
	TNF alpha Rabbit pAb	ER1919-22	Human, Mouse	WB
	TNF alpha Rabbit pAb	ER65189	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	TRAF2 Recombinant Rabbit mAb [SD205-06] - BSA and Azide free	HA750297	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
	TRAF2 Recombinant Rabbit mAb [SD205-06]	ET1612-5	Human, Mouse, Rat	WB,IHC-P,IP
Relative expression	ZAP70 Recombinant Rabbit mAb [PSH11-52]	HA723344	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	ZAP70 Recombinant Mouse mAb [A1B5-R]	HA601218	Human, Mouse, Rat	WB,IHC-P
Relative expression	ZAP70 Mouse mAb [A1B7]	EM1902-28	Human, Mouse, Rat	WB,IF-Cell,FC
Relative expression	ZAP70 Mouse mAb [A1B5]	EM1901-27	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
Relative expression	ZAP70 Recombinant Rabbit mAb [JU08-39]	ET1706-42	Human	WB,IF-Cell,IHC-P,FC,IP
Relative expression	ZAP70 Recombinant Mouse mAb [A1B5-R] - BSA and Azide free	HA610096	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层