

结肠癌研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
KD	ALIX Recombinant Rabbit mAb [JM85-31] - BSA and Azide free	HA750439	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
KD	ALIX Recombinant Rabbit mAb [JM85-31]	ET1705-74	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
KD	Annexin V Recombinant Rabbit mAb [JF50-11]	ET1702-62	Human, Mouse, Rat	WB,FC
	Calretinin Recombinant Mouse mAb [PSH08-29] - BSA and Azide free	HA610215	Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Tissue,IP
	Calretinin Recombinant Rabbit mAb [JM12-93]	ET1705-19	Human, Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Tissue
	Calretinin Recombinant Mouse mAb [PSH08-29]	HA601369	Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Tissue,IP
	iFluor™ 488 Conjugated Calretinin Recombinant Rabbit mAb [JM12-93]	HA720153F	Human	IF-Cell,FC
	CD20 Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS112MS	Human	mIHC
	APC Conjugated CD20 Recombinant Rabbit mAb [PSH08-12]	HA723838F	Human	FC
Relative expression	CD20 Recombinant Rabbit mAb [PD00-02]	HA721138	Human, Mouse	IHC-P,mIHC,WB,IF-Cell,IP,IHC-Fr
	CD20 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS113RT	Human	mIHC
	CD20 Recombinant Rabbit mAb	IRS003	Human	mIHC
Relative expression	CD20 Recombinant Rabbit mAb [PSH08-12] - BSA and Azide free	HA751207	Human	WB,IF-Cell,IHC-P,FC,IP
Relative expression	CD20 Recombinant Antibody [PD00-02] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610360	Human, Mouse	WB,IHC-P,IHC-Fr
Relative expression	CD20 Recombinant Antibody [PD00-02] - Mouse IgG1 (Chimeric)	HA601544	Human, Mouse	WB,IHC-P,IHC-Fr
	PE Conjugated CD20 Recombinant Rabbit mAb [PSH08-12]	HA723839F	Human	FC
Relative expression	CD20 Recombinant Rabbit mAb [PSH08-12]	HA722953	Human	WB,IF-Cell,IHC-P,FC,IP
Relative expression	CD20 Recombinant Rabbit mAb [SY12-01]	ET1605-33	Human, Monkey	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD20 Rabbit pAb	IRS114RB	Mouse	mIHC
Relative expression, KD	CD9 Recombinant Rabbit mAb [PSH0-95] - BSA and Azide free	HA750675	Human	WB,IHC-P
	CD9 Rabbit pAb	ER80402	Human, Mouse	WB
	CD9 Recombinant Rabbit mAb [SA35-08] - BSA and Azide free	HA750027	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	CD9 Recombinant Rabbit mAb [SA35-08]	ET1601-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
Relative expression, KD	CD9 Recombinant Rabbit mAb [PSH0-95]	HA721533	Human	WB,IHC-P
	CDH17 Mouse mAb [D10-E1]	EM1706-54	Human	WB,IHC-P
KD, Relative expression	CDX2 Recombinant Rabbit mAb [SY09-02]	ET1605-4	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,CUT&Tag-seq
KD	COX2 Mouse mAb [A3F7]	EM1902-12	Human, Rat	WB,IHC-P,FC
	COX2 Mouse mAb	RT1159	Human, Mouse, Rat	WB,IP,IF,IHC-P,ELISA
	COX2 Mouse mAb [A3F8]	EM1902-11	Human, Rat	WB,IHC-P,FC
	COX2 Recombinant Rabbit mAb [6J8]	RE6064	Human, Mouse	WB,ICC
KD	COX2/Cyclooxygenase 2 Recombinant Rabbit mAb [SC56-06]	ET1610-23	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
KD	COX2/Cyclooxygenase 2 Recombinant Rabbit mAb [SC56-06] - BSA and Azide free	HA750207	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Cytokeratin 20 Mouse mAb [A3H7]	EM1901-97	Human	IHC-P,IF-Tissue
	Cytokeratin 20 Recombinant Mouse mAb [A3H7-R]	HA601077	Human	IHC-P,IF-Tissue
	Cytokeratin 20 Recombinant mAb	IRS135RB	Human	mIHC
	Cytokeratin 20 Rabbit pAb	ER1901-76	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Cytokeratin 20 Recombinant Rabbit mAb [SA35-03]	ET1601-8	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	iFluor™ 488 Conjugated Cytokeratin 20 Recombinant Rabbit mAb [SA35-03]	HA720189F	Human	IF-Cell,IF-Tissue
	Cytokeratin 20 Recombinant Rabbit mAb [JE20-39]	ET7110-54	Human, Mouse	WB,IHC-P
	Cytokeratin 20 Rabbit pAb	ER1901-97	Human, Mouse	WB,IHC-P,FC
	Cytokeratin 20 Recombinant Rabbit mAb [SA35-03] - BSA and Azide free	HA750026	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Cytokeratin 20 Mouse mAb [A3D12]	EM1901-96	Human	WB,FC,IHC-P
KD, Relative expression	Cytokeratin 7 Recombinant Rabbit mAb [ST50-05]	ET1609-62	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,mIHC
	Cytokeratin 7 Recombinant Rabbit mAb	IRS029	Human	mIHC
KD, Relative expression	Cytokeratin 7 Recombinant Rabbit mAb [ST50-05] - BSA and Azide free	HA750185	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO	Cytokeratin 7 Mouse mAb [3-G3-D8]	EM0702	Human	WB,IHC-P,IF-Cell,IF-Tissue
	Cytokeratin 7 Recombinant Mouse mAb [3-G3-D8-R]	HA601065	Human	WB,IHC-P,IF-Cell,FC,IF-Tissue
	Cytokeratin 7 Rabbit pAb	ER1803-89	Human, Mouse, Rat	WB,IHC-P
	iFluor™ 488 Conjugated Cytokeratin 7 Recombinant Rabbit mAb [ST50-05]	HA720131F	Human, Mouse	IF-Tissue,IF-Cell,FC
KD	Cytokeratin 7 Rabbit pAb	R1309-4	Human	WB,IF-Cell,IHC-P,FC
	iFluor™ 647 Conjugated Cytokeratin 7 Recombinant Rabbit mAb [ST50-05]	HA720144F	Human, Mouse	IF-Cell,IF-Tissue,FC
	iFluor™ 594 Conjugated Cytokeratin 7 Recombinant Rabbit mAb [ST50-05]	HA720120F	Human, Mouse	IF-Cell,IF-Tissue,FC
	Cytokeratin (Pan) Mouse mAb [7H8C4]	EM1712-42	Human	WB,IHC,IF-Cell,FC
Relative expression, KD	EGFR Recombinant Rabbit mAb [PD00-75]	HA721220	Human, Monkey	WB,IHC-P,IF-Cell,FC
Relative expression, KD	EGFR Recombinant Rabbit mAb [PD00-75] - BSA and Azide free	HA750550	Human, Monkey	WB,IHC-P,IF-Cell,FC
	EGFR Rabbit pAb	R1608-3	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	EGFR Mouse mAb [B6-E5-D9]	EM1701-22	Human, Mouse, Monkey, Rat	WB,IHC-P,IF-Cell,FC
Relative expression, KO	EGFR Recombinant Rabbit mAb [SP00-86] - BSA and Azide free	HA750078	Human, Mouse, Rat, Monkey	WB,IF-Cell,IP,FC
Relative expression, KD	EGFR Mouse mAb [A8E12]	HA601073	Human, Monkey	WB,IHC-P
Relative expression, KO	EGFR Mouse mAb [7-F2-F8]	EM1901-67	Human, Monkey	WB,IHC-P,FC
	EGFR Rabbit pAb	R1511-18	Human, Monkey	WB,IF-Cell,IHC-P
	EGFR Rabbit pAb	ER1512-6	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression, KO	EGFR Recombinant Rabbit mAb [SP00-86]	ET1604-44	Human, Mouse, Rat, Monkey	WB,IF-Cell,IP,FC
	EGFR Rabbit pAb	0407-21	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P
Relative expression, KD	EGFR Recombinant Rabbit mAb [SZ40-19]	ET1603-37	Human, Mouse, Monkey, Rat	WB,IF-Cell,IP,FC,mIHC
	EGFR Recombinant Rabbit mAb	IRS046	Human	mIHC
	EpCAM Recombinant Rabbit mAb [PS01-69]	HA721644	Human	IHC-P,IF-Cell,FC
	EpCAM Mouse mAb [7E11]	EM1712-43	Human	WB,IHC-P,ELISA
	EpCAM Recombinant Rabbit mAb	IRS254RB	Human	mIHC
	EpCAM Rabbit pAb	ER1802-80	Human, Mouse, Rat	WB,IHC-P
	EpCAM Recombinant Rabbit mAb [PS01-69] - BSA and Azide free	HA750703	Human	IHC-P,IF-Cell,FC
	iFluor™ 488 Conjugated EpCAM Recombinant Rabbit mAb [PS01-69]	HA720217F	Human	IF-Cell
	iFluor™ 647 Conjugated EpCAM Recombinant Rabbit mAb [PS01-69]	HA720185F	Human	FC,IF-Cell
	EpCAM Recombinant Rabbit mAb [SU03-32]	ET1608-26	Human, Mouse, Rat	WB,IHC-P
	EpCAM Mouse mAb [3-E1-E2]	EM1111	Human, Mouse	WB,IF-Cell
	iFluor™ 555 Conjugated EpCAM Recombinant Rabbit mAb [PS01-69]	HA720218F	Human	FC
	Flotillin 1 Recombinant Rabbit mAb [JB19-45]	ET7107-82	Human, Mouse, Rat	WB,IHC-P
	HRP Conjugated GAPDH Recombinant Rabbit mAb [SA30-01]	HA721131	Human, Mouse, Rat	WB,IHC
	GAPDH Recombinant Rabbit mAb [SA30-01] - Loading control	ET1601-4	Human, Mouse, Rat, Chicken, Drosophila melanogaster, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	HRP Conjugated GAPDH Recombinant Rabbit mAb [JF81-04]	ET1702-66	Human, Mouse, Rat, Monkey, Chicken, Fish	WB
	GAPDH Rabbit pAb	R1210-1	Human, Mouse, Rat, Oryza sativa	WB,IF-Cell,IHC-P,FC
	GAPDH Rabbit pAb	ER1901-65	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	GAPDH Mouse mAb [12D6]	M1310-2	Human, Mouse, Rat, Zebrafish, Escherichia coli	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	GAPDH Recombinant Rabbit mAb [PD00-07] - BSA and Azide free	HA750478	Human, Mouse, Rat, Escherichia coli, Zebrafish, Drosophila melanogaster	WB,IHC-P,IF-Cell,FC
	GAPDH Recombinant Rabbit mAb [SA30-01] - BSA and Azide free	HA750022	Human, Mouse, Rat, Chicken, Drosophila melanogaster, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	GAPDH Mouse mAb [12D7]	EM1901-57	Human, Mouse, Rat, Zebrafish	WB,IHC-P,FC
	GAPDH Mouse mAb [5-E10]	EM1101	Human, Mouse, Rat, Zebrafish, Rabbit	WB,IF-Cell,IHC-P
	GAPDH Recombinant Rabbit mAb [PD00-07]	HA721136	Human, Mouse, Rat, Escherichia coli, Zebrafish, Drosophila melanogaster	WB,IHC-P,IF-Cell,FC
	GAPDH Rabbit pAb	ER1706-83	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	GM130 Recombinant Rabbit mAb [JE42-53]	HA721282	Human, Mouse	WB,IHC-P,IF-Cell
	GM130 (cis-Golgi Marker) Rabbit pAb	R1608-7	Human, Mouse, Rat	WB,IHC-P
	GPA33 Recombinant Rabbit mAb [JE49-54]	ET7109-93	Human	WB,IHC-P,IP,IF-Tissue
	Hexokinase 1 Recombinant Rabbit mAb [ST47-05]	ET1609-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Hexokinase 1 Recombinant Rabbit mAb [ST47-05] - BSA and Azide free	HA750168	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Hexokinase II Recombinant Rabbit mAb [PSH07-97] - BSA and Azide free	HA751199	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC,IP
Relative expression	Hexokinase II Recombinant Rabbit mAb [PSH08-81] - BSA and Azide free	HA751246	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP
Relative expression	Hexokinase II Recombinant Rabbit mAb [PSH08-81]	HA723033	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP
	Hexokinase II Rabbit pAb	HA500186	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Hexokinase II Recombinant Rabbit mAb [PSH07-97]	HA722933	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC,IP
	Hexokinase II Mouse mAb [F1-D3]	EM1710-04	Human	WB,IHC-P,FC
KD	Hsp70 Recombinant Rabbit mAb [SA0379]	ET1601-11	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Hsp70 Rabbit pAb	ER50802	Human, Mouse, Rat, Cow	WB,IF-Cell,IHC-P,FC
	Hsp70 Recombinant Rabbit mAb [PD01-65]	HA721497	Human	WB,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HSP70 Mouse mAb [3G10]	EM60062	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	HSP70 (Acetyl K246) Rabbit pAb	ER60056	Human	WB
	Biotin Conjugated Human ICAM1 Recombinant Rabbit mAb [PSH10-16] - Detector	HA723186B	Human	ELISA(Det),ELISA
	Human ICAM1 Recombinant Rabbit mAb [PSH10-15] - BSA and Azide free (Capture)	HA723184	Human	ELISA(Cap)
	Human ICAM1 Recombinant Rabbit mAb [PSH10-16] - BSA and Azide free (Detector)	HA723185	Human	ELISA(Det)
	HRP Conjugated Human ICAM1 Recombinant Rabbit mAb [PSH10-16] - Detector	HA723474H	Human	ELISA(Det),ELISA
	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-93] - Detector	HA723497B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-92] - Detector	HA723495B	Human	ELISA(Det),ELISA
	Human VEGFA Recombinant Rabbit mAb [PSH12-93] - BSA and Azide free (Detector)	HA723496	Human	ELISA(Det)
	Human VEGFA Recombinant Rabbit mAb [PSH12-92] - BSA and Azide free (Detector)	HA723494	Human	ELISA(Det)
	Human VEGFA Recombinant Rabbit mAb [PSH12-91] - BSA and Azide free (Capture)	HA723493	Human	ELISA(Cap)
Relative expression	ICAM1 Recombinant Rabbit mAb [ST0487] - BSA and Azide free	HA750176	Human	WB,IF-Tissue,IHC-P
	ICAM1 Rabbit pAb	ER131207	Human	WB,IF-Cell,IHC-P,FC
Relative expression	ICAM1 Recombinant Rabbit mAb [ST0487]	ET1609-46	Human	WB,IF-Tissue,IHC-P
	ICAM1 Mouse mAb [J0-F7]	M1511-6	Human	WB,IF-Cell,FC
KD	ICAM1 Mouse mAb [F11-A4]	M1505-7	Human	WB,IF-Cell,FC
Relative expression, KD	ICAM1 Recombinant Rabbit mAb [PSH09-31]	HA723089	Human	WB,IHC-P,IF-Tissue
	IGFBP2 Recombinant Rabbit mAb [JU37-46]	ET7106-63	Human	WB,IHC-P
	Ki67 Rabbit pAb	ER1902-75	Human	WB,IF-Cell,FC
	Ki67 Recombinant Rabbit mAb [ST50-01]	ET1609-34	Human	WB,IF-Cell,IHC-P,FC
	Ki67 Recombinant Mouse mAb [PD00-10]	HA601053	Human, Mouse	IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Ki67 Recombinant Rabbit mAb	IRS057	Mouse	mIHC
	iFluor™ 647 Conjugated Ki67 Recombinant Rabbit mAb [SR00-02]	HA720163F	Human, Mouse, Rat	IF-Cell,IF-Tissue,FC,IHC-Fr
	Ki67 Recombinant Rabbit mAb [PSH0-02]	HA721254	Human	IHC-P,FC
	Ki67 Rabbit pAb	RT1349	Human	WB,IP,IF,IHC-P
	Ki67 Recombinant Rabbit mAb [PSH0-02] - BSA and Azide free	HA750570	Human	IHC-P,FC
	Ki67 Recombinant Rabbit mAb	IRS036	Human	mIHC
	Ki67 Rabbit pAb	ER1706-46	Human	WB,IF-Cell,IHC-P,FC
	Ki67 Recombinant Rabbit mAb [SR00-02]	HA721115	Human, Mouse, Rat, Cynomolgus monkey	IHC-P,IF-Tissue,mIHC,IF-Cell,FC
	iFluor™ 488 Conjugated Ki67 Recombinant Rabbit mAb [SR00-02]	HA720162F	Human	IF-Cell,IF-Tissue,FC
KD	Lactate Dehydrogenase Recombinant Rabbit mAb [SU39-06] - BSA and Azide free	HA750155	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	Lactate Dehydrogenase Recombinant Rabbit mAb [SU39-06]	ET1608-57	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	LDHA Mouse mAb [20-179]	M1007-7	Human, Mouse	IHC-P
	LDHA Recombinant Rabbit mAb [PSH13-58] - BSA and Azide free	HA751477	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	LDHA Recombinant Rabbit mAb [PSH13-58]	HA723547	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
KD	LDHA Rabbit pAb	ER00702	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	LDHA + LDHB Recombinant Rabbit mAb [PSH12-75]	HA723486	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	LDHA + LDHB Recombinant Rabbit mAb [PSH12-75] - BSA and Azide free	HA751451	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
KD	LDHB Rabbit pAb	0807-1	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
	LI Cadherin Rabbit pAb	ER1804-02	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Relative expression	LI Cadherin Recombinant Rabbit mAb [JE60-34]	HA721895	Human, Mouse	WB,IHC-P,IF-Tissue,IF-Cell,FC
Relative expression	LI Cadherin Recombinant Rabbit mAb [PD00-49]	HA721643	Human	WB,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse VEGF Recombinant Rabbit mAb [PSH02-73] - BSA and Azide free (Detector)	HA721851	Mouse	ELISA(Det)
	Biotin Conjugated Mouse VEGF Recombinant Rabbit mAb [PSH02-73]	HA722636B	Mouse	ELISA(Det)
	Mouse VEGF Recombinant Rabbit mAb [PSH02-72] - BSA and Azide free (Capture)	HA721850	Mouse	ELISA(Cap)
	MUC2 Recombinant Rabbit mAb [JA10-92]	ET1704-06	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	MUC2 Recombinant Rabbit mAb	IRS175RB	Human	mIHC
	MUC2 Recombinant Rabbit mAb [JA10-92] - BSA and Azide free	HA750401	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	MUC2 Recombinant Rabbit mAb	IRS235RB	Mouse	mIHC
Cell treatment	p53 Recombinant Rabbit mAb [PSH03-95]	HA722074	Mouse	WB,IF-Cell,IP
KD	p53 Recombinant Mouse mAb [6B2-R] - BSA and Azide free	HA610194	Human, Monkey	WB,IF-Cell
Relative expression, KD, Cell treatment	p53 Recombinant Rabbit mAb [SA39-07] - BSA and Azide free	HA750005	Human	WB,IF-Cell,ChIP
KO	p53 Mouse mAb [7-5]	EM20603	Human	WB,IF-Cell,IHC-P,FC
KD	p53 Rabbit pAb	R1308-7	Human	WB,IF-Cell
KD	p53 Recombinant Rabbit mAb [JB42-26]	ET7107-33	Human	WB,IHC-P,FC,IP,IF-Cell,IF-Tissue
Relative expression	p53 Mouse mAb [6B2]	EM1701-91	Human, Mouse	WB,IF-Cell
Cell treatment, Relative expression, KD	p53 Recombinant Mouse mAb [6B5-R] - BSA and Azide free	HA610187	Human, Mouse, Monkey	WB,IF-Cell
Cell treatment	p53 Recombinant Rabbit mAb [PSH03-95] - BSA and Azide free	HA750920	Mouse	WB,IF-Cell,IP
	p53 Recombinant Rabbit mAb [PD01-32] - BSA and Azide free	HA750545	Human	IHC-P
Cell treatment, KD, Relative expression	p53 Recombinant Mouse mAb [6B5-R]	HA601315	Human, Mouse, Monkey	WB,IF-Cell
Relative expression, KD, Cell treatment	p53 Recombinant Rabbit mAb [SA39-07]	ET1601-13	Human	WB,IF-Cell,ChIP
Relative expression, KD	p53 Recombinant Rabbit mAb [SY010-6]	ET1605-16	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	p53 Mouse mAb [6B5]	EM1701-92	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	p53 Recombinant Rabbit mAb [PD01-32]	HA721213	Human	IHC-P
	p53 Mouse mAb [C8-A11]	M1312-2	Human	WB,IF-Cell,IHC-P
KD	p53 Recombinant Mouse mAb [6B2-R]	HA601322	Human, Monkey	WB,IF-Cell
Cell treatment	p53 (acetyl K370) Recombinant Rabbit mAb [SR40-09]	ET1602-38	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC
Cell treatment	p53 (acetyl K370) Recombinant Rabbit mAb [SR40-09] - BSA and Azide free	HA750047	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [PSH06-88] - BSA and Azide free	HA751112	Human	WB,IF-Cell
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [PSH06-87] - BSA and Azide free	HA751111	Human, Mouse	WB,IF-Cell
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [PSH06-88]	HA722763	Human	WB,IF-Cell
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [PSH06-87]	HA722762	Human, Mouse	WB,IF-Cell
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [SD0801]	ET1612-71	Human, Mouse	WB,IF-Cell,IHC-P
KD	PDHA1 Recombinant Rabbit mAb [JF996-0]	ET1702-75	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	PFKP Rabbit pAb	HA500472	Human, Mouse, Rat	WB,IHC-P
	Phospho-EGFR (S695) Recombinant Rabbit mAb [SN07-36]	ET1611-50	Human	WB,IF-Cell,IF-Tissue
Cell treatment	Phospho-EGFR (T669) Recombinant Rabbit mAb [PSH10-03] - BSA and Azide free	HA751333	Human, Mouse, Rat	WB,IP
Cell treatment	Phospho-EGFR (T669) Recombinant Rabbit mAb [PSH10-03]	HA723172	Human, Mouse, Rat	WB,IP
Cell treatment	Phospho-EGFR (Y1045) Recombinant Rabbit mAb [PSH03-28] - BSA and Azide free	HA750868	Human	WB,IF-Cell,FC
Cell treatment	Phospho-EGFR (Y1045) Recombinant Rabbit mAb [PSH03-28]	HA721970	Human	WB,IF-Cell,FC
Cell treatment	Phospho-EGFR (Y1068) Recombinant Rabbit mAb [SD2055]	ET1612-30	Human	WB,IF-Cell,IF-Tissue,IP,FC
Cell treatment	Phospho-EGFR (Y1068) Recombinant Rabbit mAb [SD2055] - BSA and Azide free	HA750290	Human	WB,IF-Cell,IF-Tissue,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-EGFR (Y1092) Recombinant Rabbit mAb [SJ0194]	ET1606-44	Human, Mouse	WB,IF-Tissue,IHC-P
	Phospho-EGFR (Y1148) Rabbit pAb	HA500580	Human	WB,IF-Cell
Cell treatment	Phospho-EGFR (Y1173) Recombinant Rabbit mAb [SC57- 04]	ET1610-4	Human, Mouse	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-EGFR (Y845) Recombinant Rabbit mAb [JE46-77]	HA721488	Human	WB
Cell treatment	Phospho-p53 (S15) Recombinant Rabbit mAb [PSH01-98]	HA721756	Human	WB,IF-Cell
Cell treatment	Phospho-p53 (S15) Recombinant Rabbit mAb [PSH01-98] - BSA and Azide free	HA750772	Human	WB,IF-Cell
Cell treatment	Phospho-p53 (S20) Recombinant Rabbit mAb [PSH06-86]	HA722761	Human	WB
Cell treatment	Phospho-p53 (S20) Recombinant Rabbit mAb [PSH06-86] - BSA and Azide free	HA751110	Human	WB
Cell treatment	Phospho-p53 (S33) Recombinant Rabbit mAb [JE41-97]	HA722703	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-p53 (S376) Recombinant Rabbit mAb [ST0440]	ET1609-14	Human, Mouse	WB
Cell treatment	Phospho-p53 (S392) Recombinant Rabbit mAb [SI17-04]	ET1606-24	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-p53 (S392) Recombinant Rabbit mAb [SI17-04] - BSA and Azide free	HA750097	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-p53 (T55) Recombinant Rabbit mAb [ST0439]	ET1609-13	Human	WB,IHC-P
Cell treatment	Phospho-PDHA1 (S293) Recombinant Rabbit mAb [JE58-27]	HA721443	Human, Mouse, Rat	WB,ICC,IHC-P
Relative expression	PKM1 Recombinant Rabbit mAb [PSH04-37]	HA722120	Human, Mouse, Rat	WB,IF-Cell,IHC- P,FC
Relative expression	PKM1 Recombinant Rabbit mAb [PSH04-37] - BSA and Azide free	HA750946	Human, Mouse, Rat	WB,IF-Cell,IHC- P,FC
	PKM1/2 Recombinant Rabbit mAb [JE07-34]	HA722524	Human, Mouse, Rat	WB,IP
Relative expression	PKM2 Recombinant Rabbit mAb [PSH11-74] - BSA and Azide free	HA751429	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC- P,FC,IP
KO	PKM2 Rabbit pAb	ER1802-70	Human	WB,IF-Cell,IHC-P
KO	PKM2 Rabbit pAb	ER1901-90	Human, Mouse	WB,FC
Relative expression	PKM2 Recombinant Rabbit mAb [PSH11-74]	HA723370	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC- P,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	PKM2 Rabbit pAb	R1603-5	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	VEGF Recombinant Rabbit mAb [SP07-01]	ET1604-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	VEGF Rabbit pAb	ER30607	Human	WB,IF-Cell,IHC-P,FC
	VEGF Recombinant Rabbit mAb [SP07-01] - BSA and Azide free	HA750074	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	VEGFA Recombinant Rabbit mAb [PSH19-17]	HA724051	Human, Mouse, Rat, Monkey	WB
	VIL1 Mouse mAb [G11-5E]	EM1706-46	Human	WB,IF-Cell,IHC-P,FC
	Villin1 Recombinant Rabbit mAb [A2] - BSA and Azide free	HA750477	Human	WB,IHC-P
Relative expression	Villin1 Recombinant Rabbit mAb [JE39-08]	ET7110-55	Human	WB,IHC-P,IP,IF-Tissue
	Villin1 Recombinant Rabbit mAb	IRS213RB	Human	mIHC
Relative expression	Villin1 Recombinant Rabbit mAb [JU34-75]	ET7106-62	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,mIHC
	Villin1 Rabbit pAb	HA500509	Human, Mouse, Rat	WB,IHC-P
	Villin1 Recombinant Rabbit mAb [A2]	HA721134	Human	WB,IHC-P
	Villin1 Recombinant Rabbit mAb [C4] - BSA and Azide free	HA750476	Human	WB,IHC-P
	Villin1 Recombinant Rabbit mAb [C4]	HA721129	Human	WB,IHC-P
	Villin1 Rabbit pAb	ER1901-31	Human, Mouse	WB,FC
	Villin1 Recombinant Rabbit mAb	IRS243RB	Mouse	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层