

2025年4月—新上市产品清单

高级验证	产品名	货号	反应种属	验证应用
	iFluor™ 647 Conjugated Alpaca anti-Rabbit IgG Fc Chimeric Goat Fc Recombinant Antibody	HA1138	Rabbit	IF-Cell,IF-Tissue
Relative expression	AMACR Recombinant Rabbit mAb [PSH13-55]	HA723544	Human, Mouse, Rat	WB,IHC-P
Relative expression	AMACR Recombinant Rabbit mAb [PSH13-54]	HA723543	Human, Mouse	WB,IHC-P,IF-Tissue
Relative expression	AMH Recombinant Rabbit mAb [PSH13-48]	HA723537	Mouse, Rat, Human	WB,IHC-Fr,IHC-P
	ASPA Recombinant Antibody [PSH08-04] - Rat IgG1 (Chimeric)	HA601430	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	BATF3 Recombinant Rabbit mAb [PSH13-52]	HA723541	Human	WB,IP,ChIP
	Bcl10 Recombinant Rabbit mAb [PSH13-49]	HA723538	Human, Mouse, Rat	WB,IF-Cell,FC,IP
	Bid Recombinant Rabbit mAb [PSH14-03]	HA723593	Human, Mouse	WB,IP
Cell treatment	Caspase-7 Recombinant Rabbit mAb [PSH13-69]	HA723558	Human, Mouse, Rat	WB,IHC-P
Relative expression	CaV2.3 / CACNA1E Recombinant Rabbit mAb [PSH13-51]	HA723540	Mouse, Rat	IHC-Fr,IHC-P,WB
	CD127 Recombinant Mouse mAb [PSH12-76]	HA601429	Human	FC
Relative expression	CD19 Recombinant Rabbit mAb [PSH13-44]	HA723533	Human	WB,IF-Cell,FC
	CD4 Recombinant Rabbit Multiclonal Antibody [PSH13-24]	HA723511	Human, Mouse	WB,IF-Cell,IHC-P
Relative expression	CD45 Recombinant Rabbit mAb [PSH13-01]	HA723506	Mouse	WB,IHC-P,IHC-Fr,FC,IF-Cell,IF-Tissue
	CD45 Rabbit pAb	HA500575	Mouse	WB,IF-Cell,FC
Relative expression	CD45 Recombinant Rabbit mAb [PSH13-02]	HA723507	Mouse, Rat	WB,IHC-Fr,IF-Tissue
Relative expression	CD45 Recombinant Rabbit Multiclonal Antibody [PSH14-19]	HA723606	Human, Mouse	WB,IHC-P,IF-Cell
	CD45RO Recombinant Mouse mAb [PSH13-94]	HA601441	Human	IHC-P,FC
Cell treatment	c-Fos Recombinant Antibody [PSH07-51] - Mouse IgG1 (Chimeric)	HA601440	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Cell,WB

	CHCHD2 Recombinant Rabbit mAb [PSH13-87]	HA723580	Human, Mouse, Rat	WB,IF-Cell,FC,IP
	Claudin 18.2 Recombinant Mouse mAb [PSH18-20]	HA601432	Human	IHC-P
	COX5A Recombinant Rabbit mAb [PSH13-57]	HA723546	Human, Mouse, Rat	WB,IHC-P
	DAG1 Recombinant Rabbit mAb [PSH13-91]	HA723584	Human, Mouse, Rat	WB
	DCC Recombinant Rabbit mAb [PSH12-73]	HA723484	Human, Mouse, Rat	WB,IHC-Fr
	DCTN1 / p150-glued Recombinant Rabbit mAb [PSH13-45]	HA723534	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Cell,FC,IF-Tissue,IP
	Decorin Recombinant Rabbit mAb [PSH12-99]	HA723504	Human, Mouse, Rat	WB,IHC-P
	Dlx2 Recombinant Rabbit mAb [PSH13-46]	HA723535	Mouse, Rat	WB
	Dlx2 Recombinant Rabbit mAb [PSH13-47]	HA723536	Mouse, Rat	IHC-P
	DR3 / LARD Recombinant Rabbit mAb [PSH13-90]	HA723583	Human, Mouse, Rat	WB,IF-Cell
Relative expression	E-Cadherin Recombinant Rabbit mAb [PSH13-75]	HA723564	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Cell,FC,IP
	FLI1 Recombinant Rabbit mAb [PSH14-17]	HA723605	Human	WB,IHC-P
	GAD65 Recombinant Antibody [PSH06-80] - Mouse IgG1 (Chimeric)	HA601455	Mouse, Rat	IHC-Fr,IHC-P,WB
	Golgin 97 Recombinant Rabbit mAb [PSH13-73]	HA723562	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	HES6 Rabbit pAb	HA500577	Human, Rat	WB,IF-Cell
Cell treatment	Histone H3 (acetyl K27) Recombinant Rabbit mAb [PSH13-93]	HA723586	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Human Apolipoprotein A-I/APOA1, Tag Free	HA211090	Human	
	Human Apolipoprotein E, Tag Free	HA211078	Human	
	Human Asprosin enzyme-linked immunoassay kit	EH0032	Human	ELISA
	Human BCMA Recombinant Rabbit mAb [PSH13-18] - BSA and Azide free (Capture)	HA725082	Human	ELISA(Cap)
	Biotin Conjugated Human BCMA Recombinant Rabbit mAb [PSH13-19] - Detector	HA725084B	Human	ELISA(Det),ELISA
	Human BCMA Recombinant Rabbit mAb [PSH13-19] - BSA and Azide free (Detector)	HA725083	Human	ELISA(Det)
	Human beta-NGF enzyme-linked immunoassay kit	EH0065	Human	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human beta-NGF, Tag Free Protein	HA211073	Human	
	Human Cathepsin D, C-His Tag Protein	HA211077	Human	
	Biotin Conjugated Human CCL17 / TARC Recombinant Rabbit mAb [PSH13-34] - Detector	HA723515B	Human	ELISA(Det),ELISA
	Human CCL17 / TARC Recombinant Rabbit mAb [PSH13-35] - BSA and Azide free (Detector)	HA723516	Human	ELISA(Det)
	Human CCL17 / TARC Recombinant Rabbit mAb [PSH13-33] - BSA and Azide free (Capture)	HA723513	Human	ELISA(Cap)
	Human CCL17 / TARC Recombinant Rabbit mAb [PSH13-34] - BSA and Azide free (Detector)	HA723514	Human	ELISA(Det)
	Biotin Conjugated Human CCL17 / TARC Recombinant Rabbit mAb [PSH13-35] - Detector	HA723517B	Human	ELISA(Det),ELISA
	Human CCL18/PARC/MIP-4 Protein	HA211080	Human	
	Human CCL20/MIP-3 alpha enzyme-linked immunoassay kit	EH0028	Human	ELISA
	Biotin Conjugated Human CCL21 Recombinant Rabbit mAb [PSH13-37] - Detector	HA723520B	Human	ELISA(Det),ELISA
	Human CCL21 Recombinant Rabbit mAb [PSH13-36] - BSA and Azide free (Capture)	HA723518	Human	ELISA(Cap)
	Human CCL21 Recombinant Rabbit mAb [PSH13-37] - BSA and Azide free (Detector)	HA723519	Human	ELISA(Det)
	Human CCL2 enzyme-linked immunoassay kit	EH0010	Human	ELISA
	HRP Conjugated Human CCL5 Recombinant Rabbit mAb [PSH13-17] - Detector	HA725057H	Human	ELISA(Det),ELISA
	Biotin Conjugated Human CCL5 Recombinant Rabbit mAb [PSH13-17] - Detector	HA725057B	Human	ELISA(Det),ELISA
	Human CCL5 Recombinant Rabbit mAb [PSH13-16] - BSA and Azide free (Capture)	HA725055	Human	ELISA(Cap)
	Human CCL5 Recombinant Rabbit mAb [PSH13-17] - BSA and Azide free (Detector)	HA725056	Human	ELISA(Det)
	Human CCL5/RANTES enzyme-linked immunoassay kit	EH0027	Human	ELISA
	Human CD147 Recombinant Rabbit mAb [PSH13-07] - BSA and Azide free (Detector)	HA725086	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human CD147 Recombinant Rabbit mAb [PSH13-07] - Detector	HA725087B	Human	ELISA(Det),ELISA
	Human CD147 Recombinant Rabbit mAb [PSH13-06] - BSA and Azide free (Capture)	HA725085	Human	ELISA(Cap)
	Human CD155/PVR, Tag Free (ECD) Protein	HA211135	Human	
	Human CD163 Recombinant Rabbit mAb [PSH13-79] - BSA and Azide free (Detector)	HA723567	Human	ELISA(Det)
	Biotin Conjugated Human CD163 Recombinant Rabbit mAb [PSH13-79] - Detector	HA723568B	Human	ELISA(Det),ELISA
	Human CD163 Recombinant Rabbit mAb [PSH13-80] - BSA and Azide free (Capture)	HA723569	Human	ELISA(Cap)
	Human CD163 Recombinant Rabbit mAb [PSH13-78] - BSA and Azide free (Capture)	HA723566	Human	ELISA(Cap)
	Human CDCP1 enzyme-linked immunoassay kit	EH0083	Human	ELISA
	Human CKMB enzyme-linked immunoassay kit	EH0072	Human	ELISA
	Human Complement C3, C-His Tag Protein	HA211092	Human	
	Human CSF-1-R Recombinant Rabbit mAb [PSH12-66] - BSA and Azide free (Capture/Detector)	HA723474	Human	ELISA(Cap),ELISA(Det)
	HRP Conjugated Human CSF-1-R Recombinant Rabbit mAb [PSH12-66] - Detector	HA723478H	Human	ELISA(Det),ELISA
	Human CSF-1-R Recombinant Rabbit mAb [PSH12-67] - BSA and Azide free (Detector)	HA723475	Human	ELISA(Det)
	Biotin Conjugated Human CSF-1-R Recombinant Rabbit mAb [PSH12-67] - Detector	HA723476B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human CSF-1-R Recombinant Rabbit mAb [PSH12-66] - Detector	HA723478B	Human	ELISA(Det),ELISA
	Human CSF-1-R Recombinant Rabbit mAb [PSH12-68] - BSA and Azide free (Capture)	HA723477	Human	ELISA(Cap)
	Human CXCL10/IP-10 enzyme-linked immunoassay kit	EH0026	Human	ELISA
	Human CXCL11/I-TAC, Tag Free Protein	HA211110	Human	
	Human CXCL1/GRO alpha enzyme-linked immunoassay kit	EH0054	Human	ELISA
	Human CXCL6/GCP-2, Tag Free Protein	HA211111	Human	

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CXCL9/MIG enzyme-linked immunoassay kit	EH0015	Human	ELISA
	Human Cytokeratin 8, Tag Free Protein	HA211101	Human	
	Human EGF enzyme-linked immunoassay kit	EH0023	Human	ELISA
	Human EGFR, C-His Tag Protein	HA211104	Human	
	Human EPHA2, C-His Tag (ECD) Protein	HA211075	Human	
	Human Factor B, C-His Tag Protein	HA211099	Human	
	Human FGF-19 Recombinant Rabbit mAb [PSH13-29] - BSA and Azide free (Capture)	HA725100	Human	ELISA(Cap)
	Biotin Conjugated Human FGF-19 Recombinant Rabbit mAb [PSH13-30] - Detector	HA725102B	Human	ELISA(Det),ELISA
	Human FGF-19 Recombinant Rabbit mAb [PSH13-30] - BSA and Azide free (Detector)	HA725101	Human	ELISA(Det)
	Human FGF1, C-His Tag Protein	HA211079	Human	
	Human FGF basic /FGF2/bFGF enzyme-linked immunoassay kit	EH0024	Human	ELISA
	Human G-CSF Recombinant Rabbit mAb [PSH13-13] - BSA and Azide free (Detector)	HA725002	Human	ELISA(Det)
	Human G-CSF Recombinant Rabbit mAb [PSH13-12] - BSA and Azide free (Capture)	HA725001	Human	ELISA(Cap)
	Biotin Conjugated Human G-CSF Recombinant Rabbit mAb [PSH13-13] - Detector	HA725003B	Human	ELISA(Det),ELISA
	Human GM-CSF enzyme-linked immunoassay kit	EH0051	Human	ELISA
	Human Granzyme B enzyme-linked immunoassay kit	EH0052	Human	ELISA
	Human Granzyme B, Tag Free Protein	HA211076	Human	
	Human Growth Hormone Recombinant Rabbit mAb [PSH12-82] - BSA and Azide free (Detector)	HA725107	Human	ELISA(Det)
	Biotin Conjugated Human Growth Hormone Recombinant Rabbit mAb [PSH12-82] - Detector	HA725108B	Human	ELISA(Det),ELISA
	Human Growth Hormone Recombinant Rabbit mAb [PSH12-81] - BSA and Azide free (Capture)	HA725106	Human	ELISA(Cap)
	Human GTPase HRAS, C-His Tag Protein	HA211122	Human	
	Human Haptoglobin/HP, C-His Tag Protein	HA211100	Human	

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human HGF enzyme-linked immunoassay kit	EH0019	Human	ELISA
	Human I-309/CCL1, Tag Free Protein	HA211081	Human	
	Human ICAM-1 / CD54 enzyme-linked immunoassay kit	EH0064	Human	ELISA
	Human IFN-gamma enzyme-linked immunoassay kit	EH0039	Human	ELISA
	Human IFN-gamma, Tag Free Protein	HA211129	Human	
	Human IFN- β enzyme-linked immunoassay kit	EH0017	Human	ELISA
	Human IGFBP-1, Tag Free Protein	HA211115	Human	
	Human IGFBP2, C-His Tag Protein	HA211103	Human	
	HRP Conjugated Human IGFBP3 Recombinant Rabbit mAb [PSH12-17] - Detector	HA723565H	Human	ELISA(Det)
	Human IGFBP4, Tag Free Protein	HA211096	Human	
	Human IGFBP-6, Tag Free Protein	HA211086	Human	
	Human IGFBP7, C-His Tag Protein	HA211108	Human	
	Human IL-10 enzyme-linked immunoassay kit	EH0009	Human	ELISA
	Human IL-12p70 enzyme-linked immunoassay kit	EH0021	Human	ELISA
	Human IL-13 enzyme-linked immunoassay kit	EH0049	Human	ELISA
	Human IL-15 enzyme-linked immunoassay kit	EH0031	Human	ELISA
	Human IL-17A enzyme-linked immunoassay kit	EH0004	Human	ELISA
	Human IL-17E/IL-25 enzyme-linked immunoassay kit	EH0066	Human	ELISA
	Human IL-18 enzyme-linked immunoassay kit	EH0053	Human	ELISA
	Human IL-1 β enzyme-linked immunoassay kit	EH0006	Human	ELISA
	Human IL-21 enzyme-linked immunoassay kit	EH0018	Human	ELISA
	Human IL-25/IL-17E, Tag Free Protein	HA211098	Human	
	Human IL-2 enzyme-linked immunoassay kit	EH0013	Human	ELISA
	Human IL-31 enzyme-linked immunoassay kit	EH0097	Human	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-33 enzyme-linked immunoassay kit	EH0047	Human	ELISA
	Human IL-3 enzyme-linked immunoassay kit	EH0050	Human	ELISA
	Human IL-4 enzyme-linked immunoassay kit	EH0007	Human	ELISA
	Human IL-5 enzyme-linked immunoassay kit	EH0008	Human	ELISA
	Human IL-6 enzyme-linked immunoassay kit	EH0001	Human	ELISA
	Human IL-7 enzyme-linked immunoassay kit	EH0045	Human	ELISA
	Human IL-8/CXCL8 enzyme-linked immunoassay kit	EH0014	Human	ELISA
	Human Leptin enzyme-linked immunoassay kit	EH0056	Human	ELISA
	Human Leptin Receptor, C-His Tag (ECD) Protein	HA211148	Human	
	Biotin Conjugated Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-68] - Detector	HA725111B	Human	ELISA(Det),ELISA
	Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-67] - BSA and Azide free (Capture)	HA725109	Human	ELISA(Cap)
	Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-68] - BSA and Azide free (Detector)	HA725110	Human	ELISA(Det)
	Human Lipoprotein A Recombinant Rabbit mAb [PSH13-82] - BSA and Azide free (Detector)	HA723571	Human	ELISA(Det)
	Human Lipoprotein A Recombinant Rabbit mAb [PSH13-81] - BSA and Azide free (Capture)	HA723570	Human	ELISA(Cap)
	Biotin Conjugated Human Lipoprotein A Recombinant Rabbit mAb [PSH13-83] - Detector	HA723574B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human Lipoprotein A Recombinant Rabbit mAb [PSH13-82] - Detector	HA723572B	Human	ELISA(Det),ELISA
	Human Lipoprotein A Recombinant Rabbit mAb [PSH13-83] - BSA and Azide free (Detector)	HA723573	Human	ELISA(Det)
	Human MBL2 enzyme-linked immunoassay kit	EH0090	Human	ELISA
	Human MERTK Recombinant Rabbit mAb [PSH12-69] - BSA and Azide free (Capture)	HA723479	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HRP Conjugated Human MERTK Recombinant Rabbit mAb [PSH12-70] - Detector	HA723481H	Human	ELISA(Det),ELISA
	Human MERTK Recombinant Rabbit mAb [PSH12-70] - BSA and Azide free (Detector)	HA723480	Human	ELISA(Det)
	Biotin Conjugated Human MERTK Recombinant Rabbit mAb [PSH12-70] - Detector	HA723481B	Human	ELISA(Det),ELISA
	HRP Conjugated Human Mesothelin Recombinant Rabbit mAb [PSH0-57] - Detector	HA723289H	Human	ELISA(Det),ELISA
	Human MICA enzyme-linked immunoassay kit	EH0044	Human	ELISA
	Human MMP-9 enzyme-linked immunoassay kit	EH0020	Human	ELISA
	Human/Mouse BDNF enzyme-linked immunoassay kit	EH0069	Human, Mouse	ELISA
	Human/Mouse/Rat/Porcine/Canine TGF- β 1 enzyme-linked immunoassay kit	EH0012	Human, Mouse, Rat, Porcine, Canine	ELISA
	Human N-Cadherin, C-His Tag Protein	HA211124	Human	
	Human NFL, Tag Free Protein	HA211087	Human	
	Human Osteopontin (OPN) enzyme-linked immunoassay kit	EH0057	Human	ELISA
	Human Osteoprotegerin Recombinant Rabbit mAb [PSH13-32] - BSA and Azide free (Detector)	HA725068	Human	ELISA(Det)
	Biotin Conjugated Human Osteoprotegerin Recombinant Rabbit mAb [PSH13-32] - Detector	HA725069B	Human	ELISA(Det),ELISA
	Human Osteoprotegerin Recombinant Rabbit mAb [PSH13-31] - BSA and Azide free (Capture)	HA725067	Human	ELISA(Cap)
	Human Osteoprotegerin/TNFRSF11B enzyme-linked immunoassay kit	EH0070	Human	ELISA
	Human Oxytocin/OXT, Tag Free Protein	HA211134	Human	
	Human PCSK9 Recombinant Rabbit mAb [PSH12-85] - BSA and Azide free (Capture)	HA725124	Human	ELISA(Cap)
	Human PCSK9 Recombinant Rabbit mAb [PSH12-86] - BSA and Azide free (Detector)	HA725125	Human	ELISA(Det)
	Biotin Conjugated Human PCSK9 Recombinant Rabbit mAb [PSH12-86] - Detector	HA725126B	Human	ELISA(Det),ELISA
	Human PD1 Recombinant Rabbit mAb [PSH13-38] - BSA and Azide free (Detector)	HA723522	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human PD1 Recombinant Rabbit mAb [PSH13-39] - Detector	HA723525B	Human	ELISA(Det),ELISA
	Human PD1 Recombinant Rabbit mAb [SJ01-91] - BSA and Azide free (Capture)	HA723521	Human	ELISA(Cap)
	Human PD1 Recombinant Rabbit mAb [PSH13-39] - BSA and Azide free (Detector)	HA723524	Human	ELISA(Det)
	Biotin Conjugated Human PD1 Recombinant Rabbit mAb [PSH13-38] - Detector	HA723523B	Human	ELISA(Det),ELISA
	Human PDGF-AA Recombinant Rabbit mAb [PSH13-28] - BSA and Azide free (Detector)	HA725062	Human	ELISA(Det)
	Biotin Conjugated Human PDGF-AA Recombinant Rabbit mAb [PSH13-28] - Detector	HA725063B	Human	ELISA(Det),ELISA
	Human PDGF-AA Recombinant Rabbit mAb [PSH13-27] - BSA and Azide free (Capture)	HA725061	Human	ELISA(Cap)
	Human PDGF-AA enzyme-linked immunoassay kit	EH0058	Human	ELISA
	Human PD-L1/B7-H1 enzyme-linked immunoassay kit	EH0048	Human	ELISA
	HRP Conjugated Human Poliovirus Receptor/PVR Recombinant Rabbit mAb [PSH12-55] - Detector	HA723466H	Human	ELISA(Det),ELISA
	Human S100B, C-His Tag Protein	HA211097	Human	
	Human Serpin A1/Alpha-1-antitrypsin, C-His Tag Protein	HA211106	Human	
	Human SIGLEC9, Tag Free (ECD) Protein	HA211074	Human	
	Human SOD1/Cu-Zn SOD, C-His Tag Protein	HA211121	Human	
	Human TGF- α enzyme-linked immunoassay kit	EH0059	Human	ELISA
	Biotin Conjugated Human Thrombopoietin Recombinant Rabbit mAb [PSH13-97] - Detector	HA725120B	Human	ELISA(Det),ELISA
	Human Thrombopoietin Recombinant Rabbit mAb [PSH13-97] - BSA and Azide free (Detector)	HA725119	Human	ELISA(Det)
	Human Thrombopoietin Recombinant Rabbit mAb [PSH13-96] - BSA and Azide free (Capture)	HA725118	Human	ELISA(Cap)
	Human TIM 4 Recombinant Rabbit mAb [PSH14-00] - BSA and Azide free (Capture)	HA723589	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human TIM 4 Recombinant Rabbit mAb [PSH14-01] - BSA and Azide free (Detector)	HA723590	Human	ELISA(Det)
	Biotin Conjugated Human TIM 4 Recombinant Rabbit mAb [PSH14-01] - Detector	HA723591B	Human	ELISA(Det),ELISA
	Human TNF- α enzyme-linked immunoassay kit	EH0002	Human	ELISA
	Human tPA Recombinant Rabbit mAb [PSH13-10] - BSA and Azide free (Capture)	HA725115	Human	ELISA(Cap)
	Biotin Conjugated Human tPA Recombinant Rabbit mAb [PSH13-11] - Detector	HA725117B	Human	ELISA(Det),ELISA
	Human tPA Recombinant Rabbit mAb [PSH13-11] - BSA and Azide free (Detector)	HA725116	Human	ELISA(Det)
	Human Transthyretin/TTR, Tag Free Protein	HA211085	Human	
	Human TSLP Recombinant Rabbit mAb [PSH12-89] - BSA and Azide free (Capture)	HA723490	Human	ELISA(Cap)
	Human TSLP Recombinant Rabbit mAb [PSH12-87] - BSA and Azide free (Capture)	HA723487	Human	ELISA(Cap)
	Biotin Conjugated Human TSLP Recombinant Rabbit mAb [PSH12-88] - Detector	HA723489B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human TSLP Recombinant Rabbit mAb [PSH12-90] - Detector	HA723492B	Human	ELISA(Det),ELISA
	Human TSLP Recombinant Rabbit mAb [PSH12-90] - BSA and Azide free (Detector)	HA723491	Human	ELISA(Det)
	Human TSLP Recombinant Rabbit mAb [PSH12-88] - BSA and Azide free (Detector)	HA723488	Human	ELISA(Det)
	Human TSLP enzyme-linked immunoassay kit	EH0061	Human	ELISA
	Human UCH-L1/PGP9.5, Tag Free Protein	HA211112	Human	
	Human uPAR Mouse mAb [PSH13-76] - BSA and Azide free (Capture)	HA601437	Human	ELISA(Cap)
	Biotin Conjugated Human uPAR Mouse mAb [PSH13-77] - Detector	HA601439B	Human	ELISA(Det),ELISA
	Human uPAR Mouse mAb [PSH13-77] - BSA and Azide free (Detector)	HA601438	Human	ELISA(Det)
	Human uPAR enzyme-linked immunoassay kit	EH0077	Human	ELISA
	Human uPAR/PLAUR, Tag Free Protein	HA211102	Human	

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HRP Conjugated Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - Detector	HA725066H	Human	ELISA(Det),ELISA
	Human VCAM-1 / CD106 enzyme-linked immunoassay kit	EH0067	Human	ELISA
	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-93] - Detector	HA723497B	Human	ELISA(Det),ELISA
	Human VEGFA Recombinant Rabbit mAb [PSH12-92] - BSA and Azide free (Detector)	HA723494	Human	ELISA(Det)
	Human VEGFA Recombinant Rabbit mAb [PSH12-93] - BSA and Azide free (Detector)	HA723496	Human	ELISA(Det)
	Human VEGFA Recombinant Rabbit mAb [PSH12-91] - BSA and Azide free (Capture)	HA723493	Human	ELISA(Cap)
	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-92] - Detector	HA723495B	Human	ELISA(Det),ELISA
	Human VEGF enzyme-linked immunoassay kit	EH0075	Human	ELISA
	Human VEGFR1/FLT1 C-His Tag (ECD) Protein	HA211083	Human	
	Human VEGFR3/FLT4, C-His Tag Protein	HA211107	Human	
	Biotin Conjugated Human VEGF Receptor 2 Recombinant Rabbit mAb [PSH12-84] - Detector	HA725123B	Human	ELISA(Det),ELISA
	Human VEGF Receptor 2 Recombinant Rabbit mAb [PSH12-84] - BSA and Azide free (Detector)	HA725122	Human	ELISA(Det)
	Human VEGF Receptor 2 Recombinant Rabbit mAb [PSH12-83] - BSA and Azide free (Capture)	HA725121	Human	ELISA(Cap)
	Human YKL-40 Recombinant Rabbit mAb [PSH13-08] - BSA and Azide free (Capture)	HA725088	Human	ELISA(Cap)
	Human YKL-40 Recombinant Rabbit mAb [PSH13-09] - BSA and Azide free (Detector)	HA725089	Human	ELISA(Det)
	Biotin Conjugated Human YKL-40 Recombinant Rabbit mAb [PSH13-09] - Detector	HA725090B	Human	ELISA(Det),ELISA
	Integrin alpha 1 Recombinant Rabbit mAb [PSH14-04]	HA723594	Human	FC
	Integrin alpha 1 Recombinant Rabbit mAb [PSH14-05]	HA723595	Human, Mouse, Rat	WB,IHC-P
Relative expression	IP3 receptor Recombinant Mouse mAb [PSH14-18]	HA601456	Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	Kv4.2 / KCND2 Recombinant Mouse mAb [PSH13-21]	HA601434	Mouse, Rat	IHC-Fr,IHC-P,WB
	LDHA Recombinant Rabbit mAb [PSH13-58]	HA723547	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	LDHA + LDHB Recombinant Rabbit mAb [PSH12-75]	HA723486	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
Relative expression	Livin Recombinant Rabbit mAb [PSH12-72]	HA723483	Human	WB,IHC-P,IP
	LRG1 Rabbit pAb	HA500576	Human, Mouse	WB
	LRRK2 Chicken pAb	HA500574	Mouse, Rat	IHC-Fr
	Mannose Receptor(CD206) Recombinant Antibody [PSH07-76] - Guinea pig IgG2 (Chimeric)	HA601452	Mouse, Rat	IHC-Fr,WB
	Mannose Receptor(CD206) Recombinant Antibody [PSH07-76] - Mouse IgG1 (Chimeric)	HA601453	Mouse, Rat	IHC-Fr,IF-Tissue,IHC-P,WB,IF-Cell
	MCU Recombinant Rabbit mAb [PSH13-74]	HA723563	Human, Mouse, Rat, Monkey	WB,IHC-Fr,IHC-P,IF-Cell,IP
	MeCP2 Recombinant Rabbit mAb [PSH14-12]	HA723602	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Relative expression	MICA Recombinant Rabbit mAb [PSH13-60]	HA723549	Human	WB,FC,IP
	Mouse alpha-Synuclein, C-His Tag Protein	HA211105	Mouse	
	Mouse Angiotensinogen, C-His Tag Protein	HA211120	Mouse	
	Mouse CCL2/MCP-1 enzyme-linked immunoassay kit	EM0025	Mouse	ELISA
	Mouse CCL3/MIP-1 alpha enzyme-linked immunoassay kit	EM0031	Mouse	ELISA
	Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-15] - BSA and Azide free (Capture)	HA725046	Mouse	ELISA(Cap)
	Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-16] - BSA and Azide free (Detector)	HA725047	Mouse	ELISA(Det)
	Biotin Conjugated Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-16] - Detector	HA725048B	Mouse	ELISA(Det),ELISA
	Mouse CCL5/RANTES, C-hFc Tag Protein	HA211095	Mouse	
	Mouse CCL5 / RANTES enzyme-linked immunoassay kit	EM0032	Mouse	ELISA
	Mouse CD68/SR-D1, C-His Tag (ECD) Protein	HA211131	Mouse	

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse CD86 Recombinant Rabbit mAb [PSH13-64] - Detector	HA723554B	Mouse	ELISA(Det),ELISA
	Mouse CD86 Recombinant Rabbit mAb [PSH13-65] - BSA and Azide free (Capture)	HA723555	Mouse	ELISA(Cap)
	Mouse CD86 Recombinant Rabbit mAb [PSH13-64] - BSA and Azide free (Detector)	HA723553	Mouse	ELISA(Det)
	Mouse CD86 Recombinant Rabbit mAb [PSH13-66] - BSA and Azide free (Detector)	HA723556	Mouse	ELISA(Det)
	Mouse CD86 Recombinant Rabbit mAb [PSH13-63] - BSA and Azide free (Capture)	HA723552	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CD86 Recombinant Rabbit mAb [PSH13-66] - Detector	HA723557B	Mouse	ELISA(Det),ELISA
	Mouse Complement C3, C-His Tag Protein	HA211093	Mouse	
	Mouse Complement C5, C-His Tag Protein	HA211091	Mouse	
	Mouse CRP/C-Reactive Protein enzyme-linked immunoassay kit	EM0035	Mouse	ELISA
	Mouse CTGF/CCN2, Tag Free Protein	HA211094	Mouse	
	Mouse CXCL1/KC enzyme-linked immunoassay kit	EM0017	Mouse	ELISA
	Mouse CXCL2/MIP-2 enzyme-linked immunoassay kit	EM0018	Mouse	ELISA
	Mouse CXCL9/MIG enzyme-linked immunoassay kit	EM0015	Mouse	ELISA
	Mouse DKK1, C-His Tag Protein	HA211082	Mouse	
	HRP Conjugated Mouse Factor XII Recombinant Rabbit mAb [PSH12-22] - Detector	HA723509H	Mouse	ELISA(Det),ELISA
	Mouse GDF-15, Tag Free Protein	HA211089	Mouse	
	Mouse GLP1R, C-His Tag Protein	HA211109	Mouse	
	Mouse GM-CSF/CSF2 enzyme-linked immunoassay kit	EM0011	Mouse	ELISA
	Mouse IFN- γ enzyme-linked immunoassay kit	EM0009	Mouse	ELISA
	Mouse IL-10 enzyme-linked immunoassay kit	EM0005	Mouse	ELISA
	Mouse IL-12/IL-23 p40 enzyme-linked immunoassay kit	EM0012	Mouse	ELISA
	Mouse IL-13 enzyme-linked immunoassay kit	EM0007	Mouse	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IL-17A enzyme-linked immunoassay kit	EM0008	Mouse	ELISA
	Mouse IL-17E enzyme-linked immunoassay kit	EM0028	Mouse	ELISA
	Mouse IL-18 enzyme-linked immunoassay kit	EM0034	Mouse	ELISA
	Biotin Conjugated Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-41] - Detector	HA723528B	Mouse	ELISA(Det),ELISA
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-41] - BSA and Azide free (Detector)	HA723527	Mouse	ELISA(Det)
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-42] - BSA and Azide free (Capture)	HA723529	Mouse	ELISA(Cap)
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-40] - BSA and Azide free (Capture)	HA723526	Mouse	ELISA(Cap)
	Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-43] - BSA and Azide free (Detector)	HA723530	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-1 alpha Recombinant Rabbit mAb [PSH13-43] - Detector	HA723531B	Mouse	ELISA(Det),ELISA
	Mouse IL-21 enzyme-linked immunoassay kit	EM0020	Mouse	ELISA
	Mouse IL-22 Recombinant Rabbit mAb [PSH13-05] - BSA and Azide free (Detector)	HA725032	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-22 Recombinant Rabbit mAb [PSH13-05] - Detector	HA725033B	Mouse	ELISA(Det),ELISA
	Mouse IL-22 Recombinant Rabbit mAb [PSH13-04] - BSA and Azide free (Capture)	HA725031	Mouse	ELISA(Cap)
	Mouse IL-22 enzyme-linked immunoassay kit	EM0021	Mouse	ELISA
	Mouse IL-25/IL-17E, Tag Free Protein	HA211130	Mouse	
	Mouse IL-2 enzyme-linked immunoassay kit	EM0002	Mouse	ELISA
	Mouse IL-33 enzyme-linked immunoassay kit	EM0033	Mouse	ELISA
	Mouse IL-4 enzyme-linked immunoassay kit	EM0003	Mouse	ELISA
	Mouse IL-5 Recombinant Rabbit mAb [PSH12-77] - BSA and Azide free (Capture)	HA725034	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-5 Recombinant Rabbit mAb [PSH12-78] - Detector	HA725036B	Mouse	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IL-5 Recombinant Rabbit mAb [PSH12-78] - BSA and Azide free (Detector)	HA725035	Mouse	ELISA(Det)
	Mouse IL-5 enzyme-linked immunoassay kit	EM0023	Mouse	ELISA
	Mouse IL-6 enzyme-linked immunoassay kit	EM0004	Mouse	ELISA
	Mouse IL-9 enzyme-linked immunoassay kit	EM0022	Mouse	ELISA
	Mouse Leptin enzyme-linked immunoassay kit	EM0027	Mouse	ELISA
	Mouse PDGF-AA, Tag Free Protein	HA211136	Mouse	
	Mouse Pro-Collagen I alpha 1/COL1A1, Tag Free Protein	HA211113	Mouse	
	Biotin Conjugated Mouse/Rat CXCL1/GRO alpha Recombinant Rabbit mAb [PSH13-15] - Detector	HA725030B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat CXCL1/GRO alpha Recombinant Rabbit mAb [PSH13-15] - BSA and Azide free (Detector)	HA725029	Mouse, Rat	ELISA(Det)
	Mouse/Rat CXCL1/GRO alpha Recombinant Rabbit mAb [PSH13-14] - BSA and Azide free (Capture)	HA725028	Mouse, Rat	ELISA(Cap)
	Mouse/Rat IGF-I/IGF-1 enzyme-linked immunoassay kit	EM0024	Mouse, Rat	ELISA
	Mouse Serpin C1/Antithrombin-III, C-His Tag Protein	HA211123	Mouse	
	Mouse TNF- α enzyme-linked immunoassay kit	EM0010	Mouse	ELISA
	Mouse Transthyretin/TTR, C-His Tag Protein	HA211118	Mouse	
	Mouse VEGF enzyme-linked immunoassay kit	EM0013	Mouse	ELISA
Relative expression	MSK1 Recombinant Rabbit mAb [PSH13-72]	HA723561	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
Relative expression	MTAP Recombinant Mouse mAb [PSH13-98]	HA601445	Human, Mouse, Rat	WB,IF-Cell
Relative expression	Myelin PLP Chicken pAb	HA500573	Mouse, Rat	IHC-Fr,WB
	Myt1 / MTF1 Recombinant Rabbit mAb [PSH12-97]	HA723502	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	NDNF Rabbit pAb	HA500571	Mouse, Rat	WB
	NDUFA9 Recombinant Rabbit mAb [PSH14-07]	HA723597	Human, Mouse, Rat, Monkey	WB,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	NDUFA9 Recombinant Rabbit mAb [PSH14-08]	HA723598	Human, Mouse, Rat, Monkey	WB,IP
Relative expression	NF-L Recombinant Antibody [PS02-10] - Mouse IgG1 (Chimeric)	HA601448	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	NF-M Recombinant Antibody [JM11-20] - Mouse IgG1 (Chimeric)	HA601443	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	NF-M Recombinant Antibody [JM11-20] - Guinea pig IgG2 (Chimeric)	HA601444	Mouse, Rat	IHC-Fr,WB
	non-muscle Myosin IIA Recombinant Rabbit mAb [PSH14-10]	HA723600	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	Osteopontin Recombinant Rabbit mAb [PSH14-11]	HA723601	Mouse, Rat	IHC-P,IHC-Fr,WB,IF-Cell
	p70 S6 Kinase Recombinant Rabbit mAb [PSH13-50]	HA723539	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	p75 NGF Receptor Recombinant Rat mAb [PSH13-26]	HA601436	Mouse, Rat	IHC-Fr,WB
	p75 NGF Receptor Recombinant Antibody [PSH13-26] - Rabbit IgG (Chimeric)	HA723532	Mouse	IHC-Fr
	PD-L1 Recombinant Rabbit Multiclonal Antibody [PSH13-23]	HA723510	Human, Mouse	WB
Relative expression	PGP9.5 Recombinant Rabbit mAb [PSH13-92]	HA723585	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue
KD	PHD2 Recombinant Rabbit mAb [PSH12-74]	HA723485	Human	WB
Cell treatment	Phospho-AMPK alpha 1 (T183) +AMPK alpha 2 (T172) Recombinant Rabbit mAb [PSH14-13]	HA723603	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-BCKDHA (S293) Recombinant Rabbit mAb [PSH14-02]	HA723592	Human, Mouse, Rat, Monkey	WB
	Phospho-MYL9 (T18 + S19) Recombinant Rabbit mAb [PSH12-71]	HA723482	Mouse, Rat	WB
Cell treatment	Phospho-PDGFR beta (Y751) Recombinant Rabbit mAb [PSH13-61]	HA723550	Mouse, Rat, Human	WB
Cell treatment	Phospho-PLC-gamma-1 (S1248) Recombinant Rabbit mAb [PSH13-62]	HA723551	Human	WB,IP
	PRAK / MK5 Recombinant Rabbit mAb [PSH14-06]	HA723596	Human, Mouse, Rat, Monkey	WB,IP
	PSMA7 / HSPC Recombinant Rabbit mAb [PSH14-09]	HA723599	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	QKI-5 Recombinant Antibody [PSH10-22] - Rabbit IgG (Chimeric)	HA723587	Mouse, Rat	IHC-Fr,IHC-P,WB
	QKI-6 Recombinant Antibody [PSH10-23] - Rabbit IgG (Chimeric)	HA723588	Mouse, Rat	IHC-Fr,IHC-P,WB
	Rat CCL2/MCP-1, Tag Free Protein	HA211114	Rat	

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Rat CRP Recombinant Rabbit mAb [PSH12-79] - BSA and Azide free (Capture)	HA725049	Rat	ELISA(Cap)
	HRP Conjugated Rat CRP Recombinant Rabbit mAb [PSH12-80] - Detector	HA725051H	Rat	ELISA(Det),ELISA
	Rat CRP Recombinant Rabbit mAb [PSH12-80] - BSA and Azide free (Detector)	HA725050	Rat	ELISA(Det)
	Biotin Conjugated Rat CRP Recombinant Rabbit mAb [PSH12-80] - Detector	HA725051B	Rat	ELISA(Det),ELISA
	Rat CRP/C-Reactive Protein enzyme-linked immunoassay kit	ER0009	Rat	ELISA
	Rat IFN- γ enzyme-linked immunoassay kit	ER0006	Rat	ELISA
	Rat IL-10 enzyme-linked immunoassay kit	ER0005	Rat	ELISA
	Rat IL-1 β enzyme-linked immunoassay kit	ER0001	Rat	ELISA
	Rat IL-2 enzyme-linked immunoassay kit	ER0002	Rat	ELISA
	Rat IL-4 enzyme-linked immunoassay kit	ER0003	Rat	ELISA
	Rat IL-6 enzyme-linked immunoassay kit	ER0004	Rat	ELISA
	Rat TNF- α enzyme-linked immunoassay kit	ER0007	Rat	ELISA
	Retinoid X Receptor beta / RXRB Recombinant Rabbit mAb [PSH13-53]	HA723542	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	ROR1 Recombinant Mouse mAb [PSH13-99]	HA601447	Human	WB,FC
	Slc12a5 / KCC2 Recombinant Mouse mAb [PSH13-20]	HA601433	Mouse, Rat	IHC-Fr,WB
Relative expression	SOX9 Recombinant Rabbit mAb [PSH13-59]	HA723548	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell,IP,FC,ChIP
	TBC1D1 Recombinant Rabbit mAb [PSH13-71]	HA723560	Human, Mouse, Rat	WB,IHC-P,IP
Relative expression	TCF7 Recombinant Rabbit mAb [PSH13-89]	HA723582	Human, Mouse, Rat	WB,IHC-P,IP,mIHC
	TGF beta Receptor I Recombinant Rabbit mAb [PSH14-14]	HA723604	Human, Mouse, Rat	WB,IP
Relative expression	TIM 3 Recombinant Antibody [PSH13-95] - Mouse IgG1 (Chimeric)	HA601442	Human	WB,IF-Cell,FC
	TMEM119 Recombinant Rabbit mAb [PSH12-98]	HA723503	Human, Mouse, Rat	WB,IF-Cell,FC,IP
	TOM70 Recombinant Rabbit mAb [PSH13-00]	HA723505	Human, Mouse, Rat, Monkey	WB,IHC-P,FC,IF-Cell,IP
	TOMM20 Recombinant Antibody [ST04-72] - Mouse IgG1 (Chimeric)	HA601454	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	TRAF1 Recombinant Rabbit mAb [PSH13-88]	HA723581	Human, Rat	WB,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	Transferrin Receptor (CD71) Recombinant Rabbit mAb [PSH13-56]	HA723545	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP,IF-Tissue
	Transketolase Recombinant Rabbit mAb [PSH13-70]	HA723559	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
Relative expression	TRPV1 Recombinant Mouse mAb [PSH13-22]	HA601435	Rat	IHC-P
Relative expression	ULBP3 Recombinant Rabbit mAb [PSH13-03]	HA723508	Human	WB,IHC-P,IP
	VGLUT1 Recombinant Antibody [PSH10-53] - Mouse IgG1 (Chimeric)	HA601450	Mouse, Rat	IHC-Fr,IHC-P,WB
	VGLUT1 Recombinant Antibody [PSH10-53] - Guinea pig IgG2 (Chimeric)	HA601449	Mouse, Rat	IHC-Fr,WB
	VGLUT2 Recombinant Antibody [PSH10-54] - Mouse IgG1 (Chimeric)	HA601451	Mouse, Rat	IHC-Fr,IHC-P,WB
	VIP Recombinant Antibody [PSH10-56] - Guinea pig IgG2 (Chimeric)	HA601446	Mouse	IHC-Fr
	ZDHHC24 Rabbit pAb	HA500572	Human, Mouse, Rat	WB,IHC-P
	ZNF526 Rabbit pAb	HA500578	Human	IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层