

突触研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	14-3-3 Recombinant Rabbit mAb [JJ083-6]	ET1701-2	Human, Mouse, Rat	WB,IHC-P
	14-3-3 Recombinant Rabbit mAb [JJ084-3]	ET1701-4	Human, Mouse, Rat	WB
KD	14-3-3 alpha+beta Rabbit pAb	0407-5	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	14-3-3 alpha+beta Mouse mAb [9-26]	M0407-14	Human, Mouse, Zebrafish, Rat	WB,IF-Cell,IHC-P
KD	14-3-3 alpha+beta Recombinant Rabbit mAb [SD0837]	ET1612-99	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	14-3-3 beta+zeta Rabbit pAb	ER62524	Human, Mouse, Rat	WB,IHC-P
KD	14-3-3 epsilon Recombinant Rabbit mAb [JJ08-40]	ET1701-1	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P
	14-3-3 family protein Rabbit pAb	ER1903-06	Wheat, Oat, Apple	WB
KD	14-3-3 gamma Recombinant Rabbit mAb [SD20-65]	ET1612-9	Human, Mouse, Rat	WB,IHC-P,FC
	14-3-3 sigma Recombinant Rabbit mAb [SD2070]	ET1612-41	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	iFluor™ 488 Conjugated 14-3-3 sigma Recombinant Rabbit mAb [SD2070]	HA723707F	Human	IF-Tissue
KD	14-3-3 Theta Recombinant Rabbit mAb [SD084-02]	ET1612-97	Human, Rat, Mouse	WB,IF-Cell,IF-Tissue,FC
	14-3-3 zeta Rabbit pAb	ER62525	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell
	5HT1A Receptor Rabbit pAb	ER1903-07	Mouse, Rat	WB
	5HT2A Receptor Rabbit pAb	ER1903-09	Human, Mouse, Rat	WB
	5HT2B Receptor Rabbit pAb	ER1903-10	Human	WB,IHC-P,IF-Tissue
	5HT2C Receptor Recombinant Rabbit mAb [JA31-31]	ET1704-87	Human, Mouse, Rat	WB,IHC-P
	5HT3B receptor Rabbit pAb	ER1903-11	Human, Rat	WB,IHC-P,IF-Tissue,FC
	5HT4 Receptor Rabbit pAb	ER1903-13	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	5-HTR3 Rabbit pAb	ER1903-14	Human, Mouse, Rat	WB
	ABAT / GABA-T Recombinant Rabbit mAb [JE50-16]	ET7110-01	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	ADAR2 Rabbit pAb	ER62602	Human, Mouse, Rat	WB,IHC-P
	Adenosine A1-R Rabbit pAb	ER62609	Human, Rat	WB,IF-Cell
	Adenosine A1 Receptor Recombinant Rabbit mAb [JE31-35]	HA721393	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Adenosine A3 Receptor Rabbit pAb	ER1903-40	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	Adenosine Receptor A2a Rabbit pAb	ER1903-41	Human, Mouse, Rat	WB
	ADK Rabbit pAb	HA500050	Human, Mouse, Rat	WB,IHC-P
	ADORA2B Rabbit pAb	ER1903-44	Human, Mouse, Rat	WB
	ADRA1B Rabbit pAb	ER1903-46	Mouse	WB
	ADRB1 Rabbit pAb	ER1903-47	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	ADRB3 Rabbit pAb	ER1903-48	Human, Rat	WB,IHC-P,IF-Tissue
	alpha 1a Adrenergic Receptor Recombinant Rabbit mAb [JE62-50]	HA720096	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	alpha 2a Adrenergic Receptor Rabbit pAb	ER62703	Human, Mouse, Rat	WB
	Amphiphysin I Rabbit pAb	ER62659	Human, Mouse, Rat	WB
	AR alpha 1D Rabbit pAb	ER62745	Human	IHC-P,IF-Cell
	AR alpha 2B Rabbit pAb	ER62704	Human	IF-Cell
	AR alpha 2C Rabbit pAb	ER62705	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	AR beta 3 Rabbit pAb	ER62746	Human, Mouse, Rat	WB
	BEGAIN Rabbit pAb	ER62810	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	beta 1 Sodium Potassium ATPase Rabbit pAb	ER60078	Mouse, Rat	WB
	beta 2 Adrenergic Receptor Recombinant Rabbit mAb [JM102-06] - BSA and Azide free	HA750368	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IP,IF-Cell,FC,IHC-Fr
	beta 2 Adrenergic Receptor Recombinant Rabbit mAb [JM102-06]	ET1703-04	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IP,IF-Cell,FC,IHC-Fr
	BRSK1 Rabbit pAb	ER62832	Human, Mouse	WB,IF-Cell
	C1QA Recombinant Rabbit mAb [JE66-59]	HA721439	Human	WB,IHC-P
	C1QA Mouse mAb [G3-C7]	EM1707-50	Human	IHC-P,FC
	Calmodulin Recombinant Rabbit mAb [SJ16-09] - BSA and Azide free	HA750106	Human, Mouse, Dog, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Calmodulin Recombinant Rabbit mAb [SJ16-09]	ET1606-46	Human, Mouse, Dog, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Calpastatin Mouse mAb [A6B1]	HA600072	Human	WB,IF-Cell,IHC-P,FC
KD	Calreticulin Recombinant Rabbit mAb [SU37-03] - BSA and Azide free	HA750156	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	Calreticulin Recombinant Antibody [SU37-03] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610363	Human, Mouse, Rat, Monkey	WB,IHC-P
KD	Calreticulin Recombinant Rabbit mAb [SU37-03]	ET1608-60	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Calreticulin Recombinant Antibody [SU37-03] - Mouse IgG1 (Chimeric)	HA601547	Human, Mouse, Rat, Monkey	WB,IHC-P
	CAMK2D Rabbit pAb	R1305-2	Human	WB,IF-Cell,IHC-P
	CAMK2G Mouse mAb [10C-G1]	EM1706-16	Human, Rat	WB,IHC-P,FC
	CaMKII Recombinant Rabbit mAb [SU03-57] - BSA and Azide free	HA750150	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF- Tissue,IHC-P,IHC- Fr,FC
	CaMKII Recombinant Rabbit mAb [SU03-57]	ET1608-47	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF- Tissue,IHC-P,IHC- Fr,FC
	CaMKII alpha Recombinant Rabbit mAb [JM11-07]	ET1703-54	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF- Tissue,IHC-Fr,IF- Cell,FC,IP
	CaMKII alpha/beta/delta Rabbit pAb	ER62885	Human, Mouse, Rat	IHC-P,WB
	CaMKII alpha/delta Rabbit pAb	ER62886	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	CaMKII beta/gamma Rabbit pAb	ER62887	Human, Mouse, Rat	WB,IF-Cell
	CaMKII beta/gamma/delta Rabbit pAb	ER62888	Human, Mouse, Rat	WB,IHC-P
KD	CaMKII delta Recombinant Rabbit mAb [JA30-07]	ET1704-37	Human, Mouse, Rat, Zebrafish	WB,IHC-P
KD	CASK Mouse mAb [A1C5]	EM1902-09	Human, Mouse, Rat	WB,IHC-P
KD	CASK Recombinant Mouse mAb [A1C5-R] - BSA and Azide free	HA610132	Human, Mouse, Rat	WB,IHC-P
	CASK Rabbit pAb	ER1803-62	Human, Mouse, Rat	WB,IF-Cell,FC
KD	CASK Recombinant Mouse mAb [A1C5-R]	HA601257	Human, Mouse, Rat	WB,IHC-P
KD	CASK Mouse mAb [A1C7]	EM1902-08	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Caspr Mouse mAb [A2D2]	EM1902-10	Human, Mouse	IHC-P,FC
	Caspr Mouse mAb [A2D3]	EM1902-07	Human, Mouse	WB,IF-Cell,IHC- P,FC
Relative expression, KD	Caspr Recombinant Rabbit mAb [JE31-54]	HA722531	Human, Mouse, Rat	WB,IHC-P,IHC-Fr
	Caspr2 Recombinant Rabbit mAb [JE41-03]	HA721598	Human, Mouse, Rat	WB,IHC-P
KD	CD90 / THY1 Recombinant Rabbit mAb [JF10-09] - BSA and Azide free	HA750364	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P,IHC-Fr
	CD90 / THY1 Recombinant Rabbit mAb [SU35-07] - BSA and Azide free	HA750149	Human, Rat	WB,IHC-P
	CD90 / THY1 Recombinant Rabbit mAb [SU35-07]	ET1608-46	Human, Rat	WB,IHC-P
	CD90 / THY1 Rabbit pAb	R1107-9	Human, Mouse, Rat	WB
KD	CD90 / THY1 Recombinant Rabbit mAb [JF10-09]	ET1702-92	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CDKAL1 Rabbit pAb	ER62968	Human, Rat	WB
	Cellubrevin Recombinant Rabbit mAb [JG36-31]	ET7108-31	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	CHRM4 Rabbit pAb	ER1906-24	Human, Mouse	WB
	Complement C3 Mouse mAb [6-A10]	M1510-39	Human	WB,IHC-P,ELISA
	Complement C3 Mouse mAb [10-B3]	M1510-40	Human	WB,IHC-P,ELISA
	Complement C3 Recombinant Rabbit mAb [JU42-68]	ET7106-75	Human	WB,IHC-P,IP
	Complement C3 Recombinant Rabbit mAb [JF10-30] - BSA and Azide free	HA750367	Human, Mouse, Rat	WB,IF-Cell,FC
	Complement C3 Recombinant Rabbit mAb [PSH18-31]	HA723984	Mouse, Rat	WB,IHC-Fr,IHC-P
	Complement C3 Recombinant Rabbit mAb [JF10-30]	ET1702-99	Human, Mouse, Rat	WB,IF-Cell,FC
	Complement C3 (beta chain) Rabbit pAb	0809-8	Human	WB
	Complement C3 (beta chain) Mouse mAb [11-10]	M0809-14	Human	WB,IF-Cell,FC
	CPEB1 Rabbit pAb	ER1902-92	Human, Rat	WB,IHC-P,FC
	CPEB1 Rabbit pAb	ER1902-89	Human, Mouse, Rat	WB,IHC-P,FC
	CPLX1 Rabbit pAb	ER1906-61	Mouse, Rat	WB,IHC-P,IF-Tissue
	CSEN Rabbit pAb	ER1906-75	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	CSP Rabbit pAb	ER63138	Human, Mouse, Rat	WB
KD	CTBP2 Recombinant Rabbit mAb [JE56-00]	ET7111-27	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Dopamine D3 receptor Rabbit pAb	ER1907-71	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Dopamine D4 receptor Rabbit pAb	ER1907-72	Human, Mouse	WB
	Dopamine Receptor D3 Recombinant Rabbit mAb [JE33-33]	HA721391	Human, Mouse, Rat	WB
	Dopamine Transporter Recombinant Antibody [PSH10-49] - Mouse IgG1 (Chimeric)	HA601476	Mouse, Rat	IHC-Fr,IHC-P
	Dopamine Transporter Recombinant Rabbit mAb [PSH10-49]	HA723213	Mouse, Rat, Human, Cynomolgus monkey, Pig	WB,IHC-Fr,IHC-P
	Dopamine Transporter Recombinant Antibody [PSH10-49] - Chicken IgY (Chimeric) - BSA and Azide free	HA610308	Mouse	IHC-Fr
	Dopamine Transporter Recombinant Antibody [PSH10-49] - Chicken IgY (Chimeric)	HA601474	Mouse	IHC-Fr
	Dopamine Transporter Recombinant Rabbit mAb [PSH10-49] - BSA and Azide free	HA751352	Mouse, Rat, Human, Cynomolgus monkey, Pig	WB,IHC-Fr,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Dopamine Transporter Recombinant Antibody [PSH10-49] - Rat IgG1 (Chimeric)	HA601404	Mouse, Rat, Human, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
	Dopamine Transporter Recombinant Antibody [PSH10-49] - Guinea pig IgG2 (Chimeric)	HA601475	Mouse	IHC-Fr
	Dopamine Transporter Recombinant Antibody [PSH10-49] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610244	Mouse, Rat, Human, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
	Dopamine Transporter Rabbit pAb	HA500007	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Dopamine Transporter Rabbit pAb	ER1907-74	Rat	WB
	Dopamine Transporter Recombinant Antibody [PSH10-49] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610310	Mouse, Rat	IHC-Fr,IHC-P
	Dopamine Transporter Recombinant Antibody [PSH10-49] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610309	Mouse	IHC-Fr
	DRD1 Recombinant Rabbit mAb [JM10-93]	ET1703-45	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
	Dynamin 1 Rabbit pAb	ER1802-77	Mouse, Rat, Human	WB,IHC-P
	Dynamin 1 Mouse mAb [A1A12]	EM1901-44	Human	WB,IF-Cell,FC
	Dynamin 1 Recombinant Mouse mAb [A1B2-R]	HA601325	Human, Mouse, Rat	WB,IHC-P
	Dynamin 1 Recombinant Mouse mAb [A1B2-R] - BSA and Azide free	HA610197	Human, Mouse, Rat	WB,IHC-P
	Dynamin 1 Recombinant Mouse mAb [A1B1-R] - BSA and Azide free	HA610105	Mouse, Rat	WB
	Dynamin 1 Mouse mAb [A1B2]	EM1901-42	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Dynamin 1 Recombinant Mouse mAb [A1B1-R]	HA601230	Mouse, Rat	WB
	Dynamin 1 Mouse mAb [A1B1]	EM1901-43	Human, Mouse, Rat	IHC-P,WB,IF-Cell
	Dynamin 2 Recombinant Rabbit mAb [JM61-38]	ET1705-2	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Dynamin 3 Recombinant Rabbit mAb [JE54-86]	HA722870	Human, Mouse	WB
	Dynamin 3 Rabbit pAb	HA500032	Human, Mouse, Rat	WB,IHC-P
	Dynamin 3 Rabbit pAb	HA500037	Human, Rat	WB,IF-Cell
	ELKS Recombinant Rabbit mAb [PSH17-33] - BSA and Azide free	HA751621	Human, Mouse, Rat	WB
	ELKS Recombinant Rabbit mAb [PSH17-33]	HA723916	Human, Mouse, Rat	WB
	ELKS Rabbit pAb	HA500376	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	EphA2/3/4 Rabbit pAb	ER63372	Human, Rat	WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	EphA2 (N-term) Rabbit pAb	ER2001-31	Human, Mouse	WB,IHC-P,FC
	EphA3 Rabbit pAb	ER1908-23	Mouse, Rat	WB,IHC-P,IF-Tissue
	EphA6 Rabbit pAb	ER63373	Human, Mouse, Rat	IHC-P,IF-Cell
	EphA7 Rabbit pAb	ER1908-24	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	EphB1 Rabbit pAb	ER63374	Human, Mouse, Rat, Monkey	WB
	EphB2 Rabbit pAb	ER1908-25	Rat, Mouse, Human	WB,IHC-P,IF-Tissue
	EphB4 Rabbit pAb	ER63378	Human, Mouse	WB,IHC-P,IF-Cell
	EphB6 Rabbit pAb	ER1908-26	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Eph receptor A2 Recombinant Rabbit mAb [PSH09-21]	HA723079	Human	WB,IF-Cell,IP,FC
	Eph receptor A2 Recombinant Rabbit mAb [PSH09-21] - BSA and Azide free	HA751275	Human	WB,IF-Cell,IP,FC
	Eph receptor B1 Recombinant Rabbit mAb [JE39-46]	HA722345	Human, Mouse, Rat, Monkey	WB
Relative expression	Eph receptor B2 Recombinant Rabbit mAb [PSH14-91]	HA723693	Human	WB,IF-Cell,IP
Relative expression	Eph receptor B2 Recombinant Rabbit mAb [PSH14-91] - BSA and Azide free	HA751535	Human	WB,IF-Cell,IP
	Eph receptor B3 Recombinant Rabbit mAb [JE54-73]	ET7110-80	Human, Mouse, Rat	WB,IHC-P
	FHIT Recombinant Rabbit mAb [JE33-58]	HA722390	Human, Mouse, Rat	WB,IHC-P,FC,IF-Tissue
	FMRP Recombinant Rabbit mAb [JM91-41]	ET1703-70	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	FXR2 Rabbit pAb	ER63480	Human, Mouse	WB,IHC-P,IF-Cell
	GABAA R delta Rabbit pAb	ER63489	Human, Mouse, Rat	WB,IHC-P
	GABA A Receptor alpha 3 Rabbit pAb	ER60110	Human, Rat, Mouse	WB,IHC-P
	GABA A Receptor gamma 1 Rabbit pAb	HA500085	Human, Mouse, Rat	WB
	GABA A Receptor gamma 2 Rabbit pAb	ER60112	Human, Rat, Mouse	WB,IHC-P
	GABA A receptor pi Rabbit pAb	ER1909-21	Mouse	WB
	GABAA R epsilon Rabbit pAb	ER65419	Human	WB
	GABA(B) R2 Mouse mAb [3G3]	RT1238	Human, Mouse, Rat	WB,IP,IF
	GABA B Receptor 1 Recombinant Rabbit mAb [JE60-09]	HA721338	Human, Mouse	WB,IHC-P,IF-Tissue
Relative expression	GABA B Receptor 1 Recombinant Rabbit mAb [PSH03-87]	HA722053	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	GABA B Receptor 1 Recombinant Rabbit mAb [PSH03-87] - BSA and Azide free	HA750912	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P
Relative expression	GABA B Receptor 1 Recombinant Rabbit mAb [PSH06-21] - BSA and Azide free	HA751059	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr
Relative expression	GABA B Receptor 1 Recombinant Rabbit mAb [PSH06-21]	HA722653	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr
	GABA B Receptor 2 Recombinant Rabbit mAb [JU31-32]	ET7106-88	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	GABARAP Recombinant Rabbit mAb [JM30-30]	ET1705-53	Human, Mouse, Rat	WB,IHC-P
	GABARAPL2 Rabbit pAb	ER1901-73	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	GABA Transporter 1 Rabbit pAb	ER60113	Human, Rat, Mouse	WB,IHC-P
	GABA Transporter 1 / GAT 1 Recombinant Rabbit mAb [JE34-61]	HA721496	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	GABA Transporter 3 Rabbit pAb	ER1909-22	Human	WB,FC
Relative expression	GABRA1 Recombinant Mouse mAb [PSH10-46]	HA601383	Mouse, Rat	WB
	GABRA1 Rabbit pAb	ER1901-58	Human, Rat	IF-Cell,IHC-P,FC
	GABRA1 Rabbit pAb	ER1901-59	Human, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	GABRA1 Recombinant Mouse mAb [PSH10-46] - BSA and Azide free	HA610227	Mouse, Rat	WB
	GABRA2 Rabbit pAb	ER1909-25	Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	GABRA2 Recombinant Rabbit mAb [PSH17-28]	HA723915	Mouse, Rat	IHC-Fr,IHC-P,WB,IP
	GABRA2 Recombinant Rabbit mAb [PSH17-28] - BSA and Azide free	HA751620	Mouse, Rat	IHC-Fr,IHC-P,WB,IP
	GABRA2 Rabbit pAb	HA500288	Human, Mouse, Rat	WB,FC
	GABRA4 Rabbit pAb	ER1901-60	Human, Mouse, Rat	WB,IHC-P
	GABRA5 Rabbit pAb	ER1902-82	Human, Rat, Mouse	WB,IF-Tissue,IHC-P,FC
	GABRA5 Recombinant Rabbit mAb [JB34-19]	ET7107-08	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P
	GABRA6 Rabbit pAb	HA500473	Human, Mouse, Rat	WB,IHC-P
	GABRB2 Rabbit pAb	HA500289	Human, Mouse, Rat	WB,IHC-P,FC
	GABRB2 Rabbit pAb	ER1909-26	Mouse, Rat	WB,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	GABR B3 Rabbit pAb	ER1909-24	Mouse, Rat	WB,FC
	GABRB3 Rabbit pAb	HA500090	Human, Mouse, Rat	WB,IHC-P
	GABRD Rabbit pAb	ER1901-54	Human	WB,FC
	GABRG3 Rabbit pAb	ER1909-27	Human	WB
	GABRR1 Rabbit pAb	ER1909-28	Human, Mouse, Rat	WB
	GAD65 Recombinant Rabbit mAb [PSH06-80]	HA722750	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
	GAD65 Recombinant Rabbit mAb [PSH06-80] - BSA and Azide free	HA751102	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
	GAD-65/67 Rabbit pAb	ER63491	Human, Mouse	WB,IHC-P,IF-Cell
	GAD67 Recombinant Mouse mAb [PSH08-30] - BSA and Azide free	HA610216	Mouse, Rat	WB,IHC-P,IHC-Fr
	GAD67 Recombinant Mouse mAb [PSH08-30]	HA601370	Mouse, Rat	WB,IHC-P,IHC-Fr
	GAD67 Recombinant Rabbit mAb [JM11-11]	ET1703-71	Human, Mouse, Rat	WB,IP,IHC-P
Relative expression, KD	GAP43 Recombinant Rabbit mAb [SC60-06]	ET1610-94	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Relative expression	GAP43 Rabbit pAb	ER40201	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Gephyrin Recombinant Rabbit mAb [JG34-47]	ET7108-04	Human, Mouse, Rat	WB,IHC-P
	GluA4 Recombinant Rabbit mAb [JE36-81]	HA721582	Human, Mouse, Rat	WB,IHC-P,FC
Relative expression	GluR2 Recombinant Rabbit mAb [JU20-43]	ET1706-06	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	GLUR4 Rabbit pAb	ER1909-72	Human, Mouse	WB,IHC-P,IF-Tissue
	Glutamate receptor 1 Recombinant Rabbit mAb [SD2010]	ET1612-10	Human, Mouse, Rat	WB,IHC-P,IHC-Fr
	Glutamate receptor 1 Recombinant Rabbit mAb [SD2010] - BSA and Azide free	HA750281	Human, Mouse, Rat	WB,IHC-P,IHC-Fr
	GPR156 Rabbit pAb	ER63583	Human	WB,IF-Cell
	GRIA3 Mouse mAb [E2-D1]	EM1709-87	Human, Rat	WB,IHC-P
	GRK3 Rabbit pAb	HA500317	Human, Rat	WB,IF-Cell,FC
	GRM3 Mouse mAb [G2-G9]	EM1706-64	Human	WB,IF-Cell,IHC-P,FC
	HOME1 Rabbit pAb	ER60781	Human, Rat, Mouse	WB
	HTR3A Mouse mAb [3B11C10]	EM1710-35	Human	IHC-P,IF-Cell,FC
	HTR3B Mouse mAb [2H9H5]	EM1710-36	Human	WB,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human Complement C3 Recombinant Rabbit mAb [PSH15-08] - Detector	HA725142B	Human	ELISA(Det),ELISA
	Human Complement C3 Recombinant Rabbit mAb [PSH15-07] - BSA and Azide free (Capture)	HA725140	Human	ELISA(Cap)
	Human Complement C3 Recombinant Rabbit mAb [PSH15-08] - BSA and Azide free (Detector)	HA725141	Human	ELISA(Det)
	HRP Conjugated Human Complement C3 Recombinant Rabbit mAb [PSH15-08] - Detector	HA725142H	Human	ELISA(Det),ELISA
	Human EphA2 Recombinant Rabbit mAb [PSH15-12] - BSA and Azide free (Capture)	HA723714	Human	ELISA(Cap)
	Biotin Conjugated Human EphA2 Recombinant Rabbit mAb [PSH15-13] - Detector	HA723716B	Human	ELISA(Det),ELISA
	Human EphA2 Recombinant Rabbit mAb [PSH15-13] - BSA and Azide free (Detector)	HA723715	Human	ELISA(Det)
	Human/Mouse/Rat Synapsin I Recombinant Rabbit mAb [PSH18-41] - BSA and Azide free (Detector)	HA723994	Human, Mouse, Rat	ELISA(Det)
	Human/Mouse/Rat Synapsin I Recombinant Rabbit mAb [PSH18-40] - BSA and Azide free (Capture)	HA723993	Human, Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Human/Mouse/Rat Synapsin I Recombinant Rabbit mAb [PSH18-41] - Detector	HA723994B	Human, Mouse, Rat	ELISA(Det),ELISA
	Biotin Conjugated Human SIRP alpha Recombinant Rabbit mAb [PSH17-64] - Detector	HA725348B	Human	ELISA(Det),ELISA
	Human SIRP alpha Recombinant Rabbit mAb [PSH17-63] - BSA and Azide free (Capture)	HA725347	Human	ELISA(Cap)
	Human SIRP alpha Recombinant Rabbit mAb [PSH17-64] - BSA and Azide free (Detector)	HA725348	Human	ELISA(Det)
	Ionotropic glutamate receptor 2 + 3 Recombinant Rabbit mAb [JE33-28]	HA721405	Mouse, Rat, Human	WB
	KCNA4 Rabbit pAb	ER1706-36	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Relative expression	KCNA5 / KV1-5 Recombinant Rabbit mAb [PSH12-30]	HA723455	Human, Mouse, Rat	WB
	KCNA5 / KV1-5 Rabbit pAb	ER1803-23	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	KCNA5 / KV1-5 Recombinant Rabbit mAb [PSH12-30] - BSA and Azide free	HA751444	Human, Mouse, Rat	WB
	KCNG3 Rabbit pAb	ER63794	Human, Mouse, Rat	WB
	KCNQ2/3/4/5 Rabbit pAb	ER63796	Human, Mouse, Rat	IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KCNQ4 Rabbit pAb	ER2001-23	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KCNQ5 Rabbit pAb	ER63798	Human, Mouse	WB
Kv10.1 Rabbit pAb	ER60123	Human, Rat, Mouse	IHC-P
KV1.1 Rabbit pAb	ER1901-23	Human, Mouse, Rat	WB,IHC-P,FC
Kv1.1 potassium channel Recombinant Rabbit mAb [SN66-06]	ET1611-7	Human, Rat, Mouse	WB,IHC-P,IP
Kv1.4 Rabbit pAb	ER1902-49	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Kv3.4 Rabbit pAb	ER63826	Human, Mouse, Monkey	WB,IHC-P,IF-Cell
LAP1B Rabbit pAb	ER63840	Human	WB
mAChR M5 Rabbit pAb	ER63888	Human	WB,IF-Cell
mGluR1 Recombinant Rabbit mAb [JM11-61]	ET1703-84	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
mGluR2 Recombinant Rabbit mAb [JJ091-09]	ET1701-74	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
mGluR3 Recombinant Rabbit mAb [JE41-18]	HA722434	Human, Mouse, Rat	WB,IHC-P
mGluR5 Recombinant Rabbit mAb [ST51-00]	ET1609-36	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
HRP Conjugated Mouse/Rat Complement C3 Recombinant Rabbit mAb [PSH16-82] - Detector	HA725145H	Mouse, Rat	ELISA(Det),ELISA
Mouse/Rat Complement C3 Recombinant Rabbit mAb [PSH16-81] - BSA and Azide free (Capture)	HA725143	Mouse, Rat	ELISA(Cap)
Mouse/Rat Complement C3 Recombinant Rabbit mAb [PSH16-82] - BSA and Azide free (Detector)	HA725144	Mouse, Rat	ELISA(Det)
Biotin Conjugated Mouse/Rat Complement C3 Recombinant Rabbit mAb [PSH16-82] - Detector	HA725145B	Mouse, Rat	ELISA(Det),ELISA
Munc18-1 Recombinant Rabbit mAb [JE38-54]	HA722250	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
muscarinic Acetylcholine Receptor 1 Rabbit pAb	ER1913-39	Human, Mouse, Rat	WB
Muscarinic Acetylcholine Receptor 2 Recombinant Rabbit mAb [JE38-35]	HA721490	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Na ⁺ CP-pan Rabbit pAb	ER64102	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Na ⁺ CP type II alpha Rabbit pAb	ER64098	Human, Mouse, Rat	WB,IHC-P
Na ⁺ CP type IV alpha Rabbit pAb	ER64100	Human, Mouse, Rat	WB
Na ⁺ /K ⁺ ATPase Mouse mAb [1G1]	RT1412	Human, Mouse, Rat	WB,IP,IF,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Na ⁺ /K ⁺ -ATPase alpha 1 Rabbit pAb	ER64103	Rat	WB,IF-Cell
	NaV1.7 Mouse mAb [1D6]	EM1701-35	Human	IHC-P,IF-Cell,FC,WB
	NaV1.7 Rabbit pAb	ER1706-57	Human, Rat	IF-Cell,IHC-P
	NaV1.7 Rabbit pAb	ER1706-44	Human	IF-Cell,IHC-P,FC
	Neurexin I Rabbit pAb	ER64160	Human, Mouse, Rat	WB
	Neurexin III beta Rabbit pAb	ER64161	Human	WB
	Neuroglycan C Rabbit pAb	ER1913-72	Human	WB
	Neurologin 1 Rabbit pAb	ER1913-74	Human, Mouse, Rat	WB,FC
	NLGN3 Rabbit pAb	ER60922	Human, Mouse, Rat	WB
	NMBR Rabbit pAb	HA500253	Human, Mouse, Rat	WB,IHC-P
	NMDA epsilon 4 Rabbit pAb	ER64208	Human, Mouse, Rat, Monkey	WB
	NMDA epsilon 4 Rabbit pAb	ER64207	Human, Mouse, Rat	IHC-P
Relative expression	NMDAR1 Recombinant Rabbit mAb [JM11-26] - BSA and Azide free	HA750391	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,IHC-Fr
Relative expression	NMDAR1 Recombinant Rabbit mAb [JM11-26]	ET1703-75	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,IHC-Fr
	NMDAR2A Recombinant Rabbit mAb [JA31-20]	ET1704-80	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue
	NMDAR2B Rabbit pAb	ER1913-90	Human, Rat	WB,IHC-P,IF-Tissue
Relative expression	NMDAR2B Recombinant Rabbit mAb [PSH05-46]	HA722284	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
	NMDAR2B Rabbit pAb	HA500208	Mouse, Rat, Human	WB,IF-Cell,IHC-P,FC
Relative expression	NMDAR2B Recombinant Rabbit mAb [PSH05-46] - BSA and Azide free	HA751005	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
	NMDAR2C Rabbit pAb	ER1913-91	Mouse, Rat	WB,IHC-P,IF-Tissue
	NPTX2 Rabbit pAb	ER1914-08	Mouse	WB,IHC-P,IF-Tissue
	NR3B Rabbit pAb	HA500079	Human	WB
	NR3B Rabbit pAb	HA500134	Human	WB
	P2RY1 Mouse mAb [4G5D6]	EM1710-48	Human, Rat	WB,IHC-P,FC
	P2RY4 Mouse mAb [3B4C11]	EM1710-51	Human	WB,IHC-P
	P2RY8 Mouse mAb [1G5A11]	EM1710-43	Human	WB,IHC-P,IF-Cell,FC
	P2Y10 Rabbit pAb	ER1914-60	Human	WB,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	P2Y11 Rabbit pAb	ER64494	Human	WB
	P2Y12 Recombinant Rabbit mAb [PSH04-48]	HA722133	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr,IF-Tissue
	P2Y12 Recombinant Antibody [PSH04-48] - Rat IgG1 (Chimeric)	HA601388	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
	P2Y12 Recombinant Rabbit mAb [JE63-23]	HA721011	Human, Mouse	WB
	P2Y12 Recombinant Antibody [PSH04-48] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610231	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
	P2Y12 Recombinant Rabbit mAb [PSH04-48] - BSA and Azide free	HA750957	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr,IF-Tissue
	P2Y4 Rabbit pAb	ER64495	Human, Mouse, Rat	WB,IF-Cell
	P2Y6 Recombinant Rabbit mAb [JG37-77]	ET7108-52	Human	WB,IP,IHC-P,IF-Cell
	P2Y8 Rabbit pAb	ER1914-62	Human	WB,FC
	P2Y9 Rabbit pAb	ER1914-63	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Phospho-CaMKII (T286) Mouse mAb [2G1]	RT1092	Human, Mouse, Rat	WB,IP,IF
Cell treatment	Phospho-CaMKII (T286) Recombinant Rabbit mAb [JE51-81]	HA721794	Human, Mouse, Rat	WB
	Phospho-GluR1 (AMPA subtype) (S831) Recombinant Rabbit mAb [JE44-96]	HA721865	Human, Mouse, Rat	WB
Cell treatment	Phospho-GluR1 (S845) Recombinant Rabbit mAb [JJ2009]	ET1701-28	Mouse, Rat, Human	WB,FC
	Phospho-Synapsin I (S9) Recombinant Rabbit mAb [SN63-04]	ET1611-26	Mouse, Rat, Human	WB,IHC-P
Relative expression	Phospho-TrkB (Y817) Recombinant Rabbit mAb [SC0556]	ET1610-35	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Relative expression	Phospho-TrkB (Y817) Recombinant Rabbit mAb [SC0556] - BSA and Azide free	HA750213	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Profilin 1 Recombinant Rabbit mAb [JA30-13]	ET1704-42	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	PSD93 Recombinant Rabbit mAb [JM41-04]	ET1705-90	Human, Mouse, Rat	WB,FC,IP
	PSD95 Mouse mAb [B5-E6]	M1511-4	Human, Mouse, Rat	WB,IF-Cell,FC
	PSD95 Recombinant Rabbit mAb [SR38-09] - BSA and Azide free	HA750035	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,IP,IHC-Fr
Relative expression	PSD95 Recombinant Antibody [SR38-09] - Rat IgG1 (Chimeric)	HA601416	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	PSD95 Recombinant Antibody [SR38-09] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610256	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,WB
	PSD95 Rabbit pAb	R1508-4	Human, Mouse, Rat	WB,IHC-P,FC
	PSD95 Recombinant Rabbit mAb [SR38-09]	ET1602-20	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,IHC-Fr
	PTPRF/LAR Rabbit pAb	R1510-12	Human	IF-Cell,IHC-P,FC
	Rab4A Recombinant Rabbit mAb [SR27-08]	ET1602-15	Human, Rat, Monkey, Mouse	WB,IF-Cell,IF-Tissue,IP
KD	Rab5 Recombinant Rabbit mAb [ST47-04]	ET1609-27	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KD	Rab5 Recombinant Rabbit mAb [ST47-04] - BSA and Azide free	HA750167	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Rab5a Mouse mAb [C11-H5]	EM1708-08	Human	WB,IF-Cell,IHC-P,FC
	Rabphilin-3A Rabbit pAb	ER64738	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	RIM 2 Mouse mAb [3G2]	RT1539	Human, Rat	WB,IP
	rSec6 Recombinant Rabbit mAb [JE64-12]	HA722470	Human, Mouse, Rat	WB
	SAP97 Rabbit pAb	ER2001-07	Human, Mouse	WB,IHC-P,FC
KD	SAP97 Recombinant Rabbit mAb [JE39-74]	HA721502	Human, Mouse, Rat	WB,IHC-P,FC
	SCFD1 Rabbit pAb	ER64914	Human, Mouse, Rat	WB
	SCFD2 Rabbit pAb	ER1916-27	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	SCN1A Rabbit pAb	ER65886	Human, Mouse, Rat	IHC-P,WB
	SERCA3 Rabbit pAb	ER1916-41	Human	WB
Relative expression	Serotonin transporter Recombinant Rabbit mAb [PSH19-18]	HA724052	Mouse, Rat	IHC-Fr,IHC-P,WB
	Serotonin transporter Recombinant Rabbit mAb [JE34-53]	HA721544	Human, Mouse, Rat	WB,FC
	Serotonin transporter Rabbit pAb	ER1916-44	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Shank 2 Rabbit pAb	ER64959	Human, Mouse, Rat	WB,IHC-P
Relative expression	SHANK2 Recombinant Rabbit mAb [PSH10-61] - BSA and Azide free	HA751364	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IF-Tissue
Relative expression	SHANK2 Recombinant Rabbit mAb [PSH10-61]	HA723226	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IF-Tissue
KD	SIRP alpha Recombinant Rabbit mAb [JE36-20]	HA721316	Human, Mouse, Rat	WB,IHC-P,FC
	SIRP alpha Mouse mAb [A3E3]	EM1902-37	Human	WB,IF-Cell,IHC-P,FC
	SIRP Alpha Rabbit pAb	ER1916-61	Human, Mouse, Rat	WB,FC
	SNAP23 Recombinant Rabbit mAb [JA73-15]	ET1704-95	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	SNAP25 Recombinant Rabbit mAb [PSH17-65] - BSA and Azide free (Capture)	HA725227	Human, Mouse	ELISA(Cap)
Relative expression	SNAP25 Recombinant Rabbit mAb [JE45-86]	HA721892	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Biotin Conjugated SNAP25 Recombinant Rabbit mAb [PSH17-66] - Detector	HA725228B	Human, Mouse	ELISA(Det),ELISA
	SNAP25 Recombinant Rabbit mAb [PSH17-66] - BSA and Azide free (Detector)	HA725228	Human, Mouse	ELISA(Det)
	SNAP25 Rabbit pAb	HA500539	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	SNAP29 Recombinant Rabbit mAb [PSH01-76]	HA721728	Human	WB
	SNAP29 Recombinant Rabbit mAb [PSH01-76] - BSA and Azide free	HA750750	Human	WB
KD	Sodium Potassium ATPase Recombinant Rabbit mAb [ST0533]	ET1609-76	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
	iFluor™ 488 Conjugated Sodium Potassium ATPase Recombinant Rabbit mAb [ST0533]	HA720174F	Human, Mouse, Rat	IF-Tissue
	iFluor™ 594 Conjugated Sodium Potassium ATPase Recombinant Rabbit mAb [ST0533]	HA720175F	Human, Mouse, Rat	IF-Tissue
	iFluor™ 647 Conjugated Sodium Potassium ATPase Recombinant Rabbit mAb [ST0533]	HA720176F	Human, Mouse, Rat	IF-Tissue,FC
KD	Sodium Potassium ATPase Recombinant Rabbit mAb [ST0533] - BSA and Azide free	HA750194	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
	SR-1A Rabbit pAb	ER65015	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	SR-1B Rabbit pAb	ER65016	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	SR-1D Rabbit pAb	ER65667	Human, Mouse, Rat	WB,IHC-P
	SR-1E Rabbit pAb	ER65017	Human	WB
	SR-3D Rabbit pAb	ER65018	Human	WB,IF-Cell
	SR-5A Rabbit pAb	ER65019	Human	WB,IF-Cell
	SR-6 Rabbit pAb	ER65020	Human, Mouse, Rat	WB
	SR-7 Rabbit pAb	ER65021	Human, Mouse, Rat	IF-Cell
	STX2 Rabbit pAb	ER61291	Human, Rat, Mouse	WB
Relative expression	SV2A Recombinant Rabbit mAb [PSH10-69]	HA723233	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IP
	SV2A Rabbit pAb	ER61279	Human, Rat, Mouse	WB
Relative expression	SV2A Recombinant Rabbit mAb [PSH10-69] - BSA and Azide free	HA751371	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	SYBL1 Rabbit pAb	ER1917-39	Mouse	WB
	Synapsin II Recombinant Rabbit mAb [JE49-03]	ET7109-88	Mouse, Rat, Human	WB,IF-Tissue,IHC-P,FC
Relative expression	Synapsin I + II Recombinant Rabbit mAb [PSH10-29]	HA723196	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-Fr,IHC-P,IP,IF-Tissue
Relative expression	Synapsin I + II Recombinant Rabbit mAb [PSH10-29] - BSA and Azide free	HA751346	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-Fr,IHC-P,IP,IF-Tissue
Relative expression	Synaptophysin Recombinant Rabbit mAb [SJ26-85]	ET1606-56	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	Synaptophysin Rabbit pAb	0407-2	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Synaptophysin Mouse mAb [A1-7]	M0407-30	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Synaptophysin Recombinant Rabbit mAb [JE42-51]	HA721599	Human, Mouse, Rat	WB,IHC
Relative expression	Synaptopodin Recombinant Rabbit mAb [PSH10-45]	HA723211	Human, Mouse, Rat	WB,IHC-Fr,IHC-P
	Synaptopodin Rabbit pAb	ER1917-40	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
Relative expression	Synaptopodin Recombinant Antibody [PSH10-45] - Rat IgG1 (Chimeric)	HA601402	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Synaptopodin Recombinant Antibody [PSH10-45] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610242	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Synaptopodin Recombinant Rabbit mAb [PSH10-45] - BSA and Azide free	HA751350	Human, Mouse, Rat	WB,IHC-Fr,IHC-P
	Synaptotagmin Rabbit pAb	ER65057	Human, Mouse, Rat	WB,IHC-P
	Synaptotagmin 1 Rabbit pAb	ER1917-41	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Syndecan 1 Recombinant Rabbit mAb [JM11-21] - BSA and Azide free	HA750381	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Syndecan 1 Recombinant mAb	IRS105RB	Human	mIHC
	Syndecan 1 Recombinant Rabbit mAb [JM11-21]	ET1703-42	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	SYNDIG1 Rabbit pAb	ER2001-43	Human, Mouse, Rat	WB,IHC-P,FC
	SynGAP Recombinant Rabbit mAb [JE37-00]	HA722728	Human, Mouse, Rat	WB,IHC-P
	SynGAP Recombinant Rabbit mAb [JE57-79]	HA721860	Mouse, Rat	WB,IHC-P,IF-Tissue
	SynGAP Rabbit pAb	ER1902-57	Human, Mouse, Rat	WB,IHC-P,FC
	SYNPR/Synaptoporin Rabbit pAb	ER1917-44	Rat	WB,IHC-P,IF-Tissue,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Syntaxin 16 Recombinant Rabbit mAb [JE62-35]	HA720107	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Syntaxin 1a Recombinant Rabbit mAb [JE54-37]	ET7110-68	Human, Mouse, Rat	WB,IHC-P,FC
KD	Syntaxin 3 Recombinant Rabbit mAb [JE40-96]	ET7109-19	Human, Mouse, Rat	WB,IHC-P
KD	Syntaxin 4 Recombinant Rabbit mAb [PSH0-80] - BSA and Azide free	HA750660	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
KD	Syntaxin 4 Recombinant Rabbit mAb [JE35-92]	HA721467	Human, Mouse, Rat	WB,IHC-P
KD	Syntaxin 4 Recombinant Rabbit mAb [PSH0-80]	HA721509	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
KD	Syntenin Recombinant Rabbit mAb [JE40-72]	ET7109-14	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	SYT2 Rabbit pAb	ER60348	Human, Mouse, Rat	WB
	Thy-1 (OX7) Mouse mAb [1G2]	RT1615	Mouse, Rat	WB,IP,IF,IHC-P,FC
Relative expression	TrkA Recombinant Rabbit mAb [SU0354] - BSA and Azide free	HA750148	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	TrkA Rabbit pAb	ER1706-86	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	TrkA Recombinant Rabbit mAb [SU0354]	ET1608-44	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Relative expression	TrkA+B+C Recombinant Rabbit mAb [JJ084-04]	ET1701-16	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	TrkB Recombinant Rabbit mAb [JE58-22]	ET7111-38	Mouse, Rat	WB,IHC-P
	Trk C Rabbit pAb	ER65473	Human	WB
Relative expression	TrkC Recombinant Rabbit mAb [PSH05-40]	HA722269	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Relative expression	TrkC Recombinant Rabbit mAb [PSH05-40] - BSA and Azide free	HA750999	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
	TRPC3 Rabbit pAb	ER65609	Human, Mouse, Rat	WB,IHC-P
	TRPM2 Rabbit pAb	HA500437	Human, Mouse, Rat	WB,FC
	TRPM5 Rabbit pAb	ER1918-17	Human, Mouse	WB
	TSLC1 Rabbit pAb	ER1918-21	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Unc18-1 Rabbit pAb	ER65267	Human, Mouse, Rat	IHC-P,IF-Cell
	UTRO Rabbit pAb	ER61389	Human	IHC-P
	VAMP1 Recombinant Rabbit mAb [JG40-19]	ET7109-03	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
	VAMP-1/2/3 Rabbit pAb	ER65291	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	VAMP2 Recombinant Rabbit mAb [JM11-00]	ET1703-50	Human, Mouse, Rat	WB,IF-Cell,FC,IP,IHC-P
	VAMP-4 Rabbit pAb	ER65292	Human, Mouse	IF-Cell
	VAMP8 Recombinant Rabbit mAb [JF0963]	ET1702-10	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	VAMP8 Recombinant Rabbit mAb [JF0963] - BSA and Azide free	HA750334	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	VAPA Recombinant Rabbit mAb [JE77-48]	HA722372	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	VGLUT1 Recombinant Antibody [PSH10-53] - Guinea pig IgG2 (Chimeric)	HA601449	Mouse, Rat	IHC-Fr,WB
	VGLUT1 Mouse mAb [A7A3]	HA601014	Human, Mouse, Rat	WB
	VGLUT1 Rabbit pAb	HA500488	Human, Mouse, Rat	WB
Relative expression	VGLuT1 Recombinant Antibody [PSH10-53] - Rat IgG1 (Chimeric)	HA601392	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB,IF-Tissue
Relative expression	VGLuT1 Recombinant Antibody [PSH10-53] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610234	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB,IF-Tissue
	VGLUT1 Recombinant Antibody [PSH10-53] - Chicken IgY (Chimeric)	HA601478	Mouse	IHC-Fr
	VGLUT1 Recombinant Antibody [PSH10-53] - Mouse IgG1 (Chimeric)	HA601450	Mouse, Rat	IHC-Fr,IHC-P,WB
	VGLUT1 Rabbit pAb	HA500083	Human, Mouse, Rat	WB,IHC-P
	VGLUT1 Recombinant Mouse mAb [PSH08-34]	HA601371	Mouse, Rat	WB,IHC-P
	VGLUT1 Recombinant Antibody [PSH10-53] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610283	Mouse, Rat	IHC-Fr,WB
Relative expression	VGLUT1 Recombinant Rabbit mAb [PSH10-53]	HA723217	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IHC-P,IF-Tissue
	VGLUT1 Recombinant Antibody [PSH10-53] - Chicken IgY (Chimeric) - BSA and Azide free	HA610312	Mouse	IHC-Fr
	VGLUT1 Recombinant Antibody [PSH10-53] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610284	Mouse, Rat	IHC-Fr,IHC-P,WB
	VGLUT1 Recombinant Mouse mAb [PSH08-34] - BSA and Azide free	HA610217	Mouse, Rat	WB,IHC-P
Relative expression	VGLUT1 Recombinant Rabbit mAb [PSH10-53] - BSA and Azide free	HA751356	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IHC-P,IF-Tissue
	v-SNARE Vti1a Rabbit pAb	ER65309	Human, Mouse, Rat	WB
	YWHAZ Rabbit pAb	HA500368	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

