

抑郁症研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	5HT1A Receptor Rabbit pAb	ER1903-07	Mouse, Rat	WB
	5HT2B Receptor Rabbit pAb	ER1903-10	Human	WB,IHC-P,IF-Tissue
	5HT2C Receptor Recombinant Rabbit mAb [JA31-31]	ET1704-87	Human, Mouse, Rat	WB,IHC-P
	ABAT / GABA-T Recombinant Rabbit mAb [JE50-16]	ET7110-01	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
KD	A-RAF Recombinant Rabbit mAb [JJ08-68]	ET1701-18	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	BDNF Recombinant Rabbit mAb [SJ12-09] - BSA and Azide free	HA750104	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC
	BDNF Recombinant Rabbit mAb [JE59-35]	HA722912	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-P,IHC-Fr,FC
	BDNF Rabbit pAb	ER130915	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	BDNF Recombinant Rabbit mAb [SJ12-09]	ET1606-42	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC
KD	BRAF Recombinant Rabbit mAb [SU34-04]	ET1608-36	Human, Mouse, Rat	WB,IHC-P,FC,IP
KD	BRAF Mouse mAb [F1-F4-E10]	EM1701-32	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	CAC1A Rabbit pAb	ER62087	Human, Mouse, Rat	IHC-P
	cGKI Rabbit pAb	ER1919-37	Human, Mouse, Rat	WB,IF-Cell
	cGKII Rabbit pAb	ER63000	Human, Mouse, Rat	WB,IHC-P
	Claudin 5 Recombinant Rabbit mAb [JM11-22] - BSA and Azide free	HA750384	Human	WB,IHC-P
	Claudin 5 Recombinant Rabbit mAb [JM11-22]	ET1703-58	Human	WB,IHC-P
	cPLA2 Rabbit pAb	ER63109	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	cPLA2-delta Rabbit pAb	ER63111	Human	IHC-P
	cPLA2-epsilon Rabbit pAb	ER63112	Human	IHC-P
	cPLA2-gamma Rabbit pAb	ER63110	Human	WB
	CRF Rabbit pAb	ER63124	Human	IHC-P
	CRHR1 Rabbit pAb	ER1906-68	Mouse, Rat	WB,IHC-P,IF-Tissue
KD	ERK1 Recombinant Rabbit mAb [SP00-74]	ET1604-32	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC

KD	ERK1 Recombinant Rabbit mAb [SP05-09]	ET1604-16	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	ERK1 Mouse mAb [7-D6-E5]	M1505-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Erk1/2 Mouse mAb [2G1]	RT1484	Human, Mouse, Rat	WB,IP,IF
	ERK1/2 Rabbit pAb	ER131011	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03] - BSA and Azide free	HA750019	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03]	ET1601-29	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
KD	ERK2 Recombinant Mouse mAb [C3B12-R]	HA601201	Human, Mouse, Rat	WB
KD	ERK2 Recombinant Mouse mAb [C3B12-R] - BSA and Azide free	HA610079	Human, Mouse, Rat	WB
KD	ERK2 Mouse mAb [B10-8]	EM1901-53	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Recombinant Rabbit mAb [SZ25-01]	ET1603-23	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
KD	ERK2 Mouse mAb [A10D2]	EM1901-54	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Rabbit pAb	ER131218	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	ERK2 Recombinant Mouse mAb [5-D2-R]	HA601202	Human, Mouse, Rat	WB
KD	ERK2 Recombinant Mouse mAb [5-D2-R] - BSA and Azide free	HA610080	Human, Mouse, Rat	WB
	ERK2 Mouse mAb [C3B12]	EM1901-55	Human, Mouse, Rat	WB,IHC-P,FC
	ERK2 Mouse mAb [5-D2]	M1407-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	GABA A Receptor alpha 3 Rabbit pAb	ER60110	Human, Rat, Mouse	WB,IHC-P
	GABA A Receptor gamma 2 Rabbit pAb	ER60112	Human, Rat, Mouse	WB,IHC-P
Relative expression	GABA B Receptor 1 Recombinant Rabbit mAb [PSH03-87] - BSA and Azide free	HA750912	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P
	GABA B Receptor 1 Recombinant Rabbit mAb [JE60-09]	HA721338	Human, Mouse	WB,IHC-P,IF-Tissue
Relative expression	GABA B Receptor 1 Recombinant Rabbit mAb [PSH06-21]	HA722653	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr
Relative expression	GABA B Receptor 1 Recombinant Rabbit mAb [PSH06-21] - BSA and Azide free	HA751059	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr
Relative expression	GABA B Receptor 1 Recombinant Rabbit mAb [PSH03-87]	HA722053	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	GABRA2 Rabbit pAb	HA500288	Human, Mouse, Rat	WB,FC
	GABRA2 Rabbit pAb	ER1909-25	Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	GABRA6 Rabbit pAb	HA500473	Human, Mouse, Rat	WB,IHC-P
	G alpha i-3 Rabbit pAb	ER63656	Human, Mouse, Rat	WB
	GCS-beta-1 Rabbit pAb	ER63524	Human, Mouse, Rat	WB,IHC-P
	Glucocorticoid Receptor Recombinant Rabbit mAb [JF0952]	ET1702-11	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC
	Glucocorticoid Receptor alpha Mouse mAb [3F2]	EM1701-66	Human, Mouse	WB,IHC-P,FC,IF-Cell
	Glucocorticoid Receptor alpha Rabbit pAb	R1701-1	Human, Mouse, Rat	WB,IF-Cell
	Glucocorticoid Receptor alpha Mouse mAb [3F1]	EM1701-65	Human	WB,IF-Cell,IHC-P,FC
Relative expression	GluR2 Recombinant Rabbit mAb [JU20-43]	ET1706-06	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	GluR-delta 2 Rabbit pAb	ER63548	Human, Mouse, Rat	WB,IF-Cell
	Glutamate receptor 1 Recombinant Rabbit mAb [SD2010]	ET1612-10	Human, Mouse, Rat	WB,IHC-P,IHC-Fr
	Glutamate receptor 1 Recombinant Rabbit mAb [SD2010] - BSA and Azide free	HA750281	Human, Mouse, Rat	WB,IHC-P,IHC-Fr
	GNA11 Rabbit pAb	ER60746	Human, Rat, Mouse	WB
	GNA12 Rabbit pAb	ER60741	Human, Mouse, Rat	WB
	GNA13 Recombinant Rabbit mAb [JE63-76]	HA721007	Human, Mouse, Rat	WB,IHC-P
	GNAI1 Rabbit pAb	ER62459	Human, Mouse, Rat	WB
KD	GNAI2 Rabbit pAb	HA500293	Human, Mouse, Rat	WB,IF-Cell
	GNAO1 Rabbit pAb	HA500136	Human, Mouse, Rat	WB,IHC-P
KD	GNAQ Recombinant Rabbit mAb [JE65-59]	HA721328	Human, Mouse, Rat	WB,IHC-P
	GNAS Mouse mAb [A8F12]	HA601074	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	GNAS Recombinant Mouse mAb [A8F11-R]	HA601275	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	GNAS Recombinant Mouse mAb [A8F11-R] - BSA and Azide free	HA610150	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	GNAS Mouse mAb [A8F11]	HA601075	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	GNAS3 Rabbit pAb	ER62067	Human	WB
	GRIA3 Mouse mAb [E2-D1]	EM1709-87	Human, Rat	WB,IHC-P
	GTPase HRAS Recombinant Rabbit mAb [JE42-15]	HA721355	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Guanylyl Cyclase alpha 1 Recombinant Rabbit mAb [JE34-21]	HA722537	Human	WB
	GUCY1A3 Mouse mAb [B2-G6]	EM1708-83	Human	WB,IF-Cell,IHC-P,FC
	GUCY1B3 Recombinant Rabbit mAb [JE41-09]	HA722314	Human, Mouse, Rat	WB
	Gz-alpha Rabbit pAb	ER63654	Human, Mouse, Rat, Monkey	WB,IHC-P
	Hck Rabbit pAb	ER63669	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA721713	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free (Capture)	HA721278	Human	ELISA(Cap)
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - Detector	HA723882B	Human	ELISA(Det),ELISA
	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free (Detector)	HA721279	Human	ELISA(Det)
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA722649B	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39]	HA721712	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - BSA and Azide free (Detector)	HA723882	Human	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-87] - BSA and Azide free (Capture)	HA723883	Human	ELISA(Cap)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free	HA750745	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-85] - BSA and Azide free (Capture)	HA723881	Human	ELISA(Cap)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free	HA750744	Human	WB
	Human IL-6 Mouse mAb [A6D2]	HA601050	Human	WB,ELISA
	Human IL-6 Mouse mAb [A6D1]	HA601051	Human	WB,ELISA
	HRP Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA721293	Human	ELISA(Det)
	Human IL-6 Recombinant Rabbit mAb [PS00-21] - BSA and Azide free (Capture)	HA721167	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-6 Recombinant Rabbit mAb [PS00-15] - BSA and Azide free (Detector)	HA721168	Human	ELISA(Det)
	Biotin Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA722233B	Human	ELISA(Det)
	Biotin Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721559	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb	IRS040	Human	mIHC
	HRP Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721294	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free	HA750593	Human	IHC-P
	Human TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free (Detector)	HA721266	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40]	HA721286	Human	IHC-P,mIHC
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free (Capture)	HA721265	Human	ELISA(Cap)
Cell treatment	IGF1 Recombinant Rabbit mAb [JE30-01]	HA722527	Human	WB
	IGF1 (E peptide) Rabbit pAb	HA500107	Human	WB,IHC-P,FC
	IGF1R Mouse mAb [D5-G2]	HA600033	Human	WB,FC
	IGF1R Rabbit pAb	R1310-6	Human	WB,IHC-P
Relative expression	IGF1 Receptor beta Recombinant Rabbit mAb [PSH11-72] - BSA and Azide free	HA751427	Human, Mouse, Rat	WB,IF-Cell
Relative expression	IGF1 Receptor beta Recombinant Rabbit mAb [PSH11-72]	HA723368	Human, Mouse, Rat	WB,IF-Cell
	IGF I Rabbit pAb	ER1911-10	Human, Mouse, Rat, Dog	WB,IHC-P,IF-Tissue
	IGF-IR Rabbit pAb	ER63734	Human, Mouse, Rat	IF-Cell,WB,IHC-P
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R] - BSA and Azide free	HA610004	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Recombinant Rabbit mAb [JJ087-3]	ET1701-39	Human	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R]	HA601036	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Mouse mAb [A7A10]	HA601035	Human	WB,ELISA,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	IL-1 beta Mouse mAb [A7F9]	HA601002	Human, Mouse, Rat	WB,ELISA,IF-Cell
	IL-6 Mouse mAb	EM30301	Human	WB
	IL-6 Recombinant Rabbit mAb	IRS160RB	Human	mIHC
	IL-6 Rabbit pAb	R1412-2	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	IL-6 Mouse mAb [3-A4]	M1303-3	Mouse	WB,IF-Cell,FC,ELISA
	IL-6 Mouse mAb [1-6]	EM1701-45	Mouse, Rat	WB,IF-Cell
	IL-6 Rabbit pAb	R1203-2	Mouse	WB,ELISA
KD, Relative expression	Insulin Receptor Beta Recombinant Rabbit mAb [SN20-13] - BSA and Azide free	HA750276	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC
KD, Relative expression	Insulin Receptor Beta Recombinant Rabbit mAb [SN20-13]	ET1611-90	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC
	Ionotropic glutamate receptor 2 + 3 Recombinant Rabbit mAb [JE33-28]	HA721405	Mouse, Rat, Human	WB
Relative expression	IP3 receptor Recombinant Mouse mAb [PSH14-18]	HA601456	Mouse, Rat	WB,IHC-P
Relative expression	IP3 receptor Recombinant Mouse mAb [PSH14-18] - BSA and Azide free	HA610290	Mouse, Rat	WB,IHC-P
	IP3 Receptor Recombinant Rabbit mAb [JA11-35]	ET1704-77	Human, Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Tissue
	IP3 Receptor Rabbit pAb	ER1803-17	Human, Mouse, Rat	WB,IHC-P
	ITPR2 Rabbit pAb	ER1803-65	Human, Mouse, Rat	Dot Blot,IF-Cell,IHC-P,FC
	KCNK2 Rabbit pAb	ER1902-06	Human, Mouse, Rat	WB,IHC-P,FC
	KCNQ1 Mouse mAb [D3-E12]	EM1709-39	Human	WB,FC
	KCNQ2/3/4/5 Rabbit pAb	ER63796	Human, Mouse, Rat	IHC-P
	KCNQ4 Rabbit pAb	ER2001-23	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	KCNQ5 Rabbit pAb	ER63798	Human, Mouse	WB
	K-Ras Rabbit pAb	ER65706	Human, Mouse, Rat	WB,IHC-P
	KRAS Mouse mAb [A8E5]	HA601059	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	KRAS Recombinant Mouse mAb [A8E6]	HA601060	Human, Mouse, Rat	WB,IHC-P
	KRAS Recombinant Mouse mAb [A8E6] - BSA and Azide free	HA610008	Human, Mouse, Rat	WB,IHC-P
	KRAS Mouse mAb [A8E4]	HA601058	Human, Mouse, Rat	WB,IF-Cell,FC
Relative expression	Lck Recombinant Rabbit mAb [PSH09-25] - BSA and Azide free	HA751279	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	Lck Rabbit pAb	HA500280	Human, Rat	WB,FC
Relative expression	Lck Recombinant Rabbit mAb [PSH09-25]	HA723083	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Lck Rabbit pAb	ER63849	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Lyn Rabbit pAb	ER1912-57	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Lyn Rabbit pAb	HA500142	Human, Mouse, Rat	WB
KD	Lyn Recombinant Rabbit mAb [JE44-75]	HA721991	Human, Mouse, Rat	WB,FC
KD	MEK1 Recombinant Rabbit mAb [SZ22-01]	ET1603-20	Human, Mouse, Rat, Cow, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK1/2 Recombinant Rabbit mAb [SR13-07] - BSA and Azide free	HA750042	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK1/2 Recombinant Rabbit mAb [SR13-07]	ET1602-3	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK2 Recombinant Rabbit mAb [SD207-08]	ET1612-6	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	mGluR1 Recombinant Rabbit mAb [JM11-61]	ET1703-84	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-60] - BSA and Azide free (Detector)	HA721413	Mouse	ELISA(Det)
	HRP Conjugated Mouse IL-6 Mouse mAb [3-A4]	M1303-2	Mouse	WB,ELISA
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PS01-83] - BSA and Azide free (Detector)	HA721602	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PS01-74] - BSA and Azide free (Capture/Detector)	HA721600	Mouse	ELISA(Cap),ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PSH0-59] - Detector	HA722646B	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-59] - BSA and Azide free (Detector)	HA721412	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-58] - BSA and Azide free (Capture)	HA721411	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PSH0-60] - Detector	HA722647B	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PS01-83] - BSA and Azide free (Detector)	HA721601	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PS01-74] - Detector	HA723643B	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IL-6 Recombinant Rabbit mAb [PSH14-47] - BSA and Azide free (Capture)	HA723642	Mouse	ELISA(Cap)
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30] - BSA and Azide free (Detector)	HA721975	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30]	HA722626B	Mouse, Rat	ELISA(Det)
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-29] - BSA and Azide free (Capture)	HA721974	Mouse, Rat	ELISA(Cap)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77] - BSA and Azide free (Detector)	HA722041	Mouse	ELISA(Det)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-76] - BSA and Azide free (Capture)	HA722040	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77]	HA722623B	Mouse	ELISA(Det)
	N/H/K-Ras Rabbit pAb	ER64095	Human, Mouse, Rat	WB,IHC-P
Relative expression	NMDAR1 Recombinant Rabbit mAb [JM11-26] - BSA and Azide free	HA750391	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,IHC-Fr
Relative expression	NMDAR1 Recombinant Rabbit mAb [JM11-26]	ET1703-75	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,IHC-Fr
	NMDAR2A Recombinant Rabbit mAb [JA31-20]	ET1704-80	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue
Relative expression	NMDAR2B Recombinant Rabbit mAb [PSH05-46]	HA722284	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
	NMDAR2B Rabbit pAb	HA500208	Mouse, Rat, Human	WB,IF-Cell,IHC-P,FC
Relative expression	NMDAR2B Recombinant Rabbit mAb [PSH05-46] - BSA and Azide free	HA751005	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
	NMDAR2B Rabbit pAb	ER1913-90	Human, Rat	WB,IHC-P,IF-Tissue
	nNOS Recombinant Rabbit mAb [ST520]	ET1609-61	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,IP
	nNOS Rabbit pAb	R1510-28	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IHC-Fr
Relative expression	PDE4B Recombinant Rabbit mAb [JE33-95]	HA722323	Human, Mouse	WB,IHC-P,FC
	PDE4B/C/D Rabbit pAb	ER64537	Human, Mouse, Rat	WB,IHC-P
	PDE4C Rabbit pAb	ER64538	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PDE4D Mouse mAb [A7F3]	HA601001	Human, Mouse	WB
	PDE4D Rabbit pAb	ER1914-89	Mouse, Rat	WB,FC
Cell treatment	Phospho-BRAF (T401) Recombinant Rabbit mAb [JJ08-72]	ET1701-20	Human, Mouse, Rat	WB,IF-Cell,IP,IHC-P
Cell treatment	Phospho-Cytosolic Phospholipase A2 (S505) Recombinant Rabbit mAb [PSH07-56]	HA722849	Human, Mouse	WB,Dot Blot
Cell treatment	Phospho-Cytosolic Phospholipase A2 (S505) Recombinant Rabbit mAb [PSH07-16] - BSA and Azide free	HA751137	Human, Mouse	WB
Cell treatment	Phospho-Cytosolic Phospholipase A2 (S505) Recombinant Rabbit mAb [PSH07-16]	HA722804	Human, Mouse	WB
Cell treatment	Phospho-Cytosolic Phospholipase A2 (S505) Recombinant Rabbit mAb [PSH07-56] - BSA and Azide free	HA751176	Human, Mouse	WB,Dot Blot
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4] - BSA and Azide free	HA750061	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4]	ET1603-22	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01] - BSA and Azide free	HA750199	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01]	ET1610-13	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
	Phospho-GluR1 (AMPA subtype) (S831) Recombinant Rabbit mAb [JE44-96]	HA721865	Human, Mouse, Rat	WB
Cell treatment	Phospho-GluR1 (S845) Recombinant Rabbit mAb [JJ2009]	ET1701-28	Mouse, Rat, Human	WB,FC
	Phospholipase C beta 1 Rabbit pAb	HA500302	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Phospholipase C beta 2 Rabbit pAb	ER1915-09	Human, Mouse	WB,IHC-P,IF-Tissue,FC
	Phospholipase C beta 3/PLCB3 Recombinant Rabbit mAb [JE39-14]	HA721561	Human	WB,IHC-P
Cell treatment, Relative expression	Phospho-LYN (Y397)/LCK (Y394)/HCK (Y411)/BLK (Y389) Recombinant Rabbit mAb [PSH09-43]	HA723105	Human, Mouse, Rat	WB,IF-Cell
Cell treatment, Relative expression	Phospho-LYN (Y397)/LCK (Y394)/HCK (Y411)/BLK (Y389) Recombinant Rabbit mAb [PSH09-43] - BSA and Azide free	HA751290	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-MEK1/2 (S218 + S222) Recombinant Rabbit mAb [ST0490]	ET1609-50	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Phospho-MEK1 (S298) Recombinant Rabbit mAb [SD206-7]	ET1612-40	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-MEK1 (T292) Recombinant Rabbit mAb [SD2088]	ET1612-42	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH15-91]	HA723796	Human, Mouse, Rat	WB,IF-Cell,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH15-91] - BSA and Azide free	HA751582	Human, Mouse, Rat	WB,IF-Cell,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH16-27] - BSA and Azide free	HA751595	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH16-27]	HA723825	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,Dot Blot
Cell treatment	Phospho-PKC alpha (T497) Recombinant Rabbit mAb [JE44-81]	HA722278	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-PKC alpha (T638) Recombinant Rabbit mAb [JF0964]	ET1702-17	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,Dot Blot
Cell treatment	Phospho-PKC delta (Y311) Recombinant Rabbit mAb [JE44-20]	HA722458	Human, Mouse, Rat	WB
Cell treatment	Phospho-PKC (S729) Recombinant Rabbit mAb [JE77-54]	HA722846	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-PKC zeta (T560) Recombinant Rabbit mAb [SN0753]	ET1611-78	Human, Mouse, Rat	WB,IHC-P
	Phospho-PP2A (Y307) Recombinant Rabbit mAb [ST49-05]	ET1609-40	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-Raf1 (S259) Recombinant Rabbit mAb [SD85-07]	ET1612-87	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-Raf1 (S259) Recombinant Rabbit mAb [SD85-07] - BSA and Azide free	HA750304	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Phospho-Raf1 (S43) Recombinant Rabbit mAb [SD087-05]	ET1612-75	Human, Mouse, Rat	WB
	Phospho-Raf1 (S621) Recombinant Rabbit mAb [JJ085-05]	ET1701-3	Human, Rat, Mouse	WB,IHC-P
	PKC Rabbit pAb	ER64620	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
KD	PKC alpha Recombinant Rabbit mAb [SU31-08]	ET1608-15	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	PKC beta Rabbit pAb	ER1802-14	Human, Mouse	WB,IF-Cell,IHC-P
	PKC beta 1 Recombinant Rabbit mAb [JM43-26]	ET1705-30	Human, Mouse, Rat	WB,IHC-P,FC
	PKC beta 2 Recombinant Rabbit mAb [ST48-06]	ET1609-44	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	PKC delta Recombinant Rabbit mAb [JJ093-02]	ET1701-85	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PKC delta Recombinant Rabbit mAb [JJ093-02] - BSA and Azide free	HA750325	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	PKC epsilon Recombinant Rabbit mAb [JA38-21]	ET1704-85	Human, Mouse, Rat	WB,IHC-P
	PKC eta Recombinant Rabbit mAb [JE32-41]	HA722389	Human, Mouse, Rat, Monkey	WB
	PKC gamma Recombinant Rabbit mAb [JG69-32]	ET7108-42	Human, Mouse, Rat	WB,IHC-P
	PKC theta Rabbit pAb	ER64621	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Relative expression	PKC zeta Recombinant Rabbit mAb [PSH06-74]	HA722773	Human	WB
Relative expression	PKC zeta Recombinant Rabbit mAb [PSH06-74] - BSA and Azide free	HA751122	Human	WB
	PLCB2 Rabbit pAb	HA500399	Human	WB,IF-Cell
	PLCB3 Rabbit pAb	HA500391	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	PLC beta 3 Rabbit pAb	ER64629	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	PP2A alpha + beta Recombinant Rabbit mAb [SN70-08]	ET1611-54	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P
	PP2A alpha + beta Mouse mAb [A0-G7]	M1603-4	Human, Mouse, Rat	WB,IF-Cell
	PP2A(alpha+beta) Rabbit pAb	R1510-31	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	PPP2R1A Rabbit pAb	HA500160	Human, Mouse, Rat	WB,IHC-P
KD	PPP2R1B Recombinant Rabbit mAb [PSH02-92] - BSA and Azide free	HA750847	Human, Mouse, Rat	WB,IHC-P
KD	PPP2R1B Recombinant Rabbit mAb [PSH02-92]	HA721911	Human, Mouse, Rat	WB,IHC-P
	PPP2R1B Rabbit pAb	HA500383	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KD	Raf1 Recombinant Rabbit mAb [JJ086-06]	ET1701-21	Human, Mouse, Rat	WB
	Ras Recombinant Rabbit mAb [JF10-11]	ET1702-94	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IP,FC
	Ras Recombinant Rabbit mAb [JE44-13]	HA721883	Human, Mouse, Rat	WB,IF-Cell,FC
	Ras Recombinant Rabbit mAb [SA39-05]	ET1601-16	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC,IP
	Ras Rabbit pAb	ER40115	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	Reelin Recombinant Rabbit mAb [JA63-13]	ET1704-98	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
	RYR1 Rabbit pAb	ER1803-19	Human, Mouse	Dot Blot,IHC-P,FC
	SR-1A Rabbit pAb	ER65015	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	TNFAIP3 Rabbit pAb	ER1706-32	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	TNFAIP3 Recombinant Rabbit mAb [SN07-31] - BSA and Azide free	HA750253	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	TNFAIP3 Recombinant Rabbit mAb [SN07-31]	ET1611-40	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07]	HA722022	Human, Mouse	WB
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free	HA750897	Human, Mouse	WB
	TNF alpha Rabbit pAb	ER1919-22	Human, Mouse	WB
	TNF alpha Rabbit pAb	ER65189	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	TNF alpha Rabbit pAb	R1203-1	Mouse	WB,IF-Cell,ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层