

胆汁酸代谢研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	17 beta-HSD4 Rabbit pAb	ER65564	Human, Mouse, Rat	WB,IHC-P
	3 beta-HSD7 Rabbit pAb	ER62529	Human, Mouse, Rat	WB
	ABCC4 Mouse mAb [H3-A7]	EM1708-59	Human	IHC-P,FC
	ABCD3 Rabbit pAb	ER60507	Human, Mouse, Rat	WB
KD	ABCG2 Recombinant Rabbit mAb [JF0994]	ET1702-40	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,FC
	ABCG5 Mouse mAb [E10-B5]	EM1708-58	Human	WB,IHC-P,FC
	ABCG8 Rabbit pAb	ER1903-17	Human	WB,IHC-P,IF-Tissue
	ACOT8 Rabbit pAb	ER62580	Human	WB,IHC-P,IF-Cell
	ACOX2 Rabbit pAb	ER60439	Human, Mouse, Rat	WB
	A Cyclase IV Rabbit pAb	ER62534	Human, Mouse, Rat	WB
	A Cyclase IX Rabbit pAb	ER62535	Human	WB,IHC-P
	A Cyclase VIII Rabbit pAb	ER62537	Human, Mouse, Rat	IHC-P,IF-Cell
	A Cyclase V/VI Rabbit pAb	ER62536	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	ADCY3 Rabbit pAb	ER1903-37	Mouse	WB
	ADCY5 Rabbit pAb	ER1903-38	Mouse, Rat	WB
	ADCY7 Rabbit pAb	ER1903-39	Human, Rat	WB
	Adenylate cyclase 1 Rabbit pAb	ER1903-42	Human, Rat	WB,IHC-P,IF-Tissue
	AK1C4 Rabbit pAb	ER60464	Human	WB
	AMACR Rabbit pAb	ER1706-60	Human, Mouse	WB,IF-Cell,IHC-P,FC
Relative expression	AMACR Recombinant Rabbit mAb [PSH13-54]	HA723543	Human, Mouse	WB,IHC-P,IF-Tissue
Relative expression	AMACR Recombinant Rabbit mAb [PSH13-55] - BSA and Azide free	HA751474	Human, Mouse, Rat	WB,IHC-P
Relative expression	AMACR Recombinant Rabbit mAb [PSH13-55]	HA723544	Human, Mouse, Rat	WB,IHC-P
	AMACR Mouse mAb [6F5]	EM1701-80	Human, Mouse, Rat	WB
	AQP1 Recombinant Rabbit mAb	IRS088RB	Human	mIHC
	AQP1 Recombinant Rabbit mAb [JM10-98] - BSA and Azide free	HA750378	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P

	AQP1 Recombinant Rabbit mAb [JM10-98]	ET1703-34	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,mIHC
	AQP1 Recombinant Rabbit mAb	IRS218RB	Mouse	mIHC
	AQP8 Rabbit pAb	ER60486	Human, Mouse, Rat	WB
	AQP9 Rabbit pAb	ER1903-85	Human	FC
Relative expression	Aquaporin 4 Recombinant Antibody [PSH06-38] - Guinea pig IgG2 (Chimeric)	HA601461	Mouse	IHC-Fr,WB
Relative expression	Aquaporin 4 Recombinant Antibody [PSH06-38] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610296	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	Aquaporin 4 Rabbit pAb	ER1903-87	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
Relative expression	Aquaporin 4 Recombinant Rabbit mAb [PSH06-38] - BSA and Azide free	HA751077	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr,IP
	Aquaporin 4 Recombinant Antibody [PSH06-38] - Chicken IgY (Chimeric) - BSA and Azide free	HA610294	Mouse	IHC-Fr
Relative expression	Aquaporin 4 Recombinant Antibody [PSH06-38] - Rat IgG1 (Chimeric)	HA601462	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Aquaporin 4 Recombinant Antibody [PSH06-38] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610295	Mouse	IHC-Fr,WB
	Aquaporin 4 Recombinant Antibody [PSH06-38] - Chicken IgY (Chimeric)	HA601460	Mouse	IHC-Fr
	Aquaporin 4 Rabbit pAb	HA500125	Human, Mouse, Rat	WB,IHC-P
Relative expression	Aquaporin 4 Recombinant Rabbit mAb [JE59-76]	HA721954	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	Aquaporin 4 Recombinant Rabbit mAb [PSH06-38]	HA722672	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr,IP
	ATP1A2 Recombinant Rabbit mAb [JE34-00]	HA721422	Human, Mouse, Rat	WB,IHC-P,FC
	ATP1A3 Recombinant Rabbit mAb [JE35-24]	HA721596	Human, Mouse, Rat	WB
Relative expression	ATP1B4 Rabbit pAb	HA500311	Human, Mouse, Rat	WB
	ATPase Na ⁺ / K ⁺ beta 2 Rabbit pAb	ER1904-02	Mouse, Rat	WB
	BAAT Rabbit pAb	ER60528	Human	WB
	beta 1 Sodium Potassium ATPase Rabbit pAb	ER60078	Mouse, Rat	WB
	Bonzo Rabbit pAb	ER65568	Human	WB,IHC-P,IF-Cell
	Carbonic anhydrase 2 Rabbit pAb	ER1802-72	Human	WB,IF-Cell,IHC-P
	Carbonic anhydrase 2 Mouse mAb [11A1]	EM1801-08	Human, Mouse, Rat	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Carbonic anhydrase 2 Mouse mAb [11A2]	EM1801-14	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	Carbonic anhydrase 2 Recombinant Rabbit mAb [JU12-37]	ET1706-47	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
	Carbonic anhydrase 2 Recombinant Mouse mAb [11A2-R] - BSA and Azide free	HA610116	Human	WB,IF-Cell,IHC-P
Relative expression	Carbonic anhydrase 2 Recombinant Rabbit mAb [JU12-37] - BSA and Azide free	HA750456	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
	Carbonic anhydrase 2 Recombinant Mouse mAb [11A2-R]	HA601241	Human	WB,IF-Cell,IHC-P
	CD182 Mouse mAb [5F7A11]	EM1711-55	Human	WB,IHC-P,IF-Cell,FC
	CD183 Mouse mAb [5C10E6]	EM1710-60	Human	WB,IHC-P,FC
	CD298 Rabbit pAb	ER62941	Human	WB
	CH25H Rabbit pAb	ER1919-29	Human, Rat	WB,IHC-P,IF-Tissue
	CH25H Rabbit pAb	ER1919-26	Human	IHC-P,IF-Tissue
	CXCR1 Rabbit pAb	HA500526	Human	WB,IF-Cell
	CXCR1 Rabbit pAb	ER1906-86	Human	WB
	CXCR2 Rabbit pAb	ER1906-87	Human	WB
	CXCR3 Recombinant Rabbit mAb [JA61-33] - BSA and Azide free	HA750422	Human	WB
	CXCR3 Rabbit pAb	ER1906-88	Human, Mouse, Rat	WB
	CXCR3 Recombinant Rabbit mAb [JA61-33]	ET1704-97	Human	WB
	CXCR4 Rabbit pAb	ER1802-28	Human, Rat	WB,IF-Cell,FC
	CXCR4 Recombinant Rabbit mAb [JE03-35]	HA722304	Human	WB,IF-Cell,IHC-P
	CXCR5 Recombinant Rabbit mAb [JB11-40]	ET7107-74	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	CXCR5 Recombinant Rabbit mAb [JB11-40] - BSA and Azide free	HA750463	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	CXCR6 Recombinant Rabbit mAb [JE66-40]	HA721542	Human	WB
	CYP27A1 Recombinant Rabbit mAb [JG40-36]	ET7109-05	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	CYP39A1 Rabbit pAb	ER63176	Human, Mouse	WB,IHC-P,IF-Cell
	CYP3A4 Recombinant Rabbit mAb [JA11-33]	ET1704-88	Human	WB,IHC-P,IP
	CYP3A4/5 Rabbit pAb	ER63177	Human	WB,IHC-P,IF-Cell
	CYP7A1 Rabbit pAb	HA500207	Human, Rat	WB,IF-Cell,FC
	CYP7B1 Rabbit pAb	ER63187	Human	WB,IHC-P
	CYP8B1 Rabbit pAb	ER1907-09	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	EDG-5 Rabbit pAb	ER63315	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	FABP6 Rabbit pAb	ER62027	Human, Rat, Mouse	WB
	FATP2 Mouse mAb [A9-B3]	EM1708-57	Human	WB,IHC-P
	FGF-19 Rabbit pAb	ER65959	Human, Pig	WB
	FGFR-4 Rabbit pAb	ER65688	Human	WB
	FGFR4 Recombinant Rabbit mAb [JE59-96]	HA721578	Human	WB,FC
	FXR1 Rabbit pAb	ER1909-14	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	FXR1 Recombinant Rabbit mAb [JE40-57]	HA721499	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	FXR1 Recombinant Rabbit mAb [PSH01-67] - BSA and Azide free	HA750743	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC,IP
	FXR1 Recombinant Rabbit mAb [PSH01-67]	HA721710	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC,IP
	FXR1 Recombinant Rabbit mAb [PSH01-35] - BSA and Azide free	HA750713	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,IP
	FXR1 Recombinant Rabbit mAb [PSH01-35]	HA721669	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,IP
	FXR2 Rabbit pAb	ER63480	Human, Mouse	WB,IHC-P,IF-Cell
Relative expression	GLP-1 Recombinant Rabbit mAb [PSH19-16]	HA724050	Human, Mouse, Rat	IHC-P,IF-Tissue
	Glucagon Recombinant Rabbit mAb [JF0960]	ET1702-20	Human, Mouse, Rat	WB,IHC-P,mIHC,IHC-Fr,IF-Tissue
	Glucagon Recombinant Rabbit mAb	IRS150RB	Human	mIHC
	Glucagon Recombinant Rabbit mAb [JF33-10]	ET1702-59	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Glucagon Mouse mAb	EM1701-34	Human, Mouse, Rat	IHC,IF-Cell
	Glucagon Recombinant Rabbit mAb	IRS230RB	Mouse	mIHC
	Glucose Transporter GLUT1 Recombinant Antibody [SA0377] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610348	Human, Mouse, Rat	IHC-Fr,IHC-P,IF-Tissue,WB
	iFluor™ 647 Conjugated Glucose Transporter GLUT1 Recombinant Rabbit mAb [SA0377]	HA720170F	Human, Mouse	IF-Tissue
	Glucose Transporter GLUT1 Recombinant Rabbit mAb	IRS151RB	Human	mIHC
KD	Glucose Transporter GLUT1 Recombinant Rabbit mAb [SA0377]	ET1601-10	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	iFluor™ 488 Conjugated Glucose Transporter GLUT1 Recombinant Rabbit mAb [SA0377]	HA720154F	Human, Mouse, Rat	IF-Cell,IF-Tissue
	Glucose Transporter GLUT1 Recombinant Antibody [SA0377] - Rat IgG1 (Chimeric)	HA601521	Human, Mouse, Rat	IHC-Fr,IHC-P,IF-Tissue,WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	Glucose Transporter GLUT1 Recombinant Rabbit mAb [SA0377] - BSA and Azide free	HA750002	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Glucose Transporter GLUT1 Rabbit pAb	ER1510-11	Human, Mouse	WB,IF-Cell,IHC-P,FC
	GNAS Recombinant Mouse mAb [A8F11-R] - BSA and Azide free	HA610150	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	GNAS Mouse mAb [A8F11]	HA601075	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	GNAS Recombinant Mouse mAb [A8F11-R]	HA601275	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	GNAS Mouse mAb [A8F12]	HA601074	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	GNAS3 Rabbit pAb	ER62067	Human	WB
	HMGCR Recombinant Rabbit mAb [JF0981]	ET1702-41	Human, Mouse, Rat	WB,IP,IHC-P
	HMGCR Recombinant Rabbit mAb [JF0981] - BSA and Azide free	HA750345	Human, Mouse, Rat	WB,IP,IHC-P
	Human FGF-19 Recombinant Rabbit mAb [PSH13-29] - BSA and Azide free (Capture)	HA725100	Human	ELISA(Cap)
	Human FGF-19 Recombinant Rabbit mAb [PSH13-30] - BSA and Azide free (Detector)	HA725101	Human	ELISA(Det)
	Biotin Conjugated Human FGF-19 Recombinant Rabbit mAb [PSH13-30] - Detector	HA725102B	Human	ELISA(Det),ELISA
	Human IL-10 Recombinant Rabbit mAb [PS00-98] - BSA and Azide free (Capture)	HA721236	Human	ELISA(Cap)
	Biotin Conjugated Human IL-10 Recombinant Rabbit mAb [PS00-81]	HA722652B	Human	ELISA(Det)
	Human IL-10 Recombinant Rabbit mAb [PS00-81] - BSA and Azide free (Detector)	HA721237	Human	ELISA(Det)
	Biotin Conjugated Human IL-17A Recombinant Rabbit mAb [PS01-24] - Detector	HA721631	Human	ELISA(Det)
	Human IL-17A Recombinant Rabbit mAb [PS01-20] - BSA and Azide free (Capture)	HA721232	Human	ELISA(Cap)
	Human IL-17A Recombinant Rabbit mAb [PS01-24] - BSA and Azide free (Detector)	HA721233	Human	ELISA(Det)
	HRP Conjugated Human IL-17A Recombinant Rabbit mAb [PS01-24]	HA721296	Human	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free (Capture)	HA721278	Human	ELISA(Cap)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA721713	Human	WB
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39]	HA721712	Human	WB
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free	HA750744	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-1 beta Recombinant Rabbit mAb [PSH16-85] - BSA and Azide free (Capture)	HA723881	Human	ELISA(Cap)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-87] - BSA and Azide free (Capture)	HA723883	Human	ELISA(Cap)
	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free (Detector)	HA721279	Human	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - BSA and Azide free (Detector)	HA723882	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free	HA750745	Human	WB
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA722649B	Human	ELISA(Det)
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - Detector	HA723882B	Human	ELISA(Det),ELISA
	Human IL-2 Mouse mAb [A8B7] - BSA and Azide free (Detector)	HA601006	Human	WB,ELISA(Det)
	Human IL-2 Recombinant Rabbit mAb [PSH01-42] - BSA and Azide free (Capture)	HA721680	Human	ELISA(Cap)
	Biotin Conjugated Human IL-2 Mouse mAb [A8B7]	HA601361B	Human	ELISA(Det)
	Human IL-23A Recombinant Rabbit mAb [PSH03-03] - BSA and Azide free (Capture)	HA721933	Human	ELISA(Cap)
	Biotin Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721559	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free	HA750593	Human	IHC-P
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free (Capture)	HA721265	Human	ELISA(Cap)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40]	HA721286	Human	IHC-P,mIHC
	Human TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free (Detector)	HA721266	Human	ELISA(Det)
	HRP Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721294	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb	IRS040	Human	mIHC
	IL-10 Rabbit pAb	ER65464	Human, Mouse, Rat	WB,IHC-P
	IL-10 Rabbit pAb	ER1911-19	Human, Mouse	WB,IHC-P,IF-Tissue
Cell treatment	IL-10 Recombinant Rabbit mAb [JE59-71]	HA722032	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	IL-17 Rabbit pAb	ER1706-91	Human	WB,IHC-P,FC
	IL17 Mouse mAb [4G2]	RT1326	Human, Mouse, Rat	WB,IP,IF

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	IL-17A Rabbit pAb	IRS056	Mouse	mIHC
	IL-17A Rabbit pAb	ER1902-37	Human, Mouse, Rat	WB,IHC-P,FC,mIHC
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R]	HA601036	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Mouse mAb [A7F9]	HA601002	Human, Mouse, Rat	WB,ELISA,IF-Cell
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R] - BSA and Azide free	HA610004	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Recombinant Rabbit mAb [JJ087-3]	ET1701-39	Human	WB,IF-Cell,IF-Tissue,IHC-P
	IL-1 beta Recombinant Rabbit mAb [PSH18-00]	HA723965	Mouse, Rat	WB
Cell treatment	IL-1 beta Mouse mAb [A7A10]	HA601035	Human	WB,ELISA,IF-Cell
	IL-2 Mouse mAb [A8B8]	HA601015	Human	WB,ELISA
	IL-2 Mouse mAb [A8B9]	HA601016	Human	WB,ELISA
	IL-2 Rabbit pAb	ER63748	Human, Rat	WB,IHC-P
	IL-23 Recombinant Rabbit mAb [JE38-99]	HA721037	Human	WB,IHC-P,IF-Tissue
	KCNN2 Rabbit pAb	ER1902-41	Human, Mouse, Rat	WB,IHC-P,FC
	LDL Receptor Rabbit pAb	R1512-12	Human, Mouse	WB,IF-Cell,FC
	LDL Receptor Recombinant Rabbit mAb [SJ0197]	ET1606-47	Human, Rat	WB,IF-Cell,FC
	LDL Receptor Recombinant Rabbit mAb [SJ0197] - BSA and Azide free	HA750107	Human, Rat	WB,IF-Cell,FC
	LXR alpha Recombinant Rabbit mAb [JA20-38]	ET1704-51	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	LXR alpha Recombinant Rabbit mAb [JA20-38] - BSA and Azide free	HA750414	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	MDR1 Rabbit pAb	ER1912-91	Human, Mouse, Rat	WB,IF-Cell,FC
	MDR1 Mouse mAb [1G6]	RT1383	Human, Mouse, Rat	WB,IP,IF
	MDR3 Rabbit pAb	ER60854	Human, Rat, Mouse	WB
	Mouse IL-10 Recombinant Rabbit mAb [PS01-93] - BSA and Azide free (Detector)	HA721518	Mouse	ELISA(Det)
	Mouse IL-10 Recombinant Rabbit mAb [PS01-85] - BSA and Azide free (Capture)	HA721517	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-10 Recombinant Rabbit mAb [PS01-93] - Detector	HA721519	Mouse	ELISA(Det)
	Mouse IL-17A Recombinant Rabbit mAb [PSH05-36] - BSA and Azide free (Detector)	HA722265	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse IL-17A Recombinant Rabbit mAb [PSH05-36]	HA722600B	Mouse, Rat	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IL-17A Recombinant Rabbit mAb [PS01-92] - BSA and Azide free (Detector)	HA721607	Mouse	ELISA(Det)
	Mouse IL-17A Recombinant Rabbit mAb [PSH05-35] - BSA and Azide free (Capture)	HA722264	Mouse, Rat	ELISA(Cap)
	Mouse IL-17A Recombinant Rabbit mAb [PS01-91] - BSA and Azide free (Capture)	HA721606	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-17A Recombinant Rabbit mAb [PS01-92] - Detector	HA721608	Mouse	ELISA(Det),ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - BSA and Azide free (Detector)	HA721977	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67] - BSA and Azide free (Detector)	HA722167	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - Detector	HA722625B	Mouse	ELISA(Det),ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH09-58] - BSA and Azide free	HA723120	Mouse	ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-31] - BSA and Azide free (Capture)	HA721976	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67]	HA722608B	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH09-59] - BSA and Azide free	HA723121	Mouse	ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-66] - BSA and Azide free (Capture)	HA722166	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30]	HA722626B	Mouse, Rat	ELISA(Det)
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-29] - BSA and Azide free (Capture)	HA721974	Mouse, Rat	ELISA(Cap)
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30] - BSA and Azide free (Detector)	HA721975	Mouse, Rat	ELISA(Det)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77] - BSA and Azide free (Detector)	HA722041	Mouse	ELISA(Det)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-76] - BSA and Azide free (Capture)	HA722040	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77]	HA722623B	Mouse	ELISA(Det)
	MRGE Rabbit pAb	ER64001	Human, Mouse, Rat	WB
Relative expression, KD	MRP2 Recombinant Rabbit mAb [JA32-01]	ET1704-47	Human	WB,IF-Cell,IHC-P,FC
	MRP3 Rabbit pAb	ER64009	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Na ⁺ /K ⁺ ATPase Mouse mAb [1G1]	RT1412	Human, Mouse, Rat	WB,IP,IF,IHC-P
	Na ⁺ /K ⁺ -ATPase alpha 1 Rabbit pAb	ER64103	Rat	WB,IF-Cell
	Na ⁺ /K ⁺ -ATPase alpha 2 Rabbit pAb	ER64104	Human, Mouse, Rat, Monkey	WB
KD	NF-kB p105 / p50 Recombinant Rabbit mAb [SZ20-01]	ET1603-18	Human, Mouse, Rat	WB,IHC-P,ICC,FC,IF-Tissue
KD	NF-kB p105 / p50 Recombinant Rabbit mAb [SZ20-01] - BSA and Azide free	HA750060	Human, Mouse, Rat	WB,IHC-P,ICC,FC,IF-Tissue
	NF-kB p105 / p50 Rabbit pAb	R1309-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	NF-kB p105 / p50 Mouse mAb [A2-7]	M1505-8	Human, Mouse	WB,IF-Cell,IHC-P
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04] - BSA and Azide free	HA750058	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04]	ET1603-12	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27]	HA721307	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	NF-kB p65 Rabbit pAb	HA500402	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	NF-kB p65 Rabbit pAb	ER0815	Human, Mouse, Rat, Zebrafish	WB,IHC-P,FC
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27] - BSA and Azide free	HA750601	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	NHE-1 Rabbit pAb	ER1902-51	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	NR0B2 Rabbit pAb	ER60945	Human, Rat, Mouse	WB
	NR1H4 Rabbit pAb	ER1914-12	Human, Mouse	WB,FC
	NR5A2 Recombinant Rabbit mAb [JE30-50]	HA721379	Human, Mouse, Rat	WB
	OATP1 Rabbit pAb	ER64249	Human	WB,IHC-P
	OATP-C Rabbit pAb	ER64250	Human	WB
	P Glycoprotein Recombinant Mouse mAb [A10F2-R]	HA601314	Human, Mouse, Rat	IHC-P,IF-Tissue
	P Glycoprotein Rabbit pAb	HA500502	Human, Mouse, Rat	WB,IHC-P,FC
KD	P Glycoprotein Recombinant Rabbit mAb [SN06-42] - BSA and Azide free	HA750249	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell
	P Glycoprotein Recombinant Mouse mAb [A10E12-R]	HA601312	Human, Mouse, Rat	IHC-P,IF-Tissue
	P Glycoprotein Recombinant Mouse mAb [A10F1-R]	HA601313	Human, Mouse, Rat	IHC-P,IF-Tissue
KD	P Glycoprotein Recombinant Rabbit mAb [SN06-42]	ET1611-30	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	P Glycoprotein Rabbit pAb	HA500474	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-NF-kB p65 (S529) Recombinant Rabbit mAb [SP07-00]	ET1604-27	Human	WB,IHC-P,FC,IP
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58]	HA723223	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58] - BSA and Azide free	HA751361	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-PKA alpha/beta/gamma (catalytic subunit) (T197) Recombinant Rabbit mAb [JE45-17]	HA721864	Human, Mouse, Rat	WB
	PKA alpha/beta cat Rabbit pAb	ER64617	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	PKA alpha/beta/gamma cat Rabbit pAb	ER64618	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	PKA beta cat Rabbit pAb	ER64619	Human, Mouse, Rat	WB,IHC-P
	PKA C alpha Rabbit pAb	ER1706-65	Human, Mouse, Rat	IF-Cell,IHC-P,FC
	PKA C-alpha Recombinant Rabbit mAb [PSH08-57] - BSA and Azide free	HA751232	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC,IP
	PKA C-alpha Recombinant Rabbit mAb [PSH08-57]	HA723003	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC,IP
KD	PMP70 Recombinant Rabbit mAb [PSH03-27]	HA721969	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	PMP70 Recombinant Rabbit mAb [PSH03-27] - BSA and Azide free	HA750867	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	PMP70 Recombinant Rabbit mAb [JE39-00]	HA721476	Human, Mouse, Rat	WB,IHC-P
	PRKACG Mouse mAb [B3-E4]	EM1709-59	Human	WB,IHC-P,FC
	Rat IL-10 Recombinant Rabbit mAb [PSH02-21] - BSA and Azide free (Capture)	HA721786	Rat	ELISA(Cap)
	Biotin Conjugated Rat IL-10 Recombinant Rabbit mAb [PSH02-22]	HA722637B	Rat	ELISA(Det)
	Rat IL-10 Recombinant Rabbit mAb [PSH02-22] - BSA and Azide free (Detector)	HA721787	Rat	ELISA(Det)
	Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - BSA and Azide free (Detector)	HA721610	Rat	ELISA(Det)
	Rat IL-1 beta Recombinant Rabbit mAb [PS01-86] - BSA and Azide free (Capture)	HA721609	Rat	ELISA(Cap)
	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-59] - BSA and Azide free (Capture)	HA721836	Rat	ELISA(Cap)
	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - Detector	HA721611	Rat	ELISA(Det),ELISA
	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60]	HA722634B	Rat	ELISA(Det)
	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60] - BSA and Azide free (Detector)	HA721837	Rat	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Retinoid X Receptor alpha Recombinant Rabbit mAb [JG99-38]	ET7108-99	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	Scavenging Receptor SR-BI Recombinant Rabbit mAb [SR37-06]	ET1602-32	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	SECR Rabbit pAb	ER62444	Human, Mouse, Rat	WB
	Secretin Receptor Rabbit pAb	ER64920	Human	WB,IHC-P,IF-Cell
	SGLT1 Rabbit pAb	ER1916-51	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	SGLT1 Rabbit pAb	ER64950	Human, Mouse, Rat	WB
	SLC10A2 Rabbit pAb	ER1903-97	Mouse	WB
	SLC22A1 Recombinant Rabbit mAb [JE37-67]	HA722335	Human, Monkey	WB
	SLC22A7 Rabbit pAb	ER1916-71	Mouse	WB
	SLC27A5 Mouse mAb [C10-B11]	EM1708-95	Human, Mouse	WB,IHC-P
	SLC4A4 Rabbit pAb	ER1916-82	Mouse	WB,IHC-P,IF-Tissue
	Slco1a1 Rabbit pAb	ER1916-86	Mouse, Rat	WB,IHC-P,IF-Tissue
	Sodium hydrogen exchanger 1 Rabbit pAb	ER1916-95	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
KD	Sodium Potassium ATPase Recombinant Rabbit mAb [ST0533] - BSA and Azide free	HA750194	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
	iFluor™ 488 Conjugated Sodium Potassium ATPase Recombinant Rabbit mAb [ST0533]	HA720174F	Human, Mouse, Rat	IF-Tissue
	iFluor™ 647 Conjugated Sodium Potassium ATPase Recombinant Rabbit mAb [ST0533]	HA720176F	Human, Mouse, Rat	IF-Tissue,FC
KD	Sodium Potassium ATPase Recombinant Rabbit mAb [ST0533]	ET1609-76	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
	iFluor™ 594 Conjugated Sodium Potassium ATPase Recombinant Rabbit mAb [ST0533]	HA720175F	Human, Mouse, Rat	IF-Tissue
	Sterol carrier protein 2 Recombinant Rabbit mAb [JE62-79]	HA720099	Human, Mouse	WB,IHC-P
	SULT2A1 Recombinant Rabbit mAb [JG36-18]	ET7108-26	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	TGR5 Rabbit pAb	ER65156	Human, Mouse, Rat	WB,IHC-P
	TNF alpha Rabbit pAb	ER1919-22	Human, Mouse	WB
Cell treatment	TNF alpha Rabbit pAb	R1203-1	Mouse	WB,IF-Cell,ELISA
	TNF alpha Rabbit pAb	ER65189	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07]	HA722022	Human, Mouse	WB
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free	HA750897	Human, Mouse	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	UGT1A7 Rabbit pAb	ER1918-38	Mouse, Rat	WB,IHC-P,IF-Tissue
	UGT1A9 Rabbit pAb	ER1918-39	Human, Rat	WB
	Vitamin D3 Receptor Rabbit pAb	ER1803-18	Human, Rat	WB,IF-Cell,FC
	Vitamin D Receptor Rabbit pAb	ER1802-40	Human	WB,FC
KO	Vitamin D Receptor Recombinant Rabbit mAb [JA11-16] - BSA and Azide free	HA750402	Human, Mouse, Rat	WB,IHC-P,FC
KO	Vitamin D Receptor Recombinant Rabbit mAb [JA11-16]	ET1704-09	Human, Mouse, Rat	WB,IHC-P,FC