

CHIP级别抗体清单

高级验证	产品名	货号	反应种属	验证应用
Cell treatment	Acetylated-Lysine Recombinant Rabbit mAb [PSH09-15] - BSA and Azide free	HA751269	Species independent	WB,IHC-P,IF-Tissue,ChIP,IP
Cell treatment	Acetylated-Lysine Recombinant Rabbit Multiclonal Antibody [PSH10-04] - BSA and Azide free	HA751334	Species independent	WB,IHC-P,ChIP,IF-Tissue
Cell treatment	Acetylated-Lysine Recombinant Rabbit mAb [PSH09-15]	HA723073	Species independent	WB,IHC-P,IF-Tissue,ChIP,IP
Cell treatment	Acetylated-Lysine Recombinant Rabbit Multiclonal Antibody [PSH10-04]	HA723173	Species independent	WB,IHC-P,ChIP,IF-Tissue
Cell treatment	ATF6 Recombinant Rabbit mAb [PSH06-59] - BSA and Azide free	HA751087	Human, Mouse, Rat	WB,IHC-P,ChIP
Cell treatment	ATF6 Recombinant Rabbit mAb [PSH06-59]	HA722722	Human, Mouse, Rat	WB,IHC-P,ChIP
	BAF170 Recombinant Rabbit mAb [PSH12-00] - BSA and Azide free	HA751433	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,ChIP
	BAF170 Recombinant Rabbit mAb [PSH12-00]	HA723413	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,ChIP
Relative expression	BATF3 Recombinant Rabbit mAb [PSH13-52]	HA723541	Human	WB,IP,ChIP
Relative expression	BATF3 Recombinant Rabbit mAb [PSH13-52] - BSA and Azide free	HA751471	Human	WB,IP,ChIP
	Beta Catenin Mouse mAb [A6-F8]	EM0306	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,ChIP
	c-Fos Recombinant Rabbit mAb [PSH07-51] - BSA and Azide free	HA751172	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IF-Cell,IHC-P,ChIP,IF-Tissue,IP
	c-Fos Recombinant Rabbit mAb [PSH07-51]	HA722839	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IF-Cell,IHC-P,ChIP,IF-Tissue,IP
	CHD8 Recombinant Rabbit mAb [PSH11-71]	HA723367	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,ChIP
	CHD8 Recombinant Rabbit mAb [PSH11-71] - BSA and Azide free	HA751426	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,ChIP
	c-Myc Recombinant Rabbit mAb [PS00-11]	HA721182	Human, Mouse, Rat	WB,IHC-P,IP,ChIP
	c-Myc Recombinant Rabbit mAb [PS00-11] - BSA and Azide free	HA750519	Human, Mouse, Rat	WB,IHC-P,IP,ChIP

	CTCF Recombinant Rabbit mAb [JM10-61] - BSA and Azide free	HA750394	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,CUT&Tag-seq,ChIP
	CTCF Recombinant Rabbit mAb [JM10-61]	ET1703-90	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,CUT&Tag-seq,ChIP
	Cyclin T1 Recombinant Rabbit mAb [PSH09-66] - BSA and Azide free	HA751308	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP,ChIP
	Cyclin T1 Recombinant Rabbit mAb [PSH09-66]	HA723131	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP,ChIP
	DYKDDDDK Tag (FLAG) Recombinant Rabbit mAb [PSH07-02] - BSA and Azide free	HA751129	Species independent	WB,IHC-P,IF-Cell,FC,IP,ChIP
	DYKDDDDK Tag (FLAG) Recombinant Rabbit mAb [PSH07-02]	HA722780	Species independent	WB,IHC-P,IF-Cell,FC,IP,ChIP
Relative expression	FRA1 Recombinant Rabbit mAb [JE33-20]	HA721447	Human	WB,IHC-P,IF-Cell,FC,ChIP
	HDAC5 Recombinant Rabbit mAb [PSH08-70] - BSA and Azide free	HA751235	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,ChIP
	HDAC5 Recombinant Rabbit mAb [PSH08-70]	HA723022	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,ChIP
KD	Histone H1.2 Recombinant Rabbit mAb [JU43-48]	ET1706-26	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP
Cell treatment	Histone H2A (acetyl K9) Recombinant Rabbit mAb [PS01-45] - BSA and Azide free	HA750591	Human, Mouse, Rat	WB,IHC-P,IF-Cell,CUT&Tag-seq,ChIP
Cell treatment	Histone H2A (acetyl K9) Recombinant Rabbit mAb [SR4-15]	ET1602-34	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP,CUT&Tag-seq
Cell treatment	Histone H2A (acetyl K9) Recombinant Rabbit mAb [PS01-45]	HA721284	Human, Mouse, Rat	WB,IHC-P,IF-Cell,CUT&Tag-seq,ChIP
	Histone H2A.Z Recombinant Rabbit mAb [ST46-09]	ET1609-23	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,ChIP
	Histone H2B Recombinant Rabbit mAb [SD20-63] - BSA and Azide free	HA750287	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP,ChIP
	Histone H2B Recombinant Rabbit mAb [SD20-63]	ET1612-25	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP,ChIP
Cell treatment	Histone H2B (acetyl K20) Recombinant Rabbit mAb [SN72-09]	ET1611-48	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP,ChIP,CUT&Tag-seq
Cell treatment	Histone H2B (acetyl K20) Recombinant Rabbit mAb [PS01-46] - BSA and Azide free	HA750592	Human, Mouse, Rat	WB,IHC-P,IF-Cell,ChIP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Histone H2B (acetyl K20) Recombinant Rabbit mAb [PS01-46]	HA721285	Human, Mouse, Rat	WB,IHC-P,IF-Cell,ChIP
	Histone H3 Recombinant Mouse mAb [6-A7-R] - BSA and Azide free	HA610203	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,ChIP
	Histone H3 Mouse mAb [6-A7]	M1306-4	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,ChIP
	Histone H3 Recombinant Rabbit mAb [JJ090-07]	ET1701-64	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP,IP
	Histone H3 Recombinant Mouse mAb [6-A7-R] - Loading control	HA601335	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,ChIP
	Histone H3 Recombinant Rabbit mAb [JJ092-08]	ET1701-81	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,ChIP
	Histone H3 Recombinant Rabbit mAb [JJ090-07] - BSA and Azide free	HA750318	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP,IP
	Histone H3 Mouse mAb [3-C4]	EM30605	Human, Mouse, Rat	WB,IHC-P,mIHC,IF-Tissue,ChIP
	Histone H3.3 Recombinant Rabbit mAb [ST50-08]	ET1609-59	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,ChIP
	Histone H3.3 Recombinant Rabbit mAb [ST50-08] - BSA and Azide free	HA750184	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,ChIP
Cell treatment	Histone H3 (acetyl K14) Recombinant Rabbit mAb [JU43-26]	ET1706-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,SNAP-ChIP,CUT&Tag-seq
Cell treatment	Histone H3 (acetyl K14) Recombinant Rabbit mAb [JU43-26] - BSA and Azide free	HA750454	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,SNAP-ChIP,CUT&Tag-seq
Cell treatment	Histone H3 (acetyl K18) Mouse mAb [A7B5]	HA600090	Human, Mouse, Rat	WB,IHC-P,ChIP
Cell treatment	Histone H3 (acetyl K27) Mouse mAb [A6D6]	HA600048	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,ChIP
Cell treatment	Histone H3 (acetyl K27) Mouse mAb [A6D7]	HA600047	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,ChIP
Cell treatment	Histone H3 (acetyl K27) Rabbit pAb	HA500046	Human, Mouse, Rat	WB,IHC-P,ChIP
Cell treatment	Histone H3 (acetyl K36) Recombinant Rabbit mAb [JE77-32]	HA722415	Human, Mouse, Rat	WB,IHC-P,IF-Cell,ChIP
Cell treatment	Histone H3 (acetyl K4) Recombinant Rabbit mAb [PSH03-71]	HA722035	Human, Mouse, Rat	WB,IF-Cell,ChIP,Dot Blot
Cell treatment	Histone H3 (acetyl K4) Recombinant Rabbit mAb [PSH03-71] - BSA and Azide free	HA750905	Human, Mouse, Rat	WB,IF-Cell,ChIP,Dot Blot
Cell treatment	Histone H3 (acetyl K56) Recombinant Rabbit mAb [SU30-10]	ET1608-9	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP,CUT&Tag-seq

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Histone H3 (acetyl K56) Recombinant Rabbit mAb [SU30-10] - BSA and Azide free	HA750163	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP,CUT&Tag-seq
Cell treatment	Histone H3 (acetyl K9) Recombinant Rabbit mAb [PSH04-47] - BSA and Azide free	HA750956	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,FC,ChIP,Dot Blot,IP
Cell treatment	Histone H3 (acetyl K9) Recombinant Rabbit mAb [PSH04-47] - ChIP Grade	HA722132	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,FC,ChIP,Dot Blot,IP
Cell treatment	Histone H3 (di methyl K27) Recombinant Rabbit mAb [PSH05-88] - BSA and Azide free	HA751035	Human, Mouse, Rat	WB,IF-Cell,ChIP
Cell treatment	Histone H3 (di methyl K27) Recombinant Rabbit mAb [PSH05-88]	HA722486	Human, Mouse, Rat	WB,IF-Cell,ChIP
	Histone H3 (di methyl K4) Recombinant Rabbit mAb [JE00-98]	HA722135	Human, Mouse, Rat, Monkey	WB,IF-Cell,ChIP,Dot Blot,IHC-P
	Histone H3 (mono+di+tri methyl K36) Recombinant Rabbit mAb [PSH07-14]	HA722802	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,ChIP
	Histone H3 (mono+di+tri methyl K36) Recombinant Rabbit mAb [PSH07-14] - BSA and Azide free	HA751135	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,ChIP
	Histone H3 (mono+di+tri methyl K4) Recombinant Rabbit mAb [PSH06-30]	HA722662	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,ChIP
	Histone H3 (mono+di+tri methyl K4) Recombinant Rabbit mAb [PSH06-30] - BSA and Azide free	HA751067	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,ChIP
	Histone H3 (mono+di+tri methyl K79) Recombinant Rabbit mAb [PSH09-10] - BSA and Azide free	HA751264	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,ChIP
	Histone H3 (mono+di+tri methyl K79) Recombinant Rabbit mAb [SR42-06]	ET1602-41	Human, Mouse, Rat	WB,IHC-P,ChIP
	Histone H3 (mono+di+tri methyl K79) Recombinant Rabbit mAb [PSH09-10]	HA723068	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,ChIP
	Histone H3 (mono+di+tri methyl K79) Recombinant Rabbit mAb [SR42-06] - BSA and Azide free	HA750051	Human, Mouse, Rat	WB,IHC-P,ChIP
	Histone H3 (mono methyl K14) Recombinant Rabbit mAb [JE43-29]	HA722232	Human, Mouse, Rat	WB,IF-Cell,ChIP,Dot Blot
	Histone H3 (mono methyl K27) Recombinant Rabbit mAb [PSH05-67] - BSA and Azide free	HA751019	Human, Mouse, Rat	WB,IF-Cell,FC,ChIP
	Histone H3 (mono methyl K27) Recombinant Rabbit mAb [PSH05-67]	HA722376	Human, Mouse, Rat	WB,IF-Cell,FC,ChIP
	Histone H3 (mono methyl K36) Recombinant Rabbit mAb [SR04-20]	ET1602-48	Human, Mouse	WB,IF-Cell,IF-Tissue,ChIP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Histone H3 (mono methyl K4) Recombinant Rabbit mAb [PSH07-15] - BSA and Azide free	HA751136	Human, Mouse, Rat	WB,IF-Cell,ChIP,Dot Blot
	Histone H3 (mono methyl K4) Recombinant Rabbit mAb [PSH07-15]	HA722803	Human, Mouse, Rat	WB,IF-Cell,ChIP,Dot Blot
KD	Histone H3 (mono methyl K79) Recombinant Rabbit mAb [PSH07-66]	HA722848	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC- P,FC,ChIP
KD	Histone H3 (mono methyl K79) Recombinant Rabbit mAb [PSH07-66] - BSA and Azide free	HA751175	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC- P,FC,ChIP
	Histone H3 (mono methyl K9) Recombinant Rabbit mAb [JE43-27]	HA722350	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC- P,FC,ChIP,IF-Tissue,Dot Blot
	Histone H3 (mono methyl R2) Recombinant Rabbit mAb [ST0427]	ET1609-8	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P,ChIP
	Histone H3 (mono methyl R2) Recombinant Rabbit mAb [ST0427] - BSA and Azide free	HA750197	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P,ChIP
Cell treatment, Relative expression	Histone H3 (mutated K27 M) Recombinant Rabbit mAb [JE46-51]	HA722200	Human, Mouse	WB,IF-Cell,IHC-P,ChIP
Cell treatment	Histone H3 (tri methyl K27) Recombinant Rabbit mAb [PSH05-15] - BSA and Azide free	HA750988	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF- Tissue,ChIP
Cell treatment	Histone H3 (tri methyl K27) Recombinant Rabbit mAb [PSH05-15]	HA722231	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF- Tissue,ChIP
	Histone H3 (tri methyl K36) Recombinant Rabbit mAb [PSH15-01] - BSA and Azide free	HA751539	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,IP,Dot Blot,ChIP
	Histone H3 (tri methyl K36) Recombinant Rabbit mAb [PSH15-01]	HA723701	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,IP,Dot Blot,ChIP
	Histone H4 Recombinant Rabbit mAb [SD209-04] - BSA and Azide free	HA750294	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P,ChIP
	Histone H4 Recombinant Rabbit mAb [SD209-04]	ET1612- 43	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P,ChIP
Cell treatment	Histone H4 (acetyl K16) Recombinant Rabbit mAb [JB21-44] - BSA and Azide free	HA750464	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell,IF- Tissue,ChIP
Cell treatment	Histone H4 (acetyl K16) Recombinant Rabbit mAb [JB21-44]	ET7107- 89	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell,IF- Tissue,ChIP
Cell treatment	Histone H4 (acetyl K5) Recombinant Rabbit mAb [SR31-07]	ET1602- 40	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC- P,IP,ChIP,CUT&Tag-seq
Cell treatment	Histone H4 (acetyl K5) Recombinant Rabbit mAb [SR31-07] - BSA and Azide free	HA750050	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC- P,IP,ChIP,CUT&Tag-seq

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Histone H4 (acetyl K8) Recombinant Rabbit mAb [PS01-25] - BSA and Azide free	HA750568	Human, Mouse, Rat	WB,IF-Cell,IHC-P,ChIP
	Histone H4 (acetyl K8) Recombinant Rabbit mAb [PS01-25]	HA721249	Human, Mouse, Rat	WB,IF-Cell,IHC-P,ChIP
Cell treatment	IRF9 Recombinant Rabbit mAb [PSH08-20] - BSA and Azide free	HA751211	Human	WB,IF-Cell,IP,ChIP
Cell treatment	IRF9 Recombinant Rabbit mAb [PSH08-20]	HA722963	Human	WB,IF-Cell,IP,ChIP
KD	JMJD6 Recombinant Rabbit mAb [PSH15-03]	HA723703	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP,ChIP
KD	JMJD6 Recombinant Rabbit mAb [PSH15-03] - BSA and Azide free	HA751541	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP,ChIP
	JunD Recombinant Rabbit mAb [PSH08-85]	HA723037	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP,ChIP
	JunD Recombinant Rabbit mAb [JE59-51]	HA722913	Human, Mouse, Rat	WB,IF-Cell,IHC-P,ChIP,IP
	JunD Recombinant Rabbit mAb [PSH08-85] - BSA and Azide free	HA751250	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP,ChIP
	KAT6A Recombinant Rabbit mAb [PSH09-88]	HA723157	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,ChIP
	KAT6A Recombinant Rabbit mAb [PSH09-88] - BSA and Azide free	HA751322	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,ChIP
	KDM2B Recombinant Rabbit mAb [PSH05-38] - BSA and Azide free	HA750997	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,IP,ChIP
	KDM2B Recombinant Rabbit mAb [PSH05-38]	HA722267	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,IP,ChIP
	KMT6 / EZH2 Recombinant Rabbit mAb [JJ089-9]	ET1701-56	Human, Mouse	IHC-P,IF-Cell,IF-Tissue,ChIP
Relative expression	Nanog Recombinant Rabbit mAb [SC05-70]	ET1610-2	Human, Mouse	WB,IF-Tissue,FC,IHC-P,ChIP
Relative expression	Nanog Recombinant Rabbit mAb [SC05-70] - BSA and Azide free	HA750204	Human, Mouse	WB,IF-Tissue,FC,IHC-P,ChIP
Cell treatment	Nrf2 Recombinant Rabbit mAb [PSH11-20] - BSA and Azide free	HA751401	Human, Mouse, Rat	WB,IF-Cell,ChIP
Cell treatment	Nrf2 Recombinant Rabbit mAb [PSH11-20]	HA723302	Human, Mouse, Rat	WB,IF-Cell,ChIP
KO	Oct4 Recombinant Rabbit mAb [SD0750]	ET1612-20	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,ChIP
Relative expression, KD, Cell treatment	p53 Recombinant Rabbit mAb [SA39-07] - BSA and Azide free	HA750005	Human	WB,IF-Cell,ChIP
Relative expression, KD, Cell treatment	p53 Recombinant Rabbit mAb [SA39-07]	ET1601-13	Human	WB,IF-Cell,ChIP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Histone H3 (S10) Recombinant Rabbit mAb [SA31-01]	ET1601-30	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,ChIP
Cell treatment	Phospho-Histone H3 (S10) Recombinant Rabbit mAb [SA31-01] - BSA and Azide free	HA750020	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,ChIP
Cell treatment	Phospho-Histone H3 (S28) Recombinant Rabbit mAb [JE56-06]	HA721795	Human, Rat, Mouse	WB,IF-Cell,IHC-P,FC,ChIP
	POLR2A Mouse mAb [A5F4]	HA600061	Human, Mouse, Rat	WB,IHC-P,FC,ChIP
Cell treatment	Recombinant Mouse IgG1 Recombinant Mouse mAb [PSH04-91] - BSA and Azide free	HA610204		WB,FC,IP,ChIP
Cell treatment	Recombinant Mouse IgG1 protein [PSH04-91] - Isotype Control	HA601336		WB,FC,IP,ChIP
Cell treatment	Recombinant Rabbit IgG Recombinant Rabbit mAb [PSH04-42] - BSA and Azide free	HA750951		WB,IHC-P,IF-Cell,FC,IP,ChIP
Cell treatment	Recombinant Rabbit IgG [PSH04-42] - Isotype control	HA722127		WB,IHC-P,IF-Cell,FC,IP,ChIP
KD	SIRT1 Recombinant Rabbit mAb [SZ04-01] - BSA and Azide free	HA750064	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP
KD	SIRT1 Recombinant Rabbit mAb [SZ04-01]	ET1603-3	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP
KD	SIRT3 Recombinant Rabbit mAb [PSH05-22] - BSA and Azide free	HA750993	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue,IP,ChIP
KD	SIRT3 Recombinant Rabbit mAb [PSH05-22]	HA722251	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue,IP,ChIP
	SIRT7 Recombinant Rabbit mAb [PSH09-17] - BSA and Azide free	HA751271	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,ChIP
	SIRT7 Recombinant Rabbit mAb [PSH09-17]	HA723075	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,ChIP
Relative expression	SOX2 Recombinant Rabbit mAb [PO00-28]	HA721155	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,WB,IF-Cell,ChIP,IHC-Fr,IF-Tissue
Relative expression	SOX9 Recombinant Rabbit mAb [PSH13-59] - BSA and Azide free	HA751478	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell,IP,FC,ChIP
Relative expression	SOX9 Recombinant Rabbit mAb [PSH13-59]	HA723548	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell,IP,FC,ChIP
	SP1 Recombinant Rabbit mAb [JF0950] - BSA and Azide free	HA750330	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,ChIP,IP
	SP1 Recombinant Rabbit mAb [JF0950]	ET1702-02	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,ChIP,IP
	SUMO-1 Recombinant Rabbit mAb [SJ20-03] - BSA and Azide free	HA750110	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,ChIP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	SUMO-1 Recombinant Rabbit mAb [SJ20-03]	ET1606-53	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,ChIP
	UHRF1 Recombinant Rabbit mAb [PSH10-63] - BSA and Azide free	HA751366	Human, Mouse, Monkey	WB,IF-Cell,FC,IP,ChIP
	UHRF1 Recombinant Rabbit mAb [PSH10-63]	HA723228	Human, Mouse, Monkey	WB,IF-Cell,FC,IP,ChIP