

自闭症研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	4E-BP1 Mouse mAb [3G1]	RT1002	Human, Mouse, Rat	WB,IP,IF,IHC-P
	4E-BP1 Rabbit pAb	ER62532	Mouse, Rat, Monkey	WB,IHC-P
	ADAR Recombinant Rabbit mAb [JU99-33]	ET7106-72	Human	WB,IF-Cell,IHC-P,IF-Tissue,FC
KO	ADNP Recombinant Rabbit mAb [PSH07-36]	HA722824	Human, Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Cell,IF-Tissue,IP
KO	ADNP Recombinant Rabbit mAb [PSH07-36] - BSA and Azide free	HA751157	Human, Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Cell,IF-Tissue,IP
KD	AKT1 Recombinant Rabbit mAb [ST05-09] - BSA and Azide free	HA750177	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
KD	AKT1 Recombinant Rabbit mAb [ST05-09]	ET1609-47	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
	AKT1 Recombinant Rabbit mAb [PSH0-50] - BSA and Azide free	HA750627	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
	AKT1 Mouse mAb [D9-9-C9]	EM40507	Human, Mouse, Rat	WB,IHC-P
	AKT1 Recombinant Rabbit mAb [PSH0-50]	HA721364	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
	AKT1/2/3 Recombinant Rabbit mAb [ST48-09] - BSA and Azide free	HA750179	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75-09] - BSA and Azide free	HA750840	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75-09]	HA721870	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	AKT1/2/3 Rabbit pAb	ER1901-33	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	AKT1/2/3 Recombinant Rabbit mAb [ST48-09]	ET1609-51	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
	Akt1/3 Rabbit pAb	ER62638	Human, Mouse, Rat	WB,IHC-P
	Akt2 Mouse mAb [2A1]	RT1027	Human, Mouse, Rat	WB,IP,IF,IHC-P
	AKT2 Rabbit pAb	HA500091	Human, Mouse, Rat	WB,IHC-P
	BAF250b Rabbit pAb	ER62788	Human, Rat	WB,IHC-P,IF-Cell
	Beta Catenin Rabbit pAb	0407-16	Human, Mouse, Rat, Zebrafish	WB,IF-Cell
	Beta Catenin Recombinant Rabbit mAb	IRS043	Human	mIHC
	Beta Catenin Recombinant Mouse mAb [A6-F8-R]	HA601179	Human, Mouse, Rat	WB,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	iFluor™ 488 Conjugated Beta Catenin Recombinant Rabbit mAb [SA30-04]	HA720182F	Human, Mouse	IF,ICC
	Beta Catenin Rabbit pAb	ER0805	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell,IF-Tissue
KD	Beta Catenin Recombinant Rabbit mAb [SA30-04]	ET1601-5	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP,mIHC,IF-Cell,IHC-Fr,FC
	Beta Catenin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS278RT	Human	mIHC
	Beta Catenin Mouse mAb [10-C0-B7]	M1405-6	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Beta Catenin Recombinant Rabbit mAb	IRS061	Mouse	mIHC
KD	Beta Catenin Recombinant Rabbit mAb [SA30-04] - BSA and Azide free	HA750023	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP,IF-Cell,IHC-Fr,FC
	Beta Catenin Recombinant Mouse mAb [A6-F8-R] - BSA and Azide free	HA610062	Human, Mouse, Rat	WB,IF-Cell,FC
	Beta Catenin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS277MS	Human	mIHC
	Beta Catenin Mouse mAb [A6-F8]	EM0306	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,ChIP
	CACNA1C Rabbit pAb	ER1803-49	Human, Mouse, Rat	Dot Blot,IF-Cell,IHC-P,FC
	CAMK2D Rabbit pAb	R1305-2	Human	WB,IF-Cell,IHC-P
	CAMK2G Mouse mAb [10C-G1]	EM1706-16	Human, Rat	WB,IHC-P,FC
	CaMKII Recombinant Rabbit mAb [SU03-57] - BSA and Azide free	HA750150	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr,FC
	CaMKII Recombinant Rabbit mAb [SU03-57]	ET1608-47	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr,FC
	CaMKII alpha Recombinant Rabbit mAb [JM11-07]	ET1703-54	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr,IF-Cell,FC,IP
	CaMKII alpha/beta/delta Rabbit pAb	ER62885	Human, Mouse, Rat	IHC-P,WB
	CaMKII alpha/delta Rabbit pAb	ER62886	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	CaMKII beta/gamma Rabbit pAb	ER62887	Human, Mouse, Rat	WB,IF-Cell
	CaMKII beta/gamma/delta Rabbit pAb	ER62888	Human, Mouse, Rat	WB,IHC-P
KD	CaMKII delta Recombinant Rabbit mAb [JA30-07]	ET1704-37	Human, Mouse, Rat, Zebrafish	WB,IHC-P
	Caspr2 Recombinant Rabbit mAb [JE41-03]	HA721598	Human, Mouse, Rat	WB,IHC-P
	CD44 Recombinant Rabbit mAb	IRS123RB	Mouse	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD44 Recombinant Rabbit mAb [ST57-03] - BSA and Azide free	HA750192	Human	WB,IHC-P
Relative expression	CD44 Recombinant Rabbit mAb [ST57-03]	ET1609-74	Human	WB,IHC-P
	CD44 Recombinant mAb	IRS122RB	Human	mIHC
Relative expression, KD	CD44 Recombinant Rabbit mAb [JE64-01]	HA722845	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,FC
	CD44 Mouse mAb [A7-1]	EM010805	Human	WB,IHC-P,FC
	CD62E Rabbit pAb	0805-6	Human	WB,IF-Cell,FC
	CD62L Rabbit pAb	ER1905-95	Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	CD62P Recombinant Rabbit mAb [JM10-95]	ET1703-49	Human	WB,IHC-P
	CHD8 Recombinant Rabbit mAb [PSH11-71]	HA723367	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,ChIP
	CHD8 Rabbit pAb	ER60609	Human	IHC-P
	CHD8 Recombinant Rabbit mAb [PSH11-71] - BSA and Azide free	HA751426	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,ChIP
	Ctip1 Recombinant Rabbit mAb [JB35-68]	ET7107-38	Human, Mouse	WB,IHC-P,FC,IF-Cell,IF-Tissue
	Ctip1 Recombinant Rabbit mAb [JB35-68] - BSA and Azide free	HA750461	Human, Mouse	WB,IHC-P,FC,IF-Cell,IF-Tissue
	eIF4E Recombinant Rabbit mAb [JE60-70]	HA721192	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	eIF4EBP1 Recombinant Rabbit mAb [JJ09-25]	ET1701-83	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	eIF4EBP1 Rabbit pAb	ER1706-64	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	E-Selectin Rabbit pAb	ER1908-35	Human	WB,IHC-P,IF-Tissue
	FMRP Recombinant Rabbit mAb [JM91-41]	ET1703-70	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	GSK3 alpha Rabbit pAb	R1312-1	Human	WB,IF-Cell,IHC-P,FC
	GSK3 alpha Recombinant Rabbit mAb [JE33-27]	HA722217	Human, Mouse, Rat, Monkey	WB,IHC-P
	GSK3 alpha/beta Rabbit pAb	ER63648	Human, Mouse, Rat	WB,IHC-P
	GSK3 beta Rabbit pAb	R1510-4	Human, Mouse, Rat	IF-Cell,IHC-P,FC
KD	GSK3 beta Recombinant Rabbit mAb [SY28-03]	ET1607-71	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
KD	GSK3 beta Recombinant Rabbit mAb [SY28-03] - BSA and Azide free	HA750131	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
	Hamartin Rabbit pAb	HA500199	Human, Mouse, Rat	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Hamartin Mouse mAb [A6F1]	HA600004	Human, Mouse	WB,IHC-P
	Hamartin Recombinant Rabbit mAb [PSH14-67]	HA723665	Human, Mouse, Rat, Monkey	WB,IP
	Hamartin Recombinant Rabbit mAb [PSH14-67] - BSA and Azide free	HA751528	Human, Mouse, Rat, Monkey	WB,IP
	HRP Conjugated Human CD62E Recombinant Rabbit mAb [PSH14-42] - Detector	HA723637H	Human	ELISA(Det),ELISA
	Human CD62E Recombinant Rabbit mAb [PSH14-43] - BSA and Azide free (Detector)	HA723638	Human	ELISA(Det)
	Human CD62E Recombinant Rabbit mAb [PSH14-42] - BSA and Azide free (Capture/Detector)	HA723636	Human	ELISA(Cap),ELISA(Det)
	Human CD62E Recombinant Rabbit mAb [PSH14-41] - BSA and Azide free (Capture)	HA723635	Human	ELISA(Cap)
	Biotin Conjugated Human CD62E Recombinant Rabbit mAb [PSH14-42] - Detector	HA723637B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human CD62E Recombinant Rabbit mAb [PSH14-43] - Detector	HA723639B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human Neurogranin Recombinant Rabbit mAb [PSH16-57] - Detector	HA725192B	Human	ELISA(Det),ELISA
	Human Neurogranin Recombinant Rabbit mAb [PSH16-56] - BSA and Azide free (Capture)	HA725191	Human	ELISA(Cap)
	Human Neurogranin Recombinant Rabbit mAb [PSH16-57] - BSA and Azide free (Detector)	HA725192	Human	ELISA(Det)
Relative expression	L-Selectin / CD62L Recombinant Rabbit mAb [PSH19-78]	HA724102	Mouse	WB,IHC-P,IF-Cell
	MeCP2 Recombinant Rabbit mAb [PSH14-12]	HA723602	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	MeCP2 Recombinant Rabbit mAb [PSH14-12] - BSA and Azide free	HA751507	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	MECP2 Mouse mAb [5B-9H]	EM1706-90	Human	WB,IF-Cell,IHC-P,FC
	mGluR5 Recombinant Rabbit mAb [ST51-00]	ET1609-36	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
	Mouse L-Selectin / CD62L Recombinant Rabbit mAb [PSH16-72] - BSA and Azide free (Capture)	HA723874	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse L-Selectin / CD62L Recombinant Rabbit mAb [PSH16-73] - Detector	HA723875B	Mouse	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HRP Conjugated Mouse L-Selectin / CD62L Recombinant Rabbit mAb [PSH16-73] - Detector	HA723875H	Mouse	ELISA(Det),ELISA
	Mouse L-Selectin / CD62L Recombinant Rabbit mAb [PSH16-73] - BSA and Azide free (Detector)	HA723875	Mouse	ELISA(Det)
	Mouse L-Selectin / CD62L Recombinant Rabbit mAb [PSH16-74] - BSA and Azide free (Capture)	HA723876	Mouse	ELISA(Cap)
	mTOR Rabbit pAb	R1510-21	Human, Mouse, Rat	IF-Cell,IHC-P
	mTOR Recombinant Rabbit mAb [SU30-00] - BSA and Azide free	HA750153	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell
	mTOR Rabbit pAb	HA500126	Human, Rat, Mouse	WB,IF-Cell,IHC-P
	mTOR Recombinant Rabbit mAb [SU30-00]	ET1608-5	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell
	Na+ CP-pan Rabbit pAb	ER64102	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Na+ CP type II alpha Rabbit pAb	ER64098	Human, Mouse, Rat	WB,IHC-P
	Na+ CP type IV alpha Rabbit pAb	ER64100	Human, Mouse, Rat	WB
	NaV1.7 Rabbit pAb	ER1706-44	Human	IF-Cell,IHC-P,FC
	NaV1.7 Rabbit pAb	ER1706-57	Human, Rat	IF-Cell,IHC-P
	NaV1.7 Mouse mAb [1D6]	EM1701-35	Human	IHC-P,IF-Cell,FC,WB
	Neurexin I Rabbit pAb	ER64160	Human, Mouse, Rat	WB
	Neurofibromin Rabbit pAb	ER64168	Human, Mouse, Rat	WB,IF-Cell
	Neurogranin Recombinant Rabbit mAb [PSH02-96] - BSA and Azide free	HA750851	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
	Neurogranin Recombinant Rabbit mAb [PSH02-96]	HA721915	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
	Neurologin 1 Rabbit pAb	ER1913-74	Human, Mouse, Rat	WB,FC
	NLGN3 Rabbit pAb	ER60922	Human, Mouse, Rat	WB
Relative expression	NMDAR1 Recombinant Rabbit mAb [JM11-26] - BSA and Azide free	HA750391	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,IHC-Fr
Relative expression	NMDAR1 Recombinant Rabbit mAb [JM11-26]	ET1703-75	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,IHC-Fr
	Oxytocin R Rabbit pAb	ER1914-56	Human, Mouse, Rat	WB
	Oxytocin Receptor Recombinant Rabbit mAb [JE34-57]	HA722539	Human, Mouse, Rat	WB
	Oxytocin Receptor Rabbit pAb	HA500203	Human, Mouse	WB
	p70 S6 Kinase Rabbit pAb	ER1916-13	Human, Rat	WB,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	p70 S6 Kinase Recombinant Rabbit mAb [PSH13-50] - BSA and Azide free	HA751469	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	p70 S6 Kinase Recombinant Rabbit mAb [PSH13-50]	HA723539	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	p70 S6 Kinase Rabbit pAb	HA500095	Human, Mouse, Rat	WB,IHC-P
	p70 S6 Kinase Recombinant Rabbit mAb [JE03-20]	HA722520	Human, Mouse, Rat	WB,IHC-P,IP
	p70 S6 Kinase beta Rabbit pAb	ER31205	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Phospho-4E-BP1/2/3 (T45) Mouse mAb [3G2]	RT1004	Human, Mouse, Rat	WB,IP,IF
Cell treatment	Phospho-AKT1 (S124) Recombinant Rabbit mAb [JJ08-46]	ET1701-36	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-AKT1 (S129) Recombinant Rabbit mAb [JE59-50]	HA721677	Human, Mouse, Rat	WB
Cell treatment	Phospho-AKT1 (T450) Recombinant Rabbit mAb [SD08-12]	ET1612-73	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell,FC
Cell treatment	Phospho-AKT2 (S474) Recombinant Rabbit mAb [PSH09-46] - BSA and Azide free	HA751293	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-AKT2 (S474) Recombinant Rabbit mAb [PSH09-46]	HA723108	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-AKT2 (S474) Rabbit pAb	HA500116	Human, Mouse	WB,IHC-P
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05] - BSA and Azide free	HA750132	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44]	HA722129	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44] - BSA and Azide free	HA750953	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05]	ET1607-73	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10]	HA722951	Human	WB,IF-Cell,FC
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10] - BSA and Azide free	HA751205	Human	WB,IF-Cell,FC
Cell treatment	Phospho-Beta Catenin (S33 + S37) Recombinant Rabbit mAb [JE59-59]	HA721580	Human, Rat, Mouse	WB
Cell treatment	Phospho-Beta Catenin (S552) Recombinant Rabbit mAb [PSH08-72]	HA723024	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	Phospho-Beta Catenin (S552) Recombinant Rabbit mAb [PSH08-72] - BSA and Azide free	HA751237	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	Phospho-Beta Catenin (T41/S45) Recombinant Rabbit mAb [JE54-02]	HA722316	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-CaMKII (T286) Recombinant Rabbit mAb [JE51-81]	HA721794	Human, Mouse, Rat	WB
	Phospho-CaMKII (T286) Mouse mAb [2G1]	RT1092	Human, Mouse, Rat	WB,IP,IF
Cell treatment	Phospho-eIF4EBP1 (S65) Recombinant Rabbit mAb [PSH09-44] - BSA and Azide free	HA751291	Human	WB,IF-Cell,IP
Cell treatment	Phospho-eIF4EBP1 (S65) Recombinant Rabbit mAb [PSH09-44]	HA723106	Human	WB,IF-Cell,IP
Cell treatment	Phospho-eIF4EBP1 (T70) Recombinant Rabbit mAb [PSH09-91]	HA723222	Human, Mouse	WB,IHC-P,IP
Cell treatment	Phospho-eIF4EBP1 (T70) Recombinant Rabbit mAb [PSH09-91] - BSA and Azide free	HA751360	Human, Mouse	WB,IHC-P,IP
Cell treatment	Phospho-eIF4E (S209) Recombinant Rabbit mAb [SU0396] - BSA and Azide free	HA750160	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
Cell treatment	Phospho-eIF4E (S209) Recombinant Rabbit mAb [SU0396]	ET1608-66	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
Cell treatment	Phospho-GSK3 (alpha + beta) (Y216 + Y279) Recombinant Rabbit mAb [SY26-05]	ET1607-54	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,IF-Tissue,FC
Cell treatment	Phospho-GSK3 (alpha + beta) (Y216 + Y279) Recombinant Rabbit mAb [SY26-05] - BSA and Azide free	HA750127	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,IF-Tissue,FC
Cell treatment	Phospho-GSK3 alpha (S21) Recombinant Rabbit mAb [PSH09-11]	HA723069	Human, Mouse	WB,IHC-P,Dot Blot
Cell treatment	Phospho-GSK3 alpha (S21) Recombinant Rabbit mAb [PSH09-11] - BSA and Azide free	HA751265	Human, Mouse	WB,IHC-P,Dot Blot
Cell treatment	Phospho-GSK3 beta (S9) Recombinant Rabbit mAb [SY02-71]	ET1607-60	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-GSK3 beta (S9) Recombinant Rabbit mAb [SY02-71] - BSA and Azide free	HA750128	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-mTOR (S2448) Mouse mAb [A5D5]	HA600094	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-mTOR (S2481) Recombinant Rabbit mAb [JE59-62]	HA721632	Human	WB
	Phospho-p70 S6 Kinase (S434) Mouse mAb [2G1]	RT1456	Human, Mouse, Rat	WB,IP,IF,IHC-P
Cell treatment	Phospho-p70 S6 Kinase (T389) Recombinant Rabbit mAb [PSH02-23]	HA721803	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-p70 S6 Kinase (T389) Recombinant Rabbit mAb [PSH02-23] - BSA and Azide free	HA750791	Human, Mouse, Rat	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Phospho-PTEN (S380) Recombinant Rabbit mAb [JJ08-81]	ET1701-46	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-Tuberin (S939) Recombinant Rabbit mAb [JE54-06]	HA721762	Human, Mouse	WB
Cell treatment	Phospho-Tuberin (T1462) Recombinant Rabbit mAb [JE43-81]	HA721792	Human, Mouse, Rat	WB
	PSD95 Rabbit pAb	R1508-4	Human, Mouse, Rat	WB,IHC-P,FC
	PSD95 Recombinant Rabbit mAb [SR38-09] - BSA and Azide free	HA750035	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,IP,IHC-Fr
Relative expression	PSD95 Recombinant Antibody [SR38-09] - Rat IgG1 (Chimeric)	HA601416	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,WB
Relative expression	PSD95 Recombinant Antibody [SR38-09] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610256	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,WB
	PSD95 Mouse mAb [B5-E6]	M1511-4	Human, Mouse, Rat	WB,IF-Cell,FC
	PSD95 Recombinant Rabbit mAb [SR38-09]	ET1602-20	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,IHC-Fr
	HRP Conjugated P-Selectin / CD62P Recombinant Rabbit mAb [PSH16-71] - Detector	HA723873H	Human, Mouse, Rat	ELISA(Det),ELISA
	Biotin Conjugated P-Selectin / CD62P Recombinant Rabbit mAb [PSH16-71] - Detector	HA723873B	Human, Mouse, Rat	ELISA(Det),ELISA
	P-Selectin / CD62P Recombinant Rabbit mAb [PSH16-71] - BSA and Azide free (Detector)	HA723873	Human, Mouse, Rat	ELISA(Det)
	P-Selectin / CD62P Recombinant Rabbit mAb [PSH16-70] - BSA and Azide free (Capture)	HA723872	Human, Mouse, Rat	ELISA(Cap)
Relative expression	P-Selectin / CD62P Recombinant Rabbit mAb [PSH19-03]	HA724039	Human, Mouse, Rat	WB,IHC-P,FC,IHC-Fr,IP
	PTEN Mouse mAb [4D8]	RT1519	Human, Mouse, Rat	WB,IP,IF
Relative expression, KD	PTEN Recombinant Rabbit mAb [JE51-96]	HA721962	Human, Mouse, Rat	WB,IHC-P,FC
Relative expression, KD	PTEN Recombinant Rabbit mAb [SJ19-03] - BSA and Azide free	HA750105	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	PTEN Rabbit pAb	ER1706-51	Human, Mouse	WB,IF-Cell,IHC-P,FC
	PTEN Rabbit pAb	ER1902-27	Human, Mouse, Rat	WB,IHC-P,FC
Relative expression, KD	PTEN Recombinant Rabbit mAb [SJ19-03]	ET1606-43	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	PTEN Rabbit pAb	ER1902-33	Human, Mouse, Rat	WB,IHC-P,FC
	QM Rabbit pAb	ER64708	Human, Mouse, Rat	WB
Relative expression	Reelin Recombinant Rabbit mAb [JA63-13]	ET1704-98	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
	SCN1A Rabbit pAb	ER65886	Human, Mouse, Rat	IHC-P,WB
	Shank 2 Rabbit pAb	ER64959	Human, Mouse, Rat	WB,IHC-P
Relative expression	SHANK2 Recombinant Rabbit mAb [PSH10-61] - BSA and Azide free	HA751364	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IF-Tissue
Relative expression	SHANK2 Recombinant Rabbit mAb [PSH10-61]	HA723226	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IF-Tissue
	SLC1A5 Recombinant Rabbit mAb [PSH07-74]	HA722890	Human, Mouse	WB,IF-Cell,IHC-P,FC,IF-Tissue
	SLC1A5 Recombinant Rabbit mAb [PSH07-74] - BSA and Azide free	HA751186	Human, Mouse	WB,IF-Cell,IHC-P,FC,IF-Tissue
	SynGAP Rabbit pAb	ER1902-57	Human, Mouse, Rat	WB,IHC-P,FC
	SynGAP Recombinant Rabbit mAb [JE37-00]	HA722728	Human, Mouse, Rat	WB,IHC-P
	SynGAP Recombinant Rabbit mAb [JE57-79]	HA721860	Mouse, Rat	WB,IHC-P,IF-Tissue
	TMM74 Rabbit pAb	ER62094	Human, Mouse	WB
	TSC2/Tuberin Rabbit pAb	R1510-2	Human, Mouse	IF-Cell,IHC-P,FC
	Tuberin Recombinant Rabbit mAb [SC05-59]	ET1610-10	Human, Mouse, Rat	WB,IHC-P,FC
	UBE3A Recombinant Rabbit mAb [JE32-48]	HA721921	Human	WB
	Wnt-1 Rabbit pAb	ER65317	Human, Mouse	WB,IHC-P,IF-Cell
	WNT2 Rabbit pAb	ER1511-4	Human, Mouse	WB,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层