

脑卒中（脑中风）研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	ABO Recombinant Rabbit mAb [JE39-25]	HA722455	Human, Rat	WB,IP
	ALDH2 Recombinant Mouse mAb [E4-D10-R] - BSA and Azide free	HA610059	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	ALDH2 Rabbit pAb	R1407-3	Human, Mouse, Rat	WB,IF-Cell
	ALDH2 Mouse mAb [E4-D10]	M1509-1	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	ALDH2 Recombinant Mouse mAb [7-D8-R]	HA601277	Human, Rat	WB,IF-Cell,IHC-P
	ALDH2 Mouse mAb [7-D8]	M1501-7	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	ALDH2 Recombinant Mouse mAb [7-D8-R] - BSA and Azide free	HA610152	Human, Rat	WB,IF-Cell,IHC-P
	ALDH2 Recombinant Mouse mAb [E4-D10-R]	HA601175	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Apolipoprotein A1 Recombinant Rabbit mAb [JF0548]	ET1702-23	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Apolipoprotein A1 Mouse mAb	M0809-13	Human	WB,IF-Cell,IHC-P
	Apolipoprotein A1 Rabbit pAb	RT1046	Human, Mouse, Rat	WB,IP,IF
	Apolipoprotein A1 Rabbit pAb	ER1802-29	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Arginase1 Recombinant Rabbit mAb	IRS081RB	Human	mIHC
	Axl Mouse mAb [C4-A8]	M1310-5	Human	WB,IF-Cell
Relative expression	Axl Recombinant Rabbit mAb [PSH06-68]	HA722740	Human, Monkey	WB,IF-Cell,FC
Relative expression	Axl Rabbit pAb	R1406-3	Human, Mouse	WB,IF-Cell,IHC-P,FC
Relative expression	Axl Recombinant Rabbit mAb [PSH06-68] - BSA and Azide free	HA751097	Human, Monkey	WB,IF-Cell,FC
KD	BDNF Recombinant Rabbit mAb [SJ12-09] - BSA and Azide free	HA750104	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC
KD	BDNF Recombinant Rabbit mAb [SJ12-09]	ET1606-42	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC
	BDNF Recombinant Rabbit mAb [JE59-35]	HA722912	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-P,IHC-Fr,FC
	BDNF Rabbit pAb	ER130915	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	CCR5 Rabbit pAb	ER65863	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD195 Mouse mAb [6G11D1]	EM1710-50	Human, Rat	WB,IHC-P,FC
	CD201 Rabbit pAb	ER65529	Human, Mouse, Rat	WB,IHC-P
	CD215 Mouse mAb [5F10D9]	EM1712-22	Human	WB,IF-Cell
	CD40L Recombinant Rabbit mAb [JM11-34]	ET1703-64	Human, Mouse, Rat	WB
	CD40L Recombinant Rabbit mAb [JM11-34] - BSA and Azide free	HA750387	Human, Mouse, Rat	WB
	CKR-5 Rabbit pAb	ER63019	Human, Mouse, Rat	IHC-P
	C Reactive Protein Recombinant Rabbit mAb [JE03-72]	HA722687	Human	WB,IHC-P
	C Reactive Protein Recombinant Rabbit Multiclonal Antibody [PSH15-35]	HA723739	Human, Mouse, Rat	WB
	CRP Mouse mAb [E8-G1]	EM1709-89	Human	WB,IHC-P,FC
	CX3CL1 Rabbit pAb	HA500304	Human	WB,IHC-P
KD	DDAH1 Recombinant Rabbit mAb [JM53-18]	ET1705-71	Human	WB,IP
Relative expression	EPCR Recombinant Rabbit mAb [PSH03-50] - BSA and Azide free	HA750884	Human	WB,IF-Cell,FC
Relative expression	EPCR Recombinant Rabbit mAb [PSH03-50]	HA722009	Human	WB,IF-Cell,FC
	FABP3 Recombinant Rabbit mAb [JE53-37]	HA721859	Human, Mouse, Rat	WB,IHC-P
	FGF21 Recombinant Rabbit mAb [JA10-97] - BSA and Azide free	HA750399	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	FGF21 Recombinant Rabbit mAb [JA10-97]	ET1704-04	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	FGF5 Rabbit pAb	ER1908-83	Human, Mouse, Rat	IHC-P,IF-Tissue
	Follistatin Rabbit pAb	ER1901-69	Human, Rat	WB,IHC-P,IF-Cell,FC
	Follistatin Recombinant Rabbit mAb [JE31-78]	HA722532	Mouse, Rat, Human	WB,IHC-P
	Gas6 Rabbit pAb	ER63503	Human	WB,IHC-P
	GFAP Recombinant Mouse mAb [PSH05-17] - BSA and Azide free (Detector)	HA601342	Human, Mouse, Rat	ELISA(Det)
	iFluor™ 488 Conjugated GFAP Mouse mAb [1-D4]	HA600101F	Human, Rat	IF-Tissue,IHC-Fr
	iFluor™ 594 Conjugated GFAP Mouse mAb [1-D4]	HA600102F	Human, Rat	IF-Tissue
	iFluor™ 647 Conjugated GFAP Mouse mAb [1-D4]	HA600103F	Human, Rat	IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	GFAP Recombinant Rabbit mAb [SA03-04]	ET1601-23	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr,mIHC,IP
	Biotin Conjugated GFAP Mouse mAb [PSH05-17]	HA601355B	Human, Mouse, Rat	ELISA(Det)
	GFAP Mouse mAb [1-D4]	EM140707	Human, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression	GFAP Recombinant Mouse mAb [PSH10-24] - BSA and Azide free	HA610225	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue
	GFAP Rabbit pAb	R1308-10	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	GFAP Recombinant Rabbit mAb [PSH05-16] - BSA and Azide free (Capture)	HA722234	Human, Mouse, Rat	ELISA(Cap)
Relative expression	GFAP Recombinant Mouse mAb [PSH10-24]	HA601380	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue
	GFAP Rabbit pAb	R1308-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	GRK 5 Rabbit pAb	ER63641	Human, Mouse, Rat	WB,IHC-P
	H-FABP Recombinant Rabbit mAb [PSH16-99] - BSA and Azide free (Detector)	HA725178	Human, Mouse	ELISA(Det)
	H-FABP Recombinant Rabbit mAb [PSH16-98] - BSA and Azide free (Capture)	HA725177	Human, Mouse	ELISA(Cap)
	Biotin Conjugated H-FABP Recombinant Rabbit mAb [PSH16-99] - Detector	HA725178B	Human, Mouse	ELISA(Det),ELISA
	Biotin Conjugated Human Apolipoprotein A1 Recombinant Rabbit mAb [PSH18-90] - Detector	HA724027B	Human	ELISA(Det),ELISA
	Human Apolipoprotein A1 Recombinant Rabbit mAb [PSH18-89] - BSA and Azide free (Capture)	HA724026	Human	ELISA(Cap)
	Human Apolipoprotein A1 Recombinant Rabbit mAb [PSH18-91] - BSA and Azide free (Detector)	HA724028	Human	ELISA(Det)
	Biotin Conjugated Human Apolipoprotein A1 Recombinant Rabbit mAb [PSH18-91] - Detector	HA724028B	Human	ELISA(Det),ELISA
	Human Apolipoprotein A1 Recombinant Rabbit mAb [PSH18-90] - BSA and Azide free (Detector)	HA724027	Human	ELISA(Det)
	HRP Conjugated Human CCL5 Recombinant Rabbit mAb [PSH13-17] - Detector	HA725057H	Human	ELISA(Det),ELISA
	Biotin Conjugated Human CCL5 Recombinant Rabbit mAb [PSH13-17] - Detector	HA725057B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CCL5 Recombinant Rabbit mAb [PSH13-16] - BSA and Azide free (Capture)	HA725055	Human	ELISA(Cap)
	Human CCL5 Recombinant Rabbit mAb [PSH13-17] - BSA and Azide free (Detector)	HA725056	Human	ELISA(Det)
	Human CRP Recombinant Rabbit mAb [PSH04-25] - BSA and Azide free (Capture)	HA722108	Human	ELISA(Cap)
	Human CRP Recombinant Rabbit mAb [PSH04-26] - BSA and Azide free (Detector)	HA722109	Human	ELISA(Det)
	Biotin Conjugated Human CRP Recombinant Rabbit mAb [PSH04-26] - Detector	HA722616B	Human	ELISA(Det)
	Human CX3CL1 Recombinant Rabbit mAb [PSH16-08] - BSA and Azide free (Capture)	HA723813	Human	ELISA(Cap)
	Biotin Conjugated Human CX3CL1 Recombinant Rabbit mAb [PSH16-09] - Detector	HA723814B	Human	ELISA(Det),ELISA
	Human CX3CL1 Recombinant Rabbit mAb [PSH16-09] - BSA and Azide free (Detector)	HA723814	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free	HA750744	Human	WB
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39]	HA721712	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free (Detector)	HA721279	Human	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-85] - BSA and Azide free (Capture)	HA723881	Human	ELISA(Cap)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - BSA and Azide free (Detector)	HA723882	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free	HA750745	Human	WB
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA721713	Human	WB
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA722649B	Human	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-87] - BSA and Azide free (Capture)	HA723883	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free (Capture)	HA721278	Human	ELISA(Cap)
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - Detector	HA723882B	Human	ELISA(Det),ELISA
	Human IL-33 Recombinant Rabbit mAb [PSH03-06] - BSA and Azide free (Detector)	HA721936	Human	ELISA(Det)
	Biotin Conjugated Human IL-33 Recombinant Rabbit mAb [PSH03-06]	HA722628B	Human	ELISA(Det)
	Human IL-33 Recombinant Rabbit mAb [PSH03-05] - BSA and Azide free (Capture)	HA721935	Human	ELISA(Cap)
	Human IL-6 Mouse mAb [A6D2]	HA601050	Human	WB,ELISA
	HRP Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA721293	Human	ELISA(Det)
	Human IL-6 Recombinant Rabbit mAb [PS00-21] - BSA and Azide free (Capture)	HA721167	Human	ELISA(Cap)
	Human IL-6 Recombinant Rabbit mAb [PS00-15] - BSA and Azide free (Detector)	HA721168	Human	ELISA(Det)
	Biotin Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA722233B	Human	ELISA(Det)
	Human IL-6 Mouse mAb [A6D1]	HA601051	Human	WB,ELISA
	Biotin Conjugated Human MERTK Recombinant Rabbit mAb [PSH12-70] - Detector	HA723481B	Human	ELISA(Det),ELISA
	HRP Conjugated Human MERTK Recombinant Rabbit mAb [PSH12-70] - Detector	HA723481H	Human	ELISA(Det),ELISA
	Human MERTK Recombinant Rabbit mAb [PSH12-69] - BSA and Azide free (Capture)	HA723479	Human	ELISA(Cap)
	Human MERTK Recombinant Rabbit mAb [PSH12-70] - BSA and Azide free (Detector)	HA723480	Human	ELISA(Det)
	Human MMP-9 Recombinant Rabbit mAb [PSH04-63] - BSA and Azide free (Detector)	HA722163	Human	ELISA(Det)
	Biotin Conjugated Human MMP-9 Recombinant Rabbit mAb [PSH04-24]	HA722615B	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human MMP-9 Recombinant Rabbit mAb [PSH04-62] - BSA and Azide free (Capture)	HA722162	Human	ELISA(Cap)
	Human MMP-9 Recombinant Rabbit mAb [PSH04-24] - BSA and Azide free (Detector)	HA722107	Human	ELISA(Det)
	Human MMP-9 Recombinant Rabbit mAb [PSH04-23] - BSA and Azide free (Capture)	HA722106	Human	ELISA(Cap)
	Human Neuropilin-1 Recombinant Rabbit mAb [PSH14-73] - BSA and Azide free (Capture)	HA725112	Human	ELISA(Cap)
	Human Neuropilin-1 Recombinant Rabbit mAb [PSH14-74] - BSA and Azide free (Detector)	HA725113	Human	ELISA(Det)
	HRP Conjugated Human Neuropilin-1 Recombinant Rabbit mAb [PSH14-74] - Detector	HA725114H	Human	ELISA(Det),ELISA
	Biotin Conjugated Human PD-L1 Recombinant Rabbit mAb [PSH15-37] - Detector	HA723742B	Human	ELISA(Det),ELISA
	Human PD-L1 Recombinant Rabbit mAb [PSH15-36] - BSA and Azide free (Capture)	HA723740	Human	ELISA(Cap)
	Human PD-L1 Recombinant Rabbit mAb [PSH15-37] - BSA and Azide free (Detector)	HA723741	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free (Detector)	HA721266	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb	IRS040	Human	mIHC
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free (Capture)	HA721265	Human	ELISA(Cap)
	HRP Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721294	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40]	HA721286	Human	IHC-P,mIHC
	Biotin Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721559	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free	HA750593	Human	IHC-P
	Biotin Conjugated Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - Detector	HA725066B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - BSA and Azide free (Detector)	HA725065	Human	ELISA(Det)
	Human VCAM-1 Recombinant Rabbit mAb [PSH12-46] - BSA and Azide free (Capture)	HA725064	Human	ELISA(Cap)
	HRP Conjugated Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - Detector	HA725066H	Human	ELISA(Det),ELISA
	Human VEGFA Recombinant Rabbit mAb [PSH12-93] - BSA and Azide free (Detector)	HA723496	Human	ELISA(Det)
	Human VEGFA Recombinant Rabbit mAb [PSH12-91] - BSA and Azide free (Capture)	HA723493	Human	ELISA(Cap)
	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-92] - Detector	HA723495B	Human	ELISA(Det),ELISA
	Human VEGFA Recombinant Rabbit mAb [PSH12-92] - BSA and Azide free (Detector)	HA723494	Human	ELISA(Det)
	Biotin Conjugated Human VEGFA Recombinant Rabbit mAb [PSH12-93] - Detector	HA723497B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human Von Willebrand Factor Recombinant Rabbit mAb [PSH09-04] - BSA and Azide free (Detector)	HA723060B	Human	ELISA(Det),ELISA
	HRP Conjugated Human Von Willebrand Factor Recombinant Rabbit mAb [PSH09-04] - Detector	HA723431H	Human	ELISA(Det)
	Human Von Willebrand Factor Recombinant Rabbit mAb [PSH09-04] - BSA and Azide free (Detector)	HA723059	Human	ELISA(Det)
	Human Von Willebrand Factor Recombinant Rabbit mAb [PSH09-03] - BSA and Azide free (Capture)	HA723058	Human	ELISA(Cap)
	ICAM3 Recombinant Rabbit mAb [JE63-07]	HA722467	Human	WB,IHC-P
Cell treatment	IGF1 Recombinant Rabbit mAb [JE30-01]	HA722527	Human	WB
	IGF1 (E peptide) Rabbit pAb	HA500107	Human	WB,IHC-P,FC
	IGF I Rabbit pAb	ER1911-10	Human, Mouse, Rat, Dog	WB,IHC-P,IF-Tissue
	IL-15R alpha Rabbit pAb	ER63744	Human, Rat	WB
Cell treatment	IL-1 beta Mouse mAb [A7F9]	HA601002	Human, Mouse, Rat	WB,ELISA,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	IL-1 beta Recombinant Rabbit mAb [JJ087-3]	ET1701-39	Human	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R]	HA601036	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R] - BSA and Azide free	HA610004	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Mouse mAb [A7A10]	HA601035	Human	WB,ELISA,IF-Cell
	IL33 Rabbit pAb	ER1911-36	Human, Rat	WB,IHC-P,IF-Tissue
	IL-6 Rabbit pAb	R1412-2	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	IL-6 Mouse mAb	EM30301	Human	WB
	IL-6 Recombinant Rabbit mAb	IRS160RB	Human	mIHC
	LIFR Rabbit pAb	ER1912-33	Human, Mouse, Rat	IHC-P,IF-Tissue
Relative expression	LIFR Recombinant Rabbit mAb [PSH08-76] - BSA and Azide free	HA751241	Human	WB
	LIFR Rabbit pAb	ER65635	Human	WB
Relative expression	LIFR Recombinant Rabbit mAb [PSH08-76]	HA723028	Human	WB
Relative expression	Liver Arginase Recombinant Rabbit mAb [PSH09-78]	HA723143	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Cell treatment	Liver Arginase Recombinant Rabbit mAb [SY09-06]	ET1605-8	Human, Mouse, Rat	WB
Relative expression	Liver Arginase Recombinant Rabbit mAb [JB21-49]	ET7107-90	Human	WB,IHC-P,IP
	Liver Arginase Recombinant Rabbit mAb [PD00-05]	HA721147	Human	WB,IHC-P,IF-Tissue
	MATN3 Rabbit pAb	ER61725	Human, Mouse	WB
KD	MERTK Recombinant Rabbit mAb [SR29-07]	ET1602-21	Human	WB,IHC-P,FC,IF-Tissue
KD	MERTK Recombinant Rabbit mAb [SR29-07] - BSA and Azide free	HA750036	Human	WB,IHC-P,IP,FC,IF-Tissue
KD	MMP-12 Recombinant Rabbit mAb [SR03-23] - BSA and Azide free	HA750052	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	MMP-12 Recombinant Rabbit mAb [SR03-23]	ET1602-42	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cell treatment, KD	MMP-9 Recombinant Rabbit mAb [JA80-73] - BSA and Azide free	HA750420	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC
	MMP-9 Mouse mAb [2G3]	RT1401	Human	WB,IP,IF,IHC-P
	MMP-9 Recombinant Mouse mAb [10A1]	EM1801-22	Human	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	MMP-9 Recombinant Rabbit mAb [PSH06-33] - BSA and Azide free	HA751070	Human	WB,IHC-P
	MMP-9 Recombinant Rabbit mAb	IRS172RB	Human	mIHC
Cell treatment	MMP-9 Recombinant Rabbit mAb [PSH06-33]	HA722665	Human	WB,IHC-P
Cell treatment, KD	MMP-9 Recombinant Rabbit mAb [JA80-73]	ET1704-69	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC
	Mouse CCL5 Recombinant Rat mAb	HA601395	Human, Mouse	ELISA
	Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-15] - BSA and Azide free (Capture)	HA725046	Mouse	ELISA(Cap)
	Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-16] - BSA and Azide free (Detector)	HA725047	Mouse	ELISA(Det)
	Biotin Conjugated Mouse CCL5/RANTES Recombinant Rabbit mAb [PSH14-16] - Detector	HA725048B	Mouse	ELISA(Det),ELISA
	Mouse FGF-21 Recombinant Rabbit mAb [PSH15-27] - BSA and Azide free (Detector)	HA723732	Mouse	ELISA(Det)
	Biotin Conjugated Mouse FGF-21 Recombinant Rabbit mAb [PSH15-27] - Detector	HA723733B	Mouse	ELISA(Det),ELISA
	Mouse FGF-21 Recombinant Rabbit mAb [PSH15-26] - BSA and Azide free (Capture)	HA723731	Mouse	ELISA(Cap)
	Mouse Merck Rabbit pAb	ER1919-54	Mouse, Rat	WB
	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-79] - BSA and Azide free (Capture)	HA723784	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - Detector	HA723788B	Mouse	ELISA(Det),ELISA
	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - BSA and Azide free (Detector)	HA723787	Mouse	ELISA(Det)
	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-80] - BSA and Azide free (Capture/Detector)	HA723785	Mouse	ELISA(Det),ELISA(Cap)
	HRP Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - Detector	HA723788H	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-80] - Detector	HA723786B	Mouse	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse PD-L1 Recombinant Rabbit mAb [PSH04-84]	HA722602B	Mouse	ELISA(Det)
	Mouse PD-L1 Recombinant Rabbit mAb [PSH04-83] - BSA and Azide free (Capture)	HA722189	Mouse	ELISA(Cap)
Cell treatment	Mouse PD-L1 Recombinant Rabbit mAb [PSH05-04]	HA722209	Mouse	FC
Cell treatment	FITC Conjugated Mouse PD-L1 Recombinant Rabbit mAb [PSH05-04]	HA720213F	Mouse	FC
	Mouse PD-L1 Recombinant Rabbit mAb [PSH04-84] - BSA and Azide free (Detector)	HA722190	Mouse	ELISA(Det)
Cell treatment	Mouse PD-L1 Recombinant Rabbit mAb [PSH05-03]	HA722208	Mouse	FC
Relative expression, KD	Neuropilin 1 Recombinant Rabbit mAb [ST05-30] - BSA and Azide free	HA750189	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Neuropilin 1 Mouse mAb [A3-D1]	M0910-5	Human, Mouse	WB,IF-Cell,IHC-P
Relative expression, KD	Neuropilin 1 Recombinant Rabbit mAb [ST05-30]	ET1609-69	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
KD	NSE Recombinant Rabbit mAb [SC06-28]	ET1610-96	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P
KD	NSE Rabbit pAb	R1401-6	Human, Mouse, Rat, Zebrafish, Hybrid fish (crucian-carp)	WB,IF-Cell
	iFluor™ 488 Conjugated NSE Recombinant Rabbit mAb [SC06-28]	HA720178F	Human, Mouse, Rat	IF-Cell,FC
KD	NSE Recombinant Rabbit mAb [SC06-28] - BSA and Azide free	HA750239	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P
	OTUB2 Rabbit pAb	ER64479	Human, Mouse, Rat	WB,IHC-P
	PD-L1 Recombinant Rabbit Multiclonal Antibody [PSH13-23]	HA723510	Human, Mouse	WB
	PD-L1 Recombinant Rabbit Multiclonal Antibody [PSH13-23] - BSA and Azide free	HA751459	Human, Mouse	WB
Cell treatment, Relative expression, KD	PD-L1 Recombinant Rabbit mAb [PD01-02]	HA721176	Human	IHC-P,mIHC,WB
Cell treatment	PD-L1 Recombinant Rabbit mAb [PSH04-80]	HA722184	Mouse	WB,IHC-P,mIHC
	PD-L1 Recombinant Rabbit mAb	IRS009	Human	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PD-L1 Recombinant Rabbit mAb	IRS059	Mouse	mIHC
Cell treatment	PD-L1 Recombinant Rabbit mAb [PSH04-80] - BSA and Azide free	HA750974	Mouse	WB,IHC-P
	PD-L1 Recombinant Rabbit mAb [PD00-97]	HA722300	Human	IHC-P
	PD-L1 Recombinant Rabbit mAb [JJ08-95]	ET1701-41	Human	WB,IHC-P
	PD-L1 Recombinant Rabbit mAb [PD00-97] - BSA and Azide free	HA751015	Human	IHC-P
	PD-L1 Recombinant Rabbit mAb [JJ08-95] - BSA and Azide free	HA750309	Human	WB,IHC-P
	PEDF Recombinant Rabbit mAb [JE41-62]	HA721661	Human, Mouse, Rat	WB
	PEDF Rabbit pAb	ER1914-97	Human, Mouse	WB,IHC-P,IF-Tissue
Cell treatment	Phospho-Axl (Y702) Recombinant Rabbit mAb [PSH17-98]	HA723963	Human	WB
Cell treatment	Phospho-Rac1+Cdc42 (S71) Recombinant Rabbit mAb [JJ086-9]	ET1701-40	Human	WB,IF-Cell,IHC-P
Cell treatment	Phospho-STAT6 (Y641) Recombinant Rabbit mAb [PSH06-98] - BSA and Azide free	HA751125	Human, Mouse	WB
Cell treatment	Phospho-STAT6 (Y641) Recombinant Rabbit mAb [PSH06-98]	HA722776	Human, Mouse	WB
	RAC1 Recombinant Rabbit mAb [PSH05-09]	HA722225	Human, Rat	WB,IF-Cell
	RAC1 Recombinant Rabbit mAb [PSH05-09] - BSA and Azide free	HA750982	Human, Rat	WB,IF-Cell
	Rac1-2-3 Recombinant Rabbit mAb [JM11-29]	ET1703-80	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	RAC1/2/3 Rabbit pAb	HA500424	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	RANTES Recombinant Rabbit mAb [JM03-45]	ET1705-70	Human, Mouse	WB,IHC-P,IF-Cell,IF-Tissue
	S100B Recombinant Rabbit mAb [SC57-02] - BSA and Azide free (Capture)	HA723728	Human, Mouse, Rat	ELISA(Cap)
	Biotin Conjugated S100B Recombinant Rabbit mAb [PSH15-25] - Detector	HA723730B	Human, Mouse, Rat	ELISA(Det),ELISA
	S100B Recombinant Rabbit mAb [PSH15-25] - BSA and Azide free (Detector)	HA723729	Human, Mouse, Rat	ELISA(Det)
	S100 beta Recombinant Antibody [SC57-02] - Guinea pig IgG2 (Chimeric)	HA601550	Mouse	IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	S100 beta Recombinant Rabbit mAb [PSH17-13]	HA723900	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Tissue
	iFluor™ 594 Conjugated S100 beta Recombinant Rabbit mAb [SC57-02]	HA720125F	Human	IF-Tissue
	S100 beta Recombinant Rabbit mAb [PSH17-13] - BSA and Azide free	HA751617	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Tissue
	S100 beta Recombinant Antibody [SC57-02] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610367	Mouse, Rat	IHC-Fr,IHC-P
	S100 beta Recombinant Antibody [SC57-02] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610366	Mouse, Rat	IHC-Fr,IHC-P
	S100 beta Recombinant Rabbit mAb [SC57-02]	ET1610-3	Human, Mouse, Rat, Zebrafish, Goat, Cynomolgus monkey, Green monkey, Pig, Rabbit	WB,IF-Tissue,IHC-P,IHC-Fr,IF-Cell,FC,IP
	S100 beta Recombinant Antibody [SC57-02] - Rat IgG1 (Chimeric)	HA601552	Mouse, Rat	IHC-Fr,IHC-P
Relative expression	S100 beta Recombinant Mouse mAb [PD00-11]	HA601098	Human, Mouse, Rat	IHC-P,WB,IF-Tissue
	S100 beta Recombinant Rabbit mAb [SC57-02] - BSA and Azide free	HA750209	Human, Mouse, Rat, Zebrafish, Goat, Cynomolgus monkey, Green monkey, Pig, Rabbit	WB,IF-Tissue,IHC-P,IHC-Fr,IF-Cell,FC,IP
	S100 beta Rabbit pAb	RT1547	Human, Mouse, Rat	WB,IP,IF
	S100 beta Recombinant Antibody [SC57-02] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610365	Mouse	IHC-Fr
	S100 beta Recombinant Antibody [SC57-02] - Mouse IgG1 (Chimeric)	HA601551	Mouse, Rat	IHC-Fr,IHC-P
	SH3BGRL3 Rabbit pAb	ER65932	Human, Mouse, Rat	IHC-P
	SH3L3 Rabbit pAb	ER61207	Human, Mouse	WB
	SPARC Rabbit pAb	ER1916-99	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
Relative expression	SPARC Recombinant Rabbit mAb [PSH06-23] - BSA and Azide free	HA751061	Human	WB
	SPARC Recombinant Rabbit mAb [PD00-57]	HA721199	Human	IHC-P,IF-Tissue
Relative expression	SPARC Recombinant Rabbit mAb [PSH06-23]	HA722655	Human	WB
KD	STAT6 Rabbit pAb	ER40121	Human	WB,IF-Cell,IHC-P,FC
KD	STAT6 Recombinant Rabbit mAb [SY13-09]	ET1605-49	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,IF-Tissue
KD	STAT6 Recombinant Rabbit mAb [SY02-72]	ET1607-61	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	TMEM16A Rabbit pAb	ER1901-30	Human	WB,IF-Cell,IHC-P,FC
Relative expression	TMEM16A Recombinant Rabbit mAb [PSH06-57]	HA722719	Human, Mouse, Rat	WB,IHC-P
Relative expression	TMEM16A Recombinant Rabbit mAb [PSH06-56] - BSA and Azide free	HA751085	Human	WB,IF-Cell,FC
Relative expression	TMEM16A Recombinant Rabbit mAb [PSH06-56]	HA722720	Human	WB,IF-Cell,FC
	TMEM16A Recombinant Rabbit mAb [PO00-27]	HA721165	Human	IHC-P
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07]	HA722022	Human, Mouse	WB
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free	HA750897	Human, Mouse	WB
	TNF alpha Rabbit pAb	ER65189	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	TNF alpha Rabbit pAb	ER1919-22	Human, Mouse	WB
Relative expression, KD	VCAM1 Recombinant Rabbit mAb [SA05-04]	ET1601-18	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	VCAM1 Rabbit pAb	R1512-13	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression, KD	VCAM1 Recombinant Rabbit mAb [SA05-04] - BSA and Azide free	HA750009	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	VEGF Rabbit pAb	ER30607	Human	WB,IF-Cell,IHC-P,FC
	VEGF Recombinant Rabbit mAb [SP07-01] - BSA and Azide free	HA750074	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	VEGF Recombinant Rabbit mAb [SP07-01]	ET1604-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	VEGFA Recombinant Rabbit mAb [PSH19-17]	HA724051	Human, Mouse, Rat, Monkey	WB
Relative expression, KD, Cell treatment	VISTA Recombinant Rabbit mAb [PSH0-64]	HA721452	Human	WB,IHC-P,IF-Tissue
Relative expression, KD, Cell treatment	VISTA Recombinant Rabbit mAb [PSH0-64] - BSA and Azide free	HA750641	Human	WB,IHC-P,IF-Tissue
	iFluor™ 488 Conjugated VISTA Recombinant Rabbit mAb [PSH0-65]	HA720183F	Human	FC
	VISTA Recombinant Rabbit mAb [PSH0-65] - BSA and Azide free	HA750642	Human	WB,IHC-P,IF-Cell,FC,IF-Tissue
	VISTA Recombinant Rabbit mAb [PSH0-65]	HA721453	Human	WB,IHC-P,IF-Cell,FC,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	Von Willebrand Factor Recombinant Rabbit mAb [PSH07-45] - BSA and Azide free	HA751166	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP
Relative expression	Von Willebrand Factor Recombinant Rabbit mAb [PSH07-45]	HA722833	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP
	Von Willebrand Factor Recombinant Rabbit mAb	IRS215RB	Human	mIHC
	Von Willebrand Factor Recombinant Rabbit mAb [JJ09-34]	ET1701-87	Human, Mouse	WB