

2025年6月—新上市产品清单

高级验证	产品名	货号	反应种属	验证应用
Cell treatment	3mC Recombinant Rabbit mAb [PSH03-21] - BSA and Azide free	HA750861	Species independent	IF-Cell, Dot Blot
Cell treatment	3mC Recombinant Rabbit mAb [PSH03-22] - BSA and Azide free	HA750862	Species independent	ELISA, IF-Cell, Dot Blot
	53BP1 Recombinant Rabbit mAb [JA66-11] - BSA and Azide free	HA750400	Human, Mouse, Rat	WB, IF-Cell, IF-Tissue, IHC-P, FC
	5-HT (Serotonin) Recombinant Rabbit mAb [PSH06-04] - BSA and Azide free	HA751044	Species independent	WB, IHC-P, IF-Tissue, IHC-Fr, ELISA
	5-HT (Serotonin) Recombinant Rabbit mAb [PSH06-05] - BSA and Azide free	HA751045	Species independent	WB, IHC-P, IF-Tissue, IHC-Fr, ELISA
Relative expression	5-Hydroxymethylcytosine (5-hmC) Recombinant Mouse mAb [PSH06-08] - BSA and Azide free	HA610209	Species independent	WB, IHC-P, IF-Cell, Dot Blot, ELISA, MeDIP
Relative expression	5-methylcytosine (5-mC) Recombinant Mouse mAb [PSH06-07] - BSA and Azide free	HA610208	Species independent	WB, IHC-P, ELISA, MeDIP
	6X His tag Recombinant Mouse mAb [A5D1-R] - BSA and Azide free	HA610016	Species independent	WB, IP
	6X His tag Recombinant Rabbit mAb [PSH07-10] - BSA and Azide free	HA751132	Species independent	WB, IF-Cell, IHC-P, FC, IP, ELISA
	6X His tag (C-terminal) Recombinant Mouse mAb [8-F0-R] - BSA and Azide free	HA610015	Species independent	WB, IP
	ABL1 Recombinant Rabbit mAb [PSH04-08] - BSA and Azide free	HA750931	Human, Mouse	WB, IF-Cell, IHC-P, FC
	ACE Recombinant Rabbit mAb [JM59-32] - BSA and Azide free	HA750429	Human, Mouse, Rat	WB, IHC-P, FC
KD	ACE2 Recombinant Rabbit mAb [SN0754] - BSA and Azide free	HA750264	Human, Mouse, Rat, Hamster	WB, IF-Cell, IHC-P, IF-Tissue, IP
Cell treatment	Acetylated-Lysine Recombinant Rabbit mAb [PSH09-16] - BSA and Azide free	HA751270	Species independent	WB, IHC-P, IF-Tissue
Cell treatment	Acetylated-Lysine Recombinant Rabbit mAb [PSH09-15] - BSA and Azide free	HA751269	Species independent	WB, IHC-P, IF-Tissue, ChIP, IP
Cell treatment	Acetylated-Lysine Recombinant Rabbit Multiclonal Antibody [PSH10-04] - BSA and Azide free	HA751334	Species independent	WB, IHC-P, ChIP, IF-Tissue
Relative expression	Acetylcholinesterase Recombinant Rabbit mAb [PSH10-55] - BSA and Azide free	HA751358	Mouse, Rat	WB, IHC-Fr, IHC-P

Relative expression	Acetylcholinesterase Recombinant Antibody [PSH10-55] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610240	Mouse, Rat	WB,IHC-Fr
KD	Acetyl CoA Carboxylase 1 (ACC1) Recombinant Rabbit mAb [ST53-08] - BSA and Azide free	HA750195	Human, Mouse, Rat	WB,IHC-P
	Acid phosphatase/ACP1 Recombinant Rabbit mAb [PSH04-57] - BSA and Azide free	HA750964	Human, Green monkey	WB,IHC-P
	ACSL3 Recombinant Mouse mAb [A8F6-R] - BSA and Azide free	HA610091	Human	WB,IHC-P
KD	ACSL4 / FACL4 Recombinant Rabbit mAb [JE56-33] - BSA and Azide free	HA750469	Human, Mouse, Rat	WB,IHC-P
	ACSL5 Recombinant Mouse mAb [A9D3-R] - BSA and Azide free	HA610053	Human, Mouse, Rat	WB,IHC-P
	ACSL5 Recombinant Mouse mAb [A9D4-R] - BSA and Azide free	HA610050	Mouse, Rat	WB,IHC-P
	Actin Recombinant Antibody [JM36-62] - Human IgG1 (Chimeric) - BSA and Azide free	HA750844	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Actin Recombinant Rabbit mAb [JJ09-29] - BSA and Azide free	HA750323	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment, KD	Active+Pro Caspase-9 Recombinant Rabbit mAb [SC65-05] - BSA and Azide free	HA750238	Human, Mouse	WB,IP
	ACY1 Recombinant Mouse mAb [15G2-R] - BSA and Azide free	HA610111	Human, Mouse, Rat	WB,IHC-P
	ADAM17 Recombinant Rabbit mAb [JM10-35] - BSA and Azide free	HA750370	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,FC,IP,IHC-Fr,IHC-P
	ADAMTS5 Recombinant Rabbit mAb [PSH03-52] - BSA and Azide free	HA750886	Human, Mouse, Rat	WB,IF-Cell,FC
KD	ADFP / PLIN2 Recombinant Rabbit mAb [JA31-81] - BSA and Azide free	HA750407	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC,IHC-P
	Adiponectin Recombinant Rabbit mAb [PSH05-48] - BSA and Azide free	HA751007	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Adiponectin receptor protein 1 Recombinant Rabbit mAb [SC69-04] - BSA and Azide free	HA750233	Human, Mouse, Rat, Zebrafish	IF-Cell,IF-Tissue,IHC-P,FC
KO	ADNP Recombinant Rabbit mAb [PSH07-36] - BSA and Azide free	HA751157	Human, Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Cell,IF-Tissue,IP
	ADRM1 Recombinant Rabbit mAb [PSH05-41] - BSA and Azide free	HA751000	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
Relative expression	AFP Recombinant Rabbit mAb [PSH11-75] - BSA and Azide free	HA751430	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	AGL/Alpha-glucosidase Recombinant Rabbit mAb [PSH04-56] - BSA and Azide free	HA750963	Human, Mouse, Rat	WB,IHC-P
	AHNAK Recombinant Rabbit mAb [PSH01-11] - BSA and Azide free	HA750687	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	AHNAK Recombinant Rabbit mAb [PSH01-10] - BSA and Azide free	HA750686	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	AIDA Recombinant Rabbit mAb [PSH05-81] - BSA and Azide free	HA751029	Human, Mouse, Rat	WB
KO	AIF Recombinant Rabbit mAb [SZ05-01] - BSA and Azide free	HA750066	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	AIMP2 Recombinant Rabbit mAb [PSH0-75] - BSA and Azide free	HA750654	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
KD	AIP Recombinant Rabbit mAb [PSH03-17] - BSA and Azide free	HA750858	Human, Mouse, Rat	WB
	AKR1C1 Recombinant Mouse mAb [5-F10-R] - BSA and Azide free	HA610058	Human	WB,IF-Cell
	AKT1 Recombinant Rabbit mAb [PSH0-50] - BSA and Azide free	HA750627	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
KD	AKT1 Recombinant Rabbit mAb [ST05-09] - BSA and Azide free	HA750177	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75-09] - BSA and Azide free	HA750840	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	AKT1/2/3 Recombinant Rabbit mAb [ST48-09] - BSA and Azide free	HA750179	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
	Albumin Recombinant Rabbit mAb [JF32-10] - BSA and Azide free	HA750352	Human, Mouse, Rat, Cow	WB,IHC-P,IF-Tissue
	Aldehyde Oxidase Recombinant Rabbit mAb [PSH0-17] - BSA and Azide free	HA750589	Human, Mouse, Rat	WB,IHC-P,FC
KD	ALDH1A1 Recombinant Mouse mAb [8F1-R] - BSA and Azide free	HA610093	Human	WB,IF-Cell,FC
	ALDH1L1 Recombinant Rabbit mAb [JU53-54] - BSA and Azide free	HA750457	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
	ALDH2 Recombinant Mouse mAb [E4-D10-R] - BSA and Azide free	HA610059	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KD	ALDH2 Recombinant Mouse mAb [7-D8-R] - BSA and Azide free	HA610152	Human, Rat	WB,IF-Cell,IHC-P
	ALDH3A1 Recombinant Rabbit mAb [PD00-73] - BSA and Azide free	HA750533	Human	WB,IHC-P
	ALDH4A1 Recombinant Mouse mAb [7E2-R] - BSA and Azide free	HA610192	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD, Relative expression	Aldolase A/B/C Recombinant Rabbit mAb [PSH01-81] - BSA and Azide free	HA750755	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	Aldolase A/B/C Recombinant Antibody [PSH01-81] - Mouse IgG1 (Chimeric)	HA601511	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Aldolase A/B/C Recombinant Antibody [PSH01-81] - Rat IgG1 (Chimeric)	HA601512	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Aldolase C Recombinant Rabbit mAb [PSH01-82] - BSA and Azide free	HA750756	Mouse, Rat, Human	WB,IHC-P
KD	ALIX Recombinant Rabbit mAb [JM85-31] - BSA and Azide free	HA750439	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	Alkaline Phosphatase Recombinant Rabbit mAb [SA40-00] - BSA and Azide free	HA750013	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	ALKBH5 Recombinant Rabbit mAb [PSH01-16] - BSA and Azide free	HA750692	Human, Mouse, Rat	WB,IHC-P
Relative expression	ALP Recombinant Rabbit mAb [PSH14-35] - BSA and Azide free	HA751517	Human, Mouse, Rat	WB,IP
	Alpha 1 Acid Glycoprotein Recombinant Mouse mAb [A4-G2-A3-R] - BSA and Azide free	HA610101	Human	WB
	alpha Actinin Recombinant Mouse mAb [A1G2-R] - BSA and Azide free	HA610074	Human, Mouse, Rat	WB,IHC-P
	Alpha-cardiac actin Recombinant Mouse mAb [33-32-R] - BSA and Azide free	HA610086	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
KD, Relative expression	alpha smooth muscle Actin Recombinant Rabbit mAb [SY02-64] - BSA and Azide free	HA750126	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD, Relative expression	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610267	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment, Relative expression	Alpha-Synuclein Recombinant Rabbit mAb [PSH08-83] - BSA and Azide free	HA751248	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	alpha Tubulin Recombinant Rabbit mAb [PSH02-95] - BSA and Azide free	HA750850	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	alpha Tubulin Recombinant Rabbit mAb [PSH02-94] - BSA and Azide free	HA750849	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	alpha Tubulin Recombinant Rabbit mAb [PSH0-14] - BSA and Azide free	HA750584	Human, Mouse, Rat	WB,IHC-P
Cell treatment	alpha Tubulin (acetyl K40) Recombinant Rabbit mAb [PSH0-13] - BSA and Azide free	HA750583	Human, Mouse, Rat	WB,IF-Cell
Relative expression	AMACR Recombinant Rabbit mAb [PSH13-55] - BSA and Azide free	HA751474	Human, Mouse, Rat	WB,IHC-P
Relative expression	AMH Recombinant Rabbit mAb [PSH13-48] - BSA and Azide free	HA751467	Mouse, Rat, Human	WB,IHC-Fr,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KO	AMPK alpha 1 Recombinant Rabbit mAb [SU03-48] - BSA and Azide free	HA750147	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC,IP,IHC-Fr,IHC-P
	AMPK alpha 1/2 Recombinant Rabbit mAb [PSH08-86] - BSA and Azide free	HA751251	Human, Mouse, Rat	WB,IHC-P
Relative expression	Angiomotin Recombinant Rabbit mAb [PSH07-46] - BSA and Azide free	HA751167	Human	WB,IHC-P,IF-Cell,IP
	Angiopoietin 2 Recombinant Rabbit mAb [JM71-34] - BSA and Azide free	HA750436	Human, Mouse, Rat	WB,IHC-P,FC
KD	ANGPTL4 Recombinant Rabbit mAb [PSH07-35] - BSA and Azide free	HA751156	Human, Monkey	WB
	Annexin A3 Recombinant Mouse mAb [A1F7-R] - BSA and Azide free	HA610176	Human, Mouse, Rat	WB
Relative expression	Antithrombin III / ATIII Recombinant Rabbit mAb [PSH12-13] - BSA and Azide free	HA751442	Human, Mouse, Rat	WB,IHC-P
	ANXA6 Recombinant Rabbit mAb [PSH05-82] - BSA and Azide free	HA751030	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP
KD	APE1 Recombinant Mouse mAb [12H1-R] - BSA and Azide free	HA610089	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	APE1 Recombinant Mouse mAb [12H2-R] - BSA and Azide free	HA610090	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment, KD	APOE4 Recombinant Rabbit mAb [PSH02-85] - BSA and Azide free	HA750838	Human	WB
Relative expression	Apolipoprotein E Recombinant Antibody [PSH10-52] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610233	Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Apolipoprotein E Recombinant Rabbit mAb [PSH10-52] - BSA and Azide free	HA751355	Mouse, Rat	WB,IHC-Fr,IHC-P
	Apolipoprotein E Recombinant Rabbit mAb [SC0536] - BSA and Azide free	HA750206	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment, Relative expression	APP/ β -Amyloid Recombinant Rabbit mAb [PSH04-90] - BSA and Azide free	HA750977	Human, Mouse	WB,IF-Cell,IHC-P
	AQP1 Recombinant Rabbit mAb [JM10-98] - BSA and Azide free	HA750378	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P
	Aquaporin 2 Recombinant Rabbit mAb [PSH09-19] - BSA and Azide free	HA751273	Human, Mouse, Rat	IHC-P,IHC-Fr,IF-Tissue
	Aquaporin 3 Recombinant Rabbit mAb [PSH06-52] - BSA and Azide free	HA751081	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Aquaporin 4 Recombinant Antibody [PSH06-38] - Chicken IgY (Chimeric) - BSA and Azide free	HA610294	Mouse	IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	Aquaporin 4 Recombinant Antibody [PSH06-38] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610295	Mouse	IHC-Fr,WB
Relative expression	Aquaporin 4 Recombinant Rabbit mAb [PSH06-38] - BSA and Azide free	HA751077	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr,IP
Relative expression	Aquaporin 4 Recombinant Antibody [PSH06-38] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610296	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Arc / Arg3.1 Recombinant Rabbit mAb [PSH10-67] - BSA and Azide free	HA751369	Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
KO	Argonaute 2 Recombinant Rabbit mAb [JF0992] - BSA and Azide free	HA750344	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	ARPC5 / p16 ARC Recombinant Rabbit mAb [SR34-02] - BSA and Azide free	HA750055	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	AS160 Recombinant Rabbit mAb [PSH10-95] - BSA and Azide free	HA751385	Human, Rat	WB,IHC-P,IP
	ASAH1 Recombinant Rabbit mAb [PSH07-71] - BSA and Azide free	HA751183	Human	WB,IHC-P
	ASC1 Recombinant Rabbit mAb [PSH0-26] - BSA and Azide free	HA750600	Human, Rat	WB,IHC-P
Relative expression	ASGPR1 Recombinant Rabbit mAb [PSH06-01] - BSA and Azide free	HA751041	Human, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,IP
Cell treatment	ASPA Recombinant Rabbit mAb [PSH08-04] - BSA and Azide free	HA751204	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
	ASPA Recombinant Antibody [PSH08-04] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610269	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	ASPA Recombinant Antibody [PSH08-04] - Mouse IgG1 (Chimeric)	HA601498	Human, Mouse, Rat	WB,IHC-P
KD	ASS1 Recombinant Mouse mAb [11F1-RA] - BSA and Azide free	HA610025	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Asymmetric Di-Methyl Arginine (ADMA) Recombinant Rabbit mAb [PSH09-55] - BSA and Azide free	HA751302	Species independent	ELISA,WB
Cell treatment	ATF-4 Recombinant Rabbit mAb [PSH07-08] - BSA and Azide free	HA751130	Human, Mouse, Rat	WB,IP
	ATF4 Recombinant Rabbit mAb [SD20-92] - BSA and Azide free	HA750292	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	ATF5 Recombinant Rabbit mAb [SD2099] - BSA and Azide free	HA750293	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cell treatment	ATF6 Recombinant Rabbit mAb [PSH06-59] - BSA and Azide free	HA751087	Human, Mouse, Rat	WB,IHC-P,ChIP
KD	ATF6 Recombinant Mouse mAb [8D3-R] - BSA and Azide free	HA610193	Human, Mouse, Rat	WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	ATG4B Recombinant Rabbit mAb [PSH02-53] - BSA and Azide free	HA750816	Human, Mouse, Rat, Monkey	WB
	ATG4C Recombinant Rabbit mAb [PSH05-23] - BSA and Azide free	HA750994	Human, Mouse, Rat	WB,IP
KO	ATG5 Recombinant Rabbit mAb [SN73-07] - BSA and Azide free	HA750251	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	ATG7 Recombinant Rabbit mAb [SC06-30] - BSA and Azide free	HA750224	Human, Rat	WB,IF-Tissue,IHC-P,IP,FC
	ATL3 Recombinant Rabbit mAb [PSH02-47] - BSA and Azide free	HA750812	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP
KO	ATM Recombinant Rabbit mAb [SI70-01] - BSA and Azide free	HA750096	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	ATOX1 Recombinant Rabbit mAb [PSH08-50] - BSA and Azide free	HA751225	Human	WB,IF-Cell,IHC-P,FC
	ATP13A2 Recombinant Rabbit mAb [PSH06-69] - BSA and Azide free	HA751098	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IHC-Fr
	ATP5F1 Recombinant Mouse mAb [A10H10-R] - BSA and Azide free	HA610210	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	ATP5O Recombinant Rabbit mAb [PSH15-04]	HA723704	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,FC,IP
KO	ATP citrate lyase Recombinant Rabbit mAb [ST51-07] - BSA and Azide free	HA750171	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	ATRX/Daxx Antibody Sampler Kit	HAK21138	Human	WB
KD	Aurora B Recombinant Rabbit mAb [SC55-08] - BSA and Azide free	HA750208	Human	WB,IHC-P,IP,IF-Cell
KD	Axin1 Recombinant Rabbit mAb [PSH03-00] - BSA and Azide free	HA750855	Human	WB,IHC-P
Relative expression	Axl Recombinant Rabbit mAb [PSH06-68] - BSA and Azide free	HA751097	Human, Monkey	WB,IF-Cell,FC
Relative expression	B7H4 Recombinant Rabbit mAb [PSH03-65] - BSA and Azide free	HA750899	Human, Mouse, Rat	WB,IF-Cell
Relative expression	B7-H6 Recombinant Rabbit mAb [PSH10-89] - BSA and Azide free	HA751381	Human	WB,IF-Cell,FC,IP
	BACE1 Recombinant Rabbit mAb [PSH05-19] - BSA and Azide free	HA750990	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	Bad Recombinant Rabbit mAb [SA32-01] - BSA and Azide free	HA750025	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	BAF170 Recombinant Rabbit mAb [PSH12-00] - BSA and Azide free	HA751433	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,ChIP
KD	BAF180 Recombinant Rabbit mAb [PSH02-34] - BSA and Azide free	HA750799	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
Relative expression	BATF3 Recombinant Rabbit mAb [PSH13-52] - BSA and Azide free	HA751471	Human	WB,IP,ChIP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KO	Bax Recombinant Rabbit mAb [SZ3-07] - BSA and Azide free	HA750065	Human, Mouse, Rat, Goat	WB,IF-Cell,FC
	BAX Recombinant Mouse mAb [PSH0-20] - BSA and Azide free	HA610038	Human	WB,IHC-P
	Bcl10 Recombinant Rabbit mAb [PSH13-49] - BSA and Azide free	HA751468	Human, Mouse, Rat	WB,IF-Cell,FC,IP
KD	Bcl-2 Recombinant Rabbit mAb [SZ10-03] - BSA and Azide free	HA750057	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Bcl-2 Recombinant Rabbit mAb [PSH09-49] - BSA and Azide free	HA751296	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
KD	Bcl-2 Recombinant Rabbit mAb [JF104-8] - BSA and Azide free	HA750350	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Bcl-2 Recombinant Mouse mAb [9F3-R] - BSA and Azide free	HA610131	Human	WB,IHC-P,FC
	BCL2A1 Recombinant Rabbit mAb [SC05-42] - BSA and Azide free	HA750205	Human, Mouse	WB,IHC-P
	BCL6 Recombinant Mouse mAb [A6C7-R] - BSA and Azide free	HA610020	Human	IHC-P,IF-Cell
	Bcl-XL Recombinant Rabbit mAb [SZ3-03] - BSA and Azide free	HA750063	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	BDH1 Recombinant Rabbit mAb [PSH14-69]	HA723667	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P
KD	BDNF Recombinant Rabbit mAb [SJ12-09] - BSA and Azide free	HA750104	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC
	Beclin 1 Recombinant Rabbit mAb [JE59-31] - BSA and Azide free	HA750546	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	beta 2 Adrenergic Receptor Recombinant Rabbit mAb [JM102-06] - BSA and Azide free	HA750368	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IP,IF-Cell,FC,IHC-Fr
Relative expression	beta 2 Microglobulin Recombinant Rabbit mAb [PSH12-14] - BSA and Azide free	HA751443	Human	WB,IF-Cell,IHC-P,FC,IP
	beta Actin Recombinant Rabbit mAb [PSH03-63] - BSA and Azide free	HA750898	Human, Mouse, Rat, Rabbit, Goat, Zebrafish, Dog, Monkey, Pig, Chicken, Cow, Hamster, Green monkey	WB,IF-Cell,FC,IP
	beta Actin Recombinant Mouse mAb [A2-F6-R] - BSA and Azide free	HA610019	Human, Mouse, Rat	WB
	beta Amyloid Recombinant Rabbit mAb [PSH08-18] - BSA and Azide free	HA751209	Human, Mouse	WB,Dot Blot
Relative expression	beta Amyloid Recombinant Rabbit mAb [PSH08-19] - BSA and Azide free	HA751210	Human, Mouse	WB,IHC-P,IF-Tissue,Dot Blot
Relative expression	beta Amyloid 1-40 Recombinant Rabbit mAb [PSH02-82] - BSA and Azide free	HA750789	Human, Mouse, Cynomolgus monkey, Pig	IHC-P,WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	beta Amyloid 1-40 Recombinant Rabbit mAb [PSH05-86] - BSA and Azide free	HA751034	Human, Mouse	WB,IHC-P,IF-Tissue
Relative expression	beta Amyloid 1-42 Recombinant Rabbit mAb [PSH02-83] - BSA and Azide free	HA750790	Human, Mouse, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue
Relative expression	beta Amyloid 1-42 Recombinant Rabbit mAb [PSH06-81] - BSA and Azide free	HA751103	Human, Mouse	WB,IHC-P,IF-Tissue
KD	Beta Catenin Recombinant Rabbit mAb [SA30-04] - BSA and Azide free	HA750023	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP,IF-Cell,IHC-Fr,FC
	Beta Catenin Recombinant Mouse mAb [A6-F8-R] - BSA and Azide free	HA610062	Human, Mouse, Rat	WB,IF-Cell,FC
	Beta III Tubulin Recombinant Rabbit mAb [SP06-00] - BSA and Azide free	HA750071	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
	beta Tubulin Recombinant Mouse mAb [1-B11-R] - BSA and Azide free	HA610049	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	beta Tubulin Recombinant Rabbit mAb [SR25-04] - BSA and Azide free	HA750049	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	beta Tubulin Recombinant Mouse mAb [A1-A4-R] - BSA and Azide free	HA610067	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IHC-P
	BHMT Recombinant Mouse mAb [9E1-R] - BSA and Azide free	HA610119	Human, Mouse, Rat	WB,IHC-P
	Bid Recombinant Rabbit mAb [PSH14-03] - BSA and Azide free	HA751499	Human, Mouse	WB,IP
Relative expression	Bik Recombinant Rabbit mAb [PSH11-51] - BSA and Azide free	HA751413	Human	WB,IF-Cell,IP
	Bim Recombinant Rabbit mAb [SU0318] - BSA and Azide free	HA750137	Human, Mouse	WB,IF-Cell,IHC-P,IP,FC
Relative expression	BLBP Recombinant Rabbit mAb [PSH10-68] - BSA and Azide free	HA751370	Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
Relative expression	BLBP Recombinant Rabbit mAb [PSH0-89] - BSA and Azide free	HA750669	Human	WB,IHC-P
Relative expression	BLBP Recombinant Rabbit mAb [PSH0-88] - BSA and Azide free	HA750668	Human	WB
Relative expression	BLBP Recombinant Antibody [PSH10-68] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610251	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
Cell treatment	BOLF1 Recombinant Rabbit mAb [PSH08-88] - BSA and Azide free	HA751253	Epstein-barr virus	WB,IF-Cell,IHC-P
Cell treatment	BOLF1 Recombinant Rabbit mAb [PSH08-87] - BSA and Azide free	HA751252	Epstein-barr virus	WB,IF-Cell,IHC-P,FC
	Bone sialoprotein 2 Recombinant Rabbit mAb [PSH09-64] - BSA and Azide free	HA751306	Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	BRCA1 Recombinant Mouse mAb [A7E3-R] - BSA and Azide free	HA610177	Human, Mouse, Rat	IHC-P,IF-Cell
	BRCA2 Recombinant Rabbit mAb [PSH07-55] - BSA and Azide free	HA751177	Human	WB,FC
	BRCC36 Recombinant Rabbit mAb [PSH05-10] - BSA and Azide free	HA750983	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	BRCC45/BRE Recombinant Rabbit mAb [PSH06-15] - BSA and Azide free	HA751053	Human, Mouse, Rat	WB,IF-Cell,FC
	BrdU Recombinant Rabbit mAb [PSH0-18] - BSA and Azide free	HA750590	Species independent	IHC-P,IF-Cell,IF-Tissue,FC
	BrdU Recombinant Antibody [PSH0-18] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610099	Species independent	IHC-P,IF-Cell,IF-Tissue
	BrdU Recombinant Mouse mAb [A9C7-R] - BSA and Azide free	HA610183	Species independent	IHC-P,IF-Tissue,IF-Cell,FC
Relative expression	BRG1 Recombinant Mouse mAb [A5C9-R] - BSA and Azide free	HA610205	Human, Mouse, Rat	WB
	BRN3A Recombinant Antibody [PSH11-65] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610252	Mouse	IHC-Fr,IHC-P
	BRN3A Recombinant Rabbit mAb [PSH11-65] - BSA and Azide free	HA751421	Mouse, Rat	IHC-Fr,IHC-P,WB
KO, Relative expression	BTK Recombinant Rabbit mAb [SN69-07] - BSA and Azide free	HA750252	Human, Mouse	WB,IHC-P,IP
	C10orf58 Recombinant Mouse mAb [7-B7-R] - BSA and Azide free	HA610088	Human, Rat	WB,IF-Cell,FC
Cell treatment	c17orf74 Recombinant Rabbit mAb [PSH02-66] - BSA and Azide free	HA750823	Human, Mouse, Rat	WB,IHC-P
	C9orf72 Recombinant Rabbit mAb [PSH11-11] - BSA and Azide free	HA751395	Human, Mouse, Rat	WB
	CA12 Recombinant Rabbit mAb [PSH05-49] - BSA and Azide free	HA751008	Human	WB,IF-Cell,FC
	CAB39 Recombinant Rabbit mAb [PSH15-16]	HA723720	Human, Mouse, Rat, Monkey	WB,IHC-P,IP
Relative expression	Calbindin Recombinant Rabbit mAb [JF05-01] - BSA and Azide free	HA750351	Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,IHC-Fr
	Calmodulin Recombinant Rabbit mAb [SJ16-09] - BSA and Azide free	HA750106	Human, Mouse, Dog, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
KD	Calnexin Recombinant Rabbit mAb [SN20-54] - BSA and Azide free	HA750275	Human, Mouse, Rat	WB,IP,IHC-P,IF-Tissue
	Calponin Recombinant Rabbit mAb [SI67-01] - BSA and Azide free	HA750094	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
KD	Calreticulin Recombinant Rabbit mAb [SU37-03] - BSA and Azide free	HA750156	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Calretinin Recombinant Mouse mAb [PSH08-29] - BSA and Azide free	HA610215	Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Tissue,IP
Relative expression	Calsequestrin 2 Recombinant Rabbit mAb [PSH03-16] - BSA and Azide free	HA750857	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	CaMKII Recombinant Rabbit mAb [SU03-57] - BSA and Azide free	HA750150	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr,FC
	CAMKK2 Recombinant Rabbit mAb [PSH07-28] - BSA and Azide free	HA751149	Human	WB
Relative expression	Cannabinoid receptor 1 Recombinant Rabbit mAb [PSH02-43] - BSA and Azide free	HA750808	Human, Mouse, Rat	WB,IHC-P
	CAPZA1 Recombinant Mouse mAb [A2E7-R] - BSA and Azide free	HA610155	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	Carbonic anhydrase 2 Recombinant Rabbit mAb [JU12-37] - BSA and Azide free	HA750456	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
	Carbonic anhydrase 2 Recombinant Mouse mAb [11A2-R] - BSA and Azide free	HA610116	Human	WB,IF-Cell,IHC-P
	Carcino Embryonic Antigen CEA Recombinant Rabbit mAb [JM93-28] - BSA and Azide free	HA750434	Human	WB,IF-Cell,IF-Tissue,IHC-P
	CARD9 Recombinant Rabbit mAb [PSH0-41] - BSA and Azide free	HA750619	Human, Mouse, Rat	WB,IF-Cell,FC
	Cardiac Troponin T Recombinant Rabbit mAb [SC64-03] - BSA and Azide free	HA750222	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	CARM1 Recombinant Rabbit mAb [PSH05-08] - BSA and Azide free	HA750981	Human, Mouse, Rat	WB,IHC-P,IP
	Casein Kinase 1 alpha Recombinant Rabbit mAb [PSH04-46] - BSA and Azide free	HA750955	Human, Mouse, Rat	WB,IF-Cell,FC,IP
KD	CASK Recombinant Mouse mAb [A1C5-R] - BSA and Azide free	HA610132	Human, Mouse, Rat	WB,IHC-P
	Caspase-1 Recombinant Rabbit mAb [SU40-07] - BSA and Azide free	HA750161	Human, Mouse, Rat	WB,FC,IF-Cell
Relative expression	Caspase-12 Recombinant Rabbit mAb [PSH07-49] - BSA and Azide free	HA751170	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Caspase-7 Recombinant Rabbit mAb [PSH13-69] - BSA and Azide free	HA751482	Human, Mouse, Rat	WB,IHC-P
Cell treatment, KD	Caspase-8 Recombinant Rabbit mAb [PSH04-77] - BSA and Azide free	HA750971	Human	WB,IF-Cell,IHC-P,FC,IP
Cell treatment	Caspase-8 (p18) Recombinant Rabbit mAb [PSH05-83] - BSA and Azide free	HA751031	Mouse	WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Caspase-9 Recombinant Rabbit mAb [SZ29-01] - BSA and Azide free	HA750062	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Relative expression	CaSR Recombinant Rabbit mAb [PSH03-47] - BSA and Azide free	HA750881	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
KD	Catalase Recombinant Rabbit mAb [JM11-12] - BSA and Azide free	HA750377	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IF-Cell
	Cathepsin D Recombinant Mouse mAb [13F4-R] - BSA and Azide free	HA610002	Human	WB,IHC-P
	Cathepsin D Recombinant Rabbit mAb [SU0360] - BSA and Azide free	HA750152	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	Cathepsin D Recombinant Mouse mAb [13F3-R] - BSA and Azide free	HA610001	Human	WB,IHC-P
KD	Cathepsin L Recombinant Rabbit mAb [PSH03-91] - BSA and Azide free	HA750916	Human, Mouse, Rat, Monkey	WB,IHC-P,FC
KD	Cathepsin L Recombinant Rabbit mAb [PSH03-92] - BSA and Azide free	HA750917	Human, Mouse, Rat, Monkey	WB,IHC-P
Relative expression	CaV2.3 / CACNA1E Recombinant Rabbit mAb [PSH13-51] - BSA and Azide free	HA751470	Mouse, Rat	IHC-Fr,IHC-P,WB
	CBLB Recombinant Rabbit mAb [PSH03-49] - BSA and Azide free	HA750883	Human, Mouse, Rat	WB,IHC-P
	Cbx8 Recombinant Rabbit mAb [PSH06-61] - BSA and Azide free	HA751091	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-44]	HA723691	Mouse	WB,IF-Cell,FC
Cell treatment, Relative expression	CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-27] - BSA and Azide free	HA751406	Human	WB,IF-Cell,FC
Cell treatment, Relative expression	CCL3 / MIP-1 alpha Recombinant Rabbit mAb [PSH11-28] - BSA and Azide free	HA751407	Human, Mouse	WB,IF-Cell,IHC-P
	CCR2 Recombinant Rabbit mAb [SN707] - BSA and Azide free	HA750266	Human, Mouse	WB,IHC-P,IP,FC
KD, Relative expression	CCR7 Recombinant Rabbit mAb [SR36-04] - BSA and Azide free	HA750037	Human, Mouse, Rat	WB,IF-Cell,IP
	CCS Recombinant Rabbit mAb [PSH09-50] - BSA and Azide free	HA751297	Human, Mouse, Rat	WB,IHC-P,IP
Relative expression	CD10 Recombinant Mouse mAb [A1G3-R] - BSA and Azide free	HA610121	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	CD10 Recombinant Rabbit mAb [SN75-07] - BSA and Azide free	HA750273	Human, Mouse, Rat	WB,FC,IHC-P,IF-Cell,IF-Tissue
Relative expression	CD11b Recombinant Rabbit mAb [PSH03-96] - BSA and Azide free	HA750921	Human, Mouse, Rat	WB,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD11b Recombinant Rabbit mAb [JU93-81] - BSA and Azide free	HA750451	Human	WB,IP,IHC-P,IF-Cell,IF-Tissue
	CD11c Recombinant Antibody [SI19-06] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610258		
Relative expression	CD11c Recombinant Rabbit mAb [PSH07-42] - BSA and Azide free	HA751163	Mouse	WB,IHC-P,IHC-Fr,IF-Tissue
Relative expression	CD13 Recombinant Rabbit mAb [PSH08-48] - BSA and Azide free	HA751223	Human	WB,IF-Cell,IHC-P,FC
Cell treatment	CD133 Recombinant Rabbit mAb [PSH06-48] - BSA and Azide free	HA751072	Mouse	WB,IHC-P,IHC-Fr
Relative expression	CD133 Recombinant Mouse mAb [A8C3-R] - BSA and Azide free	HA610009	Human	WB,IHC-P
Relative expression	CD134 Recombinant Rabbit mAb [PSH10-14] - BSA and Azide free	HA751343	Human	WB,IF-Cell,FC,IP
KD	CD14 Recombinant Rabbit mAb [SC69-02] - BSA and Azide free	HA750232	Human, Mouse	WB,IF-Cell,IHC-P,FC
KD	CD14 Recombinant Antibody [SC69-02] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610260	0	
	CD15 Recombinant Rabbit mAb [PD00-42] - BSA and Azide free	HA750565	Human	IHC-P,FC,IF-Cell
	CD163 Recombinant Rabbit mAb [JA51-30] - BSA and Azide free	HA750412	Human	WB,IHC-P,IF-Tissue,IP
Relative expression	CD163 Recombinant Rabbit mAb [PSH15-06]	HA723706	Human	IHC-P,WB
KD	CD168 Recombinant Rabbit mAb [PSH02-16] - BSA and Azide free	HA750786	Human	WB,IF-Cell,IHC-P,FC
Cell treatment	CD169 Recombinant Rabbit mAb [PSH05-21] - BSA and Azide free	HA750992	Human, Mouse	WB,FC
Relative expression	CD19 Recombinant Rabbit mAb [PSH03-94] - BSA and Azide free	HA750919	Mouse	WB,IHC-P,IF-Tissue,IP,IHC-Fr,FC,IF-Cell
Relative expression	CD19 Recombinant Rabbit mAb [PSH13-44] - BSA and Azide free	HA751463	Human	WB,IF-Cell,FC
Relative expression	CD19 Recombinant Rabbit mAb [JF100-06] - BSA and Azide free	HA750365	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	CD2 Recombinant Rabbit mAb [PSH05-85] - BSA and Azide free	HA751033	Human	WB,IHC-P,IF-Cell,FC,IF-Tissue
Relative expression	CD20 Recombinant Rabbit mAb [PSH08-12] - BSA and Azide free	HA751207	Human	WB,IF-Cell,IHC-P,FC,IP
Cell treatment	CD204 Recombinant Mouse mAb [A8G2-R] - BSA and Azide free	HA610052	Human, Rat	WB,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment, Relative expression	CD24 Recombinant Rabbit mAb [PSH10-85] - BSA and Azide free	HA751377	Human	IF-Cell,FC,WB
Relative expression	CD27 Recombinant Rabbit mAb [PSH08-90] - BSA and Azide free	HA751255	Human	WB,IF-Cell,IHC-P,FC,IP
Relative expression, KD	CD276 Recombinant Rabbit mAb [PSH02-46] - BSA and Azide free	HA750811	Human	WB,IF-Cell,FC
	CD31 Recombinant Rabbit mAb [SU03-59] - BSA and Azide free	HA750151	Human	WB,IHC-P,FC,IP,IF-Cell,IF-Tissue
Relative expression	CD31 Recombinant Mouse mAb [7-A1-R] - BSA and Azide free	HA610115	Human	WB,IF-Cell,IHC-P
	CD31 Recombinant Antibody [7-A1-R] - Rat IgG1 (Chimeric)	HA601522	Human	IHC-P
Cell treatment	CD32 Recombinant Rabbit mAb [JM10-70] - BSA and Azide free	HA750382	Human	WB,IHC-P,IF-Cell,FC,IP
Relative expression, Cell treatment	CD33 Recombinant Rabbit mAb [PSH02-50] - BSA and Azide free	HA750813	Human	WB,IF-Cell,FC
Relative expression	CD33 Recombinant Rabbit mAb [PSH02-51] - BSA and Azide free	HA750814	Human	WB,IF-Cell,FC
	CD34 Recombinant Rabbit mAb [SI16-01] - BSA and Azide free	HA750092	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Relative expression	CD35 Recombinant Mouse mAb [A4G1-R] - BSA and Azide free	HA610125	Human	WB,IHC-P,IF-Tissue
Relative expression	CD35 Recombinant Mouse mAb [A4F3-R] - BSA and Azide free	HA610095	Human	WB,IF-Tissue,IHC-P
Cell treatment	CD36 Recombinant Rabbit mAb [JJ2005] - BSA and Azide free	HA750307	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Relative expression	CD38 Recombinant Antibody [PD01-49] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610257		
	CD39 Recombinant Rabbit mAb [JA90-36] - BSA and Azide free	HA750421	Human, Mouse, Rat	WB,IHC-P
	CD3 epsilon Recombinant Rabbit mAb [SY0239] - BSA and Azide free	HA750115	Human	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
	CD4 Recombinant Rabbit Multiclonal Antibody [PSH13-24] - BSA and Azide free	HA751460	Human, Mouse	WB,IF-Cell,IHC-P
Relative expression	CD4 Recombinant Rabbit mAb [PSH08-22] - BSA and Azide free	HA751213	Mouse	WB,IF-Cell,IHC-P,IHC-Fr
	CD40 Recombinant Rabbit mAb [PSH10-01] - BSA and Azide free	HA751331	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD40 Recombinant Rabbit mAb [PSH05-90] - BSA and Azide free	HA751037	Human, Cynomolgus monkey, Pig	WB,IF-Cell,FC
	CD40L Recombinant Rabbit mAb [JM11-34] - BSA and Azide free	HA750387	Human, Mouse, Rat	WB
	CD41 Recombinant Rabbit mAb [PSH08-03] - BSA and Azide free	HA751203	Human, Mouse, Rat	WB,IHC-P
	CD41 Recombinant Rabbit mAb [PSH08-02] - BSA and Azide free	HA751202	Mouse	IHC-P,FC,IF-Cell,Depletion
	CD41 Recombinant Rabbit mAb [PSH08-01] - BSA and Azide free	HA751201	Mouse	IHC-P,FC,IF-Cell,Depletion
	CD42a Recombinant Mouse mAb [A5E6-R] - BSA and Azide free	HA610137	Human, Mouse, Rat	WB,IHC-P
	CD42b Recombinant Mouse mAb [13G2-R] - BSA and Azide free	HA610202	Human	IHC-P
Relative expression	CD44 Recombinant Rabbit mAb [ST57-03] - BSA and Azide free	HA750192	Human	WB,IHC-P
Relative expression	CD45 Recombinant Rabbit mAb [PSH13-02] - BSA and Azide free	HA751457	Mouse, Rat	WB,IHC-Fr,IF-Tissue
Relative expression	CD45 Recombinant Rabbit Multiclonal Antibody [PSH14-19] - BSA and Azide free	HA751510	Human, Mouse	WB,IHC-P,IF-Cell
Relative expression	CD45 Recombinant Rabbit mAb [PSH13-01] - BSA and Azide free	HA751456	Mouse	WB,IHC-P,IHC-Fr,FC,IF-Cell,IF-Tissue
Relative expression	CD47 Recombinant Rabbit mAb [PSH10-31] - BSA and Azide free	HA751348	Human	WB,IF-Cell,FC
Relative expression	CD48 Recombinant Rabbit mAb [PSH10-65] - BSA and Azide free	HA751368	Human	WB,FC,IP,IHC-P
	CD58 Recombinant Rabbit mAb [PSH02-74] - BSA and Azide free	HA750829	Human	WB,FC
KD	CD63 Recombinant Rabbit mAb [SY21-02] - BSA and Azide free	HA750112	Human	WB,IHC-P
Relative expression	CD63 Recombinant Rabbit mAb [PSH06-64] - BSA and Azide free	HA751088	Human	WB,IF-Cell,FC,IP,IHC-P
	CD68 Recombinant Antibody [PDM0-13] - Rabbit IgG (Chimeric)	HA723708	Human	IF-Tissue,IHC-P
Relative expression, Cell treatment	CD68 Recombinant Rabbit mAb [PSH06-65] - BSA and Azide free	HA751094	Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,IHC-Fr
Relative expression	CD68 Recombinant Mouse mAb [A3C2-R] - BSA and Azide free	HA610164	Human	WB,IF-Cell,IHC-P
Cell treatment	CD68 Recombinant Rabbit mAb [PSH05-47] - BSA and Azide free	HA751006	Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD68 Recombinant Mouse mAb [A3C4-R] - BSA and Azide free	HA610165	Human	WB,IF-Cell,IHC-P
Relative expression	CD68 Recombinant Mouse mAb [PDM0-13] - BSA and Azide free	HA610034	Human	IHC-P,IF-Cell,FC,IF-Tissue
	CD68 Recombinant Antibody [PDM0-13] - Rat IgG1 (Chimeric)	HA601520	Human	IF-Tissue,IHC-P
Relative expression	CD73 Recombinant Mouse mAb [A6G2-R] - BSA and Azide free	HA610172	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Relative expression	CD73 Recombinant Rabbit mAb [PSH09-57] - BSA and Azide free	HA751304	Human	WB,IF-Cell,FC
Relative expression	CD80 Recombinant Rabbit mAb [PSH09-79] - BSA and Azide free	HA751321	Human	WB,IF-Cell,FC
	CD80 Recombinant Mouse mAb [8-E5-R] - BSA and Azide free	HA610122	Human	WB
	CD86 Recombinant Rabbit mAb [SJ20-00] - BSA and Azide free	HA750109	Human, Mouse, Rat	WB,IF-Cell,IP
	CD8 alpha Recombinant Rabbit mAb [SI18-01] - BSA and Azide free	HA750101	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	CD8 beta Recombinant Rabbit mAb [PSH08-98] - BSA and Azide free	HA751259	Human	WB,IF-Cell,IHC-P,FC
Relative expression, KD	CD9 Recombinant Rabbit mAb [PSH0-95] - BSA and Azide free	HA750675	Human	WB,IHC-P
	CD9 Recombinant Rabbit mAb [SA35-08] - BSA and Azide free	HA750027	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
KD	CD90 / THY1 Recombinant Rabbit mAb [JF10-09] - BSA and Azide free	HA750364	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
	CD90 / THY1 Recombinant Rabbit mAb [SU35-07] - BSA and Azide free	HA750149	Human, Rat	WB,IHC-P
	CDC42 Recombinant Rabbit mAb [JJ086-04] - BSA and Azide free	HA750320	Human, Mouse, Rat	WB,FC,IF-Cell
	CDK1 Recombinant Rabbit mAb [SY26-02] - BSA and Azide free	HA750125	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
KD	Cdk2 Recombinant Rabbit mAb [SR01-01] - BSA and Azide free	HA750054	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Cdk4 Recombinant Rabbit mAb [SD205-1] - BSA and Azide free	HA750286	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Cdk4 Recombinant Mouse mAb [5-C3-G3-R] - BSA and Azide free	HA610102	Human, Mouse, Rat	WB
	Cdk6 Recombinant Rabbit mAb [SD20-50] - BSA and Azide free	HA750289	Human	IF-Cell,IF-Tissue,IHC-P,FC,WB
Relative expression	CEACAM1 Recombinant Rabbit mAb [PSH0-78] - BSA and Azide free	HA750658	Human	WB,IHC-P
	C/EBP Beta Recombinant Rabbit mAb [SC55-05] - BSA and Azide free	HA750235	Human, Rat, Mouse	WB,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CELF1 Recombinant Mouse mAb [7A1-R] - BSA and Azide free	HA610199	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P
	Cell Cycle/Checkpoint Antibody Sampler Kit	HAK21137	Human	WB
	CENPA Recombinant Rabbit mAb [PSH0-09] - BSA and Azide free	HA750578	Human	WB,IF-Cell
	CENPE Recombinant Rabbit mAb [PSH03-54] - BSA and Azide free	HA750888	Human	WB
	CENPF Recombinant Rabbit mAb [PSH0-79] - BSA and Azide free	HA750659	Human, Mouse, Rat	WB,IHC-P
Cell treatment	c-Fos Recombinant Antibody [PSH05-77] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610219	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Tissue,IF-Cell,WB
Cell treatment	c-Fos Recombinant Rabbit mAb [PSH05-77] - BSA and Azide free	HA751071	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,IHC-Fr,IF-Tissue
	c-Fos Recombinant Rabbit mAb [PSH07-51] - BSA and Azide free	HA751172	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IF-Cell,IHC-P,ChIP,IF-Tissue,IP
Cell treatment	c-Fos Recombinant Antibody [PSH05-77] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610220	Mouse, Rat, Human, Cynomolgus monkey, Pig	IHC-Fr,WB
Cell treatment	c-Fos Recombinant Antibody [PSH07-51] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610275	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Cell,WB
Cell treatment	c-Fos Recombinant Rabbit mAb [JJ0938] - BSA and Azide free	HA750329	Human, Rat	WB,IF-Cell,IHC-P
Relative expression	CFTR Recombinant Rabbit mAb [PSH09-99] - BSA and Azide free	HA751330	Human	WB,IP
Relative expression	cGAS Recombinant Rabbit mAb [PSH06-78] - BSA and Azide free	HA751093	Human	WB
	CGBP Recombinant Rabbit mAb [PSH01-54] - BSA and Azide free	HA750730	Human, Mouse, Rat	WB,IHC-P
Relative expression	CGRP Recombinant Antibody [PSH08-45] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610263	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression, Cell treatment	CGRP Recombinant Rabbit mAb [PSH08-45] - BSA and Azide free	HA751220	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IHC-Fr
	CHCHD2 Recombinant Rabbit mAb [PSH13-87] - BSA and Azide free	HA751489	Human, Mouse, Rat	WB,IF-Cell,FC,IP
	CHD8 Recombinant Rabbit mAb [PSH11-71] - BSA and Azide free	HA751426	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,ChIP
	Chk1 Recombinant Rabbit mAb [ST57-09] - BSA and Azide free	HA750191	Human	WB,IF-Cell
KD	Chk2 Recombinant Rabbit mAb [SC604] - BSA and Azide free	HA750223	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Choline Acetyltransferase Recombinant Rabbit mAb [PSH15-05]	HA723705	Human, Mouse, Rat	WB,IHC-Fr
	Choline Acetyltransferase Recombinant Rabbit mAb [JA67-11] - BSA and Azide free	HA750406	Human, Mouse, Rat	WB,IP,IF-Cell,FC
	cIAP1 Recombinant Rabbit mAb [PSH06-72] - BSA and Azide free	HA751101	Human, Mouse, Rat	WB,IF-Cell,FC
KD	cIAP1 Recombinant Rabbit mAb [PSH02-10] - BSA and Azide free	HA750780	Human	WB
KD	cIAP1 Recombinant Rabbit mAb [PSH06-71] - BSA and Azide free	HA751100	Human, Mouse, Rat	WB
	cIAP2 Recombinant Rabbit mAb [PSH06-92] - BSA and Azide free	HA751116	Human	WB
Relative expression	Cingulin Recombinant Rabbit mAb [PSH06-19] - BSA and Azide free	HA751057	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,FC
	c-Jun Recombinant Rabbit mAb [SY0298] - BSA and Azide free	HA750143	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
KD	Claudin 1 Recombinant Rabbit mAb [PS01-43] - BSA and Azide free	HA750567	Human	WB,IHC-P,IF-Tissue,IF-Cell
Relative expression	Claudin 18.2 Recombinant Rabbit mAb [PSH11-63] - BSA and Azide free	HA751419	Human, Mouse, Rat	IHC-P
	Claudin 5 Recombinant Rabbit mAb [JM11-22] - BSA and Azide free	HA750384	Human	WB,IHC-P
	Claudin 7 Recombinant Mouse mAb [6-G3-R] - BSA and Azide free	HA610206	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,IP
Cell treatment	Cleaved Caspase-3 Recombinant Rabbit mAb [SR01-02] - BSA and Azide free	HA750053	Human, Pig	WB,IF-Cell,IHC-P
Cell treatment, KD	Cleaved PARP Recombinant Rabbit mAb [SU0314] - BSA and Azide free	HA750136	Human	WB,IF-Cell,IP
Cell treatment, KO	Cleaved+pro Caspase-3 Recombinant Rabbit mAb [SU38-04] - BSA and Azide free	HA750159	Human, Mouse, Rat	WB,IHC-P,IP
Relative expression	CLP Recombinant Rabbit mAb [PSH01-73] - BSA and Azide free	HA750748	Human, Mouse, Rat	IHC,IF-Tissue
Relative expression	CLP Recombinant Rabbit mAb [PSH01-72] - BSA and Azide free	HA750747	Human, Mouse, Rat	IHC,IF-Tissue
Relative expression	c-Met Recombinant Mouse mAb [A9A4-R] - BSA and Azide free	HA610037	Human	WB,IHC-P
Relative expression	c-Met Recombinant Rabbit mAb [PSH0-19] - BSA and Azide free	HA750594	Human	WB,IHC-P,FC
	CMTM6 Recombinant Rabbit mAb [PSH15-02]	HA723702	Human, Monkey	WB,IP
	c-Myc Recombinant Rabbit mAb [PS00-11] - BSA and Azide free	HA750519	Human, Mouse, Rat	WB,IHC-P,IP,ChIP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CNOT7 Recombinant Rabbit mAb [PSH02-81] - BSA and Azide free	HA750836	Human, Mouse, Rat	WB,IF-Cell,FC
	Cofilin Recombinant Rabbit mAb [PSH01-44] - BSA and Azide free	HA750721	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	Cofilin Recombinant Rabbit mAb [PSH0-33] - BSA and Azide free	HA750612	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
	Cofilin Recombinant Rabbit mAb [PSH0-34] - BSA and Azide free	HA750613	Human, Mouse, Rat, Monkey	WB,FC
	Cofilin Recombinant Rabbit mAb [PSH0-38] - BSA and Azide free	HA750616	Human, Mouse, Rat, Monkey	WB,FC
	COL1A1 Recombinant Rabbit mAb [ST58-04] - BSA and Azide free	HA750188	Human, Cow	WB,IHC-P
	COL1A1 Recombinant Rabbit mAb [PSH06-20] - BSA and Azide free	HA751058	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP,IHC-Fr
Relative expression	COL1A2 Recombinant Rabbit mAb [PSH07-21] - BSA and Azide free	HA751142	Human, Rat	WB,IF-Cell,FC
	COL5A1 Recombinant Rabbit mAb [PSH10-08] - BSA and Azide free	HA751338	Human	WB,IF-Cell,FC,IP
Relative expression	Collagen II Recombinant Rabbit mAb [PSH06-62] - BSA and Azide free	HA751090	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
	Collagen VI Recombinant Rabbit mAb [SD83-03] - BSA and Azide free	HA750305	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Collagen X Recombinant Rabbit mAb [JF0961] - BSA and Azide free	HA750336	Rat, Human	WB
KD	Collagen XVII Recombinant Rabbit mAb [SR46-05] - BSA and Azide free	HA750030	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
	COMMD10 Recombinant Rabbit mAb [PSH09-13] - BSA and Azide free	HA751267	Human	WB,IHC-P,IP
	Complement C3 Recombinant Rabbit mAb [JF10-30] - BSA and Azide free	HA750367	Human, Mouse, Rat	WB,IF-Cell,FC
	Complement C5 Recombinant Rabbit mAb [PSH04-01] - BSA and Azide free	HA750924	Human	WB,IF-Cell,IHC-P,FC
	Connexin 43 Recombinant Rabbit mAb [PSH0-32] - BSA and Azide free	HA750606	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
KD	COX1/Cyclooxygenase 1 Recombinant Rabbit mAb [SC68-05] - BSA and Azide free	HA750241	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
KD	COX2/Cyclooxygenase 2 Recombinant Rabbit mAb [SC56-06] - BSA and Azide free	HA750207	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	COX5A Recombinant Rabbit mAb [PSH13-57] - BSA and Azide free	HA751476	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	COX IV Recombinant Rabbit mAb [JJ09-05] - BSA and Azide free	HA750317	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	CPT1A Recombinant Rabbit mAb [PSH10-94] - BSA and Azide free	HA751384	Mouse, Rat	WB
	CPT1A Recombinant Rabbit mAb [PSH04-00] - BSA and Azide free	HA750923	Human, Mouse	WB,IHC-P,IP
KD	CPT2 Recombinant Rabbit mAb [SN06-70] - BSA and Azide free	HA750265	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	CRALBP Recombinant Rabbit mAb [PSH0-74] - BSA and Azide free	HA750653	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
KD	CRBN Recombinant Rabbit mAb [PSH0-69] - BSA and Azide free	HA750646	Human, Mouse, Rat	WB,IHC-P,IP
KD	CRBN Recombinant Rabbit mAb [PSH0-70] - BSA and Azide free	HA750647	Human, Mouse, Rat	WB,IHC-P
	Creatine kinase B type Recombinant Rabbit mAb [PSH09-32] - BSA and Azide free	HA751286	Human, Mouse, Rat	WB,IF-Cell,FC
	CREB Recombinant Rabbit mAb [SA04-04] - BSA and Azide free	HA750007	Human, Mouse, Zebrafish, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	CREBBP Recombinant Rabbit mAb [PS01-16] - BSA and Azide free	HA751014	Human, Mouse	IHC-P
Cell treatment	Cre recombinase Recombinant Rabbit mAb [PSH01-93] - BSA and Azide free	HA750767	Species independent	WB,IHC-P,IF-Tissue,IF-Cell,FC
Cell treatment	CRISPR-Cas9 SP Recombinant Rabbit mAb [JM11-55] - BSA and Azide free	HA750393	Streptococcus pyogenes	WB,IF-Cell,IF-Tissue,IHC-P,FC
	CRLF3 Recombinant Rabbit mAb [PSH03-88] - BSA and Azide free	HA750913	Human, Mouse, Rat	WB,IF-Cell
Relative expression	CSF-1-R Recombinant Rabbit mAb [PSH06-70] - BSA and Azide free	HA751099	Mouse	WB,IF-Cell
Relative expression	CSF-1-R Recombinant Rabbit mAb [PSH09-77] - BSA and Azide free	HA751319	Human, Mouse, Rat	WB,IHC-P
	CSN7b Recombinant Rabbit mAb [PSH06-84] - BSA and Azide free	HA751107	Human, Mouse, Rat	WB
	CSNK1A1 + CSNK1A1L Recombinant Rabbit mAb [PSH04-58] - BSA and Azide free	HA750965	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	CSNK2A1 Recombinant Rabbit mAb [PSH04-54] - BSA and Azide free	HA750961	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	CSPP1 Recombinant Rabbit mAb [PSH06-31] - BSA and Azide free	HA751068	Human	WB,IF-Cell
	CSPP1 Recombinant Rabbit mAb [PSH06-32] - BSA and Azide free	HA751069	Human, Mouse, Rat	WB
	CTCF Recombinant Rabbit mAb [JM10-61] - BSA and Azide free	HA750394	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,CUT&Tag-seq,ChIP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Ctip1 Recombinant Rabbit mAb [JB35-68] - BSA and Azide free	HA750461	Human, Mouse	WB,IHC-P,FC,IF-Cell,IF-Tissue
	Ctip2 Recombinant Rabbit mAb [JF09-90] - BSA and Azide free	HA750358	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell
	CTLA-4 Recombinant Rabbit mAb [PS01-34] - BSA and Azide free	HA750581	Human, Mouse	IHC-P
	CTU2 Recombinant Rabbit mAb [PSH04-05] - BSA and Azide free	HA750928	Human, Mouse, Monkey	WB,IHC-P
	Cx40 / GJA5 Recombinant Rabbit mAb [PSH10-82] - BSA and Azide free	HA751374	Human, Mouse, Rat, Monkey	WB,IP
Cell treatment, Relative expression	CXCL13 Recombinant Rabbit mAb [PSH04-34] - BSA and Azide free	HA750943	Human	WB,IF-Cell,IHC-P
Cell treatment	CXCL1 / GRO alpha Recombinant Rabbit mAb [PSH07-53] - BSA and Azide free	HA751174	Human	WB
Cell treatment	CXCL2 Recombinant Rabbit mAb [PSH14-44] - BSA and Azide free	HA751523	Mouse	WB,IF-Cell,IHC-P,FC
Cell treatment	CXCL9 Recombinant Rabbit mAb [PSH12-34] - BSA and Azide free	HA751512	Human	WB,IF-Cell
	CXCR3 Recombinant Rabbit mAb [JA61-33] - BSA and Azide free	HA750422	Human	WB
	CXCR5 Recombinant Rabbit mAb [JB11-40] - BSA and Azide free	HA750463	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
Cell treatment	Cyclin A2 Recombinant Rabbit mAb [SD2052] - BSA and Azide free	HA750288	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
KD	Cyclin B1 Recombinant Rabbit mAb [SU33-03] - BSA and Azide free	HA750141	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Cyclin E2 Recombinant Mouse mAb [40-89-R] - BSA and Azide free	HA610048	Human	WB,IHC-P,IF-Cell
	Cyclin T1 Recombinant Rabbit mAb [PSH09-66] - BSA and Azide free	HA751308	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP,ChIP
	Cyclophilin 40 Recombinant Rabbit mAb [PSH02-71] - BSA and Azide free	HA750828	Human, Monkey	WB,IHC-P,FC,IF-Cell
	CYP11A1 Recombinant Rabbit mAb [PSH07-26] - BSA and Azide free	HA751147	Human, Mouse, Rat	WB
Relative expression, Cell treatment	CYR61 / CCN1 Recombinant Rabbit mAb [PSH09-12] - BSA and Azide free	HA751266	Human	WB,IF-Cell,IHC-P,IF-Tissue
KD	Cystathionase Recombinant Rabbit mAb [JE35-91] - BSA and Azide free	HA750649	Human, Mouse, Rat	WB,IP
Relative expression	Cystatin F Recombinant Rabbit mAb [PSH09-00] - BSA and Azide free	HA751261	Human, Mouse	WB,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Cytochrome C Recombinant Rabbit mAb [SC59-01] - BSA and Azide free	HA750202	Human, Mouse, Rat	WB,IHC-P
	Cytochrome C Recombinant Rabbit mAb [SC69-08] - BSA and Azide free	HA750226	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP,IF-Cell
	Cytokeratin 1 Recombinant Rabbit mAb [SN72-08] - BSA and Azide free	HA750256	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
KD	Cytokeratin 10 Recombinant Rabbit mAb [ST05-43] - BSA and Azide free	HA750193	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Relative expression	Cytokeratin 13 Recombinant Rabbit mAb [SN71-09] - BSA and Azide free	HA750261	Human, Mouse, Zebrafish	WB,IF-Cell,IHC-P,FC
Relative expression	Cytokeratin 14 Recombinant Rabbit mAb [SC65-06] - BSA and Azide free	HA750216	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Cytokeratin 15 Recombinant Rabbit mAb [ST04-85] - BSA and Azide free	HA750182	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell
	Cytokeratin 16 Recombinant Rabbit mAb [SC52-09] - BSA and Azide free	HA750203	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Cytokeratin 17 Recombinant Rabbit mAb [SR45-06] - BSA and Azide free	HA750031	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Relative expression	Cytokeratin 17 Recombinant Mouse mAb [A2B10-R] - BSA and Azide free	HA610110	Human	WB,IF-Cell,IHC-P
KD	Cytokeratin 18 Recombinant Rabbit mAb [SZ80-07] - BSA and Azide free	HA750069	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Cytokeratin 18 Recombinant Mouse mAb [A4B3-R] - BSA and Azide free	HA610135	Human, Mouse, Rat	WB,IHC-P
	Cytokeratin 18 Recombinant Mouse mAb [6-19-R] - BSA and Azide free	HA610144	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	Cytokeratin 19 Recombinant Mouse mAb [A3D1-R] - BSA and Azide free	HA610133	Human	WB,IF-Cell,IHC-P,FC
	Cytokeratin 20 Recombinant Rabbit mAb [SA35-03] - BSA and Azide free	HA750026	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	Cytokeratin 4 Recombinant Rabbit mAb [SN74-03] - BSA and Azide free	HA750271	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	Cytokeratin 5+6 Recombinant Mouse mAb [A2B2-R] - BSA and Azide free	HA610108	Human, Rat	WB,IHC-P
	Cytokeratin 5+6 Recombinant Rabbit mAb [PSH0-68] - BSA and Azide free	HA750645	Human	WB,IF-Cell,FC
	Cytokeratin 6 Recombinant Rabbit mAb [SN71-07] - BSA and Azide free	HA750270	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD, Relative expression	Cytokeratin 7 Recombinant Rabbit mAb [ST50-05] - BSA and Azide free	HA750185	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	Cytokeratin 8 Recombinant Rabbit mAb [SU0338] - BSA and Azide free	HA750145	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	Cytokeratin 9 Recombinant Rabbit mAb [SD201-09] - BSA and Azide free	HA750302	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	DAG1 Recombinant Rabbit mAb [PSH13-91] - BSA and Azide free	HA751493	Human, Mouse, Rat	WB
	DAP Kinase 1 Recombinant Mouse mAb [12B2-R] - BSA and Azide free	HA610051	Human	WB,IHC-P
	DARPP32 Recombinant Rabbit mAb [SU0329] - BSA and Azide free	HA750139	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	DCAMKL1 Recombinant Rabbit mAb [JA11-03] - BSA and Azide free	HA750403	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	DCC Recombinant Rabbit mAb [PSH12-73] - BSA and Azide free	HA751449	Human, Mouse, Rat	WB,IHC-Fr
	DCI / Eci1 Rabbit pAb	HA500581	Mouse, Rat	WB,IF-Cell,FC
	DCTN1 / p150-glued Recombinant Rabbit mAb [PSH13-45] - BSA and Azide free	HA751464	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Cell,FC,IF-Tissue,IP
Cell treatment	DDIT3 Recombinant Rabbit mAb [PSH07-65] - BSA and Azide free	HA751178	Human, Mouse, Rat	WB,IF-Cell,FC
KD	DDIT3 Recombinant Rabbit mAb [JM10-31] - BSA and Azide free	HA750369	Human, Mouse, Rat	WB,FC
Cell treatment	DDIT4 Recombinant Rabbit mAb [PS00-10] - BSA and Azide free	HA750969	Human	WB
	ddx24 Recombinant Mouse mAb [A3B10-R] - BSA and Azide free	HA610066	Zebrafish	IHC-P,IF-Tissue
KD	DDX39A Recombinant Rabbit mAb [PSH05-57] - BSA and Azide free	HA751020	Human, Mouse	WB,IF-Cell,IHC-P,IP,IF-Tissue
Relative expression	DDX4 Recombinant Rabbit mAb [PSH08-53] - BSA and Azide free	HA751228	Human, Mouse	WB,IHC-P,IHC-Fr
Relative expression	DDX4 Recombinant Rabbit mAb [PSH08-54] - BSA and Azide free	HA751229	Human, Mouse, Rat	WB,IHC-P,IHC-Fr
	Decorin Recombinant Rabbit mAb [PSH12-99] - BSA and Azide free	HA751454	Human, Mouse, Rat	WB,IHC-P
	Desmin Recombinant Rabbit mAb [SI18-00] - BSA and Azide free	HA750100	Human, Mouse, Rat, Zebrafish, Chicken	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment, Relative expression	DFNA5 / GSDME Recombinant Rabbit mAb [PSH10-83] - BSA and Azide free	HA751375	Human	WB
Cell treatment, Relative expression	DFNA5 / GSDME Recombinant Rabbit mAb [PSH10-84] - BSA and Azide free	HA751376	Human, Mouse	WB,IP
	DFNA5 / GSDME Recombinant Rabbit mAb [PSH07-48] - BSA and Azide free	HA751169	Human, Mouse, Rat	WB
Cell treatment, Relative expression	DFNA5 / GSDME Recombinant Rabbit mAb [PSH09-51] - BSA and Azide free	HA751298	Human, Mouse	WB,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	DHODH Recombinant Rabbit mAb [PSH09-69] - BSA and Azide free	HA751311	Human, Mouse, Rat	WB,IHC-P
Relative expression	Digoxigenin Recombinant Rabbit mAb [PSH04-39] - BSA and Azide free	HA750948	Species independent	WB,IHC-P,IHC-Fr,Dot Blot,ELISA,IF-Tissue,IP
Relative expression	Digoxin Recombinant Rabbit mAb [PSH04-40] - BSA and Azide free	HA750949	Species independent	WB,ELISA,IHC-P,IHC-Fr,Dot Blot,IF-Tissue,IP
	Dishevelled 2 Recombinant Rabbit mAb [PSH04-07] - BSA and Azide free	HA750930	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	DKK1 Recombinant Rabbit mAb [SC06-86] - BSA and Azide free	HA750228	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	DLK-1 Recombinant Rabbit mAb [PSH04-76] - BSA and Azide free	HA750970	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	DLST Recombinant Rabbit mAb [PSH10-07] - BSA and Azide free	HA751337	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P
	Dlx2 Recombinant Rabbit mAb [PSH13-47] - BSA and Azide free	HA751466	Mouse, Rat	IHC-P
	Dlx2 Recombinant Rabbit mAb [PSH13-46] - BSA and Azide free	HA751465	Mouse, Rat	WB
	DMT1 Recombinant Rabbit mAb [PSH09-20] - BSA and Azide free	HA751274	Human	WB,IHC-P,FC
	DNAJA1 Recombinant Rabbit mAb [PSH07-73] - BSA and Azide free	HA751185	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP
Relative expression	DNAJC6 Recombinant Rabbit mAb [PSH07-95] - BSA and Azide free	HA751197	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	DNA PKcs Recombinant Rabbit mAb [SC57-08] - BSA and Azide free	HA750198	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	DNA-PKcs/PRKDC Recombinant Mouse mAb [2-B3-R] - BSA and Azide free	HA610107	Human	WB,IHC-P,IF-Cell
	DNA Polymerase beta Recombinant Rabbit mAb [JM93-12] - BSA and Azide free	HA750423	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IF-Cell,IP
	DNA Polymerase lambda Recombinant Rabbit mAb [PSH05-24] - BSA and Azide free	HA750995	Human, Mouse, Monkey	WB,IHC-P
	DNA polymerase mu Recombinant Rabbit mAb [PSH08-75] - BSA and Azide free	HA751240	Human	WB
	Dnmt1 Recombinant Rabbit mAb [JF09-89] - BSA and Azide free	HA750359	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	DNMT1 Recombinant Mouse mAb [A3A8-R] - BSA and Azide free	HA610156	Human	WB,IF-Cell,IHC-P
KD	Dnmt2 Recombinant Rabbit mAb [PSH03-43] - BSA and Azide free	HA750877	Human, Monkey	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	Dnmt3a Recombinant Rabbit mAb [ST04-78] - BSA and Azide free	HA750169	Human, Mouse, Rat	WB,IF-Tissue,IHC-P
Relative expression	Dopamine D2 Receptor Recombinant Rabbit mAb [PSH09-94] - BSA and Azide free	HA751325	Human, Mouse, Rat	IHC-Fr,IHC-P
	Dopamine D2 Receptor Recombinant Rabbit mAb [PSH09-93] - BSA and Azide free	HA751324	Human, Mouse	WB,IHC-Fr,IF-Cell,FC
	Dopamine Transporter Recombinant Antibody [PSH10-49] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610244	Mouse, Rat, Human, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
	Dopamine Transporter Recombinant Antibody [PSH10-49] - Chicken IgY (Chimeric) - BSA and Azide free	HA610308	Mouse	IHC-Fr
	Dopamine Transporter Recombinant Rabbit mAb [PSH10-49] - BSA and Azide free	HA751352	Mouse, Rat, Human, Cynomolgus monkey, Pig	WB,IHC-Fr,IHC-P
	Dopamine Transporter Recombinant Antibody [PSH10-49] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610309	Mouse	IHC-Fr
	Dopamine Transporter Recombinant Antibody [PSH10-49] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610310	Mouse, Rat	IHC-Fr,IHC-P
	DOT1L Recombinant Rabbit mAb [PSH01-57] - BSA and Azide free	HA750733	Human, Mouse, Rat	WB,IHC-P
Relative expression	Doublecortin Recombinant Antibody [JJ0959] - Chicken IgY (Chimeric)	HA601495	Human, Mouse, Rat	IHC-Fr,WB
Relative expression	Doublecortin Recombinant Antibody [JJ0959] - Mouse IgG1 (Chimeric)	HA601497	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Doublecortin Recombinant Antibody [JJ0959] - Guinea pig IgG2 (Chimeric)	HA601496	Human, Mouse, Rat	IHC-Fr,WB
Relative expression	Doublecortin Recombinant Antibody [JJ0959] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610238	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
Cell treatment	DPP4 / CD26 Recombinant Rabbit mAb [PSH12-59] - BSA and Azide free	HA751446	Human	WB,IHC-P,IF-Cell,FC,IP
KO	DPP4 / CD26 Recombinant Rabbit mAb [JM11-42] - BSA and Azide free	HA750396	Human, Mouse, Rat	WB
	DPTOR Recombinant Rabbit mAb [PH00-03] - BSA and Azide free	HA750555	Human, Mouse, Rat	WB,IHC-P
	DR3 / LARD Recombinant Rabbit mAb [PSH13-90] - BSA and Azide free	HA751492	Human, Mouse, Rat	WB,IF-Cell
	DUSP1 Recombinant Rabbit mAb [JJ0930] - BSA and Azide free	HA750324	Human, Mouse, Rat	WB,IHC-P,IP,FC,IF-Cell,IF-Tissue
	DUSP6 Recombinant Rabbit mAb [SR39-09] - BSA and Azide free	HA750032	Human, Mouse, Rat, Chicken	WB,IF-Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	DXD Recombinant Rabbit mAb [PSH0-73] - BSA and Azide free	HA750652	Species independent	ELISA
	DYKDDDDK Tag (FLAG) Recombinant Mouse mAb [PS00-M2] - BSA and Azide free	HA610017	Species independent	WB,IP,IF-Cell
	DYKDDDDK Tag (FLAG) Recombinant Rabbit mAb [PSH07-02] - BSA and Azide free	HA751129	Species independent	WB,IHC-P,IF-Cell,FC,IP,ChIP
	DYKDDDDK Tag (FLAG) Recombinant Mouse mAb [A2-A4-R] - BSA and Azide free	HA610054	Species independent	WB,IP,IF-Cell,IHC-P
	Dynamin 1 Recombinant Mouse mAb [A1B2-R] - BSA and Azide free	HA610197	Human, Mouse, Rat	WB,IHC-P
	Dynamin 1 Recombinant Mouse mAb [A1B1-R] - BSA and Azide free	HA610105	Mouse, Rat	WB
	E2F1 Recombinant Rabbit mAb [JJ092-02] - BSA and Azide free	HA750322	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	E2F4 Recombinant Rabbit mAb [PSH0-45] - BSA and Azide free	HA750625	Human	WB
Cell treatment	EAAT1 Recombinant Rabbit mAb [JA30-35] - BSA and Azide free	HA750415	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
	EAAT2 Recombinant Antibody [PS01-62] - Chicken IgY (Chimeric)	HA601503	Mouse	IHC-Fr
Relative expression	EAAT2 Recombinant Antibody [PS01-62] - Guinea pig IgG2 (Chimeric)	HA601504	Mouse, Rat	IHC-Fr,WB
Relative expression	EAAT2 Recombinant Antibody [PS01-62] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610243	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
Relative expression	EAAT2 Recombinant Rabbit mAb [PSH07-43] - BSA and Azide free	HA751164	Human, Mouse, Rat	WB,IHC-P,IHC-Fr
Relative expression	EAAT2 Recombinant Rabbit mAb [PS01-62] - BSA and Azide free	HA750609	Mouse, Rat, Human	WB,IHC-P,IHC-Fr
Relative expression, KD	E-Cadherin Recombinant Rabbit mAb [SY0287] - BSA and Azide free	HA750133	Human, Cynomolgus monkey	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	E-Cadherin Recombinant Mouse mAb [A0-G11-2-R] - BSA and Azide free	HA610047	Human, Mouse, Rat	WB,IHC
Relative expression	E-Cadherin Recombinant Rabbit mAb [PSH13-75] - BSA and Azide free	HA751488	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Cell,FC,IP
KD	EDG1 Recombinant Rabbit mAb [JM10-66] - BSA and Azide free	HA750376	Human, Mouse, Rat	WB,IP
KD	EEF2K Recombinant Rabbit mAb [PSH11-67] - BSA and Azide free	HA751423	Human, Mouse, Rat	WB,IF-Cell,FC
	EFTUD2 Recombinant Mouse mAb [1-D10-A1-R] - BSA and Azide free	HA610094	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression, KD	EGFR Recombinant Rabbit mAb [PD00-75] - BSA and Azide free	HA750550	Human, Monkey	WB,IHC-P,IF-Cell,FC
Relative expression, KO	EGFR Recombinant Rabbit mAb [SP00-86] - BSA and Azide free	HA750078	Human, Mouse, Rat, Monkey	WB,IF-Cell,IP,FC
	EIF2S1 Recombinant Rabbit mAb [PSH04-29] - BSA and Azide free	HA750938	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP
	eIF3g Recombinant Rabbit mAb [PSH09-42] - BSA and Azide free	HA751289	Human, Mouse, Rat	WB
	eIF3g Recombinant Rabbit mAb [PSH09-41] - BSA and Azide free	HA751288	Human	WB,IF-Cell,FC,IP
	EIF3H Recombinant Mouse mAb [PSH04-49] - BSA and Azide free	HA610200	Human, Mouse, Rat	WB,IHC-P
	EIF3S1 / EIF3J Recombinant Rabbit mAb [PSH06-97] - BSA and Azide free	HA751124	Human, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	eIF4A2 Recombinant Rabbit mAb [PSH01-89] - BSA and Azide free	HA750763	Human, Mouse, Rat, Monkey	WB,IP
	eIF4A2 Recombinant Rabbit mAb [PSH01-88] - BSA and Azide free	HA750762	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	eIF4G1 Recombinant Rabbit mAb [PSH06-28] - BSA and Azide free	HA751065	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	eIF-6 Recombinant Mouse mAb [A6A2-R] - BSA and Azide free	HA610167	Human, Mouse, Rat	WB
	eIF-6 Recombinant Mouse mAb [A6A3-R] - BSA and Azide free	HA610166	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Emerin Recombinant Rabbit mAb [PSH01-85] - BSA and Azide free	HA750759	Human	WB,IHC-P,IF-Cell,FC
Relative expression, KD	Emerin Recombinant Rabbit mAb [PSH01-86] - BSA and Azide free	HA750760	Human	WB,IHC-P,IF-Cell,FC
	Engrailed 1 Recombinant Rabbit mAb [PSH08-17] - BSA and Azide free	HA751208	Human, Mouse, Rat	WB,FC
Relative expression	ENPP1 / PC1 Recombinant Rabbit mAb [PSH12-11] - BSA and Azide free	HA751440	Human	WB,IHC
Cell treatment, Relative expression	ENPP3 / B10 Recombinant Rabbit mAb [PSH06-96] - BSA and Azide free	HA751123	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	EpCAM Recombinant Rabbit mAb [PS01-69] - BSA and Azide free	HA750703	Human	IHC-P,IF-Cell,FC
Relative expression	EPCR Recombinant Rabbit mAb [PSH03-50] - BSA and Azide free	HA750884	Human	WB,IF-Cell,FC
	Eph receptor A2 Recombinant Rabbit mAb [PSH09-21] - BSA and Azide free	HA751275	Human	WB,IF-Cell,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	Eph receptor B2 Recombinant Rabbit mAb [PSH14-91]	HA723693	Human	WB,IF-Cell,IP
	Ephrin A1 Recombinant Mouse mAb [A7B10-R] - BSA and Azide free	HA610173	Human	WB
	Ephrin B2 Recombinant Rabbit mAb [JM53-21] - BSA and Azide free	HA750428	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	EPO Recombinant Rabbit mAb [PSH07-33] - BSA and Azide free	HA751154	Human	WB
	ERCC1 Recombinant Rabbit mAb [JM10-07] - BSA and Azide free	HA750374	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
KD	ERCC8 Recombinant Rabbit mAb [PSH0-86] - BSA and Azide free	HA750666	Human	WB
	ERCC8 Recombinant Rabbit mAb [PSH0-87] - BSA and Azide free	HA750667	Human	WB
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03] - BSA and Azide free	HA750019	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
KD	ERK2 Recombinant Mouse mAb [C3B12-R] - BSA and Azide free	HA610079	Human, Mouse, Rat	WB
KD	ERK2 Recombinant Mouse mAb [5-D2-R] - BSA and Azide free	HA610080	Human, Mouse, Rat	WB
KD	ERK5 Recombinant Rabbit mAb [SD2084] - BSA and Azide free	HA750299	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IP,FC
	ERK8 Recombinant Rabbit mAb [PSH07-23] - BSA and Azide free	HA751144	Human, Mouse, Rat	WB,IHC-P
	ERM / Etf5 Recombinant Rabbit mAb [PSH02-64] - BSA and Azide free	HA750821	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	ERO1L Recombinant Rabbit mAb [PSH07-30] - BSA and Azide free	HA751151	Human, Mouse, Rat	WB,IHC-P
KD	ESD Recombinant Mouse mAb [J0-8-R] - BSA and Azide free	HA610178	Human, Mouse, Rat	WB
	Estrogen Related Receptor alpha Recombinant Rabbit mAb [SC57-07] - BSA and Azide free	HA750200	Human, Mouse, Rat	WB,IP,IHC-P,IF-Tissue
	EXOSC9 Recombinant Rabbit mAb [PSH08-55] - BSA and Azide free	HA751230	Human, Mouse, Rat	WB
	EZH1 Recombinant Rabbit mAb [PS01-09] - BSA and Azide free	HA750521	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Relative expression	F4/80 Recombinant Rabbit mAb [PSH0-85] - BSA and Azide free	HA750665	Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression	F4/80 Recombinant Rabbit mAb [PSH01-87] - BSA and Azide free	HA750761	Mouse	WB,IHC-P,IF-Tissue,IHC-Fr
	F4/80 Recombinant Antibody [PSH0-85] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610319	Mouse, Rat	IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	FABP4 Recombinant Rabbit mAb [JM10-99] - BSA and Azide free	HA750398	Human, Mouse, Rat	WB,IF-Cell,IHC-Fr
	FADD Recombinant Rabbit mAb [PSH01-13] - BSA and Azide free	HA750689	Human	WB,IF-Cell,FC
	FADD Recombinant Rabbit mAb [PSH07-78] - BSA and Azide free	HA751190	Human, Mouse, Rat	WB,IHC-P
KD	FAK Recombinant Rabbit mAb [SR46-04] - BSA and Azide free	HA750039	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	FAM134B Recombinant Rabbit mAb [PSH01-94] - BSA and Azide free	HA750768	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	FAM38A / PIEZO1 Recombinant Rabbit mAb [PSH11-49] - BSA and Azide free	HA751411	Mouse	WB,IF-Cell
	FAM38A / PIEZO1 Recombinant Mouse mAb [2-10-R-A] - BSA and Azide free	HA610024	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	FAM38B / PIEZO2 Recombinant Rabbit mAb [PSH11-50] - BSA and Azide free	HA751412	Mouse	WB,IF-Cell
	FAP Recombinant Rabbit mAb [JA56-11] - BSA and Azide free	HA750409	Human, Mouse, Rat	WB,IHC-P,IHC-Fr
Relative expression	Fast Myosin Skeletal Heavy chain Recombinant Rabbit mAb [PSH05-12] - BSA and Azide free	HA750985	Human, Mouse, Rat	IHC-P,IF-Tissue
KO	Fatty Acid Synthase Recombinant Rabbit mAb [JJ0939] - BSA and Azide free	HA750327	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	FBXO3 Recombinant Rabbit mAb [PSH04-09] - BSA and Azide free	HA750932	Human, Mouse	WB
	Fbxw7 Recombinant Rabbit mAb [PD00-51] - BSA and Azide free	HA750530	Human	IHC-P
	Ferritin Light Chain Recombinant Rabbit mAb [JM10-37] - BSA and Azide free	HA750435	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	FGF18 Recombinant Rabbit mAb [PSH01-24] - BSA and Azide free	HA750700	Human, Mouse	WB
Cell treatment	FGF18 Recombinant Rabbit mAb [PSH01-25] - BSA and Azide free	HA750701	Mouse	WB
	FGF21 Recombinant Rabbit mAb [JA10-97] - BSA and Azide free	HA750399	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	FGFR2/CD332 Recombinant Mouse mAb [A5-E10-G7-R] - BSA and Azide free	HA610003	Human	WB
KD	FGFR3 Recombinant Rabbit mAb [JM81-10] - BSA and Azide free	HA750379	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	FGFR3 Recombinant Rabbit mAb [JM110-33] - BSA and Azide free	HA750388	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Fibronectin Recombinant Rabbit mAb [PSH14-21] - BSA and Azide free	HA751514	Human, Mouse	WB,IHC-P,IF-Cell,FC
	Fibronectin Recombinant Rabbit mAb [PSH14-22] - BSA and Azide free	HA751515	Human, Mouse	WB,IF-Cell,FC
KD	Fibronectin Recombinant Rabbit mAb [JF0582] - BSA and Azide free	HA750340	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
	Filamin B Recombinant Rabbit mAb [PSH07-34] - BSA and Azide free	HA751155	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
KD	FLCN Recombinant Rabbit mAb [PSH07-68] - BSA and Azide free	HA751180	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	FLIP Recombinant Rabbit mAb [PSH01-95] - BSA and Azide free	HA750769	Human	WB
KD	Flotillin-2 Recombinant Rabbit mAb [PSH03-41] - BSA and Azide free	HA750875	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
Relative expression	FOLR1 Recombinant Rabbit mAb [PSH06-93] - BSA and Azide free	HA751117	Human	WB,IHC-P,FC
	FOXA1 Recombinant Mouse mAb [A2E8-R] - BSA and Azide free	HA610171	Human	WB,IF-Cell,IHC-P
	FOXA1 Recombinant Rabbit mAb [JF10-02] - BSA and Azide free	HA750363	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	FOXC1 Recombinant Rabbit mAb [PSH03-89] - BSA and Azide free	HA750914	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	FOXG1 Recombinant Rabbit mAb [PSH09-08] - BSA and Azide free	HA751263	Mouse	WB,IP
	FOXM1 Recombinant Rabbit mAb [PSH10-30] - BSA and Azide free	HA751347	Human, Mouse	WB,IP
	FOXM1 Recombinant Rabbit mAb [PSH10-93] - BSA and Azide free	HA751383	Human, Mouse	IHC-P,IF-Cell
KD	FOXO1A Recombinant Rabbit mAb [SU33-01] - BSA and Azide free	HA750140	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	FOXP3 Recombinant Rabbit mAb [JF0898] - BSA and Azide free	HA750335	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
Relative expression	FOXP3 Recombinant Rabbit mAb [PSH07-47] - BSA and Azide free	HA751168	Human, Mouse, Rat	WB,IHC-P
	FRA1 Recombinant Rabbit mAb [PSH06-95] - BSA and Azide free	HA751120	Human	WB
Relative expression	Fragilis Recombinant Rabbit mAb [JU73-02] - BSA and Azide free	HA750453	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Frataxin Recombinant Rabbit mAb [PSH15-18]	HA723722	Human, Rat	WB,IHC-P
	F-spondin Recombinant Rabbit mAb [PSH02-84] - BSA and Azide free	HA750837	Human	WB,IF-Cell,FC
	FUNDC1 Recombinant Rabbit mAb [PSH0-92] - BSA and Azide free	HA750672	Human, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	FUS/TLS Recombinant Rabbit mAb [JJ09-31] - BSA and Azide free	HA750326	Human, Mouse, Cynomolgus monkey, Pig	WB,IHC-P,FC,IF-Cell
	FXR1 Recombinant Rabbit mAb [PSH01-67] - BSA and Azide free	HA750743	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC,IP
	FXR1 Recombinant Rabbit mAb [PSH01-35] - BSA and Azide free	HA750713	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,IP
Relative expression	GABA B Receptor 1 Recombinant Rabbit mAb [PSH06-21] - BSA and Azide free	HA751059	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr
Relative expression	GABA B Receptor 1 Recombinant Rabbit mAb [PSH03-87] - BSA and Azide free	HA750912	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P
Relative expression	GABRA1 Recombinant Mouse mAb [PSH10-46] - BSA and Azide free	HA610227	Mouse, Rat	WB
	GAD65 Recombinant Antibody [PSH06-80] - Chicken IgY (Chimeric) - BSA and Azide free	HA610307	Mouse, Rat	IHC-Fr,WB
	GAD65 Recombinant Antibody [PSH06-80] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610289	Mouse, Rat	IHC-Fr,IHC-P,WB
	GAD65 Recombinant Rabbit mAb [PSH06-80] - BSA and Azide free	HA751102	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
	GAD65 Recombinant Antibody [PSH06-80] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610232	Mouse, Rat, Human, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
	GAD67 Recombinant Mouse mAb [PSH08-30] - BSA and Azide free	HA610216	Mouse, Rat	WB,IHC-P,IHC-Fr
	GADD45B Recombinant Rabbit mAb [PSH09-52] - BSA and Azide free	HA751299	Human, Mouse, Rat	WB
	GALC Recombinant Rabbit mAb [PSH03-45] - BSA and Azide free	HA750879	Human, Mouse, Rat	WB,IF-Cell
KD	GALC Recombinant Rabbit mAb [PSH04-32] - BSA and Azide free	HA750941	Human	WB,IHC-P
Relative expression	Galectin 1 Recombinant Rabbit mAb [JM13-37] - BSA and Azide free	HA750441	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Relative expression	Galectin 9 Recombinant Rabbit mAb [PSH02-44] - BSA and Azide free	HA750809	Human, Mouse	WB,FC
	Galectin 9 Recombinant Rabbit mAb [PSH02-45] - BSA and Azide free	HA750810	Human, Rat	WB
KD	GALT Recombinant Rabbit mAb [PSH02-87] - BSA and Azide free	HA750841	Human, Mouse, Rat	WB,IF-Cell,FC
KD	GALT Recombinant Rabbit mAb [PSH02-88] - BSA and Azide free	HA750842	Human, Mouse	WB
KD	gamma Adaptin Recombinant Rabbit mAb [PSH02-65] - BSA and Azide free	HA750822	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	GANAB Recombinant Rabbit mAb [PSH01-53] - BSA and Azide free	HA750729	Human	WB,IHC-P,IF-Tissue
KD	GANAB Recombinant Rabbit mAb [PSH01-52] - BSA and Azide free	HA750728	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP
	GAPDH Recombinant Rabbit mAb [SA30-01] - BSA and Azide free	HA750022	Human, Mouse, Rat, Chicken, Drosophila melanogaster, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	GAPDH Recombinant Rabbit mAb [PD00-07] - BSA and Azide free	HA750478	Human, Mouse, Rat, Escherichia coli, Zebrafish, Drosophila melanogaster	WB,IHC-P,IF-Cell,FC
Cell treatment	Gasdermin D (N terminal) Recombinant Rabbit mAb [PD00-18] - BSA and Azide free	HA750486	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Gasdermin D (N terminal) Recombinant Rabbit mAb [PSH10-87] - BSA and Azide free	HA751379	Human	WB,IP
Relative expression	GATA2 Recombinant Rabbit mAb [PSH04-53] - BSA and Azide free	HA750960	Human, Rat	WB,IF-Cell
	GATA3 Recombinant Mouse mAb [A3G3-R] - BSA and Azide free	HA610069	Human	WB,IF-Cell,IHC-P,FC
KD	GC1q R Recombinant Mouse mAb [A8F5-R] - BSA and Azide free	HA610076	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	GC1q R Recombinant Mouse mAb [A8F4-R] - BSA and Azide free	HA610057	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Relative expression	GCH1 Recombinant Rabbit mAb [PSH08-00] - BSA and Azide free	HA751200	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue
	GCSF Receptor Recombinant Rabbit mAb [PSH03-44] - BSA and Azide free	HA750878	Human, Mouse, Rat	WB
Relative expression	GDF15 Recombinant Rabbit mAb [PSH10-09] - BSA and Azide free	HA751339	Human	WB,IHC-P
	Genetics of Parkinson's Disease: Lysosomal Dysfunction Antibody Sampler Kit	HAK21127	Human, Mouse, Rat	WB,IHC-P
Relative expression	GFAP Recombinant Antibody [SA03-04] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610317	Mouse, Rat	IHC-Fr,WB
Relative expression	GFAP Recombinant Mouse mAb [PSH10-24] - BSA and Azide free	HA610225	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue
	GFP Recombinant Rabbit mAb [SY0243] - BSA and Azide free	HA750116	Species independent	WB,IF-Cell,IF-Tissue,IP,IHC-P
Relative expression	Ghrelin Recombinant Rabbit mAb [PSH08-77] - BSA and Azide free	HA751242	Human, Mouse	IHC-P,IHC-Fr
	Giantin Recombinant Rabbit mAb [PSH03-86] - BSA and Azide free	HA750911	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	GIRK2 Recombinant Rabbit mAb [PSH05-91] - BSA and Azide free	HA751027	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IP,IHC-Fr
	GITR Recombinant Rabbit mAb [PSH04-06] - BSA and Azide free	HA750929	Human	WB,IF-Cell,FC
	Gli1 Recombinant Rabbit mAb [JF09-08] - BSA and Azide free	HA750361	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Gli2 Recombinant Rabbit mAb [PSH11-23] - BSA and Azide free	HA751402	Human, Mouse, Rat	WB
	Gluconeogenesis Antibody Sampler Kit	HAK21135	Human	WB
	Glucose Transporter GLUT1 Recombinant Antibody [SA0377] - Rat IgG1 (Chimeric)	HA601521	Human, Mouse, Rat	IHC-Fr,IHC-P,IF-Tissue,WB
KD	Glucose Transporter GLUT1 Recombinant Rabbit mAb [SA0377] - BSA and Azide free	HA750002	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	Glucose Transporter GLUT3 Recombinant Rabbit mAb [JA50-31] - BSA and Azide free	HA750416	Human, Mouse, Rat	WB,FC
Relative expression	Glucose Transporter GLUT4 Recombinant Rabbit mAb [PSH10-06] - BSA and Azide free	HA751336	Mouse, Rat	WB,IHC-Fr
KD	Glucosidase 2 subunit beta Recombinant Mouse mAb [A6H9-R] - BSA and Azide free	HA610151	Human	WB
	Glutamate receptor 1 Recombinant Rabbit mAb [SD2010] - BSA and Azide free	HA750281	Human, Mouse, Rat	WB,IHC-P,IHC-Fr
	Glutamine Synthetase Recombinant Mouse mAb [A3G2-R] - BSA and Azide free	HA610100	Human, Mouse, Rat	WB,IHC-P
	Glutamine Synthetase Recombinant Rabbit mAb [PSH10-48] - BSA and Azide free	HA751351	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr
	Glutamine Synthetase Recombinant Antibody [PSH10-48] - Chicken IgY (Chimeric)	HA601516	Mouse, Rat	IHC-Fr,WB
	Glutamine Synthetase Recombinant Antibody [PSH10-48] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610268	Human, Mouse, Rat	IHC-Fr,WB
	Glutathione peroxidase 1 Recombinant Mouse mAb [C5-A10-R] - BSA and Azide free	HA610078	Human	WB,IF-Cell,IHC-P,FC
	Glutathione peroxidase 1 Recombinant Mouse mAb [A6-C0-B9-R] - BSA and Azide free	HA610077	Human	WB,IF-Cell,IHC-P,FC
	GNAS Recombinant Mouse mAb [A8F11-R] - BSA and Azide free	HA610150	Human, Mouse, Rat	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	iFluor™ 594 Conjugated Goat anti-Chicken IgY Goat pAb	HA1140	Chicken	IHC-Fr,IF-Tissue,IF-Cell
	iFluor™ 488 Conjugated Goat anti-Chicken IgY Goat pAb	HA1139	Chicken	IHC-Fr,IF-Tissue,IF-Cell
	Golgin 97 Recombinant Rabbit mAb [PSH13-73] - BSA and Azide free	HA751486	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
KD	GOLPH4/GPP130 Recombinant Rabbit mAb [PSH02-98] - BSA and Azide free	HA750853	Human	WB,IHC-P,IF-Cell,IF-Tissue,FC
KD	GOT1 Recombinant Mouse mAb [A2D6-R] - BSA and Azide free	HA610154	Human, Mouse, Rat	WB
	GPD2 Recombinant Rabbit mAb [PSH0-31] - BSA and Azide free	HA750605	Human	WB,FC
	GNMB Recombinant Rabbit mAb [PSH0-82] - BSA and Azide free	HA750662	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Relative expression	GPR56 Recombinant Rabbit mAb [PSH05-43] - BSA and Azide free	HA751002	Human	WB,IHC-P
KD	GPX4 Recombinant Rabbit mAb [JU11-31] - BSA and Azide free	HA750455	Human, Mouse, Rat, Zebrafish, Pig, Cynomolgus monkey	WB,IHC-P,IF-Cell,FC
Relative expression	Granzyme A Recombinant Rabbit mAb [PSH14-39] - BSA and Azide free	HA751521	Human	WB
	Granzyme B Recombinant Rabbit mAb [PSH14-20] - BSA and Azide free	HA751513	Mouse	WB
Relative expression	Granzyme B Recombinant Rabbit mAb [PSH14-66]	HA723664	Human, Mouse	WB,IHC-P
Relative expression	Granzyme B Recombinant Rabbit mAb [PSH14-66] - BSA and Azide free	HA751527	Human, Mouse	WB,IHC-P
	Grp75 Recombinant Mouse mAb [6-A9-R] - BSA and Azide free	HA610128	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	GRP78 / BIP Recombinant Mouse mAb [A5-C9-G6-R] - BSA and Azide free	HA610190	Human, Mouse, Rat	WB
Cell treatment	GRP78 / BIP Recombinant Mouse mAb [C9-9-R] - BSA and Azide free	HA610013	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IF-Cell
KD	GSK3 beta Recombinant Rabbit mAb [SY28-03] - BSA and Azide free	HA750131	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
	GST Recombinant Rabbit mAb [SN70-05] - BSA and Azide free	HA750257	Species independent	WB,IF-Cell,IP
KD	GST3 Recombinant Mouse mAb [11D5-R] - BSA and Azide free	HA610147	Human	WB,IF-Cell,IHC-P
	GWL Recombinant Rabbit mAb [PSH07-69] - BSA and Azide free	HA751181	Human, Monkey	WB,IF-Cell,IP
KD	HADH Recombinant Mouse mAb [D10-E7-R] - BSA and Azide free	HA610081	Human, Mouse, Rat	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HADHA Recombinant Rabbit mAb [PSH01-28] - BSA and Azide free	HA750706	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell
	HADHA Recombinant Rabbit mAb [PSH01-29] - BSA and Azide free	HA750707	Human, Mouse, Rat, Monkey	WB
	Hamartin Recombinant Rabbit mAb [PSH14-67] - BSA and Azide free	HA751528	Human, Mouse, Rat, Monkey	WB,IP
	Hamartin Recombinant Rabbit mAb [PSH14-67]	HA723665	Human, Mouse, Rat, Monkey	WB,IP
	HAPLN1 Recombinant Rabbit mAb [SY02-17] - BSA and Azide free	HA750130	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P
	HA tag Recombinant Rabbit mAb [SN07-06] - BSA and Azide free	HA750258	Species independent	WB,IP
Independent antibody validation	HA tag Recombinant Rabbit Multiclonal Antibody [PSH01-92] - BSA and Azide free	HA750766	Species independent	WB,IF-Cell,IP
	HA tag Recombinant Mouse mAb [2-G10] - BSA and Azide free	HA610320	Species independent	WB,IHC-P,IF-Cell
Relative expression	HCLS1 Recombinant Rabbit mAb [PSH01-12] - BSA and Azide free	HA750688	Human	WB,IHC-P,IF-Cell,FC,IF-Tissue
KO	HDAC1 Recombinant Rabbit mAb [SY12-04] - BSA and Azide free	HA750083	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,IF-Tissue
KO	HDAC2 Recombinant Rabbit mAb [SY29-02] - BSA and Azide free	HA750134	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	HDAC3 Recombinant Rabbit mAb [SC55-02] - BSA and Azide free	HA750221	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	HDAC5 Recombinant Rabbit mAb [PSH08-70] - BSA and Azide free	HA751235	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,ChIP
	HE4 Recombinant Rabbit mAb [JF62-09] - BSA and Azide free	HA750355	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	Heme Oxygenase 1 (HO-1) Recombinant Rabbit mAb [PSH02-76] - BSA and Azide free	HA750831	Human, Mouse	WB,IF-Cell,IHC-P
	Heme Oxygenase 1 (HO-1) Recombinant Rabbit mAb [SP08-07] - BSA and Azide free	HA750079	Human	WB,IHC-P,IP,FC,IF-Tissue
Relative expression	Hemoglobin subunit gamma 1and2 Recombinant Rabbit mAb [JM84-10] - BSA and Azide free	HA750383	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Hes1 Recombinant Rabbit mAb [SC06-21] - BSA and Azide free	HA750240	Human, Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC
	HEXB Recombinant Rabbit mAb [PSH06-83] - BSA and Azide free	HA751106	Human	WB,IHC-P,IF-Cell,FC,IF-Tissue
	Hexokinase 1 Recombinant Rabbit mAb [ST47-05] - BSA and Azide free	HA750168	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Hexokinase II Recombinant Rabbit mAb [PSH07-97] - BSA and Azide free	HA751199	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC,IP
Relative expression	Hexokinase II Recombinant Rabbit mAb [PSH08-81] - BSA and Azide free	HA751246	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP
Cell treatment	HIF-1 alpha Recombinant Rabbit mAb [PSH07-00] - BSA and Azide free	HA751127	Human, Mouse	WB,IF-Cell,IHC-P
KD	HIF-1 beta Recombinant Rabbit mAb [ST05-25] - BSA and Azide free	HA750187	Human, Mouse, Rat	WB,IHC-P
	HIP2 / LIG Recombinant Rabbit mAb [PSH11-24] - BSA and Azide free	HA751403	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
Cell treatment	Histone H2A (acetyl K9) Recombinant Rabbit mAb [PS01-45] - BSA and Azide free	HA750591	Human, Mouse, Rat	WB,IHC-P,IF-Cell,CUT&Tag-seq,ChIP
	Histone H2B Recombinant Rabbit mAb [SD20-63] - BSA and Azide free	HA750287	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP,ChIP
Cell treatment	Histone H2B (acetyl K20) Recombinant Rabbit mAb [PS01-46] - BSA and Azide free	HA750592	Human, Mouse, Rat	WB,IHC-P,IF-Cell,ChIP
	Histone H3 Recombinant Mouse mAb [6-A7-R] - BSA and Azide free	HA610203	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,ChIP
	Histone H3 Recombinant Rabbit mAb [JJ090-07] - BSA and Azide free	HA750318	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP,IP
	Histone H3 Recombinant Mouse mAb [A11-D7-R] - BSA and Azide free	HA610153	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Histone H3.3 Recombinant Rabbit mAb [ST50-08] - BSA and Azide free	HA750184	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,ChIP
Cell treatment	Histone H3 (acetyl K14) Recombinant Rabbit mAb [JU43-26] - BSA and Azide free	HA750454	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,SNAP-ChIP,CUT&Tag-seq
Cell treatment	Histone H3 (acetyl K27) Recombinant Rabbit mAb [PSH13-93] - BSA and Azide free	HA751495	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Histone H3 (acetyl K4) Recombinant Rabbit mAb [PSH03-71] - BSA and Azide free	HA750905	Human, Mouse, Rat	WB,IF-Cell,ChIP,Dot Blot
Cell treatment	Histone H3 (acetyl K56) Recombinant Rabbit mAb [SU30-10] - BSA and Azide free	HA750163	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP,CUT&Tag-seq
Cell treatment	Histone H3 (acetyl K9) Recombinant Rabbit mAb [PSH04-47] - BSA and Azide free	HA750956	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,FC,ChIP,Dot Blot,IP
Cell treatment	Histone H3 (di methyl K27) Recombinant Rabbit mAb [PSH05-88] - BSA and Azide free	HA751035	Human, Mouse, Rat	WB,IF-Cell,ChIP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Histone H3 (di methyl K9) Recombinant Rabbit mAb [SN07-30] - BSA and Azide free	HA750259	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P
	Histone H3 (mono+di+tri methyl K36) Recombinant Rabbit mAb [PSH07-14] - BSA and Azide free	HA751135	Human, Mouse, Rat	WB,IHC-P,IF- Cell,FC,ChIP
	Histone H3 (mono+di+tri methyl K4) Recombinant Rabbit mAb [PSH06-30] - BSA and Azide free	HA751067	Human, Mouse, Rat	WB,IHC-P,IF- Cell,FC,ChIP
	Histone H3 (mono+di+tri methyl K79) Recombinant Rabbit mAb [PSH09-10] - BSA and Azide free	HA751264	Human, Mouse, Rat, Monkey	WB,IHC-P,IF- Cell,ChIP
	Histone H3 (mono+di+tri methyl K79) Recombinant Rabbit mAb [SR42-06] - BSA and Azide free	HA750051	Human, Mouse, Rat	WB,IHC-P,ChIP
	Histone H3 (mono methyl K18) Recombinant Rabbit mAb [SA42-07] - BSA and Azide free	HA750006	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P
	Histone H3 (mono methyl K27) Recombinant Rabbit mAb [PSH05-67] - BSA and Azide free	HA751019	Human, Mouse, Rat	WB,IF-Cell,FC,ChIP
	Histone H3 (mono methyl K4) Recombinant Rabbit mAb [PSH07-15] - BSA and Azide free	HA751136	Human, Mouse, Rat	WB,IF-Cell,ChIP,Dot Blot
KD	Histone H3 (mono methyl K79) Recombinant Rabbit mAb [PSH07-66] - BSA and Azide free	HA751175	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC- P,FC,ChIP
	Histone H3 (mono methyl R2) Recombinant Rabbit mAb [ST0427] - BSA and Azide free	HA750197	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P,ChIP
Cell treatment	Histone H3 (tri methyl K27) Recombinant Rabbit mAb [PSH05-15] - BSA and Azide free	HA750988	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF- Tissue,ChIP
	Histone H3 (tri methyl K36) Recombinant Rabbit mAb [PSH15-01]	HA723701	Human, Mouse, Rat, Monkey	WB,IHC-P,IF- Cell,IP,Dot Blot,ChIP
	Histone H4 Recombinant Rabbit mAb [SD209-04] - BSA and Azide free	HA750294	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P,ChIP
Cell treatment	Histone H4 (acetyl K16) Recombinant Rabbit mAb [JB21-44] - BSA and Azide free	HA750464	Human, Mouse, Rat	WB,IHC-P,FC,IF- Cell,IF-Tissue,ChIP
Cell treatment	Histone H4 (acetyl K5) Recombinant Rabbit mAb [SR31-07] - BSA and Azide free	HA750050	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC- P,IP,ChIP,CUT&Tag- seq
	Histone H4 (acetyl K8) Recombinant Rabbit mAb [PS01-25] - BSA and Azide free	HA750568	Human, Mouse, Rat	WB,IF-Cell,IHC- P,ChIP
	HLA-ABC (MHC I) Recombinant Rabbit mAb [PSH12-02] - BSA and Azide free	HA751435	Human, Mouse, Rat	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	HLA Class I Recombinant Mouse mAb [PSH12-15] - BSA and Azide free	HA610264	Human	IF-Cell,FC
	HLA-DR Recombinant Rabbit mAb [SC06-78] - BSA and Azide free	HA750230	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	HLA-DR Recombinant Mouse mAb [10-D8-R] - BSA and Azide free	HA610127	Human	WB,IF-Cell,IHC-P,IF-Tissue
	HMGA2 Recombinant Rabbit mAb [PSH05-53] - BSA and Azide free	HA751012	Human, Mouse	WB,IF-Cell,IHC-P,FC,IF-Tissue
KO	HMGB1 Recombinant Rabbit mAb [SA39-03] - BSA and Azide free	HA750011	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
	HMGCR Recombinant Rabbit mAb [JF0981] - BSA and Azide free	HA750345	Human, Mouse, Rat	WB,IP,IHC-P
Cell treatment	HMGCS1 Recombinant Rabbit mAb [PSH08-52] - BSA and Azide free	HA751227	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP
	HMGCS1 Recombinant Rabbit mAb [PSH08-51] - BSA and Azide free	HA751226	Human, Mouse, Rat, Monkey	WB,IP
	HNF-4-alpha Recombinant Rabbit mAb [SN72-03] - BSA and Azide free	HA750254	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	HNRNPA0 Recombinant Rabbit mAb [PSH02-77] - BSA and Azide free	HA750832	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	hnRNP A2B1 Recombinant Mouse mAb [PSH06-47] - BSA and Azide free	HA610211	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue
	HNRPAB Recombinant Rabbit mAb [PSH05-37] - BSA and Azide free	HA750996	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue
	Homeobox protein cut-like 1 Recombinant Rabbit mAb [PSH14-70] - BSA and Azide free	HA751530	Mouse, Rat	IHC-Fr,IHC-P
	Homeobox protein cut-like 1 Recombinant Rabbit mAb [PSH14-70]	HA723668	Mouse, Rat	IHC-Fr,IHC-P
Relative expression	Homer1 Recombinant Antibody [JE65-19] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610255	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
Relative expression	Homer1 Recombinant Antibody [JE65-19] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610300	Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Homer1 Recombinant Antibody [JE65-19] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610299	Mouse, Rat	IHC-Fr,WB
Relative expression	HOPX / HOD Recombinant Antibody [PSH10-32] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610318	Mouse, Rat	IHC-Fr,WB
Relative expression	HOPX / HOD Recombinant Rabbit mAb [PSH10-32] - BSA and Azide free	HA751349	Mouse, Rat	WB,IHC-Fr,IHC-P
Relative expression	Hormone sensitive lipase / HSL Recombinant Rabbit mAb [PSH07-22] - BSA and Azide free	HA751143	Human, Mouse, Rat	WB,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HPRT Recombinant Rabbit mAb [JU03-26] - BSA and Azide free	HA750452	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IP
	Hsc70 Recombinant Rabbit mAb [SR39-04] - BSA and Azide free	HA750044	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	hSET1/SET1 Recombinant Rabbit mAb [PH00-01] - BSA and Azide free	HA750554	Human, Mouse, Rat	WB,IHC-P
	hSET1/SET1 Recombinant Rabbit mAb [PH00-02] - BSA and Azide free	HA750556	Human, Mouse, Rat	WB,FC
	Hsp27 Recombinant Rabbit mAb [JJ09-13] - BSA and Azide free	HA750321	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Hsp40 Recombinant Rabbit mAb [PSH12-10] - BSA and Azide free	HA751439	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	HSP60 Recombinant Rabbit mAb [ST48-04] - BSA and Azide free	HA750175	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP
	Hsp90 alpha Recombinant Rabbit mAb [SY14-06] - BSA and Azide free	HA750091	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P
KD	Hsp90 beta Recombinant Rabbit mAb [SY46-01] - BSA and Azide free	HA750090	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	HSP90 Beta Recombinant Mouse mAb [2-1-G3-R] - BSA and Azide free	HA610060	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	HSPA1L Recombinant Mouse mAb [7D3-R] - BSA and Azide free	HA610120	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Human ACE, Tag Free Protein	HA211209	Human	
	Human CCL17/TARC enzyme-linked immunoassay kit	EH0164	Human	ELISA
	Human CCL19 / MIP-3 beta, Tag Free Protein	HA211211	Human	
	Human CCL3/MIP-1 enzyme-linked immunoassay kit	EH0146	Human	ELISA
	Human CCL4/MIP-1, Tag Free Protein	HA211213	Human	
	Human CCL5/RANTES enzyme-linked immunoassay kit(one step)	EHY0027	Human	ELISA
	Human CD155/PVR(Poliovirus Receptor) enzyme-linked immunoassay kit(one step)	EHY0155	Human	ELISA
	Human CD163, Tag Free (ECD) Protein	HA211143	Human	
	Human CD328/Siglec-7 enzyme-linked immunoassay kit	EH0163	Human	ELISA
	Human CD45 Recombinant Mouse mAb [PSH03-64] - BSA and Azide free	HA610158	Human	IF-Cell,FC
	Human CD83, C-His Tag Protein	HA211203	Human	
	Human Chitotriosidase/CHIT1, C-His Tag Protein	HA211215	Human	

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CNTF enzyme-linked immunoassay kit	EH0147	Human	ELISA
	Human Complement C3 Recombinant Rabbit mAb [PSH15-08] - BSA and Azide free (Detector)	HA725141	Human	ELISA(Det)
	Biotin Conjugated Human Complement C3 Recombinant Rabbit mAb [PSH15-08] - Detector	HA725142B	Human	ELISA(Det),ELISA
	Human Complement C3 Recombinant Rabbit mAb [PSH15-07] - BSA and Azide free (Capture)	HA725140	Human	ELISA(Cap)
	Human CTGF/CCN2 enzyme-linked immunoassay kit(one step)	EHY0156	Human	ELISA
	Human CXADR/CAR, C-His Tag Protein	HA211206	Human	
	Human CXCL6/GCP-2 enzyme-linked immunoassay kit(one step)	EH0165	Human	ELISA
	Human Endoglin / CD105 Recombinant Rabbit mAb [PSH14-59] - BSA and Azide free (Detector)	HA723655	Human	ELISA(Det)
	Human Endoglin / CD105 Recombinant Rabbit mAb [PSH14-60] - BSA and Azide free (Capture)	HA723657	Human	ELISA(Cap)
	Human Endoglin / CD105 Recombinant Rabbit mAb [PSH14-58] - BSA and Azide free (Capture)	HA723654	Human	ELISA(Cap)
	Biotin Conjugated Human Endoglin / CD105 Recombinant Rabbit mAb [PSH14-59] - Detector	HA723656B	Human	ELISA(Det),ELISA
	Human EphA2 Recombinant Rabbit mAb [PSH15-13] - BSA and Azide free (Detector)	HA723715	Human	ELISA(Det)
	Human EphA2 Recombinant Rabbit mAb [PSH15-12] - BSA and Azide free (Capture)	HA723714	Human	ELISA(Cap)
	Biotin Conjugated Human EphA2 Recombinant Rabbit mAb [PSH15-13] - Detector	HA723716B	Human	ELISA(Det),ELISA
	Human Fibrillin-1/FBN1, Tag Free Protein	HA211184	Human	
	Human G-CSF enzyme-linked immunoassay kit	EH0149	Human	ELISA
	Human IGFBP3/IBP3 enzyme-linked immunoassay kit(one step)	EHY0157	Human	ELISA
	Human IL-10RB, C-His Tag Protein	HA211205	Human	
Relative expression	Human IL-18 Recombinant Rabbit mAb [PSH02-09] - BSA and Azide free	HA750779	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-18 Recombinant Rabbit mAb [PSH02-08] - BSA and Azide free	HA750778	Human	WB,IHC-P,IF-Tissue
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free	HA750745	Human	WB
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free	HA750744	Human	WB
	Human IL-21, Tag Free Protein	HA211070	Human	
	Human IL-25 Recombinant Rabbit mAb [PSH15-23] - BSA and Azide free (Capture)	HA725016	Human	ELISA(Cap)
	Human IL-25 Recombinant Rabbit mAb [PSH15-24] - BSA and Azide free (Detector)	HA725017	Human	ELISA(Det)
	Biotin Conjugated Human IL-25 Recombinant Rabbit mAb [PSH15-24] - Detector	HA725018B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human IL-29 Recombinant Rabbit mAb [PSH15-11] - Detector	HA723713B	Human	ELISA(Det),ELISA
	Human IL-29 Recombinant Rabbit mAb [PSH15-11] - BSA and Azide free (Detector)	HA723712	Human	ELISA(Det)
	Human IL-29 Recombinant Rabbit mAb [PSH15-10] - BSA and Azide free (Detector)	HA723710	Human	ELISA(Det)
	Human IL-29 Recombinant Rabbit mAb [PSH15-09] - BSA and Azide free (Capture)	HA723709	Human	ELISA(Cap)
	Biotin Conjugated Human IL-29 Recombinant Rabbit mAb [PSH15-10] - Detector	HA723711B	Human	ELISA(Det),ELISA
	Human IL-35 (linked heterodimer), C-hFc Tag Protein	HA211212	Human	
	Human IL-4R enzyme-linked immunoassay kit	EH0145	Human	ELISA
	Human Insulin, Tag Free Protein	HA211183	Human	
	Human KLK3, Tag Free Protein	HA211208	Human	
	Human LAMC2 Recombinant Rabbit mAb [PSH15-15] - BSA and Azide free (Detector)	HA723718	Human	ELISA(Det)
	Human LAMC2 Recombinant Rabbit mAb [PSH15-14] - BSA and Azide free (Capture)	HA723717	Human	ELISA(Cap)
	Biotin Conjugated Human LAMC2 Recombinant Rabbit mAb [PSH15-15] - Detector	HA723719B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Leptin Receptor Recombinant Rabbit mAb [PSH14-71] - BSA and Azide free (Capture)	HA725058	Human	ELISA(Cap)
	Human Leptin Receptor Recombinant Rabbit mAb [PSH14-72] - BSA and Azide free (Detector)	HA725059	Human	ELISA(Det)
	Biotin Conjugated Human Leptin Receptor Recombinant Rabbit mAb [PSH14-72] - Detector	HA725060B	Human	ELISA(Det),ELISA
	HRP Conjugated Human Lipocalin-2 Recombinant Rabbit mAb [PSH13-68] - Detector	HA725111H	Human	ELISA(Det),ELISA
	Human MDH1, Tag Free Protein	HA211217	Human	
	Human Mesothelin enzyme-linked immunoassay kit(one step)	EHY0152	Human	ELISA
	Human MIP-5/CCL15, C-His Tag Protein	HA211207	Human	
	Human Myeloperoxidase Recombinant Rabbit mAb [PSH14-79] - BSA and Azide free (Detector)	HA723674	Human	ELISA(Det)
	Biotin Conjugated Human Myeloperoxidase Recombinant Rabbit mAb [PSH14-79] - Detector	HA723675B	Human	ELISA(Det),ELISA
	Human Myeloperoxidase Recombinant Rabbit mAb [PSH14-78] - BSA and Azide free (Capture)	HA723673	Human	ELISA(Cap)
	APC Conjugated Human NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53]	HA723757F	Human	FC
	Human Nectin-4 enzyme-linked immunoassay kit	EH0151	Human	ELISA
	Human Neuropeptide Y, Tag Free Protein	HA211142	Human	
	HRP Conjugated Human Neuropilin-1 Recombinant Rabbit mAb [PSH14-74] - Detector	HA725114H	Human	ELISA(Det),ELISA
	Human Neuropilin-1 Recombinant Rabbit mAb [PSH14-73] - BSA and Azide free (Capture)	HA725112	Human	ELISA(Cap)
	Human Neuropilin-1 Recombinant Rabbit mAb [PSH14-74] - BSA and Azide free (Detector)	HA725113	Human	ELISA(Det)
	Human NG2/MCSP/CSPG4 enzyme-linked immunoassay kit(one step)	EHY0162	Human	ELISA
	Human NPTXR, C-His Tag Protein	HA211216	Human	
	Human Oxytocin-neurophysin 1/OXT enzyme-linked immunoassay kit	EH0153	Human	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human PAI-1/Serpin E1 enzyme-linked immunoassay kit(one step)	EHY0160	Human	ELISA
	Human PD-1 enzyme-linked immunoassay kit	EH0062	Human	ELISA
	Human RBP4 Recombinant Rabbit mAb [PSH14-88] - BSA and Azide free (Capture)	HA723688	Human	ELISA(Cap)
	Human RBP4 Recombinant Rabbit mAb [PSH14-89] - BSA and Azide free (Detector)	HA723689	Human	ELISA(Det)
	Biotin Conjugated Human RBP4 Recombinant Rabbit mAb [PSH14-89] - Detector	HA723690B	Human	ELISA(Det),ELISA
	Human Resistin Recombinant Rabbit mAb [PSH14-61] - BSA and Azide free (Capture)	HA723658	Human	ELISA(Cap)
	Biotin Conjugated Human Resistin Recombinant Rabbit mAb [PSH14-62] - Detector	HA723660B	Human	ELISA(Det),ELISA
	Human Resistin Recombinant Rabbit mAb [PSH14-63] - BSA and Azide free (Capture)	HA723661	Human	ELISA(Cap)
	Human Resistin Recombinant Rabbit mAb [PSH14-62] - BSA and Azide free (Detector)	HA723659	Human	ELISA(Det)
	Human sFRP-1, Tag Free Protein	HA211218	Human	
	Human SIRPA, Tag Free Protein	HA211204	Human	
	Human SLIT2, Tag Free Protein	HA211214	Human	
	Human Thioredoxin/TXN, C-His Tag Protein	HA211189	Human	
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free	HA750593	Human	IHC-P
	Human TNF RII/TNFRSF1B enzyme-linked immunoassay kit	EH0148	Human	ELISA
	Human TNF RI/TNFRSF1A enzyme-linked immunoassay kit (one step)	EHY0128	Human	ELISA
	Human tPA enzyme-linked immunoassay kit(one step)	EHY0158	Human	ELISA
	Human TREM2 Activity Antibody Sampler Kit	HAK21130	Human	WB
	Human ULBP1/NKG2DL1 enzyme-linked immunoassay kit	EH0150	Human	ELISA
	Human Uromodulin/UMOD, Tag Free Protein	HA211187	Human	

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human VCAM-1/CD106 enzyme-linked immunoassay kit(one step)	EHY0067	Human	ELISA
	Human YKL-40/CHI3L1 enzyme-linked immunoassay kit(one step)	EHY0159	Human	ELISA
KD	Huntingtin Recombinant Rabbit mAb [JB89-34] - BSA and Azide free	HA750462	Human, Mouse, Rat, Pig	WB,IHC-P,FC,IHC-Fr,IF-Cell,IF-Tissue
Relative expression	HUWE1 Recombinant Rabbit mAb [PSH14-37] - BSA and Azide free	HA751519	Human, Mouse, Rat	WB,IHC-P
	HUWE1 Recombinant Rabbit mAb [PSH02-38] - BSA and Azide free	HA750803	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Iba1 Recombinant Mouse mAb [8G8] - BSA and Azide free	HA610321	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,FC,IP,IHC-Fr,IF-Cell,IF-Tissue
Relative expression	Iba1 Recombinant Rabbit mAb [JM36-62] - BSA and Azide free	HA750440	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC,IP
	Iba1 Recombinant Antibody [JM36-62] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610213	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
Relative expression	Iba1 Recombinant Antibody [JM36-62] - Chicken IgY (Chimeric) - BSA and Azide free	HA610221	Human, Mouse, Rat, Pig, Cynomolgus monkey	IHC-Fr,WB,FC
	Iba1 Recombinant Antibody [JM36-62] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610214	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,WB
Relative expression	ICAM1 Recombinant Rabbit mAb [ST0487] - BSA and Azide free	HA750176	Human	WB,IF-Tissue,IHC-P
	Id1 Recombinant Rabbit mAb [JM92-13] - BSA and Azide free	HA750431	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment, Relative expression	IDO Recombinant Rabbit mAb [PSH12-31] - BSA and Azide free	HA751445	Human	WB,IF-Cell,IHC-P,FC,IP
	IFNAR1 Recombinant Rabbit mAb [SR45-08] - BSA and Azide free	HA750046	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Relative expression, KD	IFNGR1 Recombinant Rabbit mAb [PSH06-11] - BSA and Azide free	HA751049	Human	WB,IF-Cell
	IFT88 Recombinant Rabbit mAb [PSH11-61] - BSA and Azide free	HA751417	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell
Relative expression	IGF1 Receptor beta Recombinant Rabbit mAb [PSH11-72] - BSA and Azide free	HA751427	Human, Mouse, Rat	WB,IF-Cell
Relative expression	IGF2BP1 Recombinant Rabbit mAb [PSH06-37] - BSA and Azide free	HA751076	Human	WB,IF-Cell,FC
	IKB alpha Recombinant Rabbit mAb [SZ00-07] - BSA and Azide free	HA750068	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	IKK alpha Recombinant Rabbit mAb [SN06-99] - BSA and Azide free	HA750246	Human, Mouse, Rat	WB,IHC-P,IP,FC
KD	IKK alpha + IKK beta Recombinant Rabbit mAb [SN63-02] - BSA and Azide free	HA750248	Human, Mouse, Rat	WB,IF-Cell,IP
Relative expression	IL-13R alpha 2 Recombinant Rabbit mAb [PSH03-85] - BSA and Azide free	HA750910	Human	WB,IF-Cell,IHC-P,FC
Relative expression	IL-18 Recombinant Rabbit mAb [PSH12-04] - BSA and Azide free	HA751511	Mouse	WB,IF-Cell,FC
	IL-18 Recombinant Rabbit mAb [PS00-41] - BSA and Azide free	HA750676	Human	WB,IHC-P
	IL-18 Recombinant Rabbit mAb [PS01-06] - BSA and Azide free	HA750677	Human	WB,FC
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R] - BSA and Azide free	HA610004	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-27A p28 Recombinant Rabbit mAb [PSH14-93]	HA723695	Mouse	WB
	IL-4R Recombinant Rabbit mAb [PSH10-86] - BSA and Azide free	HA751378	Human, Mouse, Rat	WB,IP,IF-Cell
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-96] - BSA and Azide free	HA751386	Human	WB,FC
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-97] - BSA and Azide free	HA751387	Human	WB
	ILKAP Recombinant Rabbit mAb [PSH06-85] - BSA and Azide free	HA751108	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
KD	IMP3 Recombinant Rabbit mAb [PSH05-54] - BSA and Azide free	HA751013	Human, Mouse	WB,IHC-P,IF-Tissue
	Inhibin alpha Recombinant Rabbit mAb [PD01-29] - BSA and Azide free	HA751016	Human, Mouse, Rat	IHC-P,WB
Cell treatment	iNOS Recombinant Rabbit mAb [PSH03-69] - BSA and Azide free	HA750903	Human, Mouse	WB,IHC-P,IF-Cell
Relative expression, KD	INPP4A Recombinant Rabbit mAb [PSH01-27] - BSA and Azide free	HA750705	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Relative expression	Insulin Recombinant Rabbit mAb [SA0410] - BSA and Azide free	HA750004	Human, Mouse, Rat	IF-Cell,IF-Tissue,IHC-P,IHC-Fr,WB
KD, Relative expression	Insulin Receptor Beta Recombinant Rabbit mAb [SN20-13] - BSA and Azide free	HA750276	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC
	Integrin alpha 1 Recombinant Rabbit mAb [PSH14-05] - BSA and Azide free	HA751500	Human, Mouse, Rat	WB,IHC-P
KD	Integrin alpha 2 Recombinant Rabbit mAb [SN0752] - BSA and Azide free	HA750263	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	Integrin alpha 3 Recombinant Mouse mAb [A7H2-R] - BSA and Azide free	HA610148	Human	WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Integrin alpha 6 Recombinant Rabbit mAb [SR45-00] - BSA and Azide free	HA750028	Human, Mouse, Rat	WB,IHC-P
KD	Integrin alpha V Recombinant Rabbit mAb [SC56-07] - BSA and Azide free	HA750201	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
KD	Integrin beta 1 Recombinant Rabbit mAb [SA40-08] - BSA and Azide free	HA750008	Human, Mouse, Rat, Zebrafish	WB,IHC-P
	Integrin beta 5 Recombinant Rabbit mAb [PSH11-53] - BSA and Azide free	HA751415	Human	IF-Cell,FC
	Integrin beta 5 Recombinant Rabbit mAb [PSH10-13] - BSA and Azide free	HA751342	Human, Rat, Monkey	WB,IHC-P
Cell treatment	Interferon alpha 2 Recombinant Rabbit mAb [PSH08-23] - BSA and Azide free	HA751214	Human	WB,IF-Cell,IHC-P
	Interferon gamma Recombinant Rabbit mAb [JM10-10] - BSA and Azide free	HA750373	Human	WB,FC
Cell treatment	IP10 Recombinant Rabbit mAb [JA10-82] - BSA and Azide free	HA750411	Human	WB,IF-Cell,IF-Tissue,IHC-P
Relative expression	IP3 receptor Recombinant Mouse mAb [PSH14-18] - BSA and Azide free	HA610290	Mouse, Rat	WB,IHC-P
	IRE1 Recombinant Rabbit mAb [PSH10-60] - BSA and Azide free	HA751363	Human, Mouse, Rat, Monkey	WB
	IREB2 Recombinant Mouse mAb [9D1-R] - BSA and Azide free	HA610146	Human, Mouse, Rat	WB
Cell treatment	IRF1 Recombinant Rabbit mAb [SR44-08] - BSA and Azide free	HA750041	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	IRF3 Recombinant Rabbit mAb [SD2062] - BSA and Azide free	HA750283	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	IRF5 Recombinant Rabbit mAb [SN201-05] - BSA and Azide free	HA750277	Human	WB,IHC-P,FC
	IRF7 Recombinant Rabbit mAb [SC0617] - BSA and Azide free	HA750234	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Relative expression	IRF8 Recombinant Rabbit mAb [PSH10-05] - BSA and Azide free	HA751335	Human, Mouse	WB,IF-Cell,FC,IP
Cell treatment	IRF9 Recombinant Rabbit mAb [PSH07-87] - BSA and Azide free	HA751193	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
Cell treatment	IRF9 Recombinant Rabbit mAb [PSH08-20] - BSA and Azide free	HA751211	Human	WB,IF-Cell,IP,ChIP
Cell treatment	IRG1 Recombinant Rabbit mAb [PSH07-31] - BSA and Azide free	HA751152	Mouse	WB,IF-Cell
	Islet 1 Recombinant Rabbit mAb [SC05-64] - BSA and Azide free	HA750211	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-P,FC,IP,IHC-Fr
	ISYNA1 Recombinant Rabbit mAb [PSH07-25] - BSA and Azide free	HA751146	Human, Mouse, Rat	WB,IHC-P
	JAK1 Recombinant Rabbit mAb [JM75-03] - BSA and Azide free	HA750442	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	JAK2 Recombinant Rabbit mAb [SY0245] - BSA and Azide free	HA750118	Human, Mouse, Rat	WB,IF-Cell
	Jarid2 Recombinant Rabbit mAb [PSH01-32] - BSA and Azide free	HA750710	Human, Mouse	IHC-P,IF-Cell
KD	JMJD6 Recombinant Rabbit mAb [PSH15-03]	HA723703	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP,ChIP
	JNK1+JNK2+JNK3 Recombinant Rabbit mAb [SA43-06] - BSA and Azide free	HA750018	Human, Mouse, Rat, Cow, Monkey, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Junctophilin-2 Recombinant Mouse mAb [3-2-R] - BSA and Azide free	HA610170	Human, Mouse, Rat	WB,IHC-P
	JunD Recombinant Rabbit mAb [PSH08-85] - BSA and Azide free	HA751250	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP,ChIP
	Kappa Light Chain Recombinant Mouse mAb [A8H10-R] - BSA and Azide free	HA610141	Human	WB,IHC-P
	KAT1 / HAT1 Recombinant Rabbit mAb [PSH06-13] - BSA and Azide free	HA751051	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue
	KAT3B / p300 Recombinant Rabbit mAb [PSH15-31]	HA723737	Human, Mouse, Rat, Monkey	WB,IF-Cell
	KAT3B / p300 Recombinant Rabbit mAb [PSH15-30]	HA723736	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	KAT6A Recombinant Rabbit mAb [PSH09-88] - BSA and Azide free	HA751322	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,ChIP
Relative expression	KCNA5 / KV1-5 Recombinant Rabbit mAb [PSH12-30] - BSA and Azide free	HA751444	Human, Mouse, Rat	WB
	KDM2B Recombinant Rabbit mAb [PSH05-38] - BSA and Azide free	HA750997	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,IP,ChIP
	KDM3A / JHDM2A Recombinant Rabbit mAb [PSH01-77] - BSA and Azide free	HA750751	Human, Mouse	WB,IHC
	KDM3B Recombinant Rabbit mAb [1A5] - BSA and Azide free	HA750508	Human, Mouse	WB,IHC-P,IF-Cell
	KDM3B Recombinant Rabbit mAb [1A8] - BSA and Azide free	HA750509	Human, Mouse	WB,IHC-P,IF-Cell,FC
	KDM3B Recombinant Rabbit mAb [PSH0-11] - BSA and Azide free	HA750569	Human, Mouse	WB,IF-Cell
	KDM4A Recombinant Rabbit mAb [PSH03-67] - BSA and Azide free	HA750901	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	KDM5B Recombinant Rabbit mAb [PSH03-97] - BSA and Azide free	HA750922	Human	WB,IP
KD	Keap1 Recombinant Rabbit mAb [PSH0-90] - BSA and Azide free	HA750670	Human, Mouse, Rat	WB,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	ketohexokinase Recombinant Rabbit mAb [PSH03-18] - BSA and Azide free	HA750859	Human, Mouse	WB,IF-Cell,IHC-P,FC
KD	ketohexokinase Recombinant Rabbit mAb [PSH03-19] - BSA and Azide free	HA750860	Human, Mouse, Rat	WB,IF-Cell,FC
	Ki67 Recombinant Rabbit mAb [PSH0-02] - BSA and Azide free	HA750570	Human	IHC-P,FC
	KIAA0100 Recombinant Mouse mAb [A0-F3-R] - BSA and Azide free	HA610189	Human	WB,IF-Cell,IHC-P
	KLF4 Recombinant Rabbit mAb [JF98-08] - BSA and Azide free	HA750356	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Klotho Recombinant Rabbit mAb [JM93-76] - BSA and Azide free	HA750445	Human, Mouse	WB,IHC-P
	KMT2C / MLL3 Recombinant Rabbit mAb [PSH02-11] - BSA and Azide free	HA750781	Human, Rat, Mouse	WB
	KMT3B Recombinant Rabbit mAb [1B1] - BSA and Azide free	HA750507	Human	IF-Cell
	KMT6/EZH2 Recombinant Rabbit mAb [PSH04-14] - BSA and Azide free	HA750937	Human, Mouse, Rat, Monkey	WB,IF-Cell
	KPNB1 Recombinant Rabbit mAb [PSH01-59] - BSA and Azide free	HA750735	Human, Mouse, Rat	WB,IP
	KPNB1 Recombinant Rabbit mAb [PSH01-60] - BSA and Azide free	HA750736	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP
	KRAS Recombinant Mouse mAb [A8E6] - BSA and Azide free	HA610008	Human, Mouse, Rat	WB,IHC-P
Relative expression	Kv4.2 / KCND2 Recombinant Mouse mAb [PSH13-21] - BSA and Azide free	HA610272	Mouse, Rat	IHC-Fr,IHC-P,WB
KD	Lactate Dehydrogenase Recombinant Rabbit mAb [SU39-06] - BSA and Azide free	HA750155	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-42] - BSA and Azide free	HA750620	Mouse	WB,IF-Cell
Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-43] - BSA and Azide free	HA750621	Human, Mouse	WB,IHC-P
Cell treatment	LAG-3 Recombinant Rabbit mAb [PSH0-44] - BSA and Azide free	HA750622	Human, Mouse	WB,IHC-P,IF-Cell
Relative expression	LAIR1 Recombinant Mouse mAb [PSH04-50] - BSA and Azide free	HA610201	Human	WB,IHC-P
	Lambda Light Chain Recombinant Mouse mAb [A8G5-R] - BSA and Azide free	HA610118	Human	WB,IF-Cell,IHC-P,FC
	Lambda Light Chain Recombinant Mouse mAb [A8G3-R] - BSA and Azide free	HA610142	Human	WB,IF-Cell,IHC-P,FC
Relative expression	LAMC2 Recombinant Rabbit mAb [PSH07-27] - BSA and Azide free	HA751148	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	Lamin A + Lamin C Recombinant Rabbit mAb [PSH05-20] - BSA and Azide free	HA750991	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	Lamin A + Lamin C Recombinant Mouse mAb [A6F5-R] - BSA and Azide free	HA610149	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KO	Lamin B1 Recombinant Rabbit mAb [SI17-07] - BSA and Azide free	HA750098	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,CUT&Tag-seq
	Lamin B2 Recombinant Mouse mAb [J9-C4-R] - BSA and Azide free	HA610162	Human, Monkey	WB,IF-Cell,IHC-P,IF-Tissue
KD	Lamin B Receptor / LBR Recombinant Rabbit mAb [PSH02-93] - BSA and Azide free	HA750848	Human	WB,IHC-P,IF-Cell,FC
Relative expression	Laminin 2 alpha Recombinant Rabbit mAb [PSH0-94] - BSA and Azide free	HA750674	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
KD	LAMP1 Recombinant Rabbit mAb [JJ0940] - BSA and Azide free	HA750328	Human	WB,IHC-P,IF-Tissue
Cell treatment	LAMP1 Recombinant Rabbit mAb [PSH07-39] - BSA and Azide free	HA751160	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr,IP
Cell treatment	LAMP1 Recombinant Rabbit mAb [PSH05-55] - BSA and Azide free	HA751017	Mouse	WB,IF-Cell,IHC-P,FC,IP,IF-Tissue,IHC-Fr
KD	LAMP2 Recombinant Mouse mAb [D3-D8-D3-R] - BSA and Azide free	HA610072	Human	WB,IHC-P
	LAMP2a Recombinant Rabbit mAb [SA46-01] - BSA and Azide free	HA750016	Human, Mouse, Rat	WB,IHC-P,IP
KD	LAP3 Recombinant Rabbit mAb [PSH01-34] - BSA and Azide free	HA750712	Human, Mouse, Rat	WB,IHC-P
KD	LAP3 Recombinant Rabbit mAb [PSH01-26] - BSA and Azide free	HA750704	Human, Mouse, Rat	WB,IHC-P
	LARP1 Recombinant Rabbit mAb [PSH12-09] - BSA and Azide free	HA751438	Human, Monkey	WB,IF-Cell
Relative expression	LAT Recombinant Rabbit mAb [PSH09-76] - BSA and Azide free	HA751318	Human	WB,IF-Cell,FC
	Late-Onset Alzheimer's Disease Risk Gene Antibody Sampler Kit	HAK21060	Human	WB
	LATS1 / WARTS Recombinant Rabbit mAb [PSH07-20] - BSA and Azide free	HA751141	Human, Mouse, Monkey	WB
Cell treatment, KD	LC3B Recombinant Rabbit mAb [JJ090-6] - BSA and Azide free	HA750319	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,IP
Relative expression	Lck Recombinant Rabbit mAb [PSH09-25] - BSA and Azide free	HA751279	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
	LDHA Recombinant Rabbit mAb [PSH13-58] - BSA and Azide free	HA751477	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	LDHA + LDHB Recombinant Rabbit mAb [PSH12-75] - BSA and Azide free	HA751451	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	LDL Receptor Recombinant Rabbit mAb [SJ0197] - BSA and Azide free	HA750107	Human, Rat	WB,IF-Cell,FC
	Leptin Receptor Recombinant Rabbit mAb [JA73-01] - BSA and Azide free	HA750413	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	LIAS Recombinant Rabbit mAb [PSH04-04] - BSA and Azide free	HA750927	Human, Mouse, Rat	WB,IF-Cell
Relative expression	LIFR Recombinant Rabbit mAb [PSH08-76] - BSA and Azide free	HA751241	Human	WB
Relative expression	LILRB1 Recombinant Rabbit mAb [PSH04-59] - BSA and Azide free	HA750966	Human	WB
	LIMK2 Recombinant Rabbit mAb [PSH08-82] - BSA and Azide free	HA751247	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,IP
Relative expression	LIM Kinase 1 Recombinant Rabbit mAb [PSH08-68] - BSA and Azide free	HA751233	Human, Mouse, Rat, Monkey	WB,IHC-P,FC
KD	LIMPII Recombinant Rabbit mAb [PSH03-25] - BSA and Azide free	HA750865	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell
Relative expression	Lin28A Recombinant Rabbit mAb [JA10-63] - BSA and Azide free	HA750410	Human, Mouse, Rat	WB,IHC-P
	Lipoamide Dehydrogenase Recombinant Rabbit mAb [PSH0-83] - BSA and Azide free	HA750663	Human, Mouse	WB,IHC-P,IF-Cell
Relative expression	Lipoprotein a Recombinant Rabbit mAb [PSH11-29] - BSA and Azide free	HA751408	Human	WB,IHC-P,IF-Tissue
	LIPT1 Recombinant Rabbit mAb [PSH06-54] - BSA and Azide free	HA751082	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	Livin Recombinant Rabbit mAb [PSH12-72] - BSA and Azide free	HA751448	Human	WB,IHC-P,IP
Relative expression	LKB1 Recombinant Rabbit mAb [PSH09-29] - BSA and Azide free	HA751283	Human, Mouse, Rat, Monkey	WB,IHC-P,IP
	LNP Recombinant Mouse mAb [11-1-R] - BSA and Azide free	HA610188	Human, Mouse, Rat	WB,IHC-P
	LPCAT3 Recombinant Rabbit mAb [PSH10-02] - BSA and Azide free	HA751332	Human, Monkey	WB,IF-Cell
Relative expression	Lrp2 / Megalin Recombinant Rabbit mAb [PSH04-52] - BSA and Azide free	HA750959	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	Lrp2 / Megalin Recombinant Rabbit mAb [PSH02-91] - BSA and Azide free	HA750846	Human	IHC-P,IF-Tissue,WB
	LRP4 Recombinant Mouse mAb [PSH10-25] - BSA and Azide free	HA610226	Mouse, Rat	WB
	LRRK2 Recombinant Rabbit mAb [PSH07-96] - BSA and Azide free	HA751198	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IF-Tissue,IHC-Fr
	LSD2 / AOF1 Recombinant Rabbit mAb [PSH06-18] - BSA and Azide free	HA751056	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	LXR alpha Recombinant Rabbit mAb [JA20-38] - BSA and Azide free	HA750414	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	LXR beta Recombinant Rabbit mAb [PSH0-56] - BSA and Azide free	HA750633	Human, Mouse, Rat	WB,IF-Cell
Relative expression	Ly6c Recombinant Rabbit mAb [PSH04-51] - BSA and Azide free	HA750958	Mouse	WB
	LYPD3 Recombinant Rabbit mAb [PS00-59] - BSA and Azide free	HA750615	Human	WB,IHC-P
	LYRIC Recombinant Mouse mAb [C6-C10-R] - BSA and Azide free	HA610068	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Lysozyme Recombinant Mouse mAb [A5A8-R] - BSA and Azide free	HA610104	Human	WB,IHC-P
	Lysozyme Recombinant Rabbit mAb [ST50-02] - BSA and Azide free	HA750170	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP
Relative expression	LYVE1 Recombinant Rabbit mAb [PSH07-67] - BSA and Azide free	HA751179	Mouse	WB,IHC-P
KD	LYVE1 Recombinant Rabbit mAb [JF0979] - BSA and Azide free	HA750343	Human	WB
	m6A Recombinant Rabbit mAb [1G1] - BSA and Azide free	HA750493	Species independent	Dot Blot,ELISA
	M6PR (cation dependent) Recombinant Rabbit mAb [PSH01-06] - BSA and Azide free	HA750681	Human, Mouse, Rat	WB,IHC-P
Relative expression	MAG Recombinant Antibody [PSH02-41] - Guinea pig IgG2 (Chimeric)	HA601510	Mouse, Rat	IHC-Fr,WB
Relative expression	MAG Recombinant Rabbit mAb [PSH02-41] - BSA and Azide free	HA750806	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IP,IHC-Fr
Relative expression	MAG Recombinant Rabbit mAb [PSH02-42] - BSA and Azide free	HA750807	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	MAG Recombinant Antibody [PSH02-41] - Chicken IgY (Chimeric)	HA601509	Mouse, Rat	IHC-Fr,WB
Relative expression	MAG Recombinant Antibody [PSH02-41] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610254	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB,IF-Tissue
	MAGI1 Recombinant Rabbit mAb [PSH10-91] - BSA and Azide free	HA751382	Human, Mouse, Rat	WB
KD	MALT1 Recombinant Rabbit mAb [PSH02-31] - BSA and Azide free	HA750563	Human, Mouse, Rat	WB,IHC-P
	MALT1 Recombinant Rabbit mAb [PSH02-32] - BSA and Azide free	HA750722	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Relative expression	Mammaglobin A Recombinant Mouse mAb [A9A2-R] - BSA and Azide free	HA610045	Human	IHC-P
	Mammal NF-L enzyme-linked immunoassay kit	EH0154	Mammal	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	MAN2A1 Recombinant Rabbit mAb [PSH02-67] - BSA and Azide free	HA750824	Human, Mouse, Rat	WB,IHC-P,FC
	Mannose Receptor(CD206) Recombinant Antibody [PSH07-76] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610286	Mouse, Rat	IHC-Fr,WB
	Mannose Receptor(CD206) Recombinant Antibody [PSH07-76] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610287	Mouse, Rat	IHC-Fr,IF-Tissue,IHC-P,WB,IF-Cell
KD	Mannose Receptor(CD206) Recombinant Rabbit mAb [JF0953] - BSA and Azide free	HA750332	Human, Mouse	WB,IF-Cell,FC
Cell treatment	Mannose Receptor(CD206) Recombinant Rabbit mAb [PSH07-76] - BSA and Azide free	HA751188	Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr,IF-Cell,IP
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610293	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610292	Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell
Relative expression	MAP2 Recombinant Rabbit mAb [PSH08-73] - BSA and Azide free	HA751238	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IHC-Fr,IF-Cell
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610291	Mouse, Rat	IHC-Fr,WB
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Chicken IgY (Chimeric) - BSA and Azide free	HA610301	Mouse, Rat	IHC-Fr,WB
	MAP3K1 Recombinant Rabbit mAb [PSH06-10] - BSA and Azide free	HA751048	Human, Mouse, Rat	WB
	MAP4K1 / HPK1 Recombinant Rabbit mAb [PSH07-24] - BSA and Azide free	HA751145	Human	WB,IHC-P,IP
KD	MAPK6 / ERK3 Recombinant Rabbit mAb [PSH03-39] - BSA and Azide free	HA750873	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	MARCO Recombinant Rabbit mAb [PSH04-31] - BSA and Azide free	HA750940	Human, Mouse	WB,IHC-P
	Mast Cell Tryptase Recombinant Antibody [SC68-07] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610259	0	
	Mast Cell Tryptase Recombinant Rabbit mAb [SC68-07] - BSA and Azide free	HA750229	Human, Mouse	WB,IHC-P,IP
Cell treatment	Matrix Remodeling Antibody Sampler Kit	HAK21072	Human	WB
	MAVS Recombinant Rabbit mAb [PSH01-65] - BSA and Azide free	HA750741	Human	WB,IHC-P,IF-Cell,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	MAVS Recombinant Rabbit mAb [PSH0-30] - BSA and Azide free	HA750604	Human, Mouse	WB,IHC-P,IF-Cell,IP
	MBD1 Recombinant Rabbit mAb [PSH01-78] - BSA and Azide free	HA750752	Human, Mouse, Rat	WB
	MBD1 Recombinant Rabbit mAb [PSH01-79] - BSA and Azide free	HA750753	Mouse, Rat	WB,IF-Cell,FC
KD	MCAK Recombinant Rabbit mAb [PSH01-51] - BSA and Azide free	HA750727	Human, Mouse, Rat	WB,IHC-P
KD	MCL1 Recombinant Rabbit mAb [SI16-04] - BSA and Azide free	HA750093	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
	MCM7 Recombinant Mouse mAb [15E1-R] - BSA and Azide free	HA610161	Human, Mouse, Rat, Monkey	WB
	MCU Recombinant Rabbit mAb [PSH13-74] - BSA and Azide free	HA751487	Human, Mouse, Rat, Monkey	WB,IHC-Fr,IHC-P,IF-Cell,IP
	MD2 Recombinant Rabbit mAb [PSH07-72] - BSA and Azide free	HA751184	Human, Mouse, Rat	WB
	MDA5 Recombinant Rabbit mAb [PSH03-56] - BSA and Azide free	HA750890	Human, Mouse, Rat	WB
	MeCP2 Recombinant Rabbit mAb [PSH14-12] - BSA and Azide free	HA751507	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	MED12 Recombinant Rabbit mAb [PSH06-67] - BSA and Azide free	HA751096	Human, Mouse, Rat	WB
	MEK1/2 Recombinant Rabbit mAb [SR13-07] - BSA and Azide free	HA750042	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK4 / MKK4 Recombinant Rabbit mAb [PSH11-99] - BSA and Azide free	HA751432	Human, Mouse, Rat, Monkey	WB,IHC-P,FC
	MEK6 Recombinant Rabbit mAb [PSH07-86] - BSA and Azide free	HA751194	Human, Mouse, Rat	WB,IP
KD	MELK Recombinant Rabbit mAb [PSH03-55] - BSA and Azide free	HA750889	Human, Mouse	WB,IHC-P
Relative expression, KD	Merlin Recombinant Rabbit mAb [PSH02-52] - BSA and Azide free	HA750815	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
KD	MERTK Recombinant Rabbit mAb [SR29-07] - BSA and Azide free	HA750036	Human	WB,IHC-P,IP,FC,IF-Tissue
Relative expression	Mesothelin Recombinant Rabbit mAb [PSH0-57] - BSA and Azide free	HA750634	Human	WB,IHC-P,IF-Tissue
KD	METTL14 Recombinant Rabbit mAb [PSH07-85] - BSA and Azide free	HA751192	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP
	METTL3 Recombinant Rabbit mAb [JE58-20] - BSA and Azide free	HA750470	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	MFF Recombinant Rabbit mAb [PSH14-36] - BSA and Azide free	HA751518	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	MFGE8 Recombinant Rabbit mAb [PSH05-75] - BSA and Azide free	HA751025	Human	WB
	MGST1 Recombinant Mouse mAb [14H2-R] - BSA and Azide free	HA610130	Human, Rat	WB,IHC-P
Relative expression, KD	MHC class I Recombinant Rabbit mAb [JF10-38] - BSA and Azide free	HA750347	Human	WB,IF-Cell,IHC-P,FC,IF-Tissue
	MHC Class II Recombinant Rabbit mAb [JA10-94] - BSA and Azide free	HA750404	Human	WB,IP,IF-Cell,IF-Tissue,IHC-P
Relative expression	MICA Recombinant Rabbit mAb [PSH13-60] - BSA and Azide free	HA751479	Human	WB,FC,IP
Relative expression	MICA + MICB Recombinant Rabbit mAb [PSH10-88] - BSA and Azide free	HA751380	Human	WB,FC,IP
	MICU2 Recombinant Rabbit mAb [PSH01-84] - BSA and Azide free	HA750758	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
KD	MIF Recombinant Rabbit mAb [PSH06-49] - BSA and Azide free	HA751079	Human, Mouse, Rat	WB,IF-Cell,FC
KD	MiTF Recombinant Rabbit mAb [JF100-01] - BSA and Azide free	HA750362	Human, Mouse, Rat	WB,IF-Cell,FC
	Mitofilin Recombinant Rabbit mAb [PSH01-18] - BSA and Azide free	HA750694	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
KD	Mitofilin Recombinant Rabbit mAb [PSH01-19] - BSA and Azide free	HA750695	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
	MKLP1 Recombinant Rabbit mAb [SY02-74] - BSA and Azide free	HA750129	Human, Mouse, Rat	WB,IF-Cell
	MLD Recombinant Rabbit mAb [PSH02-54] - BSA and Azide free	HA750817	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	MLH1 Recombinant Rabbit mAb [PSH01-74] - BSA and Azide free	HA750749	Human	IHC-P,IF-Tissue,IF-Cell
KO, Relative expression	MLH1 Recombinant Rabbit mAb [SP08-04] - BSA and Azide free	HA750077	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Relative expression, KD	MLKL Recombinant Rabbit mAb [SA40-04] - BSA and Azide free	HA750017	Human, Mouse	WB,IHC-P,IF-Cell,FC
Relative expression	MMP-1 Recombinant Rabbit mAb [PSH09-56] - BSA and Azide free	HA751303	Human	WB,IHC-P
	MMP-11 Recombinant Rabbit mAb [SN74-08] - BSA and Azide free	HA750250	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
KD	MMP-12 Recombinant Rabbit mAb [SR03-23] - BSA and Azide free	HA750052	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MMP-14 Recombinant Rabbit mAb [SJ18-09] - BSA and Azide free	HA750108	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Relative expression	MMP-2 Recombinant Rabbit mAb [PSH11-26] - BSA and Azide free	HA751405	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	MMP-7 Recombinant Rabbit mAb [PSH08-91] - BSA and Azide free	HA751256	Human	WB,IF-Cell,FC,IHC-P
Cell treatment	MMP-9 Recombinant Rabbit mAb [PSH06-33] - BSA and Azide free	HA751070	Human	WB,IHC-P
Cell treatment, KD	MMP-9 Recombinant Rabbit mAb [JA80-73] - BSA and Azide free	HA750420	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC
	Monoamine Oxidase A / MAO-A Recombinant Mouse mAb [5D1-R] - BSA and Azide free	HA610191	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Mono-Methyl Arginine (R*GG) Recombinant Rabbit mAb [PSH07-52] - BSA and Azide free	HA751173	Species independent	WB
	Mouse Adiponectin Recombinant Rabbit mAb [PSH14-83] - BSA and Azide free (Detector)	HA723680	Mouse	ELISA(Det)
	Mouse Adiponectin Recombinant Rabbit mAb [PSH14-82] - BSA and Azide free (Capture)	HA723679	Mouse	ELISA(Cap)
	Mouse Adiponectin Recombinant Rabbit mAb [PSH14-84] - BSA and Azide free (Capture)	HA723682	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse Adiponectin Recombinant Rabbit mAb [PSH14-83] - Detector	HA723681B	Mouse	ELISA(Det),ELISA
	Mouse ApoE, Tag Free Protein	HA211210	Mouse	
	Mouse B7-2/CD86, C-His, Flag Tag Protein	HA210879	Mouse	
	Mouse CD25/IL-2R alpha, Tag Free Protein	HA211117	Mouse	
	Mouse CD45R Recombinant Rat mAb [PSH10-80] - BSA and Azide free	HA610229	Mouse	IHC-P,FC,IF-Tissue
	Mouse CD68/SR-D1 enzyme-linked immunoassay kit	EM0042	Mouse	ELISA
	Mouse Complement Factor D, C-His Tag Protein	HA211116	Mouse	
	Mouse EGF, Tag Free Protein	HA211220	Mouse	
	Mouse FGF-21 Recombinant Rabbit mAb [PSH15-26] - BSA and Azide free (Capture)	HA723731	Mouse	ELISA(Cap)
	Mouse FGF-21 Recombinant Rabbit mAb [PSH15-27] - BSA and Azide free (Detector)	HA723732	Mouse	ELISA(Det)
	Biotin Conjugated Mouse FGF-21 Recombinant Rabbit mAb [PSH15-27] - Detector	HA723733B	Mouse	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse IFN gamma Recombinant Rabbit mAb [PSH14-81] - Detector	HA723678B	Mouse	ELISA(Det),ELISA
	Mouse IFN gamma Recombinant Rabbit mAb [PSH14-81] - BSA and Azide free (Detector)	HA723677	Mouse	ELISA(Det)
	Mouse IFN gamma Recombinant Rabbit mAb [PSH14-80] - BSA and Azide free (Capture)	HA723676	Mouse	ELISA(Cap)
	Mouse IL-17E/IL-25 Recombinant Rabbit mAb [PSH15-34] - BSA and Azide free (Detector)	HA725044	Mouse	ELISA(Det)
	Mouse IL-17E/IL-25 Recombinant Rabbit mAb [PSH15-33] - BSA and Azide free (Capture)	HA725043	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-17E/IL-25 Recombinant Rabbit mAb [PSH15-34] - Detector	HA725045B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse IL-1RA Recombinant Rabbit mAb [PSH14-99] - Detector	HA723699B	Mouse	ELISA(Det)
	Mouse IL-1RA Recombinant Rabbit mAb [PSH15-00] - BSA and Azide free (Capture)	HA723700	Mouse	ELISA(Cap)
	Mouse IL-1RA Recombinant Rabbit mAb [PSH14-98] - BSA and Azide free (Capture)	HA723697	Mouse	ELISA(Cap)
	Mouse IL-1RA Recombinant Rabbit mAb [PSH14-99] - BSA and Azide free (Detector)	HA723698	Mouse	ELISA(Det)
	Mouse IL-21 Recombinant Rabbit mAb [PSH14-57] - BSA and Azide free (Detector)	HA723652	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-21 Recombinant Rabbit mAb [PSH14-56] - Detector	HA723651B	Mouse	ELISA(Det),ELISA
	Mouse IL-21 Recombinant Rabbit mAb [PSH14-56] - BSA and Azide free (Detector)	HA723650	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-21 Recombinant Rabbit mAb [PSH14-57] - Detector	HA723653B	Mouse	ELISA(Det),ELISA
	Mouse IL-21 Recombinant Rabbit mAb [PSH04-75] - BSA and Azide free (Capture)	HA723649	Mouse	ELISA(Cap)
	Mouse IL-3 Recombinant Rabbit mAb [PSH14-85] - BSA and Azide free (Capture)	HA723683	Mouse	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse IL-3 Recombinant Rabbit mAb [PSH14-86] - Detector	HA723685B	Mouse	ELISA(Det),ELISA
	Mouse IL-3 Recombinant Rabbit mAb [PSH14-86] - BSA and Azide free (Detector)	HA723684	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-3 Recombinant Rabbit mAb [PSH14-87] - Detector	HA723687B	Mouse	ELISA(Det),ELISA
	Mouse IL-3 Recombinant Rabbit mAb [PSH14-87] - BSA and Azide free (Detector)	HA723686	Mouse	ELISA(Det)
	Mouse Interferon alpha 1 Recombinant Rabbit mAb [PSH14-76] - BSA and Azide free (Detector)	HA723670	Mouse	ELISA(Det)
	Mouse Interferon alpha 1 Recombinant Rabbit mAb [PSH14-77] - BSA and Azide free (Capture)	HA723672	Mouse	ELISA(Cap)
	Mouse Interferon alpha 1 Recombinant Rabbit mAb [PSH14-75] - BSA and Azide free (Capture)	HA723669	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse Interferon alpha 1 Recombinant Rabbit mAb [PSH14-76] - Detector	HA723671B	Mouse	ELISA(Det),ELISA
	Mouse MCSF/CSF-1, C-His Tag (33-187) Protein	HA211182	Mouse	
	HRP Conjugated Mouse Pro-Collagen I alpha 1 Recombinant Rabbit mAb [PSH07-99] - Detector	HA722936H	Mouse	ELISA(Det)
	HRP Conjugated Mouse/Rat IGF1 Recombinant Rabbit mAb [PSH12-43] - Detector	HA725039H	Mouse, Rat	ELISA(Det),ELISA
	Mouse TSLP, Tag Free Protein	HA211088	Mouse	
KD	MPG / AAG Recombinant Rabbit mAb [PSH06-34] - BSA and Azide free	HA751073	Human	WB,IHC-P,IF-Cell,FC
	MPP1 Recombinant Rabbit mAb [PSH02-69] - BSA and Azide free	HA750826	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Ms4a10 Recombinant Rabbit mAb [PSH0-71] - BSA and Azide free	HA750650	Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	MSH2 Recombinant Mouse mAb [10G2-R] - BSA and Azide free	HA610179	Human, Mouse, Rat	WB,IHC-P
Relative expression	MSH3 Recombinant Rabbit mAb [PSH07-41] - BSA and Azide free	HA751162	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	MSH6 Recombinant Mouse mAb [A4C4-R] - BSA and Azide free	HA610061	Human, Mouse	WB,IHC-P
	MSI2 Recombinant Rabbit mAb [PSH01-17] - BSA and Azide free	HA750693	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	MSK1 Recombinant Rabbit mAb [PSH13-72] - BSA and Azide free	HA751485	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	MTCO1 Recombinant Rabbit mAb [PSH07-50] - BSA and Azide free	HA751171	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
	MTERFD1 Recombinant Mouse mAb [A8-A10-R] - BSA and Azide free	HA610112	Human, Mouse	WB
	MTF1 Recombinant Rabbit mAb [PSH06-26] - BSA and Azide free	HA751063	Human	WB,IF-Cell,IHC-P,FC
KD	MTH1 Recombinant Rabbit mAb [PSH02-40] - BSA and Azide free	HA750805	Human, Mouse, Rat	WB,IHC-P,FC
	mTOR Recombinant Rabbit mAb [SU30-00] - BSA and Azide free	HA750153	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell
Relative expression, KD	MUC1 Recombinant Rabbit mAb [SN06-80] - BSA and Azide free	HA750245	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	MUC2 Recombinant Rabbit mAb [JA10-92] - BSA and Azide free	HA750401	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	MUC5AC Recombinant Mouse mAb [A10G5-R] - BSA and Azide free	HA610064	Human, Mouse	IHC-P
	MUC5AC Recombinant Mouse mAb [A10G6-R] - BSA and Azide free	HA610065	Human	IHC-P
Relative expression	Mucolipin-1 Recombinant Rabbit mAb [PSH09-01] - BSA and Azide free	HA751262	Human, Mouse, Rat	WB,IHC-P
	MUL1 Recombinant Rabbit mAb [PSH15-17]	HA723721	Human, Mouse, Rat	WB,IHC-P
	MUSK Recombinant Rabbit mAb [PSH02-29] - BSA and Azide free	HA750797	Human, Mouse, Rat	WB,IHC-P
	MUSK Recombinant Rabbit Multiclonal Antibody [PSH03-51] - BSA and Azide free	HA750885	Human, Mouse, Rat	WB,IHC-P
	MUSK Recombinant Rabbit mAb [PSH02-28] - BSA and Azide free	HA750796	Human, Mouse, Rat	WB,IHC-P
	MUSK Recombinant Rabbit mAb [PSH05-68] - BSA and Azide free	HA751022	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Cell treatment	MX1 Recombinant Rabbit mAb [PSH06-03] - BSA and Azide free	HA751043	Human	WB,IF-Cell,IHC-P
	MYCBP2 Recombinant Rabbit mAb [PSH03-68] - BSA and Azide free	HA750902	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Myc tag Recombinant Mouse mAb [A3C8-R-2] - BSA and Azide free	HA610018	Species independent	ELISA,WB,IP,IF-Cell
	Myc tag Recombinant Mouse mAb [PSH14-94]	HA601513	Species independent	WB,IF-Cell,IP
	MyD88 Recombinant Rabbit mAb [PSH09-74] - BSA and Azide free	HA751316	Human, Mouse, Rat, Monkey	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	MyD88 Recombinant Rabbit mAb [PSH05-58] - BSA and Azide free	HA751021	Human, Monkey	WB,IF-Cell,IHC-P,FC,IP
	MyD88 Recombinant Rabbit mAb [SC65-04] - BSA and Azide free	HA750231	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	Myelin Basic Protein Recombinant Rabbit mAb [JF0943] - BSA and Azide free	HA750337	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression	Myelin PLP Recombinant Antibody [PSH08-21] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610303	Mouse, Rat	IHC-Fr,WB
Relative expression	Myelin PLP Recombinant Antibody [PSH08-21] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610247	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Myelin PLP Recombinant Antibody [PSH08-21] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610304	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Myelin PLP Recombinant Rabbit mAb [PSH08-21] - BSA and Azide free	HA751212	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
	Myeloperoxidase Recombinant Mouse mAb [A1F2-R] - BSA and Azide free	HA610124	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
	MYL9 Recombinant Rabbit mAb [PSH02-07] - BSA and Azide free	HA750777	Human, Mouse, Rat	WB,FC
	MyoD1 Recombinant Rabbit mAb [JM10-72] - BSA and Azide free	HA750395	Human	WB,IP,FC
	Myosin Light Chain 2 Recombinant Rabbit mAb [SN67-09] - BSA and Azide free	HA750244	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IP
	MYST Antibody Sampler Kit	HAK21134	Human	WB
	MyT1L Recombinant Rabbit mAb [PSH09-71] - BSA and Azide free	HA751313	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IHC-Fr
	Myt1 / MTF1 Recombinant Rabbit mAb [PSH12-97] - BSA and Azide free	HA751452	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	N6-mA Methyltransferase Antibody Sampler Kit	HAK21142	Human	WB
KD	NADPH oxidase 4/NOX4 Recombinant Rabbit mAb [SY0214] - BSA and Azide free	HA750122	Human, Rat, Mouse	WB,IF-Cell,IHC-P,IP,FC
	NAK / TBK1 Recombinant Rabbit mAb [PSH11-69] - BSA and Azide free	HA751425	Human, Mouse, Rat	WB,IF-Cell,FC,IP
Relative expression	Nanog Recombinant Rabbit mAb [SC05-70] - BSA and Azide free	HA750204	Human, Mouse	WB,IF-Tissue,FC,IHC-P,ChIP
Cell treatment	Nanos3 Recombinant Rabbit mAb [PSH09-65] - BSA and Azide free	HA751307	Human	WB,IF-Cell,IHC-P
Relative expression	NAPSIN A Recombinant Mouse mAb [A9B12-R] - BSA and Azide free	HA610117	Human	IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	NAPSIN A Recombinant Mouse mAb [A9C1-R] - BSA and Azide free	HA610175	Human	IHC-P
Relative expression, KD	N Cadherin Recombinant Rabbit mAb [SY02-46] - BSA and Azide free	HA750119	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	N Cadherin Recombinant Antibody [SY02-46] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610298	Human, Mouse, Rat	IHC-Fr,IHC-P
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53] - BSA and Azide free	HA751105	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	NCAM1 / CD56 Recombinant Rabbit mAb [JF1021] - BSA and Azide free	HA750346	Human, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-82] - BSA and Azide free	HA751104	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue
	NCOA4 Recombinant Rabbit mAb [PSH05-45] - BSA and Azide free	HA751004	Human, Mouse, Rat	IHC-P,IF-Tissue
KD	NCOA4 Recombinant Rabbit mAb [PSH05-44] - BSA and Azide free	HA751003	Human	WB
	NDP52 Recombinant Rabbit mAb [PSH02-33] - BSA and Azide free	HA750623	Human, Mouse	WB,FC
	NDUFA9 Recombinant Rabbit mAb [PSH14-08] - BSA and Azide free	HA751503	Human, Mouse, Rat, Monkey	WB,IP
	NDUFA9 Recombinant Rabbit mAb [PSH14-07] - BSA and Azide free	HA751502	Human, Mouse, Rat, Monkey	WB,IHC-P,IP
	NDUFS1 Recombinant Rabbit mAb [PSH14-45] - BSA and Azide free	HA751524	Human, Mouse, Rat, Monkey	WB,IHC-P
	NDUFS2 Recombinant Rabbit mAb [PSH03-26] - BSA and Azide free	HA750866	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
KD	NDUFV2 Recombinant Rabbit mAb [PSH02-86] - BSA and Azide free	HA750839	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Relative expression	Nectin 2 Recombinant Rabbit mAb [PSH11-54] - BSA and Azide free	HA751416	Human	WB,IF-Cell,FC
Relative expression	Nectin-4 Recombinant Rabbit mAb [PSH09-73] - BSA and Azide free	HA751315	Human	WB,IHC-P
Cell treatment	NEMP1 Recombinant Rabbit mAb [PSH06-66] - BSA and Azide free	HA751095	Mouse, Rat	WB,IF-Cell
Relative expression, Cell treatment	Nephrin Recombinant Rabbit Multiclonal Antibody [PSH07-70] - BSA and Azide free	HA751182	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression, KD	Nestin Recombinant Rabbit mAb [PSH07-89] - BSA and Azide free	HA751196	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IF-Tissue,IP
	Nestin Recombinant Antibody [PSH07-89] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610246	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	Nestin Recombinant Antibody [PSH07-89] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610302	Human, Mouse, Rat	IHC-Fr,WB
	NeuN Recombinant Antibody [SR45-07] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610237	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
Relative expression	NeuN Recombinant Antibody [SR45-07] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610315	Mouse, Rat	IHC-Fr,WB
	NeuN Recombinant Antibody [SR45-07] - Chicken IgY (Chimeric) - BSA and Azide free	HA610314	Mouse	IHC-Fr
	NeuN Recombinant Mouse mAb [PD01-45] - BSA and Azide free	HA610031	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IHC-Fr
Relative expression	NeuN Recombinant Antibody [SR45-07] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610316	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	Neurogranin Recombinant Rabbit mAb [PSH02-96] - BSA and Azide free	HA750851	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression, KD	Neuropilin 1 Recombinant Rabbit mAb [ST05-30] - BSA and Azide free	HA750189	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Neutrophil Elastase Recombinant Rabbit mAb [JF098-6] - BSA and Azide free	HA750360	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
Relative expression	NFAT1 Recombinant Rabbit mAb [JA11-08] - BSA and Azide free	HA750405	Human	WB,IF-Tissue,IHC-P
Relative expression	NF-H Recombinant Rabbit mAb [JF99-07] - BSA and Azide free	HA750357	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
KD	NF-kB p105 / p50 Recombinant Rabbit mAb [SZ20-01] - BSA and Azide free	HA750060	Human, Mouse, Rat	WB,IHC-P,ICC,FC,IF-Tissue
KO, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [PSH0-27] - BSA and Azide free	HA750601	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
KO, KD, Cell treatment	NF-kB p65 Recombinant Rabbit mAb [SZ10-04] - BSA and Azide free	HA750058	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Relative expression	NF-L Recombinant Rabbit mAb [PS02-10] - BSA and Azide free	HA750678	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,IF-Tissue,IHC-Fr
Relative expression	NF-L Recombinant Antibody [PS02-10] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610249	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
Relative expression	NF-L Recombinant Antibody [PS02-10] - Chicken IgY (Chimeric) - BSA and Azide free	HA610306	Mouse, Rat	IHC-Fr,WB
Relative expression	NF-L Recombinant Antibody [PS02-10] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610282	Human, Mouse, Rat	IHC-Fr,IHC-P,WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	NF-L Recombinant Rabbit mAb [PSH0-77] - BSA and Azide free	HA750657	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
	NF-L Recombinant Antibody [PS02-10] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610297	Mouse, Rat	IHC-Fr
	NF-M Recombinant Antibody [JM11-20] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610278	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	NF-M Recombinant Antibody [JM11-20] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610279	Mouse, Rat	IHC-Fr,WB
	NF-M Recombinant Antibody [JM11-20] - Chicken IgY (Chimeric)	HA601507	Mouse	IHC-Fr
	NF-M Recombinant Antibody [JM11-20] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610248	Human, Mouse, Rat, Cynomolgus monkey	IHC-Fr,IHC-P,WB
	NF-κB Family Antibody Sampler Kit	HAK21141	Human	WB
Relative expression	NG2 Recombinant Rabbit mAb [JM10-13] - BSA and Azide free	HA750372	Human, Mouse, Rat	WB,IHC-P
Relative expression	NG2 Recombinant Rabbit mAb [PSH09-98] - BSA and Azide free	HA751329	Human	WB,IF-Cell,FC
	NGF Recombinant Rabbit mAb [SI79-01] - BSA and Azide free	HA750099	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
	Niemann Pick C2 Recombinant Rabbit mAb [PSH05-52] - BSA and Azide free	HA751011	Human, Mouse, Rat	WB,IHC-P
	Ninein Recombinant Rabbit mAb [PSH03-23] - BSA and Azide free	HA750863	Human, Mouse	WB
	Ninein Recombinant Rabbit mAb [PSH03-46] - BSA and Azide free	HA750880	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	NK1.1 Recombinant Rabbit mAb [PSH08-47] - BSA and Azide free	HA751222	Mouse	WB,IHC-P
	NK1.1 Recombinant Rabbit mAb [PSH08-46] - BSA and Azide free	HA751221	Mouse	FC
Cell treatment	NLRP3 Recombinant Rabbit mAb [SC06-23] - BSA and Azide free	HA750236	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	NM23 Recombinant Mouse mAb [12A1-R] - BSA and Azide free	HA610070	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	NM23 Recombinant Mouse mAb [13C2-R] - BSA and Azide free	HA610180	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KD	NM23 Recombinant Mouse mAb [13C1-R] - BSA and Azide free	HA610113	Human, Mouse, Rat	WB,IF-Cell
Relative expression	NMDAR1 Recombinant Rabbit mAb [JM11-26] - BSA and Azide free	HA750391	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,IHC-Fr
Relative expression	NMDAR2B Recombinant Rabbit mAb [PSH05-46] - BSA and Azide free	HA751005	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Noggin Recombinant Rabbit mAb [PSH10-62] - BSA and Azide free	HA751365	Human, Mouse, Rat	WB,IHC-P
	non-muscle Myosin IIA Recombinant Rabbit mAb [PSH14-10] - BSA and Azide free	HA751505	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
KD	Notch 1 Recombinant Rabbit mAb [SJ205] - BSA and Azide free	HA750111	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Notch2 Recombinant Rabbit mAb [PSH0-55] - BSA and Azide free	HA750632	Human, Mouse, Rat	WB
	Notch2 Recombinant Rabbit mAb [PSH0-54] - BSA and Azide free	HA750631	Human, Mouse, Rat	WB
	NOX2 / gp91phox Recombinant Rabbit mAb [SN07-16] - BSA and Azide free	HA750255	Human, Mouse, Rat	WB
Relative expression	NOX2 / gp91phox Recombinant Rabbit mAb [PSH10-59] - BSA and Azide free	HA751362	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IP
	NPHS2 Recombinant Rabbit mAb [JB51-33] - BSA and Azide free	HA750460	Human, Mouse, Rat	WB,IHC-P,FC,IF-Tissue
Relative expression, KD	NQO1 Recombinant Mouse mAb [13A1-R] - BSA and Azide free	HA610087	Human	WB
KD	NQO1 Recombinant Rabbit mAb [JF440-1] - BSA and Azide free	HA750349	Human, Mouse, Rat	WB,IP,FC,IHC-P
	NRF1 Recombinant Rabbit mAb [JM89-63] - BSA and Azide free	HA750444	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Cell treatment	Nrf2 Recombinant Rabbit mAb [PSH11-20] - BSA and Azide free	HA751401	Human, Mouse, Rat	WB,IF-Cell,ChIP
KD	NSDHL Recombinant Rabbit mAb [PSH02-37] - BSA and Azide free	HA750802	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
KD	NSDHL Recombinant Rabbit mAb [PSH02-36] - BSA and Azide free	HA750801	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
KD	NSE Recombinant Rabbit mAb [SC06-28] - BSA and Azide free	HA750239	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P
	NSUN2 Recombinant Rabbit mAb [PSH05-05] - BSA and Azide free	HA750978	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,IF-Tissue,IP
KD	NSUN3 Recombinant Rabbit mAb [PSH01-97] - BSA and Azide free	HA750771	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	NSUN4 Recombinant Mouse mAb [4G3-R] - BSA and Azide free	HA610129	Human	WB,IF-Cell
	NSUN5 Recombinant Rabbit mAb [PSH05-07] - BSA and Azide free	HA750980	Human, Mouse, Rat	WB
	NSUN6 Recombinant Rabbit mAb [PSH02-80] - BSA and Azide free	HA750835	Human, Mouse, Rat	WB
	NUR77 Recombinant Rabbit mAb [JM59-11] - BSA and Azide free	HA750397	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Nutm1 Recombinant Rabbit mAb [PSH02-78] - BSA and Azide free	HA750833	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
Cell treatment, Relative expression	Nutm1 Recombinant Rabbit mAb [PSH03-36] - BSA and Azide free	HA750870	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Nutm1 Recombinant Rabbit mAb [PSH09-09] - BSA and Azide free	HA750723	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
KD	Oct4 Recombinant Mouse mAb [A3-C11-R] - BSA and Azide free	HA610085	Human	WB,IF-Cell
	OGT Recombinant Rabbit mAb [JB44-39] - BSA and Azide free	HA750459	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	OLFM4 Recombinant Rabbit mAb [PSH01-37] - BSA and Azide free	HA750715	Human	WB,IHC-P
	OLFM4 Recombinant Rabbit mAb [PSH01-36] - BSA and Azide free	HA750714	Human	WB,IHC-P
	OPA1 Recombinant Rabbit mAb [JM81-35] - BSA and Azide free	HA750446	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC
Cell treatment	OPA1 Recombinant Rabbit mAb [PSH06-39] - BSA and Azide free	HA751078	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	Osteopontin Recombinant Rabbit mAb [PSH09-23] - BSA and Azide free	HA751277	Human	IHC-Fr,IHC-P,IF-Tissue
Cell treatment	Osteopontin Recombinant Rabbit mAb [PSH09-24] - BSA and Azide free	HA751278	Human, Mouse	WB,IF-Cell,FC
	Osteopontin Recombinant Rabbit mAb [PSH14-11] - BSA and Azide free	HA751506	Mouse, Rat	IHC-P,IHC-Fr,WB,IF-Cell
Cell treatment	Otx1 + Otx2 Recombinant Rabbit mAb [PSH01-62] - BSA and Azide free	HA750738	Human, Rat, Mouse	WB,IHC-P,IF-Cell,FC
	Otx2 Recombinant Antibody [PSH01-61] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610239	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
Cell treatment	Otx2 Recombinant Rabbit mAb [PSH01-61] - BSA and Azide free	HA750737	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC
Relative expression	Oxytocin-neurophysin 1 / OXT Recombinant Rabbit mAb [PSH08-24] - BSA and Azide free	HA751215	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	Oxytocin-neurophysin 1 / OXT Recombinant Antibody [PSH08-24] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610265	Mouse, Rat, Human, Cynomolgus monkey, Pig	IHC-Fr,IF-Tissue
	Oxytocin-neurophysin 1 / OXT Recombinant Antibody [PSH08-24] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610266	Mouse, Rat, Human, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Tissue
KD	p16INK4a Recombinant Rabbit mAb [SU0702] - BSA and Azide free	HA750158	Human	WB,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	p16INK4a Recombinant Rabbit mAb [PSH15-40]	HA723746	Mouse	WB,IHC-P,IP
	p21 Recombinant Mouse mAb [A8C11-R] - BSA and Azide free	HA610098	Human	IHC-P
	p21 Recombinant Rabbit mAb [PSH03-93] - BSA and Azide free	HA750918	Mouse	WB,IF-Cell,IHC-P,FC
	p27 KIP 1 Recombinant Rabbit mAb [SU37-04] - BSA and Azide free	HA750157	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	P2Y12 Recombinant Rabbit mAb [PSH04-48] - BSA and Azide free	HA750957	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr,IF-Tissue
	P2Y12 Recombinant Antibody [PSH04-48] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610231	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
KD, Cell treatment	p38 alpha / MAPK14 Recombinant Rabbit mAb [SR43-04] - BSA and Azide free	HA750040	Human, Mouse, Rat	WB,IF-Cell,FC
	p38 alpha/MAPK14 + p38 beta/MAPK11 + p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-60] - BSA and Azide free	HA750894	Human, Mouse, Rat, Monkey	WB
	p38 beta/MAPK11 + p38 gamma/MAPK12 + p38 delta/MAPK13 Recombinant Rabbit mAb [PSH03-59] - BSA and Azide free	HA750893	Human, Mouse, Rat, Monkey	WB
	p38 gamma/MAPK12 Recombinant Rabbit mAb [PSH03-58] - BSA and Azide free	HA750892	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	p38 MAPK Recombinant Rabbit Multiclonal Antibody [PSH03-72] - BSA and Azide free	HA750906	Human, Mouse, Rat, Monkey	WB
	p47-phox Recombinant Mouse mAb [A6B10-R] - BSA and Azide free	HA610139	Human	WB
Cell treatment, Relative expression, KD	p53 Recombinant Mouse mAb [6B5-R] - BSA and Azide free	HA610187	Human, Mouse, Monkey	WB,IF-Cell
	p53 Recombinant Rabbit mAb [PD01-32] - BSA and Azide free	HA750545	Human	IHC-P
Relative expression, KD, Cell treatment	p53 Recombinant Rabbit mAb [SA39-07] - BSA and Azide free	HA750005	Human	WB,IF-Cell,ChIP
KD	p53 Recombinant Mouse mAb [6B2-R] - BSA and Azide free	HA610194	Human, Monkey	WB,IF-Cell
Cell treatment	p53 Recombinant Rabbit mAb [PSH03-95] - BSA and Azide free	HA750920	Mouse	WB,IF-Cell,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	p53 (acetyl K370) Recombinant Rabbit mAb [SR40-09] - BSA and Azide free	HA750047	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [PSH06-88] - BSA and Azide free	HA751112	Human	WB,IF-Cell
Cell treatment	p53 (acetyl K382) Recombinant Rabbit mAb [PSH06-87] - BSA and Azide free	HA751111	Human, Mouse	WB,IF-Cell
	p62/KEAP1/NRF2 Pathway Antibody Sampler Kit	HAK21140	Human	WB
Relative expression	p63 Recombinant Rabbit mAb [SC06-31] - BSA and Azide free	HA750218	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	p70 S6 Kinase Recombinant Rabbit mAb [PSH13-50] - BSA and Azide free	HA751469	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	p75 NGF Receptor Recombinant Antibody [PSH13-26] - Rabbit IgG (Chimeric) - BSA and Azide free	HA751462	Mouse	IHC-Fr
KD	p75 NGF Receptor Recombinant Rabbit mAb [SA39-02] - BSA and Azide free	HA750014	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	p75 NGF Receptor Recombinant Rat mAb [PSH13-26] - BSA and Azide free	HA610274	Mouse, Rat	IHC-Fr,WB
	PABP Recombinant Rabbit mAb [PSH05-51] - BSA and Azide free	HA751010	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC,IF-Tissue
	PACT (PKR activating protein) / PRKRA Recombinant Rabbit mAb [PSH02-15] - BSA and Azide free	HA750785	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
KD	PACT (PKR activating protein) / PRKRA Recombinant Rabbit mAb [PSH02-14] - BSA and Azide free	HA750784	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
Cell treatment	PADI4 Recombinant Rabbit mAb [PSH01-33] - BSA and Azide free	HA750711	Human, Mouse, Rat	WB,IHC-P
	Pan-Actin Recombinant Mouse mAb [A2E1-R] - BSA and Azide free	HA610160	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	Pan-Cadherin Recombinant Rabbit mAb [ST54-01] - BSA and Azide free	HA750190	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,IHC-Fr
Relative expression	Pancreatic Polypeptide Recombinant Rabbit mAb [PSH0-84] - BSA and Azide free	HA750664	Human, Mouse, Rat	IHC-P,IHC-Fr,IF-Tissue
	pan Cytokeratin Recombinant Antibody [PDH09-10] - Rat IgG1 (Chimeric)	HA601519	Human	IHC-P
Cell treatment	Pan-Lactyl-lysine Recombinant Rabbit mAb [PSH03-74] - BSA and Azide free	HA750908	Species independent	WB,IHC-P,IF-Cell,Dot Blot
Cell treatment	Pan-Lactyl-lysine Recombinant Rabbit mAb [PSH03-73] - BSA and Azide free	HA750907	Species independent	WB,IHC-P,IF-Cell,Dot Blot
	Parkin Recombinant Rabbit mAb [PSH08-11] - BSA and Azide free	HA751206	Human, Mouse, Rat, Monkey, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Parkin Recombinant Rabbit mAb [JF82-09] - BSA and Azide free	HA750354	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cell treatment, KD	PARP1 Recombinant Rabbit mAb [SU03-68] - BSA and Azide free	HA750154	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	Parvalbumin Recombinant Rabbit mAb [JM100-08] - BSA and Azide free	HA750371	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IP
Cell treatment	Patched / PTCH1 Recombinant Rabbit mAb [PSH03-53] - BSA and Azide free	HA750887	Human	WB
	Pax2 Recombinant Rabbit mAb [JJ082-08] - BSA and Azide free	HA750306	Human, Rat	WB,IHC-P
Relative expression	PAX8 Recombinant Mouse mAb [6G5-R] - BSA and Azide free	HA610145	Human, Mouse, Monkey	WB,IHC-P,IF-Cell
	Paxillin Recombinant Rabbit mAb [SY23-02] - BSA and Azide free	HA750114	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
KD	PBR Recombinant Rabbit mAb [SA90-03] - BSA and Azide free	HA750010	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	PBX1 Recombinant Rabbit mAb [PSH0-24] - BSA and Azide free	HA750598	Human, Mouse, Rat	WB,IHC
	PCB Recombinant Rabbit mAb [PSH01-56] - BSA and Azide free	HA750732	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
KD	PCB Recombinant Rabbit mAb [PSH01-55] - BSA and Azide free	HA750731	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
KD	PCNA Recombinant Rabbit mAb [SY12-07] - BSA and Azide free	HA750084	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	PCNA Recombinant Mouse mAb [A6-G11-R] - BSA and Azide free	HA610056	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,IP
	PD1 Recombinant Rabbit mAb [SJ01-91] - BSA and Azide free	HA750103	Human	WB,IHC-P
	PDCD6 Recombinant Mouse mAb [A6A7-R] - BSA and Azide free	HA610097	Human	WB,IF-Cell,FC
	PDE10A Recombinant Rabbit mAb [PSH03-61] - BSA and Azide free	HA750895	Human, Mouse, Rat	WB,IHC-P,FC
	PDGFR beta Recombinant Rabbit mAb [SY10-08] - BSA and Azide free	HA750082	Human, Mouse, Rat	WB,IHC-P,IP,FC,IF-Tissue,IF-Cell,IHC-Fr
KD	PDIA3 / ERp57 Recombinant Rabbit mAb [PSH08-79] - BSA and Azide free	HA751244	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
KD	PDK1 Recombinant Rabbit mAb [JA67-30] - BSA and Azide free	HA750419	Human, Mouse, Rat	WB,IHC-P,IP,FC,IF-Tissue
	PD-L1 Recombinant Rabbit mAb [JJ08-95] - BSA and Azide free	HA750309	Human	WB,IHC-P
Cell treatment	PD-L1 Recombinant Rabbit mAb [PSH08-56] - BSA and Azide free	HA751231	Cat	WB,IF-Cell,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PD-L1 Recombinant Rabbit mAb [PD00-97] - BSA and Azide free	HA751015	Human	IHC-P
	PD-L1 Recombinant Rabbit Multiclonal Antibody [PSH13-23] - BSA and Azide free	HA751459	Human, Mouse	WB
Cell treatment	PD-L1 Recombinant Rabbit mAb [PSH04-80] - BSA and Azide free	HA750974	Mouse	WB,IHC-P
Cell treatment	PD-L1 Recombinant Rabbit mAb [PSH08-27] - BSA and Azide free	HA751218	Dog	WB,IF-Cell,IHC-P
	Pericentrin Recombinant Rabbit mAb [PSH0-66] - BSA and Azide free	HA750643	Human	WB
	Pericentrin Recombinant Rabbit mAb [PSH01-66] - BSA and Azide free	HA750742	Human	WB,IF-Cell
	Peripherin Recombinant Antibody [PSH10-70] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610245	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB,IF-Cell
	Peripherin Recombinant Rabbit mAb [PSH10-70] - BSA and Azide free	HA751372	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Cell
	PERK Recombinant Rabbit mAb [PSH0-81] - BSA and Azide free	HA750661	Human	WB
KD	Peroxisredoxin 2 Recombinant Mouse mAb [7F2-R] - BSA and Azide free	HA610123	Human	WB
	Peroxisredoxin 2 Recombinant Mouse mAb [7F5-R] - BSA and Azide free	HA610082	Human	WB,IF-Cell,IHC-P
KD	Peroxisredoxin 4 Recombinant Rabbit mAb [PSH03-40] - BSA and Azide free	HA750874	Human	WB,IF-Cell,IHC-P
	Peroxisredoxin 5 Recombinant Rabbit mAb [PSH03-90] - BSA and Azide free	HA750915	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
KD	Peroxisredoxin 6 Recombinant Mouse mAb [7G1-R] - BSA and Azide free	HA610174	Human, Mouse, Rat	WB
Cell treatment	Pet1 Recombinant Rabbit mAb [PSH07-77] - BSA and Azide free	HA751189	Human	WB,IF-Cell,IHC-P
	PGAM5 Recombinant Rabbit mAb [PSH11-68] - BSA and Azide free	HA751424	Human, Mouse, Rat	WB
KD	P Glycoprotein Recombinant Rabbit mAb [SN06-42] - BSA and Azide free	HA750249	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell
Relative expression, KD	PGP9.5 Recombinant Rabbit mAb [JM10-59] - BSA and Azide free	HA750375	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
Relative expression	PGP9.5 Recombinant Antibody [JM10-59] - Rat IgG1 (Chimeric)	HA601515	Human, Mouse, Rat	IHC-Fr,IF-Tissue,IHC-P,WB,IF-Cell
Relative expression	PGP9.5 Recombinant Antibody [JM10-59] - Mouse IgG1 (Chimeric)	HA601514	Human, Mouse, Rat	IHC-Fr,IF-Tissue,IHC-P,WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PHACTR3 Recombinant Rabbit mAb [PSH06-09] - BSA and Azide free	HA751047	Human	WB,IF-Cell
KD	PHD2 Recombinant Rabbit mAb [PSH12-74] - BSA and Azide free	HA751450	Human	WB
KD	PHF6 Recombinant Rabbit mAb [PSH02-35] - BSA and Azide free	HA750800	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell
Relative expression	PHGDH Recombinant Mouse mAb [A10H9-R] - BSA and Azide free	HA610207	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	PHOS / PDC Recombinant Rabbit mAb [PSH09-54] - BSA and Azide free	HA751301	Human, Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Tissue
Relative expression	PHOS / PDC Recombinant Rabbit mAb [PSH09-53] - BSA and Azide free	HA751300	Human, Mouse, Rat	WB
Cell treatment	Phospho-Acetyl Coenzyme A Carboxylase (S79) Recombinant Rabbit mAb [PSH07-13] - BSA and Azide free	HA751134	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-AKT2 (S474) Recombinant Rabbit mAb [PSH09-46] - BSA and Azide free	HA751293	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44] - BSA and Azide free	HA750953	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05] - BSA and Azide free	HA750132	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10] - BSA and Azide free	HA751205	Human	WB,IF-Cell,FC
Cell treatment	Phospho-AMPK alpha 1 (S496) Recombinant Rabbit mAb [SD0810] - BSA and Azide free	HA750300	Human, Mouse	WB,IP,IHC-Fr
Cell treatment	Phospho-AMPK alpha 1 (T183) +AMPK alpha 2 (T172) Recombinant Rabbit mAb [PSH14-13] - BSA and Azide free	HA751508	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-ATM (S1981) Recombinant Rabbit mAb [JM93-25] - BSA and Azide free	HA750432	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cell treatment	Phospho-ATR (T1989) Recombinant Rabbit mAb [PS01-18] - BSA and Azide free	HA750524	Human	WB,IF-Cell
Cell treatment	Phospho-BCKDHA (S293) Recombinant Rabbit mAb [PSH14-02] - BSA and Azide free	HA751498	Human, Mouse, Rat, Monkey	WB
Cell treatment	Phospho-Bcl-2 (S70) Recombinant Rabbit mAb [PSH11-02] - BSA and Azide free	HA751392	Human	WB,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Beta Catenin (S552) Recombinant Rabbit mAb [PSH08-72] - BSA and Azide free	HA751237	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	Phospho-CDK1 (T161) Recombinant Rabbit mAb [PSH03-42] - BSA and Azide free	HA750876	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	Phospho-Cdk2 (Y15) Recombinant Rabbit mAb [SN72-04] - BSA and Azide free	HA750260	Human, Mouse, Rat	WB,IHC-P,IP,IF-Tissue
Cell treatment	Phospho-CDK (Y15) Recombinant Rabbit mAb [PSH09-45] - BSA and Azide free	HA751292	Human, Mouse	WB,IHC-P
Cell treatment	Phospho-Chk1 (S317) Recombinant Rabbit mAb [PSH03-84] - BSA and Azide free	HA750909	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-Chk1 (S345) Recombinant Rabbit mAb [PS01-17] - BSA and Azide free	HA750523	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-c-Jun (S63) Recombinant Rabbit mAb [SY0297] - BSA and Azide free	HA750146	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit mAb [PSH07-12] - BSA and Azide free	HA751133	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-c-Myc (S62) Recombinant Rabbit mAb [ST49-08] - BSA and Azide free	HA750186	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-Cofilin (S3) Recombinant Rabbit mAb [PSH07-18] - BSA and Azide free	HA751139	Human, Rat, Monkey	WB
Cell treatment	Phospho-Creb (S133) Recombinant Rabbit mAb [JB25-40] - BSA and Azide free	HA750465	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,IF-Tissue
Cell treatment	Phospho-Cytosolic Phospholipase A2 (S505) Recombinant Rabbit mAb [PSH07-16] - BSA and Azide free	HA751137	Human, Mouse	WB
Cell treatment	Phospho-Cytosolic Phospholipase A2 (S505) Recombinant Rabbit mAb [PSH07-56] - BSA and Azide free	HA751176	Human, Mouse	WB,Dot Blot
Cell treatment	Phospho-DRP1 (S616) Recombinant Rabbit mAb [PSH06-77] - BSA and Azide free	HA751109	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC
Cell treatment	Phospho-EEF2 (T56) Recombinant Rabbit mAb [PSH06-36] - BSA and Azide free	HA751075	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Phospho-EGFR (T669) Recombinant Rabbit mAb [PSH10-03] - BSA and Azide free	HA751333	Human, Mouse, Rat	WB,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-EGFR (Y1045) Recombinant Rabbit mAb [PSH03-28] - BSA and Azide free	HA750868	Human	WB,IF-Cell,FC
Cell treatment	Phospho-EGFR (Y1068) Recombinant Rabbit mAb [SD2055] - BSA and Azide free	HA750290	Human	WB,IF-Cell,IF-Tissue,IP,FC
	Phospho-EGFR (Y1148) Rabbit pAb	HA500580	Human	WB,IF-Cell
Cell treatment, Relative expression	Phospho-EIF2S1 (S51) Recombinant Rabbit mAb [SZ01-06] - BSA and Azide free	HA750059	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cell treatment	Phospho-eIF4EBP1 (S65) Recombinant Rabbit mAb [PSH09-44] - BSA and Azide free	HA751291	Human	WB,IF-Cell,IP
Cell treatment	Phospho-eIF4EBP1 (T70) Recombinant Rabbit mAb [PSH09-91] - BSA and Azide free	HA751360	Human, Mouse	WB,IHC-P,IP
Cell treatment	Phospho-eIF4E (S209) Recombinant Rabbit mAb [SU0396] - BSA and Azide free	HA750160	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
Cell treatment	Phospho-eIF4G (S1108) Recombinant Rabbit mAb [PS01-08] - BSA and Azide free	HA750588	Human	WB,IHC-P
Cell treatment	Phospho-ELK1 (S383) Recombinant Rabbit mAb [PSH14-65] - BSA and Azide free	HA751526	Human, Mouse	WB,IHC-P,IP
Cell treatment	Phospho-ELK1 (S383) Recombinant Rabbit mAb [PSH14-65]	HA723663	Human, Mouse	WB,IHC-P,IP
Cell treatment	Phospho-eNOS (S1177) Recombinant Rabbit mAb [PS01-10] - BSA and Azide free	HA750522	Human, Dog, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4] - BSA and Azide free	HA750061	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01] - BSA and Azide free	HA750199	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-Ezrin (T567)/Radixin (T564)/Moesin (T558) Recombinant Rabbit mAb [PSH05-73] - BSA and Azide free	HA751023	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-FAK (Y397) Recombinant Rabbit mAb [SC54-07] - BSA and Azide free	HA750212	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-FOXO1A (T24)+ FOXO3A (T32)+ FOXO4 (T28)+ FOXO6 (T26) Recombinant Rabbit mAb [PSH06-00] - BSA and Azide free	HA751040	Human, Mouse	WB,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-FOXO3a (S253) Recombinant Rabbit mAb [ST49-01] - BSA and Azide free	HA750178	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-GSK3 (alpha + beta) (Y216 + Y279) Recombinant Rabbit mAb [SY26-05] - BSA and Azide free	HA750127	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,IF-Tissue,FC
Cell treatment	Phospho-GSK3 alpha (S21) Recombinant Rabbit mAb [PSH09-11] - BSA and Azide free	HA751265	Human, Mouse	WB,IHC-P,Dot Blot
Cell treatment	Phospho-GSK3 beta (S9) Recombinant Rabbit mAb [SY02-71] - BSA and Azide free	HA750128	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-HER2 / ErbB2 (Y1248) Recombinant Rabbit mAb [PSH04-03] - BSA and Azide free	HA750926	Human	WB,IF-Cell
Cell treatment	Phospho-Histone H2A.X (S139) Recombinant Rabbit mAb [SR33-09] - BSA and Azide free	HA750034	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IF-Cell,IF-Tissue
Cell treatment	Phospho-Histone H3 (S10) Recombinant Rabbit mAb [SA31-01] - BSA and Azide free	HA750020	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,ChIP
Cell treatment	Phospho-Hormone sensitive lipase (S660) Recombinant Rabbit mAb [PSH10-28] - BSA and Azide free	HA751345	Mouse, Rat, Human	WB
Cell treatment	Phospho-IKB alpha (S32) Recombinant Rabbit mAb [ST53-05] - BSA and Azide free	HA750196	Human, Mouse	WB,IF-Cell,IP
Cell treatment	Phospho-IKB alpha (S32/S36) Recombinant Rabbit mAb [PSH06-76] - BSA and Azide free	HA751119	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Phospho-IRE1 (S724) Recombinant Rabbit mAb [PSH03-35] - BSA and Azide free	HA750869	Human, Mouse	WB,IF-Cell,IHC-P,FC
Cell treatment	Phospho-IRF3 (S386) Recombinant Rabbit mAb [SU03-28] - BSA and Azide free	HA750138	Human	WB,IF-Cell,IF-Tissue
Cell treatment	Phospho-IRF3 (S396) Recombinant Rabbit mAb [PSH06-75] - BSA and Azide free	HA751121	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-JAK2 (Y1007 + Y1008) Recombinant Rabbit mAb [SY24-03] - BSA and Azide free	HA750117	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-JNK1/2/3 (T183 + T183 + T221) Recombinant Rabbit mAb [ST500] - BSA and Azide free	HA750173	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
Cell treatment	Phospho-Lamin A + Lamin C (S22) Recombinant Rabbit mAb [PSH11-47] - BSA and Azide free	HA751409	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Lamin A + Lamin C (S22) Recombinant Rabbit mAb [PSH11-48] - BSA and Azide free	HA751410	Human, Mouse, Rat	WB
Cell treatment, Relative expression	Phospho-LYN (Y397)/LCK (Y394)/HCK (Y411)/BLK (Y389) Recombinant Rabbit mAb [PSH09-43] - BSA and Azide free	HA751290	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-MEK3 (S189)+MEK6 (S207) Recombinant Rabbit mAb [PSH11-17] - BSA and Azide free	HA751398	Human, Mouse	WB
Cell treatment	Phospho-Met (Y1234/Y1235) Recombinant Rabbit mAb [PSH06-90] - BSA and Azide free	HA751114	Human	WB
Cell treatment	Phospho-MLKL (S345) Recombinant Rabbit mAb [JM92-37] - BSA and Azide free	HA750433	Mouse	WB,IHC-P,IHC-Fr
	Phospho-MYL9 (T18 + S19) Recombinant Rabbit mAb [PSH12-71] - BSA and Azide free	HA751447	Mouse, Rat	WB
Cell treatment	Phospho-MYPT1 (T696) Recombinant Rabbit mAb [PSH08-96] - BSA and Azide free	HA751257	Human, Mouse, Rat	WB
Cell treatment	Phospho-MYPT1 (T853) Recombinant Rabbit mAb [PSH08-71] - BSA and Azide free	HA751236	Human, Mouse, Rat	WB
Cell treatment	Phospho-NDRG1 (T346) Recombinant Rabbit mAb [PSH07-17] - BSA and Azide free	HA751138	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58] - BSA and Azide free	HA751361	Human, Mouse, Rat	WB,IF-Cell
Cell treatment, Relative expression	Phospho-Nrf2 (S40) Recombinant Rabbit mAb [SU0334] - BSA and Azide free	HA750142	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Cell treatment	Phospho-p38 alpha (T180) Recombinant Rabbit mAb [PSH06-25] - BSA and Azide free	HA751062	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180 + Y182) Recombinant Rabbit mAb [PSH06-35] - BSA and Azide free	HA751074	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-p53 (S15) Recombinant Rabbit mAb [PSH01-98] - BSA and Azide free	HA750772	Human	WB,IF-Cell
Cell treatment	Phospho-p53 (S20) Recombinant Rabbit mAb [PSH06-86] - BSA and Azide free	HA751110	Human	WB
Cell treatment	Phospho-p53 (S392) Recombinant Rabbit mAb [SI17-04] - BSA and Azide free	HA750097	Human, Mouse, Rat	WB,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-p70 S6 Kinase (T389) Recombinant Rabbit mAb [PSH02-23] - BSA and Azide free	HA750791	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-PDGFR beta (Y751) Recombinant Rabbit mAb [PSH13-61] - BSA and Azide free	HA751480	Mouse, Rat, Human	WB
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-39] - BSA and Azide free	HA750717	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-38] - BSA and Azide free	HA750716	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PKR (T446) Recombinant Rabbit mAb [SY230] - BSA and Azide free	HA750113	Human	WB,IHC-P,IP,IF-Cell
Cell treatment	Phospho-PLC-gamma-1 (S1248) Recombinant Rabbit mAb [PSH13-62] - BSA and Azide free	HA751481	Human	WB,IP
Cell treatment	Phospho-PLC-gamma-1 (Y783) Recombinant Rabbit mAb [PSH06-29] - BSA and Azide free	HA751066	Human, Mouse, Rat	WB
Cell treatment	Phospho-PPP1CA (T320) Recombinant Rabbit mAb [PSH10-27] - BSA and Azide free	HA751344	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-Raf1 (S259) Recombinant Rabbit mAb [SD85-07] - BSA and Azide free	HA750304	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-Rb (S608) Recombinant Rabbit mAb [PSH11-04] - BSA and Azide free	HA751394	Human, Monkey	WB
Cell treatment	Phospho-Rb (S795) Recombinant Rabbit mAb [PSH06-89] - BSA and Azide free	HA751113	Human	WB,IF-Cell
Cell treatment	Phospho-Rb (S807) Recombinant Rabbit mAb [SR3-82] - BSA and Azide free	HA750045	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-Rb (S807/S811) Recombinant Rabbit mAb [PSH11-01] - BSA and Azide free	HA751391	Human, Mouse	WB,IHC-P,IF-Cell,FC,IP
Cell treatment	Phospho-RICTOR (T1135) Recombinant Rabbit mAb [PSH09-92] - BSA and Azide free	HA751323	Human, Mouse	WB
	Phospho-RPS6 (S235 + S236) Recombinant Rabbit mAb [PS01-04] - BSA and Azide free	HA750587	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Phospho-RPS6 (S240 + S244) Recombinant Rabbit mAb [PS00-65] - BSA and Azide free	HA750684	Human, Mouse, Rat	WB,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Phosphoserine aminotransferase Recombinant Rabbit mAb [PSH08-99] - BSA and Azide free	HA751260	Human, Mouse, Rat	WB,IHC-P,FC,IP
Cell treatment	Phospho-Smad2 (S250) Recombinant Rabbit mAb [SD207-1] - BSA and Azide free	HA750291	Human, Mouse, Rat	WB,FC,IP,IHC-P
Cell treatment	Phospho-SMAD2 (S465+467) Recombinant Rabbit mAb [PSH05-76] - BSA and Azide free	HA751026	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-Smad3 (S423 + S425) Recombinant Rabbit mAb [ST0493] - BSA and Azide free	HA750172	Human, Mouse	WB,IF-Cell,IF- Tissue,IHC-P
Cell treatment	Phospho-Smad5 (S463 + S465) Recombinant Rabbit mAb [SY09-03] - BSA and Azide free	HA750089	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-SQSTM1 / p62 (S349) Recombinant Rabbit mAb [PSH08-97] - BSA and Azide free	HA751258	Human, Mouse	WB,IF-Cell,IHC-P
Cell treatment	Phospho-Src (S17) Recombinant Rabbit mAb [PSH11-03] - BSA and Azide free	HA751393	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC- P,FC,IP
Cell treatment	Phospho-Src (Y419) Recombinant Rabbit mAb [ST0800] - BSA and Azide free	HA750164	Human, Mouse, Rat	WB,FC,IF-Cell
Cell treatment, Relative expression	Phospho-STAT1 (S727) Recombinant Rabbit mAb [SN67-04] - BSA and Azide free	HA750247	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-STAT1 (Y701) Recombinant Rabbit mAb [PSH04-02] - BSA and Azide free	HA750925	Human	WB,IF-Cell
Cell treatment, KD	Phospho-STAT3 (S727) Recombinant Rabbit mAb [SY24-09] - BSA and Azide free	HA750121	Human, Rat, Mouse	WB,IF-Cell,IF- Tissue,IHC-P,IP
Cell treatment	Phospho-STAT3 (Y705) Recombinant Rabbit mAb [SZ43-01] - BSA and Azide free	HA750067	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-STAT5 (Y694) Recombinant Rabbit mAb [SC05-31] - BSA and Azide free	HA750220	Human, Rat, Mouse	WB,IHC-P,IF-Cell
Cell treatment	Phospho-STAT6 (Y641) Recombinant Rabbit mAb [PSH06-98] - BSA and Azide free	HA751125	Human, Mouse	WB
Cell treatment	Phospho-STING (S366) Recombinant Rabbit mAb [PSH09-72] - BSA and Azide free	HA751314	Human	WB,Dot Blot,IP
Cell treatment	Phospho-Syk (Y352) Recombinant Rabbit mAb [PSH11-66] - BSA and Azide free	HA751422	Human	WB,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Syk (Y525 + Y526) Recombinant Rabbit mAb [PSH03-70] - BSA and Azide free	HA750904	Human	WB,IF-Cell
Cell treatment	Phospho-TAK1 (S412) Recombinant Rabbit mAb [PSH14-64] - BSA and Azide free	HA751525	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-TAK1 (S412) Recombinant Rabbit mAb [PSH14-64]	HA723662	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-Tau (S396) Recombinant Rabbit mAb [SN62-09] - BSA and Azide free	HA750268	Mouse, Rat, Human, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
	Phospho-Tau (Ser214/T217) Signaling Antibody Sampler Kit	HAK21128	Human	WB
Cell treatment	Phospho-Tau (T181) Recombinant Rabbit mAb [PSH05-42] - BSA and Azide free	HA751001	Human, Mouse, Cynomolgus monkey	WB,IHC-P
Cell treatment	Phospho-Tau (T217) Recombinant Rabbit mAb [PSH09-33] - BSA and Azide free	HA751287	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr
	Phospho-Tau (T231) Recombinant Rabbit mAb [PSH01-05] - BSA and Azide free	HA750679	Human, Mouse, Rat	WB,IHC-P
	Phospho-Tau (T231) Recombinant Rabbit mAb [PSH01-04] - BSA and Azide free	HA750680	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Cell treatment	Phospho-Tau (T231) Recombinant Rabbit mAb [SC58-08] - BSA and Azide free	HA750210	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IP,IF-Tissue,IHC-Fr
Cell treatment	Phosphothreonine Recombinant Rabbit mAb [PSH14-34] - BSA and Azide free	HA751516	Human, Mouse, Rat, Species independent	WB,IHC-P,IP,ELISA
Relative expression	Phospho-TrkB (Y817) Recombinant Rabbit mAb [SC0556] - BSA and Azide free	HA750213	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-YAP1 (S127) Recombinant Rabbit mAb [SN0718] - BSA and Azide free	HA750269	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	PI 3 Kinase p85 alpha Recombinant Rabbit mAb [SU04-07] - BSA and Azide free	HA750162	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	PI 3 Kinase p85 alpha Recombinant Mouse mAb [A3-D0-R] - BSA and Azide free	HA610084	Human	WB,IF-Cell
	PICALM Recombinant Rabbit mAb [PSH01-63] - BSA and Azide free	HA750739	Human, Mouse, Rat	WB,IHC-P
	PICALM Recombinant Rabbit mAb [PSH01-64] - BSA and Azide free	HA750740	Human, Mouse, Rat	WB
	PIM1 Recombinant Rabbit mAb [ST0513] - BSA and Azide free	HA750183	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Pin1 Recombinant Rabbit mAb [PSH10-81] - BSA and Azide free	HA751373	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	PINK1 Recombinant Rabbit mAb [PSH08-69] - BSA and Azide free	HA751234	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
KD	PIP5K1C Recombinant Rabbit mAb [PSH01-43] - BSA and Azide free	HA750720	Human	WB
	PKA C-alpha Recombinant Rabbit mAb [PSH08-57] - BSA and Azide free	HA751232	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC,IP
	PKC delta Recombinant Rabbit mAb [JJ093-02] - BSA and Azide free	HA750325	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Relative expression	PKC zeta Recombinant Rabbit mAb [PSH06-74] - BSA and Azide free	HA751122	Human	WB
Relative expression	PKM1 Recombinant Rabbit mAb [PSH04-37] - BSA and Azide free	HA750946	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	PKM2 Recombinant Rabbit mAb [PSH11-74] - BSA and Azide free	HA751429	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP
	Plexin B2 Recombinant Rabbit mAb [PSH03-48] - BSA and Azide free	HA750882	Human	WB,IF-Cell,IHC-P
	PLK1 Recombinant Rabbit mAb [PSH04-35] - BSA and Azide free	HA750944	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P
KD	PML Protein Recombinant Rabbit mAb [PSH02-89] - BSA and Azide free	HA750843	Human	WB,IHC-P,IF-Cell,FC
KD	PML Protein Recombinant Rabbit mAb [PSH02-68] - BSA and Azide free	HA750825	Human, Rat	WB,IHC-P,IF-Cell
	PMP22 Recombinant Rabbit mAb [JM52-30] - BSA and Azide free	HA750425	Human, Mouse, Rat	WB
KD	PMP70 Recombinant Rabbit mAb [PSH03-27] - BSA and Azide free	HA750867	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	PNLDC1 Recombinant Rabbit mAb [PSH09-60] - BSA and Azide free	HA751305	Human, Mouse, Rat	WB
	POLDIP2 Recombinant Rabbit mAb [PSH01-20] - BSA and Azide free	HA750696	Human, Rat	WB
KD	PPP2R1B Recombinant Rabbit mAb [PSH02-92] - BSA and Azide free	HA750847	Human, Mouse, Rat	WB,IHC-P
	PPP2R5D Recombinant Rabbit mAb [PSH01-22] - BSA and Azide free	HA750698	Human, Mouse, Rat	WB,IHC,FC
	PPP2R5D Recombinant Rabbit mAb [PSH01-21] - BSA and Azide free	HA750697	Human	WB,IF-Cell,FC
	PPP3CB Recombinant Rabbit mAb [PSH0-21] - BSA and Azide free	HA750595	Human, Mouse, Rat	WB,IHC-P,ICC,FC
	PRA1 Recombinant Rabbit mAb [PSH01-96] - BSA and Azide free	HA750770	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell
	PRAK / MK5 Recombinant Rabbit mAb [PSH14-06] - BSA and Azide free	HA751501	Human, Mouse, Rat, Monkey	WB,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Prealbumin Recombinant Rabbit mAb [JM11-43] - BSA and Azide free	HA750392	Human	WB,IP,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Presenilin 1 / PS-1 Recombinant Rabbit mAb [PSH08-49] - BSA and Azide free	HA751224	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr,IF-Tissue
	Presenilin 1 / PS-1 Recombinant Rabbit mAb [PSH11-19] - BSA and Azide free	HA751400	Human, Mouse, Rat	WB,IP
	Presenilin 2 / AD5 Recombinant Rabbit mAb [PSH11-18] - BSA and Azide free	HA751399	Human, Mouse, Rat	WB,IF-Cell,FC,IP
	Presenilin 2 / AD5 Recombinant Rabbit mAb [PSH05-18] - BSA and Azide free	HA750989	Human, Mouse, Rat	WB,IP
	PRMT5 Recombinant Rabbit mAb [ST51-06] - BSA and Azide free	HA750174	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cell treatment, Relative expression, KO	Pro Caspase-3 Recombinant Rabbit mAb [SR03-01] - BSA and Azide free	HA750048	Human, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Prohibitin Recombinant Rabbit mAb [SR46-02] - BSA and Azide free	HA750043	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,IP,FC
Relative expression, Cell treatment	Prosurfactant Protein C Recombinant Rabbit mAb [PSH15-28]	HA723734	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue
	PSD95 Recombinant Rabbit mAb [SR38-09] - BSA and Azide free	HA750035	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Tissue,IHC-P,IP,IHC-Fr
Relative expression	PSD95 Recombinant Antibody [SR38-09] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610256	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,WB
Relative expression	PSIP1 Recombinant Rabbit mAb [PSH02-55] - BSA and Azide free	HA750818	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	PSMA7 / HSPC Recombinant Rabbit mAb [PSH14-09] - BSA and Azide free	HA751504	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
Cell treatment	PSMB10 Recombinant Rabbit mAb [PSH01-58] - BSA and Azide free	HA750734	Human, Mouse	WB,IHC-P
	PSMB2 Recombinant Rabbit mAb [PSH0-28] - BSA and Azide free	HA750602	Human	WB
Cell treatment	PSMB5 / MB1 Recombinant Rabbit mAb [PSH15-20]	HA723724	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	PSMC5 Recombinant Rabbit mAb [PSH0-22] - BSA and Azide free	HA750596	Human, Mouse, Rat	WB,IHC-P
	PSMC6 Recombinant Rabbit mAb [PSH0-12] - BSA and Azide free	HA750582	Human, Mouse, Rat	WB,IHC-P,FC
	PSMD13 Recombinant Rabbit mAb [PSH0-52] - BSA and Azide free	HA750629	Human, Mouse, Rat	WB,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PSMD7 Recombinant Rabbit mAb [PSH01-83] - BSA and Azide free	HA750757	Human, Mouse	IF-Cell
KD	PSME2 Recombinant Rabbit mAb [PSH0-10] - BSA and Azide free	HA750579	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	PTBP1 Recombinant Mouse mAb [5C1-R] - BSA and Azide free	HA610075	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression, KD	PTEN Recombinant Rabbit mAb [SJ19-03] - BSA and Azide free	HA750105	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
KD	PTGR1 Recombinant Rabbit mAb [PSH02-12] - BSA and Azide free	HA750782	Human, Rat, Mouse	WB,IHC-P,FC
	PUMA Recombinant Rabbit mAb [SR42-09] - BSA and Azide free	HA750038	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Pyrin Recombinant Rabbit mAb [PSH06-12] - BSA and Azide free	HA751050	Human	WB,FC
	QKI Recombinant Mouse mAb [PSH10-21] - BSA and Azide free	HA610222	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr,IP
	QKI-5 Recombinant Antibody [PSH10-22] - Rabbit IgG (Chimeric) - BSA and Azide free	HA751496	Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	QKI-5 Recombinant Antibody [PSH10-22] - Guinea pig IgG2 (Chimeric)	HA601501	Mouse, Rat	IHC-Fr,WB
	QKI-5 Recombinant Mouse mAb [PSH10-22] - BSA and Azide free	HA610223	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IP
	QKI-6 Recombinant Antibody [PSH10-23] - Rabbit IgG (Chimeric) - BSA and Azide free	HA751497	Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	QKI-6 Recombinant Antibody [PSH10-23] - Guinea pig IgG2 (Chimeric)	HA601502	Mouse, Rat	IHC-Fr,WB
	QKI-6 Recombinant Mouse mAb [PSH10-23] - BSA and Azide free	HA610224	Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IHC-Fr
	QKI-7 Recombinant Mouse mAb [PSH10-92] - BSA and Azide free	HA610230	Mouse, Rat, Cynomolgus monkey, Pig	WB
KD	QPRT Recombinant Rabbit mAb [PSH01-31] - BSA and Azide free	HA750709	Human	WB,IHC-P
KD	RAB10 Recombinant Rabbit mAb [PSH03-24] - BSA and Azide free	HA750864	Human, Mouse, Rat	WB
Relative expression, KD	Rab1A Recombinant Rabbit mAb [PSH03-62] - BSA and Azide free	HA750896	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Relative expression	RAB27B Recombinant Rabbit mAb [PSH14-68]	HA723666	Human, Mouse, Rat	WB,IHC-P
Relative expression	RAB27B Recombinant Rabbit mAb [PSH14-68] - BSA and Azide free	HA751529	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression, KD	Rab3A Recombinant Rabbit mAb [PSH01-08] - BSA and Azide free	HA750683	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
KD	Rab5 Recombinant Rabbit mAb [ST47-04] - BSA and Azide free	HA750167	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KD	RAB7 Recombinant Rabbit mAb [SN202-03] - BSA and Azide free	HA750279	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	RAC1 Recombinant Rabbit mAb [PSH05-09] - BSA and Azide free	HA750982	Human, Rat	WB,IF-Cell
	Rad51 Recombinant Rabbit mAb [JM54-26] - BSA and Azide free	HA750450	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	RAD54 Recombinant Mouse mAb [A10H5-R] - BSA and Azide free	HA610212	Human, Mouse	WB
Relative expression	RAGE Recombinant Rabbit mAb [JF0975] - BSA and Azide free	HA750341	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Rat CRP/C-Reactive Protein enzyme-linked immunoassay kit(one step)	ERY0009	Rat	ELISA
	Rb Recombinant Rabbit mAb [SY63-03] - BSA and Azide free	HA750135	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP
KD	RbBP5 Recombinant Rabbit mAb [PSH0-15] - BSA and Azide free	HA750585	Human, Mouse, Rat	WB
	RbBP5 Recombinant Rabbit mAb [PSH0-16] - BSA and Azide free	HA750586	Human, Mouse	IHC-P
KD	RBMX2 Recombinant Rabbit mAb [PSH09-47] - BSA and Azide free	HA751294	Human	WB
	RBX1 Recombinant Rabbit mAb [PSH07-01] - BSA and Azide free	HA751128	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,IF-Tissue
	RecA Recombinant Rabbit mAb [PSH09-96] - BSA and Azide free	HA751327	Escherichia coli	WB,IF-Cell
Cell treatment	Recombinant Mouse IgG1 Recombinant Mouse mAb [PSH04-91] - BSA and Azide free	HA610204		WB,FC,IP,ChIP
Cell treatment	Recombinant Rabbit IgG Recombinant Rabbit mAb [PSH04-42] - BSA and Azide free	HA750951		WB,IHC-P,IF-Cell,FC,IP,ChIP
	REIIBP/NSD2 Recombinant Rabbit mAb [PSH02-39] - BSA and Azide free	HA750804	Human, Mouse, Rat	WB,IHC-P
	Ret Recombinant Rabbit mAb [SN20-28] - BSA and Azide free	HA750280	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Retinoic Acid Receptor alpha/beta/gamma Recombinant Rabbit mAb [PSH08-74] - BSA and Azide free	HA751239	Human	WB
	Retinoic Acid Receptor gamma Recombinant Rabbit mAb [PSH07-19] - BSA and Azide free	HA751140	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Retinoid X Receptor beta / RXRB Recombinant Rabbit mAb [PSH13-53] - BSA and Azide free	HA751472	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	RFP Recombinant Rabbit mAb [JA10-52] - BSA and Azide free	HA750408	Species independent	WB
	RhoA Recombinant Rabbit mAb [SN0612] - BSA and Azide free	HA750242	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	RHOG Recombinant Rabbit mAb [PSH02-70] - BSA and Azide free	HA750827	Human, Monkey	WB,FC
	RICTOR Recombinant Rabbit mAb [1A5] - BSA and Azide free	HA750484	Human, Mouse	WB,IHC-P,IF-Cell,FC
	RICTOR Recombinant Rabbit mAb [1D5] - BSA and Azide free	HA750483	Human	WB,IF-Cell,FC
Relative expression	RIP3 Recombinant Rabbit mAb [PSH04-78] - BSA and Azide free	HA750972	Human	WB,IHC-P
Relative expression	RIP3 Recombinant Rabbit mAb [PSH04-79] - BSA and Azide free	HA750973	Human	WB,IHC-P
KD	RNF114 Recombinant Rabbit mAb [PSH01-99] - BSA and Azide free	HA750773	Human	WB
	RNF146 Recombinant Rabbit mAb [PSH0-29] - BSA and Azide free	HA750603	Human, Mouse	WB,IHC-P,IF-Cell,FC
	RNF168 Recombinant Rabbit mAb [PSH07-38] - BSA and Azide free	HA751159	Human, Mouse, Rat	IHC-P,IF-Cell
	RNF168 Recombinant Rabbit mAb [PSH07-37] - BSA and Azide free	HA751158	Human	WB
	ROBO1 Recombinant Mouse mAb [10E2-R] - BSA and Azide free	HA610007	Human, Mouse, Rat	WB,IHC-P
	ROCK1 Recombinant Rabbit mAb [ST05-19] - BSA and Azide free	HA750181	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
Relative expression	ROR1 Recombinant Mouse mAb [PSH13-99] - BSA and Azide free	HA610281	Human	WB,FC
Relative expression	ROR alpha / RORA Recombinant Rabbit mAb [PSH11-00] - BSA and Azide free	HA751390	Mouse, Rat	IHC-P,WB
	ROR gamma T Recombinant Rabbit mAb [PSH04-38] - BSA and Azide free	HA750947	Human, Rat	WB,IHC-P
	ROS1 Recombinant Rabbit mAb [PD01-27] - BSA and Azide free	HA750637	Human	IHC-P
	RPS6 Recombinant Mouse mAb [A6B8-R] - BSA and Azide free	HA610196	Human, Mouse, Rat	WB
	RPS6 Recombinant Mouse mAb [A6B9-R] - BSA and Azide free	HA610092	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	RSK4 Recombinant Rabbit mAb [PSH02-04] - BSA and Azide free	HA750774	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	RTN3 Recombinant Rabbit mAb [PSH0-40] - BSA and Azide free	HA750618	Human, Rat	WB,IF-Cell,FC
	Rubicon / Baron Recombinant Rabbit mAb [PSH08-25] - BSA and Azide free	HA751216	Human, Mouse, Rat, Monkey	WB
	RUNX1+RUNX2+RUNX3 Recombinant Rabbit mAb [SD0803] - BSA and Azide free	HA750296	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Relative expression, KD	RUNX2 Recombinant Rabbit mAb [SD208-0] - BSA and Azide free	HA750295	Human, Mouse, Rat	IF-Cell,IF-Tissue,IHC-P,WB,FC,IP
	RUVB2 Recombinant Mouse mAb [A1E12-R] - BSA and Azide free	HA610114	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P
Cell treatment	S100A4 Recombinant Rabbit mAb [SD200-08] - BSA and Azide free	HA750282	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Relative expression	S100A8 Recombinant Rabbit mAb [PSH09-70] - BSA and Azide free	HA751312	Human	WB,IHC-P,IF-Cell,FC,IF-Tissue,IP
	S100 alpha Recombinant Rabbit mAb [JM51-31] - BSA and Azide free	HA750389	Human, Mouse, Rat	WB,IP,IHC-P,FC,IF-Cell
KD	S100 alpha 6 Recombinant Rabbit mAb [JF0976] - BSA and Azide free	HA750342	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	S100B Recombinant Rabbit mAb [PSH15-25] - BSA and Azide free (Detector)	HA723729	Human, Mouse, Rat	ELISA(Det)
	Biotin Conjugated S100B Recombinant Rabbit mAb [PSH15-25] - Detector	HA723730B	Human, Mouse, Rat	ELISA(Det),ELISA
	S100B Recombinant Rabbit mAb [SC57-02] - BSA and Azide free (Capture)	HA723728	Human, Mouse, Rat	ELISA(Cap)
	S100 beta Recombinant Rabbit mAb [SC57-02] - BSA and Azide free	HA750209	Human, Mouse, Rat, Zebrafish, Goat, Cynomolgus monkey, Green monkey, Pig, Rabbit	WB,IF-Tissue,IHC-P,IHC-Fr,IF-Cell,FC,IP
	SAGE1 Recombinant Rabbit mAb [PSH05-89] - BSA and Azide free	HA751036	Human	WB,IHC-P,IF-Cell,IP
KD	SAMD9 Recombinant Rabbit mAb [PSH02-75] - BSA and Azide free	HA750830	Human	WB,IF-Cell,IHC-P,FC
	SATB2 Recombinant Rabbit mAb [JM52-44] - BSA and Azide free	HA750449	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC
	SATB2 Recombinant Antibody [JM52-44] - Guinea pig IgG2 (Chimeric)	HA601508	Human, Mouse, Rat	IHC-Fr,WB
Relative expression	SAV1 Recombinant Rabbit mAb [PSH04-45] - BSA and Azide free	HA750954	Human, Mouse, Rat	WB,IHC-P,FC
	SCD1 Recombinant Mouse mAb [A9A11-R] - BSA and Azide free	HA610063	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	SCGB1A1 Recombinant Rabbit mAb [JU34-03] - BSA and Azide free	HA750458	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
	SCO2 Recombinant Rabbit mAb [PSH09-97] - BSA and Azide free	HA751328	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP
	SCP3 Recombinant Rabbit mAb [PSH02-24] - BSA and Azide free	HA750792	Human	WB,IHC-P
Cell treatment	SCXA Recombinant Rabbit mAb [PSH02-25] - BSA and Azide free	HA750793	Human	WB,IF-Cell
	SEC24A Recombinant Rabbit mAb [PSH04-36] - BSA and Azide free	HA750945	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	SEC61B Recombinant Rabbit mAb [PSH09-48] - BSA and Azide free	HA751295	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,IF-Tissue
	SEC62 Recombinant Rabbit mAb [PSH0-25] - BSA and Azide free	HA750599	Human, Mouse, Rat	WB
	SEC62 Recombinant Rabbit mAb [PSH0-91] - BSA and Azide free	HA750671	Human, Mouse, Rat	WB
	SESN2 / Sestrin-2 Recombinant Rabbit mAb [PSH05-78] - BSA and Azide free	HA751038	Human, Mouse, Rat	WB
KD	SET7 Recombinant Rabbit mAb [PSH02-90] - BSA and Azide free	HA750845	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
Relative expression	SHANK2 Recombinant Rabbit mAb [PSH10-61] - BSA and Azide free	HA751364	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IF-Tissue
	SHARP2 Recombinant Rabbit mAb [PSH01-90] - BSA and Azide free	HA750764	Human	WB,IHC-P,IF-Cell
	SHARP2 Recombinant Rabbit mAb [PSH01-91] - BSA and Azide free	HA750765	Human	WB
	SHP1 Recombinant Rabbit mAb [SR41-02] - BSA and Azide free	HA750033	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	SIAH1 Recombinant Rabbit mAb [PSH01-80] - BSA and Azide free	HA750754	Human, Mouse, Rat	WB,IF-Cell
KD	SIAH1 Recombinant Rabbit mAb [PSH0-46] - BSA and Azide free	HA750626	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
KD	SIRT1 Recombinant Rabbit mAb [SZ04-01] - BSA and Azide free	HA750064	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,ChIP
	SIRT2 Recombinant Antibody [SN70-04] - Mouse IgG1 (Chimeric)	HA601505	Mouse, Rat	IHC-Fr,IHC-P,WB
	SIRT2 Recombinant Antibody [SN70-04] - Rat IgG1 (Chimeric)	HA601506	Mouse, Rat	IHC-Fr,IHC-P,WB
KD	SIRT3 Recombinant Rabbit mAb [PSH05-22] - BSA and Azide free	HA750993	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue,IP,ChIP
	SIRT4 Recombinant Rabbit mAb [PSH15-46]	HA723752	Human, Mouse, Rat, Monkey	WB,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KO	SIRT6 Recombinant Rabbit mAb [PSH10-99] - BSA and Azide free	HA751389	Human, Mouse, Rat, Monkey	WB,IHC-P
	SIRT6 Recombinant Rabbit mAb [PSH10-98] - BSA and Azide free	HA751388	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P
	SIRT7 Recombinant Rabbit mAb [PSH09-17] - BSA and Azide free	HA751271	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,ChIP
Cell treatment	Sirtuin Antibody Sampler Kit	HAK21076	Human	WB
	SKP2 Recombinant Rabbit mAb [PSH06-22] - BSA and Azide free	HA751060	Human, Monkey	WB
Relative expression	SLAM / CD150 Recombinant Rabbit mAb [PSH12-08] - BSA and Azide free	HA751437	Human	WB,IP
	Slc12a5 / KCC2 Recombinant Mouse mAb [PSH13-20] - BSA and Azide free	HA610271	Mouse, Rat	IHC-Fr,WB
Relative expression	SLC16A3 / MCT 4 Recombinant Rabbit mAb [PSH06-91] - BSA and Azide free	HA751115	Human, Monkey	WB,IF-Cell,IHC-P,FC
	SLC1A5 Recombinant Rabbit mAb [PSH07-74] - BSA and Azide free	HA751186	Human, Mouse	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	SLC25A38 Recombinant Rabbit mAb [PSH02-79] - BSA and Azide free	HA750834	Human, Mouse	WB,IF-Cell,IHC-P
	SLC32A1 / VGAT Recombinant Antibody [PSH10-47] - Chicken IgY (Chimeric) - BSA and Azide free	HA610305	Mouse	IHC-Fr
	SLC32A1 / VGAT Recombinant Mouse mAb [PSH10-47] - BSA and Azide free	HA610228	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
	SLC32A1 / VGAT Recombinant Antibody [PSH10-47] - Rabbit IgG (Chimeric) - BSA and Azide free	HA751397	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P
	SLC32A1 / VGAT Recombinant Mouse mAb [PSH08-35] - BSA and Azide free	HA610218	Mouse, Rat	IHC-P
Relative expression	SLC7A5 / LAT1 Recombinant Rabbit mAb [PSH09-28] - BSA and Azide free	HA751282	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Relative expression	Slow Skeletal Myosin Heavy chain Recombinant Rabbit mAb [PSH07-54] - BSA and Azide free	HA751092	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression	Slow Skeletal Myosin Heavy chain Recombinant Rabbit mAb [PSH06-58] - BSA and Azide free	HA751086	Human, Mouse, Rat	WB,IHC-P
Relative expression	SLP76 Recombinant Rabbit mAb [PSH09-22] - BSA and Azide free	HA751276	Human	WB,IHC-P,IF-Cell,IP
Relative expression	SLUG Recombinant Rabbit mAb [PSH07-40] - BSA and Azide free	HA751161	Human	WB,IF-Cell
KO	Smad1 Recombinant Rabbit mAb [SY0254] - BSA and Azide free	HA750124	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
KO	Smad2 Recombinant Rabbit mAb [SP06-05] - BSA and Azide free	HA750073	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KO	Smad3 Recombinant Rabbit mAb [SY25-01] - BSA and Azide free	HA750123	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KO	Smad4 Recombinant Rabbit mAb [SP05-05] - BSA and Azide free	HA750070	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	SMARCA2 Recombinant Rabbit mAb [PSH01-15] - BSA and Azide free	HA750691	Human, Monkey	WB,IHC-P,FC
	SMC2 Recombinant Rabbit mAb [PSH0-63] - BSA and Azide free	HA750640	Human, Monkey	WB,IHC-P,IF-Cell
	SMC2 Recombinant Rabbit mAb [PSH0-62] - BSA and Azide free	HA750639	Human, Monkey	WB,IHC-P,IF-Cell
	SMC5 Recombinant Rabbit mAb [PSH15-29]	HA723735	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P
	SMOC2 Recombinant Rabbit mAb [PSH11-25] - BSA and Azide free	HA751404	Human, Mouse, Rat	WB,IHC-P
	smooth muscle Myosin heavy chain 11 Recombinant Rabbit mAb [JA03-35] - BSA and Azide free	HA750417	Human, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
	SMURF1 Recombinant Rabbit mAb [PSH10-12] - BSA and Azide free	HA751341	Human	WB
	SNAP29 Recombinant Rabbit mAb [PSH01-76] - BSA and Azide free	HA750750	Human	WB
Relative expression	SOCS1 Recombinant Rabbit mAb [PSH09-68] - BSA and Azide free	HA751310	Human, Mouse, Rat, Monkey	WB
	SOD2 Recombinant Rabbit mAb [JJ089-02] - BSA and Azide free	HA750314	Human, Mouse, Rat	WB,IHC-P
KD	Sodium Potassium ATPase Recombinant Rabbit mAb [ST0533] - BSA and Azide free	HA750194	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
	Sortilin / NT3 Recombinant Rabbit mAb [PSH06-16] - BSA and Azide free	HA751054	Human, Mouse, Rat	WB,IHC-P
	SOX1 Recombinant Rabbit mAb [JJ20-40] - BSA and Azide free	HA750308	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	SOX1 Recombinant Rabbit mAb [PSH10-51] - BSA and Azide free	HA751354	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IHC-P
	SOX10 Recombinant Antibody [PDH0-03] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610253	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Tissue
	SOX2 Recombinant Antibody [PO00-28] - Chicken IgY (Chimeric)	HA601499	Human, Mouse	IHC-Fr,WB
	SOX2 Recombinant Antibody [PO00-28] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610236	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
	SOX2 Recombinant Antibody [PO00-28] - Guinea pig IgG2 (Chimeric)	HA601500	Human, Mouse	IHC-Fr,WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression, KD	SOX9 Recombinant Rabbit mAb [SN74-09] - BSA and Azide free	HA750262	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-P,IP,FC,CUT&Tag-seq
Relative expression	SOX9 Recombinant Rabbit mAb [PSH13-59] - BSA and Azide free	HA751478	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell,IP,FC,ChIP
	SP1 Recombinant Rabbit mAb [JF0950] - BSA and Azide free	HA750330	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,ChIP,IP
	SP110 Recombinant Rabbit mAb [PSH03-66] - BSA and Azide free	HA750900	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Sp7 / Osterix Recombinant Rabbit mAb [PSH07-29] - BSA and Azide free	HA751150	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression	SPARC Recombinant Rabbit mAb [PSH06-23] - BSA and Azide free	HA751061	Human	WB
	SPINK1 Recombinant Rabbit mAb [PD00-93] - BSA and Azide free	HA750551	Human	IHC-P
	SPR Recombinant Rabbit mAb [PSH06-14] - BSA and Azide free	HA751052	Human	WB,IHC-P,IF-Cell,FC
Cell treatment, KD	SQSTM1 / p62 Recombinant Rabbit mAb [PS00-61] - BSA and Azide free	HA750510	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,FC,IF-Tissue,IHC-Fr
KO	Src Recombinant Rabbit mAb [JF0947] - BSA and Azide free	HA750331	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	SREBP1 Recombinant Rabbit mAb [PSH04-61] - BSA and Azide free	HA750968	Human	WB,IHC-P,IF-Cell
	SRF Recombinant Rabbit mAb [PSH05-06] - BSA and Azide free	HA750979	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	SRPK1 Recombinant Rabbit mAb [PSH08-26] - BSA and Azide free	HA751217	Human, Mouse	WB,IF-Cell,IHC-P,IF-Tissue,FC,IP
	SSRP1 Recombinant Rabbit mAb [PSH09-75] - BSA and Azide free	HA751317	Human, Mouse, Rat, Monkey	WB,IHC-P
	SST/Somatostatin 28 Recombinant Rabbit mAb [JE43-52] - BSA and Azide free	HA750648	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-P,IHC-Fr,IF-Tissue
	S tag Recombinant Mouse mAb [G10-C8-R] - BSA and Azide free	HA610198	Species independent	WB
Cell treatment	Stanniocalcin 2 / STC-2 Recombinant Rabbit mAb [PSH12-01] - BSA and Azide free	HA751434	Human	WB,IP
	STAT1 Recombinant Mouse mAb [G3-B11-R] - BSA and Azide free	HA610181	Human	WB,IHC-P
KD	STAT1 alpha Recombinant Rabbit mAb [SJ01-89] - BSA and Azide free	HA750102	Human, Mouse	WB,IP,FC,IHC-P,IF-Cell,IF-Tissue
KD	STAT3 Recombinant Rabbit mAb [SY34-01] - BSA and Azide free	HA750087	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	STAT3 Recombinant Rabbit mAb [SY24-08] - BSA and Azide free	HA750120	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	Stat3/Stat5 Signaling Antibody Sampler Kit	HAK21126	Human	WB
KD	STAT5a Recombinant Rabbit mAb [SC66-09] - BSA and Azide free	HA750225	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	STAT 5A+B Recombinant Rabbit mAb [JJ08-78] - BSA and Azide free	HA750310	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Relative expression	Steroid sulfatase Recombinant Rabbit mAb [PSH05-74] - BSA and Azide free	HA751024	Human	WB,IHC-P
	STIL / SIL Recombinant Rabbit mAb [PSH05-79] - BSA and Azide free	HA751039	Human	WB
	STING Recombinant Rabbit mAb [JM03-47] - BSA and Azide free	HA750437	Human, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue,IP
	STING Recombinant Rabbit mAb [PSH07-44] - BSA and Azide free	HA751165	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,IP
	Stra8 Recombinant Rabbit mAb [PSH04-55] - BSA and Azide free	HA750962	Human	WB,IF-Cell,FC
	Strep-Tag II Recombinant Rabbit mAb [PSH05-14] - BSA and Azide free	HA750987	Species independent	WB,IF-Cell,IP
	Strep-Tag II Recombinant Rabbit mAb [PSH05-13] - BSA and Azide free	HA750986	Species independent	WB,IF-Cell,IP
	SUMO-1 Recombinant Rabbit mAb [SJ20-03] - BSA and Azide free	HA750110	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,ChIP
	Survivin Recombinant Rabbit mAb [PSH11-73] - BSA and Azide free	HA751428	Human, Mouse	WB,IHC-P
	Survivin Recombinant Rabbit mAb [SP07-06] - BSA and Azide free	HA750076	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	SV2A Recombinant Rabbit mAb [PSH10-69] - BSA and Azide free	HA751371	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IP
Relative expression	SV2A Recombinant Antibody [PSH10-69] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610261	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
KD	Syk Recombinant Rabbit mAb [PSH0-76] - BSA and Azide free	HA750655	Human	WB,IF-Cell
Relative expression	Synapsin I + II Recombinant Rabbit mAb [PSH10-29] - BSA and Azide free	HA751346	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-Fr,IHC-P,IP,IF-Tissue
Relative expression	Synaptopodin Recombinant Antibody [PSH10-45] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610242	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Synaptopodin Recombinant Rabbit mAb [PSH10-45] - BSA and Azide free	HA751350	Human, Mouse, Rat	WB,IHC-Fr,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Syndecan 1 Recombinant Rabbit mAb [JM11-21] - BSA and Azide free	HA750381	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Syndecan 4 Recombinant Rabbit mAb [PSH0-72] - BSA and Azide free	HA750651	Human	WB
KD	Syntaxin 4 Recombinant Rabbit mAb [PSH0-80] - BSA and Azide free	HA750660	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	T7 tag Recombinant Mouse mAb [4-C2-R] - BSA and Azide free	HA610073	Species independent	WB
KD	TAK1 Recombinant Rabbit mAb [JM73-19] - BSA and Azide free	HA750424	Human, Mouse, Rat, Monkey	WB
Relative expression	Talin 1 Recombinant Rabbit mAb [PSH09-95] - BSA and Azide free	HA751326	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment, KD	Tankyrase 2 Recombinant Rabbit mAb [PSH0-53] - BSA and Azide free	HA750630	Human, Rat	WB,IF-Cell
Relative expression	Tau Recombinant Rabbit mAb [PSH05-50] - BSA and Azide free	HA751009	Mouse, Rat	WB,IHC-P,IF-Tissue,IP
	TBC1D1 Recombinant Rabbit mAb [PSH13-71] - BSA and Azide free	HA751484	Human, Mouse, Rat	WB,IHC-P,IP
	TBR1 Recombinant Rabbit mAb [JF10-00] - BSA and Azide free	HA750366	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression	TBR1 Recombinant Antibody [JF10-00] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610262	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
Relative expression	TBR2 Recombinant Rabbit mAb [PSH07-88] - BSA and Azide free	HA751195	Human, Mouse, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
Relative expression	TCF-4 Recombinant Rabbit mAb [PSH15-19]	HA723723	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	TCF7 Recombinant Rabbit mAb [PSH13-89] - BSA and Azide free	HA751491	Human, Mouse, Rat	WB,IHC-P,IP
	TCF7 Recombinant Rabbit mAb [PSH10-10] - BSA and Azide free	HA751340	Human	WB,IHC-P,IP
	TDP1 Recombinant Rabbit mAb [PSH0-93] - BSA and Azide free	HA750673	Human, Mouse, Rat	WB,IHC-P
	TDP43 Recombinant Rabbit mAb [PSH09-14] - BSA and Azide free	HA751268	Human, Mouse, Rat, Zebrafish, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-P,IHC-Fr,IF-Tissue,FC,IP
	TDP43 Recombinant Rabbit mAb [JM51-10] - BSA and Azide free	HA750390	Human, Mouse, Rat, Zebrafish, Cynomolgus monkey, Pig	WB,IP,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
	TEAD2 Recombinant Rabbit mAb [PSH07-32] - BSA and Azide free	HA751153	Human, Mouse	WB
	TFEB Recombinant Rabbit mAb [PSH12-05] - BSA and Azide free	HA751436	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - Detector	HA723727B	Human, Mouse, Rat	ELISA(Det),ELISA
	TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - BSA and Azide free (Detector)	HA723726	Human, Mouse, Rat	ELISA(Det)
	TGF beta 1 Recombinant Rabbit mAb [PSH15-21] - BSA and Azide free (Capture)	HA723725	Human, Mouse, Rat	ELISA(Cap)
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17] - BSA and Azide free	HA750485	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	TGF beta 3 Recombinant Rabbit mAb [PSH09-26] - BSA and Azide free	HA751280	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	TGF beta Receptor I Recombinant Rabbit mAb [PSH14-14] - BSA and Azide free	HA751509	Human, Mouse, Rat	WB,IP
	TGF beta Receptor II Recombinant Rabbit mAb [PSH01-50] - BSA and Azide free	HA750726	Human	WB
	Thrombospondin 1 Recombinant Rabbit mAb [PSH02-97] - BSA and Azide free	HA750852	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
	Thymidine Phosphorylase Recombinant Mouse mAb [A1A10-R] - BSA and Azide free	HA610106	Human, Mouse, Rat	WB,IHC-P
	Thymidylate Synthase Recombinant Rabbit mAb [JM33-20] - BSA and Azide free	HA750427	Human, Mouse, Rat	WB,IHC-P,FC,IP
	TIA1 Recombinant Rabbit mAb [JM42-11] - BSA and Azide free	HA750385	Human, Mouse	WB,IF-Tissue,IHC-P,IP
Relative expression	TIE2 Recombinant Rabbit mAb [PSH06-63] - BSA and Azide free	HA751089	Human	WB,FC
Relative expression	TIM 3 Recombinant Rabbit mAb [PSH07-75] - BSA and Azide free	HA751187	Human	WB
Relative expression	TIM 3 Recombinant Antibody [PSH13-95] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610277	Human	WB,IF-Cell,FC
Relative expression	TIMP2 Recombinant Rabbit mAb [PSH11-62] - BSA and Azide free	HA751418	Human, Monkey	WB
Relative expression	TIMP3 Recombinant Rabbit mAb [PSH04-41] - BSA and Azide free	HA750950	Human, Rat	WB,IF-Cell
	Tissue Factor Recombinant Rabbit mAb [SN20-16] - BSA and Azide free	HA750278	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Tissue type plasminogen activator Recombinant Rabbit mAb [JF0958] - BSA and Azide free	HA750338	Human, Mouse, Rat	WB,IHC-P
KO	TLR2 Recombinant Rabbit mAb [JM22-41] - BSA and Azide free	HA750447	Human, Mouse, Rat	WB,FC,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	TMC1 Recombinant Rabbit mAb [PSH11-12] - BSA and Azide free	HA751396	Human, Mouse	WB,IF-Cell
	TMED9 Recombinant Mouse mAb [A9A8-R] - BSA and Azide free	HA610143	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	TMED9 Recombinant Mouse mAb [A9A7-R] - BSA and Azide free	HA610040	Human, Mouse, Rat	WB,IHC-P
	TMEM119 Recombinant Rabbit mAb [PSH12-98] - BSA and Azide free	HA751453	Human, Mouse, Rat	WB,IF-Cell,FC,IP
Relative expression	TMEM16A Recombinant Rabbit mAb [PSH06-56] - BSA and Azide free	HA751085	Human	WB,IF-Cell,FC
KD	TMEM192 Recombinant Rabbit mAb [JE64-89] - BSA and Azide free	HA750474	Human, Mouse, Rat	WB,IHC-P
	TMP21 Recombinant Rabbit mAb [PSH01-40] - BSA and Azide free	HA750718	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	TMP21 Recombinant Rabbit mAb [PSH01-41] - BSA and Azide free	HA750719	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC
	TNFAIP3 Recombinant Rabbit mAb [SN07-31] - BSA and Azide free	HA750253	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free	HA750897	Human, Mouse	WB
Relative expression	TNF Receptor II Recombinant Rabbit mAb [PSH12-12] - BSA and Azide free	HA751441	Human	FC
	TNF Receptor II Recombinant Rabbit mAb [JM113-01] - BSA and Azide free	HA750386	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP,IHC-Fr
Cell treatment, Relative expression	TNFRSF14 Recombinant Rabbit mAb [PSH06-02] - BSA and Azide free	HA751042	Human	WB,IF-Cell,FC
Cell treatment, Relative expression	TNFRSF14 Recombinant Rabbit mAb [PSH05-56] - BSA and Azide free	HA751018	Human	WB
KD	TNKS1 Recombinant Rabbit mAb [PSH0-51] - BSA and Azide free	HA750628	Human, Mouse, Rat	WB
	TNS4 Recombinant Rabbit mAb [PD00-95] - BSA and Azide free	HA750540	Human	IHC-P
KD	TOM1 Recombinant Rabbit mAb [PSH03-38] - BSA and Azide free	HA750872	Human, Mouse, Rat	WB,IF-Tissue,IHC-P
	TOM70 Recombinant Rabbit mAb [PSH13-00] - BSA and Azide free	HA751455	Human, Mouse, Rat, Monkey	WB,IHC-P,FC,IF-Cell,IP
Relative expression, KD	TOMM20 Recombinant Rabbit mAb [ST04-72] - BSA and Azide free	HA750166	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP,IHC-Fr
	TOMM20 Recombinant Antibody [ST04-72] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610288	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	TOMM40 Recombinant Rabbit mAb [PSH09-30] - BSA and Azide free	HA751284	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
KD	TopBP1 Recombinant Rabbit mAb [PSH01-14] - BSA and Azide free	HA750690	Human	WB,IHC-P
	Topoisomerase I Recombinant Rabbit mAb [SC69-03] - BSA and Azide free	HA750227	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Total OXPHOS Antibody Cocktail Recombinant Rabbit Multiclonal Antibody [PSH15-50]	HA723756	Human	WB,IHC-P
Relative expression	TPH2 Recombinant Antibody [PSH04-89] - Rat IgG1 (Chimeric)	HA601494	Mouse, Rat	WB,IHC-P
Relative expression	TPH2 Recombinant Antibody [PSH04-89] - Guinea pig IgG2 (Chimeric)	HA601492	Mouse, Rat	IHC-Fr,WB
Relative expression	TPH2 Recombinant Rabbit mAb [PSH04-89] - BSA and Azide free	HA750976	Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IP
Relative expression	TPH2 Recombinant Antibody [PSH04-89] - Mouse IgG1 (Chimeric)	HA601493	Mouse, Rat	WB,IHC-P
Relative expression	TRAF1 Recombinant Rabbit mAb [PSH13-88] - BSA and Azide free	HA751490	Human, Rat	WB,IP
	TRAF2 Recombinant Rabbit mAb [SD205-06] - BSA and Azide free	HA750297	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
	TRAF3 Recombinant Rabbit mAb [PSH01-30] - BSA and Azide free	HA750708	Human, Mouse, Rat	WB
Relative expression	TRAIL Recombinant Rabbit mAb [PSH14-40] - BSA and Azide free	HA751522	Human	WB,FC
	Transcription factor AP-2-alpha Recombinant Rabbit mAb [PSH02-05] - BSA and Azide free	HA750775	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Transcription factor AP-2-alpha Recombinant Rabbit mAb [PSH02-06] - BSA and Azide free	HA750776	Human	WB,IHC-P,IF-Cell
KD	Transferrin Receptor (CD71) Recombinant Rabbit mAb [JF0956] - BSA and Azide free	HA750333	Human, Mouse, Rat	WB,IP,FC
	Transketolase Recombinant Rabbit mAb [PSH13-70] - BSA and Azide free	HA751483	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	Transmembrane protein 200A Recombinant Mouse mAb [4-B3-R] - BSA and Azide free	HA610071	Human	WB
Relative expression	TREX1 Recombinant Rabbit mAb [PSH08-84] - BSA and Azide free	HA751249	Mouse	WB,FC,IP
	TRIB3 Recombinant Rabbit mAb [PSH02-13] - BSA and Azide free	HA750783	Human	WB
	TRIF Recombinant Rabbit mAb [PSH05-39] - BSA and Azide free	HA750998	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment, KD	TRIM21 Recombinant Rabbit mAb [PSH02-56] - BSA and Azide free	HA750819	Human	WB,IF-Cell,IHC-P,FC,IP
	TRIM29 Recombinant Rabbit mAb [PSH02-99] - BSA and Azide free	HA750854	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	TRIM32 Recombinant Rabbit mAb [PSH08-78] - BSA and Azide free	HA751243	Human, Mouse	WB
KD	TRIM56 Recombinant Rabbit mAb [PSH02-17] - BSA and Azide free	HA750787	Human, Mouse, Rat, Monkey	WB,FC
	TRIM72 Recombinant Mouse mAb [A1D4-R] - BSA and Azide free	HA610103	Human, Mouse, Rat	WB
	Triosephosphate isomerase Recombinant Rabbit mAb [PSH14-38] - BSA and Azide free	HA751520	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
Relative expression	TrkA Recombinant Rabbit mAb [SU0354] - BSA and Azide free	HA750148	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Relative expression	TrkC Recombinant Rabbit mAb [PSH05-40] - BSA and Azide free	HA750999	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
	TROP2 Recombinant Mouse mAb [A10G7-R] - BSA and Azide free	HA610138	Human	IHC-P,WB,IF-Cell,FC
	TROP2 Recombinant Mouse mAb [A5E2-R] - BSA and Azide free	HA610157	Human	WB,IF-Cell,IHC-P,FC
KD	Tropomyosin 2 Recombinant Rabbit mAb [PSH02-63] - BSA and Azide free	HA750820	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell
Relative expression	TRPV1 Recombinant Mouse mAb [PSH13-22] - BSA and Azide free	HA610273	Rat	IHC-P
KD	Tryptophan Hydroxylase 1 (TPH1) Recombinant Rabbit mAb [SC53-07] - BSA and Azide free	HA750214	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	TSG101 Recombinant Rabbit mAb [JJ0900] - BSA and Azide free	HA750316	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Twist1 Recombinant Rabbit mAb [PSH11-98] - BSA and Azide free	HA751431	Human, Mouse, Rat	WB
	TXNDC9 Recombinant Rabbit mAb [PSH04-43] - BSA and Azide free	HA750952	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
KO	TXNIP Recombinant Rabbit mAb [JM60-35] - BSA and Azide free	HA750438	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	TXNRD2 Recombinant Rabbit mAb [PSH02-26] - BSA and Azide free	HA750794	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	TXNRD2 Recombinant Rabbit mAb [PSH02-27] - BSA and Azide free	HA750795	Human, Mouse, Rat	WB,IHC-P
	Tyrosine Hydroxylase Recombinant Rabbit mAb [SD080-02] - BSA and Azide free	HA750298	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IF-Tissue,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Tyrosine Hydroxylase Recombinant Rabbit mAb [SN59-03] - BSA and Azide free	HA750243	Human, Mouse, Rat, Cynomolgus monkey	WB,IF-Cell,IHC-P,FC,IF-Tissue,IHC-Fr
	UBA3 Recombinant Mouse mAb [1B5-2-3-R] - BSA and Azide free	HA610041	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	UBA3 Recombinant Mouse mAb [1C10-5-5-R] - BSA and Azide free	HA610109	Human, Mouse, Rat	WB,IHC-P
	UBA3 Recombinant Mouse mAb [1B6-6-5-R] - BSA and Azide free	HA610083	Human, Mouse, Rat	WB
	Ubiquinol-Cytochrome C Reductase Core Protein I Recombinant Rabbit mAb [PSH01-23] - BSA and Azide free	HA750699	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell
Cell treatment	Ubiquitin Recombinant Rabbit mAb [ST46-03] - BSA and Azide free	HA750165	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	UFSP2 Recombinant Rabbit mAb [PSH03-57] - BSA and Azide free	HA750891	Human, Mouse, Rat	WB,IHC-P
	UHRF1 Recombinant Rabbit mAb [PSH10-63] - BSA and Azide free	HA751366	Human, Mouse, Monkey	WB,IF-Cell,FC,IP,ChIP
Relative expression	ULBP3 Recombinant Rabbit mAb [PSH13-03] - BSA and Azide free	HA751458	Human	WB,IHC-P,IP
	ULK1 Recombinant Rabbit mAb [JA58-36] - BSA and Azide free	HA750418	Human, Mouse, Rat	WB,IHC-P
KD	ULK3 Recombinant Rabbit mAb [PSH01-09] - BSA and Azide free	HA750685	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue
Cell treatment	uPA Receptor Recombinant Mouse mAb [A6E8-R] - BSA and Azide free	HA610169	Human	WB
	V5 tag Recombinant Rabbit mAb [SY30-01] - BSA and Azide free	HA750086	Species independent	WB,IF-Cell,IP
	VACHT / SLC18A3 Recombinant Antibody [PSH10-50] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610250	Mouse, Species independent, Pig	IHC-Fr
Relative expression	VACHT / SLC18A3 Recombinant Rabbit mAb [PSH10-50] - BSA and Azide free	HA751353	Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IHC-P
	VACHT / SLC18A3 Recombinant Antibody [PSH10-50] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610311	Mouse, Rat	IHC-Fr,IHC-P
	VAMP8 Recombinant Rabbit mAb [JF0963] - BSA and Azide free	HA750334	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	Vasopressin Recombinant Rabbit mAb [PSH11-64] - BSA and Azide free	HA751420	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,IF-Tissue
Relative expression, KD	VCAM1 Recombinant Rabbit mAb [SA05-04] - BSA and Azide free	HA750009	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	VDAC1 Recombinant Rabbit mAb [SA93-03] - BSA and Azide free	HA750012	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	VDAC1 / 2 / 3 Recombinant Rabbit mAb [PSH14-90]	HA723692	Human, Mouse, Rat	WB
KD	VE Cadherin Recombinant Rabbit mAb [PSH03-15] - BSA and Azide free	HA750856	Human	WB,IF-Cell,FC
	VEGF Recombinant Rabbit mAb [SP07-01] - BSA and Azide free	HA750074	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	VEGF Receptor 1 Recombinant Rabbit mAb [SY09-09] - BSA and Azide free	HA750081	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	VGLUT1 Recombinant Antibody [PSH10-53] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610283	Mouse, Rat	IHC-Fr,WB
	VGLUT1 Recombinant Mouse mAb [PSH08-34] - BSA and Azide free	HA610217	Mouse, Rat	WB,IHC-P
Relative expression	VGLuT1 Recombinant Antibody [PSH10-53] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610234	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB,IF-Tissue
Relative expression	VGLUT1 Recombinant Rabbit mAb [PSH10-53] - BSA and Azide free	HA751356	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IHC-P,IF-Tissue
	VGLUT1 Recombinant Antibody [PSH10-53] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610284	Mouse, Rat	IHC-Fr,IHC-P,WB
	VGLUT1 Recombinant Antibody [PSH10-53] - Chicken IgY (Chimeric) - BSA and Azide free	HA610312	Mouse	IHC-Fr
Relative expression	VGLUT2 Recombinant Rabbit mAb [PSH10-54] - BSA and Azide free	HA751357	Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IHC-P,IF-Tissue
	VGLUT2 Recombinant Antibody [PSH10-54] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610285	Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	VGLUT2 Recombinant Antibody [PSH10-54] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610235	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,WB
	VILIP1 Recombinant Rabbit mAb [PSH0-08] - BSA and Azide free	HA750577	Human, Mouse, Rat	WB,IHC-P
	Villin1 Recombinant Rabbit mAb [C4] - BSA and Azide free	HA750476	Human	WB,IHC-P
	Villin1 Recombinant Rabbit mAb [A2] - BSA and Azide free	HA750477	Human	WB,IHC-P
Relative expression	Vimentin Recombinant Mouse mAb [A6-C1-R] - BSA and Azide free	HA610126	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KO, Relative expression	Vimentin Recombinant Rabbit mAb [SC60-05] - BSA and Azide free	HA750215	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC,IP
KD	Vinculin Recombinant Rabbit mAb [JM42-43] - BSA and Azide free	HA750448	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	VIP Recombinant Antibody [PSH10-56] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610280	Mouse	IHC-Fr
	VIP Recombinant Antibody [PSH10-56] - Chicken IgY (Chimeric) - BSA and Azide free	HA610313	Mouse	IHC-Fr
	VIP Recombinant Rabbit mAb [PSH10- 56] - BSA and Azide free	HA751359	Human, Mouse, Rat, Cynomolgus monkey	IHC-Fr,IHC-P,IF- Tissue
	VIP Recombinant Antibody [PSH10-56] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610241	Human, Mouse, Rat, Cynomolgus monkey	IHC-Fr,IHC-P
	VIRMA Recombinant Rabbit mAb [PSH08-28] - BSA and Azide free	HA751219	Human, Rat, Monkey	WB,IF-Cell,IHC- P,FC,IP
	VISTA Recombinant Rabbit mAb [PSH0-65] - BSA and Azide free	HA750642	Human	WB,IHC-P,IF- Cell,FC,IF-Tissue
Relative expression, KD, Cell treatment	VISTA Recombinant Rabbit mAb [PSH0-64] - BSA and Azide free	HA750641	Human	WB,IHC-P,IF-Tissue
KO	Vitamin D Receptor Recombinant Rabbit mAb [JA11-16] - BSA and Azide free	HA750402	Human, Mouse, Rat	WB,IHC-P,FC
Relative expression	Von Willebrand Factor Recombinant Rabbit mAb [PSH07-45] - BSA and Azide free	HA751166	Human, Mouse, Rat	WB,IHC-P,IF- Tissue,IP
KD	VPS34 Recombinant Rabbit mAb [PSH10-64] - BSA and Azide free	HA751367	Human, Mouse, Rat	WB
	VSIG4 Recombinant Rabbit mAb [PSH06-55] - BSA and Azide free	HA751083	Human	WB,IHC-P
	VTI1B Recombinant Rabbit mAb [PSH06-40] - BSA and Azide free	HA751080	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	WDR77 Recombinant Rabbit mAb [PSH0-07] - BSA and Azide free	HA750576	Human	WB,IHC-P,FC
	Wee1 Recombinant Rabbit mAb [PSH05-11] - BSA and Azide free	HA750984	Human, Monkey	WB,IHC-P
KD	Wilms Tumor Protein Recombinant Rabbit mAb [SC06-41] - BSA and Azide free	HA750219	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P
	WNT8A Recombinant Rabbit mAb [PSH0-39] - BSA and Azide free	HA750617	Human, Mouse, Rat	WB,IF-Cell,FC
	xCT / SLC7A11 Recombinant Mouse mAb [A7C6-R] - BSA and Azide free	HA610011	Human, Mouse	WB,IHC-P,IF-Cell,FC
KD	XIAP Recombinant Rabbit mAb [PSH04-30] - BSA and Azide free	HA750939	Human, Mouse, Rat, Monkey	WB,IHC-P
	XLF Recombinant Rabbit mAb [PSH01-07] - BSA and Azide free	HA750682	Human	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	YTHDC1 Recombinant Rabbit mAb [PSH06-27] - BSA and Azide free	HA751064	Human, Mouse, Rat	WB,IF-Cell,FC
	YTHDF1 Recombinant Rabbit mAb [PSH0-23] - BSA and Azide free	HA750597	Human, Mouse, Rat	WB,IHC-P,FC
	YTHDF2 Recombinant Rabbit mAb [PSH06-17] - BSA and Azide free	HA751055	Human, Mouse, Rat, Monkey	WB,IHC-P,FC
Relative expression	ZAP70 Recombinant Mouse mAb [A1B5-R] - BSA and Azide free	HA610096	Human, Mouse, Rat	WB,IHC-P
	ZEB Recombinant Rabbit mAb [PSH0-61] - BSA and Azide free	HA750638	Human, Mouse, Rat	WB,IHC-P,FC
KD	ZIP Kinase Recombinant Rabbit mAb [PSH02-18] - BSA and Azide free	HA750788	Human	WB
	ZO1 Recombinant Rabbit mAb [PSH07-09] - BSA and Azide free	HA751131	Human, Mouse, Monkey, Cynomolgus monkey	WB,IF-Cell
	ZW10 Recombinant Rabbit mAb [PSH05-80] - BSA and Azide free	HA751028	Human, Mouse, Rat, Monkey	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层