

血液肿瘤研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	AKT1 Mouse mAb [D9-9-C9]	EM40507	Human, Mouse, Rat	WB,IHC-P
	AKT1 Recombinant Rabbit mAb [PSH0-50]	HA721364	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
KD	AKT1 Recombinant Rabbit mAb [ST05-09] - BSA and Azide free	HA750177	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
KD	AKT1 Recombinant Rabbit mAb [ST05-09]	ET1609-47	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
	AKT1 Recombinant Rabbit mAb [PSH0-50] - BSA and Azide free	HA750627	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
	AKT1/2/3 Recombinant Rabbit mAb [ST48-09] - BSA and Azide free	HA750179	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
	AKT1/2/3 Rabbit pAb	ER1901-33	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75-09] - BSA and Azide free	HA750840	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75-09]	HA721870	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	AKT1/2/3 Recombinant Rabbit mAb [ST48-09]	ET1609-51	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
	Akt1/3 Rabbit pAb	ER62638	Human, Mouse, Rat	WB,IHC-P
KD	Bcl-2 Recombinant Rabbit mAb [PD01-41]	HA721235	Human	WB,IHC-P,IF-Cell,FC
	Bcl-2 Mouse mAb [F7-1-A1]	M1206-4	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Bcl-2 Rabbit pAb	ER1802-97	Human, Mouse	WB,IF-Cell,FC
	Bcl-2 Recombinant Rabbit mAb [PSH09-49] - BSA and Azide free	HA751296	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	Bcl-2 Mouse mAb [9F3]	EM1701-83	Human	WB,IHC-P,FC
	Bcl-2 Recombinant Mouse mAb [9F3-R] - BSA and Azide free	HA610131	Human	WB,IHC-P,FC
KD	Bcl-2 Recombinant Rabbit mAb [SZ10-03] - BSA and Azide free	HA750057	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Bcl-2 Rabbit pAb	ER1706-47	Human	WB,IF-Cell,IHC-P,FC
KD	Bcl-2 Recombinant Rabbit mAb [SZ10-03]	ET1603-11	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Bcl-2 Recombinant Rabbit mAb [PSH09-49]	HA723111	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	Bcl-2 Mouse mAb [9F1]	EM1701-82	Human	WB,IHC-P
KD	Bcl-2 Recombinant Rabbit mAb [JF104-8] - BSA and Azide free	HA750350	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	Bcl-2 Recombinant Rabbit mAb [JF104-8]	ET1702-53	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	FITC Conjugated Bcl-2 Recombinant Rabbit mAb [JF104-8]	HA720210F	Human, Mouse	FC
	Bcl-2 Rabbit pAb	ER0602	Human	WB,IF-Cell,IHC-P,FC
	Bcl-2 Recombinant Mouse mAb [9F3-R]	HA601256	Human	WB,IHC-P,FC
	BCMA Recombinant Rabbit mAb [PSH16-80]	HA723880	Human	FC
Relative expression	BCMA Recombinant Rabbit Multiclonal Antibody [PSH18-02]	HA723967	Human	IF-Cell,FC
Relative expression	BCMA Recombinant Rabbit mAb [PSH17-94]	HA723959	Human	FC
	BCMA Recombinant Rabbit mAb [PSH16-80] - BSA and Azide free	HA751607	Human	FC
KO, Relative expression	BTK Recombinant Rabbit mAb [SN69-07] - BSA and Azide free	HA750252	Human, Mouse	WB,IHC-P,IP
KO, Relative expression	BTK Recombinant Rabbit mAb [SN69-07]	ET1611-4	Human, Mouse	WB,IHC-P,IP
	Casein Kinase 1 alpha Recombinant Rabbit mAb [PSH04-46]	HA722131	Human, Mouse, Rat	WB,IF-Cell,FC,IP
	Casein Kinase 1 alpha Recombinant Rabbit mAb [PSH04-46] - BSA and Azide free	HA750955	Human, Mouse, Rat	WB,IF-Cell,FC,IP
	Casein Kinase I alpha Rabbit pAb	ER62901	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Casein Kinase I alpha 1/1L Rabbit pAb	ER62900	Human, Mouse, Rat	IHC-P,IF-Cell,WB
	CCR9 Recombinant Rabbit mAb [JE00-35]	HA722444	Human, Mouse, Rat	WB
	CD109 Rabbit pAb	ER65770	Human	IHC-P
Relative expression	CD13 Recombinant Rabbit mAb [PSH08-48]	HA722994	Human	WB,IF-Cell,IHC-P,FC
Relative expression, KD	CD13 Recombinant Rabbit mAb [PD00-39]	HA721160	Human, Mouse	IHC-P,IF-Tissue,WB
	CD13 Rabbit pAb	ER1803-74	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	CD13 Recombinant Rabbit mAb [PSH08-48] - BSA and Azide free	HA751223	Human	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD, Relative expression	CD13 Recombinant Rabbit mAb [SN71-04]	ET1611-61	Human, Mouse, Rat	WB,IHC-P,IP
KD, Relative expression	CD13 Recombinant Rabbit mAb [SC70-01]	ET1610-59	Human, Mouse, Rat	WB,IHC-P,IP
	CD15 Recombinant Rabbit mAb [PD00-42] - BSA and Azide free	HA750565	Human	IHC-P,FC,IF-Cell
	FITC Conjugated CD15 Recombinant Rabbit mAb [PD00-42]	HA720195F	Human	FC
	CD15 Recombinant Rabbit mAb [PD00-42]	HA721246	Human	IHC-P,mIHC,FC,IF-Cell
	CD16 Recombinant Rabbit mAb [JE49-79]	ET7109-97	Human, Mouse, Rat	WB
	CD19 Recombinant Rabbit mAb [JF099-9]	ET1702-74	Human, Mouse, Rat	WB,FC
	PE Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723763F	Human	FC
Relative expression	CD19 Recombinant Rabbit mAb [PSH13-44] - BSA and Azide free	HA751463	Human	WB,IF-Cell,FC
Relative expression	CD19 Recombinant Rabbit mAb [JF100-06] - BSA and Azide free	HA750365	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	CD19 Recombinant Rabbit mAb [PSH13-44]	HA723533	Human	WB,IF-Cell,FC
	CD19 Recombinant mAb	IRS110RB	Mouse	mIHC
	CD19 Recombinant Rabbit mAb	IRS018	Human	mIHC
	CD19 Mouse mAb [A3-A8]	EM40308	Human, Mouse	WB,IHC-P,FC
Relative expression	APC Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723761F	Human	FC
	FITC Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723762F	Human	FC
Relative expression	CD19 Recombinant Rabbit mAb [JF100-06]	ET1702-93	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,mIHC
Relative expression	CD19 Recombinant Rabbit mAb [PSH03-94] - BSA and Azide free	HA750919	Mouse	WB,IHC-P,IF-Tissue,IP,IHC-Fr,FC,IF-Cell
	CD19 Rabbit pAb	R1107-7	Human, Mouse	WB
Relative expression	CD19 Recombinant Rabbit mAb [PSH03-94]	HA722073	Mouse	WB,IHC-P,IF-Tissue,IP,IHC-Fr,FC,IF-Cell,mIHC
	FITC Conjugated CD19 Recombinant Mouse mAb [PSH14-97]	HA601518F	Human	FC
	CD20 Rabbit pAb	IRS114RB	Mouse	mIHC
Relative expression	CD20 Recombinant Antibody [PD00-02] - Mouse IgG1 (Chimeric)	HA601544	Human, Mouse	WB,IHC-P,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD20 Recombinant Rabbit mAb [PSH08-12] - BSA and Azide free	HA751207	Human	WB,IF-Cell,IHC-P,FC,IP
Relative expression	CD20 Recombinant Rabbit mAb [PSH08-12]	HA722953	Human	WB,IF-Cell,IHC-P,FC,IP
	CD20 Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS112MS	Human	mIHC
Relative expression	CD20 Recombinant Antibody [PD00-02] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610360	Human, Mouse	WB,IHC-P,IHC-Fr
	APC Conjugated CD20 Recombinant Rabbit mAb [PSH08-12]	HA723838F	Human	FC
Relative expression	CD20 Recombinant Rabbit mAb [PD00-02]	HA721138	Human, Mouse	IHC-P,mIHC,WB,IF-Cell,IP,IHC-Fr
	CD20 Recombinant Rabbit mAb	IRS003	Human	mIHC
Relative expression	CD20 Recombinant Rabbit mAb [SY12-01]	ET1605-33	Human, Monkey	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP
	CD20 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS113RT	Human	mIHC
	PE Conjugated CD20 Recombinant Rabbit mAb [PSH08-12]	HA723839F	Human	FC
	CD22 Recombinant Rabbit mAb [PO00-31]	HA721157	Human	WB,IHC-P,IF-Tissue
	CD22 Mouse mAb [A2G6]	EM1902-13	Human	WB,IHC-P,FC
	CD229 Rabbit pAb	ER65519	Human	WB,IHC-P
	CD269 Rabbit pAb	ER65683	Human	WB
	CD28 Rabbit pAb	ER2001-42	Human, Mouse, Rat	WB,IHC-P,FC
	CD3 Recombinant Rabbit mAb	IRS074	Mouse	mIHC
	CD3 Recombinant Rabbit mAb [JE80-02]	HA720082	Human	WB,IHC-P,IF-Tissue,IP,mIHC,IF-Cell
	CD3 Recombinant Mouse mAb [PD01-43]	HA601110	Human	WB,IHC-P,IF-Tissue
Relative expression	CD3 Recombinant Rabbit mAb [PSH19-70]	HA724094	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	CD3 Recombinant Rabbit mAb	IRS022	Human	mIHC
Relative expression	CD30 Recombinant Rabbit mAb [PSH04-11]	HA722092	Human	WB,IF-Cell,IHC-P
Relative expression	CD30 Recombinant Rabbit mAb [JE57-71]	HA721579	Human	WB,IHC-P
Relative expression	CD30 Recombinant Rabbit mAb [PSH04-10]	HA722091	Human	WB,IHC-P,FC,IF-Tissue
Relative expression	CD30 Recombinant Rabbit mAb [PSH04-12]	HA722093	Human	WB,IF-Cell,IHC-P,IF-Tissue
	CD327 Mouse mAb [3C11C12G1]	EM1711-43	Human	WB,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD327 Rabbit pAb	ER65778	Human	IHC-P
	CD33 Rabbit pAb	ER1905-89	Mouse	WB,IHC-P,IF-Tissue
Relative expression	CD33 Recombinant Rabbit mAb [PSH02-51]	HA721827	Human	WB,IF-Cell,FC
	CD33 Mouse mAb [B7-C6]	EM1709-20	Human	IHC-P,FC
Relative expression	CD33 Recombinant Rabbit mAb [PSH02-51] - BSA and Azide free	HA750814	Human	WB,IF-Cell,FC
Relative expression, Cell treatment	CD33 Recombinant Rabbit mAb [PSH02-50]	HA721826	Human	WB,IF-Cell,FC
Relative expression, Cell treatment	CD33 Recombinant Rabbit mAb [PSH02-50] - BSA and Azide free	HA750813	Human	WB,IF-Cell,FC
	CD34 Recombinant Rabbit mAb	IRS024	Human	mIHC
	CD34 Rabbit pAb	ER0802	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	CD34 Recombinant Rabbit mAb [SI16-01] - BSA and Azide free	HA750092	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	CD34 Mouse mAb [15H1]	EM1901-01	Human, Mouse, Rat	IHC-P,FC
	CD34 Recombinant Mouse mAb [PDM0-12]	HA601109	Human	IHC-P,IF-Tissue
	CD34 Rabbit pAb	0108-2	Human, Mouse, Zebrafish	WB,IF-Cell,IHC-P
	CD34 Recombinant Rabbit mAb [SI16-01]	ET1606-11	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,mIHC,FC
	CD34 Rabbit pAb	ER1706-66	Human, Mouse	WB,IHC-P
	CD34 Recombinant Rabbit mAb	IRS050	Mouse	mIHC
	CD38 Recombinant Antibody [PD01-49] - Mouse IgG1 (Chimeric)	HA601554	Human	WB,IHC-P,FC
	CD38 Recombinant Antibody [PD01-49] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610369	Human	WB,IHC-P,FC
	CD38 Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS120MS	Human	mIHC
	CD38 Recombinant Rabbit mAb [JE29-35]	ET7110-53	Human	WB,IHC-P,FC,IF-Tissue
	CD38 Recombinant Rabbit mAb	IRS006	Human	mIHC
	CD38 Recombinant Rabbit mAb [PSH20-33]	HA724144	Mouse	WB
	CD38 Rabbit pAb	ER1905-91	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	CD38 Recombinant Rabbit mAb [SN07-20]	ET1611-66	Human	WB,IHC-P,FC,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD38 Recombinant Antibody [PD01-49] - Rat IgG1 (Chimeric)	HA601417		
Relative expression	CD38 Recombinant Rabbit mAb [PD01-49]	HA721268	Human	WB,IHC-P,IF-Cell,FC,mIHC,IF-Tissue
	CD38 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS119RT	Human	mIHC
Relative expression	CD38 Recombinant Antibody [PD01-49] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610257		
	CD3D Recombinant Rabbit mAb [JJ08-97]	ET1701-57	Human	WB,IF-Cell,IHC-P,FC,IP,IF-Tissue
	CD3E Rabbit pAb	HA500344	Human, Mouse, Rat	WB
	CD3 epsilon Recombinant Rabbit mAb [SY0239] - BSA and Azide free	HA750115	Human	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
	CD3 epsilon Recombinant Rabbit mAb [JE17-74]	ET7110-63	Human	WB,IHC-P,FC
	CD3 epsilon Rabbit pAb	ER80501	Human, Mouse, Rat	WB
	CD3 epsilon Recombinant Rabbit mAb [SY0239]	ET1607-29	Human	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
	CD3G Recombinant Rabbit mAb [JB38-29]	ET7107-55	Human, Mouse	WB,IHC-P,FC,IP,IF-Cell,IF-Tissue
	CD3 zeta Rabbit pAb	0805-7	Human	WB,IF-Cell,FC
Relative expression	CD3 zeta Recombinant Rabbit mAb [JE09-37]	HA722692	Human, Rat	WB,IHC-P,IF-Tissue
KD	CD4 Rabbit pAb	0108-6	Human, Mouse, Rat	WB,IF-Cell,FC
Relative expression	CD4 Recombinant Antibody [PSH08-22] - Mouse IgG1 (Chimeric)	HA601533	Mouse	WB,IF-Tissue,IHC-P,IF-Cell
	CD4 Recombinant Rabbit mAb	IRS005	Human	mIHC
Relative expression	CD4 Recombinant Rabbit mAb [PSH08-22]	HA722966	Mouse	WB,IF-Cell,IHC-P,mIHC,IHC-Fr
KD	CD4 Recombinant Rabbit mAb [ST0488]	ET1609-52	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,mIHC
	CD4 Rabbit pAb	ER1803-84	Human	WB,FC
Relative expression	CD4 Recombinant Rabbit mAb [PSH08-22] - BSA and Azide free	HA751213	Mouse	WB,IF-Cell,IHC-P,IHC-Fr
	CD4 Recombinant Antibody [PSH08-22] - Rat IgG1 (Chimeric)	HA601534	Mouse	WB,IF-Tissue,IHC-P,IHC-Fr
	CD4 Recombinant Rabbit mAb [JE56-36]	ET7110-64	Human	WB,IHC-P
Relative expression	CD4 Recombinant Antibody [PSH08-22] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610353	Mouse	WB,IF-Tissue,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD4 Recombinant Antibody [PSH08-22] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610354	Mouse	WB,IF-Tissue,IHC-P,IHC-Fr
	CD4 Recombinant Rabbit Multiclonal Antibody [PSH13-24] - BSA and Azide free	HA751460	Human, Mouse	WB,IF-Cell,IHC-P
	CD4 Recombinant Rabbit mAb	IRS051	Mouse	mIHC
	CD4 Recombinant Antibody [ST0488] - Mouse IgG1 (Chimeric)	HA601557	Human	WB,IF-Cell,IHC-P
	CD4 Rabbit pAb	R1211-1	Human	WB,IF-Cell
	CD4 Recombinant Rabbit Multiclonal Antibody [PSH13-24]	HA723511	Human, Mouse	WB,IF-Cell,IHC-P
Relative expression	CD44 Recombinant Rabbit mAb [ST57-03]	ET1609-74	Human	WB,IHC-P
	CD44 Recombinant mAb	IRS122RB	Human	mIHC
	CD44 Mouse mAb [A7-1]	EM010805	Human	WB,IHC-P,FC
Relative expression	CD44 Recombinant Rabbit mAb [ST57-03] - BSA and Azide free	HA750192	Human	WB,IHC-P
	CD44 Recombinant Rabbit mAb	IRS123RB	Mouse	mIHC
Relative expression, KD	CD44 Recombinant Rabbit mAb [JE64-01]	HA722845	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,FC
	CD47 Mouse mAb [1A9]	EM1701-38	Human	IF-Cell,IHC-P,FC
	iFluor™ 488 Conjugated CD47 Mouse mAb [1B1]	HA600104F	Human	FC
	CD47 Recombinant Rabbit mAb [PD00-33]	HA721158	Human	IHC-P,FC
Relative expression	CD47 Recombinant Rabbit mAb [PSH10-31] - BSA and Azide free	HA751348	Human	WB,IF-Cell,FC
Relative expression	CD47 Recombinant Rabbit mAb [PSH10-31]	HA723198	Human	WB,IF-Cell,FC
	CD5 Recombinant Rabbit mAb [JE25-41]	ET7110-56	Human	WB,IHC-P
	CD5 Recombinant Rabbit mAb [PD00-01]	HA721137	Human	IHC-P,IF-Tissue,WB
Relative expression	CD7 Recombinant Rabbit mAb [JE57-72]	HA720023	Human	WB,IF-Cell,IHC-P
	CD70 Rabbit pAb	ER65832	Human, Mouse	IHC-P
	CD70 Rabbit pAb	ER62948	Human	WB
	CD79b Recombinant Rabbit mAb [JE58-34]	HA720043	Human	WB,IHC-P,FC
	CD8 Rabbit pAb	0108-7	Mouse	WB,IF-Cell,FC
	CD8 Rabbit pAb	IRS083	Mouse	mIHC
	CD8a Mouse mAb	IRS007	Human	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD8 alpha Recombinant Rabbit mAb [PSH19-81]	HA724105	Human, Rat	WB,IHC-P,IF-Tissue,FC
	CD8 alpha Recombinant Rabbit mAb [SI18-01] - BSA and Azide free	HA750101	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	CD8 alpha Recombinant Rabbit mAb [SI18-01]	ET1606-31	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	CD99 Recombinant Rabbit mAb [PD00-13]	HA721150	Human	IHC-P,IF-Tissue
	CD99 Rabbit pAb	ER1803-81	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	CD99 Recombinant Rabbit mAb [JF0991]	ET1702-35	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	CDC123 Recombinant Rabbit mAb [JE53-14]	ET7110-32	Human, Rat	WB
	Cdk11A/B Rabbit pAb	ER62965	Human	WB
	CEBP alpha Recombinant Rabbit mAb [SD202-09]	ET1612-46	Human, Mouse, Rat	WB,IP
Relative expression	CEBP alpha Recombinant Rabbit mAb [PSH19-55]	HA724081	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
Relative expression	CEBP alpha Recombinant Rabbit mAb [PSH20-43]	HA724154	Human, Rat	WB,IHC-P,IF-Cell,FC,IP
	c-kit Rabbit pAb	0108-4	Human	WB,IF-Cell
Relative expression	c-Kit Recombinant Rabbit mAb [PSH17-79]	HA723948	Mouse	IHC-Fr,IF-Cell,FC
Relative expression, KD	c-Kit Recombinant Rabbit mAb [PD00-24]	HA721154	Human	IHC-P,WB,mIHC
Relative expression	c-Kit Recombinant Rabbit mAb [PSH17-96]	HA723961	Mouse, Rat	WB,IHC-Fr,IHC-P
Relative expression, KD	c-Kit Recombinant Rabbit mAb [ST04-99]	ET1609-60	Human	WB,IHC-P
	c-Kit Recombinant Rabbit mAb	IRS013	Human	mIHC
	CL12A Rabbit pAb	ER61745	Human, Mouse	WB
	CLL1 Mouse mAb [4D5A8]	EM1711-46	Human	WB,IF-Cell,FC
	CLSTN1 Recombinant Rabbit mAb [JE48-81]	ET7109-84	Human, Mouse, Rat	WB
	CRM1 Recombinant Rabbit mAb [JB35-22]	ET7107-27	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
	CSNK1A1 + CSNK1A1L Recombinant Rabbit mAb [PSH04-58] - BSA and Azide free	HA750965	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	CSNK1A1 + CSNK1A1L Recombinant Rabbit mAb [PSH04-58]	HA722157	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Fli-1 Rabbit pAb	ER63441	Human, Mouse	WB,IHC-P
	FLI1 Recombinant Rabbit mAb [PSH14-17]	HA723605	Human	WB,IHC-P
	FLI1 Rabbit pAb	HA500149	Human, Mouse, Rat	IHC-P,WB
	Flt3 Rabbit pAb	ER2001-25	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	FLT3 Mouse mAb [C3-B7]	EM1709-22	Human	WB,IHC-P
	GPRC5D Rabbit pAb	ER63621	Human, Mouse	WB,IF-Cell
	GPX4 Rabbit pAb	ER1803-15	Mouse, Rat	WB
KD	GPX4 Recombinant Rabbit mAb [JU11-31]	ET1706-45	Human, Mouse, Rat, Zebrafish, Pig, Cynomolgus monkey	WB,IHC-P,IF-Cell,FC
KD	GPX4 Recombinant Rabbit mAb [JU11-31] - BSA and Azide free	HA750455	Human, Mouse, Rat, Zebrafish, Pig, Cynomolgus monkey	WB,IHC-P,IF-Cell,FC
KD	HDAC3 Recombinant Rabbit mAb [SC55-02] - BSA and Azide free	HA750221	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
KD	HDAC3 Mouse mAb [B2-G9]	M1511-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	HDAC3 Recombinant Rabbit mAb [SC55-02]	ET1610-5	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Human BCMA Recombinant Rabbit mAb [PSH13-18] - BSA and Azide free (Capture)	HA725082	Human	ELISA(Cap)
	Human BCMA Recombinant Rabbit mAb [PSH13-19] - BSA and Azide free (Detector)	HA725083	Human	ELISA(Det)
	Biotin Conjugated Human BCMA Recombinant Rabbit mAb [PSH13-19] - Detector	HA725084B	Human	ELISA(Det),ELISA
	Human CD28 Recombinant Rabbit mAb [PSH12-51] - BSA and Azide free (Detector)	HA723458	Human	ELISA(Det)
	Biotin Conjugated Human CD28 Recombinant Rabbit mAb [PSH12-51] - Detector	HA723459B	Human	ELISA(Det),ELISA
	Human CD28 Recombinant Rabbit mAb [PSH12-50] - BSA and Azide free (Capture)	HA723457	Human	ELISA(Cap)
	Human CD3 Recombinant Mouse mAb [PSH04-92]	HA601337	Human	FC
	APC Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600118F	Human	FC
	FITC Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600116F	Human	FC
	PE Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA600117F	Human	FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Percp-Cy5.5 Conjugated Human CD3 Recombinant Mouse mAb [PSH04-92]	HA601548F	Human	FC
	Human CD30 Recombinant Rabbit mAb [PSH04-81] - BSA and Azide free (Capture)	HA722187	Human	ELISA(Cap)
	Biotin Conjugated Human CD30 Recombinant Rabbit mAb [PSH04-82]	HA722601B	Human	ELISA(Det)
	Human CD30 Recombinant Rabbit mAb [PSH04-82] - BSA and Azide free (Detector)	HA722188	Human	ELISA(Det)
	Human CD33 Recombinant Rabbit mAb [PSH05-26] - BSA and Azide free (Detector)	HA722255	Human	ELISA(Det)
	FITC Conjugated Human CD33 Recombinant Rabbit mAb [PSH04-96]	HA720197F	Human	FC
	Biotin Conjugated Human CD33 Recombinant Rabbit mAb [PSH05-26]	HA722599B	Human	ELISA(Det)
	Human CD33 Recombinant Rabbit mAb [PSH05-25] - BSA and Azide free (Capture)	HA722254	Human	ELISA(Cap)
	FITC Conjugated Human CD38 Recombinant Rabbit mAb [SN07-20]	HA720196F	Human	FC
	PE Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA600120F	Human	FC
	APC Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA600121F	Human	FC
	PE-Cy7 Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA601542F	Human	FC
	FITC Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA600119F	Human	FC
	Human CD4 Recombinant Mouse mAb [PSH04-93]	HA601338	Human	FC
	Biotin Conjugated Human CD47 Recombinant Rabbit mAb [PSH11-06] - Detector	HA723274B	Human	ELISA(Det),ELISA
	Human CD47 Recombinant Rabbit mAb [PSH11-06] - BSA and Azide free (Detector)	HA723273	Human	ELISA(Det)
	Human CD47 Recombinant Rabbit mAb [PSH11-05] - BSA and Azide free (Capture)	HA723272	Human	ELISA(Cap)
Relative expression	Human CD5 Recombinant Rabbit mAb [PSH05-01]	HA722210	Human	IF-Cell,FC
	FITC Conjugated Human CD5 Recombinant Rabbit mAb [PSH05-01]	HA720211F	Human	FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	FITC Conjugated Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA600122F	Human	FC
	APC Conjugated Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA600124F	Human	FC
	PE Conjugated Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA600123F	Human	FC
	APC-Cy7 Conjugated Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA601541F	Human	FC
	Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA601339	Human	FC
	Human PD1 Recombinant Rabbit mAb [SJ01-91] - BSA and Azide free (Capture)	HA723521	Human	ELISA(Cap)
	Biotin Conjugated Human PD1 Recombinant Rabbit mAb [PSH13-38] - Detector	HA723523B	Human	ELISA(Det),ELISA
	Human PD1 Recombinant Rabbit mAb [PSH13-39] - BSA and Azide free (Detector)	HA723524	Human	ELISA(Det)
	Biotin Conjugated Human PD1 Recombinant Rabbit mAb [PSH13-39] - Detector	HA723525B	Human	ELISA(Det),ELISA
	Human PD1 Recombinant Rabbit mAb [PSH13-38] - BSA and Azide free (Detector)	HA723522	Human	ELISA(Det)
	IDH1 Mouse mAb [A11-B4]	EM40705	Human, Mouse, Rat, Monkey, Zebrafish	WB,IHC-P
KD	IDH1 Recombinant Rabbit mAb [JE34-27]	HA722696	Human, Mouse, Rat	WB,IHC-P,IP
KD	IDH2 Recombinant Rabbit mAb [JA55-31]	ET1704-93	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Ikaros Recombinant Rabbit mAb [JB50-38]	ET7107-25	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	IKZF3 Recombinant Rabbit mAb [JE62-81]	HA720110	Human	WB
	ILT-3 Rabbit pAb	ER63760	Human	IF-Cell
	Integrin alpha 6 Recombinant Rabbit mAb [SR45-00] - BSA and Azide free	HA750028	Human, Mouse, Rat	WB,IHC-P
	Integrin alpha 6 Recombinant Rabbit mAb [SR45-00]	ET1602-1	Human, Mouse, Rat	WB,IHC-P
	KAT6A Recombinant Rabbit mAb [PSH09-88] - BSA and Azide free	HA751322	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,ChIP
	KAT6A Recombinant Rabbit mAb [PSH09-88]	HA723157	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,ChIP
	KDM1 / LSD1 Recombinant Rabbit mAb [JE04-60]	HA722723	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	K-Ras Rabbit pAb	ER65706	Human, Mouse, Rat	WB,IHC-P
	KRAS Recombinant Mouse mAb [A8E6]	HA601060	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	KRAS Mouse mAb [A8E5]	HA601059	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	KRAS Recombinant Mouse mAb [A8E6] - BSA and Azide free	HA610008	Human, Mouse, Rat	WB,IHC-P
	KRAS Mouse mAb [A8E4]	HA601058	Human, Mouse, Rat	WB,IF-Cell,FC
	LSD1 Rabbit pAb	ER1802-12	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	MCL1 Recombinant Rabbit mAb [SI16-04]	ET1606-14	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
	MCL1 Rabbit pAb	R1510-26	Human, Mouse	IF-Cell,IHC-P
KD	MCL1 Recombinant Rabbit mAb [SI16-04] - BSA and Azide free	HA750093	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
KD	Menin Recombinant Rabbit mAb [JE30-97]	HA721313	Human	WB,IHC-P
	Menin Rabbit pAb	ER63947	Human, Mouse, Rat	WB,IHC-P
	MLL Rabbit pAb	ER65766	Human, Mouse, Rat	IHC-P
	MN1 Rabbit pAb	ER63982	Human	WB
	PE Conjugated Mouse CD16/32 Recombinant Rabbit mAb [PSH06-60]	HA720216F	Mouse	FC
	FITC Conjugated Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA720199F	Mouse	FC
	PE Conjugated Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA720200F	Mouse	FC
	APC Conjugated Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA720201F	Mouse	FC
	Anti-Mouse CD3 Recombinant Antibody [PSH04-98] - Rabbit IgG (Chimeric)	HA723347	Mouse	FC
	FITC Conjugated Anti-Mouse CD4 Recombinant Antibody [PSH04-99] - Rabbit IgG (Chimeric)	HA720202F	Mouse	FC
	APC Conjugated Anti-Mouse CD4 Recombinant Antibody [PSH04-99] - Rabbit IgG (Chimeric)	HA720204F	Mouse	FC
	PE Conjugated Anti-Mouse CD4 Recombinant Antibody [PSH04-99] - Rabbit IgG (Chimeric)	HA720203F	Mouse	FC
	Anti-Mouse CD8 alpha Recombinant Antibody [PSH05-00] - Rabbit IgG (Chimeric) - BSA and Azide free	HA751593	Mouse	FC,IHC-Fr
	PE Conjugated Anti-Mouse CD8 alpha Recombinant Antibody [PSH05-00] - Rabbit IgG (Chimeric)	HA720206F	Mouse	FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse CD8 alpha Recombinant Rat mAb [PSH16-78] - BSA and Azide free	HA610359	Mouse	FC
	FITC Conjugated Anti-Mouse CD8 alpha Recombinant Antibody [PSH05-00] - Rabbit IgG (Chimeric)	HA720205F	Mouse	FC
	Anti-Mouse CD8 alpha Recombinant Antibody [PSH05-00] - Rabbit IgG (Chimeric)	HA723820	Mouse	FC,IHC-Fr
	APC Conjugated Anti-Mouse CD8 alpha Recombinant Antibody [PSH05-00] - Rabbit IgG (Chimeric)	HA720207F	Mouse	FC
	Mouse CD8 alpha Recombinant Rat mAb [PSH16-78]	HA601543	Mouse	FC
	Mouse PD1 Recombinant Rabbit mAb [PSH04-18] - BSA and Azide free (Detector)	HA722101	Mouse	ELISA(Det)
	Mouse PD1 Recombinant Rabbit mAb [PSH04-17] - BSA and Azide free (Capture)	HA722100	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse PD1 Recombinant Rabbit mAb [PSH04-18]	HA722613B	Mouse	ELISA(Det)
	mTOR Rabbit pAb	R1510-21	Human, Mouse, Rat	IF-Cell,IHC-P
	mTOR Rabbit pAb	HA500126	Human, Rat, Mouse	WB,IF-Cell,IHC-P
	mTOR Recombinant Rabbit mAb [SU30-00] - BSA and Azide free	HA750153	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell
	mTOR Recombinant Rabbit mAb [SU30-00]	ET1608-5	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell
	N/H/K-Ras Rabbit pAb	ER64095	Human, Mouse, Rat	WB,IHC-P
	Nucleophosmin Rabbit pAb	ER1914-19	Human, Mouse	WB,IHC-P,IF-Tissue
KD	Nucleophosmin Recombinant Rabbit mAb [PD00-91]	HA721206	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	NUP98 Recombinant Rabbit mAb [JE04-31]	HA722523	Human, Mouse	WB,IF-Cell,IHC-P
	PD1 Mouse mAb	IRS190MS	Human	mIHC
	PD1 Recombinant Rabbit mAb [SJ01-91]	ET1606-41	Human	WB,IHC-P
	PD1 Recombinant Rabbit mAb [SJ01-91] - BSA and Azide free	HA750103	Human	WB,IHC-P
Cell treatment	PD1 Mouse mAb [1F2]	EM1707-60	Human	WB,IF-Cell,IHC-P,FC
	PD1 Mouse mAb	IRS191MS	Mouse	mIHC
	PD1 Blocking Mouse mAb [3A11]	EM1701-17	Human	ELISA,FC
	PHAPI2 / APRIL Recombinant Rabbit mAb [JE35-49]	HA722541	Human	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-AKT1 (S124) Recombinant Rabbit mAb [JJ08-46]	ET1701-36	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-AKT1 (S129) Recombinant Rabbit mAb [JE59-50]	HA721677	Human, Mouse, Rat	WB
Cell treatment	Phospho-AKT1 (T450) Recombinant Rabbit mAb [SD08-12]	ET1612-73	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44] - BSA and Azide free	HA750953	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05] - BSA and Azide free	HA750132	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44]	HA722129	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05]	ET1607-73	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10]	HA722951	Human	WB,IF-Cell,FC
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10] - BSA and Azide free	HA751205	Human	WB,IF-Cell,FC
Cell treatment	Phospho-Bcl-2 (S70) Recombinant Rabbit mAb [PSH11-02] - BSA and Azide free	HA751392	Human	WB,IF-Cell,FC
Cell treatment	Phospho-Bcl-2 (S70) Recombinant Rabbit mAb [PSH11-02]	HA723269	Human	WB,IF-Cell,FC
Cell treatment	Phospho-BTK (Y223) Recombinant Rabbit mAb [ST51-05]	ET1609-55	Human, Mouse	WB,IP
Cell treatment	Phospho-mTOR (S2448) Mouse mAb [A5D5]	HA600094	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-mTOR (S2481) Recombinant Rabbit mAb [JE59-62]	HA721632	Human	WB
	Phospho-Nucleophosmin (S125) Recombinant Rabbit mAb [JE51-14]	HA722559	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-Nucleophosmin (T199) Recombinant Rabbit mAb [JE63-16]	HA722125	Human	WB
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-38] - BSA and Azide free	HA750716	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-39]	HA721673	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-39] - BSA and Azide free	HA750717	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-38]	HA721672	Human, Mouse	WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH15-91] - BSA and Azide free	HA751582	Human, Mouse, Rat	WB,IF-Cell,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH15-91]	HA723796	Human, Mouse, Rat	WB,IF-Cell,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH16-27]	HA723825	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH16-27] - BSA and Azide free	HA751595	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,Dot Blot
Cell treatment	Phospho-PKC (S729) Recombinant Rabbit mAb [JE77-54]	HA722846	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-SHP2 (Y542) Recombinant Rabbit mAb [SN61-01]	ET1611-22	Mouse, Human	WB,IP
	PI3 Kinase p110 beta Recombinant Rabbit mAb [JE65-38]	HA722474	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	PI 3 Kinase p110 delta Recombinant Rabbit mAb [JM83-24]	ET1705-46	Human, Mouse, Rat	WB,IHC-P
KD	PI3-kinase p110 subunit beta Recombinant Rabbit mAb [SC0580]	ET1610-36	Human	WB,IP,IHC-P,FC
	PI3 Kinase p110 α Recombinant Rabbit mAb [SJ0186]	ET1606-36	Human, Mouse	WB,IF-Cell,IF-Tissue,IP,IHC-P
KD	PI 3 Kinase p85 alpha Mouse mAb [A3-D0]	M1510-2	Human	WB,IF-Cell
	PI 3 Kinase p85 alpha Rabbit pAb	R1312-6	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	PI 3 Kinase p85 alpha Recombinant Rabbit mAb [SU04-07]	ET1608-70	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	PI 3 Kinase p85 alpha Mouse mAb [L3-D6]	EM1701-62	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
	PI 3 Kinase p85 alpha Recombinant Rabbit mAb [SU04-07] - BSA and Azide free	HA750162	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	PI 3 Kinase p85 alpha Recombinant Mouse mAb [A3-D0-R] - BSA and Azide free	HA610084	Human	WB,IF-Cell
KD	PI 3 Kinase p85 alpha Recombinant Mouse mAb [A3-D0-R]	HA601206	Human	WB,IF-Cell
	PI3 Kinase p85 alpha Recombinant Rabbit mAb [JE04-11]	HA722522	Human, Mouse, Rat	WB
	PI 3 Kinase p85 alpha/gamma Rabbit pAb	ER64588	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	PITSLRE Rabbit pAb	ER64610	Human, Mouse	WB,IHC-P
	PKC Rabbit pAb	ER64620	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	PKC beta Rabbit pAb	ER1802-14	Human, Mouse	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PKC beta 1 Recombinant Rabbit mAb [JM43-26]	ET1705-30	Human, Mouse, Rat	WB,IHC-P,FC
	PKC beta 2 Recombinant Rabbit mAb [ST48-06]	ET1609-44	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Ras Recombinant Rabbit mAb [JF10-11]	ET1702-94	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IP,FC
	Ras Recombinant Rabbit mAb [SA39-05]	ET1601-16	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC,IP
	Ras Recombinant Rabbit mAb [JE44-13]	HA721883	Human, Mouse, Rat	WB,IF-Cell,FC
	Ras Rabbit pAb	ER40115	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	ROBO3 Rabbit pAb	ER61894	Human, Mouse	WB
Relative expression	ROR1 Recombinant Mouse mAb [PSH13-99] - BSA and Azide free	HA610281	Human	WB,FC
Relative expression	ROR1 Recombinant Mouse mAb [PSH13-99]	HA601447	Human	WB,FC
	SHP2 Rabbit pAb	ER1706-93	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	SHP2 Recombinant Rabbit mAb [SN72-02]	ET1611-35	Human, Mouse, Rat	WB,IP,IHC-P
KD	SHP2 Recombinant Rabbit mAb [SN07-32]	ET1611-80	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Syndecan 1 Recombinant Rabbit mAb [JM11-21]	ET1703-42	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Syndecan 1 Recombinant mAb	IRS105RB	Human	mIHC
	Syndecan 1 Recombinant Rabbit mAb [JM11-21] - BSA and Azide free	HA750381	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	TCR beta Rabbit pAb	ER65125	Human	WB
	TLR8 Rabbit pAb	ER60086	Human, Mouse, Rat	IHC-P
	TRBC1 Recombinant Mouse mAb [PSH15-51] - BSA and Azide free	HA610351	Human	FC
	TRBC1 Recombinant Mouse mAb [PSH15-51]	HA601524	Human	FC
	TXNRD1 Recombinant Rabbit mAb [JA11-32]	ET1704-81	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
KD	TXNRD2 Recombinant Rabbit mAb [PSH02-26] - BSA and Azide free	HA750794	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	TXNRD2 Recombinant Rabbit mAb [PSH02-27]	HA721807	Human, Mouse, Rat	WB,IHC-P
KD	TXNRD2 Recombinant Rabbit mAb [PSH02-26]	HA721806	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	TXNRD2 Recombinant Rabbit mAb [PSH02-27] - BSA and Azide free	HA750795	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Wee1 Recombinant Rabbit mAb [PSH05-11]	HA722227	Human, Monkey	WB,IHC-P
	Wee1 Recombinant Rabbit mAb [PSH05-11] - BSA and Azide free	HA750984	Human, Monkey	WB,IHC-P
	WEE1 Rabbit pAb	ER1918-64	Human, Rat	WB,IHC-P,IF-Tissue