

2025年7月—新上市产品清单

高级验证	产品名	货号	反应种属	验证应用
	ACE Recombinant mAb	IRS086RB	Mouse	mIHC
Cell treatment	Acid sphingomyelinase Recombinant Rabbit mAb [PSH15-61]	HA723768	Human, Mouse, Rat, Monkey	WB,IHC-P,IP
	Albumin Recombinant Rabbit mAb	IRS087RB	Mouse	mIHC
	Alpha 1 antitrypsin Recombinant Rabbit mAb [PSH16-21]	HA723819	Mouse, Rat	WB,IHC-P
	alpha smooth muscle Actin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS093MS	Human	mIHC
	alpha smooth muscle Actin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS219MS	Mouse	mIHC
	alpha smooth muscle Actin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS220RT	Mouse	mIHC
	alpha smooth muscle Actin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS094RT	Human	mIHC
	AMID Recombinant Rabbit mAb [PSH15-60]	HA723767	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	Androgen receptor Recombinant mAb	IRS092RB	Mouse	mIHC
	AQP1 Recombinant Rabbit mAb	IRS218RB	Mouse	mIHC
	AQP1 Recombinant Rabbit mAb	IRS088RB	Human	mIHC
	Aquaporin 2 Recombinant mAb	IRS089RB	Mouse	mIHC
	Aquaporin 5 Recombinant mAb	IRS090RB	Mouse	mIHC
	BCKDHA Recombinant Rabbit mAb [PSH15-62]	HA723769	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	BCL6 Recombinant Mouse mAb	IRS097MS	Human	mIHC
Relative expression	Beta III Tubulin Recombinant Rabbit mAb [PSH15-32]	HA723738	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP
	Biotin Conjugated Brain-derived tau/BD-tau Recombinant Rabbit mAb [PSH16-47] - Detector	HA723849B	Human	ELISA(Det),ELISA
	Brain-derived tau/BD-tau Recombinant Rabbit mAb [PSH16-47] - BSA and Azide free (Detector)	HA723848	Human	ELISA(Det)
	Brain-derived tau/BD-tau Recombinant Rabbit mAb [PSH16-46] - BSA and Azide free (Capture)	HA723847	Human, Mouse, Rat	ELISA(Cap)
	BRD9 Recombinant Rabbit mAb [PSH15-95]	HA723800	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Calbindin Recombinant Rabbit mAb	IRS098RB	Human	mIHC
	Calbindin Recombinant Rabbit mAb	IRS221RB	Mouse	mIHC
	Cathepsin D Recombinant mAb	IRS099RB	Human	mIHC
	CD10 Mouse mAb	IRS101MS	Human	mIHC
	CD10 Mouse mAb	IRS222MS	Mouse	mIHC
	CD11c Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS102MS	Human	mIHC
	CD11c Recombinant Antibody - Rat IgG1 (Chimeric)	IRS103RT	Human	mIHC
	PE Conjugated CD123 / IL3RA Recombinant Mouse mAb [PSH14-51]	HA601528F	Human	FC
	APC Conjugated CD127 Recombinant Mouse mAb [PSH12-76]	HA601536F	Human	FC
	CD14 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS106RT	Human	mIHC
	CD15 Recombinant mAb	IRS107RB	Human	mIHC
	CD19 Recombinant mAb	IRS110RB	Mouse	mIHC
	FITC Conjugated CD19 Recombinant Mouse mAb [PSH14-97]	HA601518F	Human	FC
	FITC Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723762F	Human	FC
Relative expression	APC Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723761F	Human	FC
	PE Conjugated CD19 Recombinant Rabbit mAb [PSH13-44]	HA723763F	Human	FC
	CD1c Recombinant mAb	IRS111RB	Human	mIHC
	CD20 Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS112MS	Human	mIHC
	CD20 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS113RT	Human	mIHC
	APC Conjugated CD25 / IL-2 Receptor alpha Recombinant Mouse mAb [PSH14-46]	HA601530F	Human	FC
	APC Conjugated CD27 Recombinant Rabbit mAb [PSH08-90]	HA723828F	Human	FC
	CD31 Recombinant Antibody - Rabbit IgG (Chimeric)	IRS116RB	Human	mIHC
	CD38 Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS120MS	Human	mIHC
	CD38 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS119RT	Human	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD4 Recombinant Antibody [PSH08-22] - Rat IgG1 (Chimeric)	HA601534	Mouse	WB,IF-Tissue,IHC-P,IHC-Fr
Relative expression	CD4 Recombinant Antibody [PSH08-22] - Mouse IgG1 (Chimeric)	HA601533	Mouse	WB,IF-Tissue,IHC-P,IF-Cell
	CD44 Recombinant mAb	IRS122RB	Human	mIHC
	PE Conjugated CD45RO Recombinant Mouse mAb [PSH13-94]	HA601531F	Human	FC
Relative expression	CD62L Recombinant Rabbit mAb [PSH15-59]	HA723766	Human	FC
	CD66b pAb	IRS127RB	Human	mIHC
	APC Conjugated CD94 Recombinant Mouse mAb [PSH14-52]	HA601529F	Human	FC
	CD98 Mouse mAb	IRS130MS	Human	mIHC
	c-Fos Recombinant Rabbit mAb	IRS131RB	Mouse	mIHC
	CGRP Recombinant Antibody [PSH08-45] - Guinea pig IgG2 (Chimeric)	HA601540	Mouse	IHC-Fr
	Choline Acetyltransferase Recombinant Rabbit mAb [PSH15-49]	HA723755	Human, Mouse, Rat	IHC-P
	Chromogranin A Mouse mAb	IRS223MS	Mouse	mIHC
	Chromogranin A Mouse mAb	IRS132MS	Human	mIHC
	COL1A1 Recombinant mAb	IRS137RB	Human	mIHC
	COL1A1 Recombinant Rabbit mAb	IRS228RB	Mouse	mIHC
	Collagen XVII Recombinant mAb	IRS138RB	Mouse	mIHC
	C Reactive Protein Recombinant Rabbit Multiclonal Antibody [PSH15-35]	HA723739	Human, Mouse, Rat	WB
	CREB Recombinant Rabbit mAb [PSH15-48]	HA723754	Human, Mouse, Rat, Monkey	WB,IHC-P,IHC-Fr,IF-Cell,FC
Relative expression	CTGF Recombinant Rabbit mAb [PSH15-64]	HA723771	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	CYP2E1 Recombinant mAb	IRS139RB	Human	mIHC
	Cytokeratin 1 Recombinant Rabbit mAb	IRS163RB	Mouse	mIHC
	Cytokeratin 13 Recombinant mAb	IRS133RB	Human	mIHC
	Cytokeratin 14 Recombinant Rabbit mAb	IRS162RB	Human	mIHC
	Cytokeratin 18 Recombinant Rabbit mAb	IRS224RB	Mouse	mIHC
	Cytokeratin 19 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS134RT	Human	mIHC
	Cytokeratin 19 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS226RT	Mouse	mIHC
	Cytokeratin 19 Recombinant Rabbit mAb	IRS225RB	Mouse	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Cytokeratin 20 Recombinant mAb	IRS135RB	Human	mIHC
	Cytokeratin 5 Recombinant Rabbit mAb	IRS249RB	Human	mIHC
	DCBLD2 / ESDN Rabbit pAb	HA500585	Human	WB
	Desmin Recombinant mAb	IRS140RB	Mouse	mIHC
	E-Cadherin Recombinant Mouse mAb	IRS229MS	Mouse	mIHC
	Elastin Recombinant Rabbit mAb [PSH16-13]	HA723818	Human	IHC-P,IF-Tissue
	Endothelin 1 Recombinant Rabbit mAb [PSH15-63]	HA723770	Human, Mouse, Rat	WB
	Estriol Recombinant Rabbit mAb [PSH15-86]	HA723795	Species independent	ELISA
	Estrogen Receptor alpha Recombinant mAb	IRS142RB	Human	mIHC
	F4/80 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS143RT	Mouse	mIHC
Relative expression	F4/80 Recombinant Antibody [PSH01-87] - Mouse IgG1 (Chimeric)	HA601537	Mouse	WB,IHC-P
	F4/80 Recombinant Antibody [PSH01-87] - Rat IgG1 (Chimeric)	HA601538	Mouse	WB,IHC-Fr,IHC-P
	Factor XII Rabbit pAb	HA500586	Mouse	WB
	Fast Myosin Skeletal Heavy chain Recombinant mAb	IRS145RB	Mouse	mIHC
	Fibronectin Recombinant mAb	IRS144RB	Mouse	mIHC
	FOXA2 Recombinant mAb	IRS146RB	Mouse	mIHC
	GFAP Recombinant Antibody - Rat IgG1 (Chimeric)	IRS148RT	Mouse	mIHC
	Ghrelin Recombinant Rabbit mAb	IRS149RB	Human	mIHC
	Glucagon Recombinant Rabbit mAb	IRS230RB	Mouse	mIHC
	Glucagon Recombinant Rabbit mAb	IRS150RB	Human	mIHC
	Glucose Transporter GLUT1 Recombinant Rabbit mAb	IRS151RB	Human	mIHC
	Glutamine Synthetase Recombinant Mouse mAb	IRS231MS	Mouse	mIHC
	iFluor™ 555 Conjugated Goat anti-Chicken IgY Goat pAb	HA1141	Chicken	IHC-Fr,IF-Tissue,IF-Cell
	iFluor™ 647 Conjugated Goat anti-Chicken IgY Goat pAb	HA1142	Chicken	IHC-Fr,IF-Tissue,IF-Cell
	iFluor™ 594 Conjugated Goat anti-Guinea pig IgG Goat pAb	HA1143	Guinea pig	IHC-Fr,IF-Tissue,IF-Cell
Relative expression	Granzyme K Recombinant Rabbit mAb [PSH16-28]	HA723826	Human	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Heme Oxygenase 1 (HO-1) Recombinant Rabbit mAb	IRS232RB	Mouse	mIHC
	Heme Oxygenase 1 (HO-1) Recombinant Rabbit mAb	IRS156RB	Human	mIHC
	HLA-DR Recombinant Mouse mAb	IRS154MS	Human	mIHC
	HNF-4-alpha Recombinant Rabbit mAb	IRS155RB	Human	mIHC
	Hormone sensitive lipase / HSL Recombinant Rabbit mAb	IRS157RB	Human	mIHC
	HRP Conjugated Human Angiogenin Recombinant Rabbit mAb [PSH12-65] - Detector	HA725078H	Human	ELISA(Det),ELISA
	Human Angiotensinogen, C-His Tag Protein	HA211165	Human	
	HRP Conjugated Human Antithrombin III / ATIII Recombinant Rabbit mAb [PSH14-26] - Detector	HA723617H	Human	ELISA(Det),ELISA
	Biotin Conjugated Human BACE1 Recombinant Rabbit mAb [PSH16-23] - Detector	HA725172B	Human	ELISA(Det),ELISA
	Human BACE1 Recombinant Rabbit mAb [PSH16-23] - BSA and Azide free (Detector)	HA725172	Human	ELISA(Det)
	Human BACE1 Recombinant Rabbit mAb [PSH16-22] - BSA and Azide free (Capture)	HA725171	Human	ELISA(Cap)
	Human BST2/Tetherin, Tag Free Protein	HA211225	Human	
	Biotin Conjugated Human CCL11/Eotaxin Recombinant Rabbit mAb [PSH16-25] - Detector	HA725198B	Human	ELISA(Det),ELISA
	Human CCL11/Eotaxin Recombinant Rabbit mAb [PSH16-24] - BSA and Azide free (Capture)	HA725197	Human	ELISA(Cap)
	Human CCL11/Eotaxin Recombinant Rabbit mAb [PSH16-25] - BSA and Azide free (Detector)	HA725198	Human	ELISA(Det)
	Human CCL27, Tag Free Protein	HA211133	Human	
	Human CD115/M-CSF R/CSF-1-R enzyme-linked immunoassay kit (one step)	EHY0168	Human	ELISA
	PE Conjugated Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA601525F	Human	FC
	Human CD46, Tag Free Protein	HA211227	Human	
	HRP Conjugated Human CD62E Recombinant Rabbit mAb [PSH14-42] - Detector	HA723637H	Human	ELISA(Det),ELISA
	Human CD63 enzyme-linked immunoassay kit	EH0167	Human	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human CHIT1 Recombinant Rabbit mAb [PSH15-90] - Detector	HA725206B	Human	ELISA(Det),ELISA
	Human CHIT1 Recombinant Rabbit mAb [PSH15-90] - BSA and Azide free (Detector)	HA725206	Human	ELISA(Det)
	Human CHIT1 Recombinant Rabbit mAb [PSH15-89] - BSA and Azide free (Capture)	HA725205	Human	ELISA(Cap)
	Human Clusterin Recombinant Rabbit mAb [PSH15-55] - BSA and Azide free (Capture)	HA725091	Human	ELISA(Cap)
	Human Clusterin Recombinant Rabbit mAb [PSH15-56] - BSA and Azide free (Detector)	HA725092	Human	ELISA(Det)
	Biotin Conjugated Human Clusterin Recombinant Rabbit mAb [PSH15-56] - Detector	HA725093B	Human	ELISA(Det),ELISA
	HRP Conjugated Human Clusterin Recombinant Rabbit mAb [PSH15-56] - Detector	HA725093H	Human	ELISA(Det),ELISA
	HRP Conjugated Human Complement C3 Recombinant Rabbit mAb [PSH15-08] - Detector	HA725142H	Human	ELISA(Det),ELISA
	Human Complement Factor D/CFD, C-His Tag Protein	HA211128	Human	
	Human CX3CL1 Recombinant Rabbit mAb [PSH16-09] - BSA and Azide free (Detector)	HA723814	Human	ELISA(Det)
	Biotin Conjugated Human CX3CL1 Recombinant Rabbit mAb [PSH16-09] - Detector	HA723814B	Human	ELISA(Det),ELISA
	Human CX3CL1 Recombinant Rabbit mAb [PSH16-08] - BSA and Azide free (Capture)	HA723813	Human	ELISA(Cap)
	Human CXCL13/BCA-1 enzyme-linked immunoassay kit	EH0172	Human	ELISA
	Human DICR/CLEC4A, Tag Free Protein	HA211226	Human	
	Human Endorepellin/Perlecan, C-His Tag Protein	HA211127	Human	
	Human Endostatin Recombinant Rabbit mAb [PSH15-82] - BSA and Azide free (Capture)	HA723789	Human	ELISA(Cap)
	Human Endostatin Recombinant Rabbit mAb [PSH15-83] - BSA and Azide free (Detector)	HA723790	Human	ELISA(Det)
	Biotin Conjugated Human Endostatin Recombinant Rabbit mAb [PSH15-83] - Detector	HA723791B	Human	ELISA(Det),ELISA
	Human Fetuin A/AHSG, Tag Free Protein	HA211146	Human	

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Flt3 Ligand/FLT3LG, Tag Free Protein	HA211222	Human	
	Human GDF-15 Recombinant Rabbit mAb [PSH16-01] - BSA and Azide free (Detector)	HA723806	Human	ELISA(Det)
	Human GDF-15 Recombinant Rabbit mAb [PSH16-00] - BSA and Azide free (Capture)	HA723805	Human	ELISA(Cap)
	Human GDF-15 Recombinant Rabbit mAb [PSH16-02] - BSA and Azide free (Detector)	HA723807	Human	ELISA(Det)
	Biotin Conjugated Human GDF-15 Recombinant Rabbit mAb [PSH16-02] - Detector	HA723807B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human GDF-15 Recombinant Rabbit mAb [PSH16-01] - Detector	HA723806B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human Granzyme B Recombinant Rabbit mAb [PSH15-99] - Detector	HA723804B	Human	ELISA(Det),ELISA
	Human Granzyme B Recombinant Rabbit mAb [PSH15-98] - BSA and Azide free (Capture)	HA723803	Human	ELISA(Cap)
	Human Granzyme B Recombinant Rabbit mAb [PSH15-99] - BSA and Azide free (Detector)	HA723804	Human	ELISA(Det)
	Human ICAM-3/CD50, C-His Tag (ECD) Protein	HA211145	Human	
	FITC Conjugated Human IgD Recombinant Mouse mAb [PSH14-48]	HA601526F	Human	FC
	PE Conjugated Human IgE Recombinant Mouse mAb [PSH14-50]	HA601527F	Human	FC
	Human IL-2 Receptor alpha enzyme-linked immunoassay kit	EH0170	Human	ELISA
	Human IL-6R Recombinant Rabbit mAb [PSH15-38] - BSA and Azide free (Capture)	HA723743	Human	ELISA(Cap)
	Biotin Conjugated Human IL-6R Recombinant Rabbit mAb [PSH15-39] - Detector	HA723745B	Human	ELISA(Det),ELISA
	Human IL-6R Recombinant Rabbit mAb [PSH15-39] - BSA and Azide free (Detector)	HA723744	Human	ELISA(Det)
	HRP Conjugated Human IL-6R Recombinant Rabbit mAb [PSH15-39] - Detector	HA723745H	Human	ELISA(Det),ELISA
	Human Lipocalin-2/NGAL enzyme-linked immunoassay kit (one step)	EHY0166	Human	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HRP Conjugated Human Lipoprotein A Recombinant Rabbit mAb [PSH13-83] - Detector	HA723574H	Human	ELISA(Det),ELISA
	Human Lipoprotein A/LPA enzyme-linked immunoassay kit	EH0171	Human	ELISA
	Human MCP-2/CCL8, Tag Free Protein	HA211228	Human	
	Human MCP4/CCL13 Recombinant Rabbit mAb [PSH15-70] - BSA and Azide free (Detector)	HA725200	Human	ELISA(Det)
	Human MCP4/CCL13 Recombinant Rabbit mAb [PSH15-69] - BSA and Azide free (Capture)	HA725199	Human	ELISA(Cap)
	Biotin Conjugated Human MCP4/CCL13 Recombinant Rabbit mAb [PSH15-70] - Detector	HA725200B	Human	ELISA(Det),ELISA
	Human MMP-3 enzyme-linked immunoassay kit	EH0169	Human	ELISA
	Human Oncostatin M/OSM, C-His Tag Protein	HA211164	Human	
	Human PARK7 Recombinant Mouse mAb [PSH15-67] - BSA and Azide free (Capture)	HA620001	Human	ELISA(Cap)
	Biotin Conjugated Human PARK7 Recombinant Mouse mAb [PSH15-68] - Detector	HA620002B	Human	ELISA(Det),ELISA
	Human PARK7 Recombinant Mouse mAb [PSH15-68] - BSA and Azide free (Detector)	HA620002	Human	ELISA(Det)
	HRP Conjugated Human PCSK9 Recombinant Rabbit mAb [PSH12-86] - Detector	HA725126H	Human	ELISA(Det),ELISA
	Biotin Conjugated Human PD-L1 Recombinant Rabbit mAb [PSH15-37] - Detector	HA723742B	Human	ELISA(Det),ELISA
	Human PD-L1 Recombinant Rabbit mAb [PSH15-37] - BSA and Azide free (Detector)	HA723741	Human	ELISA(Det)
	Human PD-L1 Recombinant Rabbit mAb [PSH15-36] - BSA and Azide free (Capture)	HA723740	Human	ELISA(Cap)
	HRP Conjugated Human RBP4 Recombinant Rabbit mAb [PSH14-89] - Detector	HA723690H	Human	ELISA(Det),ELISA
	Human ROR1, Tag Free Protein	HA211231	Human	
	Human SCF/KITLG, Tag Free Protein	HA211223	Human	
	Human Siglec 9 enzyme-linked immunoassay kit	EH0161	Human	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Tau/Tau758, C-His (isoform PNS-tau) Protein	HA211246	Human	
	HRP Conjugated Human TIM 4 Recombinant Rabbit mAb [PSH14-01] - Detector	HA723591H	Human	ELISA(Det),ELISA
	Human VIPR2, C-hFc tag Protein	HA211224	Human	
	iFluor™ 594 Conjugated Iba1 Recombinant Rabbit mAb [JM36-62]	HA723821F	Mouse	IHC-Fr
	Iba1 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS158RT	Mouse	mIHC
	Id1 Recombinant Rabbit mAb	IRS159RB	Human	mIHC
	IL-12A Recombinant Rabbit mAb [PSH16-10] - BSA and Azide free	HA723815	Human	ELISA
	IL-1 alpha Recombinant Rabbit mAb [PSH15-66]	HA723773	Mouse	IF-Cell,FC
Cell treatment	IL-1 alpha Recombinant Rabbit mAb [PSH15-65]	HA723772	Mouse	WB
	IL3RA / CD123 Recombinant mAb	IRS104RB	Human	mIHC
	Insulin Recombinant Rabbit mAb	IRS233RB	Mouse	mIHC
	Insulin Recombinant Rabbit mAb	IRS161RB	Human	mIHC
	KLF5 Recombinant Rabbit mAb [PSH15-42]	HA723748	Human, Mouse, Rat	IF-Cell,FC
	Laminin beta 1 Recombinant Rabbit mAb	IRS164RB	Human	mIHC
	Laminin beta 1 Recombinant Rabbit mAb	IRS234RB	Mouse	mIHC
Cell treatment	LC3B Recombinant Antibody [JJ090-6] - Mouse IgG1 (Chimeric)	HA601532	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	LIPT2 Recombinant Rabbit mAb [PSH15-52]	HA723758	Human	WB,IHC-P
	Loricrin Recombinant Rabbit mAb	IRS165RB	Mouse	mIHC
	LOX Recombinant Rabbit mAb [PSH15-41]	HA723747	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Luteinizing Hormone beta Recombinant Rabbit mAb [PD00-78]	HA723822	Human	IHC-P
	LY6A Recombinant Rabbit mAb	IRS166RB	Mouse	mIHC
	Lysozyme Recombinant Rabbit mAb	IRS168RB	Mouse	mIHC
	MAGED2 Recombinant Rabbit mAb [PSH15-53]	HA723759	Human	WB
	Mast Cell Tryptase Recombinant Antibody - Rat IgG1 (Chimeric)	IRS209RT	Human	mIHC
	Mast Cell Tryptase Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS210MS	Human	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	MCM2 Recombinant Antibody [PD00-90] - Mouse IgG1 (Chimeric)	HA601523	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P
	MHC class I Recombinant Rabbit mAb	IRS171RB	Human	mIHC
	MMP-9 Recombinant Rabbit mAb	IRS172RB	Human	mIHC
	Monkey Albumin/ALB, C-His Tag Protein	HA211229	Monkey	
	Mouse ACE, Tag Free Protein	HA211125	Mouse	
	Mouse Alpha 1 antitrypsin Recombinant Rabbit mAb [PSH16-21] - BSA and Azide free (Detector)	HA725147	Mouse	ELISA(Det)
	Mouse Alpha 1 antitrypsin Recombinant Rabbit mAb [PSH16-20] - BSA and Azide free (Capture)	HA725146	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse Alpha 1 antitrypsin Recombinant Rabbit mAb [PSH16-21] - Detector	HA725148B	Mouse	ELISA(Det),ELISA
	HRP Conjugated Mouse Alpha 1 antitrypsin Recombinant Rabbit mAb [PSH16-21] - Detector	HA725148H	Mouse	ELISA(Det),ELISA
	Mouse APOC3, Tag Free Protein	HA211230	Mouse	
	Anti-Mouse CD8 alpha Recombinant Antibody [PSH05-00] - Rabbit IgG (Chimeric)	HA723820	Mouse	FC,IHC-Fr
	HRP Conjugated Mouse CRP Recombinant Rabbit mAb [PSH12-49] - Detector	HA725075H	Mouse	ELISA(Det),ELISA
	Mouse Fas/CD95/TNFRSF6, Tag Free (ECD) Protein	HA211151	Mouse	
	Mouse G-CSF Recombinant Rabbit mAb [PSH15-74] - BSA and Azide free (Capture)	HA723777	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse G-CSF Recombinant Rabbit mAb [PSH15-75] - Detector	HA723779B	Mouse	ELISA(Det),ELISA
	Mouse G-CSF Recombinant Rabbit mAb [PSH15-75] - BSA and Azide free (Detector)	HA723778	Mouse	ELISA(Det)
	Mouse IL-20, Tag Free Protein	HA211139	Mouse	
	Mouse L-selectin/CD62L, C-His Tag Protein	HA211126	Mouse	
	Mouse M-CSF/CSF1 Recombinant Rabbit mAb [PSH15-84] - BSA and Azide free (Capture)	HA723792	Mouse	ELISA(Cap)
	Mouse M-CSF/CSF1 Recombinant Rabbit mAb [PSH15-85] - BSA and Azide free (Detector)	HA723793	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse M-CSF/CSF1 Recombinant Rabbit mAb [PSH15-85] - Detector	HA723794B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse MDC/CCL22 Recombinant Rabbit mAb [PSH15-77] - Detector	HA723782B	Mouse	ELISA(Det),ELISA
	Mouse MDC/CCL22 Recombinant Rabbit mAb [PSH15-76] - BSA and Azide free (Capture)	HA723780	Mouse	ELISA(Cap)
	Mouse MDC/CCL22 Recombinant Rabbit mAb [PSH15-78] - BSA and Azide free (Capture)	HA723783	Mouse	ELISA(Cap)
	Mouse MDC/CCL22 Recombinant Rabbit mAb [PSH15-77] - BSA and Azide free (Detector)	HA723781	Mouse	ELISA(Det)
	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-79] - BSA and Azide free (Capture)	HA723784	Mouse	ELISA(Cap)
	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - BSA and Azide free (Detector)	HA723787	Mouse	ELISA(Det)
	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-80] - BSA and Azide free (Capture/Detector)	HA723785	Mouse	ELISA(Det),ELISA(Cap)
	Biotin Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - Detector	HA723788B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-80] - Detector	HA723786B	Mouse	ELISA(Det),ELISA
	Mouse Pro-collagen I alpha 1/COL1A1 enzyme-linked immunoassay kit (one step)	EMY0045	Mouse	ELISA
	Mouse P-Selectin/CD62P, Tag Free (ECD) Protein	HA211140	Mouse	
	Mouse RANKL/TRANCE/TNFSF11 enzyme-linked immunoassay kit	EM0044	Mouse	ELISA
	Mouse/Rat IGF1 enzyme-linked immunoassay kit (one step)	EMY0024	Mouse, Rat	ELISA
	Mouse/Rat PDGF-BB Recombinant Rabbit mAb [PSH16-07] - BSA and Azide free (Capture)	HA723812	Mouse, Rat	ELISA(Cap)
	Mouse/Rat PDGF-BB Recombinant Rabbit mAb [PSH16-06] - BSA and Azide free (Detector)	HA723811	Mouse, Rat	ELISA(Det)
	Mouse/Rat PDGF-BB Recombinant Rabbit mAb [PSH16-05] - BSA and Azide free (Capture)	HA723810	Mouse, Rat	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse/Rat PDGF-BB Recombinant Rabbit mAb [PSH16-06] - Detector	HA723811B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat Prealbumin Recombinant Rabbit mAb [PSH15-88] - BSA and Azide free (Detector)	HA725129	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse/Rat Prealbumin Recombinant Rabbit mAb [PSH15-88] - Detector	HA725129B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat Prealbumin Recombinant Rabbit mAb [PSH15-87] - BSA and Azide free (Capture)	HA725128	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse Resistin Recombinant Rabbit mAb [PSH16-04] - Detector	HA723809B	Mouse	ELISA(Det),ELISA
	Mouse Resistin Recombinant Rabbit mAb [PSH16-03] - BSA and Azide free (Capture)	HA723808	Mouse	ELISA(Cap)
	Mouse Resistin Recombinant Rabbit mAb [PSH16-04] - BSA and Azide free (Detector)	HA723809	Mouse	ELISA(Det)
	Mouse Von Willebrand Factor/vWF Protein	HA211245	Mouse	
	MUC1 Mouse mAb	IRS174MS	Human	mIHC
	MUC2 Recombinant Rabbit mAb	IRS175RB	Human	mIHC
	MUC2 Recombinant Rabbit mAb	IRS235RB	Mouse	mIHC
	Myeloperoxidase Recombinant Mouse mAb	IRS173MS	Human	mIHC
	Myogenin Rabbit pAb	IRS177RB	Mouse	mIHC
	NeuN Recombinant Antibody - Rat IgG1 (Chimeric)	IRS179RT	Mouse	mIHC
	NeuN Recombinant Rabbit mAb	IRS178RB	Mouse	mIHC
	non-muscle Myosin IIB / MYH10 Recombinant Rabbit mAb [PSH15-72]	HA723775	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Non-Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-44] - BSA and Azide free (Capture)	HA723844	Human	ELISA(Cap)
	Non-Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-45] - BSA and Azide free (Detector)	HA723845	Human	ELISA(Det)
	Biotin Conjugated Non-Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-45] - Detector	HA723846B	Human	ELISA(Det),ELISA
	Non-Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-43] - BSA and Azide free (Capture)	HA723843	Human, Mouse, Rat	ELISA(Cap)
	NPHS2 Recombinant Rabbit mAb	IRS237RB	Mouse	mIHC
	NPHS2 Recombinant Rabbit mAb	IRS182RB	Human	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Occludin Recombinant Rabbit mAb	IRS183RB	Mouse	mIHC
	p21 Recombinant Rabbit mAb	IRS260RB	Mouse	mIHC
	p38 delta / MAPK13 Rabbit pAb	HA500582	Human	WB
	Pan-Cadherin Recombinant Rabbit mAb	IRS185RB	Mouse	mIHC
	Pancreatic Polypeptide Recombinant Rabbit mAb	IRS189RB	Mouse	mIHC
	pan Cytokeratin Recombinant Antibody - Rabbit IgG (Chimeric)	IRS186RB	Human	mIHC
	pan Cytokeratin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS187RT	Human	mIHC
	Parathyroid Hormone Recombinant Rabbit mAb [PD00-92]	HA723823	Human	IHC-P
Relative expression	Perforin-1 Recombinant Rabbit mAb [PSH15-92]	HA723797	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH16-27]	HA723825	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH15-91]	HA723796	Human, Mouse, Rat	WB,IF-Cell,Dot Blot
	Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-41] - BSA and Azide free (Capture)	HA723841	Human, Mouse	ELISA(Cap)
	Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-42] - BSA and Azide free (Capture)	HA723842	Human	ELISA(Cap)
	Podoplanin Recombinant Rabbit mAb	IRS194RB	Human	mIHC
	RAGE Recombinant Rabbit mAb	IRS195RB	Human	mIHC
	Rat CXCL2/MIP-2, Tag Free Protein	HA211132	Rat	
	Rat GM-CSF/CSF2, Tag Free Protein	HA211147	Rat	
	Rat MMP-9, Tag Free Protein	HA211141	Rat	
	Rat TIMP1, Tag Free Protein	HA211163	Rat	
	RUNX2 Recombinant Rabbit mAb	IRS238RB	Mouse	mIHC
Relative expression	Ryanodine receptor 2 / RYR-2 Recombinant Rabbit mAb [PSH15-54]	HA723760	Mouse, Rat	WB,IHC-P
	SCD1 Recombinant Rabbit mAb [PSH15-47]	HA723753	Mouse, Rat	WB,IHC-Fr,IHC-P,FC
	SCNN1G Recombinant Rabbit mAb [PSH16-29]	HA723827	Human, Mouse, Rat, Monkey	WB,IF-Cell
	SIRT5 Recombinant Rabbit mAb [PSH15-71]	HA723774	Human, Mouse, Rat	WB,IHC-P
	SLC12A1 / NKCC2 Recombinant Rabbit mAb	IRS236RB	Mouse	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	SLC12A1 / NKCC2 Recombinant Rabbit mAb	IRS181RB	Human	mIHC
	SLC31A1 / CTR1 Recombinant Rabbit mAb [PSH09-67]	HA723132	Human, Mouse, Rat	WB
	Slow Skeletal Myosin Heavy chain Recombinant Rabbit mAb	IRS201RB	Mouse	mIHC
	smooth muscle Myosin heavy chain 11 Recombinant Rabbit mAb	IRS176RB	Human	mIHC
	SOX2 Recombinant Rabbit mAb	IRS240RB	Mouse	mIHC
	SOX9 Recombinant Rabbit mAb	IRS199RB	Human	mIHC
	SOX9 Recombinant Rabbit mAb	IRS241RB	Mouse	mIHC
	SST/Somatostatin 28 Recombinant Rabbit mAb	IRS197RB	Human	mIHC
	SST/Somatostatin 28 Recombinant Rabbit mAb	IRS239RB	Mouse	mIHC
	Syndecan 1 Recombinant mAb	IRS105RB	Human	mIHC
	TAGLN Recombinant Rabbit mAb	IRS202RB	Human	mIHC
	Biotin Conjugated Tau Recombinant Rabbit mAb [PSH10-42] - Detector	HA723209B	Human, Mouse, Rat	ELISA(Det),ELISA
	Biotin Conjugated Tau Recombinant Rabbit mAb [PSH16-49] - Detector	HA723852B	Human	ELISA(Det),ELISA
	Tau Recombinant Rabbit mAb [PSH16-49] - BSA and Azide free (Detector)	HA723851	Human	ELISA(Det)
	Tau Recombinant Rabbit mAb [PSH16-48] - BSA and Azide free (Capture)	HA723850	Human, Mouse, Rat	ELISA(Cap)
	Tau Recombinant Rabbit mAb [PSH10-42] - BSA and Azide free (Detector)	HA723209	Human, Mouse, Rat	ELISA(Det)
	TAZ / WWTR1 Recombinant Rabbit mAb [PSH15-58]	HA723765	Human, Mouse, Rat	WB,IHC-P,IP
	TAZ / WWTR1 Recombinant Rabbit mAb [PSH15-57]	HA723764	Human, Mouse, Rat, Monkey	WB,IP
KO	TDG Recombinant Rabbit mAb [PSH10-11]	HA723180	Human, Mouse, Rat	WB,IHC-P
Relative expression	TFPI2 Recombinant Rabbit mAb [PSH15-45]	HA723751	Human	WB,IHC-P,IP
	TRBC1 Recombinant Mouse mAb [PSH15-51]	HA601524	Human	FC
	TRP1 Recombinant Rabbit mAb	IRS211RB	Human	mIHC
Relative expression	UCP1 Recombinant Rabbit mAb [PSH15-73]	HA723776	Mouse, Rat	WB,IHC-Fr,IHC-P,IP
	USP15 Recombinant Rabbit mAb [PSH15-94]	HA723799	Human, Mouse, Rat, Monkey	WB,IHC-P,IP
	Villin1 Recombinant Rabbit mAb	IRS213RB	Human	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Villin1 Recombinant Rabbit mAb	IRS243RB	Mouse	mIHC
	Vimentin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS214RT	Human	mIHC
	Von Willebrand Factor Recombinant Rabbit mAb	IRS244RB	Mouse	mIHC
	Von Willebrand Factor Recombinant Rabbit mAb	IRS215RB	Human	mIHC
	WDR47 Rabbit pAb	HA500584	Mouse, Rat	WB
	WWOX Recombinant Rabbit mAb [PSH15-44]	HA723750	Human, Mouse, Rat, Monkey	WB,IP
	WWOX Recombinant Rabbit mAb [PSH15-43]	HA723749	Human, Mouse, Rat, Monkey	WB,IF-Cell,IP
Relative expression	YKL-40 / CHI3L1 Recombinant Rabbit mAb [PSH16-26]	HA723824	Human	WB,IF-Cell
	ZEB Recombinant Rabbit mAb	IRS216RB	Human	mIHC
	ZEB Recombinant Rabbit mAb	IRS245RB	Mouse	mIHC
	ZO1 Recombinant Rabbit mAb	IRS217RB	Mouse	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层