

神经营养研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	Amyloid Beta A4 Precursor (APP) Rabbit pAb	1007-5	Human, Mouse, Rat	WB,IF-Cell
Cell treatment, Relative expression	APP/beta Amyloid Recombinant Rabbit mAb [PSH04-90]	HA722207	Human, Mouse	WB,IF-Cell,IHC-P
Relative expression	APP/beta Amyloid Recombinant Rabbit mAb [JE30-54]	HA722139	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
Cell treatment, Relative expression	APP/ β -Amyloid Recombinant Rabbit mAb [PSH04-90] - BSA and Azide free	HA750977	Human, Mouse	WB,IF-Cell,IHC-P
	AP repeat Rabbit pAb	HA500112	Human	WB
	Artn Rabbit pAb	ER65795	Human	IHC-P
	BDNF Recombinant Rabbit mAb [JE59-35]	HA722912	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-P,IHC-Fr,FC
KD	BDNF Recombinant Rabbit mAb [SJ12-09] - BSA and Azide free	HA750104	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC
KD	BDNF Recombinant Rabbit mAb [SJ12-09]	ET1606-42	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC
	BDNF Rabbit pAb	ER130915	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	beta Amyloid Recombinant Rabbit mAb [PSH08-19] - BSA and Azide free	HA751210	Human, Mouse	WB,IHC-P,IF-Tissue,Dot Blot
	beta Amyloid Recombinant Rabbit mAb	IRS251RB	Mouse	mIHC
	beta Amyloid Recombinant Rabbit mAb [PSH08-18] - BSA and Azide free	HA751209	Human, Mouse	WB,Dot Blot
Relative expression	beta Amyloid Recombinant Rabbit mAb [PSH08-19]	HA722962	Human, Mouse	WB,IHC-P,IF-Tissue,Dot Blot
	beta Amyloid Recombinant Rabbit mAb [PSH08-18]	HA722961	Human, Mouse	WB,Dot Blot
	Biotin Conjugated beta Amyloid 1-40 Recombinant Rabbit mAb [PSH03-98]	HA722618B	Human, Mouse	ELISA(Det)
Relative expression	beta Amyloid 1-40 Recombinant Rabbit mAb [PSH02-82]	HA721788	Human, Mouse, Cynomolgus monkey, Pig	IHC-P,WB
Relative expression	beta Amyloid 1-40 Recombinant Rabbit mAb [PSH05-86]	HA722485	Human, Mouse	WB,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	beta Amyloid 1-40 Recombinant Rabbit mAb [PSH02-82] - BSA and Azide free	HA750789	Human, Mouse, Cynomolgus monkey, Pig	IHC-P,WB
Relative expression	beta Amyloid 1-40 Recombinant Rabbit mAb [PSH05-86] - BSA and Azide free	HA751034	Human, Mouse	WB,IHC-P,IF-Tissue
	beta Amyloid 1-40 Recombinant Rabbit mAb [PSH03-98] - BSA and Azide free (Detector)	HA722078	Human, Mouse	ELISA(Det)
	beta Amyloid 1-40 Recombinant Rabbit mAb [PSH02-82] - BSA and Azide free (Capture)	HA722077	Human, Mouse	ELISA(Cap)
Relative expression	beta Amyloid 1-42 Recombinant Rabbit mAb [PSH06-81]	HA722754	Human, Mouse	WB,IHC-P,IF-Tissue
Relative expression	beta Amyloid 1-42 Recombinant Rabbit mAb [PSH02-83] - BSA and Azide free	HA750790	Human, Mouse, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue
	Biotin Conjugated beta Amyloid 1-42 Recombinant Rabbit mAb [PSH03-99]	HA722293	Human	ELISA(Det)
	beta Amyloid 1-42 Recombinant Rabbit mAb [PSH03-99] - BSA and Azide free (Detector)	HA722080	Human	ELISA(Det)
Relative expression	beta Amyloid 1-42 Recombinant Rabbit mAb [PSH02-83]	HA721789	Human, Mouse, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue
	beta Amyloid 1-42 Recombinant Rabbit mAb	IRS247RB	Mouse	mIHC
	beta Amyloid 1-42 Recombinant Rabbit mAb [PSH02-83] - BSA and Azide free (Capture)	HA722079	Human	ELISA(Cap)
Relative expression	beta Amyloid 1-42 Recombinant Rabbit mAb [PSH06-81] - BSA and Azide free	HA751103	Human, Mouse	WB,IHC-P,IF-Tissue
	Biotin Conjugated Beta Nerve Growth Factor Recombinant Rabbit mAb [PSH07-81]	HA722905B	Human	ELISA(Det)
	bFGF Rabbit pAb	ER1904-22	Human, Mouse	WB,FC
	CD271 Rabbit pAb	1003-9	Mouse, Rat	WB,IF-Cell,IHC-P
	CDNF Rabbit pAb	ER65797	Human	IHC-P
	CNTF Rabbit pAb	ER60106	Human, Rat, Mouse	WB
	CT A Recombinant Rabbit mAb [JE48-23]	HA722459	Human, Monkey	WB,IF-Cell,FC
	EGF Recombinant Rabbit mAb [JA14-15]	ET1704-33	Human, Mouse, Rat	WB,IHC-P,IP
	FGF2 Recombinant Rabbit mAb [JM28-10]	ET1703-18	Human, Mouse	WB,IP,FC,IHC-P
	Flk-1/Flt-4 Rabbit pAb	ER63443	Human, Mouse, Rat	IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Flt-4 Rabbit pAb	ER65750	Human, Mouse, Rat	WB,IHC-P
	F-spondin Recombinant Rabbit mAb [PSH02-84]	HA721863	Human	WB,IF-Cell,FC
	F-spondin Recombinant Rabbit mAb [PSH02-84] - BSA and Azide free	HA750837	Human	WB,IF-Cell,FC
	GDNF Recombinant Rabbit mAb [JA93-10]	ET1704-46	Human, Mouse, Rat	WB,IF-Cell
	GDNF Receptor alpha 2 Rabbit pAb	ER1909-48	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	GFRA3 Rabbit pAb	ER61905	Human, Mouse	WB
	GFR alpha-1 Rabbit pAb	ER63527	Human, Mouse, Rat	IHC-P,WB
	GFR alpha-4 Rabbit pAb	ER63528	Human	WB
Relative expression	GLP-1 Recombinant Rabbit mAb [PSH19-16]	HA724050	Human, Mouse, Rat	IHC-P,IF-Tissue
	Glucagon Recombinant Rabbit mAb [JF0960]	ET1702-20	Human, Mouse, Rat	WB,IHC-P,mIHC,IHC-Fr,IF-Tissue
	Glucagon Recombinant Rabbit mAb [JF33-10]	ET1702-59	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Glucagon Mouse mAb	EM1701-34	Human, Mouse, Rat	IHC,IF-Cell
	Glucagon Recombinant Rabbit mAb	IRS150RB	Human	mIHC
KD	Glucose 6 phosphate isomerase Recombinant Rabbit mAb [JG33-86]	ET7108-01	Human	WB,IHC-P,FC,IF-Cell,IF-Tissue
	HRG1 Rabbit pAb	ER1910-81	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IF-Tissue
	HRG beta 1 Rabbit pAb	ER1910-80	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
	Human Beta Nerve Growth Factor Recombinant Rabbit mAb [PSH07-81] - BSA and Azide free (Detector)	HA722904	Human	ELISA(Det)
	Human Beta Nerve Growth Factor Recombinant Rabbit mAb [PSH07-80] - BSA and Azide free (Capture)	HA722903	Human	ELISA(Cap)
	Human Beta Nerve Growth Factor Recombinant Rabbit mAb [PSH07-79] - BSA and Azide free (Capture)	HA722902	Human	ELISA(Cap)
	Human CNTF Recombinant Rabbit mAb [PSH11-80] - BSA and Azide free (Detector)	HA723382	Human, Rat	ELISA(Det)
	Human CNTF Recombinant Rabbit mAb [PSH11-81] - BSA and Azide free (Capture)	HA723384	Human, Rat	ELISA(Cap)
	Biotin Conjugated Human CNTF Recombinant Rabbit mAb [PSH11-80] - Detector	HA723383B	Human, Rat	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human CNTF Recombinant Rabbit mAb [PSH11-82] - Detector	HA723386B	Human, Rat	ELISA(Det),ELISA
	Human CNTF Recombinant Rabbit mAb [PSH11-79] - BSA and Azide free (Capture)	HA723381	Human, Rat	ELISA(Cap)
	Human CNTF Recombinant Rabbit mAb [PSH11-82] - BSA and Azide free (Detector)	HA723385	Human, Rat	ELISA(Det)
	Biotin Conjugated Human EGF Recombinant Rabbit mAb [PSH02-20] - Detector	HA722638B	Human	ELISA(Det)
	Human EGF Recombinant Rabbit mAb [PSH02-20] - BSA and Azide free (Detector)	HA721783	Human	ELISA(Det)
	Human EGF Recombinant Rabbit mAb [PSH02-19] - BSA and Azide free (Capture)	HA721782	Human	ELISA(Cap)
	Human FGF-2 Recombinant Rabbit mAb [PSH02-08] - BSA and Azide free (Capture)	HA721780	Human	ELISA(Cap)
	Human FGF-2 Recombinant Rabbit mAb [PSH02-09] - BSA and Azide free (Detector)	HA721781	Human	ELISA(Det)
	Biotin Conjugated Human FGF-2 Recombinant Rabbit mAb [PSH02-09]	HA722639B	Human	ELISA(Det)
	Human NCAM1 Recombinant Rabbit mAb [PSH08-07] - BSA and Azide free (Capture)	HA722946	Human	ELISA(Cap)
	Human NCAM1 Recombinant Rabbit mAb [PSH08-08] - BSA and Azide free (Capture/Detector)	HA722947	Human	ELISA(Det),ELISA(Cap)
	Human NCAM1 Recombinant Rabbit mAb [PSH08-09] - BSA and Azide free (Detector)	HA722949	Human	ELISA(Det)
	Biotin Conjugated Human NCAM1 Recombinant Rabbit mAb [PSH08-08] - Detector	HA722948B	Human	ELISA(Det)
	Biotin Conjugated Human NCAM1 Recombinant Rabbit mAb [PSH08-09] - Detector	HA722950B	Human	ELISA(Det)
	PE Conjugated Human NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53]	HA720215F	Human	FC
	APC Conjugated Human NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53]	HA723757F	Human	FC
	Human VEGF Receptor 1 Recombinant Rabbit mAb [PSH19-22] - BSA and Azide free (Capture)	HA725365	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human VEGF Receptor 1 Recombinant Rabbit mAb [PSH19-23] - BSA and Azide free (Detector)	HA725366B	Human	ELISA(Det),ELISA
	Human VEGF Receptor 1 Recombinant Rabbit mAb [PSH19-23] - BSA and Azide free (Detector)	HA725366	Human	ELISA(Det)
	Biotin Conjugated Human VEGF Receptor 2 Recombinant Rabbit mAb [PSH12-84] - Detector	HA725123B	Human	ELISA(Det),ELISA
	Human VEGF Receptor 2 Recombinant Rabbit mAb [PSH12-83] - BSA and Azide free (Capture)	HA725121	Human	ELISA(Cap)
	Human VEGF Receptor 2 Recombinant Rabbit mAb [PSH12-84] - BSA and Azide free (Detector)	HA725122	Human	ELISA(Det)
Cell treatment	IGF1 Recombinant Rabbit mAb [JE30-01]	HA722527	Human	WB
	IGF1 (E peptide) Rabbit pAb	HA500107	Human	WB,IHC-P,FC
	IGF2 Recombinant Rabbit mAb [JJ092-3]	ET1701-88	Human	WB,IF-Cell,IF-Tissue,IHC-P
	IGF I Rabbit pAb	ER1911-10	Human, Mouse, Rat, Dog	WB,IHC-P,IF-Tissue
	KDR / VEGFR-2 Mouse mAb [1G1]	RT1347	Human, Mouse, Rat	WB,IP,IF,IHC-P
	Midkine Recombinant Rabbit mAb [JF096-5]	ET1702-64	Human	WB,IHC-P,IP
	HRP Conjugated Mouse/Rat IGF1 Recombinant Rabbit mAb [PSH12-43] - Detector	HA725039H	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat IGF1 Recombinant Rabbit mAb [PSH12-43] - BSA and Azide free (Detector)	HA725038	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse/Rat IGF1 Recombinant Rabbit mAb [PSH12-43] - Detector	HA725039B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat IGF1 Recombinant Rabbit mAb [PSH12-42] - BSA and Azide free (Capture)	HA725037	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat VEGFR3 / FLT4 Recombinant Rabbit mAb [PSH19-87] - Detector	HA724111B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat VEGFR3 / FLT4 Recombinant Rabbit mAb [PSH19-86] - BSA and Azide free (Capture)	HA724110	Mouse, Rat	ELISA(Cap)
	Mouse/Rat VEGFR3 / FLT4 Recombinant Rabbit mAb [PSH19-87] - BSA and Azide free (Detector)	HA724111	Mouse, Rat	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	NCAM1 / CD56 Recombinant Rabbit mAb	IRS026	Human	mIHC
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-82]	HA722755	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH08-89]	HA723041	Human	WB,IHC-P,IF-Tissue
	NCAM1 / CD56 Recombinant Rabbit mAb [JF1021] - BSA and Azide free	HA750346	Human, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53] - BSA and Azide free	HA751105	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-82] - BSA and Azide free	HA751104	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue
Relative expression	NCAM1 / CD56 Recombinant Mouse mAb [PDH0-07]	HA601126	Human	IHC-P,WB
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53]	HA722756	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	NCAM1 / CD56 Recombinant Rabbit mAb [JF1021]	ET1702-43	Human, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP,mIHC
	NCAM1 / CD56 Rabbit pAb	R1204-1	Human, Mouse, Rat	WB,IF-Cell
	Neuregulin-1 SMDF Rabbit pAb	ER64158	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Neurotrophin 4 Rabbit pAb	ER1913-76	Human, Rat	WB,IHC-P,IF-Tissue
	Neurturin Rabbit pAb	ER1913-77	Human, Mouse	WB,IHC-P,IF-Tissue
	NGF Rabbit pAb	ER1802-59	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	NGF Recombinant Rabbit mAb [SI79-01]	ET1606-29	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
	NGF Recombinant Rabbit mAb [SI79-01] - BSA and Azide free	HA750099	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
	NGF beta Mouse mAb [1G4]	RT1428	Human, Mouse, Rat	WB,IP,IF
	NGF Receptr Mouse mAb [2G2]	RT1429	Human, Mouse, Rat	WB,IP,IF
	NT-3 Rabbit pAb	ER65762	Human, Mouse, Rat	IHC-P
	p75 NGF Receptor Recombinant Rat mAb [PSH13-26]	HA601436	Mouse, Rat	IHC-Fr,WB
	p75 NGF Receptor Recombinant Antibody [PSH13-26] - Rabbit IgG (Chimeric) - BSA and Azide free	HA751462	Mouse	IHC-Fr
KD	p75 NGF Receptor Recombinant Rabbit mAb [SA39-02]	ET1601-22	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
KD	p75 NGF Receptor Recombinant Rabbit mAb [SA39-02] - BSA and Azide free	HA750014	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	iFluor™ 488 Conjugated p75 NGF Receptor Recombinant Rabbit mAb [SA39-02]	HA720190F	Human, Mouse, Rat	IF-Cell
	p75 NGF Receptor Recombinant Antibody [PSH13-26] - Rabbit IgG (Chimeric)	HA723532	Mouse	IHC-Fr
	p75 NGF Receptor Recombinant Rat mAb [PSH13-26] - BSA and Azide free	HA610274	Mouse, Rat	IHC-Fr, WB
Relative expression	Phospho-TrkB (Y817) Recombinant Rabbit mAb [SC0556] - BSA and Azide free	HA750213	Human, Mouse, Rat	WB, IF-Cell, IF-Tissue, IHC-P, IP
Relative expression	Phospho-TrkB (Y817) Recombinant Rabbit mAb [SC0556]	ET1610-35	Human, Mouse, Rat	WB, IF-Cell, IF-Tissue, IHC-P, IP
	Pro-neuregulin-1, membrane-bound isoform Rabbit pAb	ER1915-45	Human, Mouse, Rat	WB, IHC-P, IF-Tissue, FC
	PSP Rabbit pAb	ER65534	Human	WB
	PTN Recombinant Rabbit mAb [JE47-40]	HA722055	Human, Mouse, Rat	WB, IHC-P
	PTN Rabbit pAb	ER65608	Human, Mouse, Rat	WB, IHC-P
	TGF beta 1 Mouse mAb [B2-F6]	EM1709-93	Human	WB, IHC-P, FC
	TGF beta 1 Rabbit pAb	HA500496	Human	WB
	TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - BSA and Azide free (Detector)	HA723726	Human, Mouse, Rat	ELISA(Det)
	TGF beta 1 Recombinant Rabbit mAb	IRS205RB	Human	mIHC
	TGF beta 1 Recombinant Rabbit mAb [PSH15-21] - BSA and Azide free (Capture)	HA723725	Human, Mouse, Rat	ELISA(Cap)
	Biotin Conjugated TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - Detector	HA723727B	Human, Mouse, Rat	ELISA(Det), ELISA
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17]	HA721143	Human, Mouse, Rat, Monkey	WB, IHC-P, IF-Tissue
	TGF beta 1 Recombinant Rabbit mAb [PSH18-59]	HA724012	Human, Mouse, Monkey	WB, IF-Cell, IHC-P, FC
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17] - BSA and Azide free	HA750485	Human, Mouse, Rat, Monkey	WB, IHC-P, IF-Tissue
	TGF beta 2 Rabbit pAb	ER1917-62	Human, Mouse, Rat	WB
	TGF beta 2 Rabbit pAb	ER1917-63	Human, Mouse	WB, IHC-P, IF-Tissue
	TGF beta 3 Rabbit pAb	ER1917-64	Human, Mouse, Rat, Pig	WB, IHC-P, IF-Tissue, FC
	TGF beta 3 Recombinant Rabbit mAb [PSH09-26]	HA723084	Human, Mouse, Rat	WB, IF-Cell, IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	TGF beta 3 Recombinant Rabbit mAb [PSH09-26] - BSA and Azide free	HA751280	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	TrkA Recombinant Rabbit mAb [SU0354]	ET1608-44	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	TrkA Rabbit pAb	ER1706-86	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	TrkA Recombinant Rabbit mAb [SU0354] - BSA and Azide free	HA750148	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Relative expression	TrkA+B+C Recombinant Rabbit mAb [JJ084-04]	ET1701-16	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	TrkB Recombinant Rabbit mAb [JE58-22]	ET7111-38	Mouse, Rat	WB,IHC-P
	Trk C Rabbit pAb	ER65473	Human	WB
Relative expression	TrkC Recombinant Rabbit mAb [PSH05-40] - BSA and Azide free	HA750999	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Relative expression	TrkC Recombinant Rabbit mAb [PSH05-40]	HA722269	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
KD	TRK fused gene Recombinant Rabbit mAb [SD082-0]	ET1612-52	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	VEGF Receptor 1 Rabbit pAb	RT1651	Human, Mouse, Rat	WB,IP,IF,IHC-P
KD	VEGF Receptor 1 Recombinant Rabbit mAb [SY09-09]	ET1605-11	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	VEGF Receptor 1 Rabbit pAb	ER50803	Human, Mouse, Rat	WB,IHC-P
KD	VEGF Receptor 1 Recombinant Rabbit mAb [SY09-09] - BSA and Azide free	HA750081	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Relative expression	VEGF Receptor 2 Recombinant Rabbit mAb [SU03-42]	ET1608-33	Human, Rat	WB,IF-Cell,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层