

人源M1型巨噬细胞极化相关标志物

高级验证	产品名	货号	反应种属	验证应用
	beta Glucan Receptor Rabbit pAb	ER1904-20	Human, Mouse, Rat	WB,FC
	CD16 Mouse mAb [Q32]	EM60067	Human	WB,IHC-P
Relative expression	CD16 Recombinant Rabbit mAb [PD00-12]	HA721149	Human	IHC-P,IF-Tissue,IF-Cell,mIHC,WB
	CD16 Recombinant Rabbit mAb	IRS017	Human	mIHC
	CD16 Recombinant Rabbit mAb [JE49-79]	ET7109-97	Human, Mouse, Rat	WB
Cell treatment	CD32 Recombinant Rabbit mAb [JM10-70] - BSA and Azide free	HA750382	Human	WB,IHC-P,IF-Cell,FC,IP
Cell treatment	CD32 Recombinant Rabbit mAb [JM10-70]	ET1703-43	Human	WB,IHC-P,IF-Cell,FC,IP
	CD32A + CD32B + CD32C Recombinant Rabbit mAb [JE64-17]	HA721930	Human	WB,IHC-P
	CD64 Mouse mAb [5F5F8]	EM1710-99	Human	WB,IHC-P,FC
	CD64 Mouse mAb [A7E4]	HA601003	Human	WB
Relative expression	CD68 Recombinant Mouse mAb [A3C2-R] - BSA and Azide free	HA610164	Human	WB,IF-Cell,IHC-P
	CD68 Recombinant Mouse mAb	IRS004	Human	mIHC
	CD68 Recombinant Antibody [PDM0-13] - Rabbit IgG (Chimeric)	HA723708	Human	IF-Tissue,IHC-P
	CD68 Recombinant Antibody [PDM0-13] - Rabbit IgG (Chimeric) - BSA and Azide free	HA751544	Human	IF-Tissue,IHC-P
Relative expression	CD68 Mouse mAb [A3C4]	EM1901-95	Human	WB,IHC-P,FC,mIHC
	CD68 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS264RT	Human	mIHC
Relative expression	CD68 Recombinant Mouse mAb [A3C4-R] - BSA and Azide free	HA610165	Human	WB,IF-Cell,IHC-P
	CD68 Mouse mAb [F3-D7]	EM1706-11	Human	WB,IHC-P,IF-Cell,FC
Relative expression	CD68 Recombinant Mouse mAb [PDM0-13]	HA601115	Human	IHC-P,IF-Cell,FC,mIHC,IF-Tissue
	CD68 Recombinant Antibody [PDM0-13] - Rat IgG1 (Chimeric)	HA601520	Human	IF-Tissue,IHC-P

Relative expression	CD68 Recombinant Mouse mAb [A3C2-R]	HA601291	Human	WB,IF-Cell,IHC-P
Relative expression	CD68 Recombinant Mouse mAb [A3C4-R]	HA601292	Human	WB,IF-Cell,IHC-P
	CD68 Recombinant Antibody - Rabbit IgG (Chimeric)	IRS266RB	Human	mIHC
	CD68 Recombinant Antibody [PDM0-13] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610347	Human	IF-Tissue,IHC-P
Relative expression	CD68 Mouse mAb [A3C3]	EM1901-94	Human	WB,IHC-P,FC
Relative expression	CD68 Recombinant Mouse mAb [PDM0-13] - BSA and Azide free	HA610034	Human	IHC-P,IF-Cell,FC,IF-Tissue
	CD68 Rabbit pAb	ER1901-32	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	CD68 Recombinant Mouse mAb [A3C3-R]	HA601290	Human	WB,IF-Cell,IHC-P,FC
Relative expression	CD68 Mouse mAb [A3C2]	EM1901-93	Human	WB,IHC-P,FC
	CD68 Recombinant Rabbit mAb [SN07-27]	ET1611-53	Human	WB
	CD80 Rabbit pAb	ER1905-99	Human, Mouse, Rat	WB
Relative expression	CD80 Recombinant Rabbit mAb [PSH09-79] - BSA and Azide free	HA751321	Human	WB,IF-Cell,FC
Relative expression	CD80 Mouse mAb [8-E5]	M1007-10	Human, Mouse	WB,IHC-P,FC
	CD80 Rabbit pAb	1007-8	Human	WB,IF-Cell
	CD80 Recombinant Mouse mAb [8-E5-R]	HA601247	Human	WB
Relative expression	CD80 Recombinant Rabbit mAb [PSH09-79]	HA723144	Human	WB,IF-Cell,FC
Relative expression	CD80 Recombinant Rabbit mAb [JF100-4]	ET1702-95	Human	WB,IHC-P
	CD80 Recombinant Mouse mAb [8-E5-R] - BSA and Azide free	HA610122	Human	WB
	CD86 Recombinant Rabbit mAb [SJ20-00]	ET1606-50	Human, Mouse, Rat	WB,IF-Cell,IP
	CD86 Recombinant Rabbit mAb [SJ20-00] - BSA and Azide free	HA750109	Human, Mouse, Rat	WB,IF-Cell,IP
	CD86 Rabbit pAb	ER1906-01	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	FCGR1A/CD64 Rabbit pAb	ER1908-78	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
	Human CD16 Recombinant Mouse mAb [PSH04-95]	HA601340	Human	FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PE Conjugated Human CD16 Recombinant Mouse mAb [PSH04-95]	HA600129F	Human	FC
	FITC Conjugated Human CD68 Recombinant Mouse mAb [A3C3-R]	HA600132F	Human	FC
	Human CD80 Recombinant Rabbit mAb [PSH09-61] - BSA and Azide free (Capture)	HA723123	Human	ELISA(Cap)
	Biotin Conjugated Human CD80 Recombinant Rabbit mAb [PSH09-63] - Detector	HA723127B	Human	ELISA(Det),ELISA
	Human CD80 Recombinant Rabbit mAb [PSH09-63] - BSA and Azide free (Detector)	HA723125	Human	ELISA(Det)
	Biotin Conjugated Human CD80 Recombinant Rabbit mAb [PSH09-62] - Detector	HA723126B	Human	ELISA(Det),ELISA
	Human CD80 Recombinant Rabbit mAb [PSH09-62] - BSA and Azide free (Detector)	HA723124	Human	ELISA(Det)
	Biotin Conjugated Human IL-12B p40 Recombinant Rabbit mAb [PSH02-58]	HA722633B	Human	ELISA(Det)
	Biotin Conjugated Human IL-12B p40 Recombinant Rabbit mAb [PSH03-04]	HA722629B	Human	ELISA(Det)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH02-58] - BSA and Azide free (Detector)	HA721835	Human	ELISA(Det)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH03-04] - BSA and Azide free (Detector)	HA721934	Human	ELISA(Det)
	Human IL-12B p40 Recombinant Rabbit mAb [PSH02-57] - BSA and Azide free (Capture)	HA721834	Human	ELISA(Cap)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-85] - BSA and Azide free (Capture)	HA723881	Human	ELISA(Cap)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-87] - BSA and Azide free (Capture)	HA723883	Human	ELISA(Cap)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - BSA and Azide free (Detector)	HA723882	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39]	HA721712	Human	WB
	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free (Detector)	HA721279	Human	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free (Capture)	HA721278	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - Detector	HA723882B	Human	ELISA(Det),ELISA
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free	HA750745	Human	WB
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA722649B	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA721713	Human	WB
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free	HA750744	Human	WB
	Human IL-23A Recombinant Rabbit mAb [PSH03-03] - BSA and Azide free (Capture)	HA721933	Human	ELISA(Cap)
	HRP Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA721293	Human	ELISA(Det)
	Human IL-6 Recombinant Rabbit mAb [PS00-21] - BSA and Azide free (Capture)	HA721167	Human	ELISA(Cap)
	Biotin Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA722233B	Human	ELISA(Det)
	Human IL-6 Mouse mAb [A6D2]	HA601050	Human	WB,ELISA
	Human IL-6 Mouse mAb [A6D1]	HA601051	Human	WB,ELISA
	Human IL-6 Recombinant Rabbit mAb [PS00-15] - BSA and Azide free (Detector)	HA721168	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-47] - BSA and Azide free (Detector)	HA721685	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PS01-75] - BSA and Azide free (Capture)	HA721574	Human	ELISA(Cap)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-45] - BSA and Azide free (Capture)	HA721683	Human	ELISA(Cap)
	Biotin Conjugated Human Interferon gamma Recombinant Rabbit mAb [PSH01-47]	HA722645B	Human	ELISA(Det)
	Human Interferon gamma Recombinant Rabbit mAb [PSH01-46] - BSA and Azide free (Capture)	HA721684	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Interferon gamma Recombinant Rabbit mAb [PS01-76] - BSA and Azide free (Detector)	HA721575	Human	ELISA(Det)
	Biotin Conjugated Human Interferon gamma Recombinant Rabbit mAb [PS01-76] - Detector	HA721576	Human	ELISA(Det),ELISA
	Human MERTK Recombinant Rabbit mAb [PSH12-70] - BSA and Azide free (Detector)	HA723480	Human	ELISA(Det)
	Human MERTK Recombinant Rabbit mAb [PSH12-69] - BSA and Azide free (Capture)	HA723479	Human	ELISA(Cap)
	HRP Conjugated Human MERTK Recombinant Rabbit mAb [PSH12-70] - Detector	HA723481H	Human	ELISA(Det),ELISA
	Biotin Conjugated Human MERTK Recombinant Rabbit mAb [PSH12-70] - Detector	HA723481B	Human	ELISA(Det),ELISA
	Human TNF alpha Recombinant Rabbit mAb [PS00-40]	HA721286	Human	IHC-P,mIHC
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free (Capture)	HA721265	Human	ELISA(Cap)
	Human TNF alpha Recombinant Rabbit mAb	IRS040	Human	mIHC
	HRP Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721294	Human	ELISA(Det)
	Biotin Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721559	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free (Detector)	HA721266	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free	HA750593	Human	IHC-P
	IL12B Rabbit pAb	ER1803-31	Human, Mouse	WB
	IL-1 alpha Rabbit pAb	ER1802-36	Human, Mouse	IHC-P,FC,WB
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R]	HA601036	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Recombinant Rabbit mAb [JJ087-3]	ET1701-39	Human	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	IL-1 beta Mouse mAb [A7F9]	HA601002	Human, Mouse, Rat	WB,ELISA,IF-Cell
Cell treatment	IL-1 beta Mouse mAb [A7A10]	HA601035	Human	WB,ELISA,IF-Cell
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R] - BSA and Azide free	HA610004	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	IL-23 Recombinant Rabbit mAb [JE38-99]	HA721037	Human	WB,IHC-P,IF-Tissue
	IL-6 Rabbit pAb	R1412-2	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	IL-6 Recombinant Rabbit mAb	IRS160RB	Human	mIHC
	IL-6 Mouse mAb	EM30301	Human	WB
	Interferon gamma Recombinant Rabbit mAb [JM10-10] - BSA and Azide free	HA750373	Human	WB,FC
	Interferon gamma Recombinant Rabbit mAb [JM10-10]	ET1703-17	Human	WB,FC
	IRF5 Recombinant Rabbit mAb [SN201-05] - BSA and Azide free	HA750277	Human	WB,IHC-P,FC
	IRF5 Recombinant Rabbit mAb [SD203-07]	ET1612-15	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
	IRF5 Recombinant Rabbit mAb [SN201-05]	ET1611-94	Human	WB,IHC-P,FC
KD	MERTK Recombinant Rabbit mAb [SR29-07]	ET1602-21	Human	WB,IHC-P,FC,IF-Tissue
KD	MERTK Recombinant Rabbit mAb [SR29-07] - BSA and Azide free	HA750036	Human	WB,IHC-P,IP,FC,IF-Tissue
	MHC Class II Recombinant Rabbit mAb [JA10-94]	ET1704-13	Human	WB,IP,IF-Cell,IF-Tissue,IHC-P,mIHC
Cell treatment, Relative expression	Phospho-STAT1 (S727) Recombinant Rabbit mAb [SN67-04] - BSA and Azide free	HA750247	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment, Relative expression	Phospho-STAT1 (S727) Recombinant Rabbit mAb [SN67-04]	ET1611-20	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-STAT1 (Y701) Recombinant Rabbit mAb [PSH04-02] - BSA and Azide free	HA750925	Human	WB,IF-Cell
Cell treatment	Phospho-STAT1 (Y701) Recombinant Rabbit mAb [PSH04-02]	HA722083	Human	WB,IF-Cell
KD	STAT1 Recombinant Rabbit mAb [SD20-75]	ET1612-22	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	STAT1 Recombinant Mouse mAb [G3-B11-R] - BSA and Azide free	HA610181	Human	WB,IHC-P
	STAT1 Recombinant Mouse mAb [G3-B11-R]	HA601309	Human	WB,IHC-P
	STAT1 Rabbit pAb	R1408-2	Human, Mouse, Rat, Zebrafish	WB,IF-Cell
	STAT1 Mouse mAb [G3-B11]	M1407-1	Human	WB,IF-Cell,IHC-P,IF-Tissue
KD	STAT1 alpha Recombinant Rabbit mAb [SJ01-89]	ET1606-39	Human, Mouse	WB,IP,FC,IHC-P,IF-Cell,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	STAT1 alpha Recombinant Rabbit mAb [SJ01-89] - BSA and Azide free	HA750102	Human, Mouse	WB,IP,FC,IHC-P,IF-Cell,IF-Tissue
	STAT1 alpha/beta Rabbit pAb	ER31221	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	TNF alpha Rabbit pAb	ER65189	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	TNF alpha Rabbit pAb	ER1919-22	Human, Mouse	WB
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07]	HA722022	Human, Mouse	WB
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free	HA750897	Human, Mouse	WB