

磷酸化抗体-产品清单

高级验证	产品名	货号	反应种属	验证应用
	Phospho-4E-BP1/2/3 (T45) Mouse mAb [3G2]	RT1004	Human, Mouse, Rat	WB,IP,IF
Cell treatment	Phospho-Acetyl Coenzyme A Carboxylase (S79) Recombinant Rabbit mAb [PSH07-13]	HA722801	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-Acetyl Coenzyme A Carboxylase (S79) Recombinant Rabbit mAb [JE63-95]	HA721714	Human, Mouse, Rat	WB
Cell treatment	Phospho-Acetyl Coenzyme A Carboxylase (S79) Recombinant Rabbit mAb [JE55-54]	HA721516	Human, Mouse, Rat	WB
Cell treatment	Phospho-AKT1 (S124) Recombinant Rabbit mAb [JJ08-46]	ET1701-36	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-AKT1 (S129) Recombinant Rabbit mAb [JE59-50]	HA721677	Human, Mouse, Rat	WB
Cell treatment	Phospho-AKT1 (T450) Recombinant Rabbit mAb [SD08-12]	ET1612-73	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell,FC
Cell treatment	Phospho-AKT2 (S474) Recombinant Rabbit mAb [PSH09-46]	HA723108	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-AKT2 (S474) Rabbit pAb	HA500116	Human, Mouse	WB,IHC-P
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05]	ET1607-73	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44]	HA722129	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10]	HA722951	Human	WB,IF-Cell,FC
	Phospho-alpha Synuclein (S129) Recombinant Rabbit mAb [JB22-44]	ET7107-30	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-AMPK alpha 1 (S496) Recombinant Rabbit mAb [SD0810]	ET1612-72	Human, Mouse	WB,IP,IHC-Fr
Cell treatment	Phospho-AMPK alpha 1 (T183) +AMPK alpha 2 (T172) Recombinant Rabbit mAb [PSH14-13]	HA723603	Human, Mouse, Rat	WB,IHC-P
	Phospho-AMPK alpha 2 (S345) Recombinant Rabbit mAb [JJ08-19]	ET1701-37	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-AMPK alpha 2 (S491) Recombinant Rabbit mAb [JE59-49]	HA722567	Human	WB
	Phospho-AP2M1 (Thr156) Rabbit pAb	ER65961	Human, Mouse, Rat	WB
Cell treatment	Phospho-AS160 (T642) Recombinant Rabbit mAb [JE60-97]	HA721489	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-ATF2 (T71) Recombinant Rabbit mAb [SC05-90]	ET1610-30	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IP,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-ATM (S1981) Recombinant Rabbit mAb [JM93-25]	ET1705-50	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cell treatment	Phospho-ATR (T1989) Recombinant Rabbit mAb [PS01-18]	HA721190	Human	WB,IF-Cell
Cell treatment	Phospho-ATR (T1989) Rabbit pAb	HA500405	Human	WB
Cell treatment	Phospho-BCAR1 (Y410) Recombinant Rabbit mAb [JE47-97]	HA722124	Human	WB
Cell treatment	Phospho-BCKDHA (S293) Recombinant Rabbit mAb [PSH14-02]	HA723592	Human, Mouse, Rat, Monkey	WB
Cell treatment	Phospho-Bcl-2 (S70) Recombinant Rabbit mAb [PSH11-02]	HA723269	Human	WB,IF-Cell,FC
Cell treatment	Phospho-Beta Catenin (S33 + S37) Recombinant Rabbit mAb [JE59-59]	HA721580	Human, Rat, Mouse	WB
Cell treatment	Phospho-Beta Catenin (S552) Recombinant Rabbit mAb [PSH08-72]	HA723024	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	Phospho-Beta Catenin (T41/S45) Recombinant Rabbit mAb [JE54-02]	HA722316	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-BRAF (T401) Recombinant Rabbit mAb [JJ08-72]	ET1701-20	Human, Mouse, Rat	WB,IF-Cell,IP,IHC-P
Cell treatment	Phospho-BTK (Y223) Recombinant Rabbit mAb [ST51-05]	ET1609-55	Human, Mouse	WB,IP
	Phospho-CaMKII (T286) Mouse mAb [2G1]	RT1092	Human, Mouse, Rat	WB,IP,IF
Cell treatment	Phospho-CaMKII (T286) Recombinant Rabbit mAb [JE51-81]	HA721794	Human, Mouse, Rat	WB
Cell treatment	Phospho-CBL (Y774) Recombinant Rabbit mAb [JE62-93]	HA722786	Human	WB
Cell treatment	Phospho-CDC37 (S13) Recombinant Rabbit mAb [SD08-50]	ET1612-74	Human, Mouse, Rat	WB,IP
Cell treatment	Phospho-Cdc6 (S54) Recombinant Rabbit mAb [SD08-47]	ET1612-96	Human, Mouse, Rat	WB,IF-Tissue,IHC-P
Cell treatment	Phospho-Cdk1/2/3 (T14) Recombinant Rabbit mAb [JE04-37]	HA722272	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-CDK1 (T161) Recombinant Rabbit mAb [PSH03-42]	HA721987	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	Phospho-CDK1 (T161) + CDK2/CDK3 (T160) Recombinant Rabbit mAb [JE75-37]	HA722582	Human	WB,Dot Blot
Cell treatment	Phospho-Cdk2 (Y15) Recombinant Rabbit mAb [SN72-04]	ET1611-52	Human, Mouse, Rat	WB,IHC-P,IP,IF-Tissue
Cell treatment	Phospho-CDK (Y15) Recombinant Rabbit mAb [PSH09-45]	HA723107	Human, Mouse	WB,IHC-P
Cell treatment	Phospho-Chk1 (S280) Recombinant Rabbit mAb [JE47-66]	HA721471	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Chk1 (S296) Recombinant Rabbit mAb [SN06-50]	ET1611-76	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Phospho-Chk1 (S317) Recombinant Rabbit mAb [PSH03-84]	HA722050	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-Chk1 (S345) Recombinant Rabbit mAb [PS01-17]	HA721189	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-Chk2 (T68) Recombinant Rabbit mAb [JE43-12]	HA721633	Human, Mouse, Rat	WB
Cell treatment	Phospho-c-Jun (S63) Recombinant Rabbit mAb [SY0297]	ET1608-4	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit mAb [JE65-46]	HA722475	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-c-Jun (S73)+JunD (S100) Recombinant Rabbit mAb [PSH07-12]	HA722800	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-c-Jun (T91) Recombinant Rabbit mAb [JJ080-9]	ET1701-32	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-c-Myc (S62) Recombinant Rabbit mAb [ST49-08]	ET1609-64	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-c-Myc (T58) Recombinant Rabbit mAb [SN60-01]	ET1611-24	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC,IHC-P
Cell treatment	Phospho-c-Myc (T58 + S62) Recombinant Rabbit mAb [SZ02-06]	ET1603-24	Human, Mouse	WB,IP,FC
Cell treatment	Phospho-Cofilin (S3) Recombinant Rabbit mAb [PSH07-18]	HA722806	Human, Rat, Monkey	WB
Cell treatment	Phospho-Creb (S133) Recombinant Rabbit mAb [JB25-40]	ET7107-93	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,IF-Tissue
Cell treatment	Phospho-CrkL (Y207) Recombinant Rabbit mAb [JE55-58]	HA721675	Human, Mouse, Rat	WB
Cell treatment	Phospho-CRMP2 (T514) Recombinant Rabbit mAb [JE57-25]	HA721674	Human, Mouse, Rat	WB
	Phospho-Cyclin E1 (T77) Recombinant Rabbit mAb [SD2025]	ET1612-31	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-Cytosolic Phospholipase A2 (S505) Recombinant Rabbit mAb [PSH07-16]	HA722804	Human, Mouse	WB
Cell treatment	Phospho-Cytosolic Phospholipase A2 (S505) Recombinant Rabbit mAb [PSH07-56]	HA722849	Human, Mouse	WB,Dot Blot
Cell treatment	Phospho-DNA PKcs (S2056) Recombinant Rabbit mAb [JE41-92]	HA722066	Human	WB,IF-Cell
Cell treatment	Phospho-DRP1 (S616) Recombinant Rabbit mAb [PSH06-77]	HA722760	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC
Cell treatment	Phospho-EEF2 (T56) Recombinant Rabbit mAb [PSH06-36]	HA722670	Human, Mouse, Rat	WB,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Phospho-EGFR (S695) Recombinant Rabbit mAb [SN07-36]	ET1611-50	Human	WB,IF-Cell,IF-Tissue
Cell treatment	Phospho-EGFR (T669) Recombinant Rabbit mAb [PSH10-03]	HA723172	Human, Mouse, Rat	WB,IP
Cell treatment	Phospho-EGFR (Y1045) Recombinant Rabbit mAb [PSH03-28]	HA721970	Human	WB,IF-Cell,FC
Cell treatment	Phospho-EGFR (Y1068) Recombinant Rabbit mAb [SD2055]	ET1612-30	Human	WB,IF-Cell,IF-Tissue,IP,FC
Cell treatment	Phospho-EGFR (Y1092) Recombinant Rabbit mAb [SJ0194]	ET1606-44	Human, Mouse	WB,IF-Tissue,IHC-P
	Phospho-EGFR (Y1148) Rabbit pAb	HA500580	Human	WB,IF-Cell
Cell treatment	Phospho-EGFR (Y1173) Recombinant Rabbit mAb [SC57-04]	ET1610-4	Human, Mouse	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-EGFR (Y845) Recombinant Rabbit mAb [JE46-77]	HA721488	Human	WB
Cell treatment, Relative expression	Phospho-EIF2S1 (S51) Recombinant Rabbit mAb [SZ01-06]	ET1603-14	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Phospho-EIF2S1 (S52) Rabbit pAb	RT1198	Human, Mouse, Rat	WB,IP,IF,IHC-P
Cell treatment	Phospho-eIF4B (S406) Recombinant Rabbit mAb [JE58-32]	HA721676	Human, Mouse, Rat	WB
Cell treatment	Phospho-eIF4B (S406) Recombinant Rabbit mAb [JE63-80]	HA722149	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	Phospho-eIF4EBP1 (S65) Recombinant Rabbit mAb [PSH09-44]	HA723106	Human	WB,IF-Cell,IP
Cell treatment	Phospho-eIF4EBP1 (T70) Recombinant Rabbit mAb [PSH09-91]	HA723222	Human, Mouse	WB,IHC-P,IP
Cell treatment	Phospho-eIF4E (S209) Recombinant Rabbit mAb [SU0396]	ET1608-66	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
Cell treatment	Phospho-eIF4G (S1108) Recombinant Rabbit mAb [PS01-08]	HA721276	Human	WB,IHC-P
	Phospho-ELK1 (S383) Rabbit pAb	ER65967	Human, Mouse, Rat	WB,IHC
Cell treatment	Phospho-ELK1 (S383) Recombinant Rabbit mAb [PSH14-65]	HA723663	Human, Mouse	WB,IHC-P,IP
Cell treatment	Phospho-eNOS (S1177) Recombinant Rabbit mAb [PS01-10]	HA721188	Human, Dog, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-ErbB3/HER3 (Y1289) Recombinant Rabbit mAb [JE51-25]	HA722178	Human	WB,IF-Cell
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4]	ET1603-22	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01]	ET1610-13	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Estrogen Receptor alpha (S118) Recombinant Rabbit mAb [SC05-87]	ET1610-32	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-Ezrin (T567)/Radixin (T564)/Moesin (T558) Recombinant Rabbit mAb [PSH05-73]	HA722440	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Phospho-FAK (S722) Mouse mAb [4G1]	RT1216	Human, Mouse, Rat	WB,IP,IF,IHC-P
Cell treatment	Phospho-FAK (Y397) Recombinant Rabbit mAb [SC54-07]	ET1610-34	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-FOXO1A (T24)+ FOXO3A (T32)+ FOXO4 (T28)+ FOXO6 (T26) Recombinant Rabbit mAb [PSH06-00]	HA722499	Human, Mouse	WB,IF-Cell,FC
Cell treatment	Phospho-FOXO3a (S253) Recombinant Rabbit mAb [ST49-01]	ET1609-49	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-GATA3 (S308) Recombinant Rabbit mAb [ST44-09]	ET1609-17	Human, Mouse, Rat	WB,IHC-P,IP,IF-Tissue
Cell treatment	Phospho-GCN2 (T899) Recombinant Rabbit mAb [JE00-51]	HA722249	Human, Mouse, Rat	WB
	Phospho-Glucose Transporter GLUT4 (S488) Recombinant Rabbit mAb [JE43-23]	HA721866	Mouse, Human, Rat	WB
	Phospho-GluR1 (AMPA subtype) (S831) Recombinant Rabbit mAb [JE44-96]	HA721865	Human, Mouse, Rat	WB
Cell treatment	Phospho-GluR1 (S845) Recombinant Rabbit mAb [JJ2009]	ET1701-28	Mouse, Rat, Human	WB,FC
Cell treatment	Phospho-Glycogen synthase 1 (S641) Recombinant Rabbit mAb [SR46-06]	ET1602-13	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-GSK3 (alpha + beta) (Y216 + Y279) Recombinant Rabbit mAb [SY26-05]	ET1607-54	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,IF-Tissue,FC
Cell treatment	Phospho-GSK3 alpha (S21) Recombinant Rabbit mAb [PSH09-11]	HA723069	Human, Mouse	WB,IHC-P,Dot Blot
Cell treatment	Phospho-GSK3 beta (S9) Recombinant Rabbit mAb [SY02-71]	ET1607-60	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Phospho-HER2 / ErbB2 (T686) Mouse mAb [4G3]	RT1275	Human, Mouse, Rat	WB,IP
Cell treatment	Phospho-HER2 / ErbB2 (Y1139) Recombinant Rabbit mAb [JE60-90]	HA722201	Human	WB,IF-Cell,FC
Cell treatment	Phospho-HER2 / ErbB2 (Y1221 + Y1222) Recombinant Rabbit mAb [JE44-12]	HA721433	Human	WB
Cell treatment	Phospho-HER2 / ErbB2 (Y1248) Recombinant Rabbit mAb [PSH04-03]	HA722084	Human	WB,IF-Cell
Cell treatment	Phospho-HER3/ErbB3 (Y1222) Recombinant Rabbit mAb [JE51-89]	ET7110-10	Human	WB,IHC-P
Cell treatment	Phospho-Histone H1.3 (T17)+Histone H1.4 (T17) Recombinant Rabbit mAb [SR38-03]	ET1602-11	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Histone H2A.X (S139) Recombinant Rabbit mAb [SR33-09]	ET1602-2	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IF-Cell,IF- Tissue
Cell treatment	Phospho-Histone H3.3 (S31) Recombinant Rabbit mAb [JE62-92]	HA721034	Human, Rat	WB,IF-Cell
Cell treatment	Phospho-Histone H3 (S10) Recombinant Rabbit mAb [SA31-01]	ET1601- 30	Human, Mouse, Rat	WB,IF-Cell,IHC- P,IP,FC,ChIP
Cell treatment	Phospho-Histone H3 (S28) Recombinant Rabbit mAb [JE56-06]	HA721795	Human, Rat, Mouse	WB,IF-Cell,IHC- P,FC,ChIP
Cell treatment	Phospho-Histone H3 (T3) Recombinant Rabbit mAb [JE42-48]	HA721796	Human, Mouse, Rat	WB,IHC-P,IF- Cell,FC
Cell treatment	Phospho-Hormone sensitive lipase (S660) Recombinant Rabbit mAb [PSH10-28]	HA723195	Mouse, Rat, Human	WB
Cell treatment	Phospho-Hormone sensitive lipase (S853) Recombinant Rabbit mAb [SN06-39]	ET1611- 19	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-HSF1 (S326) Recombinant Rabbit mAb [SU31-03]	ET1608- 11	Human	WB,IHC-P,IP
Cell treatment	Phospho-Hsp27 (S78) Recombinant Rabbit mAb [JJ08-70]	ET1701- 19	Human	WB,IF-Cell,IF- Tissue,IHC-P,IP
	Phospho-Hsp27 (S82) Recombinant Rabbit mAb [SN06-73]	ET1611- 16	Human, Mouse, Rat	WB,IP,IF-Cell,FC
Cell treatment	Phospho-IKB alpha (S32) Recombinant Rabbit mAb [ST53-05]	ET1609- 78	Human, Mouse	WB,IF-Cell,IP
Cell treatment	Phospho-IKB alpha (S32/S36) Recombinant Rabbit mAb [PSH06-76]	HA722770	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Phospho-IKB alpha (S36) Recombinant Rabbit mAb [JE56-65]	HA721802	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-IRE1 (S724) Recombinant Rabbit mAb [PSH03-35]	HA721980	Human, Mouse	WB,IF-Cell,IHC- P,FC
Cell treatment	Phospho-IRF3 (S386) Recombinant Rabbit mAb [SU03-28]	ET1608- 22	Human	WB,IF-Cell,IF-Tissue
Cell treatment	Phospho-IRF3 (S396) Recombinant Rabbit mAb [PSH06-75]	HA722772	Human, Mouse, Rat	WB,IHC-P,IF- Cell,FC
Cell treatment	Phospho-JAK2 (Y1007 + Y1008) Recombinant Rabbit mAb [SY24-03]	ET1607- 34	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-JNK1/2/3 (T183 + T183 + T221) Recombinant Rabbit mAb [ST500]	ET1609- 42	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC- P,IP,FC,IHC-Fr
Cell treatment	Phospho-JunD (S255) Recombinant Rabbit mAb [JJ08-21]	ET1701- 35	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P
Cell treatment	Phospho-KAP1 (S824) Recombinant Rabbit mAb [JE50-99]	ET7110- 11	Human, Mouse, Rat	WB,IHC-P,IP,Dot Blot,IF-Cell
Cell treatment	Phospho-Lamin A + Lamin C (S22) Recombinant Rabbit mAb [PSH11-47]	HA723339	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Lamin A + Lamin C (S22) Recombinant Rabbit mAb [PSH11-48]	HA723340	Human, Mouse, Rat	WB
	Phospho-LATS1/2 (Ser909) Rabbit pAb	ER65962	Human, Mouse, Rat	WB,IHC
Cell treatment	Phospho-LRRK2 (S935) Recombinant Rabbit mAb [JE42-11]	HA722315	Mouse, Human	WB,IF-Cell,FC
Cell treatment, Relative expression	Phospho-LYN (Y397)/LCK (Y394)/HCK (Y411)/BLK (Y389) Recombinant Rabbit mAb [PSH09-43]	HA723105	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-MAP2 (S136) Recombinant Rabbit mAb [JE44-25]	HA722351	Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Tissue
	Phospho-MCM2 (S108) Recombinant Rabbit mAb [JE63-52]	HA721021	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-MCM2 (S27) Recombinant Rabbit mAb [JE63-40]	HA722574	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-MCM2 (S40) Recombinant Rabbit mAb [JE35-01]	HA722216	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-MCM2 (S41) Recombinant Rabbit mAb [JE50-04]	HA722279	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
Cell treatment	Phospho-MDM2 (S166) Recombinant Rabbit mAb [JE64-20]	HA721791	Human, Mouse, Rat	WB
Cell treatment	Phospho-MEK1/2 (S218 + S222) Recombinant Rabbit mAb [ST0490]	ET1609-50	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Phospho-MEK1 (S298) Recombinant Rabbit mAb [SD206-7]	ET1612-40	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-MEK1 (T292) Recombinant Rabbit mAb [SD2088]	ET1612-42	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-MEK3 (S189)+MEK6 (S207) Recombinant Rabbit mAb [PSH11-17]	HA723299	Human, Mouse	WB
Cell treatment	Phospho-Met (Y1234/Y1235) Recombinant Rabbit mAb [PSH06-90]	HA722765	Human	WB
Cell treatment	Phospho-Met (Y1349) Recombinant Rabbit mAb [JE51-68]	ET7110-14	Human, Rat	WB,IHC-P
Cell treatment	Phospho-MLKL (S345) Recombinant Rabbit mAb [JM92-37]	ET1705-51	Mouse	WB,IHC-P,IHC-Fr
Cell treatment	Phospho-MSK1 (S376) Recombinant Rabbit mAb [SN63-05]	ET1611-17	Human, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-MST1/MST2 (T180 + T183) Recombinant Rabbit mAb [JE46-64]	HA721737	Human, Mouse, Rat	WB
Cell treatment	Phospho-mTOR (S2448) Mouse mAb [A5D5]	HA600094	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-mTOR (S2481) Recombinant Rabbit mAb [JE59-62]	HA721632	Human	WB
	Phospho-MYL9 (T18 + S19) Recombinant Rabbit mAb [PSH12-71]	HA723482	Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Phospho-Myosin light chain (S20) Rabbit pAb	ER1913-46	Human, Mouse	WB,IHC-P,IF-Tissue
Cell treatment	Phospho-MYPT1 (T853) Recombinant Rabbit mAb [PSH08-71]	HA723023	Human, Mouse, Rat	WB
Cell treatment	Phospho-NDRG1 (T346) Recombinant Rabbit mAb [PSH07-17]	HA722805	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue
Cell treatment	Phospho-NF-kB p65 (S529) Recombinant Rabbit mAb [SP07-00]	ET1604-27	Human	WB,IHC-P,FC,IP
Cell treatment	Phospho-NF-kB p65 (S536) Recombinant Rabbit mAb [PSH10-58]	HA723223	Human, Mouse, Rat	WB,IF-Cell
	Phospho-NLRC4 (S533) Rat mAb [B7-B7]	EM1707-24	Human, Mouse	WB,IF-Cell
Cell treatment, Relative expression	Phospho-Nrf2 (S40) Recombinant Rabbit mAb [SU0334]	ET1608-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Cell treatment	Phospho-Nucleolin (T76) Recombinant Rabbit mAb [JE64-03]	HA722942	Human	WB,IHC-P
Cell treatment	Phospho-Nucleolin (T84) Recombinant Rabbit mAb [JE64-18]	HA722920	Human	WB,IHC-P
	Phospho-Nucleophosmin (S125) Recombinant Rabbit mAb [JE51-14]	HA722559	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-Nucleophosmin (T199) Recombinant Rabbit mAb [JE63-16]	HA722125	Human	WB
	Phospho-p130 (T986) Recombinant Rabbit mAb [JE64-95]	HA722472	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-p38 alpha (T180) Recombinant Rabbit mAb [PSH06-25]	HA722657	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-p38 alpha (T180 + Y182) Recombinant Rabbit mAb [PSH06-35]	HA722669	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-p38 (T180) Rabbit pAb	ER2001-51	Human, Mouse, Rat	WB
Cell treatment	Phospho-p38 (T180) Recombinant Rabbit mAb [JE77-37]	HA722150	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-p38 (T180 + Y182) Rabbit pAb	ER2001-52	Human, Mouse, Rat	WB
	Phospho-p38 (T180 + Y182) Rabbit pAb	ER1903-01	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
Cell treatment	Phospho-p53 (S15) Recombinant Rabbit mAb [PSH01-98]	HA721756	Human	WB,IF-Cell
Cell treatment	Phospho-p53 (S20) Recombinant Rabbit mAb [PSH06-86]	HA722761	Human	WB
Cell treatment	Phospho-p53 (S33) Recombinant Rabbit mAb [JE41-97]	HA722703	Human, Mouse, Rat	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-p53 (S376) Recombinant Rabbit mAb [ST0440]	ET1609-14	Human, Mouse	WB
Cell treatment	Phospho-p53 (S392) Recombinant Rabbit mAb [SI17-04]	ET1606-24	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-p53 (T55) Recombinant Rabbit mAb [ST0439]	ET1609-13	Human	WB,IHC-P
	Phospho-p70 S6 Kinase (S434) Mouse mAb [2G1]	RT1456	Human, Mouse, Rat	WB,IP,IF,IHC-P
Cell treatment	Phospho-p70 S6 Kinase (T389) Recombinant Rabbit mAb [PSH02-23]	HA721803	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-p95/NBS1 (S343) Recombinant Rabbit mAb [SY0215]	ET1607-5	Human, Mouse	WB,IF-Tissue,IHC-P,IP
	Phospho-PAK1 (S144)+PAK2 (S141)+PAK3 (S139) Recombinant Rabbit mAb [SU03-70]	ET1608-58	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Cell treatment	Phospho-PAK4+PAK5+PAK6 (S474 + S560 + S602) Recombinant Rabbit mAb [JE57-42]	HA722146	Human	WB
Cell treatment	Phospho-Paxillin (Y118) Recombinant Rabbit mAb [JE09-10]	HA722526	Human	WB
Cell treatment	Phospho-PDGFR beta (Y1009) Recombinant Rabbit mAb [JE46-72]	HA721736	Human, Mouse, Rat	WB
Cell treatment	Phospho-PDGFR beta (Y1021) Recombinant Rabbit mAb [JE48-54]	HA721833	Mouse, Rat, Human	WB
Cell treatment	Phospho-PDGFR beta (Y740) Recombinant Rabbit mAb [SN0646]	ET1611-29	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue
Cell treatment	Phospho-PDGFR beta (Y751) Recombinant Rabbit mAb [PSH13-61]	HA723550	Mouse, Rat, Human	WB
Cell treatment	Phospho-PDHA1 (S293) Recombinant Rabbit mAb [JE58-27]	HA721443	Human, Mouse, Rat	WB,ICC,IHC-P
Cell treatment	Phospho-Phospholamban (S16) Recombinant Rabbit mAb [JE47-45]	HA722291	Mouse, Rat	WB
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-39]	HA721673	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-38]	HA721672	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PKA alpha/beta/gamma (catalytic subunit) (T197) Recombinant Rabbit mAb [JE45-17]	HA721864	Human, Mouse, Rat	WB
Cell treatment	Phospho-PKA R2 (S99) Recombinant Rabbit mAb [SC54-04]	ET1610-29	Human, Mouse, Rat, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH16-27]	HA723825	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH15-91]	HA723796	Human, Mouse, Rat	WB,IF-Cell,Dot Blot

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-PKC alpha (T497) Recombinant Rabbit mAb [JE44-81]	HA722278	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-PKC alpha (T638) Recombinant Rabbit mAb [JF0964]	ET1702-17	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,Dot Blot
Cell treatment	Phospho-PKC delta (Y311) Recombinant Rabbit mAb [JE44-20]	HA722458	Human, Mouse, Rat	WB
Cell treatment	Phospho-PKC (S729) Recombinant Rabbit mAb [JE77-54]	HA722846	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-PKC zeta (T560) Recombinant Rabbit mAb [SN0753]	ET1611-78	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-PKR (T446) Recombinant Rabbit mAb [SY230]	ET1607-20	Human	WB,IHC-P,IP,IF-Cell
Cell treatment	Phospho-PKR (T451) Recombinant Rabbit mAb [JE46-70]	HA722067	Human	WB,FC,IF-Cell
Cell treatment	Phospho-PLC-gamma-1 (S1248) Recombinant Rabbit mAb [PSH13-62]	HA723551	Human	WB,IP
Cell treatment	Phospho-PLC-gamma-1 (Y783) Recombinant Rabbit mAb [PSH06-29]	HA722661	Human, Mouse, Rat	WB
	Phospho-PLK1 (T210) Recombinant Rabbit mAb [JJ20-58]	ET1701-33	Human, Mouse	WB,IF-Cell,IF-Tissue,FC,IHC-P
Cell treatment	Phospho-POLR2A (S2) Recombinant Rabbit mAb [JM11-51]	ET1703-86	Human, Mouse, Rat	WB,IP,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-POLR2A (S5) Recombinant Rabbit mAb [JM51-21]	ET1703-87	Human, Mouse, Rat	WB,IP,IF-Cell,IF-Tissue,IHC-P,FC
	Phospho-PP2A (Y307) Recombinant Rabbit mAb [ST49-05]	ET1609-40	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-PPP1CA (T320) Recombinant Rabbit mAb [PSH10-27]	HA723194	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-PRAS40 (T246) Recombinant Rabbit mAb [JE53-94]	HA721790	Human, Mouse, Rat	WB
	Phospho-PTEN (S380) Recombinant Rabbit mAb [JJ08-81]	ET1701-46	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-Rac1+Cdc42 (S71) Recombinant Rabbit mAb [JJ086-9]	ET1701-40	Human	WB,IF-Cell,IHC-P
Cell treatment	Phospho-Rad17 (S656) Recombinant Rabbit mAb [JE63-30]	HA722148	Human, Mouse	WB,IF-Cell,FC
Cell treatment	Phospho-Raf1 (S259) Recombinant Rabbit mAb [SD85-07]	ET1612-87	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Phospho-Raf1 (S43) Recombinant Rabbit mAb [SD087-05]	ET1612-75	Human, Mouse, Rat	WB
	Phospho-Raf1 (S621) Recombinant Rabbit mAb [JJ085-05]	ET1701-3	Human, Rat, Mouse	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Rb (S608) Recombinant Rabbit mAb [PSH11-04]	HA723271	Human, Monkey	WB
Cell treatment	Phospho-Rb (S780) Recombinant Rabbit mAb [JE44-61]	HA721513	Human, Mouse, Rat	WB
Cell treatment	Phospho-Rb (S795) Recombinant Rabbit mAb [PSH06-89]	HA722764	Human	WB,IF-Cell
Cell treatment	Phospho-Rb (S807) Recombinant Rabbit mAb [SR3-82]	ET1602-36	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-Rb (S807/S811) Recombinant Rabbit mAb [PSH11-01]	HA723268	Human, Mouse	WB,IHC-P,IF-Cell,FC,IP
Cell treatment	Phospho-RICTOR (T1135) Recombinant Rabbit mAb [PSH09-92]	HA723161	Human, Mouse	WB
Cell treatment	Phospho-RIP3 (S232) Recombinant Rabbit mAb [JE44-62]	HA721428	Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-RPA32/RPA2 (T21) Recombinant Rabbit mAb [SN06-36]	ET1611-18	Human, Mouse, Rat	WB,IP,IF-Cell
	Phospho-RPS6 (S235 + S236) Recombinant Rabbit mAb [PS01-04]	HA721275	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Phospho-RPS6 (S240 + S244) Recombinant Rabbit mAb [PS00-65]	HA721589	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Phospho-RSK1 p90 (T359 + S363) Recombinant Rabbit mAb [SU03-65]	ET1608-53	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-RSK1 p90 (T573) Recombinant Rabbit mAb [JE50-12]	HA721793	Human, Mouse, Rat	WB
Cell treatment	Phospho-RSK1 (S380) Recombinant Rabbit mAb [SC05-32]	ET1610-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IP
	Phospho-SHIP (Y1020) Recombinant Rabbit mAb [JJ082-03]	ET1701-5	Human	WB,IF-Cell,IP
Cell treatment	Phospho-SHP2 (Y542) Recombinant Rabbit mAb [SN61-01]	ET1611-22	Mouse, Human	WB,IP
Cell treatment	Phospho-SIRT1 (S47) Recombinant Rabbit mAb [SN06-40]	ET1611-31	Human	WB,IF-Cell,IHC-P
	Phospho-SIRT1 (T530) Recombinant Rabbit mAb [JJ206-6]	ET1701-27	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-Smad1/5/9 (S463/S465/S467) Recombinant Rabbit mAb [JE59-46]	HA722566	Human, Mouse, Rat	WB,IHC,IF-Cell
Cell treatment	Phospho-Smad2 (S250) Recombinant Rabbit mAb [SD207-1]	ET1612-32	Human, Mouse, Rat	WB,FC,IP,IHC-P
Cell treatment	Phospho-Smad2 (S255) Recombinant Rabbit mAb [JF0882]	ET1702-34	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-SMAD2 (S465+467) Recombinant Rabbit mAb [PSH05-76]	HA722443	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-Smad3 (S423 + S425) Recombinant Rabbit mAb [ST0493]	ET1609-41	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Smad5 (S463 + S465) Recombinant Rabbit mAb [SY09-03]	ET1605-5	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-SMC1 (S957) Recombinant Rabbit mAb [SC65-03]	ET1610-55	Human	WB,IHC-P,IF-Cell
Cell treatment	Phospho-Src (S17) Recombinant Rabbit mAb [PSH11-03]	HA723270	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP
Cell treatment	Phospho-Src (Y419) Recombinant Rabbit mAb [ST0800]	ET1609-15	Human, Mouse, Rat	WB,FC,IF-Cell
Cell treatment, Relative expression	Phospho-STAT1 (S727) Recombinant Rabbit mAb [SN67-04]	ET1611-20	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Cell treatment	Phospho-STAT1 (Y701) Recombinant Rabbit mAb [PSH04-02]	HA722083	Human	WB,IF-Cell
Cell treatment, KD	Phospho-STAT3 (S727) Recombinant Rabbit mAb [SY24-09]	ET1607-39	Human, Rat, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-STAT3 (Y705) Recombinant Rabbit mAb [SZ43-01]	ET1603-40	Human, Mouse, Rat	WB,IHC-P,IP
	Phospho-STAT4 (Y693) Rabbit pAb	ER1919-55	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-STAT5 (Y694) Recombinant Rabbit mAb [SC05-31]	ET1610-48	Human, Rat, Mouse	WB,IHC-P,IF-Cell
Cell treatment	Phospho-STAT6 (Y641) Recombinant Rabbit mAb [PSH06-98]	HA722776	Human, Mouse	WB
Cell treatment	Phospho-STING (S366) Recombinant Rabbit mAb [PSH09-72]	HA723137	Human	WB,Dot Blot,IP
Cell treatment	Phospho-Syk (Y352) Recombinant Rabbit mAb [PSH11-66]	HA723362	Human	WB,IF-Cell,FC
Cell treatment	Phospho-Syk (Y525 + Y526) Recombinant Rabbit mAb [PSH03-70]	HA722034	Human	WB,IF-Cell
	Phospho-Synapsin I (S9) Recombinant Rabbit mAb [SN63-04]	ET1611-26	Mouse, Rat, Human	WB,IHC-P
Cell treatment	Phospho-TAK1 (S412) Recombinant Rabbit mAb [PSH14-64]	HA723662	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-TAK1 (S439) Recombinant Rabbit mAb [JE44-93]	HA722242	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S198) Recombinant Rabbit mAb [JE66-79]	HA721668	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S199) Recombinant Rabbit mAb [JE66-80]	HA721667	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S199) Recombinant Rabbit mAb [JE46-95]	HA721800	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S202) Recombinant Rabbit mAb [JE66-81]	HA721666	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Tau (S214) Recombinant Rabbit mAb [JE66-83]	HA721797	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S324) Recombinant Rabbit mAb [JE30-73]	HA721798	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S356) Recombinant Rabbit mAb [JE46-86]	HA721799	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S396) Recombinant Rabbit mAb [SN62-09]	ET1611-68	Mouse, Rat, Human, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-Tau (S404) Recombinant Rabbit mAb [JE47-88]	HA721801	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-Tau (S404) Recombinant Rabbit mAb [JE66-82]	HA721665	Human, Mouse, Rat	WB,IHC
Cell treatment	Phospho-Tau (T181) Recombinant Rabbit mAb [PSH05-42]	HA722271	Human, Mouse, Cynomolgus monkey	WB,IHC-P
Cell treatment	Phospho-Tau (T217) Recombinant Rabbit mAb [PSH09-33]	HA723091	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr
	Phospho-Tau (T231) Recombinant Rabbit mAb [PSH01-04]	HA721570	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Cell treatment	Phospho-Tau (T231) Recombinant Rabbit mAb [SC58-08]	ET1610-31	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IP,IF-Tissue,IHC-Fr
	Phospho-Tau (T231) Recombinant Rabbit mAb [PSH01-05]	HA721569	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phosphothreonine Recombinant Rabbit mAb [PSH14-34]	HA723628	Human, Mouse, Rat, Species independent	WB,IHC-P,IF-Cell,IP,ELISA
Relative expression	Phospho-TrkB (Y817) Recombinant Rabbit mAb [SC0556]	ET1610-35	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Cell treatment	Phospho-Tuberin (S939) Recombinant Rabbit mAb [JE54-06]	HA721762	Human, Mouse	WB
Cell treatment	Phospho-Tuberin (T1462) Recombinant Rabbit mAb [JE43-81]	HA721792	Human, Mouse, Rat	WB
Cell treatment	Phospho-UBF1 (S484) Recombinant Rabbit mAb [JE77-43]	HA722477	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-ULK1 (S556) Recombinant Rabbit mAb [JE59-28]	HA721445	Human, Mouse	WB
Cell treatment	Phospho-VASP (S156) Recombinant Rabbit mAb [JE45-91]	HA722145	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-Vimentin (S39) Recombinant Rabbit mAb [JE43-26]	HA721442	Human, Mouse, Rat	WB,IF-Cell
	Phospho-V-Myb+C-Myb (S11) Recombinant Rabbit mAb [SZ04-81]	ET1603-41	Human	IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-YAP1 (S127) Recombinant Rabbit mAb [SN0718]	ET1611-69	Human, Mouse, Rat	WB,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

