

2025年8月—新上市产品清单

高级验证	产品名	货号	反应种属	验证应用
	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Mouse IgG1 (Chimeric)	HA601546	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Arginase1 Recombinant Rabbit mAb	IRS081RB	Human	mIHC
	Arginase1 Recombinant Rabbit mAb	IRS082RB	Mouse	mIHC
	ATPase Inhibitory Factor 1 / IF1 Recombinant Rabbit mAb [PSH16-94]	HA723890	Human, Monkey	WB,IF-Cell,IHC-P,FC,IP
	BAHD1 Rabbit pAb	HA500589	Human, Mouse, Rat	WB,IHC-P
	BCMA Recombinant Rabbit mAb [PSH16-80]	HA723880	Human	FC
	beta Amyloid Recombinant Rabbit mAb	IRS251RB	Mouse	mIHC
	beta Amyloid 1-42 Recombinant Rabbit mAb	IRS247RB	Mouse	mIHC
	beta COP Rabbit pAb	HA500588	Human, Mouse, Rat	WB
	Cadherin-Catenin Antibody Sampler Kit	HAK21157	Human	WB,IHC-P
	Calreticulin Recombinant Antibody [SU37-03] - Mouse IgG1 (Chimeric)	HA601547	Human, Mouse, Rat, Monkey	WB,IHC-P
	Caveolin-1 Recombinant Rabbit mAb	IRS100RB	Mouse	mIHC
Relative expression	CD105 Recombinant Rabbit mAb [PSH15-96]	HA723801	Human	FC
Relative expression	CD105 Recombinant Rabbit mAb [PSH15-97]	HA723802	Human	FC
	CD11c Recombinant Antibody [SI19-06] - Mouse IgG1 (Chimeric)	HA601545	Human	IF-Tissue,IHC-P
	CD123 / IL3RA Recombinant Mouse mAb [PSH14-51]	HA601490	Human	FC
	CD163 Rabbit pAb	IRS109RB	Mouse	mIHC
	PE Conjugated CD20 Recombinant Rabbit mAb [PSH08-12]	HA723839F	Human	FC
Relative expression	CD20 Recombinant Antibody [PD00-02] - Mouse IgG1 (Chimeric)	HA601544	Human, Mouse	WB,IHC-P,IHC-Fr
	APC Conjugated CD20 Recombinant Rabbit mAb [PSH08-12]	HA723838F	Human	FC
	CD20 Rabbit pAb	IRS114RB	Mouse	mIHC
	CD25 / IL-2 Receptor alpha Recombinant Mouse mAb [PSH14-46]	HA601486	Human, Mouse	FC,IHC-P,IF-Tissue
	CD27 Recombinant Rabbit mAb	IRS115RB	Human	mIHC
	CD31 Recombinant Rabbit mAb	IRS117RB	Mouse	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD36 Recombinant Rabbit mAb	IRS118RB	Human	mIHC
	CD44 Recombinant Rabbit mAb	IRS123RB	Mouse	mIHC
	CD45RO Recombinant Mouse mAb	IRS124MS	Human	mIHC
Relative expression	CD62L Recombinant Rabbit mAb [PSH14-92]	HA723694	Human	FC,IF-Cell
	CD8 Recombinant Rabbit mAb	IRS129RB	Human	mIHC
	CD94 Recombinant Mouse mAb [PSH14-52]	HA601491	Human	FC
	PE Conjugated CD94 Recombinant Mouse mAb [PSH14-52]	HA601539F	Human	FC
	CIP2A Recombinant Rabbit mAb [PSH16-12]	HA723817	Human, Mouse, Rat	WB,IP
	Cleaved Caspase-3 Rabbit pAb	IRS227RB	Mouse	mIHC
	Cleaved Caspase-3 Recombinant Rabbit mAb	IRS136RB	Human	mIHC
	COL4A3 Rabbit pAb	HA500591	Rat	WB,IF-Cell
	CX3CR1 Recombinant Rabbit mAb	IRS250RB	Mouse	mIHC
	Cyclin D2 Recombinant Rabbit mAb [PSH16-53]	HA723857	Human, Mouse, Rat	WB,IHC-P,FC
	Dishevelled 3 Recombinant Rabbit mAb [PSH15-93]	HA723798	Human, Mouse, Rat, Monkey	WB,IP
	Doublecortin Rabbit pAb	IRS252RB	Mouse	mIHC
	EMA Mouse mAb	IRS253MS	Human	mIHC
	Endomucin Rabbit pAb	IRS141RB	Mouse	mIHC
Relative expression	eNOS Recombinant Rabbit mAb [PSH16-96]	HA723892	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	EpCAM Recombinant Rabbit mAb	IRS254RB	Human	mIHC
	FAM38A / PIEZO1 Recombinant Mouse mAb	IRS257MS	Mouse	mIHC
	FOXP3 Rabbit pAb	IRS147RB	Mouse	mIHC
	iFluor™ 750 Conjugated Goat anti-Chicken IgY Goat pAb	HA1145	Chicken	FC,IHC-Fr,IF-Tissue,IF-Cell
	iFluor™ 488 Conjugated Goat anti-Guinea pig IgG Goat pAb	HA1144	Guinea pig	IHC-Fr,IF-Tissue,IF-Cell
	Goat anti-Rat IgG(H+L) Goat pAb	HA1034	Rat	WB,ELISA
	HIV-1 lentiviral BCA protein assay kit	PQ01	Human	ELISA
	Human AGRP, Tag Free Protein	HA211302	Human	
	Human C1q R1 / CD93, Tag Free Protein	HA211294	Human	
	Human Cardiotrophin-1/CT-1, Tag Free Protein	HA211250	Human	
	Human CCL16/LEC, Tag Free Protein	HA211292	Human	

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CCL22 / MDC, C-hFc Tag Protein	HA211312	Human	
	FITC Conjugated Human CD10 Recombinant Rabbit mAb [SN75-07]	HA720209F	Human	FC
	Human CD200 enzyme-linked immunoassay kit	EH0173	Human	ELISA
	FITC Conjugated Human CD200R1 Recombinant Rabbit mAb [PSH05-02]	HA720212F	Human	FC
	Human CD30 Ligand / TNFSF8, Tag Free Protein	HA211293	Human	
	FITC Conjugated Human CD38 Recombinant Rabbit mAb [SN07-20]	HA720196F	Human	FC
	PE-Cy7 Conjugated Human CD4 Recombinant Mouse mAb [PSH04-93]	HA601542F	Human	FC
	Human CD45RA Recombinant Mouse mAb [PSH04-97]	HA601341	Human	FC,IF-Cell,IHC-P
	Human CD73 enzyme-linked immunoassay kit	EH0093	Human	ELISA
	APC-Cy7 Conjugated Human CD8 alpha Recombinant Mouse mAb [PSH04-94]	HA601541F	Human	FC
	Human CEACAM-5/CD66e, Tag Free Protein	HA211260	Human	
	Human cIAP1/BIRC2, Tag Free Protein	HA211301	Human	
	Human CSF3R/G-CSFR/CD114, C-His Tag Protein	HA211262	Human	
	Human CTLA-4 enzyme-linked immunoassay kit	EH0179	Human	ELISA
	Human CXCL16, Tag Free Protein	HA211138	Human	
	Human DR3/TNFRSF25 Recombinant Rabbit mAb [PSH16-69] - BSA and Azide free (Detector)	HA723871	Human	ELISA(Det)
	Biotin Conjugated Human DR3/TNFRSF25 Recombinant Rabbit mAb [PSH16-68] - Detector	HA723870B	Human	ELISA(Det),ELISA
	Human DR3/TNFRSF25 Recombinant Rabbit mAb [PSH16-68] - BSA and Azide free (Detector)	HA723870	Human	ELISA(Det)
	Human DR3/TNFRSF25 Recombinant Rabbit mAb [PSH16-67] - BSA and Azide free (Capture)	HA723869	Human	ELISA(Cap)
	Biotin Conjugated Human DR3/TNFRSF25 Recombinant Rabbit mAb [PSH16-69] - Detector	HA723871B	Human	ELISA(Det),ELISA
	Human FAM171A2, Tag Free Protein	HA211284	Human	
	Biotin Conjugated Human Fas / CD95 Recombinant Rabbit mAb [PSH16-38] - Detector	HA723836B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Fas / CD95 Recombinant Rabbit mAb [PSH16-38] - BSA and Azide free (Detector)	HA723836	Human	ELISA(Det)
	Human Fas / CD95 Recombinant Rabbit mAb [PSH16-39] - BSA and Azide free (Capture)	HA723837	Human	ELISA(Cap)
	Human Fas / CD95 Recombinant Rabbit mAb [PSH16-37] - BSA and Azide free (Capture)	HA723835	Human	ELISA(Cap)
	Human Galectin-9 / LGALS9, Tag Free Protein	HA211310	Human	
	Human GITR Ligand/TNFSF18, Tag Free Protein	HA211265	Human	
	Human GRO gamma/CXCL3, Tag Free Protein	HA211249	Human	
	Human IgD Recombinant Mouse mAb [PSH14-48]	HA601487	Human	FC
	Human IgE Recombinant Mouse mAb [PSH14-50]	HA601489	Human	FC
	Human IgM Recombinant Mouse mAb [PSH14-49]	HA601488	Human	FC
	Human IL-13R alpha 2 enzyme-linked immunoassay kit	EH0182	Human	ELISA
	Human IL17B, Tag Free Protein	HA211270	Human	
	Human IL-18, Tag Free (HEK293) Protein	HA211285	Human	
	Human IL-19, Tag Free Protein	HA211269	Human	
	Human IL1R1, C-His Tag Protein	HA211261	Human	
	Human IL-1RA/IL-1F3 enzyme-linked immunoassay kit	EH0177	Human	ELISA
	Human IL-20, C-hFc Tag Protein	HA211268	Human	
	Human IL-22 enzyme-linked immunoassay kit	EH0086	Human	ELISA
	Human IL-34, Tag Free Protein	HA211267	Human	
	Human IL-7 R alpha / CD127, Tag Free Protein	HA211291	Human	
	Human Interferon omega-1, Tag Free Protein	HA211271	Human	
	Human LAMC2, Tag Free Protein	HA211219	Human	
	Human LAMP3/CD208, Tag Free Protein	HA211259	Human	
	Human LDLR/LDL Receptor, Tag Free Protein	HA211287	Human	
	Biotin Conjugated Human Neurogranin Recombinant Rabbit mAb [PSH16-57] - Detector	HA725192B	Human	ELISA(Det),ELISA
	Human Neurogranin Recombinant Rabbit mAb [PSH16-57] - BSA and Azide free (Detector)	HA725192	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Neurogranin Recombinant Rabbit mAb [PSH16-56] - BSA and Azide free (Capture)	HA725191	Human	ELISA(Cap)
	Human Neuropilin-1/NRP1 enzyme-linked immunoassay kit (one step)	EHY0175	Human	ELISA
	Human NKG2D, Tag Free Protein	HA211290	Human	
	Human NKp46/NCR1, Tag Free Protein	HA211272	Human	
	Human Progranulin/GRN, Tag Free Protein	HA211273	Human	
	Human Pro-neuropeptide Y (29-64 peptide) Protein	HA211303	Human	
	Human Proprotein Convertase 9/PCSK9 enzyme-linked immunoassay kit	EH0071	Human	ELISA
	Human Proprotein Convertase 9/PCSK9 enzyme-linked immunoassay kit (one step)	EHY0071	Human	ELISA
	Human RAGE, C-His Tag Protein	HA211263	Human	
	Human RANKL Recombinant Rabbit mAb [PSH16-60] - BSA and Azide free (Capture)	HA723862	Human	ELISA(Cap)
	Human RANKL Recombinant Rabbit mAb [PSH16-58] - BSA and Azide free (Capture)	HA723860	Human	ELISA(Cap)
	Biotin Conjugated Human RANKL Recombinant Rabbit mAb [PSH16-59] - Detector	HA723861B	Human	ELISA(Det),ELISA
	Human RANKL Recombinant Rabbit mAb [PSH16-59] - BSA and Azide free (Detector)	HA723861	Human	ELISA(Det)
	Human REST, Tag Free Protein	HA211258	Human	
	Human S100B, Tag Free Protein	HA211252	Human	
	Human SCF/KITLG enzyme-linked immunoassay kit	EH0174	Human	ELISA
	Human SCF/KITLG, Tag Free (26-189) Protein	HA211286	Human	
	Human SMN1/SMN2, Tag Free Protein	HA211276	Human	
	Human Sonic Hedgehog/SHH (C24II), Tag Free Protein	HA211279	Human	
	Human Sonic Hedgehog/SHH, Tag Free Protein	HA211278	Human	
	Human TACI/TNFRSF13B, Tag Free Protein	HA211248	Human	
	Human TL1A/TNFSF15, Tag Free Protein	HA211266	Human	
	Human TRAF1, Tag Free Protein	HA211289	Human	
	Human TRAIL/TNFSF10 enzyme-linked immunoassay kit	EH0176	Human	ELISA
	IAP Family Antibody Sampler Kit	HAK21148	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	IKK Isoform Antibody Sampler Kit	HAK21155	Human	WB
	IL-6 Recombinant Rabbit mAb	IRS160RB	Human	mIHC
	IRAK Isoform Antibody Sampler Kit	HAK21147	Human	WB
	Jak Isoform Antibody Sampler Kit	HAK21158	Human	WB,IHC-P
	PE Conjugated Ly6g Recombinant Rabbit mAb [PSH14-95]	HA723853F	Mouse	FC
	Ly6g Recombinant Rabbit mAb	IRS167RB	Mouse	mIHC
	Ly6g Recombinant Rabbit mAb [PSH14-95]	HA723696	Mouse	FC,IHC-P
	LYVE1 Recombinant Rabbit mAb	IRS169RB	Mouse	mIHC
	LYVE1 Recombinant Rabbit mAb	IRS256RB	Human	mIHC
	Mammal GFAP enzyme-linked immunoassay kit	EH0178	Mammal	ELISA
	Mammal Insulin enzyme-linked immunoassay kit	EH0038	Mammal	ELISA
	mCherry Mouse mAb	IRS170MS	Species-independent	mIHC
	Mouse CCL19/MIP-3b, Tag Free Protein	HA211277	Mouse	
	Mouse CD86 enzyme-linked immunoassay kit	EM0046	Mouse	ELISA
	Mouse CD8 alpha Recombinant Rat mAb [PSH16-78]	HA601543	Mouse	FC
	Mouse CXCL12 / SDF-1, Tag Free Protein	HA211308	Mouse	
	Mouse FGF basic/FGF2/bFGF, C-His Tag Protein	HA211283	Mouse	
	Mouse Flt-3 Ligand/FLT3LG, Tag Free Protein	HA211282	Mouse	
	Mouse ICAM1 Recombinant Rabbit mAb [PSH16-34] - BSA and Azide free (Capture)	HA723832	Mouse	ELISA(Cap)
	Mouse ICAM1 Recombinant Rabbit mAb [PSH16-35] - BSA and Azide free (Detector)	HA723833	Mouse	ELISA(Det)
	Mouse ICAM1 Recombinant Rabbit mAb [PSH16-36] - BSA and Azide free (Capture)	HA723834	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse ICAM1 Recombinant Rabbit mAb [PSH16-35] - Detector	HA723833B	Mouse	ELISA(Det),ELISA
	Mouse IL-1 alpha/IL-1 F1 enzyme-linked immunoassay kit	EM0047	Mouse	ELISA
	Mouse MMP-3, Tag Free Protein	HA211280	Mouse	
	HRP Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - Detector	HA723788H	Mouse	ELISA(Det),ELISA
	Mouse Noggin/NOG, Tag Free Protein	HA211281	Mouse	
	Mouse/Rat Complement C3 Recombinant Rabbit mAb [PSH16-82] - BSA and Azide free (Detector)	HA725144	Mouse, Rat	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse/Rat Complement C3 Recombinant Rabbit mAb [PSH16-82] - Detector	HA725145B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat Complement C3 Recombinant Rabbit mAb [PSH16-81] - BSA and Azide free (Capture)	HA725143	Mouse, Rat	ELISA(Cap)
	Mouse/Rat GDF15 Recombinant Rabbit mAb [PSH16-64] - BSA and Azide free (Capture)	HA723866	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat GDF15 Recombinant Rabbit mAb [PSH16-65] - Detector	HA723867B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat GDF15 Recombinant Rabbit mAb [PSH16-65] - BSA and Azide free (Detector)	HA723867	Mouse, Rat	ELISA(Det)
	Mouse/Rat GDF15 Recombinant Rabbit mAb [PSH16-66] - BSA and Azide free (Capture)	HA723868	Mouse, Rat	ELISA(Cap)
	HRP Conjugated Mouse Resistin Recombinant Rabbit mAb [PSH16-04] - Detector	HA723809H	Mouse	ELISA(Det),ELISA
	Mouse R-spondin-1/RSPO1, Tag Free Protein	HA211306	Mouse	
	Mouse Thrombopoietin / TPO, C-His Tag Protein	HA211311	Mouse	
	Mouse VCAM1 Recombinant Rabbit mAb [PSH16-32] - BSA and Azide free (Detector)	HA723830	Mouse	ELISA(Det)
	Biotin Conjugated Mouse VCAM1 Recombinant Rabbit mAb [PSH16-32] - Detector	HA723830B	Mouse	ELISA(Det),ELISA
	Mouse VCAM1 Recombinant Rabbit mAb [PSH16-33] - BSA and Azide free (Capture)	HA723831	Mouse	ELISA(Cap)
	Mouse VCAM1 Recombinant Rabbit mAb [PSH16-31] - BSA and Azide free (Capture)	HA723829	Mouse	ELISA(Cap)
Cell treatment	Neuropeptide Y Recombinant Rabbit mAb [PSH16-52]	HA723856	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	Orexin Recombinant Rabbit mAb [PSH16-77]	HA723879	Mouse, Rat	IHC-Fr,IHC-P
	Osteopontin Recombinant Rabbit mAb	IRS200RB	Human	mIHC
	OTULIN Rabbit pAb	HA500587	Human, Mouse, Rat, Monkey	WB
	OTULIN Recombinant Rabbit mAb [PSH16-95]	HA723891	Human	WB,FC
	pan Cytokeratin Recombinant Rabbit mAb	IRS188RB	Mouse	mIHC
	PD1 Mouse mAb	IRS191MS	Mouse	mIHC
	PD1 Mouse mAb	IRS190MS	Human	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PD-L1 Recombinant Antibody - Rat IgG2a (Chimeric)	IRS193RT	Human	mIHC
	Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-75] - BSA and Azide free (Capture)	HA723877	Human	ELISA(Cap)
	PKC Isoform Antibody Sampler Kit	HAK21154	Human	WB
	Pro-Apoptosis Bcl-2 Family Antibody Sampler Kit	HAK21159	Human	WB
	Protein Folding and Stability Antibody Sampler Kit	HAK21156	Human	WB
	Biotin Conjugated P-Selectin / CD62P Recombinant Rabbit mAb [PSH16-71] - Detector	HA723873B	Human, Mouse, Rat	ELISA(Det),ELISA
	P-Selectin / CD62P Recombinant Rabbit mAb [PSH16-70] - BSA and Azide free (Capture)	HA723872	Human, Mouse, Rat	ELISA(Cap)
	P-Selectin / CD62P Recombinant Rabbit mAb [PSH16-71] - BSA and Azide free (Detector)	HA723873	Human, Mouse, Rat	ELISA(Det)
	Rab Family Antibody Sampler Kit	HAK21152	Human	WB
	RBP4 Recombinant Rabbit mAb [PSH16-97]	HA723893	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	ROR alpha / RORA Recombinant Rabbit mAb	IRS258RB	Mouse	mIHC
	S100A8+S100A9 Recombinant Rabbit Multiclonal Antibody	IRS196RB	Mouse	mIHC
	FITC Conjugated Siglec 8 Recombinant Mouse mAb [PSH14-96]	HA601517F	Human	FC
	SOX2 Rabbit pAb	IRS198RB	Human	mIHC
Cell treatment	TAP2 Recombinant Rabbit mAb [PSH16-55]	HA723859	Human, Mouse, Rat	WB,IHC-P,IP
	Biotin Conjugated Tau Recombinant Rabbit mAb [PSH17-11] - Detector	HA723898B	Human	ELISA(Det),ELISA
	Tau Recombinant Rabbit mAb [PSH17-27] - BSA and Azide free	HA723914	Human	ELISA
	Tau Recombinant Rabbit mAb [PSH17-11] - BSA and Azide free (Detector)	HA723898	Human	ELISA(Det)
	Tau Recombinant Rabbit mAb [PSH17-26] - BSA and Azide free	HA723913	Human	ELISA
	TCF7 Recombinant Rabbit mAb	IRS204RB	Human	mIHC
	TCR beta chain Recombinant Mouse mAb [PSH16-30]	HA601535	Mouse	FC
	TGF beta 1 Recombinant Rabbit mAb	IRS205RB	Human	mIHC
Relative expression	TIM 3 Recombinant Rabbit mAb [PSH16-50]	HA723854	Human, Mouse, Rat	WB,IHC-P,IP
	TROP2 Recombinant Rabbit mAb	IRS208RB	Human	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	TRPV4 Recombinant Rabbit mAb [PSH16-51]	HA723855	Human, Mouse	WB,IHC-Fr,IHC-P,IP
	Type II keratin Antibody Sampler Kit Antibody Sampler Kit	HAK21151	Human	WB
	Type I keratin Antibody Sampler Kit Antibody Sampler Kit	HAK21150	Human	WB
	Tyrosine Hydroxylase Recombinant Rabbit mAb	IRS206RB	Mouse	mIHC
	VE Cadherin Mouse Antibody	IRS212MS	Mouse	mIHC
KD	VPS35 Recombinant Rabbit mAb [PSH16-54]	HA723858	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IHC-P,IP