

主要免疫细胞分泌因子-ELISA Kits

产品名	货号	反应种属	灵敏度
Human CCL2 enzyme-linked immunoassay kit	EH0010	Human	0.04 ng/ml(50 µl) ; 0.03 ng/ml(10 µl)
Human CCL3/MIP-1 enzyme-linked immunoassay kit	EH0146	Human	0.05 pg/ml (50 µl) ; 0.52 pg/ml (10 µl)
Human CCL5/RANTES enzyme-linked immunoassay kit	EH0027	Human	0.7 pg/ml (50 µl) ; 4.65 pg/ml (10 µl)
Human CCL5/RANTES enzyme-linked immunoassay kit(one step)	EHY0027	Human	0.7 pg/ml (50 µl) ; 4.99 pg/ml (10 µl) ;
Human EGF enzyme-linked immunoassay kit	EH0023	Human	0.12 pg/ml (50 µl) ; 0.59 pg/ml (10 µl)
Human FGF-19 enzyme-linked immunoassay kit	EH0140	Human	0.48 pg/ml (50 µl) ; 4.11 pg/ml (10 µl)
Human FGF basic /FGF2/bFGF enzyme-linked immunoassay kit	EH0024	Human	1.51 pg/ml (50 µl) ; 9.93 pg/ml (10 µl)
Human GM-CSF enzyme-linked immunoassay kit	EH0051	Human	3.57 pg/ml (50 µl) ; 5.45 pg/ml (10 µl)
Human Granzyme A/GZMA enzyme-linked immunoassay kit	EH0129	Human	0.46 pg/ml (50 µl) ; 1.46 pg/ml (10 µl)
Human Granzyme B enzyme-linked immunoassay kit	EH0052	Human	1.12 pg/ml (50 µl) ; 5.96 pg/ml (10 µl)
Human IFN-gamma enzyme-linked immunoassay kit	EH0039	Human	0.36 pg/ml (50 µl) ; 0.85 pg/ml (10 µl)
Human IL-10 enzyme-linked immunoassay kit	EH0009	Human	0.89 pg/ml (10 µl)
Human IL-12B/IL-12 p40 enzyme-linked immunoassay kit	EH0098	Human	0.29 pg/ml (50 µl) ; 1.84 pg/ml (10 µl)
Human IL-12p70 enzyme-linked immunoassay kit	EH0021	Human	3.42 pg/ml (50 µl) ; 11.63 pg/ml (10 µl)
Human IL-13 enzyme-linked immunoassay kit	EH0049	Human	38.56 pg/ml (50 µl) ; 179.32 pg/ml (10 µl)
Human IL-15 enzyme-linked immunoassay kit	EH0031	Human	0.14 pg/ml (50 µl) ; 0.79 pg/ml (10 µl)
Human IL-17A enzyme-linked immunoassay kit	EH0004	Human	0.26 pg/ml (50 µl) ; 0.65 pg/ml (10 µl)
Human IL-17E/IL-25 enzyme-linked immunoassay kit	EH0066	Human	0.69 pg/ml (50 µl) ; 2.74 pg/ml (10 µl)
Human IL-18 enzyme-linked immunoassay kit	EH0053	Human	1.38 pg/ml (50 µl) ; 2.13 pg/ml (10 µl)

Human IL-1 β enzyme-linked immunoassay kit	EH0006	Human	0.18 pg/ml (10 μ l) ; 0.09 pg/ml (50 μ l) ;
Human IL-21 enzyme-linked immunoassay kit	EH0018	Human	21.7 pg/ml (50 μ l) ; 67.9 pg/ml (10 μ l)
Human IL-22 enzyme-linked immunoassay kit	EH0086	Human	0.03 pg/ml (50 μ l) ; 0.21 pg/ml (10 μ l)
Human IL-2 enzyme-linked immunoassay kit	EH0013	Human	1.67 pg/ml (50 μ l) ; 3.08 pg/ml (10 μ l)
Human IL-31 enzyme-linked immunoassay kit	EH0097	Human	0.4 pg/ml (50 μ l) ; 2.68 pg/ml (10 μ l)
Human IL-3 enzyme-linked immunoassay kit	EH0050	Human	3.74 pg/ml (50 μ l) ; 20.64 pg/ml (10 μ l)
Human IL-4 enzyme-linked immunoassay kit	EH0007	Human	0.03 pg/ml (50 μ l) ; 0.28 pg/ml (10 μ l)
Human IL-5 enzyme-linked immunoassay kit	EH0008	Human	0.15 pg/ml (50 μ l) ; 0.73 pg/ml (10 μ l)
Human IL-6 enzyme-linked immunoassay kit	EH0001	Human	0.26 pg/ml (50 μ l) ; 0.43 pg/ml (10 μ l)
Human IL-8/CXCL8 enzyme-linked immunoassay kit	EH0014	Human	1.45 pg/ml (50 μ l) ; 3.63 pg/ml (10 μ l)
Human/Mouse/Rat/Porcine/Canine TGF- β 1 enzyme-linked immunoassay kit	EH0012	Human, Mouse, Rat, Porcine, Canine	33.32 pg/ml (50 μ l)
Human Osteopontin (OPN) enzyme-linked immunoassay kit	EH0057	Human	0.04 ng/ml (50 μ l) ; 0.12 ng/ml (10 μ l)
Human PDGF-AA enzyme-linked immunoassay kit	EH0058	Human	0.10 ng/ml (50 μ l) ; 0.17 ng/ml (10 μ l)
Human PDGF-BB enzyme-linked immunoassay kit	EH0068	Human	7.56 pg/ml (50 μ l) ; 2.79 pg/ml (10 μ l)
Human TNF- α enzyme-linked immunoassay kit	EH0002	Human	0.21 pg/ml (50 μ l) ; 0.26 pg/ml (10 μ l)
Human TRAIL/TNFSF10 enzyme-linked immunoassay kit	EH0176	Human	0.90 pg/ml (50 μ l) ; 11.03 pg/ml (10 μ l) ;
Human VEGF enzyme-linked immunoassay kit	EH0075	Human	1.04 pg/ml (50 μ l) ; 3.21 pg/ml (10 μ l)
Mammal VEGF-C enzyme-linked immunoassay kit	EH0113	Mammal	1.82 pg/ml (50 μ l)
Mouse CCL2/MCP-1 enzyme-linked immunoassay kit	EM0025	Mouse	1.78 pg/ml (10 μ l)
Mouse CCL3/MIP-1 alpha enzyme-linked immunoassay kit	EM0031	Mouse	0.51 pg/ml (10 μ l)
Mouse CCL5 / RANTES enzyme-linked immunoassay kit	EM0032	Mouse	1.86 pg/ml (10 μ l)
Mouse GM-CSF/CSF2 enzyme-linked immunoassay kit	EM0011	Mouse	0.19 pg/ml (10 μ l)
Mouse IFN- γ enzyme-linked immunoassay kit	EM0009	Mouse	3.24 pg/ml (10 μ l)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Mouse IL-10 enzyme-linked immunoassay kit	EM0005	Mouse	5.75 pg/ml (10 µl)
Mouse IL-12/IL-23 p40 enzyme-linked immunoassay kit	EM0012	Mouse	1.65 pg/ml (10 µl)
Mouse IL-13 enzyme-linked immunoassay kit	EM0007	Mouse	0.37 pg/ml (10 µl)
Mouse IL-17A enzyme-linked immunoassay kit	EM0008	Mouse	0.25 pg/ml (10 µl)
Mouse IL-17E enzyme-linked immunoassay kit	EM0028	Mouse	12.33 pg/ml (10 µl)
Mouse IL-18 enzyme-linked immunoassay kit	EM0034	Mouse	1.46 pg/ml (10 µl)
Mouse IL-1 alpha/IL-1 F1 enzyme-linked immunoassay kit	EM0047	Mouse	0.05 pg/ml (10 µl)
Mouse IL-1 beta Super Sensitivity enzyme-linked immunoassay kit	EMSS0029	Mouse	0.12 pg/ml (50 µl) ; 0.36 pg/ml (10 µl)
Mouse IL-1β enzyme-linked immunoassay kit	EM0029	Mouse	0.91 pg/ml (10 µl)
Mouse IL-21 enzyme-linked immunoassay kit	EM0020	Mouse	3.22 pg/ml (10 µl)
Mouse IL-22 enzyme-linked immunoassay kit	EM0021	Mouse	11.2 pg/ml (10 µl)
Mouse IL-2 enzyme-linked immunoassay kit	EM0002	Mouse	0.16 pg/ml (10 µl)
Mouse IL-4 enzyme-linked immunoassay kit	EM0003	Mouse	0.35 pg/ml (10 µl)
Mouse IL-5 enzyme-linked immunoassay kit	EM0023	Mouse	0.24 pg/ml (10 µl)
Mouse IL-6 enzyme-linked immunoassay kit	EM0004	Mouse	0.36 pg/ml (10 µl)
Mouse IL-9 enzyme-linked immunoassay kit	EM0022	Mouse	0.10 pg/ml (10 µl)
Mouse TNF-α enzyme-linked immunoassay kit	EM0010	Mouse	0.097 pg/ml (10 µl)
Mouse VEGF enzyme-linked immunoassay kit	EM0013	Mouse	0.57 pg/ml (10 µl)
Rat IFN-γ enzyme-linked immunoassay kit	ER0006	Rat	7.39 pg/ml (50 µl)
Rat IL-10 enzyme-linked immunoassay kit	ER0005	Rat	5.76 pg/ml (50 µl)
Rat IL-2 enzyme-linked immunoassay kit	ER0002	Rat	8.13 pg/ml (50 µl)
Rat IL-4 enzyme-linked immunoassay kit	ER0003	Rat	17.29 pg/ml (50 µl)
Rat IL-6 enzyme-linked immunoassay kit	ER0004	Rat	7.05 pg/ml (50 µl)
Rat TNF-α enzyme-linked immunoassay kit	ER0007	Rat	0.08 pg/ml (50 µl) ; 0.16 pg/ml (10 µl)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层