

2025年9月—新上市产品清单

高级验证	产品名	货号	反应种属	验证应用
	3×Flag Peptide	APK051		IP Elution
	Acetylation and deacetylation dynamics Antibody Sampler Kit	HAK21160	Human	WB
Cell treatment	Acid sphingomyelinase Recombinant Rabbit mAb [PSH15-61] - BSA and Azide free	HA751572	Human, Mouse, Rat, Monkey	WB,IHC-P,IP
	AFF4 Rabbit pAb	HA500590	Human, Mouse, Rat	WB,IF-Cell,FC
Relative expression	Aldolase A/B/C Recombinant Antibody [PSH01-81] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610341	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Aldolase A/B/C Recombinant Antibody [PSH01-81] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610342	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	Alpha 1 antitrypsin Recombinant Rabbit mAb [PSH16-21] - BSA and Azide free	HA751592	Mouse, Rat	WB,IHC-P
	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610362	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	AMID Recombinant Rabbit mAb [PSH15-60] - BSA and Azide free	HA751571	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	ASPA Recombinant Antibody [PSH08-04] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610328	Human, Mouse, Rat	WB,IHC-P
Relative expression, Cell treatment	ATF3 Recombinant Rabbit mAb [PSH17-99]	HA723964	Human, Mouse	WB
	ATP5O Recombinant Rabbit mAb [PSH15-04] - BSA and Azide free	HA751542	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,FC,IP
	ATPase Inhibitory Factor 1 / IF1 Recombinant Rabbit mAb [PSH16-94] - BSA and Azide free	HA751608	Human, Monkey	WB,IF-Cell,IHC-P,FC,IP
	BCKDHA Recombinant Rabbit mAb [PSH15-62] - BSA and Azide free	HA751573	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	Bcl2-L-13 Recombinant Rabbit mAb [PSH17-12] - BSA and Azide free	HA751616	Human, Mouse, Rat, Monkey	WB,IHC-P,IP
	Bcl2-L-13 Recombinant Rabbit mAb [PSH17-12]	HA723899	Human, Mouse, Rat, Monkey	WB,IHC-P,IP
	BCMA Recombinant Rabbit mAb [PSH16-80] - BSA and Azide free	HA751607	Human	FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	BCMA Recombinant Rabbit Multiclonal Antibody [PSH18-02]	HA723967	Human	IF-Cell,FC
	BDH1 Recombinant Rabbit mAb [PSH14-69] - BSA and Azide free	HA751532	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P
	beta Amyloid 1-40 Protein	HA210498	Human	
	beta Amyloid 1-42 Protein	HA210499	Human	
	Beta Catenin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS278RT	Human	mIHC
	Beta Catenin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS277MS	Human	mIHC
Relative expression	Beta III Tubulin Recombinant Rabbit mAb [PSH15-32] - BSA and Azide free	HA751554	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP
Relative expression	beta IV Tubulin Recombinant Rabbit mAb [PSH17-78]	HA723947	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	Binding Control Magnetic beads	100-100-200	Species independent	IP,ChIP,MS,Purification
	Binding Control Nanoselector Agarose	100-100-100	Species independent	IP,ChIP,MS,Purification
	Biotin Nanoselector Magnetic beads	078-101-003	Species independent	IP,ChIP,MS
	BRD9 Recombinant Rabbit mAb [PSH15-95] - BSA and Azide free	HA751586	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,IP
	CAB39 Recombinant Rabbit mAb [PSH15-16] - BSA and Azide free	HA751545	Human, Mouse, Rat, Monkey	WB,IHC-P,IP
	Calreticulin Recombinant Antibody [SU37-03] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610363	Human, Mouse, Rat, Monkey	WB,IHC-P
	CBL Recombinant Rabbit mAb [PSH17-61]	HA723936	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	CCL2 / MCP-1 Recombinant Rabbit mAb [PSH12-44] - BSA and Azide free	HA751533	Mouse	WB,IF-Cell,FC
Relative expression	CD105 Recombinant Rabbit mAb [PSH15-97] - BSA and Azide free	HA751588	Human	FC
Relative expression	CD105 Recombinant Rabbit mAb [PSH15-96] - BSA and Azide free	HA751587	Human	FC
	CD11c Recombinant Antibody [SI19-06] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610361	Human	IF-Tissue,IHC-P
	CD163 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS276RT	Human	mIHC
	CD163 Recombinant Antibody [JA51-30] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610364	Human	IHC-P
	CD163 Recombinant Antibody [JA51-30] - Rat IgG1 (Chimeric)	HA601549	Human	IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	CD20 Recombinant Antibody [PD00-02] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610360	Human, Mouse	WB,IHC-P,IHC-Fr
Relative expression	CD200 / OX2 Recombinant Rabbit mAb [PSH17-62]	HA723937	Mouse	IHC-Fr,IHC-P,WB
	CD21 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS265RT	Human	mIHC
Relative expression	CD31 Recombinant Antibody [7-A1-R] - Rabbit IgG (Chimeric)	HA723897	Human	WB,IF-Cell,IHC-P,FC
Relative expression	CD31 Recombinant Antibody [7-A1-R] - Rabbit IgG (Chimeric) - BSA and Azide free	HA751615	Human	WB,IF-Cell,IHC-P,FC
	CD31 Recombinant Antibody [7-A1-R] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610349	Human	IHC-P
	CD38 Recombinant Antibody [PD01-49] - Mouse IgG1 (Chimeric)	HA601554	Human	WB,IHC-P,FC
	CD38 Recombinant Antibody [PD01-49] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610369	Human	WB,IHC-P,FC
	CD4 Recombinant Antibody [ST0488] - Mouse IgG1 (Chimeric)	HA601557	Human	WB,IF-Cell,IHC-P
	CD4 Recombinant Antibody [PSH08-22] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610354	Mouse	WB,IF-Tissue,IHC-P,IHC-Fr
Relative expression	CD4 Recombinant Antibody [PSH08-22] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610353	Mouse	WB,IF-Tissue,IHC-P,IF-Cell
Relative expression	CD62L Recombinant Rabbit mAb [PSH15-59] - BSA and Azide free	HA751570	Human	FC
Relative expression	CD62L Recombinant Rabbit mAb [PSH14-92] - BSA and Azide free	HA751536	Human	FC,IF-Cell
	CD68 Recombinant Antibody [PDM0-13] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610347	Human	IF-Tissue,IHC-P
	CD68 Recombinant Antibody [PDM0-13] - Rabbit IgG (Chimeric) - BSA and Azide free	HA751544	Human	IF-Tissue,IHC-P
	CD68 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS264RT	Human	mIHC
	CD68 Recombinant Antibody - Rabbit IgG (Chimeric)	IRS266RB	Human	mIHC
	CGRP Recombinant Antibody [PSH08-45] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610358	Mouse	IHC-Fr
	Choline Acetyltransferase Recombinant Rabbit mAb [PSH15-49] - BSA and Azide free	HA751564	Human, Mouse, Rat	IHC-P
	Choline Acetyltransferase Recombinant Rabbit mAb [PSH15-05] - BSA and Azide free	HA751543	Human, Mouse, Rat	WB,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CIP2A Recombinant Rabbit mAb [PSH16-12] - BSA and Azide free	HA751590	Human, Mouse, Rat	WB,IP
Relative expression	c-Kit Recombinant Rabbit mAb [PSH17-79]	HA723948	Mouse	IHC-Fr,IF-Cell,FC
	CMTM6 Recombinant Rabbit mAb [PSH15-02] - BSA and Azide free	HA751540	Human, Monkey	WB,IP
	c-Myc Peptide	APK052		IP Elution
	CREB Recombinant Rabbit mAb [PSH15-48] - BSA and Azide free	HA751563	Human, Mouse, Rat, Monkey	WB,IHC-P,IHC-Fr,IF-Cell,FC
Relative expression	CTGF Recombinant Rabbit mAb [PSH15-64] - BSA and Azide free	HA751575	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Cyclin D2 Recombinant Rabbit mAb [PSH16-53] - BSA and Azide free	HA751602	Human, Mouse, Rat	WB,IHC-P,FC
	Cytokeratin 19 Recombinant Antibody [SA30-06] - Rat IgG1 (Chimeric)	HA601555	Human, Rat	WB,IF-Cell,IHC-P,FC
	Cytokeratin 19 Recombinant Antibody [SA30-06] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610370	Human, Rat	WB,IF-Cell,IHC-P,FC
	Dishevelled 3 Recombinant Rabbit mAb [PSH15-93] - BSA and Azide free	HA751584	Human, Mouse, Rat, Monkey	WB,IP
Relative expression	Doublecortin Recombinant Antibody [JJ0959] - Chicken IgY (Chimeric) - BSA and Azide free	HA610325	Human, Mouse, Rat	IHC-Fr,WB
Relative expression	Doublecortin Recombinant Antibody [JJ0959] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610327	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Doublecortin Recombinant Antibody [JJ0959] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610326	Human, Mouse, Rat	IHC-Fr,WB
	DYKDDDDK tag Nanoselector Agarose	016-101-002	Species independent	IP,ChIP,MS,Purification
	DYKDDDDK tag Nanoselector Magnetic beads	016-101-003	Species independent	IP,ChIP,MS,Purification
	EAAT2 Recombinant Antibody [PS01-62] - Chicken IgY (Chimeric) - BSA and Azide free	HA610333	Mouse	IHC-Fr
Relative expression	EAAT2 Recombinant Antibody [PS01-62] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610334	Mouse, Rat	IHC-Fr,WB
	E-Cadherin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS284RT	Mouse	mIHC
	Elastin Recombinant Rabbit mAb [PSH16-13] - BSA and Azide free	HA751591	Human	IHC-P,IF-Tissue
	ELKS Recombinant Rabbit mAb [PSH17-33] - BSA and Azide free	HA751621	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	ELKS Recombinant Rabbit mAb [PSH17-33]	HA723916	Human, Mouse, Rat	WB
	Endothelin 1 Recombinant Rabbit mAb [PSH15-63] - BSA and Azide free	HA751574	Human, Mouse, Rat	WB
Relative expression	eNOS Recombinant Rabbit mAb [PSH16-96] - BSA and Azide free	HA751610	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
Relative expression	Eph receptor B2 Recombinant Rabbit mAb [PSH14-91] - BSA and Azide free	HA751535	Human	WB,IF-Cell,IP
	Estriol Recombinant Rabbit mAb [PSH15-86] - BSA and Azide free	HA751581	Species independent	ELISA
Relative expression	F4/80 Recombinant Antibody [PSH01-87] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610356	Mouse	WB,IHC-P
	F4/80 Recombinant Antibody [PSH01-87] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610357	Mouse	WB,IHC-Fr,IHC-P
	FAP Recombinant Rabbit mAb	IRS289RB	Human	mIHC
	Frataxin Recombinant Rabbit mAb [PSH15-18] - BSA and Azide free	HA751547	Human, Rat	WB,IHC-P
	GABRA2 Recombinant Rabbit mAb [PSH17-28] - BSA and Azide free	HA751620	Mouse, Rat	IHC-Fr,IHC-P,WB,IP
	GABRA2 Recombinant Rabbit mAb [PSH17-28]	HA723915	Mouse, Rat	IHC-Fr,IHC-P,WB,IP
	Gankyrin Recombinant Rabbit mAb [PSH17-80]	HA723949	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IF-Tissue
	GFP Nanoselector Magnetic beads	019-101-003	Species independent	IP,ChIP,MS,Purification
	Glucose Transporter GLUT1 Recombinant Antibody [SA0377] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610348	Human, Mouse, Rat	IHC-Fr,IHC-P,IF-Tissue,WB
	Glutamine Synthetase Recombinant Antibody [PSH10-48] - Chicken IgY (Chimeric) - BSA and Azide free	HA610346	Mouse, Rat	IHC-Fr,WB
Relative expression	Glutathione Peroxidase 3 / GPx-3 Recombinant Rabbit mAb [PSH17-06] - BSA and Azide free	HA751612	Mouse, Rat	WB,IHC-P
Relative expression	Glutathione Peroxidase 3 / GPx-3 Recombinant Rabbit mAb [PSH17-06]	HA723894	Mouse, Rat	WB,IHC-P
	Glycophorin A / CD235a Recombinant Mouse mAb [PSH17-08] - BSA and Azide free	HA610368	Human	FC
	Glycophorin A / CD235a Recombinant Mouse mAb [PSH17-08]	HA601553	Human	FC
	HRP Conjugated Goat anti-Chicken IgY Goat pAb	HA1146	Chicken	WB,ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	iFluor™ 680se Conjugated Goat anti-mouse IgG polyclonal Antibody(For Odyssey system)	HA1130	Mouse	WB
Relative expression	Granzyme K Recombinant Rabbit mAb [PSH16-28] - BSA and Azide free	HA751596	Human	WB,IHC-P
	GST Nanoselector Agarose	010-101-002	Species independent	IP,ChIP,MS,Purification
	GST Nanoselector Magnetic beads	010-101-003	Species independent	IP,ChIP,MS,Purification
	Halo Nanoselector Agarose	012-101-002	Species independent	IP,ChIP,MS,Purification
	HA Peptide	APK053		IP Elution
	HA tag Nanoselector Magnetic beads	003-101-003	Species independent	IP,ChIP,MS,Purification
	HERC6 Rabbit pAb	HA500594	Human, Mouse, Rat	WB
	H-FABP Recombinant Rabbit mAb [PSH16-98] - BSA and Azide free (Capture)	HA725177	Human, Mouse	ELISA(Cap)
	H-FABP Recombinant Rabbit mAb [PSH16-99] - BSA and Azide free (Detector)	HA725178	Human, Mouse	ELISA(Det)
	Biotin Conjugated H-FABP Recombinant Rabbit mAb [PSH16-99] - Detector	HA725178B	Human, Mouse	ELISA(Det),ELISA
	HGF Recombinant Rabbit mAb [PSH03-79]	HA723896	Human	WB
	HGF Recombinant Rabbit mAb [PSH03-79] - BSA and Azide free	HA751614	Human	WB
	His tag Nanoselector Agarose	004-101-002	Species independent	IP,ChIP,MS,Purification
	His tag Nanoselector Magnetic beads	004-101-003	Species independent	IP,ChIP,MS,Purification
	Histone H3 (tri methyl K36) Recombinant Rabbit mAb [PSH15-01] - BSA and Azide free	HA751539	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,IP,Dot Blot,ChIP
	Biotin Conjugated Human ACE2 Recombinant Rabbit mAb [PSH17-18] - Detector	HA723905B	Human	ELISA(Det),ELISA
	Human ACE2 Recombinant Rabbit mAb [PSH17-17] - BSA and Azide free (Capture)	HA723904	Human	ELISA(Cap)
	Human ACE2 Recombinant Rabbit mAb [PSH17-18] - BSA and Azide free (Detector)	HA723905	Human	ELISA(Det)
	Human Acetylcholinesterase Recombinant Rabbit mAb [PSH17-02] - BSA and Azide free (Capture)	HA725223	Human	ELISA(Cap)
	Biotin Conjugated Human Acetylcholinesterase Recombinant Rabbit mAb [PSH17-03] - Detector	HA725224B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Acetylcholinesterase Recombinant Rabbit mAb [PSH17-03] - BSA and Azide free (Detector)	HA725224	Human	ELISA(Det)
	Biotin Conjugated Human ANGPTL3 Recombinant Rabbit mAb [PSH17-40] - Detector	HA725081B	Human	ELISA(Det),ELISA
	Human ANGPTL3 Recombinant Rabbit mAb [PSH17-40] - BSA and Azide free (Detector)	HA725080	Human	ELISA(Det)
	Human ANGPTL3 Recombinant Rabbit mAb [PSH17-39] - BSA and Azide free (Capture)	HA725079	Human	ELISA(Cap)
	Human APOC3, Tag Free Protein	HA211274	Human	
	Human CA125/MUC16, C-His Tag Protein	HA211195	Human	
	Human Caspase-9 Recombinant Rabbit mAb [PSH17-19] - BSA and Azide free (Capture)	HA723906	Human	ELISA(Cap)
	Human Caspase-9 Recombinant Rabbit mAb [PSH17-20] - BSA and Azide free (Detector)	HA723907	Human	ELISA(Det)
	Biotin Conjugated Human Caspase-9 Recombinant Rabbit mAb [PSH17-20] - Detector	HA723907B	Human	ELISA(Det),ELISA
	Human CCL21 enzyme-linked immunoassay kit	EH0184	Human	ELISA
	Human CD163 enzyme-linked immunoassay kit	EH0185	Human	ELISA
	Human CD28 enzyme-linked immunoassay kit	EH0183	Human	ELISA
	Human CFH, Tag Free Protein	HA211200	Human	
	Human Coagulation factor VIII, Tag Free Protein	HA211190	Human	
	Human Corticoliberin/CRH, C-hFc Tag Protein	HA211316	Human	
	Human CXADR Recombinant Rabbit mAb [PSH17-47] - BSA and Azide free (Detector)	HA725294	Human	ELISA(Det)
	Biotin Conjugated Human CXADR Recombinant Rabbit mAb [PSH17-47] - Detector	HA725294B	Human	ELISA(Det),ELISA
	Human CXADR Recombinant Rabbit mAb [PSH17-46] - BSA and Azide free (Capture)	HA725293	Human	ELISA(Cap)
	Human CXCL2/MIP-2 Recombinant Rabbit mAb [PSH17-88] - BSA and Azide free (Capture)	HA725247	Human	ELISA(Cap)
	Biotin Conjugated Human CXCL2/MIP-2 Recombinant Rabbit mAb [PSH17-89] - Detector	HA725248B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CXCL2/MIP-2 Recombinant Rabbit mAb [PSH17-89] - BSA and Azide free (Detector)	HA725248	Human	ELISA(Det)
	Human DR3/TNFRSF25, Tag Free Protein	HA211322	Human	
	HRP Conjugated Human Endostatin Recombinant Rabbit mAb [PSH15-83] - Detector	HA723791H	Human	ELISA(Det),ELISA
	Human GDNF, Tag Free Protein	HA211323	Human	
	Human HBEGF, Tag Free Protein	HA211320	Human	
	Human ICAM-2/CD102, C-His Tag Protein	HA211314	Human	
	Biotin Conjugated Human IFN beta Recombinant Rabbit mAb [PSH17-30] - Detector	HA725054B	Human	ELISA(Det),ELISA
	Human IFN beta Recombinant Rabbit mAb [PSH17-29] - BSA and Azide free (Capture)	HA725052	Human	ELISA(Cap)
	Human IFN beta Recombinant Rabbit mAb [PSH17-30] - BSA and Azide free (Detector)	HA725053	Human	ELISA(Det)
	Biotin Conjugated Human IGFBP1 Recombinant Rabbit mAb [PSH17-57] - Detector	HA723932B	Human	ELISA(Det),ELISA
	Human IGFBP1 Recombinant Rabbit mAb [PSH17-56] - BSA and Azide free (Capture)	HA723931	Human	ELISA(Cap)
	Human IGFBP1 Recombinant Rabbit mAb [PSH17-58] - BSA and Azide free (Detector)	HA723933	Human	ELISA(Det)
	Human IGFBP1 Recombinant Rabbit mAb [PSH17-57] - BSA and Azide free (Detector)	HA723932	Human	ELISA(Det)
	Biotin Conjugated Human IGFBP1 Recombinant Rabbit mAb [PSH17-58] - Detector	HA723933B	Human	ELISA(Det),ELISA
	Human IGF-II/IGF2, Tag Free Protein	HA211319	Human	
	Human IgG Nanoselector Magnetic beads	023-101-003	Species independent	IP,ChIP,Purification
	Human IL-10RB Recombinant Rabbit mAb [PSH17-45] - BSA and Azide free (Detector)	HA725296	Human	ELISA(Det)
	Human IL-10RB Recombinant Rabbit mAb [PSH17-44] - BSA and Azide free (Capture)	HA725295	Human	ELISA(Cap)
	Biotin Conjugated Human IL-10RB Recombinant Rabbit mAb [PSH17-45] - Detector	HA725296B	Human	ELISA(Det),ELISA
	Human IL-17C, C-hFc Tag Protein	HA211315	Human	
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - BSA and Azide free (Detector)	HA723882	Human	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-85] - BSA and Azide free (Capture)	HA723881	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - Detector	HA723882B	Human	ELISA(Det),ELISA
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-87] - BSA and Azide free (Capture)	HA723883	Human	ELISA(Cap)
	Human LAP3, C-His Tag (Isoform 2) Protein	HA211317	Human	
	Human LDHB, C-His Tag Protein	HA211198	Human	
	Human Livin/BIRC7, C-His Tag Protein	HA211325	Human	
	Human Mucin-6/MUC6, Tag Free Protein	HA211305	Human	
	HRP Conjugated Human Myeloperoxidase Recombinant Rabbit mAb [PSH14-79] - Detector	HA723675H	Human	ELISA(Det),ELISA
	Human NCR1/NKp46 Recombinant Rabbit mAb [PSH17-31] - BSA and Azide free (Capture)	HA725321	Human	ELISA(Cap)
	Human NCR1/NKp46 Recombinant Rabbit mAb [PSH17-32] - BSA and Azide free (Detector)	HA725322	Human	ELISA(Det)
	Biotin Conjugated Human NCR1/NKp46 Recombinant Rabbit mAb [PSH17-32] - Detector	HA725322B	Human	ELISA(Det),ELISA
	Human NKG2D Recombinant Rabbit mAb [PSH17-04] - BSA and Azide free (Capture)	HA725299	Human	ELISA(Cap)
	Human NKG2D Recombinant Rabbit mAb [PSH17-05] - BSA and Azide free (Detector)	HA725300	Human	ELISA(Det)
	Biotin Conjugated Human NKG2D Recombinant Rabbit mAb [PSH17-05] - Detector	HA725300B	Human	ELISA(Det),ELISA
	HRP Conjugated Human PARK7 Recombinant Mouse mAb [PSH15-68] - Detector	HA620002H	Human	ELISA(Det),ELISA
	Human Peroxiredoxin 6/PRDX6, C-His Tag Protein	HA211191	Human	
	Biotin Conjugated Human Prealbumin Recombinant Rabbit mAb [PSH17-43] - Detector	HA725260B	Human	ELISA(Det),ELISA
	Human Prealbumin Recombinant Rabbit mAb [PSH17-42] - BSA and Azide free (Capture)	HA725259	Human	ELISA(Cap)
	Human Prealbumin Recombinant Rabbit mAb [PSH17-43] - BSA and Azide free (Detector)	HA725260	Human	ELISA(Det)
	HRP Conjugated Human Prealbumin Recombinant Rabbit mAb [PSH17-43] - Detector	HA725260H	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human Renin Recombinant Rabbit mAb [PSH17-01] - Detector	HA725170B	Human	ELISA(Det),ELISA
	Human Renin Recombinant Rabbit mAb [PSH17-01] - BSA and Azide free (Detector)	HA725170	Human	ELISA(Det)
	Human Renin Recombinant Rabbit mAb [PSH17-00] - BSA and Azide free (Capture)	HA725169	Human	ELISA(Cap)
	Human S100A12 Recombinant Rabbit mAb [PSH17-68] - BSA and Azide free (Detector)	HA725216	Human	ELISA(Det)
	Human S100A12 Recombinant Rabbit mAb [PSH17-67] - BSA and Azide free (Capture)	HA725215	Human	ELISA(Cap)
	Biotin Conjugated Human S100A12 Recombinant Rabbit mAb [PSH17-68] - Detector	HA725216B	Human	ELISA(Det),ELISA
	HRP Conjugated Human S100A12 Recombinant Rabbit mAb [PSH17-68] - Detector	HA725216H	Human	ELISA(Det),ELISA
	Human SIRP alpha Recombinant Rabbit mAb [PSH17-63] - BSA and Azide free (Capture)	HA725347	Human	ELISA(Cap)
	Human SIRP alpha Recombinant Rabbit mAb [PSH17-64] - BSA and Azide free (Detector)	HA725348	Human	ELISA(Det)
	Biotin Conjugated Human SIRP alpha Recombinant Rabbit mAb [PSH17-64] - Detector	HA725348B	Human	ELISA(Det),ELISA
	Human TFAF5/FAM19A5, Tag Free Protein	HA211192	Human	
	Human TAG1/Contactin 2/CNTN2, Tag Free Protein	HA211194	Human	
	Human TGF-alpha, Tag Free Protein	HA211318	Human	
	Human TIM-1/KIM-1/HAVCR , C-Flag Tag Protein	HA211313	Human	
	Human TNF beta/Lymphotoxin alpha enzyme-linked immunoassay kit	EH0181	Human	ELISA
	Human TREM1 Recombinant Rabbit mAb [PSH17-86] - BSA and Azide free (Detector)	HA725234	Human	ELISA(Det)
	Human TREM1 Recombinant Rabbit mAb [PSH17-85] - BSA and Azide free (Capture)	HA725233	Human	ELISA(Cap)
	Biotin Conjugated Human TREM1 Recombinant Rabbit mAb [PSH17-86] - Detector	HA725234B	Human	ELISA(Det),ELISA
	ICAM1 Recombinant Rabbit mAb [PSH17-38]	HA723921	Mouse	IHC-Fr,IHC-P,WB,IF-Cell,FC,IP
	ICAM1 Recombinant Rabbit mAb [PSH17-38] - BSA and Azide free	HA751626	Mouse	IHC-Fr,IHC-P,WB,IF-Cell,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	ICAM1 Recombinant Rabbit mAb [PSH16-76] - BSA and Azide free	HA751605	Mouse, Rat	WB,IHC-P
Cell treatment	ICAM1 Recombinant Rabbit mAb [PSH16-76]	HA723878	Mouse, Rat	WB,IHC-P
	IL-1 alpha Recombinant Rabbit mAb [PSH15-66] - BSA and Azide free	HA751577	Mouse	IF-Cell,FC
Cell treatment	IL-1 alpha Recombinant Rabbit mAb [PSH15-65] - BSA and Azide free	HA751576	Mouse	WB
	IL-1 beta Recombinant Rabbit mAb [PSH18-00]	HA723965	Mouse, Rat	WB
Cell treatment	IL-27A p28 Recombinant Rabbit mAb [PSH14-93] - BSA and Azide free	HA751537	Mouse	WB
Relative expression, Cell treatment	IL-33 Recombinant Rabbit mAb [PSH17-37]	HA723920	Mouse, Rat	WB,IHC-P
Relative expression, Cell treatment	IL-33 Recombinant Rabbit mAb [PSH17-37] - BSA and Azide free	HA751625	Mouse, Rat	WB,IHC-P
	intestinal alkaline phosphatase Recombinant Rabbit mAb [PSH17-83]	HA723952	Human, Mouse, Rat	IHC-P
	intestinal alkaline phosphatase Recombinant Rabbit mAb [PSH17-82]	HA723951	Human, Mouse	WB
KD	JMJD6 Recombinant Rabbit mAb [PSH15-03] - BSA and Azide free	HA751541	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP,ChIP
	KAT3B / p300 Recombinant Rabbit mAb [PSH15-31] - BSA and Azide free	HA751553	Human, Mouse, Rat, Monkey	WB,IF-Cell
	KAT3B / p300 Recombinant Rabbit mAb [PSH15-30] - BSA and Azide free	HA751552	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	KLF5 Recombinant Rabbit mAb [PSH15-42] - BSA and Azide free	HA751557	Human, Mouse, Rat	IF-Cell,FC
Cell treatment	LC3B Recombinant Antibody [JJ090-6] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610352	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	LIPT2 Recombinant Rabbit mAb [PSH15-52] - BSA and Azide free	HA751565	Human	WB,IHC-P
	LOX Recombinant Rabbit mAb [PSH15-41] - BSA and Azide free	HA751556	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Ly6g Recombinant Rabbit mAb [PSH14-95] - BSA and Azide free	HA751538	Mouse	FC,IHC-P
Relative expression	MAG Recombinant Antibody [PSH02-41] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610340	Mouse, Rat	IHC-Fr,WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	MAG Recombinant Antibody [PSH02-41] - Chicken IgY (Chimeric) - BSA and Azide free	HA610339	Mouse, Rat	IHC-Fr,WB
	MAGED2 Recombinant Rabbit mAb [PSH15-53] - BSA and Azide free	HA751566	Human	WB
	MBP Nanoselector Magnetic beads	015-101-003	Species independent	IP,ChIP,MS,Purification
	MCM2 Recombinant Antibody [PD00-90] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610350	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P
	MESA Human Immune Checkpoint–Antigen Presentation Panel	IRSK6005		mIHC,CmTSA
	MESA Human Immune Checkpoint–Proliferation Panel	IRSK6004		mIHC,CmTSA
	MESA Human Macrophage-Dendritic-Epithelial Panel	IRSK3002		mIHC,CmTSA
	MESA Human Macrophage-Tcell-Dendritic Panel	IRSK3003		mIHC,CmTSA
	MESA Human Pan-Immune Microenvironment 21 Marker Panel	IRSK21001		mIHC,CmTSA
	MESA Human PD1–PDL1 Axis Panel	IRSK6003		mIHC,CmTSA
	MESA Human PD-L1 Cytotoxic-Epithelial Panel	IRSK3004		mIHC,CmTSA
	MESA Human Stromal-Epithelial-Myofibroblast Panel	IRSK3008		mIHC,CmTSA
	MESA Human T/B-Epithelial Panel1	IRSK3005		mIHC,CmTSA
	MESA Human T/B-Epithelial Panel2	IRSK3006		mIHC,CmTSA
	MESA Human T/B Lymphocyte Panel	IRSK3001		mIHC,CmTSA
	MESA Human T/B Lymphocyte Panel	IRSK6002		mIHC,CmTSA
	MESA Human T/B Lymphocyte Panel	IRSK6001		mIHC,CmTSA
	MESA Human Vascular-Bcell-Epithelial Panel	IRSK3009		mIHC,CmTSA
	MESA Human Vascular-Macrophage Panel	IRSK3007		mIHC,CmTSA
	mNeongreen Nanoselector Agarose	013-101-002	Species independent	IP,ChIP,MS,Purification
	MOB1 Recombinant Rabbit mAb [PSH17-36]	HA723919	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP
	MOB1 Recombinant Rabbit mAb [PSH17-36] - BSA and Azide free	HA751624	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP
	HRP Conjugated Mouse Angiotensinogen Recombinant Rabbit mAb [PSH17-10] - Detector	HA725133H	Mouse	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse Angiotensinogen Recombinant Rabbit mAb [PSH17-10] - Detector	HA725133B	Mouse	ELISA(Det),ELISA
	Mouse Angiotensinogen/AGT Recombinant Rabbit mAb [PSH17-10] - BSA and Azide free (Detector)	HA725132	Mouse	ELISA(Det)
	Mouse Angiotensinogen/AGT Recombinant Rabbit mAb [PSH17-09] - BSA and Azide free (Capture)	HA725131	Mouse	ELISA(Cap)
	Mouse Apolipoprotein B/Apo B Recombinant Rabbit mAb [PSH17-77] - BSA and Azide free (Capture)	HA723946	Mouse	ELISA(Cap)
	Mouse Apolipoprotein B/Apo B Recombinant Rabbit mAb [PSH17-76] - BSA and Azide free (Detector)	HA723945	Mouse	ELISA(Det)
	Biotin Conjugated Mouse Apolipoprotein B/Apo B Recombinant Rabbit mAb [PSH17-76] - Detector	HA723945B	Mouse	ELISA(Det),ELISA
	Mouse Apolipoprotein B/Apo B Recombinant Rabbit mAb [PSH17-75] - BSA and Azide free (Capture)	HA723944	Mouse	ELISA(Cap)
	Mouse Apolipoprotein E Recombinant Rabbit mAb [PSH16-84] - BSA and Azide free (Detector)	HA725222	Mouse	ELISA(Det)
	Mouse Apolipoprotein E Recombinant Rabbit mAb [PSH16-83] - BSA and Azide free (Capture)	HA725221	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse Apolipoprotein E Recombinant Rabbit mAb [PSH16-84] - Detector	HA725222B	Mouse	ELISA(Det),ELISA
	Mouse Cadherin-3/CDH3, Tag Free Protein	HA211188	Mouse	
	Mouse CD30 Recombinant Rabbit mAb [PSH17-55] - BSA and Azide free (Detector)	HA723930	Mouse	ELISA(Det)
	Mouse CD30 Recombinant Rabbit mAb [PSH17-54] - BSA and Azide free (Capture)	HA723929	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CD30 Recombinant Rabbit mAb [PSH17-55] - Detector	HA723930B	Mouse	ELISA(Det),ELISA
	Anti-Mouse CD8 alpha Recombinant Antibody [PSH05-00] - Rabbit IgG (Chimeric) - BSA and Azide free	HA751593	Mouse	FC,IHC-Fr
	Mouse CD8 alpha Recombinant Rat mAb [PSH16-78] - BSA and Azide free	HA610359	Mouse	FC
	Biotin Conjugated Mouse CXCL16 Recombinant Rabbit mAb [PSH17-72] - Detector	HA723941B	Mouse	ELISA(Det),ELISA
	Mouse CXCL16 Recombinant Rabbit mAb [PSH17-71] - BSA and Azide free (Capture)	HA723940	Mouse	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse CXCL16 Recombinant Rabbit mAb [PSH17-72] - BSA and Azide free (Detector)	HA723941	Mouse	ELISA(Det)
	Mouse CXCL16 Recombinant Rabbit mAb [PSH17-73] - BSA and Azide free (Capture)	HA723942	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CXCL16 Recombinant Rabbit mAb [PSH17-74] - Detector	HA723943B	Mouse	ELISA(Det),ELISA
	Mouse CXCL16 Recombinant Rabbit mAb [PSH17-74] - BSA and Azide free (Detector)	HA723943	Mouse	ELISA(Det)
	Mouse EPO/Erythropoietin Recombinant Rabbit mAb [PSH17-23] - BSA and Azide free (Capture)	HA723910	Mouse	ELISA(Cap)
	Mouse EPO/Erythropoietin Recombinant Rabbit mAb [PSH17-22] - BSA and Azide free (Detector)	HA723909	Mouse	ELISA(Det)
	Mouse EPO/Erythropoietin Recombinant Rabbit mAb [PSH17-21] - BSA and Azide free (Capture)	HA723908	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse EPO/Erythropoietin Recombinant Rabbit mAb [PSH17-22] - Detector	HA723909B	Mouse	ELISA(Det)
	Mouse gp130/IL6ST, C-His Tag Protein	HA211196	Mouse	
	HRP Conjugated Mouse ICAM1 Recombinant Rabbit mAb [PSH16-35] - Detector	HA723833H	Mouse	ELISA(Det),ELISA
	Mouse IGFBP-2, Tag Free Protein	HA211197	Mouse	
	Mouse IGFBP-4, Tag Free Protein	HA211240	Mouse	
	Mouse IGFBP5, Tag Free Protein	HA211201	Mouse	
	Mouse IgG Nanoselector Agarose	001-101-002	Species independent	IP,ChIP,MS,Purification
	Mouse IgG Nanoselector Magnetic beads	001-101-003	Species independent	IP,ChIP,Purification
	Mouse IgG&Rabbit IgG Nanoselector Agarose	001-300-002	Species independent	IP,ChIP,Purification
	Mouse IgG&Rabbit IgG Nanoselector Magnetic beads	001-300-003	Species independent	IP,ChIP,Purification
	Mouse IL-15, Tag Free Protein	HA211321	Mouse	
	Biotin Conjugated Mouse IL-17F Recombinant Rabbit mAb [PSH17-15] - Detector	HA723902B	Mouse	ELISA(Det)
	Mouse IL-17F Recombinant Rabbit mAb [PSH17-15] - BSA and Azide free (Capture/Detector)	HA723902	Mouse	ELISA(Det),ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse IL-17F Recombinant Rabbit mAb [PSH17-16] - Detector	HA723903B	Mouse	ELISA(Det)
	Mouse IL-17F Recombinant Rabbit mAb [PSH17-14] - BSA and Azide free (Capture)	HA723901	Mouse	ELISA(Cap)
	Mouse IL-17F Recombinant Rabbit mAb [PSH17-16] - BSA and Azide free (Detector)	HA723903	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-33 Recombinant Rabbit mAb [PSH16-62] - Detector	HA723864B	Mouse	ELISA(Det),ELISA
	Mouse IL-33 Recombinant Rabbit mAb [PSH16-63] - BSA and Azide free (Capture)	HA723865	Mouse	ELISA(Cap)
	Mouse IL-33 Recombinant Rabbit mAb [PSH16-61] - BSA and Azide free (Capture)	HA723863	Mouse	ELISA(Cap)
	Mouse IL-33 Recombinant Rabbit mAb [PSH16-62] - BSA and Azide free (Detector)	HA723864	Mouse	ELISA(Det)
	Biotin Conjugated Mouse L-Selectin / CD62L Recombinant Rabbit mAb [PSH16-73] - Detector	HA723875B	Mouse	ELISA(Det),ELISA
	Mouse L-Selectin / CD62L Recombinant Rabbit mAb [PSH16-73] - BSA and Azide free (Detector)	HA723875	Mouse	ELISA(Det)
	Mouse L-Selectin / CD62L Recombinant Rabbit mAb [PSH16-74] - BSA and Azide free (Capture)	HA723876	Mouse	ELISA(Cap)
	Mouse L-Selectin / CD62L Recombinant Rabbit mAb [PSH16-72] - BSA and Azide free (Capture)	HA723874	Mouse	ELISA(Cap)
	Mouse Osteopontin Recombinant Rabbit mAb [PSH17-49] - BSA and Azide free (Capture)	HA723924	Mouse	ELISA(Cap)
	Mouse Osteopontin Recombinant Rabbit mAb [PSH17-50] - BSA and Azide free (Detector)	HA723925	Mouse	ELISA(Det)
	Biotin Conjugated Mouse Osteopontin Recombinant Rabbit mAb [PSH17-50] - Detector	HA723925B	Mouse	ELISA(Det),ELISA
	Mouse PCSK9, Tag Free Protein	HA211202	Mouse	
	HRP Conjugated Mouse/Rat Complement C3 Recombinant Rabbit mAb [PSH16-82] - Detector	HA725145H	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-24] - BSA and Azide free (Capture)	HA723911	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-25] - Detector	HA723912B	Mouse, Rat	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-25] - BSA and Azide free (Detector)	HA723912	Mouse, Rat	ELISA(Det)
	Mouse/Rat SCF Recombinant Rabbit mAb [PSH16-93] - BSA and Azide free (Capture)	HA723889	Mouse, Rat	ELISA(Cap)
	Mouse/Rat SCF Recombinant Rabbit mAb [PSH16-91] - BSA and Azide free (Capture)	HA723887	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat SCF Recombinant Rabbit mAb [PSH16-92] - Detector	HA723888B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat SCF Recombinant Rabbit mAb [PSH16-92] - BSA and Azide free (Detector)	HA723888	Mouse, Rat	ELISA(Det)
	Mouse ST2/IL-1 RL1, C-His Tag Protein	HA211193	Mouse	
	HRP Conjugated Mouse VCAM1 Recombinant Rabbit mAb [PSH16-32] - Detector	HA723830H	Mouse	ELISA(Det),ELISA
	mScarlet Nanoselector Agarose	028-101-002	Species independent	IP,ChIP,MS
	mScarlet Nanoselector Magnetic beads	028-101-003	Species independent	IP,ChIP,MS
	MUL1 Recombinant Rabbit mAb [PSH15-17] - BSA and Azide free	HA751546	Human, Mouse, Rat	WB,IHC-P
	mWasabi Nanoselector Agarose	049-101-002	Species independent	IP,ChIP,MS,Purification
	Myc tag Recombinant Rabbit mAb [PSH17-81]	HA723950	Species independent	WB,IF-Cell,IHC-P,FC,IP
	Myc tag Recombinant Mouse mAb [PSH14-94] - BSA and Azide free	HA610343	Species independent	WB,IF-Cell,IP
	Myc tag Nanoselector Agarose	002-101-002	Species independent	IP,ChIP,MS,Purification
	Myc tag Nanoselector Magnetic beads	002-101-003	Species independent	IP,CHIP,MS,Purification
	Neuropeptide Y Recombinant Rabbit mAb [PSH16-88] - BSA and Azide free (Capture)	HA723884	Human, Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Neuropeptide Y Recombinant Rabbit mAb [PSH16-89] - Detector	HA723885B	Human, Mouse, Rat	ELISA(Det),ELISA
	Biotin Conjugated Neuropeptide Y Recombinant Rabbit mAb [PSH16-90] - Detector	HA723886B	Human, Mouse, Rat	ELISA(Det),ELISA
	Neuropeptide Y Recombinant Rabbit mAb [PSH16-89] - BSA and Azide free (Detector)	HA723885	Human, Mouse, Rat	ELISA(Det)
Cell treatment	Neuropeptide Y Recombinant Rabbit mAb [PSH16-52] - BSA and Azide free	HA751601	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	Neuropeptide Y Recombinant Rabbit mAb [PSH16-90] - BSA and Azide free (Detector)	HA723886	Human, Mouse, Rat	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	NF-M Recombinant Antibody [JM11-20] - Chicken IgY (Chimeric) - BSA and Azide free	HA610337	Mouse	IHC-Fr
	non-muscle Myosin IIB / MYH10 Recombinant Rabbit mAb [PSH15-72] - BSA and Azide free	HA751579	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
Relative expression	NP-I / NPTX1 Recombinant Rabbit mAb [PSH17-07] - BSA and Azide free	HA751613	Human, Mouse, Rat	WB,IHC-P,FC
Relative expression	NP-I / NPTX1 Recombinant Rabbit mAb [PSH17-07]	HA723895	Human, Mouse, Rat	WB,IHC-P,FC
	Orexin Recombinant Rabbit mAb [PSH16-77] - BSA and Azide free	HA751606	Mouse, Rat	IHC-Fr,IHC-P
Relative expression	OSMR Recombinant Rabbit mAb [PSH17-60]	HA723935	Human	WB
	OTULIN Recombinant Rabbit mAb [PSH16-95] - BSA and Azide free	HA751609	Human	WB,FC
Relative expression	p16INK4a Recombinant Rabbit mAb [PSH15-40] - BSA and Azide free	HA751555	Mouse	WB,IHC-P,IP
	Pan-Cadherin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS287RT	Human	mIHC
	Pan-Cadherin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS288RT	Mouse	mIHC
	Pan-Cadherin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS286MS	Mouse	mIHC
	Pan-Cadherin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS285MS	Human	mIHC
	PDE1B Recombinant Rabbit mAb [PSH17-59]	HA723934	Mouse, Rat	WB,IHC-P,IHC-Fr
	PDGFR beta Recombinant Antibody - Rat IgG1 (Chimeric)	IRS283RT	Mouse	mIHC
Relative expression	Perforin-1 Recombinant Rabbit mAb [PSH15-92] - BSA and Azide free	HA751583	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	PGP9.5 Recombinant Antibody [JM10-59] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610344	Human, Mouse, Rat	IHC-Fr,IF-Tissue,IHC-P,WB,IF-Cell
Relative expression	PGP9.5 Recombinant Antibody [JM10-59] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610345	Human, Mouse, Rat	IHC-Fr,IF-Tissue,IHC-P,WB,IF-Cell
Cell treatment	Phospho-Beclin 1 (S30) Recombinant Rabbit mAb [PSH16-40] - BSA and Azide free	HA751598	Human, Mouse, Rat	WB
Cell treatment	Phospho-Beclin 1 (S30) Recombinant Rabbit mAb [PSH16-40]	HA723840	Human, Mouse, Rat	WB
Cell treatment	Phospho-MNK1 (T250 + T255) Recombinant Rabbit mAb [PSH17-95]	HA723960	Human, Mouse	WB,IHC-P,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH16-27] - BSA and Azide free	HA751595	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH15-91] - BSA and Azide free	HA751582	Human, Mouse, Rat	WB,IF-Cell,Dot Blot
	Phospho-Tau (T217) Recombinant Rabbit mAb [PSH17-41] - BSA and Azide free (Capture)	HA723922	Human, Mouse	ELISA(Cap)
Relative expression	PIGR Rabbit pAb	HA500593	Human	WB,IHC-P
Relative expression	PIGR Recombinant Rabbit mAb [PSH17-84]	HA723953	Human, Mouse, Rat	WB,IP
Relative expression, Cell treatment	Prosurfactant Protein C Recombinant Rabbit mAb [PSH15-28] - BSA and Azide free	HA751550	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue
Cell treatment	PSMB5 / MB1 Recombinant Rabbit mAb [PSH15-20] - BSA and Azide free	HA751549	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
Relative expression	QKI-5 Recombinant Antibody [PSH10-22] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610331	Mouse, Rat	IHC-Fr,WB
Relative expression	QKI-6 Recombinant Antibody [PSH10-23] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610332	Mouse, Rat	IHC-Fr,WB
	Rabbit IgG Nanoselector Agarose	025-101-002	Species independent	IP,ChIP,MS,Purification
	Rabbit IgG Nanoselector Magnetic beads	025-101-003	Species independent	IP,ChIP,Purification
	Biotin Conjugated Rat ICAM1 Recombinant Rabbit mAb [PSH17-52] - Detector	HA723927B	Rat	ELISA(Det),ELISA
	Rat ICAM1 Recombinant Rabbit mAb [PSH17-52] - BSA and Azide free (Detector)	HA723927	Rat	ELISA(Det)
	Rat ICAM1 Recombinant Rabbit mAb [PSH17-51] - BSA and Azide free (Capture)	HA723926	Rat	ELISA(Cap)
	Rat ICAM1 Recombinant Rabbit mAb [PSH17-53] - BSA and Azide free (Capture)	HA723928	Rat	ELISA(Cap)
	Rat IL-18, Tag Free Protein	HA211199	Rat	
	RBP4 Recombinant Rabbit mAb [PSH16-97] - BSA and Azide free	HA751611	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	RFP Nanoselector Agarose	020-101-002	Species independent	IP,ChIP,MS,Purification
	RFP Nanoselector Magnetic beads	020-101-003	Species independent	IP,ChIP,MS,Purification
Relative expression	Ryanodine receptor 2 / RYR-2 Recombinant Rabbit mAb [PSH15-54] - BSA and Azide free	HA751567	Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	S100 beta Recombinant Antibody [SC57-02] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610366	Mouse, Rat	IHC-Fr,IHC-P
	S100 beta Recombinant Antibody [SC57-02] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610367	Mouse, Rat	IHC-Fr,IHC-P
	S100 beta Recombinant Antibody [SC57-02] - Rat IgG1 (Chimeric)	HA601552	Mouse, Rat	IHC-Fr,IHC-P
	S100 beta Recombinant Antibody [SC57-02] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610365	Mouse	IHC-Fr
Relative expression	S100 beta Recombinant Rabbit mAb [PSH17-13]	HA723900	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Tissue
	S100 beta Recombinant Antibody [SC57-02] - Guinea pig IgG2 (Chimeric)	HA601550	Mouse	IHC-Fr
	S100 beta Recombinant Antibody [SC57-02] - Mouse IgG1 (Chimeric)	HA601551	Mouse, Rat	IHC-Fr,IHC-P
	S100 beta Recombinant Rabbit mAb [PSH17-13] - BSA and Azide free	HA751617	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Tissue
	SATB2 Recombinant Antibody [JM52-44] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610338	Human, Mouse, Rat	IHC-Fr,WB
	SCD1 Recombinant Rabbit mAb [PSH15-47] - BSA and Azide free	HA751562	Mouse, Rat	WB,IHC-Fr,IHC-P,FC
	SCNN1G Recombinant Rabbit mAb [PSH16-29] - BSA and Azide free	HA751597	Human, Mouse, Rat, Monkey	WB,IF-Cell
	SIRT2 Recombinant Antibody [SN70-04] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610335	Mouse, Rat	IHC-Fr,IHC-P,WB
	SIRT2 Recombinant Antibody [SN70-04] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610336	Mouse, Rat	IHC-Fr,IHC-P,WB
	SIRT2 Recombinant Rabbit mAb [PSH17-48]	HA723923	Human, Mouse, Rat, Cynomolgus monkey	IHC-Fr,IHC-P,WB
	SIRT4 Recombinant Rabbit mAb [PSH15-46] - BSA and Azide free	HA751561	Human, Mouse, Rat, Monkey	WB,IHC-P,IP
	SIRT5 Recombinant Rabbit mAb [PSH15-71] - BSA and Azide free	HA751578	Human, Mouse, Rat	WB,IHC-P
	SMC5 Recombinant Rabbit mAb [PSH15-29] - BSA and Azide free	HA751551	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P
	SNAP25 Recombinant Rabbit mAb [PSH17-66] - BSA and Azide free (Detector)	HA725228	Human, Mouse	ELISA(Det)
	Biotin Conjugated SNAP25 Recombinant Rabbit mAb [PSH17-66] - Detector	HA725228B	Human, Mouse	ELISA(Det),ELISA
	SNAP25 Recombinant Rabbit mAb [PSH17-65] - BSA and Azide free (Capture)	HA725227	Human, Mouse	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	SNAP tag Nanoselector Agarose	011-101-002	Species independent	IP,ChIP,MS,Purification
	SOX2 Recombinant Antibody [PO00-28] - Chicken IgY (Chimeric) - BSA and Azide free	HA610329	Human, Mouse	IHC-Fr,WB
	SOX2 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS273RT	Human	mIHC
	SOX2 Recombinant Antibody [PO00-28] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610330	Human, Mouse	IHC-Fr,WB
	SOX2 Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS281MS	Human	mIHC
	SOX2 Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS282MS	Mouse	mIHC
	SOX2 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS274RT	Mouse	mIHC
	SOX9 Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS279MS	Human	mIHC
	SOX9 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS280RT	Human	mIHC
	SpyCatcher Magnetic beads	086-101-003	Species independent	IP,ChIP,MS,Purification
	SpyCatcher Nanoselector Agarose	086-101-002	Species independent	IP,ChIP,MS,Purification
	StayGold/mBaojin Nanoselector Agarose	082-101-002	Species independent	IP,ChIP,MS,Purification
	StayGold/mBaojin Nanoselector Magnetic beads	082-101-003	Species independent	IP,ChIP,MS,Purification
	Streptavidin Magnetic beads	067-101-003	Species independent	IP,ChIP,MS,Purification
	TagFP Nanoselector Agarose	017-101-002	Species independent	IP,ChIP,MS,Purification
Cell treatment	TAP2 Recombinant Rabbit mAb [PSH16-55] - BSA and Azide free	HA751604	Human, Mouse, Rat	WB,IHC-P,IP
	TAZ / WWTR1 Recombinant Rabbit mAb [PSH15-58] - BSA and Azide free	HA751569	Human, Mouse, Rat	WB,IHC-P,IP
	TAZ / WWTR1 Recombinant Rabbit mAb [PSH15-57] - BSA and Azide free	HA751568	Human, Mouse, Rat, Monkey	WB,IP
Relative expression	TCF-4 Recombinant Rabbit mAb [PSH15-19] - BSA and Azide free	HA751548	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	TCR beta chain Recombinant Mouse mAb [PSH16-30] - BSA and Azide free	HA610355	Mouse	FC
KO	TDG Recombinant Rabbit mAb [PSH10-11] - BSA and Azide free	HA751531	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	tdTomato Nanoselector Agarose	026-101-002	Species independent	IP,ChIP,MS
	tdTomato Nanoselector Magnetic beads	026-101-003	Species independent	IP,ChIP,MS
Relative expression	TFPI2 Recombinant Rabbit mAb [PSH15-45] - BSA and Azide free	HA751560	Human	WB,IHC-P,IP
Relative expression	TIM 3 Recombinant Rabbit mAb [PSH16-50] - BSA and Azide free	HA751599	Human, Mouse, Rat	WB,IHC-P,IP
	TMEM184A Rabbit pAb	HA500595	Human, Mouse	WB
Relative expression	TMS1 / ASC Recombinant Rabbit mAb [PSH17-35]	HA723918	Mouse, Rat	WB
Relative expression	TMS1 / ASC Recombinant Rabbit mAb [PSH17-34]	HA723917	Human, Mouse	WB,IHC-Fr,IHC-P
Relative expression	TMS1 / ASC Recombinant Rabbit mAb [PSH17-34] - BSA and Azide free	HA751622	Human, Mouse	WB,IHC-Fr,IHC-P
Relative expression	TMS1 / ASC Recombinant Rabbit mAb [PSH17-35] - BSA and Azide free	HA751623	Mouse, Rat	WB
Relative expression	TPH2 Recombinant Antibody [PSH04-89] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610322	Mouse, Rat	IHC-Fr,WB
Relative expression	TPH2 Recombinant Antibody [PSH04-89] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610324	Mouse, Rat	WB,IHC-P
Relative expression	TPH2 Recombinant Antibody [PSH04-89] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610323	Mouse, Rat	WB,IHC-P
	TRBC1 Recombinant Mouse mAb [PSH15-51] - BSA and Azide free	HA610351	Human	FC
Relative expression	TRPV4 Recombinant Rabbit mAb [PSH16-51] - BSA and Azide free	HA751600	Human, Mouse	WB,IHC-Fr,IHC-P,IP
	TurboGFP Nanoselector Agarose	014-101-002	Species independent	IP,ChIP,MS,Purification
Relative expression	UCP1 Recombinant Rabbit mAb [PSH15-73] - BSA and Azide free	HA751580	Mouse, Rat	WB,IHC-Fr,IHC-P,IP
	USP15 Recombinant Rabbit mAb [PSH15-94] - BSA and Azide free	HA751585	Human, Mouse, Rat, Monkey	WB,IHC-P,IP
	V5 tag Nanoselector Agarose	064-101-002	Species independent	IP,ChIP,MS,Purification
	V5 tag Nanoselector Magnetic beads	064-101-003	Species independent	IP,ChIP,MS,Purification
	VDAC1 / 2 / 3 Recombinant Rabbit mAb [PSH14-90] - BSA and Azide free	HA751534	Human, Mouse, Rat	WB
	Vimentin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS263RT	Mouse	mIHC
KD	VPS35 Recombinant Rabbit mAb [PSH16-54] - BSA and Azide free	HA751603	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	WVOX Rabbit pAb	HA500592	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
	WVOX Recombinant Rabbit mAb [PSH15-43] - BSA and Azide free	HA751558	Human, Mouse, Rat, Monkey	WB,IF-Cell,IP
	WVOX Recombinant Rabbit mAb [PSH15-44] - BSA and Azide free	HA751559	Human, Mouse, Rat, Monkey	WB,IP
Relative expression	YKL-40 / CHI3L1 Recombinant Rabbit mAb [PSH16-26] - BSA and Azide free	HA751594	Human	WB,IF-Cell
	荧光淬灭试剂盒 IRISKit®HyperView Quench Buffer	900102		mIHC
	虹和CmTSA多重免疫荧光染色试剂盒(20T-4C , 1Rabbit+1Mouse+1Rat)IRISKitCmTSA Kit	900801		mIHC
	虹和CmTSA多重免疫荧光染色试剂盒(套装 20T-4C , 1Mouse+2Rat)IRISKitCmTSA Kit	900807		mIHC
	虹和CmTSA多重免疫荧光染色试剂盒(套装 20T-4C , 1Rabbit+2Mouse)IRISKitCmTSA Kit	900803		mIHC
	虹和CmTSA多重免疫荧光染色试剂盒(套装 20T-4C , 1 Rabbit+2Rat)IRISKitCmTSA Kit	900805		mIHC
	虹和CmTSA多重免疫荧光染色试剂盒(套装 20T-4C , 2Mouse+1Rat)IRISKitCmTSA Kit	900806		mIHC
	虹和CmTSA多重免疫荧光染色试剂盒(套装 20T-4C , 2Rabbit+1Mouse)IRISKitCmTSA Kit	900802		mIHC
	虹和CmTSA多重免疫荧光染色试剂盒(套装 20T-4C , 2Rabbit+1Rat)IRISKitCmTSA Kit	900804		mIHC
	虹和CmTSA多重免疫荧光染色试剂盒(套装 20T-4C , 3Mouse)IRISKitCmTSA Kit	900809		mIHC
	虹和CmTSA多重免疫荧光染色试剂盒(套装 20T-4C , 3Rabbit)IRISKitCmTSA Kit	900808		mIHC
	虹和CmTSA多重免疫荧光染色试剂盒(套装 20T-4C , 3Rat)IRISKitCmTSA Kit	900810		mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层