

细胞外基质（ECM）研究-抗体&ELISA kit产品清单

高级验证	产品名	货号	反应种属	验证应用
	4.1R Rabbit pAb	ER62531	Human, Mouse	WB
	67kDa Laminin Receptor Recombinant Rabbit mAb [JF0955]	ET1702-05	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	ACTG Rabbit pAb	ER61700	Human, Rat, Mouse	WB
	Actin Recombinant Rabbit mAb [JJ09-29]	ET1701-80	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P
	Actin Recombinant Antibody [JM36-62] - Human IgG1 (Chimeric)	HA721878	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	iFluor™ 488 Conjugated Actin Recombinant Rabbit mAb [JJ09-29]	HA720160F	Human, Mouse, Rat	IF-Cell,IF-Tissue,FC,WB
	Actin Recombinant Rabbit mAb [JF47-01]	ET1702-52	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IF-Cell,FC
	Actin Recombinant Antibody [JM36-62] - Human IgG1 (Chimeric) - BSA and Azide free	HA750844	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	HRP Conjugated Actin Recombinant Rabbit mAb [JJ09-29]	HA721186	Human, Mouse, Rat	WB
	Actin Recombinant Rabbit mAb [JJ09-29] - BSA and Azide free	HA750323	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P
	Actin Rabbit pAb	ER62585	Human, Mouse, Rat	WB,IHC-P
	Actin alpha 3 Rabbit pAb	ER62586	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Actinin alpha 2/3 Rabbit pAb	ER62588	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Actinin alpha 3 Rabbit pAb	ER62589	Human, Mouse, Rat	WB,IF-Cell
	Activin A Receptor Type IB Recombinant Rabbit mAb [JE38-48]	HA721550	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Activin Receptor Type IA Rabbit pAb	HA500117	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Activin Receptor Type IIB / ACVR2B Recombinant Rabbit mAb [JE33-52]	HA722725	Human, Mouse, Rat	WB
	ACTL6A Recombinant Rabbit mAb [JU40-32]	ET7106-60	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP
	ACTR1B Recombinant Rabbit mAb [JE65-44]	HA721050	Human	WB,IHC-P,FC
	ACTR-IB Rabbit pAb	ER62592	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	ACVR2A Mouse mAb [A7C8]	HA601091	Human	WB
	ACVR2A Rabbit pAb	HA500084	Human, Mouse, Rat	WB

	ACVR2B Mouse mAb [A5H7]	HA600067	Human, Mouse, Rat, Zebrafish	WB
	ACVR2B Mouse mAb [A5H6]	HA600066	Human, Mouse, Rat, Zebrafish	WB
KD	ADAM10 Recombinant Rabbit mAb [JM32-11]	ET1703-60	Human, Mouse, Rat	WB,IHC-P,IP
	ADAM10 Rabbit pAb	ER1706-97	Human	WB,IHC-P
	ADAM11 Rabbit pAb	ER1903-24	Human, Mouse	WB
	ADAM12 Rabbit pAb	ER1903-25	Human, Mouse	WB,FC
	ADAM15 Recombinant Rabbit mAb [JE39-01]	HA722402	Human, Mouse, Rat	WB,IHC-P
	ADAM15 Rabbit pAb	ER1903-26	Human	WB
	ADAM17 Recombinant Rabbit mAb [JM10-35]	ET1703-06	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,FC,IP,IHC-Fr,IHC-P
	ADAM17 Recombinant Rabbit mAb [JM10-35] - BSA and Azide free	HA750370	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,FC,IP,IHC-Fr,IHC-P
	ADAM19 Rabbit pAb	ER1903-28	Human	WB
	ADAM2 Rabbit pAb	ER1903-29	Human	WB
	ADAM22 Rabbit pAb	ER2001-24	Human, Mouse	WB
	ADAM23 Rabbit pAb	ER1903-30	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	ADAM9 Rabbit pAb	ER60448	Human	WB
Cell treatment	ADAMTS1 Recombinant Rabbit mAb [PSH18-51]	HA724004	Human, Mouse, Rat	WB
	ADAMTS-12 Rabbit pAb	ER62597	Human	IHC-P
	ADAMTS13 Recombinant Rabbit mAb [JG39-21]	ET7108-86	Human, Rat	WB,IHC-P
	ADAMTS-17 Rabbit pAb	ER62599	Human, Mouse	IHC-P
	ADAMTS-19 Rabbit pAb	ER62600	Human, Mouse	WB
	ADAMTS4 Rabbit pAb	ER1903-31	Human, Mouse, Rat	WB
	ADAMTS5 Recombinant Rabbit mAb [PSH03-52]	HA722011	Human, Mouse, Rat	WB,IF-Cell,FC
	ADAMTS5 Rabbit pAb	ER1903-32	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	ADAMTS5 Recombinant Rabbit mAb [PSH03-52] - BSA and Azide free	HA750886	Human, Mouse, Rat	WB,IF-Cell,FC
	ADAMTS-7 Rabbit pAb	ER62601	Human	IHC-P
	ADAMTS7(NT) Rabbit pAb	ER1903-33	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	ADAMTSL4 Rabbit pAb	ER1903-35	Mouse	WB
	Adducin alpha/beta Rabbit pAb	ER62606	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Adducin beta Rabbit pAb	ER62607	Human, Mouse, Rat	WB,IF-Cell
	Adducin gamma Rabbit pAb	ER62608	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	ADRM1 Recombinant Rabbit mAb [PSH05-41]	HA722270	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	ADRM1 Recombinant Rabbit mAb [PSH05-41] - BSA and Azide free	HA751000	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	Afadin/I Afadin Rabbit pAb	ER1903-50	Mouse, Rat	IHC-P,IF-Tissue
	Aggrecan Recombinant Rabbit mAb [JA40-23]	ET1704-57	Human, Rat	IHC-P
	AGRIN Rabbit pAb	ER60457	Human, Mouse, Rat	WB
Relative expression	ALP Recombinant Rabbit mAb [PSH14-35] - BSA and Azide free	HA751517	Human, Mouse, Rat	WB,IP
Relative expression	ALP Recombinant Rabbit mAb [PSH14-35]	HA723629	Human, Mouse, Rat	WB,IP
	alpha 1 Catenin Recombinant Rabbit mAb [SC0619]	ET1610-46	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	alpha 1 Spectrin Rabbit pAb	HA500357	Human, Mouse	WB
	alpha Actinin Mouse mAb [A1G2]	EM1901-52	Human, Rat	WB,IF-Cell,IHC-P,FC
	alpha Actinin Recombinant Mouse mAb [A1G2-R]	HA601194	Human, Mouse, Rat	WB,IHC-P
	alpha Actinin Recombinant Mouse mAb [A1G2-R] - BSA and Azide free	HA610074	Human, Mouse, Rat	WB,IHC-P
	alpha Actinin Rabbit pAb	ER1803-60	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	Alpha Actinin 2 Recombinant Rabbit mAb [JF0980]	ET1702-30	Human, Mouse, Rat	WB,IHC-P,IP
KD	alpha Actinin 4 Recombinant Rabbit mAb [JU20-23]	ET1706-05	Human, Mouse, Rat	WB,IF-Cell,IP,IHC-P,FC
	Alpha-cardiac actin Mouse mAb [33-32]	M1206-1	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,IF-Tissue,IHC-Fr
	Alpha-cardiac actin Recombinant Mouse mAb [33-32-R]	HA601208	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Alpha-cardiac actin Recombinant Mouse mAb [33-32-R] - BSA and Azide free	HA610086	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Alpha-cardiac actin Rabbit pAb	0407-3	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P
	Alpha-dystroglycan Rabbit pAb	ER1901-72	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Alpha Fodrin Recombinant Rabbit mAb [JU32-09]	ET1706-44	Human, Mouse, Rat	WB,IHC-P
	alpha sarcoglycan Recombinant Rabbit mAb [JA51-81]	ET1704-25	Human, Mouse, Rat	WB,IP,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD, Relative expression	alpha smooth muscle Actin Recombinant Rabbit mAb [SY02-64] - BSA and Azide free	HA750126	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P,FC
KD, Relative expression	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Rat IgG1 (Chimeric)	HA601427	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD, Relative expression	alpha smooth muscle Actin Recombinant Rabbit mAb [SY02-64]	ET1607-53	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P,FC,mIHC
	alpha smooth muscle Actin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS094RT	Human	mIHC
	alpha smooth muscle Actin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS220RT	Mouse	mIHC
	alpha smooth muscle Actin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS093MS	Human	mIHC
	alpha smooth muscle Actin Recombinant Rabbit mAb	IRS048	Mouse	mIHC
	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610362	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	alpha smooth muscle Actin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS219MS	Mouse	mIHC
KD, Relative expression	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610267	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	alpha smooth muscle Actin Recombinant Rabbit mAb [SY25-03]	ET1607-43	Human, Mouse, Rat	WB,IHC-P,FC,IP
	alpha smooth muscle Actin Recombinant Rabbit mAb	IRS012	Human	mIHC
	alpha smooth muscle Actin Recombinant Antibody [SY02-64] - Mouse IgG1 (Chimeric)	HA601546	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	alpha Tubulin Rabbit pAb	ER130905	Human, Mouse, Rat, Zebrafish	WB,IHC-P,FC,IF- Cell,IF-Tissue
	alpha Tubulin Recombinant Rabbit mAb [PSH02-95]	HA721914	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	alpha Tubulin Rabbit pAb	HA500366	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	alpha Tubulin Recombinant Rabbit mAb [PSH0-14] - BSA and Azide free	HA750584	Human, Mouse, Rat	WB,IHC-P
	alpha Tubulin Recombinant Rabbit mAb [PSH02-94]	HA721913	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	alpha Tubulin Recombinant Rabbit mAb [PSH02-94] - BSA and Azide free	HA750849	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	alpha Tubulin Mouse mAb [A8-6]	M1501-1	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC
	alpha Tubulin Recombinant Rabbit mAb [PSH0-14]	HA721272	Human, Mouse, Rat	WB,IHC-P
	alpha Tubulin Rabbit pAb	R1308-14	Human, Rat, Mouse	WB,IF-Cell
	alpha Tubulin Recombinant Rabbit mAb [PSH02-95] - BSA and Azide free	HA750850	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
KD	alpha Tubulin 4A Recombinant Rabbit mAb [JM73-24]	ET1705-31	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	alpha Tubulin (Acetyl K112) Rabbit pAb	ER60054	Human, Mouse, Rat	WB
Cell treatment	alpha Tubulin (acetyl K40) Recombinant Rabbit mAb [PSH0-13] - BSA and Azide free	HA750583	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	alpha Tubulin (acetyl K40) Recombinant Rabbit mAb [PSH0-13]	HA721271	Human, Mouse, Rat	WB,IF-Cell
	Ameloblastin Rabbit pAb	ER62654	Human	WB
	ANG I Rabbit pAb	ER65619	Human	WB,IHC-P
Relative expression	Angiomotin Recombinant Rabbit mAb [PSH07-46]	HA722834	Human	WB,IHC-P,IF-Cell,IP
Relative expression	Angiomotin Recombinant Rabbit mAb [PSH07-46] - BSA and Azide free	HA751167	Human	WB,IHC-P,IF-Cell,IP
	Angiomotin-L1 Rabbit pAb	ER62666	Human, Mouse	IHC-P
	Annexin A1 Rabbit pAb	ER1803-70	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Annexin A1 Recombinant Rabbit mAb [PSH04-13]	HA722094	Human	IHC-P,IF-Tissue
	Annexin A10 Recombinant Rabbit mAb [SC06-02]	ET1610-57	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Annexin A13 Rabbit pAb	ER1903-72	Human	WB,IHC-P,IF-Tissue
KD	Annexin A2 Recombinant Rabbit mAb [JA42-30]	ET1704-49	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
	Annexin IV Recombinant Rabbit mAb [JG37-83]	ET7108-56	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	ANXA6 Recombinant Rabbit mAb [PSH05-82]	HA722481	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP
	ANXA6 Recombinant Rabbit mAb [PSH05-82] - BSA and Azide free	HA751030	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP
KD	ANXA7 Recombinant Rabbit mAb [JE65-60]	HA721068	Human	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	ANXA7 Recombinant Rabbit mAb [JE65-36]	HA721078	Human, Rat	WB
	AQP0 Rabbit pAb	ER62698	Human, Mouse, Rat	WB,IHC-P
	ARHGAP17 Rabbit pAb	ER62715	Human, Mouse, Rat	WB
	ARHGDI A Mouse mAb [G8-C3]	EM1710-11	Human, Mouse, Monkey	WB,IF-Cell,FC
	Arp2 Recombinant Rabbit mAb [JE40-62]	HA721348	Human, Rat, Mouse	WB,IHC-P,IF-Cell
	ARPC1A Rabbit pAb	ER1903-94	Human, Rat	IHC-P,IF-Cell,IF-Tissue
	ARPC2 Recombinant Rabbit mAb [JE58-79]	HA720030	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	ARPC5 / p16 ARC Recombinant Rabbit mAb [SR34-02] - BSA and Azide free	HA750055	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	ARPC5 / p16 ARC Recombinant Rabbit mAb [SR34-02]	ET1602-9	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	ASAP1 Rabbit pAb	ER1903-96	Human, Mouse, Rat	WB
	ASAP2 Rabbit pAb	ER60163	Human, Mouse	WB
	ATG4A Rabbit pAb	ER1901-38	Human, Mouse, Rat	WB,IHC-P
	ATG4A Recombinant Rabbit mAb [JE32-31]	HA721010	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	ATG4B Recombinant Rabbit mAb [JE59-97]	HA721334	Human, Mouse, Rat, Monkey	WB
KD	ATG4B Recombinant Rabbit mAb [PSH02-53]	HA721829	Human, Mouse, Rat, Monkey	WB
KD	ATG4B Recombinant Rabbit mAb [PSH02-53] - BSA and Azide free	HA750816	Human, Mouse, Rat, Monkey	WB
	ATG4B Rabbit pAb	ER62759	Human, Mouse, Rat, Monkey	WB,IHC-P
	Bestrophin-1 Rabbit pAb	ER60094	Human, Rat, Mouse	WB,IHC-P
	beta-1,4-GalNAc-T Rabbit pAb	ER65392	Human, Mouse, Rat	WB,IHC-P
	beta Actin Recombinant Rabbit mAb [PSH03-63] - BSA and Azide free	HA750898	Human, Mouse, Rat, Rabbit, Goat, Zebrafish, Dog, Monkey, Pig, Chicken, Cow, Hamster, Green monkey	WB,IF-Cell,FC,IP
	beta Actin Rabbit pAb	R1207-1	Human, Mouse, Rat, Zebrafish, Bamboo	WB,IF-Cell,IHC-P,FC
	beta Actin Recombinant Rabbit mAb [PSH03-63]	HA722023	Human, Mouse, Rat, Rabbit, Goat, Zebrafish, Dog, Monkey, Pig, Chicken, Cow, Hamster, Green monkey	WB,IF-Cell,FC,IP
	beta Actin Mouse mAb [A2-F6]	EM21002	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	beta Actin Recombinant Mouse mAb [A2-F6-R]	HA601082	Human, Mouse, Rat	WB
	HRP Conjugated beta Actin Recombinant Mouse mAb [A2-F6-R]	HA601037	Human, Mouse, Rat	WB
	beta Actin Mouse mAb [B4-B2]	M1210-2	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	HRP Conjugated beta Actin Mouse mAb [A2-F6]	M1210-5	Human, Mouse, Rat	WB
	beta Actin Recombinant Mouse mAb [A2-F6-R] - BSA and Azide free	HA610019	Human, Mouse, Rat	WB
	HRP Conjugated beta Actin Recombinant Rabbit mAb [JF53-10]	ET1702-67	Human, Mouse, Rat, Zebrafish, Monkey, Hamster, Arabidopsis thaliana, Rice	WB
	beta Actin Rabbit pAb	R1102-1	Human, Mouse, Rat	WB,IHC-P
	Beta Catenin Recombinant Rabbit mAb	IRS043	Human	mIHC
	Beta Catenin Recombinant Mouse mAb [A6-F8-R] - BSA and Azide free	HA610062	Human, Mouse, Rat	WB,IF-Cell,FC
	Beta Catenin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS277MS	Human	mIHC
KD	Beta Catenin Recombinant Rabbit mAb [SA30-04] - BSA and Azide free	HA750023	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP,IF-Cell,IHC-Fr,FC
	Beta Catenin Rabbit pAb	ER0805	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell,IF-Tissue
	Beta Catenin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS278RT	Human	mIHC
KD	Beta Catenin Recombinant Rabbit mAb [SA30-04]	ET1601-5	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IP,mIHC,IF-Cell,IHC-Fr,FC
	Beta Catenin Recombinant Rabbit mAb	IRS061	Mouse	mIHC
	Beta Catenin Mouse mAb [10-C0-B7]	M1405-6	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	iFluor™ 488 Conjugated Beta Catenin Recombinant Rabbit mAb [SA30-04]	HA720182F	Human, Mouse	IF,ICC
	Beta Catenin Rabbit pAb	0407-16	Human, Mouse, Rat, Zebrafish	WB,IF-Cell
	Beta Catenin Mouse mAb [A6-F8]	EM0306	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,ChIP
	Beta Catenin Recombinant Mouse mAb [A6-F8-R]	HA601179	Human, Mouse, Rat	WB,IF-Cell,FC
	Beta III Tubulin Rabbit pAb	0805-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	iFluor™ 594 Conjugated Beta III Tubulin Mouse mAb [A8-D10]	HA600133F	Human, Rat	IF,ICC,FC
Relative expression	Beta III Tubulin Recombinant Rabbit mAb [PSH15-32]	HA723738	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP
Relative expression	Beta III Tubulin Recombinant Rabbit mAb [PSH15-32] - BSA and Azide free	HA751554	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP
	iFluor™ 488 Conjugated Beta III Tubulin Mouse mAb [A8-D10]	HA600106F	Human, Mouse, Rat	IF-Cell,FC,IF-Tissue
	Beta III Tubulin Recombinant Rabbit mAb [SP06-00]	ET1604-17	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
	Beta III Tubulin Recombinant Rabbit mAb [SP06-00] - BSA and Azide free	HA750071	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
	Beta III Tubulin Rabbit pAb	R1401-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Beta III Tubulin Mouse mAb [A8-D10]	M0805-8	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
	Beta III Tubulin Rabbit pAb	0807-2	Human, Mouse, Rat	WB,IF-Cell,FC
	beta II Tubulin Recombinant Rabbit mAb [ST52-04]	ET1609-48	Mouse, Rat, Human	WB,IF-Tissue,IHC-P
	beta I Tubulin Recombinant Rabbit mAb [JE01-23]	HA722290	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IF-Tissue
Relative expression	beta IV Tubulin Recombinant Rabbit mAb [PSH17-78]	HA723947	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	beta Tubulin Recombinant Rabbit mAb [SR25-04]	ET1602-4	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	iFluor™ 594 Conjugated beta Tubulin Recombinant Mouse mAb [A1-A4-R]	HA600110F	Human, Mouse	IF-Cell,IF-Tissue,FC
	beta Tubulin Recombinant Rabbit mAb [SR25-04] - BSA and Azide free	HA750049	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	beta Tubulin Recombinant Mouse mAb [A1-A4-R] - BSA and Azide free	HA610067	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IHC-P
	HRP Conjugated beta Tubulin Recombinant Rabbit mAb [JF41-50]	ET1702-68	Human, Mouse, Rat	WB
	beta Tubulin Mouse mAb [1-B11]	EM0103	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC,IF-Tissue
	beta Tubulin Recombinant Mouse mAb [1-B11-R] - BSA and Azide free	HA610049	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	beta Tubulin Recombinant Mouse mAb [1-B11-R]	HA601147	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	beta Tubulin Mouse mAb [A1-A4]	M1305-2	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	HRP Conjugated beta Tubulin Mouse mAb [1-B11]	EM1701-59	Human, Mouse, Rat, Zebrafish	WB,IHC-P
	beta Tubulin Recombinant Mouse mAb [A1-A4-R]	HA601187	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biglycan Recombinant Rabbit mAb [JB71-31]	ET7107-48	Human, Mouse	WB,IHC-P,FC,IP
	Bmf Recombinant Rabbit mAb [JE33-57]	HA721487	Human, Mouse, Rat	WB
	BMP-1 Rabbit pAb	ER65804	Human, Mouse, Rat	IHC-P
	BMP2 Rabbit pAb	ER80602	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	BMP2 Recombinant Rabbit mAb [JE10-29]	HA722693	Human, Mouse, Rat	WB
	BMP-3 Rabbit pAb	ER62817	Human, Mouse, Rat	WB,IHC-P
	BMP-5 Rabbit pAb	ER62818	Human, Mouse	WB,IHC-P
	BMP7 Rabbit pAb	ER1904-26	Human, Mouse, Rat, Rabbit	WB,IHC-P,IF-Tissue
	BMP-8 Rabbit pAb	ER62819	Human	IHC-P
	BMP-8B Rabbit pAb	ER62820	Human	WB
	BMP-9 Rabbit pAb	ER65787	Human, Mouse	IHC-P
	BMPR1B Rabbit pAb	R1311-4	Human, Mouse	WB,IF-Cell,IHC-P
	BMPR2 Mouse mAb [F8-G6]	EM1709-92	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P
	Bone Sialoprotein Rabbit pAb	ER1904-28	Mouse	WB
	Bone sialoprotein 2 Recombinant Rabbit mAb [PSH09-64] - BSA and Azide free	HA751306	Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Tissue
	Bone sialoprotein 2 Recombinant Rabbit mAb [PSH09-64]	HA723129	Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Tissue
	Brain-derived tau/BD-tau Recombinant Rabbit mAb [PSH16-47] - BSA and Azide free (Detector)	HA723848	Human	ELISA(Det)
	Biotin Conjugated Brain-derived tau/BD-tau Recombinant Rabbit mAb [PSH16-47] - Detector	HA723849B	Human	ELISA(Det),ELISA
	Brain-derived tau/BD-tau Recombinant Rabbit mAb [PSH16-46] - BSA and Azide free (Capture)	HA723847	Human, Mouse, Rat	ELISA(Cap)
	Cadherin 10 Rabbit pAb	HA500059	Human, Mouse	WB
	Cadherin 16 Recombinant Rabbit mAb [JE60-37]	HA722408	Human, Mouse	WB,IHC-P,IF-Tissue
	Cadherin-19 Rabbit pAb	ER62867	Human	WB,IHC-P,IF-Cell
	Cadherin-23 Rabbit pAb	ER62870	Human, Mouse, Rat	IF-Cell
	Cadherin-7 Rabbit pAb	ER62873	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Cadherin-9 Rabbit pAb	ER62875	Human, Mouse	WB,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Cadherin EGF LAG seven-pass G-type receptor 2 Rabbit pAb	0905-6	Human, Mouse	WB,IF-Cell
	Caldesmon Recombinant Rabbit mAb [SU03-19]	ET1608-16	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	Calponin Recombinant Rabbit mAb [SI67-01] - BSA and Azide free	HA750094	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Calponin Recombinant Rabbit mAb [SI67-01]	ET1606-17	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Calponin-1 Rabbit pAb	R1510-7	Human, Mouse, Rat	IF-Cell,IHC-P,FC
	Calponin 2 Rabbit pAb	ER62884	Human, Mouse	WB,IF-Cell
	CAP1 Recombinant Rabbit mAb [JE61-08]	HA720069	Human, Rat	WB,IHC-P,IF-Cell
	CAPG Recombinant Rabbit mAb [JE34-75]	HA722236	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	CAPZA1 Recombinant Mouse mAb [A2E7-R] - BSA and Azide free	HA610155	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	CAPZA1 Recombinant Mouse mAb [A2E7-R]	HA601280	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	CAPZA1 Mouse mAb [A2E7]	EM1901-82	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	CAPZB Rabbit pAb	HA500486	Human, Mouse	WB
	CAPZB Rabbit pAb	HA500003	Human	WB
	CAR Rabbit pAb	ER62892	Human, Mouse, Rat	WB
Relative expression	Cardiac Troponin I Recombinant Rabbit mAb [JF0550]	ET1702-37	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Cardiac Troponin T Mouse mAb [5-C7]	M1302-1	Human, Mouse	WB,IF-Cell,IHC-P,FC
	Cardiac Troponin T Rabbit pAb	ER21102	Human, Mouse, Rat	WB,IHC-P
	Cardiac Troponin T Recombinant Rabbit mAb [SC64-03] - BSA and Azide free	HA750222	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Cardiac Troponin T Recombinant Rabbit mAb [SC64-03]	ET1610-51	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Cardiac Troponin T Mouse mAb [3-G6]	EM1701-39	Human, Mouse	WB,IHC-P,IF-Cell,IF-Tissue
Relative expression	CaSR Recombinant Rabbit mAb [PSH03-47]	HA722006	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	CaSR Recombinant Rabbit mAb [PSH03-47] - BSA and Azide free	HA750881	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	CaSR Rabbit pAb	ER62910	Human, Mouse, Rat	WB,IHC-P
	Cathepsin B Rabbit pAb	ER1905-48	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	Cathepsin B Mouse mAb [J11-A11]	M1506-1	Human, Zebrafish, Mouse, Rat	WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	Cathepsin B Recombinant Rabbit mAb [JA11-02]	ET1704-03	Human	WB,IF-Cell,IF-Tissue,IHC-P
Relative expression	Cathepsin B Rabbit pAb	R1407-1	Human	WB
	Cathepsin D Mouse mAb [13F4]	EM1901-17	Human	WB,IHC-P
	Cathepsin D Mouse mAb [13F2]	EM1901-15	Human	WB,IHC-P
	Cathepsin D Recombinant Mouse mAb [13F4-R]	EM1901-03	Human	WB,IHC-P
	Cathepsin D Recombinant Rabbit mAb [SU0360]	ET1608-49	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	Cathepsin D Rabbit pAb	R1601-2	Human, Mouse, Rat	WB,IHC-P,FC
	Cathepsin D Recombinant Rabbit mAb [JE24-31]	ET7110-58	Human	IF-Cell,IHC-P
	Cathepsin D Mouse mAb [13F3]	EM1901-16	Human	WB,IHC-P
	Cathepsin D Recombinant Mouse mAb [13F3-R]	EM1901-02	Human	WB,IHC-P
	Cathepsin D Rabbit pAb	R1512-9	Human, Mouse	WB,IF-Cell,IHC-P,FC
	Cathepsin D Recombinant mAb	IRS099RB	Human	mIHC
	Cathepsin D Recombinant Rabbit mAb [SU0360] - BSA and Azide free	HA750152	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	Cathepsin D Recombinant Mouse mAb [13F4-R] - BSA and Azide free	HA610002	Human	WB,IHC-P
	Cathepsin D Recombinant Mouse mAb [13F3-R] - BSA and Azide free	HA610001	Human	WB,IHC-P
	Cathepsin G Rabbit pAb	ER62914	Human	WB
	Cathepsin K Rabbit pAb	HA500469	Human	WB,IF-Cell
KD	Cathepsin L Recombinant Rabbit mAb [PSH03-92] - BSA and Azide free	HA750917	Human, Mouse, Rat, Monkey	WB,IHC-P
KD	Cathepsin L Recombinant Rabbit mAb [PSH03-92]	HA722064	Human, Mouse, Rat, Monkey	WB,IHC-P
KD	Cathepsin L Recombinant Rabbit mAb [PSH03-91]	HA722063	Human, Mouse, Rat, Monkey	WB,IHC-P,FC
KD	Cathepsin L Recombinant Rabbit mAb [PSH03-91] - BSA and Azide free	HA750916	Human, Mouse, Rat, Monkey	WB,IHC-P,FC
	Cathepsin L Rabbit pAb	ER1905-49	Human, Rat	WB,IHC-P,IF-Tissue,FC
KD	Cathepsin L/V/K/H Recombinant Rabbit mAb [JM10-78]	ET1703-44	Human, Mouse, Rat	WB,IF-Tissue,IHC-P
	CCDC39 Rabbit pAb	ER1905-59	Human	WB
	CCP2 Rabbit pAb	ER62928	Human, Mouse	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CCP3 Rabbit pAb	ER62929	Human	IHC-P
Relative expression	CD11a Recombinant Rabbit mAb [SC56-01]	ET1610-1	Human	WB,IF-Tissue,IHC-P,IP
Relative expression	CD11b Recombinant Rabbit mAb [PSH03-96]	HA722075	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	CD11b Recombinant Rabbit mAb [JU93-81] - BSA and Azide free	HA750451	Human	WB,IP,IHC-P,IF-Cell,IF-Tissue
Relative expression	CD11b Recombinant Rabbit mAb [PSH03-96] - BSA and Azide free	HA750921	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	CD11b Recombinant Rabbit mAb	IRS014	Human	mIHC
Relative expression	CD11b Mouse mAb [A2-9]	M1211-6	Human	WB,IHC-P,FC,IF-Tissue
Relative expression	CD11b Recombinant Rabbit mAb [JU93-81]	ET1706-04	Human	WB,IP,IHC-P,IF-Cell,IF-Tissue,mIHC
	CD11b Mouse mAb [4-6]	EM1701-42	Human	IHC-P,IF-Cell,FC,WB,IF-Tissue
	CD11b Mouse mAb [A2-B5-D2]	EM1902-35	Human, Mouse	WB
Relative expression	CD11b Recombinant Rabbit mAb [JE01-17]	HA722856	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue
	CD11b/ITAM Mouse mAb [4G2]	RT1107	Human, Mouse	WB,IP,IF,FC
	CD11c Recombinant Antibody - Rat IgG1 (Chimeric)	IRS103RT	Human	mIHC
	CD11c Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS102MS	Human	mIHC
Relative expression	CD11c Recombinant Rabbit mAb [PSH07-42]	HA722830	Mouse	WB,IHC-P,IHC-Fr,IF-Tissue
	CD11c Recombinant Rabbit mAb	IRS084	Mouse	mIHC
Relative expression	CD11c Recombinant Rabbit mAb [PSH07-42] - BSA and Azide free	HA751163	Mouse	WB,IHC-P,IHC-Fr,IF-Tissue
	CD11c Recombinant Rabbit mAb	IRS001	Human	mIHC
	CD11c Recombinant Antibody [SI19-06] - Mouse IgG1 (Chimeric)	HA601545	Human	IF-Tissue,IHC-P
	CD11c Recombinant Antibody [SI19-06] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610258		
	CD11c Recombinant Antibody [SI19-06] - Rat IgG1 (Chimeric)	HA601418		
	CD11c Recombinant Antibody [SI19-06] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610361	Human	IF-Tissue,IHC-P
	CD11c Recombinant Rabbit mAb [SI19-06]	ET1606-19	Human	WB,IHC-P,IP,mIHC,IF-Tissue
	CD11c / ITGAX Mouse mAb [3G1]	RT1108	Human	WB,IP,IF,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD146 Rabbit pAb	0806-5	Human, Mouse, Rat	WB,IP
KD	CD146 Recombinant Rabbit mAb [JB42-35]	ET7107-91	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	CD147 Rabbit pAb	1003-10	Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	CD147 Recombinant Rabbit mAb [JF1-045]	ET1702-58	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Relative expression	CD151 Mouse mAb [A6G11]	HA600059	Human, Mouse	WB,IF-Cell,IHC-P
	CD151 Mouse mAb [3G1]	RT1111	Human	WB,IP,IF,IHC-P
	CD156 Mouse mAb [2C1F3]	EM1711-73	Human	IHC-P,IF-Cell,FC
	CD164 Rabbit pAb	ER62936	Human, Mouse, Rat	IHC-P
KD	CD168 Recombinant Rabbit mAb [PSH02-16]	HA721775	Human	WB,IF-Cell,IHC-P,FC
KD	CD168 Recombinant Rabbit mAb [PSH02-16] - BSA and Azide free	HA750786	Human	WB,IF-Cell,IHC-P,FC
	CD18 Recombinant Rabbit mAb [JM71-31]	ET1703-94	Human	WB,IP
	CD299 Mouse mAb [8A1B3]	EM1711-85	Human	WB,FC
	CD299 Rabbit pAb	ER65752	Human, Mouse, Rat	WB,IHC-P
	CD31 Recombinant Rabbit mAb [SU03-59] - BSA and Azide free	HA750151	Human	WB,IHC-P,FC,IP,IF-Cell,IF-Tissue
	CD31 Recombinant Antibody [7-A1-R] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610349	Human	IHC-P
Relative expression	CD31 Recombinant Mouse mAb [7-A1-R]	HA601240	Human	WB,IF-Cell,IHC-P
Relative expression	CD31 Rabbit pAb	ER31219	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
	iFluor™ 488 Conjugated CD31 Recombinant Rabbit mAb [SU03-59]	HA720220F	Human	IF-Tissue
	CD31 Mouse mAb	IRS023	Human	mIHC
	CD31 Recombinant Rabbit mAb	IRS117RB	Mouse	mIHC
Relative expression	CD31 Mouse mAb [7-A1]	M1511-8	Human, Mouse	WB,IF-Cell,IHC-P,mIHC
Relative expression	CD31 Recombinant Mouse mAb [7-A1-R] - BSA and Azide free	HA610115	Human	WB,IF-Cell,IHC-P
	CD31 Recombinant Antibody - Rabbit IgG (Chimeric)	IRS116RB	Human	mIHC
Relative expression	CD31 Recombinant Antibody [7-A1-R] - Rabbit IgG (Chimeric) - BSA and Azide free	HA751615	Human	WB,IF-Cell,IHC-P,FC
	CD31 Recombinant Rabbit mAb	IRS075	Mouse	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD31 Recombinant Rabbit mAb [SU03-59]	ET1608-48	Human	WB,IHC-P,FC,IP,IF-Cell,IF-Tissue
Relative expression	CD31 Recombinant Antibody [7-A1-R] - Rabbit IgG (Chimeric)	HA723897	Human	WB,IF-Cell,IHC-P,FC
	CD31 Recombinant Antibody [7-A1-R] - Rat IgG1 (Chimeric)	HA601522	Human	IHC-P
	CD322 Mouse mAb [7E4C5]	EM1711-67	Human	WB,IF-Cell,FC
	CD329 Rabbit pAb	ER65523	Human	WB
	CD41 Recombinant Rabbit mAb [PSH08-01] - BSA and Azide free	HA751201	Mouse	IHC-P,FC,IF-Cell,Depletion
	CD41 Recombinant Rabbit mAb [PSH08-03]	HA722940	Human, Mouse, Rat	WB,IHC-P
	CD41 Recombinant Rabbit mAb [PSH08-03] - BSA and Azide free	HA751203	Human, Mouse, Rat	WB,IHC-P
	CD41 Recombinant Rabbit mAb [PSH08-02]	HA722939	Mouse	IHC-P,FC,IF-Cell,Depletion
	CD41 Recombinant Rabbit mAb [PSH08-02] - BSA and Azide free	HA751202	Mouse	IHC-P,FC,IF-Cell,Depletion
	CD41 Recombinant Rabbit mAb [JM11-19]	ET1703-37	Human	WB,IHC-P
	CD41 Recombinant Rabbit mAb [PSH08-01]	HA722938	Mouse	IHC-P,FC,IF-Cell,Depletion
Relative expression, KD	CD44 Recombinant Rabbit mAb [JE64-01]	HA722845	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,FC
	CD44 Recombinant Rabbit mAb	IRS123RB	Mouse	mIHC
Relative expression	CD44 Recombinant Rabbit mAb [ST57-03]	ET1609-74	Human	WB,IHC-P
	CD44 Recombinant mAb	IRS122RB	Human	mIHC
	CD44 Mouse mAb [A7-1]	EM010805	Human	WB,IHC-P,FC
Relative expression	CD44 Recombinant Rabbit mAb [ST57-03] - BSA and Azide free	HA750192	Human	WB,IHC-P
Relative expression	CD63 Recombinant Rabbit mAb [PSH06-64]	HA722731	Human	WB,IF-Cell,FC,IP,IHC-P
KD	CD63 Recombinant Rabbit mAb [SY21-02] - BSA and Azide free	HA750112	Human	WB,IHC-P
KD	CD63 Recombinant Rabbit mAb [SY21-02]	ET1607-2	Human	WB,IHC-P
Relative expression	CD63 Recombinant Rabbit mAb [PSH06-64] - BSA and Azide free	HA751088	Human	WB,IF-Cell,FC,IP,IHC-P
	CD99 Recombinant Rabbit mAb [PD00-13]	HA721150	Human	IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CD99 Rabbit pAb	ER1803-81	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	CD99 Recombinant Rabbit mAb [JF0991]	ET1702-35	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	CD99L2 Rabbit pAb	ER1906-03	Human	WB,IHC-P,IF-Tissue
	CDC42 Recombinant Rabbit mAb [JJ086-04] - BSA and Azide free	HA750320	Human, Mouse, Rat	WB,FC,IF-Cell
	CDC42 Recombinant Rabbit mAb [JJ086-04]	ET1701-7	Human, Mouse, Rat	WB,FC,IF-Cell
	CDC42EP1 Recombinant Rabbit mAb [SD08-54]	ET1612-95	Human	WB
	Cdc42EP5 Rabbit pAb	ER62955	Human, Mouse	WB,IHC-P
Cell treatment	CDCP1 Mouse mAb [A10F10]	EM1701-84	Human, Mouse	WB
	CDH11 Mouse mAb [A6-C12]	EM1707-03	Human	WB,IF-Cell,IHC-P,FC
	CDH17 Mouse mAb [D10-E1]	EM1706-54	Human	WB,IHC-P
	CDHF9 Rabbit pAb	ER62962	Human, Mouse, Rat	IF-Cell
	CDK5 Activator-binding C48 Rabbit pAb	ER62967	Human	IHC-P
	CENPR Rabbit pAb	ER61447	Human	WB
	ChemR23 Rabbit pAb	ER63003	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	CHST3 Rabbit pAb	ER1906-28	Human, Mouse	WB,IHC-P,IF-Tissue
	CHST7 Rabbit pAb	ER1906-30	Mouse	WB,IHC-P,IF-Tissue
	CILP Rabbit pAb	ER1906-32	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	Cingulin Recombinant Rabbit mAb [PSH06-19] - BSA and Azide free	HA751057	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,FC
Relative expression	Cingulin Recombinant Rabbit mAb [PSH06-19]	HA722516	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IP,FC
KD	Cip4 Recombinant Rabbit mAb [JE63-45]	HA721029	Human, Mouse, Rat	WB
KD	Claudin 1 Recombinant Rabbit mAb [PS01-43]	HA721248	Human	WB,IHC-P,IF-Tissue,IF-Cell
	Claudin 1 Mouse mAb [8G1]	RT1141	Human, Mouse, Rat	WB,IP,IF,IHC-P
	Claudin 1 Rabbit pAb	ER1906-37	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression, KD	Claudin 1 Recombinant Rabbit mAb [JE30-33]	HA721999	Human, Mouse, Rat	WB,IHC-P
KD	Claudin 1 Recombinant Rabbit mAb [PS01-43] - BSA and Azide free	HA750567	Human	WB,IHC-P,IF-Tissue,IF-Cell
	Claudin 10 Rabbit pAb	ER63022	Human, Mouse	WB,IHC-P,IF-Cell
	Claudin 15 Rabbit pAb	ER1906-39	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Claudin 17 Rabbit pAb	ER65463	Human	WB
	Claudin18 Rabbit pAb	HA500131	Human	WB,IHC-P
	Claudin 18.1 Rabbit pAb	ER1902-11	Human, Mouse	WB,FC
Relative expression	Claudin 18.2 Recombinant Rabbit mAb [PSH11-63]	HA723359	Human, Mouse, Rat	IHC-P
	Claudin 18.2 Rabbit pAb	ER1902-86	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	Claudin 18.2 Recombinant Rabbit mAb [PSH11-63] - BSA and Azide free	HA751419	Human, Mouse, Rat	IHC-P
	Claudin 18.2 Recombinant Mouse mAb [PSH18-20]	HA601432	Human	IHC-P
	Claudin 2 Mouse mAb [8H11]	EM60133	Human, Mouse, Rat	WB,IHC-P
	Claudin 3 Rabbit pAb	ER63024	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Claudin 4 Rabbit pAb	ER63025	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Claudin 5 Recombinant Rabbit mAb [JM11-22] - BSA and Azide free	HA750384	Human	WB,IHC-P
	Claudin 5 Recombinant Rabbit mAb [JM11-22]	ET1703-58	Human	WB,IHC-P
	Claudin 6 Rabbit pAb	ER63026	Human, Mouse	IHC-P
	Claudin 7 Recombinant Mouse mAb [6-G3-R]	HA601346	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,IP
	Claudin 7 Mouse mAb [6-G3]	M1412-2	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	Claudin 7 Recombinant Mouse mAb [6-G3-R] - BSA and Azide free	HA610206	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,IP
	Claudin 8 Rabbit pAb	ER63027	Human	WB
	CLD14 Rabbit pAb	ER60614	Human, Mouse	WB
	CLD20 Rabbit pAb	ER60615	Human	WB
	CLD9 Rabbit pAb	ER60618	Human, Mouse	WB
	Cleaved-Cathepsin G (I21) Rabbit pAb	ER60016	Human, Mouse, Rat	WB
	Cleaved-Cathepsin L2 (L114) Rabbit pAb	ER60017	Human, Monkey	WB
	Cleaved-COL4A3 (L1425) Rabbit pAb	ER60021	Human	WB
	Cleaved-MMP-10 (F99) Rabbit pAb	ER60028	Human	WB
	Cleaved-MMP-23 (Y79) Rabbit pAb	ER60029	Human, Mouse, Rat	WB
	Cleaved-MMP-27 (Y99) Rabbit pAb	ER60030	Human, Monkey	WB
	Cleaved-Tumstatin (P1426) Rabbit pAb	ER60036	Human, Monkey	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	CLECSF6 Rabbit pAb	ER65733	Human	WB
	CLIP170 Recombinant Rabbit mAb [JE63-82]	HA721004	Human	WB,IF-Cell,IHC-P
	CNN3 Rabbit pAb	ER60562	Human, Mouse, Rat	WB
	Cofilin Recombinant Rabbit mAb [PSH0-33] - BSA and Azide free	HA750612	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
	Cofilin Recombinant Rabbit mAb [PSH0-33]	HA721326	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
	Cofilin Recombinant Rabbit mAb [PSH01-44] - BSA and Azide free	HA750721	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	Cofilin Recombinant Rabbit mAb [PSH01-44]	HA721682	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	Cofilin Recombinant Rabbit mAb [PSH0-34]	HA721327	Human, Mouse, Rat, Monkey	WB,FC
	Cofilin Rabbit pAb	ER63062	Human, Mouse, Rat	WB,IHC-P
	Cofilin Rabbit pAb	HA500523	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
	Cofilin Recombinant Rabbit mAb [PSH0-38] - BSA and Azide free	HA750616	Human, Mouse, Rat, Monkey	WB,FC
	Cofilin Recombinant Rabbit mAb [PSH0-34] - BSA and Azide free	HA750613	Human, Mouse, Rat, Monkey	WB,FC
	Cofilin Recombinant Rabbit mAb [PSH0-38]	HA721340	Human, Mouse, Rat, Monkey	WB,FC
	COL11A1 Rabbit pAb	ER63063	Human, Mouse	WB,IHC-P,IF-Cell
	COL11A1 Recombinant Rabbit mAb [JE40-84]	HA721718	Human	WB
	COL11A2 Rabbit pAb	ER63064	Human, Mouse	WB,IHC-P
	COL12A1 Rabbit pAb	ER63065	Human	IHC-P,IF-Cell
	COL13A1 Rabbit pAb	ER65862	Human, Mouse	WB
	COL16A1 Rabbit pAb	ER63066	Human, Mouse	IHC-P
	COL19A1 Rabbit pAb	ER63067	Human	WB,IHC-P,IF-Cell
	COL1A1 Recombinant Rabbit mAb [ST58-04] - BSA and Azide free	HA750188	Human, Cow	WB,IHC-P
	COL1A1 Recombinant Rabbit mAb [ST58-04]	ET1609-68	Human, Cow	WB,IHC-P
	COL1A1 Recombinant Rabbit mAb [PSH06-20]	HA722517	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP,mIHC,IHC-Fr
	COL1A1 Recombinant mAb	IRS137RB	Human	mIHC
	COL1A1 Recombinant Rabbit mAb	IRS228RB	Mouse	mIHC
	COL1A1 Recombinant Rabbit mAb [PSH06-20] - BSA and Azide free	HA751058	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC,IP,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	COL1A2 Recombinant Rabbit mAb [PSH07-21]	HA722809	Human, Rat	WB,IF-Cell,FC
Relative expression	COL1A2 Recombinant Rabbit mAb [PSH07-21] - BSA and Azide free	HA751142	Human, Rat	WB,IF-Cell,FC
	COL23A1 Rabbit pAb	ER63068	Human, Mouse, Rat	IF-Cell
	COL4A1 Rabbit pAb	ER63070	Human, Mouse, Rat	WB
	COL4A2 Rabbit pAb	ER63071	Human, Mouse	WB,IHC-P,IF-Cell
	COL4A3 Rabbit pAb	ER63072	Human	IHC-P,IF-Cell
	COL4A5 Rabbit pAb	ER63074	Human, Mouse	IHC-P,IF-Cell,WB
	COL4A6 Rabbit pAb	ER63075	Human	WB,IHC-P,IF-Cell
	COL5A1 Recombinant Rabbit mAb [PSH10-08]	HA723177	Human	WB,IF-Cell,FC,IP
	COL5A1 Rabbit pAb	ER63076	Human	WB,IHC-P,IF-Cell
	COL5A1 Recombinant Rabbit mAb [PSH10-08] - BSA and Azide free	HA751338	Human	WB,IF-Cell,FC,IP
	COL5A2 Rabbit pAb	ER65613	Human	WB,IHC-P
	COL6A2 Rabbit pAb	ER63079	Human, Mouse, Monkey	WB,IHC-P,IF-Cell
	COL6A3 Rabbit pAb	ER63080	Human, Mouse	IHC-P,IF-Cell
	COL7A1 Rabbit pAb	ER63081	Human, Mouse	IHC-P,IF-Cell
	COL9A3 Rabbit pAb	ER63083	Human	IHC-P,IF-Cell
	Collagen alpha-1(XXVIII)pAb Rabbit Polyclonal	ER65722	Human, Mouse	WB
	Collagen II Rabbit pAb	ER1906-48	Human, Rabbit	WB,IHC-P,IF-Tissue
Relative expression	Collagen II Recombinant Rabbit mAb [PSH06-62]	HA722733	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
	Collagen II Rabbit pAb	ER1906-49	Human, Mouse	WB,IHC-P,IF-Tissue
Relative expression	Collagen II Recombinant Rabbit mAb [PSH06-62] - BSA and Azide free	HA751090	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
	Collagen III Rabbit pAb	ER1906-50	Human	WB,FC,IF-Cell
	Collagen III Recombinant Rabbit mAb [JE60-79]	HA720050	Human, Mouse, Rat	WB,IF-Cell,FC
	Collagen IV Rabbit pAb	HA500197	Human	WB,IF-Cell
	Collagen VI Recombinant Rabbit mAb [SD83-03]	ET1612-91	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Collagen VI Recombinant Rabbit mAb [SD83-03] - BSA and Azide free	HA750305	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	iFluor™ 488 Conjugated Collagen VI Recombinant Rabbit mAb [SD83-03]	HA720155F	Human, Mouse, Rat	IF-Cell,IF-Tissue
	Collagen X Recombinant Rabbit mAb [JF0961]	ET1702-13	Rat, Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Collagen X Recombinant Rabbit mAb [JF0961] - BSA and Azide free	HA750336	Rat, Human	WB
	Collagen XVII Recombinant mAb	IRS138RB	Mouse	mIHC
KD	Collagen XVII Recombinant Rabbit mAb [SR46-05] - BSA and Azide free	HA750030	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
KD	Collagen XVII Recombinant Rabbit mAb [SR46-05]	ET1602-14	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
	COLQ Rabbit pAb	ER65874	Human, Mouse, Rat	IHC-P, WB
Relative expression	COMP Mouse mAb [A10H2]	HA601354	Human, Mouse, Rat	WB, IHC-P
	COMP Rabbit pAb	ER1906-53	Human	WB
	Connexin-26 Rabbit pAb	ER1902-42	Human, Mouse, Rat	WB, IHC-P, FC
	Connexin-40 Rabbit pAb	ER1906-57	Human, Mouse	WB, IHC-P, IF-Tissue
	Connexin 43 Recombinant Rabbit mAb [PSH0-32] - BSA and Azide free	HA750606	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB, IHC-P, IHC-Fr
	Connexin 43 Recombinant Rabbit mAb [PSH0-32]	HA721312	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB, IHC-P, IHC-Fr
	Connexin 45 / GJA7 / Cx45 Rabbit pAb	ER1803-93	Human, Mouse, Rat	WB, IF-Cell, IHC-P
	Connexin 46 Rabbit pAb	ER63086	Human, Mouse, Rat	IHC-P
	Connexin 47 Rabbit pAb	ER63087	Human	WB, IF-Cell
	Contactin 1 Rabbit pAb	ER65455	Human, Mouse, Rat	WB
	Coronin 1A Rabbit pAb	ER63096	Human, Mouse, Rat	WB, IHC-P
KD	Contactin Recombinant Rabbit mAb [SC61-08]	ET1610-87	Human, Mouse, Rat	WB, IF-Cell, IF-Tissue, IHC-P, IP, FC
	CRDL2 Rabbit pAb	ER61936	Human, Mouse	WB
	C Reactive Protein Recombinant Rabbit mAb [JE03-72]	HA722687	Human	WB, IHC-P
	C Reactive Protein Recombinant Rabbit Multiclonal Antibody [PSH15-35]	HA723739	Human, Mouse, Rat	WB
	CRP Mouse mAb [E8-G1]	EM1709-89	Human	WB, IHC-P, FC
	CTHRC1 Mouse mAb [5H1]	EM1701-36	Human, Rat, Mouse	WB, IHC-P
	CTNNA3 Rabbit pAb	HA500192	Human, Mouse	WB, IHC-P
	CTNNA3 Recombinant Rabbit mAb [JE55-32]	ET7111-07	Human, Mouse, Rat	WB, IHC-P
	CTXL Rabbit pAb	ER65844	Human	WB
	CX3CR1 Recombinant Rabbit mAb	IRS250RB	Mouse	mIHC
	CX3CR1 Mouse mAb [A5C7]	HA600074	Human	IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Cx40 / GJA5 Recombinant Rabbit mAb [PSH10-82] - BSA and Azide free	HA751374	Human, Mouse, Rat, Monkey	WB,IP
	Cx40 / GJA5 Recombinant Rabbit mAb [PSH10-82]	HA723249	Human, Mouse, Rat, Monkey	WB,IP
Relative expression, Cell treatment	CYR61 / CCN1 Recombinant Rabbit mAb [PSH09-12]	HA723070	Human	WB,IF-Cell,IHC-P,IF-Tissue
	CYR61 / CCN1 Rabbit pAb	ER1905-70	Human, Mouse	WB,IHC-P,IF-Tissue
Relative expression, Cell treatment	CYR61 / CCN1 Recombinant Rabbit mAb [PSH09-12] - BSA and Azide free	HA751266	Human	WB,IF-Cell,IHC-P,IF-Tissue
	Cystatin A Rabbit pAb	ER63192	Human	WB,IHC-P,IF-Cell
	Cytokeratin 1 Recombinant Rabbit mAb [SN72-08]	ET1611-46	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	iFluor™ 594 Conjugated Cytokeratin 1 Recombinant Rabbit mAb [SN72-08]	HA720124F	Human, Mouse, Rat	IF-Tissue
	Cytokeratin 1 Recombinant Rabbit mAb	IRS163RB	Mouse	mIHC
	iFluor™ 488 Conjugated Cytokeratin 1 Recombinant Rabbit mAb [SN72-08]	HA720127F	Human, Mouse, Rat	IF-Tissue
	Cytokeratin 1 Recombinant Rabbit mAb [SN72-08] - BSA and Azide free	HA750256	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	iFluor™ 647 Conjugated Cytokeratin 1 Recombinant Rabbit mAb [SN72-08]	HA720141F	Human, Mouse, Rat	IF-Tissue
	iFluor™ 488 Conjugated Cytokeratin 10 Recombinant Rabbit mAb [ST05-43]	HA720134F	Human, Rat	IF-Tissue
KD	Cytokeratin 10 Recombinant Rabbit mAb [ST05-43] - BSA and Azide free	HA750193	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	iFluor™ 594 Conjugated Cytokeratin 10 Recombinant Rabbit mAb [ST05-43]	HA720119F	Human, Rat	IF-Tissue
	iFluor™ 647 Conjugated Cytokeratin 10 Recombinant Rabbit mAb [ST05-43]	HA720146F	Human, Rat, Mouse	IF-Tissue
KD	Cytokeratin 10 Recombinant Rabbit mAb [ST05-43]	ET1609-75	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Cytokeratin 13 Recombinant mAb	IRS133RB	Human	mIHC
Relative expression	Cytokeratin 13 Recombinant Rabbit mAb [SN71-09]	ET1611-55	Human, Mouse, Zebrafish	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	Cytokeratin 13 Recombinant Rabbit mAb [SN71-09] - BSA and Azide free	HA750261	Human, Mouse, Zebrafish	WB,IF-Cell,IHC-P,FC
	Cytokeratin 13 Rabbit pAb	R1601-8	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	iFluor™ 488 Conjugated Cytokeratin 14 Recombinant Rabbit mAb [SC65-06]	HA720135F	Mouse, Rat, Human	IF-Tissue,IF-Cell
Relative expression	Cytokeratin 14 Recombinant Mouse mAb [A2C11-R]	HA601259	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	Cytokeratin 14 Recombinant Mouse mAb [A2C10-R]	HA601323	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	Cytokeratin 14 Recombinant Rabbit mAb [SC65-06] - BSA and Azide free	HA750216	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	Cytokeratin 14 Mouse mAb [A2C11]	EM1901-32	Human	WB,IHC-P,IF-Cell
	Cytokeratin 14 Mouse mAb [A2C12]	EM1901-31	Human	WB,IHC-P,FC
	Cytokeratin 14 Recombinant Rabbit mAb	IRS162RB	Human	mIHC
Relative expression	Cytokeratin 14 Mouse mAb [A2C10]	EM1901-33	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
	iFluor™ 594 Conjugated Cytokeratin 14 Recombinant Rabbit mAb [SC65-06]	HA720117F	Rat, Human, Mouse	IF-Tissue
	iFluor™ 647 Conjugated Cytokeratin 14 Recombinant Rabbit mAb [SC65-06]	HA720147F	Human, Mouse, Rat	IF-Cell,IF-Tissue
	Cytokeratin 14 Rabbit pAb	ER1901-08	Human, Rat	WB,IHC-P
Relative expression	Cytokeratin 14 Recombinant Rabbit mAb [SC65-06]	ET1610-42	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	iFluor™ 488 Conjugated Cytokeratin 15 Recombinant Rabbit mAb [ST04-85]	HA720136F	Human, Mouse, Rat	IF-Tissue,IF-Cell,FC
	Cytokeratin 15 Recombinant Rabbit mAb [ST04-85]	ET1609-54	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell
	iFluor™ 647 Conjugated Cytokeratin 15 Recombinant Rabbit mAb [ST04-85]	HA720148F	Human, Mouse, Rat	IF-Tissue
	Cytokeratin 15 Recombinant Rabbit mAb [ST04-85] - BSA and Azide free	HA750182	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell
	iFluor™ 594 Conjugated Cytokeratin 15 Recombinant Rabbit mAb [ST04-85]	HA720122F	Human, Mouse, Rat	IF-Tissue
	iFluor™ 594 Conjugated Cytokeratin 16 Recombinant Rabbit mAb [SC52-09]	HA720114F	Human	IF-Cell,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	iFluor™ 488 Conjugated Cytokeratin 16 Recombinant Rabbit mAb [SC52-09]	HA720137F	Human	IF-Tissue,IF-Cell,FC
	Cytokeratin 16 Recombinant Rabbit mAb [SC52-09] - BSA and Azide free	HA750203	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Cytokeratin 16 Recombinant Rabbit mAb [SC52-09]	ET1610-17	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	iFluor™ 647 Conjugated Cytokeratin 16 Recombinant Rabbit mAb [SC52-09]	HA720149F	Human	IF-Tissue
	Cytokeratin 17 Rabbit pAb	0407-4	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Cytokeratin 17 Mouse mAb [83-13]	M0407-13	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Cytokeratin 17 Recombinant Rabbit mAb [SR45-06]	ET1602-16	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Cytokeratin 17 Recombinant Rabbit mAb [SR45-06] - BSA and Azide free	HA750031	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Cytokeratin 17 Rabbit pAb	ER1901-18	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	iFluor™ 488 Conjugated Cytokeratin 17 Recombinant Rabbit mAb [SR45-06]	HA720138F	Human, Mouse	IF-Tissue,IF-Cell
Relative expression	Cytokeratin 17 Mouse mAb [A2B9]	EM1901-28	Human	WB,IHC-P,IF-Cell
Relative expression	Cytokeratin 17 Recombinant Mouse mAb [A2B10-R] - BSA and Azide free	HA610110	Human	WB,IF-Cell,IHC-P
Relative expression	Cytokeratin 17 Mouse mAb [A2B10]	EM1901-29	Human	WB,IHC-P,FC,IF-Cell
Relative expression	Cytokeratin 17 Recombinant Mouse mAb [A2B10-R]	HA601235	Human	WB,IF-Cell,IHC-P
	Cytokeratin 17 Mouse mAb [A2B11]	EM1901-30	Human	WB,IHC-P,FC
	iFluor™ 594 Conjugated Cytokeratin 17 Recombinant Rabbit mAb [SR45-06]	HA720116F	Human, Mouse	IF-Tissue,IF-Cell
	iFluor™ 647 Conjugated Cytokeratin 17 Recombinant Rabbit mAb [SR45-06]	HA720150F	Human, Mouse, Rat	IF-Cell,IF-Tissue
	Cytokeratin 17 Rabbit pAb	ER1901-01	Human	WB
	Cytokeratin 17 Rabbit pAb	ER1803-94	Human, Mouse, Rat	WB,IHC-P,FC
	Cytokeratin 18 Recombinant Mouse mAb [6-19-R] - BSA and Azide free	HA610144	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Cytokeratin 18 Rabbit pAb	0407-1	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	Cytokeratin 18 Recombinant Rabbit mAb [SZ80-07] - BSA and Azide free	HA750069	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	Cytokeratin 18 Recombinant Rabbit mAb [SZ80-07]	ET1603-8	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,mIHC
	Cytokeratin 18 Mouse mAb [A4B3]	HA600021	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	Cytokeratin 18 Mouse mAb [6-19]	M0407-19	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Cytokeratin 18 Recombinant Rabbit mAb	IRS224RB	Mouse	mIHC
	iFluor™ 488 Conjugated Cytokeratin 18 Recombinant Rabbit mAb [SZ80-07]	HA720139F	Human, Mouse, Rat	IF-Tissue,IF-Cell
	Cytokeratin 18 Recombinant Rabbit mAb	IRS028	Human	mIHC
	Cytokeratin 18 Mouse mAb [A4B4]	HA600020	Human	WB,IHC-P,FC
	iFluor™ 647 Conjugated Cytokeratin 18 Recombinant Rabbit mAb [SZ80-07]	HA720151F	Human, Mouse, Rat	IF-Tissue
	Cytokeratin 18 Recombinant Mouse mAb [A4B3-R]	HA601260	Human, Mouse, Rat	WB,IHC-P
	Cytokeratin 18 Recombinant Mouse mAb [6-19-R]	HA601269	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Cytokeratin 18 Recombinant Mouse mAb [A4B3-R] - BSA and Azide free	HA610135	Human, Mouse, Rat	WB,IHC-P
	Cytokeratin 19 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS226RT	Mouse	mIHC
	Cytokeratin 19 Recombinant Antibody [SA30-06] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610370	Human, Rat	WB,IF-Cell,IHC-P,FC
	Cytokeratin 19 Rabbit pAb	ER1803-79	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	Cytokeratin 19 Recombinant Mouse mAb [A3D1-R] - BSA and Azide free	HA610133	Human	WB,IF-Cell,IHC-P,FC
KD	Cytokeratin 19 Recombinant Mouse mAb [A3D1-R]	HA601258	Human	WB,IF-Cell,IHC-P,FC
KD	Cytokeratin 19 Mouse mAb [A3D1]	EM1901-75	Human	WB,IHC-P,IF-Cell,FC
KD	Cytokeratin 19 Mouse mAb [A3D5]	EM1901-76	Human	WB,IHC-P,IF-Cell,FC
KO	Cytokeratin 19 Recombinant Rabbit mAb [SA30-06]	ET1601-6	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,mIHC
	iFluor™ 594 Conjugated Cytokeratin 19 Recombinant Rabbit mAb [SA30-06]	HA720115F	Human, Mouse	IF-Cell,IF-Tissue,FC
	Cytokeratin 19 Recombinant Antibody [SA30-06] - Rat IgG1 (Chimeric)	HA601555	Human, Rat	WB,IF-Cell,IHC-P,FC
	Cytokeratin 19 Recombinant Rabbit mAb	IRS225RB	Mouse	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Cytokeratin 19 Recombinant Antibody - Rat IgG1 (Chimeric)	IRS134RT	Human	mIHC
	iFluor™ 488 Conjugated Cytokeratin 19 Recombinant Rabbit mAb [SA30-06]	HA720140F	Human, Mouse	IF-Cell,IF-Tissue,FC
	Cytokeratin 19 Recombinant Rabbit mAb	IRS011	Human	mIHC
	iFluor™ 647 Conjugated Cytokeratin 19 Recombinant Rabbit mAb [SA30-06]	HA720152F	Human, Mouse	IF-Cell,IF-Tissue,FC
	Cytokeratin 20 Recombinant Rabbit mAb [SA35-03]	ET1601-8	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	Cytokeratin 20 Recombinant Rabbit mAb [JE20-39]	ET7110-54	Human, Mouse	WB,IHC-P
	Cytokeratin 20 Rabbit pAb	ER1901-97	Human, Mouse	WB,IHC-P,FC
	Cytokeratin 20 Mouse mAb [A3H7]	EM1901-97	Human	IHC-P,IF-Tissue
	Cytokeratin 20 Recombinant Rabbit mAb [SA35-03] - BSA and Azide free	HA750026	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	Cytokeratin 20 Rabbit pAb	ER1901-76	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Cytokeratin 20 Recombinant Mouse mAb [A3H7-R]	HA601077	Human	IHC-P,IF-Tissue
	Cytokeratin 20 Recombinant mAb	IRS135RB	Human	mIHC
	iFluor™ 488 Conjugated Cytokeratin 20 Recombinant Rabbit mAb [SA35-03]	HA720189F	Human	IF-Cell,IF-Tissue
	Cytokeratin 20 Mouse mAb [A3D12]	EM1901-96	Human	WB,FC,IHC-P
	Cytokeratin 2e Recombinant Rabbit mAb [JJ201-8]	ET1701-13	Human	WB,IHC-P
	iFluor™ 594 Conjugated Cytokeratin 4 Recombinant Rabbit mAb [SN74-03]	HA720123F	Human	IF-Tissue
	iFluor™ 488 Conjugated Cytokeratin 4 Recombinant Rabbit mAb [SN74-03]	HA720128F	Human	IF-Tissue,IF-Cell
	Cytokeratin 4 Recombinant Rabbit mAb [SN74-03]	ET1611-71	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	iFluor™ 647 Conjugated Cytokeratin 4 Recombinant Rabbit mAb [SN74-03]	HA720142F	Human	IF-Cell,IF-Tissue
	Cytokeratin 4 Recombinant Rabbit mAb [SN74-03] - BSA and Azide free	HA750271	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	Cytokeratin 4 Recombinant Rabbit mAb [SN2001]	ET1611-91	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	iFluor™ 488 Conjugated Cytokeratin 5 Recombinant Rabbit mAb [SC62-04]	HA720129F	Human	IF-Tissue,IF-Cell
	iFluor™ 594 Conjugated Cytokeratin 5 Recombinant Rabbit mAb [SC62-04]	HA720121F	Human	IF-Tissue
	Cytokeratin 5 Recombinant Rabbit mAb	IRS249RB	Human	mIHC
	Cytokeratin 5 Recombinant Rabbit mAb [SY88-19]	ET1605-17	Human, Mouse	WB,IHC-P
Relative expression, KD	Cytokeratin 5 Recombinant Rabbit mAb [SC62-04]	ET1610-43	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,FC
	Cytokeratin 5 Mouse mAb [6D-4D]	EM1706-36	Human, Mouse, Monkey	WB,IHC-P,IF-Cell,FC
	iFluor™ 647 Conjugated Cytokeratin 5 Recombinant Rabbit mAb [SC62-04]	HA720126F	Human	IF-Cell,IF-Tissue
	Cytokeratin 5+6 Mouse mAb [A2A12]	EM1901-07	Human, Rat	WB,IHC-P,FC
	Cytokeratin 5+6 Rabbit pAb	ER1901-03	Human, Mouse, Rat	WB,IHC-P
	Cytokeratin 5+6 Recombinant Rabbit mAb [PSH0-68] - BSA and Azide free	HA750645	Human	WB,IF-Cell,FC
	Cytokeratin 5+6 Mouse mAb [A2B2]	EM1901-08	Human, Rat	WB,IHC-P,FC
	Cytokeratin 5+6 Recombinant Mouse mAb [A2B2-R]	HA601233	Human, Rat	WB,IHC-P
	Cytokeratin 5+6 Rabbit pAb	HA500508	Human	WB,IHC-P
	Cytokeratin 5+6 Recombinant Mouse mAb [A2B2-R] - BSA and Azide free	HA610108	Human, Rat	WB,IHC-P
	Cytokeratin 5+6 Recombinant Rabbit mAb [PSH0-67]	HA721455	Human, Mouse, Rat	IHC-P,IF-Tissue,WB
	Cytokeratin 5+6 Rabbit pAb	ER1901-86	Human, Mouse, Rat	WB,IHC-P,FC
	Cytokeratin 5+6 Recombinant Rabbit mAb [PSH0-68]	HA721456	Human	WB,IF-Cell,FC
	Cytokeratin 6 Recombinant Rabbit mAb [SN71-07]	ET1611-70	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	iFluor™ 488 Conjugated Cytokeratin 6 Recombinant Rabbit mAb [SN71-07]	HA720130F	Human, Rat	IF-Cell
	Cytokeratin 6 Recombinant Rabbit mAb [SN71-07] - BSA and Azide free	HA750270	Human, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	iFluor™ 647 Conjugated Cytokeratin 6 Recombinant Rabbit mAb [SN71-07]	HA720143F	Human, Rat	IF-Cell,IF-Tissue
	iFluor™ 594 Conjugated Cytokeratin 7 Recombinant Rabbit mAb [ST50-05]	HA720120F	Human, Mouse	IF-Cell,IF-Tissue,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Cytokeratin 7 Recombinant Mouse mAb [3-G3-D8-R]	HA601065	Human	WB,IHC-P,IF-Cell,FC,IF-Tissue
	Cytokeratin 7 Rabbit pAb	ER1803-89	Human, Mouse, Rat	WB,IHC-P
KD	Cytokeratin 7 Rabbit pAb	R1309-4	Human	WB,IF-Cell,IHC-P,FC
KO	Cytokeratin 7 Mouse mAb [3-G3-D8]	EM0702	Human	WB,IHC-P,IF-Cell,IF-Tissue
KD, Relative expression	Cytokeratin 7 Recombinant Rabbit mAb [ST50-05]	ET1609-62	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,mIHC
	iFluor™ 647 Conjugated Cytokeratin 7 Recombinant Rabbit mAb [ST50-05]	HA720144F	Human, Mouse	IF-Cell,IF-Tissue,FC
	iFluor™ 488 Conjugated Cytokeratin 7 Recombinant Rabbit mAb [ST50-05]	HA720131F	Human, Mouse	IF-Tissue,IF-Cell,FC
	Cytokeratin 7 Recombinant Rabbit mAb	IRS029	Human	mIHC
KD, Relative expression	Cytokeratin 7 Recombinant Rabbit mAb [ST50-05] - BSA and Azide free	HA750185	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Cytokeratin 8 Rabbit pAb	ER1802-43	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	Cytokeratin 8 Recombinant Rabbit mAb [SU0338]	ET1608-32	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Cytokeratin 8 Rabbit pAb	ER1803-87	Human	WB,IF-Cell,IHC-P,FC
	iFluor™ 647 Conjugated Cytokeratin 8 Recombinant Rabbit mAb [SU0338]	HA720145F	Human, Mouse	IF-Cell,IF-Tissue
	iFluor™ 488 Conjugated Cytokeratin 8 Recombinant Rabbit mAb [SU0338]	HA720132F	Human	IF-Tissue,IF-Cell
KD	Cytokeratin 8 Recombinant Rabbit mAb [SU0338] - BSA and Azide free	HA750145	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Cytokeratin 8 Mouse mAb [A1-B11]	M1603-2	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Cytokeratin 8 Rabbit pAb	ER1706-92	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	iFluor™ 594 Conjugated Cytokeratin 8 Recombinant Rabbit mAb [SU0338]	HA720118F	Human, Mouse	IF-Cell,IF-Tissue,FC
	Cytokeratin 80 Rabbit pAb	ER1911-82	Mouse	WB
KD	Cytokeratin 9 Recombinant Rabbit mAb [SD201-09]	ET1612-77	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,IP
	iFluor™ 488 Conjugated Cytokeratin 9 Recombinant Rabbit mAb [SD201-09]	HA720133F	Human, Mouse	IF-Tissue
KD	Cytokeratin 9 Recombinant Rabbit mAb [SD201-09] - BSA and Azide free	HA750302	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Cytokeratin (Pan) Mouse mAb [7H8C4]	EM1712-42	Human	WB,IHC,IF-Cell,FC
KD	D4 GDI Rabbit pAb	HA500217	Human, Mouse, Rat	WB,IF-Cell,FC
	DAAM1 Rabbit pAb	ER65757	Human, Mouse, Rat	WB
	DAG1 Recombinant Rabbit mAb [PSH13-91] - BSA and Azide free	HA751493	Human, Mouse, Rat	WB
	DAG1 Recombinant Rabbit mAb [PSH13-91]	HA723584	Human, Mouse, Rat	WB
	DCBLD2 / ESDN Rabbit pAb	HA500585	Human	WB
	DC-SIGN Rabbit pAb	ER1907-23	Human	WB,IHC-P,IF-Tissue,FC
	DC-SIGN Recombinant Rabbit mAb [JE31-24]	HA721382	Human	WB
	DCTN1 / p150-glued Recombinant Rabbit mAb [PSH13-45] - BSA and Azide free	HA751464	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Cell,FC,IF-Tissue,IP
	DCTN1 / p150-glued Recombinant Rabbit mAb [PSH13-45]	HA723534	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Cell,FC,IF-Tissue,IP
	DCTN4 Mouse mAb [D7-G9]	EM1708-82	Human	WB,IF-Cell,IHC-P,FC
	DCTN5 Recombinant Rabbit mAb [JE58-30]	HA722706	Human, Mouse, Rat	WB,IHC-P,IP
	Decorin Recombinant Rabbit mAb [PSH12-99] - BSA and Azide free	HA751454	Human, Mouse, Rat	WB,IHC-P
	Decorin Recombinant Rabbit mAb [PSH12-99]	HA723504	Human, Mouse, Rat	WB,IHC-P
	delta 1 Catenin/CAS Rabbit pAb	ER1901-07	Human, Mouse	WB,IF-Cell,IHC-P
	delta 1 Catenin/CAS Recombinant Rabbit mAb [PD00-43]	HA721193	Human	WB,IHC-P,IF-Tissue
	delta 1 Catenin/CAS Rabbit pAb	ER1901-80	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	delta 1 Catenin/CAS Rabbit pAb	ER1803-71	Human, Rat	WB,IF-Cell
KD	delta 1 Catenin/CAS Rabbit pAb	ER1901-92	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Delta 1 Tubulin Rabbit pAb	ER1907-27	Mouse	WB
Relative expression	Delta-catenin Recombinant Rabbit mAb [JE75-34]	HA722438	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IP,IHC-Fr
	Delta-catenin Rabbit pAb	HA500428	Human, Mouse, Rat	WB
	delta Sarcoglycan Recombinant Rabbit mAb [JM61-10]	ET1703-88	Human	WB,IP,IHC-P
	Desmin Recombinant Rabbit mAb [SI18-00]	ET1606-30	Human, Mouse, Rat, Zebrafish, Chicken	WB,IF-Cell,IF-Tissue,IHC-P,FC,mIHC
KO	Desmin Mouse mAb [3-F7]	M1501-10	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Desmin Recombinant Rabbit mAb [SI18-00] - BSA and Azide free	HA750100	Human, Mouse, Rat, Zebrafish, Chicken	WB,IF-Cell,IF-Tissue,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Desmin Recombinant mAb	IRS140RB	Mouse	mIHC
	Desmin Rabbit pAb	ER40113	Human, Mouse, Zebrafish	WB,IF-Cell,IHC-P
	Desmocollin 1 Recombinant Rabbit mAb [JE66-57]	HA721418	Mouse, Rat	WB
Relative expression	Desmocollin 3 Recombinant Rabbit mAb [JE55-86]	ET7111-26	Human, Mouse, Rat	WB
	Desmoglein 2 Recombinant Rabbit mAb [JB96-37]	ET7107-65	Human, Mouse, Rat	WB,IHC-P,FC
	Desmoglein 3 Rabbit pAb	ER1901-81	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Desmoglein 3 Recombinant Rabbit mAb [JE40-30]	HA722930	Human, Mouse, Rat	WB,IHC-P
	Desmoplakin Recombinant Rabbit mAb [JE51-43]	HA722560	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	DIAP2 Rabbit pAb	ER65875	Human, Mouse	IHC-P,WB
	DIAPH1 Recombinant Rabbit mAb [JE40-59]	ET7109-11	Human	WB,IF-Cell,IF-Tissue,IHC-P
Relative expression	DIAPH3 Rabbit pAb	ER1802-54	Human	WB,IHC-P
	DLRB1 Rabbit pAb	ER60676	Human, Mouse, Rat	WB
	DNJB6 Rabbit pAb	ER61577	Human, Mouse, Rat	WB
	DOCK 2 Rabbit pAb	ER63277	Human, Mouse, Rat	IHC-P
Relative expression	DOCK2 Recombinant Rabbit mAb [JE53-62]	HA721993	Human	WB
	DPP6 Rabbit pAb	ER1902-07	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Drebrin Rabbit pAb	ER1902-55	Human	WB,IHC-P,FC
Relative expression	Drebrin Recombinant Rabbit mAb [JE30-20]	HA722062	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	DSC2 Rabbit pAb	ER62009	Human, Mouse	WB
	DSG1 Recombinant Rabbit mAb [JB19-46]	ET7107-83	Human, Mouse	WB,IHC-P
	DSPP Rabbit pAb	ER1907-81	Human	WB,IHC-P,IF-Tissue
	DSTN Recombinant Rabbit mAb [JE36-06]	HA722397	Human, Mouse, Rat	WB,IHC-P
	DYLT1 Rabbit pAb	ER60677	Human	WB
	DYLT3 Rabbit pAb	ER60678	Human	WB
	Dynein IC2 Rabbit pAb	ER63298	Human, Mouse, Rat	WB
	Dynein intermediate chain 1 Recombinant Rabbit mAb [JM11-38]	ET1705-8	Mouse, Human, Rat	WB,IHC-P
	Dynein LC 2 Rabbit pAb	ER63300	Human, Mouse, Rat	IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	DYNLL1 Recombinant Rabbit mAb [SD08-04]	ET1612-64	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
Relative expression	Dystrophin Recombinant Rabbit mAb [JF1-022]	ET1702-98	Human, Mouse, Rat	WB,IHC-P,IHC-Fr
	EB2 Rabbit pAb	ER1907-89	Human, Mouse	WB,IHC-P,IF-Tissue
	EB3 Recombinant Rabbit mAb [JU38-30]	ET7106-99	Human, Mouse, Rat	WB,IP
	EBP50 Recombinant Rabbit mAb [JJ0946]	ET1701-96	Human, Rat	WB,IHC-P
	E-Cadherin Recombinant Rabbit mAb	IRS030	Human	mIHC
	E-Cadherin Mouse mAb [3-F9]	M1405-3	Human, Mouse	IHC-P
	E-Cadherin Recombinant Mouse mAb	IRS031	Human	mIHC
	E-Cadherin Recombinant Mouse mAb	IRS229MS	Mouse	mIHC
	iFluor™ 594 Conjugated E-Cadherin Recombinant Rabbit mAb [SY0287]	HA720164F	Human	IF-Tissue
Relative expression	E-Cadherin Rabbit pAb	ER63312	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IHC-Fr,FC
Relative expression	E-Cadherin Recombinant Rabbit mAb [PSH13-75]	HA723564	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Cell,FC,IP
Relative expression, KD	E-Cadherin Recombinant Rabbit mAb [SY0287]	ET1607-75	Human, Cynomolgus monkey	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP,mIHC
	E-Cadherin Recombinant Mouse mAb [A0-G11-2-R] - BSA and Azide free	HA610047	Human, Mouse, Rat	WB,IHC
	E-Cadherin Recombinant Mouse mAb [A0-G11-2-R]	HA601143	Human, Mouse, Rat	WB,IHC,mIHC
	iFluor™ 488 Conjugated E-Cadherin Recombinant Rabbit mAb [SY0287]	HA720159F	Human	IF-Cell,FC
	E-Cadherin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS284RT	Mouse	mIHC
KD	E-Cadherin Mouse mAb [A0-G11-2]	EM0502	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	E-Cadherin Rabbit pAb	0407-25	Human, Mouse	WB,IF-Cell,FC
Relative expression	E-Cadherin Recombinant Rabbit mAb [PSH13-75] - BSA and Azide free	HA751488	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Cell,FC,IP
Relative expression, KD	E-Cadherin Recombinant Rabbit mAb [SY0287] - BSA and Azide free	HA750133	Human, Cynomolgus monkey	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	EDG1 Recombinant Rabbit mAb [JM10-66] - BSA and Azide free	HA750376	Human, Mouse, Rat	WB,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	EDG1 Recombinant Rabbit mAb [JM10-66]	ET1703-27	Human, Mouse, Rat	WB,IP
	Eg5 Recombinant Rabbit mAb [JU83-16]	ET1706-03	Human	WB,IP
	Eg5 Rabbit pAb	ER1803-25	Human, Mouse, Rat	IF-Cell,IHC-P,FC
	Elastin Rabbit pAb	ER1908-02	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	Elastin Recombinant Rabbit mAb [PSH16-13] - BSA and Azide free	HA751591	Human	IHC-P,IF-Tissue
	Elastin Recombinant Rabbit mAb [PSH16-13]	HA723818	Human	IHC-P,IF-Tissue
	ELMO1 Recombinant Rabbit mAb [JG34-59]	ET7108-06	Human, Mouse, Rat	WB,IHC-P,IP
	ELMO2 Rabbit pAb	ER1919-56	Human, Mouse	WB,IHC-P
Relative expression, KD	Emerin Recombinant Rabbit mAb [PSH01-86]	HA721744	Human	WB,IHC-P,IF-Cell,FC
	Emerin Recombinant Rabbit mAb [PSH01-85] - BSA and Azide free	HA750759	Human	WB,IHC-P,IF-Cell,FC
Relative expression, KD	Emerin Recombinant Rabbit mAb [PSH01-86] - BSA and Azide free	HA750760	Human	WB,IHC-P,IF-Cell,FC
	Emerin Recombinant Rabbit mAb [PSH01-85]	HA721743	Human	WB,IHC-P,IF-Cell,FC
	ENAH Rabbit pAb	HA500168	Human, Mouse	WB,IHC-P,FC
	ENAH Rabbit pAb	HA500026	Human, Mouse, Rat	WB,IF-Cell,FC
	Endomucin Rabbit pAb	IRS141RB	Mouse	mIHC
	Endomucin Rabbit pAb	ER1908-09	Human, Mouse	IHC-P,IF-Tissue
	Endostatin Rabbit pAb	ER1908-10	Human, Mouse	WB,IHC-P,IF-Tissue
	ENPP1 Rabbit pAb	ER1908-13	Human	WB
Relative expression	ENPP1 / PC1 Recombinant Rabbit mAb [PSH12-11] - BSA and Azide free	HA751440	Human	WB,IHC
Relative expression	ENPP1 / PC1 Recombinant Rabbit mAb [PSH12-11]	HA723427	Human	WB,IHC
	Eps8L1 Rabbit pAb	ER63385	Human, Mouse	WB,IHC-P,IF-Cell
	Eps8L2 Rabbit pAb	ER63386	Human, Mouse	WB,IHC-P,IF-Cell
	EPYC Rabbit pAb	ER61950	Human, Mouse	WB
KD	ERK1 Recombinant Rabbit mAb [SP05-09]	ET1604-16	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	ERK1 Mouse mAb [7-D6-E5]	M1505-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	ERK1 Recombinant Rabbit mAb [SP00-74]	ET1604-32	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	ERK1/2 Rabbit pAb	ER131011	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03]	ET1601-29	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03] - BSA and Azide free	HA750019	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
	Ermin Rabbit pAb	ER1908-31	Mouse, Rat	WB,IHC-P,IF-Tissue
	EZR Mouse mAb [F1-9A]	EM1706-9	Human, Mouse, Monkey	WB,IF-Cell,FC
	Ezrin Recombinant Rabbit mAb [JM10-65]	ET1703-25	Human	WB,IHC-P
	FAK Rabbit pAb	RT1215	Human, Mouse, Rat	WB,IP,ICC,IHC-P
KD	FAK Recombinant Rabbit mAb [SR46-04] - BSA and Azide free	HA750039	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	FAK Recombinant Rabbit mAb [SR46-04]	ET1602-25	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	FAM48A Rabbit pAb	ER65743	Human, Mouse, Rat	WB
	FAM83H Rabbit pAb	ER1908-67	Human, Mouse, Rat	IHC-P,IF-Tissue
KD	Fascin Recombinant Rabbit mAb [JM12-53]	ET1705-18	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	Fascin Rabbit pAb	ER1901-71	Human, Mouse, Rat	WB,IHC-P,FC
Relative expression	Fast Myosin Skeletal Heavy chain Recombinant Rabbit mAb [PSH05-12]	HA722228	Human, Mouse, Rat	IHC-P,IF-Tissue
Relative expression	Fast Myosin Skeletal Heavy chain Recombinant Rabbit mAb [PSH05-12] - BSA and Azide free	HA750985	Human, Mouse, Rat	IHC-P,IF-Tissue
	Fast Myosin Skeletal Heavy chain Recombinant mAb	IRS145RB	Mouse	mIHC
	FAT1 Rabbit pAb	0905-4	Human	WB,IHC-P
	FAT tumor suppressor homolog 4 Rabbit pAb	0905-1	Human	WB,IF-Cell
	FBN1 Rabbit pAb	ER63417	Human, Mouse, Rat	IHC-P,WB
	Fetuin A Recombinant Rabbit mAb [JE54-96]	HA721436	Human	WB,IHC-P,FC
	Fetuin A/AHSG Rabbit pAb	ER1908-80	Human, Mouse	WB
Relative expression	FGF23 Recombinant Rabbit mAb [PSH19-00]	HA724036	Human, Mouse, Rat	WB
	FGF23 Rabbit pAb	ER1908-82	Human, Rat	WB
	FHOD1 Rabbit pAb	ER1908-85	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Fibrinogen Rabbit pAb	ER1802-26	Human, Mouse	WB,IF-Cell,FC,IHC-P
	Fibrinogen Mouse mAb [1-A7]	M1510-13	Human	WB,IHC-P,ELISA
	Fibrinogen Mouse mAb [D5-B3]	M1510-11	Human	IHC-P,WB
	Fibrinogen Mouse mAb [11-B3-D8]	M1510-12	Human	WB,IHC-P,ELISA
	Fibrinogen alpha chain Recombinant Rabbit mAb [JE37-07]	HA722860	Human, Mouse, Rat	WB,IHC-P
	Fibronectin Mouse mAb [3G4]	RT1224	Human, Mouse, Rat	WB,IP,IF,IHC-P
	Fibronectin Recombinant Rabbit mAb [PSH14-21]	HA723610	Human, Mouse	WB,IHC-P,IF-Cell,FC
	Fibronectin Recombinant mAb	IRS144RB	Mouse	mIHC
KD	Fibronectin Recombinant Rabbit mAb [JF0582] - BSA and Azide free	HA750340	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
	Fibronectin Recombinant Rabbit mAb [PSH14-21] - BSA and Azide free	HA751514	Human, Mouse	WB,IHC-P,IF-Cell,FC
KD	Fibronectin Recombinant Rabbit mAb [JF0582]	ET1702-25	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
	Fibronectin Recombinant Rabbit mAb [PSH14-22]	HA723611	Human, Mouse	WB,IF-Cell,FC
	Fibronectin Recombinant Rabbit mAb [PSH14-22] - BSA and Azide free	HA751515	Human, Mouse	WB,IF-Cell,FC
	Fibulin-2 Rabbit pAb	ER63433	Human	WB,IHC-P,IF-Cell
	Fibulin-4 Rabbit pAb	ER63435	Human	WB,IHC-P,IF-Cell
	Fibulin 5 Mouse mAb [A1C12]	HA600044	Human	WB,IF-Cell,FC
	Fibulin 5 Rabbit pAb	ER1802-95	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	FILA Rabbit pAb	ER62428	Human	IHC-P
KD	Filamin A Recombinant Rabbit mAb [SA30-08]	ET1601-3	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Filamin A Rabbit pAb	R1601-3	Human, Mouse, Rat	WB
	Filamin B Recombinant Rabbit mAb [PSH07-34] - BSA and Azide free	HA751155	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Filamin B Recombinant Rabbit mAb [PSH07-34]	HA722822	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	FKBP52 Recombinant Rabbit mAb [JG96-30]	ET7108-93	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Flightless 1 Rabbit pAb	HA500167	Human, Mouse, Rat	WB,IHC-P
	FLNB Rabbit pAb	ER1908-91	Human	WB
	Flotillin 1 Recombinant Rabbit mAb [JB19-45]	ET7107-82	Human, Mouse, Rat	WB,IHC-P
	Formin 2 Rabbit pAb	ER63448	Human, Mouse	IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD, Relative expression	Furin Recombinant Rabbit mAb [JB35-53]	ET7107-37	Human, Mouse, Rat	WB,IHC-P
	GAL4 Rabbit pAb	ER1902-63	Human, Mouse, Rat	WB,IHC-P,FC
	Galectin 8 Recombinant Rabbit mAb [JB85-35]	ET7107-56	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
	gamma Catenin Recombinant Rabbit mAb [JF0973]	ET1702-33	Human, Mouse, Rat	WB,IF-Cell,IF- Tissue,IHC-P,FC
	gamma Sarcoglycan Recombinant Rabbit mAb [JA11-57]	ET1704-24	Human, Mouse	WB,IP,IHC-P
	gamma tubulin Mouse mAb [10-A4]	M1701-13	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	gamma Tubulin Recombinant Rabbit mAb [JF0972]	ET1702-32	Human, Mouse, Rat, Monkey, Chicken, Hamster, Cow, Dog, Fish, Xenopus tropicalis	WB,IF-Cell,IF- Tissue,IHC-P,IP,FC
	Gamma tubulin Rabbit pAb	R1411-2	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	GAS2 Recombinant Rabbit mAb [JE48-38]	HA721647	Human, Mouse, Rat	WB,IHC-P
	Gas8 Rabbit pAb	ER63504	Human, Mouse, Rat	IHC-P,WB
	GCDFP 15 Recombinant Rabbit mAb [JE48-53]	ET7109-80	Human	WB,IF-Cell,IHC-P
	GCNT3 Rabbit pAb	ER63517	Human, Mouse, Rat	IHC-P,IF-Cell
	GCP4 Rabbit pAb	ER63521	Human, Mouse	WB,IHC-P
	GCP5 Rabbit pAb	ER63522	Human, Mouse, Rat	WB,IHC-P
Relative expression	GDF15 Recombinant Rabbit mAb [PSH10-09] - BSA and Azide free	HA751339	Human	WB,IHC-P
	GDF15 Recombinant Rabbit mAb [JE39-90]	HA722551	Human	WB
	GDF15 Rabbit pAb	ER1909-46	Mouse, Rat	WB
Relative expression	GDF15 Recombinant Rabbit mAb [PSH10-09]	HA723178	Human	WB,IHC-P
	GEF H1 Recombinant Rabbit mAb [JG36-46]	ET7108-34	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	Gelsolin Recombinant Rabbit mAb [JB36-68]	ET7107-45	Human, Mouse, Monkey	WB,IF-Tissue,IHC-P
	GFAP Rabbit pAb	R1308-10	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Biotin Conjugated GFAP Mouse mAb [PSH05-17]	HA601355B	Human, Mouse, Rat	ELISA(Det)
	iFluor™ 647 Conjugated GFAP Mouse mAb [1-D4]	HA600103F	Human, Rat	IF-Tissue
	GFAP Recombinant Antibody - Rat IgG1 (Chimeric)	IRS148RT	Mouse	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	GFAP Recombinant Rabbit mAb	IRS069	Mouse	mIHC
	GFAP Mouse mAb [1-D4]	EM140707	Human, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
	iFluor™ 594 Conjugated GFAP Mouse mAb [1-D4]	HA600102F	Human, Rat	IF-Tissue
Relative expression	GFAP Recombinant Antibody [SA03-04] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610317	Mouse, Rat	IHC-Fr,WB
Relative expression	GFAP Recombinant Mouse mAb [PSH10-24] - BSA and Azide free	HA610225	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue
	iFluor™ 488 Conjugated GFAP Mouse mAb [1-D4]	HA600101F	Human, Rat	IF-Tissue,IHC-Fr
	GFAP Rabbit pAb	R1308-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	GFAP Recombinant Mouse mAb [PSH05-17] - BSA and Azide free (Detector)	HA601342	Human, Mouse, Rat	ELISA(Det)
Relative expression	GFAP Recombinant Antibody [SA03-04] - Rat IgG1 (Chimeric)	HA601483	Mouse, Rat	IHC-Fr,WB
	GFAP Recombinant Rabbit mAb [PSH05-16] - BSA and Azide free (Capture)	HA722234	Human, Mouse, Rat	ELISA(Cap)
	GFAP Recombinant Rabbit mAb [SA03-04]	ET1601-23	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr,mIHC,IP
Relative expression	GFAP Recombinant Mouse mAb [PSH10-24]	HA601380	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue
	Girdin Rabbit pAb	ER63536	Human, Mouse	IHC-P
	Glypican 4 Rabbit pAb	HA500432	Human, Mouse, Rat	WB,IF-Cell
	GPC1 Rabbit pAb	R1512-15	Human	WB,IF-Cell,IHC-P,FC
	GPC3 Mouse mAb [G3-C2]	EM1709-60	Human, Mouse, Rat	WB,IHC-P,FC
	GRAF Recombinant Rabbit mAb [JE62-23]	HA722409	Human, Mouse, Rat	WB,IHC-P
	Gremlin Rabbit pAb	ER1910-23	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	GWL Recombinant Rabbit mAb [PSH07-69]	HA722884	Human, Monkey	WB,IF-Cell,IP
	GWL Recombinant Rabbit mAb [PSH07-69] - BSA and Azide free	HA751181	Human, Monkey	WB,IF-Cell,IP
	HAPLN1 Recombinant Rabbit mAb [SY02-17] - BSA and Azide free	HA750130	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P
	HAPLN1 Recombinant Rabbit mAb [SY02-17]	ET1607-7	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P
	HAS2 Mouse mAb [C4-E7]	EM1709-26	Human	WB,IF-Cell,IHC-P
	HAS3 Mouse mAb [B2-C9]	EM1710-05	Human	IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	H Cadherin Recombinant Rabbit mAb [JB14-48]	ET7107-76	Human, Mouse, Rat	IHC-P, WB
Relative expression	heavy chain Myosin / MYH3 Recombinant Rabbit mAb [PSH18-12]	HA723975	Human, Mouse, Rat	WB, IHC-P, IP
	Hemopexin Recombinant Rabbit mAb [JU83-33]	ET1706-10	Human	WB, IF-Cell, IF-Tissue, IHC-P, IP
	HEPACAM Rabbit pAb	ER1910-45	Human	WB, IHC-P, IF-Tissue
	HGFA Rabbit pAb	ER65420	Human, Mouse, Rat	WB
	HPA1 Recombinant Rabbit mAb	ER1910-48	Human, Mouse, Rat	WB, IHC-P
	HPLN4 Rabbit pAb	ER62000	Human, Mouse	WB
	Hsp27 Recombinant Rabbit mAb [JJ09-13] - BSA and Azide free	HA750321	Human, Mouse, Rat	WB, IF-Cell, IF-Tissue, IHC-P, IP, FC
	Hsp27 Mouse mAb [10-8]	M1505-3	Human, Mouse, Rat	WB, IHC-P, IF-Cell
	Hsp27 Recombinant Rabbit mAb [JJ09-13]	ET1701-70	Human, Mouse, Rat	WB, IF-Cell, IF-Tissue, IHC-P, IP, FC
	Hsp27 Rabbit pAb	ER31216	Human	WB, IF-Cell, IHC-P
	Hsp90 Recombinant Rabbit mAb [JE05-06]	HA722689	Human, Mouse, Rat, Monkey	WB, IF-Cell, IHC-P, FC, IP
	Hsp90 alpha Recombinant Rabbit mAb [SY14-06]	ET1605-57	Human, Mouse, Rat, Monkey	WB, IF-Cell, IF-Tissue, IHC-P
	Hsp90 alpha Mouse mAb [A10-A5-D2]	M1603-3	Human, Mouse, Rat	WB, IF-Cell, IHC-P, FC
	Hsp90 alpha Recombinant Rabbit mAb [SY14-06] - BSA and Azide free	HA750091	Human, Mouse, Rat, Monkey	WB, IF-Cell, IF-Tissue, IHC-P
	Hsp90 alpha Rabbit pAb	R1510-29	Human, Mouse, Rat	WB, IF-Cell, IHC-P, FC
	HSPB7 Rabbit pAb	ER61587	Human, Mouse	WB
	Human Alpha-Synuclein enzyme-linked immunoassay kit	EH0191	Human	ELISA
	HRP Conjugated Human Angiogenin Recombinant Rabbit mAb [PSH12-65] - Detector	HA725078H	Human	ELISA(Det), ELISA
	Human Angiogenin Recombinant Rabbit mAb [PSH12-65] - BSA and Azide free (Detector)	HA725077	Human	ELISA(Det)
	Biotin Conjugated Human Angiogenin Recombinant Rabbit mAb [PSH12-65] - Detector	HA725078B	Human	ELISA(Det), ELISA
	Human Angiogenin Recombinant Rabbit mAb [PSH12-64] - BSA and Azide free (Capture)	HA725076	Human	ELISA(Cap)
	Human Angiogenin enzyme-linked immunoassay kit (one step)	EHY0190	Human	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human B7-H6/NCR3LG1 enzyme-linked immunoassay kit	EH0116	Human	ELISA
	Human BCMA/TNFRSF17 enzyme-linked immunoassay kit	EH0143	Human	ELISA
	Human beta-NGF enzyme-linked immunoassay kit	EH0065	Human	ELISA
	Human CCL13/MCP-4 enzyme-linked immunoassay kit	EH0199	Human	ELISA
	Human CCL17/TARC enzyme-linked immunoassay kit	EH0164	Human	ELISA
	Human CCL20/MIP-3 alpha enzyme-linked immunoassay kit	EH0028	Human	ELISA
	Human CCL21 enzyme-linked immunoassay kit	EH0184	Human	ELISA
	Human CCL23/MPIF-1 enzyme-linked immunoassay kit	EH0130	Human	ELISA
	Human CCL25 enzyme-linked immunoassay kit	EH0117	Human	ELISA
	Human CCL2 enzyme-linked immunoassay kit	EH0010	Human	ELISA
	Human CCL3/MIP-1 enzyme-linked immunoassay kit	EH0146	Human	ELISA
	Human CCL5/RANTES enzyme-linked immunoassay kit	EH0027	Human	ELISA
	Human CCL5/RANTES enzyme-linked immunoassay kit(one step)	EHY0027	Human	ELISA
	FITC Conjugated Human CD11b Recombinant Mouse mAb [A2-9-R]	HA600130F	Human	FC
	Biotin Conjugated Human CD147 Recombinant Rabbit mAb [PSH13-07] - Detector	HA725087B	Human	ELISA(Det),ELISA
	Human CD147 Recombinant Rabbit mAb [PSH13-07] - BSA and Azide free (Detector)	HA725086	Human	ELISA(Det)
	Human CD147 Recombinant Rabbit mAb [PSH13-06] - BSA and Azide free (Capture)	HA725085	Human	ELISA(Cap)
	Biotin Conjugated Human CD31 Recombinant Rabbit mAb [PSH07-83] - Detector	HA722908B	Human	ELISA(Det)
	Human CD31 Recombinant Rabbit mAb [PSH07-82] - BSA and Azide free (Capture)	HA722906	Human	ELISA(Cap)
	Human CD31 Recombinant Rabbit mAb [PSH07-83] - BSA and Azide free (Detector)	HA722907	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human CD63 Recombinant Rabbit mAb [PSH08-44] - Detector	HA722990B	Human	ELISA(Det),ELISA
	Human CD63 Recombinant Rabbit mAb [PSH08-44] - BSA and Azide free (Detector)	HA722989	Human	ELISA(Det)
	Human CD63 Recombinant Rabbit mAb [PSH08-43] - BSA and Azide free (Capture)	HA722988	Human	ELISA(Cap)
	Human cIAP2 enzyme-linked immunoassay kit	EH0123	Human	ELISA
	Human Complement C5 enzyme-linked immunoassay kit	EH0108	Human	ELISA
	Human Complement C5 enzyme-linked immunoassay kit (one step)	EHY0108	Human	ELISA
	Biotin Conjugated Human CRP Recombinant Rabbit mAb [PSH04-26] - Detector	HA722616B	Human	ELISA(Det)
	Human CRP Recombinant Rabbit mAb [PSH04-25] - BSA and Azide free (Capture)	HA722108	Human	ELISA(Cap)
	Human CRP Recombinant Rabbit mAb [PSH04-26] - BSA and Azide free (Detector)	HA722109	Human	ELISA(Det)
	Human CRP enzyme-linked immunoassay kit	EH0076	Human	ELISA
	Human CSPG4 Recombinant Rabbit mAb [PSH10-78] - BSA and Azide free (Capture)	HA723245	Human	ELISA(Cap)
	HRP Conjugated Human CSPG4 Recombinant Rabbit mAb [PSH10-79] - Detector	HA723376H	Human	ELISA(Det)
	Biotin Conjugated Human CSPG4 Recombinant Rabbit mAb [PSH10-79] - Detector	HA723247B	Human	ELISA(Det),ELISA
	Human CSPG4 Recombinant Rabbit mAb [PSH10-79] - BSA and Azide free (Detector)	HA723246	Human	ELISA(Det)
	Human CSPG4 Recombinant Rabbit mAb [PSH10-77] - BSA and Azide free (Capture)	HA723244	Human	ELISA(Cap)
	Human CTLA-4 enzyme-linked immunoassay kit	EH0179	Human	ELISA
	HRP Conjugated Human CXADR Recombinant Rabbit mAb [PSH17-47] - Detector	HA725294H	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human CXADR Recombinant Rabbit mAb [PSH17-47] - BSA and Azide free (Detector)	HA725294	Human	ELISA(Det)
	Biotin Conjugated Human CXADR Recombinant Rabbit mAb [PSH17-47] - Detector	HA725294B	Human	ELISA(Det),ELISA
	Human CXADR Recombinant Rabbit mAb [PSH17-46] - BSA and Azide free (Capture)	HA725293	Human	ELISA(Cap)
	Human CXCL10/IP-10 enzyme-linked immunoassay kit	EH0026	Human	ELISA
	Human CXCL13/BCA-1 enzyme-linked immunoassay kit	EH0172	Human	ELISA
	Human CXCL1/GRO alpha enzyme-linked immunoassay kit	EH0054	Human	ELISA
	Human CXCL6/GCP-2 enzyme-linked immunoassay kit(one step)	EH0165	Human	ELISA
	Human CXCL9/MIG enzyme-linked immunoassay kit	EH0015	Human	ELISA
	Biotin Conjugated Human E-Cadherin Recombinant Rabbit mAb [SY0287] - BSA and Azide free (Detector)	HA723057B	Human	ELISA(Det),ELISA
	Human E-Cadherin Recombinant Rabbit mAb [PSH09-02] - BSA and Azide free (Capture)	HA723055	Human	ELISA(Cap)
	Human E-Cadherin Recombinant Rabbit mAb [SY0287] - BSA and Azide free (Detector)	HA723056	Human	ELISA(Det)
	Human E-Cadherin enzyme-linked immunoassay kit	EH0115	Human	ELISA
	Human EGF enzyme-linked immunoassay kit	EH0023	Human	ELISA
	HRP Conjugated Human Endostatin Recombinant Rabbit mAb [PSH15-83] - Detector	HA723791H	Human	ELISA(Det),ELISA
	Human Endostatin Recombinant Rabbit mAb [PSH15-83] - BSA and Azide free (Detector)	HA723790	Human	ELISA(Det)
	Biotin Conjugated Human Endostatin Recombinant Rabbit mAb [PSH15-83] - Detector	HA723791B	Human	ELISA(Det),ELISA
	Human Endostatin Recombinant Rabbit mAb [PSH15-82] - BSA and Azide free (Capture)	HA723789	Human	ELISA(Cap)
	Human FGF-19 enzyme-linked immunoassay kit	EH0140	Human	ELISA
	Human FGF basic /FGF2/bFGF enzyme-linked immunoassay kit	EH0024	Human	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human Fibromodulin Recombinant Rabbit mAb [PSH18-73] - Detector	HA724016B	Human	ELISA(Det),ELISA
	Human Fibromodulin Recombinant Rabbit mAb [PSH18-72] - BSA and Azide free (Capture)	HA724015	Human	ELISA(Cap)
	Human Fibromodulin Recombinant Rabbit mAb [PSH18-73] - BSA and Azide free (Detector)	HA724016	Human	ELISA(Det)
	Human Fibromodulin Recombinant Rabbit mAb [PSH18-74] - BSA and Azide free (Capture)	HA724017	Human	ELISA(Cap)
	Human FOLR1 enzyme-linked immunoassay kit	EH0137	Human	ELISA
	Human G-CSF enzyme-linked immunoassay kit	EH0149	Human	ELISA
	Biotin Conjugated Human GDF-15 Recombinant Rabbit mAb [PSH16-01] - Detector	HA723806B	Human	ELISA(Det),ELISA
	Human GDF-15 Recombinant Rabbit mAb [PSH16-02] - BSA and Azide free (Detector)	HA723807	Human	ELISA(Det)
	Human GDF-15 Recombinant Rabbit mAb [PSH16-01] - BSA and Azide free (Detector)	HA723806	Human	ELISA(Det)
	Human GDF-15 Recombinant Rabbit mAb [PSH16-00] - BSA and Azide free (Capture)	HA723805	Human	ELISA(Cap)
	Biotin Conjugated Human GDF-15 Recombinant Rabbit mAb [PSH16-02] - Detector	HA723807B	Human	ELISA(Det),ELISA
	Human GITR/TNFRSF18 enzyme-linked immunoassay kit	EH0095	Human	ELISA
	Human GM-CSF enzyme-linked immunoassay kit	EH0051	Human	ELISA
	Human Granzyme A/GZMA enzyme-linked immunoassay kit	EH0129	Human	ELISA
	Human Granzyme B enzyme-linked immunoassay kit	EH0052	Human	ELISA
	Human Growth Hormone enzyme-linked immunoassay kit	EH0142	Human	ELISA
	Human HGF enzyme-linked immunoassay kit	EH0019	Human	ELISA
	Human HVEM/TNFRSF14 enzyme-linked immunoassay kit	EH0107	Human	ELISA
	Human IDO1 enzyme-linked immunoassay kit	EH0118	Human	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IFN-gamma enzyme-linked immunoassay kit	EH0039	Human	ELISA
	Human IFN-gamma R1 enzyme-linked immunoassay kit	EH0104	Human	ELISA
	Human IFN- α enzyme-linked immunoassay kit	EH0085	Human	ELISA
	Human IFN- β enzyme-linked immunoassay kit	EH0017	Human	ELISA
	Biotin Conjugated Human IL-10 Recombinant Rabbit mAb [PS00-81]	HA722652B	Human	ELISA(Det)
	Human IL-10 Recombinant Rabbit mAb [PS00-98] - BSA and Azide free (Capture)	HA721236	Human	ELISA(Cap)
	Human IL-10 Recombinant Rabbit mAb [PS00-81] - BSA and Azide free (Detector)	HA721237	Human	ELISA(Det)
	Human IL-10 enzyme-linked immunoassay kit	EH0009	Human	ELISA
	Human IL-12B/IL-12 p40 enzyme-linked immunoassay kit	EH0098	Human	ELISA
	Human IL-12p70 enzyme-linked immunoassay kit	EH0021	Human	ELISA
	Human IL-13 enzyme-linked immunoassay kit	EH0049	Human	ELISA
	Human IL-13R alpha 2 enzyme-linked immunoassay kit	EH0182	Human	ELISA
	Biotin Conjugated Human IL-15 Recombinant Rabbit mAb [PSH04-20]	HA722612B	Human	ELISA(Det)
	Human IL-15 Recombinant Rabbit mAb [PSH04-20] - BSA and Azide free (Detector)	HA722103	Human	ELISA(Det)
	Human IL-15 Recombinant Rabbit mAb [PSH04-19] - BSA and Azide free (Capture)	HA722102	Human	ELISA(Cap)
	Human IL-15 enzyme-linked immunoassay kit	EH0031	Human	ELISA
	Human IL-17A enzyme-linked immunoassay kit	EH0004	Human	ELISA
	Human IL-17E/IL-25 enzyme-linked immunoassay kit	EH0066	Human	ELISA
	Human IL-18BP enzyme-linked immunoassay kit	EH0131	Human	ELISA
	Human IL-18 enzyme-linked immunoassay kit	EH0053	Human	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free (Capture)	HA721278	Human	ELISA(Cap)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA721713	Human	WB
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39]	HA721712	Human	WB
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free	HA750745	Human	WB
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - Detector	HA723882B	Human	ELISA(Det),ELISA
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-87] - BSA and Azide free (Capture)	HA723883	Human	ELISA(Cap)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-85] - BSA and Azide free (Capture)	HA723881	Human	ELISA(Cap)
	Human IL-1 beta Recombinant Rabbit mAb [PSH16-86] - BSA and Azide free (Detector)	HA723882	Human	ELISA(Det)
	Human IL-1 beta Recombinant Rabbit mAb [PS00-82] - BSA and Azide free (Detector)	HA721279	Human	ELISA(Det)
Cell treatment	Human IL-1 beta Recombinant Rabbit mAb [PS01-39] - BSA and Azide free	HA750744	Human	WB
	Biotin Conjugated Human IL-1 beta Recombinant Rabbit mAb [PS00-82]	HA722649B	Human	ELISA(Det)
	Human IL-1RA/IL-1F3 enzyme-linked immunoassay kit	EH0177	Human	ELISA
	Human IL-1 β enzyme-linked immunoassay kit	EH0006	Human	ELISA
	Human IL-2 Recombinant Rabbit mAb [PSH01-42] - BSA and Azide free (Capture)	HA721680	Human	ELISA(Cap)
	Human IL-2 Mouse mAb [A8B7] - BSA and Azide free (Detector)	HA601006	Human	WB,ELISA(Det)
	Biotin Conjugated Human IL-2 Mouse mAb [A8B7]	HA601361B	Human	ELISA(Det)
	Human IL-21 enzyme-linked immunoassay kit	EH0018	Human	ELISA
	Human IL-22 enzyme-linked immunoassay kit	EH0086	Human	ELISA
	Human IL-23 p19/IL-12 p40 enzyme-linked immunoassay kit	EH0109	Human	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-27 enzyme-linked immunoassay kit	EH0094	Human	ELISA
	Human IL-2 enzyme-linked immunoassay kit	EH0013	Human	ELISA
	Human IL-2 Receptor alpha enzyme-linked immunoassay kit	EH0170	Human	ELISA
	Human IL-31 enzyme-linked immunoassay kit	EH0097	Human	ELISA
	Human IL-33 enzyme-linked immunoassay kit	EH0047	Human	ELISA
	Human IL-3 enzyme-linked immunoassay kit	EH0050	Human	ELISA
	Human IL-4 enzyme-linked immunoassay kit	EH0007	Human	ELISA
	Human IL-4R enzyme-linked immunoassay kit	EH0145	Human	ELISA
	Human IL-5 enzyme-linked immunoassay kit	EH0008	Human	ELISA
	Human IL-6 Recombinant Rabbit mAb [PS00-21] - BSA and Azide free (Capture)	HA721167	Human	ELISA(Cap)
	Human IL-6 Recombinant Rabbit mAb [PS00-15] - BSA and Azide free (Detector)	HA721168	Human	ELISA(Det)
	HRP Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA721293	Human	ELISA(Det)
	Human IL-6 Mouse mAb [A6D1]	HA601051	Human	WB,ELISA
	Human IL-6 Mouse mAb [A6D2]	HA601050	Human	WB,ELISA
	Biotin Conjugated Human IL-6 Recombinant Rabbit mAb [PS00-15] - Detector	HA722233B	Human	ELISA(Det)
	Human IL-6 enzyme-linked immunoassay kit	EH0001	Human	ELISA
	Human IL-7 enzyme-linked immunoassay kit	EH0045	Human	ELISA
	Human IL-8 Recombinant Rabbit mAb [PS01-22] - BSA and Azide free (Detector)	HA722099	Human	ELISA(Det)
	Human IL-8 Recombinant Rabbit mAb [PS00-83] - BSA and Azide free (Capture)	HA722098	Human	ELISA(Cap)
	Biotin Conjugated Human IL-8 Recombinant Rabbit mAb [PS01-22]	HA722611B	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-8/CXCL8 enzyme-linked immunoassay kit	EH0014	Human	ELISA
	Biotin Conjugated Human LAMC2 Recombinant Rabbit mAb [PSH15-15] - Detector	HA723719B	Human	ELISA(Det),ELISA
	Human LAMC2 Recombinant Rabbit mAb [PSH15-15] - BSA and Azide free (Detector)	HA723718	Human	ELISA(Det)
	Human LAMC2 Recombinant Rabbit mAb [PSH15-14] - BSA and Azide free (Capture)	HA723717	Human	ELISA(Cap)
	Human Leptin enzyme-linked immunoassay kit	EH0056	Human	ELISA
	Human LIF enzyme-linked immunoassay kit	EH0096	Human	ELISA
	Human Lipoprotein A/LPA enzyme-linked immunoassay kit	EH0171	Human	ELISA
	Human MBL2 enzyme-linked immunoassay kit	EH0090	Human	ELISA
	Human M-CSF enzyme-linked immunoassay kit	EH0126	Human	ELISA
	Human MICA enzyme-linked immunoassay kit	EH0044	Human	ELISA
	Human MICB enzyme-linked immunoassay kit	EH0127	Human	ELISA
	Human MMP-1 Recombinant Rabbit mAb [PSH07-59] - BSA and Azide free (Capture)	HA722874	Human	ELISA(Cap)
	Biotin Conjugated Human MMP-1 Recombinant Rabbit mAb [PSH07-60]	HA722876B	Human	ELISA(Det)
	Human MMP-1 Recombinant Rabbit mAb [PSH07-60] - BSA and Azide free (Detector)	HA722875	Human	ELISA(Det)
	Biotin Conjugated Human MMP13 Recombinant Rabbit mAb [PSH11-14] - Detector	HA723295B	Human	ELISA(Det),ELISA
	Human MMP13 Recombinant Rabbit mAb [PSH11-13] - BSA and Azide free (Capture)	HA723293	Human	ELISA(Cap)
	Human MMP13 Recombinant Rabbit mAb [PSH11-14] - BSA and Azide free (Detector)	HA723294	Human	ELISA(Det)
	Human MMP-13 enzyme-linked immunoassay kit	EH0135	Human	ELISA
	Human MMP-3 Recombinant Rabbit mAb [PSH07-61] - BSA and Azide free (Capture)	HA722877	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human MMP-3 Recombinant Rabbit mAb [PSH07-62] - BSA and Azide free (Detector)	HA722878	Human	ELISA(Det)
	Biotin Conjugated Human MMP-3 Recombinant Rabbit mAb [PSH07-62]	HA722879B	Human	ELISA(Det)
	Human MMP-3 enzyme-linked immunoassay kit	EH0169	Human	ELISA
	Human MMP-7 Recombinant Rabbit mAb [PSH07-94] - BSA and Azide free (Detector)	HA722927	Human	ELISA(Det)
	Biotin Conjugated Human MMP-7 Recombinant Rabbit mAb [PSH07-94] - Detector	HA722928B	Human	ELISA(Det)
	Human MMP-7 Recombinant Rabbit mAb [PSH07-93] - BSA and Azide free (Capture)	HA722926	Human	ELISA(Cap)
	Human MMP-7 enzyme-linked immunoassay kit	EH0125	Human	ELISA
	Human MMP-8 Recombinant Rabbit mAb [PSH03-34] - BSA and Azide free (Detector)	HA721979	Human	ELISA(Det)
	Biotin Conjugated Human MMP-8 Recombinant Rabbit mAb [PSH03-34]	HA722624B	Human	ELISA(Det)
	Human MMP-8 Recombinant Rabbit mAb [PSH03-33] - BSA and Azide free (Capture)	HA721978	Human	ELISA(Cap)
	Human MMP-8 enzyme-linked immunoassay kit	EH0103	Human	ELISA
	Biotin Conjugated Human MMP-9 Recombinant Rabbit mAb [PSH04-24]	HA722615B	Human	ELISA(Det)
	Human MMP-9 Recombinant Rabbit mAb [PSH04-23] - BSA and Azide free (Capture)	HA722106	Human	ELISA(Cap)
	Human MMP-9 Recombinant Rabbit mAb [PSH04-62] - BSA and Azide free (Capture)	HA722162	Human	ELISA(Cap)
	Human MMP-9 Recombinant Rabbit mAb [PSH04-63] - BSA and Azide free (Detector)	HA722163	Human	ELISA(Det)
	Human MMP-9 Recombinant Rabbit mAb [PSH04-24] - BSA and Azide free (Detector)	HA722107	Human	ELISA(Det)
	Human MMP-9 enzyme-linked immunoassay kit	EH0020	Human	ELISA
	Human/Mouse/Rat Irisin enzyme-linked immunoassay kit	EH0081	Human, Mouse, Rat	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-42] - BSA and Azide free (Detector)	HA723331	Human, Mouse, Rat	ELISA(Det)
	Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-43] - BSA and Azide free (Capture)	HA723333	Human, Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-42] - Detector	HA723332B	Human, Mouse, Rat	ELISA(Det),ELISA
	Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-44] - BSA and Azide free (Detector)	HA723334	Human, Mouse, Rat	ELISA(Det)
	Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-41] - BSA and Azide free (Capture)	HA723330	Human, Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Human/Mouse/Rat NF-L Recombinant Rabbit mAb [PSH11-44] - Detector	HA723335B	Human, Mouse, Rat	ELISA(Det),ELISA
	Human/Mouse/Rat/Porcine/Canine TGF- β 1 enzyme-linked immunoassay kit	EH0012	Human, Mouse, Rat, Porcine, Canine	ELISA
	Biotin Conjugated Human NCAM1 Recombinant Rabbit mAb [PSH08-08] - Detector	HA722948B	Human	ELISA(Det)
	Human NCAM1 Recombinant Rabbit mAb [PSH08-09] - BSA and Azide free (Detector)	HA722949	Human	ELISA(Det)
	Biotin Conjugated Human NCAM1 Recombinant Rabbit mAb [PSH08-09] - Detector	HA722950B	Human	ELISA(Det)
	Human NCAM1 Recombinant Rabbit mAb [PSH08-08] - BSA and Azide free (Capture/Detector)	HA722947	Human	ELISA(Det),ELISA(Cap)
	Human NCAM1 Recombinant Rabbit mAb [PSH08-07] - BSA and Azide free (Capture)	HA722946	Human	ELISA(Cap)
	APC Conjugated Human NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53]	HA723757F	Human	FC
	PE Conjugated Human NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53]	HA720215F	Human	FC
	Biotin Conjugated Human Osteopontin Recombinant Rabbit mAb [PSH08-60] - Detector	HA723008B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Osteopontin Recombinant Rabbit mAb [PSH08-59] - BSA and Azide free (Detector)	HA723005	Human	ELISA(Det)
	Human Osteopontin Recombinant Rabbit mAb [PSH08-60] - BSA and Azide free (Detector)	HA723006	Human	ELISA(Det)
	Biotin Conjugated Human Osteopontin Recombinant Rabbit mAb [PSH08-59] - Detector	HA723007B	Human	ELISA(Det),ELISA
	Human Osteopontin Recombinant Rabbit mAb [PSH08-58] - BSA and Azide free (Capture)	HA723004	Human	ELISA(Cap)
	Human Osteopontin (OPN) enzyme-linked immunoassay kit	EH0057	Human	ELISA
	Human Osteoprotegerin Recombinant Rabbit mAb [PSH13-32] - BSA and Azide free (Detector)	HA725068	Human	ELISA(Det)
	Human Osteoprotegerin Recombinant Rabbit mAb [PSH13-31] - BSA and Azide free (Capture)	HA725067	Human	ELISA(Cap)
	Biotin Conjugated Human Osteoprotegerin Recombinant Rabbit mAb [PSH13-32] - Detector	HA725069B	Human	ELISA(Det),ELISA
	Human Osteoprotegerin/TNFRSF11B enzyme-linked immunoassay kit	EH0070	Human	ELISA
	Human OX40/TNFRSF4 enzyme-linked immunoassay kit	EH0120	Human	ELISA
	Human PD-1 enzyme-linked immunoassay kit	EH0062	Human	ELISA
	Human PDGF-AA enzyme-linked immunoassay kit	EH0058	Human	ELISA
	Human PDGF-BB enzyme-linked immunoassay kit	EH0068	Human	ELISA
	Human PD-L1/B7-H1 enzyme-linked immunoassay kit	EH0048	Human	ELISA
	Human PD-L2 Recombinant Rabbit mAb [PSH09-81] - BSA and Azide free (Capture)	HA723146	Human	ELISA(Cap)
	HRP Conjugated Human PD-L2 Recombinant Rabbit mAb [PSH09-82] - Detector	HA723373H	Human	ELISA(Det)
	Human PD-L2 Recombinant Rabbit mAb [PSH09-82] - BSA and Azide free (Detector)	HA723147	Human	ELISA(Det)
	Human PD-L2 Recombinant Rabbit mAb [PSH09-80] - BSA and Azide free (Capture)	HA723145	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human PD-L2 Recombinant Rabbit mAb [PSH09-82] - Detector	HA723148B	Human	ELISA(Det),ELISA
	Human Pentraxin 3/PTX3 enzyme-linked immunoassay kit	EH0141	Human	ELISA
	Human Prolactin enzyme-linked immunoassay kit	EH0089	Human	ELISA
	Human RANKL Recombinant Rabbit mAb [PSH16-58] - BSA and Azide free (Capture)	HA723860	Human	ELISA(Cap)
	Human RANKL Recombinant Rabbit mAb [PSH16-60] - BSA and Azide free (Capture)	HA723862	Human	ELISA(Cap)
	Human RANKL Recombinant Rabbit mAb [PSH16-59] - BSA and Azide free (Detector)	HA723861	Human	ELISA(Det)
	Biotin Conjugated Human RANKL Recombinant Rabbit mAb [PSH16-59] - Detector	HA723861B	Human	ELISA(Det),ELISA
	Human Siglec 9 enzyme-linked immunoassay kit	EH0161	Human	ELISA
	Human SIRP alpha Recombinant Rabbit mAb [PSH17-63] - BSA and Azide free (Capture)	HA725347	Human	ELISA(Cap)
	Biotin Conjugated Human SIRP alpha Recombinant Rabbit mAb [PSH17-64] - Detector	HA725348B	Human	ELISA(Det),ELISA
	Human SIRP alpha Recombinant Rabbit mAb [PSH17-64] - BSA and Azide free (Detector)	HA725348	Human	ELISA(Det)
	Human TGF- α enzyme-linked immunoassay kit	EH0059	Human	ELISA
	Human TIMP1 Recombinant Rabbit mAb [PSH18-79] - BSA and Azide free (Detector)	HA725368	Human	ELISA(Det)
	Human TIMP1 Recombinant Rabbit mAb [PSH18-78] - BSA and Azide free (Capture)	HA725367	Human	ELISA(Cap)
	Biotin Conjugated Human TIMP1 Recombinant Rabbit mAb [PSH18-79] - Detector	HA725368B	Human	ELISA(Det),ELISA
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free (Capture)	HA721265	Human	ELISA(Cap)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40]	HA721286	Human	IHC-P, mIHC
	Biotin Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721559	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free	HA750593	Human	IHC-P
	HRP Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721294	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb	IRS040	Human	mIHC
	Human TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free (Detector)	HA721266	Human	ELISA(Det)
	Human TNF beta/Lymphotoxin alpha enzyme-linked immunoassay kit	EH0181	Human	ELISA
	Human TNF RII/TNFRSF1B enzyme-linked immunoassay kit	EH0148	Human	ELISA
	Human TNF RI/TNFRSF1A enzyme-linked immunoassay kit	EH0128	Human	ELISA
	Human TNF RI/TNFRSF1A enzyme-linked immunoassay kit (one step)	EHY0128	Human	ELISA
	Human TNF- α enzyme-linked immunoassay kit	EH0002	Human	ELISA
	Human Total HSP60 enzyme-linked immunoassay kit	EH0060	Human	ELISA
	Human TRAIL/TNFSF10 enzyme-linked immunoassay kit	EH0176	Human	ELISA
	Human TREM2 enzyme-linked immunoassay kit (one step)	EHY0122	Human	ELISA
	Human TSLP enzyme-linked immunoassay kit	EH0061	Human	ELISA
	Human uPAR Mouse mAb [PSH13-76] - BSA and Azide free (Capture)	HA601437	Human	ELISA(Cap)
	Biotin Conjugated Human uPAR Mouse mAb [PSH13-77] - Detector	HA601439B	Human	ELISA(Det),ELISA
	Human uPAR Mouse mAb [PSH13-77] - BSA and Azide free (Detector)	HA601438	Human	ELISA(Det)
	Human uPAR enzyme-linked immunoassay kit	EH0077	Human	ELISA
	HRP Conjugated Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - Detector	HA725066H	Human	ELISA(Det),ELISA
	Human VCAM-1 Recombinant Rabbit mAb [PSH12-46] - BSA and Azide free (Capture)	HA725064	Human	ELISA(Cap)
	Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - BSA and Azide free (Detector)	HA725065	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human VCAM-1 Recombinant Rabbit mAb [PSH12-47] - Detector	HA725066B	Human	ELISA(Det),ELISA
	Human VE Cadherin Recombinant Rabbit mAb [PSH02-61] - BSA and Azide free (Capture)	HA721838	Human	ELISA(Cap)
	Human VE Cadherin Recombinant Rabbit mAb [PSH02-62] - BSA and Azide free (Detector)	HA721839	Human	ELISA(Det)
	Biotin Conjugated Human VE Cadherin Recombinant Rabbit mAb [PSH02-62] - Detector	HA722635B	Human	ELISA(Det)
	Human VE Cadherin enzyme-linked immunoassay kit	EH0091	Human	ELISA
	Human VEGF enzyme-linked immunoassay kit	EH0075	Human	ELISA
	Human VEGFR2/KDR enzyme-linked immunoassay kit	EH0200	Human	ELISA
	Human von Willebrand Factor enzyme-linked immunoassay kit (one step)	EHY0138	Human	ELISA
	HRP Conjugated Human YKL-40 Recombinant Rabbit mAb [PSH13-09] - Detector	HA725090H	Human	ELISA(Det),ELISA
	Human YKL-40 Recombinant Rabbit mAb [PSH13-09] - BSA and Azide free (Detector)	HA725089	Human	ELISA(Det)
	Biotin Conjugated Human YKL-40 Recombinant Rabbit mAb [PSH13-09] - Detector	HA725090B	Human	ELISA(Det),ELISA
	Human YKL-40 Recombinant Rabbit mAb [PSH13-08] - BSA and Azide free (Capture)	HA725088	Human	ELISA(Cap)
	Hyaluronidase2 Rabbit pAb	ER1910-95	Human, Mouse	WB,IHC-P,IF-Tissue
	Hyaluronidase3 Rabbit pAb	ER1910-96	Human, Mouse, Rat	WB,FC
	ICAM4 Rabbit pAb	ER1911-01	Human, Mouse	WB,FC
	IGF1R Rabbit pAb	R1310-6	Human	WB,IHC-P
	IGF1R Mouse mAb [D5-G2]	HA600033	Human	WB,FC
Relative expression	IGF1 Receptor beta Recombinant Rabbit mAb [PSH11-72]	HA723368	Human, Mouse, Rat	WB,IF-Cell
Relative expression	IGF1 Receptor beta Recombinant Rabbit mAb [PSH11-72] - BSA and Azide free	HA751427	Human, Mouse, Rat	WB,IF-Cell
	IGFBP2 Recombinant Rabbit mAb [JU37-46]	ET7106-63	Human	WB,IHC-P
	IGFBP2 Rabbit pAb	ER1911-11	Human	WB,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	IGFBP4 Rabbit pAb	ER1911-13	Mouse, Rat	WB,IHC-P,IF-Tissue
	IGF-IR Rabbit pAb	ER63734	Human, Mouse, Rat	IF-Cell,WB,IHC-P
Cell treatment	IL-10 Recombinant Rabbit mAb [JE59-71]	HA722032	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	IL-10 Rabbit pAb	ER65464	Human, Mouse, Rat	WB,IHC-P
	IL-10 Rabbit pAb	ER1911-19	Human, Mouse	WB,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Mouse mAb [A7F9]	HA601002	Human, Mouse, Rat	WB,ELISA,IF-Cell
	IL-1 beta Recombinant Rabbit mAb [PSH18-00]	HA723965	Mouse, Rat	WB
Cell treatment	IL-1 beta Recombinant Rabbit mAb [JJ087-3]	ET1701-39	Human	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R]	HA601036	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Recombinant Mouse mAb [A7F7-R] - BSA and Azide free	HA610004	Human, Mouse, Rat	WB,ELISA,IHC-P,IF-Tissue
Cell treatment	IL-1 beta Mouse mAb [A7A10]	HA601035	Human	WB,ELISA,IF-Cell
	IL-2 Mouse mAb [A8B8]	HA601015	Human	WB,ELISA
	IL-2 Mouse mAb [A8B9]	HA601016	Human	WB,ELISA
	IL-2 Rabbit pAb	ER63748	Human, Rat	WB,IHC-P
	IL-6 Mouse mAb [1-6]	EM1701-45	Mouse, Rat	WB,IF-Cell
	IL-6 Rabbit pAb	R1203-2	Mouse	WB,ELISA
	IL-6 Mouse mAb [3-A4]	M1303-3	Mouse	WB,IF-Cell,FC,ELISA
	IL-6 Recombinant Rabbit mAb	IRS160RB	Human	mIHC
	IL-6 Rabbit pAb	R1412-2	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	IL-6 Mouse mAb	EM30301	Human	WB
Cell treatment	IL-8 Rabbit pAb	ER1901-61	Human	WB,IHC-P,IF-Cell,FC
	IL-8 Rabbit pAb	ER1706-67	Human	ELISA,IF-Cell,IHC-P
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-96]	HA723263	Human	WB,FC
Cell treatment	IL-8 Rabbit pAb	R1511-15	Human	WB,IHC-P,IF-Cell,FC
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-97]	HA723264	Human	WB
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-97] - BSA and Azide free	HA751387	Human	WB
Cell treatment	IL-8 Recombinant Rabbit mAb [PSH10-96] - BSA and Azide free	HA751386	Human	WB,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Integrin alpha 1 Recombinant Rabbit mAb [PSH14-04]	HA723594	Human	FC
	Integrin alpha 1 Rabbit pAb	ER1911-50	Human, Rat	WB,FC
	Integrin alpha 1 Recombinant Rabbit mAb [PSH14-05]	HA723595	Human, Mouse, Rat	WB,IHC-P
	Integrin alpha 1 Recombinant Rabbit mAb [PSH14-05] - BSA and Azide free	HA751500	Human, Mouse, Rat	WB,IHC-P
KD	Integrin alpha 2 Recombinant Rabbit mAb [SN0752]	ET1611-57	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
KD	Integrin alpha 2 Recombinant Rabbit mAb [SN0752] - BSA and Azide free	HA750263	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	Integrin alpha 3 Recombinant Mouse mAb [A7H2-R]	HA601273	Human	WB,IF-Cell
	Integrin alpha 3 Recombinant Mouse mAb [A7H2-R] - BSA and Azide free	HA610148	Human	WB,IF-Cell
	Integrin alpha 3 Rabbit pAb	HA500111	Human, Rat	WB,IHC-P
	Integrin alpha 3 Mouse mAb [A7H2]	HA600100	Human	WB,IF-Cell
	Integrin alpha 4 Rabbit pAb	ER65508	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Integrin alpha 4 Recombinant Rabbit mAb [JA09-36]	ET1704-31	Human, Mouse, Rat	WB,IHC-P,FC
KD	Integrin alpha 5 Recombinant Rabbit mAb [JJ08-94]	ET1701-58	Human, Mouse, Monkey, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Integrin alpha 6 Recombinant Rabbit mAb [SR45-00] - BSA and Azide free	HA750028	Human, Mouse, Rat	WB,IHC-P
	Integrin alpha 6 Recombinant Rabbit mAb [SR45-00]	ET1602-1	Human, Mouse, Rat	WB,IHC-P
	Integrin alpha 9 Recombinant Rabbit mAb [JE41-58]	HA722406	Human, Mouse, Rat	WB
	Integrin alpha D Rabbit pAb	ER65586	Human	WB,IHC-P
	Integrin alpha E Recombinant Rabbit mAb [SN62-06]	ET1611-27	Human	IHC-P
KD	Integrin alpha V Recombinant Rabbit mAb [SC56-07]	ET1610-15	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
KD	Integrin alpha V Recombinant Rabbit mAb [SC56-07] - BSA and Azide free	HA750201	Human, Mouse, Rat	WB,IF-Cell,FC,IHC-P
	Integrin Alpha V Rabbit pAb	ER1911-52	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,IP
KD	Integrin beta 1 Recombinant Rabbit mAb [SA40-08] - BSA and Azide free	HA750008	Human, Mouse, Rat, Zebrafish	WB,IHC-P
KD	Integrin beta 1 Recombinant Rabbit mAb [SA40-08]	ET1601-17	Human, Mouse, Rat, Zebrafish	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Integrin beta 1 Mouse mAb [B10-A5]	EM1901-77	Human, Mouse	WB,IHC-P,FC
KD	Integrin beta 1 Rabbit pAb	ER31001	Human, Mouse	WB,IHC-P,FC
KD	Integrin beta 1 Recombinant Rabbit mAb [SR30-03]	ET1602-10	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Integrin beta 3 Recombinant Rabbit mAb [JE22-64]	ET7110-61	Human, Mouse	WB,IHC-P
Relative expression	Integrin beta 3 Recombinant Rabbit mAb [PSH09-27]	HA723085	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
KD	Integrin beta 3 Recombinant Rabbit mAb [SJ19-09]	ET1606-49	Human, Mouse	WB,IHC-P,IF-Tissue
	Integrin beta 3 Rabbit pAb	ER1601-4	Human, Mouse, Rat	WB,IHC-P
	Integrin beta 4 Recombinant Rabbit mAb [JM11-06]	ET1703-52	Human, Mouse, Rat	WB,IHC-P
	Integrin beta 5 Recombinant Rabbit mAb [PSH11-53] - BSA and Azide free	HA751415	Human	IF-Cell,FC
	Integrin beta 5 Recombinant Rabbit mAb [PSH10-13] - BSA and Azide free	HA751342	Human, Rat, Monkey	WB,IHC-P
	Integrin beta 5 Rabbit pAb	ER63768	Human, Mouse, Rat, Monkey	WB
	Integrin beta 5 Recombinant Rabbit mAb [PSH11-53]	HA723345	Human	IF-Cell,FC
	Integrin beta 5 Recombinant Rabbit mAb [PSH10-13]	HA723182	Human, Rat, Monkey	WB,IHC-P
	Integrin beta 7 Rabbit pAb	ER65852	Human, Mouse	WB
	Inversin Rabbit pAb	HA500442	Human, Rat, Mouse	WB
	Involucrin Recombinant Rabbit mAb [JU30-27]	ET1706-32	Human, Rat	WB,IHC-P,FC
	IRTKS Rabbit pAb	ER63780	Human, Mouse, Rat	WB,IHC-P
	ITBP2 Rabbit pAb	ER60304	Human, Mouse	WB
	ITGA7 Rabbit pAb	ER1911-61	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	JAM-A Recombinant Rabbit mAb [SC60-07]	ET1610-90	Human	WB,IHC-P
	JAM-B Rabbit pAb	ER65451	Human, Mouse, Rat	WB,IHC-P
	JAMC Rabbit pAb	ER1911-63	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	JAML1 Rabbit pAb	ER61607	Human	WB
	K1C23 Rabbit pAb	ER60307	Human	WB
	Kallikrein 11 Rabbit pAb	HA500216	Human, Mouse	WB
	Kallikrein 5 Recombinant Rabbit mAb [JE54-94]	ET7110-90	Human	WB,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Kallikrein 8 Recombinant Rabbit mAb [JE62-46]	HA720101	Human	WB,IHC-P
	K-cadherin Rabbit pAb	ER65448	Human, Mouse, Rat	WB,IHC-P
	K Cadherin Recombinant Rabbit mAb [JE40-27]	HA722277	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	Keratin 12 / K12 Recombinant Rabbit mAb [JE55-35]	ET7111-24	Human, Mouse, Rat	WB,IHC-P
	Keratin 77 Rabbit pAb	ER65895	Human	IHC-P,WB
	Keratocan Rabbit pAb	ER1911-83	Mouse	WB
	KID Mouse mAb [D1-E3]	EM1709-64	Human	WB,IHC-P,FC
	KIF13B Rabbit pAb	ER1911-88	Mouse, Rat	WB,IHC-P,IF-Tissue
	KIF1C Rabbit pAb	ER60141	Human, Mouse, Rat	WB
	KIF20A Rabbit pAb	ER63804	Human	WB,IHC-P,IF-Cell
	KIF3A Recombinant Rabbit mAb [JE31-39]	HA721478	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	KIF3A Rabbit pAb	ER1803-40	Human, Mouse, Rat	WB,IHC-P,FC
	KIF4A Rabbit pAb	ER63806	Human	WB
KD	KIF4A Recombinant Rabbit mAb [JE38-91]	HA721591	Human, Mouse	WB,IHC-P
	KIF5B Recombinant Rabbit mAb [JE64-16]	HA721045	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Kif 7 Mouse mAb [3F8]	EM60063	Human, Mouse, Rat	IHC-P,IF-Cell
	KIFC1 Recombinant Rabbit mAb [JE00-46]	HA722197	Human	WB,IF-Cell,IHC-P,FC
	Kindlin 2 Mouse mAb [A10G12]	HA601343	Human, Mouse, Rat	WB,IF-Cell,IP
	Kininogen 1 Recombinant Rabbit mAb [JE43-72]	ET7109-34	Human, Mouse	WB,IHC-P
	KLC1 Recombinant Rabbit mAb [JE54-47]	ET7110-72	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	KLC2 Rabbit pAb	ER1911-96	Human, Mouse	WB
	KLC3 Recombinant Rabbit mAb [JE77-61]	HA722375	Human, Mouse, Rat	WB,IF-Cell,FC
	KLK1 Rabbit pAb	ER65582	Human	WB
	KLKB1 Recombinant Rabbit mAb [JE31-29]	HA721397	Human	WB,IHC-P
	KPTN Rabbit pAb	ER2001-17	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	KRT34 Rabbit pAb	ER60310	Human	WB
	KRT81 Rabbit pAb	ER1912-12	Mouse	WB
	KT33A Rabbit pAb	ER61609	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	KT33B Rabbit pAb	ER60313	Human	WB
Relative expression	L1CAM Recombinant Rabbit mAb [JM11-05]	ET1703-51	Human	WB,IHC-P,IF-Tissue
Relative expression	L1CAM Rabbit pAb	0805-5	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	LAMB3 Rabbit pAb	HA500480	Human, Mouse, Rat	WB,IHC-P
Relative expression	LAMC2 Recombinant Rabbit mAb [PSH07-27] - BSA and Azide free	HA751148	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,IP
Relative expression	LAMC2 Recombinant Rabbit mAb [PSH07-27]	HA722815	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,IP
KD	Lamin A + Lamin C Recombinant Mouse mAb [A6F5-R]	HA601274	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	Lamin A + Lamin C Recombinant Rabbit mAb [JE51-60]	ET7110-12	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	Lamin A + Lamin C Recombinant Rabbit mAb [PSH05-20]	HA722245	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	Lamin A + Lamin C Recombinant Rabbit mAb [PSH05-20] - BSA and Azide free	HA750991	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	Lamin A + Lamin C Recombinant Mouse mAb [A6F5-R] - BSA and Azide free	HA610149	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	Lamin A + Lamin C Mouse mAb [A6F5]	HA600093	Human, Mouse, Rat	WB,IHC-P
	Lamin B1 Rabbit pAb	RT1355	Human, Mouse, Rat	WB,IP,IF,IHC-P
	iFluor™ 647 Conjugated Lamin B1 Recombinant Rabbit mAb [SI17-07]	HA720172F	Human, Mouse, Rat	IF-Tissue
	iFluor™ 488 Conjugated Lamin B1 Recombinant Rabbit mAb [SI17-07]	HA720171F	Human, Mouse, Rat	IF-Tissue
KD	Lamin B1 Rabbit pAb	R1508-1	Human, Rat, Mouse	WB
KO	Lamin B1 Recombinant Rabbit mAb [SI17-07] - BSA and Azide free	HA750098	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,CUT&Tag-seq
KO	Lamin B1 Recombinant Rabbit mAb [SI17-07]	ET1606-27	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,CUT&Tag-seq
	Lamin B2 Rabbit pAb	R1408-4	Human	WB
	Lamin B2 Recombinant Mouse mAb [J9-C4-R]	HA601289	Human, Monkey	WB,IF-Cell,IHC-P,IF-Tissue
	Lamin B2 Mouse mAb [D5-C2]	M1505-11	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P
	Lamin B2 Recombinant Mouse mAb [J9-C4-R] - BSA and Azide free	HA610162	Human, Monkey	WB,IF-Cell,IHC-P,IF-Tissue
	Lamin B2 Mouse mAb [J9-C4]	M1505-10	Human	WB,IHC-P
	Lamin B2 Rabbit pAb	ER0403	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Lamin B Receptor Rabbit pAb	ER1912-18	Human, Mouse	WB,IF-Cell,FC
KD	Lamin B Receptor / LBR Recombinant Rabbit mAb [PSH02-93]	HA721912	Human	WB,IHC-P,IF-Cell,FC
KD	Lamin B Receptor / LBR Recombinant Rabbit mAb [PSH02-93] - BSA and Azide free	HA750848	Human	WB,IHC-P,IF-Cell,FC
Relative expression	Laminin 2 alpha Recombinant Rabbit mAb [PSH0-94] - BSA and Azide free	HA750674	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	Laminin 2 alpha Recombinant Rabbit mAb [PSH0-94]	HA721532	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Laminin 5 alpha 3 Recombinant Rabbit mAb [JJ0957]	ET1701-99	Human	WB,IF-Cell
	Laminin alpha-1 Rabbit pAb	ER63832	Human, Mouse	WB,IHC
	Laminin alpha-2 Rabbit pAb	ER63833	Human, Mouse	IHC-P,IF-Cell
	Laminin alpha-4 Rabbit pAb	ER63834	Human, Mouse	WB,IHC-P,IF-Cell
	Laminin alpha-5 Rabbit pAb	ER63835	Human, Mouse	WB,IHC-P
	Laminin beta 1 Recombinant Rabbit mAb	IRS234RB	Mouse	mIHC
	Laminin beta 1 Recombinant Rabbit mAb [JM099-03]	ET1703-14	Human, Mouse, Rat	WB,IHC-P,FC,mIHC
	Laminin beta 1 Recombinant Rabbit mAb	IRS164RB	Human	mIHC
	Laminin beta-2 Rabbit pAb	ER63836	Human, Mouse, Rat	WB,IHC-P
KD	Laminin gamma 1 Recombinant Rabbit mAb [JE40-75]	HA721662	Human	WB
	Laminin gamma-1 Rabbit pAb	ER63838	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell
	Laminin gamma-2 Rabbit pAb	ER65562	Human, Rat	WB,IHC-P
	Laminin gamma-3 Rabbit pAb	ER63839	Human, Mouse	WB,IHC-P,IF-Cell
	LAP1B Rabbit pAb	ER63840	Human	WB
	LASP-1 Rabbit pAb	ER65502	Human, Mouse, Rat	WB
	Layilin Rabbit pAb	ER65494	Human, Mouse, Rat	WB,IHC-P
Cell treatment, KD	LC3B Recombinant Rabbit mAb [JJ090-6] - BSA and Azide free	HA750319	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,IP
Cell treatment	LC3B Recombinant Antibody [JJ090-6] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610352	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment, KD	LC3B Recombinant Rabbit mAb [JJ090-6]	ET1701-65	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,IP,mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	LC3B Recombinant Antibody [JJ090-6] - Mouse IgG1 (Chimeric)	HA601532	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	LECT1 Rabbit pAb	ER62052	Human, Rat, Mouse	WB
	LECT2 Rabbit pAb	ER61907	Human, Mouse	WB
	Legumain Recombinant Rabbit mAb [JE35-56]	HA721645	Human, Mouse, Rat	WB,IHC-P
	Legumain Rabbit pAb	ER1912-26	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	LI Cadherin Rabbit pAb	ER1804-02	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Relative expression	LI Cadherin Recombinant Rabbit mAb [PD00-49]	HA721643	Human	WB,IHC-P,IF-Tissue
Relative expression	LI Cadherin Recombinant Rabbit mAb [JE60-34]	HA721895	Human, Mouse	WB,IHC-P,IF-Tissue,IF-Cell,FC
	LIMA1 Rabbit pAb	ER62400	Human	WB
	LIMK1 Rabbit pAb	ER1912-36	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	LIMK2 Recombinant Rabbit mAb [PSH08-82]	HA723034	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,IP
	LIMK2 Recombinant Rabbit mAb [PSH08-82] - BSA and Azide free	HA751247	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,IP
Relative expression	LIM Kinase 1 Recombinant Rabbit mAb [PSH08-68]	HA723020	Human, Mouse, Rat, Monkey	WB,IHC-P,FC
Relative expression	LIM Kinase 1 Recombinant Rabbit mAb [PSH08-68] - BSA and Azide free	HA751233	Human, Mouse, Rat, Monkey	WB,IHC-P,FC
	LIM Kinase 1 Rabbit pAb	HA500189	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	LIN7A Rabbit pAb	ER60802	Human, Rat, Mouse	WB
	Liprin alpha 2 Rabbit pAb	ER1912-38	Mouse, Rat	WB,IHC-P,IF-Tissue
	L-Lactyl-alpha Tubulin (K40) Rabbit pAb	HA500554	Human, Mouse, Rat, Zebrafish	WB,IHC-P
	Loricrin Recombinant Rabbit mAb [JE40-03]	HA721547	Human, Mouse, Rat	IHC-P,IF-Tissue
	Loricrin Recombinant Rabbit mAb	IRS165RB	Mouse	mIHC
	LOX Recombinant Rabbit mAb [PSH15-41]	HA723747	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	LOX Recombinant Rabbit mAb [JU30-23]	ET1706-31	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IHC-P,FC,IF-Tissue
	LOX Recombinant Rabbit mAb [PSH15-41] - BSA and Azide free	HA751556	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	LOX-1 Mouse mAb [2G5]	RT1367	Human, Mouse, Rat	WB,IP,IF,IHC-P
	LOXL2 Rabbit pAb	ER1803-41	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	LOXL2 Recombinant Rabbit mAb [JU32-53]	ET1706-11	Mouse, Human, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	LOXL4 Rabbit pAb	ER1912-43	Human	WB
	LPP Recombinant Rabbit mAb [JE55-29]	ET7111-02	Human, Mouse, Rat	WB,IHC-P
	LRG1 Rabbit pAb	HA500576	Human, Mouse	WB
	LRG1 Recombinant Rabbit mAb [JB34-31]	ET7107-11	Human	WB,IHC-P,IP
	LTBP1 Rabbit pAb	ER61833	Human, Mouse, Rat	WB
	Lubricin Rabbit pAb	ER1912-53	Human, Mouse, Rat	WB
	Lumican Recombinant Rabbit mAb [JE11-45]	ET7109-24	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	LYRIC Mouse mAb [6-F4]	EM1901-61	Human	WB,FC
	LYRIC Mouse mAb [C6-C10]	EM1901-60	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	LYRIC Recombinant Mouse mAb [C6-C10-R] - BSA and Azide free	HA610068	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	LYRIC Recombinant Mouse mAb [C6-C10-R]	HA601188	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	LYVE1 Recombinant Rabbit mAb [PSH18-29]	HA723982	Human	IHC-P
KD	LYVE1 Recombinant Rabbit mAb [JF0979]	ET1702-29	Human	WB
	LYVE1 Recombinant Rabbit mAb	IRS169RB	Mouse	mIHC
Relative expression	LYVE1 Recombinant Rabbit mAb [PSH07-67]	HA722867	Mouse	WB,IHC-P
KD	LYVE1 Recombinant Rabbit mAb [JF0979] - BSA and Azide free	HA750343	Human	WB
Relative expression	LYVE1 Recombinant Rabbit mAb [PSH07-67] - BSA and Azide free	HA751179	Mouse	WB,IHC-P
	LYVE1 Recombinant Rabbit mAb	IRS256RB	Human	mIHC
	MAdCAM-1 Rabbit pAb	ER1912-60	Human, Mouse	WB,IHC-P,IF-Tissue
	MAdCAM-1 Rat mAb [IgG2a]	RT1375	Human, Mouse, Rat	WB,IP,IF,FC
	MAGI1 Recombinant Rabbit mAb [PSH10-91]	HA723258	Human, Mouse, Rat	WB
	MAGI1 Recombinant Rabbit mAb [PSH10-91] - BSA and Azide free	HA751382	Human, Mouse, Rat	WB
	MAGI3 Rabbit pAb	ER1912-66	Mouse	WB
	Mammal VEGF-C enzyme-linked immunoassay kit	EH0113	Mammal	ELISA
	MAP1A Rabbit pAb	ER60830	Human, Mouse, Rat	IHC-P
	MAP1B Rabbit pAb	ER60831	Human, Rat, Mouse	IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	MAP1LC3A Recombinant Rabbit mAb [SN07-40]	ET1611-37	Human, Mouse, Rat	WB,IHC-P
Cell treatment	MAP1LC3A Recombinant Rabbit mAb [ST47-03]	ET1609-26	Human, Mouse, Rat	WB,IF-Cell,FC,IP
Relative expression	MAP2 Recombinant Rabbit mAb [PSH08-73] - BSA and Azide free	HA751238	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IHC-Fr,IF-Cell
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Mouse IgG1 (Chimeric)	HA601458	Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610293	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Guinea pig IgG2 (Chimeric)	HA601457	Mouse, Rat	IHC-Fr,WB
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Rat IgG1 (Chimeric)	HA601459	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell
	MAP2 Rabbit pAb	HA500177	Human, Mouse, Rat	WB,IF-Cell,IHC-P,mIHC
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Chicken IgY (Chimeric)	HA601467	Mouse, Rat	IHC-Fr,WB
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610291	Mouse, Rat	IHC-Fr,WB
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Chicken IgY (Chimeric) - BSA and Azide free	HA610301	Mouse, Rat	IHC-Fr,WB
Relative expression	MAP2 Recombinant Antibody [PSH08-73] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610292	Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Cell
Relative expression	MAP2 Recombinant Rabbit mAb [PSH08-73]	HA723025	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IHC-Fr,IF-Cell
	MAP2 Recombinant Rabbit mAb [PSH08-73]	IRS071	Mouse	mIHC
	MAP4 Rabbit pAb	ER1912-71	Human, Mouse	WB
KD	MAPKAP Kinase 2 Recombinant Rabbit mAb [JU30-31]	ET7106-83	Human, Mouse	WB,IP
	MARK1/2/3/4 Rabbit pAb	ER63911	Human, Mouse, Rat	WB,IHC-P,IF-Cell
KD	MARK2 Recombinant Rabbit mAb [JE40-97]	ET7109-20	Human, Mouse, Rat	WB,IHC-P,FC
	MARK3 Mouse mAb [H12-G2]	EM1709-06	Human, Mouse	WB,IHC-P,FC
	MARK3 Recombinant Rabbit mAb [JE44-84]	ET7109-39	Human, Mouse, Rat	WB,IP,IF-Cell,IF-Tissue,IHC-P,FC
	MARK4 Rabbit pAb	ER63912	Human, Mouse	IF-Cell
	Mast Cell Chymase Recombinant Rabbit mAb [JB74-32]	ET7107-19	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	MATN3 Rabbit pAb	ER61725	Human, Mouse	WB
	Matriptase Rabbit pAb	ER1912-79	Human	IF-Cell
	MBPHL Rabbit pAb	ER65968	Human, Mouse, Rat	WB
	M Cadherin Rabbit pAb	ER1912-59	Human	WB
KD	MCAK Recombinant Rabbit mAb [PSH01-51]	HA721694	Human, Mouse, Rat	WB,IHC-P
KD	MCAK Recombinant Rabbit mAb [PSH01-51] - BSA and Azide free	HA750727	Human, Mouse, Rat	WB,IHC-P
	Metadherin / LYRIC Rabbit pAb	R1511-10	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	MFAP5 Rabbit pAb	ER62023	Human, Mouse	WB
	MGT5A Rabbit pAb	ER62106	Human, Mouse, Rat	WB
	Mindin Rabbit pAb	ER63958	Human, Mouse, Rat, Monkey	WB
	MIP-T3 Rabbit pAb	ER63963	Human	WB,IHC-P,IF-Cell
	MKLP1 Recombinant Rabbit mAb [SY02-74] - BSA and Azide free	HA750129	Human, Mouse, Rat	WB,IF-Cell
	MKLP1 Recombinant Rabbit mAb [SY02-74]	ET1607-63	Human, Mouse, Rat	WB,IF-Cell
	MLC-2 Rabbit pAb	ER63968	Human, Mouse, Rat	IHC-P,IF-Cell
	MLP3C Rabbit pAb	ER62516	Human	WB
	MMP-1 Rabbit pAb	ER31211	Human, Mouse, Monkey	WB,IHC-P
Relative expression	MMP-1 Recombinant Rabbit mAb [PSH09-56] - BSA and Azide free	HA751303	Human	WB,IHC-P
Relative expression	MMP-1 Recombinant Rabbit mAb [PSH09-56]	HA723118	Human	WB,IHC-P
	MMP-10 Rabbit pAb	ER63978	Human, Mouse, Rat	IHC-P,IF-Cell,WB
	MMP-11 Recombinant Rabbit mAb [SN74-08] - BSA and Azide free	HA750250	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	MMP-11 Recombinant Rabbit mAb [SN74-08]	ET1611-33	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
KD	MMP-12 Recombinant Rabbit mAb [SR03-23] - BSA and Azide free	HA750052	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
KD	MMP-12 Recombinant Rabbit mAb [SR03-23]	ET1602-42	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MMP-13 Recombinant Rabbit mAb [JF0893]	ET1702-14	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	MMP-14 Recombinant Rabbit mAb [SJ18-09]	ET1606-48	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	MMP-14 Recombinant Rabbit mAb [SJ18-09] - BSA and Azide free	HA750108	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	MMP-14 Mouse mAb [11A-1D]	EM1706-30	Human	IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	MMP-15 Rabbit pAb	ER1913-06	Mouse	WB
	MMP-16 Rabbit pAb	ER64067	Human, Mouse, Rat	WB,IHC-P
	MMP-17 Recombinant Rabbit mAb [JU33-03]	ET7106-84	Human, Mouse, Rat	IF-Cell,IF-Tissue,IHC-P,WB
	MMP-19 Rabbit pAb	ER1913-07	Human, Mouse	WB
	MMP-2 Recombinant Rabbit mAb [SI15-04]	ET1606-4	Human, Mouse, Rat	WB,IHC-P,IP,IF-Tissue
Relative expression	MMP-2 Recombinant Rabbit mAb [PSH11-26] - BSA and Azide free	HA751405	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	MMP-2 Mouse mAb [B11-E6]	EM1706-99	Human, Mouse	WB,FC
Relative expression	MMP-2 Recombinant Rabbit mAb [PSH11-26]	HA723308	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	MMP-23 Rabbit pAb	ER63979	Human, Mouse, Rat	WB
	MMP-3 Rabbit pAb	ER1706-77	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	MMP-3 Recombinant Rabbit mAb [JM46-22]	ET1705-98	Human, Mouse, Rat	WB,IHC-P
	MMP-7 Recombinant Rabbit mAb [PSH08-91]	HA723043	Human	WB,IF-Cell,FC,IHC-P
	MMP-7 Rabbit pAb	HA500305	Human	WB,IF-Cell,IHC-P,FC
	MMP-7 Rabbit pAb	ER1913-08	Human, Mouse	WB,IHC-P,IF-Tissue
	MMP-7 Recombinant Rabbit mAb [PSH08-91] - BSA and Azide free	HA751256	Human	WB,IF-Cell,FC,IHC-P
	MMP-8 Recombinant Rabbit mAb [JE59-38]	HA721686	Human, Mouse, Rat	IHC-P
	MMP-8 Rabbit pAb	ER63980	Human, Mouse, Rat	WB,IHC-P
Cell treatment, KD	MMP-9 Recombinant Rabbit mAb [JA80-73]	ET1704-69	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC
	MMP-9 Mouse mAb [2G3]	RT1401	Human	WB,IP,IF,IHC-P
Cell treatment	MMP-9 Recombinant Rabbit mAb [PSH06-33] - BSA and Azide free	HA751070	Human	WB,IHC-P
	MMP-9 Recombinant Rabbit mAb	IRS172RB	Human	mIHC
Cell treatment, KD	MMP-9 Recombinant Rabbit mAb [JA80-73] - BSA and Azide free	HA750420	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC
Cell treatment	MMP-9 Recombinant Rabbit mAb [PSH06-33]	HA722665	Human	WB,IHC-P
	MMP-9 Recombinant Mouse mAb [10A1]	EM1801-22	Human	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	Moesin Recombinant Rabbit mAb [SC69-01]	ET1610-65	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Mouse alpha 1 antitrypsin enzyme-linked immunoassay kit (one step)	EMY0052	Mouse	ELISA
	Mouse CCL2/MCP-1 enzyme-linked immunoassay kit	EM0025	Mouse	ELISA
	Mouse CCL3/MIP-1 alpha enzyme-linked immunoassay kit	EM0031	Mouse	ELISA
	Mouse CCL5 / RANTES enzyme-linked immunoassay kit	EM0032	Mouse	ELISA
	Mouse Coagulation Factor XII enzyme-linked immunoassay kit	EM0041	Mouse	ELISA
	Mouse Coagulation Factor XII enzyme-linked immunoassay kit(one step)	EMY0001	Mouse	ELISA
	HRP Conjugated Mouse CRP Recombinant Rabbit mAb [PSH12-49] - Detector	HA725075H	Mouse	ELISA(Det),ELISA
	Mouse CRP Recombinant Rabbit mAb [PSH12-49] - BSA and Azide free (Detector)	HA725074	Mouse	ELISA(Det)
	Mouse CRP Recombinant Rabbit mAb [PSH12-48] - BSA and Azide free (Capture)	HA725073	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CRP Recombinant Rabbit mAb [PSH12-49] - Detector	HA725075B	Mouse	ELISA(Det),ELISA
	Mouse CRP/C-Reactive Protein enzyme-linked immunoassay kit	EM0035	Mouse	ELISA
	Mouse CXCL1/KC enzyme-linked immunoassay kit	EM0017	Mouse	ELISA
	Mouse CXCL2/MIP-2 enzyme-linked immunoassay kit	EM0018	Mouse	ELISA
	Mouse CXCL9/MIG enzyme-linked immunoassay kit	EM0015	Mouse	ELISA
	Mouse GM-CSF/CSF2 enzyme-linked immunoassay kit	EM0011	Mouse	ELISA
	Mouse IFN alpha-1 enzyme-linked immunoassay kit	EM0049	Mouse	ELISA
	Mouse IFN- β enzyme-linked immunoassay kit	EM0040	Mouse	ELISA
	Mouse IFN- γ enzyme-linked immunoassay kit	EM0009	Mouse	ELISA
	Mouse IL-10 Recombinant Rabbit mAb [PS01-93] - BSA and Azide free (Detector)	HA721518	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IL-10 Recombinant Rabbit mAb [PS01-85] - BSA and Azide free (Capture)	HA721517	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse IL-10 Recombinant Rabbit mAb [PS01-93] - Detector	HA721519	Mouse	ELISA(Det)
	Mouse IL-10 enzyme-linked immunoassay kit	EM0005	Mouse	ELISA
	Mouse IL-12/IL-23 p40 enzyme-linked immunoassay kit	EM0012	Mouse	ELISA
	Mouse IL-13 enzyme-linked immunoassay kit	EM0007	Mouse	ELISA
	Mouse IL-17A enzyme-linked immunoassay kit	EM0008	Mouse	ELISA
	Mouse IL-17E enzyme-linked immunoassay kit	EM0028	Mouse	ELISA
	Mouse IL-18 enzyme-linked immunoassay kit	EM0034	Mouse	ELISA
	Mouse IL-1 alpha/IL-1 F1 enzyme-linked immunoassay kit	EM0047	Mouse	ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - BSA and Azide free (Detector)	HA721977	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-32] - Detector	HA722625B	Mouse	ELISA(Det),ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH03-31] - BSA and Azide free (Capture)	HA721976	Mouse	ELISA(Cap)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67] - BSA and Azide free (Detector)	HA722167	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-67]	HA722608B	Mouse	ELISA(Det)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH04-66] - BSA and Azide free (Capture)	HA722166	Mouse	ELISA(Cap)
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH09-58] - BSA and Azide free	HA723120	Mouse	ELISA
	Mouse IL-1 beta Recombinant Rabbit mAb [PSH09-59] - BSA and Azide free	HA723121	Mouse	ELISA
	Mouse IL-1 beta Super Sensitivity enzyme-linked immunoassay kit	EMSS0029	Mouse	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse IL-1 β enzyme-linked immunoassay kit	EM0029	Mouse	ELISA
	Biotin Conjugated Mouse IL-2 Recombinant Rabbit mAb [PSH04-28]	HA722617B	Mouse	ELISA(Det)
	Mouse IL-2 Recombinant Rabbit mAb [PSH04-27] - BSA and Azide free (Capture)	HA722110	Mouse	ELISA(Cap)
	Mouse IL-2 Recombinant Rabbit mAb [PSH04-28] - BSA and Azide free (Detector)	HA722111	Mouse	ELISA(Det)
	Mouse IL-21 enzyme-linked immunoassay kit	EM0020	Mouse	ELISA
	Mouse IL-22 enzyme-linked immunoassay kit	EM0021	Mouse	ELISA
	Mouse IL-2 enzyme-linked immunoassay kit	EM0002	Mouse	ELISA
	Mouse IL-33 enzyme-linked immunoassay kit	EM0033	Mouse	ELISA
	Mouse IL-3 enzyme-linked immunoassay kit	EM0048	Mouse	ELISA
	Mouse IL-4 enzyme-linked immunoassay kit	EM0003	Mouse	ELISA
	Mouse IL-5 enzyme-linked immunoassay kit	EM0023	Mouse	ELISA
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-58] - BSA and Azide free (Capture)	HA721411	Mouse	ELISA(Cap)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-60] - BSA and Azide free (Detector)	HA721413	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PSH0-60] - Detector	HA722647B	Mouse	ELISA(Det)
	HRP Conjugated Mouse IL-6 Mouse mAb [3-A4]	M1303-2	Mouse	WB,ELISA
	Mouse IL-6 Recombinant Rabbit mAb [PS01-83] - BSA and Azide free (Detector)	HA721601	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PS01-83] - BSA and Azide free (Detector)	HA721602	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PSH14-47] - BSA and Azide free (Capture)	HA723642	Mouse	ELISA(Cap)
	Mouse IL-6 Recombinant Rabbit mAb [PSH0-59] - BSA and Azide free (Detector)	HA721412	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PS01-74] - Detector	HA723643B	Mouse	ELISA(Det)
	Mouse IL-6 Recombinant Rabbit mAb [PS01-74] - BSA and Azide free (Capture/Detector)	HA721600	Mouse	ELISA(Cap),ELISA(Det)
	Biotin Conjugated Mouse IL-6 Recombinant Rabbit mAb [PSH0-59] - Detector	HA722646B	Mouse	ELISA(Det)
	Mouse IL-6 enzyme-linked immunoassay kit	EM0004	Mouse	ELISA
	Mouse IL-9 enzyme-linked immunoassay kit	EM0022	Mouse	ELISA
	Mouse Leptin enzyme-linked immunoassay kit	EM0027	Mouse	ELISA
	Mouse LIF enzyme-linked immunoassay kit	EM0038	Mouse	ELISA
	Biotin Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-80] - Detector	HA723786B	Mouse	ELISA(Det),ELISA
	HRP Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - Detector	HA723788H	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - Detector	HA723788B	Mouse	ELISA(Det),ELISA
	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-81] - BSA and Azide free (Detector)	HA723787	Mouse	ELISA(Det)
	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-79] - BSA and Azide free (Capture)	HA723784	Mouse	ELISA(Cap)
	Mouse MMP-9 Recombinant Rabbit mAb [PSH15-80] - BSA and Azide free (Capture/Detector)	HA723785	Mouse	ELISA(Det),ELISA(Cap)
	HRP Conjugated Mouse Osteopontin Recombinant Rabbit mAb [PSH17-50] - Detector	HA723925H	Mouse	ELISA(Det),ELISA
	Mouse Osteopontin Recombinant Rabbit mAb [PSH17-49] - BSA and Azide free (Capture)	HA723924	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse Osteopontin Recombinant Rabbit mAb [PSH17-50] - Detector	HA723925B	Mouse	ELISA(Det),ELISA
	Mouse Osteopontin Recombinant Rabbit mAb [PSH17-50] - BSA and Azide free (Detector)	HA723925	Mouse	ELISA(Det)
	Mouse PD-1 enzyme-linked immunoassay kit	EM0037	Mouse	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse PD-L1 enzyme-linked immunoassay kit	EM0036	Mouse	ELISA
	Biotin Conjugated Mouse Pro-Collagen I alpha 1 Recombinant Rabbit mAb [PSH07-99]	HA722936B	Mouse	ELISA(Det)
	Mouse Pro-Collagen I alpha 1 Recombinant Rabbit mAb [PSH07-99] - BSA and Azide free (Detector)	HA722935	Mouse	ELISA(Det)
	HRP Conjugated Mouse Pro-Collagen I alpha 1 Recombinant Rabbit mAb [PSH07-99] - Detector	HA722936H	Mouse	ELISA(Det)
	Mouse Pro-Collagen I alpha 1 Recombinant Rabbit mAb [PSH07-98] - BSA and Azide free (Capture)	HA722934	Mouse	ELISA(Cap)
	Mouse Pro-collagen I alpha 1/COL1A1 enzyme-linked immunoassay kit (one step)	EMY0045	Mouse	ELISA
	Mouse RANKL/TRANCE/TNFSF11 enzyme-linked immunoassay kit	EM0044	Mouse	ELISA
	Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-24] - BSA and Azide free (Capture)	HA723911	Mouse, Rat	ELISA(Cap)
	Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-25] - BSA and Azide free (Detector)	HA723912	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse/Rat E-Cadherin Recombinant Rabbit mAb [PSH17-25] - Detector	HA723912B	Mouse, Rat	ELISA(Det),ELISA
	Biotin Conjugated Mouse/Rat GDF15 Recombinant Rabbit mAb [PSH16-65] - Detector	HA723867B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat GDF15 Recombinant Rabbit mAb [PSH16-66] - BSA and Azide free (Capture)	HA723868	Mouse, Rat	ELISA(Cap)
	Mouse/Rat GDF15 Recombinant Rabbit mAb [PSH16-64] - BSA and Azide free (Capture)	HA723866	Mouse, Rat	ELISA(Cap)
	Mouse/Rat GDF15 Recombinant Rabbit mAb [PSH16-65] - BSA and Azide free (Detector)	HA723867	Mouse, Rat	ELISA(Det)
	Mouse/Rat IGF1 enzyme-linked immunoassay kit (one step)	EMY0024	Mouse, Rat	ELISA
	Mouse/Rat IGF-I/IGF-1 enzyme-linked immunoassay kit	EM0024	Mouse, Rat	ELISA
	Biotin Conjugated Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30]	HA722626B	Mouse, Rat	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-29] - BSA and Azide free (Capture)	HA721974	Mouse, Rat	ELISA(Cap)
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30] - BSA and Azide free (Detector)	HA721975	Mouse, Rat	ELISA(Det)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77] - BSA and Azide free (Detector)	HA722041	Mouse	ELISA(Det)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-76] - BSA and Azide free (Capture)	HA722040	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77]	HA722623B	Mouse	ELISA(Det)
	Mouse TNF- α enzyme-linked immunoassay kit	EM0010	Mouse	ELISA
	Mouse TREM2 enzyme-linked immunoassay kit	EM0039	Mouse	ELISA
	Biotin Conjugated Mouse VCAM1 Recombinant Rabbit mAb [PSH16-32] - Detector	HA723830B	Mouse	ELISA(Det),ELISA
	Mouse VCAM1 Recombinant Rabbit mAb [PSH16-31] - BSA and Azide free (Capture)	HA723829	Mouse	ELISA(Cap)
	Mouse VCAM1 Recombinant Rabbit mAb [PSH16-33] - BSA and Azide free (Capture)	HA723831	Mouse	ELISA(Cap)
	HRP Conjugated Mouse VCAM1 Recombinant Rabbit mAb [PSH16-32] - Detector	HA723830H	Mouse	ELISA(Det),ELISA
	Mouse VCAM1 Recombinant Rabbit mAb [PSH16-32] - BSA and Azide free (Detector)	HA723830	Mouse	ELISA(Det)
	Mouse VEGF enzyme-linked immunoassay kit	EM0013	Mouse	ELISA
	MRas Recombinant Rabbit mAb [JE64-57]	HA722579	Human, Mouse, Rat	WB
	MRCK alpha Rabbit pAb	ER63998	Human, Mouse, Rat	WB
	M-RIP Rabbit pAb	ER64007	Human, Mouse, Rat	IHC-P,WB
	MTAP Rabbit pAb	ER1901-47	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	MTSS1 Rabbit pAb	ER1913-35	Human, Mouse	WB
	Muc1 Rabbit pAb	ER1902-10	Human	WB,IF-Cell,IHC-P
Relative expression	MUC1 Mouse mAb [PSH0-47]	HA601140	Human	IHC-P,WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	iFluor™ 647 Conjugated MUC1 Mouse mAb [PSH0-49]	HA600108F	Human	FC
	MUC1 Mouse mAb	IRS174MS	Human	mIHC
Relative expression	MUC1 Mouse mAb [PSH0-48]	HA601141	Human	IHC-P,FC,IF-Cell
Relative expression	MUC1 Mouse mAb [PSH0-49]	HA601142	Human	IHC-P,FC,IF-Cell
	iFluor™ 647 Conjugated MUC1 Mouse mAb [PSH0-48]	HA600107F	Human	FC
Relative expression	MUC1 Mouse mAb [A4D9]	EM1902-32	Human	IHC-P,FC,IF-Cell
Relative expression	MUC1 Mouse mAb [A9G4]	HA601139	Human	IHC-P,IF-Cell
Relative expression, KD	MUC1 Recombinant Rabbit mAb [SN06-80] - BSA and Azide free	HA750245	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Relative expression, KD	MUC1 Recombinant Rabbit mAb [SN06-80]	ET1611-14	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	MUC13 Rabbit pAb	ER1913-37	Human, Mouse	WB,IHC-P,IF-Tissue
	MUC15 Rabbit pAb	ER60851	Human, Mouse	WB
Relative expression	MUC16 Recombinant Rabbit mAb [PD01-13]	HA721217	Human	IHC-P,WB,IF-Tissue
Relative expression	MUC16 Recombinant Rabbit mAb [SN0715]	ET1611-73	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	MUC16 Recombinant Rabbit mAb [SN69-06]	ET1611-74	Human	WB,IF-Cell,IF-Tissue,IHC-P,FC
	MUC2 Recombinant Rabbit mAb	IRS175RB	Human	mIHC
	MUC2 Recombinant Rabbit mAb [JA10-92] - BSA and Azide free	HA750401	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	MUC2 Recombinant Rabbit mAb	IRS235RB	Mouse	mIHC
	MUC2 Recombinant Rabbit mAb [JA10-92]	ET1704-06	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	MUC3A Rabbit pAb	ER62460	Human	IHC-P
	MUC4 Recombinant Rabbit mAb [JM19-36]	ET1705-13	Human, Rat	WB,IHC-P,FC
	MUC5AC Mouse mAb [A3H3]	EM1902-02	Human	IHC-P,FC
	MUC5AC Recombinant Mouse mAb [A10G6-R] - BSA and Azide free	HA610065	Human	IHC-P
	MUC5AC Recombinant Mouse mAb [A10G5-R]	HA601181	Human, Mouse	IHC-P
	MUC5AC Recombinant Mouse mAb [A10G5-R] - BSA and Azide free	HA610064	Human, Mouse	IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	MUC5AC Mouse mAb [G6-H8]	EM1707-36	Human	IHC-P,FC
	MUC5AC Recombinant Mouse mAb [A10G6-R]	HA601182	Human	IHC-P
	MUC5B Rabbit pAb	ER60852	Human	IHC-P
	MUC5B Mouse mAb [7F5A1]	EM1711-75	Human	IHC-P,FC
	MUC7 Rabbit pAb	ER60853	Human	WB
	MUSK Recombinant Rabbit mAb [PSH02-29]	HA721809	Human, Mouse, Rat	WB,IHC-P
	MUSK Recombinant Rabbit mAb [PSH05-68]	HA722420	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	MUSK Recombinant Rabbit Multiclonal Antibody [PSH03-51]	HA722010	Human, Mouse, Rat	WB,IHC-P
	MUSK Recombinant Rabbit mAb [PSH05-68] - BSA and Azide free	HA751022	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	MUSK Recombinant Rabbit mAb [PSH02-28] - BSA and Azide free	HA750796	Human, Mouse, Rat	WB,IHC-P
	MUSK Recombinant Rabbit mAb [PSH02-28]	HA721808	Human, Mouse, Rat	WB,IHC-P
	MUSK Recombinant Rabbit Multiclonal Antibody [PSH03-51] - BSA and Azide free	HA750885	Human, Mouse, Rat	WB,IHC-P
	MUSK Recombinant Rabbit mAb [PSH02-29] - BSA and Azide free	HA750797	Human, Mouse, Rat	WB,IHC-P
	MY18B Rabbit pAb	ER60878	Human	IHC-P
	MYBPC1 Rabbit pAb	ER64071	Human, Mouse, Rat	WB
	MYH14 Rabbit pAb	ER64075	Human, Mouse	WB,IHC-P
Relative expression	MYH2 Recombinant Rabbit mAb [JE53-30]	HA721862	Human, Mouse, Rat	WB,IHC-P
	MYH3 Rabbit pAb	ER60864	Human, Rat, Mouse	WB
	MYH4 Rabbit pAb	ER64077	Human	IHC-P
	MYH7B Recombinant Rabbit mAb [JE54-44]	ET7110-69	Rat, Mouse, Human	WB
	MYH-pan Rabbit pAb	ER65575	Human, Mouse, Rat	WB,IHC-P
	MYL3 Rabbit pAb	ER64078	Human	WB
	MYL4 Rabbit pAb	HA500471	Human, Mouse, Rat	WB
	MYL6B Rabbit pAb	ER60329	Human	WB
	MYL7 Rabbit pAb	ER64080	Human, Mouse, Rat, Monkey	WB
	MYL9 Recombinant Rabbit mAb [PSH02-07]	HA721766	Human, Mouse, Rat	WB,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	MYL9 Recombinant Rabbit mAb [JE34-63]	HA721361	Human, Mouse, Rat	WB,IHC-P
	MYL9 Recombinant Rabbit mAb [PSH02-07] - BSA and Azide free	HA750777	Human, Mouse, Rat	WB,FC
	MYLIP Rabbit pAb	ER64081	Human, Mouse, Rat	WB,IHC-P
	MYLK2 Rabbit pAb	ER60330	Human, Mouse, Rat	WB
	MYLPF Rabbit pAb	ER1913-44	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	MYO10 Rabbit pAb	ER60875	Human, Rat, Mouse	WB
	MYO19 Recombinant Rabbit mAb [JE64-30]	HA721051	Human, Mouse, Rat	WB,IF-Cell
	MYO1A Rabbit pAb	ER60868	Human, Rat, Mouse	WB
	MYO5B Rabbit pAb	ER60872	Human, Rat, Mouse	WB
	MYO7A Rabbit pAb	ER60874	Human, Mouse	WB
	Myomegalin Rabbit pAb	ER2001-16	Mouse	WB
	Myomesin-1 Rabbit pAb	ER64083	Human	IHC-P
	Myomesin-2 Rabbit pAb	ER64084	Human, Mouse, Monkey	WB
	Myosin heavy chain Recombinant Rabbit mAb [JF097-7]	ET1702-88	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Myosin Id Rabbit pAb	ER64085	Human, Mouse, Rat	IHC-P
	Myosin IXb Rabbit pAb	ER64086	Human, Mouse, Rat	WB,IHC-P
	Myosin Light Chain 2 Recombinant Rabbit mAb [SN67-09] - BSA and Azide free	HA750244	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IP
	Myosin Light Chain 2 Recombinant Rabbit mAb [SN67-09]	ET1611-13	Human, Mouse, Rat, Zebrafish	WB,IHC-P,IP
	Myosin Light Chain 2 Rabbit pAb	R1512-19	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Myosin light chain kinase / MLCK Rabbit pAb	R1510-18	Human, Mouse	IF-Cell,IHC-P,FC
	Myosin light chain kinase / MLCK Rabbit pAb	R1511-26	Human, Mouse	IF-Cell,IHC-P,FC,WB
	Myosin light chain kinase / MLCK Recombinant Rabbit mAb [SU40-06]	ET1608-68	Human, Mouse, Rat	WB,IF-Tissue,IHC-P
	Myosin VA Rabbit pAb	ER64087	Human, Mouse, Rat	WB,IHC-P
	Myosin VI Rabbit pAb	ER65422	Human, Mouse, Rat	WB,IHC-P
	Myosin XVI Rabbit pAb	ER64088	Human	WB,IHC-P
	MYOZ2 Rabbit pAb	ER65742	Human, Mouse, Rat	WB
	MYPT1 Rabbit pAb	ER64090	Human, Mouse, Rat	WB
	MYPT3 Rabbit pAb	ER64091	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	NAT10 Recombinant Rabbit mAb [JE55-38]	ET7111-23	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	N Cadherin Recombinant Antibody [SY02-46] - Mouse IgG1 (Chimeric)	HA601464	Human, Mouse, Rat	IHC-Fr,IHC-P
	N Cadherin Recombinant Antibody [SY02-46] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610298	Human, Mouse, Rat	IHC-Fr,IHC-P
Relative expression, KD	N Cadherin Recombinant Rabbit mAb [SY02-46]	ET1607-37	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
Relative expression, KD	N Cadherin Recombinant Rabbit mAb [SY02-46] - BSA and Azide free	HA750119	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	N-Cadherin Rabbit pAb	ER0503	Human, Mouse, Rat	WB,IHC,IF
	N-Cadherin Mouse mAb [A10-4]	M1304-1	Human, Mouse	WB,IHC-P
	NCAM1 / CD56 Recombinant Rabbit mAb	IRS026	Human	mIHC
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-82] - BSA and Azide free	HA751104	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH08-89]	HA723041	Human	WB,IHC-P,IF-Tissue
	NCAM1 / CD56 Recombinant Rabbit mAb [JF1021]	ET1702-43	Human, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP,mIHC
Relative expression	NCAM1 / CD56 Recombinant Mouse mAb [PDH0-07]	HA601126	Human	IHC-P,WB
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53]	HA722756	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
	NCAM1 / CD56 Rabbit pAb	R1204-1	Human, Mouse, Rat	WB,IF-Cell
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-53] - BSA and Azide free	HA751105	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IP
Relative expression	NCAM1 / CD56 Recombinant Rabbit mAb [PSH06-82]	HA722755	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Cell,FC,IF-Tissue
	NCAM1 / CD56 Recombinant Rabbit mAb [JF1021] - BSA and Azide free	HA750346	Human, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	NCK1 Recombinant Rabbit mAb [JU41-80]	ET7106-57	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Nck-2 Rabbit pAb	ER64117	Human, Mouse	IHC-P
	Nck-2 Recombinant Rabbit mAb [JE63-36]	HA722573	Human	WB,IF-Cell,IHC-P,IF-Tissue
	Nephrin Rabbit pAb	ER1803-67	Human, Rat	WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression, Cell treatment	Nephrin Recombinant Rabbit Multiclonal Antibody [PSH07-70] - BSA and Azide free	HA751182	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression, Cell treatment	Nephrin Recombinant Rabbit Multiclonal Antibody [PSH07-70]	HA722886	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
	NES/Kallikrein 10 Rabbit pAb	ER1913-63	Human	WB
	Nesprin 3 Rabbit pAb	ER1913-64	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Nesprin3 Recombinant Rabbit mAb [JE65-22]	HA722581	Human	WB
	Netrin 1 Recombinant Rabbit mAb [JE38-87]	HA721545	Human, Mouse	WB
	Netrin 1 Rabbit pAb	ER1913-65	Human, Rat	IHC-P,IF-Cell,IF-Tissue
	Neurabin 2 Rabbit pAb	ER1913-68	Mouse, Rat	WB,IHC-P,IF-Tissue
	Neurocan Recombinant Rabbit mAb [JE54-15]	ET7110-59	Human	IHC-P
	Neurofilament Rabbit pAb	R1706-17	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Neutrophil Elastase Recombinant Rabbit mAb [JF098-6]	ET1702-78	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
	Neutrophil Elastase Recombinant Rabbit mAb [JF098-6] - BSA and Azide free	HA750360	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P
	NF-H Mouse mAb [2G1]	RT1424	Human, Mouse, Rat	WB,IP,IF,IHC-P
Relative expression	NF-H Recombinant Rabbit mAb [JF99-07] - BSA and Azide free	HA750357	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
Relative expression	NF-H Recombinant Rabbit mAb [JF99-07]	ET1702-72	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
Relative expression	NF-H Recombinant Rabbit mAb [PSH18-14]	HA723977	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	NF-L Recombinant Rabbit mAb [PS02-10]	HA721538	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,IF-Tissue,IHC-Fr
	NF-L Recombinant Antibody [PS02-10] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610297	Mouse, Rat	IHC-Fr
Relative expression	NF-L Recombinant Antibody [PS02-10] - Chicken IgY (Chimeric) - BSA and Azide free	HA610306	Mouse, Rat	IHC-Fr,WB
Relative expression	NF-L Recombinant Antibody [PS02-10] - Rat IgG1 (Chimeric)	HA601409	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
Relative expression	NF-L Recombinant Rabbit mAb [PS02-10] - BSA and Azide free	HA750678	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Cell,IF-Tissue,IHC-Fr

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	NF-L Recombinant Antibody [PS02-10] - Guinea pig IgG2 (Chimeric)	HA601463	Mouse, Rat	IHC-Fr
Relative expression	NF-L Recombinant Rabbit mAb [PSH0-77] - BSA and Azide free	HA750657	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
Relative expression	NF-L Recombinant Antibody [PS02-10] - Chicken IgY (Chimeric)	HA601472	Mouse, Rat	IHC-Fr,WB
Relative expression	NF-L Recombinant Antibody [PS02-10] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610249	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB
Relative expression	NF-L Recombinant Antibody [PS02-10] - Mouse IgG1 (Chimeric)	HA601448	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	NF-L Recombinant Antibody [PS02-10] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610282	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	NF-L Recombinant Rabbit mAb [PSH0-77]	HA721498	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
	NF-L Rabbit pAb	ER65439	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	NF-M Recombinant Antibody [JM11-20] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610279	Mouse, Rat	IHC-Fr,WB
	NF-M Recombinant Antibody [JM11-20] - Guinea pig IgG2 (Chimeric)	HA601444	Mouse, Rat	IHC-Fr,WB
	NF-M Recombinant Antibody [JM11-20] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610278	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	NF-M Recombinant Antibody [JM11-20] - Chicken IgY (Chimeric)	HA601507	Mouse	IHC-Fr
	NF-M Recombinant Antibody [JM11-20] - Mouse IgG1 (Chimeric)	HA601443	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	NF-M Recombinant Antibody [JM11-20] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610248	Human, Mouse, Rat, Cynomolgus monkey	IHC-Fr,IHC-P,WB
	NF-M Recombinant Antibody [JM11-20] - Chicken IgY (Chimeric) - BSA and Azide free	HA610337	Mouse	IHC-Fr
	NF-M Recombinant Antibody [JM11-20] - Rat IgG1 (Chimeric)	HA601408	Human, Mouse, Rat, Cynomolgus monkey	IHC-Fr,IHC-P,WB
Relative expression	NF-M Recombinant Rabbit mAb [JM11-20]	ET1703-57	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IHC-Fr,IF-Tissue
Relative expression	NG2 Recombinant Rabbit mAb [JM10-13] - BSA and Azide free	HA750372	Human, Mouse, Rat	WB,IHC-P
Relative expression	NG2 Recombinant Rabbit mAb [JM10-13]	ET1703-16	Human, Mouse, Rat	WB,IHC-P
Relative expression	NG2 Recombinant Rabbit mAb [PSH09-98] - BSA and Azide free	HA751329	Human	WB,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	NG2 Recombinant Rabbit mAb [PSH09-98]	HA723167	Human	WB,IF-Cell,FC
	NHERF-2 Rabbit pAb	ER65675	Human, Mouse, Rat	WB
	Nidogen Rabbit pAb	ER64192	Human, Mouse, Rat	WB
	non-muscle Myosin IIA Recombinant Rabbit mAb [PSH14-10] - BSA and Azide free	HA751505	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	non-muscle Myosin IIA Recombinant Rabbit mAb [PSH14-10]	HA723600	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	non-muscle Myosin IIA Rabbit pAb	ER62480	Human, Mouse, Rat	WB,IF-Cell,FC
	non-muscle Myosin IIB / MYH10 Recombinant Rabbit mAb [PSH15-72] - BSA and Azide free	HA751579	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	non-muscle Myosin IIB / MYH10 Recombinant Rabbit mAb [PSH15-72]	HA723775	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Non-Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-43] - BSA and Azide free (Capture)	HA723843	Human, Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Non-Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-45] - Detector	HA723846B	Human	ELISA(Det),ELISA
	Non-Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-45] - BSA and Azide free (Detector)	HA723845	Human	ELISA(Det)
	Non-Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-44] - BSA and Azide free (Capture)	HA723844	Human	ELISA(Cap)
	NPR-C Recombinant Rabbit mAb [JE34-50]	HA722538	Human, Mouse, Rat	WB,FC
	NPRC Rabbit pAb	ER1914-07	Human	WB,IHC-P,IF-Tissue
	NUAK2 Rabbit pAb	ER61680	Human	WB
	N WASP Recombinant Rabbit mAb [JU07-31]	ET1706-39	Human, Mouse	WB,IF-Cell,IHC-P,FC
KD	Occludin Recombinant Rabbit mAb [JJ091-08]	ET1701-76	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IP,FC,IF-Cell,IHC-P,IHC-Fr,IF-Tissue
	Occludin Recombinant Rabbit mAb	IRS183RB	Mouse	mIHC
	Occludin Rabbit pAb	R1510-33	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	OMD Rabbit pAb	ER61731	Human, Rat, Mouse	WB
	OSTCN Rabbit pAb	ER61846	Human, Rat, Mouse	WB
	Osteocalcin Rabbit pAb	ER1919-20	Human, Mouse	WB,IHC-P,IF-Tissue
	Osteoglycin Rabbit pAb	ER65533	Human, Mouse	WB
	Osteopontin Rabbit pAb	ER1802-15	Human, Mouse	IHC-P,ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Osteopontin Rabbit pAb	ER1802-16	Human, Mouse	IF-Cell,IHC-P,ELISA
Cell treatment	Osteopontin Recombinant Rabbit mAb [PSH09-24] - BSA and Azide free	HA751278	Human, Mouse	WB,IF-Cell,FC
	Osteopontin Recombinant Rabbit mAb	IRS200RB	Human	mIHC
Cell treatment	Osteopontin Recombinant Rabbit mAb [PSH09-24]	HA723082	Human, Mouse	WB,IF-Cell,FC
	Osteopontin Recombinant Rabbit mAb [PSH14-11] - BSA and Azide free	HA751506	Mouse, Rat	IHC-P,IHC-Fr,WB,IF-Cell
	Osteopontin Rabbit pAb	HA500464	Human, Mouse	WB,FC
	Osteopontin Recombinant Rabbit mAb [PSH09-23] - BSA and Azide free	HA751277	Human	IHC-Fr,IHC-P,IF-Tissue
	Osteopontin Recombinant Rabbit mAb [PSH14-11]	HA723601	Mouse, Rat	IHC-P,IHC-Fr,WB,IF-Cell
Cell treatment	Osteopontin Rabbit pAb	0806-6	Mouse, Human, Rat	WB
	Osteopontin Rabbit pAb	0806-8	Human, Zebrafish	WB,IHC-P
	Osteopontin Recombinant Rabbit mAb [PSH09-23]	HA723081	Human	IHC-Fr,IHC-P,IF-Tissue
	Osteoprotegerin Mouse mAb [9C1]	EM1701-98	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Osteoprotegerin Rabbit pAb	R1608-4	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Osteoprotegerin Mouse mAb [9C2]	EM1701-99	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	p21 ARC Rabbit pAb	HA500005	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	p21-ARC Rabbit pAb	ER1914-57	Mouse	WB,IHC-P,IF-Tissue
KD	P4HB Recombinant Rabbit mAb [JE54-98]	ET7110-92	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	p50 dynamitin Rabbit pAb	HA500029	Human, Mouse	WB,IF-Cell,FC
	PAK1 + PAK2 + PAK3 Rabbit pAb	ER64516	Human, Mouse, Rat	WB,IHC-P
	PAK3 Recombinant Rabbit mAb [JE57-39]	ET7111-35	Human, Mouse, Rat	WB
	PAK3 Mouse mAb [F11-G8]	EM1708-48	Human, Monkey	WB
	PAK5 Rabbit pAb	ER64513	Human, Mouse, Rat	WB,IHC-P
	PAK5/6 Rabbit pAb	ER64514	Human, Mouse	WB
	Palladin Recombinant Rabbit mAb [JE66-56]	HA721440	Human	WB
	Pan-Actin Recombinant Mouse mAb [A2E1-R] - BSA and Azide free	HA610160	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Pan-Actin Mouse mAb [A2E1]	HA600032	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Pan-Actin Recombinant Mouse mAb [A2E1-R]	HA601287	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	Pan Cadherin Recombinant Rabbit mAb	RE6045	Human, Mouse, Rat	WB,ICC,IHC-P,FC,IP
	Pan-Cadherin Recombinant Rabbit mAb [ST54-01]	ET1609-70	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,IHC-Fr
	Pan-Cadherin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS286MS	Mouse	mIHC
	Pan-Cadherin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS287RT	Human	mIHC
	Pan-Cadherin Recombinant Rabbit mAb	IRS185RB	Mouse	mIHC
	Pan-Cadherin Recombinant Rabbit mAb [ST54-01] - BSA and Azide free	HA750190	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,IHC-Fr
	Pan-Cadherin Recombinant Antibody - Mouse IgG1 (Chimeric)	IRS285MS	Human	mIHC
	Pan-Cadherin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS288RT	Mouse	mIHC
	pan Cytokeratin Recombinant Mouse mAb [PD00-15]	HA601094	Human	IHC-P,IF-Cell,FC,mIHC,IF-Tissue
	pan Cytokeratin Recombinant Mouse mAb [PDH09-10]	HA601138	Human	IHC-P,mIHC,IF-Cell,IF-Tissue,WB,FC
	iFluor™ 647 Conjugated pan Cytokeratin Recombinant Mouse mAb [PD00-15]	HA600105F	Human	IF-Tissue,FC
	pan Cytokeratin Recombinant Mouse mAb	IRS010	Human	mIHC
	iFluor™ 555 Conjugated pan Cytokeratin Recombinant Mouse mAb [PDH09-10]	HA600134F	Human	IF-Cell
	pan Cytokeratin Recombinant Antibody - Rabbit IgG (Chimeric)	IRS186RB	Human	mIHC
	iFluor™ 488 Conjugated pan Cytokeratin Recombinant Mouse mAb [PDH0-10]	HA600111F	Human	IF-Tissue
	pan Cytokeratin Recombinant Rabbit mAb	IRS188RB	Mouse	mIHC
	pan Cytokeratin Recombinant Antibody [PDH09-10] - Rat IgG1 (Chimeric)	HA601519	Human	IHC-P
	pan Cytokeratin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS187RT	Human	mIHC
	iFluor™ 488 Conjugated pan Cytokeratin Recombinant Mouse mAb [PDH09-10]	HA600135F	Human	IF-Cell,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Pan Cytokeratin Rabbit pAb	ER1914-74	Human, Mouse, Rat, Cow	WB,IHC-P,IF-Cell,IF-Tissue,FC
	PARD3A Rabbit pAb	ER64522	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	PARD6B Rabbit pAb	ER64521	Human, Mouse	WB,IHC-P
	Paxillin Rabbit pAb	RT1466	Human, Mouse, Rat	WB,IP,IF
	Paxillin Recombinant Rabbit mAb [SY23-02] - BSA and Azide free	HA750114	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Paxillin Recombinant Rabbit mAb [SY23-02]	ET1607-22	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	PCDH20 Rabbit pAb	ER1914-83	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	PCDH7 Rabbit pAb	ER1914-84	Human, Mouse	WB,IHC-P,IF-Tissue
KD	PDK1 Recombinant Rabbit mAb [JA67-30] - BSA and Azide free	HA750419	Human, Mouse, Rat	WB,IHC-P,IP,FC,IF-Tissue
KD	PDK1 Recombinant Rabbit mAb [JA67-30]	ET1704-66	Human, Mouse, Rat	WB,IHC-P,IP,FC,IF-Tissue
Relative expression	PD-L2 Mouse mAb [A10D8]	HA601331	Human	WB,IF-Cell
	PDLIM5 Rabbit pAb	ER64540	Human	WB
	PDPK1 Recombinant Rabbit mAb [SD204-9]	ET1612-27	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Periostin Recombinant Rabbit mAb [JA63-40]	ET1704-65	Human	WB
	Peripherin Recombinant Rabbit mAb [PSH10-70] - BSA and Azide free	HA751372	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Cell
	Peripherin Recombinant Antibody [PSH10-70] - Rat IgG1 (Chimeric)	HA601405	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB,IF-Cell
	Peripherin Recombinant Rabbit mAb [PSH10-70]	HA723234	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Cell
	Peripherin Recombinant Antibody [PSH10-70] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610245	Human, Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB,IF-Cell
	PERP Recombinant Rabbit mAb [JU04-51]	ET7106-59	Human, Mouse	WB,IHC-P
	PHACTR4 Rabbit pAb	ER64574	Human	WB
	Phospho1 Rabbit pAb	ER1915-07	Human, Mouse	WB,IHC-P,IF-Tissue
Cell treatment	Phospho-Beta Catenin (S33 + S37) Recombinant Rabbit mAb [JE59-59]	HA721580	Human, Rat, Mouse	WB
Cell treatment	Phospho-Beta Catenin (S552) Recombinant Rabbit mAb [PSH08-72]	HA723024	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Beta Catenin (S552) Recombinant Rabbit mAb [PSH08-72] - BSA and Azide free	HA751237	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Cell treatment	Phospho-Beta Catenin (T41/S45) Recombinant Rabbit mAb [JE54-02]	HA722316	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-Cofilin (S3) Recombinant Rabbit mAb [PSH07-18] - BSA and Azide free	HA751139	Human, Rat, Monkey	WB
Cell treatment	Phospho-Cofilin (S3) Recombinant Rabbit mAb [PSH07-18]	HA722806	Human, Rat, Monkey	WB
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4]	ET1603-22	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4] - BSA and Azide free	HA750061	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01] - BSA and Azide free	HA750199	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01]	ET1610-13	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-Ezrin (T567)/Radixin (T564)/Moesin (T558) Recombinant Rabbit mAb [PSH05-73]	HA722440	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-Ezrin (T567)/Radixin (T564)/Moesin (T558) Recombinant Rabbit mAb [PSH05-73] - BSA and Azide free	HA751023	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Phospho-FAK (S722) Mouse mAb [4G1]	RT1216	Human, Mouse, Rat	WB,IP,IF,IHC-P
Cell treatment	Phospho-FAK (Y397) Recombinant Rabbit mAb [SC54-07]	ET1610-34	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-FAK (Y397) Recombinant Rabbit mAb [SC54-07] - BSA and Azide free	HA750212	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-Hsp27 (S78) Recombinant Rabbit mAb [JJ08-70]	ET1701-19	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Phospho-Hsp27 (S82) Recombinant Rabbit mAb [SN06-73]	ET1611-16	Human, Mouse, Rat	WB,IP,IF-Cell,FC
Cell treatment	Phospho-Lamin A + Lamin C (S22) Recombinant Rabbit mAb [PSH11-47]	HA723339	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-Lamin A + Lamin C (S22) Recombinant Rabbit mAb [PSH11-48] - BSA and Azide free	HA751410	Human, Mouse, Rat	WB
Cell treatment	Phospho-Lamin A + Lamin C (S22) Recombinant Rabbit mAb [PSH11-47] - BSA and Azide free	HA751409	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Lamin A + Lamin C (S22) Recombinant Rabbit mAb [PSH11-48]	HA723340	Human, Mouse, Rat	WB
Cell treatment	Phospho-MAP2 (S136) Recombinant Rabbit mAb [JE44-25]	HA722351	Mouse, Rat	WB,IHC-P,IHC-Fr,IF-Tissue
	Phospho-MYL9 (T18 + S19) Recombinant Rabbit mAb [PSH12-71] - BSA and Azide free	HA751447	Mouse, Rat	WB
	Phospho-MYL9 (T18 + S19) Recombinant Rabbit mAb [PSH12-71]	HA723482	Mouse, Rat	WB
	Phospho-Myosin light chain (S20) Rabbit pAb	ER1913-46	Human, Mouse	WB,IHC-P,IF-Tissue
Cell treatment	Phospho-MYPT1 (T696) Recombinant Rabbit mAb [PSH08-96] - BSA and Azide free	HA751257	Human, Mouse, Rat	WB
Cell treatment	Phospho-MYPT1 (T696) Recombinant Rabbit mAb [PSH08-96]	HA723049	Human, Mouse, Rat	WB
Cell treatment	Phospho-MYPT1 (T853) Recombinant Rabbit mAb [PSH08-71] - BSA and Azide free	HA751236	Human, Mouse, Rat	WB
Cell treatment	Phospho-MYPT1 (T853) Recombinant Rabbit mAb [PSH08-71]	HA723023	Human, Mouse, Rat	WB
	Phospho-PAK1 (S144)+PAK2 (S141)+PAK3 (S139) Recombinant Rabbit mAb [SU03-70]	ET1608-58	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Cell treatment	Phospho-PAK4+PAK5+PAK6 (S474 + S560 + S602) Recombinant Rabbit mAb [JE57-42]	HA722146	Human	WB
Cell treatment	Phospho-Paxillin (Y118) Recombinant Rabbit mAb [JE09-10]	HA722526	Human	WB
Cell treatment	Phospho-Rac1+Cdc42 (S71) Recombinant Rabbit mAb [JJ086-9]	ET1701-40	Human	WB,IF-Cell,IHC-P
Cell treatment	Phospho-Tau (S198) Recombinant Rabbit mAb [JE66-79]	HA721668	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S199) Recombinant Rabbit mAb [JE46-95]	HA721800	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S199) Recombinant Rabbit mAb [JE66-80]	HA721667	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S202) Recombinant Rabbit mAb [JE66-81]	HA721666	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S214) Recombinant Rabbit mAb [JE66-83]	HA721797	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S324) Recombinant Rabbit mAb [JE30-73]	HA721798	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Tau (S356) Recombinant Rabbit mAb [JE46-86]	HA721799	Human, Mouse, Rat	WB
Cell treatment	Phospho-Tau (S396) Recombinant Rabbit mAb [SN62-09] - BSA and Azide free	HA750268	Mouse, Rat, Human, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-Tau (S396) Recombinant Rabbit mAb [SN62-09]	ET1611-68	Mouse, Rat, Human, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-Tau (S404) Recombinant Rabbit mAb [JE47-88]	HA721801	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-Tau (S404) Recombinant Rabbit mAb [JE66-82]	HA721665	Human, Mouse, Rat	WB,IHC
Cell treatment	Phospho-Tau (T181) Recombinant Rabbit mAb [PSH05-42] - BSA and Azide free	HA751001	Human, Mouse, Cynomolgus monkey	WB,IHC-P
Cell treatment	Phospho-Tau (T181) Recombinant Rabbit mAb [PSH05-42]	HA722271	Human, Mouse, Cynomolgus monkey	WB,IHC-P
	Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-42] - BSA and Azide free (Capture)	HA723842	Human	ELISA(Cap)
	Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-75] - BSA and Azide free (Capture)	HA723877	Human	ELISA(Cap)
Cell treatment	Phospho-Tau (T217) Recombinant Rabbit mAb [PSH09-33]	HA723091	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr
	Phospho-Tau (T217) Recombinant Rabbit mAb [PSH16-41] - BSA and Azide free (Capture)	HA723841	Human, Mouse	ELISA(Cap)
Cell treatment	Phospho-Tau (T217) Recombinant Rabbit mAb [PSH09-33] - BSA and Azide free	HA751287	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr
	Phospho-Tau (T217) Recombinant Rabbit mAb [PSH17-41] - BSA and Azide free (Capture)	HA723922	Human, Mouse	ELISA(Cap)
	Phospho-Tau (T231) Recombinant Rabbit mAb [PSH01-04]	HA721570	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Cell treatment	Phospho-Tau (T231) Recombinant Rabbit mAb [SC58-08]	ET1610-31	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IP,IF-Tissue,IHC-Fr
	Phospho-Tau (T231) Recombinant Rabbit mAb [PSH01-05]	HA721569	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-Tau (T231) Recombinant Rabbit mAb [SC58-08] - BSA and Azide free	HA750210	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IP,IF-Tissue,IHC-Fr
	Phospho-Tau (T231) Recombinant Rabbit mAb [PSH01-05] - BSA and Azide free	HA750679	Human, Mouse, Rat	WB,IHC-P
	Phospho-Tau (T231) Recombinant Rabbit mAb [PSH01-04] - BSA and Azide free	HA750680	Human, Mouse, Rat	WB,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-VASP (S156) Recombinant Rabbit mAb [JE45-91]	HA722145	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-Vimentin (S39) Recombinant Rabbit mAb [JE43-26]	HA721442	Human, Mouse, Rat	WB,IF-Cell
	Photomedin-1 Rabbit pAb	ER64584	Human, Mouse	WB,IF-Cell,IHC
	PINCH Mouse mAb [B9-G7]	EM1709-96	Human	WB,IF-Cell,FC
	PIP5K2 gamma Rabbit pAb	HA500139	Human, Mouse, Rat	WB
	PIX beta Rabbit pAb	ER64612	Human, Mouse, Rat	WB,IHC-P
	Plakophilin 2 Rabbit pAb	ER64628	Human, Rat	WB
	Plakophilin 3 Rabbit pAb	HA500034	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Plakophilin 4 Rabbit pAb	ER1915-12	Human	WB
	PLASTIN3 Rabbit pAb	ER1915-14	Mouse	WB,IHC-P,IF-Tissue
KD	Plastin L Recombinant Rabbit mAb [JE48-95]	ET7109-87	Human	WB,IHC-P,FC
KD	Plectin Recombinant Rabbit mAb [SY29-04]	ET1607-80	Human, Mouse, Rat	WB,IHC-P,IP,IF-Tissue,IF-Cell
	PLEKHA7 Rabbit pAb	ER1915-17	Human	WB,IHC-P,IF-Tissue
	PLOD2 Rabbit pAb	ER62020	Human, Rat, Mouse	WB
	PODN Rabbit pAb	ER1915-21	Mouse	WB
	PODXL Rabbit pAb	ER1915-24	Mouse	WB
	PODXL Mouse mAb [2-A4]	M1403-5	Human	IF-Cell,IHC-P
	PODXL Mouse mAb [D0-E8]	M1403-4	Human	WB
	POPD1 Rabbit pAb	ER60530	Human, Mouse, Rat	WB
	PPM1E Recombinant Rabbit mAb [JE61-81]	HA720063	Human	WB
	PRDC Rabbit pAb	HA500225	Human	WB,FC
	PROF2 Rabbit pAb	ER62405	Human, Mouse, Rat	WB
	Profilin 1 Recombinant Rabbit mAb [JA30-13]	ET1704-42	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Prostate Specific Antigen Recombinant Rabbit mAb [JE19-84]	ET7110-62	Human	IHC-P
	Prostate Specific Antigen Recombinant Rabbit mAb [SD07-45]	ET1612-19	Human	WB,IF-Cell,IF-Tissue,IHC-P
	Prostate-Specific Antigen (PSA) Rabbit pAb	ER1706-31	Human	IHC-P
	Protocadherin-11 Rabbit pAb	ER64678	Human	IHC-P,IF-Cell
	Protocadherin-16 Rabbit pAb	0905-2	Human	WB
	PRS23 Rabbit pAb	ER61169	Human, Rat, Mouse	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PRTN3 Rabbit pAb	ER1706-25	Human, Mouse	WB,IHC-P
	PTN Rabbit pAb	ER65608	Human, Mouse, Rat	WB,IHC-P
	PTN Recombinant Rabbit mAb [JE47-40]	HA722055	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Pyrin Recombinant Rabbit mAb [PSH06-12] - BSA and Azide free	HA751050	Human	WB,FC
Cell treatment	Pyrin Recombinant Rabbit mAb [PSH06-12]	HA722509	Human	WB,FC
	Rab13 Mouse mAb	EM1707-79	Human	IF-Cell,IHC-P,FC
	RAB21 Rabbit pAb	ER61099	Human, Rat, Mouse	WB
	Radixin Recombinant Rabbit mAb [SC06-22]	ET1610-41	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	RANK Rabbit pAb	ER1915-70	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	RANKL Rabbit pAb	HA500369	Human, Mouse, Rat	WB,IF-Cell
	Rat CRP Recombinant Rabbit mAb [PSH12-80] - BSA and Azide free (Detector)	HA725050	Rat	ELISA(Det)
	Rat CRP Recombinant Rabbit mAb [PSH12-79] - BSA and Azide free (Capture)	HA725049	Rat	ELISA(Cap)
	Biotin Conjugated Rat CRP Recombinant Rabbit mAb [PSH12-80] - Detector	HA725051B	Rat	ELISA(Det),ELISA
	HRP Conjugated Rat CRP Recombinant Rabbit mAb [PSH12-80] - Detector	HA725051H	Rat	ELISA(Det),ELISA
	Rat CRP/C-Reactive Protein enzyme-linked immunoassay kit	ER0009	Rat	ELISA
	Rat CRP/C-Reactive Protein enzyme-linked immunoassay kit(one step)	ERY0009	Rat	ELISA
	Rat IFN- γ enzyme-linked immunoassay kit	ER0006	Rat	ELISA
	Rat IL-10 Recombinant Rabbit mAb [PSH02-21] - BSA and Azide free (Capture)	HA721786	Rat	ELISA(Cap)
	Rat IL-10 Recombinant Rabbit mAb [PSH02-22] - BSA and Azide free (Detector)	HA721787	Rat	ELISA(Det)
	Biotin Conjugated Rat IL-10 Recombinant Rabbit mAb [PSH02-22]	HA722637B	Rat	ELISA(Det)
	Rat IL-10 enzyme-linked immunoassay kit	ER0005	Rat	ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-59] - BSA and Azide free (Capture)	HA721836	Rat	ELISA(Cap)
	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - Detector	HA721611	Rat	ELISA(Det),ELISA
	Rat IL-1 beta Recombinant Rabbit mAb [PS01-86] - BSA and Azide free (Capture)	HA721609	Rat	ELISA(Cap)
	Biotin Conjugated Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60]	HA722634B	Rat	ELISA(Det)
	Rat IL-1 beta Recombinant Rabbit mAb [PS01-87] - BSA and Azide free (Detector)	HA721610	Rat	ELISA(Det)
	Rat IL-1 beta Recombinant Rabbit mAb [PSH02-60] - BSA and Azide free (Detector)	HA721837	Rat	ELISA(Det)
	Rat IL-1 β enzyme-linked immunoassay kit	ER0001	Rat	ELISA
	Rat IL-2 Recombinant Rabbit mAb [PSH05-72] - BSA and Azide free (Detector)	HA722424	Rat	ELISA(Det)
	Rat IL-2 Recombinant Rabbit mAb [PS01-88] - BSA and Azide free (Capture)	HA721612	Rat	ELISA(Cap)
	Biotin Conjugated Rat IL-2 Recombinant Rabbit mAb [PS01-89] - Detector	HA721614	Rat	ELISA(Det)
	Rat IL-2 Recombinant Rabbit mAb [PS01-89] - BSA and Azide free (Detector)	HA721613	Rat	ELISA(Det)
	Rat IL-2 Recombinant Rabbit mAb [PSH05-71] - BSA and Azide free (Capture)	HA722423	Rat	ELISA(Cap)
	Biotin Conjugated Rat IL-2 Recombinant Rabbit mAb [PSH05-72]	HA722590B	Rat	ELISA(Det)
	Rat IL-2 enzyme-linked immunoassay kit	ER0002	Rat	ELISA
	Rat IL-4 enzyme-linked immunoassay kit	ER0003	Rat	ELISA
	Rat IL-6 Recombinant Rabbit mAb [PS01-90] - BSA and Azide free (Capture)	HA721615	Rat	ELISA(Cap)
	Rat IL-6 Recombinant Rabbit mAb [PS01-94] - BSA and Azide free (Detector)	HA721616	Rat	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Rat IL-6 Recombinant Rabbit mAb [PS01-94] - BSA and Azide free (Detector)	HA721617	Rat	ELISA(Det)
	Rat IL-6 enzyme-linked immunoassay kit	ER0004	Rat	ELISA
	Rat TNF alpha Recombinant Rabbit mAb [PSH04-71] - BSA and Azide free (Detector)	HA722171	Rat	ELISA(Det)
	Rat TNF alpha Recombinant Rabbit mAb [PSH04-70] - BSA and Azide free (Capture)	HA722170	Rat	ELISA(Cap)
	Biotin Conjugated Rat TNF alpha Recombinant Rabbit mAb [PSH04-71]	HA722606B	Rat	ELISA(Det)
	Rat TNF- α enzyme-linked immunoassay kit	ER0007	Rat	ELISA
	R-cadherin Rabbit pAb	ER64773	Human, Mouse	IHC-P,IF-Cell,WB
Relative expression	Reelin Recombinant Rabbit mAb [JA63-13]	ET1704-98	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr
	RhoA Recombinant Rabbit mAb [SN0612]	ET1611-10	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	RhoA Recombinant Rabbit mAb [SN0612] - BSA and Azide free	HA750242	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Rho A+B+C Recombinant Rabbit mAb [SR38-00]	ET1602-23	Human, Mouse, Rat	WB,IF-Cell,FC
	Rho B Rabbit pAb	ER64810	Human, Mouse, Rat	WB
	Rho D Rabbit pAb	ER64811	Human, Mouse	WB,IHC-P
KD	RhoGDI Rabbit pAb	HA500067	Human, Mouse, Rat	WB,IHC-P
	RND1 Rabbit pAb	ER61148	Human, Mouse	WB
	ROBO4 Rabbit pAb	ER61895	Human, Rat, Mouse	WB
	ROCK1 Recombinant Rabbit mAb [ST05-19] - BSA and Azide free	HA750181	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
KD	ROCK1 Rabbit pAb	ER40122	Human, Mouse, Rat	WB,IHC-P
	ROCK1 Recombinant Rabbit mAb [ST05-19]	ET1609-53	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,IP
	ROCK2 Recombinant Rabbit mAb [JE31-36]	HA722176	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	ROCK2 Rabbit pAb	ER1706-48	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	ROR2 Recombinant Rabbit mAb [JE37-72]	HA721417	Human	WB,IHC-P
	RSPH4A Rabbit pAb	ER1916-15	Human, Mouse	WB
KD	S100A10 Recombinant Rabbit mAb [JF0987]	ET1702-38	Human, Mouse, Rat	WB,IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	S100A4 Recombinant Rabbit mAb [SD200-08] - BSA and Azide free	HA750282	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Cell treatment	S100A4 Recombinant Rabbit mAb [SD200-08]	ET1612-13	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	S100 alpha Recombinant Rabbit mAb [JM51-31]	ET1703-72	Human, Mouse, Rat	WB,IP,IHC-P,FC,IF-Cell
	S100 alpha Recombinant Rabbit mAb [JM51-31] - BSA and Azide free	HA750389	Human, Mouse, Rat	WB,IP,IHC-P,FC,IF-Cell
	Biotin Conjugated S100B Recombinant Rabbit mAb [PSH15-25] - Detector	HA723730B	Human, Mouse, Rat	ELISA(Det),ELISA
	S100B Recombinant Rabbit mAb [SC57-02] - BSA and Azide free (Capture)	HA723728	Human, Mouse, Rat	ELISA(Cap)
	S100B Recombinant Rabbit mAb [PSH15-25] - BSA and Azide free (Detector)	HA723729	Human, Mouse, Rat	ELISA(Det)
	S100 beta Recombinant Antibody [SC57-02] - Mouse IgG1 (Chimeric)	HA601551	Mouse, Rat	IHC-Fr,IHC-P
	S100 beta Rabbit pAb	RT1547	Human, Mouse, Rat	WB,IP,IF
	S100 beta Recombinant Antibody [SC57-02] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610366	Mouse, Rat	IHC-Fr,IHC-P
Relative expression	S100 beta Recombinant Mouse mAb [PD00-11]	HA601098	Human, Mouse, Rat	IHC-P,WB,IF-Tissue
	S100 beta Recombinant Antibody [SC57-02] - Guinea pig IgG2 (Chimeric)	HA601550	Mouse	IHC-Fr
	S100 beta Recombinant Antibody [SC57-02] - Rat IgG1 (Chimeric)	HA601552	Mouse, Rat	IHC-Fr,IHC-P
Relative expression	S100 beta Recombinant Rabbit mAb [PSH17-13]	HA723900	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Tissue
	S100 beta Recombinant Rabbit mAb [SC57-02] - BSA and Azide free	HA750209	Human, Mouse, Rat, Zebrafish, Goat, Cynomolgus monkey, Green monkey, Pig, Rabbit	WB,IF-Tissue,IHC-P,IHC-Fr,IF-Cell,FC,IP
	S100 beta Recombinant Antibody [SC57-02] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610367	Mouse, Rat	IHC-Fr,IHC-P
	S100 beta Recombinant Antibody [SC57-02] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610365	Mouse	IHC-Fr
	S100 beta Recombinant Rabbit mAb [PSH17-13] - BSA and Azide free	HA751617	Human, Mouse, Rat	IHC-Fr,IHC-P,WB,IF-Tissue
	iFluor™ 594 Conjugated S100 beta Recombinant Rabbit mAb [SC57-02]	HA720125F	Human	IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	S100 beta Recombinant Rabbit mAb [SC57-02]	ET1610-3	Human, Mouse, Rat, Zebrafish, Goat, Cynomolgus monkey, Green monkey, Pig, Rabbit	WB,IF-Tissue,IHC-P,IHC-Fr,IF-Cell,FC,IP
	Sarcoglycan beta Rabbit pAb	ER64909	Human, Mouse	WB,IHC-P
	Sarcomeric Alpha Actinin Recombinant Rabbit mAb [JE45-22]	ET7109-42	Human, Mouse, Rat	WB,IP,IHC-P
	SCIN Recombinant Rabbit mAb [JE62-89]	HA722465	Human, Mouse, Rat	WB
	SFRP3 Rabbit pAb	ER61167	Human, Mouse	WB
	SgK269 Rabbit pAb	ER64948	Human	WB
	SH2D4A Rabbit pAb	ER64954	Human	IHC-P
	SH3PX1 Rabbit pAb	HA500109	Human, Mouse, Rat	WB,IF-Cell,FC
	SHARP2 Recombinant Rabbit mAb [PSH01-91]	HA721749	Human	WB
	SHARP2 Recombinant Rabbit mAb [PSH01-91] - BSA and Azide free	HA750765	Human	WB
	SHARP2 Recombinant Rabbit mAb [PSH01-90] - BSA and Azide free	HA750764	Human	WB,IHC-P,IF-Cell
	SHARP2 Recombinant Rabbit mAb [PSH01-90]	HA721748	Human	WB,IHC-P,IF-Cell
	SHIP-2 Rabbit pAb	ER65496	Human, Mouse, Rat	WB,IHC-P
	SIAL Rabbit pAb	ER62043	Human, Rat, Mouse	WB
	FITC Conjugated Siglec 8 Recombinant Mouse mAb [PSH14-96]	HA601517F	Human	FC
	SIRP alpha Mouse mAb [A3E3]	EM1902-37	Human	WB,IF-Cell,IHC-P,FC
KD	SIRP alpha Recombinant Rabbit mAb [JE36-20]	HA721316	Human, Mouse, Rat	WB,IHC-P,FC
	SIRP Alpha Rabbit pAb	ER1916-61	Human, Mouse, Rat	WB,FC
	SIRT2 Recombinant Rabbit mAb [SN70-04]	ET1611-72	Human, Cynomolgus monkey	WB,IP,IHC-P
KD	SIRT2 Recombinant Rabbit mAb [JJ08-83]	ET1701-43	Human	WB
	SIRT2 Recombinant Antibody [SN70-04] - Mouse IgG1 (Chimeric)	HA601505	Mouse, Rat	IHC-Fr,IHC-P,WB
	SIRT2 Recombinant Antibody [SN70-04] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610335	Mouse, Rat	IHC-Fr,IHC-P,WB
	SIRT2 Recombinant Rabbit mAb [PSH17-48]	HA723923	Human, Mouse, Rat, Cynomolgus monkey	IHC-Fr,IHC-P,WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	SIRT2 Recombinant Antibody [SN70-04] - Rat IgG1 (Chimeric)	HA601506	Mouse, Rat	IHC-Fr,IHC-P,WB
	SIRT2 Recombinant Antibody [SN70-04] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610336	Mouse, Rat	IHC-Fr,IHC-P,WB
	SLINGSHOT-1L Mouse mAb [G9-C11]	EM1707-22	Human	IHC-P,FC
Relative expression	Slow Skeletal Myosin Heavy chain Recombinant Rabbit mAb [PSH07-54]	HA722735	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
	Slow Skeletal Myosin Heavy chain Recombinant Rabbit mAb	IRS201RB	Mouse	mIHC
Relative expression	Slow Skeletal Myosin Heavy chain Recombinant Rabbit mAb [PSH06-58] - BSA and Azide free	HA751086	Human, Mouse, Rat	WB,IHC-P
Relative expression	Slow Skeletal Myosin Heavy chain Recombinant Rabbit mAb [PSH07-54] - BSA and Azide free	HA751092	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression	Slow Skeletal Myosin Heavy chain Recombinant Rabbit mAb [PSH06-58]	HA722721	Human, Mouse, Rat	WB,IHC-P
	Smoothelin Rabbit pAb	ER2001-40	Human	WB,IHC-P
	smooth muscle Myosin heavy chain 11 Recombinant Rabbit mAb [JE60-21]	HA722863	Human, Mouse, Rat	WB,IHC-P
	smooth muscle Myosin heavy chain 11 Recombinant Rabbit mAb [JA03-35]	ET1704-61	Human, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
	smooth muscle Myosin heavy chain 11 Recombinant Rabbit mAb	IRS176RB	Human	mIHC
	smooth muscle Myosin heavy chain 11 Recombinant Rabbit mAb [JA03-35] - BSA and Azide free	HA750417	Human, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC
	SP-A Rabbit pAb	ER65482	Human, Mouse, Rat	WB,IHC-P
	SPARC Recombinant Rabbit mAb [PD00-57]	HA721199	Human	IHC-P,IF-Tissue
	SPARC Rabbit pAb	ER1916-99	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
Relative expression	SPARC Recombinant Rabbit mAb [PSH06-23]	HA722655	Human	WB
Relative expression	SPARC Recombinant Rabbit mAb [PSH06-23] - BSA and Azide free	HA751061	Human	WB
	Spastin Rabbit pAb	ER1917-02	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Spectrin beta V Rabbit pAb	ER65008	Human	IHC-P
	SPOCK1 Rabbit pAb	ER1917-13	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	SPOCK2 Rabbit pAb	ER1917-14	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
	SPTBN1 Rabbit pAb	HA500014	Human, Mouse, Rat	WB,FC,IF-Cell
	SPTN4 Rabbit pAb	ER61257	Human	IHC-P
	SSH3 Rabbit pAb	ER65033	Human, Mouse, Rat	WB,IHC-P
	SSX2IP Recombinant Rabbit mAb [JE65-50]	HA721065	Human	WB
	STAB2 Rabbit pAb	ER61939	Human, Rat, Mouse	IHC-P
Relative expression	Stathmin Recombinant Rabbit mAb [PD00-54]	HA721179	Human, Mouse, Rat	IHC-P,WB,IF-Cell,IF-Tissue,FC
	Stathmin 1 Mouse mAb [11E1]	EM1801-19	Human	WB,IF-Cell,IHC-P,FC
KD	Stathmin 1 Recombinant Rabbit mAb [SS0453]	ET1609-20	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
KD	Staufen Recombinant Rabbit mAb [JU30-85]	ET1706-41	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC
	SUN2 Recombinant Rabbit mAb [JG39-52]	ET7108-92	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell,IF-Tissue
	Supervillin Rabbit pAb	ER1917-34	Mouse	IHC-P,IF-Tissue
	Syndecan 1 Recombinant Rabbit mAb [JM11-21]	ET1703-42	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Syndecan 1 Recombinant mAb	IRS105RB	Human	mIHC
	Syndecan 1 Recombinant Rabbit mAb [JM11-21] - BSA and Azide free	HA750381	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
	Syndecan-2 Rabbit pAb	ER65062	Human	WB
	Syndecan 3 Recombinant Rabbit mAb [JE64-19]	HA722471	Human, Mouse, Rat	WB
	Syndecan 4 Recombinant Rabbit mAb [PSH0-72] - BSA and Azide free	HA750651	Human	WB
	Syndecan 4 Recombinant Rabbit mAb [PSH0-72]	HA721473	Human	WB
	Syndecan-4 Rabbit pAb	ER65063	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	SYNE2 Rabbit pAb	ER60902	Human, Mouse	IHC-P
	Syntrophin alpha 1 Recombinant Rabbit mAb [JG35-62]	ET7108-15	Human, Mouse, Rat	WB,IHC-P,FC,IF-Tissue
	Syntrophin alpha 1 Rabbit pAb	ER1901-40	Mouse, Rat, Human	WB
	SYT5 Rabbit pAb	ER61286	Human, Rat, Mouse	WB
	TAG1 Rabbit pAb	HA500013	Human, Mouse, Rat	WB,IF-Cell
	TAGLN Recombinant Rabbit mAb [JE60-72]	HA721928	Human, Mouse, Rat	WB,IHC-P,FC
	Talin 1 Rabbit pAb	ER65629	Human, Mouse, Rat	WB,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	Talin 1 Recombinant Rabbit mAb [PSH09-95]	HA723164	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	Talin 1 Recombinant Rabbit mAb [PSH09-95] - BSA and Azide free	HA751326	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Tartrate-resistant acid phosphatase type 5 Recombinant Rabbit mAb [JE35-95]	HA721897	Human	WB,IHC-P
	Tau Recombinant Rabbit mAb [PSH17-27] - BSA and Azide free	HA723914	Human	ELISA
	Biotin Conjugated Tau Recombinant Rabbit mAb [PSH16-49] - Detector	HA723852B	Human	ELISA(Det),ELISA
	Biotin Conjugated Tau Recombinant Rabbit mAb [PSH17-11] - Detector	HA723898B	Human	ELISA(Det),ELISA
	Tau Recombinant Rabbit mAb [PSH17-26] - BSA and Azide free	HA723913	Human	ELISA
Relative expression	Tau Recombinant Rabbit mAb [PSH05-50]	HA722288	Mouse, Rat	WB,IHC-P,IF-Tissue,IP
KD	Tau Recombinant Rabbit mAb [SZ03-03]	ET1603-2	Human, Mouse, Rat	WB,IHC-P,IP,IHC-Fr
	Tau Recombinant Rabbit mAb [PSH16-49] - BSA and Azide free (Detector)	HA723851	Human	ELISA(Det)
	Tau Recombinant Rabbit mAb [SD205-09]	ET1612-44	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
	Tau Recombinant Rabbit mAb [PSH16-48] - BSA and Azide free (Capture)	HA723850	Human, Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Tau Recombinant Rabbit mAb [PSH10-42] - Detector	HA723209B	Human, Mouse, Rat	ELISA(Det),ELISA
	Tau Recombinant Rabbit mAb [PSH10-42] - BSA and Azide free (Detector)	HA723209	Human, Mouse, Rat	ELISA(Det)
Relative expression	Tau Recombinant Rabbit mAb [PSH05-50] - BSA and Azide free	HA751009	Mouse, Rat	WB,IHC-P,IF-Tissue,IP
	Tau Recombinant Rabbit mAb [PSH17-11] - BSA and Azide free (Detector)	HA723898	Human	ELISA(Det)
	Tau Rabbit pAb	R1401-7	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	TBCB Rabbit pAb	ER61346	Human, Mouse	WB
	TBCE Rabbit pAb	ER1917-52	Mouse, Rat	WB,IHC-P,IF-Tissue
	TCP1 alpha / CCTA Recombinant Rabbit mAb [JE48-48]	ET7109-79	Human, Mouse	WB,IHC-P,IF-Cell,FC
	TCP-1 epsilon Rabbit pAb	ER65121	Human, Mouse, Rat	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Tenascin C Recombinant Rabbit mAb [SU36-01]	ET1608-50	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	Tenascin-R Rabbit pAb	ER65764	Human, Mouse, Rat	IHC-P
	Tenascin-X Rabbit pAb	ER65136	Human	IHC-P
	Tensin-1 Rabbit pAb	ER65137	Human, Mouse	WB
	Tensin3 Rabbit pAb	ER65138	Human	IHC-P,IF-Cell
	TESK2 Rabbit pAb	ER65140	Human, Mouse, Rat	WB
	Tetranectin Recombinant Rabbit mAb [JG37-16]	ET7108-45	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	TFF3 Rabbit pAb	ER1917-60	Human, Mouse, Rat	IHC-P,IF-Tissue
Relative expression	TFF3/Trefoil Factor 3 Recombinant Rabbit mAb [JE32-01]	HA721541	Human	WB,IHC-P,FC
	Biotin Conjugated TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - Detector	HA723727B	Human, Mouse, Rat	ELISA(Det),ELISA
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17] - BSA and Azide free	HA750485	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	TGF beta 1 Recombinant Rabbit mAb [PSH15-21] - BSA and Azide free (Capture)	HA723725	Human, Mouse, Rat	ELISA(Cap)
	TGF beta 1 Recombinant Rabbit mAb [PSH15-22] - BSA and Azide free (Detector)	HA723726	Human, Mouse, Rat	ELISA(Det)
	TGF beta 1 Recombinant Rabbit mAb [PSH18-59]	HA724012	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC
	TGF beta 1 Mouse mAb [B2-F6]	EM1709-93	Human	WB,IHC-P,FC
	TGF beta 1 Rabbit pAb	HA500496	Human	WB
Relative expression	TGF beta 1 Recombinant Rabbit mAb [PD00-17]	HA721143	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	TGF beta 1 Recombinant Rabbit mAb	IRS205RB	Human	mIHC
	TGFBI Recombinant Rabbit mAb [JM24-53]	ET1705-29	Human, Mouse, Rat	WB,IHC-P
	Thrombin Activation peptide fragment 1 Rabbit pAb	ER1917-70	Mouse, Rat	WB
	Thrombin light chain Rabbit pAb	ER1917-71	Mouse	ELISA
	Thrombospondin 4 Rabbit pAb	ER65162	Human	WB,IHC-P,IF-Cell
	TIAF1 Rabbit pAb	ER60340	Human	WB
	TIMP-1 Mouse mAb [4C2]	EM2001-03	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	TIMP1 Recombinant Rabbit mAb [JE59-72]	HA721741	Human	WB,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	TIMP2 Recombinant Rabbit mAb [PSH11-62]	HA723358	Human, Monkey	WB
Relative expression	TIMP2 Recombinant Rabbit mAb [PSH11-62] - BSA and Azide free	HA751418	Human, Monkey	WB
Relative expression	TIMP3 Recombinant Rabbit mAb [PSH04-41] - BSA and Azide free	HA750950	Human, Rat	WB,IF-Cell
Relative expression	TIMP3 Recombinant Rabbit mAb [PSH04-41]	HA722126	Human, Rat	WB,IF-Cell
	TIMP4 Rabbit pAb	ER65170	Human, Mouse	WB
	TLN2 Rabbit pAb	ER61892	Human, Mouse	IHC-P
	TMPO Rabbit pAb	ER65745	Human, Mouse, Rat	WB,IHC-P
Cell treatment	TNF alpha Rabbit pAb	R1203-1	Mouse	WB,IF-Cell,ELISA
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free	HA750897	Human, Mouse	WB
	TNF alpha Rabbit pAb	ER1919-22	Human, Mouse	WB
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07]	HA722022	Human, Mouse	WB
	TNF alpha Rabbit pAb	ER65189	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	TNFSF11 Mouse mAb [8A7B9]	EM1710-42	Human, Monkey	WB,IHC-P,IF-Cell,FC
	TNNC1 Rabbit pAb	HA500407	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	TNNC1 Rabbit pAb	ER1902-25	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	TNNT1 Rabbit pAb	ER61345	Human, Mouse, Rat	WB
	TNNT3 Rabbit pAb	ER61344	Human	WB
	TNS4 Recombinant Rabbit mAb [PD00-95]	HA721208	Human	IHC-P
	TNS4 Recombinant Rabbit mAb [PD00-95] - BSA and Azide free	HA750540	Human	IHC-P
	TPPP Rabbit pAb	ER60161	Human, Mouse	WB
	TPPP Recombinant Rabbit mAb [JE47-35]	ET7109-56	Human, Mouse, Rat	IHC-P,WB
	TRAK1 Rabbit pAb	ER61308	Human, Mouse	WB
	Tropomodulin 2 Recombinant Rabbit mAb [JE53-06]	ET7110-27	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Tropomyosin 1 Recombinant Rabbit mAb [JG38-66]	ET7108-74	Human, Mouse, Rat	WB,IHC-P,FC
KD	Tropomyosin 2 Recombinant Rabbit mAb [PSH02-63] - BSA and Azide free	HA750820	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell
KD	Tropomyosin 2 Recombinant Rabbit mAb [PSH02-63]	HA721841	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Troponin I-C Rabbit pAb	ER65220	Mouse, Rat	WB,IHC-P,IF-Cell
Cell treatment	Troponin I fast skeletal muscle Recombinant Rabbit mAb [JE35-67]	HA722542	Human, Mouse, Rat	WB,IHC
	TSG-6 Rabbit pAb	ER65230	Human, Mouse	WB,IHC-P,IF-Cell
	TSLC1 Rabbit pAb	ER1918-21	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	TSPAN5 Rabbit pAb	ER1918-25	Mouse, Rat	WB,IHC-P,IF-Tissue
	TTN Rabbit pAb	ER65816	Human, Mouse, Rat	IHC-P
	TUBA8 Mouse mAb [D7-D6]	EM1709-38	Human, Rat	WB,IF-Cell,IHC-P,FC
	TUBE1 Mouse mAb [B7-F3]	EM1706-94	Human	WB,IHC-P,FC
	Tubulin alpha 1/3/4 Rabbit pAb	ER65239	Human, Mouse, Rat	WB
	TWF1 Mouse mAb [A2-C11]	EM1710-07	Human	WB,IF-Cell,IHC-P,FC
Cell treatment	uPA Receptor Mouse mAb [A6E7]	HA601040	Human	WB
Cell treatment	uPA Receptor Recombinant Mouse mAb [A6E8-R] - BSA and Azide free	HA610169	Human	WB
Cell treatment	uPA Receptor Mouse mAb [A6E8]	HA601039	Human	WB
	uPA Receptor Rabbit pAb	ER1918-43	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Cell treatment	uPA Receptor Recombinant Mouse mAb [A6E8-R]	HA601296	Human	WB
	UTRO Rabbit pAb	ER61389	Human	IHC-P
	VAPA Recombinant Rabbit mAb [JE77-48]	HA722372	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	VASP Recombinant Rabbit mAb [SY02-16]	ET1607-6	Human	WB,IF-Cell,IHC-P
	VCAM1 Rabbit pAb	R1512-13	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression, KD	VCAM1 Recombinant Rabbit mAb [SA05-04]	ET1601-18	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
Relative expression, KD	VCAM1 Recombinant Rabbit mAb [SA05-04] - BSA and Azide free	HA750009	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	VE Cadherin Recombinant Rabbit mAb [JJ0941]	ET1701-93	Mouse	WB,IP
KD	VE Cadherin Recombinant Rabbit mAb [PSH03-15] - BSA and Azide free	HA750856	Human	WB,IF-Cell,FC
KD	VE Cadherin Recombinant Rabbit mAb [PSH03-15]	HA721943	Human	WB,IF-Cell,FC
	Versican Recombinant Rabbit mAb [JB42-32]	ET7107-09	Human, Mouse, Rat	IHC-P,WB,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	VGLUT1 Recombinant Mouse mAb [PSH08-34]	HA601371	Mouse, Rat	WB,IHC-P
	VGLUT1 Recombinant Antibody [PSH10-53] - Mouse IgG1 (Chimeric)	HA601450	Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	VGLUT1 Recombinant Antibody [PSH10-53] - Rat IgG1 (Chimeric)	HA601392	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB,IF-Tissue
	VGLUT1 Recombinant Antibody [PSH10-53] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610283	Mouse, Rat	IHC-Fr,WB
Relative expression	VGLUT1 Recombinant Rabbit mAb [PSH10-53] - BSA and Azide free	HA751356	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IHC-P,IF-Tissue
	VGLUT1 Rabbit pAb	HA500488	Human, Mouse, Rat	WB
Relative expression	VGLUT1 Recombinant Antibody [PSH10-53] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610234	Mouse, Rat, Cynomolgus monkey, Pig	IHC-Fr,IHC-P,WB,IF-Tissue
	VGLUT1 Recombinant Mouse mAb [PSH08-34] - BSA and Azide free	HA610217	Mouse, Rat	WB,IHC-P
	VGLUT1 Recombinant Antibody [PSH10-53] - Guinea pig IgG2 (Chimeric)	HA601449	Mouse, Rat	IHC-Fr,WB
	VGLUT1 Rabbit pAb	HA500083	Human, Mouse, Rat	WB,IHC-P
	VGLUT1 Recombinant Antibody [PSH10-53] - Chicken IgY (Chimeric)	HA601478	Mouse	IHC-Fr
	VGLUT1 Recombinant Antibody [PSH10-53] - Chicken IgY (Chimeric) - BSA and Azide free	HA610312	Mouse	IHC-Fr
	VGLUT1 Recombinant Antibody [PSH10-53] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610284	Mouse, Rat	IHC-Fr,IHC-P,WB
	VGLUT1 Mouse mAb [A7A3]	HA601014	Human, Mouse, Rat	WB
Relative expression	VGLUT1 Recombinant Rabbit mAb [PSH10-53]	HA723217	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-Fr,IHC-P,IF-Tissue
	VIL1 Mouse mAb [G11-5E]	EM1706-46	Human	WB,IF-Cell,IHC-P,FC
	Villin1 Recombinant Rabbit mAb	IRS243RB	Mouse	mIHC
	Villin1 Rabbit pAb	ER1901-31	Human, Mouse	WB,FC
	Villin1 Rabbit pAb	HA500509	Human, Mouse, Rat	WB,IHC-P
	Villin1 Recombinant Rabbit mAb [C4] - BSA and Azide free	HA750476	Human	WB,IHC-P
	Villin1 Recombinant Rabbit mAb [A2]	HA721134	Human	WB,IHC-P
Relative expression	Villin1 Recombinant Rabbit mAb [JE39-08]	ET7110-55	Human	WB,IHC-P,IP,IF-Tissue
Relative expression	Villin1 Recombinant Rabbit mAb [JU34-75]	ET7106-62	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Villin1 Recombinant Rabbit mAb [C4]	HA721129	Human	WB,IHC-P
	Villin1 Recombinant Rabbit mAb [A2] - BSA and Azide free	HA750477	Human	WB,IHC-P
	Villin1 Recombinant Rabbit mAb	IRS213RB	Human	mIHC
	Vimentin Recombinant Rabbit mAb	IRS044	Human	mIHC
	Vimentin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS214RT	Human	mIHC
	iFluor™ 488 Conjugated Vimentin Recombinant Rabbit mAb [SC60-05]	HA720166F	Human, Mouse, Rat	IF-Cell,IF-Tissue,FC
KO	Vimentin Mouse mAb [D4-B11]	EM0401	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue,IHC-Fr,IP
Relative expression	Vimentin Recombinant Mouse mAb [A6-C1-R] - BSA and Azide free	HA610126	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	iFluor™ 647 Conjugated Vimentin Recombinant Rabbit mAb [SC60-05]	HA720165F	Human, Mouse, Rat	IF-Cell,IF-Tissue,FC
	Vimentin Recombinant Rabbit mAb	IRS060	Mouse	mIHC
	Vimentin Recombinant Rabbit mAb [PDH0-01]	HA721174	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,IF-Tissue,mIHC
KO, Relative expression	Vimentin Recombinant Rabbit mAb [SC60-05] - BSA and Azide free	HA750215	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC,IP
KO, Relative expression	Vimentin Recombinant Rabbit mAb [SC60-05]	ET1610-39	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue,IF-Cell,FC,IP
Relative expression, KO	Vimentin Mouse mAb [A6-C1]	M1412-1	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Relative expression	Vimentin Recombinant Mouse mAb [A6-C1-R]	HA601251	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KO	Vimentin Rabbit pAb	R1308-6	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Vimentin Recombinant Antibody - Rat IgG1 (Chimeric)	IRS263RT	Mouse	mIHC
KD	Vinculin Recombinant Rabbit mAb [JM42-43] - BSA and Azide free	HA750448	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
KD	Vinculin Recombinant Rabbit mAb [JM42-43]	ET1705-94	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	WASF3 Rabbit pAb	ER1918-61	Mouse, Rat	WB,IHC-P,IF-Tissue
KD	WASL Rabbit pAb	ER1803-16	Human, Mouse, Rat	WB,IF-Cell,FC
	WAVE1 Rabbit pAb	ER65311	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	WDR1 Recombinant Rabbit mAb [JE64-32]	HA721053	Human, Mouse, Rat	WB,IHC-P
	X11 gamma Rabbit pAb	ER65318	Human	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	YKL-40 / CHI3L1 Mouse mAb [A3G10]	EM1902-31	Human	WB,IHC-P
	YKL-40 / CHI3L1 Mouse mAb [A3G13]	HA600075	Human	WB
	YKL-40 / CHI3L1 Mouse mAb [10C1]	EM1902-14	Human, Mouse	WB,IHC-P,IF-Cell,FC,IF-Tissue
Relative expression	YKL-40 / CHI3L1 Recombinant Rabbit mAb [PSH16-26]	HA723824	Human	WB,IF-Cell
Relative expression	YKL-40 / CHI3L1 Recombinant Rabbit mAb [PSH16-26] - BSA and Azide free	HA751594	Human	WB,IF-Cell
	ZO2 Rabbit pAb	R1402-2	Human, Mouse	WB,IF-Cell
	ZO3 Rabbit pAb	ER61305	Human, Mouse, Monkey	WB,IF-Cell
	Zyxin Recombinant Rabbit mAb [JU83-12]	ET1706-18	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层