

TSA染色试剂盒清单

高级验证	产品名	货号	反应种属	验证应用
	抗体洗脱试剂盒 IRISKit® HyperView Advanced Ab Stripping Kit	900402		mIHC
	荧光淬灭试剂盒 IRISKit®HyperView Quench Buffer	900102		mIHC
	荧光淬灭试剂盒 IRISKit®HyperView Quench Buffer	900103		mIHC
	虹和CmTSA多重免疫荧光染色试剂盒(20T-4C, 1Rabbit+1Mouse+1Rat)IRISKitCmTSA Kit	900801		mIHC
	虹和CmTSA多重免疫荧光染色试剂盒(套装20T-4C, 1Mouse+2Rat)IRISKitCmTSA Kit	900807		mIHC
	虹和CmTSA多重免疫荧光染色试剂盒(套装20T-4C, 1Rabbit+2Mouse)IRISKitCmTSA Kit	900803		mIHC
	虹和CmTSA多重免疫荧光染色试剂盒(套装20T-4C, 1Rabbit+2Rat)IRISKitCmTSA Kit	900805		mIHC
	虹和CmTSA多重免疫荧光染色试剂盒(套装20T-4C, 2Mouse+1Rat)IRISKitCmTSA Kit	900806		mIHC
	虹和CmTSA多重免疫荧光染色试剂盒(套装20T-4C, 2Rabbit+1Mouse)IRISKitCmTSA Kit	900802		mIHC
	虹和CmTSA多重免疫荧光染色试剂盒(套装20T-4C, 2Rabbit+1Rat)IRISKitCmTSA Kit	900804		mIHC
	虹和CmTSA多重免疫荧光染色试剂盒(套装20T-4C, 3Mouse)IRISKitCmTSA Kit	900809		mIHC
	虹和CmTSA多重免疫荧光染色试剂盒(套装20T-4C, 3Rabbit)IRISKitCmTSA Kit	900808		mIHC
	虹和CmTSA多重免疫荧光染色试剂盒(套装20T-4C, 3Rat)IRISKitCmTSA Kit	900810		mIHC