

泛修饰抗体-产品清单

高级验证	产品名	货号	反应种属	验证应用
Cell treatment	Acetylated-Lysine Recombinant Rabbit mAb [PSH09-16]	HA723074	Species independent	WB,IHC-P,IF-Tissue
Cell treatment	Acetylated-Lysine Recombinant Rabbit Multiclonal Antibody [PSH10-04]	HA723173	Species independent	WB,IHC-P,ChIP,IF-Tissue
Cell treatment	Acetylated-Lysine Recombinant Rabbit mAb [PSH09-15]	HA723073	Species independent	WB,IHC-P,IF-Tissue,ChIP,IP
Cell treatment	Mono-Methyl Arginine (R*GG) Recombinant Rabbit mAb [PSH07-52]	HA722840	Species independent	WB
Cell treatment	Pan-Lactyl-lysine Recombinant Rabbit mAb [PSH03-73]	HA722037	Species independent	WB,IHC-P,IF-Cell,Dot Blot
Cell treatment	Pan-Lactyl-lysine Recombinant Rabbit mAb [PSH03-74]	HA722038	Species independent	WB,IHC-P,IF-Cell,Dot Blot
	Phosphoserine Mouse mAb [5B12]	ER60006	Species independent	WB,IHC-P
Cell treatment	Phosphothreonine Recombinant Rabbit mAb [PSH14-34]	HA723628	Human, Mouse, Rat, Species independent	WB,IHC-P,IF-Cell,IP,ELISA
Cell treatment	Phosphotyrosine Recombinant Rabbit mAb [JA10-49]	ET1704-20	Species independent	WB,IF-Cell,IF-Tissue,IP,FC