

mTOR信号通路研究-产品清单

高级验证	产品名	货号	反应种属	验证应用
	4E-BP1 Mouse mAb [3G1]	RT1002	Human, Mouse, Rat	WB,IP,IF,IHC-P
	4E-BP1 Rabbit pAb	ER62532	Mouse, Rat, Monkey	WB,IHC-P
KD	AKT1 Recombinant Rabbit mAb [ST05-09]	ET1609-47	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
	AKT1 Mouse mAb [D9-9-C9]	EM40507	Human, Mouse, Rat	WB,IHC-P
	AKT1 Recombinant Rabbit mAb [PSH0-50] - BSA and Azide free	HA750627	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
KD	AKT1 Recombinant Rabbit mAb [ST05-09] - BSA and Azide free	HA750177	Human, Mouse, Rat	WB,IHC-P,IP,FC,IHC-Fr,IF-Cell
	AKT1 Recombinant Rabbit mAb [PSH0-50]	HA721364	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
	AKT1/2/3 Recombinant Rabbit mAb [ST48-09]	ET1609-51	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75-09]	HA721870	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	AKT1/2/3 Recombinant Rabbit mAb [ST48-09] - BSA and Azide free	HA750179	Human, Mouse, Rat, Monkey	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
	AKT1/2/3 Rabbit pAb	ER1901-33	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	AKT1/2/3 Recombinant Rabbit mAb [JE75-09] - BSA and Azide free	HA750840	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	Akt1/3 Rabbit pAb	ER62638	Human, Mouse, Rat	WB,IHC-P
	AKT1S1 Mouse mAb [6H1C6]	EM1712-37	Human	WB,IHC-P,IF-Cell
	Akt2 Mouse mAb [2A1]	RT1027	Human, Mouse, Rat	WB,IP,IF,IHC-P
	AKT2 Rabbit pAb	HA500091	Human, Mouse, Rat	WB,IHC-P
KO	AMPK alpha 1 Recombinant Rabbit mAb [SU03-48]	ET1608-40	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC,IP,IHC-Fr,IHC-P
	AMPK alpha 1 Mouse mAb [B1-E10]	M1409-1	Human	WB,IF-Cell
KO	AMPK alpha 1 Recombinant Rabbit mAb [SU03-48] - BSA and Azide free	HA750147	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,FC,IP,IHC-Fr,IHC-P
	AMPK alpha 1/2 Recombinant Rabbit mAb [PSH08-86] - BSA and Azide free	HA751251	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	AMPK alpha 1/2 Recombinant Rabbit mAb [PSH08-86]	HA723038	Human, Mouse, Rat	WB,IHC-P
	AMPK alpha 1/2 Rabbit pAb	ER62660	Human, Mouse, Rat, Monkey	WB,IHC-P
	AMPK alpha 2 Mouse mAb [A6A10]	HA600078	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	AMPK alpha 2 Mouse mAb [A6A12]	HA600080	Human, Mouse, Rat	WB
	AMPK alpha 2 Mouse mAb [A6A11]	HA600079	Human, Mouse, Rat	IF-Cell,IHC-P
KD	BRAF Mouse mAb [F1-F4-E10]	EM1701-32	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
KD	BRAF Recombinant Rabbit mAb [SU34-04]	ET1608-36	Human, Mouse, Rat	WB,IHC-P,FC,IP
	CAB39 Recombinant Rabbit mAb [PSH15-16] - BSA and Azide free	HA751545	Human, Mouse, Rat, Monkey	WB,IHC-P,IP
	CAB39 Recombinant Rabbit mAb [PSH15-16]	HA723720	Human, Mouse, Rat, Monkey	WB,IHC-P,IP
	CAB39 Rabbit pAb	ER60546	Human, Mouse, Rat, Monkey	WB
	CB39L Rabbit pAb	ER60547	Human, Mouse	WB
	CD98 Mouse mAb [A9D8]	HA601124	Human	WB,IHC-P
	CD98 Mouse mAb [A9D7]	HA601125	Human	WB,IHC-P
	CD98 Mouse mAb	IRS130MS	Human	mIHC
	CD98 Rabbit pAb	ER65634	Human, Mouse, Rat	WB
	CLIP170 Recombinant Rabbit mAb [JE63-82]	HA721004	Human	WB,IF-Cell,IHC-P
Cell treatment	DDIT4 Recombinant Rabbit mAb [PS00-10]	HA722161	Human	WB
Cell treatment	DDIT4 Recombinant Rabbit mAb [PS00-10] - BSA and Azide free	HA750969	Human	WB
	Dishevelled 2 Recombinant Rabbit mAb [PSH04-07]	HA722088	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Dishevelled 2 Recombinant Rabbit mAb [PSH04-07] - BSA and Azide free	HA750930	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Dishevelled 3 Recombinant Rabbit mAb [PSH15-93]	HA723798	Human, Mouse, Rat, Monkey	WB,IP
	Dishevelled 3 Recombinant Rabbit mAb [PSH15-93] - BSA and Azide free	HA751584	Human, Mouse, Rat, Monkey	WB,IP
	DPTOR Recombinant Rabbit mAb [PH00-03] - BSA and Azide free	HA750555	Human, Mouse, Rat	WB,IHC-P
	DPTOR Recombinant Rabbit mAb [PH00-03]	HA721225	Human, Mouse, Rat	WB,IHC-P
	DVL1 Rabbit pAb	ER1907-85	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Dvl-2 Rabbit pAb	ER63295	Human, Mouse	WB,IHC-P
	Dvl-3 Rabbit pAb	ER63296	Human, Mouse	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	EIF4B Mouse mAb [D7-F5]	EM1709-90	Human, Mouse	WB,IF-Cell,IHC-P,FC
	eIF4E Recombinant Rabbit mAb [JE60-70]	HA721192	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	eIF4EBP1 Recombinant Rabbit mAb [JJ09-25]	ET1701-83	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	eIF4EBP1 Rabbit pAb	ER1706-64	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	ERK1 Recombinant Rabbit mAb [SP05-09]	ET1604-16	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	ERK1 Recombinant Rabbit mAb [SP00-74]	ET1604-32	Human	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	ERK1 Mouse mAb [7-D6-E5]	M1505-9	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Erk1/2 Mouse mAb [2G1]	RT1484	Human, Mouse, Rat	WB,IP,IF
	ERK1/2 Rabbit pAb	ER131011	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03]	ET1601-29	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
KD	ERK1/2 Recombinant Rabbit mAb [SA43-03] - BSA and Azide free	HA750019	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IP,FC,IHC-P,IHC-Fr
KD	ERK2 Recombinant Mouse mAb [5-D2-R] - BSA and Azide free	HA610080	Human, Mouse, Rat	WB
KD	ERK2 Recombinant Mouse mAb [5-D2-R]	HA601202	Human, Mouse, Rat	WB
	ERK2 Mouse mAb [C3B12]	EM1901-55	Human, Mouse, Rat	WB,IHC-P,FC
	ERK2 Mouse mAb [5-D2]	M1407-3	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	ERK2 Mouse mAb [A10D2]	EM1901-54	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Mouse mAb [B10-8]	EM1901-53	Human, Mouse, Rat	WB,IHC-P,FC
KD	ERK2 Recombinant Mouse mAb [C3B12-R] - BSA and Azide free	HA610079	Human, Mouse, Rat	WB
KD	ERK2 Recombinant Mouse mAb [C3B12-R]	HA601201	Human, Mouse, Rat	WB
KD	ERK2 Rabbit pAb	ER131218	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	ERK2 Recombinant Rabbit mAb [SZ25-01]	ET1603-23	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC,IHC-Fr
KD	FLCN Recombinant Rabbit mAb [PSH07-68] - BSA and Azide free	HA751180	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
KD	FLCN Recombinant Rabbit mAb [PSH07-68]	HA722883	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Frizzled-1 Rabbit pAb	ER63469	Human	WB,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Frizzled 10/CD350 Rabbit pAb	ER1909-02	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	Frizzled-2 Rabbit pAb	ER63470	Human, Mouse, Rat	WB,IF-Cell
	Frizzled-4 Rabbit pAb	ER65833	Human, Mouse, Rat	IHC-P
	Frizzled-5 Rabbit pAb	0905-3	Human	WB,IF-Cell,FC
	Frizzled-6 Rabbit pAb	ER63471	Human	WB,IHC-P,IF-Cell
	Frizzled 7 Rabbit pAb	HA500433	Human	WB
	Frizzled 7 Rabbit pAb	ER1909-03	Human, Mouse, Rat	WB,FC
	Frizzled 8 Recombinant Rabbit mAb [JB40-16]	ET7107-68	Human, Mouse	WB,IHC-P
	Frizzled-9 Rabbit pAb	ER63472	Human, Mouse, Monkey	WB,IF-Cell
	FZD2 Rabbit pAb	HA500406	Human, Mouse	WB,IHC-P
	FZD3 Rabbit pAb	ER1909-16	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	G21 protein/NPR2L Rabbit pAb	ER1909-19	Human, Rat	WB,IHC-P,IF-Tissue
	GRB10 Rabbit pAb	ER63631	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	GRB2 Rabbit pAb	ER31203	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	GRB2 Recombinant Rabbit mAb [JE49-73]	HA721923	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
	GRB2 Mouse mAb [11-B11]	EM1701-41	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	GSK3 alpha Recombinant Rabbit mAb [JE33-27]	HA722217	Human, Mouse, Rat, Monkey	WB,IHC-P
KD	GSK3 alpha Rabbit pAb	R1312-1	Human	WB,IF-Cell,IHC-P,FC
	GSK3 alpha/beta Rabbit pAb	ER63648	Human, Mouse, Rat	WB,IHC-P
KD	GSK3 beta Recombinant Rabbit mAb [SY28-03]	ET1607-71	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
KD	GSK3 beta Recombinant Rabbit mAb [SY28-03] - BSA and Azide free	HA750131	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
	GSK3 beta Rabbit pAb	R1510-4	Human, Mouse, Rat	IF-Cell,IHC-P,FC
	GTPase HRAS Recombinant Rabbit mAb [JE42-15]	HA721355	Human, Mouse, Rat	WB,IHC-P
	Hamartin Recombinant Rabbit mAb [PSH14-67] - BSA and Azide free	HA751528	Human, Mouse, Rat, Monkey	WB,IP
	Hamartin Rabbit pAb	HA500199	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Hamartin Mouse mAb [A6F1]	HA600004	Human, Mouse	WB,IHC-P
	Hamartin Recombinant Rabbit mAb [PSH14-67]	HA723665	Human, Mouse, Rat, Monkey	WB,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	HBXIP Rabbit pAb	ER60748	Human	WB
Cell treatment	HIF-1 alpha Rabbit pAb	ER1802-41	Human	WB,IHC-P
Cell treatment	HIF-1 alpha Recombinant Rabbit mAb [PSH07-00]	HA722778	Human, Mouse	WB,IF-Cell,IHC-P
Cell treatment	HIF-1 alpha Recombinant Rabbit mAb [PSH07-00] - BSA and Azide free	HA751127	Human, Mouse	WB,IF-Cell,IHC-P
Cell treatment	HIF-1 alpha Recombinant Rabbit mAb [JE75-33]	HA721997	Human, Mouse, Rat	WB,IF-Cell
	HIF-1 alpha Mouse mAb [5H7]	RT1278	Human, Mouse, Rat	WB,IP,ICC,IHC-P
KD	HIF-1 beta Recombinant Rabbit mAb [ST05-25] - BSA and Azide free	HA750187	Human, Mouse, Rat	WB,IHC-P
KD	HIF-1 beta Recombinant Rabbit mAb [ST05-25]	ET1609-67	Human, Mouse, Rat	WB,IHC-P
	HIF-2 alpha Recombinant Rabbit mAb [JB24-42]	ET7107-32	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	HRP Conjugated Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - Detector	HA723374H	Human	ELISA(Det)
	Biotin Conjugated Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - BSA and Azide free (Detector)	HA723160B	Human	ELISA(Det),ELISA
	Human sTNF RI Recombinant Rabbit mAb [PSH09-89] - BSA and Azide free (Capture)	HA723158	Human	ELISA(Cap)
	Human sTNF RI Recombinant Rabbit mAb [PSH09-90] - BSA and Azide free (Detector)	HA723159	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb	IRS040	Human	mIHC
	Human TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free (Detector)	HA721266	Human	ELISA(Det)
	Biotin Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721559	Human	ELISA(Det)
	HRP Conjugated Human TNF alpha Recombinant Rabbit mAb [PS01-07] - Detector	HA721294	Human	ELISA(Det)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free	HA750593	Human	IHC-P
	Human TNF alpha Recombinant Rabbit mAb [PS00-40] - BSA and Azide free (Capture)	HA721265	Human	ELISA(Cap)
	Human TNF alpha Recombinant Rabbit mAb [PS00-40]	HA721286	Human	IHC-P,mIHC
	Biotin Conjugated Human WNT7A Recombinant Rabbit mAb [PSH19-62] - Detector	HA725314B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human WNT7A Recombinant Rabbit mAb [PSH19-61] - BSA and Azide free (Capture)	HA725313	Human	ELISA(Cap)
	Human WNT7A Recombinant Rabbit mAb [PSH19-62] - BSA and Azide free (Detector)	HA725314	Human	ELISA(Det)
Cell treatment	IGF1 Recombinant Rabbit mAb [JE30-01]	HA722527	Human	WB
	IGF1 (E peptide) Rabbit pAb	HA500107	Human	WB,IHC-P,FC
	IGF1R Mouse mAb [D5-G2]	HA600033	Human	WB,FC
	IGF1R Rabbit pAb	R1310-6	Human	WB,IHC-P
Relative expression	IGF1 Receptor beta Recombinant Rabbit mAb [PSH11-72]	HA723368	Human, Mouse, Rat	WB,IF-Cell
Relative expression	IGF1 Receptor beta Recombinant Rabbit mAb [PSH11-72] - BSA and Azide free	HA751427	Human, Mouse, Rat	WB,IF-Cell
	IGF I Rabbit pAb	ER1911-10	Human, Mouse, Rat, Dog	WB,IHC-P,IF-Tissue
	IGF-IR Rabbit pAb	ER63734	Human, Mouse, Rat	IF-Cell,WB,IHC-P
	IKK alpha Recombinant Rabbit mAb [SN06-99] - BSA and Azide free	HA750246	Human, Mouse, Rat	WB,IHC-P,IP,FC
	IKK alpha Rabbit pAb	ER30911	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	IKK alpha Recombinant Rabbit mAb [SN06-99]	ET1611-15	Human, Mouse, Rat	WB,IHC-P,IP,FC
KD	IKK alpha + IKK beta Recombinant Rabbit mAb [SN63-02] - BSA and Azide free	HA750248	Human, Mouse, Rat	WB,IF-Cell,IP
KD	IKK alpha + IKK beta Recombinant Rabbit mAb [SN63-02]	ET1611-23	Human, Mouse, Rat	WB,IF-Cell,IP
KD	IKK beta Recombinant Rabbit mAb [JE32-81]	HA721715	Human	WB,FC
	IKK beta Rabbit pAb	R1706-13	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Insulin Recombinant Rabbit mAb	IRS233RB	Mouse	mIHC
	Insulin Mouse mAb [A6-6]	EM80714	Human, Mouse, Rat	IHC-P,ELISA,IF-Tissue,IHC-Fr
Relative expression	Insulin Recombinant Rabbit mAb [SA0410]	ET1601-12	Human, Mouse, Rat	IF-Cell,IF-Tissue,IHC-P,mIHC,IHC-Fr,WB
	Insulin Recombinant Rabbit mAb [PSH04-15] - BSA and Azide free (Capture)	HA722096	Human, Mouse, Rat	ELISA(Cap)
	Insulin Recombinant Rabbit mAb [PSH04-16] - BSA and Azide free (Detector)	HA722097	Human, Mouse, Rat	ELISA(Det)
Relative expression	Insulin Recombinant Rabbit mAb [SA0410] - BSA and Azide free	HA750004	Human, Mouse, Rat	IF-Cell,IF-Tissue,IHC-P,IHC-Fr,WB
	Insulin Recombinant Rabbit mAb	IRS161RB	Human	mIHC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Insulin Mouse mAb [5E1]	EM1701-33	Human	IHC-P,ELISA
	Biotin Conjugated Insulin Recombinant Rabbit mAb [PSH04-16] - Detector	HA722610B	Human, Mouse, Rat	ELISA(Det),ELISA
	Insulin B Chain Rabbit pAb	0807-11	Human, Mouse, Rat	IHC-P,ELISA,IF-Tissue
KD, Relative expression	Insulin Receptor Beta Recombinant Rabbit mAb [SN20-13]	ET1611-90	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC
KD, Relative expression	Insulin Receptor Beta Recombinant Rabbit mAb [SN20-13] - BSA and Azide free	HA750276	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC
	IRS1 Rabbit pAb	ER2001-35	Human, Rat	WB,IHC-P,IF-Cell,FC
	IRS1 Recombinant Rabbit mAb [SD2046]	ET1612-45	Human	WB,IF-Cell,IF-Tissue,IHC-P
	K-Ras Rabbit pAb	ER65706	Human, Mouse, Rat	WB,IHC-P
	KRAS Mouse mAb [A8E4]	HA601058	Human, Mouse, Rat	WB,IF-Cell,FC
	KRAS Recombinant Mouse mAb [A8E6] - BSA and Azide free	HA610008	Human, Mouse, Rat	WB,IHC-P
	KRAS Mouse mAb [A8E5]	HA601059	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	KRAS Recombinant Mouse mAb [A8E6]	HA601060	Human, Mouse, Rat	WB,IHC-P
	Lipin 1 Recombinant Rabbit mAb [JG37-69]	ET7108-51	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
Relative expression	LKB1 Recombinant Rabbit mAb [PSH09-29] - BSA and Azide free	HA751283	Human, Mouse, Rat, Monkey	WB,IHC-P,IP
	LKB1 Rabbit pAb	HA500143	Human, Mouse, Rat	WB,IHC-P
Relative expression	LKB1 Recombinant Rabbit mAb [PSH09-29]	HA723087	Human, Mouse, Rat, Monkey	WB,IHC-P,IP
	LKB1 Rabbit pAb	ER1912-41	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,FC
	LRP6 Mouse mAb [A5E7]	HA600065	Human	WB,IHC-P
	LRP6 Mouse mAb [1G8]	RT1369	Human, Mouse	WB,IP,IF,IHC-P
	LTOR2 Rabbit pAb	ER60325	Human, Mouse	WB
	LTOR3 Rabbit pAb	ER61476	Human, Mouse, Rat	WB
KD	MEK1 Recombinant Rabbit mAb [SZ22-01]	ET1603-20	Human, Mouse, Rat, Cow, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK1/2 Recombinant Rabbit mAb [SR13-07]	ET1602-3	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MEK1/2 Recombinant Rabbit mAb [SR13-07] - BSA and Azide free	HA750042	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	MEK2 Recombinant Rabbit mAb [SD207-08]	ET1612-6	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	MLST8 Rabbit pAb	HA500372	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-29] - BSA and Azide free (Capture)	HA721974	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30]	HA722626B	Mouse, Rat	ELISA(Det)
	Mouse/Rat TNF alpha Recombinant Rabbit mAb [PSH03-30] - BSA and Azide free (Detector)	HA721975	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77]	HA722623B	Mouse	ELISA(Det)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-77] - BSA and Azide free (Detector)	HA722041	Mouse	ELISA(Det)
	Mouse TNF alpha Recombinant Rabbit mAb [PSH03-76] - BSA and Azide free (Capture)	HA722040	Mouse	ELISA(Cap)
	mTOR Rabbit pAb	HA500126	Human, Rat, Mouse	WB,IF-Cell,IHC-P
	mTOR Recombinant Rabbit mAb [SU30-00]	ET1608-5	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell
	mTOR Rabbit pAb	R1510-21	Human, Mouse, Rat	IF-Cell,IHC-P
	mTOR Recombinant Rabbit mAb [SU30-00] - BSA and Azide free	HA750153	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell
	N/H/K-Ras Rabbit pAb	ER64095	Human, Mouse, Rat	WB,IHC-P
	p70 S6 Kinase Rabbit pAb	HA500095	Human, Mouse, Rat	WB,IHC-P
	p70 S6 Kinase Recombinant Rabbit mAb [PSH13-50] - BSA and Azide free	HA751469	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	p70 S6 Kinase Rabbit pAb	ER1916-13	Human, Rat	WB,IHC-P,IF-Tissue
	p70 S6 Kinase Recombinant Rabbit mAb [PSH13-50]	HA723539	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
	p70 S6 Kinase Recombinant Rabbit mAb [JE03-20]	HA722520	Human, Mouse, Rat	WB,IHC-P,IP
	p70 S6 Kinase beta Rabbit pAb	ER31205	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	P70 S6 Kinase beta Recombinant Rabbit mAb [JE35-69]	HA721354	Human, Mouse	WB,IF-Cell,FC,IHC-P
	PDPK1 Recombinant Rabbit mAb [SD204-9]	ET1612-27	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Phospho-4E-BP1/2/3 (T45) Mouse mAb [3G2]	RT1004	Human, Mouse, Rat	WB,IP,IF
Cell treatment	Phospho-AKT1 (S124) Recombinant Rabbit mAb [JJ08-46]	ET1701-36	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-AKT1 (S129) Recombinant Rabbit mAb [JE59-50]	HA721677	Human, Mouse, Rat	WB
Cell treatment	Phospho-AKT1 (T450) Recombinant Rabbit mAb [SD08-12]	ET1612-73	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell,FC
Cell treatment	Phospho-AKT2 (S474) Rabbit pAb	HA500116	Human, Mouse	WB,IHC-P
Cell treatment	Phospho-AKT2 (S474) Recombinant Rabbit mAb [PSH09-46] - BSA and Azide free	HA751293	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-AKT2 (S474) Recombinant Rabbit mAb [PSH09-46]	HA723108	Human, Mouse, Rat	WB,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05] - BSA and Azide free	HA750132	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44] - BSA and Azide free	HA750953	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH04-44]	HA722129	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [SY28-05]	ET1607-73	Human, Mouse, Rat, Dog	WB,IF-Cell,IF-Tissue,IHC-P,IHC-Fr
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10] - BSA and Azide free	HA751205	Human	WB,IF-Cell,FC
Cell treatment	Phospho-AKT (T308) Recombinant Rabbit mAb [PSH08-10]	HA722951	Human	WB,IF-Cell,FC
	Phospho-AMPK alpha 1 (S485) Recombinant Rabbit mAb [PSH18-01]	HA723966	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-AMPK alpha 1 (S496) Recombinant Rabbit mAb [SD0810]	ET1612-72	Human, Mouse	WB,IP,IHC-Fr
Cell treatment	Phospho-AMPK alpha 1 (S496) Recombinant Rabbit mAb [SD0810] - BSA and Azide free	HA750300	Human, Mouse	WB,IP,IHC-Fr
Cell treatment	Phospho-AMPK alpha 1 (T183) +AMPK alpha 2 (T172) Recombinant Rabbit mAb [PSH14-13]	HA723603	Human, Mouse, Rat	WB,IHC-P
	Phospho-AMPK alpha 1 (T183) +AMPK alpha 2 (T172) Recombinant Rabbit mAb [PSH20-31]	HA724142	Human	WB
Cell treatment	Phospho-AMPK alpha 1 (T183) +AMPK alpha 2 (T172) Recombinant Rabbit mAb [PSH14-13] - BSA and Azide free	HA751508	Human, Mouse, Rat	WB,IHC-P
	Phospho-AMPK alpha 2 (S345) Recombinant Rabbit mAb [JJ08-19]	ET1701-37	Human, Mouse, Rat	WB,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-AMPK alpha 2 (S491) Recombinant Rabbit mAb [JE59-49]	HA722567	Human	WB
Cell treatment	Phospho-BRAF (T401) Recombinant Rabbit mAb [JJ08-72]	ET1701-20	Human, Mouse, Rat	WB,IF-Cell,IP,IHC-P
Cell treatment	Phospho-eIF4B (S406) Recombinant Rabbit mAb [JE63-80]	HA722149	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-eIF4B (S406) Recombinant Rabbit mAb [JE58-32]	HA721676	Human, Mouse, Rat	WB
Cell treatment	Phospho-eIF4EBP1 (S65) Recombinant Rabbit mAb [PSH09-44]	HA723106	Human	WB,IF-Cell,IP
Cell treatment	Phospho-eIF4EBP1 (S65) Recombinant Rabbit mAb [PSH09-44] - BSA and Azide free	HA751291	Human	WB,IF-Cell,IP
Cell treatment	Phospho-eIF4EBP1 (T70) Recombinant Rabbit mAb [PSH09-91] - BSA and Azide free	HA751360	Human, Mouse	WB,IHC-P,IP
Cell treatment	Phospho-eIF4EBP1 (T70) Recombinant Rabbit mAb [PSH09-91]	HA723222	Human, Mouse	WB,IHC-P,IP
Cell treatment	Phospho-eIF4E (S209) Recombinant Rabbit mAb [SU0396]	ET1608-66	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
Cell treatment	Phospho-eIF4E (S209) Recombinant Rabbit mAb [SU0396] - BSA and Azide free	HA750160	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4]	ET1603-22	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202)+Erk2 (T185) Recombinant Rabbit mAb [SZ2-4] - BSA and Azide free	HA750061	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IF-Tissue
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01]	ET1610-13	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-Erk1 (T202 + Y204) + Erk2 (T185 + Y187) Recombinant Rabbit mAb [SC58-01] - BSA and Azide free	HA750199	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue,FC,IP
Cell treatment	Phospho-GSK3 (alpha + beta) (Y216 + Y279) Recombinant Rabbit mAb [SY26-05]	ET1607-54	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,IF-Tissue,FC
Cell treatment	Phospho-GSK3 (alpha + beta) (Y216 + Y279) Recombinant Rabbit mAb [SY26-05] - BSA and Azide free	HA750127	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,IF-Tissue,FC
Cell treatment	Phospho-GSK3 alpha (S21) Recombinant Rabbit mAb [PSH09-11] - BSA and Azide free	HA751265	Human, Mouse	WB,IHC-P,Dot Blot
Cell treatment	Phospho-GSK3 alpha (S21) Recombinant Rabbit mAb [PSH09-11]	HA723069	Human, Mouse	WB,IHC-P,Dot Blot
Cell treatment	Phospho-GSK3 beta (S9) Recombinant Rabbit mAb [SY02-71] - BSA and Azide free	HA750128	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-GSK3 beta (S9) Recombinant Rabbit mAb [SY02-71]	ET1607-60	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-MEK1/2 (S218 + S222) Recombinant Rabbit mAb [ST0490]	ET1609-50	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Phospho-MEK1 (S298) Recombinant Rabbit mAb [SD206-7]	ET1612-40	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
Cell treatment	Phospho-MEK1 (T292) Recombinant Rabbit mAb [SD2088]	ET1612-42	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-mTOR (S2448) Mouse mAb [A5D5]	HA600094	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-mTOR (S2481) Recombinant Rabbit mAb [JE59-62]	HA721632	Human	WB
	Phospho-p70 S6 Kinase (S434) Mouse mAb [2G1]	RT1456	Human, Mouse, Rat	WB,IP,IF,IHC-P
Cell treatment	Phospho-p70 S6 Kinase (T389) Recombinant Rabbit mAb [PSH02-23]	HA721803	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-p70 S6 Kinase (T389) Recombinant Rabbit mAb [PSH02-23] - BSA and Azide free	HA750791	Human, Mouse, Rat	WB,IF-Cell,IHC-P
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-39] - BSA and Azide free	HA750717	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-39]	HA721673	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-38]	HA721672	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PI3K p85 (Y467) + PI3K p55 (Y199) Recombinant Rabbit mAb [PSH01-38] - BSA and Azide free	HA750716	Human, Mouse	WB,IF-Cell
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH16-27]	HA723825	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH16-27] - BSA and Azide free	HA751595	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH15-91] - BSA and Azide free	HA751582	Human, Mouse, Rat	WB,IF-Cell,Dot Blot
Cell treatment	Phospho-PKC alpha/beta II (T638/T641) Recombinant Rabbit mAb [PSH15-91]	HA723796	Human, Mouse, Rat	WB,IF-Cell,Dot Blot
Cell treatment	Phospho-PKC alpha (T497) Recombinant Rabbit mAb [JE44-81]	HA722278	Human, Mouse, Rat	WB,IF-Cell,FC
Cell treatment	Phospho-PKC alpha (T638) Recombinant Rabbit mAb [JF0964]	ET1702-17	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,Dot Blot
Cell treatment	Phospho-PKC delta (Y311) Recombinant Rabbit mAb [JE44-20]	HA722458	Human, Mouse, Rat	WB
Cell treatment	Phospho-PKC (S729) Recombinant Rabbit mAb [JE77-54]	HA722846	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-PKC zeta (T560) Recombinant Rabbit mAb [SN0753]	ET1611-78	Human, Mouse, Rat	WB,IHC-P
Cell treatment	Phospho-PRAS40 (T246) Recombinant Rabbit mAb [JE53-94]	HA721790	Human, Mouse, Rat	WB
	Phospho-PTEN (S380) Recombinant Rabbit mAb [JJ08-81]	ET1701-46	Human, Mouse, Rat	WB,IF-Cell
Cell treatment	Phospho-Raf1 (S259) Recombinant Rabbit mAb [SD85-07]	ET1612-87	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Raf1 (S259) Recombinant Rabbit mAb [SD85-07] - BSA and Azide free	HA750304	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	Phospho-Raf1 (S43) Recombinant Rabbit mAb [SD087-05]	ET1612-75	Human, Mouse, Rat	WB
	Phospho-Raf1 (S621) Recombinant Rabbit mAb [JJ085-05]	ET1701-3	Human, Rat, Mouse	WB,IHC-P
Cell treatment	Phospho-RICTOR (T1135) Recombinant Rabbit mAb [PSH09-92] - BSA and Azide free	HA751323	Human, Mouse	WB
Cell treatment	Phospho-RICTOR (T1135) Recombinant Rabbit mAb [PSH09-92]	HA723161	Human, Mouse	WB
	Phospho-RPS6 (S235 + S236) Recombinant Rabbit mAb [PS01-04] - BSA and Azide free	HA750587	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Phospho-RPS6 (S235 + S236) Recombinant Rabbit mAb [PS01-04]	HA721275	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Phospho-RPS6 (S240 + S244) Recombinant Rabbit mAb [PS00-65]	HA721589	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Phospho-RPS6 (S240 + S244) Recombinant Rabbit mAb [PS00-65] - BSA and Azide free	HA750684	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Phospho-RSK1 p90 (T359 + S363) Recombinant Rabbit mAb [SU03-65]	ET1608-53	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Cell treatment	Phospho-RSK1 p90 (T573) Recombinant Rabbit mAb [JE50-12]	HA721793	Human, Mouse, Rat	WB
Cell treatment	Phospho-RSK1 (S380) Recombinant Rabbit mAb [SC05-32]	ET1610-28	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IP
Cell treatment	Phospho-Tuberin (S939) Recombinant Rabbit mAb [JE54-06]	HA721762	Human, Mouse	WB
Cell treatment	Phospho-Tuberin (T1462) Recombinant Rabbit mAb [JE43-81]	HA721792	Human, Mouse, Rat	WB
Cell treatment	Phospho-ULK1 (S556) Recombinant Rabbit mAb [JE59-28]	HA721445	Human, Mouse	WB
	PI3 Kinase p110 beta Recombinant Rabbit mAb [JE65-38]	HA722474	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	PI 3 Kinase p110 delta Recombinant Rabbit mAb [JM83-24]	ET1705-46	Human, Mouse, Rat	WB,IHC-P
KD	PI3-kinase p110 subunit beta Recombinant Rabbit mAb [SC0580]	ET1610-36	Human	WB,IP,IHC-P,FC
	PI3 Kinase p110 α Recombinant Rabbit mAb [SJ0186]	ET1606-36	Human, Mouse	WB,IF-Cell,IF-Tissue,IP,IHC-P
	PI 3 Kinase p55 gamma Recombinant Rabbit mAb [SR44-06]	ET1602-27	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	PI 3 Kinase p85 alpha Rabbit pAb	R1312-6	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	PI 3 Kinase p85 alpha Recombinant Rabbit mAb [SU04-07] - BSA and Azide free	HA750162	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KD	PI 3 Kinase p85 alpha Mouse mAb [A3-D0]	M1510-2	Human	WB,IF-Cell
KD	PI 3 Kinase p85 alpha Mouse mAb [L3-D6]	EM1701-62	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
	PI 3 Kinase p85 alpha Recombinant Rabbit mAb [SU04-07]	ET1608-70	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
KD	PI 3 Kinase p85 alpha Recombinant Mouse mAb [A3-D0-R] - BSA and Azide free	HA610084	Human	WB,IF-Cell
KD	PI 3 Kinase p85 alpha Recombinant Mouse mAb [A3-D0-R]	HA601206	Human	WB,IF-Cell
	PI3 Kinase p85 alpha Recombinant Rabbit mAb [JE04-11]	HA722522	Human, Mouse, Rat	WB
	PI 3 Kinase p85 alpha/gamma Rabbit pAb	ER64588	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Tissue
	PI 3 Kinase p85 beta Recombinant Rabbit mAb [ST04-77]	ET1609-30	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	PKC Rabbit pAb	ER64620	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
KD	PKC alpha Recombinant Rabbit mAb [SU31-08]	ET1608-15	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IP
	PKC beta Rabbit pAb	ER1802-14	Human, Mouse	WB,IF-Cell,IHC-P
	PKC beta 1 Recombinant Rabbit mAb [JM43-26]	ET1705-30	Human, Mouse, Rat	WB,IHC-P,FC
	PKC beta 2 Recombinant Rabbit mAb [ST48-06]	ET1609-44	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	PKC delta Recombinant Rabbit mAb [JJ093-02]	ET1701-85	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	PKC delta Recombinant Rabbit mAb [JJ093-02] - BSA and Azide free	HA750325	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	PKC epsilon Recombinant Rabbit mAb [JA38-21]	ET1704-85	Human, Mouse, Rat	WB,IHC-P
	PKC eta Recombinant Rabbit mAb [JE32-41]	HA722389	Human, Mouse, Rat, Monkey	WB
	PKC gamma Recombinant Rabbit mAb [JG69-32]	ET7108-42	Human, Mouse, Rat	WB,IHC-P
	PKC theta Rabbit pAb	ER64621	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Relative expression	PKC zeta Recombinant Rabbit mAb [PSH06-74] - BSA and Azide free	HA751122	Human	WB
Relative expression	PKC zeta Recombinant Rabbit mAb [PSH06-74]	HA722773	Human	WB
	PRAS40 Recombinant Rabbit mAb [JE30-60]	HA721581	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	PRAS40 Rabbit pAb	ER64668	Human, Mouse, Rat	WB,IHC-P
	PRKAA2 Mouse mAb [H5-E11]	EM1708-49	Human, Monkey, Mouse, Rat	WB,IHC-P,FC,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PRR5 Rabbit pAb	HA500295	Human, Mouse, Rat	WB
	PTEN Mouse mAb [4D8]	RT1519	Human, Mouse, Rat	WB,IP,IF
KD	PTEN Rabbit pAb	ER1902-33	Human, Mouse, Rat	WB,IHC-P,FC
Relative expression, KD	PTEN Recombinant Rabbit mAb [SJ19-03] - BSA and Azide free	HA750105	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	PTEN Rabbit pAb	ER1706-51	Human, Mouse	WB,IF-Cell,IHC-P,FC
Relative expression, KD	PTEN Recombinant Rabbit mAb [SJ19-03]	ET1606-43	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	PTEN Rabbit pAb	ER1902-27	Human, Mouse, Rat	WB,IHC-P,FC
Relative expression, KD	PTEN Recombinant Rabbit mAb [JE51-96]	HA721962	Human, Mouse, Rat	WB,IHC-P,FC
KD	Raf1 Recombinant Rabbit mAb [JJ086-06]	ET1701-21	Human, Mouse, Rat	WB
	Rag B Rabbit pAb	ER64751	Human, Mouse, Rat	IHC-P
	Rag C Rabbit pAb	ER64752	Human, Mouse, Rat	WB
	Raptor Recombinant Rabbit mAb [JE59-39]	HA722028	Human, Mouse, Rat	WB,IHC-P
	Raptor Rabbit pAb	ER1802-57	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Ras Recombinant Rabbit mAb [JE44-13]	HA721883	Human, Mouse, Rat	WB,IF-Cell,FC
	Ras Recombinant Rabbit mAb [SA39-05]	ET1601-16	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,FC,IP
	Ras Rabbit pAb	ER40115	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	Ras Recombinant Rabbit mAb [JF10-11]	ET1702-94	Human, Mouse, Rat, Zebrafish	WB,IF-Cell,IP,FC
KD	RHEB Recombinant Rabbit mAb [JG37-12]	ET7108-44	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,FC
	RhoA Recombinant Rabbit mAb [SN0612] - BSA and Azide free	HA750242	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	RhoA Recombinant Rabbit mAb [SN0612]	ET1611-10	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	Rho A+B+C Recombinant Rabbit mAb [SR38-00]	ET1602-23	Human, Mouse, Rat	WB,IF-Cell,FC
	Rho B Rabbit pAb	ER64810	Human, Mouse, Rat	WB
	Ribosomal Protein S6 Rabbit pAb	ER64863	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	RICTOR Recombinant Rabbit mAb [1A5]	HA721142	Human, Mouse	WB,IHC-P,IF-Cell,FC
	RICTOR Mouse mAb [E9-B3]	EM1709-50	Human, Mouse, Monkey	WB,IF-Cell,IHC-P,FC

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	RICTOR Recombinant Rabbit mAb [1A5] - BSA and Azide free	HA750484	Human, Mouse	WB,IHC-P,IF-Cell,FC
	RICTOR Recombinant Rabbit mAb [1D5] - BSA and Azide free	HA750483	Human	WB,IF-Cell,FC
	RICTOR Recombinant Rabbit mAb [1D5]	HA721141	Human	WB,IF-Cell,FC
	RPS6 Mouse mAb [A6B9]	HA600084	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	RPS6 Recombinant Mouse mAb [A6B8-R]	HA601324	Human, Mouse, Rat	WB
	RPS6 Recombinant Mouse mAb [A6B8-R] - BSA and Azide free	HA610196	Human, Mouse, Rat	WB
	RPS6 Mouse mAb [A6B8]	HA600081	Human, Mouse, Rat	WB,IHC-P
	RPS6 Recombinant Mouse mAb [A6B9-R] - BSA and Azide free	HA610092	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	RPS6 Recombinant Mouse mAb [A6B9-R]	HA601214	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	RPS6KA2 Mouse mAb [8C-4C]	EM1706-43	Human, Mouse	WB,IF-Cell,IHC-P,FC
	RRAGA Rabbit pAb	HA500373	Human, Mouse, Rat	WB
	RRAGD Rabbit pAb	ER61094	Human, Mouse	WB
	RSK1 p90 Recombinant Rabbit mAb [JE48-76]	ET7109-83	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
KD	RSK2 Recombinant Rabbit mAb [JE63-19]	HA721033	Human, Mouse, Rat	WB,IF-Cell,FC
	Rsk-4 Rabbit pAb	ER64892	Human	WB,IHC-P
	RSK4 Recombinant Rabbit mAb [PSH02-04] - BSA and Azide free	HA750774	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	RSK4 Recombinant Rabbit mAb [PSH02-04]	HA721763	Human, Mouse, Rat	WB,IHC-P,IF-Cell,FC
	SESN2 / Sestrin-2 Recombinant Rabbit mAb [PSH05-78]	HA722489	Human, Mouse, Rat	WB
	SESN2 / Sestrin-2 Recombinant Rabbit mAb [PSH05-78] - BSA and Azide free	HA751038	Human, Mouse, Rat	WB
	SGK1 Recombinant Rabbit mAb [SC05-71]	ET1610-19	Human, Mouse, Rat	WB,IHC-P,IP,FC,IF-Tissue
	SIN1 Mouse mAb [A6G10]	HA601008	Human	WB
	SKP2 Recombinant Rabbit mAb [PSH06-22]	HA722654	Human, Monkey	WB
	SKP2 Rabbit pAb	HA500427	Human	WB,IHC-P,IF-Cell
	SKP2 Rabbit pAb	ER1916-62	Human, Mouse, Rat	WB,IHC-P,IF-Cell,IF-Tissue
	SKP2 Recombinant Rabbit mAb [PSH06-22] - BSA and Azide free	HA751060	Human, Monkey	WB
	SOS1 Recombinant Rabbit mAb [JE30-06]	HA720077	Human, Mouse, Rat	WB,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Sos 2 Rabbit pAb	ER64998	Human, Mouse	WB,IHC-P
	Strad Rabbit pAb	ER65049	Human, Mouse	WB,IHC-P
	TNF alpha Rabbit pAb	ER1919-22	Human, Mouse	WB
Cell treatment	TNF alpha Rabbit pAb	R1203-1	Mouse	WB,IF-Cell,ELISA
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07] - BSA and Azide free	HA750897	Human, Mouse	WB
Cell treatment	TNF alpha Recombinant Rabbit mAb [PS01-07]	HA722022	Human, Mouse	WB
	TNF alpha Rabbit pAb	ER65189	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	TNFR1 Mouse mAb [2G2]	RT1624	Human, Mouse, Rat	WB,IP,IF,IHC-P
	TNF Receptor I Rabbit pAb	HA500140	Human	WB,IHC-P,FC
	TSC2/Tuberin Rabbit pAb	R1510-2	Human, Mouse	IF-Cell,IHC-P,FC
	Tuberin Recombinant Rabbit mAb [SC05-59]	ET1610-10	Human, Mouse, Rat	WB,IHC-P,FC
	ULK1 Recombinant Rabbit mAb [JA58-36]	ET1704-63	Human, Mouse, Rat	WB,IHC-P
	ULK1 Recombinant Rabbit mAb [JA58-36] - BSA and Azide free	HA750418	Human, Mouse, Rat	WB,IHC-P
	ULK2 Mouse mAb	EM1707-69	Human, Mouse, Rat	WB,FC
	VATC1 Rabbit pAb	ER61392	Human, Rat, Mouse	WB
	VATE2 Rabbit pAb	ER61394	Human, Mouse	WB
	VATF Rabbit pAb	ER61395	Human, Rat, Mouse	WB
	VATG2 Rabbit pAb	ER61396	Human, Mouse	WB
	VATG3 Rabbit pAb	ER61397	Human	WB
	V-ATPase B1 Rabbit pAb	ER65293	Human, Mouse	WB,IHC-P
	V-ATPase C2 Rabbit pAb	ER65294	Human	WB
	V-ATPase D Rabbit pAb	ER65295	Human, Mouse, Rat	WB
	V-ATPase H Rabbit pAb	ER65296	Human, Mouse	WB,IHC-P
	Wnt-1 Rabbit pAb	ER65317	Human, Mouse	WB,IHC-P,IF-Cell
	Wnt10a Rabbit pAb	HA500356	Human, Mouse	WB
	Wnt11 Rabbit pAb	ER1918-68	Human, Mouse	WB,IHC-P,IF-Tissue
	Wnt16 Rabbit pAb	HA500145	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	WNT2 Rabbit pAb	ER1511-4	Human, Mouse	WB,FC
	Wnt2b Recombinant Rabbit mAb [JB48-34]	ET7107-20	Human, Mouse, Rat	WB,IF-Cell,FC
	Wnt3a Recombinant Rabbit mAb [PSH19-96]	HA724116	Human, Mouse	WB,IHC-P,IP
Cell treatment	Wnt3a Rabbit pAb	HA500193	Human	WB,IF-Cell,IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	WNT4 Mouse mAb [4G1]	RT1662	Human, Mouse, Rat	WB,IP,IF,IHC-P
	Wnt5a Recombinant Rabbit mAb [JU30-30]	ET1706-33	Human	WB,IHC-P
	WNT5B Recombinant Rabbit mAb [JE53-67]	ET7110-34	Human	WB,FC
	Wnt6 Rabbit pAb	ER1918-70	Human, Mouse	WB,IHC-P,IF-Tissue
	WNT7A Rabbit pAb	ER1918-71	Human, Mouse, Rat	WB
	WNT7B Rabbit pAb	ER60344	Human, Mouse	WB
	WNT8A Recombinant Rabbit mAb [PSH0-39] - BSA and Azide free	HA750617	Human, Mouse, Rat	WB,IF-Cell,FC
	WNT8A Recombinant Rabbit mAb [PSH0-39]	HA721341	Human, Mouse, Rat	WB,IF-Cell,FC
	WNT8B Rabbit pAb	ER61525	Human, Mouse	WB
	WNT9A Rabbit pAb	ER60345	Human, Mouse	WB
	WNT9B Rabbit pAb	HA500412	Human, Mouse	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: sales@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层