

2026年2月—新上市产品清单

高级验证	产品名	货号	反应种属	验证应用
	AGO4 Rabbit pAb	HA500617	Human	WB
	Angiotensinogen Recombinant Rabbit mAb [PSH21-83] - BSA and Azide free	HA751811	Human, Mouse, Rat	WB,IHC-P
	Angiotensinogen Recombinant Rabbit mAb [PSH21-83]	HA724254	Human, Mouse, Rat	WB,IHC-P
	Anti-DYKDDDDK (Flag) Affinity Gel	HAK21024		IP,Co-IP
	Anti-Human IgG(CH1 Fragment specific), AlpSdAbs® VHH(HRP)	023-105-005	Human	WB,ELISA
	Mouse IgG1 (isotype control) Magnetic Beads	HAK21190		IP
Relative expression	Apolipoprotein A1 Recombinant Rabbit mAb [PSH21-76] - BSA and Azide free	HA751809	Human	WB,IF-Cell,IHC-P,FC
	Apolipoprotein A1 Recombinant Rabbit mAb [PSH21-75] - BSA and Azide free	HA751808	Human, Mouse, Rat	WB,IF-Cell,IHC-P
	Aurora A Recombinant Rabbit mAb [PSH22-00] - BSA and Azide free	HA751816	Human	WB,IF-Cell,FC,IP
	Aurora A Recombinant Rabbit mAb [PSH22-00]	HA724269	Human	WB,IF-Cell,FC,IP
Cell treatment	BARD1 Recombinant Rabbit mAb [PSH22-30]	HA724291	Human	WB,IP
Relative expression	Bassoon / BSN Recombinant Rabbit mAb [PSH21-32] - BSA and Azide free	HA751795	Mouse, Rat	IHC-Fr,IHC-P
	BST1 Rabbit pAb	HA500620	Human	WB
	CD130 / GP130 Recombinant Rabbit mAb [PSH22-56]	HA724306	Mouse	WB
Relative expression	CD25 / IL-2 Receptor alpha Recombinant Rabbit mAb [PSH21-79] - BSA and Azide free	HA751810	Mouse	WB,IHC-P
	CD38 Rabbit pAb	HA500613	Mouse, Rat	WB,IF-Cell,FC
Cell treatment	CXCL5 Recombinant Rabbit mAb [PSH22-57]	HA724307	Human	WB,IF-Cell,FC
Cell treatment	CXCL5 Recombinant Rabbit mAb [PSH22-58]	HA724308	Human	WB,IF-Cell,FC
	CXCR2 Recombinant Rabbit mAb [PSH21-28] - BSA and Azide free	HA751793	Human, Mouse, Rat	WB
	Cynomolgus IFN-gamma, Tag Free Protein	HA211458	Monkey	

Relative expression	Cytokeratin 7 Recombinant Antibody [ST50-05] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610395	Human	WB,IF-Cell,IHC-P,FC
	Cytokeratin (CAM 5.2) Recombinant Mouse mAb [PSH22-15]	HA601586	Human	IHC-P
	EDC4 Recombinant Rabbit mAb [PSH22-01]	HA724270	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
	EDC4 Recombinant Rabbit mAb [PSH22-01] - BSA and Azide free	HA751817	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC
Relative expression	EGFR Recombinant Rabbit mAb [PSH21-66] - BSA and Azide free	HA751805	Human	WB,IF-Cell,FC
	FAM48A / p38IP Rabbit pAb	HA500624	Human, Mouse, Rat	IHC-P
	FPR2 Recombinant Rabbit mAb [PSH21-63] - BSA and Azide free	HA751802	Human, Mouse, Rat	WB,IF-Cell,FC
	HCST Rabbit pAb	HA500615	Rat	WB
	HRP Conjugated Histone H3 Recombinant Rabbit mAb [JJ090-07] - Loading control	HA724276H	Human, Mouse, Rat, Monkey	WB
	HNF-1B Rabbit pAb	HA500622	Human	WB,IF-Cell,FC
	HTR3A Rabbit pAb	HA500614	Mouse, Rat	WB
	Human 14-3-3 gamma/YWHAG, Tag Free Protein	HA211350	Human	
	Human Angiopoietin-2/ANGPT2, C-His Tag Protein	HA211354	Human	
	Human Angiotensinogen Recombinant Rabbit mAb [PSH22-02] - BSA and Azide free (Capture)	HA725245	Human	ELISA(Cap)
	Human Angiotensinogen Recombinant Rabbit mAb [PSH22-03] - BSA and Azide free (Detector)	HA725246	Human	ELISA(Det)
	Biotin Conjugated Human Angiotensinogen Recombinant Rabbit mAb [PSH22-03] - Detector	HA725246B	Human	ELISA(Det),ELISA
	Human Antithrombin III/AT III enzyme-linked immunoassay kit(one step)	EHY0196	Human	ELISA
	Human Apolipoprotein M/APOM, Tag Free	HA211307	Human	
	Human Beta ig-h3/BIGH3/TGFBI, Tag Free Protein	HA211459	Human	
	Human Caspase-1, Tag Free Protein	HA211426	Human	
	Human Cathepsin S/CTSS, Tag Free Protein	HA211465	Human	
	HRP Conjugated Human CCL16 Recombinant Rabbit mAb [PSH19-91] - Detector	HA725282H	Human	ELISA(Det),ELISA
	Human CDCP1, Tag Free Protein	HA211462	Human	

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Chemerin/RARRES2, Tag Free Protein	HA211352	Human	
	Biotin Conjugated Human Complement factor H Recombinant Rabbit mAb [PSH22-21] - Detector	HA724282B	Human	ELISA(Det),ELISA
	Human Complement factor H Recombinant Rabbit mAb [PSH22-21] - BSA and Azide free (Detector)	HA724282	Human	ELISA(Det)
	Human Complement factor H Recombinant Rabbit mAb [PSH22-20] - BSA and Azide free (Capture)	HA724281	Human	ELISA(Cap)
	Biotin Conjugated Human C-peptide Recombinant Rabbit mAb [PSH22-39] - Detector	HA725095B	Human	ELISA(Det),ELISA
	Human C-peptide Recombinant Rabbit mAb [PSH22-38] - BSA and Azide free (Capture)	HA725094	Human	ELISA(Cap)
	Human C-peptide Recombinant Rabbit mAb [PSH22-39] - BSA and Azide free (Detector)	HA725095	Human	ELISA(Det)
	Human Dopamine beta-hydroxylase/DBH, C-His Tag Protein	HA211471	Human	
	Human DR6 / TNFRSF21, C-His Tag Protein	HA211300	Human	
	Human FABP1/L-FABP, Tag Free Protein	HA211456	Human	
	Human FABP2 Recombinant Rabbit mAb [PSH22-41] - BSA and Azide free (Detector)	HA725372	Human	ELISA(Det)
	Biotin Conjugated Human FABP2 Recombinant Rabbit mAb [PSH22-41] - Detector	HA725372B	Human	ELISA(Det),ELISA
	Human FABP2 Recombinant Rabbit mAb [PSH22-40] - BSA and Azide free (Capture)	HA725371	Human	ELISA(Cap)
	Biotin Conjugated Human Fas / CD95 Recombinant Rabbit mAb [PSH16-39] - Detector	HA723837B	Human	ELISA(Det)
	Biotin Conjugated Human Fas / CD95 Recombinant Rabbit mAb [PSH16-37] - Detector	HA723835B	Human	ELISA(Det)
	Human FGF-21 Recombinant Rabbit mAb [PSH22-34] - BSA and Azide free (Capture)	HA725235	Human	ELISA(Cap)
	Biotin Conjugated Human FGF-21 Recombinant Rabbit mAb [PSH22-33] - Detector	HA725236B	Human	ELISA(Det),ELISA
	Human FGF-21 Recombinant Rabbit mAb [PSH22-33] - BSA and Azide free (Detector)	HA725236	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human GDI1, C-His Tag Protein	HA211295	Human	
	Human GFR alpha-like, Tag Free Protein	HA211457	Human	
	Human GRN Recombinant Rabbit mAb [PSH22-43] - BSA and Azide free (Detector)	HA725310	Human	ELISA(Det)
	Human GRN Recombinant Rabbit mAb [PSH22-42] - BSA and Azide free (Capture)	HA725309	Human	ELISA(Cap)
	Biotin Conjugated Human GRN Recombinant Rabbit mAb [PSH22-43] - Detector	HA725310B	Human	ELISA(Det),ELISA
	Biotin Conjugated Human IGFBP6 Recombinant Rabbit mAb [PSH22-19] - Detector	HA724280B	Human	ELISA(Det),ELISA
	Human IGFBP6 Recombinant Rabbit mAb [PSH22-18] - BSA and Azide free (Capture)	HA724279	Human	ELISA(Cap)
	Human IGFBP6 Recombinant Rabbit mAb [PSH22-19] - BSA and Azide free (Detector)	HA724280	Human	ELISA(Det)
	Human IGF-I R/IGF1R, C-His Tag Protein	HA211349	Human	
	Human IL-11, Tag Free Protein	HA211299	Human	
	Human IL-18 R1 / CD218a, C-His Tag Protein	HA211256	Human	
	Human IL-1 alpha Recombinant Rabbit mAb [PSH22-28] - BSA and Azide free (Capture)	HA724289	Human	ELISA(Cap)
	Biotin Conjugated Human IL-1 alpha Recombinant Rabbit mAb [PSH22-27] - Detector	HA724288B	Human	ELISA(Det)
	Human IL-1 alpha Recombinant Rabbit mAb [PSH22-26] - BSA and Azide free (Capture)	HA724287	Human	ELISA(Cap)
	Human IL-1 alpha Recombinant Rabbit mAb [PSH22-27] - BSA and Azide free (Detector)	HA724288	Human	ELISA(Det)
	Human Kininogen 1, C-His Tag Protein	HA211254	Human	
	Human Lipoprotein A/LPA enzyme-linked immunoassay kit(one step)	EHY0171	Human	ELISA
	HRP Conjugated Human/Monkey Albumin Recombinant Rabbit mAb [PSH20-17] - Detector	HA724128H	Human, Dog, Cynomolgus monkey	ELISA(Det),ELISA
	Human/Monkey Coagulation factor VIII Recombinant Rabbit mAb [PSH21-69] - BSA and Azide free (Capture)	HA724240	Human, Cynomolgus monkey	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human/Monkey Coagulation factor VIII Recombinant Rabbit mAb [PSH21-70] - Detector	HA724241B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human/Monkey Coagulation factor VIII Recombinant Rabbit mAb [PSH21-70] - BSA and Azide free (Detector)	HA724241	Human, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human/Monkey VEGFR3 / FLT4 Recombinant Rabbit mAb [PSH22-17] - Detector	HA724278B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human/Monkey VEGFR3 / FLT4 Recombinant Rabbit mAb [PSH22-17] - BSA and Azide free (Detector)	HA724278	Human, Cynomolgus monkey	ELISA(Det)
	Human/Monkey VEGFR3 / FLT4 Recombinant Rabbit mAb [PSH22-16] - BSA and Azide free (Capture)	HA724277	Human, Cynomolgus monkey	ELISA(Cap)
	Human/Mouse Hemoglobin subunit alpha Recombinant Rabbit mAb [PSH22-08] - BSA and Azide free (Capture)	HA725427	Human, Mouse	ELISA(Cap)
	Human/Mouse Hemoglobin subunit alpha Recombinant Rabbit mAb [PSH22-09] - BSA and Azide free (Detector)	HA725428	Human, Mouse	ELISA(Det)
	Biotin Conjugated Human/Mouse Hemoglobin subunit alpha Recombinant Rabbit mAb [PSH22-09] - Detector	HA725428B	Human, Mouse	ELISA(Det),ELISA
	Human/Mouse IL-20 Recombinant Rabbit mAb [PSH21-86] - BSA and Azide free (Detector)	HA725274	Human, Mouse	ELISA(Det)
	Biotin Conjugated Human/Mouse IL-20 Recombinant Rabbit mAb [PSH21-86] - Detector	HA725274B	Human	ELISA(Det),ELISA
	Human/Mouse IL-20 Recombinant Rabbit mAb [PSH21-85] - BSA and Azide free (Capture)	HA725273	Human, Mouse	ELISA(Cap)
	Biotin Conjugated Human/Mouse MMP-2 Recombinant Rabbit mAb [PSH22-05] - Detector	HA725380B	Human, Mouse	ELISA(Det),ELISA
	Human/Mouse MMP-2 Recombinant Rabbit mAb [PSH22-04] - BSA and Azide free (Capture)	HA725379	Human, Mouse	ELISA(Cap)
	Human/Mouse MMP-2 Recombinant Rabbit mAb [PSH22-05] - BSA and Azide free (Detector)	HA725380	Human, Mouse	ELISA(Det)
	Human Myeloperoxidase enzyme-linked immunoassay kit(one step)	EHY0205	Human	ELISA
	Human Myoglobin, Tag Free Protein	HA211453	Human	
	Human Neutrophil Elastase Recombinant Rabbit mAb [PSH22-07] - BSA and Azide free (Detector)	HA725418	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Neutrophil Elastase Recombinant Rabbit mAb [PSH22-06] - BSA and Azide free (Capture)	HA725417	Human	ELISA(Cap)
	Biotin Conjugated Human Neutrophil Elastase Recombinant Rabbit mAb [PSH22-07] - Detector	HA725418B	Human	ELISA(Det),ELISA
	Human S100A8 & S100A9 Heterodimer, N-Strep II, C-His Tag Protein	HA211445	Human	
	Human Tissue Factor Recombinant Rabbit mAb [PSH22-45] - BSA and Azide free (Capture)	HA725401	Human	ELISA(Cap)
	Biotin Conjugated Human Tissue Factor Recombinant Rabbit mAb [PSH22-46] - Detector	HA725402B	Human	ELISA(Det),ELISA
	Human Tissue Factor Recombinant Rabbit mAb [PSH22-46] - BSA and Azide free (Detector)	HA725402	Human	ELISA(Det)
	Human TLR3 Recombinant Rabbit mAb [PSH22-25] - BSA and Azide free (Detector)	HA724286	Human	ELISA(Det)
	Biotin Conjugated Human TLR3 Recombinant Rabbit mAb [PSH22-25] - Detector	HA724286B	Human	ELISA(Det),ELISA
	Human TLR3 Recombinant Rabbit mAb [PSH22-24] - BSA and Azide free (Capture)	HA724285	Human	ELISA(Cap)
	Human VSTM1, C-His Tag Protein	HA211255	Human	
	IGF1 Recombinant Rabbit mAb [PSH21-44] - BSA and Azide free	HA751799	Mouse, Rat	IF-Cell,IHC-P
	IRS2 Recombinant Rabbit mAb [PSH22-29]	HA724290	Human, Mouse, Rat	IHC-P,IF-Cell,FC
Cell treatment, Relative expression	ISG15 Mouse mAb	RT1674	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC
	METTL9 Rabbit pAb	HA500618	Human, Mouse	WB
Relative expression, Cell treatment	MMP-9 Recombinant Rabbit mAb [PSH19-80] - BSA and Azide free	HA751790	Mouse, Rat	WB,IHC-P
	Mouse Albumin Recombinant Rabbit mAb [PSH21-67] - BSA and Azide free (Capture)	HA724238	Mouse	ELISA(Cap)
	Mouse Albumin Recombinant Rabbit mAb [PSH21-68] - BSA and Azide free (Detector)	HA724239	Mouse	ELISA(Det)
	Biotin Conjugated Mouse Albumin Recombinant Rabbit mAb [PSH21-68] - Detector	HA724239B	Mouse	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Biotin Conjugated Mouse CD200 / OX2 Recombinant Rabbit mAb [PSH21-91] - Detector	HA724260B	Mouse	ELISA(Det),ELISA
Mouse CD200 / OX2 Recombinant Rabbit mAb [PSH21-91] - BSA and Azide free (Detector)	HA724260	Mouse	ELISA(Det)
Mouse CD200 / OX2 Recombinant Rabbit mAb [PSH21-90] - BSA and Azide free (Capture)	HA724259	Mouse	ELISA(Cap)
Mouse CD38 Recombinant Rabbit mAb [PSH21-80] - BSA and Azide free (Capture)	HA724251	Mouse	ELISA(Cap)
Biotin Conjugated Mouse CD80 Recombinant Rabbit mAb [PSH21-95] - Detector	HA724264B	Mouse	ELISA(Det),ELISA
Mouse CD80 Recombinant Rabbit mAb [PSH21-95] - BSA and Azide free (Detector)	HA724264	Mouse	ELISA(Det)
Mouse CD80 Recombinant Rabbit mAb [PSH21-94] - BSA and Azide free (Capture)	HA724263	Mouse	ELISA(Cap)
Mouse FGF-21 enzyme-linked immunoassay kit	EM0053	Mouse	ELISA
Biotin Conjugated Mouse/Rat CD38 Recombinant Rabbit mAb [PSH21-81] - Detector	HA724252B	Mouse, Rat	ELISA(Det),ELISA
Mouse/Rat CD38 Recombinant Rabbit mAb [PSH21-82] - BSA and Azide free (Capture)	HA724253	Mouse, Rat	ELISA(Cap)
Mouse/Rat CD38 Recombinant Rabbit mAb [PSH21-81] - BSA and Azide free (Detector)	HA724252	Mouse, Rat	ELISA(Det)
Mouse/Rat IGFBP6 Recombinant Rabbit mAb [PSH22-22] - BSA and Azide free (Capture)	HA724283	Mouse, Rat	ELISA(Cap)
Mouse/Rat IGFBP6 Recombinant Rabbit mAb [PSH22-23] - BSA and Azide free (Detector)	HA724284	Mouse, Rat	ELISA(Det)
Biotin Conjugated Mouse/Rat IGFBP6 Recombinant Rabbit mAb [PSH22-23] - Detector	HA724284B	Mouse, Rat	ELISA(Det),ELISA
Mouse/Rat NF-H Recombinant Rabbit mAb [PSH21-71] - BSA and Azide free (Capture)	HA724242	Mouse, Rat	ELISA(Cap)
Biotin Conjugated Mouse/Rat NF-H Recombinant Rabbit mAb [PSH21-72] - Detector	HA724243B	Mouse, Rat	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse/Rat NF-H Recombinant Rabbit mAb [PSH21-72] - BSA and Azide free (Detector)	HA724243	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse/Rat P-Cadherin Recombinant Rabbit mAb [PSH21-93] - Detector	HA724262B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat P-Cadherin Recombinant Rabbit mAb [PSH21-93] - BSA and Azide free (Detector)	HA724262	Mouse, Rat	ELISA(Det)
	Mouse/Rat P-Cadherin Recombinant Rabbit mAb [PSH21-92] - BSA and Azide free (Capture)	HA724261	Mouse, Rat	ELISA(Cap)
	Mouse/Rat VEGFA Recombinant Rabbit mAb [PSH21-87] - BSA and Azide free (Capture)	HA724256	Mouse, Rat	ELISA(Cap)
	Mouse/Rat VEGFA Recombinant Rabbit mAb [PSH21-89] - BSA and Azide free (Detector)	HA724258	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse/Rat VEGFA Recombinant Rabbit mAb [PSH21-88] - Detector	HA724257B	Mouse, Rat	ELISA(Det),ELISA
	Biotin Conjugated Mouse/Rat VEGFA Recombinant Rabbit mAb [PSH21-89] - Detector	HA724258B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat VEGFA Recombinant Rabbit mAb [PSH21-88] - BSA and Azide free (Detector)	HA724257	Mouse, Rat	ELISA(Det)
	MRPS27 Rabbit pAb	HA500623	Human, Mouse, Rat	WB
	MTA1 Recombinant Rabbit mAb [PSH21-25] - BSA and Azide free	HA751791	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP
Relative expression	MTBR-TAU 243 Isoform 2N4R(236-250) Recombinant Rabbit mAb [PSH22-59]	HA724309	Mouse, Rat, Pig, Human	WB,IP,FC
	Myelin Basic Protein Recombinant Antibody [JF0943] - Chicken IgY (Chimeric) - BSA and Azide free	HA610392	Human, Mouse, Rat	IHC-Fr,IF-Tissue
Relative expression, Cell treatment	Myogenin Recombinant Mouse mAb [PSH22-35]	HA601587	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr
Relative expression	Nav1.6 / SCN8A Recombinant Rabbit mAb [PSH22-36]	HA724294	Mouse, Rat	WB,IHC-P,IP
	NCR1 Recombinant Rabbit mAb [PSH21-26] - BSA and Azide free	HA751792	Mouse	IHC-P,FC
	NDUFA8 Recombinant Rabbit mAb [PSH21-96]	HA724265	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC,IP
	NDUFA8 Recombinant Rabbit mAb [PSH21-96] - BSA and Azide free	HA751812	Human, Mouse, Rat, Monkey	WB,IHC-P,IF-Cell,FC,IP
	NDUFV1 Recombinant Rabbit mAb [PSH21-97]	HA724266	Human, Mouse, Rat, Monkey, Pig	WB,IF-Cell,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	NDUFV1 Recombinant Rabbit mAb [PSH21-97] - BSA and Azide free	HA751813	Human, Mouse, Rat, Monkey, Pig	WB,IF-Cell,IHC-P,IP
	NF-H Recombinant Antibody [JF99-07] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610394	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	NF-H Recombinant Antibody [JF99-07] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610393	Human, Mouse, Rat	IHC-P
	NLRX1 Rabbit pAb	HA500621	Human, Rat	WB
Relative expression	NPTX2 Recombinant Rabbit mAb [PSH21-53] - BSA and Azide free	HA751800	Human, Mouse, Rat	WB,IHC-P
Relative expression, Cell treatment	p21 Recombinant Rabbit mAb [PSH22-37]	HA724295	Human	WB,IF-Cell,IHC-P,IP
Relative expression	P2X7 Recombinant Rabbit mAb [PSH22-44]	HA724296	Mouse	WB
Relative expression	PCSK9 Recombinant Rabbit mAb [PSH21-65] - BSA and Azide free	HA751804	Human, Mouse	WB,FC,IP
	PDGFR beta Recombinant Antibody [SY10-08] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610390	Human, Mouse, Rat	WB,IHC-Fr,IHC-P
	PDGFR beta Recombinant Antibody [SY10-08] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610391	Human, Mouse, Rat	WB,IHC-Fr,IHC-P
Cell treatment	Phospho-c-Fos (S32) Recombinant Rabbit mAb [PSH21-99] - BSA and Azide free	HA751815	Human	WB,IF-Cell
Cell treatment	Phospho-c-Fos (S32) Recombinant Rabbit mAb [PSH21-99]	HA724268	Human	WB,IF-Cell
Cell treatment	Phospho-LATS1 (T1079) Recombinant Rabbit mAb [PSH21-98]	HA724267	Human, Mouse	WB,IP
Cell treatment	Phospho-LATS1 (T1079) Recombinant Rabbit mAb [PSH21-98] - BSA and Azide free	HA751814	Human, Mouse	WB,IP
	PINK1 Recombinant Mouse mAb [PSH22-50]	HA601588	Human, Mouse, Rat	WB
	PMS2 Recombinant Rabbit mAb [PSH22-60]	HA724310	Human	IHC-P
	Rat IL-12 (linked heterodimer), C-His Tag Protein	HA211345	Rat	
	Rat IL-12 p40/IL12B, C-His Tag Protein	HA211346	Rat	
	Rat IL-3, C-His Tag Protein	HA211455	Rat	
	Rat IL-5, Tag Free Protein	HA211454	Rat	
	Rat NFkB1/p105, Tag Free Protein	HA211348	Rat	
	Rat Prolactin, Tag Free Protein	HA211351	Rat	
	Rat SOD1 / Cu-Zn SOD, Tag Free Protein	HA211309	Rat	

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	SAA1 + SAA2 Recombinant Rabbit mAb [PSH21-34] - BSA and Azide free	HA751797	Mouse	WB,IHC-P
	SAA1 + SAA2 Recombinant Rabbit mAb [PSH21-33] - BSA and Azide free	HA751796	Mouse	WB,IHC-Fr
Relative expression	SATB2 Recombinant Rabbit mAb [PSH21-84]	HA724255	Human	IHC-P
	SEL1L Recombinant Rabbit mAb [PSH22-54]	HA724304	Human, Mouse, Rat	WB,IF-Cell,FC
Relative expression	SP-D / SFTPD Recombinant Rabbit mAb [PSH21-62] - BSA and Azide free	HA751801	Mouse, Rat	WB,IHC-P,IHC-Fr,IP,IF-Tissue
	SSBP1 Recombinant Rabbit mAb [PSH21-31] - BSA and Azide free	HA751794	Human, Mouse, Rat, Monkey	WB,IHC-P
	TATA binding protein TBP Recombinant Rabbit mAb [PSH22-13]	HA724274	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC,IP
Relative expression	TIMP1 Recombinant Rabbit mAb [PSH21-35] - BSA and Azide free	HA751798	Mouse	WB
	TIMP1 Recombinant Rabbit mAb [PSH21-64] - BSA and Azide free	HA751803	Human, Monkey	WB,IHC-P
	TMEM106C Rabbit pAb	HA500625	Human	IHC-P
Relative expression	TROP2 Recombinant Rabbit mAb [PSH22-14]	HA724275	Human, Mouse, Rat	WB,IF-Cell,IHC-Fr,IHC-P,IF-Tissue
Relative expression	TRP2 / DCT Recombinant Rabbit mAb [PSH22-31]	HA724292	Human	WB,IP
	UMOD Recombinant Rabbit mAb [PSH21-74] - BSA and Azide free	HA751807	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	UMPS Rabbit pAb	HA500616	Human, Mouse	WB,IF-Cell,IHC-P,FC
Relative expression	VAV1 Recombinant Rabbit mAb [PSH21-73] - BSA and Azide free	HA751806	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	ZnT1 Rabbit pAb	HA500619	Mouse, Rat	IHC-P

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层