

一步法ELISA Kit产品清单

产品名	货号	反应种属	灵敏度
Human Angiogenin enzyme-linked immunoassay kit(one step)	EHY0190	Human	0.01 ng/ml (50 μ l) ; 0.07 ng/ml (10 μ l) ;
Human Antithrombin III/AT III enzyme-linked immunoassay kit(one step)	EHY0196	Human	0.32 ng/ml (10 μ l)
Human CD115/M-CSF R/CSF-1-R enzyme-linked immunoassay kit (one step)	EHY0168	Human	0.01 ng/ml (50 μ l) ; 0.08 ng/ml (10 μ l) ;
Human CD155/PVR(Poliiovirus Receptor) enzyme-linked immunoassay kit(one step)	EHY0155	Human	0.02 ng/ml (50 μ l) ; 0.11 ng/ml (10 μ l) ;
Human CD273/B7-DC/PD-L2 enzyme-linked immunoassay kit (one step)	EHY0132	Human	0.01 ng/ml (50 μ l) ; 0.10 ng/ml (10 μ l) ;
Human Complement C3 enzyme-linked immunoassay kit(one step)	EHY0202	Human	0.09 pg/ml (10 μ l)
Human Complement C5 enzyme-linked immunoassay kit (one step)	EHY0108	Human	0.01 ng/ml (50 μ l) ; 0.18 ng/ml (10 μ l)
Human CTGF/CCN2 enzyme-linked immunoassay kit(one step)	EHY0156	Human	0.07 ng/ml (50 μ l) ; 0.35 ng/ml (10 μ l) ;
Human ICAM-1/CD54 enzyme-linked immunoassay kit (one step)	EHY0064	Human	1.22 pg/ml (50 μ l) ; 0.05 ng/ml (10 μ l)
Human IGFBP3/IBP3 enzyme-linked immunoassay kit(one step)	EHY0157	Human	0.004 ng/ml (50 μ l) ; 13.60 ng/ml (10 μ l) ;
Human IL-6R enzyme-linked immunoassay kit(one step)	EHY0207	Human	0.05 ng/ml (50 μ l) ; 0.56 ng/ml (10 μ l) ;
Human Lipoprotein A/LPA enzyme-linked immunoassay kit(one step)	EHY0171	Human	0.08 ng/ml (10 μ l) ;
Human MERTK enzyme-linked immunoassay kit(one step)	EHY0189	Human	0.03 ng/ml (50 μ l) ; 0.41 ng/ml (10 μ l) ;
Human Mesothelin enzyme-linked immunoassay kit(one step)	EHY0152	Human	0.16 ng/ml (50 μ l) ; 0.23 ng/ml (10 μ l) ;
Human/Monkey/Canine CCL5/RANTES enzyme-linked immunoassay kit(one step)	EHY0027	Human, Monkey, Canine	0.7 pg/ml (50 μ l) ; 4.99 pg/ml (10 μ l) ;
Human/Monkey Clusterin enzyme-linked immunoassay kit(one step)	EHY0073	Human, Monkey	0.14 ng/ml (10 μ l)
Human/Monkey Lipocalin-2/NGAL enzyme-linked immunoassay kit(one step)	EHY0166	Human, Monkey	0.01 ng/ml (50 μ l) ; 0.03 ng/ml (10 μ l) ;
Human/Monkey YKL-40/CHI3L1 enzyme-linked immunoassay kit(one step)	EHY0159	Human, Monkey	0.01 ng/ml (50 μ l) ; 0.05 ng/ml (10 μ l) ;
Human Myeloperoxidase enzyme-linked immunoassay kit(one step)	EHY0205	Human	0.02 ng/ml (50 μ l) ; 0.07 ng/ml (10 μ l)

Human Neuropilin-1/NRP1 enzyme-linked immunoassay kit (one step)	EHY0175	Human	0.02 ng/ml (50 µl) ; 0.03 ng/ml (10 µl) ;
Human NG2/MCSP/CSPG4 enzyme-linked immunoassay kit(one step)	EHY0162	Human	0.004 ng/ml (50 µl) ; 0.21 ng/ml (10 µl) ;
Human PAI-1/Serpin E1 enzyme-linked immunoassay kit(one step)	EHY0160	Human	2.18 pg/ml (50 µl) ; 40.00 pg/ml (10 µl)
Human Proprotein Convertase 9/PCSK9 enzyme-linked immunoassay kit (one step)	EHY0071	Human	0.39 ng/ml (50 µl) ; 0.73 ng/ml (10 µl)
Human RAGE enzyme-linked immunoassay kit(one step)	EHY0222	Human	21.39 pg/ml (50 µl) ; 65.37 pg/ml (10 µl) ;
Human RBP4 enzyme-linked immunoassay kit(one step)	EHY0212	Human	0.002 ng/ml (50 µl) ; 0.03 ng/ml (10 µl) ;
Human TIM4 enzyme-linked immunoassay kit(one step)	EHY0194	Human	0.001 ng/ml (50 µl) ; 0.03 ng/ml (10 µl) ;
Human TNF RI/TNFRSF1A enzyme-linked immunoassay kit (one step)	EHY0128	Human	0.01 ng/ml (50 µl) ; 0.13 ng/ml (10 µl)
Human tPA enzyme-linked immunoassay kit(one step)	EHY0158	Human	0.03 ng/ml (50 µl) ; 0.10 ng/ml (10 µl) ;
Human TREM2 enzyme-linked immunoassay kit (one step)	EHY0122	Human	0.004 ng/ml (50 µl) ; 0.05 ng/ml (50 µl)
Human VCAM-1/CD106 enzyme-linked immunoassay kit(one step)	EHY0067	Human	0.04 ng/ml (50 µl) ; 0.08 ng/ml (10 µl) ;
Human von Willebrand Factor enzyme-linked immunoassay kit (one step)	EHY0138	Human	0.02 ng/ml (50 µl) ; 0.69 ng/ml (10 µl) ;
Lentivirus HIV-1 p24 enzyme-linked immunoassay kit (one step)	EHY0078	Human	0.12 ng/ml (10 µl)
Mouse alpha 1 antitrypsin enzyme-linked immunoassay kit (one step)	EMY0052	Mouse	1.57 ng/ml (10 µl)
Mouse Coagulation Factor XII enzyme-linked immunoassay kit(one step)	EMY0001	Mouse	0.22 ng/ml (10 µl)
Mouse IL-27 P28/IL-30 enzyme-linked immunoassay kit(one step)	EM0051	Mouse	2.88 pg/ml (10 µl)
Mouse Pro-collagen I alpha 1/COL1A1 enzyme-linked immunoassay kit(one step)	EMY0045	Mouse	0.01 pg/ml (10 µl)
Mouse/Rat IGF1 enzyme-linked immunoassay kit (one step)	EMY0024	Mouse, Rat	0.10 pg/ml (50 µl) ; 0.37 pg/ml (10 µl) ;
Rat CRP/C-Reactive Protein enzyme-linked immunoassay kit(one step)	ERY0009	Rat	2.9 pg/ml (50 µl) ; 42.15 pg/ml (10 µl)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层