

## 2026年3月—新上市产品清单

高级验证	产品名	货号	反应种属	验证应用
Relative expression	Adiponectin Recombinant Rabbit mAb [PSH22-55]	HA724305	Human, Mouse, Rat	WB,IHC-P
	AGRP Recombinant Rabbit mAb [PSH22-63]	HA724311	Mouse, Rat	IHC-Fr,IHC-P
	A-Myb Rabbit pAb	HA500631	Human	WB
	Argonaute Antibody Sampler Kit	HAK21191	Human	WB
	ATM Substrates Antibody Sampler Kit	HAK21187	Human	WB
Relative expression	BAP1 Recombinant Rabbit mAb [PSH22-53]	HA724303	Human, Mouse, Rat	WB,IP
	BCKDHB Rabbit pAb	HA500628	Human, Mouse, Rat	IHC-P
	BDNF Recombinant Rabbit mAb [JE59-35] - BSA and Azide free	HA751818	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-P,IHC-Fr,FC
	Brd4 Recombinant Rabbit mAb [JE04-04] - BSA and Azide free	HA751820	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Relative expression	Caspase-1 + p10 + p12 Recombinant Rabbit mAb [JE56-35] - BSA and Azide free	HA751822	Human, Mouse, Rat	WB,IF-Cell,IHC-P
KD, Relative expression	CD13 Recombinant Rabbit mAb [SN71-04] - BSA and Azide free	HA751849	Human, Mouse, Rat	WB,IHC-P,IP,IF-Tissue
Relative expression	CD1a Recombinant Rabbit mAb [JM21-33] - BSA and Azide free	HA751834	Human	WB,IHC-P,IF-Tissue
Relative expression	CD21 Recombinant Rabbit mAb [SC0681] - BSA and Azide free	HA751853	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP,FC,IF-Tissue
Relative expression	CD3 Recombinant Rabbit mAb [PSH22-52]	HA724302	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue,IF-Cell,FC,IP
Relative expression, KD	CD44 Recombinant Rabbit mAb [JE64-01] - BSA and Azide free	HA751819	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,FC
	CD45 Recombinant Rabbit mAb [JE03-05] - BSA and Azide free	HA751827	Human	WB,IHC-P,FC,IF-Tissue
Cell treatment	CD68 Recombinant Antibody [PSH06-65] - Guinea pig IgG2 (Chimeric)	HA601595	Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue
Cell treatment, Relative expression	CD68 Recombinant Antibody [PSH06-65] - Rat IgG1 (Chimeric)	HA601596	Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Cell,IF-Tissue
	Chaperone-Mediated Autophagy Antibody Sampler Kit	HAK21185	Human	WB,IHC-P

	Collagen III Recombinant Rabbit mAb [JE60-79] - BSA and Azide free	HA751826	Human, Mouse, Rat	WB,IF-Cell,FC
	Dnmt3b Recombinant Rabbit mAb [SY09-07] - BSA and Azide free	HA751851	Human, Rat	WB,IF-Cell,ChIP
	Biotin Conjugated Dxd Recombinant Mouse mAb [PSH0-37]	HA601134B	Species independent	ELISA
	EPS15 Recombinant Rabbit mAb [PSH23-03]	HA724343	Human, Mouse, Rat, Monkey	WB,IF-Cell,FC
	FAF1 Recombinant Rabbit mAb [PSH22-51]	HA724301	Human, Mouse, Rat, Monkey	WB,IF-Cell,IHC-P,FC,IP
	Fbx32 Recombinant Rabbit mAb [JE41-27] - BSA and Azide free	HA751828	Human, Mouse, Rat	WB,IF-Cell
	FGF2 Recombinant Rabbit mAb [PSH22-75]	HA724320	Human	WB,IF-Cell,FC
Relative expression	FGFR1 Recombinant Rabbit mAb [PSH22-64]	HA724312	Human, Mouse, Rat	WB,IHC-P
	FOXA2 Recombinant Rabbit mAb [JM10-64] - BSA and Azide free	HA751839	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IF-Cell,IHC-P,IHC-Fr,IF-Tissue
	FPR3 Rabbit pAb	HA500627	Human, Monkey	WB
	GABA B Receptor 2 Recombinant Rabbit mAb [JU31-32] - BSA and Azide free	HA751832	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	GABA Transporter 1 / GAT 1 Recombinant Rabbit mAb [JE34-61] - BSA and Azide free	HA751824	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	Glucagon Recombinant Rabbit mAb [JF33-10] - BSA and Azide free	HA751846	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P
	Glutathione metabolism Antibody Sampler Kit	HAK21181	Human	WB
	iFluor™ 800cwse Conjugated Goat anti-rabbit IgG polyclonal Antibody(For Odyssey system)	HA1152	Rabbit	WB
	Hemopexin Recombinant Rabbit mAb [PSH23-09]	HA724349	Human, Mouse, Rat	WB
	HIF-2 alpha Recombinant Rabbit mAb [JB24-42] - BSA and Azide free	HA751829	Human, Mouse, Rat	WB,IHC-P,FC,IF-Cell
	HRP Conjugated Human Adiponectin Recombinant Rabbit mAb [PSH05-32] - Detector	HA722595H	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human CD21/CR2, C-His Tag Protein	HA211498	Human	
	Human CD38, C-His Tag Protein	HA211487	Human	
	Human CD44, C-His Tag Protein	HA211442	Human	

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

HRP Conjugated Human Complement factor B Recombinant Rabbit mAb [PSH19-47] - Detector	HA725163H	Human	ELISA(Det),ELISA
Human C-peptide Protein	HA211473	Human	
Human Dermatotontin/DPT, Tag Free Protein	HA211485	Human	
Human Dystrophin/DMD, C-His Tag Protein	HA211478	Human	
Human EBI3/IL-27B, C-hFc Tag Protein	HA211495	Human	
Human FABP4/A-FABP, C-His Tag Protein	HA211489	Human	
Human FGF R3 (IIIc)/FGF R3C/CD333, C-His Tag Protein	HA211488	Human	
Human FLT3, Tag Free Protein	HA211437	Human	
Human HE4/WFDC2, Tag Free Protein	HA211484	Human	
Biotin Conjugated Human Hemopexin Recombinant Rabbit mAb [PSH22-71] - Detector	HA724315B	Human	ELISA(Det),ELISA
Human Hemopexin Recombinant Rabbit mAb [PSH22-71] - BSA and Azide free (Detector)	HA724315	Human	ELISA(Det)
Biotin Conjugated Human Hemopexin Recombinant Rabbit mAb [PSH22-72] - Detector	HA724316B	Human	ELISA(Det),ELISA
Human Hemopexin Recombinant Rabbit mAb [PSH22-70] - BSA and Azide free (Capture)	HA724314	Human	ELISA(Cap)
Human Hemopexin Recombinant Rabbit mAb [PSH22-72] - BSA and Azide free (Detector)	HA724316	Human	ELISA(Det)
Human IL-17C Recombinant Rabbit mAb [PSH22-62] - BSA and Azide free (Detector)	HA725270	Human	ELISA(Det)
Biotin Conjugated Human IL-17C Recombinant Rabbit mAb [PSH22-62] - Detector	HA725270B	Human	ELISA(Det),ELISA
Human IL-17C Recombinant Rabbit mAb [PSH22-61] - BSA and Azide free (Capture)	HA725269	Human	ELISA(Cap)
Biotin Conjugated Human IL-18R1 Recombinant Rabbit mAb [PSH22-86] - Detector	HA724331B	Human	ELISA(Det),ELISA
Human IL-18R1 Recombinant Rabbit mAb [PSH22-86] - BSA and Azide free (Detector)	HA724331	Human	ELISA(Det)

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-18R1 Recombinant Rabbit mAb [PSH22-85] - BSA and Azide free (Capture)	HA724330	Human	ELISA(Cap)
	Human IL-36 beta/IL-1 F8, Tag Free Protein	HA211356	Human	
	Human LECT2, Tag Free Protein	HA211492	Human	
	Human LRG1, C-His Tag Protein	HA211480	Human	
	HRP Conjugated Human Lumican Recombinant Rabbit mAb [PSH19-05] - Detector	HA725318H	Human	ELISA(Det),ELISA
	Human MASP2 Protein	HA210282	Human	
	HRP Conjugated Human Neutrophil Elastase Recombinant Rabbit mAb [PSH22-07] - Detector	HA725418H	Human	ELISA(Det),ELISA
	Human NT-proBNP, Tag Free Protein	HA211494	Human	
	Biotin Conjugated Human POMC Recombinant Rabbit mAb [PSH22-66] - Detector	HA725396B	Human	ELISA(Det),ELISA
	Human POMC Recombinant Rabbit mAb [PSH22-66] - BSA and Azide free (Detector)	HA725396	Human	ELISA(Det)
	Human POMC Recombinant Rabbit mAb [PSH22-65] - BSA and Azide free (Capture)	HA725395	Human	ELISA(Cap)
	Human Procalcitonin/CALCA, Tag Free Protein	HA211353	Human	
	Human Pro-Collagen I alpha 1/COL1A1, Tag Free Protein	HA211375	Human	
	Human Prorelaxin H2/RLN2, Tag Free Protein	HA211491	Human	
	Human REG4, C-His Tag Protein	HA211483	Human	
	Human Serpin A12, C-His Tag Protein	HA211482	Human	
	Human Serum Amyloid P-component/APCS, Tag Free Protein	HA211481	Human	
	Human TFF2, Tag Free Protein	HA211486	Human	
	Human Thrombomodulin, C-His Tag Protein	HA211324	Human	
	Human TIGIT, C-His Tag Protein	HA211400	Human	
	Human Uteroglobin, C-His Tag Protein	HA211490	Human	
	IgD Recombinant Rabbit mAb [JE39-23] - BSA and Azide free	HA751821	Human	WB,IHC-P,IF-Tissue
Relative expression	IL-6R Recombinant Rabbit mAb [PSH22-78]	HA724323	Human	WB,FC

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated INHBE Recombinant Rabbit mAb [PSH22-84] - Detector	HA724329B	Human, Mouse, Rat, Cynomolgus monkey	ELISA(Det),ELISA
	INHBE Recombinant Rabbit mAb [PSH22-83] - BSA and Azide free (Capture)	HA724328	Human, Mouse, Rat, Cynomolgus monkey	ELISA(Cap)
	INHBE Recombinant Rabbit mAb [PSH22-84] - BSA and Azide free (Detector)	HA724329	Human, Mouse, Rat, Cynomolgus monkey	ELISA(Det)
	INSIG1 Recombinant Rabbit mAb [PSH22-32]	HA724293	Human, Mouse, Rat	WB
Relative expression	L1CAM Recombinant Rabbit mAb [JM11-05] - BSA and Azide free	HA751842	Human	WB,IHC-P,IF-Tissue
	Laminin beta 1 Recombinant Rabbit mAb [JM099-03] - BSA and Azide free	HA751844	Human, Mouse, Rat	WB,IHC-P,FC
Cell treatment	LAMP1 Recombinant Antibody [PSH07-39] - Mouse IgG1 (Chimeric)	HA601597	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	LDHB Recombinant Rabbit mAb [PSH22-91]	HA724336	Human, Mouse, Rat, Monkey, Pig	WB,IF-Cell,FC,IP
Relative expression	LI Cadherin Recombinant Rabbit mAb [JE60-34] - BSA and Azide free	HA751823	Human, Mouse	WB,IHC-P,IF-Tissue,IF-Cell,FC
	Lysosomal Marker Antibody Sampler Kit	HAK21188	Human	WB,IHC-P
	LYVE1 Recombinant Rabbit mAb [PSH23-01]	HA724340	Mouse	WB,IHC-Fr,IHC-P
	MCP-3/CCL7 Recombinant Rabbit mAb [PSH23-04]	HA724344	Human, Mouse	WB
	Metabolic-Exosome Pathway Antibody Sampler Kit	HAK21186	Human	WB
	mGluR1 Recombinant Rabbit mAb [JM11-61] - BSA and Azide free	HA751838	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IHC-Fr
Relative expression	MLH1 Recombinant Antibody [PSH20-67] - Rabbit IgG (Chimeric)	HA724342	Human	IHC-P,IF-Cell,WB
Relative expression	MOG Recombinant Rabbit mAb [JM23-06] - BSA and Azide free	HA751833	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IF-Tissue,IHC-Fr
	Mouse CD163 Recombinant Rabbit mAb [PSH22-49] - BSA and Azide free (Detector)	HA724299	Mouse	ELISA(Det)
	Mouse CD163 Recombinant Rabbit mAb [PSH22-47] - BSA and Azide free (Capture)	HA724297	Mouse	ELISA(Cap)
	Mouse CD163 Recombinant Rabbit mAb [PSH22-48] - BSA and Azide free (Capture/Detector)	HA724298	Mouse	ELISA(Cap),ELISA(Det)
	Biotin Conjugated Mouse CD163 Recombinant Rabbit mAb [PSH22-48] - Detector	HA724298B	Mouse	ELISA(Det),ELISA

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse CD163 Recombinant Rabbit mAb [PSH22-49] - Detector	HA724299B	Mouse	ELISA(Det),ELISA
	Mouse FSH alpha & beta Heterodimer, Tag Free Protein	HA211423	Mouse	
	Mouse IL-16 Recombinant Rabbit mAb [PSH22-69] - BSA and Azide free (Detector)	HA725424	Mouse	ELISA(Det)
	Biotin Conjugated Mouse IL-16 Recombinant Rabbit mAb [PSH22-69] - Detector	HA725424B	Mouse	ELISA(Det),ELISA
	Mouse IL-16 Recombinant Rabbit mAb [PSH22-68] - BSA and Azide free (Capture)	HA725423	Mouse	ELISA(Cap)
	Mouse Osteocalcin Recombinant Rabbit mAb [PSH22-74] - BSA and Azide free (Detector)	HA724318	Mouse	ELISA(Det)
	Biotin Conjugated Mouse Osteocalcin Recombinant Rabbit mAb [PSH22-74] - Detector	HA724318B	Mouse	ELISA(Det),ELISA
	Mouse Osteocalcin Recombinant Rabbit mAb [PSH22-73] - BSA and Azide free (Capture)	HA724317	Mouse	ELISA(Cap)
	Mouse Osteocalcin/BGLAP (50-95) peptide Protein	HA211500	Mouse	
	Mouse/Rat IGFBP3 Recombinant Rabbit mAb [PSH22-89] - BSA and Azide free (Capture)	HA724334	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat IGFBP3 Recombinant Rabbit mAb [PSH22-90] - Detector	HA724335B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat IGFBP3 Recombinant Rabbit mAb [PSH22-90] - BSA and Azide free (Detector)	HA724335	Mouse, Rat	ELISA(Det)
	Mouse TNFR1/CD120a/TNFRSF1A, C-His Tag Protein	HA211493	Mouse	
	Mre11 Recombinant Rabbit mAb [JM11-18] - BSA and Azide free	HA751837	Human, Mouse, Rat	WB,IHC-P
	NeuroD1 Recombinant Rabbit mAb [JM11-10] - BSA and Azide free	HA751840	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IP,IHC-P,IHC-Fr
Relative expression	NF-M Recombinant Rabbit mAb [JM11-20] - BSA and Azide free	HA751854	Human, Mouse, Rat, Cynomolgus monkey	WB,IHC-P,IHC-Fr,IF-Tissue
	NK1.1 Recombinant Rabbit mAb [PSH22-99]	HA724338	Mouse	IHC-Fr,IF-Cell,FC
	NKG2D Rabbit pAb	HA500626	Human	WB
	NKX3.1 Recombinant Rabbit mAb [JE36-27] - BSA and Azide free	HA751825	Human	WB,IHC-P

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	PD-L1 Recombinant Rabbit mAb [PSH22-67]	HA724313	Human	IHC-P
KD	PFKM Recombinant Rabbit mAb [JU53-31] - BSA and Azide free	HA751831	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Cell treatment	Phospho-NF-kB p65 (S529) Recombinant Rabbit mAb [SP07-00] - BSA and Azide free	HA751843	Human	WB,IHC-P,FC,IP
	Plasminogen Rabbit pAb	HA500630	Human, Rat	WB
KD	Podoplanin Recombinant Rabbit mAb [JM31-21] - BSA and Azide free	HA751841	Human	WB,IHC-P
KD	PPP1R1A Recombinant Rabbit mAb [SN0756] - BSA and Azide free	HA751850	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,IP
	Rat Caspase-3, C-His Tag Protein	HA211376	Rat	
	Rat EGF, Tag Free Protein	HA211479	Rat	
	Rat LIF, C-His Tag Protein	HA211373	Rat	
Relative expression	RIP3 Recombinant Rabbit mAb [PSH23-11]	HA724351	Mouse	WB,IP
	S100A9 Recombinant Rabbit mAb [JF096-8] - BSA and Azide free	HA751845	Human	WB,IF-Tissue,IHC-P
	FITC Conjugated S100 alpha 6 Recombinant Rabbit mAb [JF0976]	HA724357F	Human, Mouse, Rat	IF-Cell,IF-Tissue
	SRSF2 Recombinant Rabbit mAb [PSH22-77]	HA724322	Human, Mouse, Rat	IHC-P,IF-Cell,FC
	STEAP2 Rabbit pAb	HA500629	Human	WB
Relative expression	Synaptophysin Recombinant Rabbit mAb [SJ26-85] - BSA and Azide free	HA751835	Human, Mouse, Rat, Cynomolgus monkey, Pig	WB,IHC-P,IHC-Fr,IF-Tissue
	TAG1 Recombinant Rabbit mAb [PSH22-87] - BSA and Azide free (Capture)	HA724332	Human, Mouse, Rat	ELISA(Cap)
	Biotin Conjugated TAG1 Recombinant Rabbit mAb [PSH22-88] - Detector	HA724333B	Human, Mouse, Rat	ELISA(Det),ELISA
	TAG1 Recombinant Rabbit mAb [PSH22-88] - BSA and Azide free (Detector)	HA724333	Human, Mouse, Rat	ELISA(Det)
KD	TAPA1/CD81 Recombinant Rabbit mAb [SN206-01] - BSA and Azide free	HA751848	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	T-bet / Tbx21 Recombinant Rabbit mAb [PSH23-00]	HA724339	Human	IHC-P,IF-Tissue
	TCF-4 / TCF7L2 Recombinant Rabbit mAb [SC06-90] - BSA and Azide free	HA751852	Human, Mouse, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC
Relative expression	TCF7 Recombinant Rabbit mAb [PSH22-98]	HA724337	Human	WB,IHC-P,IF-Tissue,IF-Cell,IP

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression, Cell treatment	Thrombomodulin Recombinant Rabbit mAb [PSH23-10]	HA724350	Human	WB,IF-Cell,IHC-P,FC,IP
	TREM2-dependent mTOR Metabolic Fitness Antibody Sampler Kit	HAK21193	Human, Mouse	WB
KD	TRPA1 Recombinant Rabbit mAb [PSH23-56]	HA724319	Human, Monkey	WB
Relative expression	Tyrosinase Recombinant Rabbit mAb [JA52-11] - BSA and Azide free	HA751836	Human, Mouse	WB,IF-Cell,IF-Tissue,IHC-P,IP,FC
	Versican Recombinant Rabbit mAb [JB42-32] - BSA and Azide free	HA751830	Human, Mouse, Rat	IHC-P,WB,IF-Tissue
Relative expression	Wilms Tumor Protein Recombinant Rabbit mAb [PSH23-02]	HA724341	Human	IHC-P
	YY1 Recombinant Rabbit mAb [SY29-01] - BSA and Azide free	HA751847	Human, Mouse, Rat	WB,IHC-P,IP,IF-Cell
	荧光淬灭试剂盒 IRISKit® Quench Buffer Kit	900105		mIHC

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层