

## 2026年4月—新上市产品清单

高级验证	产品名	货号	反应种属	验证应用
Cell treatment	4-Hydroxynonenal Recombinant Mouse mAb [PSH23-22]	HA601598	Species independent	WB,IF-Cell
Cell treatment	4-Hydroxynonenal Recombinant Mouse mAb [PSH23-22] - BSA and Azide free	HA610401	Species independent	WB,IF-Cell
Relative expression	ABCC4 Recombinant Rabbit mAb [PSH23-61]	HA724387	Human, Rat	WB,IF-Cell,IHC-P,FC,IP
Relative expression	ABCC4 Recombinant Rabbit mAb [PSH23-61] - BSA and Azide free	HA751903	Human, Rat	WB,IF-Cell,IHC-P,FC,IP
Relative expression	Adiponectin Recombinant Rabbit mAb [PSH22-55] - BSA and Azide free	HA751868	Human, Mouse, Rat	WB,IHC-P
	AGRP Recombinant Rabbit mAb [PSH22-63] - BSA and Azide free	HA751873	Mouse, Rat	IHC-Fr,IHC-P
Relative expression	Ak8 Rabbit pAb	HA500638	Human, Mouse, Rat	WB
Relative expression	ALDH1A3 Recombinant Rabbit mAb [PSH23-87]	HA724405	Human	WB,IF-Cell,FC,IP
Relative expression	ALDH1A3 Recombinant Rabbit mAb [PSH23-87] - BSA and Azide free	HA751916	Human	WB,IF-Cell,FC,IP
	iFluor™ 488 Conjugated Alpaca anti-Rabbit IgG Fc Chimeric Goat Fc Recombinant Antibody	HA1137	Rabbit	IF-Cell,IHC-Fr
	AMACR Recombinant Antibody [PSH13-54] - Mouse IgG1 (Chimeric)	HA601602	Human	IHC-P
	Anti-DYKDDDDK (Flag) Magnetic Beads	HAK21011		IP,Co-IP
Cell treatment	Apelin Receptor Rabbit pAb	HA500634	Mouse, Rat	WB
Relative expression	ATRX Rabbit pAb	HA500635	Human, Mouse, Rat	WB,IF-Cell,FC
Relative expression	BAP1 Recombinant Rabbit mAb [PSH22-53] - BSA and Azide free	HA751866	Human, Mouse, Rat	WB,IP
Cell treatment	BARD1 Recombinant Rabbit mAb [PSH22-30] - BSA and Azide free	HA751858	Human	WB,IP
Relative expression	Bassoon / BSN Recombinant Mouse mAb [PSH24-11]	HA601618	Mouse, Rat	IHC-Fr,IHC-P
Relative expression	Bassoon / BSN Recombinant Mouse mAb [PSH24-11] - BSA and Azide free	HA610416	Mouse, Rat	IHC-Fr,IHC-P
	CACNA1G Recombinant Mouse mAb [PSH23-59]	HA601604	Mouse, Rat	WB,IHC-P

	CACNA1G Recombinant Mouse mAb [PSH23-59] - BSA and Azide free	HA610406	Mouse, Rat	WB,IHC-P
	Cardiac Troponin I Recombinant Rabbit mAb [PSH23-75] - BSA and Azide free	HA751910	Mouse, Rat	WB
	Cardiac Troponin I Recombinant Rabbit mAb [PSH23-75]	HA724398	Mouse, Rat	WB
	CD130 / GP130 Recombinant Rabbit mAb [PSH22-56] - BSA and Azide free	HA751869	Mouse	WB
Relative expression	CD21 Recombinant Antibody [PD00-23] - Rat IgG1 (Chimeric)	HA601605	Human	WB,IHC-P
Relative expression	CD21 Recombinant Antibody [PD00-23] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610405	Human	WB,IHC-P
Relative expression	CD27 Recombinant Rabbit mAb [PSH23-81] - BSA and Azide free	HA751914	Mouse	WB,IHC-P,IP,IHC-Fr,FC
Relative expression	CD27 Recombinant Rabbit mAb [PSH24-12]	HA724423	Human	WB,IF-Cell,IHC-P,FC(Intra)
Relative expression	CD27 Recombinant Rabbit mAb [PSH24-12] - BSA and Azide free	HA751926	Human	WB,IF-Cell,IHC-P,FC(Intra)
Relative expression	CD27 Recombinant Rabbit mAb [PSH23-81]	HA724402	Mouse	WB,IHC-P,IP,IHC-Fr,FC
Relative expression	CD3 Recombinant Rabbit mAb [PSH22-52] - BSA and Azide free	HA751865	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue,IF-Cell,FC,IP
Relative expression	CD48 Recombinant Rabbit mAb [PSH23-53]	HA724384	Mouse	WB,FC
Relative expression	CD48 Recombinant Rabbit mAb [PSH23-53] - BSA and Azide free	HA751901	Mouse	WB,FC
Cell treatment	CD68 Recombinant Antibody [PSH06-65] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610398	Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Tissue
Cell treatment, Relative expression	CD68 Recombinant Antibody [PSH06-65] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610399	Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Cell,IF-Tissue
	CEP55 Recombinant Rabbit mAb [PSH23-80]	HA724401	Human, Mouse, Rat	WB,IF-Cell,FC
	CEP55 Recombinant Rabbit mAb [PSH23-80] - BSA and Azide free	HA751913	Human, Mouse, Rat	WB,IF-Cell,FC
Relative expression	Claudin 11 Recombinant Rabbit mAb [PSH24-17]	HA724428	Mouse, Rat	WB,IP
Relative expression	Claudin 11 Recombinant Rabbit mAb [PSH24-17] - BSA and Azide free	HA751931	Mouse, Rat	WB,IP
Relative expression	Claudin 4 Recombinant Rabbit mAb [PSH23-48]	HA724379	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC(Intra),IP
Relative expression	Claudin 4 Recombinant Rabbit mAb [PSH23-48] - BSA and Azide free	HA751896	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC(Intra),IP

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	Claudin 4 Recombinant Rabbit mAb [PSH23-49]	HA724380	Human, Mouse, Rat	WB,IHC-P,FC(Intra),IP
Relative expression	Claudin 4 Recombinant Rabbit mAb [PSH23-49] - BSA and Azide free	HA751897	Human, Mouse, Rat	WB,IHC-P,FC(Intra),IP
	COX5A Recombinant Rabbit mAb [PSH23-69] - BSA and Azide free	HA751908	Human, Mouse, Rat	WB,IHC-P
	COX5A Recombinant Rabbit mAb [PSH23-69]	HA724392	Human, Mouse, Rat	WB,IHC-P
Cell treatment	CXCL5 Recombinant Rabbit mAb [PSH22-58] - BSA and Azide free	HA751871	Human	WB,IF-Cell,FC
Cell treatment	CXCL5 Recombinant Rabbit mAb [PSH22-57] - BSA and Azide free	HA751870	Human	WB,IF-Cell,FC
Relative expression	CXCR4 Recombinant Rabbit mAb [PSH22-76]	HA724321	Human, Mouse	WB,IHC-Fr
Relative expression	CXCR4 Recombinant Rabbit mAb [PSH22-76] - BSA and Azide free	HA751877	Human, Mouse	WB,IHC-Fr
	Cytokeratin 10 Recombinant Antibody [ST05-43] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610415	Human, Mouse, Rat	IHC-P,IF-Tissue
	Cytokeratin 10 Recombinant Antibody [ST05-43] - Mouse IgG1 (Chimeric)	HA601617	Human, Mouse, Rat	IHC-P,IF-Tissue
	Cytokeratin 8 Recombinant Antibody [SU0338] - Mouse IgG1 (Chimeric)	HA601599	Human	IF-Cell,IHC-P
	Cytokeratin 8 Recombinant Antibody [SU0338] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610402	Human	IF-Cell,IHC-P
Relative expression	Dopamine beta Hydroxylase Recombinant Rabbit mAb [PSH23-88]	HA724406	Mouse, Rat	WB,IHC-Fr,IHC-P,IP
Relative expression	Dopamine beta Hydroxylase Recombinant Rabbit mAb [PSH23-88] - BSA and Azide free	HA751917	Mouse, Rat	WB,IHC-Fr,IHC-P,IP
	EAAT1 Recombinant Antibody [JA30-35] - Rat IgG1 (Chimeric)	HA601614	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	EAAT1 Recombinant Antibody [JA30-35] - Guinea pig IgG2 (Chimeric)	HA601613	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	EAAT1 Recombinant Antibody [JA30-35] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610414	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
	EAAT1 Recombinant Antibody [JA30-35] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610413	Human, Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression	Endomucin Recombinant Rabbit mAb [PSH24-00]	HA724412	Mouse	WB,IHC-P,mIHC
Relative expression	Endomucin Recombinant Rabbit mAb [PSH24-00] - BSA and Azide free	HA751922	Mouse	WB,IHC-P
Relative expression	Eph receptor B2 Recombinant Rabbit mAb [PSH23-74]	HA724397	Human, Mouse, Rat	WB,IHC-Fr,IP

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	Eph receptor B2 Recombinant Rabbit mAb [PSH23-74] - BSA and Azide free	HA751909	Human, Mouse, Rat	WB,IHC-Fr,IP
	EPS15 Recombinant Rabbit mAb [PSH23-03] - BSA and Azide free	HA751885	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC
	FAF1 Recombinant Rabbit mAb [PSH22-51] - BSA and Azide free	HA751864	Human, Mouse, Rat, Green monkey	WB,IF-Cell,IHC-P,FC,IP
	FGF2 Recombinant Rabbit mAb [PSH22-75] - BSA and Azide free	HA751876	Human	WB,IF-Cell,FC
Relative expression	FGFR1 Recombinant Rabbit mAb [PSH22-64] - BSA and Azide free	HA751874	Human, Mouse, Rat	WB,IHC-P
Relative expression	Filaggrin Recombinant Rabbit mAb [PSH23-66] - BSA and Azide free	HA751907	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,FC
Relative expression	Filaggrin Recombinant Rabbit mAb [PSH23-66]	HA724391	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,FC
Relative expression	FOLR1 Recombinant Mouse mAb [PSH23-62]	HA601607	Human	WB,IF-Cell,FC
Relative expression	FOLR1 Recombinant Mouse mAb [PSH23-62] - BSA and Azide free	HA610407	Human	WB,IF-Cell,FC
Relative expression	FOXP2 Recombinant Rabbit mAb [PSH23-79]	HA724400	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	FOXP2 Recombinant Rabbit mAb [PSH23-79] - BSA and Azide free	HA751912	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
Relative expression	GAD67 Recombinant Antibody [PSH08-30] - Rabbit IgG (Chimeric) - BSA and Azide free	HA751895	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	GAD67 Recombinant Antibody [PSH08-30] - Rabbit IgG (Chimeric)	HA724362	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	GAD67 Recombinant Antibody [PSH08-30] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610408	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	GAD67 Recombinant Antibody [PSH08-30] - Guinea pig IgG2 (Chimeric)	HA601606	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	GATA4 Recombinant Rabbit mAb [PSH23-98]	HA724409	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC(Intra)
Relative expression	GATA4 Recombinant Rabbit mAb [PSH23-98] - BSA and Azide free	HA751920	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC(Intra)
	GC / Vitamin D Binding protein Recombinant Rabbit mAb [PSH24-22]	HA724429	Human, Mouse, Rat	WB
	GC / Vitamin D Binding protein Recombinant Rabbit mAb [PSH24-22] - BSA and Azide free	HA751932	Human, Mouse, Rat	WB
	Gemin 1 Recombinant Rabbit mAb [PSH23-89] - BSA and Azide free	HA751918	Human	WB,IF-Cell,FC,IP
	Gemin 1 Recombinant Rabbit mAb [PSH23-89]	HA724407	Human	WB,IF-Cell,FC,IP

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	Granzyme K Recombinant Rabbit mAb [PSH23-21]	HA724356	Human	WB,IF-Cell,IHC-P,IF-Tissue
Relative expression	Granzyme K Recombinant Rabbit mAb [PSH23-21] - BSA and Azide free	HA751890	Human	WB,IF-Cell,IHC-P,IF-Tissue
	Hemopexin Recombinant Rabbit mAb [PSH23-09] - BSA and Azide free	HA751887	Human, Mouse, Rat	WB
Cell treatment	Histone H3 (citrulline R2) Recombinant Rabbit mAb [PSH23-25]	HA724360	Human, Mouse	WB,IP
Cell treatment	Histone H3 (citrulline R2) Recombinant Rabbit mAb [PSH23-25] - BSA and Azide free	HA751893	Human, Mouse	WB,IP
	HOXB5 Rabbit pAb	HA500633	Mouse, Rat	WB
	Human Adiponectinenzyme-linked immunoassay kit ( one step )	EHY0099	Human	ELISA
	Human Alpha 2 Antiplasmin/SerpinF2, C-His Tag Protein	HA211518	Human	
	Human APLP1, C-His Tag Protein	HA211520	Human	
	Human Ataxin 3 Recombinant Rabbit mAb [PSH23-45] - BSA and Azide free (Detector)	HA724376	Human	ELISA(Det)
	Biotin Conjugated Human Ataxin 3 Recombinant Rabbit mAb [PSH23-45] - Detector	HA724376B	Human	ELISA(Det),ELISA
	Human Ataxin 3 Recombinant Rabbit mAb [PSH23-44] - BSA and Azide free (Capture)	HA724375	Human	ELISA(Cap)
	Human BAFFR/TNFRSF13C, Tag Free Protein	HA211435	Human	
	Human BST2 Recombinant Rabbit mAb [PSH23-96] - BSA and Azide free (Capture)	HA725329	Human	ELISA(Cap)
	Biotin Conjugated Human BST2 Recombinant Rabbit mAb [PSH23-97] - Detector	HA725330B	Human	ELISA(Det),ELISA
	Human BST2 Recombinant Rabbit mAb [PSH23-97] - BSA and Azide free (Detector)	HA725330	Human	ELISA(Det)
	Biotin Conjugated Human CA125/MUC16 Recombinant Rabbit mAb [PSH23-20] - Detector	HA724355B	Human	ELISA(Det),ELISA
	Human CA125/MUC16 Recombinant Rabbit mAb [PSH23-20] - BSA and Azide free (Detector)	HA724355	Human	ELISA(Det)
	Human CA125/MUC16 Recombinant Rabbit mAb [PSH23-19] - BSA and Azide free (Capture)	HA724354	Human	ELISA(Cap)
	Human CCL1 Recombinant Rabbit mAb [PSH24-05] - BSA and Azide free (Detector)	HA724417	Human	ELISA(Det)

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Human CCL1 Recombinant Rabbit mAb [PSH24-04] - BSA and Azide free (Capture)	HA724416	Human	ELISA(Cap)
Biotin Conjugated Human CCL1 Recombinant Rabbit mAb [PSH24-05] - Detector	HA724417B	Human	ELISA(Det),ELISA
Biotin Conjugated Human CCL15/MIP-5 Recombinant Rabbit mAb [PSH23-12] - Detector	HA725342B	Human	ELISA(Det),ELISA
Human CCL15/MIP-5 Recombinant Rabbit mAb [PSH23-11] - BSA and Azide free (Capture)	HA725341	Human	ELISA(Cap)
Human CCL15/MIP-5 Recombinant Rabbit mAb [PSH23-12] - BSA and Azide free (Detector)	HA725342	Human	ELISA(Det)
Human CCL24/Eotaxin 2 Recombinant Rabbit mAb [PSH23-39] - BSA and Azide free (Detector)	HA724370	Human	ELISA(Det)
Human CCL24/Eotaxin 2 Recombinant Rabbit mAb [PSH23-38] - BSA and Azide free (Capture)	HA724369	Human	ELISA(Cap)
Biotin Conjugated Human CCL24/Eotaxin 2 Recombinant Rabbit mAb [PSH23-39] - Detector	HA724370B	Human	ELISA(Det),ELISA
Biotin Conjugated Human CD40 Recombinant Rabbit mAb [PSH22-80] - Detector	HA724325B	Human	ELISA(Det),ELISA
Human CD40 Recombinant Rabbit mAb [PSH22-80] - BSA and Azide free (Detector)	HA724325	Human	ELISA(Det)
Human CD40 Recombinant Rabbit mAb [PSH22-79] - BSA and Azide free (Capture)	HA724324	Human	ELISA(Cap)
Human Coagulation factor XII/F12, Tag Free Protein	HA211374	Human	
Human Desmoglein-3, C-His Tag Protein	HA211402	Human	
Human FGF23 Recombinant Rabbit mAb [PSH23-70] - BSA and Azide free (Capture)	HA724393	Human	ELISA(Cap)
Biotin Conjugated Human FGF23 Recombinant Rabbit mAb [PSH23-71] - Detector	HA724394B	Human	ELISA(Det),ELISA
Human FGF23 Recombinant Rabbit mAb [PSH23-71] - BSA and Azide free (Detector)	HA724394	Human	ELISA(Det)
Human Ficolin 2 Recombinant Rabbit mAb [PSH23-90] - BSA and Azide free (Capture)	HA725209	Human	ELISA(Cap)
Biotin Conjugated Human Ficolin 2 Recombinant Rabbit mAb [PSH23-91] - Detector	HA725210B	Human	ELISA(Det),ELISA

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Human Ficolin 2 Recombinant Rabbit mAb [PSH23-91] - BSA and Azide free (Detector)	HA725210	Human	ELISA(Det)
Human Follistatin-like 1/FSTL1, C-His Tag Protein	HA211443	Human	
Human GDF-15, Tag Free Protein	HA211526	Human	
Human Gemin 1 Recombinant Rabbit mAb [PSH24-06] - BSA and Azide free (Capture)	HA724418	Human	ELISA(Cap)
Human Gemin 1 Recombinant Rabbit mAb [PSH24-07] - BSA and Azide free (Detector)	HA724419	Human	ELISA(Det)
Biotin Conjugated Human Gemin 1 Recombinant Rabbit mAb [PSH24-07] - Detector	HA724419B	Human	ELISA(Det),ELISA
Human GLP1R, C-His Tag Protein	HA211370	Human	
Human Glypican 3/GPC3, C-His Tag Protein	HA211441	Human	
Human HO1/HMOX1, C-His Tag Protein	HA211372	Human	
Human IFNB1/IFN beta, Tag Free Protein	HA211523	Human	
Human IFN- $\gamma$ Super Sensitivity enzyme-linked immunoassay kit	EHSS0039	Human	ELISA
Human IL-19 Recombinant Rabbit mAb [PSH23-29] - BSA and Azide free (Detector)	HA725272	Human	ELISA(Det)
Human IL-19 Recombinant Rabbit mAb [PSH23-28] - BSA and Azide free (Capture)	HA725271	Human	ELISA(Cap)
Biotin Conjugated Human IL-19 Recombinant Rabbit mAb [PSH23-29] - Detector	HA725272B	Human	ELISA(Det),ELISA
Human IL-20, Tag Free Protein	HA211525	Human	
Human IL-6 Super Sensitivity enzyme-linked immunoassay kit	EHSS0001	Human	ELISA
Human L-FABP Recombinant Rabbit mAb [PSH22-97] - BSA and Azide free (Detector)	HA725448	Human	ELISA(Det)
Human L-FABP Recombinant Rabbit mAb [PSH22-96] - BSA and Azide free (Capture)	HA725447	Human	ELISA(Cap)
Biotin Conjugated Human L-FABP Recombinant Rabbit mAb [PSH22-97] - Detector	HA725448B	Human	ELISA(Det),ELISA
Human MCP-3/CCL7 Recombinant Rabbit mAb [PSH22-95] - BSA and Azide free (Detector)	HA725360	Human	ELISA(Det)

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Biotin Conjugated Human MCP-3/CCL7 Recombinant Rabbit mAb [PSH22-95] - Detector	HA725360B	Human	ELISA(Det),ELISA
Human MCP-3/CCL7 Recombinant Rabbit mAb [PSH22-94] - BSA and Azide free (Capture)	HA725359	Human	ELISA(Cap)
Human MMP-12 Recombinant Rabbit mAb [PSH22-92] - BSA and Azide free (Capture)	HA725323	Human	ELISA(Cap)
Biotin Conjugated Human MMP-12 Recombinant Rabbit mAb [PSH22-93] - Detector	HA725324B	Human	ELISA(Det),ELISA
Human MMP-12 Recombinant Rabbit mAb [PSH22-93] - BSA and Azide free (Detector)	HA725324	Human	ELISA(Det)
Biotin Conjugated Human/Monkey ACE Recombinant Rabbit mAb [PSH22-82] - Detector	HA724327B	Human, Cynomolgus monkey	ELISA(Det),ELISA
Human/Monkey ACE Recombinant Rabbit mAb [PSH22-81] - BSA and Azide free (Capture)	HA724326	Human, Cynomolgus monkey	ELISA(Cap)
Human/Monkey ACE Recombinant Rabbit mAb [PSH22-82] - BSA and Azide free (Detector)	HA724327	Human, Cynomolgus monkey	ELISA(Det)
Human/Monkey Mannose Receptor(CD206) Recombinant Rabbit mAb [PSH23-47] - BSA and Azide free (Detector)	HA724378	Human, Cynomolgus monkey	ELISA(Det)
Human/Monkey Mannose Receptor(CD206) Recombinant Rabbit mAb [PSH23-46] - BSA and Azide free (Capture)	HA724377	Human, Cynomolgus monkey	ELISA(Cap)
Biotin Conjugated Human/Monkey Mannose Receptor(CD206) Recombinant Rabbit mAb [PSH23-47] - Detector	HA724378B	Human, Cynomolgus monkey	ELISA(Det),ELISA
Human/Monkey PCSK9 enzyme-linked immunoassay kit (one step)	EHY0240	Human	ELISA
HRP Conjugated Human/Monkey VEGFR3 / FLT4 Recombinant Rabbit mAb [PSH22-17] - Detector	HA724278H	Human, Cynomolgus monkey	ELISA(Det),ELISA
Human/Mouse GLP-1 Recombinant Rabbit mAb [PSH23-58] - BSA and Azide free (Detector)	HA725354	Human, Mouse	ELISA(Det)
Biotin Conjugated Human/Mouse GLP-1 Recombinant Rabbit mAb [PSH23-58] - Detector	HA725354B	Human, Mouse	ELISA(Det),ELISA
Human/Mouse GLP-1 Recombinant Rabbit mAb [PSH23-57] - BSA and Azide free (Capture)	HA725353	Human, Mouse	ELISA(Cap)

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

HRP Conjugated Human/Mouse MASP2 Recombinant Rabbit mAb [PSH20-52] - Detector	HA725392H	Human, Mouse, Rhesus monkey	ELISA(Det),ELISA
Human/Mouse MASP2 enzyme-linked immunoassay kit	EH0231	Human	ELISA
Biotin Conjugated Human/Mouse PNS-tau Recombinant Rabbit mAb [PSH23-08] - BSA and Azide free (Detector)	HA724348B	Human, Mouse	ELISA(Det),ELISA
Human/Mouse PNS-tau Recombinant Rabbit mAb [PSH23-07] - BSA and Azide free (Capture)	HA724347	Human, Mouse	ELISA(Cap)
Human/Mouse PNS-tau Recombinant Rabbit mAb [PSH23-08] - BSA and Azide free (Detector)	HA724348	Human, Mouse	ELISA(Det)
Human NGFR/TNFRSF16, Tag Free Protein	HA211371	Human	
Human PDGF alpha receptor/PDGFR $\alpha$ , C-His Tag Protein	HA211519	Human	
Human PDLIM5, C-His Tag Protein	HA211357	Human	
Human phosphorylated MAPK1/ERK-2 (pT185 and pY187) Protein	HA211528	Human	
Human PLGF Recombinant Rabbit mAb [PSH23-86] - BSA and Azide free (Detector)	HA725378	Human	ELISA(Det)
Human PLGF Recombinant Rabbit mAb [PSH23-85] - BSA and Azide free (Capture)	HA725377	Human	ELISA(Cap)
Biotin Conjugated Human PLGF Recombinant Rabbit mAb [PSH23-86] - Detector	HA725378B	Human	ELISA(Det),ELISA
Human proBNP, Tag Free Protein	HA211505	Human	
Human RANK/TNFRSF11A, C-His Tag Protein	HA211522	Human	
Human Renin enzyme-linked immunoassay kit	EH0211	Human	ELISA
Human SCGB1A1 Recombinant Rabbit mAb [PSH23-84] - BSA and Azide free (Detector)	HA725471	Human	ELISA(Det)
Human SCGB1A1 Recombinant Rabbit mAb [PSH23-83] - BSA and Azide free (Capture)	HA725470	Human	ELISA(Cap)
Biotin Conjugated Human SCGB1A1 Recombinant Rabbit mAb [PSH23-84] - Detector	HA725471B	Human	ELISA(Det),ELISA
Human Sclerostin Recombinant Rabbit mAb [PSH23-42] - BSA and Azide free (Capture)	HA724373	Human	ELISA(Cap)

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human Sclerostin Recombinant Rabbit mAb [PSH23-43] - Detector	HA724374B	Human	ELISA(Det),ELISA
	Human Sclerostin Recombinant Rabbit mAb [PSH23-43] - BSA and Azide free (Detector)	HA724374	Human	ELISA(Det)
	Human Serpin A12 Recombinant Rabbit mAb [PSH24-21] - BSA and Azide free (Detector)	HA725481	Human	ELISA(Det)
	Human Serpin A12 Recombinant Rabbit mAb [PSH24-20] - BSA and Azide free (Capture)	HA725480	Human	ELISA(Cap)
	Biotin Conjugated Human Serpin A12 Recombinant Rabbit mAb [PSH24-21] - Detector	HA725481B	Human	ELISA(Det),ELISA
	Human Serpin F1/PEDF, C-His Tag Protein	HA211403	Human	
	Human SERPING1/C1 inhibitor, C-His Tag Protein	HA211517	Human	
	Human SNAP25 enzyme-linked immunoassay kit	EH0228	Human	ELISA
	Biotin Conjugated Human SP-D Recombinant Rabbit mAb [PSH23-95] - Detector	HA725218B	Human	ELISA(Det),ELISA
	Human SP-D Recombinant Rabbit mAb [PSH23-95] - BSA and Azide free (Detector)	HA725218	Human	ELISA(Det)
	Human SP-D Recombinant Rabbit mAb [PSH23-94] - BSA and Azide free (Capture)	HA725217	Human	ELISA(Cap)
	Human Thrombomodulin Recombinant Rabbit mAb [PSH23-37] - BSA and Azide free (Detector)	HA724368	Human	ELISA(Det)
	Biotin Conjugated Human Thrombomodulin Recombinant Rabbit mAb [PSH23-37] - Detector	HA724368B	Human	ELISA(Det),ELISA
	Human Thrombomodulin Recombinant Rabbit mAb [PSH23-36] - BSA and Azide free (Capture)	HA724367	Human	ELISA(Cap)
	Human TNNT3, Tag Free Protein	HA211524	Human	
	Human TREM-1 enzyme-linked immunoassay kit	EH0227	Human	ELISA
	Human VSTM1 Recombinant Rabbit mAb [PSH23-40] - BSA and Azide free (Capture)	HA724371	Human	ELISA(Cap)
	Biotin Conjugated Human VSTM1 Recombinant Rabbit mAb [PSH23-41] - Detector	HA724372B	Human	ELISA(Det),ELISA

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human VSTM1 Recombinant Rabbit mAb [PSH23-41] - BSA and Azide free (Detector)	HA724372	Human	ELISA(Det)
Relative expression	IL-6R Recombinant Rabbit mAb [PSH22-78] - BSA and Azide free	HA751879	Human	WB,FC
	INSIG1 Recombinant Rabbit mAb [PSH22-32] - BSA and Azide free	HA751860	Human, Mouse, Rat	WB
Relative expression	Involucrin Recombinant Rabbit mAb [PSH23-99] - BSA and Azide free	HA751921	Human	WB,IHC-P
Relative expression	Involucrin Recombinant Rabbit mAb [PSH23-99]	HA724410	Human	WB,IHC-P
	IRS2 Recombinant Rabbit mAb [PSH22-29] - BSA and Azide free	HA751857	Human, Mouse, Rat	IHC-P,IF-Cell,FC
Cell treatment	LAMP1 Recombinant Antibody [PSH07-39] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610400	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	LDHB Recombinant Rabbit mAb [PSH22-91] - BSA and Azide free	HA751880	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,FC,IP
	LYVE1 Recombinant Rabbit mAb [PSH23-01] - BSA and Azide free	HA751884	Mouse	WB,IHC-Fr,IHC-P
Cell treatment	Mannose Receptor(CD206) Recombinant Rabbit mAb [PSH23-93] - BSA and Azide free	HA751919	Human	WB,IHC-P,IF-Tissue,IP
Cell treatment	Mannose Receptor(CD206) Recombinant Rabbit mAb [PSH23-93]	HA724408	Human	WB,IHC-P,IF-Tissue,IP
	MCP-3/CCL7 Recombinant Rabbit mAb [PSH23-04] - BSA and Azide free	HA751886	Human, Mouse	WB
Relative expression, KD	MDC1 Recombinant Rabbit mAb [PSH23-24]	HA724359	Human, Mouse, Rat	WB,IF-Cell
Relative expression, KD	MDC1 Recombinant Rabbit mAb [PSH23-24] - BSA and Azide free	HA751892	Human, Mouse, Rat	WB,IF-Cell
	HRP Conjugated Monkey PCSK9 Recombinant Mouse mAb [PSH23-55] - Detector	HA601603H	Cynomolgus monkey	ELISA(Det),ELISA
	Monkey PCSK9 Recombinant Rabbit mAb [PSH23-54] - BSA and Azide free (Capture)	HA724385	Cynomolgus monkey	ELISA(Cap)
	Monkey PCSK9 Recombinant Mouse mAb [PSH23-55] - BSA and Azide free (Detector)	HA601603	Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Monkey PCSK9 Recombinant Mouse mAb [PSH23-55] - Detector	HA601603B	Cynomolgus monkey	ELISA(Det),ELISA
	Mouse ACE2 Recombinant Rabbit mAb [PSH23-73]	HA724396	Mouse	ELISA(Det)

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Mouse ACE2 Recombinant Rabbit mAb [PSH23-72]	HA724395	Mouse	ELISA(Cap)
Biotin Conjugated Mouse ACE2 Recombinant Rabbit mAb [PSH23-73]	HA724396B	Mouse	ELISA(Det),ELISA
HRP Conjugated Mouse Adiponectin Recombinant Rabbit mAb [PSH14-83] - Detector	HA723681H	Mouse	ELISA(Det),ELISA
Mouse Adiponectin enzyme-linked immunoassay kit(one step)	EMY0064	Mouse	ELISA
HRP Conjugated Mouse Apolipoprotein E Recombinant Rabbit mAb [PSH16-84] - Detector	HA725222H	Mouse	ELISA(Det),ELISA
Biotin Conjugated Mouse CD25 / IL-2 Receptor alpha Recombinant Rabbit mAb [PSH24-25] - Detector	HA724432B	Mouse	ELISA(Det),ELISA
Mouse CD25 / IL-2 Receptor alpha Recombinant Rabbit mAb [PSH24-25] - BSA and Azide free (Detector)	HA724432	Mouse	ELISA(Det)
Mouse CD25 / IL-2 Receptor alpha Recombinant Rabbit mAb [PSH24-24] - BSA and Azide free (Capture)	HA724431	Mouse	ELISA(Cap)
Biotin Conjugated Mouse CD27 Recombinant Rabbit mAb [PSH23-18] - Detector	HA724353B	Mouse	ELISA(Det),ELISA
Mouse CD27 Recombinant Rabbit mAb [PSH23-17] - BSA and Azide free (Capture)	HA724352	Mouse	ELISA(Cap)
Mouse CD27 Recombinant Rabbit mAb [PSH23-18] - BSA and Azide free (Detector)	HA724353	Mouse	ELISA(Det)
HRP Conjugated Mouse Complement Factor D Recombinant Rabbit mAb [PSH19-08] - Detector	HA724042H	Mouse	ELISA(Det),ELISA
Biotin Conjugated Mouse CXCL10 / IP-10 Recombinant Rabbit mAb [PSH24-09] - Detector	HA724421B	Mouse	ELISA(Det),ELISA
Mouse CXCL10 / IP-10 Recombinant Rabbit mAb [PSH24-09] - BSA and Azide free (Detector)	HA724421	Mouse	ELISA(Det)
Mouse CXCL10 / IP-10 Recombinant Rabbit mAb [PSH24-08] - BSA and Azide free (Capture)	HA724420	Mouse	ELISA(Cap)
Mouse Fas/CD95/TNFRSF6, C-His Tag Protein	HA211527	Mouse	
Mouse GPT, Tag Free Protein	HA211503	Mouse	
Mouse IGFBP1, Tag Free Protein	HA211409	Mouse	
Mouse iNOS Recombinant Rabbit mAb [PSH23-67] - BSA and Azide free (Capture)	HA725431	Mouse	ELISA(Cap)

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse iNOS Recombinant Rabbit mAb [PSH23-68] - Detector	HA725432B	Mouse	ELISA(Det),ELISA
	Mouse iNOS Recombinant Rabbit mAb [PSH23-68] - BSA and Azide free (Detector)	HA725432	Mouse	ELISA(Det)
	Mouse LBP, Tag Free Protein	HA211401	Mouse	
	Mouse NOS2/iNOS, C-His Tag Protein	HA211501	Mouse	
	Mouse PCSK9 Recombinant Rabbit mAb [PSH23-33] - BSA and Azide free (Detector)	HA724364	Mouse	ELISA(Det)
	Mouse PCSK9 Recombinant Rabbit mAb [PSH23-32] - BSA and Azide free (Capture)	HA724363	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse PCSK9 Recombinant Rabbit mAb [PSH23-33] - Detector	HA724364B	Mouse	ELISA(Det),ELISA
	HRP Conjugated Mouse PCSK9 Recombinant Rabbit mAb [PSH23-33] - Detector	HA724364H	Mouse	ELISA(Det),ELISA
	Mouse/Rat Alpha-synuclein Recombinant Rabbit mAb [PSH23-06] - BSA and Azide free (Detector)	HA724346	Mouse, Rat	ELISA(Det)
	Mouse/Rat Alpha-synuclein Recombinant Rabbit mAb [PSH23-05] - BSA and Azide free (Capture)	HA724345	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat Alpha-synuclein Recombinant Rabbit mAb [PSH23-06] - Detector	HA724346B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat GDF-15 enzyme-linked immunoassay kit	EM0058	Mouse	ELISA
	Mouse/Rat MMP-3 Recombinant Rabbit mAb [PSH23-35] - BSA and Azide free (Detector)	HA724366	Mouse, Rat	ELISA(Det)
	Mouse/Rat MMP-3 Recombinant Rabbit mAb [PSH23-34] - BSA and Azide free (Capture)	HA724365	Mouse, Rat	ELISA(Cap)
	Biotin Conjugated Mouse/Rat MMP-3 Recombinant Rabbit mAb [PSH23-35] - Detector	HA724366B	Mouse, Rat	ELISA(Det),ELISA
	HRP Conjugated Mouse/Rat P-Cadherin Recombinant Rabbit mAb [PSH21-93] - Detector	HA724262H	Mouse, Rat	ELISA(Det),ELISA
	HRP Conjugated Mouse/Rat Prealbumin Recombinant Rabbit mAb [PSH15-88] - Detector	HA725129H	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat Von Willebrand Factor Recombinant Rabbit mAb [PSH23-15] - BSA and Azide free (Capture)	HA725158	Mouse, Rat	ELISA(Cap)

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Mouse/Rat Von Willebrand Factor Recombinant Rabbit mAb [PSH23-16] - Detector	HA725159B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat Von Willebrand Factor Recombinant Rabbit mAb [PSH23-16] - BSA and Azide free (Detector)	HA725159	Mouse, Rat	ELISA(Det)
	Mouse RBP4, Tag Free Protein	HA211521	Mouse	
	Biotin Conjugated Mouse sTNF RI Recombinant Rabbit mAb [PSH23-31] - Detector	HA725434B	Mouse	ELISA(Det),ELISA
	Mouse sTNF RI Recombinant Rabbit mAb [PSH23-31] - BSA and Azide free (Detector)	HA725434	Mouse	ELISA(Det)
	Mouse sTNF RI Recombinant Rabbit mAb [PSH23-30] - BSA and Azide free (Capture)	HA725433	Mouse	ELISA(Cap)
	MSH2 Recombinant Rabbit mAb [PSH24-02]	HA724414	Human, Mouse, Rat	IHC-P,IF-Cell,FC(Intra)
Relative expression, Cell treatment	MSR1 / CD204 Recombinant Rabbit mAb [PSH23-60]	HA724386	Mouse	WB
Relative expression, Cell treatment	MSR1 / CD204 Recombinant Rabbit mAb [PSH23-60] - BSA and Azide free	HA751902	Mouse	WB
Relative expression	MTBR-TAU 243 Isoform 2N4R (236-250) Recombinant Rabbit mAb [PSH22-59] - BSA and Azide free	HA751872	Mouse, Rat, Pig, Human	WB,IP,FC,IHC-P
Relative expression	MUC5B Recombinant Rabbit mAb [PSH23-64]	HA724389	Human	IHC-P,IF-Tissue
Relative expression	MUC5B Recombinant Rabbit mAb [PSH23-64] - BSA and Azide free	HA751905	Human	IHC-P,IF-Tissue
	Myelin Basic Protein Recombinant Antibody [JF0943] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610404	Mouse, Rat	IHC-Fr,IHC-P,WB
	Myelin Basic Protein Recombinant Antibody [JF0943] - Guinea pig IgG2 (Chimeric)	HA601600	Mouse, Rat	IHC-Fr,WB
	Myelin Basic Protein Recombinant Antibody [JF0943] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610403	Mouse, Rat	IHC-Fr,WB
	Myelin Basic Protein Recombinant Antibody [JF0943] - Rat IgG1 (Chimeric)	HA601601	Mouse, Rat	IHC-Fr,IHC-P,WB
Relative expression, Cell treatment	Myogenin Recombinant Mouse mAb [PSH22-35] - BSA and Azide free	HA610396	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IHC-Fr
Relative expression	Nav1.6 / SCN8A Recombinant Rabbit mAb [PSH22-36] - BSA and Azide free	HA751861	Mouse, Rat	WB,IHC-P,IP

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	NK1.1 Recombinant Rabbit mAb [PSH22-99] - BSA and Azide free	HA751882	Mouse	IHC-Fr,IF-Cell,FC
Cell treatment	NR4A3 Recombinant Rabbit mAb [PSH24-01] - BSA and Azide free	HA751923	Human	WB
Cell treatment	NR4A3 Recombinant Rabbit mAb [PSH24-01]	HA724413	Human	WB
Relative expression	NUPR1 Rabbit pAb	HA500637	Human	WB,IF-Cell,FC(Intra)
	NuRD Complex Antibody Sampler Kit	HAK21198	Human, Mouse	WB,IHC-P
Cell treatment, Relative expression	OASL Recombinant Rabbit mAb [PSH24-15]	HA724426	Human	WB,IHC-P,FC,IP
Cell treatment, Relative expression	OASL Recombinant Rabbit mAb [PSH24-15] - BSA and Azide free	HA751929	Human	WB,IHC-P,FC,IP
Relative expression, Cell treatment	p21 Recombinant Rabbit mAb [PSH22-37] - BSA and Azide free	HA751862	Human	WB,IF-Cell,IHC-P,IP
Relative expression	P2X7 Recombinant Rabbit mAb [PSH22-44] - BSA and Azide free	HA751863	Mouse	WB
	pan Cytokeratin Recombinant Antibody [PDH09-10] - Rabbit IgG (Chimeric)	HA724411F	Human	IF-Cell
	PARK7 Recombinant Mouse mAb [PSH23-76]	HA601611	Human	WB,IF-Cell,IHC-P,FC
	PARK7 Recombinant Mouse mAb [PSH23-76] - BSA and Azide free	HA610411	Human	WB,IF-Cell,IHC-P,FC
	PARK7 Recombinant Mouse mAb [PSH23-77]	HA601612	Human, Mouse, Rat	WB,IHC-P
	PARK7 Recombinant Mouse mAb [PSH23-77] - BSA and Azide free	HA610412	Human, Mouse, Rat	WB,IHC-P
	PAX5 Recombinant Antibody [PSH06-99] - Mouse IgG1 (Chimeric)	HA601616	Human, Mouse, Rat	IHC-P
Relative expression	PCK1 Recombinant Rabbit mAb [PSH23-82] - BSA and Azide free	HA751915	Mouse, Rat	WB
Relative expression	PCK1 Recombinant Rabbit mAb [PSH23-82]	HA724404	Mouse, Rat	WB
	Phospho-ACSL4 (T484) Rabbit pAb	HA500600	Human	WB
	Phospho-ALOX15 (T537) Rabbit pAb	HA500601	Human	WB
	Phospho-ALOX15 (T560) Rabbit pAb	HA500599	Human	WB
Cell treatment	Phospho-EGFR (Y992) Recombinant Rabbit mAb [PSH23-63]	HA724388	Human	WB,IF-Cell,FC

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-EGFR (Y992) Recombinant Rabbit mAb [PSH23-63] - BSA and Azide free	HA751904	Human	WB,IF-Cell,FC
Cell treatment	Phospho-Syk (Y526) Recombinant Rabbit mAb [PSH23-78]	HA724399	Human	WB,IP
Cell treatment	Phospho-Syk (Y526) Recombinant Rabbit mAb [PSH23-78] - BSA and Azide free	HA751911	Human	WB,IP
Cell treatment	Phospho-TAK1 (T184, T187) Recombinant Rabbit mAb [PSH24-03]	HA724415	Human, Mouse	WB
Cell treatment	Phospho-TAK1 (T184, T187) Recombinant Rabbit mAb [PSH24-03] - BSA and Azide free	HA751924	Human, Mouse	WB
Cell treatment	Phospho-TDP43 (S409, S410) Rabbit pAb	HA500636	Human, Mouse, Rat	WB,IHC-P
	PINK1 Recombinant Mouse mAb [PSH22-50] - BSA and Azide free	HA610397	Human, Mouse, Rat	WB
Relative expression	RAB27A Recombinant Rabbit mAb [PSH24-10]	HA724422	Human, Mouse, Rat	WB,IHC-P,FC,IP
Relative expression	RAB27A Recombinant Rabbit mAb [PSH24-10] - BSA and Azide free	HA751925	Human, Mouse, Rat	WB,IHC-P,FC,IP
	HRP Conjugated Rat ICAM1 Recombinant Rabbit mAb [PSH17-52] - Detector	HA723927H	Rat	ELISA(Det),ELISA
	Rat IL-5 Recombinant Rabbit mAb [PSH23-14] - BSA and Azide free (Detector)	HA725442	Rat	ELISA(Det)
	Biotin Conjugated Rat IL-5 Recombinant Rabbit mAb [PSH23-14] - Detector	HA725442B	Rat	ELISA(Det),ELISA
	Rat IL-5 Recombinant Rabbit mAb [PSH23-13] - BSA and Azide free (Capture)	HA725441	Rat	ELISA(Cap)
	Rat RELA/NF-kappa B p65, Tag Free Protein	HA211404	Rat	
Relative expression	RIP3 Recombinant Rabbit mAb [PSH23-11] - BSA and Azide free	HA751889	Mouse	WB,IP
Relative expression	ROR alpha Recombinant Rabbit mAb [PSH24-13] - BSA and Azide free	HA751927	Human, Mouse, Rat	WB,IHC-P
Relative expression	ROR alpha Recombinant Rabbit mAb [PSH24-13]	HA724424	Human, Mouse, Rat	WB,IHC-P
	RPS10 Recombinant Rabbit mAb [PSH23-26]	HA724361	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	RPS10 Recombinant Rabbit mAb [PSH23-26] - BSA and Azide free	HA751894	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC
	SEL1L Recombinant Rabbit mAb [PSH22-54] - BSA and Azide free	HA751867	Human, Mouse, Rat	WB,IF-Cell,FC
Relative expression	SLC12A2 / NKCC1 Recombinant Rabbit mAb [PSH23-65] - BSA and Azide free	HA751906	Human, Mouse, Rat	WB,IF-Cell,IHC-Fr,IHC-P,FC(Intra),IP

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	SLC12A2 / NKCC1 Recombinant Rabbit mAb [PSH23-65]	HA724390	Human, Mouse, Rat	WB,IF-Cell,IHC-Fr,IHC-P,FC(Intra),IP
Relative expression	SLC12A3 Recombinant Rabbit mAb [PSH24-16]	HA724427	Mouse, Rat	WB,IHC-Fr,IHC-P,IP
Relative expression	SLC12A3 Recombinant Rabbit mAb [PSH24-16] - BSA and Azide free	HA751930	Mouse, Rat	WB,IHC-Fr,IHC-P,IP
	SNX27 Rabbit pAb	HA500632	Human, Mouse, Rat	WB
	SQSTM1 / p62 Recombinant Antibody [PS00-61] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610409	Human, Mouse, Rat	WB,IHC-P
	SQSTM1 / p62 Recombinant Antibody [PS00-61] - Guinea pig IgG2 (Chimeric)	HA601609	Human, Mouse, Rat	WB,IHC-P
	SQSTM1 / p62 Recombinant Antibody [PS00-61] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610410	Human, Mouse, Rat	WB,IHC-P
	SQSTM1 / p62 Recombinant Antibody [PS00-61] - Rat IgG1 (Chimeric)	HA601610	Human, Mouse, Rat	WB,IHC-P
	SRSF2 Recombinant Rabbit mAb [PSH22-77] - BSA and Azide free	HA751878	Human, Mouse, Rat	IHC-P,IF-Cell,FC
	SST/Somatostatin 28 Recombinant Antibody [JE43-52] - Rat IgG1 (Chimeric)	HA601619	Human, Mouse, Rat	IHC-Fr,IHC-P
	SST/Somatostatin 28 Recombinant Antibody [JE43-52] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610417	Human, Mouse, Rat	IHC-Fr,IHC-P
Relative expression	STEAP1 Recombinant Rabbit mAb [PSH23-51]	HA724382	Human	WB,IF-Cell,IHC-P,FC
Relative expression	STEAP1 Recombinant Rabbit mAb [PSH23-51] - BSA and Azide free	HA751899	Human	WB,IF-Cell,IHC-P,FC
	STX17 Recombinant Rabbit mAb [PSH23-50] - BSA and Azide free	HA751898	Human, Mouse, Rat, Green monkey	WB,IHC-P,FC
	STX17 Recombinant Rabbit mAb [PSH23-50]	HA724381	Human, Mouse, Rat, Green monkey	WB,IHC-P,FC
	TATA binding protein TBP Recombinant Rabbit mAb [PSH22-13] - BSA and Azide free	HA751855	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC,IP
	T-bet / Tbx21 Recombinant Rabbit mAb [PSH23-00] - BSA and Azide free	HA751883	Human	IHC-P,IF-Tissue
Relative expression	TCF7 Recombinant Rabbit mAb [PSH22-98] - BSA and Azide free	HA751881	Human	WB,IHC-P,IF-Tissue,IF-Cell,IP
	TCPTP Recombinant Rabbit mAb [PSH24-14]	HA724425	Human, Mouse, Rat, Green monkey	WB,IHC-P
	TCPTP Recombinant Rabbit mAb [PSH24-14] - BSA and Azide free	HA751928	Human, Mouse, Rat, Green monkey	WB,IHC-P
	Thioredoxin Recombinant Rabbit mAb [PSH23-52]	HA724383	Human, Mouse, Rat, Green monkey	WB,IP

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Thioredoxin Recombinant Rabbit mAb [PSH23-52] - BSA and Azide free	HA751900	Human, Mouse, Rat, Green monkey	WB,IP
Relative expression, Cell treatment	Thrombomodulin Recombinant Rabbit mAb [PSH23-10] - BSA and Azide free	HA751888	Human	WB,IF-Cell,IHC-P,FC,IP
	Toll-like Receptor Antibody Sampler Kit	HAK21199	Human	WB
Relative expression	TROP2 Recombinant Rabbit mAb [PSH22-14] - BSA and Azide free	HA751856	Human, Mouse, Rat	WB,IF-Cell,IHC-Fr,IHC-P,IF-Tissue
Relative expression	TRP2 / DCT Recombinant Rabbit mAb [PSH22-31] - BSA and Azide free	HA751859	Human	WB,IP
KD	TRPA1 Recombinant Rabbit mAb [PSH23-56] - BSA and Azide free	HA751875	Human, Green monkey	WB
Cell treatment	XBP1 Recombinant Rabbit mAb [PSH23-23] - BSA and Azide free	HA751891	Human, Mouse, Rat	WB
Cell treatment	XBP1 Recombinant Rabbit mAb [PSH23-23]	HA724358	Human, Mouse, Rat	WB

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层