

超敏ELISA Kit-产品清单

产品名	货号	反应种属	灵敏度
Human IFN- γ Super Sensitivity enzyme-linked immunoassay kit	EHSS0039	Human	0.02 pg/ml (50 μ l) ; 0.17 pg/ml (10 μ l)
Human IL-6 Super Sensitivity enzyme-linked immunoassay kit	EHSS0001	Human	0.11 pg/ml (50 μ l) ; 0.33 pg/ml (10 μ l)
Mouse IFN- γ Super Sensitivity enzyme-linked immunoassay kit	EMSS0009	Mouse	0.38 pg/ml (10 μ l)
Mouse IL-1 beta Super Sensitivity enzyme-linked immunoassay kit	EMSS0029	Mouse	0.12 pg/ml (50 μ l) ; 0.36 pg/ml (10 μ l)