

## 2026年5月—新上市产品清单

高级验证	产品名	货号	反应种属	验证应用
	HRP Conjugated 6X His tag Recombinant Rabbit mAb [PSH07-10]	HA724498H	Species independent	WB
KO	ADNP Recombinant Rabbit mAb [PSH24-23]	HA724430	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Cell,IP,ChIP
KO	ADNP Recombinant Rabbit mAb [PSH24-23] - BSA and Azide free	HA751933	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Cell,IP,ChIP
Relative expression	Aggrecan Recombinant Rabbit mAb [PSH24-95]	HA724490	Human, Mouse, Rat	IHC-P,IF-Tissue
Relative expression	Aggrecan Recombinant Rabbit mAb [PSH24-95] - BSA and Azide free	HA751961	Human, Mouse, Rat	IHC-P,IF-Tissue
Relative expression	Albumin Recombinant Rabbit mAb [PSH24-72]	HA724468	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,FC(Intra)
Relative expression	Albumin Recombinant Rabbit mAb [PSH24-72] - BSA and Azide free	HA751950	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,FC(Intra)
	Albumin from human serum Protein	HA211532	Human	
Relative expression	ALK Recombinant Rabbit mAb [PSH24-97]	HA724492	Human	WB,IHC-P
Relative expression	ALK Recombinant Rabbit mAb [PSH24-98]	HA724493	Human	WB,IHC-P
	Apolipoprotein CIII Recombinant Rabbit mAb [PSH24-47] - BSA and Azide free	HA751943	Human	WB,IHC-P
	Apolipoprotein CIII Recombinant Rabbit mAb [PSH24-47]	HA724449	Human	WB,IHC-P
	Assay Accessory Kit	DSFZ01		ELISA
	ATG13 Rabbit pAb	HA500639	Human, Mouse, Rat	WB
	CACNA1C Recombinant Rabbit mAb [PSH24-46]	HA724448	Mouse, Rat	IHC-Fr,IHC-P
	CACNA1C Recombinant Rabbit mAb [PSH24-46] - BSA and Azide free	HA751942	Mouse, Rat	IHC-Fr,IHC-P
Relative expression	Cadherin 3 Recombinant Rabbit mAb [PSH24-84] - BSA and Azide free	HA751960	Human, Mouse, Rat	WB,IHC-P,FC
Relative expression	Cadherin 3 Recombinant Rabbit mAb [PSH24-84]	HA724479	Human, Mouse, Rat	WB,IHC-P,FC
	Calpain 1 Recombinant Rabbit mAb [PSH24-37] - BSA and Azide free	HA751938	Human, Mouse, Rat, Green monkey	WB
	Calpain 1 Recombinant Rabbit mAb [PSH24-37]	HA724444	Human, Mouse, Rat, Green monkey	WB

Relative expression	CCL5 Recombinant Rabbit mAb [PSH24-74]	HA724471	Human	WB,IF-Cell,FC
Relative expression	CCL5 Recombinant Rabbit mAb [PSH24-74] - BSA and Azide free	HA751952	Human	WB,IF-Cell,FC
Relative expression	CD25 / IL-2 Receptor alpha Recombinant Rabbit mAb [PSH21-78]	HA724249	Mouse	FC
	PE Conjugated CD27 Recombinant Rabbit mAb [PSH23-81]	HA724403F	Mouse	FC
	Cell Junction Markers Antibody Sampler Kit	HAK21200	Human	WB,IHC-P
Relative expression	CEMIP Recombinant Rabbit mAb [PSH24-78] - BSA and Azide free	HA751956	Mouse, Rat, Human	WB,IP
Relative expression	CEMIP Recombinant Rabbit mAb [PSH24-78]	HA724475	Mouse, Rat, Human	WB,IP
Relative expression	Chil3 / Ym-1 Recombinant Rabbit mAb [PSH24-49] - BSA and Azide free	HA751945	Mouse	WB,IHC-P,IF-Tissue
Relative expression	Chil3 / Ym-1 Recombinant Rabbit mAb [PSH24-49]	HA724451	Mouse	WB,IHC-P,IF-Tissue
	Choline Acetyltransferase Recombinant Antibody [PSH15-05] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610431	Mouse	IHC-Fr
	Choline Acetyltransferase Recombinant Antibody [PSH15-05] - Guinea pig IgG2 (Chimeric)	HA601633	Mouse	IHC-Fr
Relative expression, KD	Claudin 3 Recombinant Rabbit mAb [PSH24-82] - BSA and Azide free	HA751958	Human, Mouse, Rat	WB,IF-Tissue,IP
Relative expression, KD	Claudin 3 Recombinant Rabbit mAb [PSH24-82]	HA724477	Human, Mouse, Rat	WB,IF-Tissue,IP
Relative expression	CRTC1 Recombinant Rabbit mAb [PSH25-02] - BSA and Azide free	HA751965	Human, Mouse, Rat	WB,FC(Intra),IP
Relative expression	CRTC1 Recombinant Rabbit mAb [PSH25-02]	HA724497	Human, Mouse, Rat	WB,FC(Intra),IP
Relative expression	CSPG4 Recombinant Rabbit mAb [PSH24-56] - BSA and Azide free	HA751946	Human, Mouse, Rat	WB,IF-Cell,FC,IP
Relative expression	CSPG4 Recombinant Rabbit mAb [PSH24-57] - BSA and Azide free	HA751947	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,FC
Relative expression	CSPG4 Recombinant Rabbit mAb [PSH24-57]	HA724453	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue,FC
Relative expression	CSPG4 Recombinant Rabbit mAb [PSH24-56]	HA724452	Human, Mouse, Rat	WB,IF-Cell,FC,IP
	Ctip2 Recombinant Antibody [JF09-90] - Guinea pig IgG2 (Chimeric)	HA601623	Human, Mouse, Rat	WB,IHC-Fr,IHC-P

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Ctip2 Recombinant Antibody [JF09-90] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610421	Human, Mouse, Rat	WB,IHC-Fr,IHC-P
	Ctip2 Recombinant Antibody [JF09-90] - Rat IgG1 (Chimeric)	HA601624	Human, Mouse, Rat	WB,IHC-Fr,IHC-P
	Ctip2 Recombinant Antibody [JF09-90] - Chicken IgY (Chimeric) - BSA and Azide free	HA610420	Human, Mouse	WB,IHC-Fr
	Ctip2 Recombinant Antibody [JF09-90] - Chicken IgY (Chimeric)	HA601622	Human, Mouse	WB,IHC-Fr
	Ctip2 Recombinant Antibody [JF09-90] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610422	Human, Mouse, Rat	WB,IHC-Fr,IHC-P
Cell treatment	CXCL11 Recombinant Rabbit mAb [PSH24-81] - BSA and Azide free	HA751957	Human	WB,IF-Cell,IP
Cell treatment	CXCL11 Recombinant Rabbit mAb [PSH24-81]	HA724476	Human	WB,IF-Cell,IP
	Cytokeratin 19 Recombinant Antibody [SA30-06] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610428	Human, Mouse	IHC-P,IF-Cell
	Cytokeratin 19 Recombinant Antibody [SA30-06] - Mouse IgG1 (Chimeric)	HA601630	Human, Mouse	IHC-P,IF-Cell
	HRP Conjugated Cytokeratin 7 Recombinant Mouse mAb [3-G3-D8-R]	HA601065H	Human	IHC-P
	Disulfidptosis Antibody Sampler Kit	HAK21206	Human, Mouse, Rat	WB,IHC-P,IF-Cell
Relative expression	DNAJC5 Recombinant Rabbit mAb [PSH24-38] - BSA and Azide free	HA751939	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,IHC-Fr,IHC-P,FC(Intra),IP
Relative expression	DNAJC5 Recombinant Rabbit mAb [PSH24-38]	HA724445	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,IHC-Fr,IHC-P,FC(Intra),IP
Relative expression	Endomucin Recombinant Rabbit mAb [PSH24-71]	HA724467	Mouse	WB,IHC-P,IF-Tissue,mIHC
Relative expression	Endomucin Recombinant Rabbit mAb [PSH24-71] - BSA and Azide free	HA751949	Mouse	WB,IHC-P,IF-Tissue
	EPO Recombinant Rabbit mAb [PSH24-53] - BSA and Azide free (Detector)	HA725352	Human, Mouse, Rat	ELISA(Det)
	EPO Recombinant Rabbit mAb [PSH24-52] - BSA and Azide free (Capture)	HA725351	Human, Mouse, Rat	ELISA(Cap)
	Biotin Conjugated EPO Recombinant Rabbit mAb [PSH24-53] - Detector	HA725352B	Human, Mouse, Rat	ELISA(Det),ELISA
	iFluor™ 647 Conjugated Estrogen Receptor Recombinant Rabbit mAb [PD00-04]	HA723312F	Human	IF-Cell,IF-Tissue
Relative expression	FCER1G Recombinant Rabbit mAb [PSH24-96]	HA724491	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	FCER1G Recombinant Rabbit mAb [PSH24-96] - BSA and Azide free	HA751962	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	Ganglioside GD2 Recombinant Mouse mAb [PSH24-41] - BSA and Azide free	HA610426	Human	IF-Cell,FC(Intra)
	Ganglioside GD2 Recombinant Mouse mAb [PSH24-41]	HA601628	Human	IF-Cell,FC(Intra)
	Glucose Transporter GLUT1 Recombinant Antibody [SA0377] - Mouse IgG1 (Chimeric)	HA601621	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Cell
	Glucose Transporter GLUT1 Recombinant Antibody [SA0377] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610419	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IF-Cell
	Biotin Conjugated Human APCS Recombinant Rabbit mAb [PSH24-40] - Detector	HA725491B	Human	ELISA(Det),ELISA
	Human APCS Recombinant Rabbit mAb [PSH24-39] - BSA and Azide free (Capture)	HA725490	Human	ELISA(Cap)
	HRP Conjugated Human APCS Recombinant Rabbit mAb [PSH24-40] - Detector	HA725491H	Human	ELISA(Det),ELISA
	Human APCS Recombinant Rabbit mAb [PSH24-40] - BSA and Azide free (Detector)	HA725491	Human	ELISA(Det)
	Human APRIL/TNFSF13, N-hFc Tag Protein	HA211550	Human	
	Human Ataxin-3/ATXN3, C-His Tag Protein	HA211529	Human	
	Human Bikunin/AMBP, C-His Tag Protein	HA211461	Human	
	Human BRAK/CXCL14, C-hFc Tag Protein	HA211551	Human	
	Biotin Conjugated Human Caspase-1 Recombinant Rabbit mAb [PSH24-90] - Detector	HA724485B	Human	ELISA(Det),ELISA
	Human Caspase-1 Recombinant Rabbit mAb [PSH24-91] - BSA and Azide free (Capture)	HA724486	Human	ELISA(Cap)
	Human Caspase-1 Recombinant Rabbit mAb [PSH24-89] - BSA and Azide free (Capture)	HA724484	Human	ELISA(Cap)
	Human Caspase-1 Recombinant Rabbit mAb [PSH24-90] - BSA and Azide free (Detector)	HA724485	Human	ELISA(Det)
	Human Caspase-9 enzyme-linked immunoassay kit	EH0225	Human	ELISA

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human Cathepsin L / CTSL1, C-His Tag Protein	HA211475	Human	
	Human CD38 Recombinant Rabbit mAb [PSH24-45] - BSA and Azide free (Detector)	HA725477	Human	ELISA(Det)
	Biotin Conjugated Human CD38 Recombinant Rabbit mAb [PSH24-45] - Detector	HA725477B	Human	ELISA(Det),ELISA
	Human CD38 Recombinant Rabbit mAb [PSH24-44] - BSA and Azide free (Capture)	HA725476	Human	ELISA(Cap)
	iFluor™ 594 Conjugated Human CD45 Recombinant Mouse mAb [PSH03-64]	HA601608F	Human	FC
	Human CD74 / DHLAG, C-His Tag Protein	HA211536	Human	
	Human Chromogranin A, Tag Free Protein	HA211470	Human	
	Human Complement C2, C-His Tag Protein	HA211546	Human	
	HRP Conjugated Human Complement factor H Recombinant Rabbit mAb [PSH22-21] - Detector	HA724282H	Human	ELISA(Det),ELISA
	Human C-peptide enzyme-linked immunoassay kit	EH0244	Human	ELISA
	Human Cystatin B/CSTB, Tag Free Protein	HA211466	Human	
	Human DKK-3, Tag Free Protein	HA211460	Human	
	Human ESM1, Tag Free Protein	HA211552	Human	
	Human Follistatin, Tag Free Protein	HA211463	Human	
	Human Galectin-4, Tag Free Protein	HA211468	Human	
	Human HAO1, Tag Free Protein	HA211469	Human	
	Human IgE Fc, Tag Free Protein	HA211545	Human	
	Biotin Conjugated Human IL-16 Recombinant Rabbit mAb [PSH24-80] - Detector	HA725214B	Human	ELISA(Det),ELISA
	Human IL-16 Recombinant Rabbit mAb [PSH24-79] - BSA and Azide free (Capture)	HA725213	Human	ELISA(Cap)
	Human IL-16 Recombinant Rabbit mAb [PSH24-80] - BSA and Azide free (Detector)	HA725214	Human	ELISA(Det)
	Biotin Conjugated Human IL-17RA Recombinant Rabbit mAb [PSH24-66] - Detector	HA724462B	Human	ELISA(Det),ELISA

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL-17RA Recombinant Rabbit mAb [PSH24-65] - BSA and Azide free (Capture)	HA724461	Human	ELISA(Cap)
	Human IL-17RA Recombinant Rabbit mAb [PSH24-66] - BSA and Azide free (Detector)	HA724462	Human	ELISA(Det)
	Human IL-5RA Recombinant Rabbit mAb [PSH24-51] - BSA and Azide free (Detector)	HA725254	Human	ELISA(Det)
	Biotin Conjugated Human IL-5RA Recombinant Rabbit mAb [PSH24-51] - Detector	HA725254B	Human	ELISA(Det),ELISA
	Human IL-5RA Recombinant Rabbit mAb [PSH24-50] - BSA and Azide free (Capture)	HA725253	Human	ELISA(Cap)
	Biotin Conjugated Human LDL Receptor Recombinant Rabbit mAb [PSH24-60] - Detector	HA724456B	Human	ELISA(Det),ELISA
	Human LDL Receptor Recombinant Rabbit mAb [PSH24-59] - BSA and Azide free (Capture)	HA724455	Human	ELISA(Cap)
	Human LDL Receptor Recombinant Rabbit mAb [PSH24-60] - BSA and Azide free (Detector)	HA724456	Human	ELISA(Det)
	Human/Monkey Apolipoprotein CIII Recombinant Rabbit mAb [PSH24-26] - BSA and Azide free (Capture)	HA724433	Human, Rhesus monkey	ELISA(Cap)
	Biotin Conjugated Human/Monkey Apolipoprotein CIII Recombinant Rabbit mAb [PSH24-27] - Detector	HA724434B	Human, Rhesus monkey	ELISA(Det),ELISA
	Human/Monkey Apolipoprotein CIII Recombinant Rabbit mAb [PSH24-27] - BSA and Azide free (Detector)	HA724434	Human, Rhesus monkey	ELISA(Det)
	Human/Monkey Apo-M Recombinant Rabbit mAb [PSH24-28] - BSA and Azide free (Capture)	HA724435	Human, Rhesus monkey	ELISA(Cap)
	Human/Monkey Apo-M Recombinant Rabbit mAb [PSH24-29] - BSA and Azide free (Detector)	HA724436	Human, Rhesus monkey	ELISA(Det)
	Human/Monkey Apo-M Recombinant Rabbit mAb [PSH24-30] - BSA and Azide free (Capture)	HA724437	Human, Rhesus monkey	ELISA(Cap)
	Biotin Conjugated Human/Monkey Apo-M Recombinant Rabbit mAb [PSH24-29] - Detector	HA724436B	Human, Rhesus monkey	ELISA(Det),ELISA
	Human/Monkey ESM1 Recombinant Rabbit mAb [PSH24-32] - BSA and Azide free (Detector)	HA724439	Human, Rhesus monkey	ELISA(Det)

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Human/Monkey ESM1 Recombinant Rabbit mAb [PSH24-31] - BSA and Azide free (Capture)	HA724438	Human, Rhesus monkey	ELISA(Cap)
Biotin Conjugated Human/Monkey ESM1 Recombinant Rabbit mAb [PSH24-32] - Detector	HA724439B	Human, Rhesus monkey	ELISA(Det),ELISA
HRP Conjugated Human/Mouse MMP-2 Recombinant Rabbit mAb [PSH22-05] - Detector	HA725380H	Human, Mouse	ELISA(Det),ELISA
Human phosphorylated Tau181 Protein	HA211530	Human	
Human phosphorylated Tau231 Protein	HA211531	Human	
Human Pro-Collagen alpha-2 C-Terminal Propeptide Recombinant Rabbit mAb [PSH24-68] - BSA and Azide free (Detector)	HA724464	Human	ELISA(Det)
Biotin Conjugated Human Pro-Collagen alpha-2 C-Terminal Propeptide Recombinant Rabbit mAb [PSH24-68] - Detector	HA724464B	Human	ELISA(Det),ELISA
Human Pro-Collagen alpha-2 C-Terminal Propeptide Recombinant Rabbit mAb [PSH24-67] - BSA and Azide free (Capture)	HA724463	Human	ELISA(Cap)
Human Prokineticin-1/PROK1, Tag Free Protein	HA211535	Human	
Human Prorelaxin H2 / RLN2 Recombinant Rabbit mAb [PSH24-55] - BSA and Azide free (Detector)	HA725467	Human	ELISA(Det)
Human Prorelaxin H2 / RLN2 Recombinant Rabbit mAb [PSH24-54] - BSA and Azide free (Capture)	HA725466	Human	ELISA(Cap)
Biotin Conjugated Human Prorelaxin H2 / RLN2 Recombinant Rabbit mAb [PSH24-55] - Detector	HA725467B	Human	ELISA(Det),ELISA
Human Prostate specific Antigen/KLK3 enzyme-linked immunoassay kit	EH0234	Human	ELISA
Human TFF1 Recombinant Rabbit mAb [PSH24-18] - BSA and Azide free (Capture)	HA725468	Human	ELISA(Cap)
Biotin Conjugated Human TFF1 Recombinant Rabbit mAb [PSH24-19] - Detector	HA725469B	Human	ELISA(Det),ELISA
Human TFF1 Recombinant Rabbit mAb [PSH24-19] - BSA and Azide free (Detector)	HA725469	Human	ELISA(Det)
HRP Conjugated Human Thymidine Kinase 1 / TK1 Recombinant Rabbit mAb [PSH18-63] - Detector	HA725240H	Human	ELISA(Det),ELISA

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human UMOD Recombinant Rabbit mAb [PSH24-69] - BSA and Azide free (Capture)	HA724465	Human	ELISA(Cap)
	Human UMOD Recombinant Rabbit mAb [PSH24-70] - BSA and Azide free (Detector)	HA724466	Human	ELISA(Det)
	Biotin Conjugated Human UMOD Recombinant Rabbit mAb [PSH24-70] - Detector	HA724466B	Human	ELISA(Det),ELISA
	Human WIF1, Tag Free Protein	HA211467	Human	
Cell treatment	IDO Recombinant Antibody [PD00-62] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610430	Human	WB,IHC-P
Cell treatment	IDO Recombinant Antibody [PD00-62] - Mouse IgG1 (Chimeric)	HA601632	Human	WB,IHC-P
Cell treatment	iNOS Recombinant Rabbit mAb [PSH24-75]	HA724472	Human, Mouse, Rat	WB,IF-Cell,FC(Intra)
Cell treatment	iNOS Recombinant Rabbit mAb [PSH24-75] - BSA and Azide free	HA751953	Human, Mouse, Rat	WB,IF-Cell,FC(Intra)
Relative expression	ITPR3 Rabbit pAb	HA500640	Human	WB
	HRP Conjugated Ki67 Recombinant Rabbit mAb [SR00-02]	HA721115H	Mouse	IHC-P
	LAMP1 Recombinant Antibody [PSH05-55] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610418	Mouse	WB,IHC-Fr,IF-Cell
	LAMP1 Recombinant Antibody [PSH05-55] - Rat IgG1 (Chimeric)	HA601620	Mouse	WB,IHC-Fr,IF-Cell
	LRPPRC Recombinant Rabbit mAb [PSH24-58] - BSA and Azide free	HA751948	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,IHC-P,FC(Intra)
	LRPPRC Recombinant Rabbit mAb [PSH24-58]	HA724454	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,IHC-P,FC(Intra)
	iFluor™ 594 Conjugated F(ab') <sub>2</sub> Fragment Anti-Ly6g Antibody [PSH14-95]	HA724234F	Mouse	FC
	Mammalian INHBE enzyme-linked immunoassay kit (one step)	EHY0245	Mammalian	ELISA
	Mitoxyperilysis Antibody Sampler Kit	HAK21205	Human	WB
	Biotin Conjugated Mouse CD28 Recombinant Rabbit mAb [PSH24-62] - Detector	HA724458B	Mouse	ELISA(Det),ELISA
	Mouse CD28 Recombinant Rabbit mAb [PSH24-62] - BSA and Azide free (Detector)	HA724458	Mouse	ELISA(Det)

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse CD28 Recombinant Rabbit mAb [PSH24-61] - BSA and Azide free (Capture)	HA724457	Mouse	ELISA(Cap)
	Mouse EPO/Erythropoietin enzyme-linked immunoassay kit	EM0059	Mouse	ELISA
	Mouse Follitropin Subunit beta/FSHB, Tag Free Protein	HA211464	Mouse	
	Mouse Fractalkine/CX3CL1, Tag Free Protein	HA211451	Mouse	
	Mouse IFN- $\gamma$ Super Sensitivity enzyme-linked immunoassay kit	EMSS0009	Mouse	ELISA
	Mouse IL-17F enzyme-linked immunoassay kit	EM0061	Mouse	ELISA
	Mouse/Rat IL-12 p40 Recombinant Rabbit mAb [PSH24-94] - BSA and Azide free (Capture)	HA724489	Mouse, Rat	ELISA(Cap)
	Mouse/Rat IL-12 p40 Recombinant Rabbit mAb [PSH24-92] - BSA and Azide free (Capture/Detector)	HA724487	Mouse, Rat	ELISA(Cap),ELISA(Det)
	Biotin Conjugated Mouse/Rat IL-12 p40 Recombinant Rabbit mAb [PSH24-92] - Detector	HA724487B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat RAGE Recombinant Rabbit mAb [PSH24-85] - BSA and Azide free (Capture)	HA724480	Mouse, Rat	ELISA(Cap)
	Mouse/Rat RAGE Recombinant Rabbit mAb [PSH24-86] - BSA and Azide free (Detector)	HA724481	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse/Rat RAGE Recombinant Rabbit mAb [PSH24-86] - Detector	HA724481B	Mouse, Rat	ELISA(Det),ELISA
	Mouse Resistin enzyme-linked immunoassay kit ( one step )	EMY0063	Mouse	ELISA
	Mouse VE-Cadherin Recombinant Rabbit mAb [PSH24-63] - BSA and Azide free (Capture)	HA724459	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse VE-Cadherin Recombinant Rabbit mAb [PSH24-64] - Detector	HA724460B	Mouse	ELISA(Det),ELISA
	Mouse VE-Cadherin Recombinant Rabbit mAb [PSH24-64] - BSA and Azide free (Detector)	HA724460	Mouse	ELISA(Det)
Relative expression	MTBR-TAU 243 Isoform 2N4R (225-242) Recombinant Rabbit mAb [PSH24-42]	HA724446	Mouse, Rat, Pig, Human	WB,IHC-P,IP
Relative expression	MTBR-TAU 243 Isoform 2N4R (225-242) Recombinant Rabbit mAb [PSH24-42] - BSA and Azide free	HA751940	Mouse, Rat, Pig, Human	WB,IHC-P,IP

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	MTHFD2 Recombinant Rabbit mAb [PSH24-73] - BSA and Azide free	HA751951	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	MTHFD2 Recombinant Rabbit mAb [PSH24-73]	HA724469	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
Relative expression	MyoD1 Recombinant Rabbit mAb [PSH24-48] - BSA and Azide free	HA751944	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
Relative expression	MyoD1 Recombinant Rabbit mAb [PSH24-48]	HA724450	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC,IP
	Neurodegenerative disease-associated synaptic plasticity Antibody Sampler Kit	HAK21203	Human, Mouse	WB,IHC-Fr
KD	Neurofibromin Recombinant Rabbit mAb [PSH25-01] - BSA and Azide free	HA751964	Human, Mouse, Rat, Green monkey	WB
KD	Neurofibromin Recombinant Rabbit mAb [PSH25-01]	HA724496	Human, Mouse, Rat, Green monkey	WB
	NuRD Complex Antibody Sampler Kit	HAK21198	Human, Mouse	WB,IHC-P
	Olig2 Recombinant Antibody [SP07-02] - Mouse IgG1 (Chimeric)	HA601629	Human, Mouse, Rat	IHC-P,IF-Tissue
	Olig2 Recombinant Antibody [SP07-02] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610423	Human, Mouse, Rat	IHC-Fr,IHC-P
	Olig2 Recombinant Antibody [SP07-02] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610427	Human, Mouse, Rat	IHC-P,IF-Tissue
	Olig2 Recombinant Antibody [SP07-02] - Guinea pig IgG2 (Chimeric)	HA601625	Human, Mouse, Rat	IHC-Fr,IHC-P
Relative expression	p63 Recombinant Antibody [SC06-31] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610429	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
Relative expression	p63 Recombinant Antibody [SC06-31] - Mouse IgG1 (Chimeric)	HA601631	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
Relative expression, Cell treatment	PAI1 Recombinant Rabbit mAb [PSH24-33] - BSA and Azide free	HA751934	Human, Mouse, Green monkey	WB,IP
Relative expression, Cell treatment	PAI1 Recombinant Rabbit mAb [PSH24-33]	HA724440	Human, Mouse, Green monkey	WB,IP
Relative expression	Periostin Recombinant Rabbit mAb [PSH24-36]	HA724443	Mouse, Rat	WB,IHC-Fr
Relative expression	Periostin Recombinant Rabbit mAb [PSH24-35] - BSA and Azide free	HA751936	Mouse	WB,IHC-Fr,IHC-P
Relative expression	Periostin Recombinant Rabbit mAb [PSH24-36] - BSA and Azide free	HA751937	Mouse, Rat	WB,IHC-Fr
Relative expression	Periostin Recombinant Rabbit mAb [PSH24-35]	HA724442	Mouse	WB,IHC-Fr,IHC-P

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH24-43] - BSA and Azide free	HA751941	Human, Mouse, Rat	WB,IP
Cell treatment	Phospho-AKT (S473) Recombinant Rabbit mAb [PSH24-43]	HA724447	Human, Mouse, Rat	WB,IP
Cell treatment	Phospho-RPA32/RPA2 (S4, S8) Recombinant Rabbit mAb [PSH24-99]	HA724494	Human, Mouse, Rat	WB,IF-Cell,IP
Cell treatment	Phospho-RPA32/RPA2 (S4, S8) Recombinant Rabbit mAb [PSH24-99] - BSA and Azide free	HA751963	Human, Mouse, Rat	WB,IF-Cell,IP
Cell treatment	Phospho-SAPK/JNK (T183/Y185) Recombinant Rabbit mAb [PSH24-76] - BSA and Azide free	HA751954	Human, Mouse, Rat	WB,IHC-P,IP
Cell treatment	Phospho-SAPK/JNK (T183/Y185) Recombinant Rabbit mAb [PSH24-76]	HA724473	Human, Mouse, Rat	WB,IHC-P,IP
Relative expression	RAB27A+RAB27B Recombinant Rabbit mAb [PSH24-34] - BSA and Azide free	HA751935	Human, Mouse, Rat	WB,FC(Intra),IP
Relative expression	RAB27A+RAB27B Recombinant Rabbit mAb [PSH24-34]	HA724441	Human, Mouse, Rat	WB,FC(Intra),IP
	Rabbit Anti-DYKDDDDK (Flag) Magnetic Beads	HAK21021		IP
	Rabbit IgG (isotype control) Magnetic Beads	HAK21202		IP
	Rat Albumin/ALB, Tag Free Protein	HA211547	Rat	
	Biotin Conjugated Rat IL-12 p40 Recombinant Rabbit mAb [PSH24-93] - Detector	HA724488B	Rat	ELISA(Det),ELISA
	Rat IL-12 p40 Recombinant Rabbit mAb [PSH24-93] - BSA and Azide free (Detector)	HA724488	Rat	ELISA(Det)
	Rat Prolactin Recombinant Rabbit mAb [PSH24-87] - BSA and Azide free (Capture)	HA724482	Rat	ELISA(Cap)
	Rat Prolactin Recombinant Rabbit mAb [PSH24-88] - BSA and Azide free (Detector)	HA724483	Rat	ELISA(Det)
	Biotin Conjugated Rat Prolactin Recombinant Rabbit mAb [PSH24-88] - Detector	HA724483B	Rat	ELISA(Det),ELISA
	Rat Renin, C-His Tag Protein	HA211446	Rat	
	iFluor™ 488 Conjugated Recombinant Rabbit IgG [PSH04-42] - Isotype control	HA720231F		FC
Relative expression	RYR1 Recombinant Rabbit mAb [PSH24-83] - BSA and Azide free	HA751959	Mouse, Rat	WB,IHC-Fr,IHC-P
Relative expression	RYR1 Recombinant Rabbit mAb [PSH24-83]	HA724478	Mouse, Rat	WB,IHC-Fr,IHC-P

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	SATB2 Recombinant Antibody [PD01-42] - Mouse IgG1 (Chimeric)	HA601626	Human, Mouse, Rat	IHC-Fr,IHC-P
	SATB2 Recombinant Antibody [PD01-42] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610424	Human, Mouse, Rat	IHC-Fr,IHC-P
	SMCR7L / MID51 Recombinant Rabbit mAb [PSH24-77]	HA724474	Human, Mouse, Rat, Green monkey, Pig	WB,IP
	SMCR7L / MID51 Recombinant Rabbit mAb [PSH24-77] - BSA and Azide free	HA751955	Human, Mouse, Rat, Green monkey, Pig	WB,IP
	TDP43 Recombinant Antibody [PSH09-14] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610425	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell
	TDP43 Recombinant Antibody [PSH09-14] - Rat IgG1 (Chimeric)	HA601627	Human, Mouse, Rat	WB,IHC-P,IF-Tissue,IF-Cell
	Thyroid Cancer Multi - target Differential Diagnosis IHC Antibody Sampler Kit	HAK21204	Human	IHC-P,IF-Tissue
	Toll-like Receptor Antibody Sampler Kit	HAK21199	Human	WB
Relative expression	TRP2 / DCT Recombinant Rabbit mAb [PSH25-05] - BSA and Azide free	HA751966	Human, Mouse	WB,IHC-P,FC(Intra),IP
Relative expression	TRP2 / DCT Recombinant Rabbit mAb [PSH25-05]	HA724499	Human, Mouse	WB,IHC-P,FC(Intra),IP
	HRP Conjugated Vinculin Recombinant Rabbit mAb [JM42-43]	HA724470H	Human, Mouse, Rat, Green monkey, Pig	WB

## 杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: [market@huabio.cn](mailto:market@huabio.cn)

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层