

2026年6月—新上市产品清单

高级验证	产品名	货号	反应种属	验证应用
	ACTN3 Recombinant Rabbit mAb [PSH25-74]	HA724565	Human, Mouse, Rat	IHC-P
	Anti-Biotin, AlpHcAbs® Rabbit antibody	078-201-001		ICC/IF, Flow Cyt, WB, ELISA, IP
	Anti-Biotin, AlpHcAbs® Rabbit antibody (HRP)	078-201-005		WB, ELISA
	Anti-Human IgG (CH2 Fragment specific), AlpHcAbs® Mouse antibody	023-312-001	Human	WB, ELISA, IP, ICC/IF, Flow Cyt
	Anti-Human IgG (Fcγ Fragment specific), AlpHcAbs® Goat antibody (iFluor488)	023-401-007	Human	WB, ELISA, ICC/IF, Flow Cyt
	Anti-Human IgG (Fcγ Fragment specific), AlpHcAbs® Goat antibody (iFluor647)	023-401-009	Human	WB, ELISA, ICC/IF, Flow Cyt
	Anti-Human IgG (Fcγ Fragment specific), AlpSdAbs® VHH (iFluor488)	023-101-007	Human	WB, ELISA, ICC/IF, Flow Cyt
	Anti-Human IgG (Fcγ Fragment specific), AlpSdAbs® VHH (iFluor647)	023-101-009	Human	WB, ELISA, ICC/IF, Flow Cyt
	Anti-Human IgG (Fcγ Fragment specific), AlpSdAbs® VHH (PE)	023-101-010	Human	Flow Cyt
	Anti-Human IgG (H+L), AlpHcAbs® Goat antibody (iFluor488)	023-406-007	Human	WB, ELISA, ICC/IF, Flow Cyt
	Anti-Human IgG (H+L), AlpHcAbs® Goat antibody (iFluor647)	023-406-009	Human	WB, ELISA, ICC/IF, Flow Cyt
	Anti-Human IgG (H+L), AlpSdAbs® VHH (iFluor488)	023-106-007	Human	WB, ELISA, ICC/IF, Flow Cyt
	Anti-Human IgG (H+L), AlpSdAbs® VHH (iFluor594)	023-106-008	Human	WB, ELISA, ICC/IF, Flow Cyt
	Anti-Human IgG (H+L), AlpSdAbs® VHH (iFluor647)	023-106-009	Human	WB, ELISA, ICC/IF, Flow Cyt
	Anti-Rabbit IgG (Fcγ Fragment specific), AlpHcAbs® Goat antibody (iFluor647)	025-401-009	Rabbit	ELISA, WB, ICC/IF, Flow Cyt
	Anti-Rabbit IgG (H+L), AlpSdAbs® VHH (APC)	025-102-011	Rabbit	ICC/IF, Flow Cyt
	BAF Complex Antibody Sampler Kit	HAK21077	Human	WB
	Cannabinoid Receptor 1 Downstream Signaling Antibody Sampler Kit	HAK21139	Human, Mouse	WB
	CD16 Recombinant Antibody [PD00-12] - Mouse IgG1 (Chimeric)	HA601634	Human	IHC-P, IF-Tissue

	iFluor™ 647 Conjugated CD63 Recombinant Rabbit mAb [SY21-02]	HA724566F	Human	IF-Tissue
Relative expression	CD81 Recombinant Rabbit mAb [PSH10-26] - BSA and Azide free	HA751976	Human, Rat	WB,IHC-P,FC
Relative expression	CD86 Recombinant Rabbit mAb [PSH25-20]	HA724514	Mouse	WB,IF-Tissue,IHC-P
Relative expression	CD86 Recombinant Rabbit mAb [PSH25-20] - BSA and Azide free	HA751967	Mouse	WB,IF-Tissue,IHC-P
Relative expression	CerS2 Recombinant Rabbit mAb [PSH25-62] - BSA and Azide free	HA751981	Human, Mouse, Rat	WB,IF-Cell,IHC-Fr,IHC-P,IP
Relative expression	CerS2 Recombinant Rabbit mAb [PSH25-62]	HA724549	Human, Mouse, Rat	WB,IF-Cell,IHC-Fr,IHC-P,IP
Relative expression	Claudin 11 Recombinant Rabbit mAb [PSH25-21]	HA724515	Mouse, Rat	IHC-Fr,IHC-P
Relative expression	Claudin 11 Recombinant Rabbit mAb [PSH25-21] - BSA and Azide free	HA751968	Mouse, Rat	IHC-Fr,IHC-P
Relative expression	Complement factor H Recombinant Rabbit mAb [PSH25-71]	HA724558	Human, Mouse, Rat	WB,IHC-P
Relative expression	Complement factor H Recombinant Rabbit mAb [PSH25-71] - BSA and Azide free	HA751990	Human, Mouse, Rat	WB,IHC-P
Cell treatment	CXCL10 Recombinant Rabbit mAb [PSH25-64]	HA724551	Mouse	WB
Cell treatment	CXCL10 Recombinant Rabbit mAb [PSH25-64] - BSA and Azide free	HA751983	Mouse	WB
	CYP17A1 Recombinant Rabbit mAb [PSH25-58] - BSA and Azide free	HA751977	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
	CYP17A1 Recombinant Rabbit mAb [PSH25-58]	HA724532	Human, Mouse, Rat	WB,IHC-P,IF-Tissue
Relative expression	DHCR24 Recombinant Rabbit mAb [PSH25-92] - BSA and Azide free	HA752004	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC(Intra)
Relative expression	DHCR24 Recombinant Rabbit mAb [PSH25-92]	HA724584	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC(Intra)
	DNA Methylation Antibody Sampler Kit	HAK21056	Human, Rat	WB
Cell treatment	dsRNA Recombinant Mouse mAb [PSH25-25]	HA601635	Species independent	IF-Cell
Cell treatment	dsRNA Recombinant Mouse mAb [PSH25-25] - BSA and Azide free	HA610432	Species independent	IF-Cell
	iFluor™ 750 Conjugated EpCAM Recombinant Rabbit mAb [PS01-69]	HA720230F	Human	IF-Cell
	Erk1 (pT202+pY204) + Erk2 (pT185+pY187) Recombinant Rabbit mAb [PD00-38]	HA724560	Human	WB,IHC-P,IF-Cell
	Exosomal Marker Antibody Sampler Kit	HAK21118	Human	WB

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

KO	Ferritin Heavy Chain Recombinant Rabbit mAb [SC0620] - BSA and Azide free	HA751971	Human, Mouse, Zebrafish, Rat	WB,IF-Cell,IF-Tissue,IHC-P,FC,IHC-Fr
	GABAergic Neuron Identification Antibody Sampler Kit	HAK21161	Mouse, Rat	IHC-P,IHC-Fr
Relative expression	GATA3 Recombinant Rabbit mAb [PSH25-24]	HA724518	Human	WB,IF-Cell,IHC-P
	Gata6 Recombinant Rabbit mAb [JE31-94]	HA724563	Human	IHC-P,FC
	GLA Recombinant Rabbit mAb [PSH25-93]	HA724585	Human	WB,IF-Cell,IHC-P,FC
	GLA Recombinant Rabbit mAb [PSH25-93] - BSA and Azide free	HA752005	Human	WB,IF-Cell,IHC-P,FC
	GM130 Recombinant Rabbit mAb [JE42-53] - BSA and Azide free	HA751992	Human, Mouse	WB,IHC-P,IF-Cell
	iFluor™ 647 Conjugated Goat anti-Rat IgG(H+L) Goat pAb	HA1153	Rat	IF-Tissue,IHC-Fr,IF-Cell
	GPI Recombinant Rabbit mAb [PSH25-26] - BSA and Azide free	HA751972	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
	GPI Recombinant Rabbit mAb [PSH25-26]	HA724519	Human, Mouse, Rat, Green monkey	WB,IHC-P,IP
	GRP94 Recombinant Rabbit mAb [PSH25-59] - BSA and Azide free	HA751978	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,IHC-P,FC(Intra),IP
	GRP94 Recombinant Rabbit mAb [PSH25-59]	HA724533	Human, Mouse, Rat, Green monkey, Pig	WB,IF-Cell,IHC-P,FC(Intra),IP
	GSK-3 Antibody Sampler Kit	HAK21153	Human, Mouse	WB,IHC-P
Relative expression	HCN4 Recombinant Rabbit mAb [PSH25-88]	HA724580	Mouse	WB
Relative expression	HCN4 Recombinant Rabbit mAb [PSH25-88] - BSA and Azide free	HA752000	Mouse	WB
Cell treatment, Relative expression	Histone H3.3 (mutated G34R) Recombinant Rabbit mAb [PSH25-65]	HA724552	Human	WB,IHC-P
Cell treatment, Relative expression	Histone H3.3 (mutated G34R) Recombinant Rabbit mAb [PSH25-65] - BSA and Azide free	HA751984	Human	WB,IHC-P
Cell treatment, Relative expression	Histone H3.3 (mutated G34V) Recombinant Rabbit mAb [PSH25-66] - BSA and Azide free	HA751985	Human	WB,IF-Cell,IHC-P,IP
Cell treatment, Relative expression	Histone H3.3 (mutated G34V) Recombinant Rabbit mAb [PSH25-66]	HA724553	Human	WB,IF-Cell,IHC-P,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human ACE2 enzyme-linked immunoassay kit	EH0224	Human	ELISA
	Human Adenosine deaminase/ADA, Tag Free Protein	HA211502	Human	
	Human Angiotensinogen enzyme-linked immunoassay kit	EH0239	Human	ELISA
	Human Apolipoprotein B-100/APOB, Tag Free	HA211573	Human	
	Human Autotaxin/ENPP-2, C-His Tag Protein	HA211477	Human	
	Human AZGP1, C-His Tag Protein	HA211513	Human	
	Human CCL26 Recombinant Rabbit mAb [PSH25-04] - BSA and Azide free (Detector)	HA725202	Human	ELISA(Det)
	Biotin Conjugated Human CCL26 Recombinant Rabbit mAb [PSH25-04] - Detector	HA725202B	Human	ELISA(Det),ELISA
	Human CCL26 Recombinant Rabbit mAb [PSH25-03] - BSA and Azide free (Capture)	HA725201	Human	ELISA(Cap)
	Human CD14, Tag Free Protein	HA211512	Human	
	Human CEACAM-6/CD66c, C-His Tag Protein	HA211497	Human	
	Human Complement C4-A/C4A, C-His Tag Protein	HA211575	Human	
	Human Complement factor B enzyme-linked immunoassay kit(one step)	EHY0238	Human	ELISA
	Human Cystatin F/CST7, C-His Tag Protein	HA211557	Human	
	Biotin Conjugated Human DLL3 Recombinant Rabbit mAb [PSH25-30] - Detector	HA724521B	Human	ELISA(Det),ELISA
	Human DLL3 Recombinant Rabbit mAb [PSH25-29] - BSA and Azide free (Capture)	HA724520	Human	ELISA(Cap)
	Human DLL3 Recombinant Rabbit mAb [PSH25-30] - BSA and Azide free (Detector)	HA724521	Human	ELISA(Det)
	Human EWSR1, C-His Tag Protein	HA211507	Human	
	Human FABP4 Recombinant Rabbit mAb [PSH25-37] - BSA and Azide free (Capture)	HA725472	Human	ELISA(Cap)
	Human FABP4 Recombinant Rabbit mAb [PSH25-38] - BSA and Azide free (Detector)	HA725473	Human	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human FABP4 Recombinant Rabbit mAb [PSH25-38] - Detector	HA725473B	Human	ELISA(Det),ELISA
	HRP Conjugated Human FGF-21 Recombinant Rabbit mAb [PSH22-33] - Detector	HA725236H	Human	ELISA(Det),ELISA
	Biotin Conjugated Human FLT3 Recombinant Rabbit mAb [PSH25-36] - Detector	HA724527B	Human	ELISA(Det),ELISA
	Human FLT3 Recombinant Rabbit mAb [PSH25-36] - BSA and Azide free (Detector)	HA724527	Human	ELISA(Det)
	Human FLT3 Recombinant Rabbit mAb [PSH25-35] - BSA and Azide free (Capture)	HA724526	Human	ELISA(Cap)
	Biotin Conjugated Human GATA4 Recombinant Rabbit mAb [PSH25-52] - Detector	HA724543B	Human	ELISA(Det),ELISA
	Human GATA4 Recombinant Rabbit mAb [PSH25-51] - BSA and Azide free (Capture)	HA724542	Human	ELISA(Cap)
	Human GATA4 Recombinant Rabbit mAb [PSH25-52] - BSA and Azide free (Detector)	HA724543	Human	ELISA(Det)
	Human Heparanase/HPSE, Tag Free Protein	HA211559	Human	
	Human HMGB1, C-His Tag Protein	HA211570	Human	
	Human HSP60/HSPD1, Tag Free Protein	HA211499	Human	
	Human HSPA1L, C-His Tag Protein	HA211506	Human	
	Human IAPP/Amylin, Tag Free Protein	HA211569	Human	
	Human ICAM2 Recombinant Rabbit mAb [PSH25-48] - BSA and Azide free (Detector)	HA724539	Human	ELISA(Det)
	Biotin Conjugated Human ICAM2 Recombinant Rabbit mAb [PSH25-48] - Detector	HA724539B	Human	ELISA(Det),ELISA
	Human ICAM2 Recombinant Rabbit mAb [PSH25-47] - BSA and Azide free (Capture)	HA724538	Human	ELISA(Cap)
	Human IgG Recombinant Rabbit mAb [SJ30-56]	HA722858	Human	IHC-P
	Human IgG1 Fc, C-His Tag Protein	HA211534	Human	
	Human IL1R1 Recombinant Rabbit mAb [PSH25-43] - BSA and Azide free (Capture)	HA725343	Human	ELISA(Cap)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human IL1R1 Recombinant Rabbit mAb [PSH25-44] - BSA and Azide free (Detector)	HA725344	Human	ELISA(Det)
	Biotin Conjugated Human IL1R1 Recombinant Rabbit mAb [PSH25-44] - Detector	HA725344B	Human	ELISA(Det),ELISA
	Human IL1RL1/DER4/ST2 enzyme-linked immunoassay kit	EH0235	Human	ELISA
	Human Kallikrein-14, Tag Free Protein	HA211508	Human	
	Human Ki67 Recombinant Rabbit mAb [PSH25-46] - BSA and Azide free (Capture)	HA724536	Human	ELISA(Cap)
	Biotin Conjugated Human Ki67 Recombinant Rabbit mAb [SR00-02] - Detector	HA724537B	Human	ELISA(Det),ELISA
	Human Ki67 Recombinant Rabbit mAb [SR00-02] - BSA and Azide free (Detector)	HA724537	Human	ELISA(Det)
	Human/Monkey Albumin enzyme-linked immunoassay kit(one step)	EHY0251	Human	ELISA
	Human/Monkey Alpha-1-Microglobulin Recombinant Rabbit mAb [PSH25-81] - BSA and Azide free (Capture)	HA724573	Human, Cynomolgus monkey	ELISA(Cap)
	Human/Monkey Alpha-1-Microglobulin Recombinant Rabbit mAb [PSH25-82] - BSA and Azide free (Detector)	HA724574	Human, Cynomolgus monkey	ELISA(Det)
	Biotin Conjugated Human/Monkey Alpha-1-Microglobulin Recombinant Rabbit mAb [PSH25-82] - Detector	HA724574B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	HRP Conjugated Human/Monkey Apolipoprotein CIII Recombinant Rabbit mAb [PSH24-27] - Detector	HA724434H	Human, Rhesus monkey	ELISA(Det),ELISA
	Biotin Conjugated Human/Monkey IGFBP2 Recombinant Rabbit mAb [PSH25-34] - Detector	HA724525B	Human, Cynomolgus monkey	ELISA(Det),ELISA
	Human/Monkey IGFBP2 Recombinant Rabbit mAb [PSH25-33] - BSA and Azide free (Capture)	HA724524	Human, Cynomolgus monkey	ELISA(Cap)
	Human/Monkey IGFBP2 Recombinant Rabbit mAb [PSH25-34] - BSA and Azide free (Detector)	HA724525	Human, Cynomolgus monkey	ELISA(Det)
	Human/Monkey Pro-COL1A1 Recombinant Rabbit mAb [PSH25-16] - BSA and Azide free (Detector)	HA724510	Human, Rhesus monkey	ELISA(Det)
	Biotin Conjugated Human/Monkey Pro-COL1A1 Recombinant Rabbit mAb [PSH25-16] - Detector	HA724510B	Human, Rhesus monkey	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Human/Monkey Pro-COL1A1 Recombinant Rabbit mAb [PSH25-15] - BSA and Azide free (Capture)	HA724509	Human, Rhesus monkey	ELISA(Cap)
	Human/Mouse GLP-1 Super Sensitivity enzyme-linked immunoassay kit	EHSS0248	Human	ELISA
	Biotin Conjugated Human MUC6 Recombinant Rabbit mAb [PSH25-18] - Detector	HA724512B	Human	ELISA(Det),ELISA
	Human MUC6 Recombinant Rabbit mAb [PSH25-19] - BSA and Azide free (Detector)	HA724513	Human	ELISA(Det)
	Human MUC6 Recombinant Rabbit mAb [PSH25-18] - BSA and Azide free (Detector)	HA724512	Human	ELISA(Det)
	Biotin Conjugated Human MUC6 Recombinant Rabbit mAb [PSH25-19] - Detector	HA724513B	Human	ELISA(Det),ELISA
	Human MUC6 Recombinant Rabbit mAb [PSH25-17] - BSA and Azide free (Capture)	HA724511	Human	ELISA(Cap)
	Human NEFH, Tag Free Protein	HA211576	Human	
	Human NKp46/NCR1 enzyme-linked immunoassay kit	EH0221	Human	ELISA
	Human NOV/CCN3, Tag Free Protein	HA211556	Human	
	Human OX40 Ligand/TNFSF4, Tag Free Protein	HA211558	Human	
	Human PGD2 synthase/PTGDS, C-His Tag Protein	HA211514	Human	
Cell treatment	Human Reactive Inflammatory Cytokine Antibody Sampler Kit	HAK21123	Human	WB
	Human SOD2, C-His Tag Protein	HA211476	Human	
	Human SPINT1/HAI1, Tag Free Protein	HA211496	Human	
	Human Thymidine Kinase 1/TK1, Tag Free Protein	HA211577	Human	
	Human TIA1, Tag Free Protein	HA211574	Human	
	HRP Conjugated Human UMOD Recombinant Rabbit mAb [PSH24-70] - Detector	HA724466H	Human	ELISA(Det),ELISA
	Human VAMP2, Tag Free Protein	HA211509	Human	
	Human VEGFD Recombinant Rabbit mAb [PSH25-11] - BSA and Azide free (Capture)	HA724505	Human	ELISA(Cap)
	Biotin Conjugated Human VEGFD Recombinant Rabbit mAb [PSH25-12] - Detector	HA724506B	Human	ELISA(Det),ELISA

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Biotin Conjugated Human VEGFD Recombinant Rabbit mAb [PSH25-14] - Detector	HA724508B	Human	ELISA(Det),ELISA
	Human VEGFD Recombinant Rabbit mAb [PSH25-14] - BSA and Azide free (Detector)	HA724508	Human	ELISA(Det)
	Human VEGFD Recombinant Rabbit mAb [PSH25-13] - BSA and Azide free (Capture)	HA724507	Human	ELISA(Cap)
	Human VEGFD Recombinant Rabbit mAb [PSH25-12] - BSA and Azide free (Detector)	HA724506	Human	ELISA(Det)
	Human VGF, Tag Free Protein	HA211510	Human	
Relative expression	ICAM2 Recombinant Rabbit mAb [PSH25-72]	HA724559	Human, Mouse	WB,IF-Cell,FC
Relative expression	ICAM2 Recombinant Rabbit mAb [PSH25-72] - BSA and Azide free	HA751991	Human, Mouse	WB,IF-Cell,FC
Cell treatment	IL-11 Rabbit pAb	HA500641	Human	WB
	Iron Homeostasis/Ferroptosis Antibody Sampler Kit	HAK21122	Human	WB
Relative expression	MAPRE1 / EB1 Recombinant Rabbit mAb [PSH25-89] - BSA and Azide free	HA752001	Human, Mouse, Rat, Pig, Green monkey	WB,IF-Cell,IHC-P,FC(Intra),IP
Relative expression	MAPRE1 / EB1 Recombinant Rabbit mAb [PSH25-89]	HA724581	Human, Mouse, Rat, Pig, Green monkey	WB,IF-Cell,IHC-P,FC(Intra),IP
	Mesenchymal Stromal Cell Marker	HAK21111	Human	WB,IHC-P,IF-Tissue
	Methyl-Histone H3 (Lys9) Antibody Sampler Kit	HAK21069	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IF-Tissue
	Mouse Angiotensinogen/AGT enzyme-linked immunoassay kit(one step)	EMY0055	Mouse	ELISA
	Mouse Apolipoprotein B enzyme-linked immunoassay kit (one step)	EMY0069	Mouse	ELISA
	Mouse Apolipoprotein CIII Recombinant Rabbit mAb [PSH25-27] - BSA and Azide free (Capture)	HA725155	Mouse	ELISA(Cap)
	Mouse Apolipoprotein CIII Recombinant Rabbit mAb [PSH25-28] - BSA and Azide free (Detector)	HA725156	Mouse	ELISA(Det)
	Biotin Conjugated Mouse Apolipoprotein CIII Recombinant Rabbit mAb [PSH25-28] - Detector	HA725156B	Mouse	ELISA(Det),ELISA
	Mouse CCL22/MDC enzyme-linked immunoassay kit	EM0054	Mouse	ELISA
	Mouse CD48 Recombinant Rabbit mAb [PSH25-10] - BSA and Azide free (Detector)	HA724504	Mouse	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse CD48 Recombinant Rabbit mAb [PSH25-09] - BSA and Azide free (Capture)	HA724503	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse CD48 Recombinant Rabbit mAb [PSH25-10] - Detector	HA724504B	Mouse	ELISA(Det),ELISA
	Biotin Conjugated Mouse CTGF Recombinant Rabbit mAb [PSH25-95] - Detector	HA725153B	Mouse	ELISA(Det),ELISA
	Mouse CTGF Recombinant Rabbit mAb [PSH25-94] - BSA and Azide free (Capture)	HA725152	Mouse	ELISA(Cap)
	Mouse CTGF Recombinant Rabbit mAb [PSH25-95] - BSA and Azide free (Detector)	HA725153	Mouse	ELISA(Det)
	Mouse Insulin enzyme-linked immunoassay kit(one step)	EMY0067	Mouse	ELISA
	Mouse Insulin Super Sensitivity enzyme-linked immunoassay kit	EMSS0067	Mouse	ELISA
	Mouse/Rat Complement C3 enzyme-linked immunoassay kit(one step)	EMY0057	Mouse	ELISA
	Biotin Conjugated Mouse/Rat DLL4 Recombinant Rabbit mAb [PSH25-32] - Detector	HA724523B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat DLL4 Recombinant Rabbit mAb [PSH25-31] - BSA and Azide free (Capture)	HA724522	Mouse, Rat	ELISA(Cap)
	Mouse/Rat DLL4 Recombinant Rabbit mAb [PSH25-32] - BSA and Azide free (Detector)	HA724523	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse/Rat Myeloperoxidase Recombinant Rabbit mAb [PSH25-80] - Detector	HA724572B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat Myeloperoxidase Recombinant Rabbit mAb [PSH25-79] - BSA and Azide free (Capture)	HA724571	Mouse, Rat	ELISA(Cap)
	Mouse/Rat Myeloperoxidase Recombinant Rabbit mAb [PSH25-80] - BSA and Azide free (Detector)	HA724572	Mouse, Rat	ELISA(Det)
	Biotin Conjugated Mouse/Rat NFkB1 (p105/p50) Recombinant Rabbit mAb [PSH25-78] - Detector	HA724570B	Mouse, Rat	ELISA(Det),ELISA
	Mouse/Rat NFkB1 (p105/p50) Recombinant Rabbit mAb [PSH25-77] - BSA and Azide free (Capture)	HA724569	Mouse, Rat	ELISA(Cap)
	Mouse/Rat NFkB1 (p105/p50) Recombinant Rabbit mAb [PSH25-78] - BSA and Azide free (Detector)	HA724570	Mouse, Rat	ELISA(Det)

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

	Mouse/Rat Prealbumin/Transthyretin enzyme-linked immunoassay kit(one step)	EMY0066	Mouse	ELISA
	HRP Conjugated Mouse/Rat Von Willebrand Factor Recombinant Rabbit mAb [PSH23-16] - Detector	HA725159H	Mouse, Rat	ELISA(Det),ELISA
	Mouse Reactive PANoptosis Antibody Sampler Kit	HAK21133	Mouse	WB
	Mouse Reactive Senescence Associated Secretory Phenotype (SASP) Antibody Sampler Kit	HAK21136	Mouse	WB
	Mouse Reactive Senescence Marker Antibody Sampler Kit	HAK21192	Mouse	WB
	Biotin Conjugated Mouse SAA Recombinant Rabbit mAb [PSH25-07] - Detector	HA724501B	Mouse	ELISA(Det),ELISA
	Mouse SAA Recombinant Rabbit mAb [PSH25-07] - BSA and Azide free (Detector)	HA724501	Mouse	ELISA(Det)
	Mouse SAA Recombinant Rabbit mAb [PSH25-06] - BSA and Azide free (Capture)	HA724500	Mouse	ELISA(Cap)
	Mouse SAA1 Recombinant Rabbit mAb [PSH25-08] - BSA and Azide free (Capture)	HA724502	Mouse	ELISA(Cap)
	Biotin Conjugated Mouse Thrombopoietin Recombinant Rabbit mAb [PSH25-50] - Detector	HA724541B	Mouse	ELISA(Det),ELISA
	Mouse Thrombopoietin Recombinant Rabbit mAb [PSH25-49] - BSA and Azide free (Capture)	HA724540	Mouse	ELISA(Cap)
	Mouse Thrombopoietin Recombinant Rabbit mAb [PSH25-50] - BSA and Azide free (Detector)	HA724541	Mouse	ELISA(Det)
	HRP Conjugated Mouse VE-Cadherin Recombinant Rabbit mAb [PSH24-64] - Detector	HA724460H	Mouse	ELISA(Det),ELISA
	Mouse VEGFR1 Recombinant Rabbit mAb [PSH25-99] - BSA and Azide free (Detector)	HA724589	Mouse	ELISA(Det)
	Biotin Conjugated Mouse VEGFR1 Recombinant Rabbit mAb [PSH25-99] - Detector	HA724589B	Mouse	ELISA(Det),ELISA
	Mouse VEGFR1 Recombinant Rabbit mAb [PSH25-98] - BSA and Azide free (Capture)	HA724588	Mouse	ELISA(Cap)
Relative expression	MSR1 / CD204 Recombinant Rabbit mAb [PSH25-57] - BSA and Azide free	HA751975	Mouse, Rat	WB,IF-Cell,FC,IP

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	MSR1 / CD204 Recombinant Rabbit mAb [PSH25-57]	HA724531	Mouse, Rat	WB,IF-Cell,FC,IP
Relative expression	NDRG2 Recombinant Rabbit mAb [PSH25-90]	HA724582	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
Relative expression	NDRG2 Recombinant Rabbit mAb [PSH25-90] - BSA and Azide free	HA752002	Human, Mouse, Rat	WB,IF-Cell,IHC-P,IP
	Neural Progenitor Cell Marker Antibody Sampler Kit	HAK21174	Human, Mouse	WB,IHC-Fr
	Nipah virus Glycoprotein G, C-His Tag	HA211571		
	NRAS (mutated Q61R) Recombinant Rabbit mAb [PD00-64]	HA724561	Human	IHC-P
Relative expression	NUMB Recombinant Rabbit mAb [PSH25-86]	HA724578	Human, Rat, Green monkey	WB,FC(Intra),IP
Relative expression	NUMB Recombinant Rabbit mAb [PSH25-86] - BSA and Azide free	HA751998	Human, Rat, Green monkey	WB,FC(Intra),IP
Relative expression	Olig1 Recombinant Rabbit mAb [PSH25-70] - BSA and Azide free	HA751989	Mouse, Rat	IHC-Fr,IHC-P
Relative expression	Olig1 Recombinant Rabbit mAb [PSH25-70]	HA724557	Mouse, Rat	IHC-Fr,IHC-P
Relative expression	Parvalbumin Recombinant Antibody [JM100-08] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610439	Human, Mouse, Rat	WB,IHC-P
Relative expression	Parvalbumin Recombinant Antibody [JM100-08] - Guinea pig IgG2 (Chimeric)	HA601644	Human, Mouse, Rat	WB,IHC-P
Relative expression	Parvalbumin Recombinant Antibody [JM100-08] - Rat IgG1 (Chimeric)	HA601643	Human, Mouse, Rat	WB,IHC-Fr,IHC-P
Relative expression	Parvalbumin Recombinant Antibody [JM100-08] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610438	Human, Mouse, Rat	WB,IHC-Fr,IHC-P
Relative expression	Parvalbumin Recombinant Antibody [JM100-08] - Chicken IgY (Chimeric)	HA601645	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,FC(Intra)
Relative expression	Parvalbumin Recombinant Antibody [JM100-08] - Chicken IgY (Chimeric) - BSA and Azide free	HA610440	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,FC(Intra)
Cell treatment	Phospho-Beta Catenin (S675) Recombinant Rabbit mAb [PSH25-60] - BSA and Azide free	HA751979	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IP
Cell treatment	Phospho-Beta Catenin (S675) Recombinant Rabbit mAb [PSH25-60]	HA724534	Human, Mouse, Rat	WB,IHC-Fr,IHC-P,IP
Cell treatment	Phospho-FOXM1 (S35) Rabbit pAb	HA500642	Human	WB
Cell treatment	Phospho-Histone H2A.X (S139) Recombinant Antibody [SR33-09] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610435	Human, Mouse, Rat	WB,IHC-P,IF-Cell

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Cell treatment	Phospho-Histone H2A.X (S139) Recombinant Antibody [SR33-09] - Mouse IgG1 (Chimeric)	HA601639	Human, Mouse, Rat	WB,IHC-P,IF-Cell
	PKC Recombinant Rabbit mAb [JE77-73]	HA722847	Human, Mouse	WB,IF-Cell,FC,IP
Relative expression	PRAME Recombinant Rabbit mAb [PSH25-61] - BSA and Azide free	HA751980	Human	WB,IHC-P
Relative expression	PRAME Recombinant Rabbit mAb [PSH25-61]	HA724548	Human	WB,IHC-P
	PSGL-1 Recombinant Rabbit mAb [PSH25-73]	HA724564	Human	IHC-P
Relative expression	PTPRZ1 Recombinant Rabbit mAb [PSH25-63]	HA724550	Mouse, Rat	WB
Relative expression	PTPRZ1 Recombinant Rabbit mAb [PSH25-63] - BSA and Azide free	HA751982	Mouse, Rat	WB
	Anti-Rabbit IgG (Fc specific), F(ab') ₂ Fragment - BSA and Azide free	HA601615	Rabbit	ELISA
Relative expression	RIP3 Recombinant Rabbit mAb [PSH25-87]	HA724579	Mouse	WB,IF-Cell,IHC-P,FC(Intra),IP
Relative expression	RIP3 Recombinant Rabbit mAb [PSH25-87] - BSA and Azide free	HA751999	Mouse	WB,IF-Cell,IHC-P,FC(Intra),IP
	S100A1 Recombinant Rabbit mAb [SJ52-42]	HA722862	Human, Mouse, Rat	IHC-P
	SART1 Recombinant Rabbit mAb [PSH25-39]	HA724528	Human, Mouse, Rat	WB,IF-Cell
Relative expression	SHANK3 Recombinant Rabbit mAb [PSH25-85] - BSA and Azide free	HA751997	Mouse, Rat	WB,IHC-Fr,IHC-P
Relative expression	SHANK3 Recombinant Rabbit mAb [PSH25-85]	HA724577	Mouse, Rat	WB,IHC-Fr,IHC-P
	SLC6A2 Recombinant Rabbit mAb [PSH25-42] - BSA and Azide free	HA751974	Mouse, Rat	IHC-Fr
	SLC6A2 Recombinant Rabbit mAb [PSH25-42]	HA724530	Mouse, Rat	IHC-Fr
	Syntaxin 6 Recombinant Rabbit mAb [PSH25-41] - BSA and Azide free	HA751973	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC(Intra)
	Syntaxin 6 Recombinant Rabbit mAb [PSH25-41]	HA724529	Human, Mouse, Rat	WB,IF-Cell,IHC-P,FC(Intra)
	TFF1 Recombinant Rabbit mAb [PSH25-84]	HA724576	Human	IF-Cell,IHC-P,FC(Intra)
	TFF1 Recombinant Rabbit mAb [PSH25-84] - BSA and Azide free	HA751996	Human	IF-Cell,IHC-P,FC(Intra)
	iFluor™ 594 Conjugated TTF1 Recombinant Rabbit mAb [JE61-73]	HA724535F	Mouse	IHC-Fr,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层

Relative expression	TTF1 Recombinant Antibody [JE61-73] - Mouse IgG1 (Chimeric) - BSA and Azide free	HA610433	Human, Mouse, Rat	IF-Cell,IHC-P
Relative expression	TTF1 Recombinant Antibody [JE61-73] - Mouse IgG1 (Chimeric)	HA601636	Human, Mouse, Rat	IF-Cell,IHC-P
	USP25 Recombinant Rabbit mAb [PSH25-91] - BSA and Azide free	HA752003	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC(Intra),IP
	USP25 Recombinant Rabbit mAb [PSH25-91]	HA724583	Human, Mouse, Rat, Green monkey	WB,IF-Cell,FC(Intra),IP
	Vasopressin Recombinant Antibody [PSH11-64] - Rat IgG1 (Chimeric) - BSA and Azide free	HA610437	Mouse, Rat	IHC-P
	Vasopressin Recombinant Antibody [PSH11-64] - Guinea pig IgG2 (Chimeric) - BSA and Azide free	HA610436	Mouse, Rat	IHC-P
	Vasopressin Recombinant Antibody [PSH11-64] - Chicken IgY (Chimeric) - BSA and Azide free	HA610434	Mouse, Rat	IHC-Fr,IHC-P,IF-Tissue
	Vasopressin Recombinant Antibody [PSH11-64] - Chicken IgY (Chimeric)	HA601638	Mouse, Rat	IHC-Fr,IHC-P,IF-Tissue
	Vasopressin Recombinant Antibody [PSH11-64] - Guinea pig IgG2 (Chimeric)	HA601641	Mouse, Rat	IHC-P
	Vasopressin Recombinant Antibody [PSH11-64] - Rat IgG1 (Chimeric)	HA601642	Mouse, Rat	IHC-P
	iFluor™ 647 Conjugated Vimentin Mouse mAb [D4-B11]	HA601646F	Human	IF-Tissue,FC
	ZBTB20 Recombinant Rabbit mAb [PSH25-22]	HA724516	Mouse, Rat	IHC-P,IF-Tissue
	ZBTB20 Recombinant Rabbit mAb [PSH25-23] - BSA and Azide free	HA751970	Mouse	IHC-P,IF-Tissue
	ZBTB20 Recombinant Rabbit mAb [PSH25-22] - BSA and Azide free	HA751969	Mouse, Rat	IHC-P,IF-Tissue
	ZBTB20 Recombinant Rabbit mAb [PSH25-23]	HA724517	Mouse	IHC-P,IF-Tissue

杭州华安生物技术有限公司

电话: 0571-88062880 传真: 0571-88060400 邮箱: market@huabio.cn

地址: 杭州市钱塘区白杨街道6号大街452号高科技企业孵化园2号楼16层